BRAUN


WelchAllyn®

PRO 4000



Type 6021

ThermoScan




English	4
Français	16
Español	28



Internet:
www.welchallyn.com

6-021-360/00/VIII-07 USA/
CDN/MEX

English

The Braun ThermoScan thermometer has been carefully developed for accurate, safe and fast human body temperature measurements in the ear. The shape of the thermometer probe prevents it from being inserted too far into the ear canal which could perforate the tympanic membrane.

However, as with any thermometer, proper technique is critical to obtaining accurate temperatures. Therefore, please read all instructions carefully and thoroughly before using this product.

Important

- Do not expose the thermometer to temperature extremes (below -4°F / -20°C or over 122°F / 50°C) nor excessive humidity ($> 95\%$ RH).
- The probe window must be kept clean, dry, and undamaged at all times to ensure accurate measurements. To protect the probe window, always keep the thermometer in the storage cover while transporting or when not in use.
- This thermometer must only be used with genuine Braun ThermoScan probe covers (PC 800).
- To avoid inaccurate measurements, always attach a new, clean probe cover for each temperature measurement.
- If the thermometer is accidentally used without a probe cover attached, clean the lens (see section «Care and cleaning»).

Braun ThermoScan measures the infrared heat generated by the eardrum and surrounding tissues. To help ensure accurate temperature measurements, the sensor itself is warmed to a temperature close to that of the human body. When the Braun ThermoScan is placed in the ear, it continuously monitors the infrared energy until a temperature equilibrium has been reached and an accurate measurement can be taken.

Why measure in the ear?

Clinical studies have shown that the ear is an excellent site for temperature measurement because temperatures taken in the ear reflect the body's core temperature. Body temperature is regulated by the hypothalamus, which shares the same blood supply as the tympanic membrane. Changes in core body temperature are usually seen sooner at the tympanic membrane than at other sites, such as the rectum, mouth or under the arm.

Advantages of taking temperatures at the ear versus traditional sites:

- **Axillary** temperature measurements reflect skin temperature which may not reliably indicate the internal body temperature.
- **Rectal** temperatures often lag significantly behind internal body temperature changes, especially at times of rapidly changing temperatures. Also, there is a risk of cross contamination.
- **Oral** temperatures are often influenced by eating, drinking, thermometer placement, breathing through the mouth, or the inability of the person to close their mouth completely.

1. Guyton A C, Textbook of medical physiology, W.B. Saunders, Philadelphia, 1996, p 919
 2. Guyton A C, Textbook of medical physiology, W.B. Saunders, Philadelphia, 1996, pp 754-5
 3. Netter H F, Atlas of Human Anatomy, Novartis Medical Education, East Hanover, NJ, 1997.

Body temperature

Normal body temperature is a range. The following table shows that ranges of normal also vary by site. Therefore, readings from different sites, even if taken at the same time, should not be directly compared.

Normal ranges by site:


Axillary ¹ :	94.5 – 99.1 °F	34.7 – 37.3 °C
Oral ¹ :	95.9 – 99.5 °F	35.5 – 37.5 °C
Rectal ¹ :	97.9 – 100.4 °F	36.6 – 38.0 °C
ThermoScan ² :	96.4 – 100.4 °F	35.8 – 38.0 °C

- A person's normal temperature range tends to decrease with age. The following table shows normal ThermoScan ranges by age:

Normal ThermoScan ranges by age:

0 – 2 years	97.5 – 100.4 °F	36.4 – 38.0 °C	3 – 10 years	
97.0 – 100.0 °F	36.1 – 37.8 °C	11 – 65 years	96.6	
– 99.7 °F	35.9 – 37.6 °C	> 65 years	96.4 – 99.5 °F	

The range of normal varies from person to person and can be influenced by many factors such as time of day, level of activity, medications and gender.

1. Chamberlain, J.M., Terndrup, T.E., New Light on Ear Thermometer Readings, *Contemporary Pediatrics*, March 1994.
2. Chamberlain, J.M., et al., Determination of Normal Ear Temperature with an Infrared Emission Detection Thermometer, *Annals of Emergency Medicine*, January 1995.

Package components




- Braun ThermoScan thermometer
- Use and Care Manual
- 2 Duracell Alkaline batteries
- Additional items may be included as noted on outer packaging

Product description

1. Probe cover
2. Probe
3. Probe cover detector/automatic «on»
4. Probe cover ejector
5. Display
6. «I/O mem» button (On/memory function)
7. «ExacTemp» light
8. Start button
9. Battery door
10. Protective cover

How to use your Braun ThermoScan

1.
To achieve accurate measurements, make sure a new, clean probe cover (1) is in place before each measurement. Use the following temperature taking technique:
 2.
When the probe cover is in place, the thermometer turns on automatically. Wait for the ready signal beep.
 3.
Fit the probe snugly into the ear canal, then push and release the Start button (8).

If the probe has been properly inserted into the ear canal during the complete measuring process, a long beep will signal the end of the measuring process. The thermometer detects that an accurate temperature measurement has been taken. The result is shown on the display (5).

When you take a temperature measurement, the «ExacTemp» light (7) can be of help. It flashes during

the measuring process when the probe is correctly positioned, and lights up continuously when the thermometer detects that an accurate measurement has been taken.

4.

If the probe has not been positioned correctly in the ear canal or was moved during the measuring process, a sequence of short beeps will sound, the «ExacTemp» light will go out and the display will show an error message («POS» = position error).

5.


For the next measurement, eject the used probe cover (push ejector (4)) and put on a new, clean probe cover. The thermometer will turn on automatically.

Wait for the ready signal. Fit the probe snugly into the ear canal, then push and release the Start button.

The Braun ThermoScan ear thermometer turns off automatically after 60 seconds of inactivity. It can also be turned off by pressing the «I/O» button for at least three seconds. The display will shortly flash «OFF» and steadily display the word «OFF» after releasing the button.

Temperature taking hints

- A temperature measurement taken in the right ear may differ from a measurement taken in the left ear. Therefore, always take the temperature in the same ear.
- The ear must be free from obstructions or excess cerumen build-up in order to take an accurate measurement.
- External factors may influence ear temperatures, particularly when an individual has:
 - been lying on one ear or the other
 - had their ears covered
 - been exposed to very hot or very cold temperatures, or
 - been recently swimming or bathing.
 In these cases, remove the individual from the situation and wait 20 minutes prior to taking a temperature.
- For persons wearing hearing aids or ear plugs, remove the device and wait 20 minutes prior to taking a temperature.
- Use the untreated ear if ear drops or other ear medications have been placed in the ear canal.

98.0
MEM °F

°F ... SET ...

43.5

Memory mode

The last temperature taken before the Braun ThermoScan is turned off is stored in its memory and will be automatically displayed when it is turned on again. The display will show «MEM». For longer display of the stored measurement, the thermometer must be turned on and the «I/O» button must be pressed for at least 1 second. The memory mode is left automatically after pressing the «I/O» button for at least 1 second.

Changing the temperature scale

Your Braun ThermoScan is shipped with the Fahrenheit temperature scale activated. If you wish to switch to Celsius «°C», proceed as follows:

Make sure the thermometer is turned off.

Press and hold down the «I/O» button. After about 3 seconds the display will show this sequence: «°F» / «SET» / «°C» / «SET» ... Release the «I/O» button when «°C» is shown.

There will be a short beep to confirm the new setting, then the thermometer is turned off automatically. When turned on again, temperatures will be shown in Celsius, until the above procedure is repeated.

Base station for storage, charging, anti-theft function (optional)

For your Braun ThermoScan PRO 4000, there is a Base Station available. Contact your distributor for details (see page 14).

The station automatically recharges the thermometer (when using the rechargeable batteries supplied with the base station), and it has an electronic, individually adjustable anti-theft function. It also serves as a convenient storage base (can be wall mounted).

It requires returning the thermometer to the station within an individually pre-selected time. Otherwise the thermometer will be locked.

When leaving the factory, the anti-theft function of the PRO 4000 is turned off.

If the anti-theft function is enabled, and when not in the measuring process (stand-by mode) the thermometer will show the hours remaining before it is locked.

To turn on the anti-theft function, please refer to the base station manual.

If the thermometer shows «RET» (without the battery symbol) on its display, the anti-theft time has elapsed and the thermometer must be returned to the station for unlocking.

If the anti-theft function is not enabled, the display will show «off» when in stand-by mode.

Care and cleaning

Braun ThermoScan probe covers are intended for single use only!

To ensure accurate measurements, always use Braun ThermoScan disposable probe covers with the Pro 4000 thermometer. Damaged, perforated, or soiled probe covers can lead to inaccurate temperature measurements. Additionally, the accuracy of your temperature measurement can be affected by damage to the probe window, or the presence of dirt or cerumen on the probe window.

Probe Window


To assure a high level of accuracy, it is very important to keep this window clean, dry and undamaged. Fingerprints, cerumen, dust and other soiling compounds reduce transparency of the window and result in lower temperature measurements. If the thermometer is accidentally placed in the ear without a probe cover installed, the window must be cleaned immediately. Additionally, the window and/or probe tip area should never be modified, changed, or adjusted. These changes will affect calibration and accuracy of the thermometer.

To clean the window, gently wipe its surface with a cotton swab slightly moistened with alcohol and immediately wipe dry with a clean cotton swab. After cleaning, allow at least 5 minutes drying time before taking temperatures. Make sure probe window is clean and dry. Avoid touching the window except when cleaning is required. If the window is damaged, the thermometer must be returned to Welch Allyn for service.

Note: Do not use any chemical other than alcohol to clean the probe window.

Thermometer

Use a soft cloth slightly moistened with alcohol to clean the thermometer display and exterior. Do not use abrasive cleaners. Never submerge the thermometer in water or any other liquid. Store thermometer and probe covers in a dry location (the thermometer is not protected against ingress of water), free from dust and contamination and away from direct sunlight. The ambient temperature at the storage location should remain fairly constant and within the range of 50 °F to 104 °F (10 °C to 40 °C).


Trouble-shooting

Error message

Situation

Solution

No probe cover is attached.

Attach new, clean probe cover.


(flashing)

(POS = position error)

The infrared monitor cannot find a temperature equilibrium and allows no measurement.

Clear the display by pushing the «Start» button once.

Ensure that the positioning of the probe is correct and remains stable.

Ambient temperature is not within the allowed operating range 50–104 °F (10–40 °C) or changing too rapidly.

Allow the thermometer to remain for 30 minutes in a room where the temperature is between 50 and 104 °F / 10 and 40 °C.


Temperature taken is not within typical human temperature range (93.2–108 °F (34–42.2 °C)).

HI = too high

Make sure the probe tip and lens are clean and a new, clean probe cover is attached. Make sure the thermometer is properly inserted into the ear canal.

Then, take a new temperature.

LO = too low

System error

Wait 1 minute until the thermo-meter turns off automatically, then turn on again.

If error persists,

... reset the thermometer by re-moving the batteries and putting them back in.

If error still persists,


... contact Welch Allyn at 1-800-535-6663 (in Canada at 1-800-561-8797).

Trouble-shooting

Error message


Situation

Solution

Battery is low, but thermometer will still operate correctly.

Insert new batteries.



Battery is too low to take correct temperature measurements.

Insert new batteries.

Do you have any further questions?

... contact Welch Allyn at 1-800-535-6663 (in Canada at 1-800-561-8797).

Replacing the batteries

The thermometer is supplied with two 1.5 V batteries type AA (LR 6). For best performance, we recommend Duracell® alkaline batteries. Insert new batteries when the battery symbol appears on the display.

Open the battery compartment. Remove the batteries and replace with new batteries, making sure the poles are in the right direction. Slide battery door back until it snaps in place.

To protect the environment, dispose of empty batteries according to national or local regulations.

Calibration

The thermometer is initially calibrated at the time of manufacture. If this thermometer is used according to the use instructions, periodic readjustment is not required. If at any time you question the accuracy of temperature measurements, please contact Welch Allyn (see page 14).

Manufacturing date is given by the LOT number located in the battery compartment. The first number after LOT represents the last digit of the year of manufacture. The next three digits give the day of the year of manufacture.

An example: LOT 0116 xx xx – this product has been manufactured on the 116th day of the year 2000.

Special situations

Clinical studies have shown the Braun ThermoScan thermometer can obtain accurate temperature measurements on persons of all ages. However, there are special situations when an ear thermometer should not be used. These include, but may not be limited to the following situations:

- Do not use an ear thermometer if there is blood or drainage in the external ear canal.
- An ear thermometer should not be used on a patient who exhibits symptoms of an acute or chronic inflammatory condition of the external ear canal.
- Common situations like the presence of moderate amounts of cerumen in the ear canal, otitis media and tympanostomy tubes do not significantly impact temperature readings. However, complete ear canal occlusion due to cerumen can result in lower temperature readings.
- For persons using hearing aids or ear plugs, remove the device and wait 20 minutes prior to taking a temperature.
- Use the untreated ear if prescription ear drops or other ear medications have been placed in the ear canal.
- Patients who have deformities of the face and/or ear may not be able to have a temperature taken with an ear thermometer.

Fever facts

Many persons may not have an elevated temperature even when they are ill. These include, but are not limited to infants under 90 days of age, persons taking steroids, antibiotics or antipyretics (acetaminophen, ibuprofen, aspirin), persons with compromised immune systems, including the elderly or persons with some chronic illnesses.

Fever is described as an elevation of body temperature over an individual's «normal» temperature.

An elevated temperature or fever is often viewed as a danger sign. In fact, a fever can be very beneficial, and helps our immune system work more effectively. A fever should be evaluated along with other physical symptoms. With the exception of newborn infants, the presence of fever, without any other symptoms of illness, or in a child who is behaving normally may not be cause of concern.

Some people, such as the elderly, infants under 90 days of age, those with compromised immune systems, or persons who take steroids for example are often unable to build a response to illness or environmental conditions. These individuals may not have a fever when they are ill, or their elevated temperature may be lower than expected for the severity of their illness. Other medications such as anti-inflammatories and some analgesics may also mask fever.

The presence or absence of fever should not be used as the only measurement of the significance of the illness.

Febrile seizures, or convulsions, which usually occur in children six months to six years of age, are thought to occur not because a fever is present, but because of the rapid rise of the child's temperature.

Product specifications

Displayed temperature range:	68–108 °F (20–42.2 °C)
Operating ambient temperature range:	50–104 °F (10–40 °C)
Display resolution:	0.1 °F or °C
Accuracy for displayed temperature range:	±0.4 °F (95.9–107.6 °F) (35.5 °C–42 °C) ±0.5 °F (outside this temperature range)

Long term storage ranges

Temperature:	–4 °F to 122 °F (–20 °C to 50 °C)
Humidity:	95% non condensing
Shock:	withstands drop of 3 feet
Measurement time:	3–7 seconds
Automatic «power down»:	60 seconds
Battery life:	6 months/1000 measurements
Battery type:	2 x MN 1500 or 1.5 V AA (LR6)
Thermometer dimension:	6" x 1.7" x 1.3" (152 mm x 44 mm x 33 mm)
Thermometer weight:	3.6 oz (100 g) without batteries

This infrared thermometer meets requirements established in ASTM Standard E 1965-98 (for the thermometer system [thermometer with probe cover]). Full responsibility for the conformance of the product to the standard is assumed by Kaz, Inc. Hudson NY, 12534 USA.

ASTM laboratory accuracy requirements for the thermometer only (not including a probe cover) in the display range of 96.8 °F to 102.2 °F (36 °C to 39 °C) for infrared thermometers is ± 0.4 °F (± 0.2 °C), whereas for mercury-in-glass and electronic thermometers, the requirement per ASTM Standards E 667-86 and E 1112-86 is ± 0.2 °F (± 0.1 °C).

A clinical summary is available upon request.

This device conforms to the following standards:
 DIN EN 60601-1: 3/96 «Medical electrical equipment» –
 Part 1: General requirements for safety
 DIN EN 12470-5: 2003 «Clinical thermometers» –
 Part 5: Performance of infrared ear thermometers (with maximum device)

The Braun ThermoScan thermometer has been proven to be safe. To ensure the best results, please take the time to read this manual completely and keep it handy for future reference. In the unlikely event you experience any difficulty using the Braun ThermoScan thermometer, simply call Welch Allyn toll free at: 1-800-535-6663 (in Canada:1-800-561-8797).

This products conforms to the provisions of the EC directive 93/42/EEC (Medical Device Directive).

TYPE BF Attention, consult
EQUIPMENT ACCOMPANYING
DOCUMENTS

EAR THERMOMETER
WITH RESPECT TO ELECTRICAL SHOCK, FIRE
MECHANICAL HAZARDS ONLY
IN ACCORDANCE WITH UL60601 – 1/10 HK

Internally Powered Equipment
Continuous Operation
Not protected against ingress of water
Patents Pending

Service

The Braun ThermoScan thermometer has a limited three year warranty. Should service be required during or after the warranty period you must call Welch Allyn's Customer Service Department at 1-800-535-6663 (in Canada: 1-800-561-8797) to obtain a return authorization number and shipping address. Repackage the thermometer carefully in its original packaging or securely packed to avoid damage during shipping. Include the original sales slip indicating the date of purchase, a note describing the problem, your authorization number, and your return address. Send the Braun ThermoScan thermometer prepaid and insured.

Note: Shipping address varies by country.

Answers to any questions regarding operation of the Braun ThermoScan thermometer may be obtained by calling Welch Allyn toll-free at:

1-800-535-6663 (in Canada: 1-800-561-8797) Monday–Friday 8:00 A.M. to 8:00 P.M. Eastern Time www.welchallyn.com

Limited three year manufacturer's warranty

For Model Pro 4000

The PRO 4000 instrument is warranted by Kaz to be free from defects in material and workmanship under normal use and service for a period of three years from the date of delivery to the first professional user who purchases the instrument.

This warranty does not cover batteries, damage to the probe window, or damage to the instrument caused by misuse, negligence or accident, and extends to only the first purchaser of the product. Additionally this warranty becomes void, if the thermometer is operated with anything other than genuine Braun ThermoScan® Probe Covers.

Thermometer service or replacement will be carried out by Kaz subject to the terms of this warranty. Loss or damage in return shipment to Kaz shall be at purchaser's risk. Kaz reserves the right to replace a non-conforming product with the most comparable product currently available. This product has a five year minimum expected lifetime.

IN NO EVENT SHALL KAZ BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH THE PURCHASE OR USE OF ANY PRODUCT. This warranty shall not apply to, and Kaz shall not be responsible for, any loss arising in connection with the purchase or use of any product which has been serviced by anyone other than a Kaz authorized service representative or altered in any way so as, in Kaz's judgment, to affect its stability or reliability, or which has been subject to misuse or negligence or accident, or which has the serial number or lot number altered, effaced or removed, or which has been used other than in accordance with the instructions furnished by Kaz. IN NO EVENT SHALL KAZ'S LIABILITY EXCEED THE ORIGINAL COST OF THE PRODUCT.

This warranty is in lieu of all other warranties, express or implied, and of all other obligations or liabilities on Kaz's part, and Kaz neither assumes nor authorizes any representatives or other person to assume for it any other liability in connection with the sale of its products.

TO THE EXTENT PERMITTED BY LAW, KAZ DISCLAIMS ALL STATUTORY AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

Manufactured by / Fabriqué par / Fabricado por
Kaz Canada, Inc.
510 Bronte St. S.
Milton, ON L9T 2X6

Distributed by/distribué par:

USA

Welch Allyn, Inc.
4341 State Street Road, P.O. Box 220
Skaneateles Falls, NY 13153-0220

Canada

Welch Allyn Canada Ltd.
160 Matheson Blvd. East, Unit 2
Mississauga, Ontario L4Z 1V4

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

1. Guyton A C. Textbook of medical physiology, W.B. Saunders, Philadelphia, 1996, p 919
 2. Guyton A C. Textbook of medical physiology, W.B. Saunders, Philadelphia, 1996, pp 754-5
 3. Netter H F, Atlas of Human Anatomy, Novartis Medical Education, East Hanover, NJ, 1997, pp 63, 95

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1-800-535-6663 (en Canadá: 1-800-561-8797)

Lunes-viernes

De 8:00 a.m. a 8:00 p.m. Hora oriental www.welchallyn.com

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A person's normal temperature range tends to decrease with age. The following table shows normal ThermoScan ranges by age.

Normal ThermoScan temperature ranges:

0 – 2 years	97.5°F – 100.4°F	36.4°C – 38.0°C	3
– 10 years	97.0°F – 100.0°F	36.1°C – 37.8°C	11
– 65 years	96.6°F – 99.7°F	35.9°C – 37.6°C	
> 65 years	96.4°F – 99.5°F	35.8°C – 37.5°C	

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How to use Braun ThermoScan

1.

To help to ensure accurate readings, make sure the probe tip and lens are clean and a new, clean probe cover is in place before each measurement.

2.

When the probe cover is in place, the thermometer turns on automatically. Wait for the ready signal beep.

Fit the probe snugly into the ear canal, then push and release the Start button.

4.

If the probe has been fitted correctly into the ear canal during the complete measuring process, a long beep will signal the end of the measuring process. You can be assured to have taken an accurate temperature reading. The thermometer detects that an accurate temperature measurement has been taken.

If you take the temperature of another person, the «ExacTemp» light can be of help. It flashes during the measuring process as long as the probe is correctly positioned, and lights up continuously when the thermometer detects that an accurate reading has been taken.

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