# ArthroTunneler™ TunnelPro™ System

It is important to read the Instructions for Use in its entirety prior to using the ArthroTunneler™ TunnelPro™ System

Manufacturer:

Rx Only

**T.A.G. Medical Products Corporation LTD.**

Kibbutz Gaaton 25130

Israel

Tel: +972-4-9858400

Fax: +972-4-9858404

**Distributed by:  
Tornier Inc.**

Edina, MN 55435

USA: 1-888-867-6437

International: 1 281-494-7900

Authorized Representative:

**MEDNET GmbH**

Borkstrasse 10, 48163 Muenster

Germany

Fax: +49 251 3226622

P/N: IFU551tr Rev: A

2012/05

**Table of Contents**

English..... 1  
Français ..... 9  
Deutsch..... 17  
Italiano ..... 25  
Português..... 33  
Español..... 41  
Nederlands ..... 49  
日本語..... 57

# English

## **Table of Contents**

Description.....	2
Indications.....	2
Contraindications .....	2
Warnings.....	3
Precautions .....	3
Use of the ArthroTunneler™ TunnelPro™ System .....	3
Cleaning, Disinfection and Sterilization - Reusable Instruments Only.....	4
Limitations on Reprocessing .....	4
Containment and Transportation .....	4
Maintenance, Inspection and Testing .....	5
Packaging for Sterilization.....	5
Storage .....	6
How supplied .....	6
Symbols.....	7

## Description

### Single Patient Use /Sterile kit

ArthroTunneler™ TunnelPro™ System Cat. # SMB100101 consists of the following:

- |                  |  |
|------------------|--|
| Cat. # SMB100201 | TunnelPro™ (box contains 2 units) - see figure 1 below<br>The TunnelPro™ is provided preloaded on a disposable inserter. |
| Cat. # SMB000101 | Arthroscopic Transosseous Tunneling Device - see figure 2 below<br>(Hereafter ArthroTunneler™ Device).                   |

### Reusable Instrument Set / Non-Sterile (supplied separately)

- |                  |                      |
|------------------|----------------------|
| Cat. # SMB000201 | AT Drill/Punch Guide |
| Cat. # SMB000301 | AT Obturator         |
| Cat. # SMB000401 | AT 2.9 Drill - M     |
| Cat. # SMB000501 | AT 2.9 Punch - M     |
| Cat. # SMB000601 | AT 2.5 Drill - L     |
| Cat. # SMB000701 | AT Suture Inserter   |


Fig. 1: The TunnelPro™


Fig. 2: The ArthroTunneler™ device

## Indications

The ArthroTunneler™ TunnelPro™ System is intended to be used for Rotator Cuff repair procedures.

The TunnelPro™ is intended to protect the suture/bone during transosseous fixation procedures.

The ArthroTunneler™ is intended for transosseous fixation of tendons to bone in the following applications: Rotator Cuff repair.

## Contraindications

The ArthroTunneler™ TunnelPro™ System is NOT intended for use in:

- Pathological conditions of the soft tissue to be repaired or reconstructed which would adversely affect suture fixation.
- Insufficient quantity and/or quality or cortical bone integrity that would impair the ability to securely fix the TunnelPro™ implants.
- Compromised bone surface that would not allow secure fixation of the TunnelPro™.
- Patients with an active infection.
- Physical conditions that would retard healing, such as blood supply limitation and infection.
- Conditions which tend to limit the patient's ability or willingness to follow instructions during the healing period.

## **⚠ Warnings**

- Read these instructions prior to use of the ArthroTunneler™ TunnelPro™ System.
- Prior to use, inspect the device package and the product itself for signs of damage or tampering. If tampering or damage exists, do not use.
- The TunnelPro™ and the ArthroTunneler™ device are supplied sterile, and are intended for single patient use only:
  - Do not resterilize / reuse.
  - Do not use if sterile barrier is compromised. Discard any open, unused product.
  - Do not use after the expiration date.
  - Do not use if the label is missing.
  - Reprocessing may result in biocontamination, degraded performance or loss of functionality. The single use devices are not designed or validated to be cleaned, disinfected or sterilized by the user.
- The reusable instruments are supplied non-sterile. Clean, disinfect and sterilize the instruments prior to the first use and after every subsequent use.
- Handle with care: The Drills and the Drill Guide in the Set have sharp cutting edges.

## **Precautions**

- It is imperative that the surgeon and operating room staff are fully conversant with the appropriate surgical technique for using each component of the ArthroTunneler™ TunnelPro™ System.
- Inspect the instruments before use for integrity and compatibility in order to ensure proper functionality and safety. Do not use the ArthroTunneler™ TunnelPro™ System if they appear to be damaged and notify your Tornier representative.
- Test the ArthroTunneler™ device before surgical use by operating the Anvil and Loop actuators. Verify free movement of the Anvil and the Loop.
- In case repeated insertion of the implant is required, care shall be taken in order to avoid tunnel enlargement.
- Ensure TunnelPro™ implant is threaded to its inserter.
- Replace all dull instruments.

## **Use of the ArthroTunneler™ TunnelPro™ System**

1. Load a free strand of suture on the forked distal end of the AT Suture Inserter and secure the suture tails on the proximal end of the AT Suture Inserter handle.
2. Clean the area prior to drilling in order to enable clear identification of the tunnel when inserting the TunnelPro™ implant(s).
3. Use the AT Drill Guide either through a standard anterior portal (with the arm in adduction and internal rotation to provide controlled access), or from an accessory portal or stab incision adjacent to the acromion. Create a medial tunnel, using either the AT 2.9 Drill – M, or the AT 2.9 Punch – M to the depth indicated by the positive stop.
4. Mark the location of the medial hole for identification.
5. From the lateral portal, gently introduce the ArthroTunneler™ device to the superior aspect of the greater tuberosity. During insertion, keep the Anvil Actuator in the extended (normal) position.
6. Withdraw the Anvil by moving the Anvil Actuator to the proximal position and fully introduce the Hook into the medial tunnel.
  - During this step, position the Anvil on the lateral portion of the greater tuberosity.
  - Check for free movement of the Loop by sliding the Loop Actuator between the Deploy and Retrieve positions.
  - The superior beam of the Hook should be approximately parallel to the plateau of the tuberosity.
7. Deploy the Loop of the ArthroTunneler™ device by sliding the Loop Actuator to the Deploy position.
8. Drill the lateral tunnel through the Loop deployed in the medial tunnel using the AT 2.5 Drill - L:
  - Introduce the drill through the hole at the proximal end of the ArthroTunneler™ device.
  - Drill to the full depth of the drill as indicated by the positive stop.
  - Move the loop actuator toward the Retrieve position; the actuator should be moved only partially, indicating that the drill reaches the intersection.
9. Remove the AT 2.5 Drill - L, while maintaining the position of the ArthroTunneler™ device.
 

**Caution:** Changing the position of the ArthroTunneler™ device at this stage may result in losing the entry point into the lateral hole.
10. It is recommended to use tipped tail suture # 2 or # 3/4.
 

Using more than 3 strands may cause difficulty in threading the sutures through the TunnelPro™ implants and inserters.
11. Insert the loaded Suture Inserter through the hole at the proximal end of the ArthroTunneler™ device into the lateral bone tunnel until the tab of the Suture Inserter is flush with the device.
12. Move the Loop Actuator toward the Retrieve position. The Actuator should move only partially, indicating that the Suture Inserter is captured within the Loop.
13. Release the suture ends from the handle and remove the Suture Inserter only, leaving the suture in the lateral hole.
14. Slide the Loop Actuator to the Retrieve position, retract the Anvil and remove the ArthroTunneler™ device together with the suture captured from the medial hole.

15. Once the suture protrudes from the medial hole, slide the Loop Actuator to the Deploy position and release the suture from the device.  
**Note:** The loop of suture retrieved through the medial hole can also be used as a shuttle to pass multiple sutures through the tunnels simultaneously.
16. Using a suture manipulator, grasp the suture tail exiting the lateral tunnel and pull them through the cannula.
17. Load the TunnelPro™ with the end(s) of the suture(s) that exit the lateral tunnel.
18. Secure the end(s) of the suture(s) exiting the medial tunnel to avoid threading it over the TunnelPro™ implant as well as pulling out of the tunnel.
19. With slight tension on the suture tail slide the TunnelPro™ implant along the suture down the cannula to the opening of the lateral tunnel.
20. Follow the drilling direction when inserting the TunnelPro™ implant(s).
21. Use a mallet to advance the TunnelPro™ into the lateral tunnel until the bottom surface of the head of the TunnelPro™ implant is flush with the bone.
22. Disengage the inserter from the TunnelPro™ by turning 1/2 turn counterclockwise and pulling only the inserter.
23. If implant manipulation is desired reengage the TunnelPro™ inserter with the implant by placing the TunnelPro™ inserter into the implant and rotate inserter 1/2 turn clockwise.  
**Note:** Use only the TunnelPro™ inserter to reposition the TunnelPro™ implant.
24. Pass the suture(s) through the tendon and tie securely.
25. Repeat steps for additional tunnels as needed or indicated.

### **Cleaning, Disinfection and Sterilization - Reusable Instruments Only**

These reprocessing instructions are provided in accordance with ISO 17664. While they have been validated by the manufacturer of the medical device as being capable of preparing the instruments for re-use, it remains the responsibility of the processor to ensure that the reprocessing as actually performed, using equipment, materials, and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process.

#### **⚠ Warnings**

- The following cleaning, disinfection and sterilization methods have been validated. Failure to follow reprocessing instructions as provided may result in incomplete decontamination.
- Wear appropriate protective equipment when handling contaminated devices (gloves, eye protection, etc.).

#### **Limitations on Reprocessing**

1. Repeated reprocessing has minimal effect on instrument life and function. End of useful instrument life is generally determined by wear or damage from handling or surgical use.
2. Do not leave the instruments in solutions longer than necessary, as this may accelerate normal product aging.
3. Damage incurred by improper reprocessing will not be the responsibility of T.A.G. Medical Products Corporation Ltd.

#### **Containment and Transportation**

1. Reprocess the instruments as soon as reasonably practical following use.
2. Transport the instruments in an appropriate manner to prevent damage.

#### **Cleaning**

##### **⚠ Warnings**

1. Cleaning should be performed as soon as possible<sup>(1)</sup> after each procedure, before blood, saline and debris are dry.
2. Do not use metallic brushes, scrub pads, or other abrasive cleaning aids, as these are likely to damage the instruments.
3. Do not use highly alkaline solutions or harsh chemicals such as chlorine, caustic soda, and organic or ammoniated acids or solvents (such as acetone), as these are likely to damage the instruments.
4. Follow the instructions and warnings issued by the suppliers of any cleaning and disinfection agents and equipment used.
5. Complex devices, such as those with tubes, hinges, retractable features, and matte or textured surface finishes, require special attention during cleaning.

<sup>(1)</sup> 30 minutes was used for the Cleaning Validation.

#### **Manual cleaning**

1. Rinse instruments with distilled or demineralized water.
2. Prepare an enzymatic detergent<sup>(2)</sup> according to the manufacturer's instructions.
3. Submerge the instruments in the enzymatic detergent solution for 20 minutes. Use a 5 mL syringe and needle to flush enzymatic detergent through lumens, and use a soft bristled brush to scrub the lumens, crevices, and slots.
4. Rinse instruments with distilled or demineralized water, agitating for one minute; repeat rinse twice.
5. Dry all instruments thoroughly before storage.

6. Inspect the instruments after each cleaning cycle to ensure proper functionality and safety:
  - 6.1. If you suspect a problem, do not use the instrument.
  - 6.2. Inspect all components for cleanliness: activate all moving parts to their maximum positions while inspecting for blood or debris. Place a pipe cleaner down any cannulation or crevice, and then inspect the pipe cleaner for blood or debris. If blood or debris is present, repeat the above cleaning procedure.

<sup>(2)</sup> Cleaning efficacy was validated using ENZOL<sup>®</sup> at one ounce per gallon of water at 35-40°C (95-104°F).

### Automated Cleaning

1. Prepare an enzymatic detergent<sup>(3)</sup> according to the manufacturer's instructions.
2. Perform pre-cleaning by flushing an enzymatic detergent through the lumens using a 5 mL syringe and needle and using a soft bristled brush to scrub the lumens, crevices, and slots.
3. Rinse the instruments with distilled or demineralized water until no visible detergent is pre-sent.
4. Place the instruments in the washer in a manner to facilitate draining.
5. Clean in a validated automatic washer using the "INSTRUMENTS" cycle and detergents<sup>(4)</sup> intended for use with automated cleaning. Effective automated cleaning can be achieved using the following parameters:
  - 5.1. Cold water wash (minimum of 35 seconds at 27°C / 80°F)
  - 5.2. Hot water wash (minimum of 2 minutes at 60°C / 140°F)
  - 5.3. Cold water wash (minimum of 45 seconds at 27°C / 80°F)
  - 5.4. Hot water wash (minimum of 2 minutes at 60°C / 140°F)
  - 5.5. Thermal rinse (minimum of 1 minute at 82.2°C / 180°F)
  - 5.6. Drying (minimum of 6 minutes at 82.2°C / 180°F)
6. Inspect the instruments after each cleaning cycle to ensure proper functionality and safety:
  - 6.1. If you suspect a problem, do not use the instrument.
  - 6.2. Inspect all components for cleanliness: activate all moving parts to their maximum positions while inspecting for blood or debris. Place a pipe cleaner down any cannulation or crevice, and then inspect the pipe cleaner for blood or debris. If blood or debris is present, repeat the above cleaning procedure.

<sup>(3)</sup> Validation was performed using ENZOL<sup>®</sup>.

<sup>(4)</sup> Validation was performed using Renuzyme Plus<sup>®</sup> and Tec Wash III<sup>®</sup>.

### Disinfection

#### Warnings

Disinfection alone is not adequate to render the instruments sterile. Sterilization must be performed before each procedure.

1. After cleaning, disinfection should be performed by immersion in 2.4% CIDEX<sup>®(5)</sup> (or equivalent) with a minimum immersion time of 45 minutes at 25°C / 77°F following the manufacturer's instructions.
2. Rinse with distilled or demineralized water, agitating for one minute; repeat rinse twice.
3. Dry all instruments thoroughly before storage.

<sup>(5)</sup> Validation was performed using CIDEX<sup>®</sup>.

### Maintenance, Inspection and Testing

Inspect the instruments for damage. Ensure proper functionality. If a problem is observed or suspected, contact your Tornier representative for replacement or repair.

### Packaging for Sterilization

Use no more than two layers of sterilization wrap.

### Sterilization

#### Warnings

- Sterilization is not a substitute for cleaning. Instruments must be thoroughly cleaned prior to sterilization, or incomplete decontamination will result.
- The sterilization parameters provided for use with a sterilization tray might differ from those provided for the instruments alone. When sterilizing the instruments inside a tray, select the more stringent set of parameters.
- These devices were not validated for ethylene oxide (EO) or gas plasma sterilization methods. Use only steam (autoclave) sterilization as outlined in these instructions, or incomplete sterilization may result.
- Drying time depends on several variables, including altitude, humidity, type of wrap, preconditioning, size of chamber, mass of load, material of load, and placement in chamber. Users must verify that drying time set in their autoclave yields dry surgical equipment.
- The adequacy of any sterilization procedure must be tested. To substantiate the process and its reliability and reproducibility, it is critical that appropriate process parameters be validated for each facility's sterilization equipment and product/load configuration by persons who have training and expertise in sterilization processes.

- The following autoclave processes have been validated by the manufacturer; however, the sterility assurance level for any given sterilization cycle is dependent on the actual conditions selected by the user and the correct calibration of the user's sterilization equipment and instrumentation. Because these factors are not within the control of T.A.G. Medical Products Corporation Ltd., the user is responsible for the sterility assurance level.

#### Gravity - Displacement Steam Sterilization:

- ◊ Sterilizer Type: Gravity
- ◊ Minimum Temperature: 132° C (270°F)
- ◊ Minimum Exposure Time: 10 minutes
- ◊ Minimum Dry Time: 40 minutes

#### Pre-vacuum Steam Sterilization:

- ◊ Sterilizer Type: Pre-vacuum
- ◊ Minimum Temperature: 132° C (270°F)
- ◊ Minimum Exposure Time: 4 minutes
- ◊ Minimum Dry Time: 15 minutes

#### Pre-vacuum Steam Sterilization (alternate cycle):

- ◊ Sterilizer Type: Pre-vacuum
- ◊ Minimum Temperature: 134° C (273°F)
- ◊ Minimum Exposure Time: 3 minutes<sup>(6)</sup>
- ◊ Minimum Dry Time: 15 minutes

<sup>(6)</sup> Exposure time of 18 minutes is permitted (as mandatory in some European countries).

Note: This cycle has not been validated for prion inactivation in the US.

Do not exceed 140°C (284°F) during sterilization.

## **Storage**

- Between uses, store instruments in a clean, cool and dry environment.
- To prevent damage, minimize contact between instruments and with other equipment.

**Warning:** These storage parameters are not intended to preserve instrument sterility. If necessary, reprocess the instruments again prior to use.

The Tornier logo is a registered trademark of Tornier Inc.

ENZOL® and CIDEX® are registered trademarks of Advanced Sterilization Products.

Renuzyme Plus® and Tec Wash III® are registered trademarks of Getinge USA, Inc.

## **How supplied**

- The ArthroTunneler™ device is provided sterile for single patient use. Discard after use.
- All other instruments in the Set are reusable and are provided non-sterile.
- The TunnelPro™ is provided sterile for single patient use.

## Symbols

	Caution / Attention: See Instructions for Use
	Do not reuse
	Do not re-sterilize
	Non-Sterile
	Use by
	Manufacturer
	Lot number
	Authorized Representative in the European Community
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Store in a cool dry place
	Not made with natural rubber latex
	Do not use if package is damaged or open
	Catalog number
	Date of manufacture
	Number of units
	Sterilization using irradiation
	CE Mark (Class I)
	CE Mark (Class IIa)

The ArthroTunneler™ TunnelPro™ System is a trademark of Tornier Inc.

This product is covered by one or more patents pending in the US and worldwide.

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