K033437 1 q 5

Attachment V

510(k) Summary

1.General Information

Submitter:

AllMed Systems Inc. 9232 Klemetson Drive

Pleasanton CA 94588

Phone:

925-468-0433

Fax

925-399-5984

Contact Person

Peter Allen

Date Prepared

27th October 2003

2. Names

Device Name

Sphinx

Common Name

2.1micron Laser System

Classification Name

Laser Surgical Instrument and accessories

3. Predicate Device

Lumenis – Holmium VersaPulse Power Suite Holmium Laser Trimedyne - OmniPulse Max 80 Watt Holmium Laser System

4. Product Description

The Sphinx system is a surgical laser system operating at a wavelength of 2.1 micron. The purpose of the laser is the ablation, coagulation, dissection and resection of soft tissue. The laser is designed for open surgery and surgical applications in aqueous media. The laser power is delivered via standard silica laser fibers. The distal tip is guided by a handpiece or endoscopic surgical instrument.

K033437 295

It consists of:

Laser Console with Internal Computer Control Panel and Display A fiber optic delivery system Footswitch

5. Indications for Use

The Sphinx laser system and its fiber optic delivery system are intended for use in surgical procedures using open, laparoscopic and endoscopic inclsion, excision, resection, ablation, vaporization, coagulation and hemostasis of soft tissue in use in medical specialties including: Urology, Urinary Lithotripsy, Gasteroenterology, Arthroscopy, Discetomy Pulmonary, Gynecology, ENT, Dermatology, Plastic Surgery and General Surgery.

Urology

Open and endoscopic surgery (incision, excision, resection, ablation, vaporization, coagulation and hemostasis) including:

Urethral Strictures

Bladder Neck Incisions (BNI)

Ablation and resection of Bladder Tumors, Uretheral Tumors and Ureteral Tumors,

Ablation of Benign Prostatic Hypertrophy (BHP),

Transurethral incision of the prostate (TUIP)

Holmium Laser Resection of the Prostrate (HoLRP)

Holmium Laser Enuculeation of the Prostate (HoLÉP)

Holmium laser Ablation of the Prostate (HoLAP)

Condylomas

Lesions of external genitalia

Lithotripsy and Percutaneous Urinary Lithotripsy

Endoscopic fragmentation of urethral, ureteral, bladder and renal calculi including cystine, calcium oxalate, monohydrate and calcium oxalate dehydrate stones.

Endoscopic fragmentation of kidney calculi

Treatment of distal impacted fragments of steinstrasse when guide wire cannot be passed.

Gasteroenterlogy

Open and endoscopic gasteroenterlogy surgery (incision, excision, resection, ablation, vaporization, coagulation and hemostasis) including:

Appendectomy

Polyps

Biopsy

K033437

Gall Bladder calculi Biliary/Bile duct calculi

Ulcers

Gastric ulcers

Duodenal ulcers

Non Bleeding Ulcers

Pancreatitas

Hemorrhoids

Cholecystectomy

Benign and Malignant Neoplasm

Angiodysplasia

Colorectal cancer

Telangiectasias

Telangiectasias of the Osler-Weber-Renu disease

Vascular Malformation

Gastritis

Esophagitis

Esophageal ulcers

Varices

Colitis

Mallory-Weiss tear

Gastric Erosions

Arthroscopy

Arthroscopy/Orthopedic surgery (excision, ablation and coagulation of soft and cartilaginous tissue) in small and large joints of the body, excluding the spine but including:

Ligament and tendon Release

Countouring and sculpting of articular surfaces

Capsulectomy in the Knee

Chondreplasty in the Knee

Debridement of inflamed synovial tissue

Chondromalacia Ablation

Chondromalacia and tears

Plica Removal

Meniscectomy

Loose Body Debridement

Lateral retinecular release

Ablation of soft, cartilaginous and bony tissue in Minimal Invasive Spinal Surgery including

Percutaneous Laser Disc Decompression/Discectomy of the L4-5 and L5-

S1 lumbar discs, including Foraminoplasty

Percutaneous Cervical Disc Decompression/Discectomy

Percutaneous Thoracic Disc Decompression/Discectomy

Pulmonary

Open and endoscopic pulmonary surgery (incision, excision, resection, ablation, vaporization, coagulation and hemostasis of soft tissue)

Gynecology

Open and laparoscopic gynecological surgery (incision, excision, resection, ablation, vaporization, coagulation and hemostasis) of soft tissue

ENT

Endoscopic endonasal surgery (incision, exclsion, resection, ablation, vaporization, coagulation and hemostasis of soft tissue and cartilage) including:

Endonasal/sinus Surgery
Partial turbinectomy
Polypectomy
Dacryocystorhinostomy
Frontal Sinusotomy
Ethmoidectomy
Maxillary antrostomy
Functional endoscopic sinus surgery

Dematology and Plastic Surgery

Incision, excision, resection, ablation, vaporization, coagulation and hemostasis of soft,mucosal, fatty and cartilaginous tissue, in therapeutic plastic, dermatologic and aesthetic surgical procedures including:

Basal Cell Carcinomas
Lesions of skin and subcutaneous tissue
Skin tags
Plantar warts
Lesions of skin and subcutaneous tissue
Port Wine Stains
Papillomas

General Surgery

Open, laparoscopic and endoscopic surgery (incision, excision, resection, ablation, vaporization, coagulation and hemostasis) including:

Appendectomy Skin incision

Excision of external and internal lesions

Complete of partial resection of internal organs, tumors and lesions Biopsy

K033437 595

6. Rationale for Substantial Equivalence

The Sphinx laser system with fiber optic delivery devices share the same intended use, indications for use, similar design features and functional features and therefore are substantially equivalent to the Lumenis VersaPulse PowerSuite 100 watt Holmium Laser, and the Trimedyne OmniMax 80 watt Holmium Laser.

7. Conclusion

The Sphinx Laser System with fiber optic delivery devices were found to be substantially equivalent to similar currently marketed and predicate surgical laser systems and delivery devices.





Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

JAN 2 0 2004

Mr. Peter Allen AllMed Systems, Inc. 9232 Klemetson Drive Pleasanton, California 94588

Re: K033437

Trade/Device Name: Sphinx

Regulation Number: 21 CFR 878.4810

Regulation Name: Laser surgical instrument for use in general and plastic surgery and

dermatology

Regulatory Class: II Product Code: GEX Dated: October 26, 2003 Received: October 30, 2003

Dear Mr. Allen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html

Sincerely yours,

Miriam C Provost For Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative and Neurological Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

510(k) Number: K033437
Device Name: Sphinx
Indications For Use:
The Sphinx laser system and its fiber optic delivery system are intended for use in surgical procedures using open, laparoscopic and endoscopic incision, excision, resection, ablation, vaporization, coagulation and hemostasis of soft tissue in use in medical specialties including: Urology, Urinary Lithotripsy, Gasteroenterology, Arthroscopy, Discetomy Pulmonary, Gynecology, ENT, Dermatology, Plastic Surgery and General Surgery.
Urology
Open and endoscopic surgery (incision, excision, resection, ablation, vaporization, coagulation and hemostasis) including: Urethral Strictures Bladder Neck Incisions (BNI) Ablation and resection of Bladder Tumors, Uretheral Tumors and Ureteral Tumors. Ablation of Benign Prostatic Hypertrophy (BHP), Transurethral incision of the prostate (TUIP) Holmium Laser Resection of the Prostrate (HoLRP) Holmium Laser Enuculeation of the Prostate (HoLEP) Holmium laser Ablation of the Prostate (HoLAP) Condylomas Lesions of external genitalia
Prescription Use AND/OR Over-The-Counter Use (Part 21 CFR 801 Subpart D) (21 CFR 807 Subpart C) (PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

510(k) Number: K033437					
Device Name: Sphinx					
Indications For Use:					
Lithotripsy and Percutaneous Urinary Lithotripsy Endoscopic fragmentation of urethral, ureteral, bladder and renal calculi including cystine, calcium oxalate, monohydrate and calcium oxalate dehydrate stones. Endoscopic fragmentation of kidney calculi Treatment of distal impacted fragments of steinstrasse when guide wire cannot be passed.					
Gasteroenterlogy					
Open and endoscopic gasteroenterlogy surgery (incision, excision, resection, ablation, vaporization, coagulation and hemostasis) including: Appendectomy Polyps Biopsy Gall Bladder calculi Biliary/Bile duct calculi Ulcers Gastric ulcers Duodenal ulcers Non Bleeding Ulcers Pancreatitas Hemorrhoids Cholecystectomy Benign and Malignant Neoplasm					
Prescription Use AND/OR Over-The-Counter Use (Part 21 CFR 801 Subpart D) (21 CFR 807 Subpart C) (PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)					
Concurrence of CDRH, Office of Device Evaluation (ODE)					

510(k) Number: K033437					
Device Name: Sphinx					
Indications For Use:					
Angiodysplasia Colorectal cancer Telangiectasias Telangiectasias of the Osler-Weber-Renu disease Vascular Malformation Gastritis Esophagitis Esophageal ulcers Varices Colitis Mallory-Weiss tear Gastric Erosions					
Arthroscopy					
Arthroscopy/Orthopedic surgery (excision, ablation and coagulation of soft and cartilaginous tissue) in small and large joints of the body, excluding the spine but including: Ligament and tendon Release Countouring and sculpting of articular surfaces Capsulectomy in the Knee Chondreplasty in the Knee Debridement of inflamed synovial tissue Chondromalacia Ablation Chondromalacia and tears					
Prescription Use AND/OR Over-The-Counter Use (Part 21 CFR 801 Subpart D) (21 CFR 807 Subpart C) (PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)					

Concurrence of CDRH, Office of Device Evaluation (ODE)

510(k) Number: K033437

Device Name:	Sphinx			
Indications Fo	· Use:			
N L	rlica Removal feniscectomy oose Body Debridement ateral retinecular release			
includinç F lı F	of soft, cartilaginous and bony tissue in Minimal Invasive Spinal Surgery Percutaneous Laser Disc Decompression/Discectomy of the L4-5 and L5-S1 Jumbar discs, including Foraminoplasty Percutaneous Cervical Disc Decompression/Discectomy Percutaneous Thoracic Disc Decompression/Discectomy			
Pulmonary				
Open and endoscopic pulmonary surgery (incision, excision, resection, ablation, vaporization, coagulation and hemostasis of soft tissue)				
Gynecology				
Open and laparoscopic gynecological surgery (incision, excision, resection, ablation, vaporization, coagulation and hemostasis) of soft tissue				
Prescription U (Part 21 CFR 80				
(PLEASE DO NEEDED)	NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF			
	Concurrence of CDRH, Office of Device Evaluation (ODE)			

510(k) Number: K033437
Device Name: Sphinx
Indications For Use:
ENT
Endoscopic endonasal surgery (incision, excision, resection, ablation, vaporization, coagulation and hemostasis of soft tissue and cartilage) including:
Endonasal/sinus Surgery Partial turbinectomy Polypectomy Dacryocystorhinostomy Frontal Sinusotomy Ethmoidectomy Maxillary antrostomy Functional endoscopic sinus surgery
Dermatology and Plastic Surgery
Incision, excision, resection, ablation, vaporization, coagulation and hemostasis of soft, mucosal, fatty and cartilaginous tissue, in therapeutic plastic, dermatologic and aesthetic surgical procedures including:
Basal Cell Carcinomas Lesions of skin and subcutaneous tissue Skin tags Plantar warts Lesions of skin and subcutaneous tissue
Prescription Use AND/OR Over-The-Counter Use (Part 21 CFR 801 Subpart D) (21 CFR 807 Subpart C) (PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

510(k) N umber: 1		idications to	
Device Name:	Sphinx		
Indications For U	se:		
	t Wine Stains illomas		
General Surgery			
vaporizatio App Skir Exc	n, coagulation a pendectomy n incision dision of extemal mplete of partial	nd hemostasis) incli and internal lesions	
Prescription Use (Part 21 CFR 801 S		AND/OR	Over-The-Counter Use (21 CFR 807 Subpart C)
(PLEASE DO N NEEDED)	NOT WRITE BE	ELOW THIS LINE-	CONTINUE ON ANOTHER PAGE
C	oncurrence of	CDRH Office of D	evice Evaluation (ODE)

Myram C Provost

Thereo Gign Off)

Mision of Conend, Restorative
and Neurological Devices

1910 - Ko33437

Page 6 of 6