PAXgene[®]
Blood ccfDNA Tube

EN – English

I. Intended Use

The PAXgene Blood ccfDNA Tube is intended for the collection, storage, and transport of human blood and stabilization of DNA in a closed tube for preparation of circulating, cell-free DNA (ccfDNA) from plasma and genomic DNA (gDNA) from nucleated cellular fraction for use with molecular test methods that require ccfDNA or gDNA.

The performance characteristics for the PAXgene Blood ccfDNA Tube have been established with 18S ribosomal and DYS14 Y-chromosomal ccfDNA fragments in plasma. Users must validate use of product for their specific molecular diagnostic assay.

Product Features

PAXgene Blood ccfDNA Tube • 768165

CE-IVD • Sterile • 100 tubes/case

10.0 ml draw volume

16 × 100 mm tube • 1.5 ml liquid additive

Safety-engineered BD Hemogard™ closure

Blue stopper • Pearlescent shield

- 2D barcode with serial number to uniquely identify the individual tube for sample tracking.
- Demonstrated ccfDNA stability and gDNA yield and purity: blood-filled PAXgene Blood ccfDNA Tubes can be stored for up to 10 days at temperatures up to 25°C, 7 days at temperatures up to 30°C, or 3 days at temperatures up to 37°C prior to processing.

Note: Do not store blood-filled tubes below 2°C.

- Compatible with magnetic bead- and silica membrane-based ccfDNA and gDNA isolation and purification methods that are automated or manual.
- Purified ccfDNA is compatible with PCR and next-generation sequencing (NGS) based molecular test methods.
- For in vitro diagnostic use.

II. Summary and Explanation

The PAXgene Blood ccfDNA Tube provides a means for the collection of whole blood for the isolation of ccfDNA from plasma and genomic DNA from the buffy coat or nucleated cellular fraction. Blood is collected under a standard phlebotomy protocol into an evacuated tube that contains a proprietary ccfDNA stabilization additive without formaldehyde or formaldehyde releasing substances. After collection, the whole blood sample can be stored for up to 10 days for temperatures up to 25°C, up to 7 days for temperatures up to 30°C, or for up to 3 days at temperatures up to 37°C prior to processing. After tube centrifugation, ccfDNA isolation from plasma is carried out using manual or automated methods that are based on magnetic bead or silica membrane technologies for ccfDNA isolation. The remaining buffy coat, a leukocyte and platelet rich layer at the interface of plasma and red blood cell fractions, or the remaining mixture of the buffy coat and red blood cells (cellular fraction) can be harvested and used to isolate gDNA using manual or automated methods that are based on magnetic bead or silica membrane technologies for DNA isolation.

This product was developed in accordance with the European Committee for Standardization Standard ISO 20186-3, *Molecular in vitro diagnostic examinations — Specifications for pre-examination processes for venous whole blood — Part 3: Isolated circulating cell free DNA from plasma*.

III. Specimen Collection and Processing

A. Required Blood Collection Accessories (not included with the PAXgene Blood ccfDNA Tube)

1. Blood collection set with tubing such as the BD Vacutainer® UltraTouch™ Push Button Blood Collection Set (see Section XI. Ordering Information).
2. BD Vacutainer Tube Holder (see Section XI. Ordering Information).
3. Labels for positive patient specimen identification.
4. Alcohol swab for cleansing site.
5. Dry, clean disposable gauze.
6. Tourniquet.
7. Biohazard container for used needle or needle/holder combination.

B. Recommended Order of Draw

When using single venipuncture to collect multiple tubes with different additives for clinical laboratory testing, follow the recommended tube order of draw per CLSI GP41, 7th ed. Collection of Diagnostic Venous Blood Specimens. PAXgene Blood ccfDNA Tubes should be collected last to avoid possible test result error due to additive carryover.


1. Tubes for sterile samples.
2. Tubes for coagulation studies (e.g., citrate).
3. Serum tubes with or without gel and clot activator.
4. Tubes with heparin additives.
5. Tubes with other additives (e.g., EDTA, fluoride).
6. PAXgene Blood ccfDNA Tube.

C. Prevention of Backflow

Since PAXgene Blood ccfDNA Tubes contain a liquid additive, it is important to avoid possible backflow from the tube, with the possibility of adverse subject reactions. To guard against backflow, observe the following precautions:

1. Place donor's arm in a downward position.
2. Hold tube with the stopper uppermost.
3. Release tourniquet as soon as blood starts to flow into tube.
4. Make sure tube additives do not touch stopper or non-patient end of the needle during venipuncture.

D. Procedure for Specimen Collection

1. Ensure that the PAXgene Blood ccfDNA Tube is at room temperature (15–25°C) prior to use.
2. Collect blood into the PAXgene Blood ccfDNA Tube using your institution's recommended procedure for standard venipuncture technique. In addition, follow instructions for Prevention of Backflow in Section III.C.
3. Fill the tube to its maximum stated draw volume. Continue blood collection until vacuum is exhausted and blood has stopped flowing into the tube before removing the tube from the holder.
4. After blood collection, gently invert the PAXgene Blood ccfDNA Tube 8 times to mix the blood with the additive. One complete inversion is defined as turning the filled tube upside-down, and then returning it to an upright position.
5. Store the filled PAXgene Blood ccfDNA Tube at room temperature (15–25°C). See Performance Characteristics for ccfDNA stability and gDNA yield and purity in blood samples at room temperature (15–25°C), 30°C or 37°C.

E. Plasma Preparation for Isolation of ccfDNA

1. Centrifuge the PAXgene Blood ccfDNA Tube at room temperature (15–25°C) for 15 minutes at 1900 × g using a balanced swing out bucket centrifuge. If braking is preferred, it is recommended to use medium level braking, but should be validated for your specific workflow.
2. Pipet the plasma into a 15 ml conical bottom centrifugation tube, making sure to not disturb the buffy coat and the cellular fraction.
Optional second centrifugation: For applications that require further purification of the plasma, centrifuge for 10 minutes at room temperature (15–25°C) and 1900 × g using a balanced centrifuge. Pipet the plasma into a 15 ml conical bottom centrifugation tube, making sure to not disturb the residual blood cell pellet, if present.
Note: Do not exceed the secondary tube manufacturer's maximum recommended centrifugation speed.
3. Process the plasma sample in accordance with the instructions provided with the ccfDNA sample preparation kit, or freeze the plasma (see Section III.F. Freezing and Thawing Plasma).
Note: For maximum ccfDNA yield, process the maximum volume of plasma available.
Note: When comparing ccfDNA yield to other blood collection tubes, take into account the dilution of the plasma by the PAXgene Blood ccfDNA Tube additive. If similar input volumes of blood are used, ccfDNA yields will be similar to plasma from a spray-dried K₂EDTA tube separated directly after blood draw. For more information, contact Technical Support.

F. Freezing and Thawing Plasma

1. For long-term storage, aliquot and freeze plasma at –20°C or –70°C / –80°C in cryogenic tubes.
2. Thaw aliquot tubes at room temperature (15–25°C).
Note: Do not thaw at lower temperatures (e.g., 4°C).
3. If cryoprecipitates form in the plasma, vortex the tube for 30 seconds after thawing and continue according to the instructions provided with the ccfDNA isolation kit.
Note: It is not recommended to centrifuge the plasma to remove cryoprecipitates because they may contain ccfDNA.
Note: To avoid the formation of cryoprecipitates, tubes can be thawed for 30 minutes at 30°C instead of at room temperature.

G. Isolation of gDNA

After transferring plasma for ccfDNA isolation, the remaining buffy coat or nucleated cellular fraction can be used to isolate gDNA, if desired. Buffy coat samples can be aspirated and processed in accordance with the buffy coat processing instructions provided with the DNA sample preparation kit. Kit examples are provided in Section XI. Ordering Information.

H. Freezing and Thawing the Nucleated Cellular Fraction in the PAXgene Blood ccfDNA Tube

1. Stand the PAXgene Blood ccfDNA Tube upright in a wire rack. Do not freeze tubes upright in an expanded polystyrene (EPS) tray as this may cause the tubes to crack.
2. The PAXgene Blood ccfDNA Tubes can be stored at -20°C and below. If tubes are to be kept at temperatures below -20°C , freeze them first at -20°C for 24 hours, then transfer them to -70°C or -80°C .
3. Thaw the PAXgene Blood ccfDNA Tubes in a wire rack at room temperature ($15\text{--}25^{\circ}\text{C}$).
4. Carefully invert the thawed PAXgene Blood ccfDNA Tubes 10 times.

Note: The frozen PAXgene Blood ccfDNA Tubes are subject to breakage upon impact. To reduce the risk of breakage during shipment, frozen tubes should be treated in the same manner as glass tubes. Users must validate their own freezing and shipping protocol for the PAXgene Blood ccfDNA Tubes.


IV. Performance Characteristics

A. ccfDNA Performance:

The performance characteristics of the PAXgene Blood ccfDNA Tube for ccfDNA stability were established using two validated qPCR assays (targeting 18S rDNA and DYS14 Y-chromosomal DNA sequences).

After blood is collected into the tube, the ccfDNA remains stable in whole blood for up to 10 days at temperatures up to 25°C . Based on ccfDNA stability studies, blood can be stored in the tube for up to 10 days at temperatures up to 25°C , 7 days at temperatures up to 30°C , or 3 days at temperatures up to 37°C (Figure 1).

Figure 1. ccfDNA Relative Yield for Blood Samples Stored in the PAXgene Blood ccfDNA Tube

Change in plasma ccfDNA yield after whole blood sample storage in comparison to plasma separated within 2 hours of blood collection (Day 0). Blood was drawn from a donor pool of approximately 200 consented, apparently healthy adult subjects and stored at various temperatures for the indicated number of days followed by tube centrifugation and ccfDNA purification from plasma using the QIASymphony PAXgene Blood ccfDNA Kit on the QIAGEN QIASymphony instrument. The relative ccfDNA yield was calculated as the ratio of the 18S rDNA C_T value after sample storage compared to the C_T value at Day 0. Medians and the 25th and 75th percentiles are denoted with box plots.

Table 1: Performance testing summary (automated, magnetic bead-based ccfDNA purification)


	Relative Yield Compared to EDTA Plasma at Day 0	Relative Yield (In Situ Stability) Compared to Day 0		
		Day 10, 25°C	Day 7, 30°C	Day 3, 37°C
n	120	120	120	120
Mean \pm SD	0.92 ± 0.26	1.30 ± 0.54	1.37 ± 0.59	1.52 ± 0.81
Median	0.93	1.21	1.24	1.27
Interquartile range	0.74–1.05	0.97–1.52	1.00–1.62	1.00–1.78
Range	0.26–1.87	0.32–3.15	0.31–3.27	0.48–4.26
95% of samples	0.33–1.50	0.51–2.86	0.50–3.15	0.59–4.82

B. gDNA Performance

The performance characteristics of the PAXgene Blood ccfDNA Tube for gDNA were established measuring integrity, concentration, and purity.

After blood is collected into the tube, the gDNA remains stable in whole blood for up to 10 days at temperatures up to 25°C . Based on gDNA stability studies, blood can be stored in the tube for up to 10 days at temperatures up to 25°C , 7 days at temperatures up to 30°C , or 3 days at temperatures up to 37°C (Figure 2).

Figure 2. gDNA Concentration for Blood Samples Stored in the PAXgene Blood ccfDNA Tube


gDNA Concentration, Yield, and Purity using an Automated, Magnetic Bead-Based DNA Purification System: Blood was drawn from a donor pool of approximately 200 consented, apparently healthy adult subjects into PAXgene Blood ccfDNA Tubes. Tubes were centrifuged within 2 hours of blood collection, and an aliquot was extracted from 400 µl of nucleated cellular fraction for processing. The remaining sample in each tube were stored at 25°C, 30°C, or 37°C for the indicated number of days. DNA was purified from 180 specimens using the QIASymphony DSP DNA Mini Kit (elution volume: 400 µl) on the QIAGEN QIASymphony instrument. Medians and the 25th and 75th percentiles are denoted with box plots.

Table 2: Performance testing summary (automated, magnetic bead-based DNA purification)

	Yield (µg DNA / 400 µl input sample volume)				Purity (A_{260}/A_{280})			
	Day 0	Day 10, 25°C	Day 7, 30°C	Day 3, 37°C	Day 0	Day 10, 25°C	Day 7, 30°C	Day 3, 37°C
n	120	120	120	120	120	120	120	120
Mean ± SD	16.61 ± 5.54	17.80 ± 4.95	17.03 ± 4.95	16.80 ± 4.76	1.83 ± 0.07	1.84 ± 0.06	1.84 ± 0.07	1.85 ± 0.06
Median	16.30	17.62	17.04	16.82	1.83	1.84	1.82	1.84
Interquartile range	12.73–19.69	14.08–20.12	12.93–20.34	13.06–19.95	1.79–1.86	1.80–1.88	1.79–1.87	1.81–1.88
Range	5.74–36.68	8.08–32.76	8.02–28.82	8.18–34.34	1.58–2.02	1.74–2.05	1.68–2.09	1.71–2.05
95% of samples	≥8.27	≥9.89	≥9.19	≥9.34	1.68–1.98	1.71–1.98	1.69–1.98	1.71–1.98

Figure 3. gDNA Integrity after Storage

Blood was drawn from 20 consented, apparently healthy adult subjects into PAXgene Blood ccfDNA Tubes. Total gDNA was purified from 400 µl of nucleated cellular fraction using the QIASymphony DSP DNA Mini Kit (elution volume: 400 µl) on the QIAGEN QIASymphony instrument. For each subject, 400 ng gDNA was separated by agarose gel electrophoresis (lanes 2–21). In lane 1 and 22, a Lambda x Hind III marker was loaded. The upper band of this marker represents a DNA fragment of 23 kb. The image shows gDNA isolation after storage at 30°C for 7 days. The data for all other time and temperature processing conditions are similar.

V. Limitations

1. The quantity of blood drawn should be approximately 10.0 ml per PAXgene Blood ccfDNA Tube, but this volume may vary depending on various factors such as altitude, ambient temperature, barometric pressure, tube age, venous pressure, and filling technique.
2. The PAXgene Blood ccfDNA Tube is not designed for use with open blood collection systems (manual filling of tube with the BD Hemogard closure removed) due to the increased risk of blood exposure and risk of causing incorrect blood-to-additive ratio that may affect product performance. Blood should be collected directly into the tube, or a transfer device should be used if blood is collected and transferred from a syringe.
3. ccfDNA and gDNA yields depend on the patient, the quality of the specimen, and the method used for isolation.

VI. Warnings and Precautions

A. Precautions:

1. Examine tubes prior to use. Do not use tubes if foreign matter is present inside the tube.
2. Since this PAXgene Blood ccfDNA Tube contains a liquid additive, precautions should be taken to prevent possible backflow from the tube during blood draw (see Section III.C. Prevention of Backflow).
3. Do not shake vigorously since this may cause hemolysis.
4. If tubes are not mixed 8 times immediately after collection, clotting of blood or fibrin formation may occur. This may also reduce effectiveness of blood stabilization that requires mixing with the additive immediately upon collection for best performance.
5. Remove BD Hemogard closures with a twist and pull motion. Removal by rolling with the thumb is not recommended (for detailed instructions see Section VIII. Removal of BD Hemogard Closure).
6. After venipuncture, the top of the stopper may contain residual blood. Take proper precautions when handling tubes to avoid contact with this blood.
7. Under-filling of PAXgene Blood ccfDNA Tubes will result in an incorrect blood-to-additive ratio and may lead to incorrect analytical results or poor product performance.
8. Endotoxin is not controlled. Blood and blood components collected and processed in the tube are not intended for infusion or introduction into the human body.

B. Warnings:


1. The liquid additive in this tube is a mild irritant and may cause irritation to the eyes and skin upon direct contact:
 - After skin contact, wash skin with soap and water. Get medical attention if irritation persists after washing.
 - After eye contact, flush eyes with water as a precaution. If irritation occurs, get medical assistance.
 - After swallowing, rinse mouth. Get medical attention if any discomfort occurs.
2. Practice universal precautions. Use gloves, gowns, eye protection, other personal protective equipment, and engineering controls to protect from blood splatter, blood leakage, and potential exposure to blood borne pathogens.
3. Handle all blood samples and blood collection devices (needles, luer adapters, and blood collection sets) according to the policies and procedures of your facility. Obtain appropriate medical attention in the event of any blood exposure (for example, through a needlestick injury), since there is a risk of transmission of viral hepatitis, HIV, or other infectious diseases. Utilize any built-in needle protector after blood collection if the blood collection device provides one. PreAnalytiX does not recommend re-shielding used needles; however, the policies and procedures of your facility may differ and must always be followed.
4. A blood collection set must be used with the PAXgene Blood ccfDNA Tube (see Section XI. Ordering Information).
5. Excessive centrifugation speed (over 10,000 x g) may cause PAXgene Blood ccfDNA Tube breakage, exposure to blood and possible injury.
6. Transferring a blood sample from a syringe with a needle into the PAXgene Blood ccfDNA Tube is not recommended due to the increased risk for both needlestick injury and incorrect blood-to-additive ratio.
7. If blood is collected from an intravenous (I.V.) line, ensure that line has been cleared of I.V. solution before beginning to fill blood collection tubes. This is critical to avoid erroneous laboratory test results due to I.V. fluid contamination.
8. After use, discard all blood collection tubes and accessories in biohazard containers approved for their disposal.
9. Do not re-use the PAXgene Blood ccfDNA Tubes.
10. Do not use the PAXgene Blood ccfDNA Tubes after the expiration date printed on the tube label.

VII. Storage

Store the unused PAXgene Blood ccfDNA Tubes at 4–25°C. Brief temperature excursions from –10°C to 50°C are permitted.


The PAXgene Blood ccfDNA Tube additive may have a slightly yellow appearance; this does not affect the performance of the additive. Do not use tubes after their expiration date

VIII. Removal of BD Hemogard Closure

1. Grasp the PAXgene Blood ccfDNA Tube with one hand, placing the thumb under the BD Hemogard closure. (For added stability, place arm on solid surface). With the other hand, twist the BD Hemogard closure while simultaneously pushing up with the thumb of the other hand **ONLY UNTIL THE TUBE STOPPER IS LOOSENED**.
2. Move thumb away before lifting closure. **DO NOT** use thumb to push closure off tube. If the tube contains blood, an exposure hazard exists.
3. Lift closure off tube. In the unlikely event of the plastic shield separating from the rubber stopper, **DO NOT REASSEMBLE CLOSURE**. Carefully remove rubber stopper from tube.

IX. Reinsertion of BD Hemogard Closure

1. Replace closure over tube.
2. Twist and push down firmly until stopper is fully resealed. Complete reinsertion of the stopper is necessary for the closure to remain securely on the tube during handling.

X. Technical Assistance

If you have any questions regarding the PAXgene Blood ccfDNA Tube, contact one of the BD Technical Services Departments listed in the BD - Customer Service section.

XI. Ordering Information

PAXgene Products

PAXgene Blood ccfDNA Tube / 768165

Content: 100 blood collection tubes.

10.0 ml • 16 × 100 mm	1.5 ml additive	Sterile
BD Hemogard closure	100 tubes/case	CE-IVD

The PAXgene Blood ccfDNA Tube (CE-IVD) is available in Europe and some parts of the world outside the United States. For more information please visit www.preanalytix.com.

To order PAXgene Blood ccfDNA Tubes or Kits:

Go to: www.preanalytix.com

BD Products*

BD Vacutainer UltraTouch Push Button Blood Collection Set / 367364

Content: 23G 3/4 inch (0.8 × 19 mm) needle, 12 inch (305 mm) tubing with luer adapter. 50/box, 200/case

BD Vacutainer Safety-Lok™ Blood Collection Set / 367281 US / 367286 CE

Content: 21G 3/4 inch (0.8 × 19 mm) needle, 12 inch (305 mm) tubing with luer adapter. 50/box, 200/case

BD Vacutainer One Use Holder / 364815

Content: Case only for 13 mm and 16 mm diameter. 1000/case

* These catalog numbers represent typical products that can be used with the PAXgene Blood ccfDNA Tube.

For more information on BD Blood Collection Accessories:

Go to: www.bd.com

QIAGEN Products*

To order QIAGEN ccfDNA Preparation Kits:

Go to: www.qiagen.com/shop

For Genomic DNA Isolation:

QIASymphony DSP DNA Mini Kit (192) / 937236

2 reagent cartridges and enzyme racks and accessories.

* These catalog numbers represent typical sample preparation kits that can be used with the PAXgene Blood ccfDNA Tube.

To order QIAGEN gDNA Preparation Kits:

Go to: www.qiagen.com/shop

BD – Customer Service

Argentina, Uruguay and Paraguay

Orders: 0800.444.5523
E-mail: crc_argentina@bd.com

Australia

Orders: 1.800.656.100
Fax: 1.800.656.110
E-mail: bd_anz@bd.com

Austria

Orders: 43.1.7063660
Fax: 43.1.706366011
E-mail: customercare.at@bd.com

Belgium

Orders: 32.53.720.556
Fax: 32.53.720.549
E-mail: orders.be@bd.com

Brazil

Orders: 0800.055.56.54
E-mail: consultoria_vacutainer@bd.com

Canada

Technical support: 1.800.631.0174
Orders: 1.866.979.9408
Fax: 1.800.565.0897
E-mail: customer.service.canada@bd.com

Central and Eastern Europe

Orders: 48.22.377.11.11
Fax: 48.22.377.11.02
Bulgaria orders: info_bulgaria@bd.com
Czech Republic orders: info_czech@bd.com
Croatia orders: info_croatia@bd.com
Hungary orders: info_hungary@bd.com
Poland orders: info_poland@bd.com
Romania orders: info_romania@bd.com
Southeast Europe orders: info_balkan@bd.com
Serbia orders: info_serbia@bd.com
Slovakia orders: info_slovakia@bd.com
Slovenia orders: info_slovenia@bd.com

Denmark

Orders: 45.43.43.45.66
Fax: 45.43.96.56.76
Orders: ordre.dk@bd.com
Technical support: bddenmark@bd.com

Finland

Orders: 358.9.88.70.780
Fax: 358.9.88.70.7816
Orders: tilaukset.fi@bd.com
E-mail: bdsuomi@bd.com

France

Orders: 33.476.68.36.36
Fax: 33.476.68.36.93
E-mail: serviceclientbdf@bd.com
Orders: commandesfr@bd.com
Technical support: vacutainerfr@bd.com

Germany

Orders: 49.6221.3050
Fax: 49.6221.305.216
E-mail: customercare.de@bd.com

India

Orders: 91.124.3949390
Orders: bd_india@bd.com

Ireland (Aquilant Specialist Healthcare Services)

Customer support: 353.1.404.8350
Fax: 353.1.404.8352
E-mail: contactus@aquilantscientific.ie

Israel (Lapidot Medical)

Customer Support: 972.700.70.90.22
E-mail: cs@lapidot.com

Italy

Orders: 39.02.48240.500
Fax: 39.02.48240.775
Technical support: 39.3450655140
E-mail: ordini.it@bd.com

Middle East & Africa

Orders: 971.45.592.555
Fax: 971.45.592.599
E-mail: EMA_PAS@bd.com

The Netherlands

Orders: 31.20.582.94.20
Fax: 31.20.582.94.21
Orders: orders.nl@bd.com

New Zealand

Orders: 0800.572.468
Fax: 0800.572.469
E-mail: nz_customerservice@bd.com

Norway

Customer Support: 64.00.99.00
E-mail: bdnorge@bd.com
Orders: ordre.no@bd.com

Southeast Asia

E-mail: PAS.SEA@bd.com
Indonesia orders: 622.1577.1920
Malaysia orders: 603.2093.8788
Philippines orders: 63.2478.8881
Singapore orders: 65.6861.0633
Thailand orders: 662.646.1800
Vietnam orders: 848.3822.7409

South Korea

Orders: 02.3404.3706

Fax: 02.3404.3785

Technical: 02.3404.3706

Technical support: Korea_PAS@bd.com

Spain, Portugal and Andorra

Orders: 34.91.848.8174

Customer support: 34.902.27.17.27

Fax: 34.91.848.8115

E-mail: info.spain@bd.com

Sweden

Orders: 46.8.775.51.00

Fax: 46.8.645.08.08

Orders: order.se@bd.com

Technical support: bds sweden@bd.com

Switzerland

Orders: 41.61.485.22.24

Fax: 41.61.485.22.00

E-mail: infoch@bd.com

UK

Orders: 0800.917.8776

E-mail: bduk_customerservice@bd.com

USA


Customer support: 800.631.0174

E-mail: productcomplaints@bd.com

QIAGEN – Customer Service

Ordering www.qiagen.com/shop | Technical Support support.qiagen.com | Website www.qiagen.com

Symbol and Mark Key

	Batch Code
	In Vitro Diagnostic Medical Device
	Sterilized Using Irradiation
	Catalog Number
	Keep Away from Sunlight
	Temperature Limit
	Do Not Reuse
	Manufacturer
	Use-By Date
	Recyclable
	This End Up
	Fragile, Handle With Care
	Contains sufficient for <n> tests
	Telephone Instructions for Use
	Online Instructions for Use
	CE Marking

The PAXgene Blood ccfDNA Tube is made in the UK by BD for PreAnalytiX GmbH.

PAXgene and PreAnalytiX are trademarks of PreAnalytiX GmbH; QIAGEN and QIASymphony are trademarks of QIAGEN NV; all other trademarks are property of Becton, Dickinson and Company.

© 2019 PreAnalytiX GmbH. Unless otherwise noted, PreAnalytiX, the PreAnalytiX Logo and all other trademarks are property of PreAnalytiX GmbH, Hombrechtikon, CH.

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
05/2020
VDP40336-WEB-04

Production Information
This document is intended for web viewing only

Code	Language	Page
EN	English	01 – 10
BG	Български	11 – 21
CZ	čeština	22 – 31
DA	Dansk	32 – 41
DE	Deutsch	42 – 51
EL	Ελληνικά	52 – 62
ES	Español	63 – 73
ET	Eesti	74 – 83
FI	Suomi	84 – 93
FR	Français	94 – 103
HR	Hrvatski	104 – 113
HU	Magyar	114 – 123
IS	Íslenska	125 – 133
IT	Italiano	134 – 143
KK	Қазақша	144 – 154
LT	Lietuvių k.	155 – 165
LV	Latviski	166 – 176
NL	Nederlands	177 – 187
NO	Norsk	188 – 197
PL	Polski	198 – 207
PT	Português	208 – 218
RO	Română	219 – 228
RU	Русский	229 – 239
SK	Slovenský (jazyk)	240 – 249
SL	Slovenski (jezik)	250 – 259
SV	Svenska	260 – 269
TR	Türkçe	270 – 279
UK	Українська	280 – 290

BD - Preanalytical Systems • Franklin Lakes, NJ

REV LEVEL	ALOA #	ECR/O#	DESCRIPTION	DESIGNER	DATE	SAP
01	00996-02	5-103790	Release new web product circular for PAXgene Blood ccfDNA Tube.	Premanand M	02/FEB/2018	N/A
02	01027-01	5-120313	Updated Table 1 information & Figure 1 for all translations. Revision date has been updated.	Premanand M	08/Mar/2018	N/A
03	01028-01	5-122091	Updated the correct translations for Bulgarian (BG) language and few other languages. Updated the formatting changes. Updated USA technical service email address.	Premanand M	24/Apr/2018	N/A
04	01159-01	5-173933	This revision is to update the shelf life of the PAXgene Blood ccfDNA Tube (Catalog Number 768165) from 13 months to 15 months in order to meet market demand for a longer shelf life. The following changes were made in the EN version of this document and its respective translations: In page 2, updated "Intended Use", "Product Feature" and "Summary and Explanation" sections. In page 3, updated "Procedure for Specimen Collection", "Plasma Preparation for Isolation of ccfDNA" and "Freezing and Thawing Plasma" sections. In page 4, added "Freezing and Thawing the Nucleated Cellular Fraction in the PAXgene Blood ccfDNA Tube" topic under section IV; updated "ccfDNA Performance", "gDNA Performance" section, Figure1 and Table1. In Page 5, updated Figure 2, Table 2 and image/text under Figure 3. In Page 6, updated "Warnings" and "Storage" section. In page 7, updated " Ordering Information" section. In page 8 and 9, updated "BD – Customer Service" information. In page 10, updated Copyright year, document design date and revision number.	Chaitanya R	20/May/2020	N/A
PACKAGING LEVEL			Product Circular for Web only	DRAWING #	VDP40336-WEB-04	