Codman®

CODMAN® Cranial Access Kits

Integra LifeSciences Production Corporation
11 Cabot Boulevard
Mansfield, MA 02048 USA

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The CE mark on this insert pertains
only to Catalog No. 82-6617


LCN 197142-001 Rev R 2020-09 1301379-2
©2020 Integra LifeSciences Corporation.
All Rights Reserved.

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

STERILE EO

Indications

The CODMAN® Cranial Access Kit is indicated when a craniotomy is required for placement of an intracranial pressure (ICP) monitoring device and/or cerebrospinal fluid drainage procedures.

Contraindications

This device is not designed, sold, or intended for any use except as indicated.

Description

The CODMAN Cranial Access Kit is a disposable intracranial procedural kit which contains all items used during each step of the cranial access procedure. Catalog Nos. 82-6614 and 82-6616 contain medical devices with sharp injury prevention features to help minimize occupational exposure to bloodborne pathogens.

The kits contain the following components (see Figures 1–3):

Preparation and Closure Kit:

- A. Double edge razor (1) (non-sterile in Catalog No. 82-6614; sterile in Catalog Nos. 82-6616 and 82-6617)
- B. DuraPrep® Surgical Solution, 6 mL (1) (sterile)*
- C. Ruler (1) (non-sterile)*
- D. Marking pen (1) (non-sterile)*
- E. Sodium chloride, 0.9%, 10 mL flip top vial (2) (sterile)*
- F. Bone wax (1) (sterile)
- G. 2.0 silk suture (1) (sterile)
- H. 3.0 nylon suture (1) (sterile)
- I. XYLOCAINE® 1% Lidocaine with epinephrine, 30 mL flip top vial

* Catalog No. 82-6614 only

Cranial Access Kit (Sterile Content):

Prep for Incision:

- J. Ruler, 6" (1)
- K. Marking pen (1)
- L. Syringe (2)
- M. Safety Needle, 18 gauge x 1 ½" (2)*

N. Safety Needle, 25 gauge x 5/8" (1)*

* Catalog No. 82-6617 includes Needle

18 gauge x 1 1/2" (2), and Needle 25 gauge x 5/8" (1)

Burr Hole Instruments:

- O. Spinal needle, 18 gauge x 3 ½" (1)
- P. Ventricular needle, 12 gauge x 3 ½" (1)
- Q. Safety scalpel with retractable cover, no. 11 (1)*
- R. Safety scalpel with retractable cover, no. 15 (1)*
- S. Retractor, blunt (1)
- T. Hand crank drill (1)
- U. 2.7 mm drill bit w/stop (1)
- V. 5.8 mm drill bit w/stop (1)
- W. Allen wrench (1)

* Catalog No. 82-6617 includes scalpels no. 11 and 15 without retractable covers

Closure Instruments:

- X. Mosquito forceps, curved (2)
- Y. Needleholder, serrated (1)
- Z. Adson forceps w/teeth (1)
- AA. Adson forceps (1)
- BB. Suture scissors (1)

Also included in kits but not shown in Figures 1–3 are the following sterile components:

- Fenestrated drape w/barrier (1)
- Gauze sponge 4" x 4", 12-ply (10)
- Absorbent towel (3)
- Culture tube w/screw cap
- Medicine cups (2)

Adverse Events

Bleeding may occur at the site of the drill hole, originating from the scalp, bone, dural, or cerebral areas. Due to the possibility of cerebral or extra-cerebral hemorrhage or any other complication of this procedure, drill hole placement is the sole responsibility of the attending neurosurgeon.

WARNINGS

The drill guide will not stop the drill. The guide is designed only to provide the neurosurgeon with a marker for drilling depth. The guide must be adjusted to the proper position prior to drilling.

The drilling procedure must be performed with the drill held within 10° of the perpendicular position to the incision site. When used with the CODMAN Skull Bolt, installing at an angle may result in a fracture of the device.

Exercise caution when perforating the dura to avoid damage to the underlying structures.

Precautions

Inspect the sterile package carefully. Do not use if:

- the package or seal appears damaged,
- contents appear damaged, or
- the expiry date has passed.

Do not reuse. This product is for **single use only**.

It is essential to maintain strict aseptic technique during craniotomy procedures.

Care should be exercised when applying sutures to ventricular catheters or the Codman Microsensor®. Overtightening could result in catheter occlusion or breakage.

Several items have sharp tips. Exercise caution when removing items from trays and when removing protective tips.

The following items have individual expiration dates marked on their packaging: XYLOCAINE®, DuraPrep® Surgical Solution, sodium chloride, bone wax, silk suture, and nylon suture. Verify product has not expired prior to use.

All instruments are designed for single use only. Do not resterilize or reuse.

The XYLOCAINE should not be used if its color is pinkish or darker than slightly yellow or if it contains a precipitate.

The sodium chloride should not be used if particulate matter and/or discoloration appear.

Instructions for Use

1. After prepping and draping the patient, perform the incision and retraction techniques required to expose the skull.
2. Select the appropriate drill bit. Use the 5.8 mm drill bit for ventriculostomy procedures (82-6633) and plastic bolt procedures (82-6632). Use the 2.7 mm drill bit for subdural, intraparenchymal (82-6631), and bolt procedures (82-6638). Refer to the instructions included with the intended drainage or monitoring device to verify the correct bit size.
3. Place the drill bit into the chuck.
4. While holding the drill handle in place, tighten the drill bit by turning the chuck counterclockwise.
Note: To remove the drill bit, repeat this step turning the chuck clockwise.
5. Loosen the drill guide using the appropriate hex wrench.
6. Carefully slide the drill guide towards the tip of the drill bit until the determined skull depth is reached.
7. Tighten the drill guide in place with the hex wrench.
8. Begin drilling procedure.

9. Using appropriate technique, close and dress the surgical site as indicated by the procedure. Silk suture is included to secure the Codman Microsensor or ventricular catheter to the scalp.

Monitoring/Drainage

Placement of a monitoring device and/or drainage system should be completed as indicated by the procedure.

How Supplied

The CODMAN Cranial Access Kit is intended for SINGLE USE ONLY; DO NOT RESTERILIZE. Integra single-use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or resterilization, after a single patient use. Reuse can potentially compromise device performance and any usage beyond the design intent of this single-use device may result in unpredictable loss of functionality.

Integra will not be responsible for any product that is resterilized, nor accept for credit or exchange any product that has been opened, but not used.

The following components have been tested and were determined to be nonpyrogenic:

Ventricular needle

Accessories

The following accessories are available directly from Integra:

Catalog Number	Description
82-1705	CODMAN® External Drainage Ventricular Catheter Set
82-1720	CODMAN® External Drainage System II
82-1724	CODMAN® External Drainage System II with THE GHAJAR GUIDE® Ventriculostomy Kit
82-1730	CODMAN® EDS 3™ CSF External Drainage System
82-1735	CODMAN® EDS 3™ Clear CSF Catheter
82-1739	CODMAN® EVD Catheter Set, Large Lumen
82-1745	CODMAN® BACTISEAL® EVD Catheter Set
82-1749	CODMAN® EVD BACTISEAL Catheter Set, Large Lumen
82-6724	CODMAN® Double Lumen Bolt (5.8 mm)
H00001600	QFLOW 500 Probe
NS9011	CODMAN® NS9011 Skull Bolt (2.7 mm)

PRODUCT INFORMATION DISCLOSURE

INTEGRA LIFESCIENCES CORPORATION (“INTEGRA”) HAS EXERCISED REASONABLE CARE IN THE SELECTION OF MATERIALS AND THE MANUFACTURE OF THESE PRODUCTS. INTEGRA WARRANTS THAT THESE PRODUCTS SHALL CONFORM TO THE PRODUCT LIMITED WARRANTY AS PROVIDED IN THE PRODUCT LABELING OR APPLICABLE PRODUCT CATALOG. THIS WARRANTY IS EXCLUSIVE, AND INTEGRA DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. INTEGRA SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM USE OF THESE PRODUCTS. INTEGRA NEITHER ASSUMES NOR AUTHORIZES ANY PERSON TO ASSUME ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THESE PRODUCTS.

Codman, Codman Microsensor, and Bactiseal are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries. EDS 3 is a trademark of Integra LifeSciences Corporation or its subsidiaries. THE GHJAR GUIDE is a trademark of Neurodynamics, Inc. XYLOCAINE is a trademark of Fresenius Kabi USA, LLC. DuraPrep is a trademark of 3M Company.

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