


Online Neurocognitive
Concussion Testing:
Interpreting Reports and
Case Study Examples



Concussion
Vital Signs®

CONCUSSION WEBINAR SERIES

Dr. Peter C. Entwistle
Dr. Anne-Marie Kimbell

PEARSON

Concussion Vital Signs Product Training

AGENDA:

- Overview Concussion Vital Signs
 - Baseline & follow up testing
 - Medical history
 - Sideline app (SCAT2)
 - Symptom scale
 - Memory Function
 - Balance testing
 - Neurocognitive testing
 - Clinician's portal.
- Interpretation



2 | Copyright © 2013. All rights reserved. PEARSON

Concussion Vital Signs- overview

Concussion Vital Signs Product Training



Concussion Management Example

Pre-Participation or Pre-Season Screen to Athlete

Involuntary Concussion Assessment
e.g. Sideline

Evaluation & Management of Concussed Athlete

Multidisciplinary Multimedical Multitask

How Concussion Vital Signs can help...

- ✓ Athlete Education
- ✓ Neurocognitive Testing
- ✓ Athlete Information & Medical History
- ✓ Concussion Symptom Scale
- Balance Testing e.g. BESS
- ✓ Athletic Contest Evaluation with Pocket SCAT 2 – iPad, iPhone, Android, XBOX, Blackberry.
- ✓ Informant Reports e.g. Coach, Parent, ATC, Nurse, etc.
- Thorough History & Physical: e.g. Head and Neck Exam, Neurological Exam, Gait and Balance Assessment
- Previous History of Head Injuries (from baseline)
- ✓ Concussion Symptom Scale (baseline and Post-Injury)
- ✓ Neurocognitive Test: baseline and Post-Injury
- Balance Testing e.g. BESS, NeuroCom, Babix, etc.
- Possibly Neuroimaging

3 | Copyright © 2013. All rights reserved. PEARSON

Concussion Vital Signs Product Training

Symptoms and History

Pocket SCAT2

Sideline Baseline / Post-Injury

Brain Function

Computerized Neuropsychological Tests
e.g. Concussion Vital Signs
Paper and Pencil

Balance Testing

Brain imaging

Lab Values

4 | Copyright © 2013. All rights reserved. PEARSON

Concussion Vital Signs Product Training

All components are electronic: tests, sideline assessment, history and symptoms scales
Helps facilitate documentation, saves time

Pre-Season Baseline Testing

- Neurocognitive Evaluation
- Concussion History
- Concussion Symptom Scale
- Balance Testing

Key Advantage:
Parents can Login and do the Concussion History.

Sideline Assessment

- Pocket SCAT 2
- Emergency Procedures e.g., Call 911

Key Advantage:
Complete Sideline and Access Info Remotely

iPad, iPhone
Xbox, Droid
Blackberry Enabled

Post-Injury Assessment

- Rest Brain Until Symptom Free Using Concussion Symptom Scale
- After Symptom Free with Graduated Exercise then Re-administer the Neurocognitive Test
- Update Concussion History
- Balance Testing

5 | Copyright © 2013. All rights reserved. PEARSON

Concussion Vital Signs Product Training

Why use Baseline Neurocognitive Tests?

Invisible effects of concussion.

Broglio, et al, 2007
N=21 D1 college athletes, (16 men, 5 women).

- Neurocognitive decrements may persist when athletes no longer report concussion-related symptoms.
- Reliance on athlete-reported, post concussion symptoms when making return-to-play decisions may expose athletes to subsequent injury if complete recovery has not occurred.
- A multifaceted approach to concussion assessment that includes evaluation of a myriad of functions is warranted.
- Risk for **Second Impact Syndrome**.
- A VALID AND RELIABLE TEST IS REQUIRED.
- The exclusive use of symptom reports in making a return-to-play decision is not advised.*

6 | Copyright © 2013. All rights reserved. PEARSON

OPTIMIZED for Easier Management
Concussion Vital Signs Product Training



First step for athlete is to complete
Baseline testing
Concussion Vital Signs Product Training

Conditions that create highly valid test results:

- Computers with internet access.
- Time: 30-40 minutes.
- Rested, Unhurried athletes.
- Use a Setting with limited distractions.
- Use Similar settings across administrations.



Set the frame clearly for athletes:

"Your Brain is Your Life. Take this seriously."

- Do not distract each other.
- If you don't listen to the instructions your report will be invalid.
- **You will take it again until it is valid.**
- Until your test is valid you don't practice for your sport....

Note to those administering the test: Take it yourself a few times so you will know how to orient athletes to the test. Also, read:
Resources: Test Administration Guide.

Concussion Vital Signs Product Training

- Peer Reviewed Valid & Reliable Tests
- Tests are electronic versions of widely used venerable Paper & Pencil tests used in mTBI assessment
- Tests are being used worldwide by over 6000 clinical users, in over 2000 clinical investigator research sites, in 52 countries, and in major TBI research projects
- Enables an improved longitudinal collection of important clinical endpoints... baseline, sideline, post-injury
- Neurocognitive tests are available in 50+ languages
- OPTIMIZED to meet Consensus Concussion Guidelines
- OPTIMIZED to help identify athlete sandbagging
- OPTIMIZED for easier management
- OPTIMIZED for Life Span Testing norms from ages 8-90

Concussion Sideline Assessment: Pocket SCAT

Athlete Reference/ID: Test Date/Time:

Full Name: Test Administrator Name/Position:

Concussion should be suspected in the presence of any one or more of the following: symptoms (such as the physical signs (such as unsteadiness), or impaired brain function (e.g. confusion) or abnormal behavior.

1. Symptoms: The presence of any of the following signs and symptoms may suggest a concussion.
Check the SYMPTOMS exhibited by the athlete.

		Yes	No		
1	Loss of consciousness	<input type="checkbox"/>	<input type="checkbox"/>	13	Feeling slowed down
2	Seizure or convulsion	<input type="checkbox"/>	<input type="checkbox"/>	14	Feeling like "in a fog"
3	Amnesia	<input type="checkbox"/>	<input type="checkbox"/>	15	"Don't feel right"
4	Headache	<input type="checkbox"/>	<input type="checkbox"/>	16	Difficulty concentrating
5	"Pressure in head"	<input type="checkbox"/>	<input type="checkbox"/>	17	Difficulty remembering
6	Neck Pain	<input type="checkbox"/>	<input type="checkbox"/>	18	Fatigue or low energy
7	Nausea or vomiting	<input type="checkbox"/>	<input type="checkbox"/>	19	Confusion
8	Dizziness	<input type="checkbox"/>	<input type="checkbox"/>	20	Drowsiness
9	Blurred vision	<input type="checkbox"/>	<input type="checkbox"/>	21	More emotional

2. Memory function: Failure to answer all questions correctly may suggest a concussion.

	Incorrect	Correct	Additional comments:
At what venue are we at today?	<input type="checkbox"/>	<input type="checkbox"/>	
Which half is it now?	<input type="checkbox"/>	<input type="checkbox"/>	
Who scored last in this game? Incorrect	<input type="checkbox"/>	<input type="checkbox"/>	
What team did you play last week/game?	<input type="checkbox"/>	<input type="checkbox"/>	
Did your team win the last game?	<input type="checkbox"/>	<input type="checkbox"/>	

3. Balance testing: Instructions for tandem stance "Now stand heel-to-toe with your non-dominant foot in front of your dominant foot. Your weight should be evenly distributed across both feet. You should try to maintain stability for 20 seconds with your hips and your eyes closed. I will be counting the number of times you move out of this position. If you move out of this position, open your eyes and return to the start position and continue balancing. I will start timing when you have closed your eyes."

Athlete was observed for 20 seconds. If they make more than 5 errors (such as lift their hands off their hips; open their eyes; lift their forefoot or heel; step, stumble, or fall; or remain out of the start position for more than 5 seconds) then this may suggest a concussion.

Number of Errors:

Any athlete with a suspected concussion should be IMMEDIATELY REMOVED from play and should be urgently assessed medically, should not be left alone and should not drive a motor vehicle.

Sideline or abbreviated testing is designed to assist with the immediate assessment or screening of sports related concussion (e.g. SIDELINE) and should not replace computerized or comprehensive neuropsychological testing. The Concussion Sideline Assessment should not be used as a stand alone management tool. The Concussion Sideline Assessment is designed to be a support for recognizing sports related concussions and to document a concussion. It may assist a qualified health professional in their return-to-play decision making.

What are concussion danger signs?

In rare cases, a dangerous blood clot may form on the brain in a person with a concussion and crowd the brain against the skull. The student should be taken to an emergency department right away if s/he exhibits any of the following danger signs after a bump, blow, or jolt to the head or body:

- One pupil larger than the other
- Is drowsy or cannot be awakened
- A headache that gets worse and does not go away
- Weakness, numbness, or decreased coordination
- Repeated vomiting or nausea
- Slurred speech
- Convulsions or seizures
- Cannot recognize people or places
- Becomes increasingly confused, restless, or agitated
- Has unusual behavior
- Loses consciousness (even a brief loss of consciousness should be taken seriously)

For more information and tool kits for youth sports coaches and high school coaches, visit www.cdc.gov/Concussion.

PEARSON

Concussion Vital Signs Product Training
 Clinician's Portal maintains confidentiality and allows freedom to choose your own clinician

- Required information is provided to Physician by Parent
- Physician administers post-injury
- Pre-pay for an account to speed testing
- Purchase bundle for "Concussion Clinic"
- Do not have to have Impact approved doctor
- Improves communications between parents, trainers, physicians
- Allows physician to see full concussion history if parents provide access
- Billing codes: 96118, 96120

25 | Copyright © 2013. All rights reserved. **PEARSON**

Concussion Vital Signs Product Training
 Who should interpret the CVS neurocognitive test results?

- State legislation and medical guidelines generally require a medical or physician release for athletes to return-to-play.
- Interpretation of the Concussion Vital Signs neurocognitive test results should be done by a qualified health professional.

26 | Copyright © 2013. All rights reserved. **PEARSON**

Concussion Vital Signs Product Training
 How do you interpret Concussion Vital signs data?

1. Is it valid?
2. Do you have baseline data for comparison?
3. Is there any evidence of lack of effort?
4. What is the neuro-cognitive index?
5. Is there evidence of a possible concussion? Signs & symptoms?
6. Is there a decline in functioning?
7. Is this within 5% of the baseline?
8. Has repeat testing been done once the symptoms have resolved?

INTERPRETATION

27 | Copyright © 2013. All rights reserved. **PEARSON**

Concussion Vital Signs

Concussion Vital Signs Product Training



- Contains seven venerable computerized neuropsychological tests and the clinical domains, scored from the tests, measure the speed and accuracy of an athlete's neurocognitive function.
- The Concussion Vital Signs report auto-scores the athlete's performance using:
 - Subject Scores
 - Percentile Scores
 - Valid Scores

28 | Copyright © 2013. All rights reserved.

PEARSON

What does Neurocognitive Index (NCI) mean?

Concussion Vital Signs Product Training



- The Neurocognitive Index – NCI, reflects the overall neurocognitive functioning of the athlete test taker. It is an average of all the domains into a global summary score.
- Because many concussions are complex and diagnosis is difficult; clinicians should take a multidimensional approach to their assessment.
- The NCI and the other neurocognitive domain scores should be taken in context with the symptom scores, history and physical, as well as other tests and relevant clinical endpoints.

29 | Copyright © 2013. All rights reserved.

PEARSON

Concussion Vital Signs Test Report

Concussion Vital Signs Product Training



- Concussion Vital Signs helps collect information on an athlete's **BASELINE** and **POST-INJURY** status:
 - ***Athlete-Reported Concussion Symptoms***
 - ***Neurocognitive Function***
 - ***Athletes-Reported Concussion History***
- Re-testing at intervals enables a longitudinal view of an athlete's condition.

30 | Copyright © 2013. All rights reserved.

PEARSON

Clinical Domains Measured

Clinical Domains	Clinical Domain Score Calculations	Clinical Domain Description
Neurocognitive Index (NIC)	The average of the Composite Memory, Psychomotor Speed, Cognitive Flexibility, Reaction Time, and Complex Attention Domains.	Measure: An average score derived from the domain scores or a general assessment of the overall neurocognitive status of the patient. Relevance: Summary views tend to be most informative when evaluating a population, a condition category, and outcomes.
Verbal Memory	Verbal Memory is the score for the Verbal Memory Test. $VIM\ Correct\ Hits\ Immediate + VIM\ Correct\ Passes\ Immediate + VIM\ Correct\ Hits\ Delay + VIM\ Correct\ Passes\ Delay$	Measure: How well subject can recognize, remember, and retrieve words. Relevance: Remembering a scheduled test, recalling an appointment, taking medications, and attending class.
Visual Memory	Visual Memory is the score for the Visual Memory Test. $VIM\ Correct\ Hits\ Immediate + VIM\ Correct\ Passes\ Immediate + VIM\ Correct\ Hits\ Delay + VIM\ Correct\ Passes\ Delay$	Measure: How well subject can recognize, remember and retrieve geometric figures. Relevance: Remembering graphic instructions, navigating, operating machines, recalling images, and/or remember a calendar of events.
Psychomotor Speed	Psychomotor Speed is the combined score for both the Finger Tapping and the Symbol Digit Coding Test. $FTT\ Right\ Taps\ Average + FTT\ Left\ Taps\ Average + SDC\ Correct\ Responses$	Measure: How well a subject recognizes and processes information i.e., perceiving, attending/responding to incoming information, motor speed, fine motor coordination, and visual-perceptual ability. Relevance: Distractibility, fitness-to-drive, occupation issues, obsessive concern with accuracy and detail.

PEARSON

Clinical Domains Measured cont.

Clinical Domains	Clinical Domain Score Calculations	Clinical Domain Description
Executive Functioning	Executive Function reflects performance on the Shifting Attention Test. $SAT\ Correct\ Responses - SAT\ Errors$	Measure: How well a subject recognizes set shifting and manages multiple tasks simultaneously. Relevance: Ability to sequence tasks and manage multiple tasks simultaneously as well as tracking and responding to a set of simple instructions.
Cognitive Flexibility	Cognitive Flexibility reflects performance on the Shifting Attention and Stroop Tests. $SAT\ Correct\ Responses - SAT\ Errors - Stroop\ Commission\ Errors$	Measure: How well subject is able to adapt to rapidly changing and increasingly complex set of directions and/or to manipulate the information. Relevance: Reasoning, switching tasks, decision-making, impulse control, strategy formation, attending to conversation.
CPT Correct Responses	CPT Correct Responses is the number of correct responses on the Continuous Performance Test.	Measure: Ability to track and respond to information over lengthy periods of time and/or perform mental tasks requiring vigilance quickly and accurately. Relevance: Self-regulation and behavioral control.
Reaction Time*	Reaction Time* is the average reaction time on parts 2 and 1 of the Stroop Tests. $(ST\ Complex\ Reaction\ Time\ Correct + Stroop\ Reaction\ Time\ Correct) / 2$	Measure: How quickly the subject can react, in milliseconds, to a simple and increasingly complex direction set. Relevance: Driving a car, attending to conversation, tracking and responding to a set of simple instructions, taking longer to decide what response to make.

PEARSON

Reaction Time Detail

Simple Reaction Time*	Simple Reaction Time* is the average reaction time on part 1 of the Stroop Tests. Time required to press the spacebar from the time a word first appears on the display. $Average\ Reaction\ Time\ on\ Part\ 1\ of\ the\ Stroop\ Test$
Choice Reaction Time Correct*	Choice Reaction Time Correct* is the average correct reaction time on the Continuous Performance Test. Time required to press the spacebar from the time a B first appears on the display.
Shifting Attention Correct RT*	Shifting Attention Correct RT* is the average correct reaction time on the Shifting Attention Test.

An * denotes that "lower is better" in the Subject Score column, otherwise higher scores are better. With Percentile scores, higher is always better.

PEARSON

Concussion Vital Signs Concussion Symptom Scale Post-Injury Example

Post-Injury Concussion Symptom Scale	
Athlete Reference/ID: Symptom Scale Example	Test Date GMT: January 28, 2011 17:40:29
Full Name:	Age: 16
Administrator: Athletic Trainer	Language: English (United States)
Total Test Time: 2:14 (min:secs) for all tests in this report	Test Date Local: January 28, 2011 12:40:29
Testing Supervision: Unsupervised Testing	Testing Environment: Alone

[This scale was administered using Concussion Vital Signs](#)

PEARSON

Concussion Symptom Scale

CSI - Symptoms*	Baseline (Jan 18, 2011)		Post-Injury	
	Absent	Present	Absent	Present
1 Headache	0			3
3 Nausea	0			1
5 Poor balance	0			2
6 Dizziness	0		0	
7 Fatigue or loss of energy	0			3
9 Drowsiness or feeling sleepy	0		0	
14 Feeling like "in-a-fog"	0		0	
15 Difficulty concentrating		3		2
16 Difficulty remembering		1		2
10 Sensitivity to light	0		0	
11 Sensitivity to noise	0		0	
17 Blurred vision		1		1
24 Feeling slowed down	0			0

PEARSON


Additional Concussion Symptoms**

	Baseline (Jan 18, 2011)		Post-Injury	
	Absent	Present	Absent	Present
8 Difficulty falling or staying asleep	0		0	
12 Irritability, easily annoyed or frustrated	0			2
13 Sadness		2		1
2 Feeling numbness or tingling	0			4
18 Ringing in the ear	0		0	
19 Neck pain	0		0	
20 More Emotional	0		0	
21 Feeling Nervous	0		0	
22 Feeling anxious or tense		1	0	
23 Feeling Confused	0			1
4 Vomiting	0		0	

Do Symptoms get worse with Physical Activity: No
Do Symptoms get worse with Mental or Academic Activity: No

PEARSON


If a student athlete does not have a baseline, can he/she be given a post-injury test? Concussion Vital Signs Product Training



- **Baseline** testing can serve as a valuable “premorbid” (state prior to condition) point of comparison for the testing that is conducted after the concussion injury.
- However, even if baseline neuropsychological testing has not been performed, post-injury neurocognitive testing can still be a very useful source of information about the effects of the concussion.
- Using standardized PERCENTILE scores can help clinicians identify poor cognitive function performance which can be an important indicator that the brain is not working normally.
- However, there are many reasons test performance can be abnormal, including concussion.

52 | Copyright © 2013. All rights reserved. PEARSON

What combinations of what test scores should cause school personnel/clinicians to pause and look for some underlying conditions? Concussion Vital Signs Product Training



- Every student athlete is different; there is no “one-size fits all” answer to assessing concussion.
- Neurocognitive domain score performance may vary depending on a number of factors that include testing effort, type of blow to the head, location or site of the blow, and the patient’s individual history.
- The Consensus statement on concussion in sport held in Zurich, November 2008 states “...the assessment of cognitive function should be an important component in any return to play protocol.”
- It must be emphasized, however, that NP assessment should not be the sole basis of management decisions; rather it should be seen as an aid to the clinical decision-making process in conjunction with a range of clinical domains and investigational results.”

53 | Copyright © 2013. All rights reserved. PEARSON

Computerized Neuropsychological Testing in the Management of Sports-Related Concussions; Concussion Vital Signs Product Training
Athletic Training & Sports Health Care | Vol. 4 No. 1 2012

**WHEN DO YOU RE-TEST?
CONCUSSION VITAL SIGNS IS OPTIMIZED FOR THE ZURICH GUIDELINES**

According to NATA’s position statement... self-reported symptom resolution should be used as an indicator to begin neurocognitive testing.

When a concussed athlete is asymptomatic with activity of daily living, **he or she should progress through the Zurich Consensus Statement’s graduated return-to play protocol, and before starting full-contact practice a CNT should be administered to obtain objective data to guide the health care provider regarding return-to-play decisions. This allows for minimal testing and allows physical exercise to elicit symptoms prior to CNT administration.**

54 | Copyright © 2013. All rights reserved. PEARSON

Concussion Considerations

Concussion Vital Signs Product Training



- A player with a diagnosed concussion should not be allowed to return to play **on the day of injury... An important consideration in return-to-play is that athletes should not only be symptom free but also should not be taking any medications that may mask or modify the symptoms of concussion.**
- The cornerstone of concussion management is physical and cognitive rest until symptoms resolve and then a **graded program of exertion prior to medical clearance and return to play... If any one or more of these components is present, a concussion should be suspected and the appropriate management strategy instituted:**
 1. Symptoms: Somatic (headache), cognitive (feeling in a fog) and/or emotional symptoms.
 2. Physical signs (loss of consciousness, amnesia)
 3. Behavioral changes (irritability)
 4. Cognitive impairment (slowed reaction)

McCrary et al, Consensus Statement on Concussion in Sport: The 3rd International Conference on Concussion in Sport Held in Zurich, November 2008. Journal of Athletic Training, Vol. 44, No. 4, August 2009.

55 | Copyright © 2013. All rights reserved.

PEARSON

Return to play protocol

Concussion Vital Signs Product Training



- Should follow a stepwise progression as follows:
 - Step 1 No activity until asymptomatic at rest and with exertion
 - Step 2 Light aerobic exercise (walking, stationary bike)
 - Step 3 Sport specific, non-contact activities (running drills)
 - Step 4 Non-contact training drills (passing drills, begin weight lifting)
 - Step 5 Full contact practice (following medical clearance)
 - Step 6 Return-to-play (normal game play)
- Each step above should take 24 hours. If any symptoms occur, the athlete should drop back to the previous level and try to progress again after 24 hours of rest has passed.

McCrary et al, Consensus Statement on Concussion in Sport: The 3rd International Conference on Concussion in Sport Held in Zurich, November 2008. Journal of Athletic Training, Vol. 44, No. 4, August 2009.

56 | Copyright © 2013. All rights reserved.

PEARSON

What is "Executive Functioning"?

Concussion Vital Signs Product Training

- Executive Functioning, sometimes called executive control system, is generally considered a frontal lobe (see blue section of picture) cognitive system that controls and manages other cognitive processes.
- It is considered a higher-order brain function which includes attention, behavioral planning and response inhibition, and the manipulation of information in problem-solving tasks.
- Sometimes referred to as the "command and control" function (frontal lobe), the executive function can be viewed as the "conductor" of many cognitive skills.
- The SAT - Shifting Attention Test (rules, categories, rapid decision-making) results are used to calculate this frontal lobe domain.



57 | Copyright © 2013. All rights reserved.

PEARSON

PEARSON
1-800-627-7271 (USA)
1-866-335-8418 (Canada)

www.pearsonassessments.com

Webinar questions?

Peter.Entwistle@Pearson.com

Anne-Marie.Kimbell@Pearson.com

Concussion Vital Signs®

Product questions?

CVS Assessment Representative

Kent.Willette@Pearson.com

800-627-7271, ext. 267147
952-236-9885

PEARSON
