

# diasend<sup>®</sup> Transmitter

Instructions for use – for diasend<sup>®</sup> users



# Contents

Find your device .....	4-6
About diasend® .....	7
How to start using diasend® Clinic.....	8
Plug diasend® Transmitter into your power outlet.....	8
The diasend® Transmitter is ready to use when.....	8
Connect patient device to diasend® Transmitter.....	9
Transfer data.....	9
Log in to your diasend® Clinic account.....	9
Viewing the data.....	10
Creating a patient profile and saving patient devices.....	10
Displaying the data from the patient devices.....	11
Sharing the data.....	12
Invite patients to upload from home .....	12
Creating reports and printing as PDF .....	12
Explanation of the diasend® Transmitter indicators.....	13
Transmitter behavior when everything is working properly.....	13
List of compatible devices .....	14

# Find your device



## Abbott

FreeStyle	14
FreeStyle Flash	14
FreeStyle Freedom	14
FreeStyle Freedom Lite	14
FreeStyle InsuLinx	14
FreeStyle Lite	14
Freestyle Mini	14
FreeStyle Optium	14
FreeStyle Optium H	14
FreeStyle Optium Neo	14
FreeStyle Papillon	14
FreeStyle Papillon InsuLinx	14
FreeStyle Papillon Lite	14
FreeStyle Papillon Mini	14
FreeStyle Papillon Vision	14
FreeStyle Precision	14
FreeStyle Precision H	14
FreeStyle Xido	14
Optium Mini	14
Optium Xceed	14
Optium Xido	14
Precision Xceed	14
Precision Xtra	14
FreeStyle Libre	14
FreeStyle Navigator	23
FreeStyle Navigator II	14

## Acon Laboratories

On Call® Advanced	15
On Call® Express	15
On Call® Express II	15
On Call® Express Voice	15
On Call® Plus	15
On Call® Resolve	15
On Call® Sure	15
On Call® Sure Sync	15
On Call® Vivid	15

## AgaMatrix

Wavesense Jazz™	15
-----------------	----

## AktivMed

Glucocheck Gold	15
Glucocheck XL	15

## Animas

IR1200	25
IR1250	25
2020	25
OneTouch Ping	25
OneTouch Vibe Plus	25
Vibe	25

## Apex Bio

GlucoSure Extra	15
GlucoSure Max	15

## Arkray

GlucoCard 01	15
GlucoCard Expression	15
GlucoCard Shine	15
GlucoCard Shine Connex	15
GlucoCard Shine Express	15
GlucoCard Shine XL	15
GlucoCard Vital	15
GlucoCard X-meter	15
GlucoCard X-mini plus	15
Relion Confirm	15
Relion Premier Blu	15
Relion Premier Classic	15
Relion Premier Voice	15
Relion Prime	15

## Ascensia (Bayer)

Bayer - Contour next	16
Breeze 2	16
Contour	16
Contour care	16
Contour Fit	16
Contour link	16
Contour next	16

Contour next link	16
Contour next link 2.4	16
Contour next one	16
Contour next USB	16
Contour Plus	16
Contour plus one	16
Contour TS	16
Contour XT	16
Contour USB	16

## Beurer

GL50 evo	16
----------	----

## B. Braun

Omnitest 3	16
Omnitest 5	16

## Cambridge Sensors

Microdot+	16
-----------	----

## CARE-CHEK

TD-4206	16
TD-4255	16


## Dexcom

G4 Platinum	16
G5 receiver	16
Touchscreen receiver G5/G6	16
SEVEN PLUS	16

## Diatesse

XPER	17
------	----





<b>Ekoweb</b>	
 GlucoX (TD-4183) .....	17
 GlucoX Pro (TD-4183).....	17







<b>Everymed</b>	
 EveryX3.....	17

<b>Equil/Wellion</b>	
 MICRO-pump .....	25






<b>FaCare</b>	
 TD-4277.....	17

<b>ForaCare</b>	
 FORA Diamond MINI.....	23
 FORA Diamond Prima.....	17

<b>GE</b>	
 GE100.....	17




<b>GlucoRx</b>	
 HCT .....	17
 Nexus .....	17
 Nexus Mini .....	17
 Nexus Mini Ultra.....	17
 Nexus Voice.....	17
 Q.....	17


<b>GlucoX</b>	
 TD-4183.....	17
 TD-4257.....	17
 TD-4285.....	17



<b>i-SENS</b>	
 alphacheck professional NFC/ACURA PLUS..	17
 alphacheck professional.....	17
 COOL/Bravo .....	17
 CareSens BGMS/TEE2.....	17
 CareSens Dual .....	17


CareSens N.....	17
CareSens N POP.....	17
CareSens N Premier.....	17
CareSens N Voice.....	17
CareSens Pro.....	17
NoCoding1 Plus BLE.....	17
Tee2 Plus.....	17










<b>iCare</b>	
 PalmDoc I .....	18
 PalmDoc II .....	18

<b>Infopia</b>	
 Element.....	18
 Finetest Lite .....	18
 GlucoLab .....	18
 GluNEO.....	18


<b>Inside Biometrics</b>	
 KEYA Smart.....	18


<b>Insulet</b>	
 OmniPod® System.....	24
 OmniPod® DASH™ System.....	24

















<b>Intuity</b>	
 POGO® Automatic™ .....	18

<b>LifeScan</b>	
 OneTouch Ping .....	18
 OneTouch Select.....	18
 OneTouch Select Plus .....	18
 OneTouch Select Plus Flex.....	18
 OneTouch Ultra .....	18
 OneTouch Ultra2.....	18
 OneTouch UltraEasy .....	18
 OneTouch UltraMini .....	18
 OneTouch UltraSmart.....	18



OneTouch Ultra Plus Reflect .....	18
OneTouch Verio .....	18
OneTouch Verio (old) .....	18
OneTouch Verio Flex .....	18
OneTouch Verio IQ.....	18
OneTouch Verio Pro .....	18
OneTouch Verio Reflect .....	18
OneTouch Verio Sync.....	18
OneTouch Vita.....	18

<b>Lifesmart</b>	
 LS-946 (TD-4279) .....	19


<b>Medcore</b>	
 Care .....	19
 GlucoSure Extra.....	19
 GlucoSure Max.....	19

<b>Menarini</b>	
 GLUCOCARD MX.....	19
 GLUCOCARD SM.....	19
 GLUCOFIX iD .....	19
 GLUCOFIX miò .....	19
 GLUCOFIX miò Plus.....	19
 GLUCOFIX Premium.....	19
 GLUCOFIX TECH .....	19
 GlucoMen areo .....	19
 GlucoMen areo 2K .....	19
 GlucoMen G+ .....	19
 GlucoMen Gm .....	19
 GlucoMen LX.....	19
 GlucoMen LX2.....	19
 GlucoMen LX Plus .....	19
 GlucoMen READY .....	19
 GlucoMen Visio .....	19


<b>Mendor</b>	
 Discreet .....	20


<b>Nipro</b>	
 4SURE Smart .....	20
 4SURE Smart Duo.....	20


	<b>Nova</b>	
	Max Link.....	20



	<b>Pharma Supply</b>	
	Advocate Redi-Code+.....	20
	Max Plus.....	20


	<b>Pendiq</b>	
	pendiq.....	27


	<b>Pic Solution</b>	
	BeGlic.....	20
	GoGlic.....	20
	Pic Gluco Test.....	20
	Pic Gluco Test Diary.....	20
	Pic Safe Gluco Monitort.....	20


	<b>Prodigy</b>	
	AutoCode.....	20


	<b>Roche</b>	
	Accu-Chek Active.....	21
	Accu-Chek Aviva.....	24
	Accu-Chek Aviva Nano.....	24
	Accu-Chek Aviva Combo.....	24
	Accu-Chek Aviva Connect.....	21
	Accu-Chek Aviva Expert.....	24
	Accu-Chek Compact.....	24
	Accu-Chek Compact Plus.....	24
	Accu-Chek Compact Plus GT.....	24
	Accu-Chek Guide.....	21
	Accu-Chek Guide Me.....	21
	Accu-Chek Instant.....	21
	Accu-Chek Instant S.....	21
	Accu-Chek Mobile (IR).....	24
	Accu-Chek Mobile (USB).....	21
	Accu-Chek Nano.....	24
	Accu-Chek Performa.....	24
	Accu-Chek Performa Nano.....	24
	Accu-Chek Performa Connect.....	21


	Accu-Chek Insight.....	26
	Accu-Chek Spirit.....	26
	Accu-Chek Spirit Combo.....	26


	<b>Sanofi</b>	
	BGStar.....	21
	MyStar Extra.....	21

	<b>Senseonics</b>	
	Eversense CGM Transmitter.....	21


	<b>Smiths Medical</b>	
	CoZmo 1700.....	25

	<b>SOOIL</b>	
	DANA Diabecare R.....	27

	<b>TaiDoc</b>	
	TD-4183.....	21
	TD-4141.....	21
	TD-4206.....	21
	TD-4255.....	21
	TD-4257 (v1).....	21
	TD-4257 (v2).....	21
	TD-4277.....	21
	TD-4279.....	21
	TD-4280.....	21
	TD-4285.....	21
	TD-4289.....	21

	<b>Tandem</b>	
	t:flex.....	25
	t:slim.....	25
	t:slim X2.....	25
	t:slim G4.....	25

	<b>Terumo</b>	
	Finetouch.....	21


	<b>Trividia</b>	
	TRUE METRIX.....	21
	TRUE METRIX AIR.....	21
	TRUResult.....	21
	TRUEyou.....	21
	TRUEyou mini.....	21


	<b>URIGHT</b>	
	TD-4255.....	22
	TD-4277, PLUS One (TD-4277).....	22
	TD-4280.....	22


	<b>ViCentra</b>	
	Kaleido.....	25


	<b>VivaChek</b>	
	Fad.....	22

	<b>VTRUST</b>	
	TD-4279.....	22

	<b>Wellion</b>	
	Calla Dialog.....	22
	Calla Light.....	22
	Calla Mini.....	22
	Calla Premium.....	22
	Galileo Compact.....	22
	Galileo GLU/CHOL.....	22
	Galileo GLU/KET.....	22
	Leonardo Compact.....	22
	Leonardo GLU/CHOL.....	22
	Leonardo GLU/KET.....	22
	LUNA Duo.....	22
	Linus.....	22

	<b>XPER Technology</b>	
	E1 (TD-4141).....	22
	A1 (TD-4183).....	22
	Procheck Advance (TD-4206).....	22
	I1 (TD-4289).....	22

	<b>Ypsomed</b>	
	mylife Pura.....	22
	mylife Unio.....	22

	mylife OmniPod.....	25
---	---------------------	----

# About diasend®

*diasend® offers health care providers an online solution that collects and stores all their diabetes patients' data centrally, without worrying about any software installation. No matter if the device transmits using a cable, Bluetooth, infrared or a mobile app, all uploaded information will immediately be shown online in a secure diasend.com account. All data from multiple devices will be consolidated and presented in one report.*

## Intended use

diasend® is indicated for use by individuals or health care professionals in the home or health care facilities for transmission of data from home monitoring devices such as glucose meters and insulin pumps to a server database to support diabetes management. diasend® is indicated for professional use and over-the-counter sales.

## Warnings


All patient medical diagnoses and treatments are to be performed under supervision of qualified healthcare professionals! diasend® is not intended to provide automated treatment decisions or to be used as a substitute for professional healthcare judgment. diasend® is not intended for emergency calls or for transmission or indication of any real-time alarms or time-critical data! diasend® is not intended as a substitute for direct medical supervision or emergency intervention.

# How to start using diasend® Clinic

- 1 Plug diasend® Transmitter into your power outlet



- 2 The diasend® Transmitter is ready to use when:



*i* Note: if you have a USB hub connected to the diasend® Transmitter, please note that the USB hub must be powered by its own power supply.

*i* Note: Please see page 11 for a thorough explanation of the indicators on the diasend® Transmitter.

### 3 Connect patient device to the diasend® Transmitter



Each patient device connects in one of the following three manners:

- USB cable (USB port or USB hub)
- Bluetooth (requires specific Bluetooth device)
- Infrared (built-in IR receiver or for some devices a SmartPix cable).

*i* Please see page 11 of this guide for specific instructions on how to upload all compatible patient devices.

### 4 Transfer data



1 Connect the patient device - the blue send indicator will flash rapidly whilst uploading the data.

2 Wait until the blue send indicator stops flashing rapidly and the green ok indicator is lit.

3 Once the green ok indicator is lit, you can log in and view the data at [www.diasend.com](http://www.diasend.com) or continue uploading more patient devices. Click Reload page if you were already logged in. Please note that since diasend® is cloud-based, you can review the uploaded data at any time.

### 5 Log in to your diasend® Clinic account



Go to [www.diasend.com](http://www.diasend.com) to log in to diasend® and retrieve the uploaded data.

## 6 Viewing the data

**Clinic**

← In the clinic tab it's possible to view the data without creating a patient profile.

Recent uploads at clinic

In this view you see all the devices that you have uploaded today at the clinic.

Make sure that you select the right device for the current patient.  
Check that the serial number is the same as the one on the device.

Reload page ← Click here after you have performed a transmission with your diasend® Transmitter.

Show selected devices Create patient profile for selected devices Save device(s) to existing patient profile

Time	Patient	Serial number	Last Glucose	Last Ketones	Select
09:42	Pierce, Annette	7910778 (Bayer Contour)	13/03/2011	-	<input type="checkbox"/>
10:16	Powell, Arlene	8229504 (Bayer Contour)	14/03/2011	-	<input type="checkbox"/>
10:40	Patterson, Dave	4188077 (Bayer Contour)	02/08/2011	-	<input type="checkbox"/>
16:25	Garner, Dennis	9345638 (Bayer Contour)	19/07/2011	-	<input type="checkbox"/>
09:09	Gilbert, Charles	3854310 (Bayer Contour)	13/06/2011	-	<input type="checkbox"/>

1 Tick the box of the upload you wish to view.

2 Click on **Show selected devices** to view the upload.

*i* Please note that all uploads are removed from the Clinic view after midnight the same day, but no data is erased.

## 7 Creating a patient profile and saving patient devices

There are two ways to create a patient profile. In the Clinic tab and in the patient tab.

**Clinic**

← How to create a patient profile in the clinic tab:

Recent uploads at clinic

In this view you see all the devices that you have uploaded today at the clinic.

Make sure that you select the right device for the current patient.  
Check that the serial number is the same as the one on the device.

Reload page ← Click here after you have performed a transmission with your diasend® Transmitter.

Show selected devices Create patient profile for selected devices Save device(s) to existing patient profile

Time	Patient	Serial number	Last Glucose	Last Ketones	Select
09:42	Pierce, Annette	7910778 (Bayer Contour)	13/03/2011	-	<input type="checkbox"/>
10:16	Powell, Arlene	8229504 (Bayer Contour)	14/03/2011	-	<input type="checkbox"/>
10:40	Patterson, Dave	4188077 (Bayer Contour)	02/08/2011	-	<input type="checkbox"/>

1 Tick the box of the upload you wish to create a patient profile for.

2 Click here to create a patient profile for the selected device.

*i* We recommend this method of registering devices to a patient, in order to avoid registering the patient device manually in the patient tab.

## Patients ← How to create a patient profile in the patient tab:

1 Click on **New patient**.

2 The mandatory fields are first and last names. You will find further options under **Expand registration form**.

3 To allow a patient to view their data from home, tick the box: "Allow patient to log in to the diasend® system and view reports"

- If you do this in the registration process you will have to create a username and a password for the patient.
- If you tick the box after you have created a patient account, you will notice that diasend® has created a username for you. In this case, you only need to create a password.

*Please note that this only allows the patient to log in and view previously uploaded data. It does however not mean that a patient can log in and upload data from home. If patients want to upload data from home they need to create their own diasend® Personal account and share their data with the clinic.*

4 When you're done click on: **"Save patient information"**.


## 8 Displaying the data from the patients devices

Summary   Glucose   CGM   Insulin   Comparison

The data is displayed in the five tabs above, where you can access the glucose, CGM and insulin data. Please refer to the diasend® Report reference guide on our website for further details on our reports.


## 9 Sharing the data

Patients with their own diasend® Personal account can upload their patient devices at home, and they can share their uploaded data with up to two different healthcare facilities. To enable patients to share their data with you, they will require your unique identifier in diasend® called “Clinic ID”. You will find it under “Account profile” > “Patient links”.




## 10 Invite patients to upload from home

To make it easier for clinics and patients to share data we've added an invite button in the web application. This means clinics can invite patients to upload data from home.




## 11 Creating reports and printing as PDF

There are two ways to generate a report as a pdf. You can either click the “Print to PDF” icon to immediately generate a report of the current view. Or you can click the PDF Wizard button to generate a report with preset variables. You can adjust this tool under “Settings”.





## Explanation of the diasend® Transmitter indicators:



**ok**

- SLOW FLASHING  
The ok indicator flashes slowly when the transmitter is connecting to the GSM network.
- RAPID FLASHING  
The Transmitter is updating itself. The ok indicator flashes rapidly and during a sustained period of time. You can start using the transmitter once the blue send indicator flashes slowly again.
- STEADY LIGHT ON  
The indicator signals a successful upload.


**send**

- SLOW FLASHING  
The transmitter is ready for use.
- RAPID FLASHING  
Data is being uploaded.

**error**

- STEADY LIGHT ON  
The indicator signals that the latest upload was not successful. This is most likely due to temporary issues related to the wireless GSM network. Please try to upload the patient device again.
- RAPID FLASHING  
The indicator flashes alongside the ok indicator when the Transmitter cannot find a network. Please note that you can check the strength and quality of your signal under "Account Profile" on your account.

## diasend® Transmitter behavior when everything is working properly:




- 1** The Transmitter is connected to a power outlet.
- 2** The green **ok** indicator flashes slowly for approx 2 mins. when starting up.
- 3** The blue **send** indicator starts flashing slowly. The transmitter is now in stand-by mode, i.e ready for use.
- 4** A device is connected to the transmitter.
- 5** The blue **send** indicator starts flashing rapidly as the data is being uploaded.
- 6** If the upload was successful, the green **ok** indicator has a steady light on and the blue **send** indicator starts flashing slowly again. The transmitter is now ready for another transmission.

# List of compatible devices


Below is a list of compatible devices with diasend®. Unless explicit upload instructions are provided for the actual device, you simply need to connect it to the transmitter or USB hub with its specific USB cable and it will upload automatically. (Visit diasend.com to see a current list of compatible devices and known limitations.)

*Please note that some of the following patient devices may not be available in your market.*



## Glucose meters and CGMs - connected with USB cable

*Please make sure that the USB hub is connected to a power outlet. And that the USB hub is connected to the diasend® Transmitter.*



### Abbott

FreeStyle FreeStyle Flash FreeStyle Freedom FreeStyle Freedom Lite FreeStyle Lite FreeStyle Mini FreeStyle Optium	FreeStyle Optium H FreeStyle Papillon FreeStyle Papillon Lite FreeStyle Papillon Mini FreeStyle Papillon Vision FreeStyle Xido Optium Mini
---	--



FreeStyle InsuLinx FreeStyle Libre*	FreeStyle Optium Neo FreeStyle Papillon InsuLinx
--	---



\*Compatibility for FreeStyle Libre is restricted to Sweden and Norway by Abbott.


FreeStyle Navigator II
------------------------




FreeStyle Precision FreeStyle Precision H Optium Xceed	Optium Xido Precision Xceed Precision Xtra
--	--




# Glucose meters and CGMs - connected with USB cable




# Glucose meters and CGMs - connected with USB cable




# Glucose meters and CGMs - connected with USB cable




# Glucose meters and CGMs - connected with USB cable



# Glucose meters and CGMs - connected with USB cable



## Glucose meters and CGMs - connected with USB cable






# Glucose meters and CGMs - connected with USB cable

**Roche Accu-Chek**


Active Mobile (USB)	Guide Guide Me
Aviva Connect	Instant
Performa Connect	Instant S

Micro USB 

Once connected, upload should start within 20 seconds after connecting the cable. If not, go to the main menu on the meter. Go to **My data > Data Transfer > Reports**, and connect the meter to the USB cable again. (Please do not hold down the two arrow buttons.)

**Sanofi**

BGStar  
MyStar Extra

Sanofi USB 

marked with: ★

**Senseonics**


Eversense CGM Transmitter

**NOTE!**  
The Eversense Smart Transmitter (CGM) can not be uploaded via diasend® Transmitter.

The data from the Eversense CGM can only be uploaded to diasend® via **"Connect app"** in the patient's diasend® account. Please see further instructions in the user guide for patients at: <https://support.diasend.com/hc/sections/360000427618>

**TaiDoc**

TD-4141	TD-4257(v2)
TD-4183	TD-4279
TD-4206	

TaiDoc Strip Connector USB 

TD-4255	TD-4280
TD-4257(v1)	TD-4285

GlucRx Nexus Voice USB 

TD-4289 

Micro USB

TD-4277 


Mini USB

**Terumo**

Finetouch

Terumo USB 

**Trividia**







TRUEresult	Nipro TRUEresult USB docking station 	Connect and place the meter in the docking station.
TRUE METRIX		
TRUE METRIX AIR		
TRUEyou		

TRUEyou mini

Nipro TRUEyou mini USB docking station 

Connect and place the meter in the docking station.

# Glucose meters and CGMs - connected with USB cable



## Bluetooth glucose meters

To be able to transmit data from Bluetooth devices, you need a Bluetooth USB adapter. Bluetooth as a transmission method requires the device to be placed within a close proximity to the diasend® Transmitter.



### Abbott

FreeStyle Navigator



On the meter, go to **Menu > System > Upload**.  
Select and confirm to activate.

### FORACARE

FORA Diamond MINI



1. To enable Bluetooth, turn the meter on by holding down its power button.
2. Turn the meter off by holding down its power button. The LED indicator in the lower right-hand corner of the meter flashes blue to indicate that Bluetooth is activated. It then turns solid to indicate that the data is being uploaded.

Note: The supported version is the one with integrated Bluetooth.



# Infrared glucose meters

Set the meter in upload mode by following the instructions below. Place the glucose meter's top end directed towards the IR eye on the diasend® Transmitter.

*Please make sure to keep the device completely still during the upload. Do not remove the device until the green ok indicator is lit.*



## Roche

Accu-Chek Aviva  
Accu-Chek Aviva Nano  
Accu-Chek Nano  
Accu-Chek Performa  
Accu-Chek Performa Nano



Activate IR transmission by pressing and holding the two arrow buttons until two arrows appear in the display.

Accu-Chek Compact



Activate IR transmission by pressing and holding the meter's two topmost buttons until two arrows appear in the display.

Accu-Chek Compact Plus  
Accu-Chek Compact Plus GT



Activate IR transmission by pressing and holding the buttons S and M below the display until two arrows appear in the display.

Accu-Chek Aviva Combo  
Accu-Chek Aviva Expert



In the main menu, go to **My data** and select **Data transfer**.

Accu-Chek Mobile (IR)



1. Press and hold the two arrow buttons.  
2. In the **data transfer** menu, select the first alternative (**for analysis**).  
NOTE: If the meter's serial number begins with U1, please refer to the instructions for Accu-Chek Mobile (USB) on page 19.

## Insulet

OmniPod® System

Mini USB



OmniPod® DASH™ System

Micro USB



1. Plug DASH™ into diasend® Transmitter with the Micro USB cable.  
2. Tap Export on DASH™.  
3. Wait for diasend® Transmitter to confirm the upload was successful.

# Insulin pumps - connected with standard USB cable

## Equil/Wellion

MICRO-pump

Micro USB



Please note that the serial number used to connect the pump to the diasend® system is the serial number on the patch pump.

## Tandem

t:flex  
t:slim

t:slim X2  
t:slim G4

Micro USB



## ViCentra

Kaleido

Micro USB



Go to the yellow menu and choose data upload.

## Ypsomed

mylife OmniPod

Mini USB



## Infrared insulin pumps

For the Animas pumps: Please make sure that the pump is suspended and the screen is illuminated when transmission begins (the screen does however not need to be illuminated during the entire transmission). Keep the pump completely still and do not remove the pump until the green ok indicator is lit.



## Animas

IR1200  
IR1250  
2020  
OneTouch Ping  
OneTouch Vibe Plus  
Vibe



1. Suspend the pump. (**MAIN MENU**, scroll to **Suspnd/Resum** and press **OK**; **Suspend** is highlighted).
2. Press **OK**.
3. Place the pump with its back towards the IR eye on the Transmitter.
4. When the transmission is finished: Resume pump (**MAIN MENU**, scroll to **Suspnd/Resum** and press **OK**; **Resume** is highlighted). Press **OK**.

## Smiths Medical

CoZmo 1700



1. Press the down arrow button in order to activate the pump.
2. Place the pump with its back towards the centre of the front panel on the diasend® Transmitter.

Note: The battery in the CoZmo insulin pump must be well-charged in order to transmit successfully!



# Insulin pumps - connected with infrared cables

## Compatible infrared cables


Roche Accu-Chek Smart Pix cable



Roche Accu-Chek Smart Pix 2 cable



Roche  
Accu-Chek Insight



Smart Pix 2  
Note! Connect the handset of your pump with a micro USB cable to the Smart Pix 2.

1. Connect the Smart Pix 2 cable to the Transmitter and wait for the blue indicator on the Smart Pix 2 to start flashing slowly.
2. Connect the handset to the Smart Pix 2 with a micro USB cable.
3. In the handset menu select **Connect to PC** and press **OK**.
4. Ensure that the pump is within reach of the handset.
5. The blue indicator on the Smart Pix 2 is on while the data is being transmitted from the device to the Smart Pix 2.
6. The blue indicator turns itself off and both Smart Pix 2 and device emit a sound to signal that the data is being transmitted from the Smart Pix 2 to the Transmitter. Do not disconnect the Smart Pix 2 yet!
7. Wait for the Transmitter to signal the upload was successful (green ok indicator has a steady light).

Roche  
Accu-Chek Spirit  
Accu-Chek Spirit Combo



Smart Pix or Smart Pix 2  
Note! Smart Pix software version 3.02 or higher is required when using Smart Pix.

1. Connect the Smart Pix to the Transmitter and wait for the blue indicator on the Smart Pix to start flashing slowly.
2. Put the pump in **Data Transfer** mode (stop the pump and enter the **Data Transfer** menu).
3. Place the pump with its IR eye towards the Smart Pix. The blue indicator signals the data is being transmitted to the Smart Pix.
4. Wait for the blue indicator to start flashing slowly again or stop flashing altogether.
5. The data is now being transmitted from the Smart Pix to the Transmitter. Do not disconnect the Smart Pix yet!
6. Wait for the Transmitter to signal the upload was successful (green ok indicator has a steady light).

## Bluetooth insulin pumps

To be able to transmit data from Bluetooth devices, you need a Bluetooth USB adapter. Bluetooth as a transmission method requires the device to be placed within a close proximity to the diasend® Transmitter.



Bluetooth USB adapter

SOOIL


DANA Diabecare R



1. Select **Main menu**.
2. Go to **Setting**.
3. Select **Discovery**.
4. The send indicator should begin to flash rapidly. If it doesn't and the pump turns itself off, try again.
5. Wait until the send indicator stops flashing rapidly and the ok indicator is lit.

## Insulin pens - connected with standard USB cable

To be able to transmit data from an insulin pen, please follow the instructions below.

 Please make sure that the USB hub is connected to a power outlet. And that the USB hub is connected to the diasend® Transmitter.



Pendiq

pendiq

Mini USB



Make sure that the cable is properly inserted. The insulin pen's control lamp will start flashing, and the upload will be initiated after 10-15 seconds.



Glooko AB  
+46 (0)31-762 08 88  
Nellickevägen 20  
412 63 Göteborg  
SWEDEN

Support  
+44 (0)20 7795 8191  
+46 (0)31-762 08 88  
uk@diasend.com  
www.diasend.com

A2032v2WIFU\_GB\_diasend\_Transmitter (2022-02-01)