# EasyMax™ 102

QuickStart Guide

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Démarrer EasyMax™ 102

– notice succincte

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## EasyMax™ 102 QuickStart Guide

Experiment Setup, Data Visualization and Export

METTLER TOLEDO


# EasyMax™ 102 Quick Start Guide

This document will guide you through the steps of experiment setup, changing parameters in real time, data visualization and data export on the EasyMax™ 102.

For in-depth information on installation and use of the EasyMax™ 102, please refer to the installation guide and the operating instructions.

## 1. Experiment Setup

- Insert the reactor into the thermostat.
- You have the choice of magnetic or mechanical stirring. If you select magnetic stirring, place the magnetic bar into the reactor.


- Fill in starting materials, solvents, reagents, etc.

- Secure the cover of the 100-mL reactor with the chain clamp. Make sure the Tr (temperature of the reactor contents) sensor is immersed at least 1.5 cm into the reaction mass to provide correct data. When working with one zone only, you may still want to plug in the Tr sensor for the unused zone to avoid an error message for that zone. Further information can be found in the Operating Instructions.

- Switch on the EasyMax™ 102 by pressing the push button on the front panel.
- After startup, the touchscreen displays the **System Overview**.


In the **System Overview** you can monitor and change for both zones:


- Tr** – temperature of the reactor contents  
**Tr-Tj** – heat flow trending  
**Tj** – jacket temperature  
**R** – stirrer speed

## 2. Running an Experiment

- Starting an experiment is done in Single Zone View. Switch to **Single Zone View** by touching either the **Zone 1 View** or **Zone 2 View** buttons (see above). The **Experiment** button lets you name, start, and stop an EasyMax™ 102 experiment.

- Stirrer
- Touching the **R** button opens a numeric keypad for entering the desired stirrer speed. Mechanical stirring is the default whenever a stirrer is plugged in, otherwise magnetic stirring is activated.

- Heating in Tr Mode

Touching the **Tr** button opens a similar numeric keypad for entering the desired temperature of the reactor contents. While in Tr mode, the jacket temperature Tj is also monitored and recorded.

- The **Advanced** button opens the menu to enter heating or cooling ramps. The desired end temperature of the reactor contents can be approached over a predefined length of time or at a given rate. The start of a ramp may be postponed.

- Heating in Tj Mode

Identical features for controlling the jacket temperature can be accessed by touching the **Tj** button. While in Tj mode, the temperature of the reactor contents Tr is also monitored and recorded.

- Reflux or Distillation Mode


Heating in reflux mode maintains Tj above Tr by the margin  $\Delta T$  and up to the maximum temperature Tend.

### 3. Show Graph

- The graph is also accessed from Single Zone View.

**Show graph** opens this page. You can scale the x-axis with the button in the lower left hand corner.

**Advanced** allows you to customize both y-axes as well as trend colors.

### 4. Data Export

- Touching **Export** in Single Zone View allows you to export experiments with a USB storage device. Data is exported as a \*.csv file, which can be opened by MS Excel® or similar programs. You may limit the export to a specific area within an experiment.

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©06/2008 Mettler-Toledo AG
## Mettler-Toledo AutoChem Inc.
## Mettler-Toledo AG, AutoChem