EnCor ENSPIRE™

BREAST BIOPSY SYSTEM

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ENCOR ENSPIRE™ Breast Biopsy System Instructions for Use

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

A. DEVICE DESCRIPTION

The ENCOR ENSPIRE™ Breast Biopsy System provides control operations for specialized biopsy instruments intended to acquire tissue samples of suspected breast abnormalities. The ENCOR ENSPIRE™ Biopsy System is designed to recognize individual biopsy drivers and provide a user interface and to accept inputs when used for specific procedures. Each biopsy driver type performs unique functions.

The ENCOR ENSPIRE™ Breast Biopsy System may be used with ENCOR®, ENCOR® MRI and ENCOR® 360 drivers, foot pedals and probes. ENCOR® Biopsy Probes, Tubing Cassettes, Canisters, adapters and needle guides are available separately. Reference **Figure 1 – Figure 4** for photographs of the above mentioned accessories.

Figure 1 – ENCOR® and ENCOR® 360 Biopsy Probes

Figure 2 – ENCOR® 360 Driver and Probe

Figure 3 – ENCOR® Driver and Probe

Figure 4 – ENCOR® Foot Pedal

B. INDICATIONS FOR USE

The ENCOR ENSPIRE™ Breast Biopsy System is indicated to provide breast tissue samples for diagnostic sampling of breast abnormalities.

- It is intended to provide breast tissue for histologic examination with partial or complete removal of the imaged abnormality.
- It is intended to provide breast tissue for histologic examination with partial removal of a palpable abnormality.

The extent of a histologic abnormality cannot always be readily determined from palpation or imaged appearance. Therefore, the extent of removal of the palpated or imaged evidence of an abnormality does not predict the extent of removal of a histologic abnormality, e.g., malignancy. When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedures.

In instances when a patient presents with a palpable abnormality that has been classified as benign through clinical and/or radiological criteria (e.g. fibroadenoma, fibrocystic lesion), the ENCOR ENSPIRE™ Breast Biopsy System may also be used to partially remove such palpable lesions. Whenever breast tissue is removed, histological evaluation of the tissue is the standard of care. When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedures.

C. CONTRAINDICATIONS

1. This device is not intended for use except as indicated.
2. The ENCOR ENSPIRE™ Breast Biopsy System is contraindicated for those patients where, in the physician's judgment, there is an increased risk of complications associated with percutaneous removal of tissue samples.

D. WARNINGS

1. The ENCOR ENSPIRE™ Breast Biopsy System must be properly grounded to ensure patient safety. The system is supplied with a medical grade power cord with AC plug. Do not connect the included power cord to extension cords or three-prong to two-prong adapters. To avoid the risk of electric shock, this equipment must only be connected to supply mains with protective earth.
2. To minimize interference with other equipment, cables should be positioned in such a manner to prevent contact with other cables.
3. Use of accessories not compatible with the ENCOR ENSPIRE™ Breast Biopsy System may create potentially hazardous conditions.
4. Only use ENCOR® and ENCOR® MRI drivers with script version 1.19 or greater, or ENCOR® 360 drivers with script version 1.05 or greater with the ENCOR ENSPIRE™ Breast Biopsy System. The system is not compatible with earlier driver scripts. The script version is identified on the touch screen display during system initialization.
5. The ENCOR ENSPIRE™ Breast Biopsy System console may not be placed in an MRI suite. Place the console outside of the MRI suite and use the appropriate ENCOR® MRI accessories when performing a biopsy under MRI guidance.
6. No modification of this equipment is allowed. Do not remove the ENCOR ENSPIRE™ Breast Biopsy System housing. Removal of the housing may cause electrical shock. Contact Bard for service.
7. The ENCOR ENSPIRE™ Breast Biopsy System is not classified as an AP or APG classified device. The system is not suitable for use in the presence of flammable anesthetic.

E. PRECAUTIONS

1. This equipment should only be used by a physician trained in its indicated use, limitations, and possible complications of percutaneous needle techniques.
2. Locate the ENCOR ENSPIRE™ Breast Biopsy System as far as possible from other electronic equipment to minimize interference.
3. Inspect accessories and cords for breaks, cracks, nicks, or other damage before every use. If damaged, do not use. Failure to observe this precaution may result in injury or electrical shock to the patient or to the operator.
4. Inspect tubing connections to the vacuum canister and the vacuum tubing cassette to ensure proper vacuum levels are achieved and maintained during use.
5. Inspect the vacuum canister to ensure the lid is secure and that no damage has occurred during shipping or installation. A heavily scratched canister can break during use.
6. Do not leave the ENCOR ENSPIRE™ Breast Biopsy System powered on overnight. Damage may occur to the vacuum or vacuum and rinse tubing cassette.
7. Connect the power cord to a hospital grade wall outlet having the correct voltage or product damage may result.
8. Patients who may have a bleeding disorder or who are receiving anticoagulant therapy may be at increased risk of complications.
9. As with any biopsy instrument, there is a potential for infection.
10. All breast biopsies should be performed under imaging guidance to confirm the probe position relative to the target region to be sampled and to help mitigate the occurrence of a false negative biopsy.
11. When performing a biopsy with ENCOR® and ENCOR® MRI probes, the orientation of the sample notch is dictated by the image guidance selected. Prior to initiating the procedure, confirm that the sample notch orientation is correct for the image guidance being used.
12. Ensure that the ENCOR ENSPIRE™ Breast Biopsy System is positioned in such a way that the power cord and retainer are accessible. In the event that the system power switch is inoperable, release retainer and remove cord to shut off system power.

F. POTENTIAL COMPLICATIONS

1. Potential complications are those associated with any percutaneous removal/biopsy technique for tissue collection. Potential complications are limited to the region surrounding the biopsy site and include hematoma, hemorrhage, infection, a non-healing wound, pain and tissue adherence to the biopsy probe while removing it from the breast.

G. EQUIPMENT REQUIRED

The following equipment is required for a biopsy procedure:

- Appropriate imaging modality and accessories
- ENCOR®, ENCOR® MRI or ENCOR® 360 Driver
- ENCOR®, ENCOR® MRI or ENCOR® 360 Probe
- ENCOR ENSPIRE™ Vacuum Tubing Cassette or ENCOR ENSPIRE™ Vacuum and Rinse Tubing Cassette
- Vacuum Canister
- Saline (optional)
- Tissue Marker (optional)
- Surgical gloves and drapes
- Local anesthetic
- Scalpel
- Other equipment as necessary

H. DIRECTIONS FOR USE

Unpacking Instructions

Thoroughly inspect the packaging materials and the ENCOR ENSPIRE™ Breast Biopsy System for any signs of transit damage. Do not use the unit if it appears to have any damage. Contact Bard Customer Service for instructions on how to report shipping damage.

The ENCOR ENSPIRE™ Breast Biopsy System is delivered assembled with the exception of the touch screen which is detached and packaged separately. Carefully remove the main body and touch screen from the packaging. Unlock the release knob on the touch screen as illustrated in **Figure 5**. Align the connecting track on the back of the display screen and slide in to place on the arm of the main body as illustrated in **Figure 6**. Attach the power, ground and communication cables from the main body to the touch screen. Attach the ground cable to the ground stud and hand tighten the supplied nut. Use the supplied flat head screwdriver to tighten the screws on the communication cable. Reference **Figure 7** for an illustration of proper communication and power cable installation locations.


Figure 5 – Release knob on the touch screen

Figure 6 – Attaching the touch screen to the main body

Figure 7 – Connecting the touch screen cables
Note: Actual connector style may differ from illustration

Procedure Setup

Prior to transporting the system, unlock the system casters and lower the tray and monitor to their lowest positions. Align the tray and monitor to face the front of the system. The monitor and tray should be positioned as illustrated on the cover of this IFU. When transporting the system, push or pull using the handles located on the front and back of the tray. Do not transport the system by pushing or pulling sideways. To transport the system outside of normal hospital conditions (e.g. outdoors or to another facility) use Bard supplied packaging.

1. Position the ENCOR ENSPIRE™ Breast Biopsy System as required for safe use and plug the power cord (located on the rear of the console) into a hospital grade wall outlet.
 - a. The system includes four casters with locking brakes on the front casters. After the system is positioned, lock the casters as illustrated in **Figure 8**.
 - b. The touch screen height may be adjusted by releasing the locking lever on the main post and pressing the release button on the top of main body pole as illustrated in **Figure 9**. The touch screen tilt may be adjusted by loosening the lock behind the touch screen and moving the screen to the desired position as illustrated in **Figure 10**.
 - c. The adjustable tray can accommodate two drivers as well as driver cables and probe tubing. The position of the tray may be adjusted by releasing the locks on the tray and rotating, lifting or lowering the tray to the desired height as illustrated in **Figure 11**.

Figure 8 – Locking the system casters

Figure 9 – Adjusting the touch screen height

Figure 10 – Adjusting the touch screen tilt

Figure 11 – Adjusting the tray position

2. Install the vacuum canister and tubing cassette per the directions provided in the ENCOR ENSPIRE™ Vacuum and Rinse Tubing Cassette Instructions for Use.

- Turn on the main power via the power switch on the rear of the console as illustrated in **Figure 12**.

Figure 12 – Main power switch

- Turn on the standby power via the switch on the front of the system as illustrated in **Figure 13**.

Figure 13 – Standby power switch

- If desired, connect the appropriate driver to the rectangular connector on the rear of the console illustrated in **Figure 14**. The active driver is connected by a green LED light above the connector. Reference the applicable driver Instructions for Use for further information of driver setup. If desired, connect the appropriate foot pedal to the connector on the rear of the console as illustrated in **Figure 14**.

Figure 14 – Driver and foot pedal connectors

- Follow the touch screen prompts. Install and calibrate the biopsy probe. Reference the applicable probe Instructions for Use for further information on probe setup.

Operating the system via the touch screen

- Complete the biopsy procedure per the “Directions for Use” section of the Instructions for Use packaged with the driver being used.
- The ENCOR ENSPIRE™ Breast Biopsy System status is displayed on the touch screen (**Figure 15**). The sample acquisition options may be selected via this screen.
 - To change the sampling pattern, place a finger on desired clock position within the circle on the touch screen as illustrated in **Figure 15**.
 - To change the sampling from full to half, touch the “half sample” button as illustrated in **Figure 15**. Touch the “full sample” button to return to full sampling.
 - To change the sample mode from normal tissue to dense tissue, touch the “dense tissue” button as illustrated in **Figure 15**. Touch the “normal tissue button” to return to normal tissue sampling.

- During the procedure, the ENCOR ENSPIRE™ Breast Biopsy System has an optional rinse mode when using the vacuum and rinse tubing cassette and saline. This mode allows the operator to apply a saline rinse to samples collected in the sample trap. Initiate the rinse by pressing the “rinse” button on the touch screen.

- If needed, the ENCOR ENSPIRE™ Breast Biopsy System allows the operator to administer additional anesthesia through the ENCOR® Biopsy probe. To enter Anesthetic Mode, touch the “Anesthetic Mode” tab as illustrated in **Figures 15 and 16**. Select the desired anesthetic option to sweep the sample notch through the desired pattern while administering anesthesia. Leave Anesthetic Mode by pressing the “sample” button on either the handswitch or footswitch.

- Following the procedure, the ENCOR ENSPIRE™ Breast Biopsy System allows the operator to deploy a marker through the probe. Enter Marker Mode by pressing the “Marker Mode” tab on the touch screen as illustrated in **Figure 17**. Deploy the marker per the marker Instructions for Use. Rotate the sample notch 180° by pressing the “180°” button on the touch screen as illustrated in **Figure 17**. Exit Marker Mode by pressing the “sample” button on either the handswitch or footswitch.


Figure 15 – The touch screen – Main Screen

Figure 16 – The touch screen – Anesthetic Mode

Marker Mode tab

180° button

After Use

- Remove disposables and discard per your facility's procedures. Clean and disinfect the system hardware per the instructions located under the cleaning and maintenance section.

Note: The system may be left on if additional procedures are to be performed; it is not recommended to leave the main power on overnight. It is not recommended to unplug the drivers after each use. The drivers and foot pedal may be left connected to the ENCOR ENSPIRE™ Breast Biopsy system if all procedures are to be performed with the same drivers.

- To set the system to standby, press and hold the switch on the front of the console (illustrated in **Figure 13**) until the system shuts down. To turn off the main power, press the main switch located on the rear of the console (illustrated in **Figure 12**).


Figure 17 – The touch screen – Marker Mode

I. CLEANING AND MAINTENANCE

To clean the ENCOR ENSPIRE™ Breast Biopsy System, turn off power and disconnect the power cord from the AC outlet. Leave all other cables connected. Clean all exposed surfaces of the ENCOR ENSPIRE™ Breast Biopsy System, cables, and foot pedal with a soft, lint-free cloth dampened with Dispatch or a similar cleaning solution or disinfectant. Follow the procedures approved by your facility or use a validated infection control procedure. Always follow the instructions provided by the manufacturer of the cleaner or disinfectant used.

Do not use abrasive cleaners or spray any fluids directly onto any part of the ENCOR ENSPIRE™ Breast Biopsy System or cables. Spraying the system may cause it to malfunction and will void the warranty.

Do not submerge any component of the ENCOR ENSPIRE™ Breast Biopsy System in fluids. Submerging in fluids may cause the system to malfunction and will void the warranty. Sterilization and exposure to fluids may damage the electrical components of the device. If the system is inappropriately cleaned it may cause the system to malfunction and will void the warranty. Do not autoclave any components of the ENCOR ENSPIRE™ Breast Biopsy System. Do not heat over 54°C (129°F).

Disposal of Equipment

After following the cleaning recommendations above, there are no biohazard risks involved with the disposal of the ENCOR ENSPIRE™ Breast Biopsy System or any of the reusable accessories at the end of the ENCOR ENSPIRE™ Breast Biopsy System's useable life per Waste of Electrical and Electronic Equipment (WEEE) directive [Directive 2002/96/EC]. Dispose of any single use device per the Instructions for Use for the disposable device.

Do not dispose of the ENCOR ENSPIRE™ Breast Biopsy System by placing into trash receptacles.

Routine Maintenance

No calibrations or adjustments are necessary for the ENCOR ENSPIRE™ Breast Biopsy System outside of normal use. In addition, the following routine maintenance activities may be completed:

Electrical safety testing should be performed at intervals no greater than 12 months, using a standard medical safety analyzer.

Regularly inspect all cables and connectors for possible wear or damage.

Odor Control Filter Replacement

Turn the Power Switch to off at the rear of the console. Remove the old Odor Control Filter by rotating it counter-clockwise and discard. Install the new Odor Control Filter (F3000) by rotating clockwise until snug. The location of the Odor Control Filter is identified in **Figure 18**.

Figure 18 – Odor Control Filter

Fuse Replacement

The fuse holder is accessible after the power cord is removed. Unplug the ENCOR ENSPIRE™ Breast Biopsy System from the wall outlet, and remove the power cord from the power inlet connector on the rear panel of the console. Using a small slotted screwdriver, press the release tab located near the center of the fuse holder. Pull the fuse holder assembly outward to remove it. Use Fuse Kit Model Number FK4000. Install the fuse holder, as illustrated in **Figure 19**. Remove and replace any blown fuses in the fuse holder assembly, and snap it back into place. The fuse holder orientation determines the voltage. **Figure 19** illustrates the 115V orientation. Rotate the fuse 180 degrees for 230V orientation. Replace the fuse holder assembly, and reconnect the power cord to the power inlet connector.

Figure 19 – Replacement Fuse Orientation: Orientation for 115V shown

J. TROUBLESHOOTING

The ENCOR ENSPIRE™ Breast Biopsy System is programmed to attempt to determine a reason for any alert condition that occurs. If the system can identify the cause of the alert, the information is displayed on the touch screen.

Depending on the cause of the alert, different actions may be required to clear the alert. Certain alert conditions, when corrected, allow the system to recover automatically, and return to a ready state. Other alerts require pressing a button on the touch screen to clear. If a driver is connected and is the cause of the alert, follow the specific Instructions for Use for that driver to clear the alert.

When no other actions are indicated, clear the alert condition by turning off main power using the switch on the rear of the system (illustrated in **Figure 12**). Wait 20 seconds, and turn the main power back on, and restart the ENCOR ENSPIRE™ Breast Biopsy System with the standby button on the front of the console (illustrated in **Figure 13**).

If any alert condition persists after correcting known problems and clearing the alert condition from the screen, contact Bard for service.

Error Screen Text	Error Number(s)	Troubleshooting
System Error	001 – 002, 005 – 071, 073, 103 – 118, 120 – 121, 123 – 124, 126 – 128, 136, 139	A system error occurred. Contact your Bard representative for assistance.
System Error	003 – 004, 072, 074 – 102, 119, 122, 125, 129, 131 – 135, 137 – 138, 140 – 150, 153 – 163	A system error occurred. Restart the system. If the problem persists, contact your Bard representative.
Failure Reading Handpiece Memory	130	An error with the driver/handpiece occurred. Restart the system. If the problem persists, contact your Bard representative.
System Error	151, 152	A button on the footpedal is stuck. Check the footpedal and restart the system. If problem persists, contact your Bard representative.
Control Module API Initialization Error	1001, 2001	An error occurred during initialization. Restart the system. If the problem persists, contact your Bard representative.
Failed Vacuum – Initialization Error	1002, 2002	An error occurred during initialization. Restart the system. If the problem persists, contact your Bard representative.
Invalid Gearbox I.D./Driver Type	1003	An error with the driver/handpiece occurred. Restart the system. If the problem persists, contact your Bard representative.
Probe Not Removed	1004, 2003	The probe is not removed. Remove the probe to continue.
Driver Button or Contact Stuck during Initialization	1005, 2004	A button on the driver is stuck. Check the driver and restart the system. If problem persists, contact your Bard representative.
Driver Error – Initialization Failed	1006 – 1007, 2005 – 2006	An error with the driver occurred. Disconnect the driver to continue. Connect Driver to retry. If the problem persists, contact your Bard representative.
Driver Error – Control Module Software Failure	1008	The system software is out of date. Contact your Bard representative to update the system software.
Driver Not Connected	1009, 2007	Driver is not connected. Connect Driver to continue.
System Power Down Detected	1010, 2008	The system is powering down. Please wait.
Calibration Failed	1011, 2009 – 2010	Calibration failed. Re-calibrate probe or remove probe to continue.
Probe Error - Wrong Probe inserted	1012	MRI Probe must be used with MRI Driver. Non-MRI Probe must be used with Non-MRI Driver. Remove Probe to continue.
Restoring Vacuum	1013, 2012	Restoring Vacuum. Please wait.
Vacuum Failure Detected	1014, 2013	Vacuum failure detected. Check secure tubing connections. Check for pinched tubing. Ensure Tubing Cartridge is loaded properly. Press Reset to retry.
Probe Error - Hold Current Out of Range	1015	A probe error occurred. Remove probe to continue.
Probe Removed	1016	Probe Removed. Please wait.
Driver Error – Unrecoverable Error during Stroke	1017	Driver Error Occurred - Remove Probe or Disconnect Driver to continue.
Indexing Error	1018	Indexing Error occurred - press reset to continue.
Driver Error – Tool State Integrity Failure	1019	An error with the driver occurred. Disconnect the driver to continue. Connect Driver to retry. If the problem persists, contact your Bard representative.
Vacuum Failure Detected	1020	Vacuum failure detected. Restart System. If problem persists, connect a different Driver.
Error Opening Cutter	1021	An error with the probe occurred. Remove the probe or disconnect the driver to continue.
Cutter Not Centered	1022	An error with the probe occurred. Remove the probe or disconnect the driver to continue.
Probe Unlatched	2011	The probe is unlatched. Insert the probe and press Reset.
Vacuum Failure Detected	2014, 2017	A vacuum failure occurred. Restart the system. If problem persists, contact your Bard representative.
Driver Failure– Sample/Open Failed	2015	An error with the driver occurred. Remove the probe to continue.
Driver Error – Tool State Integrity Failure	2016	An error with the driver occurred. Disconnect the driver to continue.
UI Application Fault Detected	3000, 3001, 3004	An error occurred. Restart the system. If problem persists, contact your Bard representative.
System Communication Fault Detected	3002	A communication error occurred. Restart the system. If problem persists, contact your Bard representative.
Incompatible Driver Script Version	3003	The driver Script must be updated. Follow the onscreen instructions to initiate the update process, or contact Bard representative for assistance.
System detected UI Communications Fault	3005	A communication error occurred. Press Reset to continue. If problem persists, contact your Bard representative.
Driver Script Update Recommended	3006	System detected a driver script update is available. User may update via on-screen prompts, or skip.

K. SPECIFICATIONS

Environmental Conditions	Operate the device in ambient conditions (between 15°C to 25°C, 30% to 75% relative humidity and 700 to 1060 hPa atmospheric pressure). Transport and store in a cool dry place between -20°C to +60°C, 10% to 90% relative humidity and 500 to 1060 hPa atmospheric pressure.
Ingress of Water	Not protected against the ingress of water. IPX0.
Flammability	Equipment not suitable for use in the presence of flammable anesthetic.
Dimensions	16 inches Wide x 26 inches Deep x 46 inches (retracted) to 57 inches (extended) High (41 cm x 66 cm x 117 cm retracted, 145 cm extended)
Weight	105 lb (48 Kg)

SPECIFICATIONS	
System Power Ratings	Maximum: 100-120 VAC, 600VA, 50/60 Hz, 3-wire, grounded with removable power cord 220-240 VAC, 600VA, 50/60 Hz, 3-wire, grounded with removable power cord
Electrical Conformity	This medical equipment has passed all required testing for electric shock, fire and mechanical hazards in accordance with UL60601-1, IEC/EN 60601-1, CAN/CSA C22.2 No 601-1.
Classification	Class I, Type BF equipment.

Nominal Ratings	Standby
Voltage/Frequency	115 VAC ~ 50/60 Hz 230 VAC ~ 50/60 Hz
Current	4
VA (Calculated)	600
Watts	450
Power factor	0.75

Nominal Ratings	Operating
Voltage/Frequency	115 VAC ~ 50/60 Hz 230 VAC ~ 50/60 Hz
Current	4
VA (Calculated)	600
Watts	450
Power factor	0.75

Electromagnetic Emissions and Immunity

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided below.

Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment.

Guidance and Manufacturer's Declaration – Emissions All Equipment and Systems

Guidance and Manufacturer's Declaration - Emissions		
The ENCOR ENSPIRE™ Breast Biopsy System is intended for use in the electromagnetic environment specified below. The customer or user of the ENCOR ENSPIRE™ Breast Biopsy System should ensure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 1, Class A	The ENCOR ENSPIRE™ Breast Biopsy System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Harmonics IEC 61000-3-2	Class A - Complies	The ENCOR ENSPIRE™ Breast Biopsy System is suitable for use in all establishments, other than domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Flicker IEC 61000-3-3	Complies	

Guidance and Manufacturer's Declaration – Immunity All Equipment and Systems

Guidance and Manufacturer's Declaration – Immunity			
The ENCOR ENSPIRE™ Breast Biopsy System is intended for use in the electromagnetic environment specified below. The customer or user of the ENCOR ENSPIRE™ Breast Biopsy System should ensure that it is used in such an environment.			
Immunity Test	EN/IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
ESD EN/IEC 61000-4-2	±2kV, ±4kV, ±6kV Contact discharge ±2kV, ±4kV, ±6kV, ±8kV Air discharge	Pass	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the r/h should be at least 30%.
EFT EN/IEC 61000-4-4	±2kV Mains ±1kV I/Os	Pass	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN/IEC 61000-4-5	±1kV Differential ±2kV Common	Pass	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips/Dropout EN/IEC 61000-4-11	>95% Dip for 0.5 Cycle	Pass	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ENCOR ENSPIRE™ Breast Biopsy System requires continued operation during power mains interruptions, it is recommended that the ENCOR ENSPIRE™ Breast Biopsy System be powered from an uninterruptible power supply or battery.
	60% Dip for 5 Cycles	Pass	
	30% Dip for 25 Cycles	Pass	
	>95% Dip for 5 Seconds	Pass	
Power Frequency 50/60Hz Magnetic Field EN/IEC 61000-4-8	3A/m	Pass	Power frequency magnetic fields should be that of a typical commercial or hospital environment.

Guidance and Manufacturer's Declaration – Emissions Equipment and Systems that are NOT Life-supporting

Guidance and Manufacturer's Declaration – Emissions
The ENCOR ENSPIRE™ Breast Biopsy System is intended for use in the electromagnetic environment specified below. The customer or user of the ENCOR ENSPIRE™ Breast Biopsy System should ensure that it is used in such an environment.

Immunity Test	EN/IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted RF EN/IEC 61000-4-6 Radiated RF EN/IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	Pass Pass	Portable and mobile communications equipment should be separated from the ENCOR ENSPIRE™ Breast Biopsy System by no less than the distances calculated/listed below: D=(3.5/V1)(Sqrt P) D=(3.5/E1)(Sqrt P) 80 to 800 MHz D=(7/E1)(Sqrt P) 800 MHz to 2.5 GHz where P is the max power in watts and D is the recommended separation distance in meters. Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (V1 and E1). Interference may occur in the vicinity of equipment containing a transmitter.

Recommended Separation Distances between portable and mobile RF Communications equipment and the ENCOR ENSPIRE™ Breast Biopsy System Equipment and Systems that are NOT Life-supporting

Recommended Separations Distances for the ENCOR ENSPIRE™ Breast Biopsy System
The ENCOR ENSPIRE™ Breast Biopsy System is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of the ENCOR ENSPIRE™ Breast Biopsy System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications Equipment and the ENCOR ENSPIRE™ Breast Biopsy System as recommended below, according to the maximum output power of the communications equipment.

Max Output Power (Watts)	Separation (m) 150kHz to 80MHz D=(3.5/V1)(Sqrt P)	Separation (m) 80 to 800MHz D=(3.5/E1)(Sqrt P)	Separation (m) 800MHz to 2.5GHz D=(7/E1)(Sqrt P)
0.01	.1166	.1166	.2333
0.1	.3689	.3689	.7378
1	1.1666	1.1666	2.3333
10	3.6893	3.6893	7.3786
100	11.6666	11.6666	23.3333

Cable Length

Cable	Length (feet)
AC Power Cord (PC3100/PC3200/PC3500/PC3600)	12 ft (US) / 2.5m (International)

Equipotential Post:

The purpose of the Industry Standard Equipotential Post on the rear of the ENCOR ENSPIRE™ Breast Biopsy System is to provide a connection point to reduce the possibility of voltage potentials between touchable conductive parts of all the different Medical Equipment Systems within the Patient and Operator area.

To connect the ENCOR ENSPIRE™ Breast Biopsy System to an Equipotential Network, snap the Equipotential Network Socket over the Equipotential Post on the rear of the ENCOR ENSPIRE™ Breast Biopsy System base unit. To remove The ENCOR ENSPIRE™ Breast Biopsy System from the Equipotential Network, pull the equipotential Network Socket off the Equipotential Post.

L. HOW SUPPLIED

The ENCOR ENSPIRE™ Breast Biopsy System is supplied non-sterile and should be cleaned prior to first use per the Cleaning and Maintenance instructions.

- ENCOR®, ENCOR® MRI and ENCOR® 360 drivers are sold separately. The drivers are supplied non-sterile. Reference the Instructions for Use included with each driver for Cleaning and Maintenance instructions.
- The Vacuum Canister, Vacuum Tubing Cassette, Vacuum and Rinse Tubing Cassette and ENCOR®, ENCOR® MRI and ENCOR® 360 probes are sold separately. The vacuum canister, probes and tubing cassettes are supplied sterile for single patient use only.

M. REFERENCES

PK1237100 (ENCOR ENSPIRE™ Tubing Cassette IFU), IU0055 (ENCOR® 360 driver and probe IFU), IU0072 (ENCOR® driver and probe IFU) and IU0073 (ENCOR® MRI driver and probe IFU)

N. WARRANTY

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

Extended warranty service agreements are available. Consult a Bard representative for details regarding terms and conditions.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

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**Quality Assurance Seal at Intertek**

Label d'assurance qualité chez Intertek
QS-Siegel von Intertek
Sigillo di assicurazione della qualità a Intertek
Sello de garantía de calidad Intertek
Zegel van kwaliteitswaarborging bij Intertek
Selo de Garantia de Qualidade Intertek
Σφραγίδα διασφάλισης ποιότητας στην Intertek
Kvalitetssikringsforsegling ved Intertek
Bekräftad kvalitetssäkring på Intertek
Intertekin laatumerkintä
Kvalitetssikringsmerke i Intertek
Pieczęć kontroli jakości Intertek
Minőségbiztosítás az Interteknél
Pečeť jakosti Intertek
Intertek Kalite Güvence Damgası
Intertek 品質保證封條
Intertek의 품질 보증 인종
Печать контроля качества Intertek

**Date Of Manufacture**

Date de fabrication
Herstellungsdatum
Data di produzione
Fecha de fabricación
Productiedatum
Data de fabrico
Ημερομηνία κατασκευής
Fremstillingsdato
Tillverkningsdatum
Valmistuspäivämäärä
Produksjonsdato
Data produkcji
Gyártás dátuma
Datum výroby
Üretim Tarihi
製造日期
제조 연월일
Дата производства

SN**Serial Number**

Numéro de série
Seriennummer
Numero di serie
Número de serie
Seriennummer
Número de série
Σειριακός αριθμός
Seriennummer
Seriennummer
Sarjanumero
Seriennummer
Numer seryjny
Sorozatszám
Sériové číslo
Ser Numarası
序號
일련 번호
Серийный номер

**Non-Sterile**

Non stérile
Nicht steril
Non sterile
Sin esterilizar
Niet steriel
Não esterilizado
Μη στέρπο
Usteril
Icke-steril
Steriloimaton
Ikke-steril
Niesterylno
Nem steril
Nesterilni
Steril Degildir
非无菌
비살균
Нестерильно

**Do Not Re-Sterilize**

Ne pas restériliser
Nicht resterilisieren
Non risterilizzare
No reesterilizar
Niet opnieuw steriliseren
Não reesterilizar
Μην επαναστεριώνετε
Må ikke resteriliseres
Får ej resteriliseras
Ei saa steriloida uudelleen
Må ikke resteriliseres
Nie steryliżować powtórnie
Ne sterilizálja újra
Neprovádajte resterilizaci
Tekrar Sterilize Etmeyin
請勿重新消毒
다시 멸균하지 마십시오
Повторная стерилизация запрещается

**Do Not Immerse Device In Fluids**

Ne pas plonger le dispositif dans des liquides
Produkt nicht in Flüssigkeiten eintauchen
Non immergere il dispositivo in alcun fluido
No sumergir el dispositivo en líquidos
Hulpmiddel niet in vloeistoffen onderdempelen
Não mergulhar o dispositivo em líquidos
Μην εμβυθίζετε τη συσκευή σε υγρά
Anordningen må ikke nedsænkes i væsker
Enheten får inte nedsänkas i vätska
Ei saa upottaa nesteisiin
Utstyret må ikke legges i væske
Nie zanurzać urządzenia w płynach
Az eszköz folyadékba merítése tilos
Neponojujte zařízení do kapalín
Cihazı Sıvılara Batırmayın
請勿將裝置浸泡在液體中
장치를 액체에 담그지 마십시오
Запрещается погружать устройство в жидкости

**Wipe With Damp Cloth**

Frotter avec un chiffon humide
Mit feuchtem Tuch abwischen
Strofinare con un panno umido
Limpiar con un paño húmedo
Afneamen met vochtige doek
Limpar com um pano húmido
Σκουπίστε με υγρό τσάβι
Åttöras med en fuktig klud
Torkas av med en fuktig duk
Pyyhi kostealla kangaspyyhkeellä
Tørk av med en fuktig klut
Wycierać wilgotną szmatką
Nedves ruhával törölnöd
Otrějte vlhkou utěrkou
Nemli Bezle Sliin
以濕布擦拭
젖은 천으로 닦으십시오
Протирать влажной тканью

**Do Not Spray Fluids Directly Into The Device**

Ne pas vaporiser de liquides directement dans le dispositif
Flüssigkeiten nicht direkt in das Produkt sprühen
Non spruzzare fluidi direttamente nel dispositivo
No rociar líquidos directamente en el dispositivo
Vloeistoffen niet rechtstreeks in het hulpmiddel spuiten
Não pulverizar líquidos directamente no dispositivo
Μην ψεκάζετε υγρά απευθείας επάνω στη συσκευή
Der må ikke sprøjtes væsker direkte ind i anordningen
Vätska får inte sprayas direkt in i enheten
Älä suihkuta laitteeseen nesteitä
Ikke spray væske direkte inn i utstyret
Nie rozpryskiwać płynów bezpośrednio na urządzenie
Folyadékot közvetlenül az eszközbe permetezni tilos
Nestříkejte kapaliny přímo na zařízení
Cihazı Doğrudan Sıvı Püskürtmeyin
請勿直接將液體噴在裝置上
장치에 액체를 직접 분사하지 마십시오
Запрещается опрыскивать устройство какими-либо жидкостями

**Fuse Inside.**


Fusible interne.
Fusible interne.
Mit Sicherung
Parte interna del fusibile.
Fusible interno.
Bevat zekering.
Fusivel no interior.
Περιέχει ασφάλεια.
Sikring inden i.
Säkning inuti.
Sisällä on varoke.
Innendig sikring.
Bezpieczniki wewnątrz.
A készülék belsejében biztosítékok találhatók.
Pojistka uvnitř.
Sigorta içindedir.
內部保險絲。
내부 퓨즈.
Установлен предохранитель.

**Authorised Representative in the European Community**

Représentant autorisé au sein de la Communauté européenne
Berechtigter Vertreter in der Europäischen Gemeinschaft
Rappresentante autorizzato nella Comunità Europea
Representante autorizado en la Comunidad Europea
Gemachtigde vertegenwoordiger binnen de Europese Gemeenschap
Representante autorizado na Comunidade Europeia
Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα
Autoriseret repræsentant i Det Europæiske Fællesskab
Auktoriserad representant i EG
Valtuutettu edustaja Euroopan yhteisössä
Autoriseret representant i EU
Autorizovaný predstaviteľ na terenie Unii Európskejkej
Bejegyztett képviselő az Európai Közösségben
Autorizovaný zástupce pro Evropskou unii
Avrupa Topluluğu Yetkili Temsilcisi
歐洲共同體的授權代表
유럽 공동체 내의 인증받은 대리업체
Уполномоченный представитель в Европейском сообществе

**Do Not Use if the Product Sterilization Barrier or its Packaging is Compromised**

Ne pas utiliser si la barrière de stérilisation ou l'emballage du produit sont endommagés
Bei beeinträchtigtem Sterilschutz oder beschädigter Verpackung das Produkt nicht verwenden
Non utilizzzare se la barriera di sterilizzazione del prodotto o la confezione sono compromessi
No utilizar si la barrera de esterilización del producto o su envase están dañados
Niet gebruiken wanneer de sterilisatiebarrière van het product of de verpakking is aangetast
Não utilizar se a barreira de esterilização do produto ou respectiva embalagem estiverem comprometidas
Μη χρησιμοποιείτε το προϊόν εάν έχει παραβιαστεί ο φραγμός αποστείρωσης ή η συσκευασία του
Må ikke anvendes, hvis produktets steriliseringsbarriere eller emballagen er kompromitteret
Använd inte om produktens sterila barriär eller dess förpackning är skadad
Ei saa käyttää, jos tuotteen steriloitinsuojaus tai pakkaus on vaurioitunut
Bruk ikke produktet hvis den sterile barrieren eller emballasjen er brutt
Nie używać, jeśli naruszono sterylność produktu lub jego opakowanie
Ne használja, ha a termékét védő steril zár vagy a csomagolás sérült
Pokud je narušena sterilní ochrana nebo poškozen obal, výrobek nepoužívejte
Úrün Steril Bariyeri veya Ambalajı zarar görmüşse kullanmayın
如果產品消毒屏障或包裝受損，請勿使用
제품의 멸균 장벽이나 포장에 손상이 된 경우 사용하지 마십시오
Запрещается применять изделие, если стерильная упаковка или внешняя упаковка повреждена

**Atmospheric Pressure Limitation**

Limite de pression atmosphérique
Luftdruckbereich
Limite di pressione atmosferica
Limitación de presión atmosférica
Limiet atmosferische druk
Limite de pressão atmosférica
Περιορισμός ατμοσφαιρικής πίεσης
Atmosfærisk trykbeogrænsning
Lufttrycksbeogrænsning
Ilmanpaineen rajoitus
Grense for atmosfærisk trykk
Ograniczenia dotyczące ciśnienia atmosferycznego
Légköri nyomás-korlátozás
Rozmezi atmosférického tlaku
Atmosfer Basinci Sınırlaması
大氣壓力限制
대기압 제한
Ограничение атмосферного давления

**Consult Instructions for Use**

Consulter le mode d'emploi
Gebrauchsanweisung beachten
Leggere le istruzioni per l'uso
Consulte las Instrucciones de uso
Raadpleeg de gebruiksaanwijzing
Consultar as instruções de utilização
Συμβουλευτείτε τις οδηγίες χρήσης
Se bruksanvisningen
Se bruksanvisning
Lue käyttöohjeet
Se bruksanvisning
Zapoznać się z instrukcją obsługi
Lásd a használati útmutatót
Řiďte se návodem k použití
Kullanma Talimatlarına Başvurun
請參考使用說明
사용 지침 참조
Обратиться к инструкциям по применению

**Humidity Limitation**

Limite d'humidité
Feuchtigkeitsbereich
Limite di umidità
Limitación de humedad
Vochtigheidslimiet
Limite de humidade
Περιορισμός υγρασίας
Fugtighedsbegrænsning
Luftfuktighetsbegrænsning
Ilmankosteuden rajoitus
Fuktighetsgrense
Ograniczenia dotyczące wilgotności
Páratartalom-korlátozás
Rozmezi vlhkosti
Nemlilik Sınırlaması
濕度限制
습도 제한
Ограничение влажности

**Serial Port**

Port série
Serieller Anschluss
Porta seriale
Puerto serie
Seriële poort
Porta de série
Σειριακή θύρα
Seriel port
Serieport
Sarjalitántá
Serieport
Port szeregowy
Soros port
Sériový port
Seri Bağlantı Noktası
序列連接埠
직렬 포트
Последовательный порт

**Temperature Limitation**

Limite de température
Temperaturbereich
Limite di temperatura
Limitación de temperatura
Temperatuurlimiet
Limite de temperatura
Περιορισμός θερμοκρασίας
Temperaturbegrænsning
Temperaturbegrænsning
Lämpötilan rajoitus
Temperaturgrense
Ograniczenia dotyczące temperatury
Hőmérséklet-korlátozás
Rozmezi teplot
Sıcaklık Sınırlaması
溫度限制
온도 제한
Ограничение температуры

**VGA Port**

Port VGA
VGA-Anschluss
Porta VGA
Puerto VGA
VGA-poort
Porta VGA
Θύρα VGA
VGA-port
VGA-port
VGA-liitäntä
VGA-port
Port VGA
VGA port
Připojka VGA
VGA Bağlantı Noktası
VGA 連接埠
VGA 포트
Port VGA

**Dangerous Voltage**

Tension dangereuse
Gefährliche Spannung
Tensione pericolosa
Tensión peligrosa
Gevaarlijke spanning
Tensão perigosa
Επικίνδυνη τάση
Farlig spænding
Farlig spänning
Vaarallinen jännite
Farlig spenning
Niebezpieczne napięcie
Veszélyes feszültség
Nebezpečné napětí
Tehlikeli Gerilim
危險電壓
위험 전압
Опасное напряжение

**Equipotential**

Équipotentialité
Spannungsausgleich
Equipotenziale
Equipotencial
Equipotentiaal
Equipotential
Ισοδυναμικό
Ækviopotientiel
Ekvipotential
Tasapotentiaalinen
Ekvipotensial
Ekvipotencialność
Ekvipotenciális
Ekvipotenciální hladina
Espotansiyel
等電位
등전위
Эквипотенциальный

**USB Port**

Port USB
USB-Anschluss
Porta USB
Puerto USB
USB-poort
Porta USB
Θύρα USB
USB-port
USB-port
USB-port
USB-liitäntä
USB-port
Port USB
USB-port
Připojka USB
USB Bağlantı Noktası
USB 連接埠
USB 포트
Port USB

**Ethernet Port**


Port Ethernet
Ethernet-Anschluss
Porta Ethernet
Puerto Ethernet
Ethernetpoort
Porta Ethernet
Θύρα Ethernet
Ethernet-port
Ethernet-port
Ethernet-liitäntä
Ethernet-port
Port Ethernet
Ethernet port
Připojka ethernet
Ethernet Bağlantı Noktası
乙太網路連接埠
이더넷 포트
Port Ethernet

**Foot pedal**

Pédale de commande au pied
Fußpedal
Pedaliera
Pedal
Voetpedaal
Pedal
Ποδομοχλός
Fodpedal
Fotpedal
Poljin
Pedal
Pedal nożny
Lábpedál
Pedál
Ayak pedali
踏板
풋 페달
Ножная педаль

**AC Power Port (Per PK1238200)**

Port d'alimentation secteur (conformément à PK1238200)
Netzanschluss (PK1238200)
Porta di alimentazione CA (Per PK1238200)
Puerto de CA (PK1238200)
AC-voedingspoort (conform PK1238200)
Porta de alimentação CA (de acordo com PK1238200)
Θύρα εναλλασσόμενου ρεύματος (σύμφωνα με PK1238200)
Vekselstrømsport (iht. PK1238600)
Strömport (enligt PK1238200)
Verkkovirtaliitäntä (PK1238200)
Strömport (iht. PK1238200)
Port zasilania AC (na PK1238200)
AC Port (PK1238200 szerint)
Připojka napájení stř. proudem (PK1238200)
AC Güç Bağlantı Noktası (PK1238200'a göre)
交流電連接埠 (Per PK1238200)
AC 전원 포트(PK1238200에 따름)
Разъем питания переменного тока (для PK1238200)

**AC Power Port (Per PK1238600)**

Port d'alimentation secteur (conformément à PK1238600)
Netzanschluss (PK1238600)
Porta di alimentazione CA (Per PK1238600)
Puerto de CA (PK1238600)
AC-voedingspoort (conform PK1238600)
Porta de alimentação CA (de acordo com PK1238600)
Θύρα εναλλασσόμενου ρεύματος (σύμφωνα με PK1238600)
Vekselstrømsport (iht. PK1238600)
Strömport (enligt PK1238600)
Verkkovirtaliitäntä (PK1238600)
Strömport (iht. PK1238600)
Port zasilania AC (na PK1238600)
AC Port (PK1238600 szerint)
Připojka napájení stř. proudem (PK1238600)
AC Güç Bağlantı Noktası (PK1238600'a göre)
交流電連接埠 (Per PK1238600)
AC 전원 포트(PK1238600에 따름)
Разъем питания переменного тока (для PK1238600)

Copyright © 2013 C.R. Bard, Inc. All rights reserved.

Copyright © 2013 C.R. Bard, Inc. Tous droits réservés.
 Copyright © 2013 C.R. Bard, Inc. Alle Rechte vorbehalten.
 Copyright © 2013 C.R. Bard, Inc. Tutti i diritti riservati.
 Copyright © 2013 C.R. Bard, Inc. Todos los derechos reservados.
 Copyright © 2013 C.R. Bard, Inc. Alle rechten voorbehouden.
 Copyright © 2013 C.R. Bard, Inc. Todos os direitos reservados.
 Πνευματικά δικαιώματα © 2013 C.R. Bard, Inc. Με την επιφύλαξη παντός δικαιώματος.
 Copyright © 2013 C.R. Bard, Inc. Alle rettigheder forbeholdt.
 Copyright © 2013 C.R. Bard, Inc. Eftertryck förbjudes.
 Copyright © 2013 C.R. Bard, Inc. Kaikki oikeudet pidätetään.
 Copyright © 2013 C.R. Bard, Inc. Med enerrett.
 Copyright © 2013 C.R. Bard, Inc. Wszelkie prawa zastrzeżone.
 Copyright © 2013 C.R. Bard, Inc. Minden jog fenntartva.
 Copyright © 2013 C.R. Bard, Inc. Všechna práva vyhrazena.
 Telif Hakkı © 2013 C.R. Bard, Inc. Tüm hakları saklıdır.
 Copyright © 2013 C.R. Bard, Inc. 保留所有權利。
 저작권 © 2013 C.R. Bard, Inc. All rights reserved.
 Авторские права © 2013 Компания C. R. Bard, Inc. Все права защищены.

Only For Use Indoors/Do Not Use Outdoors

Réservé à une utilisation en intérieur/Ne pas utiliser à l'extérieur
 Nur zur Verwendung in Innenräumen/Nicht im Freien verwenden
 Solo per uso al coperto/non usare all'aperto
 Solo para uso en interiores/No utilizar en exteriores
 Uitsluitend voor gebruik binnenshuis/niet buiten gebruiken
 Apenas para utilização em interiores/Não utilizar em exteriores
 Αποκλειστικά για χρήση σε εσωτερικό χώρο/Μην το χρησιμοποιείτε σε εξωτερικό χώρο
 Kun til indendørsbrug/Må ikke bruges udendørs
 Endast för inomhusbruk/Får ej användas utomhus
 Saa käyttää vain sisätiloissa / Ei saa käyttää ulkona
 Kun til innendørs bruk / Må ikke brukes utendørs
 Do stosowania wyłącznie wewnątrz pomieszczeń / Nie stosować na zewnątrz
 Kizárólag beltéri használatra / Szabadban nem használható
 Pouze pro použití v uzavřených prostorech, nikoliv pro venkovní použití
 Sadece İç Mekanda Kullanılır/Dış Mekanda Kullanılmayın
 仅限室內使用 / 請勿在戶外使用
 실내에서만 사용/야외에서 사용하지 말 것
 Только для эксплуатации в помещении/эксплуатация вне помещения запрещена

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician

MISE EN GARDE : Selon la loi fédérale américaine, ce dispositif ne peut être vendu que par un médecin ou sur prescription médicale.
 ACHTUNG: Gemäß der Bundesgesetzgebung der USA darf dieses Produkt nur von einem Arzt bzw. auf ärztliche Anordnung verkauft werden.
 ATTENZIONE! La legge federale statunitense limita la vendita di questo dispositivo ai medici o su presentazione di prescrizione medica.
 ATENCIÓN: Las leyes federales (EE. UU.) limitan la venta de este dispositivo a médicos o bajo prescripción facultativa.
 LET OP: Volgens de Amerikaanse wetgeving mag dit hulpmiddel uitsluitend door of op voorschrift van een arts worden gekocht.
 ATENÇÃO: A lei federal (EUA) limita a venda deste dispositivo a um médico ou mediante receita médica.
 ΠΡΟΣΟΧΗ: Η νοσοπονοδιακή νομοθεσία (των Η.Π.Α.) περιορίζει την πώληση της συσκευής αυτής σε ιατρούς ή κατόπιν εντολής ιατρού.
 FORSGTIG: Amerikansk lovgivning begrænser denne anordning til salg af eller på ordning af en læge.
 OBS! Enligt federal (USA) lag får denna produkt endast säljas av eller på ordination av läkare.
 VAROITUS: Yhdysvaltain lain mukaan tämän laitteen saa myydä ainoastaan lääkäri tai lääkäarin määrääksestä.
 FORSIKTIG! Ifølge amerikansk federal lovgivning kan dette utstyret selges bare av lege eller på resept.
 PRZESTROGA: Na mocy prawa federalnego (USA) sprzedaż niniejszego urządzenia może być prowadzona wyłącznie przez lekarza lub na jego zlecenie.
 FIGYELEM: Az Egyesült Államok szövetségi törvénye alapján az eszköz értékesítése csak orvos által vagy orvosi utasításra történhet.
 UPOZORNĚNÍ: Federální zákony USA povolují prodej tohoto přístroje pouze lékařům nebo na lékařský předpis.
 DIKKAT: A.B.D. federal kanunlarına göre bu cihazın satışı yalnızca bir doktor tarafından ya da doktor talimatıyla yapılabilir.
 警告：美國聯邦法律規定本裝置僅能由醫師銷售或憑醫師囑銷售。
 주의: 미연방법은 이 기구를 의사에게 직접 또는 의사가 주문한 경우에만 판매하도록 제한합니다.
 ПРЕДУПРЕЖДЕНИЕ: В соответствии с федеральным законодательством (США) это устройство может быть продано только врачом или по его указанию.

Do not dispose of the EnCor EnSPIRE™ Breast Biopsy System by placing into trash receptacles

Ne pas jeter le système à biopsie mammaire EnCOR EnSPIRE™ en le plaçant dans des récipients à ordures
 Das EnCOR EnSPIRE™-Brustbiopsie-System NICHT im normalen Müll entsorgen
 Non smaltire il sistema per biopsia della mammella EnCOR EnSPIRE™ in recipienti per rifiuti
 No deseche el sistema de biopsia mamaria EnCOR EnSPIRE™ en un contenedor de basura
 Het EnCOR EnSPIRE™ borstbiopsiesysteem niet weggoeien bij het normale huisvuil
 Não eliminar o Sistema de Biopsia Mamária EnCOR EnSPIRE™ colocando-o em caixotes destinados ao lixo comum
 Μην απορρίπτετε το σύστημα βιοψίας μαστού EnCOR EnSPIRE™ σε κάδους απορριμμάτων
 EnCOR EnSPIRE™ brystbiopsisystem må ikke bortskaffes ved anbringelse i affaldsbeholdere
 Kasserer inte EnCOR EnSPIRE™ bröstbiopsisystem genom att kasta det i en soptunna
 EnCOR EnSPIRE™ -rintabiopsiajärjestelmä ei saa hävittää tavonomaisten jätteen mukana
 Ikke kasser EnCOR EnSPIRE™-brystbiopsisystemet ved å legge det i søppelcontainere
 Nie należy utylizować systemu do biopsji piersi EnCOR EnSPIRE™, wyrzucając go do pojemnika z odpadami
 EnCOR EnSPIRE™ emlőbiopsziás rendszeret ne helyezze hulladékgyűjtőbe
 Systém pro biopsii prsní tkáně EnCOR EnSPIRE™ nesmí být likvidován jako komunální odpad
 EnCOR EnSPIRE™ Meme Biyopsi Sistemini çöp kutularına atmayın
 棄置 EnCOR EnSPIRE™ 乳房切片系統時，請勿放在垃圾桶
 EnCOR EnSPIRE™ 유방 생검 시스템을 쓰레기통에 버리지 마십시오
 Запрещается выбрасывать систему для биопсии молочной железы EnCOR EnSPIRE™ в мусорные контейнеры

Attention: Ground Reliability can only be achieved when equipment is connected to a receptacle marked "HOSPITAL ONLY" or "HOSPITAL GRADE". For Canada and USA only

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Bard, SenoRx, EnCor and EnCor EnSPIRE are trademarks and/or registered trademarks of C.R. Bard, Inc. or an Affiliate

Bard, SenoRx, EnCor et EnCor EnSPIRE sont des marques commerciales et/ou marques déposées de C.R. Bard, Inc. ou d'une de ses filiales
 Bard, SenoRx, EnCor und EnCor EnSPIRE sind Marken und/oder eingetragene Marken von C.R. Bard, Inc. oder einer Tochtergesellschaft
 Bard, SenoRx, EnCor e EnCor EnSPIRE sono marchi e/o marchi registrati di C.R. Bard, Inc. o di una sua affiliata
 Bard, SenoRx, EnCor y EnCor EnSPIRE son marcas comerciales o marcas comerciales registradas de C.R. Bard, Inc. o de una filial
 Bard, SenoRx, EnCor en EnCor EnSPIRE zijn handelsmerken en/of gedeponeerde handelsmerken van C. R. Bard, Inc. of een geaffilieerd bedrijf
 Bard, SenoRx, EnCor e EnCor EnSPIRE são marcas comerciais e/ou marcas comerciais registradas da C.R. Bard, Inc. ou de uma filial
 Τα Bard, SenoRx, EnCor και EnCor EnSPIRE αποτελούν εμπορικά σήματα και/ή σήματα κατατεθέντα της C.R. Bard, Inc. ή κάποιας συνδεδεμένης με αυτήν εταιρίας
 Bard, SenoRx, EnCor og EnCor EnSPIRE er varemærker og/eller registrerede varemærker tilhørende C.R. Bard, Inc. eller et tilknyttet selskab
 Bard, SenoRx, EnCor and EnCor EnSPIRE är varumärken och/eller inregistrerade varumärken som tillhör C.R. Bard, Inc. eller en till företagert anknuten tillverkare
 Bard, SenoRx, EnCor ja EnCor EnSPIRE ovat C.R. Bard, Inc:n tai sen tytäryhtiön tavaramerkkejä tai rekisteröityjä tavaramerkkejä
 Bard, SenoRx, EnCor og EnCor EnSPIRE er varemærker og/eller registrerte varemærker som tilhører C.R. Bard, Inc. eller et tilknyttet selskap
 Bard, SenoRx, EnCor i EnCor EnSPIRE są znakami towarowymi i/lub zarejestrowanymi znakami towarowymi firmy of C.R. Bard, Inc. lub spółek zależnych
 A Bard, SenoRx, EnCor és EnCor EnSPIRE a C.R. Bard, Inc. vagy valamilyen leányvállalata védjegyei és/vagy bejegyzett védjegyei
 Bard, SenoRx, EnCor a EnCor EnSPIRE jsou zapsané ochranné známky společnosti C.R. Bard, Inc. nebo její pobočky
 Bard, SenoRx, EnCor ve EnCor EnSPIRE, C.R. Bard, Inc. veya bir Bağlı Kuruluşu ait ticari markalar ve/veya tescilli ticari markalardır
 Bard · SenoRx · EnCor 和 EnCor EnSPIRE 是 C.R. Bard, Inc. 或其附屬機構的商標和/或註冊商標
 Bard, SenoRx, EnCor 및 EnCor EnSPIRE는 C.R. Bard, Inc. 또는 계열사의 상표 또는 등록 상표입니다
 Bard, SenoRx, EnCor и EnCor EnSPIRE — товарные знаки и/или зарегистрированные товарные знаки компании C. R. Bard, Inc. или ее дочерних компаний

Manufacturer:

SenoRx, Inc.

1625 West 3rd Street

Tempe, AZ 85281

USA

TEL: 1-480-894-9515

1-800-321-4254

FAX: 1-480-966-7062

1-800-440-5376

www.bardbiopsy.com

**Authorised Representative in
the European Community**

Bard Limited

Forest House

Tilgate Forest Business Park

Brighton Road, Crawley

West Sussex

RH11 9BP, UK

BARID

BIOPSY SYSTEMS