## EnVision™ FLEX, High pH (Dako Omnis)

Code GV800

2nd edition/ 2ème édition / 2. Auflage

The kit contains reagents sufficient for 600 tests. For Dako Omnis instruments.

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Page/ Page/ Seite

### ENGLISH

Intended Use .....	4
Summary and Explanation.....	4
Reagents .....	5
A. Materials provided .....	5
B. Materials required but not provided.....	5
C. Optional reagents .....	6
Precautions .....	7
Storage.....	8
Reagent Preparation.....	9
Specimen Collection and Preparation .....	9
Procedure.....	10
Staining procedure.....	10
Quality Control.....	11
Interpretation of Results.....	11
Limitations .....	12
Explanation of symbols.....	32

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### **Intended Use**

For in vitro diagnostic use.

EnVision™ FLEX, High pH (Dako Omnis) visualization system is intended for use in immunohistochemistry together with Dako Omnis instruments. The system detects primary mouse and rabbit antibodies and the reaction is visualized by EnVision™ FLEX DAB+ Chromogen. If used with EnVision™ FLEX+ Rabbit LINKER (Code GV809) or EnVision™ FLEX+ Mouse LINKER (Code GV821) signal amplification of primary rabbit antibodies or primary mouse antibodies, respectively, can be achieved. The EnVision™ FLEX reagents are intended for use on formalin-fixed, paraffin-embedded tissue sections.

### **Summary and Explanation**

EnVision™ FLEX visualization system is a flexible system giving optimal staining on Dako Omnis instruments, when using the protocol recommended in this package insert.

Prior to staining formalin-fixed, paraffin-embedded tissue sections are subjected to onboard deparaffinization and hydration followed by heat-induced epitope retrieval (HIER) using the target retrieval reagent specified in the package insert for the primary antibody. This kit includes EnVision™ FLEX Target Retrieval Solution, High pH (50x) (DM848). Alternatively, EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code GV804) or EnVision™ FLEX Target Retrieval Solution, Low pH (50x) (Code GV805) is recommended. Some primary antibodies require enzymatic pre-treatment of tissue for optimal staining instead of HIER.

Endogenous peroxidase should be blocked with EnVision™ FLEX Peroxidase-Blocking Reagent (DM841) included in the kit. Due to an effective washing procedure and the presence of stabilizing proteins in the Dako reagents, extra blocking steps to reduce non-specific background staining are unnecessary.

Wash Buffer (20x) (GC807) available as optional reagent is recommended.

Primary antibodies are not provided with the kit. We recommend the use of Dako FLEX Ready-to-Use Primary Antibodies or Dako concentrated Primary Antibodies (see respective package inserts for guidelines).

EnVision™ FLEX+ Rabbit LINKER (Code GV809) may be applied for an optional signal amplification of rabbit primary antibodies. For signal amplification when using mouse primary antibodies, EnVision™ FLEX+ Mouse LINKER (Code GV821) may be applied. As a guideline: EnVision™ FLEX+ Rabbit LINKER provides 2-3 fold signal amplification and EnVision™ FLEX+ Mouse LINKER provides 4-5 fold signal amplification.

Dako EnVision™ FLEX /HRP reagent (DM842) in the kit consists of a dextran backbone to which a large number of horseradish peroxidase (HRP) molecules and secondary antibody molecules have been coupled. A unique chemistry is used for the coupling reaction, which permits the binding of up to 100 HRP molecules and up to 20 antibody molecules per backbone.

The substrate system in the kit consists of two components: EnVision™ FLEX DAB+ Chromogen (DM847), a concentrated diaminobenzidine (DAB) solution, and EnVision™ FLEX Substrate Buffer (DM843) containing hydrogen peroxide. Before use EnVision™ FLEX DAB+ Chromogen must be diluted in EnVision™ FLEX Substrate Buffer. The substrate system produces a crisp brown end product at the site of the target antigen. Timely mixing of chromogen and buffer is automatically performed onboard Dako Omnis.

Hematoxylin (Code GC808) is recommended for counterstaining. The reagent provides a clear blue, nuclear staining.

The stained tissue sections should be mounted with organic-solvent-based mounting medium.

## Reagents

### A. Materials provided

DM841

EnVision™ FLEX  
PEROXIDASE-BLOCKING  
REAGENT  
(Dako Omnis)

**EnVision™ FLEX Peroxidase-Blocking Reagent (Dako Omnis)**

8 x 22.5 mL, ready-to-use

Phosphate buffer containing hydrogen peroxide, 15 mmol/L NaN<sub>3</sub> and detergent.

DM842

EnVision™ FLEX  
/HRP  
(Dako Omnis)

**EnVision™ FLEX /HRP (Dako Omnis)**

8 x 22.5 mL, ready-to-use

Dextran coupled with peroxidase molecules and goat secondary antibody molecules against rabbit and mouse immunoglobulins. In buffered solution containing stabilizing protein and preservative.

DM847

EnVision™ FLEX  
DAB+ CHROMOGEN  
(Dako Omnis)

**EnVision™ FLEX DAB+ Chromogen (Dako Omnis)**

8 x 1 mL

3,3'-diaminobenzidine tetrahydrochloride in organic solvent.

The color of this reagent may vary from strong violet to colorless without having any influence on the performance of the kit.

DM843

EnVision™ FLEX  
SUBSTRATE BUFFER  
(Dako Omnis)

**EnVision™ FLEX Substrate Buffer (Dako Omnis)**

16 x 26 mL

Buffered solution containing hydrogen peroxide and preservative.

DM848

EnVision™ FLEX  
TARGET RETRIEVAL SOLUTION  
HIGH pH (50x)  
(Dako Omnis)

**EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Dako Omnis)**

9 x 68 mL, 50x concentrated

Tris/EDTA buffer, pH 9

### B. Materials required but not provided

Antibody diluent, e.g. Dako Antibody Diluent (Code S0809)

Clarify Clearing Agent (Code GC810)

Dako FLEX Ready-to-Use Primary Antibodies or suitably diluted concentrated primary rabbit or mouse antibodies from Dako

Dako Omnis instrument

Dako Omnis instrumentation utensils

Dako Proteinase K, Ready-to-Use (if necessary) (Code S3020)

FLEX Universal Negative Control, Rabbit, Ready-to-Use (Code GA600), FLEX Universal Negative Control, Mouse, Ready-to-Use (Code GA750), or suitable negative control reagent for the primary antibody

Hematoxylin (Dako Omnis) (Code GC808) or suitable hematoxylin reagent for counterstaining

Microscope slides, e.g. Dako FLEX IHC Microscope Slides (Code K8020)

Mounting medium (organic-solvent-based) and coverslips

Wash Buffer (20x) (Dako Omnis) (Code GC807)

0.3 M Sulfuric Acid (Code GC203)

### C. Optional reagents

**GV804**  
DM848

EnVision™ FLEX  
TARGET RETRIEVAL SOLUTION  
HIGH pH (50x)  
(Dako Omnis)

**EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Dako Omnis)**  
3 x 68 mL, 50x concentrated  
Tris/EDTA buffer, pH 9

**GV805**  
DM849

EnVision™ FLEX  
TARGET RETRIEVAL SOLUTION  
LOW pH (50x)  
(Dako Omnis)

**EnVision™ FLEX Target Retrieval Solution, Low pH (50x) (Dako Omnis)**  
3 x 68 mL, 50x concentrated  
Citrate buffer, pH 6.1

**GC807**  
DM851

WASH BUFFER (20x)  
(Dako Omnis)

**Wash Buffer (20x) (Dako Omnis)**  
20 x 175 mL, 20x concentrated  
Tris-buffered saline solution containing Tween 20, pH 7.6 (±0.1).

**GC808**  
DM846

HEMATOXYLIN  
(Dako Omnis)

**Hematoxylin (Dako Omnis)**  
8 x 22.5 mL, ready-to-use  
Aqueous solution of hematoxylin.

**GV809**  
DM845

EnVision™ FLEX+  
RABBIT LINKER  
(Dako Omnis)

**EnVision™ FLEX+ Rabbit LINKER (Dako Omnis)**  
1 x 22.5 mL, ready-to-use  
Buffered solution containing stabilizing protein and an antimicrobial agent.

**GV821**  
DM844

EnVision™ FLEX+  
MOUSE LINKER  
(Dako Omnis)

**EnVision™ FLEX+ Mouse LINKER (Dako Omnis)**  
1 x 22.5 mL, ready-to-use  
Buffered solution containing stabilizing protein and an antimicrobial agent.

**GV825**  
DM847 +  
DM843

**EnVision™ FLEX DAB+ Substrate Chromogen System (Dako Omnis)**

EnVision™ FLEX  
DAB+ CHROMOGEN  
(Dako Omnis)

**EnVision™ FLEX DAB+ Chromogen (Dako Omnis)**  
2 x 1 mL  
3,3'-diaminobenzidine tetrahydrochloride in organic solvent.  
The color of this reagent may vary from strong violet to colorless without having any influence on the performance of the kit.

EnVision™ FLEX  
SUBSTRATE BUFFER  
(Dako Omnis)

**EnVision™ FLEX Substrate Buffer (Dako Omnis)**  
4 x 26 mL  
Buffered solution containing hydrogen peroxide and preservative.

## Precautions

1. For professional users.
2. EnVision™ FLEX /HRP contains material of animal origin and it cannot be excluded that trace amounts of human material could be present due to manufacturing procedures. As with any product derived from biological sources, proper handling should be used.
3. Do not expose EnVision™ FLEX Peroxidase-Blocking Reagent, EnVision™ FLEX /HRP, EnVision™ FLEX Substrate Buffer, EnVision™ FLEX DAB+ Chromogen or EnVision™ FLEX Substrate Working Solution to strong light.
4. EnVision™ FLEX Target Retrieval Solution, High pH (50x) contains 5-10% trometamol, 1-5% nonylphenol, ethoxylated, and 1-5% acetic acid, (ethylenedinitrilo)tetra-, disodium salt, dihydrate. EnVision™ FLEX Target Retrieval Solution, High pH (50x) is labeled:

### Warning

H319	Causes serious eye irritation.
H411	Toxic to aquatic life with long lasting effects.
P280	Wear eye or face protection.
P273	Avoid release to the environment.
P264	Wash hands thoroughly after handling.
P305 + P351+ P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337 + P313	If eye irritation persists: Get medical attention.
P501	Dispose of contents and container in accordance with all local, regional, national and international regulations.

5. For EnVision™ FLEX Peroxidase-Blocking Reagent, a safety data sheet is available for professional users on request.
6. EnVision™ FLEX DAB+ Chromogen contains 1-5% biphenyl-3,3',4,4'-tetrayltetraammonium tetrachloride, and is labeled:

### Danger

H350	May cause cancer.
H341	Suspected of causing genetic defects.
P201	Obtain special instructions before use.
P202	Do not handle until all safety precautions have been read and understood.
P280	Wear protective gloves. Wear eye or face protection. Wear protective clothing.
P308 + P313	IF exposed or concerned: Get medical attention.
P405	Store locked up.
P501	Dispose of contents and container in accordance with all local, regional, national and international regulations.

As a general rule, persons under 18 years of age are not allowed to work with this product. Users must be carefully instructed in the proper work procedure, the dangerous properties of the product and the necessary safety instructions. Please refer to the Safety Data Sheet for additional information.

7. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
8. Unused solution should be disposed of according to local, State and Federal regulations.

## Storage

Store EnVision™ FLEX Target Retrieval Solution, High pH (50x) at 2-8 °C. Store EnVision™ FLEX Peroxidase-Blocking Reagent, EnVision™ FLEX /HRP, EnVision™ FLEX DAB+ Chromogen and EnVision™ FLEX Substrate Buffer in the dark at 2-8 °C. During storage the cap should be closed.

Onboard stability of EnVision™ FLEX Peroxidase-Blocking Reagent, EnVision™ FLEX /HRP, EnVision™ FLEX DAB+ Chromogen and EnVision™ FLEX Substrate Buffer is 75 hours. Onboard stability of diluted EnVision FLEX™ Target Retrieval Solution, High pH and Low pH is 7 days. Remaining onboard stability is tracked by the Dako Omnis software.

Do not use after expiration date stamped on the vials. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the product is suspected, contact Dako Technical Support.

## Reagent Preparation

### **A.1 EnVision™ FLEX Target Retrieval Solution, High pH (50x) (DM848, Code GV804) or EnVision™ FLEX Target Retrieval Solution, Low pH (50x) (DM849, Code GV805)**

Prepare a working solution by diluting the EnVision™ FLEX Target Retrieval Solution, High pH (50x) or Low pH (50x) concentrate 1:50 in distilled or deionized water. The following procedure is recommended:

1. Fill a Dako Omnis bulk bottle to the fill line with distilled or deionized water (3.325 L). Ensure to place the bulk bottle on a horizontal surface before filling.
2. Empty one 68 mL EnVision™ FLEX Target Retrieval Solution, High pH (50x) or Low pH (50x) (Dako Omnis) concentrate bottle into the bulk bottle.
3. Transfer the detachable label from the concentrate bottle and attach it to the bulk bottle.
4. Fasten the lid and gently invert the bulk bottle 2-3 times.

Unused diluted solution may be stored at 2-8 °C for one month. The diluted solution can be used within 7 days when stored onboard the Dako Omnis instrument. Discard solution if cloudy in appearance.

### **A.2 Wash Buffer (20x) (DM851, Code GC807)**

Prepare a working solution by diluting the Wash Buffer (20x) (Dako Omnis) concentrate 1:20 in Dako Omnis bulk bottle using distilled or deionized water. To minimize foaming the following procedure is recommended:

1. Fill a Dako Omnis bulk bottle to the fill line with distilled or deionized water (3.325 L). Ensure to place the bulk bottle on a horizontal surface before filling.
2. Empty one 175 mL Wash Buffer (20x) (Dako Omnis) concentrate bottle into the bulk bottle.
3. Transfer the detachable label from the concentrate bottle and attach it to the bulk bottle.
4. Fasten the lid and gently invert the bulk bottle 2-3 times.

Unused diluted solution may be stored at 2-8 °C for one month. The diluted solution can be used within 7 days when stored onboard the Dako Omnis instrument. Remaining onboard stability is tracked by the Dako Omnis software. Discard solution if cloudy in appearance.

### **A.3 EnVision™ FLEX Substrate Working Solution**

The DAB-containing EnVision™ FLEX Substrate Working Solution is prepared onboard the Dako Omnis instrument. The automated process mixes one part of EnVision™ FLEX DAB+ Chromogen (DM847) with 50 parts of EnVision™ FLEX Substrate Buffer (DM843). Substrate Working Solution will always be prepared for one rack (containing 1-5 slides) at a time. The volumes supplied are sufficient for loading of 200 racks.

## Specimen Collection and Preparation

The specimens should be formalin-fixed, paraffin-embedded tissue sections.

Fixation time is dependent on fixative and tissue type/thickness. For example, tissue blocks with a thickness of 3-4 mm should be fixed in neutral-buffered formalin for 18-24 hours.

The optimal thickness of paraffin-embedded sections is 4 µm.

The specimens should be mounted on microscope slides. For greater adherence of tissue sections to glass slides, the use of FLEX IHC Microscope Slides (Code K8020) is recommended.

The sections should be mounted on the slides as flat and wrinkle-free as possible. Too many wrinkles will have an impact on the staining results.

**NOTE:** The microscope slides must have a width suitable for the Dako Omnis instrument. Please refer to Dako Omnis Basic User Guide for definition of usable microscope slides.

Paraffin sections should be mounted from a pre-heated water bath containing distilled or deionized water. The water bath should contain no additives (such as gelatin, poly-L-lysine etc.). Sections should be dried by heating, generally at a temperature not above 60 °C for up to 60 minutes. To ensure proper adherence of sections to slides it is important to drain the water from beneath the sections prior to the oven drying process.

## Procedure

EnVision™ FLEX visualization system has been developed for use on the Dako Omnis instrument in combination with Dako primary antibodies.

The automated procedure on Dako Omnis includes deparaffinization of tissue sections, heat-induced epitope retrieval, immunohistochemical staining and counterstaining of formalin-fixed, paraffin-embedded tissue sections. The protocol steps for the automated procedures can be displayed, printed and edited according to the Dako Omnis Basic and Advanced User Guides.

The pretreatment procedure and staining protocol depends on the the antibody. Please refer to the package insert for the individual Dako primary antibody.

Before running protocols on your Dako Omnis instrument, please read carefully the Dako Omnis Basic User Guide.

The EnVision™ FLEX protocol is pre-programmed in the Dako Omnis software. Please refer to Dako Omnis Basic User Guide for detailed instructions on preparing the instrument for staining and Dako Omnis Advanced User Guide for detailed instructions on protocols.

## Staining procedure

When using FLEX RTU Primary Antibodies for Dako Omnis the optimal staining protocol is programmed into the Dako Omnis software.

Please refer to the package insert for the individual FLEX RTU Primary Antibody or Dako concentrated Primary Antibody for optimal pre-treatment method and incubation times. For easy setup of new protocols, it is recommended to use one of the protocol templates described below.

From the Dako Omnis Workstation, chose one of three protocol templates containing the recommended staining steps and incubation times for the EnVision™ FLEX visualization system (refer to Table 2 for overview):

FLEX  
FLEX+ Mouse  
FLEX+ Rabbit

The protocol templates can be used for the generation of custom protocols for specific primary antibodies. The following protocol steps can be customized:

- Target retrieval reagent (for heat-induced epitope retrieval), temperature and incubation time
- Primary antibody and incubation time
- Labeled polymer incubation time (reduction only)
- Counterstain reagent and incubation time

The user must verify the applied protocol.

If the protocols are not available in the Dako Omnis system, please contact Dako Technical Services.

Table 2. Overview of EnVision™ FLEX/FLEX+ protocol templates.

Phase	Step	Reagent	Staining protocol template		
			FLEX	FLEX+ Mouse	FLEX+ Rabbit
Dewax	Dewax	Clarify Clearing Agent	yes	yes	yes
Target Retrieval	Target Retrieval	EnVision™ FLEX Target Retrieval Solution, High pH or EnVision™ FLEX Target Retrieval Solution, Low pH or Custom reagent	97±2 °C* 30 min*	97±2 °C* 30 min*	97±2 °C* 30 min*
Staining	<i>Wash</i>	Wash Buffer (Dako Omnis)	wash	wash	wash
	Primary Antibody	20 min*	20 min*	20 min*	20 min*
	<i>Wash</i>	Wash Buffer (Dako Omnis)	wash	wash	wash
	Endogenous Enzyme Block	EnVision™ FLEX Peroxidase-Blocking Reagent	3 min	3 min	3 min
	<i>Wash</i>	Wash Buffer (Dako Omnis)	wash	wash	wash
	Secondary Reagent	EnVision™ FLEX+ Mouse or Rabbit LINKER	-	Mouse LINKER 10 min	Rabbit LINKER 10 min
	<i>Wash</i>	Wash Buffer (Dako Omnis)	-	wash	wash
	Labeled polymer	EnVision™ FLEX/HRP	20 min*	20 min*	20 min*
	<i>Wash</i>	Wash Buffer (Dako Omnis)	wash	wash	wash
	<i>Wash</i>	Wash Buffer (Dako Omnis)	wash	wash	wash
	Substrate Chromogen	EnVision™ FLEX Substrate Working Solution	5 min	5 min	5 min
	<i>Wash</i>	Wash Buffer (Dako Omnis)	wash	wash	wash
Counter-staining	Counterstain	Hematoxylin (Dako Omnis)*	3 min*	3 min*	3 min*
	<i>Wash</i>	Deionized Water	wash	wash	wash
	<i>Wash</i>	Wash Buffer (Dako Omnis)	wash	wash	wash

\* Customization is possible.

### Post-staining procedure

1. After completion of the run, remove the slide rack from the Unloading Station and transfer racks to fresh deionized water.
2. Dehydrate, clear, mount with permanent mounting media and coverslip the slides.

### Quality Control

Each staining run should include a known positive control specimen to ascertain a proper performance of all the applied reagents. If the positive control specimen fails to demonstrate positive staining, labeling of test specimens should be considered invalid.

A negative control reagent should be used with each specimen to identify any non-specific staining. If non-specific staining cannot be clearly differentiated from the specific staining, the labeling of the test specimen should be considered invalid.

### Interpretation of Results

The diaminobenzidine-containing Substrate Working Solution gives a brown color at the site of the target antigen recognized by the primary antibody. The brown color should be present on the positive control specimen at the expected localization of the target antigen. If non-specific staining is present, this will be recognized as a rather diffuse, brown staining on the slides treated with the negative control reagent. Nuclei will be stained blue by the hematoxylin counterstain.

## **Limitations**

1. For automated use on Dako Omnis instruments only. When performing staining using Dako concentrated primary antibodies, the optimal procedure must be determined and verified by the user.
2. For use on formalin-fixed, paraffin-embedded tissue sections only. When using frozen sections or cell smears, the optimal procedure must be determined and verified by the user, unless specified in the Package Insert for the specific Dako Primary Antibody. Changes in the procedure may have impact on the number of tests provided with the kit.

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Dako Denmark A/S