



BLitz
User Guide

fortéBIO

Table of Contents

| | |
|-------------------------------------------------------|----------------------------------------------------|
| Chapter 1: | Chapter 4: |
| Welcome 1 | Working with Experiments31 |
| <i>Introduction</i>2 | <i>Introduction</i>32 |
| <i>BLItz System Package</i>2 | <i>Accessing Experiment Modules</i>32 |
| <i>Conventions Used in This Guide</i>4 | <i>Hydrating Biosensors</i>33 |
| <i>BLItz Support Options</i>5 | <i>Running a Quick Yes/No Experiment</i>35 |
| <i>Technical Support</i>5 | <i>Creating a Standard Curve</i>41 |
| <i>Depot Repair</i>5 | <i>Running a Quantitate Sample Experiment</i> 48 |
| Chapter 2: | <i>Running a Basic Kinetics Experiment</i>54 |
| BLItz System Specifications 7 | <i>Running an Advanced Kinetics Experiment</i> 61 |
| <i>BLItz System Specifications</i>8 | Appendix A: |
| <i>Basic Functions of the BLItz System</i> 10 | Running the BLItz Pro Software on |
| <i>Biosensors</i> 12 | Computers Running Octet Systems ..67 |
| <i>Dip and Read Assays</i> 12 | <i>Running Only the BLItz System on the</i> |
| <i>Biosensor Types</i> 13 | <i>Computer</i>68 |
| <i>BLI Technology</i> 15 | <i>Running Experiments on Octet and BLItz</i> |
| Chapter 3: | <i>Systems on the Same Computer</i>68 |
| BLItz Installation 17 | Appendix B: |
| <i>Installing the BLItz System</i> 18 | Troubleshooting71 |
| <i>Using the BLItz Pro Software User Interface</i> 20 | <i>Troubleshooting BLItz System Errors with an</i> |
| <i>Toolbars</i> 21 | <i>Octet System</i>72 |
| <i>Menu Commands</i> 22 | <i>Scenario 1—The BLItz System Cannot</i> |
| <i>Icons</i> 28 | <i>Initialize Successfully</i>72 |
| <i>Uninstalling the BLItz Pro Software</i> 30 | <i>Scenario 2—The BLItz System Is Not</i> |
| | <i>Recognized</i>72 |
| | Index 1 |

CHAPTER 1:
Welcome

Introduction 2

BLitz System Package 2

Conventions Used in This Guide 4

BLitz Support Options 5

 Technical Support 5

 Depot Repair 5

INTRODUCTION

Welcome to the *ForteBio BLItz User Guide*. This guide explains how to:

- Configure, install, and operate your BLItz system and BLItz Pro software.
- Set up and run kinetics experiments on the BLItz instrument.

The BLItz system and BLItz Pro software user guide contains information on materials required for using the BLItz system, installation procedure, setting up and running experiments and performing data analysis.

If you would like further information on any topic, contact technical support using information provided in “BLItz Support Options” on page 5.

The guide provides a set of specific instructions in the “Working with Experiments” on page 31 based on the Starter kit included with your BLItz system. We highly recommend to new users that they first perform the tests illustrated in this guide to familiarize themselves with BLItz system operations before embarking on their own research.

We wish you the best in your research!

BLITZ SYSTEM PACKAGE

The BLItz system enables real-time quantitation of solution-based analytes or kinetic characterization of molecular interactions.

The BLItz system package consists of the BLItz instrument, the BLItz Pro™ software CD, power cord, USB cable, *BLItz Quickstart Guide*, and a starter kit (Figure 1-1).



Figure 1-1: *BLItz System Package Contents*

Table 1-1 lists the BLItz system starter kit contents.



IMPORTANT: Immediately store the sample diluents, human IgG, mouse IgG, and hydrochloric acid reagents at 2–8°C.

| Reagents | Accessories | Materials Required (Not Provided) |
|------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| One (1) tray of 32 Protein A biosensors | 500 μ L of 0.5M Hydrochloric acid for deep cleaning drop holder | Pipettors: 2–20 μ L 20–200 μ L 100–1000 μ L |
| 200 μ L of 10 mg/mL human IgG stock | 20 Eppendorf tubes (0.5 mL) | Appropriate pipette tips |
| 100 μ L of 10 mg/mL mouse IgG in carrier protein containing buffer | 96-well, black, flat-bottom polypropylene (Greiner Bio-One, #655209) for hydrating biosensors | Gloves |
| 50 mL of sample diluent | Two (2) drop holders | Laptop or desktop computer—Ensure that your computer meets the following minimum requirements: <ul style="list-style-type: none"> • 2.5 GHz Dual Core CPU • 2 GB RAM • 32-bit or 64-bit Windows XP and Windows 7 • One (1) USB 2.0 port • One (1) monitor (minimum 1,024 X 768) |
| | Swab, pack of 20, to wipe drop holder | |

Table 1-1: BLItz System Kit Contents

For information about BLItz system specifications and basic functions, see “BLItz System Specifications” on page 7.

For information about installing and configuring the BLItz Pro software to use with the BLItz system, see “Software Overview” on page 9.

CONVENTIONS USED IN THIS GUIDE



NOTE: A note presents pertinent details on a topic; for example, general information about tips or alternate options.






IMPORTANT: An important message for instances where the experiment or procedure will not work if not properly followed.



WARNING: A warning warns you that your actions could cause irreversible consequences or damage.

Table 1-2: BLItz Instrument Labels

| Symbol | Definition |
|-------------------------------------------------------------------------------------|-------------------|
|  | Electrical hazard |
|  | Heat/hot |
|  | Fuse |

BLITZ SUPPORT OPTIONS

ForteBio wants to ensure that you are completely satisfied with the BLitz system and software, and will address any concerns in a timely manner—it's our way of supporting your research. We are also extremely interested in your product feedback and application needs, and would be happy to discuss either with you.

Technical Support

For technical questions or to speak with one of our support staff, email technical support via our online form which also lets you attach files. You can contact BLitz technical support at:

ForteBio, Inc.
1360 Willow Road, Suite 201
Menlo Park, CA 94025
USA
Tel: +1-855-BLITZ ME
Fax: +1-650-322-1370
E-mail: blitzsupport@fortebio.com
http://www.blitzmenow.com/email_support.html

Depot Repair

ForteBio offers a depot repair service. If your BLitz system needs service, please contact technical support (http://www.blitzmenow.com/email_support.html) for more information.

CHAPTER 2:

BLitz System Specifications

| | |
|-------------------------------------------|----|
| BLitz System Specifications | 8 |
| Basic Functions of the BLitz System | 10 |
| Biosensors | 12 |
| Dip and Read Assays | 12 |
| Biosensor Types | 13 |
| BLI Technology | 15 |

BLITZ SYSTEM SPECIFICATIONS

Figure 2-1 illustrates the BLItz system components (front and back).




Figure 2-1: BLItz System—Front and Back

Table 2-1 lists and defines the BLItz system specifications.

Table 2-1: BLItz System Specifications

| Specification | Descriptions |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Equipment Classifications | <ul style="list-style-type: none"> Product Classification: Class 1: Detachable power cord Installation/Overvoltage Category: Category II Pollution Degree: Degree 2 EMC Classification: Group I, Class A, ISM Equipment (EN55011, emissions), {EN61326, immunity} |
| Environmental | <ul style="list-style-type: none"> Storage Temperature: –20 to 70 °C Optimum Operating Temperature: 22 ± 4 °C Safe Operating Temperature: 15 to 30 °C Humidity: Non-condensing 10 to 80% Relative Humidity Indoor Use Only Operating Altitude: 0 to 2,000 meters |
| Compliance | CE |
| Capabilities | <ul style="list-style-type: none"> Protein quantitation Kinetic and affinity analyses (k_a, k_d, K_D) Binding specificity and cooperativity Kinetic screening of proteins, peptides, and other bio-molecules |
| Sampling Format | Drop holder and Eppendorf tube |

Table 2-1: BLItz System Specifications (Continued)

| <i>Specification</i> | <i>Descriptions</i> |
|-------------------------|------------------------------------------------------------------------------------------|
| Sampling Volume | 4 μ L in drop holder; 250 μ L in tube |
| Sample Types | Purified samples, common culture media, crude lysates |
| Biosensor Type | Disposable fiber-optic biosensors |
| Optics and Mechanics | 1-channel biosensor manifold, one (1) spectrometer |
| Shake Rate | Static or 1000–2600 rpm (default is 2200) |
| Dimensions | 6.8" H x 6.0" W x 8.7" D (17.4 cm H x 15.3 cm W x 22.2 cm D) |
| Weight | 7.2 lb (3.3 Kg) |
| Electrical Requirements | Power consumption: 8w (18W peak) Mains: AC 100–240 V, 0.2–0.1A, 50/60Hz, single phase |
| Connections | Power input, USB data output port |
| Software | BLItz Pro™ software for data acquisition and data analysis |

BASIC FUNCTIONS OF THE BLITZ SYSTEM

Table 2-2 lists the basic functions of the BLItz system.







| <i>Function</i> | |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | <p>Hydrating Biosensors—Always hydrate biosensors for at least ten minutes before use. The buffer your sample is in is the best buffer to use for biosensor hydration. Do the hydration on the bench in the blue biosensor tray using 200 μL of buffer per well. For more information, refer to “Hydrating Biosensors” on page 33.</p> |
|  | <p>Setting Up an Experiment—Always start with the BLItz Pro software and follow the onscreen instructions.</p> <ol style="list-style-type: none"> 1. Choose an Experiment type from the Experiments menu. 2. Input experiment information and click NEXT. 3. Follow the prompts to load the sample and biosensor on the BLItz system. |
|  | <p>Placing a Drop of Sample in the Drop Holder Pipette a 4 μL drop of sample inside the dimple in the drop holder. Ensure that there are no air bubbles in the liquid.</p> |
|  | <p>Placing and Removing Drop Holder—The drop holder is magnetic and allows easy snapping into place next to the tube holder. When required, simply pick up the drop holder from its place. The drop holder can be replaced with a new one when damaged or when extensive contamination is suspected.</p> |
|  | <p>Wiping sample with a swab (or other lint-free lab wipe)—As soon as an experiment is complete, wipe the sample from the drop holder with a swab provided in the starter kit. Rinse the drop holder three (3) times with buffer, blotting each rinse with the swab. A kimwipe may be used instead of the swab. If the drop dries on the drop holder, perform deep cleaning by repeating the above actions with 0.5N HCl instead of buffer. Follow up with rinse and blot with buffer to remove all HCl.</p> |
|  | <p>What is the use of the tube holder? The drop holder is the only sample location you should use in Quick Yes/No, Create Standard Curve, and Quantitate Sample experiments. In Basic Kinetics and Advanced Kinetics experiments, the tube holder is recommend for all baseline and dissociation steps, and for association steps longer than five (5) minutes. The volume of reagent to use in the tube is 250 μL.</p> |

Table 2-2: Basic Functions of the BLItz System



| <i>Function</i> | |
|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Selecting Drop Holder and Tube Holder with Slider —Moving the slider to the right brings the drop holder in line with the black triangular mark and thus, to the read position. Moving the slider to the left brings the tube holder in line with the black triangular mark and thus, to the read position. |
|  | Mounting a Biosensor —Apply a hydrated biosensor on to the mount and give a final quarter-turn-upward twist to ensure a snug fit. Take care not to touch the tip of the biosensor on solid surfaces. Do not let the biosensor dry. This will result in loss of performance. Keep the biosensor hydrated until ready to use. After mounting a biosensor on the BLItz system, perform the experiment immediately. |

Table 2-2: Basic Functions of the BLItz System (Continued)

BIOSENSORS

Dip and Read Assays

BLItz uses ForteBio's Dip and Read™ label-free assays. These direct binding assays take place on a disposable biosensor made from a biocompatible matrix that is uniform, non-denaturing and minimizes non-specific binding (Figure 2-2). Only molecules that bind directly to the biosensor surface are detected, providing exceptional specificity for individual applications, even in crude media.




Figure 2-2: *Dip and Read Assays*

Biosensor Types

A variety of disposable, off-the-shelf biosensors are available for the BLItz system, making it easy to run a wide range of assays on the same system. The BLItz system uses the same Dip and Read biosensors that are used on BLItz's Octet platform, and over four million assays have been run to date. You can be confident that the assays you run on the BLItz system will provide the same high level of accuracy and precision. Our biosensors are manufactured in an ISO 9001:2008 certified facility. To find a current list of biosensors and to place an order, go to <http://www.blitzmenow.com/biosensors.html>.

Table 2-3: Biosensor Types

| Type | Required Capture Molecule | Analyte Measured | Description |
|----------------------|---------------------------|------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anti-hIgG-Fc | None | hIgG, human Fc-fused proteins | Binds specifically to the Fc portion of human IgGs and other proteins containing a human Fc region. Applications include quantitation of human IgG and proteins containing the human Fc region in cell line development, clone selection, process optimization, and production monitoring. |
| Anti-mIgG-Fv | None | mIgG, rat IgG, mouse-Fab, rat Fab | Binds specifically to the Fv(ab') ₂ portion of mouse and rat IgGs. Applications include quantitation of mouse and rat IgG in cell line development, clone selection, process optimization, and production monitoring. |
| Protein A | None | Many human and other IgG types | Binds with high affinity to the Fc region of human IgGs. Binds with moderate affinity to many subtypes of mouse and rabbit IgG. Applications include quantitation of IgG in cell line development, clone selection, process optimization, and production monitoring. |
| Protein G | None | Many murine and other IgG types | Binds with high affinity to the Fc region of murine, rat, goat, and bovine IgGs. Binds with moderate affinity to many subtypes of human IgG. Applications include quantitation of IgG in cell line development, clone selection, process optimization, and production monitoring. |
| Protein L | None | Most mouse, rat, human IgG and Fab | Binds with high affinity to most mouse, rat and human immunoglobulins containing a kappa light chain. Does not bind goat, bovine, rabbit or sheep IgG. Applications include quantitation of FAb fragments and of IgG in serum-based culture. |
| Anti-Penta-HIS (HIS) | None | His-tagged proteins, peptides | Uses the Qiagen Penta-HIS antibody to bind with high affinity HIS tagged recombinant proteins. Applications include quantitation of HIS-tagged proteins. |

Table 2-3: Biosensor Types (Continued)

| Type | Required Capture Molecule | Analyte Measured | Description |
|----------------------------|------------------------------------------|----------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anti-GST (GST) | None | GST-tagged proteins, peptides | High affinity anti-GST antibody on biosensor binds GST-tagged proteins. Applications include quantitation of GST-tagged proteins and kinetics of proteins and peptides binding to GST-tagged proteins captured on biosensor. |
| Streptavidin (SA) | Biotin-tagged peptides, oligos, proteins | Proteins, peptides, oligos | Streptavidin-coated biosensor. Immobilizes biotinylated antibodies, proteins and nucleic acids to form a stable surface. Quantitation applications include antibody and protein quantitation in cell line development, clone selection, process optimization and production monitoring. Kinetic applications include protein and antibody kinetic and affinity analysis (k_a , k_d , K_D), and epitope binning. |
| Amine Reactive (AR2G) | Proteins, peptides, oligos | Proteins, antibody fragments | Second generation carboxylate functionalized surface allows covalent coupling of proteins via EDC/s-NHS mediated amide bond formation. Kinetic applications include protein and antibody kinetic and affinity characterization (k_a , k_d , K_D). Requires AR2G Assay Kit (P/N 18-5095). |
| Anti-hIgG-Fc Capture (AHC) | hIgG, human fc fusion protein | Proteins, peptides, antibody fragments | Immobilization of human IgG or other human Fc-containing proteins by binding to the human Fc region. Kinetic applications include protein and antibody kinetic and affinity characterization (k_a , k_d , K_D) and epitope binning. |
| Anti-mIgG-Fc Capture (AMC) | mIgG, mouse fc fusion protein | Proteins, peptides, antibody fragments | Binds the Fc portion of IgG1, IgG2a and IgG2b for capture-based immobilization. Applications include kinetic analysis of antibody-antigen interactions (k_a , k_d , K_D) and off-rate screening. IgG3 should be evaluated on a case-by-case basis. |
| Aminopropylsilane (APS) | Proteins, peptides | Proteins, peptides | Adsorption of proteins and membrane fractions through hydrophobic moieties. Kinetic applications include adsorption of proteins and membrane fractions through hydrophobic moieties for kinetic and affinity analysis (k_a , k_d , K_D). |

BLI Technology

The BLI (bio-layer interferometry) technology (Figure 2-3) used by the BLItz system provides real-time data on protein interactions. The BLItz system emits white light down the biosensor, and then collects any light reflected back. Reflected wavelengths are affected by the thickness of the coating on the optical layer. Some wavelengths show constructive interference (blue), others show destructive interference (red).




Figure 2-3: BLI Technology—Displaying Constructive and Destructive Interference

This interference is captured by a spectrometer as a unique spectral signature and is reported in relative intensity units (nm). Any change in the number of molecules bound to the biosensor causes a shift in the interference pattern that is measured in real time. This wavelength shift is a direct measure of the change in optical thickness (nm) of the biological layer.

Want to learn more about BLI technology? Call us at 1-855-BLITZ-ME or email at blitzsupport@fortebio.com.

CHAPTER 3:

BLitz Installation

| | |
|---------------------------------------------------|----|
| Installing the BLitz System | 18 |
| Using the BLitz Pro Software User Interface | 20 |
| Toolbars | 21 |
| Menu Commands | 22 |
| Icons | 28 |
| Uninstalling the BLitz Pro Software | 30 |

INSTALLING THE BLITZ SYSTEM

To install and configure the BLItz system and BLItz Pro™ software in your lab:

1. Remove the BLItz system, power cord, USB cable, and BLItz Pro software CD from the box. Place the BLItz system on a dry, level surface with minimal vibration, and away from direct sunlight.
2. Insert the BLItz Pro Software CD in a laptop or desktop computer and run the **BLItzInstallCD.exe** file. For minimum computer specifications, see Table 1-1 on page 3.
3. Click **Install BLItz Pro Software** (Figure 3-1) to install the BLItz Pro software.
4. Click **Install BLItz Driver** (Figure 3-1) to install the BLItz driver.




Figure 3-1: Installing the BLItz Pro Software

5. Click **Exit** to complete software installation.
6. Connect the power cord to the BLItz system and plug the other side into an outlet (Figure 3-2).
7. Connect one side of the data cord to the BLItz system and the other to a USB port on your computer (Figure 3-2).



Figure 3-2: BLItz System (Back)

- Place a drop holder in its slot on the BLItz system (Figure 3-3). Store the second drop holder provided with the system for later use.





Figure 3-3: Placing a drop holder on the BLItz system. The drop holder will be held intact on the BLItz system by a magnet.

- Power on the BLItz system and wait for the Windows plug-and-play manager to assign drivers.
- Turn on the BLItz system. It may take up to five (5) minutes for the hardware to initialize and start communicating with the BLItz Pro software.



IMPORTANT: Turn on the BLItz system for one (1) hour prior to acquiring data.

Once you have installed and configured the BLItz Pro software, you can open (launch) the software. On the desktop, double-click the **BLItz Pro 1.0**  icon.



IMPORTANT: To install the BLItz system and BLItz Pro software, you should first install the software on your computer and go through the steps outlined below in the order specified.



IMPORTANT: Octet customers may wish to install the BLItz system and BLItz Pro software on the same computer that runs their Octet system. Such customers should be aware that ForteBio has released new versions of Octet software (Versions 6.4.1.3 and 7.0.1.3) that are compatible with BLItz Pro software. Please contact BLItz's technical support team by calling 1-855-BLITZ-ME or email blitzsupport@fortebio.com to get a copy of the new Octet software versions.

Until you update your Octet software, you can still install and run BLItz system and BLItz Pro software on the same computer that runs Octet system and software. See Appendix A, "Running the BLItz Pro Software on Computers Running Octet Systems" on page 67 for instructions.

USING THE BLITZ PRO SOFTWARE USER INTERFACE

This section describes the BLItz Pro software user interface (Figure 3-4) elements.







Figure 3-4: BLItz Pro Software Welcome Page

Toolbars

Table 3-1 describes the BLItz Pro software toolbars.

Table 3-1: BLItz Pro Software Toolbars

| Experiment | Description |
|----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Title Bar | <p>The Title Bar (Figure 3-6) is located at the top of the GUI window page and refers to the file name and BLItz Pro software version. Until the user saves the file, the file name is “Untitled”.</p>  <p><i>Figure 3-5: BLItz Pro Title Bar</i></p> |
| Main Menu Bar | <p>The Main Menu Bar (Figure 3-6) is located at the top of the GUI window, and displays the main menus available in the BLItz Pro software. See “Menu Commands” on page 22. for a definition of main menu bar commands.</p>  <p><i>Figure 3-6: BLItz Pro Main Menu Bar</i></p> |
| Status Bar | <p>The Status Bar (Figure 3-7) is located at the bottom of the GUI window page and displays current system and experiment status. See Table 3-6 for a definition of the status bar icons.</p>  <p><i>Figure 3-7: Status Bar</i></p> |

Menu Commands

File Menu

The **File** menu (Figure 3-8) allows users to open and save method files, view experiments, print files, and set system and software options. Table 3-2 describes the menu commands and their respective functions.




Figure 3-8: File Menu

Table 3-2: File Menu Commands and Functions


| Menu Command | Function |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| File > New | Opens a new file after clearing out the existing experiment information that may exist on the software screen from a previous run. |
|  NOTE: This menu option performs the same action as clicking the New Experiment icon. | |
| File > Open | Opens a previously saved experiment. |
| File > Save or File > Save As | Saves an experiment data file. |
| File > Export to CSV | Saves the results table to a .csv file that can be opened in a spreadsheet application. |
| File > Create Report | Creates a report of the experiment in PDF file format. |

Table 3-2: File Menu Commands and Functions

| Menu Command | Function |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| File > Options | <p>Defines the system options available. To view the BLItz system options:</p> <ol style="list-style-type: none"> 1. On the main menu, click File > Options. 2. On the Options dialog box (see Figure 3-9), confirm the default settings or enter new settings, and click OK. |



Figure 3-9: Options Dialog Box

| | |
|-------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Data Options— Significant digits | Specifies the number of decimal places (from two to nine) for the computed data. The default setting is 4 . |
| Shaker— Shaker Speed (rpm) | Shaker speed (rpm)—Sample platform orbital shaking speed (rotations per minute). The default setting is 2200 . The shaker can be set to speeds in the range of 1000–2600 rpm. |
| | <p>NOTE: The default of 2200 rpm should be used in all kinetics experiments. The default speed also works well for the Quick Yes/No, Create Standard Curve and Quantitate Sample experiments. In these experiments, to work with high concentration samples, you may reduce the shake speed to lower levels.</p> |
| Notifica- tion Sounds | When enabled (checked), plays sounds through computer speakers at the end of the hydration process, the end of a step, or at the end of a run. Click Reset Sounds to revert back to the default sound settings. |

Table 3-2: File Menu Commands and Functions

| <i>Menu Command</i> | <i>Function</i> |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------|
| | Click the Always open in simulation mode if not connected to an instrument check box if you anticipate to work offline. |
| File > Exit | Closes the software after prompting you to save any changes. |

Experiment Menu

The **Experiment** menu (Figure 3-10) displays the experiment modules available to run. Once you launch the BLItz Pro software, you may begin your experiment using one of the experiment modules available in the left pane (Figure 3-4: on page 20).

Table 3-3 describes the menu commands and their respective functions.




Figure 3-10: Experiment Menu

Table 3-3: Experiment Menu Commands

| Menu Command | Function |
|----------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Experiment > Welcome | Displays the Welcome page upon launching the BLItz Pro software. Provides a high-level overview of the BLItz system and BLItz Pro software capabilities. |
| Experiment > Quick Yes/No | Opens the Quick Yes/No module for experimentation. Determines the presence of an analyte with a simple dip and read. Can include positive/negative control samples to compare the unknown results. |
| Experiment > Create Standard Curve | Opens the Create Standard Curve module for experimentation. Runs multiple standard concentrations to create a standard curve. For best results in quantitation, proteins used as standards should be identical to unknown and should be measured in the same sample milieu. |
| Experiment > Quantitate Sample | Opens the Quantitate Sample module for experimentation. Quantitates unknown samples against a standard curve. |
| Experiment > Basic Kinetics | Opens the Basic Kinetics module for experimentation. Generates binding curves for proteins of known. Use this module to run kinetics experiments when the ligand is already bound to the biosensor, either as supplied by ForteBio, or, through an offline incubation of biosensor with ligand on the bench top. |

Table 3-3: Experiment Menu Commands (Continued)

| Menu Command | Function |
|---------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Experiment > Advanced Kinetics | Opens the Advanced Kinetics module for experimentation. Enables ligand loading, followed by the generation of binding curves, and measures kinetic constants (k_a , k_d , K_D) for proteins of known concentrations. |
| Experiment > Convert to Create Standard Curve | Switches the current module to the Create Standard Curve module. Use if the data generated in a Quantitate Sample experiment needs to be transferred to the Create Standard Curve module for treatment as standard samples. |
| Experiment > Convert to Quantitate Sample | Switches the current module to the Quantitate Sample module. Use if the data generated in a Create Standard Curve experiment needs to be transferred to the Quantitate Sample module for quantitation. |

Instrument Menu

The **Instrument** menu (Figure 3-11) resets or stops the BLItz Pro system during an experiment. Table 3-4 describes the menu commands and their respective functions.




Figure 3-11: Instrument Menu

Table 3-4: Instrument Menu Commands

| Menu Command | Function |
|------------------------------|-----------------------------------|
| Instrument > Reset | Resets the system. |
| Instrument > Stop! | Stops the current experiment run. |

Help Menu

The **Help** menu (Figure 3-12) displays the links to BLItz user guides, the BLItz Pro software version, and license information. Table 3-5 describes the menu command and its respective function.




Figure 3-12: Help Menu

Table 3-5: Help Menu Commands

| <i>Menu Command</i> | <i>Function</i> |
|-----------------------------------|-----------------------------------------------------------------------------------------------------------------|
| Help > Quickstart Guide | Opens the latest PDF of the <i>BLItz System Quickstart Guide</i> . |
| Help > User Guide | Opens the latest PDF of the <i>BLItz System and BLItz Pro Software User Guide</i> . |
| Help > BLItz Web Site | Opens your browser to the BLItz web site (http://www.blitzmenow.com). |
| Help > About BLItz | Displays the BLItz properties and license information. |

Icons

Table 3-6 defines the BLItz Pro software icons and status markers.

Table 3-6: BLItz Pro Software Icons



















| Icon | Function |
|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | BLItz Pro 1.0 —Represents the BLItz Pro 1.0 software. |
|  | Ready —Indicates that the BLItz Pro software is ready to use. |
|  | Drop Holder —Indicates that the drop holder is in the read position. |
|  | Tube —Indicates that the tube position is in the read position. |
|  | Hydrate Biosensor —A timer to use for the biosensor hydration step. Set up the hydration time and click the Hydrate Biosensor icon to start. You may stop the hydration timer by clicking the Stop icon. Otherwise, after hydration, the timer runs down to zero (0) and the icon changes to Hydration Done . The window shakes until the Hydration Done icon is clicked. |
|  | Stop —Stops the hydration or experiment. |
|  | Hydration Done —Displays when the hydration process is complete. Is preceded by a ringtone to indicate completion. See the Table 3-2 on page 22 about notification sounds. The window shakes until the Hydration Done icon is clicked. |
|  | New Experiment —Creates a new experiment, or clears your existing experiment for a new experiment. |
|  | NOTE: This menu option performs the same action as clicking the File > New command sequence. |
|  | Next —Guides you to the next step of the experiment run process. |
|  | Initial Baseline —Mandatory step type used in Advanced Kinetics experiments. |
|  | Loading —Step type used in Advanced Kinetics experiments to load ligand on biosensors. |

Table 3-6: BLItz Pro Software Icons

| Icon | Function |
|------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  Baseline | Baseline —Step type used in Basic and Advanced Kinetics experiments for the Baseline immediately preceding the Association step.  NOTE: Baseline and Dissociation steps should be performed in the tube. |
|  Association | Association —Step type used in Basic and Advanced Kinetics experiments to monitor the binding of an analyte to a ligand on the biosensor. Association should be monitored in tube for incubation times more than five (5) minutes. Up to five (5) minutes of association can be performed in the drop holder or tube. |
|  Dissociation | Dissociation —Step type used in Basic and Advanced Kinetics experiments to monitor the dissociation of an analyte from a ligand. Dissociation should be monitored in the tube only. |

UNINSTALLING THE BLITZ PRO SOFTWARE

To uninstall the BLItz Pro software:

1. From the Windows **Start** menu, click the  icon (or the  button) and click **All Programs > ForteBio > Uninstall BLItz Pro 1.0**.

You are prompted to confirm the request remove the BLItz Pro 1.0 software and all of its components (Figure 3-13).




Figure 3-13: Uninstall Confirmation Request

2. Click **Yes** to confirm the removal of the BLItz Pro software and all of its contents. Another prompt displays confirming the removal (Figure 3-14).




Figure 3-14: Confirmation of BLItz Pro Software Removal

3. Click **OK** to exit.

CHAPTER 4:

Working with Experiments

| | |
|-----------------------------------------------|----|
| Introduction | 32 |
| Accessing Experiment Modules | 32 |
| Hydrating Biosensors | 33 |
| Running a Quick Yes/No Experiment | 35 |
| Creating a Standard Curve | 41 |
| Running a Quantitate Sample Experiment | 48 |
| Running a Basic Kinetics Experiment | 54 |
| Running an Advanced Kinetics Experiment | 61 |

INTRODUCTION

The section provides a set of specific instructions based on the Starter kit included with your BLItz system. We highly recommend to new users that you first perform the tests illustrated in this guide to familiarize yourselves with BLItz system operations before embarking on your own experiments. The specific experiment instructions should help you in designing and running your own experiments.

ACCESSING EXPERIMENT MODULES

To access the experiment modules:

1. Launch the BLItz Pro™ software: on the desktop, double-click the **BLItz Pro 1.0** icon.
The **BLItz Pro Welcome** page (Figure 4-1) displays with the **Experiments** listed in the left pane.
2. You can access any Experiment type by clicking the icons in the left pane, or by clicking the **Experiment** menu.

For a description of the **Experiment** menu commands and their respective functions and experiment type definitions, see Table 3-3 on page 25.




Figure 4-1: BLItz Pro Software Welcome Page

HYDRATING BIOSENSORS



IMPORTANT: Biosensors need to be hydrated before use in any experiment.

* Hydrate only those sensors to be used in an experiment.

* Hydrate just before running an experiment.

* Hydrate for at least 10 minutes.


* Do not leave the biosensors to dry on the biosensor mount.

ForteBio's biosensors need to be dipped in buffer to hydrate the bio-layer at the tip and prepare it for use in a run. Hydration is performed on the bench top in the blue tray provided. A 96-well microplate is needed for the hydration.

In the BLItz Pro software, each experiment module reminds you to perform this hydration step. The timer included in the software may be used to time your hydration. Alternatively, you may use a lab stop watch to time your hydration. You do not have to hydrate exactly for 10 minutes; a minimum of 10 minutes is all that is needed.

Use the following instructions to set up biosensor hydration.

To hydrate biosensors:

1. Add 200 μ L of sample diluent or buffer only in wells of a 96-well microplate corresponding to the biosensors that will be used in the experiment; for example, if using biosensors A1 to A4 from the biosensor tray, add buffer only to wells A1 to A4 in the 96-well microplate.
2. Remove the clear lid of the Protein A biosensor tray.
3. Lift the green biosensor rack carefully out of the blue biosensor tray holder avoiding touching of the biosensor tips on any solid surface.
4. Place the hydration microplate inside the blue biosensor tray holder.
5. Carefully lower the green biosensor rack on top of the hydration microplate in the biosensor tray holder to begin hydrating the biosensors.
6. Select a hydration time of **10 minutes** (the default) or more, and click the **Hydrate Biosensor** icon  (Figure 4-2).




Figure 4-2: Hydrate Biosensors

- Biosensors should be hydrated in a solution buffer that is as similar to the sample matrix as possible.
- The hydration timer in the BLItz Pro software may be used to time the biosensor hydration process.
- Users may also use other lab stop watch methods to monitor the hydration time.


The **Hydrate Biosensors** icon changes, indicating that the hydration process begins. The **Hydration time (min)** field timer counts down the 10 minutes. You can click the **Stop Hydration** icon  at any time to stop the hydration process (Figure 4-3).



Figure 4-3: Stop Hydration Process

After hydration, the timer runs down to zero (0) and the icon changes to **Hydration Done** (Figure 4-4). The window shakes until the **Hydration Done** icon is clicked.




Figure 4-4: Hydration Done



NOTE: After hydrating, leaving a biosensor on the biosensor mount or in air to dry will affect performance of the biosensor. If you do not plan to perform an experiment immediately, leave the biosensor tip to hydrate in buffer.

7. Access any Experiment type by clicking the icons in the left pane, or by clicking the **Experiment** menu (Figure 3-11: on page 21).

RUNNING A QUICK YES/NO EXPERIMENT

When using the Starter kit to perform a **Quick Yes/No** experiment, perform the following steps for human IgG samples binding Protein A biosensors.

To develop a **Quick Yes/No** experiment for your samples on the BLItz system, consider these guidelines:


- Test a negative control and positive control sample along with your unknown.
- Use an appropriate biosensor type from ForteBio's list of off-the-shelf biosensors. If one is not readily available, you may prepare a custom biosensor by loading your ligand on any one of the biosensors available.
- Use the drop holder for all **Quick Yes/No** measurements. Use 4 μL of sample diluent.
- Start with the default Run time value of **120** seconds to test your positive and negative controls. You may then optimize the Run time as needed. Typically, longer run times allow more binding of analyte to the biosensor, thus providing greater signal.
- Keep the Shaker in **Enable** mode to allow good mixing of sample. You can tweak the shaker rate via **File > Options** menu to a value between 1000 and 2600 rpm to suit your needs. Typically, increasing shaker rate leads to faster binding rates.
- Always input experiment info first and click **Next**. Follow the onscreen prompts to load sample and biosensor on the BLItz system.


In this **Quick Yes/No** experiment, hIgG is the analyte you will use to prepare three (3) samples: positive control (1000 $\mu\text{g}/\text{mL}$), "unknown" (100 $\mu\text{g}/\text{mL}$), and negative control (sample diluent). The 100 $\mu\text{g}/\text{mL}$ hIgG sample will be considered a mock "unknown" sample to show how BLItz can evaluate the presence of your protein in a sample. The negative control (sample diluent) will also be used as the sample blank (reference).

Everything you need:


- Protein A biosensors
- hIgG stock (10 mg/mL)
- sample diluent
- 96-well, black flat-bottom microplate


To run a **Quick Yes/No** experiment:

1. Click the **Quick Yes/No** icon  to display the **Quick Yes/No** experiment module (Figure 4-5).



Quick Yes/No

Have you hydrated your biosensor?
Hydrate biosensors prior to use for a minimum of 10 minutes. Hydration time (min): 10 

Enter your Experiment Settings 
Experiment Name:
Description:


Enter your Run Settings
Biosensor type: Protein A Run time (s): 120
Sample ID: Shaker: Enable
Concentration: µg/ml 

Figure 4-5: Quick Yes/No Experiment Module

Preparation

2. Hydrate three (3) biosensors, as described in “Hydrating Biosensors” on page 33.
3. Prepare the hlgG test samples as follows:
 - 1000 µg/mL hlgG (positive control): 20 µL of stock hlgG + 180 µL of sample diluent
 - 100 µg/mL hlgG (unknown sample): 10 µL of Stock A + 90 µL of sample diluent



NOTE: hlgG samples should be run in this order:

Run 1: sample blank (negative control)

Run 2: 100 µg/mL (unknown)

Run 3: 1000 µg/mL (positive control)

Experiment

- Input **Run Settings** for **Negative Control**, as shown in Figure 4-6.

The screenshot shows the BLitz software interface. The top section is titled "Enter your Experiment Settings" and contains two input fields: "Experiment Name" with the value "Quick Yes/No Demo" and "Description" with the value "Check for presence of human IgG in sample using Protein A biosensors". A "New Experiment" button is visible on the right. The bottom section is titled "Enter your Run Settings" and contains several fields: "Biosensor type" set to "Protein A", "Run time (s)" set to "30", "Sample ID" set to "Negative Control", "Shaker" set to "Enable", and "Concentration" set to "0" with a unit of "µg/ml". A "Run 1" button is visible on the right.

Figure 4-6: Entering Run Settings for Negative Control


- Enter a Run time of **30 seconds**.
 - Click .
- Follow the onscreen instructions (Figure 4-7).



Figure 4-7: User Action Required—Experiment Operations

- Add 4 µL of the sample diluent for the first run into the drop holder (Figure 4-8).



Figure 4-8: Adding 4 µL of Sample into Drop Holder

- Load the Protein A biosensor on to the biosensor mount; give it a final quarter-turn upward twist to ensure a snug fit (Figure 4-9).



Figure 4-9: Mounting a Biosensor

8. Move the slider to the right to move the drop holder to the read position (Figure 4-10).




Figure 4-10: Moving Drop Holder to Read Position

9. Close the BLItz system cover to start data acquisition each time the BLItz Pro software prompts you (Figure 4-11).



Figure 4-11: Closing the BLItz System Cover

Data is collected in real-time on screen. The sensorgram (Figure 4-12) displays the results for Negative Control.




Figure 4-12: Data Collected in Real-Time—Binding Signal Measured in nanometers (nm) as a Function of Time (seconds)

After the run, a **Save As** dialog box displays and prompts you to save the experiment file (Figure 4-13).




Figure 4-13: Saving Dataset File

10. Enter a file name and **Save** the dataset.
11. Lift the cover, discard the used biosensor, wipe the drop holder with a swab or Kimwipe (Figure 4-14), and rinse 3X with sample diluent—wiping each time with a swab or Kimwipe.



Figure 4-14: Wiping the Drop Holder with a Kimwipe

12. Perform subsequent sample runs by repeating steps 4 to 11 for each sample:
 - Run 2: 100 $\mu\text{g}/\text{mL}$ (unknown)
 - Run 3: 1000 $\mu\text{g}/\text{mL}$ (positive control)

Data

13. Inspect the **Run** data.

The sensorgram (Figure 4-15) displays the results after all runs. The unknown sample displays binding signal, indicating the presence of human IgG.




Figure 4-15: Inspecting the Run Data—Binding Signal Measured in nanometers (nm) as a Function of Time (seconds)

- Check the Run List table for the Binding Rate value for the unknown sample. A higher value than that of the Negative Control should mean that protein is present in the unknown sample (Figure 4-16).

| Run List | | | | | | | | | | Report point (s): 0 | Add | Remove All |
|----------|-------------------------------------|-------------------------------------|-------------------------------------|------------------|---------------|---------------------|--------------|--------------------|----------------|---------------------|-----|------------|
| Index | Show | Analyze | Ref. | Sample ID | Conc. (µg/ml) | Binding rate (nm/s) | Run time (s) | Shaker speed (rpm) | Biosensor Type | Information | | |
| 1 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | Negative Control | 0 | | 30 | 2200 | Protein A | | | |
| 2 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Unknown | | 0.5029 | 30 | 2200 | Protein A | | | |
| 3 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Positive Control | 1000 | 3.393 | 30 | 2200 | Protein A | | | |

Export to CSV... Create Report...

Figure 4-16: Saving Your Experiment

- Report point (s)**—Specifies a time point at which to report signal for all samples in the graph.
 - Add**—Adds a report point to the table.
 - Remove All**—Removes all the report points from the table.
- Export to CSV**—Exports the data in the graph to .csv files, one for each sample.
- Create Report**—Creates a report of the experiment in .PDF file format.
- Ref**—Select one of the samples as reference to subtract from all sample data. Selecting more than one sample as reference takes an average of the references to subtract from all data.
- Binding rate (nm/s)**—Rate of change of binding signal with time; slope of the early part of the binding curve for each sample.
- Shaker speed (rpm)**—Shaking speed of sample; this value is set in **File > Options**.

Result

The “unknown” sample, which we know is a 100 µg/mL hlgG solution, displayed a binding signal that indicates the presence of protein. In just a few minutes, you have successfully tested a sample for presence of protein. Bravo!

CREATING A STANDARD CURVE

When using the starter kit to perform a Create Standard Curve experiment, follow the instructions provided below for human IgG samples binding Protein A biosensors.

To create a Standard Curve for your samples on the BLItz system, consider these guidelines:

- Use standards that are representative of the analyte/unknown in your Quantitate Sample experiment.
- Use a sample milieu that is representative of that in which the analyte/unknown will be present in your Quantitate Sample experiment.
- Start with a broad concentration range of standard samples that will cover the analyte/unknown samples that you wish to quantitate.
- When optimizing experiment conditions, use fewer standard samples distributed over the concentration range. After optimizing experiment conditions, populate more standard concentrations to maximize reliability of the standard curve.
- Run standards from low to high concentrations to minimize interferences from carry-over between samples. Cleaning the drop holder according to recommendations should typically prevent such carryover.
- Use an appropriate biosensor type from ForteBio's list of off-the-shelf biosensors. If one is not readily available, you may prepare a custom biosensor by loading your ligand on any one of the biosensors available.
- Use the drop holder for all Create Standard Curve runs. Use 4 μ L of standard samples.
- Start with the default Run time value of 120 seconds. You may then optimize the Run time as needed. Typically, longer run times allow more binding of analyte to the biosensor, thus providing greater signal.
- Keep the Shaker in **Enable** mode to allow good mixing of sample. You may tweak shaker rate in **File > Options** to a value between 1000 and 2600 rpm to suit your needs. Typically, increasing shaker rate leads to faster binding rates.
- Always input experiment info first and click **Next**. Follow the onscreen prompts to load sample and biosensor on the BLItz system.
- To fit a curve to the data, choose from the available options. When in doubt about which curve fit is best, try all three options.

The purpose of this experiment is for using samples of known concentration to determine the signal the system shows for each concentration.

This experiment should display a signal vs. concentration graph, which is the standard graph. This is essential for using with the Quantitate Sample experiment. When you perform a Quantitate Sample experiment, you should have already performed the standard curve experiment. If you do not do the Quantitate Sample experiment first, the system will generate some data, but you will not know the details of the data, which is necessary to use for the Quantitate Sample experiment. In this way, you can make assumptions from the concentration provided.


In this **Standard Curve** experiment, hlgG is the analyte you will use to prepare four (4) standards ranging from 15.6–1000 $\mu\text{g}/\text{mL}$. With sample diluent as the reference, there will be a total of five (5) samples measured on the BLItz system.

Everything you need:


- Protein A biosensors
- hlgG stock (10 mg/mL)
- sample diluent
- 96-well, black flat-bottom microplate


To run a **Create a Standard Curve** experiment:

1. Click the **Create Standard Curve** icon  to display the **Create Standard Curve** module (Figure 4-17).



Create Standard Curve

Have you hydrated your biosensor?
Hydrate biosensors prior to use for a minimum of 10 minutes. Hydration time (min): 10 

Enter your Experiment Settings
Experiment Name:
Description: 


Enter your Run Settings
Biosensor type: Protein A Run time (s): 120
Sample ID: Shaker: Enable
Concentration: $\mu\text{g}/\text{mL}$ 

Figure 4-17: Create Standard Curve Module

Preparation

2. Hydrate five (5) biosensors, as described in “Hydrating Biosensors” on page 33.
3. Prepare the hlgG test samples as follows:
 - 1000 $\mu\text{g}/\text{mL}$ hlgG: 20 μL of stock hlgG + 180 μL of sample diluent. Label this as “**Stock A**”.
 - Using **Stock A**, perform the following dilutions:
 - Step 1—4-fold dilution: 10 μL Stock A + 30 μL sample diluent for a [hlgG] = 40 μL of 250 $\mu\text{g}/\text{mL}$
 - Step 2—4-fold dilution: 10 μL 250 $\mu\text{g}/\text{mL}$ + 30 μL sample diluent for a [hlgG] = 40 μL of 62.5 $\mu\text{g}/\text{mL}$
 - Step 3—4-fold dilution: 10 μL 62.5 $\mu\text{g}/\text{mL}$ + 30 μL sample diluent for a [hlgG] = 40 μL of 15.6 $\mu\text{g}/\text{mL}$



NOTE: hIgG samples should be run in this order:

Run 1: 0 $\mu\text{g}/\text{mL}$ (Reference)

Run 2: 15.6 $\mu\text{g}/\text{mL}$

Run 3: 62.5 $\mu\text{g}/\text{mL}$

Run 4: 250 $\mu\text{g}/\text{mL}$

Run 5: 1000 $\mu\text{g}/\text{mL}$ (Stock A)

Experiment

4. Input **Run Settings** per Figure 4-18 and click  .




Figure 4-18: Entering Run Settings

5. Follow the onscreen instructions (Figure 4-19).



Figure 4-19: User Action Required—Experiment Operations

6. Add 4 μL of the sample diluent for the first run into the drop holder (Figure 4-20).



Figure 4-20: Adding 4 μL of Sample Diluent into Drop Holder

7. Load the Protein A biosensor on to the biosensor mount; give it a final quarter-turn upward twist to ensure a snug fit (Figure 4-21).



Figure 4-21: Mounting a Biosensor

8. Move the slider to the right to move the drop holder to the read position (Figure 4-22).



Figure 4-22: Moving Drop Holder to Read Position

9. Close the BLItz system cover to start data acquisition each time the BLItz Pro software prompts you (Figure 4-23).



Figure 4-23: Closing the BLItz System Cover

Data is collected in real-time on screen (Figure 4-24). The sensorgram shows the results after the first run.




Figure 4-24: Data for Run 1 Collected in Real-Time—Binding Signal Measured in nanometers (nm) as a Function of Time (seconds)

After the run, a **Save As** dialog box displays and prompts you to save the experiment file (Figure 4-25).




Figure 4-25: Saving Dataset File

10. Enter a file name and **Save** the dataset.
11. Lift the cover, discard the used biosensor, wipe the drop holder with a swab or Kimwipe (Figure 4-26), and rinse 3X with sample diluent—wiping each time with a swab or Kimwipe.



Figure 4-26: Wiping the Drop Holder with a Kimwipe

12. Perform subsequent sample runs as follows, repeating steps 4 to 11 for each sample:
 - Run 2: 15.6 µg/mL
 - Run 3: 62.5 µg/mL
 - Run 4: 250 µg/mL
 - Run 5: 1000 µg/mL

Data

13. Inspect the Run data. The sensorgram (Figure 4-27) displays the results after five (5) runs.




Figure 4-27: Inspecting the Run Data—Binding Signal Measured in nanometers (nm) as a Function of Time (seconds)

14. Click **Save** to save your experiment.

Run List Report point (s): 0

| Index | Show | Analyze | Ref. | Sample ID | Conc. (µg/ml) | Binding rate (nm/s) | Run time (s) | Shaker speed (rpm) | Biosensor Type | Information |
|-------|-------------------------------------|-------------------------------------|-------------------------------------|------------------|---------------|---------------------|--------------|--------------------|----------------|-------------|
| 1 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | Negative Control | 0 | | 30 | 2200 | Protein A | |
| 2 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Unknown | | 0.5029 | 30 | 2200 | Protein A | |
| 3 | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Positive Control | 1000 | 3.393 | 30 | 2200 | Protein A | |

Figure 4-28: Saving Your Experiment

- **Report point (s)**—Specifies a time point at which to report signal for all samples in the graph.
 - **Add**—Adds a report point to the table.
 - **Remove All**—Removes all the report points from the table.
- **Export to CSV**—Exports the data in the graph to .csv files, one for each sample.
- **Create Report**—Creates a report of the experiment in .PDF file format.
- **Ref**—Select one of the samples as reference to subtract from all sample data. Selecting more than one sample as reference takes an average of the references to subtract from all data.
- **Binding rate (nm/s)**—Rate of change of binding signal with time; slope of the early part of the binding curve for each sample.
- **Shaker speed (rpm)**—Shaking speed of sample; this value is set in **File > Options**.

15. Click **Save Standard Curve** (Figure 4-29) to save the Standard Curve.



NOTE: This action saves the generated Standard Curve, rather than the experiment file itself.




Figure 4-29: Save Standard Curve— Binding Rate Measured in nanometers per second (nm/s) as a Function of Sample Concentration

- The drop-down menu provides curve fit options of (a) Linear point-to-point:
 - (b) Linear
 - (c) 5-PL unweighted
- **R² value** reports goodness of fit.
- **Save Standards**—Saves the standard curve as a .fsc file.

Result

Your standard curve has been generated. You are now ready to perform quantitation experiments.

RUNNING A QUANTITATE SAMPLE EXPERIMENT

When using the starter kit to perform a Quantitate Sample experiment, follow the instructions provided below for human IgG samples binding Protein A biosensors.

To measure concentration of proteins in your samples on the BLItz system, consider these guidelines:

- Create a Standard Curve that is representative of your unknown samples first before testing your unknowns. You may choose to perform the standard curve experiments after the Quantitate Sample experiments, but you will not be able to obtain concentrations of your unknowns until the standard curve is available.
- Ensure that the sample milieu in the unknowns and the standards are the same.
- Use the same experiment parameters such as Run time and Shaker speed for the Quantitate Sample experiment as that used in the creating the Standard Curve.
- Use an appropriate biosensor type from ForteBio's list of off-the-shelf biosensors. If one is not readily available, you may prepare a custom biosensor by loading your ligand on any one of the biosensors available.
- Use the drop holder for all Quantitate Sample measurements. Use 4 μ L of sample.
- Always input experiment info first and click **Next**. Follow the onscreen prompts to load sample and biosensor on the BLItz system.

This module allows you to quantitate unknown samples against a standard curve. This experiment assumes that you have already created a Standard Curve (see "Creating a Standard Curve" on page 41).

In this module, you will prepare two (2) samples of hIgG of known concentrations, and use the BLItz system to calculate their concentrations and confirm that the system reliably measures protein concentrations.

Everything you need:

- Protein A biosensors
- hIgG stock (10 mg/mL)
- sample diluent
- 96-well, black flat-bottom microplate

To run a Quantitate Sample experiment:

1. Click the **Quantitate Sample** icon  to display the **Quantitate Sample** module (Figure 4-30).

Quantitate Sample

Have you hydrated your biosensor?
Hydrate biosensors prior to use for a minimum of 10 minutes. Hydration time (min): 10

Enter your Experiment Settings

Experiment Name:
Description:

Enter your Run Settings

Biosensor type: Protein A Run time (s): 120
Sample ID: Shaker: Enable
Concentration units: µg/ml

Figure 4-30: Quantitate Sample Module

Preparation

- Hydrate two (2) biosensors, as described in “Hydrating Biosensors” on page 33.
- Prepare hIgG “unknowns” as follows:

Using the Stock A (1000 µg/mL) that you created to perform the standard curve experiment, perform the following dilutions:

- 5-fold dilution: 8 µL Stock A + 32 µL sample diluent for a [hIgG] = 200 µg/mL
- 4-fold dilution: 10 µL 200 µg/mL + 30 µL sample diluent for a [hIgG] = 50 µg/mL



NOTE: hIgG samples should be run in this order:

Run 1: 50 µg/mL

Run 2: 200 µg/mL

Quantitate Samples

- Input **Run Settings** per Figure 4-31 and click .

Enter your Experiment Settings

Experiment Name: Quantitate Samples Demo
Description: Quantitation of mock "unknown" human IgG samples

Enter your Run Settings

Sensor type: Protein A Run time (s): 30 Use tube
Sample ID: Unknown1 Shaker: Enable
Concentration units: µg/ml

Figure 4-31: Entering Run Settings

5. Follow the onscreen instructions (Figure 4-32).




Figure 4-32: User Action Required—Experiment Operation

6. Add 4 μL of the sample diluent for the first run into the drop holder (Figure 4-33).



Figure 4-33: Adding 4 μL of Sample Diluent into Drop Holder

7. Load the Protein A biosensor on to the biosensor mount; give it a final quarter-turn upward twist to ensure a snug fit (Figure 4-34).



Figure 4-34: Mounting a Biosensor

8. Move the slider to the right to move the drop holder to the read position (Figure 4-35).



Figure 4-35: Moving Drop Holder to Read Position

9. Close the BLItz system cover to start data acquisition each time the BLItz Pro software prompts you (Figure 4-36).



Figure 4-36: Closing the BLItz System Cover

Data is collected in real-time on screen. The sensorgram (Figure 4-37) displays the results after the first run.




Figure 4-37: Data Collected in Real-Time for Run 1—Binding signal measured in nanometers (nm) as a Function of Time (seconds)

After the run, a **Save As** dialog box displays and prompts you to save the experiment file (Figure 4-38).




Figure 4-38: Saving Dataset File

10. Enter a file name and **Save** the dataset.
11. Lift the cover, discard the used biosensor, wipe the drop holder with a swab or Kimwipe (Figure 4-39), and rinse 3X with sample diluent—wiping each time with a swab or Kimwipe.




Figure 4-39: Wiping the Drop Holder with a Kimwipe

12. Perform subsequent “unknown” runs:

Run 2: 200 $\mu\text{g}/\text{mL}$

Select Standard Curve

13. Click the **Load Standards** button on the Standard Curve graph (Figure 4-40) and select your standard curve (.fsc file) to load. You must load the standard curve generated using the protocol described in this guide.




Figure 4-40: Loading Standard Curve—Binding Signal Measured in nanometers per second (nm/s) as a Function of Concentration

- The drop-down menu provides curve fit options of (a) Linear point-to-point:
 - (b) Linear
 - (c) 5-PL unweighted
- **R² value** reports goodness of fit.
- **Load Standards**—Loads an existing standard curve as a .fsc file.
- **Save Standards**—Saves the standard curve as a .fsc file.
- **Clear**—Clears the loaded standard curve from the graph.

Data

14. Inspect the interpolated concentrations. Unknown1 should be about 50 $\mu\text{g}/\text{mL}$ and Unknown2 should be about 200 $\mu\text{g}/\text{mL}$.

Figure 4-41 displays the interpolated unknown concentrations.

15. Click **Save** to save your experiment.

| | | Report point (s): 0 | | Add | Remove All | | | |
|--------------------------|------|---------------------|---------------|---------------------|--------------|--------------------|----------------|-------------|
| ze | Ref. | Sample ID | Conc. (µg/ml) | Binding rate (nm/s) | Run time (s) | Shaker speed (rpm) | Biosensor Type | Information |
| <input type="checkbox"/> | | Unknown1 | 55.4 | 0.2374 | 30 | 2200 | Protein A | |
| <input type="checkbox"/> | | Unknown2 | 220.5 | 0.8324 | 30 | 2200 | Protein A | |

Create Report...

Export to CSV... Create Report...

Figure 4-41: Interpolated Unknown Concentrations in Table

- **Report point (s)**—Specifies a time point at which to report signal for all samples in the graph.
 - **Add**—Adds a report point to the table.
 - **Remove All**—Removes all the report points from the table.
- **Export to CSV**—Exports the data in the graph to .csv files, one for each sample
- **Create Report**—Creates a report of the experiment in .PDF file format
- **Ref**—Select one of the samples as reference to subtract from all sample data. Selecting more than one sample as reference takes an average of the references to subtract from all data.
- **Binding rate (nm/s)**—Rate of binding signal change with time; slope of the early part of the binding curve for each sample.
- **Shaker speed (rpm)**—Shaking speed of sample; this value is set in **File > Options**.

Result

In just a few minutes, you have measured the specific concentration of protein in your unknown samples.

RUNNING A BASIC KINETICS EXPERIMENT

When using the starter kit to perform a Basic Kinetics experiment, follow the instructions provided below for mouse IgG samples binding Protein A biosensors.

To measure kinetic parameters for your ligand-analyte interactions on the BLItz system consider these guidelines:

- Ligand indicates the binding partner loaded on the biosensor. Analyte refers to the binding partner present in solution.
- The Basic Kinetics module is appropriate for use when the ligand has previously been loaded on the biosensor, either as provided by ForteBio, or loaded offline on the bench top, or loaded on BLItz system in a previous experiment.
- Perform baseline and dissociation steps in an eppendorf tube. Perform association steps in either drop holder or tube for Duration of 5 minutes or less. When association step duration is more than five minutes, perform step in eppendorf tube.
- Use 250 μL of reagent in eppendorf tube and 4 μL of reagent in drop holder.
- When optimizing experiment conditions for measuring kinetic constants for a binding interaction, use fewer analyte concentrations distributed over the range of $100\times K_D$ to $0.01\times K_D$. After optimizing experiment conditions, populate more concentrations to maximize reliability of the measured kinetic constants.
- Input concentration and molecular weight values, or, the molar concentration values for the analyte samples used. The BLItz Pro software needs these values to calculate kinetic constants.
- Keep the shaker in **Enable** mode for all kinetics experiments. In **File > Options**, leave the shaker speed at the default setting of 2200 rpm for all kinetics experiments.
- Use an appropriate biosensor type from ForteBio's list of off-the-shelf biosensors. If one is not readily available, you can prepare a custom biosensor by loading your ligand on any one of the biosensors available.
- For the baseline step, duration of 30 to 60 seconds is typically recommended. You may modify this time as desired to achieve a stable Baseline before moving to the Association step.
- The duration for Association and Dissociation steps should be optimized based on the concentration of your analyte samples and the expected affinity of the interaction. Typically, high affinity interactions can require longer dissociation times.
- Always include a reference sample in your experiment to correct for background signal.
- Always input experiment info first and click **Next**. Follow the onscreen prompts to add reagents to the appropriate sample locations and to mount biosensor on the BLItz system.


To fit a curve to the data, choose from the available options. To learn more about the differences between local and global fitting, consult available literature or contact BLItz technical support.


In this Basic Kinetics experiment, you will develop a series of real-time data graphs for multiple concentrations of mouse IgG analyte binding to pre-loaded Protein A ligand on the Protein A biosensors. You will analyze the collected data to measure kinetic constants (k_a , k_d , K_D). You will use 4 μ L of sample on the drop holder for the association step and 250 μ L of sample diluent in the tube for baseline and dissociation.

Everything you need:


- Protein A biosensors
- mIgG stock (10 mg/mL)
- Sample diluent
- 96-well, black flat-bottom microplate
- Eppendorf tubes


To run a Basic Kinetics experiment:

1. Click the **Basic Kinetics** icon  to display the default **Basic Kinetics** module (Figure 4-42).



Basic Kinetics

Have you hydrated your biosensor?
Hydrate biosensors prior to use for a minimum of 10 minutes. Hydration time (min): 10 

Enter your Experiment Settings
Experiment Name:
Description: 

Enter your Run Settings
Sample ID: Concentration: μ g/ml
Biosensor type: SA (Streptavidin) Molecular Weight: kD
Shaker: Enable Molar Concentration: nM

| | Step Type | Duration (s) | Position |
|---|--------------|--------------|----------|
| 1 | Baseline | 30 | Tube |
| 2 | Association | 120 | Drop |
| 3 | Dissociation | 120 | Tube |


Run duration: > 04:30 

Figure 4-42: Basic Kinetics Module

Preparation

2. Hydrate five (5) biosensors, as described in “Hydrating Biosensors” on page 33.

3. Prepare mIgG sample dilutions:

Using the mIgG stock (10 mg/mL), perform the following dilutions:

- Step 1—20-fold dilution: 4 μ L mIgG stock + 76 μ L sample diluent for a [mIgG] = 0.5 mg/mL. Label this as “**Stock A**”.
- Step 2—50-fold dilution: 4 μ L Stock A + 196 μ L sample diluent for a [mIgG] = 10 μ g/mL
- Step 3— 2-fold dilution: 50 μ L of 10 μ g/mL + 50 μ L sample diluent for a [mIgG] = 5 μ g/mL
- Step 4—2-fold dilution: 50 μ L of 5 μ g/mL + 50 μ L sample diluent for a [mIgG] = 2.5 μ g/mL
- Step 5—2-fold dilution: 50 μ L of 2.5 μ g/mL + 50 μ L sample diluent for a [mIgG] = 1.25 μ g/mL



NOTE: mIgG samples will be run in this order:

Run 1: sample blank

Run 2: 1.25 μ g/mL

Run 3: 2.5 μ g/mL


Run 4: 5 μ g/mL

Run 5: 10 μ g/mL

Run Kinetics Experiment

4. Input **Run Settings** per Figure 4-43 and click .

Enter your Experiment Settings

Experiment Name: 

Description:

Enter your Run Settings

Sample ID: Concentration: μ g/ml

Biosensor type: Molecular Weight: kD

Shaker: Molar Concentration: nM

| Step | Step Type | Duration (s) | Position |
|------|--------------|-------------------------------------------------------------------|----------|
| 1 | Baseline | 30 | Tube |
| 2 | Association | <input style="border: 1px solid black;" type="text" value="300"/> | Drop |
| 3 | Dissociation | 300 | Tube |


Run duration: > 10:30  Run 1

Figure 4-43: Entering Run Settings

5. Follow the onscreen instructions (Figure 4-44).



Figure 4-44: User Action Required—Experiment Operation

6. Ensure that the sample diluent (4 μL in the drop holder) or (250 μL in the tube) are available for reading (Figure 4-45).




Figure 4-45: Perform Baseline and Dissociation Measurements in Tube and Association Measurement in Drop Holder

7. Load the Protein A biosensor on to the biosensor mount; give it a final quarter-turn upward twist to ensure a snug fit (Figure 4-46).



Figure 4-46: Mounting a Biosensor

8. Move the slider to the appropriate position to bring either the tube or drop holder to the read position, as appropriate for each step (Figure 4-47):
 - a. left position for **Baseline**
 - b. right position for **Association**
 - c. left position for **Dissociation**



Figure 4-47: Moving Drop Holder or Tube to Read Position

9. Close the BLItz system cover to start data acquisition each time the BLItz Pro software prompts you (Figure 4-48).



Figure 4-48: Closing the BLItz System Cover

Data is collected in real-time on screen. The sensorgram (Figure 4-49) displays the results after the first run.




Figure 4-49: Data for Run 1 Collected in Real-Time

After the run, a **Save As** dialog box displays and prompts you to save the experiment file (Figure 4-50).




Figure 4-50: Saving Dataset File

10. Enter a file name and **Save** the dataset.

11. Lift the cover, discard the used biosensor, wipe the drop holder with a swab or Kimwipe (Figure 4-51), and rinse 3X with sample diluent—wiping each time with a swab or Kimwipe.




Figure 4-51: Wiping the Drop Holder with a Kimwipe

12. Perform subsequent sample runs, repeating steps 4 to 11 for each concentration:
 - Run 2: 1.25 $\mu\text{g}/\text{mL}$
 - Run 3: 2.5 $\mu\text{g}/\text{mL}$
 - Run 4: 5 $\mu\text{g}/\text{mL}$
 - Run 5: 10 $\mu\text{g}/\text{mL}$

Data

13. Observe the Run data and save.
The sensorgram displays the results after five (5) runs (Figure 4-52).




Figure 4-52: Observing Run Data

14. Select the settings in Figure 4-53 and click **Analyze** to display the results in Figure 4-54.




Figure 4-53: Analysis Settings

The step correction feature corrects misalignment between two steps due to system artifacts. The association step can be aligned to the dissociation step or to the baseline.

- **Start of association**—Moves the association step on the Y axis to align the beginning of the association step with the end of the adjacent baseline step.
- **Start of dissociation**—Moves the dissociation step on the Y axis to align the end of the adjacent association step with the beginning of the dissociation step.




Figure 4-54: Analysis Results—Red Line Overlaying Data Are Curve Fits

You have just generated binding curves for mIgG against Protein A and a table of kinetic constants (Figure 4-55).

| Index | Sample ID | Conc. (nM) | Information | KD (M) | ka (1/Ms) | ka Error | kd (1/s) | kd Error | Rmax | Rmax Error | R equilibrium |
|-------|------------|------------|-------------|---------|-----------|----------|----------|----------|-------|------------|---------------|
| 1 | blank | 0 | | | | | | | | | |
| 2 | 2.5 ug/mL | 16.67 | | 3.56e-9 | 1.235e5 | 3.467e2 | 4.398e-4 | 2.539e-6 | 4.46 | 0.01916 | 3.675 |
| 3 | 5 ug/mL | 33.33 | | 3.56e-9 | 1.235e5 | 3.467e2 | 4.398e-4 | 2.539e-6 | 4.577 | 0.01420 | 4.135 |
| 4 | 10 ug/mL | 66.67 | | 3.56e-9 | 1.235e5 | 3.467e2 | 4.398e-4 | 2.539e-6 | 4.622 | 0.007895 | 4.387 |
| 5 | 1.25 ug/mL | 8.333 | | 3.56e-9 | 1.235e5 | 3.467e2 | 4.398e-4 | 2.539e-6 | 4.128 | 0.02068 | 2.892 |

Figure 4-55: Binding Curve Results for mIgG against Protein A and Kinetic Constants

- **Conc.**—Concentration of analyte
- **K_D (M)**—Measured affinity of interaction; affinity constant in Molar
- **k_a (1/Ms)**—Association rate constant
- **k_a Error**—Calculated standard error in the measured association rate constant
- **k_d (1/s)**—Dissociation rate constant
- **k_d (Error)**—Calculated error in the measured dissociation rate constant
- **Rmax**—Maximum response determined from the fit of the binding data
- **Rmax Error**—Calculated standard error in Rmax
- **Requilibrium**—Calculated response at equilibrium that is determined from a fit of the binding data.

Result

You have now generated rate constants and affinity constant for the mouse IgG2-Protein A binding interaction. Congratulations!

With this, your experiment is complete. We hope that BLItz's powerful experiment features combined with its simplicity is what you are looking for in your work.

RUNNING AN ADVANCED KINETICS EXPERIMENT

The Advanced Kinetics module can be used to perform kinetic analysis of ligand-analyte interactions where the ligand is first loaded on biosensors on the BLItz system followed by monitoring the association and dissociation of analyte to the ligand. In addition, the Advanced Kinetics module can be used to set up custom binding assays involving multiple binding steps, such as the formation of a sandwich complex (capture antibody-antigen-secondary antibody).

The **Initial Baseline** (Figure 4-56) can typically be of 30 or 60 seconds duration, or until stable baseline is reached in the medium you choose for your experiment. The initial baseline is performed in the tube. The loading step allows you to load your ligand on the biosensor. The loading step can be performed on different biosensors listed in the Biosensor type field.




Figure 4-56: Advanced Kinetics—Step Type List

Table 4-1 lists the type of ligand required for loading to the biosensor.

Table 4-1: Required Ligands

| Biosensor Type | Required Ligand |
|---------------------------------------------|------------------------------------------|
| Streptavidin | Biotin-tagged proteins, oligos, peptides |
| Super Streptavidin | Biotin-tagged proteins, oligos, peptides |
| Anti-GST | GST-fused proteins |
| Anti-hIgG Fc Capture (AHC) | Human IgG, human Fc fusion proteins |
| Anti-mIgG Fc Capture (AMC) | Mouse IgG, mouse Fc fusion proteins |
| Amine reactive 2nd generation (AR2G) | Proteins, oligos, peptides |
| Amine reactive (AR) | Proteins, oligos, peptides |
| Aminopropylsilane (APS) | Proteins, peptides |

To prepare biotin-tagged proteins, follow the minimum biotinylation ratio method described in *Technical Note 28, Biotinylation of Protein for Immobilization onto Streptavidin Biosensors*, available at www.blitzmenow.com/literature.html.

To load ligand on AR and AR2G biosensors, the Amine Coupling Second Generation Reagent Kit is required.

A Baseline step follows the loading step to establish a new baseline before association and dissociation steps with the analyte. In the Association step, the analyte in solution binds to the ligand on the biosensor. In the Dissociation step, the analyte bound to the biosensor dissociates into solution.

Table 4-2 lists the applicable rules when performing an Advanced Kinetics experiment.

Table 4-2: General Step Rules for an Advanced Kinetics Experiment

| <i>Step Type</i> | <i>Position</i> |
|-------------------------|-----------------------------------------------------------------------------------------------|
| Initial Baseline | Tube |
| Loading | Drop holder or Tube if Duration is less than 300 seconds; only Tube if more than 300 seconds. |
| Baseline | Tube |
| Association | Drop holder or Tube if Duration is less than 300 seconds; only Tube if more than 300 seconds |
| Dissociation | Tube |
| Custom | Drop holder or Tube if Duration is less than 300 seconds; only Tube if more than 300 seconds. |



NOTE: In Advanced Kinetics, to create a successful experiment, there must (at a minimum) be **Initial Baseline, Association, and Dissociation** step, included in that order.

To run an Advanced Kinetics experiment, you will need:


- Biosensors (choice of biosensor made by user)
- Ligand for loading on biosensor
- Analyte or sample (one or more concentrations, depending on type of experiment)
- Sample diluent
- 96-well, black, polypropylene flat-bottom microplate
- Eppendorf tubes


To run an Advanced Kinetics experiment:

1. With the software launched, on the left pane, click the **Advanced Kinetics** icon to display the default **Advanced Kinetics** module (Figure 4-57).



Advanced Kinetics

Have you hydrated your biosensor?
Hydrate biosensors prior to use for a minimum of 10 minutes. Hydration time (min):  Hydrate Biosensor

Enter your Experiment Settings  New Experiment

Experiment Name:

Description:

Enter your Run Settings

Sample ID: Concentration: $\mu\text{g/ml}$

Biosensor type: SA (Streptavidin) Molecular Weight: kD

Shaker: Molar Concentration: nM

| | Step Type | Duration (s) | Position |
|---|------------------|--------------|----------|
| 1 | Initial Baseline | 30 | Tube |
| 2 | Loading | 120 | Drop |
| 3 | Baseline | 30 | Tube |
| 4 | Association | 120 | Drop |
| 5 | Dissociation | 120 | Tube |


Run duration: > 07:00 

Figure 4-57: Advanced Kinetics Module

Preparation

- Hydrate biosensors, as described in the “Hydrating Biosensors” on page 33.
- Prepare ligand solution for loading step. Perform dilution as required to bring the ligand solution to a concentration of 10–50 $\mu\text{g/mL}$.



NOTE: When performing kinetics experiments on AR and AR2G biosensors, two separate Advanced Kinetics experiments are needed. The first Advanced Kinetics experiment should be performed according to the following protocol to load protein on the biosensors:

- Set up an experiment with the steps in Figure 4-58. In Step 2, activation of the hydrated biosensor with EDC/Sulfo NHS mixture is performed in a tube. In Step 3, the ligand is loaded on the activated biosensor either in the drop holder or the tube. In Step 5, the biosensor is quenched with quenching solution.

| | Step Type | Duration (s) | Position |
|---|------------------|--------------|----------|
| 1 | Initial Baseline | 15 | Tube |
| 2 | Custom | 300 | Tube |
| 3 | Loading | 300 | Tube |
| 4 | Custom | 300 | Tube |
| 5 | Baseline | 60 | Tube |


Figure 4-58: Step Types for This Experiment

All biosensors needed for the experiment should be loaded one after another using this procedure and kept hydrated in solution until ready for the binding experiment with analyte.


- b. The analyte binding kinetics experiment can now be performed using the Basic Kinetics module.
4. Prepare analyte solutions.


Run Experiment

5. Input **Run Settings** per Figure 4-59 and click .



Advanced Kinetics

Have you hydrated your biosensor?
Hydrate biosensors prior to use for a minimum of 10 minutes. Hydration time (min): 

Enter your Experiment Settings 
Experiment Name:
Description:

Enter your Run Settings

Sample ID: Concentration:
Biosensor type: Molecular Weight:
Shaker: Molar Concentration:

| Step | Step Type | Duration (s) | Position |
|------|------------------|----------------------------------|----------|
| 1 | Initial Baseline | 30 | Tube |
| 2 | Loading | 300 | Drop |
| 3 | Baseline | 30 | Tube |
| 4 | Association | 300 | Drop |
| 5 | Dissociation | <input type="text" value="600"/> | Tube |


Run duration: > 21:00  Run 1

Figure 4-59: Entering Run Settings

6. Follow the onscreen instructions. For instructions on performing the functions prompted by the software, refer to “Basic Functions of the BLItz System” on page 10.
7. To clean the drop holder during the experiment, perform the following actions: Lift the BLItz system cover, discard the used biosensor, wipe the drop holder with a swab or Kimwipe (Figure 4-40), and rinse 3X with sample diluent—wiping each time with a swab or Kimwipe.
8. Repeat the experiment for all samples.
9. Save the experiment.

Data Analysis




Figure 4-60: Observing Run Data from an Advanced Kinetics Experiment

10. Select **Step Correction and Fitting (1:1)** options and click **Analyze** to display curve fits and results.

That is the general procedure to use to run Advanced Kinetics experiments.

APPENDIX A:

Running the BLitz Pro Software on Computers Running Octet Systems

| | |
|---------------------------------------------------------------------------|----|
| Running Only the BLitz System on the Computer | 68 |
| Running Experiments on Octet and BLitz Systems on the Same Computer | 68 |

For Octet™ customers who also purchase a BLItz™ system, this section addresses the procedures with having the two systems share a computer. Those Octet customers include:


- Those that have the Octet v4.x software
- Those that have a relatively new v5.x. and v6.x software
- Those that have the latest v7.0 software



NOTE: To have BLItz and Octet systems share a computer seamlessly, update your Octet software to v6.4.1.3 or v7.0.1.3. If you do not wish to update your Octet software, earlier versions of Octet software can still run on the same computer alongside BLItz Pro software. The following sections describe how to perform the procedure.

RUNNING ONLY THE BLITZ SYSTEM ON THE COMPUTER

To run ONLY the BLItz system on the computer:

1. Ensure that the Octet software is closed. If the Octet software is running, close it.
2. Connect the BLItz system power and plug its USB cable into the computer.
3. Turn on the BLItz system.
4. Double-click the **BLItz Pro 1.0** desktop icon .

The BLItz Pro software is launched, the BLItz system will initialize, and an “Initializing” prompt will display. The BLItz system’s front light will blink; after the initialization process completes successfully, the “Initializing” prompt will close automatically, and the light will stop blinking.

If no error message appears, you can start an experiment on the BLItz system; otherwise, refer to “Troubleshooting” on page 71.

RUNNING EXPERIMENTS ON OCTET AND BLITZ SYSTEMS ON THE SAME COMPUTER

To run experiments on both an Octet system and a BLItz system on the same computer:

1. Ensure that the Octet software is closed, that the BLItz system is powered off, and that the BLItz Pro software is closed.
2. Double-click the Octet software desktop icon to launch the Octet software.

The Octet software is launched, after the Octet system initializes properly. A “Ready. Date, Time” statement displays on the Instrument Status window (Figure A-1).






Figure A-1: Instrument Status Window—Octet System

3. Turn on the BLItz system.
4. Double-click the **BLItz Pro 1.0** desktop icon .

The BLItz Pro software is launched, the BLItz system will initialize, and an “Initializing” prompt will display. The BLItz system’s front light will blink; after the initialization process completes successfully, the “Initializing” prompt will close automatically, and the light will stop blinking.

If no error message appears, you can start an experiment on the BLItz system; otherwise, refer to “Troubleshooting” on page 71.



CAUTION: While running an experiment on an Octet system, do not click the Stop icon  to abort the experiment. If this occurs: 1) close the Octet software and the BLItz Pro software, 2) power off the BLItz system, and 3) then re-launch the Octet software.



CAUTION: On the Octet software, do not click **Instrument > Reset** menu command to reset the Octet system. If an experiment was aborted on the Octet software: 1) close the Octet software and the BLItz Pro software, 2) power off the BLItz system, and 3) then re-launch the Octet software.

APPENDIX B:

Troubleshooting

| | |
|------------------------------------------------------------------|----|
| Troubleshooting BLitz System Errors with an Octet System | 72 |
| Scenario 1—The BLitz System Cannot Initialize Successfully | 72 |
| Scenario 2—The BLitz System Is Not Recognized | 72 |

This appendix describes how to troubleshoot some scenarios of BLitz system errors.

TROUBLESHOOTING BLITZ SYSTEM ERRORS WITH AN OCTET SYSTEM

Scenario 1—The BLItz System Cannot Initialize Successfully

When you launch the BLItz Pro™ software, the BLItz™ system cannot initialize successfully (the BLItz system light keeps blinking), and either of the following error messages display (Figure B-1).




Figure B-1: Error Message Possibilities

How to Recover

1. Close the BLItz Pro software.
2. Switch off the BLItz system.
3. Ensure that the USB cable is connected properly between the BLItz system and the computer.
4. Switch on the BLItz system. Wait for five (5) minutes.
5. Turn on the BLItz Pro software.



NOTE: If the symptom still exists, restart steps 1 to 5 above.

Scenario 2—The BLItz System Is Not Recognized

When you launch BLItz Pro software, the following prompt (Figure B-2) displays.




Figure B-2: Simulation Mode Prompt

How to Recover

1. Close the BLItz Pro software.
2. Ensure that the BLItz system is powered on (the light is on) and that the power and USB cable are properly connected to the BLItz system and computer.
3. Re-launch the BLItz Pro software until the system initialization completes successfully (no error message displays).

Index

Numerics

96-well, black, flat-bottom polypropylene 3

A

About BLItz (menu) 27

accessories, part of the BLItz system kit 3

active method file, closing 22

Advanced Kinetics (menu) 26

Advanced Kinetics experiment, running 62

Advanced Kinetics icon 62

Aminopropylsilane (APS) biosensor type 14

analytes measured, biosensor types 13

Anti-GST (GST) biosensor type 14

Anti-hlgG-Fc biosensor type 13

Anti-hlgG-Fc Capture (AHC) biosensor type 14

Anti-mlgG-Fc biosensor type 13

Anti-mlgG-Fc Capture (AMC) biosensor type 14

Anti-Penta-HIS (HIS) biosensor type 13

Association

described 29

icon (figure) 29

B

Baseline

described 29

icon (figure) 29

basic functions of the BLItz system 10

Basic Kinetics (menu) 25

Basic Kinetics experiment, running 55

Basic Kinetics icon 55

binding assays 12

Binding rate (nm/s) menu 40

Biosensor Type, specification for the BLItz system 9

biosensor types

Amine Reactive (AR2G) Amine Reactive (AR2G) biosensor type 14

Aminopropylsilane (APS) 14

analytes measured, listed 13

Anti-GST (GST) 14

Anti-hlgG-Fc 13

Anti-hlgG-Fc Capture (AHC) 14

Anti-mlgG-Fc 13

Anti-mlgG-Fc Capture (AMC) 14

Anti-Penta-HIS (HIS) 13

listed (table) 13

Protein A 13

Protein G 13

Protein L 13

- required capture molecules, listed 13
 - Streptavidin (SA) 14
- biosensors, hydrating 33
- BLI (bio-layer interferometry) technology 15
- BLItz driver, installing 18
- BLItz Pro 1.0
 - described— 28
 - icon (figure) 28
- BLItz Pro software
 - closing 24
 - configuring 18
 - displaying properties and license information 27
 - installing 18
 - installing on Octet instruments 68
 - launching 19
 - Main Menu bar 21
 - Title bar 21
 - toolbars 21
 - uninstalling 30
- BLItz Pro Welcome screen (figure) 32
- BLItz support options
 - described 5
 - technical support 5
- BLItz system
 - basic functions 10
 - Biosensor Type specification 9
 - configuring 18
 - installing 18
 - package contents, listed (figure) 2
 - resetting 26
 - running only on the computer 68
 - specifications
 - Capabilities 8
 - Compliance 8
 - Connections 9
 - Dimensions 9
 - Electrical Requirements 9
 - Environmental 8
 - Equipment Classifications 8
 - Mechanics 9
 - Optics 9
 - Sample Types specification 9
 - Sampling Format 8
 - Sampling volume specification 9
 - Shake Rate 9
 - Software 9
 - Weight 9
 - starter kit contents
 - materials required 3
 - starter kit contents, accessories 3
 - system options, viewing 23

BLItz system components, front and back (figure) 8

BLItz Web Site (menu) 27

BLItzInstallCD.exe file 18

C

Capabilities, specification for the BLItz system 8

Close Method File (menu) 22

closing

- active method file 22
- software 24

Compliance, specification for the BLItz system 8

configuring

- BLItz Pro software 18
- BLItz system 18

connecting

- one side of the BLItz system to a USB port on your computer 18, 19
- one side of the data cord to the BLItz system 18, 19
- power cord to the BLItz system 18

Connections, specification for the BLItz system 9

contacting ForteBio technical support 5

Convert to Create Standard Curve (menu) 26

Convert to Quantitate Sample (menu) 26

Create Report (menu) 22

Create Report menu 40

Create Standard Curve (menu) 25

Create Standard Curve icon 42

Creating a Standard Curve experiment, running 42

D

data display, significant digits 23

destructive interference 15

Dimensions, specification for the BLItz system 9

Dip and Read biosensors 13

Dip and Read label-free assays 12

displaying

- BLItz software properties 27
- license information 27

Dissociation

- described 29
- icon (figure) 29

Drop Holder

- described 28
- icon (figure) 28

E

- Electrical hazard symbol 4
- Electrical Requirements, specification for the BLItz system 9
- environmental features, of the BLItz system 8
- Eppendorf tubes 3
- Equipment Classifications specification, for the BLItz system 8
- Exit (menu) 24
- Experiment menu
 - described 25
 - figure 25
 - list of menu commands 25
- experiment method files, saving 22
- experiments, setting up 10
- Export to CSV menu 40

F

- File menu
 - commands, listed (table) 22
 - described 22
- ForteBio technical support, contacting 5
- Fuse symbol 4

H

- Heat/hot symbol 4
- Help menu
 - commands, listed (table) 27
 - described 27
 - figure 27
- Hydrate Biosensor
 - described 28
 - icon (figure) 28
- hydrating biosensors 33
- Hydration Done icon (figure) 28

I

- Initial Baseline
 - described 28
 - icon (figure) 28
- Install BLItz Driver (menu) 18
- Install BLItz Pro Software (menu) 18
- installing
 - BLItz driver 18
 - BLItz Pro software in your lab 18
 - BLItz Pro software on Octet instruments 68
 - BLItz system in your lab 18
- Instrument menu
 - commands 26, 27

- commands, listed (table) 26
- described 26
- figure 26

Instrument Status window

- figure 69

Instrument Status window, logging 26

L

launching the BLItz Pro software 19

license information, displaying 27

lint-free lab wipe, wiping sample 10

load ligand on AR and AR2G biosensors 62

Load Standards button 52

Loading

- described 28
- icon (figure) 28

loading standard curve (figure) 52

logging in the Instrument Status window 26

M

Main Menu bar

- described 21
- figure 21

materials required, part of the BLItz system starter kit 3

Mechanics, specification for the BLItz system 9

method files

- closing 22
- saving 22

minimum requirements, for laptop or computer 3

mounting a biosensor (figure) 11

N

New Experiment

- described 28
- icon (figure) 28

Next

- described 28
- icon (figure) 28

O

Octet instruments, installing BLItz Pro software 68

Open Method File (menu) 22

optical thickness 15

Optics, specification for the BLItz system 9

Options (menu) 23

P

- pipettors, used for the experiments 3
- placing a drop of sample in the drop holder 10
- plugging the BLItz system into an outlet 18
- Protein A biosensor type 13
- Protein G biosensor type 13
- Protein L biosensor type 13

Q

- Quantitate Sample (menu) 25
- Quantitate Sample experiment, running 48
- Quantitate Sample icon 48
- Quick Yes/No (menu) 25
- Quick Yes/No experiment, running 36
- Quickstart Guide (menu) 27

R

- Ready
 - described 28
 - status (figure) 28
- Ref menu 40
- removing sample drop holder (figure) 10
- replacing sample drop holder (figure) 10
- Report point (s) 40

- required capture molecules, biosensor types 13
- requirements, minimum, for laptop or computer 3
- Reset (menu) 26
- resetting BLItz system 26
- running
 - Advanced Kinetics experiment 62
 - Basic Kinetics experiment 55
 - Creating a Standard Curve experiment 42
 - experiments on both an Octet system and a BLItz system on the same computer 68
 - only the BLItz system on the computer 68
 - Quantitate Sample experiment 48
 - Quick Yes/No experiment 36

S

- sample drop holder, removing and replacing 10
- sample plate shaker, stopper 26
- Sample Types specification, for the BLItz system 9
- Sampling Format, specification for the BLItz system 8
- Sampling Volume specification, BLItz system 9

- Save As dialog box 38
 - Save Method File (menu) 22
 - Save Method File As (menu) 22
 - saving experiment method files 22
 - setting up an experiment 10
 - Shake Rate, specification for the BLItz system 9
 - Shaker speed (rpm) menu 40
 - shifting position with slider 11
 - significant digits, data display 23
 - slider 11
 - Software, specification for the BLItz system 9
 - starter kit contents, listed 2
 - Status Bar
 - described 21
 - figure 21
 - Steptavidin (SA) biosensor type 14
 - Stop icon
 - described 28
 - Stop Shaker (menu) 26
 - stopping sample plate shaker 26
 - symbols
 - electrical hazard 4
 - fuse 4
 - heat/hot 4
 - used in this guide 4, 5
 - system options, viewing 23
- ## T
- technical support, contacting 5
 - Title bar
 - described 21
 - Title bar, BLItz Pro software 21
 - toolbars
 - available in the BLItz Pro software 21
 - troubleshooting BLItz system errors with an Octet system 71
 - Tube
 - described 28
 - icon (figure) 28
 - tubes, Eppendorf 3
- ## U
- uninstalling the BLItz Pro software 30
 - User Guide (menu) 27
- ## V
- viewing BLItz system options 23

W

wavelength shift 15

weight, specification for the BLItz system 9

Welcome (menu) 25

wiping sample with Kimwipe (figure) 10

