



# Table of Contents

- I. Introduction: Genomic Technologies for Cancer Research 3
- II. Approaches for Detecting Somatic Mutations 4
  - Targeted Sequencing Solutions for Somatic Mutation Detection 4
    - Exome Sequencing 4
    - Focused Sequencing Panels 4
    - Custom Targeted Sequencing 4
  - Whole-Genome Sequencing Solutions 4
  - Data Analysis Tools for Somatic Variant Detection 5
- III. Evaluating Germline Mutations in Cancer 6
  - Targeted Sequencing to Detect Common Germline Mutations 7
  - Microarray-Based Approaches 7
- IV. Structural Variant Detection in Cancer 7
  - DNA and RNA Sequencing for Translocation Detection 8
  - Copy Number Variation Arrays 8
- V. Investigating Gene Regulation in Cancer 8
  - DNA-Protein Interactions 8
  - DNA Methylation 9
  - RNA Sequencing 9
  - Targeted RNA Sequencing 9
  - Small RNA Sequencing 10
  - Data Analysis Tools for the Study of Gene Regulation 11
- VI. Summary 11









## Whole-Genome Sequencing Solutions

Through WGS, researchers can compare T/N sample pairs to identify somatic mutations in coding and noncoding regions across the entire genome. As a hypothesis-free approach, WGS is well suited for the discovery of novel driver mutations.

For researchers who prefer a T/N WGS service, the Illumina Genome Network™ (IGN) Cancer Analysis Service performs medium- to large-cohort T/N sequencing studies for researchers looking for a cost-effective solution. To learn more about IGN Cancer Analysis Services, visit [www.illumina.com/ign](http://www.illumina.com/ign).

## Data Analysis Tools for Somatic Variant Detection

Whether performing whole-genome, exome, or targeted sequencing, Illumina offers seamless workflow solutions that include library preparation kits, sequencing platforms, and data analysis software packages (Figure 2).



**Figure 2: Illumina Seamless Workflow Solutions**—Illumina sequencing solutions are fully integrated, DNA-to-data solutions from library preparation to final data analysis. Prepare T/N paired libraries with optimized kits for whole-genome, exome, transcriptome, or targeted sequencing. Perform sequencing using the MiniSeq™, MiSeq®, NextSeq®, or HiSeq® sequencing systems. Access Illumina data analysis tools for alignment, variant calling, T/N comparative analysis, or expression profiling.

Data analysis, the final step in the workflow, is essential to experimental success and must be tailored to the research question at hand. For example, it is important that the variant calling method properly model the complexities of multiple cancer subclones versus normal sample contamination. Illumina primary data analysis and somatic variant detection software solutions are easily accessible through on-instrument analysis software, such as MiSeq Reporter Software and Local Run Manager, or through BaseSpace® Sequence Hub, the Illumina cloud-based genomic computing environment. Data are streamed from the MiniSeq, MiSeq, NextSeq, or HiSeq Systems directly and seamlessly into BaseSpace Sequence Hub, which offers a suite of apps tailored to various data analysis needs. These tools are packaged into a user interface designed to be accessible to any user, regardless of bioinformatics experience. Simple prompts guide users through the entire process, starting from selecting the files generated by the sequencer to data filtration and analysis.

For somatic variant detection, the BaseSpace Tumor Normal App can be used to report single-nucleotide polymorphisms (SNPs), indels, copy number variants (CNVs), and structural variations found only within the tumor sample. For amplicon-based sequencing panels, MiSeq Reporter with Somatic Variant Caller, or the TruSeq® Amplicon BaseSpace App can be used. Depending on the panel used, Illumina also offers VariantStudio, an annotation and filtering software that allows customers to create customized reports of variant data. For a comprehensive view of Illumina informatics solutions for variant detection, visit [www.illumina.com/informatics.html](http://www.illumina.com/informatics.html).

For Research Use Only. Not for use in diagnostic procedures.

GGAATGATAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCCACTGATCAACGTACCCTAACGAAACGTATCAATTGAGACTAAATATTAACGTACCAATTAAGAGCTACCGTGAACGACGAAAGAATGATAACAGTAAACACACTTCTGTAAACCTH  
TCAAAGTACCGTAAACGAAACGTATCAATTGAGATTACTTGATCCCACTGATTCACAGTACCGTAAACGAAACGTATCAATTGAGACTAAATATTAACGTACCAATTAAGAGCTACCGTGAACGACGAAAGAATGATAACAGTAAACACACTTCTGTAAACCTH  
CGACGAAAAGAATGATAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCCACTGATTCACAGTACCGTAAAGATTACTTGATCCCACTGATTCACAGTACCGTAAACGAAACGTATCAATTGAGACTAAATATTAACGTACCAATTAAGAGCTACCG  
TACCGTACCAATTAAGAGCTAACCGTGAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCCACTGATTCACAGTACCGTAAACGAAACGTATCAATTGAGACTAAATATTAACGTACCAATTAAGAGCTACCGTGAACGACGAAAGAATGATAA  
GGAATGATAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCCACTGATTCACAGTACCGTAAACGAAACGTATCAATTGAGACTAAATATTAAGGATCAATTAGAGCTACCGTCTTTGTGTAAACCTTAAGATTACTTGATCCCACTGATTCACAG  
GATTACTTGATCCCACTGATTCACAGTAAAGATTACTTGATCCCACTGATTCACAGTACCGTAAACGAAACGTATCAATTGAGACTAAATATTAACGTACCAATTAAGAGCTACCGTGAACGAAAGAATGATAACAG  
CGTATCAATTGAGACTAAATATTAACGTACCAATTAAGATTCTGTAAACCTTAAGATTACTTGATCCCACTGATTCACAGTACCGTAAACGAAACGTATCAATTGAGACTAAATATTAACGTACCAATTAAGAGCTACCGTGAACGAAAGAATGATAACAG










## V. Investigating Gene Regulation in Cancer

NGS and microarray technology can be used to monitor changes in the transcriptome and epigenome (Figure 3). Protein–DNA interactions that play a role in cancer-related gene regulation can be assayed at the whole-genome level. Both sequencing and microarray options are available for researchers to detect these regulatory mechanisms (Table 4).

### Protein–DNA Interactions

ChIP–Seq provides hypothesis-free information on the regulation of gene expression and gives a complete snapshot of DNA-associated protein activity across the genome. Whole-genome ChIP–Seq offers a simple, cost-effective solution for obtaining visibility into the mechanics of protein-mediated gene regulation. Samples can be multiplexed for high-throughput processing, and deep sequencing enables detection of lower abundance protein–DNA interactions, such as those associated with transcription factor studies.



**Figure 3: NGS Transcriptomics Application**— Transcriptomics includes a spectrum of methods from total RNA–Seq to small RNA–Seq. Illumina provides library preparation kits for both coding and noncoding RNA sequencing applications such as mRNA–Seq, small RNA–Seq, and total RNA–Seq. Illumina also offers targeted RNA sequencing panels for apoptosis, NFκB pathway, p53 pathway, stem cells, the Wnt pathway, and more. Targeted sequencing panels can be custom designed for specific regions of interest.

### DNA Methylation

Cancer progression is often caused or influenced by epigenetic changes that alter gene expression. Assessing DNA methylation status can provide insight into the regulatory drivers of gene expression, affecting cancer development and progression. Sequencing-based methylation analysis enables researchers to identify and track methylation profiles across the entire genome, detecting variations in methylation signatures at single-base resolution. With methylation microarrays, researchers can quantitatively interrogate expert-selected methylation sites while balancing throughput capacity and affordability.

AAGATGATAACAGTAACACACTCTGTTAACCTTAAGATTACTTGATCCAAGTGATCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAGAGCTACCGTGCAACGACGAAAAGAATGATAACAGTAACACACTCTGTTAACCTTA  
TCAACGTACCGTAACGAAACGTATCAATTAGATTACTTGATCCAAGTGATCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAGAGCTACCGTGCAACGACGAAAAGAATGATAACAGTAACACACTCTGTTAACCTTA  
CGACGAAAGAATGATAACAGTAACACACTCTGTTAACCTTAAGATTACTTGATCCAAGTGATCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAGAGCTACCGTGCAACGACGAAAAGAATGATAAC  
ACGTACCATTAAGAGCTACCGTGCAACAGTAACACACTCTGTTAACCTTAAGATTACTTGATCCAAGTGATCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAGAGCTACCGTGCAACGACGAAAAGAATGATAAC  
AGATTGATAACAGTAACACACTCTGTTAACCTTAAGATTACTTGATCCAAGTGATCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAGAGCTACCGTGCAACGACGAAAAGAATGATAACAGT  
GATTACTTGATCCAAGTAACAGTTAAGATTACTTGATCCAAGTGATCAACGTACCGTAACGAACGTATCAATTGAGCTTCTGTTAACCTTAAGATTACTTGATCCAAGTGATCAACGTACCGTAACGAACGTATCAATTGAGACTAGCAACGACG  
CGTATCAATTGAGACTAAATATTAACGTACCATTAAGAGCTCTGTTAACCTTAAGATTACTTGATCCAAGTGATCAACGTACCGTAACGAACGTATCAATTGAGACTAAATATTAACGTACCATTAAGAGCTACCGTGCAACGAAAAGAATGATAACAG





