

HCS 2.4/3.0. The countersinkable compression screw.

Surgical Technique



This publication is not intended for distribution in the USA.

Instruments and implants approved by the AO Foundation.

Table of Contents

Introduction	Features and Benefits	2
	Functional Principle	3
	Indications	4
Surgical Technique	Hand – Scaphoid	5
	Foot – Chevron Osteotomy for Hallux Valgus	14
	Using the Optional Drill Guide with Stop	21
	Screw Extraction	22
Product Information	Implants HCS 2.4	23
	Implants HCS 3.0	25
	Instruments HCS 2.4 and 3.0	27
	Optional Instruments for HCS 2.4 and 3.0	29
	Setlists HCS 2.4 and 3.0	30
MRI Information		31



Image intensifier control

This description alone does not provide sufficient background for direct use of DePuy Synthes products. Instruction by a surgeon experienced in handling these products is highly recommended.

Processing, Reprocessing, Care and Maintenance

For general guidelines, function control and dismantling of multi-part instruments, as well as processing guidelines for implants, please contact your local sales representative or refer to:

<http://emea.depuyssynthes.com/hcp/reprocessing-care-maintenance>


For general information about reprocessing, care and maintenance of Synthes reusable devices, instrument trays and cases, as well as processing of Synthes non-sterile implants, please consult the Important Information leaflet (SE_023827) or refer to:

<http://emea.depuyssynthes.com/hcp/reprocessing-care-maintenance>

Features and Benefits


Tip with self-drilling and self-tapping flutes Simplified surgical technique

Shorter surgery due to simplified surgical technique



Cannulation

For minimally invasive technique and guided insertion




Head with self-tapping flutes

Facilitates countersinking of screw



Two different thread lengths of the shaft

The optimal implant for every case due to threads available in different lengths



Identical pitch of head and shaft threads

For controlled closure and compression of the fracture gap



Available in steel and titanium


All Headless Compression screws from Synthes are available both in stainless implant-grade steel and high-quality bio-compatible titanium alloy (TAN)

Lag Screw Technique with Compression Sleeve



Step 1: Screw insertion

Insertion of the screw into the bone with the compression sleeve.



Step 2: Closure of gap and compression

Once the tip of the compression sleeve lies on the bone, the fracture gap is closed and compressed by further turning of the sleeve.



Step 3: Countersinking

Once the desired degree of compression is reached, the screw is countersunk into the bone with the screwdriver while the compression sleeve is held stationary. During countersinking no additional compression is generated.

Indications

-
- Fixation of intra-articular and extra-articular fractures and non-unions of small bones and small bone fragments
 - Arthrodeses of small joints
 - Unionectomies and osteotomies

Examples include, but are not limited to scaphoid and other carpal bones, metacarpals, tarsals, metatarsals, patella, ulnar styloid, capitellum, radial head and radial styloid.

The following simplified surgical technique for a scaphoid fracture serves as example for the use of the HCS 2.4 or 3.0 in the hand.

1


Insert the guide wire

Instruments

292.623	Guide Wire Ø 1.1 mm with trocar tip, length 150 mm, Stainless Steel
or	
292.622	Guide Wire Ø 1.1 mm with threaded tip, with trocar, length 150 mm, Stainless Steel
312.151	Double Drill Guide 2.0/1.1

- While monitoring with the image intensifier, advance the guide wire through the drill guide into the bone until the thread tip is anchored in the far cortex.

Precaution: Do not forcefully insert the guide wire. This may cause it to bend.



2

Option: ream the trapezium

Instruments



03.226.003	Trapezium Burr, cannulated, for HCS – Headless Compression Screw Ø 2.4/3.0 mm
03.226.005	Protection Sleeve for HCS – Headless Compression Screw Ø 2.4/3.0 mm, for Trapezium Burr
311.430	Handle with Quick Coupling

To facilitate screw insertion, the flank of the trapezium can be removed with the trapezium burr.

Slide the trapezium burr with the protection sleeve over the guide wire and carefully ream the trapezium.

Ensure that the trapezium burr does not damage the scaphoid by using image intensification.

Note: Do not forcefully insert the trapezium burr since this may damage the guide wire.



3

Determine screw and thread length

Instruments


03.226.002	Direct Measuring Device for HCS – Headless Compression Screw Ø 2.4/3.0 mm
292.623	Guide Wire Ø 1.1 mm with trocar tip, length 150 mm, Stainless Steel
or	
292.622	Guide Wire Ø 1.1 mm with threaded tip, with trocar, length 150 mm, Stainless Steel

Slide the narrow end of the measuring device over the guide wire to the bone.

The measurement on the measuring device shows the depth of the guide wire in the bone in millimeters and directly the appropriate screw length.

If the screw is to be countersunk below the surface of the bone, subtract the appropriate screw length. If a large fracture gap needs to be closed or if the screw is inserted at an angle to the bone surface, subtract more.


Note: Only use the guide wire in its original length to ensure correct measurement.



The position of the fracture line determines the thread length


Correctly selected thread length

The shaft thread lies completely within the proximal fragment during compression. Fragments can hence be compressed.



Incorrect thread length

The shaft thread lies over the fracture gap. Fragments cannot be compressed.



Note: If there is no good bone quality in the distal part of the bone, the distal screw thread can strip-out if too much compression is applied.

4

Predrilling

Instruments

310.221	Drill Bit Ø 2.0/1.15 mm, cannulated, length 150/50 mm, 3-flute, for Quick Coupling
312.151	Double Drill Guide 2.0/1.1

Predrilling makes it substantially easier to insert the screw in dense bone.

Slide the drill guide with the drill bit over the guide wire and pre-drill to the desired depth.

- Remove the drill guide and verify the effective drilling depth with the image intensifier.

Precaution: Avoid forcefully advancing or bending the drill bit as this may cause the drill bit to break. Do not advance the drill bit beyond the wire. When drilling is complete, slowly pull the drill bit straight out while running the power tool in “forward mode” to ensure the guide wire stays in place.



5

Pick up screw

Instruments


For HCS 2.4

03.226.016	Compression Sleeve for HCS – Headless Compression Screw Ø 2.4 mm
------------	---

For HCS 3.0

03.226.000	Compression Sleeve for HCS – Headless Compression Screw Ø 3.0 mm
------------	---

Twist the compression sleeve over the head thread of the screw to remove the screw from the screw rack.



6

Insert screw and compress fragment

Instruments

03.226.006	Handle for Compression Sleeve, for HCS – Headless Compression Screw Ø 1.5–3.0 mm
------------	--

For HCS 2.4

03.226.016	Compression Sleeve for HCS – Headless Compression Screw Ø 2.4 mm
------------	---

For HCS 3.0

03.226.000	Compression Sleeve for HCS – Headless Compression Screw Ø 3.0 mm
------------	---

Slide the handle into the selected compression sleeve. Insert the screw into the bone until the fracture gap or the osteotomy is closed and compressed.

Notes


- Verify the correct position of the shaft thread in the proximal fragment using the image intensifier. If the thread lies over the fracture gap or the osteotomy, the gap cannot be compressed.

Precaution

- Carefully tighten the screw with the compression sleeve. Forceful tightening could cause stripping of the shaft thread.
- If the thread strips, some or all of the compression will be lost. If the screw is then countersunk correctly, the thread will regain purchase, thereby reducing the danger of post-operative screw loosening.

Note

- If loss of compression makes screw extraction necessary, follow the instructions on screw extraction on page 22.
-



7

Countersink screw

Instruments

03.226.004	Screwdriver Shaft, cannulated, Stardrive, SD8, with coloured marking, for HCS – Headless Compression Screw Ø 2.4 / 3.0 mm
------------	---

311.430	Handle with Quick Coupling
---------	----------------------------

For HCS 2.4

03.226.016	Compression Sleeve for HCS – Headless Compression Screw Ø 2.4 mm
------------	--


For HCS 3.0

03.226.000	Compression Sleeve for HCS – Headless Compression Screw Ø 3.0 mm
------------	--

Remove the compression sleeve handle and slide the cannulated screwdriver through the compression sleeve.


Countersink the screw by turning the screwdriver shaft while simultaneously holding the compression sleeve stationary.

Verify the screw position with the image intensifier. Ensure that the screw tip does not penetrate the proximal cortex. Remove and dispose of the guide wire.




Color markings

The color markings on the screwdriver shaft show the position of the screwdriver tip and the head thread of the screw.




Green mark at the top end of the compression sleeve

The screwdriver tip is seated correctly in the Stardrive recess of the screw.



Yellow mark at the top end of the compression sleeve

The top end of the head thread is even with the bone surface.



Red mark at the top end of the compression sleeve

The top end of the head thread is approximately 2 mm below the bone surface.

Note: If the screw is inserted at an angle, it must be countersunk further than the yellow mark so that it does not project from the surface.

Foot – Chevron Osteotomy for Hallux Valgus

The following simplified surgical technique for a chevron osteotomy serves as example for the use of the HCS 2.4 or 3.0 in the foot.

1

Remove bunion and perform V-shaped osteotomy

- Remove the bunion on the medial side of the first metatarsal with a saw blade under image intensification.

Perform a V-shaped osteotomy (inner angle approx. 55°), with the peak approximately 2 mm distal from the center of the head of the first metatarsal.



2

Move the distal fragment laterally

Move the distal fragment in a lateral direction to correct the alignment.

3


Insert the guide wire

Instruments

292.623	Guide Wire Ø 1.1 mm with trocar tip, length 150 mm, Stainless Steel
or	
292.622	Guide Wire Ø 1.1 mm with threaded tip, with trocar, length 150 mm, Stainless Steel
312.151	Double Drill Guide 2.0/1.1

- While monitoring with the image intensifier, advance the guide wire through the double drill guide from proximal dorsal to distal plantar through the osteotomy into the bone until the thread tip is anchored in the far cortex.

Precaution: Do not forcefully insert the guide wire. This may cause it to bend.



4

Determine screw and thread length

Instrument


03.226.002	Direct Measuring Device for HCS – Headless Compression Screw Ø 2.4/3.0 mm
------------	---

Slide the narrow end of the measuring device over the guide wire to the bone.

The measurement on the measuring device shows the depth of the guide wire in the bone in millimeters.

If the screw is to be countersunk below the surface of the bone, subtract the appropriate screw length. If a large fracture gap needs to be closed or if the screw is inserted at an angle to the bone surface, subtract more.

Note: The position of the osteotomy line determines the thread length (see page 8).



5

Predrilling

Instruments


310.221	Drill Bit Ø 2.0/1.15 mm, cannulated, length 150/50 mm, 3-flute, for Quick Coupling
312.151	Double Drill Guide 2.0/1.1

Predrilling makes it substantially easier to insert the screw in dense bone.

Slide the double drill guide with the drill bit over the guide wire and predrill to the desired depth.

- Verify the effective drilling depth with the image intensifier.

Precaution: Avoid forcefully advancing or bending the drill bit as this may cause the drill bit to break. Do not advance the drill bit beyond the wire. When drilling is complete, slowly pull the drill bit straight out while running the power tool in “forward mode” to ensure the guide wire stays in place.



6

Insert screw and compress osteotomy

Instruments

03.226.000	Compression Sleeve for HCS – Headless Compression Screw Ø 3.0 mm
03.226.006	Handle for Compression Sleeve, for HCS – Headless Compression Screw Ø 1.5-3.0 mm

For HCS 2.4

03.226.016	Compression Sleeve for HCS – Headless Compression Screw Ø 2.4 mm
------------	---

For HCS 3.0

03.226.000	Compression Sleeve for HCS – Headless Compression Screw Ø 3.0 mm
------------	---



Twist the selected compression sleeve over the head thread of the screw to remove the screw from the screw rack.

Slide the handle into the compression sleeve. Insert the screw into the bone until the osteotomy is closed and compressed.

Notes

- Verify the correct position of the shaft thread in the distal fragment using the image intensifier. If the thread lies over the osteotomy, the gap cannot be compressed.

Precaution

- Carefully tighten the screw with the compression sleeve. Forceful tightening could cause stripping of the shaft thread.
- If the thread strips, some or all of the compression will be lost. If the screw is then countersunk correctly, the thread will regain purchase, thereby reducing the danger of post-operative screw loosening.

Note

- If loss of compression makes screw extraction necessary, follow the instructions on screw extraction on page 22.

7

Countersink screw

Instruments

03.226.000	Compression Sleeve for HCS – Headless Compression Screw Ø 3.0 mm
03.226.004	Screwdriver Shaft, cannulated, Stardrive, SD8, with coloured marking, for HCS – Headless Compression Screw Ø 2.4/3.0 mm
311.430	Handle with Quick Coupling

For HCS 2.4

03.226.016	Compression Sleeve for HCS – Headless Compression Screw Ø 2.4 mm
------------	---

For HCS 3.0


03.226.000	Compression Sleeve for HCS – Headless Compression Screw Ø 3.0 mm
------------	---

Remove the compression sleeve handle and slide the cannulated screwdriver through the compression sleeve.

Countersink the screw by turning the screwdriver shaft while simultaneously holding the compression sleeve stationary.

- ⓘ Verify the screw position with the image intensifier. Ensure that the screw tip does not penetrate the distal cortex. Remove and dispose of the guide wire.

Note: The color markings on the screwdriver shaft show the position of the screwdriver tip and head thread of the screw (see page 13).



8

Remove protruding bone

Remove the protruding bone of the proximal fragment.



Using the Optional Drill Guide with Stop

The drill guide with stop allows control of the drilling depth and can be used for drilling instead of the double drill guide 2.0/1.1 (312.151).


Instruments

310.221	Drill Bit Ø 2.0/1.15 mm, cannulated, length 150/50 mm, 3-flute, for Quick Coupling
03.226.007	Drill Guide with Stop for Drill Bits Ø 2.0/1.15 mm No. 310.221
03.226.008	Direct Measuring Device for Drill Guide with Stop No. 03.226.007

To set the drilling depth, insert the drill bit in the drill guide with stop, and slide the measuring device over the drill bit until the retaining device is engaged.

Release the locking ring, and set the drilling depth by rotating the tip of the drill guide. The measurement on the measuring device indicates the set drilling depth in millimeters.

Tighten the locking ring to fix the drilling depth.



Screw Extraction

Instruments

314.467	Screwdriver Shaft, Stardrive, SD8, self-holding
or	
03.226.004	Screwdriver Shaft, cannulated, Stardrive, SD8, with coloured marking, for HCS – Headless Compression Screw Ø 2.4/3.0 mm

311.430	Handle with Quick Coupling
---------	----------------------------

For HCS 2.4

03.226.016	Compression Sleeve for HCS – Headless Compression Screw Ø 2.4 mm
------------	--


For HCS 3.0

03.226.000	Compression Sleeve for HCS – Headless Compression Screw Ø 3.0 mm
------------	--

For the extraction of the HCS use a Stardrive screwdriver or Screwdriver Shaft in combination with the handle.

Note: If the screw strips, use the following procedure:
Twist the compression sleeve over the head thread and insert the screwdriver through the compression sleeve into the Stardrive recess of the screw.
Remove the screw by simultaneously pulling on the compression sleeve and turning both the screwdriver and the compression sleeve in counterclockwise direction.

Note: If necessary, expose the recess and part of the head thread with a hollow reamer (e.g. 309.035) or preferred method.



Implants HCS 2.4

HCS 2.4 mm – Headless Compression Screw, short thread

Art. No.	Screw length (mm) L	Shaft thread length (mm) S
0X.226.209	9	4
0X.226.210	10	4
0X.226.211	11	4
0X.226.212	12	4
0X.226.213	13	4
0X.226.214	14	4
0X.226.215	15	4
0X.226.216	16	4
0X.226.217	17	4
0X.226.218	18	4
0X.226.219	19	4
0X.226.220	20	4
0X.226.221	21	4
0X.226.222	22	4
0X.226.223	23	4
0X.226.224	24	5
0X.226.225	25	5
0X.226.226	26	5
0X.226.227	27	6
0X.226.228	28	6
0X.226.229	29	6
0X.226.230	30	7
0X.226.232	32	7
0X.226.234	34	8
0X.226.236	36	9
0X.226.238	38	9
0X.226.240	40	10




X = 2: Stainless Steel

X = 4: Titanium Alloy (TAN)

All implants are also available sterile packed. Add suffix "S" to article number.

HCS 2.4 mm – Headless Compression Screw,
long thread

Art. No.	Screw length (mm) L	Shaft thread length (mm) S
0X.226.316	16	5
0X.226.317	17	6
0X.226.318	18	6
0X.226.319	19	7
0X.226.320	20	7
0X.226.321	21	8
0X.226.322	22	8
0X.226.323	23	8
0X.226.324	24	8
0X.226.325	25	8
0X.226.326	26	10
0X.226.327	27	10
0X.226.328	28	10
0X.226.329	29	10
0X.226.330	30	12
0X.226.332	32	12
0X.226.334	34	14
0X.226.336	36	14
0X.226.338	38	16
0X.226.340	40	16



X = 2: Stainless Steel
X = 4: Titanium Alloy (TAN)

All implants are also available sterile packed. Add suffix “S”
to article number.

Implants HCS 3.0

HCS 3.0 mm – Headless Compression Screw, short thread

Art. No.	Screw length (mm) L	Shaft thread length (mm) S
0X.226.010	10	4
0X.226.011	11	4
0X.226.012	12	4
0X.226.013	13	4
0X.226.014	14	4
0X.226.015	15	4
0X.226.016	16	4
0X.226.017	17	4
0X.226.018	18	4
0X.226.019	19	4
0X.226.020	20	4
0X.226.021	21	4
0X.226.022	22	4
0X.226.023	23	4
0X.226.024	24	5
0X.226.025	25	5
0X.226.026	26	5
0X.226.027	27	6
0X.226.028	28	6
0X.226.029	29	6
0X.226.030	30	7
0X.226.032	32	7
0X.226.034	34	8
0X.226.036	36	9
0X.226.038	38	9
0X.226.040	40	10




X = 2: Stainless Steel

X = 4: Titanium Alloy (TAN)

All implants are also available sterile packed. Add suffix "S" to article number.

HCS 3.0 mm – Headless Compression Screw,
long thread

Art. No.	Screw length (mm)	Shaft thread length (mm)
	L	S
0X.226.116	16	5
0X.226.117	17	6
0X.226.118	18	6
0X.226.119	19	7
0X.226.120	20	7
0X.226.121	21	8
0X.226.122	22	8
0X.226.123	23	8
0X.226.124	24	8
0X.226.125	25	8
0X.226.126	26	10
0X.226.127	27	10
0X.226.128	28	10
0X.226.129	29	10
0X.226.130	30	12
0X.226.132	32	12
0X.226.134	34	14
0X.226.136	36	14
0X.226.138	38	16
0X.226.140	40	16



X = 2: Stainless Steel
X = 4: Titanium Alloy (TAN)

All implants are also available sterile packed. Add suffix "S"
to article number.

Instruments HCS 2.4 and 3.0

292.623 Guide Wire Ø 1.1 mm with trocar tip
length 150 mm, Stainless Steel



312.151 Double Drill Guide 2.0/1.1
For protecting soft tissue during insertion
of guide wires and predrilling



03.226.002 Direct Measuring Device for HCS –
Headless Compression Screw
Ø 2.4/3.0 mm



310.221 Drill Bit Ø 2.0/1.15 mm, cannulated,
length 150/50 mm, 3-flute,
for Quick Coupling
For predrilling



03.226.006 Handle for Compression Sleeve, for HCS –
Headless Compression Screw
Ø 1.5–3.0 mm



311.430 Handle with Quick Coupling
For Stardrive SD8 Screwdriver Shafts
(03.226.004 and 314.467)



03.226.004 Screwdriver Shaft, cannulated, Stardrive,
SD8, with coloured marking,
for HCS – Headless Compression Screw
Ø 2.4/3.0 mm
For countersinking the screw; with color
markings to control countersink depth



314.467 Screwdriver Shaft, Stardrive SD8,
self-holding
For screw extraction; with self-retaining tip



319.970 Screw Forceps, self-holding, length 85 mm



319.292 Cleaning Stylet Ø 1.1 mm,
for Cannulated Instruments
For cleaning cannulated instruments
during surgery









For HCS 2.4
03.226.016 Compression Sleeve for HCS – Headless
Compression Screw Ø 2.4 mm
For closing the fracture gap and
compressing the bone fragments



For HCS 3.0
03.226.000 Compression Sleeve for HCS – Headless
Compression Screw Ø 3.0 mm
For closing the fracture gap and
compressing the bone fragments



Optional Instruments for HCS 2.4 and 3.0

292.622	Guide Wire Ø 1.1 mm with threaded tip with trocar, length 150 mm, Stainless Steel	
03.226.003	Trapezium Burr, cannulated For freeing the palmar approach to the distal pole of the scaphoid	
03.226.005	Protection Sleeve for HCS – Headless Compression Screw Ø 2.4/3.0 mm for Trapezium Burr For protecting soft tissue during use of the trapezium burr	
03.226.007	Drill Guide with Stop for Drill Bits Ø 2.0/1.15 mm No. 310.221 For controlled drilling	
03.226.008	Direct Measuring Device for Drill Guide with Stop No. 03.227.007 For determining the drilling depth	
398.408	Periosteal Elevator, slightly curved blade, round tip, width 5 mm For manipulating small bones and bone fragments	
398.409	Sharp Reduction Hook, graded For levering up carpal bones	

Setlists HCS 2.4 and 3.0

HCS 2.4

01.226.012	Instrument and Implant Set for HCS – Headless Compression Screw Ø 2.4 mm (Stainless Steel) for Vario Case
------------	---

01.226.014	Instrument and Implant Set for HCS – Headless Compression Screw Ø 2.4 mm (Titanium Alloy) for Vario Case
------------	--

To adapt for use with the HCS 3.0, the following articles must additionally be ordered:

03.226.000	Compression Sleeve for HCS – Headless Compression Screw Ø 3.0 mm
------------	--

68.111.443	Insert for Screw Rack Module, for HCS – Headless Compression Screw Ø 3.0 mm
------------	---

HCS 3.0

01.226.002	Set for Instruments and Implants for HCS – Headless Compression Screw Ø 3.0 mm (Stainless Steel), for Vario Case
------------	--

01.226.004	Set for Instruments and Implants for HCS – Headless Compression Screw Ø 3.0 mm (Titanium Alloy), for Vario Case
------------	---

To adapt for use with the HCS 2.4, the following articles must additionally be ordered:

03.226.016	Compression Sleeve for HCS – Headless Compression Screw Ø 2.4 mm
------------	--

68.111.446	Insert for Screw Rack Module, for HCS – Headless Compression Screw Ø 2.4 mm
------------	---

Torque, Displacement and Image Artifacts according to ASTM F 2213-06, ASTM F 2052-06e1 and ASTM F2119-07

Non-clinical testing of worst case scenario in a 3 T MRI system did not reveal any relevant torque or displacement of the construct for an experimentally measured local spatial gradient of the magnetic field of 3.69 T/m. The largest image artifact extended approximately 169 mm from the construct when scanned using the Gradient Echo (GE). Testing was conducted on a 3 T MRI system.

Radio-Frequency-(RF-)induced heating according to ASTM F2182-11a

Non-clinical electromagnetic and thermal testing of worst case scenario lead to peak temperature rise of 9.5 °C with an average temperature rise of 6.6 °C (1.5 T) and a peak temperature rise of 5.9 °C (3 T) under MRI Conditions using RF Coils (whole body averaged specific absorption rate [SAR] of 2 W/kg for 6 minutes [1.5 T] and for 15 minutes [3 T]).

Precautions: The above mentioned test relies on non-clinical testing. The actual temperature rise in the patient will depend on a variety of factors beyond the SAR and time of RF application. Thus, it is recommended to pay particular attention to the following points:

- It is recommended to thoroughly monitor patients undergoing MR scanning for perceived temperature and/or pain sensations.
 - Patients with impaired thermoregulation or temperature sensation should be excluded from MR scanning procedures.
 - Generally, it is recommended to use a MR system with low field strength in the presence of conductive implants. The employed specific absorption rate (SAR) should be reduced as far as possible.
 - Using the ventilation system may further contribute to reduce temperature increase in the body.
-

