

Instruments for Measurement of Dissolved CO₂ and O₂ in Beverages

CboxQC | **CboxQC At-line**
CarboQC | **CarboQC At-line**
OxyQC | **OxyQC Wide Range**
PFD filling device (Plus)



Product portfolio



CarboQC
CarboQC At-line



OxyQC
OxyQC Wide Range

CboxQC
CboxQC At-line



CarboQC ME
plus Option O₂

Product portfolio

System



Modularity for measurement setup

CO₂
Air Index
N₂ Index
O₂ Index



CarboQC ME



CO₂
O₂
Air Index
N₂ Index

CarboQC ME + Option O₂

O₂
O₂ WR

Laboratory



Most precise measurement



CarboQC




CboxQC



OxyQC | OxyQC Wide Range

At-line



Compact, easy robust measurement



CarboQC At-line



CboxQC At-line



OxyQC | OxyQC Wide Range

CboxQC and **CboxQC At-line** are **combined CO₂ and O₂ meters** to determine **the true amount of dissolved CO₂ and O₂** in beverages for use

- in the lab in combination with the PFD filling device
- and at-line.



Product portfolio

System



Modularity for measurement setup

CO₂
Air Index
N₂ Index
O₂ Index



CarboQC ME

CO₂
O₂
Air Index
N₂ Index



CarboQC ME + Option O₂

O₂
O₂ WR

Laboratory



Most precise measurement



CarboQC




CboxQC



OxyQC | OxyQC Wide Range

At-line



Compact, easy robust measurement



CarboQC At-line



CboxQC At-line



OxyQC | OxyQC Wide Range

CarboQC and CarboQC At-line are **CO₂ meters** to determine the **true amount of dissolved CO₂** in beverages for use

- in the lab in combination with the PFD filling device
- and at-line.



Product portfolio

System



Modularity for measurement setup



CarboQC ME

CO₂
Air Index
N₂ Index
O₂ Index



CarboQC ME + Option O₂

CO₂
O₂
Air Index
N₂ Index

O₂
O₂ WR

Laboratory



Most precise measurement



CarboQC




CboxQC



OxyQC | OxyQC Wide Range

At-line



Compact, easy robust measurement



CarboQC At-line



CboxQC At-line



OxyQC | OxyQC Wide Range

CarboQC ME (plus Option O₂) is a CO₂ and O₂ measuring module to determine the true amount of dissolved CO₂ and O₂ in beverages

— designed to be used with Anton Paar's Generation M instruments



Product portfolio

System



Modularity for measurement setup



CarboQC ME

CO₂
Air Index
N₂ Index
O₂ Index



CarboQC ME + Option O₂

CO₂
O₂
Air Index
N₂ Index

O₂
O₂ WR

Laboratory



Most precise measurement



CarboQC




CboxQC



OxyQC|OxyQC Wide Range

At-line



Compact, easy robust measurement



CarboQC At-line



CboxQC At-line



OxyQC|OxyQC Wide Range

OxyQC and OxyQC Wide Range are **O₂ meters** to determine **the true amount of dissolved O₂** in beverages for use

- in the lab in combination with the PFD filling device
- as well as **at-line**

Note: the instruments are designed to be used for both lab and at-line applications and are equipped with a rubber protection.



Product portfolio

System



Modularity for measurement setup

CO₂
Air Index
N₂ Index
O₂ Index



CarboQC ME



CO₂
O₂
Air Index
N₂ Index

CarboQC ME + Option O₂

O₂
O₂ WR

Laboratory



Most precise measurement



CarboQC




CboxQC



OxyQC | OxyQC Wide Range

At-line



Compact, easy robust measurement



CarboQC At-line



CboxQC At-line



OxyQC | OxyQC Wide Range

Whether directly
at the production line,
in the **lab** or as part of
an **analyzing system,**
*...the **best instrument***
for **every application...**



Why is dissolved CO₂ measured?

The CO₂ content

- ▶ strongly **influences** the **fresh taste** of beverages
- and
- ▶ is a **considerable cost factor**

CO₂ measurement ensures


- ▶ **consistent taste**
- ▶ **cost-efficient dosing** due to production control



Where is dissolved CO₂ measured?

CO₂ is measured in...

- ▶ Soft drinks
- ▶ Beer, mixed beer, cider
- ▶ Mineral water
- ▶ Wine and sparkling wine



Why is dissolved O₂ measured?

The O₂ content

- ▶ strongly **reduces** the **shelf life** of beverages

and

- ▶ can **cause corrosion** of cans

O₂ measurement ensures

- ▶ **product safety**
- ▶ consistent **beverage quality**



Where is dissolved O₂ measured?

O₂ is measured in...

- ▶ Beer, mixed beer, cider
- ▶ Wine and sparkling wine
- ▶ Juices
- ▶ Some soft drinks



growing demand due to trend towards natural flavorings, addition of vitamins and avoiding of preservatives



Basics of CO₂ measurement

Fundamental gas laws in beverage CO₂ measurement

CO₂
measurement

Henry's law



William Henry [1775 – 1836]

$$c_{\text{CO}_2} = \xi_{\text{CO}_2} \cdot p_{\text{CO}_2}$$


Dalton's law



John Dalton [1766 – 1844]

$$p_{\text{tot}} = p_{\text{CO}_2} + p_{\text{Air}} + p_{\text{H}_2\text{O}}$$

'Real gas law'



$$m = \frac{(p \cdot M \cdot V)}{(R \cdot T \cdot Z)}$$

c ... CO₂ concentration
p ... CO₂ partial pressure
ξ ... absorption coefficient

m ... CO₂ mass
M ... molar mass
Z ... compressibility factor

V ... volume
R ... gas constant
T ... temperature

Basics of CO₂ measurement

Henry's law

$$c_{\text{CO}_2} = \xi_{\text{CO}_2} * p_{\text{CO}_2}$$

c_{CO₂} ... Concentration of CO₂ in a liquid

ξ_{CO₂} ... Absorption coefficient/solubility

p_{CO₂} ... Partial pressure of CO₂


Absorption coefficient ξ is influenced by:



Temperature



Beverage composition
(sugar, alcohol, flavors, ...)



Other influences
(CO₂ pressure, dissolved air ...)



The CO₂ solubility in beverages is up to 20% lower than it is in pure water!

Basics of CO₂ measurement

Measurement units

- ▶ Gram per liter [g/L] Grams of CO₂ per liter of liquid
- ▶ Volumes [vol.] Volumes of CO₂ at 0 °C and 1.01326 bar
(std. atmosphere) dissolved per liquid volume

▶ 1 vol. = 1.976 g/L Conversion

Anton Paar's CO₂ measuring principle

What makes the Multiple Volume Expansion method that outstanding?

Anton Paar's CO₂ meter **selectively** determines the **concentration of dissolved CO₂** in beverages **not influenced** by other residual gases such as **air** or **nitrogen**.

Anton Paar's CO₂ measuring principle

How does the
Multiple Volume
Expansion
Method work?



The **measuring chamber is completely filled** with sample and **closed by valves**.



The **volume of the measuring chamber is expanded**.
Pressure and temperature equilibrium is generated.
Equilibrium pressure and temperature are measured.



The measuring chamber **volume is further expanded**,
equilibrium is generated and **pressure and temperature are measured again**.

The **two pressures/temperatures** are used for
CO₂ determination and **dissolved air compensation**.



Anton Paar's CO₂ measuring principle

Why does the selective measurement work?


The Multiple Volume Expansion method makes use of the fact that the **solubility of air** in beverages **is much lower** than that of **CO₂**!

How does this work? ...see the next slides!



Anton Paar's CO₂ measuring principle

Why is there no influence of other residual gases in the CO₂ result?



- The CO₂ pressure **decreases gently** due to the **large solubility of CO₂** in beverages
- The **air** pressure **decreases excessively** due to the very **low solubility of air** in beverages

Anton Paar's CO₂ measuring principle

MVE method – selective measurement

Additional output of Anton Paar's CO₂ meter:

- ▶ Air Index (sum of dissolved N₂ + O₂)
- ▶ O₂ Index (dissolved O₂) → if beverage is filled under O₂
- ▶ N₂ Index (dissolved N₂) → if beverage is filled under N₂
- ▶ Package pressure

***Results just for
orientation***



Influences of CO₂ measurement

Why do CO₂ results differ depending on the measuring method used?

▶ **Influence of altitude**

No altitude correction vs. integrated absolute pressure sensor of Anton Paar's CO₂ meters for accurate correction

▶ **Influence of sniffing**

Traditional p/T methods require sniffing by releasing the overpressure in the headspace to the ambience in order to correct the influence of other gases; CO₂ also gets released

▶ **Influence of methods and tables**

Traditional p/T methods require tables which compensate the effect of sniffing


▶ **Influence of dissolved air in the sample**

No pressure differentiation vs. partial pressures of air and CO₂ within the MVE method; the MVE method measures CO₂ independent of other dissolved gases




Measurement influence: sniffing with p/T methods

Large headspace - High temperature:




Small headspace - Low temperature:



- ▶ Sniffing effect depends on:
 - extent of pressure
 - equilibrium
 - headspace percentage
 - duration of sniffing
 - temperature
 - CO₂ content, etc.
- ▶ Thus variable amounts of CO₂, O₂, and N₂ are removed from the package prior to the equilibration and p/T measurement



Measurement influence: gas distribution in a package



Headspace:
Cans: 3 – 5 %
Bottles: 5 – 12%

Non-equilibrium state at bottle/can filling requires vigorous shaking of bottles/cans after filling!



Measurement influence: temperature dependency

Bottle with e.g. 10% headspace

Equilibrium state:
at 5 °C (41 °F)



Equilibrium state:
at 25 °C (77 °F)



**Measurement in a small
range of package
temperatures is
necessary!
Temperature-
compensated measure-
ment is possible.**

Basics of O₂ measurement

Optochemical O₂ measurement


The new series of oxygen meter with CboxQC, OxyQC and OxyQC Wide Range use an **optochemical sensor** for **fast** and **highly accurate** determination of dissolved oxygen.

Oxygen is measured based on the **effect of dynamic luminescence quenching** and is **measured while liquids are flowing through the cell** – it's **not a stationary** measurement.




Physical principle of O₂ measurement

- ▶ **Optochemical O₂ sensor**
- ▶ Inside the sensor, a light-emitting diode illuminates the light-sensitive layer.



O₂ measuring principle


Luminescence process in absence of oxygen



- ▶ Dye molecule absorbs light
- ▶ Transfer of electrons to higher energy level
- ▶ Return to the original energy level, emission of light (luminescence)

O₂ measuring principle


In the presence of oxygen: the quenching process



- ▶ Dye molecule absorbs light → Transfer to a higher energy level
- ▶ Collision with oxygen: deactivation of the dye molecule
- ▶ Return to the original energy level, no emission of light

Fast response time and ideal temperature behavior


- ▶ The sensor characteristic shows a **response time of less than 20 s** when changing from **water to beer** at **4 °C to 7 °C**
- ▶ Excellent measurement **stability**
- ▶ Even with **cold samples**



High repeatability of the O₂ sensor

Measurement conditions:

- ▶ 0.5 L beer
- ▶ 6 °C to 11 °C
- ▶ Rinsing volume before start: **240 mL**
- ▶ Rinsing volume between each measurement: **50 mL**



Basics of O₂ measurement

Measurement units

- ▶ [mg/L] Milligrams of O₂ per liter of liquid
- ▶ [ppm] Parts per million
- ▶ [ppb] Parts per billion
- ▶ [% Air-sat.] Percent saturation in Air
- ▶ [% O₂-sat.] Precent saturation in Oxygen

1 mg/L = 1 ppm = 1000 ppb


Basics of TPO measurement

What is TPO?

- ▶ Total Package Oxygen is the sum of dissolved and headspace oxygen

Why determine TPO?

- ▶ Dissolved and headspace oxygen in beverages reacts with certain beverage components.
- ▶ Oxidation products may cause changes in aroma and taste, color changes, loss of nutritional value.



Basics of TPO measurement

- ▶ **Measure dissolved oxygen, calculate TPO!**
 - ▶ Reliable information on Total Package Oxygen in beverage packages can be obtained by measuring dissolved oxygen and calculating TPO under the following conditions:
 - ▶ equilibrium established.
 - ▶ correct determination of dissolved oxygen.
 - ▶ correct temperature measurement.
 - ▶ correct brimful volume and actual volume information.
 - ▶ correct empty weight and total weight information.



Basics of TPO measurement

- ▶ **The affordable way: measure dissolved oxygen, calculate TPO!**

CboxQC (At-line) | OxyQC | OxyQC Wide Range:

AP-SoftPrint (from V1.5) provides calculation via volume and weight information.

O ₂ ppm	temp. C	Brim volume ml	Liquid volume ml	TPO mg	TPO/L mg/L	Air Content ml	Total Air mg
2,665	24	110	100	1,1	11,2	4,1	4,4
5,641	26,6	3000	2700	71,5	26,5	263,1	284,9
5,641	26,6	130	120	2,6	21,3	9,4	10
5,641	26,6	500	450	11,9	26,5	43,8	47,5
1,76	20	500	450	3,5	7,7	12,4	13,6

AP-SoftPrint

Option O₂, CboxQC in combination with Generation M series: Master instrument provides calculation.

(available from SW-release 2.40)



DMA Generation M software