


iTero STL export instructions for Restorative Cases

1. Please log into www.myaligntech.com
2. Make sure to install the iTero online software, prior the export, you can download the program from http://www.itero.com/en-us/resources/itero_element/software_downloads
3. Upon logging in, there are two ways to review case details. There is a “Quick Start” to search patient details or a complete list of patients can be viewed by clicking “Orders” at the top left menu and choosing “Orders List”.



6. Once you have located your patient’s information page, click “Export”; which requires Itero Viewer 4.1 or higher to begin the export process.

Actions: [Open File](#) [Export \(iTero 4.1 or higher\)](#) [Open RX](#)

7. Choose a specific CAD/CAM System or generic option



8. After choosing a CAD/CAM System or generic export option there is an option to name the export folder located within the default export location. Each patient file export can have its own folder. The default name is the iTero case ID but it can be changed to any desired name (ie. Patient name, chart number, etc...).



9. Press "Export" to begin the STL export process. Once completed, the following window will appear; indicating the export was completed.



How to export the STL file from My iTero website.


1. Please visit <https://myitero.com/login>
2. Login with your iTero user.




3. Go to orders tap.



4. Choose the file you would like to export.



5. Then choose the export specifications; Export type, Data format. The file type will always be "STL (No color)"



- 6. Please make sure to allow the pop-ups for this site. You can find the icon next to the URL bar;



Your files should now be saved to your computer and are ready to upload and transfer to our lab.



Want to transfer your digital scan STL files to us?

Once you have exported and saved your STL files, you will need to transfer and upload them to our lab.

We recommend using the free website WeTransfer.com for their security and ease of use.


Here's how to upload your files to our lab:

1. Go to [WeTransfer.com](https://www.wetransfer.com).
2. Choose the "Take Me to Free" option.
3. Choose "I Agree" to the Terms and Condition
4. Click "+ Add your files" and select the files from your computer that you want to transfer.
5. In the "Email to" field, enter: support@mabeldental.com
6. Enter your email address in the "Your email" field
7. Enter an optional message
8. Click the "Transfer" button

Voila! The file(s) has been sent!

We recommend that you fax/email a Prescription Rx.

You'll receive a confirmation email/phone call within 1 business day.

A screenshot of the WeTransfer website interface. At the top, there is a blue circular button with a white plus sign and the text "Add your files". Below this are three input fields: "Email to", "Your email", and "Message". At the bottom, there is a grey button labeled "Transfer" and a small circular icon with three dots to its left.

Email to

Your email

Message



Transfer

*please note this may not reflect the exact look of the screen