


# Lab 2:


## Metric Measurement and Microscopy



# The Metric System

- “SI” Metric is the standard system of measurement used in the sciences.
- Why use metric?
  - Official system of almost every country.
  - Less confusing if scientists use the same measurement system.
  - Based on units of ten, which makes converting units easy!

# How To Read A Metric Ruler



# Metric Basics

- Most measurements have two parts, a **BASE UNIT** and **PREFIX**.
- The prefix tells you **how to modify the base unit** – either larger or smaller.
- Start with the base unit and then examine the prefix to determine what to do next.

# Know Your Prefixes!

PREFIX	PREFIX MEANING	EXPONENT
Kilo-	?	$10^3$
Centi-	?	$10^{-2}$
Milli-	?	$10^{-3}$
Micro-	?	$10^{-6}$
Nano-	?	$10^{-9}$

# Know Your Prefixes!

PREFIX	PREFIX MEANING	EXPONENT
Kilo-	<b>Thousand</b> (1,000.0)	$10^3$
Centi-	<b>Hundredth</b> 0.01 or 1/100	$10^{-2}$
Milli-	<b>Thousandth</b> 0.001 or 1/1,000	$10^{-3}$
Micro-	<b>Millionth</b> 1/1,000,000	$10^{-6}$
Nano-	<b>Billionth</b> (1/1,000,000,000)	$10^{-9}$

# Quick Conversions

5.0 m  $\rightarrow$  \_\_\_\_ mm

- Find out what your target unit's exponent is:

$$1.0 \text{ mm} = 10^{-3} \text{ m}$$

- Which unit is larger? What direction do you need to move the decimal point?

Since the unit you are converting to (mm) is smaller, then move decimal to the right by 3.

# Conversion Factors

1. Take what you have... Know what you need.

$$5 \text{ m} \rightarrow \text{___ mm}$$

2. Figure out your conversion factor.

$$\text{Since } 1,000 \text{ mm} = 1 \text{ m}$$

Your conversion factor would be:

$$\frac{1,000 \text{ mm}}{1 \text{ m}}$$

The unit you want to  
get rid of ALWAYS goes  
on the bottom



# Conversion Factors

3. Multiply by your conversion factor.

$$5 \text{ m} \times \frac{1,000 \text{ mm}}{1 \text{ m}} = ?$$

$$5 \cancel{\text{ m}} \times \frac{1,000 \text{ mm}}{1 \cancel{\text{ m}}} = 5,000 \text{ mm}$$

# More Difficult Conversion Factors

1. Sometimes you may want to use more than one conversion factor...

68 mm  $\rightarrow$  \_\_\_\_  $\mu\text{m}$  (micrometers)

$$1 \text{ m} = 1,000 \text{ mm}$$

$$1 \text{ m} = 1,000,000 \mu\text{m}$$

2. Create your conversion factors:

$$\frac{1 \text{ m}}{1,000 \text{ mm}}$$

$$\frac{1,000,000 \mu\text{m}}{1 \text{ m}}$$

# More Difficult Conversion Factors


## 3. Multiply

$$68 \text{ mm} \times \frac{1 \text{ m}}{1000 \text{ mm}} \times \frac{1,000,000 \text{ } \mu\text{m}}{1 \text{ m}} = ? \text{ } \mu\text{m}$$

$$68 \text{ ~~mm~~} \times \frac{1 \text{ ~~m~~}}{1000 \text{ ~~mm~~}} \times \frac{1,000,000 \text{ } \mu\text{m}}{1 \text{ ~~m~~}} = ? \text{ } \mu\text{m}$$

$$68 \times 1 \times 1,000 \text{ } \mu\text{m} = 68,000 \text{ } \mu\text{m}$$

# How To Calculate Solid Volume



Volume of a rectangular prism:

$$V = \text{Length} \times \text{Width} \times \text{Depth}$$

$$V = \underline{\hspace{2cm}} \text{ cm}^3 \text{ (cubic cm)}$$

# How To Calculate Liquid Volume: Use the Meniscus!

- Sometimes liquid molecules are more attracted to surface of container than each other.
- The **MENISCUS** is the bottom of the curve – that is the true volume!




## Common Units:

- Liters (L)
- Milliliters (mL)

# How To Read Volume: Use the Meniscus!

- Make sure to read graduated cylinder at **EYE LEVEL**.



- If you read at an angle you will either under- or overestimate.

# Types of Microscopes

## Optical Microscopes

- Examine specimen using your light and optical lenses
- Lower magnification
- Can view live organisms
- May require dyes to see detail
- Lower resolution

## Electron Microscopes

- Examine specimen with electrons
- Higher magnification
- Can view only dead organisms
- Requires specimen to be coated with heavy metals
- Higher resolution

# Resolution

- **Resolution** = the minimum distance between two adjacent objects required so they can be distinguished.
- Compound light microscope – 200 nm
- Transmission electron microscope – 0.1 nm



# Optical Microscopes

## Stereomicroscope

- Allows you to view the surface of an opaque 3D specimen




## Compound Light Microscope

- Allows you to view flat, translucent specimens




# The Compound Optical Microscope



- Light enters the microscope from the bottom.
- Travels through multiple lenses (hence the name compound).
- The lenses magnify & bend the light to your eye.
- **Note:** The image may be inverted if there is no projector lens.
- **Parfocal** = if focused in low power, it will be focused in higher power

# “Anatomy” of the Compound Optical Microscope



# Compound Optical Microscope: Objective Lenses

- Most of these microscopes have four lenses:
  - **Scanning** objective
  - **Low power** objective
  - **High power** objective
  - **Oil immersion** objective
- Eyepiece is also a lens (**ocular** lens)
- **Total magnification = ocular x objective**


# Electron Microscopes

## Transmission Electron Microscope

- Allows you to view flat, translucent specimens

## Scanning Electron Microscope

- Allows you to view the surface of a 3D specimen




# Methylene Blue / Iodine Solution

- Common staining solutions
- **Careful!** May stain clothing, hands, equip.
- It is an irritant - don't taste or inhale
- Avoid contact with eyes and skin.

# Today's Specimens

- Human epithelial cells
  - Epithelial tissue lines all the inner and outer surfaces of the body
- Onion epidermal cells
- Euglena
  - Unicellular protist
  - Capable of consuming other organisms (phagocytosis) for food or creating food by photosynthesis

# What Can You Identify in Euglena?



- **Nucleus** contains the cell's DNA
- **Eye spot** detects light
- **Contractile vacuoles** are multipurpose storage containers
- **Mitochondria** are the “power plants”
- **Flagella** rotate in a corkscrew to move the cell



# Handling Microscopes



**Hold microscopes with two hands.**  
One holding the arm and the other underneath the base.