Analytes	Siemens Multistix [®] 10SG Urinalysis Reagent Strips							
	Visual		Siemens Clinitek [®] Status [®] /Advantus [®] Urine Analyzers					
	LOT	Level 1 2101193K	LOT	Level 2 2102189B				
Glucose		Negative		100 - 2000 mg/dL or more		Negative		100 - ≥1000 mg/dL (Tr - 3+) or 5.6 - ≥56 mmol/L
Bilirubin		Negative		Small - Large (1+ - 3+)		Negative		Small - Large (1+ - 3+)
Ketones		Negative		5 - 160 mg/dL (Tr - Large)		Negative		5 - ≥160 mg/dL (Tr - Large) (Tr - 4+) or 0.5 - ≥16.0 mmol/L
Specific Gravity		1.020 - 1.030		1.000 - 1.020		1.015 - ≥1.030		≤1.005 - 1.020
Blood		Negative		Trace - Large (Tr - 3+)		Negative		Trace - Large (Tr - 3+) or ≤5 - Ca. 200 Ery/μL
pH		5.0 - 6.0		6.5 - 8.5		5.0 - 6.5		7.0 - ≥9.0
Protein		Negative		Trace - 2000 mg/dL or more (Tr - 4+)		Negative		Trace - ≥2000 mg/dL (Tr - 4+) or ≤0.5-20.0 g/L
Urobilinogen		Normal (0.2 - 1 mg/dL/EU/dL)		2 - 8 mg/dL		Normal (0.2 mg/dL)		1 - ≥8 mg/dL or 16 - ≥128 μmol/L
Nitrite		Negative		Positive		Negative		Positive
Leukocytes		Negative		Trace - Large (Tr - 3+)		Negative		Trace - Large (Tr - 3+) or ca 15-500 cells Leu/μL

Analytes	Siemens Clinitek [®] Microalbumin Reagent Strips							
	Visual		Siemens Clinitek [®] Status [®] Urine Analyzers					
	LOT	Level 1 2101193K	LOT	Level 2 2102189B				
Microalbumin		NA		NA		0 - 10 mg/L		30 - 150 mg/L
Creatinine		NA		NA		10 - 200 mg/dL or 4.4-17.7 mmol/L		100 - 300 mg/dL or 8.84 - ≥26.5 mmol/L

Analytes	Siemens Multistix PRO [®] 10LS Urinalysis Reagent Strips							
	Visual		Siemens Clinitek [®] Status [®] /Advantus [®] Urine Analyzers					
	LOT	Level 1 2101193K	LOT	Level 2 2102189B				
Glucose		Negative		100 - 2000 mg/dL or more		Negative		100 - 2000 mg/dL or more (5.5 - 111 mmol/L)
Protein - High		Negative		30 - 2000 mg/dL (1+ - 4+)		NA		NA
Protein - Low		Negative		15 mg/dL (Low)		NA		NA
Ketones		Negative		5 - 160 mg/dL (Tr - Large)		Negative		5 - ≥160 mg/dL (Tr - Large) (Tr - 4+) or 0.5 - ≥16.0 mmol/L
Specific Gravity		1.020 - 1.030		1.000 - 1.020		1.020 - ≥1.030		1.000 - 1.020
Creatinine		10 - 200 mg/dL		100 - 300 mg/dL		10 - 200 mg/dL		100 - 300 mg/dL or 8.84 - ≥26.5 mmol/L
pH		5.0 - 6.5		6.5 - 8.5		5.0 - 6.5		6.5 - ≥9.0
Blood		Negative		Trace - Large (Tr - 3+)		Negative - Trace		Trace - Large (Tr - 3+) or ≤5 - Ca. 200 Ery/μL
Nitrite		Negative		Positive		Negative		Positive
Leukocytes		Negative		Trace - Large (Tr - 3+)		Negative		Trace - Large (Tr - 3+) or ca 15-500 cells Leu/μL
Protein		NA		NA		Negative		15 - 2000 mg/dL or more dL (Tr - 4+) or ≤0.5-20.0 g/L
Protein:Creatinine Ratio		NA		NA		Normal		Abnormal

Analytes	Roche ChemStrip [®] Urinalysis Reagent Strips							
	Visual: 10SG, 10UA, 10MD		Roche Urisys 1100 [®] Urine Analyzers					
	LOT	Level 1 2101193K	LOT	Level 2 2102189B				
Specific Gravity		1.015 - 1.030		1.000 - 1.020		1.015 - 1.030		1.000 - 1.020
pH		5.0 - 6.0		7.0 - 9.0		5.0 - 6.5		7.0 - 9.0
Leukocytes		Negative		Trace - 2+		Negative		25 - 500 Leu/μL (Tr - 2+)
Nitrite		Negative		Positive		Negative		Positive
Protein		Negative		Trace - 500 mg/dL (Tr - 3+)		Negative		15 - 500 mg/dL (Tr - 3+) or .15-5.0 g/L
Glucose		Normal		50 - 1000 mg/dL		Normal		50 - ≥1000 mg/dL (Tr - 4+) or 2.8 - ≥56 mmol/L
Ketones		Negative		Small - Large (1+ - 3+)		Negative		5 - 150 mg/dL (Tr - 3+) or 0.5 - ≥15.0 mmol/L
Urobilinogen		Normal		1 - 12 mg/dL		Normal (0.2 mg/dL)		1 - ≥8 mg/dL (1+ - 4+) or 16 - ≥128 μmol/L
Bilirubin		Negative		1+ - 3+		Negative		1 - 6 mg/dL (1+ - 3+) or 17-100 μmol/L
Blood		Negative		Trace - about 250 Ery/μL		Negative		5 - 250 Ery/μL (Tr - 4+)

Confirmatory and Other Tests			
Test	LOT	Level 1 2101193K	Level 2 2102189B
Ictotest® Reagent Tablets		Negative	Positive
Specific Gravity (refractometer)		1.020 - 1.030	1.000 - 1.020
hCG (Consult*, Genzyme OSOM®, Clinitest®)		Negative (~0 mIU/mL)	Positive (~1000 mIU/mL)
pH paper		4.0 - 6.0	7.0 - 9.0
Total Protein (Sulfosalicylic Acid)		Negative	Positive

* hCG Dipstick, hCG Urine Cassette, hCG Combo Cassette

Microscopic Evaluation of Urine Sediment					
Method	Level 1			Level 2	
	Red Blood Cells	White Blood Cells	Casts	Crystals	Bacteria
	LOT 2101193K			2023-01-31	
Kova® Glasstic™ Slide 10 with grids	0 - 10 µL	0 - 10 µL	None	None	May Be Present
Non-grid slides -1.0 mL resuspended sediment	0 - 2 hpf	0 - 2 hpf	None	None	May Be Present
Slide & Coverslip -1.0 mL resuspended sediment	0 - 5 hpf	0 - 5 hpf	None	None	May Be Present
	LOT 2102189B			2023-03-31	
Kova® Glasstic™ Slide 10 with grids	30 - 100 µL	30 - 100 µL	None	Present	May Be Present
Non-grid slides -1.0 mL resuspended sediment	3 - 100 hpf	3 - 100 hpf	None	Present	May Be Present
Slide & Coverslip -1.0 mL resuspended sediment	3 - 100 hpf	3 - 100 hpf	None	Present	May Be Present

32448C_19 OTC 3/11/21

For in vitro diagnostic use


See package insert for proper use

Made in
the
USA

CLINIQA CORPORATION
495 Enterprise St.
San Marcos, CA 92078
USA

**SALES AND
TECHNICAL SUPPORT**
P: 800 728 5205
+1 760 744 1900
F: +1 760 571 5198
www.cliniga.com

**FOR ORDERS AND
CUSTOMER SERVICE**
P: 800 728 5205
+1 760 744 1900
F: +1 760 571 5197
sma_csgroup@bio-techne.com

CEpartner4U
Esdoornlaan 13
3951 DB Maarn, The Netherlands
www.cepartner4u.eu

INTENDED USE

FOR IN VITRO DIAGNOSTIC USE

CLINIQA Liquid QC DipScopics Urinalysis Control is intended for use as an assayed quality control material for various urinalysis reagent strips, qualitative hCG methods as well as confirmatory tests such as Ictotest[®]. This control is also intended to evaluate microscopic examination test procedures.

SUMMARY

The use of quality control materials to objectively monitor the accuracy and precision of procedures is well established in clinical laboratories. CLINIQA Liquid QC DipScopics Urinalysis Control is provided at two levels to assist in the monitoring of analytical systems within the clinical range.

PRODUCT DESCRIPTION

CLINIQA Liquid QC DipScopics Urinalysis Control is human urine-based containing constituents of human and animal origin as well as purified chemicals. Preservatives and stabilizers have been added to maintain product integrity. CLINIQA Liquid QC DipScopics Urinalysis Control is a ready-to-use liquid control requiring no reconstitution.

STORAGE AND STABILITY

CLINIQA Liquid QC DipScopics Urinalysis Control is stable until the expiration date on the label when stored unopened at 2-8°C. Do not freeze.

Once opened, bottles of control are stable for 18 months when stored tightly capped at 2-8°C or for 30 days at 20-25°C.

PROCEDURE

CLINIQA Liquid QC DipScopics Urinalysis Control should be treated the same as patient specimens and run in accordance with the instructions accompanying the test system (instrument, kit, or reagent) being used. Allow the product to reach room temperature prior to use.

QC materials should be used in accordance with local, state, and/or federal regulations or accreditation requirements.

INSTRUCTIONS

1. Replace the shipping cap on the control bottle with the spout cap supplied in the kit.
2. Gently invert the bottle several times before sampling to ensure homogeneity. This is important to obtain reproducible results.
3. Pour 12 mL of the controls into standard 15 mL centrifuge tubes.
4. **For samples to be utilized for dipstick testing, microscopic evaluation and hCG testing:** remove the sample to be used for the hCG and confirmatory tests after centrifugation before discarding the supernatant and without disturbing the sediment.
For urinalysis, microalbumin and creatinine testing: dip the reagent strip in the control-filled centrifuge tube as if it were a patient specimen.
 In accordance with the manufacturer's instructions, read the reagent strips either visually or with an instrument reader.

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For microscopic evaluation of urine sediment: use a standardized system in order to obtain reproducible results which will allow the reporting of abnormal sediment elements per unit of volume¹.

For hCG and confirmatory testing: in accordance with the hCG test kit manufacturer's instructions, use controls as if they were patient specimens.

For the control of confirmatory pregnancy test kits, be sure to pour a small amount of control into a separate receptacle and use the transfer pipette included in the kit to deliver the correct sample amount to the test device.

5. Close the spout cap immediately after use and store appropriately.

LIMITATIONS OF PROCEDURE

CLINIQA Liquid QC DipScopics Urinalysis Control should not be used past the expiration date on the bottle label. If there is evidence of microbial contamination or excessive turbidity in the product, discard the bottle.

CLINIQA Liquid QC DipScopics Urinalysis Control is a stabilized liquid product. To obtain consistent assay values, CLINIQA Liquid QC DipScopics Urinalysis Control requires storage and handling as detailed in STORAGE and STABILITY.

Accurate and reproducible results are dependent upon properly functioning instruments and reagents. The ranges given are intended only as a guideline. Each laboratory should establish their own acceptable ranges and tolerance limits based on their test system.

EXPECTED VALUES

The expected values printed in this insert were derived from replicate analyses of representative samples of the product and are specific to this lot of CLINIQA Liquid QC DipScopics Urinalysis Control testing data used to establish the expected values were derived from multiple instruments. All values have been assigned with manufacturer's reagents available at the time of assay. Subsequent instrument or reagent modifications may invalidate these expected values.

Refer to www.cliniqa.com/products for lot specific value assignment. To request a faxed or emailed copy of the value assignment, contact your local distributor or sma_csgr@bio-techne.com.

CAUTION

Human source material. Treat as potentially infectious.

Each donor unit used in the manufacture of this product has been tested by FDA accepted methods and found non-reactive for the presence of HBsAg and antibody to HIV-1/2, HCV and HIV-1 Ag. While these methods are highly accurate, they do not guarantee that all infected units will be detected. Because no known test method can offer complete assurance the hepatitis B virus, hepatitis C virus, human immunodeficiency virus (HIV) or other infectious agents are absent, all products containing human source material should be considered potentially infectious and handled with the same precautions used with patient specimens.

This product contains 0.09% sodium azide as a preservative. Sodium azide may react with lead and copper plumbing to form potentially explosive compounds. Flush with excess water upon disposal.

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