

---

**ZOLL**<sup>®</sup>

# M Series<sup>®</sup> Service Manual

---

**M**Series

---

The issue date for the M Series Service Manual (**REF** 9650-0450-01 Rev. R) is **March, 2013**.

If more than 3 years have elapsed since the issue date, contact ZOLL Medical Corporation to determine if additional product information updates are available.

Copyright © 2013 ZOLL Medical Corporation. All rights reserved. Rectilinear Biphasic, M Series, and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries.

Masimo is a trademark or registered trademark of Masimo Corporation in the United States and/or other countries.



**ZOLL Medical Corporation**

269 Mill Road  
Chelmsford, MA USA  
01824-4105



**ZOLL International Holding B.V.**

Newtonweg 18  
6662 PV ELST  
The Netherlands



**0123**

# TABLE OF CONTENTS

<b>PREFACE</b> .....	<b>V</b>
Overview .....	v
Safety Considerations .....	v
Additional Reference Material .....	vi
Conventions .....	vii
Service Policy Warranty .....	vii
Technical Service .....	vii
Technical Service for International Customers .....	viii
<b>CHAPTER 1 MAINTENANCE TESTS</b> .....	<b>1</b>
Overview .....	1
Before You Begin the Maintenance Tests .....	2
Equipment that You Need to Perform the M Series Maintenance Tests .....	3
Equipment You Need for the M Series Options Maintenance Tests .....	4
1.0 Physical Inspection of the Unit.....	5
2.0 Front Panel Button Test.....	7
3.0 3, 5, and 12 Leads Test.....	10
4.0 Power Supply Test (Optional) .....	11
5.0 Leakage Current Test.....	15
6.0 Paddles Test (If applicable).....	16
7.0 Heart Rate Display Test.....	18
8.0 Calibrating Pulses on Strip Chart Test .....	19
9.0 Notch Filter Test.....	20
10.0 Heart Rate Alarm Test.....	21
11.0 Defibrillator Self Test.....	23
12.0 Synchronized Cardioversion Test.....	25
13.0 Shock Test .....	26
14.0 Summary Report Test.....	29
15.0 Advisory Message Test .....	30
16.0 Pacer Test .....	31

17.0 SpO2 Monitor Test for SpO2 Option .....	33
18.0 EtCO2 Monitor Test (for EtCO2 Option).....	36
19.0 Temperature Test.....	38
20.0 IBP Monitoring Test.....	39
21.0 NIBP Transducer Calibration Test .....	44
22.0 NIBP Monitor Test.....	47
23.0 NIBP Volume Leak Test with Fluke Biomedical NIBP Analyzer.....	49
24.0 Bluetooth Test.....	51
<b>CHAPTER 2 TROUBLESHOOTING .....</b>	<b>53</b>
Overview .....	53
Troubleshooting .....	54
Zoll M Series Error Messages .....	58
<b>CHAPTER 3 DISASSEMBLY PROCEDURES .....</b>	<b>75</b>
Overview .....	75
Required Equipment .....	76
Parts That May Need Replacing After Disassembly .....	76
Safety Precautions .....	77
Overview of Modules .....	78
1. Removing the ZIF Keeper .....	81
2. Removing the Front Panel .....	82
2A.Removing the Display .....	83
2B.Removing the Control Board .....	84
3. Removing the Upper Housing Assembly .....	86
4. Removing the System Board Assembly .....	87
5. Removing the Battery Interconnect Board Assembly .....	90
6. Removing the High Voltage/Charger Assembly .....	91
7. Removing the High Voltage Module Assembly .....	94
8. Removing the High Voltage Capacitor Assembly .....	96
9. Removing the System Interconnect Board .....	98
10.Removing the Printer/Recorder Motor .....	100
11.Removing the Lower Housing Assembly .....	102
12.Removing the Print Head Assembly .....	103
13.Removing the PCMCIA Card Slot Assembly .....	104

14.Removing the Paddle Release Latch .....	105
<b>CHAPTER 4 REPLACEMENT PARTS .....</b>	<b>107</b>
Overview .....	107
Replacement Parts .....	108
Field Replacement Parts .....	112
<b>CHAPTER 5 FUNCTIONAL DESCRIPTION .....</b>	<b>115</b>
Overview .....	115
Main System Board .....	116
Main System Board Functions .....	118
Power Supply .....	120
ECG Front End .....	121
Multifunction Electrode (MFE)/PADS (System Board and High Voltage Module) .....	121
CPU and EPU .....	121
High Voltage Module .....	122
Defibrillator Charging and Discharging .....	123
High Voltage Capacitor Monitor .....	124
Pacer/Defibrillator Control Signals .....	125
Internal Discharge Resistor Module .....	127
AC/DC Charger Module .....	127
System Interconnect Module .....	127
Stripchart Recorder .....	128
PCMCIA Slots .....	128
Front Panel and Controls PWBA .....	128
M Series Options .....	128
Isolated Power Supply Module .....	129
12 Lead Option .....	129
Pulse Oximetry (SpO2) .....	130
End Tidal Carbon Dioxide (EtCO2) .....	130
Biphasic Waveform .....	131
<b>APPENDIX .....</b>	<b>135</b>
Overview .....	135
Interconnect Diagram for the M Series Monophasic Unit .....	136
Interconnect Diagram for the M Series Biphasic Unit .....	137
Interconnect Diagram for the M Series CCT Biphasic Unit .....	138

**INDEX .....139**

# Preface

## Overview

ZOLL® Medical Corporation's M Series® Service Manual is intended for the service technician whose responsibility is to identify malfunctions and/or make repairs at the subassembly level. The Zoll M Series Service Manual has five main sections and one appendix.

**Preface**—Contains safety warnings and an overview of the manual's contents. Be sure to review this section thoroughly before attempting to use or service the M Series unit.

**Chapter 1—Maintenance Tests** explains how to check the defibrillator's performance using a series of recommended checkout procedures to be conducted every six months.

**Chapter 2—Troubleshooting** provides a listing of the procedures and error messages to help the service technician detect faults and repair them.

**Chapter 3—Disassembly Procedures** describes step-by-step procedures for removing subassemblies from the M Series unit.

**Chapter 4—Replacement Parts List** displays a complete list of ZOLL part numbers for field replaceable parts available for the M Series unit, allowing the service person to identify and order replacement parts from ZOLL.

**Chapter 5—Functional Description** provides technical descriptions for the M Series major subassembly modules.

**Appendix A**—M Series Operator's Manual.

## Safety Considerations

The following section describes general warnings and safety considerations for operators and patients. Service technicians should review the safety considerations prior to servicing any equipment and read the manual carefully before attempting to disassemble the unit. Only qualified personnel should service the M Series unit.

Federal (U.S.A.) law restricts this unit for use by or on the order of a physician.

Safety and effectiveness data submitted by ZOLL Medical Corporation to the Food and Drug Administration (FDA) under section 510(K) of the Medical Device Act to obtain approval to market is based upon the use of ZOLL accessories such as disposable electrodes, patient cables and batteries. The use of external pacing/defibrillation

electrodes and adapter units from sources other than ZOLL is not recommended. ZOLL makes no representations or warranties regarding the performance or effectiveness of its products when used in conjunction with pacing/defibrillation electrodes and adapter units from other sources. If unit failure is attributable to pacing/defibrillation electrodes or adapter units not manufactured by ZOLL, this may void ZOLL's warranty.

Only qualified personnel should disassemble the M Series unit.

---

### **WARNING!**

**This unit can generate up to 4500 volts with sufficient current to cause lethal shocks.**

---

All persons near the equipment must be warned to “STAND CLEAR” prior to discharging the defibrillator.

Do not discharge the unit's internal energy more than three times in one minute or damage to the unit may result.

Do not discharge a battery pack except in a Base PowerCharger<sup>4x4</sup> or compatible ZOLL Battery Charging/Testing unit.

Do not use the M Series in the presence of flammable agents (such as gasoline), oxygen-rich atmospheres, or flammable anesthetics. Using the unit near the site of a gasoline spill may cause an explosion.

Do not use the unit near or within puddles of water.

### **NOTE**

The M Series is protected against interference from radio frequency emissions typical of two-way radios and cellular phones (digital and analog) used in emergency service/public safety activities. Users of the M Series should assess the unit's performance in their typical environment of use for the possibility of radio frequency interference from high-power sources. Radio Frequency Interference (RFI) may be observed as shifts in monitor baseline, trace compression, or transient spikes on the display.

## Additional Reference Material

In addition to this guide, there are several other components to the Zoll M Series documentation. They include:

- Operator's Guide - A comprehensive reference work that describes all the user tasks needed to operate the M Series.
- Configuration Guide - Describes the M Series features and functions whose operation can be customized by authorized users.



## Conventions

---

**WARNING!**

Warning statements describe conditions or actions that can result in personal injury or death.

---

**CAUTION**

Caution statements describe conditions or actions that can result in damage to the unit.

---

**NOTE**

Notes contain additional information on using the defibrillator.

## Service Policy Warranty

In North America: Consult your purchasing agreement for terms and conditions associated with your warranty. Outside of North America, consult ZOLL authorized representative.

In order to maintain this warranty, the instructions and procedures contained in this manual must be strictly followed. For additional information, please call the ZOLL Technical Service Department 1-800-348-9011 in North America.

## Technical Service

If the ZOLL M Series unit requires service, contact the ZOLL Technical Service Department:

Telephone: 1-978-421-9655; 1-800-348-9011

Fax 1-978-421-0010

Have the following information available for the Technical Service representative:

- Unit serial number.
- Description of the problem.
- Department where equipment is used.
- Sample chart recorder strips documenting the problem, if applicable.
- Purchase Order to allow tracking of loan equipment.
- Purchase Order for a unit with an expired warranty.

If the unit needs to be sent to ZOLL Medical Corporation, obtain a service order request number from the Technical Service representative. Return the unit in its original container to:

ZOLL Medical Corporation

269 Mill Road

Chelmsford, Massachusetts 01824-4105

Attn: Technical Service Department

Telephone: 1-800-348-9011; 1-978-421-9655 FAX: 978-421-0010

## Technical Service for International Customers

International customers should return the unit in its original container to the nearest authorized ZOLL Medical Corporation Service Center. To locate an authorized service center, contact the International Sales Department at ZOLL Medical at the above address.

Units are available on loan while your unit is being repaired.

# Chapter 1

## Maintenance Tests

### Overview

The M Series has two checkout procedures: the operator's shift checklist and the extensive six-month maintenance tests checkout procedures.

Because the M Series units must be maintained ready for immediate use, it is important for users to conduct the Operator's Shift Checklist procedure at the beginning of every shift. This procedure can be completed in a few minutes and requires no additional test equipment. (See the ZOLL M Series Operator's Guide for the Operator's Shift Checklist.)

A qualified biomedical technician must perform a more thorough maintenance test checkout every six months to ensure that the functions of the M Series unit work properly. This chapter describes the step by step procedures for performing the six month maintenance test checkout. Use the checklist at the back of this document (ZOLL M Series Maintenance Tests Checklist) to record your results of the maintenance tests.

This chapter describes the following M Series maintenance tests:

1. Physical Inspection of the Unit
2. Front Panel Button
3. 3, 5, and 12 Leads
4. Power Supply
5. Leakage Current
6. Paddles
7. Heart Rate Display
8. Calibrating Pulses on Strip Chart
9. Notch Filter
10. Heart Rate Alarm

11. Defibrillator Self Test
12. Synchronized Cardioversion
13. Shock
14. Summary Report
15. Advisory Message
16. Pacer
17. SpO<sub>2</sub> Monitor
18. EtCO<sub>2</sub> Monitor
19. Temperature Test
20. IBP Monitoring Test
21. NIBP Transducer Calibration Test
22. NIBP Monitor Test
23. NIBP Volume Leak Test with Fluke Biomedical NIBP Analyzer
24. Bluetooth Test

## Before You Begin the Maintenance Tests

- Assemble the tools or specialized testing equipment listed in the “Equipment You Need to Perform the Maintenance Tests” section shown below.
- Keep an extra fully charged ZOLL M Series compatible battery available.
- Schedule an hour to conduct the entire maintenance test.
- Photocopy the checklist at the back of this document and use the copy to record your results. As you conduct each step of a procedure, mark the Pass/Fail/NA check boxes on your checklist and then save it for your maintenance file.
- Perform the tests in the order presented.
- Perform all the steps of each test procedure.
- Complete all the steps of the procedure before evaluating the test results.

## Equipment that You Need to Perform the M Series Maintenance Tests

This section lists equipment that we use to perform the maintenance tests that we describe in this chapter. You can substitute an equivalent device for a listed device; however, *not all simulators and analyzers will produce the same results*. Be sure to follow the manufacturer's recommendations for conducting the maintenance tests.

We recommend the use of the following equipment when performing M Series Maintenance Tests:

- ZOLL Medical Electrode Adapter from Fluke Biomedical (ZOLL Part Number 3010-0378).
- QED 6 Defibrillator Analyzer
- MFC Test Port Connector, 1004-0053-99
- Fluke Biomedical 601 Pro Series International Safety Analyzer.
- ECG Simulator; 12 Lead Simulator for 12 Lead test (e.g., Symbio CS1201).
- Stop watch.
- Standard series II PC flash memory cards.
- 1 red miniature alligator to miniature alligator test lead.
- 1 black miniature alligator to miniature alligator test lead.
- DC power supply (15 amp minimum).
- 0.1 ohm 1% resistor (1/4 watt or greater).
- 1000 ohm 1% resistor (1/4 watt or greater).
- Fluke 75 Multimeter or equivalent.
- Printer Paper.
- Battery.
- AC line cord.
- 3 lead, 5 lead and 12 lead ECG cables. (12 lead cable needed if 12 lead option is installed).

## Equipment You Need for the M Series Options Maintenance Tests

- Fluke Biomedical Index 2 SpO<sub>2</sub> Simulator or equivalent. (For SpO<sub>2</sub> units only.)
- SpO<sub>2</sub> cable and sensor (if option is installed).
- EtCO<sub>2</sub> cable, and CAPNOSTAT 3 Mainstream cable with airway adapter, or CAPNOSTAT 3 Sidestream cable with cannula (if option is installed).
- Novametrix Capnostat Simulator Yb 1265/7100
- Paddles (if used).
- )
- Fluke Biomedical BP Pump 2 NIBP Monitor Analyzer (For NIBP units only) with NIBP cable and cuff (if NIBP option is installed), or
- Fluke Biomedical Cufflink Analyzer (if NIBP option is installed)
- IBP Temperature Simulator 9100-0402-TF
- DNI Nevada 214B Patient Simulator, with M Series interface cable

**NOTE** The Fluke Biomedical BP Pump NIBP Monitor Analyzer and the Fluke Biomedical BP Pump NIBP Analyzer use different technologies for testing NIBP monitors and therefore, the manual provides two different procedures for performing the NIBP Volume Leak test with each of these types of test equipment.

## 1.0 Physical Inspection of the Unit

Tools Needed

Battery.

Test Setup

None.

Observe this...		Pass/Fail/NA
1.1	<b>Housing</b> Is the unit clean and undamaged?	0 0
1.2	Does the unit show signs of excessive wear?	0 0
1.3	Does the handle work properly?	0 0
1.4	Does the recorder door open and close properly?	0 0
1.5	Are input connectors clean and undamaged?	0 0
1.6	Are there any cracks in the housing?	0 0
1.7	Do the front panel or selector switches have any damage or cracks?	0 0
1.8	Are there any loose housing parts?	0 0
1.9	<b>Cables</b> Are all cables free of cracks, cuts, exposed or broken wires?	0 0
1.10	Are all cable bend/strain reliefs undamaged and free of excessive cable wear?	0 0
1.11	<b>Paddles</b> Do the adult and pedi plates have major scratches or show signs of damage?	0 0 0
1.12	Do the adult shoes slide on and off easily to expose the covered pedi plates?	0 0 0
1.13	Are the paddles clean (e.g., free of gel) and undamaged? (if applicable)	0 0 0

Observe this...		Pass/Fail/NA
1.14	<b>Battery</b> Place battery in battery well.	0 0
1.15	Is the battery seated in the battery well correctly?	0 0
Record your results on the Maintenance Test Checklist.		



## 2.0 Front Panel Button Test

Tools Needed

QED 6 Defibrillator Analyzer

Test Setup

- Install strip chart paper into the recorder tray.
- Install the battery in the unit or connect the A/C power cord to the unit and then plug the cord into an electrical outlet.
- Connect the universal cable and ECG cable (3 lead, 5 lead, or 12 lead) to the QED 6 Analyzer (or equivalent).

	Do this...	Verify that...	Pass/Fail/NA
2.1	Turn the selector switch to <b>MONITOR</b> . (For AED units, turn the selector switch to <b>ON</b> and select Manual mode.)	The unit sounds 4 beep tones. PADS and MONITOR display on the monitor. <b>Note:</b> PADS is a factory default setting.	o o
2.2	Set the QED 6 analyzer to NSR of 120 BPM. To check the size of the ECG waveform, press the <b>SIZE</b> button.	As you press the SIZE button five times (0.5, 1.0, 1.5, 2.0, 3.0), the size of the ECG waveform appropriately changes on the display.	o o
2.3	Press the <b>ALARM SUSPEND</b> button.	Alarm symbol changes from disabled to enabled. If the alarm sounds, press the ALARM SUSPEND button to turn it off. The alarm will only be suspended for 90 seconds at this point. Press and hold the <b>ALARM SUSPEND</b> button for 3 seconds to disable alarms.	o o o
2.4	Press the <b>RECORDER</b> button.	The strip chart paper moves out of the unit from the paper tray. Check that the correct time, date, ECG lead annotation and waveform are recorded on the paper. (Set Time and Date, if necessary.)	o o o
2.5	Open the paper compartment door, then press the <b>RECORDER</b> button.	The <i>CHECK RECORDER</i> message appears on the monitor.	o o o
2.6	Close the paper compartment door, then press the <b>RECORDER</b> button.	The strip chart paper flows out of paper tray and the monitor no longer displays the <i>CHECK RECORDER</i> message.	o o o

	Do this...	Verify that...	Pass/Fail/NA
2.7	Press the <b>RECORDER</b> button.	The strip chart paper stops flowing out of the paper tray.	0 0 0
2.8	Press the <b>VOLUME</b> softkey.	The volume bar graph displays. <b>Note:</b> The QRS tone is on or off. There is no gradual change in volume. If equipped, voice prompts are gradual. Note: The voice volume has 5 settings. Setting 3 is in the mid-range.	0 0
2.9	Press the <b>INC</b> softkey.	The bar graph appears on the display, indicating an increase in volume. <b>Note:</b> This action does not increase the volume of the unit's audio prompts.	0 0
2.10	Press the <b>DEC</b> softkey.	The bar graph appears on the display, indicating a decrease in volume. Verify that pressing the <b>DEC</b> softkey until only the last bar appears silences the unit's audio prompts, otherwise the volume of the audio prompts does not decrease.	0 0
2.11	Press the <b>CONTRAST</b> button.	The monitor screen displays the Contrast Menu.	0 0
2.12	Press the <b>INC</b> softkey.	The screen contrast and bar graph increases on the display.	0 0 0
2.13	Press the <b>DEC</b> softkey.	The screen contrast and bar graph decreases on the display.	0 0
2.14	Press the <b>SUMMARY</b> button.	The monitor screen displays the Summary menu, showing the summary report options.	0 0 0
2.15	Press the <b>CODEMARKER</b> button.	The monitor screen displays the Code marker menu.	0 0 0
2.16	Connect to AC or DC current and install the battery. Turn the unit off.	The CHARGER ON indicator lights (either the Green or Amber indicator lights).	0 0 0

	Do this...	Verify that...	Pass/Fail/NA
2.17	Remove the battery.	The charger light alternately flashes green and amber.	0 0
2.18	Replace the battery and the turn unit on.	The amber charger light illuminates.	0 0
2.19	Press the <b>ANALYZE</b> button.	The <i>SELECT DEFIB MODE</i> message appears on the monitor (for manual devices).	0 0 0
2.20	Connect the universal cable to the QED 6. Move the selector switch to DEFIB. Select 2J. Press the <b>CHARGE</b> button.	The display shows that the unit is charging. The <b>SHOCK</b> button lights when the unit is charged, and the Ready tone for DEFIB sounds.	0 0 0
2.21	Press and hold the <b>ENERGY SELECT</b> down arrow.	The unit discharges internally and selected energy decrements to 1J.	0 0 0
2.22	Press and release the <b>ENERGY SELECT</b> up arrow 19 times.	The following selections appear: Verify the following energy amounts display incrementally: Biphasic:2-10, 15, 20,30, 50, 75, 100,120, 150, 200J Monophasic(DSW): 2-10, 15,20,30,50,75, 100, 150, 200, 300, 360J	0 0 0
2.23	Press the <b>CHARGE</b> button.	The display shows the unit charged up to 200J (360J monophasic) and the <b>SHOCK</b> button lights.	0 0 0
2.24	Press the <b>SHOCK</b> button.	The unit discharges and the <b>SHOCK</b> button is no longer illuminated. A 15 second strip chart automatically prints, displaying the number of joules delivered (if equipped with recorder and configured to print event).	0 0 0
Record your results on the Maintenance Test Checklist.			

### 3.0 3, 5, and 12 Leads Test

Tools Needed: 3 lead, 5 lead, and 12 lead cables. Test each cable separately. QED 6 Defibrillator Analyzer.

Test Setup: The M Series unit must be configured to display *ECG LEAD OFF* message.  
Connect the lead wires appropriate to each test to the QED 6 Defibrillator Analyzer.

	Do this...	Observe this...	Pass/Fail/NA
3.1	Turn the selector switch to <b>MONITOR</b> . Select leads.	The monitor displays the <i>NO ECG LEADS OFF</i> message.	0 0 0
3.2	Disconnect one lead from the simulator.	The monitor displays the <i>ECG LEAD OFF</i> message within 3 seconds (if configured).	0 0 0
3.3	Reconnect the lead. Repeat step 3.2 with the remaining leads.	The <i>ECG LEAD OFF</i> message appears when the lead is disconnected and clears the lead is reconnected.	0 0 0
3.4	Repeat 3.2 and 3.3 for 5 lead and 12 lead cables.	<p><b>Note:</b> If heart rate alarm sounds, press and hold the <b>ALARM SUSPEND</b> button for 4 seconds to disable the alarms.</p> <p><b>Note:</b> When testing the 12 lead cable, the <i>ECG LEAD OFF</i> message displays when you pull off a limb lead. When you pull off a V lead, the <i>ECG VX LEAD OFF</i> message displays, where "X" is the number between 1 and 6.</p>	0 0 0
Record your results on the Maintenance Tests Checklist.			

## 4.0 Power Supply Test (Optional)

Tools Needed:

- 1 red miniature alligator to miniature alligator test lead.
- 1 black miniature alligator to miniature alligator test lead.
- DC power supply (15 Amp minimum).
- 0.1 ohm 1% resistor (1/4W or greater).
- 1000 ohm 1% resistor (1/4W or greater).
- Fluke 75 multimeter or equivalent.

Test Setup:

Make sure the unit and power supply are turned off.

Connect one end of the black lead to the “-” terminal in the battery well.

Connect the other end of the black lead to the “-” terminal of the power supply.

Connect the red lead to “+” terminal socket of the battery well. Use the middle pin with the plastic guard around it.


Connect the other end of the red lead to the “+” terminal of the power supply.

Set the power supply voltage to 7V.

---

**CAUTION** Be sure to connect the power supply properly to the M Series battery well terminals or damage to the unit may result. Do NOT raise the power supply voltage above 12V.


---



	Do this...	Verify that...	Pass/ Fail
4.1	Turn the selector switch to <b>MONITOR</b> . (For AED units, turn the selector switch to <b>ON</b> and select Manual mode.)	The unit does not turn on.	0 0
4.2	Turn the selector switch to the <b>OFF</b> position.		
4.3	Adjust the power supply voltage to 10.3V and turn the selector switch to <b>MONITOR</b> (for AED units, turn the selector switch to the <b>ON</b> position).	The unit turns on.	0 0
4.4	Low Battery Test Set voltage to 9.8V.	The unit does not display the <i>LOW BATTERY</i> message.	0 0
4.5	Set voltage to 9.4V.	The unit displays the <i>LOW BATTERY</i> message within 30 seconds.	0 0
4.6	Shut Down Voltage Test Set voltage to 8.5V.	The unit shuts off within 30 seconds.	0 0
Record your results on the Maintenance Tests Checklist.			

Test Setup

Remove red lead from power supply and connect to 0.1ohm resistor.  
 Connect other end of resistor to “+” terminal of power supply using a second red label.  
 Connect multimeter across the resistor.  
 Set voltage scale (if DVM is not autoranging) to 220mV.



	Do this....	Observe this...	Pass/ Fail/
4.7	System Current Test Set power supply to 10.3V.		

	Do this....	Observe this...	Pass/ Fail/
4.8	Turn the selector switch to MONITOR. (For AED units, turn the selector switch to ON and select Manual mode.	Voltage across resistor should be 80 mV or less (<800 mA of ON current). Note: Without optional parameters. a) With green screen or LCD and no options <80mV. b) With yellow screen and no options <81mV. c) With yellow screen and SpO2 <104mV. d) With yellow screen and voice recording <91mV. e) With yellow screen, voice recording, and SpO2 <114mV. f) All devices with EtCO2 <121mV.	<input type="radio"/> <input type="radio"/>
4.9	Turn unit off.		<input type="radio"/> <input type="radio"/>
Record your results on the Maintenance Test Checklist			

Off Current Test

Test Setup

- Remove 0.1 ohm resistor and replace with 1K
- Connect DMM across resistor.
- Set voltage scale to DCV.
- Measure voltage across resistor.

4.10	Measure across resistor with unit turned off.	Voltage should be less than 450mV (<450uA of current).	<input type="radio"/> <input type="radio"/>
Record your results on the Maintenance Test Checklist			



## 5.0 Leakage Current Test

Tools Needed:

See the manufacturer's instructions or supplied specifications for the leakage tester you use.

Test Setup:

See the manufacturer's instructions or supplied specifications for the leakage tester you use. Repeat leakage test with accessories: MFC, external paddles, and anterior/posterior paddles.

**Table 1:**

Maximum Leakage Acceptance Limits		
	Normal Condition	Single Fault Condition*
ECG	10 $\mu$ A	50 $\mu$ A
Universal	100 $\mu$ A	100 $\mu$ A
Earth	500 $\mu$ A	1000 $\mu$ A
*Single fault considered AC mains on applied part.		

## 6.0 Paddles Test (If applicable)

Tools Needed:

QED 6 Defibrillator Analyzer.

Test Setup:

If applicable, connect the universal cable to the paddles and place the paddles in paddle wells.

	Do this...	Verify that...	Pass/Fail/NA
6.1	Turn the selector switch to <b>DEFIB (For AED units, turn the selector switch to ON and select Manual mode.)</b> Press and hold the <b>ENERGY DOWN</b> button on the sternum paddle.	The energy selection decreases to 1J.	0 0
6.2	Press and release the <b>ENERGY UP</b> button on the sternum paddle for each setting.	The joules selection increases as follows: Biphasic 1-10, 15, 20, 30, 50, 75, 100, 120, 150, 200J. Monophasic (DSW): 1-10, 15, 20, 30, 50, 75, 100, 150, 200, 300, 360.	0 0
6.3	Press and release the <b>RECORDER</b> button on the sternum paddle.	The recorder starts printing.	0 0 0
6.4	Press and release the <b>RECORDER</b> button again.	The recorder stops printing.	0 0
6.5	Select 30J using the paddle <b>ENERGY</b> button. Press the <b>CHARGE</b> button on the Apex paddle.	The unit charges to 30J, then the red LED charge indicator illuminates and the charge tone sounds. <b>Note:</b> The front panel shock button does not illuminate.	0 0
6.6	Press and release the <b>APEX SHOCK</b> button.	The unit does not discharge.	0 0
6.7	Press and release the <b>STERNUM SHOCK</b> button.	The unit does not discharge.	0 0
6.8	Press and hold both paddles <b>SHOCK</b> buttons.	The unit discharges. The <i>TEST OK</i> message displays and the red LED turns off. If configured, the recorder prints a strip chart.	0 0

	Do this...	Verify that...	Pass/Fail/ NA
Record your results on the Maintenance Tests Checklist.			

## 7.0 Heart Rate Display Test

Tools Needed:

- ECG Cable (3, 5 or 12 leads).
- QED 6 Defibrillator Analyzer.

Test Setup:

Turn the selector switch to MONITOR (For AED units, turn the selector switch to ON and select Manual override.)  
Press **LEAD** button until “T” displays.

Connect the ECG leads to the QED 6 Defibrillator Analyzer.

Connect the ECG cable to the unit.

	Do this...	Verify that...	Pass/Fail/NA
7.1	Set the QED 6 Defibrillator Analyzer to 120BPM.	The Heart Rate displays as 120 +/- 2 bpm.	<input type="radio"/> <input type="radio"/> <input type="radio"/>
Record your results on the Maintenance Tests Checklist.			

## 8.0 Calibrating Pulses on Strip Chart Test

Tools Needed: None

Test Setup: None.

	Do this...	Verify that...	Pass/Fail/NA
8.1	Press the <b>RECORDER</b> button.		
8.2	Press and hold <b>SIZE</b> button to activate the calibration signal.	The strip chart displays a signal of 300 ppm with an amplitude of 10 mm +/- 1 mm. The signal also appears on the video display.	o o o
Record your results on the Maintenance Tests Checklist.			

## 9.0 Notch Filter Test

Tools Needed:

QED 6 Defibrillator Analyzer.

Test Setup:

Connect the ECG cable to the QED 6 Defibrillator Analyzer.

Connect the ECG cable to the unit.

	Do this...	Verify that...	Pass/Fail/NA
9.1	Turn the selector switch to <b>MONITOR</b> mode. (For AED units, turn the selector switch to <b>ON</b> and select Manual mode.)		
9.2	Select lead I, size 3x. Select 60Hz (or 50 Hz for a 50Hz unit) on the QED 6 Defibrillator Analyzer.		
9.3	Press <b>RECORDER</b> button.	Verify that the waveform amplitude on the strip chart is less than 1.5 mm.	○ ○ ○
9.4	Press <b>RECORDER</b> button to stop printing.		
Record your results on the Maintenance Tests Checklist.			

## 10.0 Heart Rate Alarm Test

Tools Needed:

QED 6 Defibrillator Analyzer

	Do this...	Verify that...	Pass Fail/ NA
10.1	Turn the selector switch to <b>MONITOR</b> mode. (For AED units, turn the selector switch to <b>ON</b> and select Manual mode.) Connect the ECG leads to the QED 6 Defibrillator Analyzer. Set the QED 6 Analyzer to 120 BPM and the defibrillator to lead II.	The monitor screen displays the <i>Lead II</i> message and the waveform for NSR ECG at 120 BPM +/- 2 .	0 0 0
10.2	Press <b>ALARMS</b> button.	The monitor displays the alarm menu.	0 0 0
10.3	Press <b>SELECT PARAM</b> softkey until ECG HR displays.	The cursor scrolls through parameters.	0 0 0
10.4	Press <b>INC&gt;</b> for state.	The cursor scrolls through ENABLE, AUTO and DISABLE.	0 0 0
10.5	Press <b>DEC&gt;</b> for state.	The cursor scrolls through ENABLE, DISABLE, AND AUTO.	0 0 0
10.6	Press <b>INC&gt;</b> until ENABLE displays.	ENABLE appears.	0 0 0
10.7	Press <b>ENTER</b> softkey, then press the <b>NEXT FIELD</b> softkey to select the heart rate limit.	The cursor scrolls to Low field. Press RETURN softkey to exit out of Alarm Menu.	0 0 0
10.8	Press <b>ALARM SUSPEND</b> button.	No alarm sounds.	0 0 0

	Do this...	Verify that...	Pass Fail/ NA
10.9	Remove an EEG lead wire from the QED 6 analyzer.	The alarm symbol flashes and the heart symbol stops flashing. The ECG LEAD OFF alarm tone sounds. Recorder prints a stripchart showing a low heart rate, if enabled.	0 0 0
10.10	Reattach ECG Lead wire to the QED 6 Analyzer and hold the <b>ALARM SUSPEND</b> button on unit for 4 seconds.	The unit displays the alarm symbol with an "X" through it. The heart symbol flashes with each QRS wave.	0 0 0
10.11	Press the <b>ALARM SUSPEND</b> button to enable the alarms.	The alarm symbol does not have an "X" through it.	0 0 0
10.12	Set QED 6 Analyzer to 160 BPM or higher.	The Heart Rate Value is highlighted, the alarm tone sounds, and the alarm and the heart symbol both flash.	0 0 0
10.13	Press the <b>ALARM SUSPEND</b> button on the M Series unit.	The alarm is suspended for 90 seconds. The unit displays the alarm symbol with an "X" through it, and the heart symbol flashes with each QRS wave. Note: Software version 35.00-36.XX, the alarm is suspended until the alarm button is press again, or the fault is cleared.	0 0 0
10.14	Press and hold <b>ALARM SUSPEND</b> for 4 seconds to disable alarms.		0 0 0
Record your results on the Maintenance Tests Checklist.			



## 11.0 Defibrillator Self Test

### **WARNING!** **SHOCK HAZARD!**



**TAKE THE NECESSARY PRECAUTIONS TO GUARD AGAINST SHOCK OR INJURY BEFORE YOU START CONDUCTING THE DEFIBRILLATOR TESTS.**

**Keep hands and all other objects clear of the multi-function cable connections and defibrillator analyzer when discharging the defibrillator.**

**Before you discharge the defibrillator, warn everyone near the equipment to STAND CLEAR.**

---

#### **CAUTION**

Do NOT internally discharge the unit more than 3 times in 1 minute. Note that the rapid repetition of internal discharges at more than 30 Joules may damage the unit.

---

#### Tools Needed:

- MFC Test Port Connector 1004-0053-99 with universal cable.
- MFC Test Adaptor Connector (Fluke Biomedical Part Number 3010-0378 or equivalent).
- QED 6 Defibrillator Analyzer.
- ECG Cable.
- Stop watch.

#### Test Setup:

Ensure the unit is turned off

Ensure the ECG cable is connected to the unit and analyzer.

**NOTE:** The universal cable should not be connected to any equipment at the beginning of this test.

	Do this...	Verify that...	Pass/ Fail
11.1	Turn the selector switch to <b>DEFIB</b> mode. (For AED units, turn the selector switch to <b>ON</b> and select Manual mode.) Set leads to PADS.	<i>The unit displays the CHECK PADS/POOR PAD CONTACT message.</i>	O O
11.2	Connect the universal cable to the MFC test port.	<i>DEFIB PAD SHORT</i> message displays.	O O
11.3	Select energy level of 100J and press the <b>CHARGE</b> button.	The charge time is greater than 2 seconds and less than 10 seconds and the unit displays the message, <i>SELECT 30J FOR TEST.</i>	O O
11.4	Press the <b>SHOCK</b> button.	The unit does not discharge. The unit displays the message, <i>DEFIB PAD SHORT.</i>	O O
11.5	Set energy level to 30J.	The unit discharges internally.	O O
11.6	Press the <b>CHARGE</b> button.	The unit charges to 30J and displays message, <i>DEFIB 30J READY.</i> The charge ready tone sounds.	O O
11.7	Press and hold <b>SHOCK</b> button.	The unit discharges and displays the number of joules delivered and the <i>TEST OK</i> message. For example, the message at the top of the printed strip chart would read as follows: 30 JOULES TEST OK. TEST_CUR=10-14A DEFIB_IMPED=0.  The impedance value may range from 0 to 5 ohms. If configured, the recorder prints a strip chart.	O O
Record your results on the Maintenance Tests Checklist.			

## 12.0 Synchronized Cardioversion Test

Tools Needed:

QED 6 Defibrillator Analyzer.

Test Setup:

Connect the universal cable via the adapter (D.N.I #3010-0378) to the defibrillator analyzer.  
 Select Sync on analyzer and set ECG on analyzer to 60 -120 BPM.

	Do this...	Verify that...	Pass/ Fail
12.1	Press <b>LEAD</b> button to select <b>PADS</b> and Size X1.		
12.2	Press the <b>SYNC</b> softkey on the defibrillator. Enter synchronized cardioversion timing test mode on the defibrillator analyzer.	Sync markers appear on display. <b>Note:</b> Sync markers display on the monitor. The sync marker appears as a down arrow over the ECG R-wave peaks on strip chart and display.	o o
12.3	Select 200J (360J for Monophasic units), and then press the <b>CHARGE</b> button. When the <b>SHOCK</b> button lights, press and hold the <b>SHOCK</b> button.	On the analyzer display, the R-wave to shock delay is less than 60 milliseconds.	o o
Record your results on the Maintenance Tests Checklist.			

## 13.0 Shock Test

Tools Needed:

QED 6 Defibrillator Analyzer, Stop Watch.

Test Setup:

Connect the universal cable via the adapter (D.N.I #3010-0378) to the defibrillator analyzer.

Ensure that a fully charged battery is installed in the unit.

**NOTE** If your M Series AED does not have manual override capability, do not perform this test.

	Do this...	Verify that...	Pass/Fail/NA
13.1	Turn the selector switch to <b>DEFIB</b> mode. (For AED units, turn the selector switch to <b>ON</b> and select Manual mode.)		
13.2	Press the <b>ENERGY SELECT</b> down arrow until 5J displays.	The unit displays <i>DEFIB 5J SEL.</i>	0 0 0
13.3	Press the <b>CHARGE</b> button. Wait for the <b>SHOCK</b> button to illuminate.	The unit displays <i>DEFIB 5J RDY.</i>	0 0 0
13.4	Press the <b>SHOCK</b> button.	The unit discharges 3J-5J into the simulator. <b>Note:</b> The displayed rhythm may change shape for 30 seconds before it returns to an original rhythm. This is caused by the operation of the adaptive bandwidth defibrillator recovery circuit.)	0 0 0
13.5	Press the <b>ENERGY SELECT</b> up arrow until the unit displays 50J.	The unit displays the message, <i>DEFIB 50J SEL.</i>	0 0 0
13.6	Press the <b>CHARGE</b> button. Wait for the <b>SHOCK</b> button to illuminate.	The unit displays the message, <i>DEFIB 50J RDY.</i>	0 0 0
13.7	Press the <b>SHOCK</b> button.	The unit discharges 41J-55J (monophasic 43J-57J) into the simulator.	0 0 0

	Do this...	Verify that...	Pass/Fail/NA
13.8	Press the <b>ENERGY SELECT</b> up arrow until the unit displays <i>100J</i> .	The unit displays the message, <i>DEFIB 100J SEL.</i>	0 0 0
13.9	Press the <b>CHARGE</b> button. Wait for the <b>SHOCK</b> button to illuminate.	The unit displays the message, <i>DEFIB 100J RDY.</i>	0 0 0
13.10	Press the <b>SHOCK</b> button.	The unit discharges 83J-122J (monophasic 85J-115J) into the simulator.	0 0 0
13.11	Press the <b>ENERGY SELECT</b> up arrow until the unit displays <i>200J</i> .	The unit displays the message, <i>DEFIB 200J SEL.</i>	0 0 0
13.12	Press the <b>CHARGE</b> button. Wait for the <b>SHOCK</b> button to illuminate.	The unit displays the message, <i>DEFIB 200J RDY.</i>	0 0 0
13.13	Press the <b>SHOCK</b> button.	The unit discharges 182J-247J (monophasic 170J-230J) into the simulator.	0 0 0
13.14	(Monophasic Units) Press the <b>ENERGY SELECT</b> up arrow unit the unit displays 300J.	The unit displays the message, <i>DEFIB 300J SEL.</i>	0 0 0
13.15	Press the <b>CHARGE</b> button. Wait for the <b>SHOCK</b> button to illuminate.	The unit displays the message, <i>DEFIB 300J RDY.</i>	0 0 0
13.16	Press the <b>SHOCK</b> button.	The unit discharges 255J-345J into the simulator.	0 0 0
13.17	(Monophasic Units) Press the <b>ENERGY SELECT</b> up arrow unit the unit displays 360J.	The unit displays the message, <i>DEFIB 360J SEL.</i>	0 0 0
13.18	Press the <b>CHARGE</b> button. Wait for the <b>SHOCK</b> button to illuminate.	The unit displays the message, <i>DEFIB 360J RDY.</i>	0 0 0
13.19	Press the <b>SHOCK</b> button.	The unit discharges 306J-414J into the simulator.	0 0 0

	Do this...	Verify that...	Pass/Fail/NA
13.20	(Biphasic Units) Press the <b>CHARGE</b> button and start timing with a stopwatch. Stop timing when the <b>SHOCK</b> button lights.	The charge time is between 3-6 seconds.	0 0 0
13.21	Press the <b>SHOCK</b> button.	On the strip chart, the Patient Current is between 20-24 Amps, and the Defib Impedance is between 46-54 Ohms. Note: This will occur if configured to Auto-Generate Strips on event.	0 0 0
13.22	(Monophasic Units) Press the <b>CHARGE</b> button and start timing with a stopwatch. Stop timing when the <b>SHOCK</b> button illuminates.	The charge time is between 4-8 seconds.	0 0 0
13.23	Press the <b>SHOCK</b> button.	On the strip chart, the Patient Current is between 20-24 Amps, and the Defib Impedance is between 46-54 Ohms. Note: This will occur if configured to Auto-Generate Strips on event.	0 0 0
13.24	(AED units only.) Disconnect the universal cable from the analyzer.	The unit will display <i>CHECK PADS</i> , and give an audio prompt if configured for audio.	0 0 0
Record your results on the Maintenance Tests Checklist.			

## 14.0 Summary Report Test

Tools Needed:

QED 6 Defibrillator Analyzer.

Test Setup:

Connect the universal cable to the QED 6 analyzer.

If you are using paddles, place the paddles on the analyzer's discharge plates.

	Do this...	Verify that...	Pass/ Fail
14.1	Press and hold the <b>SUMMARY</b> softkey for 4 to 8 seconds to display erase options. Press the <b>ERASE SUMMARY</b> softkey to erase any previously stored data.	The unit displays the message, <i>ERASING REPORT</i> .	0 0
14.2	Set selector switch to <b>DEFIB</b> . Select 200J (360J for Monophasic) using the <b>ENERGY SELECT</b> button, and press the <b>CHARGE</b> button. When charged, press the <b>SHOCK</b> button to discharge into the defibrillator analyzer.	The unit successfully discharges and prints a strip chart.	0 0
14.3	Wait 18 seconds, then press the <b>Code Marker</b> softkey. Press the <b>CPR</b> softkey.	The Code Markers appear.	0 0
14.4	Turn the unit off. Wait 10 seconds and then turn the unit on. Press the <b>SUMMARY</b> softkey, then press the <b>PRINT CHART</b> softkey, then press PRINT ALL.	The summary report prints. The report displays the correct date, time, 1 shock delivered and Code Marker event.	0 0
Record your results on the Maintenance Tests Checklist.			

## 15.0 Advisory Message Test

Tools Needed:

QED 6 Defibrillator Analyzer.

Test Setup:

Connect the universal cable via the adapter (D.N.I #3010-0378).


Attach the M Series to the QED 6 Defibrillator Analyzer.

	Do this...	Verify that...	Pass/ Fail
15.1	Turn the selector switch to <b>DEFIB</b> mode. (For AED units, turn the selector switch to <b>ON</b> .)		
15.2	Select VF (ventricular fibrillation) on the QED 6, then press the <b>ANALYZE</b> button.	The unit displays the following messages: <i>ANALYZING ECG.</i> <i>STAND CLEAR.*</i> <i>SHOCK ADVISED.*</i> <i>PRESS SHOCK.*+</i> *AED's audio prompts are standard. Advisory audio prompts are user configurable. +If configured for auto charge.	O O
15.3	Press the <b>SHOCK</b> button.	The unit discharges.	O O
15.4	Select the <b>NSR</b> (normal sinus rhythm) on the simulator, then press the <b>ANALYZE</b> button.	The unit displays the following messages: <i>ANALYZING ECG.</i> <i>STAND CLEAR.*</i> <i>NO SHOCK ADVISED.*</i> *AED's audio prompts are standard.	O O
Record your results on the Maintenance Tests Checklist.			



## 16.0 Pacer Test

Tools Needed:



QED 6 Defibrillator Analyzer.

**NOTE** The following tests are to be performed only on M Series units equipped with the optional pacing function.

The pacer output can be measured using an oscilloscope set to DC coupling connected across a load resistor. (See diagram in column for universal cable connector polarity.) The load resistor is a 100 ohm, 5 watt or greater. The pacer output is a positive going pulse, 40 +/- 2 ms duration with an amplitude of 0.1 volt per milliamp of selected output (e.g., 40 milliamps of selected output has an amplitude of 4 +/- 0.5 volts the specified tolerance displayed on the oscilloscope).

If an external non-invasive pacer analyzer is being used, then follow the manufacturer’s guidelines for measuring the frequency (ppm), output (mA) and the pulse width measured in milliseconds. Note that the analyzer pace load resistor must be less than 250 ohms.

Test Setup:

Connect the universal cable from the M Series to the QED 6 Defibrillator Analyzer.

	Do this...	Verify that...	Pass/ Fail
16.1	Turn the Main Selector knob of the unit to Pacer Mode. Set the <b>PACER OUTPUT</b> to 14 mA and disconnect MFC connector from the QED 6 analyzer.	The unit displays the <i>CHECK PADS</i> and <i>POOR PAD CONTACT</i> messages displays, and the pace alarm is active.	0 0
16.2	Reconnect the universal cable to the QED 6 analyzer. Press <b>Clear Pace Alarm</b> softkey.	<i>CHECK PADS</i> and <i>POOR PAD CONTACT</i> message disappears. The pace alarm is cleared.	0 0
16.3	Set rate to 180 ppm; output to 0mA.	No output appears on the QED 6 analyzer.	0 0
16.4	Increase the output to 40mA.	Output on the QED 6 analyzer is 40mA +/- 5mA.	0 0
16.5	Increase the output to 120mA.	Output on the QED 6 analyzer is 120mA +/- 6mA.	0 0

	Do this...	Verify that...	Pass/ Fail
16.6	Increase the output to 140mA.	Output on the QED 6 analyzer is 140mA +/- 7 mA.	<input type="radio"/> <input type="radio"/>
16.7	Check the pulse width.	Pulse width is 40mS +/-2mS.	<input type="radio"/> <input type="radio"/>
16.8	Decrease the output to 60mA. Decrease the rate to 30ppm.	Pacer rate on QED 6 analyzer is 29-31 ppm.	<input type="radio"/> <input type="radio"/>
16.9	Increase the rate to 180ppm.	Pacer rate on QED 6 analyzer is 177-183 ppm.	<input type="radio"/> <input type="radio"/>
16.10	Connect the ECG cable to the M Series and QED 6 analyzer. Select the ECG at 60 BPM on the QED 6 analyzer.	The unit displays ECG at 60 BPM with no stimulus markers.	<input type="radio"/> <input type="radio"/>
16.11	Press the <b>Async</b> Pace softkey.	The unit displays ECG at 60 BPM with stimulus markers, and displays the <i>Async pace</i> message.	<input type="radio"/> <input type="radio"/>
16.12	Turn off the QED 6 analyzer. Set Pacer Rate to 100ppm. Press the <b>RECORDER ON</b> button.	The unit displays pace stimulus markers every 15mm +/- 1mm.	<input type="radio"/> <input type="radio"/>
16.13	Press and hold 4:1 button.	The unit displays pace stimulus markers every 60 mm +/- 1.5 mm.	<input type="radio"/> <input type="radio"/>
Record your results on the Maintenance Tests Checklist.			

## 17.0 SpO<sub>2</sub> Monitor Test for SpO<sub>2</sub> Option

Tools Needed:

- Masimo® Reusable Sensor.
- Masimo Patient Cable
- Fluke Index Series 2 SpO<sub>2</sub> Simulator.

Test Setup:

Connect the universal cable to the MFC test plug.

DO NOT connect the ECG cable to the simulator.

Connect the Masimo Patient Cable and attach the Masimo sensor to the patient cable.

Connect the Masimo sensor to the finger simulation post.

Place a fully charged battery into the battery well or connect to ac power (dc power, if equipped).

Ensure that the SpO<sub>2</sub> Simulator is off.

	Do this...	Verify that...	Pass/ Fail
17.1	Turn the selector switch to <b>MONITOR</b> . (For AED units, turn the selector switch to <b>ON</b> and select Manual mode.)	The SpO <sub>2</sub> saturation percentage appears as a dashed line on the monitor.	o o
17.2	Wait ten seconds. Turn on the SpO <sub>2</sub> simulator. Press the <b>SIM</b> softkey on the Index SpO <sub>2</sub> Simulator, then press the <b>MAN</b> softkey.	The unit displays the SpO <sub>2</sub> <i>PULSE SEARCH</i> message.	o o
17.3	Press the <b>02+</b> or <b>02-</b> softkey of the simulator until the SpO <sub>2</sub> output is at 98%.	The SpO <sub>2</sub> reading of 98 +/- 1% appears on the E Series monitor. <b>Note:</b> You may need to wait up to 2 minutes for the information to appear on the ZOLL display.	o o

	Do this...	Verify that...	Pass/ Fail
17.4	Using the Index SpO <sub>2</sub> Simulator, press the <b>BPM+</b> or <b>BPM-</b> softkey until the heart rate is 230 BPM.	The SpO <sub>2</sub> simulator screen displays an SpO <sub>2</sub> rate of 230 BPM. Note that you may need to wait up to 2 minutes for the information to appear on the ZOLL display. The M Series monitor displays an SpO <sub>2</sub> saturation of 96-100%. The M Series monitor displays a heart rate of 226-234 BPM.	0 0
17.5	Using the Index SpO <sub>2</sub> Simulator, press the <b>BPM-</b> softkey until the heart rate is 50 BPM.	The SpO <sub>2</sub> simulator screen displays an SpO <sub>2</sub> saturation of 96-100%. The M Series monitor displays a heart rate of 46-54 BPM.	0 0 0 0
17.6	Using the Index SpO <sub>2</sub> Simulator, press the <b>02+</b> softkey until the SpO <sub>2</sub> output is at 72%.	The SpO <sub>2</sub> simulator screen displays an SpO <sub>2</sub> saturation of 70-74%. The M Series monitor displays a heart rate of 46-54 BPM.	0 0 0 0
17.7	Press <b>Wave 2</b> softkey. Select the SpO <sub>2</sub> waveform. <u>For CCT units:</u> Press the <b>Traces</b> softkey, and then select <b>SET TRACE 2</b> , and press <b>Enter</b> . <b>Select SPO2, press Enter.</b>	Plethysmographic waveform appears on the ZOLL display.	0 0
17.8	Press <b>RECORDER</b> .	The plethysmographic waveform prints on the strip chart paper.	0 0
17.9	Using the Index SpO <sub>2</sub> Simulator, press the <b>BPM-</b> softkey until the heart rate is at 230 BPM.	The SpO <sub>2</sub> simulator screen displays an SpO <sub>2</sub> saturation rate of 70-74%. The M Series monitor displays a heart rate in the heart position of 226-234 BPM.	0 0
17.10	Remove the Masimo® patient cable.		

	Do this...	Verify that...	Pass/ Fail
Record your results on the Maintenance Tests Checklist.			

## 18.0 EtCO<sub>2</sub> Monitor Test (for EtCO<sub>2</sub> Option)

Tools Needed: Novametrix Capnostat Simulator Tb 1265/7100.

Test Setup:

- Install the battery
- On the Novatrix Simulator, set the folling
- Set inspired CO<sub>2</sub> to OFF
- Set% CO<sub>2</sub> to 0
- Set Sensor Location to ZERO CELL
- Set Source Current to NORMAL
- Set CO<sub>2</sub> mode to CONTINUOUS. Set Temperature to NORMAL
- Attach simulator to the unit.

	Do this...	Verify that...	Pass/ Fail
18.1	Turn the selector switch to MONITOR mode (For AED units, turn the selector switch to ON, and select Manual mode.) Attach the EtCO <sub>2</sub> simulator to the M Series input connection.	CO <sub>2</sub> SENSO WARMUP message displays. Note: You may need to wait up to 5 minutes for the warm-up message to disappear.	0 0
18.2	Perfom manual zeroing. Press PARAM softkey, select EtCO <sub>2</sub> , press ENTER, and then press ZERO softkey, then Enter.	Verify "ZEROING CO <sub>2</sub> " Adapter" is displayed, then changes to "ZERO DONE" when complete.	0 0
18.3	On the Novametrix Simulator, set SENSOR LOCATION to AA CELL, set% CO <sub>2</sub> to 10% CO <sub>2</sub> mode to RESPIRATION.	Verify the EtCO <sub>2</sub> reading of 74-84mmHg appears on the M Series display. Note that you may need to wait up to 10 seconds for the unit to stabilize.	0 0

	Do this...	Verify that...	Pass/ Fail
18.4	On the Novamatrix simulator, set% CO <sub>2</sub> to 5%.	The EtCO <sub>2</sub> reading of 34-42mmHg displays on the M Series monitor. The Respiration Rate (RR) of 22-24 displays on the M Series monitor. Note that you may need to wait up to 10 seconds for the unit to stabilize.	0 0
18.5	Press WAVE 2 softkey <u>For CCT units:</u> Press the <b>Trace</b> softkey, select <b>SET TRACE 2</b> , and press <b>Enter</b> . Select EtCO <sub>2</sub> , press <b>Enter</b> .	Verify the EtCO <sub>2</sub> waveform displays.	0 0
18.6	Press the RECORDER button.	Verify the EtCO <sub>2</sub> waveform prints Note that the CO <sub>2</sub> waveform is displayed and printed at 12.5 millimeters per second scale.	0 0
Record your results on the Maintenance Tests Checklist.			

## 19.0 Temperature Test

Tools Needed: IBP Temp Simulator 9100-0402-TF (or equivalent<sup>1</sup>)

Test Setup: None

	Do this...	Verify that...	Pass/Fail/NA
19.1	Connect the 2-channel Y adapter to the temperature connector (T1/T2) on the M Series CCT Unit.		
19.2	Connect the Temperature Probe Simulator cables to the 2-channel Y adapter.		
19.3	On the Temperature Simulator, set the T1 temperature channel to 59.2 degrees F (15 degrees C).	Verify that the value displayed for the T1 temperature channel between 59.0 degrees to 59.4 degrees F (14.9 to 15.1 degrees C).	O O O
19.4	On the Temperature Simulator, set the T1 temperature channel to 112.8 degrees F (45 degrees C).	Verify that the value displayed for the T1 temperature channel between 112.6 degrees to 112.8 degrees F (44.9 to 45.1 degrees C).	O O O
19.5	On the Temperature Simulator, set the T2 temperature channel to 59.2 degrees F (15 degrees C).	Verify that the value displayed for the T2 temperature channel between 59.0 degrees to 59.4 degrees F (14.9 to 15.1 degrees C).	O O O
19.6	On the Temperature Simulator, set the T2 temperature channel to 112.8 degrees F (45 degrees C).	Verify that the value displayed for the T2 temperature channel between 112.6 degrees to 112.8 degrees F (44.9 to 45.1 degrees C).	O O O
Record your results on the Maintenance Test Checklist.			

1. If you are using a Temperature Simulator that does not have the same setting above, use the closest value available. Verify the display temperature value is within 0.2 F or 0.1 C of the simulator setting



## 20.0 IBP Monitoring Test

Tools Needed:

DNI Nevada 214B Patient Simulator with M Series interface cables, FOGG Transducer Simulator

Test Setup:

Set the BP output sensitivity scale to 5 uV/V/mmHg on the Patient Simulator  
 Connect the Patient Simulator to the IBP connector on the M Series CCT unit

	Do this...	Verify that...	Pass/Fail/NA
20.1	Set the selector switch on the M Series CCT unit to MONITOR and press the PARAM softkey on the front panel to display the Parameters menu.	Verify the parameter menu displayed.	0 0 0
20.2	Select the IBP option and press the Enter softkey.	Verify the unit displays the IBP menu.	0 0 0
20.3	Press the Label softkey to display the IBP Label menu.	Verify the unit displays the IBP label menu.	0 0 0
20.4	Select the Label P1 option and press the Enter softkey.	Verify the M Series CCT unit displays the Label Options menu.	0 0 0
20.5	Select the P1 option and press the Enter softkey.	Verify the unit assigns the P1 label to the first pressure channel (P1) and returns to the IBP menu.	0 0 0
20.6	Press the Label softkey to display the IBP Label menu.	Verify the unit displays the IBP Label menu.	0 0 0
20.7	Select the Label P2 option and press the Enter softkey.	Verify the M Series CCT unit display the Label Options menu.	0 0 0
20.8	Select the ART option and press the Enter sofkey.	Verify the unit assigns the ART label to the second pressure channel (P2) and return to the IBP menu. Verify that the M Series CCT unit displays advisory messages to zero the pressure channels.	0 0 0 0 0 0

	Do this...	Verify that...	Pass/Fail/NA
20.9	On the Patient Simulator, set BP channel 1 to an atmospheric pressure of 0 mmHg (waveform 10) and press EXECUTE. Set BP channel 2 to duplicate BP channel 1 (waveform 22) and press EXECUTE.		
20.10	On the unit go to the IBP menu (press the Param softkey and select the IBP option), then press the Zero softkey to display the Zero menu, select the Zero Both option and press the Enter softkey.	<p>Verify the unit zeroes both pressure channels and displays the P1 ZEROED and P2 ZEROED messages within 15 seconds.</p> <p>Verify that the values displayed for the P1 and ART pressure channels are accurate within <math>\pm 3</math> mmHg of atmospheric pressure (that is, 0 mmHg):</p> <p>systolic pressure: -3 to 3 mmHg                      diastolic pressure: -3 to 3 mmHg                      mean pressure: -3 to 3 mmHg</p> <p>Note: The unit displays pressure values for each pressure channel within 15 seconds.</p>	<p>O O O</p> <p>O O O</p>
20.11	On the Patient Simulator, set to an arterial pressure of 120/80 mmHg (waveform 00) and press EXECUTE.	<p>Verify the values displayed for the P1 and ART pressure channels:</p> <p>systolic pressure: 116 - 124 mmHg                      diastolic pressure: 77 - 83 mmHg                      mean pressure: 90 - 98 mmHg</p> <p>Note: The unit displays pressure values for each pressure channel within 15 seconds.</p>	O O O

	Do this...	Verify that...	Pass/Fail/NA
20.12	Zero the IBP channels (see step 20.9).		
20.13	Press the Traces softkey.	Verify the unit displays the Traces menu.	0 0 0
20.14	Select the Set Trace 2 option and press the Enter softkey.	Verify the M Series CCT unit displays the Trace Options menu.	0 0 0
20.15	Verify the M Series CCT unit displays the Trace Options menu.	Verify the unit assigns the P1 waveform to the second trace display area and returns to the main menu.	0 0 0
20.16	Press the Traces softkey.	Verify the unit displays the Traces menu.	0 0 0
20.17	Select the Set Trace 3 option and press the Enter softkey.	Verify the unit displays the Trace Options menu.	0 0 0
20.18	Select the ART option and press the Enter softkey.	Verify the unit assigns the ART waveform to the third trace display area and returns to the main menu.	0 0 0
20.19	Go to the IBP menu (press the Param softkey and select the IBP option).	Verify the unit displays the IBP menu.	0 0 0
20.20	Press the Range softkey.	Verify the unit displays the IBP Range menu.	0 0 0
20.21	Select the P1 Range option and press the Enter softkey.	Verify the unit displays the Range Options menu.	0 0 0
20.22	Select the 0-30 mmHg option and press the Enter softkey.	Verify the unit displays the P1 waveform on a 0-30 mmHg scale and returns to the IBP menu.	0 0 0
20.23	Press the Range softkey.	Verify the unit displays the IBP Range menu.	0 0 0
20.24	Select the ART option and press the Enter softkey.	Verify the unit displays the Range Options menu.	0 0 0
20.25	Select the 0-30 mmHg option and press the Enter softkey.	Verify the unit displays the ART waveform on a 0-30 mmHg scale and returns to the IBP menu.	0 0 0

	Do this...	Verify that...	Pass/Fail/NA
20.26	On the Patient Simulator, set BP channel 1 to a static pressure of 20 mmHg (waveform 11) and press EXECUTE.		
20.27	Press the RECORDER button.	Verify that there is no deviation from the baseline greater than 1 mm for the waveforms printed on stripchart and the unit displays the waveforms for 20mmHg on the monitor screen, and prints a stripchart of the traces.	0 0 0
20.28	Set the input and output resistance to 350 ohms on the Transducer Simulator, and set the gain to x10. Connect the Transducer Simulator cable to the P1 connector on the unit and on the Pressure Transducer Simulator, set the pressure polarity to positive, and select a pressure of 200 mmHg.		
20.29	Zero pressure channel 1.	Verify the unit calibrates pressure channel 1 against a positive pressure of 200 mmHg and displays the P1 ZEROED message.	0 0 0
20.30	On the Pressure Transducer Simulator, set the pressure polarity to positive, and select a pressure of 500 mmHg.	Verify that the value displayed for the P1 pressure channel is in the range of 291 - 300 mmHg. (Note: value will be displayed within 15 seconds.)	0 0 0
20.31	Connect the simulator cable to the P2 connector on the unit. On the Pressure Transducer Simulator, set the pressure polarity to positive, and select a pressure of 200 mmHg.		
20.32	Zero ART channel.	Verify the unit calibrates pressure channel 1 against a positive pressure of 200 mmHg and displays the ART ZEROED message.	0 0 0

	Do this...	Verify that...	Pass/Fail/NA
20.33	On the Pressure Transducer Simulator, set the pressure polarity to positive, and select a pressure of 500 mmHg.	On the unit, verify that the value displayed for the ART P2 pressure channel is in the range of 291 - 300 mmHg. (Note: value will be displayed within 15 seconds.)	0 0 0
20.34	Connect the Pressure Transducer Simulator to the first IBP channel connector (P1). On the Pressure Transducer Simulator, set the pressure polarity to negative, and select a pressure of -200 mmHg.		
20.35	Zero the pressure channel.	The unit calibrates the pressure channel against a negative pressure of -200 mmHg and displays the P1 ZEROED message.	0 0 0
20.36	On the Pressure Transducer Simulator, select a negative pressure of -250 mmHg.	On the unit, verify that the value displayed for the P1 pressure channel is in the range of -47 to -50 mmHg. (Note: value will be displayed within 15 seconds.)	0 0 0
20.37	Connect the Pressure Transducer Simulator to the second IBP channel connector (P2). On the Pressure Transducer Simulator, set the pressure polarity to negative, and select a pressure of -200 mmHg.		
20.38	Zero the pressure channel.	The M Series CCT unit calibrates the pressure channel against a negative pressure of -200 mmHg and displays the P2 ZEROED message.	0 0 0
20.39	On the Pressure Transducer Simulator, select a negative pressure of -250 mmHg.	Verify that the value displayed for the P2 pressure channel is in the range of -47 to -50 mmHg. Note: Value will be displayed within 15 seconds.	0 0 0

## 21.0 NIBP Transducer Calibration Test

The NIBP module’s pressure transducers are factory-calibrated prior to shipment. However, you can perform a two-point calibration procedure periodically to ensure accurate pressure measurements.

This procedure is optional at 6 months, but should be performed annually or every 10,000 readings, whichever comes first.

Tools Needed:

Fluke NIBP simulator (the values and procedure that we provide are specific to the BP Pump 2).

Test Setup:

Connect the Analyzer’s simulator hose to the NIBP connector on the M Series unit.

Configure the NIBP Analyzer to simulate cuff pressure. For example, on the Fluke Biomedical BP Pump:

Press the MODE button three (3) times to go into Tests mode.

Press the SELECT button once to access the Pressure Simulator screen.

These instructions apply to the Fluke Biomedical BP Pump; for equivalent devices, follow the manufacturer’s instructions.

Make sure the ECG cable is not connected to the M Series unit.

If the SpO2 option is installed, make sure that the Masimo cable is NOT connected to the M Series unit.

	Do this...	Verify that...	Pass/ Fail
21.1	Turn the Selector Switch to <b>OFF</b> . After 10 seconds, press and hold the fourth softkey from the left and turn the Selector Switch to <b>MONITOR</b> (for AED units, turn the Selector switch to <b>ON</b> ).	The M Series unit powers on in the NIBP Service Mode.	○ ○
21.2	Press the <b>NIBP Calib</b> softkey.	The M Series unit displays the NIBP Transducer Calibration Screen.	○ ○
21.3	<i>On the NIBP simulator</i> set the pressure parameter to 0 mmHg.	The NIBP Simulator displays a pressure reading of 0 mmHg.	○ ○

	Do this...	Verify that...	Pass/ Fail
21.4	On the M Series unit, press the <b>Set Low</b> softkey to calibrate the transducer to a 0 mmHg pressure reading.	The NIBP pressure transducer registers its voltage output at a known pressure of 0 mmHg. The field adjacent to the 0 mmHg value changes to PASS. <b>Note:</b> If the M Series unit displays a FAIL reading, verify the NIBP simulator's pressure setting and connection to the M Series and repeat the step.	0 0
21.5	On the NIBP Simulator, set the pressure parameter to 250 mmHg, and then press <b>Start</b> on the NIBP simulator.	The NIBP Simulator displays a pressure reading of 250 mmHg.	0 0
21.6	On the M Series unit, press the <b>Set High</b> softkey to calibrate the transducer to a 250 mmHg pressure reading.	The NIBP pressure transducer registers its voltage output at a known pressure of 250 mmHg. The field adjacent to the 250 mmHg value changes to PASS. <b>Note:</b> If the M Series unit displays a FAIL reading, verify the NIBP Analyzer's pressure setting and connection to the M Series and repeat the step.	0 0
21.7	On the NIBP simulator, set the pressure parameter to stimulate a different cuff pressure (for example, 205 mmHg).	The NIBP simulator displays the specified pressure reading.	0 0
21.8	On the M Series unit, press the <b>Read Cuff</b> softkey when the simulator's Measured screen reaches 205.	The value that the M Series unit displays is accurate within +/- 3mmHg of the pressure parameter value set on the NIBP simulator.	0 0
21.9	On the M Series unit, press the <b>EXIT</b> softkey twice.	The M Series returns to the main NIBP Service Mode screen, then to normal Monitor mode operation.	0 0
Record your results on the Maintenance Tests Checklist.			

---

**Warning!** NIBP transducer calibration can affect clinical readings of the NIBP parameter. Ensure that the NIBP Transducer Calibration procedure is performed correctly, followed by an NIBP Monitor Test for verify proper operation.

---



## 22.0 NIBP Monitor Test

The NIBP monitor test verifies the repeatability of the systolic, diastolic, and mean blood pressure measurements, as well as the patient pulse rate calculation.

Tools Needed:

Fluke NIBP simulator (the values and procedure that we provide are specific to the BP Pump 2).

**Note:** The primary purpose of an NIBP simulator is to reproduce a pressure profile similar to a live patient to be used for testing the repeatability and functionality of the system. There are many different NIBP simulators on the market, each manufacturer uses a different method to develop their algorithm. Consequently, readings from different simulators may vary. To test for repeatability, you should first establish the offset<sup>1</sup> of your simulator. The offset value should then be used to determine the expected values. NIBP simulators cannot be used as a source for testing the accuracy of the non-invasive blood pressure measurements of devices such as the ZOLL M Series monitor/defibrillator.

Test Setup:

Connect the Analyzer's simulator hose to the NIBP connector on the M Series unit.

Set the following parameters on the NIBP Analyzer:<sup>2</sup>

**Table 1-1:**

Parameter	Value
Systolic pressure	120 mmHg
Diastolic pressure	80 mmHg
Mean pressure	93 mmHg <sup>3</sup>
Heart pressure	80 bpm

Make sure the ECG cable is not connected to the M Series unit.

If the SpO2 option is installed, make sure that the Masimo cable is NOT connected to the M Series unit.

<sup>1</sup> NIBP Simulators may produce a reading on the NIBP monitor that is shifted from the simulator's setting. The offset value must be established based on a statistical sample of monitors and readings. Please contact ZOLL Technical Support if you require assistance establishing the offset of the simulator and test set-up that you are utilizing.

<sup>2</sup> If you are using the Fluke® BiomedicalCuftlink, you must change the shift value of the Blood Pressure Envelope to +3 on the Pressure Curve Adjust Menu.

<sup>3</sup> Not all simulators have a setting of 93mmHg. Check the simulators user's manual for recommendations.

	Do this...	Verify that...	Pass/ Fail
22.1	Turn the selector switch to <b>MONITOR</b> mode. (For AED units, turn the selector switch to <b>ON</b> and select Manual mode.)		
22.2	Ensure that the LEADS parameters is set to PADS (default). If necessary, press the <b>LEADS</b> button to cycle through the values to select PADS.	The M Series unit displays <i>PADS</i> in the Lead selection field on the monitor.	O O
22.3	Press the <b>NIBP</b> button on the E Series front panel.	The M Series displays the following measurements: <ul style="list-style-type: none"> <li>• Systolic pressure (115 - 125 mmHg)</li> <li>• Diastolic pressure (75 - 85 mmHg)</li> <li>• Mean pressure (88 - 98 mmHg).</li> </ul>	O O
22.4	Press the <b>SUMMARY</b> softkey, then press the <b>TREND</b> softkey, followed by the <b>NIBP Trend</b> softkey.	The M Series unit displays a summary of the NIBP measurements, including the pulse rate reading (in the range of 77 - 83 bpm).	O O
Record your results on the Maintenance Tests Checklist.			

<sup>4</sup> These values only apply for test set-ups utilizing the BPPump 2 Simulator. Variations of the test set-up or different simulators may produce readings outside the provided values and will require end-user facility to establish the appropriate offset and tolerances. Please contact ZOLL Technical Support if you require assistance establishing the offset of your simulator and test set-up.

## 23.0 NIBP Volume Leak Test with Fluke Biomedical NIBP Analyzer

The volume leak test verifies the integrity of the pneumatic system on the NIBP module. This test is optional at 6 months, but should be performed annually or every 10,000 readings, whichever comes first.

Tools Needed:

Fluke NIBP simulator (the values and procedure that we provide are specific to the BP Pump 2).

Test Setup:

Connect the Analyzer’s simulator hose to the NIBP connector on the M Series unit.

Configure the NIBP Analyzer for the volume leak test. For example, on the Fluke Biomedical BP Pump:

- Press the MODE button three times to go into Tests mode.
- Press the SELECT button twice to access the volume leak test.

Make sure the ECG cable is not connected to the M Series unit.

If the SpO2 option is installed, make sure that the Masimo cable is NOT connected to the M Series unit.

	Do this...	Verify that...	Pass/ Fail
23.1	Turn the Selector Switch to <b>OFF</b> . After 10 seconds, press and hold the fourth softkey from the left and turn the Selector Switch to <b>MONITOR</b> (For AED units, turn the selector switch to <b>ON</b> ).	The M Series unit powers on in the NIBP Service Mode.	0 0
23.2	Press the <b>Leak Test</b> softkey.	The M Series unit displays the NIBP Leak Test Screen.	0 0
23.3	<i>On the NIBP simulator</i> , set the pressure parameter to 200 mmHg.	The NIBP simulator displays a pressure reading of 200 mmHg.	0 0
23.4	<i>On the M Series unit</i> , press the <b>Close Valves</b> softkey.	The Valves status changes from <b>OPEN</b> to <b>CLOSED</b> .	0 0

	Do this...	Verify that...	Pass/ Fail
23.5	<p>Press the <b>START TEST</b> softkey within 30 seconds after closing the valve.</p> <p><b>Note:</b> This test takes approximately 3 minutes.</p>	<p>After approximately 1 minutes, a number appears in the upper middle area of the NIBP simulator display.</p> <p>If the simulator:</p> <ul style="list-style-type: none"> <li>• Displays a Volume Leak reading &lt;5,<sup>5</sup> then the M Series unit has passed the test.</li> <li>• Displays a Volume Leak reading &gt;5,<sup>6</sup> then the M Series unit has failed the test.</li> <li>• Displays no Volume Leak reading, but maintains a stable pressure reading at or above 200 mmHg, then the M Series unit has passed the test; there is no volume leak.</li> </ul> <p>In addition, the M Series unit displays the simulator's pressure reading in the "Cuff Pressure" field.</p> <p>After approximately 3 minutes, the valves open on the M Series unit.</p>	<p>○ ○</p>
23.6	<p>On the NIBP Analyzer, press the <b>STOP TEST</b> softkey.</p>	<p>The NIBP simulator terminates the Volume Leak Test.</p>	<p>○ ○</p>
23.7	<p>On the M Series unit, press the <b>EXIT</b> softkey twice.</p>	<p>The M Series unit returns to the main NIBP Service Mode screen, then to normal Monitor mode operation.</p>	<p>○ ○</p>
<p>Record your results on the Maintenance Tests Checklist.</p>			

<sup>5</sup>If you are using the Fluke® BiomedicalCuffLink Simulator, the volume leak reading should be <10.

<sup>6</sup>If you are using the Fluke® BiomedicalCuffLink Simulator, the volume leak reading for a failure should be >10.

## 24.0 Bluetooth Test

Tools Needed: USB Bluetooth dongle with driver CD, Bluetooth application software BlueSoleil.

Test Setup: Verify that the M Series unit is configured as follows: Bluetooth INSTALLED, and Baud Rate set to 115200 or 38400.

	Do this...	Verify that...	Pass/ Fail
24.1	Power up the M Series unit while pressing in and holding the left-most softkey.	The M Series unit displays the System Utilities Menu.	o o
24.2	Insert a PCMCIA card containing data into the card reader.		
24.3	Press the <b>UPLOAD CARD</b> softkey, then press the softkey <b>UPLOAD CARD</b> .	Verify that the Bluetooth LED is on and displays a green light.	o o
24.4	Press the <b>SEND</b> softkey.		
24.5	Double click on the HyperTerminal shortcut (BT115200) from the desktop.  <b>Note:</b> If setting up Hyperterminal for the first time set COM PORT programmed for Bluetooth adapter, 115200bps, 8, N, 1, N		
24.6	Start Bluetooth Utility software (Blue Soleil) on the computer.		
24.7	In the IVT Corporation Blue Soleil window under "View", select refresh devices.		

	Do this...	Verify that...	Pass/ Fail
24.8	Find the unit under test Serial Number and double click on it.		
24.9	Click on the Serial Port Icon that is highlighted at the bottom of the window	Verify the following: At the bottom of the window it states Connected, The file transfer window pops up and the file upload progress bar goes to 100%. The unit displays the file upload message, <i>Card Uploaded.</i>	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>
Record your results on the Maintenance Tests Checklist.			

# Chapter 2

## Troubleshooting

- 

### Overview

This chapter describes the most common technical problems that biomedical technicians experience when checking the M Series during routine maintenance or when there is a malfunction of the unit. It also contains a list of error messages that users may see if the unit is not operating properly.

This chapter contains the following:

- Troubleshooting tables for ECG Leads Off Messages and Monitor Displays
- Zoll M Series Error Messages

If the problems you encounter are not listed below, call ZOLL Medical Corporation's Technical Service Department for further assistance. (See page vii for contact information.)

## Troubleshooting

The following tables show the most common troubleshooting issues and their solutions.

First, attempt to solve the problem with “Recommended User Action.” If these steps do not solve the problem, follow the steps listed in the “Recommended Technical Action” column.

Reported Problem	Recommended User Action	Recommended Technical Action
<p>ECG LEAD OFF message displays. (3, 5, 12 lead cable)</p>	<ul style="list-style-type: none"> <li>• Check preparation of ECG electrode site by cleaning the site, lightly abrading the patient’s skin and/or clipping the patient’s hair at the electrode site.</li> <li>• If electrode gels are dry, replace electrodes with new ones from a freshly opened package.</li> <li>• Verify that all leads are attached.</li> <li>• Set monitor to another lead.</li> <li>• Verify that the electrodes have not exceeded their expiration date.</li> </ul>	<ul style="list-style-type: none"> <li>• Try to reproduce the problem using a simulator.</li> <li>• Inspect the ECG cables looking for corrosion or broken connector pins.</li> <li>• Check the cable for intermittent connections by flexing the cable at the yoke and snap connectors.</li> <li>• Check the cable connection to the defibrillator.</li> <li>• Inspect the ECG input connector and its pins. Replace it, if necessary.</li> <li>• Inspect the ECG cable connection to the system board.</li> <li>• Inspect the system board ECG shielding.</li> <li>• Remove and replace the system board.</li> </ul>



Reported Problem	Recommended User Action	Recommended Technical Action
<p>V LEADS OFF message displays.</p>	<ul style="list-style-type: none"> <li>• If the user is not using V leads, attach V lead connector terminator plug to the cable's V lead connector.</li> <li>• If a V1 lead wire metal snap comes in contact with the patient's skin, then the system will show all V leads as OFF.</li> <li>• Remove V1 leads and others away from the patient. Turn off the unit and wait ten seconds before turning it back on.</li> </ul>	
<p>CHECK PADS/POOR PAD CONTACT message displays.</p>	<ul style="list-style-type: none"> <li>• Remove and reinsert PADS connector into the universal cable.</li> <li>• Check for damaged defibrillator pads, wires and or connector.</li> <li>• Check for dried out or expired defibrillator pads.</li> <li>• Clip (not shave) the patient's hair and wipe pad contact area dry.</li> <li>• Connect the cable to the test plug. The DEFIB PAD SHORT message displays to indicate that the cable is functioning properly.</li> <li>• If the DEFIB PAD SHORT message displays, then check the connections of the pads to the patient and to the defibrillator cable.</li> <li>• If the DEFIB PAD SHORT message does not display, remove the defibrillator from service.</li> </ul>	<ul style="list-style-type: none"> <li>• Connect universal cable to the shorting plug. The DEFIB PAD SHORT message should display, when you SELECT PADS. If the message does not display, then:</li> <li>• Try another universal cable.</li> <li>• Check the cable from the universal cable connector to the High Voltage Module.</li> <li>• Check the cable from the High Voltage Module to the system board.</li> <li>• Remove and replace the High Voltage Module.</li> <li>• Remove and replace the system board.</li> <li>• Call ZOLL Technical Support for assistance.</li> </ul>

Reported Problem	Recommended User Action	Recommended Technical Action
Flash or arcing under defibrillator pad.	<ul style="list-style-type: none"> <li>• Avoiding using alcohol and betadine in and around the treatment area because these skin preparations may lead to increased conductivity and/or bonding between the electrode's adhesive and skin.</li> <li>• Check for gel droop. If the gel has leaked out of the gel treatment area, replace the electrode.</li> <li>• Ensure pads are coupling to the patient's skin and connected to the universal cable.</li> <li>• Check for dried out gel on the defibrillator pad.</li> <li>• Clip patient's excessive hair. Do not shave hair.</li> <li>• Check expiration date. Replace pad if date has expired.</li> <li>• Do not conduct chest compression through the pads because the pads could be damaged leading to the possibility of arcing and skin burns.</li> <li>• Apply the back electrode first. If the front electrode is already in place when the patient is being maneuvered for placement on the back, the front may become partially lifted, possibly causing arching and skin burns.</li> </ul>	<ul style="list-style-type: none"> <li>• Ensure that wet gel pads are stored flat.</li> </ul>
Displayed HR not accurate. No artifact present.	<p>Verify heart rate flashes with each QRS on display.</p> <ul style="list-style-type: none"> <li>• Change lead selection.</li> <li>• Change ECG size.</li> <li>• Reposition ECG electrodes.</li> </ul>	

Reported Problem	Recommended User Action	Recommended Technical Action
<p>Displayed HR not accurate; artifact present.</p>	<ul style="list-style-type: none"> <li>• Reduce or eliminate ECG artifact due to electrode or patient cable movement. Route cables so that they don't pull on electrodes or swing excessively.</li> <li>• Ensure patient is motionless.</li> <li>• Check for possible excessive radio frequency interference.</li> <li>• Verify a good connection of electrodes to the patient.</li> <li>• Prepare the patient's skin prior to the electrode attachment.</li> <li>• Move patient cables away from other electrical equipment, especially any RFI source.</li> <li>• Ensure ECG cable fits snugly in unit.</li> <li>• Change ECG cable.</li> <li>• Replace/reposition ECG electrodes.</li> </ul>	<ul style="list-style-type: none"> <li>• Check for contamination on snaps. Ensure springs are intact.</li> <li>• Check for intermittent ECG patient cable or connector wiring.</li> <li>• Replace ECG input connector.</li> <li>• Replace ECG connector to the system board cable.</li> <li>• Replace system board.</li> </ul>
<p>Wandering baseline.</p>	<p>See "Displayed HR not accurate." above. Note that in 90% of electrode issues, size and lead changes don't help.</p>	<p>Same as above example.</p>
<p>Electronic interference.</p>	<p>Check for possible excessive radio frequency interference. Move patient cables away from other electrical equipment.</p>	<ul style="list-style-type: none"> <li>• Turn off sources of excessive RFI.</li> <li>• Move M Series unit away from RFI source.</li> </ul>

## Zoll M Series Error Messages

The following is a list of Zoll M Series error messages that may appear on your display. The “User Advisory” column informs you about an action in progress or provides feedback on a user correctable situation that typically does not require further technical support. The “Technical Action” column describes what you as a technician can do to correct the situation. Note that these messages will sometimes overlap part of the waveform display.

First, attempt to clear the message by turning the Selector Switch to OFF for ten seconds, then back to the desired operating mode. If the fault persists, call ZOLL Technical Service.

Error Message	Explanation	User Advisory	Technical Action
200J MAX BIPHASIC	User attempted to set defibrillation energy >200J on Biphasic Unit. No higher energy is available.	✓	
50J MAX	Energy $\leq$ 50J for internal paddles. No higher energy is available.	✓	
ADJUST ECG	Unit is in sync mode and heart rate is < 20 BPM. Or, QRS size set too small for proper synchronization.	✓	
ANALYSIS HALTED	<ul style="list-style-type: none"> <li>• ECG analysis halted due to user interaction such as:</li> <li>• Lead/size change</li> <li>• Analyze button was pressed again</li> <li>• Impedance fault</li> <li>• Charging error detected in auto defib mode</li> </ul>	✓	
ANALYSIS RESTARTED	This is a user prompt issued simultaneously with ECG TOO LARGE or ECG TOO SMALL. Device detected ECG signals out of range, automatically adjusted ECG size and is now restarting its shockable rhythm analysis sequence.	✓	

Error Message	Explanation	User Advisory	Technical Action
AUDIO FAULT 136	Audio DSP hardware error.		Replace audio board. Replace system board.Turn unit off and back on again.
AUDIO NOT RECORDING	Audio is not recording.		Install PCMCIA card. Replace system board.
AUDIO QUEUE FULL	Indicates that the audio output queue is full. Additional voice prompts can't be queued at this time.		None.
BATT HIGH CURRENT	Battery is charged and battery current is >.1 A or: Battery is not charged and battery current is > 1.6 A.		Unplug from A/C. Remove the battery for 20 seconds. Reconnect all above. If the problem persists, replace battery and or charger.
BATT HIGH VOLTAGE	Battery voltage > 15.5 v.		Replace battery and or charger
BATT LOW CURRENT	Battery is not charged and battery current is <.35 A.		Replace battery and or charger.
BATT LOW VOLTAGE	Battery voltage < 9.5 v.		Replace battery and or charger.
BATT OVERCHARGE	Charger on for > 4 hours.		Replace battery and or charger.

Error Message	Explanation	User Advisory	Technical Action
BRIDGE SHORT	Current higher than expected was detected during the Biphasic bridge test or immediately following a discharge.		<p>Ensure pads/paddles are used properly.</p> <p>Attempt to clear the message by turning the Selector switch to off then back to the desired operating modes.</p> <p>Replace bridge or high voltage module.</p>
BRIDGE TEST FAILED	Biphasic module not operating properly while charging.		<p>Charge again.</p> <p>Attempt to clear the message by turning the Selector switch to OFF, then back to the desired operating mode.</p> <p>Replace bridge or high voltage module.</p>
CABLE FAULT	(Auto defib mode only.) Incorrect A/D reading for paddle ID (similar to PADDLE FAULT).	✓	Replace paddle set, universal cable and/or system board.
CANNOT CHARGE	Cannot charge when charge button pressed.		Replace high voltage module or capacitor.
REPLACE CARD	Write errors during manual or semi-automated modes.	✓	May have configuration card installed or write protection on.

Error Message	Explanation	User Advisory	Technical Action
CARD FULL	Memory Card Full.	✓	
CHECK CO <sub>2</sub> SENSOR	EtCO <sub>2</sub> Sensor is unplugged or defective.	✓	Check that sensor cable is plugged in and seated properly. Check that sensor is not exposed to excessive heat. If problem persists, replace the sensor.
CHECK CO <sub>2</sub> ADAPTER	Airway adapter is removed, occluded or adapter zeroing needs to be performed or was performed incorrectly.	✓	Replace/Clean airway adapter. Zeroing performed automatically.
CHECK MEMORY CARD	No card detected during manual or semi-automated modes.	✓	
CHECK PADS	Message displayed in conjunction with either POOR PAD CONTACT or DEFIB PAD SHORT.		Ensure pads are coupled to patient. Check /replace pads and universal cable. Replace system board.
CHECK PATIENT	Background ECG analysis detects shockable rhythm.	✓	
CHECK PULSE	Alternate message for NO SHOCK ADVISED message. Message also shown after delivering third shock when auto analyze 3 times option is enabled.	✓	
CHECK RECORDER	Produced when paper tray is empty, paper jams or recorder door is opened.		Replace paper sensor board, system interconnect board, and/or system board.

Error Message	Explanation	User Advisory	Technical Action
CHECK SPO <sub>2</sub> SITE	Low or no perfusion in monitored finger or toe.	✓	
CHECK SPO <sub>2</sub> SENSOR	Reposition SpO <sub>2</sub> sensor on patient.		
CLOCK FAULT 11	Real time clock oscillator failure.		Replace system board.
CLOCK FAULT 12	Real time clock back-up power supply failure. Found oscillator stopped at power-up, but oscillator now running when the system is running. (Oscillator only runs when main power is applied).		Replace system board.
CLOCK FAULT 13	One of the set time units (seconds, minutes, year, etc.) is out of range.		Replace system board.
CO <sub>2</sub> COMM ERROR	No or invalid communication from the EtCO <sub>2</sub> module.		Replace EtCO <sub>2</sub> module and or system board.
CO <sub>2</sub> SENSOR WARM UP	EtCO <sub>2</sub> Sensor warming up.	✓	Wait for sensor to warm up. This process takes up to approximately one minute.
CONFIRM MANUAL MODE	Displayed when manual mode is entered. Alerts user to confirm that manual mode is desired.	✓	
DEFIB DISABLED	User prompt issued simultaneously with other faults if defib is disabled.		Possible configuration problem. Replace high voltage module. Call ZOLL Technical Support.
DEFIB FAULT 71	More than 50 internal dumps occurred in less than 20 minutes.		Turn the unit to OFF and back on. If fault persists, replace high voltage module.



Error Message	Explanation	User Advisory	Technical Action
DEFIB FAULT 72	General defib error.		Turn the unit to OFF and back on. If fault persists, replace high voltage module.
DEFIB FAULT 76	Capacitor voltage too high for selected energy.		Replace high voltage module or capacitor.
DEFIB FAULT 77	Capacitor voltage > than absolute rated max.		Replace high voltage module or capacitor.
DEFIB FAULT 78	Unable to charge defib cap.		Replace high voltage module or capacitor.
DEFIB FAULT 79	Defibrillator charging too slowly.		Replace high voltage module or capacitor.
DEFIB FAULT 80	4 defibrillator faults detected within 20 second period.		Replace high voltage module or capacitor.
DEFIB FAULT 81	Discharge switch in undefined state.		Replace high voltage module or capacitor.
DEFIB FAULT 84	“Upper” discharge transistor shorted (measured via applicable A/D channel).		Replace high voltage module.
DEFIB FAULT 85	“Lower” discharge transistor shorted (measured via applicable A/D channel).		Replace high voltage module.
DEFIB FAULT 86	One discharge switch closed during power up test.		Replace paddles, control board or system board.
DEFIB FAULT 87	Both discharge switches closed during power up test.		Replace paddles, control board or system board.

Error Message	Explanation	User Advisory	Technical Action
DEFIB FAULT 94	Processor fault causing safety monitor port to be non-functional.		Replace system board, high voltage module or capacitor.
DEFIB FAULT 95	Safe or shutdown line is not functional.		Replace high voltage module.
DEFIB FAULT 96	XPATREL or XPAT_ENABLE is faulted or one of the discharge transistors has shorted.		Replace high voltage module.
DEFIB FAULT 108	VMON voltage is less than the target energy during charging.		Replace high voltage module or capacitor.
DEFIB FAULT 109	Defib capacitor voltage is greater than selected energy when defibrillator is charging or ready.		Replace high voltage module or capacitor.
DEFIB FAULT 111	Defib capacitor voltage has exceeded the absolute maximum acceptable voltage.		Replace high voltage module, capacitor, and or system board.
DEFIB NOT CHARGED	Discharge button is pressed but the unit is not charged.	✓	
DEFIB PAD SHORT	Measured impedance between high voltage leads of MFC.	✓	Ensure pads are coupled to patient. Check /replace pads or universal cable. Replace system board.
DISABLE SYNC	Sync mode active when analyze pressed in defib.	✓	
DISCHARGE FAULT	Defib capacitor voltage is not decreasing.		Replace high voltage module, capacitor, and/or system board.

Error Message	Explanation	User Advisory	Technical Action
ECG FAULT 4	Communication fault between ECG processor and main processor.		Turn off unit and then turn on to reset. If fault persists, replace system board.
ECG FAULT 5	ECU RAM test failure, or ROM checksum test failure.		Turn off unit and then turn on to reset. If fault persists, replace system board.
ECG LEAD OFF	One or more ECG leads are not properly connected when leads are selected as input.	✓	Check cable and patient connection. Change electrodes. Prepare patient's skin.
ECG TOO LARGE	ECG signal too large for accurate shockable rhythm analysis.	✓	Reduce ECG size.
ECG TOO SMALL	ECG signal too small for accurate shockable rhythm analysis.	✓	Increase ECG size.
ECG V1 LEAD OFF	Chest lead V1 is not properly attached to patient.	✓	Reattach V lead. Check cable.
ECG V2 LEAD OFF	Chest lead V2 is not properly attached to patient.	✓	Reattach V lead. Check cable.
ECG V3 LEAD OFF	Chest lead V3 is not properly attached to patient.	✓	Reattach V lead. Check cable.
ECG V4 LEAD OFF	Chest lead V4 is not properly attached to patient.	✓	Reattach V lead. Check cable.
ECG V5 LEAD OFF	Chest lead V5 is not properly attached to patient.	✓	Reattach V lead. Check cable.

Error Message	Explanation	User Advisory	Technical Action
ECG V6 LEAD OFF	Chest lead V6 is not properly attached to patient.	✓	Reattach V lead. Check cable.
ENTER ACCESS CODE	Manual mode access code needed.	✓	Enter access code to enter manual mode with AED.
ERASING REPORT	Summary report being erased.	✓	
ECU CRC FAULT	Invalid ECG samples detected over a one second period.		Turn off unit and then turn on to reset. If fault persists, replace system board.
EtCO <sub>2</sub> COM ERROR	No or invalid communication from EtCO <sub>2</sub> module.		Return unit for service to ZOLL Technical Service Department.
FAX DIALING	Preparation for sending fax.	✓	
FAX DONE	Transmission complete.	✓	
FAX PREPARING	Preparing fax for transmission.	✓	
FAX SENDING	Transmitting fax.	✓	
INSERT CARD Check memory card	No card installed in unit during manual or semi-automated modes.	✓	
LOW BATTERY	Low battery.	✓	Replace battery or plug into AC power. Replace charger.
NO QRS DETECT	Unit is in sync mode and heart rate is < 20 BPM or QRS amplitude is too low for proper synchronization.	✓	Increase ECG size and/or change lead.

Error Message	Explanation	User Advisory	Technical Action
NO SHOCK ADV	No shock advised. Advisory message when analysis finds non-shockable rhythm.	✓	
NOISY ECG	Number of noisy analysis intervals exceeds threshold.	✓	Stop all patient movement. Check connections. Press Analyze button again.
OPEN AIR DISCHARGE	Cap voltage too high after discharge attempt, e.g., full energy discharge did not occur.		Replace paddles, and, or high voltage module and system board.
PACER DISABLED	User prompt issued simultaneously with other pace faults if pacing is disabled.		Replace high voltage module or system board.
PACER FAULT 115	Flyback pulse width control circuit is not under proper control of the processor and gate array.		Replace high voltage module, capacitor, or system board.
PACER FAULT 116	Failure to detect XPACE_ON.		Replace high voltage module, capacitor, and/or system board.
PACER FAULT 117	Pace relay is stuck closed.		Replace high voltage module, capacitor, and/or system board.
PACER FAULT 121	During pace, the pace pulse width <30ms or >50ms.		Replace high voltage module, or system board.
PACER FAULT 122	Pace current is more than 15mA above and below selected value.		Replace high voltage module, or system board.

Error Message	Explanation	User Advisory	Technical Action
PACER FAULT 123	Measured pace rate is too fast compared to selected rate.		Replace high voltage module or system board.
PACER FAULT 126	Issued in conjunction with message 122. Pace current is more than 15mA and below selected value.		Replace high voltage module.
PADDLE FAULT	Cannot detect type of accessory attached to the universal cable.		Replace paddles, internal paddles, system board, high voltage module and/or universal cable.
PERFORM CPR	Advisory message in AED auto defib mode.	✓	
PLACE ON ZERO CELL	EtCO <sub>2</sub> sensor cable plugged into unit for the first time. Zeroing error or probe drift error detected.	✓	Zero sensor. Replace sensor. Return to ZOLL for service.
POOR LEAD CONTACT	One or more ECG leads are poorly connected or not connected to patient. (User configurable.)		Check electrode attachment to patient, cable connector to electrode, cable to unit connector. Broken unit.
POOR PAD CONTACT	Electrode impedance exceeds threshold.		Ensure pads are coupled to patient. Check /replace pads or universal cable. Check impedance circuit calibration. Replace system board.
PRESS ANALYZE	Alternate message for check patient prompt.	✓	

Error Message	Explanation	User Advisory	Technical Action
PRESS CHARGE	Advisory message in conjunction with shock advised.	✓	
PRESS SHOCK	Prompt issued in AED auto defib mode when defib is charged (ready).	✓	
RECORDER FAULT 142	Strip chart system error.		Check paper tray and paper path. Replace the print head, system interconnect board and or the system board.
RECORDER FAULT 143	Strip chart failed power-up echo test. Communications error.		Check paper tray and paper path. Replace the system interconnect board and/or the system board. Turn unit off and back on again.
RECORDER FAULT 147	Strip chart printhead over temperature.		Check paper tray and paper path. Replace the print head, system interconnect board and/or the system board.
RELEASE BUTTONS	Simultaneous external paddle button presses detected before unit reached full defib charge (ready state).	✓	Release buttons.
RELEASE SHOCK	Discharge switch(es) closed when pressing charge button. Discharge button pressed before defib reached ready state.	✓	Release shock button. Check paddles. Replace controls board.
REPLACE BATTERY	Battery voltage is less than absolute minimum. Shutdown imminent.	✓	Replace with charged battery.

Error Message	Explanation	User Advisory	Technical Action
REPLACE EtCO <sub>2</sub> SENSOR	EtCO <sub>2</sub> SENSOR WARM UP message displays for more than five minutes. Sensor defective.	✓	Replace sensor cable.
REPORT FULL	Summary report memory full.	✓	Erase summary report.
REPORT HALTED	Summary report stops printing unexpectedly.		Turn unit off and then back on again. Print Summary again. If fault persists, replace system board.
RESERVED 1	The watchdog timer is not functional in the unit.		Turn off unit and then turn on to reset. If fault persists, replace system board.
RETRY ANALYSIS	Advisory message in conjunction with noisy ECG. Analysis halted.	✓	
SELECT 30J FOR TEST	Attempt to run a self test at an energy other than 30J.	✓	
SELECT DEFIB MODE	Analyze button pressed in pace or monitor mode.	✓	
SELECT LIMB LEADS	Paddles or augmented ECG leads selected when continuous analysis active or started.	✓	Select limb leads I, II, III or MFE
SELECT PADS	Lead I, II, or III selected when analyze pressed.	✓	
SET CLOCK	Real time clock failure: invalid date or time.		Set date and time and/or replace system board.
SET PACE mA	Multiple copy errors are the product of intended software or memory errors. If error reoccurs other than on entering pace the first time or after more than 10 minutes in other mode, the unit could be broken.	✓	Set pace current. If broken, replace system board.



Error Message	Explanation	User Advisory	Technical Action
SET PACE RATE	Multiple copy errors are the product of intended software or memory errors. Multiple copies of pace rate don't match. If error persists, unit could be broken	✓	Set pace rate. If broken, replace system board.
SHOCK ADVISED	Advisory message when analysis finds a shockable rhythm.	✓	
SpO <sub>2</sub> AMBIENT LIGHT	Ambient light is too bright.		Shield sensor from ambient light. Replace SpO <sub>2</sub> sensor. Replace SpO <sub>2</sub> module
SpO <sub>2</sub> COMM ERR	No transmissions from SpO <sub>2</sub> unit received. Communication error or no communication from SpO <sub>2</sub> module.		Replace SpO <sub>2</sub> module and/or system board.
SpO <sub>2</sub> PULSE SEARCH	Pulse search in progress.	✓	
SpO <sub>2</sub> SENSOR FAULT 1	Defective sensor.		Replace SpO <sub>2</sub> sensor.
SpO <sub>2</sub> SENSOR FAULT 5	Unrecognized sensor.		Replace SpO <sub>2</sub> sensor.
STAND CLEAR	(Auto defib mode only.) Single analysis mode just turned on and defib idle. Patient rhythm is being analyzed.	✓	
SYSTEM FAULT 2	MCU ROM checksum test failure or MCU RAM test failure.		Turn off unit and then turn on to reset. If fault persists, replace system board.

Error Message	Explanation	User Advisory	Technical Action
SYSTEM FAULT 5	ECU RAM test failure or ROM checksum test failure.		Turn off unit and then turn on to reset. If fault persists, replace system board.
SYSTEM FAULT 6	No communications received from ECU for 4 seconds.		Turn off unit and then turn on to reset. If fault persists, replace system board.
SYSTEM FAULT 7	The A/D converter is not performing conversions in a timely manner.		Replace system board.
SYSTEM FAULT 36	P <sub>1</sub> MON is less than 412 counts or greater than 612 A/D counts. Pace/defib is disabled as long as condition exists.		Replace system board.
SYSTEM FAULT 37	Disable pace/defib and MFE monitoring.		Replace system board.
SYSTEM FAULT 38	Failure to shutdown after “shutdown order” is written to the RTC.		Replace system board.
TEST FAILED	MCU performed izeak test (defib peak current) and unit failed during 30J self test.	✓	Replace universal cable, paddles or high voltage module, capacitor, or system board.
TEST OK	MCU performed izeak test (defib peak current) and unit passed 30J self test.	✓	
USE PADDLE DISCHG	Front Panel discharge button is pressed when either external paddles or internal spoons with discharge buttons are connected.	✓	

Error Message	Explanation	User Advisory	Technical Action
USE PADS	(AUTO DEFIB MODE ONLY.) Attempt to defib with paddles in auto defib (AED) mode. Defib only allowed using PADS in AED modes.	✓	
USE PADS TO PACE	External paddles detected in pace mode.	✓	
USE ROOM AIR ADAPTER	Adapter zeroing started with EtCO <sub>2</sub> in the adapter or the adapter is on the REF or "0" cell.	✓	Place CO <sub>2</sub> sensor on adapter in room air.
USER SETUP REQ	Both copies of stored cal/config data are bad or have never been programmed.	✓	Perform configuration setup.
VF ALARMS OFF	VF alarms disabled in pace mode or when paddles are selected as leads.	✓	
VX LEADS OFF	V lead not properly attached to patient. "X" denotes lead number.	✓	Reattach V lead.
ZERO CO <sub>2</sub> SENSOR	New EtCO <sub>2</sub> sensor needs to be zero calibrated.	✓	Zero EtCO <sub>2</sub> sensor.
ZERO CO <sub>2</sub> ADAPTER	New EtCO <sub>2</sub> airway adapter needs to be zero calibrated.	✓	Zero EtCO <sub>2</sub> adapter.



# Chapter 3

## Disassembly Procedures

### Overview

This chapter provides instructions on how to disassemble and reassemble the M Series unit, and includes the following sections:

- Required Equipment
- Parts That May Need Replacing After Disassembly
- Safety Precautions
- Overview of Modules
- 1. Removing the ZIF Keeper
- 2. Removing the Front Panel
- 2A. Removing the Display
- 2B. Removing the Control Board
- 3. Removing the Upper Housing Assembly
- 4. Removing the System Board Assembly
- 5. Removing the Battery Interconnect Board Assembly
- 6. Removing the High Voltage/Charger Assembly
- 7. Removing the High Voltage Module Assembly
- 8. Removing the High Voltage Capacitor Assembly
- 9. Removing the System Interconnect Board
- 10. Removing the Printer/Recorder Motor
- 11. Removing the Lower Housing Assembly
- 12. Removing the Print Head Assembly

- 13. Removing the PCMCIA Card Slot Assembly
- 14. Removing the Paddle Release Latch

### Required Equipment

- No. 1 Phillips screwdriver.
- No. 2 Phillips screwdriver.
- Exacto-knife.
- Orange wooden stick. (Available from H.A. Stiles: 1-800-447-8537)
- 90° dental pick.
- Needle nose pliers.
- Kapton tape.
- 3M copper adhesive tape, or equivalent.
- 1/2" nut driver.
- Large diagonal cutters.
- Strong glue, such as Loctite 420 or equivalent.
- Loctite needle tip dispenser.

### Parts That May Need Replacing After Disassembly

If you are removing the Control Board from the Front Panel, you need to have:

- Main Selector knob replacement (ZOLL Part Number 9310-0521)
- Pacer/Output/Rate knob replacement (ZOLL Part Number 9310-0520)

If you are removing the Battery Interconnect Board, you may need to replace it with a new one, using ZOLL Part Number 9301-0302, if connectors have been UV welded as in older M Series models.

## Safety Precautions

---

**WARNING!**

**SHOCK HAZARD!**

---



---

**CAUTION** TAKE THE NECESSARY PRECAUTIONS TO GUARD AGAINST SHOCK OR INJURY BEFORE YOU CONDUCT DEFIBRILLATOR TESTS OR REPAIRS.

---


- Only properly trained technicians should service the unit.
- The unit can contain deadly voltages even if the unit is turned off.
- Make sure to discharge the unit before working with it.
- Make sure you take the necessary precautions when working with static sensitive units. For example, you must wear a conductive wrist strap (which touches your skin) connected to a grounding mat and to the earth ground. You must remove the wrist strap when you discharge high voltage or when you are working on energized equipment.

## Overview of Modules

The M Series unit contains 14 modules, as shown below.



Isolated Power Supply with EtCO<sub>2</sub>




SpO<sub>2</sub> module with bracket for EtCO<sub>2</sub>



System Board Assembly



System Interconnect Board



Battery Interconnect Board Assembly with 3 Battery Pin Gaskets



AC Charger Assembly





Control Board (from Front Panel)



High Voltage Module Assembly




High Voltage Capacitor Assembly



PCMCIA Card Slot



Recorder Motor



SpO<sub>2</sub> Module (without bracket)



Isolated Power Supply for SpO<sub>2</sub> Module



Biphasic Capacitor and Bridge Assembly



# 1. Removing the ZIF Keeper

## Tools Required


- Orange stick

Note: It is important to know this procedure before you start disassembling the unit. Removing the ZIF (Zero Insertion Force) Keeper incorrectly can damage the unit's system board.


## To reinstall the ZIF Keeper:

1. Place the ZIF Keeper over the laminate cable and insert the flex cable into the connector. Latch the connector.
2. Lower the left end of the ZIF Keeper over the connector end until it touches the printed wire assembly (PWBA). The other end of the connector should be angled.
3. Press the end of the ZIF Keeper down over the end of the connector. Be careful that the ZIF Keeper snaps over the end of the connector.


**Step 1:** The connector must be facing you as shown in the diagram.



**Step 2:** Angle and lift up the right end of the ZIF Keeper from the connector and the board. Slide the Keeper approximately 1 mm to the left, then gently lift the left side to clear the connector.



**Step 3:** Gently pull the ZIF (Zero Insertion Force) Keeper over the laminate cable and rotate the latch upwards.



## 2. Removing the Front Panel

### Tools Required

- No. 2 Phillips screwdriver
- Orange stick

### To reinstall the Front Panel:

1. Reinstall the laminate cable first with the black band facing up and towards the system circuit board.
2. Reconnect the multi-wire cable from the display.
3. Reverse steps 1, 2, and 3 above to reinstall the front panel.

**Step 1:** Remove the battery from the battery well and place it in front of the unit.




**Step 2:** Rotate unit on to its back side. Remove the two Phillips head screws located on the left and right sides on the bottom of the unit.



**Step 3:** Place your thumbs in the main selector switch cup and push up on the front panel to release the panel from the unit. After the front panel is removed, use the battery as a support for the panel. **Do not use the main selector switch as point of leverage.**



**Step 4:** Disconnect the multi-wire cable from the system board by gently pulling the beige connector by its sides towards the front of the unit. Remove the ZIF keeper from the laminated ribbon cable and then disconnect it. Lift the right side first with the unit facing you.



# 2A. Removing the Display

*Tools Required*

- No. 2 Phillips Screwdriver
- Exacto-Knife
- 3M Copper Adhesive Tape

*To reinstall the Display*

Reverse steps 1 through 3.


**Step 1:** Remove the grounding copper tape from the outer display shield.



**Step 2:** Remove the Video Display Assembly by rotating the display upwards from the lower portion of the display panel assembly.



**Step 3:** Remove the Video Display Assembly.



# 2B. Removing the Control Board

### Tools Required

- No. 2 Phillips Screwdriver
- Orange (Wooden) Stick
- 1/2" Nut Driver
- Large Diagonal Cutters
- Strong Glue, such as Loctite 420 or equivalent


Note: If you are removing the Control Board from the Front Panel you may need to replace the following parts:

- Main Selector knob replacement (ZOLL Part Number 9310-0521)
- Pacer Output/Rate knob replacement (ZOLL Part Number 9310-0520).

**Step 1:** (Caution: The knob will be damaged during this step.) Gently insert the cutters at the edge of the main selector knob and pry outward until the knob is removed. Then carefully remove the 1/2" nuts without damaging the Selector Switch.




**Step 2:** Gently insert the cutters at the edge of the pacer knobs and pry outward until the knobs are removed. Then carefully remove the 1/2" nuts without damaging the Pacer Switches and green pace cups.



**Step 3:** Rock the foam packaging back and forth from the control panel board to remove the foam.



**Step 4:** Remove the panel overlay circuit cable from the control board by lifting the sides of the lock lever located under the front of the control board. (See 1. Removing the ZIF Keeper.)




## 2B. Removing the Control Board (Continued)

**Step 5:** Disconnect the speaker microphone cable. Remove the Control Panel Board from the Front Panel Assembly. **Important:** Note the position of the dip switches because they must be in the same position for reinstallation.



**Step 6:** Remove the 50 pin laminate cable from the control board. See *1.0 Removing the ZIF Keeper*.



### *To reinstall the Control Board*

1. Set the dip switches and attach the SHOCK button LED actuator.
2. Place a small amount of strong glue on the end of the mode selector switch. Slide the replacement knob on. After the glue dries, rotate the selector knob to ensure that it is properly glued in place.
3. Replace the Pacer Output/Rate knobs, if applicable.

# 3. Removing the Upper Housing Assembly

## Tool Required

- No. 1 and 2 Phillips screwdriver

## Set-up

Remove Front Panel (See Step 2.)

1. Remove the two screws securing the universal cable.
2. Do not lose the O-ring when removing the universal cable.
3. Remove the three Phillips screws in the front of the Upper Housing.

To reinstall the Upper Housing Assembly, reverse the above steps.

Ensure the battery pin gaskets are properly set.


**Step 1:** Remove two screws from the back side of the Upper Housing Assembly and three screws from the front.



**Step 2:** Secure lower housing to the table by pressing downward in the paddle well. Using the carry handle, lift the Upper Housing upward.



**Step 3:** Make sure that the rubber gaskets are covering the battery contact pins. If the gaskets are still seated in the housing, remove and place them onto the contact pins. Before installing the Upper Housing ensure that the battery pin gaskets are set properly.





# 4. Removing the System Board Assembly

2. Remove the Front Panel Assembly.
3. Remove the Upper Housing Assembly.

*To reinstall the System Board Assembly, reverse the steps.*

## ***Tools Required***

- #1 and #2 Phillips Screwdriver
- Grounding Mat
- Grounding Wrist Strap
- Needle Nose Pliers
- Orange Stick


**WARNING!** You can damage the hardware of the unit. You must use ESD grounding before you handle any printed circuit boards on the unit.

## ***Setup***


Before you begin this procedure, make sure you are grounded.

1. Know or review *1.0 Removing the ZIF Keeper* procedure.


**Step 1:** Using needle nose pliers, remove left and right multi-wire cables on the back side of the system board. To avoid damage to the cable, do not pull the wires. Hold system board securely with one hand. **DO NOT** let it fall forward to prevent excess tension on the unit's wires.



**Step 2:** Remove the two wire speaker cable (if applicable).



**Step 3:** Lower the system board with one hand and remove the battery interconnect laminate cable from the center of the system board. (See *1. Removing the ZIF Keeper.*)



**Step 4:** Gently roll the system board forward and rest it on the battery.



**Step 5:** Remove the ZIF Keeper. (See 1. Removing the ZIF Keeper.)



**Step 6:** For Biphasic M Series units only: Remove the 20 pin power cable by lifting the slide locking tab upwards.



**Step 7:** For M Series units with the SpO<sub>2</sub> and/or SpO<sub>2</sub> with EtCO<sub>2</sub> only: Remove the 20 pin power cable by lifting the slide locking tab upwards.



# 5. Removing the Battery Interconnect Board Assembly

## Tools Required

- Orange stick

## Setup

1. Remove the Front Panel Assembly.
2. Remove the Upper Housing Assembly
3. Remove the System Board Battery Cable.


To reinstall the Battery Interconnect Board, reverse the steps.

**NOTE** For Step 4, remember to carefully remove the cable from the Charger Assembly when disassembling and reinstalling it.


**Step 1:** Identify the Battery Interconnect Board.




**Step 2:** Remove the wide laminate cable from the high voltage module connector by lifting the cable vertically.



**Step 3:** Remove the push pin and small insert collar.



**Step 4:** Rotate the unit around so that the rear of the unit faces you. Rotate the Battery Interconnect Board upwards toward the front of the device. Hold the board vertically while removing the 10-pin laminate cable from connectors by lifting the connector lock.



# 6. Removing the High Voltage/ Charger Assembly

DO NOT SHORT THE TERMINAL ENDS OF THE CAPACITOR. Refer to Section 8, Step 2.

## *Tools Required*

- #2 Phillips screwdriver
- Orange stick
- Small needle nose pliers


## *Setup*

1. Remove the Front Panel Assembly.
2. Remove the Upper Housing Assembly.
3. Remove the System Board.
4. Remove the Battery Interconnect Board.


*To reinstall the High Voltage Charger Assembly, reverse the steps.*

**WARNING!** This unit may contain lethal voltages. You **MUST** completely discharge the high voltage capacitor before removing from unit.


**Step 1:** Remove the High Voltage/Charger Assembly from the main housing by lifting the High Voltage/Charger Assembly upwards and rotating it towards the back of the unit.




**Step 2:** Set the Assembly on the table. Pull apart the High Voltage Module from the Charger Assembly. These three components are referred to as the “Module Cluster”.



**Step 3:** Depress the lever on the High Voltage Cable to remove it from the mounting panel. Cut the tie wrap attached to the mounting bracket that secures the cables.




**Step 4:** (SpO<sub>2</sub> units only) To remove the SpO<sub>2</sub> module from the Isolated Power Supply, remove the Mylar® tape. EtCO<sub>2</sub> module and power supply are one assembly at the same location.




**Step 5:** Lift the SpO<sub>2</sub> module straight up from the foam and disconnect the ribbon cable from the Isolated Power Supply.




**Step 6:** (Biphasic units only) Disconnect the biphasic defibrillator capacitor from the rear panel connector.



**Step 7:** (Biphasic units only) Disconnect the Molex connector from the biphasic module.



**Step 8:** (EtCO<sub>2</sub> units only) Remove the EtCO<sub>2</sub> connector by gently lifting the connector away from the back panel. To remove the EtCO<sub>2</sub> and isolated power supply, refer to Step 5.



# 7. Removing the High Voltage Module Assembly

5. Remove the High Voltage/Charger Assembly.  
*To reinstall the High Voltage Module Assembly, reverse the steps.*

## *Tools Required*

- #2 Phillips screwdriver
- Exacto-Knife

**WARNING!** This unit may contain lethal voltages. You must completely discharge the high voltage capacitor by changing the energy selection during the charge cycle of the defibrillator. Power off unit. Wait three minutes before disassembly.


**DO NOT SHORT THE TERMINAL ENDS OF THE CAPACITOR.** Refer to Section 8, Step 2.

## *Setup*

1. Remove the Front Panel Assembly.
2. Remove the Upper Housing Assembly.
3. Remove the System Board.
4. Remove the Battery Interconnect Board.



**Step 1:** Remove the signal cable from the High Voltage module.



**Step 2:** Remove the Kapton tape on the bottom of the foam surrounding the High Voltage Module.



**Step 3:** Separate the High Voltage Module from the foam.



# 8. Removing the High Voltage Capacitor Assembly

## Tools Required

- #2 Phillips screwdriver
- Exacto-Knife
- Mylar<sup>®</sup> tape


## Setup

1. Remove the Front Panel Assembly.
2. Remove the Upper Housing Assembly.
3. Remove the System Board.
4. Remove the Battery Interconnect Board.
5. Remove the High Voltage Charger Assembly.
6. Lift the Capacitor Assembly upwards from the chassis. (The Capacitor Assembly is still connected to the High Voltage Module Assembly.)


*To reinstall the High Voltage Capacitor Assembly, reverse the steps.*

**WARNING!** This unit may contain lethal voltages. You MUST completely discharge the high voltage capacitor before removing from unit. DO NOT SHORT THE TERMINAL ENDS OF THE CAPACITOR.


**Step 1:** Open the High Voltage Capacitor plastic isolator by lifting the Mylar tape.



**Step 2:** Bleed the excess voltage using a resistor with values of approximately 5 kohms, 25 watts for 10-20 seconds. Measure the voltage on the capacitor terminals.



**Step 3:** When the capacitor voltage is at 0 VDC, disconnect the inductor wires and capacitor terminal wires.



# 9. Removing the System Interconn ect Board

3. Install two screws.
4. Reconnect all the cables and reverse the steps.

## *Tools Required*

- #2 Phillips screwdriver

## *Setup*

1. Remove the Front Panel Assembly.
2. Remove the Upper Housing Assembly.
3. Remove the System Board.
4. Remove the Battery Interconnect Board.
5. Remove the High Voltage/Charger module.
6. Remove the clear plastic print head isolater.
7. Remove the motor cable, cable print head and paper sensor cable.

## *To reinstall the System Interconnect Board:*

1. Install the laminate cables and keepers.
2. Ensure that the PCMCIA slots with gaskets are located to the right rear of the board and seated properly.


**Step 1:** Remove the clear plastic head isolater.



**Step 2:** Remove two screws from each side of the System Interconnect Board.



**Step 3:** Disconnect the cables and lift the System Interconnect Board out of the unit.



# 10. Removing To reinstall the Printer/Recorder Motor, reverse the steps. the Printer/ Recorder Motor


## *Tools Required*

- #2 Phillips screwdriver
- Exacto-Knife
- Kapton tape


## *Setup*

1. Remove the paper tray, pull out and press up on the locking tab at the rear of the tray.
2. Remove the Front Panel Assembly.
3. Remove the Upper Housing Assembly.
4. Remove the System Board.
5. Remove the Battery Interconnect Board.
6. Remove the High Voltage Module Assembly.
7. Remove the sensor, printhead and motor cables.
8. Remove the System Interconnect Board.


**Step 1:** Using the needle nose pliers, carefully lift the spacer block upwards.



**Step 2:** Remove the motor bracket mounting screw.



**Step 3:** Lift the motor upward by prying on the motor bracket with the screwdriver.



# 11. Removing the Lower Housing Assembly

## Tools Required

- #2 Phillips screwdriver

## Setup

- Remove all power sources, such as the battery and power cord.

*To reinstall the Lower House Assembly, reverse the steps.*

**Step 1:** Remove the paper tray by pulling it out and pressing the locking tab upwards at the rear of the tray.



**Step 2:** Remove the screw at the bottom of the unit to remove the PCMCIA FAX/Modem card's plastic protector.



**Step 3:** Remove the remaining screws from the bottom of the unit.



**Step 4:** Lift the lower housing assembly straight up and out from the unit





# 12. Removing the Print Head Assembly

## Tools Required


- #1 Phillips screwdriver

## Setup


- Remove the Lower Housing Assembly.

To reinstall the Print Head Assembly, reverse the steps.


**Step 1:** Remove two screws from the paper tray guide and remove the paper tray guide.



**Step 2:** Disconnect the ribbon cable from the print head.




**Step 3:** Remove the pivot pin for re-use.



# 13. Removing the PCMCIA Card Slot Assembly

**Step 1:** Remove the two screws holding the card slot retainer. Lift up card slot assembly.



## *Tools Required*

- #2 Phillips screwdriver


## *Setup*

- Remove the Lower Housing.


*To reinstall the PCMCIA card slot, reverse the step.*

# 14. Removing the Paddle Release Latch


**Step 1:** Using the Exacto Knife, remove the adhesive seal on the unit's chassis.



**Step 2:** Insert the screwdriver into the opening under the adhesive seal.



**Step 3:** Gently push the screwdriver until the Paddle Release Latch is dislodged.



**Step 4:** Pull Paddle Release button away from the unit.



## Tools Required

- #2 Phillips screwdriver
- Exacto-Knife
- Upper Latch Seal (9330-0304)

## Setup

1. Remove the Front Panel Assembly.
2. Remove the Upper Housing Assembly.

*To reinstall the Paddle Release Latch, reverse the steps.*



# Chapter 4

## Replacement Parts

### Overview

This section contains a listing of the replacement parts available for the ZOLL M Series devices.

Replacement parts may be ordered through an authorized ZOLL distributor or directly from ZOLL Medical Corporation. The prices for parts are available from ZOLL Medical Corporation's Technical Service Department.

When ordering parts, please provide the following information:

- ZOLL M Series device model and serial number
- Field Replaceable unit part number
- Description of the replacement part

To order by mail from ZOLL Medical Corporation, address your request to:

ZOLL Medical Corporation

269 Mill Road

Chelmsford, MA. 01824-4105

Attention: Technical Service Department

1-978-421-9655; 1-800-348-9011; Fax: 1-978-421-0010

Note

ZOLL reserves the right to substitute different parts to reflect modifications and improvements in ZOLL M Series circuitry and design.

## Replacement Parts

<b>Description</b>	<b>Part Number</b>
12 Lead ECG Cable	1001-0031-01
3M 12 Conductor Cable - LCD	0500-6000
4-40 x 1/4, System Interconnect Board, Mounting Screw	0163-0626
4-40 x 3/8 Modem Bezel Screw	0163-0153
6/32 x 3/8 Screws, Housings	0163-0912
AC Power cord with ferrite	1001-0195-01
AC Receptical Assembly	1001-0114
Access/Detect Cable	9500-0500
Adhesive Barrier, Speaker	9330-0324
AED Non-Pacing Membrane	1001-0187-01
Assembly Cable, EL (Yellow) Display	1001-0171
Assembly Cable, ECG out, Cable Sense	9500-0500
Assembly Capacitor & Choke (DSW)	1001-0134
Assembly Paddle Test Harness	1001-0102
B/W LCD Cable Assembly	1001-0173
Barrier, Isolation, Chassis, Print Head Cable	9330-0332
Barrier, Moisture, Speaker	9330-0305
Bezel, Modem, Center Open	9310-0512
Bracket, Speaker Mount	9320-0409

<b>Description</b>	<b>Part Number</b>
Brush, Static Dissipation	9340-0102
Bumper, Chassis, High Voltage, Capacitor	9330-0205
Cable Flex CRT. 312 TO System Board	9500-0519
Cable Laminate, 50 Pin	9500-0501
Cable, Laminate, ECG Connector to System Board	9500-0506
Cable, Recorder Print Head	9500-0400
Configuration Manual	9650-0201-01
Connector Assembly SPO2 (din connector)	1001-0160
ECG Connector 12 Lead	1001-0232
ECG Input Connector Assembly	1001-0132
ECG Laminated Cable	9500-0506
EL (Yellow) Display only	9355-0505
EL Display	1011-0029
FED Display	0208-0011
Front Panel Membrane Switch Assembly	1001-0135-01
Front Panel Screw, 6-32x 1 3/4	0163-0211
Gasket, Battery Pin	9330-0317
Gasket, Die Cut, Printer Motor	9340-0101
Gasket, MPPM Port and Motor Support, Chassis	9330-0303
Gasket, PCMCIA Connector, System Interconnect Board	9330-0312

<b>Description</b>	<b>Part Number</b>
Gasket, Print Head Cable	9330-0314
Handle (Part One)	9310-0545
Handle Insert (Part Two)	9310-0546
Isolator, High Voltage Module, Folded	1001-0146v
Isolator, Left Side Chassis	9310-0552
Keeper, 50 Pin Laminate Cable	9310-0573
Knob, Main Selector	9310-0521
Knob, Pace Control	9310-0520
Latch Paddle Release	9310-1514
Latch Spring	0190-0100
LCD Display	9355-0510-99
LCD Frame	9310-0595
Leafspring, Print Head	9320-0300
Lower Housing	1001-0124
Lower Latch, Paddle	9310-0511
Lower Latch, Pivot	9320-0065
Main Chassis	9310-2502-90
MFC Signal Cable	9500-0517
Motor Support (Teflon Piece)	9330-0355
M Series Main Label Set	9305-0527-01



<b>Description</b>	<b>Part Number</b>
Noncharger Substitute Board	9301-0311
O-ring, MFC, M Series	0160-6950
Paddle Shoes	1001-0148
Paper Tray Assembly	1001-0103
PCMCIA Conn., Assembly	3001-0101
PCMCIA Stabilizer Board	9320-0304
Pivot Rod, Printhead	9320-0401
Print Head Assembly	1001-0101
Print Head Support/Paper Tray Guide	9310-0510
Recorder Motor Assembly	1001-0104
Retainer, Latch Paddle Release	9310-1515
Rubber Foot x 4	0310-0311
Seal, Pace and Main Selector Knobs	9310-0548
Snap Rivet, Battery Interconnect Board	0163-1709
Spacer, Lower Latch, Paddle well	9320-0315
Speaker Assembly	1001-0115
Spring, Recorder Tray	0190-0101
Tape, Copper Adhesive, Roll (54')	0550-0037
Tape, Kapton, Roll (108')	0550-0003
Tape, Yellow Mylar, Roll (108')	0550-0125

<b>Description</b>	<b>Part Number</b>
Universal cable	1001-0196-01
Upper Housing	1001-0126
Upper Latch Seal	9330-0304

## Field Replacement Parts

<b>Description</b>	<b>Part Number</b>
Battery Interconnect for units EtCO <sub>2</sub>	9301-0303
HV Assembly Module (DSW) Defib Only	1001-0105-03
System Board 50 MHz 9301-0300-01	1001-0130-01
System Board 50 MHz (LCD only) 9301-0300-02	1001-0130-02
System Board 3/5 lead with audio 9301-0337-02	1001-0130-14
System Board 12 lead with audio 9301-0337-02 and 9301-0304	1001-0130-15
LCD Interface PWB for use with New Inverter	9301-0313
Biphasic Bridge Cap Assembly	1001-0182
Biphasic HV Module	1001-0181
AC Charger with EtCO <sub>2</sub> Heatsink	1001-0108-02
HV Assembly Module (DSW) Pace/Defib	1001- 0105
Battery Interconnect for units EtCO <sub>2</sub>	9301-0303

Description	Part Number
HV Assembly Module (DSW) Defib Only	1001-0105-03
Battery Interconnect (pre EtCO <sub>2</sub> )	9301-0302
Paper Sensor Board	1001-0131
System Interconnect	9301-0306
AC Charger Board with Heatsink	1001-0108
Control Board with Pace	9301-0312
LCD Display Interface Board	9301-0308
SpO <sub>2</sub> Board (Masimo)	1001-0158
Power Supply, Isolated, SpO <sub>2</sub>	1001-0159
DC Charger Board with Heatsink	1001-0107
100 Mhz AED/SpO <sub>2</sub> non 12 Lead (901-0300-04)	1001-0130-04
System PCB. (Biphasic only 9301-0300-03)	1001-0130-03
100 Mhz 12 Lead/SpO <sub>2</sub> (9301-0307)	1001-0130-05
System Interconnect 12 Lead Faxing	9301-030602
PCB Control with Pace Non-AED	9301-0312-05
Assy PCB Connector EtCO <sub>2</sub>	9301-0325
EtCO <sub>2</sub> Isolated power supply	1001-0015
EtCO <sub>2</sub> Isolated power supply (with EtCO <sub>2</sub> PCB)	3001-0103
PCB control no pace	9301-0312-02
PCB control AED with pace	9301-0312-03

<b>Description</b>	<b>Part Number</b>
PCB control AED no pace	9301-0312-04
12 Lead Based (9301-0337-01)	1001-0130-09
3/5 Lead With Audio 9301-03 37-02 and 9301-0304	1001-0130-10
12 Lead (9301-0307-01)	1001-0130-06
3/5 Lead/Biphasic (9301-0307-02)	1001-0130-07
System Interconnect Board with Faxing Capability	9301-0306-03
PCB control without pace non AED	9301-0312-06
PCB control with pace AED	9301-0312-07
PCB control without pace	9301-0312-08

# Chapter 5

## Functional Description

### Overview

This chapter provides functional descriptions of the components contained in the ZOLL M Series and the M Series options. Refer to the interconnect diagram that delineates the different components of the defibrillator.

This chapter includes:

#### **Main System Board**

- Main System Board Functions
- Power Supply
- User Interface
- Audio I/O Module
- ECG Front End
- Multifunction (MFE) Paddles
- Main CPU and EPU

#### **High Voltage Module**

- Defibrillator Charging and Discharging
- High Voltage Capacitor Module
- Pacer and Defibrillator Control Signals
- Internal Discharge Resistor Module

#### **AC/DC Charger**

- AC/DC Charger Module

### **System Interconnect Board**

Stripchart Recorder  
PCMCIA Slots  
Front Panel and Display

### **M Series Options**

12 Lead Monitoring  
Pulse Oximetry (SpO<sub>2</sub>)  
End Tidal Carbon Dioxide (EtCO<sub>2</sub>)  
Biphasic Module

## **Main System Board**

The M Series electrical circuitry consists of several functional modules. Each module is physically located on one or more of the printed wiring board assemblies (PWBA). In some cases, a functional module is distributed across several assemblies within the unit. The main components of the M Series include:

- Display
- Main System Board
- High Voltage Module
- AC/DC Charger
- Battery Interconnect Module
- System Interconnect Module
- High Voltage Capacitor

Some units are equipped with M Series options. These options include:

- 12 Lead ECG
- Pulse Oximetry (SpO<sub>2</sub>)
- End Tidal Carbon Dioxide (EtCO<sub>2</sub>)
- Biphasic Defibrillation Waveform

Refer to the M Series Interconnect diagram to identify unit components described in this manual.

MODULE	LOCATION
Main Central Processing Unit (CPU) and ECU	Main System PWBA
Pacer/Defib Charging and High Voltage Control (Defib/Pace)	High Voltage Module Biphasic Bridge Module
Internal Discharge Resistor	HV Module
AC or DC Power Supply/Battery Charger	AC Charger PWBA DC Charger PWBA Main System PWBA
ECG Front End Signal Acquisition	Main System PWBA
Battery Pack or Smart Battery™	Battery Interconnect PWBA
User Interface and Controls	Controls PWBA
PCMCIA Interface	System Interconnect, Main System
Stripchart Recorder	System Interconnect, Main System
Audio I/O (optional)	Main System PWBA Audio Display PWBA Controls PWBA
SpO <sub>2</sub> (optional)	Pulse Oximetry PWBA Isolated Power Supply PWBA
EtCO <sub>2</sub> (optional)	EtCO <sub>2</sub> PWBA Isolated Power Supply PWBA
12 Lead (optional)	Main System PWBA

## Main System Board Functions

The Main System Board contains the major computing and control elements for the M Series unit. The printed wiring board assembly (PWBA) receives signals from the front panel control switches, ECG input connectors and functional modules, such as the AC/DC charger, pacer/defibrillator modules, stripchart recorder and PCMCIA card interfaces, and if applicable, the SpO<sub>2</sub> and EtCO<sub>2</sub> modules. The Main System Board monitors and processes these input signals to produce other signals that: 1) control the operation of other modules within the system; 2) drive the unit's front panel display and audio outputs and; 3) store data for retrieval via Summary Reports, PCMCIA cards and/ or a modem.

The electronic circuitry and software contained on the main system board performs the following major M Series functions:

1. Main CPU and memory.
2. ECG signal acquisition and processing for ECG leads including,
  - A/D conversion.
  - ECG signal filtering.
  - QRS detection.
  - Implanted pacemaker detection.
  - Heart Rate counting.
  - Shockable ECG rhythm analysis.
3. Data communications with and control over the SpO<sub>2</sub> and EtCO<sub>2</sub> modules.
4. Control over and safety monitoring of pacer and defibrillator functions performed by the High Voltage module and Biphasic Bridge module.
5. Physiological alarm processing.
6. Control switch monitoring for the front panel, accessory connection monitoring and control over the system response to switch activation or accessory connection to the M Series.
7. Format and updating of the front panel display.
8. Primary power supplies for the unit.
9. Audio output generation and control (e.g., alarms, voice prompts, warning tones).
10. Audio signal processing, data compression and storage of voice recording data.
11. Real time clock and other time keeping functions.
12. Summary Report, 12 Lead ECG Reports data storage.



- 13.* Monitoring of battery status and control over Battery Charging functions performed by the AC or DC Charger PWBA.
- 14.* Data transmission to and control over the System Interconnect PCB functions, including the stripchart recorder and the PCMCIA slot functions.

## Power Supply

The power supply converts DC power from a removable battery or the AC/DC Battery Charger module to voltages required by the M Series hardware.

The power supply circuit converts the raw battery or the Charger PWBA output voltages of +8.5 VDC to +16 VDC into the voltages shown in the table below, including load and line regulation.

VOLTAGE	DESCRIPTION	VOLTAGE VDC (Nominal)	COMMENTS
FUSE_PWR	Fused Input Power from Battery/Charger	12	
SW_PWR	Switched Input Power after Power Switch	12	
3VDD	+ 3.3 VDC Power for Digital circuits	3.3	Switching @ 300 kHz
5VDD	+ 5.0 VDC Power for Digital circuits	5.0	Switching @ 300 kHz
-5VSS	- 5.0 VDC Power for Analog circuits	-5.0	Linear
12VEE	+ 12.0 VDC Power for Analog circuits	12.0	Linear
15VDD	+ 15.0 VDC Power for VPP and 12VEE	15.0	Switching @ 300 kHz
LCD_BS	LCD BIAS Power for LCD display	- 18	Switching @ 100 kHz
3_3REF	+ 3.3 VDC Reference	3.3	Linear

## ECG Front End

The ECG front end provides an electrically isolated serial interface between the main system board functions and patient interface ECG connectors. It performs the following:

- Analog ECG amplification and signal conditioning.
- Pacemaker pulse detection.
- ECG signals acquisition and analog to digital conversion.
- ECG 3/5/12-lead detection.
- ECG leads off detection.
- Front-end defibrillator protection.
- Isolated power conditioning.
- Patient impedance measurement via MFE or paddles.

## Multifunction Electrode (MFE)/PADS (System Board and High Voltage Module)

Selected for optimal performance for the application, a dedicated ECG amplifier with a limited bandwidth processes the signal. It is then chopper modulated and coupled to the system side via an isolation transformer. On the system side, the signal is synchronously demodulated, converted by a 10 bit A/D at 250 samples per second and digitally processed by the main control unit of the system board.

To measure thoracic impedance, a high frequency (HF) measuring current passes through the patient's chest and measures the resulting voltage across the electrodes. After amplification, the impedance signal is synchronously demodulated. It is then converted to a stream of pulses with frequency proportional to the measured impedance.

## CPU and EPU

The Main System Board contains two microprocessors. A Motorola HC-11 single chip microprocessor is used to acquire, convert and process ECG signals, (ECU). An Hitachi SH-3 RISC microprocessor acts as the system's main CPU. The SH-3 CPU and has an integrated on-chip multiplier, a cache memory, a memory management unit as well

as data protection and virtual memory functions. It also has a timer, a real time clock, an interrupt controller, a serial communication interface (SCI), and other peripheral functions necessary for the system operation. The memory circuitry includes Flash ROM, internal flash non-volatile memory and DRAM.

The EPU acquires ECG data and runs the A/D convertor that sends data in the form of a serial stream to the CPU.

## High Voltage Module

The High Voltage (HV) module includes the high voltage circuitry required for pacing and defibrillation, including the defib charge circuitry, solid state patient relay, safety relay, defib capacitor, defib choke and front end protection circuitry for the MFC ECG. There are two different types of HV Modules: a damped sinusoidal waveform HV module and a Biphasic (HV) module.

The following table describes the high voltage board components:

<b>Component</b>	<b>Function</b>
Solid State Patient Relay	Controls the delivery of therapeutic energy to patient.
Safety Relay	Discharges Defib capacitor into the internal discharge resistor when defibrillator is not in use.
Defibrillator Capacitor	Stores energy for therapy.
Defibrillator Choke	Conditions waveform delivered to the patient.(DSW)
Front End Protection Circuitry for the MFC ECG	Protects ECG front end against defibrillator pulses.
Monophasic HV	Provides damped sinusoidal waveform therapeutic energy. (Monophasic units only.)
Biphasic HV	Provides biphasic waveform therapeutic energy. (Biphasic units only.)

## Defibrillator Charging and Discharging

The defibrillator charges and discharges high voltage capacitor energy. A user can initiate a charge in three ways by 1) pressing the charge button on the front panel; or (2) pressing the charge button on the paddles; or (3) configuring unit to charge automatically when it detects a shockable rhythm following an ECG analysis. To initiate a discharge, a user depresses both shock buttons on the paddles or depresses a single shock button on the front panel.

The defibrillator circuit charges the high voltage capacitor to the energy level the user specifies. This circuit also provides feedback to the main system board on the high voltage capacitor's voltage level and discharges the high voltage capacitor energy through paddles or the universal cable. The defibrillator portion of the high voltage circuitry is active only when the front panel selector switch is set to DEFIB.

### Charging

The charging process starts when the Main System Board detects a charge request. The defibrillator circuits begin charging the high voltage capacitor to the target voltage or energy that the user selects on the front panel display. The Main System Board continuously monitors the capacitor voltage signal to ensure that the high voltage capacitor charges at the proper rate. When the target voltage is reached, the Main System Board initiates a continuous beeper tone to indicate that defibrillator is ready to discharge. The target energy level displays on the display screen.

The defibrillator holds the energy for 60 seconds for manual units and 15 seconds for AED units, refreshing the energy level as necessary. An intermittent beep tone sounds during the last ten seconds (five seconds for AED unit) of the hold period. After the 60 second period, if the defibrillator has not been discharged, the energy is dissipated into the internal discharge resistor by closing the safety relay (XSAFREL). The unit discharges internally and displays a warning message if it is not functioning properly.

Unlike previous ZOLL designs that isolated the patient from defib circuitry via an electromechanical patient relay, the M Series utilizes a bank of silicon-controlled rectifiers (SCRs). As the defibrillator capacitor is charged, the voltage is monitored via R1 - R4, which drive differential amplifiers referred to the system ground. These resistor dividers split the capacitor voltage more or less equally above and below ground in order that the positive capacitor terminal is approximately 2500 volts above ground, and the negative capacitor terminal is approximately 2500 volts below ground (at 360J setting). The voltage at the patient electrodes is set by the divider RN1 and RN2. These networks are each 5X 25 M (125 M total) whose total resistance is specified to be 125 M  $\pm$  1%. As a result, the patient is nominally at ground and the hot switch bank is split into a 'positive' side and a 'negative' side.

### Discharging

Initiating a discharge provides voltage to the solid state patient relay and notification to the Main System Board through the PADMON signal. The Main System Board then controls activation of the solid state patient relay (for DSW only). Energy delivered to the patient goes through a wave shaping inductor to create a defibrillation waveform compliant with AAMI Standards. When the patient discharge SCRs are deactivated, the safety relay closes to internally dissipate any remaining energy.

If the M Series is in the self test mode, the energy is delivered internally. The microprocessor calculates the actual delivered energy from the current waveform and displays a TEST OK message on the display, if the self test meets the appropriate criteria. If the criteria are not met, a TEST FAILED message displays.

## High Voltage Capacitor Monitor

Before charging the defibrillator, the High Voltage Capacitor monitor runs a self test to check the pace relay. The pace relay controls the high voltage circuitry configuration either for generating pace pulses or for charging the high voltage defibrillation capacitor.

The defibrillator capacitor is shunted for safety reasons with a resistor and relay to internally dissipate any energy remaining. When the Main System Board initiates a charge, this relay opens by providing a low level on signal XSAFREL. The safety relay is a biased reed switch. The relay is driven by Q318.

The pace relay driver is a grounded source switch Q308 that is biased on by R593. It is held off by Q330 when XPACEREL is '1' false. When XPACEREL comes true, Q330 is turned off, and Q308 is no longer clamped off.

The high voltage capacitor is charged by converting the system battery voltage to a pulsed high voltage by way of transformer T1. The basic operating frequency signal that is used to switch transistor Q1 providing current in the primary windings of the transformer T1 originates in the system board's gate array.

When the high voltage capacitor is charging, the Main System Board independently monitors the capacitor voltage through signal VMON. If the Main System Board detects an improper level, it halts operation by setting SAFE high. This disables the SCR discharge circuitry and flyback transformer drive.

The solid state patient relay discharges via the signal PATREL\_DRV generated by XPATREL and Q304, Q323, and Q322. PATREL\_DRV is disabled when XPACE\_SEL is at a logic low.

When the solid state patient relay activation completes, the Main System Board releases the XPATREL signal. Several hundred milliseconds later, the safety relay closes to ensure the high voltage capacitor energy is completely dissipated.

The Pacer circuit produces and delivers user-controllable pace pulses to the pacing electrodes. To initiate pacing, the front panel switch is turned to PACER and the OUTPUT and RATE controls are set. Pacing current amplitude is constant during the pulse and is determined by the position of the front panel PACER OUTPUT dial. Pacing pulse rate is determined by the position of the front panel PACER RATE dial. The pacing pulse duration is fixed at 40 milliseconds.

## Pacer/Defibrillator Control Signals

The Pacer/Defibrillator Control charges the high voltage capacitor to a voltage requested by the main system board in response to user energy selections. It delivers defibrillator energy to the patient through the patient connector to the paddles and pacer electrodes or multi-function electrodes (PADS). This control also generates pacing pulses at rates and amplitudes requested by the main system board in response to user selections, controls damped sinusoidal waveform and biphasic waveform defibrillation, and measures pace current and high voltage capacitor voltage by two independent channels.

The following signals control the operation of the Pacer/Defibrillator subsystem:

ANALOG VOLTAGE	OPERATION	COMMENT
VCAP	Analog voltage spanning 0 - 2.5 V for 0 - 5000V capacitor voltage.	Used by the defib charging controller.
VMON	Analog voltage spanning 0 - 2.5 V for 0 - 5000V capacitor voltage.	Used by the defib monitor.
VSENS	Pace duty cycle voltage, scaled as 0 - 2.5 V for 0 - 100% duty cycle.	When multiplied by the pulse width (as read from PW_READ) battery voltage is proportional to the actual pace current.
VCTL	Analog control voltage scaled 0 - 2.5 V for pace current of 0 - 140 ma.	Only active in pace mode.

<b>ANALOG VOLTAGE</b>	<b>OPERATION</b>	<b>COMMENT</b>
FET_MON	Analog voltage monitors the condition of the discharge transistors.	Provides a signature voltage in case of a fault.
PAT_CUR	Bound on the range of 0 - 2.5 V and accommodates defib currents of -50 to +100 A	Analog signal representing the patient current during a defib discharge.

<b>LOGIC CONTROL SIGNAL</b>	<b>OPERATION</b>
XPWR_ENABLE:	This logic signal from the gate array enables the charging circuit when true, and inhibits the charging circuit when false.
SAFE:	This logic control signal is generated by the Main System Board to halt the pace/defib function in the event of a detected fault.
SHUTDOWN:	This logic signal is true during reset and fault conditions. (VCC error, watchdog error, etc.) and halts operation of the PD generator.
XPACE_ON:	This logic signal is generated by an optocoupler, and indicates that the pace output circuit is active. It is '0' true when pace current is flowing.
XSAFREL:	Logic signal from the GA that operates the safety relay when '0' true.
XPACEREL:	Logic signal from the GA that operates the pace relay when '0' true. XPATREL: Logic signal from the GA that operates the solid state patient relay when '0' true.
XPAT_ENABLE:	Logic signal from the processor controlled by the monitor that grants operation of the solid state patient relay. It is false during pacing.



LOGIC CONTROL SIGNAL	OPERATION
XPACE_SEL:	Hardware only signal from the front panel switch that is at '0' during pace. Used as an additional safety interlock on the solid state patient relay so that operation of the relay during pacing is additionally disabled.

## Internal Discharge Resistor Module

The Internal Discharge Resistor Module contains the internal discharge resistor, and a means for dissipating the heat generated by the internal discharge.

## AC/DC Charger Module

The AC/DC Power/Battery Charger provides a universal (IEC 320) connection to the AC mains or to a DC source, input line filtering and double-pole fusing (for a mains input), AC-DC and DC-DC conversion and isolation barrier between the M Series and power sources. This module also provides the power necessary to run the M Series in any mode of operation, as well as providing additional charging current to the battery. When the M Series is turned off but connected to an external AC or DC source, the charger module controls battery charging currents and voltages needed to charge the M Series battery. These voltages and currents are controlled in response to the main CPU signals that manage the battery charging process.

## System Interconnect Module

The system interconnect PWBA receives signals from the Main System Board and in turn controls operation of the stripchart recorder and PCMCIA functions.

## Stripchart Recorder

The Stripchart Recorder module includes a microprocessor, serial interface to the main system board and circuitry which drive the stripchart recorder's motor and printhead in response to the main CPU signals. Based upon signals sent by the main CPU, the recorder's main processor drives the recorder stripchart motor, formats data for printing on the chart and drives the printhead. It also detects when the sensor drawer is not properly fitted into the unit, when the paper supply is out and needs to be refilled and the print head temperature.

## PCMCIA Slots

The PCMCIA interface module supports two PCMCIA slots which accept Type I and/or Type II PCMCIA cards. These cards may be read or written to. Data sent by the main CPU is passed to the installed PCMCIA card via the system interconnect PWBA.

## Front Panel and Controls PWBA

User Interface Module provides several functions that enable the user to operate the unit. The user interface has a display monitor and three rotary selector switches. One selector switch is for three modes: pacer, monitor and defibrillation. The two other knobs are for pacer output and pacer rate. The unit interface also has specific buttons for defibrillation, including the ENERGY SELECT button, the CHARGE button, the ANALYZE button and the SHOCK button. The five softkeys underneath the display provide specific operations depending on the unit's configuration. The other push buttons (from left to right) are used for volume control, monitor illumination, summary report, and code markers. The CHARGER ON indicator displays the status of the unit's power supply.

This input module on the front panel and the Main System Board provides a beeper for the AC/DC Power/Battery Charger.

The Controls PWBA is physically located in the front panel assembly. Units that are equipped for voice recording include a microphone and audio signal conditioning circuitry on the Controls PWBA.

## M Series Options

This section describes four options for the M Series unit.

## Isolated Power Supply Module

The Isolated Power Supply Module provides electrically isolated power to the EtCO<sub>2</sub> and SpO<sub>2</sub> modules. It also provides the electrically isolated serial communications and isolated control signals between the EtCO<sub>2</sub> and SpO<sub>2</sub> modules and the main system PWBA.

## 12 Lead Option

The ZOLL M Series 12 lead option is used to acquire ECG data needed to assist in the diagnosis of myocardial infarction (“heart attack”), often caused by a coronary artery occlusion. The 12 lead ECG can be viewed on the display one lead at a time in monitoring and diagnostic bandwidths and printed in the standard 4x3 format with 12 simultaneously acquired leads.

The 12 lead option provides for the recording, printing and automated analysis of 12 lead ECG using GE Marquette 12SL™ Analysis and supports the transmission of these reports by fax to a remote location, such as a hospital. In the pre-hospital environment, the 12 lead reports can be faxed to a physician as the patient is en route to the Emergency Department. As a result, the physician can initiate hospital accommodations immediately, such as activating the staff of the cardiac catheterization lab, prior to the patient’s arrival and subsequent treatment. Or the patient may be treated in the pre-hospital environment with thrombolytic agents.

The 12 lead cable is required to produce 12 lead reports. M Series unit must have the 12 lead option installed. All limb leads and at least one V-lead must be connected to initiate a 12 lead acquisition. Printed 12 Lead bandwidth is user configurable to be either 0.05-150 Hz (per AAMI EC11) or 0.05-40 Hz. The 0.05-40 Hz bandwidth selection is used to reduce noise artifact in the high end of the diagnostic frequency range. Reports can be printed in a standard 4x3 or Cabrera format. Faxed reports can be configured in a 2x6 format in addition to 4x3 and Cabrera formats.

The GE Marquette 12SL™ Analysis algorithm provides measurements of the 12 lead waveforms along with interpretive statements. The algorithm is interpretive, not “diagnostic.” (A physician should always confirm interpretive statements. A diagnosis requires a complete clinical assessment including other modalities, such as a physical examination.) 12SL™ produces global waveform measurements as well as a measurement matrix containing measurements on each lead. Both the interpretive statements and measurement matrix are configurable to be printed or not printed.

The acquired 12 lead with 12SL™ may be faxed to a remote location using landline or cellular phone technology. Several PCMCIA fax modems are supported and the modem determines the specific phone compatibility. Cellular phone support includes analog AMPS phones (in the U.S.) and GSM phones (internationally). The M-Series supports Group 3 facsimile, Class 1 and Class 2. 12 lead reports may be re-printed or re-transmitted using the Patient Records capability. Individual patient records may be selected based on patient ID, date, and time.

### Pulse Oximetry (SpO<sub>2</sub>)

The ZOLL M Series pulse oximetry option enables the user to continuously, noninvasively, and painlessly monitor the percentage of oxygen saturation of arterial hemoglobin at a peripheral measurement site (i.e. foot, toe or finger.)

The oximetry sensor contains two light emitting diodes, or LEDs, that transmit red and infrared light through the body's extremities. A photodetector receives the transmitted light. Oxygen saturated blood absorbs light differently than unsaturated blood. Thus the amount of red and infrared light absorbed by the blood flowing through a suitable peripheral area of the body, such as the finger in adults and the foot in neonates, can be used to calculate the ratio of oxygenated hemoglobin to total hemoglobin in the arterial blood. The monitor displays this ratio as percent SpO<sub>2</sub>. Normal values typically range from 95% to 100% at sea level.

The M Series uses a Masimo® Pulse Oximetry Circuit Board which features a fundamentally distinct method of acquiring, processing and reporting arterial oxygen saturation and pulse rate. The M Series' SpO<sub>2</sub> module (Masimo® Circuit Board) connects to the Masimo sensors and reports monitoring results (oxygen saturation, pulse rate, pulse waveform, etc.) via a serial digital interface to the M Series system board. The M Series system provides isolated DC power and serial communication to the SpO<sub>2</sub> Board via the Isolated Power Supply board.

### End Tidal Carbon Dioxide (EtCO<sub>2</sub>)

The ZOLL M Series EtCO<sub>2</sub> Novamatrix technology and Capnostat® sensor option continually and noninvasively monitors the patient's carbon dioxide in respiratory gases and from these measurements computes End Tidal CO<sub>2</sub> and respiration rate. The unit can display and print a recording of EtCO<sub>2</sub> readings, respiration rates, and capnograph waveforms. In addition, the unit can configure an alarm to sound when the unit detects EtCO<sub>2</sub> values and respiration rates that are above or below acceptable ranges as set by the user. This option is intended for use in all critical monitoring environments including ventilator support, patient transport, and anesthesia and is intended for monitoring all patient types, including adult, pediatric, and neonatal.

The EtCO<sub>2</sub> option incorporates, without modification, the Novamatrix respiratory carbon dioxide technology, including a printed circuit board, and a patented and proprietary mainstream CO<sub>2</sub> sensor and reusable/disposable airway adapters. Carbon dioxide measurements are monitored by using a solid-state infrared Capnostat<sup>®</sup> sensor which works on the principles of infrared absorption.

The M Series system board software transmits commands to the Novamatrix board to set parameters (averaging mode, oxygen, and nitrous oxide compensation) and retrieves CO<sub>2</sub> waveform, EtCO<sub>2</sub>, and respiration rate data. The M Series software formats the data for output to the display and strip chart recorder. In addition, the M Series software performs range checking of the Novamatrix data to determine the presence of low or high EtCO<sub>2</sub> and respiration rate alarm limit violations. If the limits have been violated, the software issues audio and visual alarms.

The Capnostat<sup>®</sup> mainstream sensor is attached to an airway adapter that connects to an endotracheal tube, similar airway, or disposable mouthpiece. The sensor generates infrared light and beams it through the airway adapter to a detector on the opposite side of the airway. As a result of respiration, CO<sub>2</sub> flowing through the airway adapter absorbs some of this infrared energy. The monitor relates the amount of detected energy to the amount of CO<sub>2</sub> in the airway adapter. The CO<sub>2</sub> display reflects the maximum concentration of CO<sub>2</sub> detected during expiration. End-Tidal Carbon Dioxide (EtCO<sub>2</sub>) displays as a numerical value in millimeters of mercury (mmHg), percent (%), or kilopascals (kPa) on the unit. In addition, a capnogram waveform may display underneath the ECG waveform.

The EtCO<sub>2</sub> option uses Non-Dispersive Infrared (NDIR) absorption and dual wavelength ratiometric-True Single Beam Optics. The Capnostat<sup>®</sup> sensor contains the NDIR light source. Carbon dioxide, flowing in the airway adapter as a result of respiration, absorbs some of this light energy.

## Biphasic Waveform

The M Series Biphasic Waveform Defibrillator produces a Rectilinear Biphasic™ waveform similar to the chart shown at the end of this section. The electrical energy is delivered in two successive current phases of opposite polarity. As compared to the monophasic waveform, the biphasic waveform typically defibrillates with substantially less current than the earlier monophasic waveform used by most defibrillators.

The M Series Biphasic Defibrillator system consists of circuitry and software located on the following assemblies:

1. Main System PWBA.
2. High Voltage Module.
- 3) Biphasic Bridge PWBA Module.


The High Voltage module controls the monitoring, charging, and internal discharges of the defibrillator capacitor. The Bridge module controls capacitor discharge, including waveform polarity control, patient current and voltage monitoring during discharge and real-time control over defibrillator internal impedance to create the rectilinear waveshape. The main CPU controls the waveform timing and resistor switching performed by the Bridge Module based upon measured patient transthoracic impedance.

The Biphasic Bridge Assembly is constructed as two printed wire board assemblies (PWBA) containing the power and isolation circuits required to deliver a biphasic defibrillation pulse. These boards produce a defibrillation pulse consisting of a positive current pulse followed by a negative current pulse. The positive portion of the pulse is shaped to be rectilinear by switching resistors in series with the patient to compensate for droop in the capacitor voltage as it delivers energy to the patient load.

The ZOLL Biphasic units produce a rectilinear waveform whose shape remains essentially constant from patient to patient. The rectilinear biphasic waveform consists of a 6 millisecond, essentially constant current first phase followed by a 4 millisecond, truncated exponential second phase. The first and second phases of the defibrillation waveform are of opposite polarity and their amplitudes vary based on the user selected therapeutic energy level. The initial amplitude of this waveform's second phase is approximately equal to the first phase's final amplitude. This waveform has an integrated patient impedance measurement sensing pulse at the beginning of the waveform. The positive and negative phases are separated by 100  $\mu$ sec.

Electronics and software control the shape of the waveform's first phase and compensate for different transthoracic impedances to maintain an essentially constant current throughout the first phase. When the highest energy setting is selected and patient impedance exceeds 85 ohms, the first phase of the waveform will droop. All other waveform parameters (phase duration, inter-phase delay and integrated impedance measurement sensing pulse) remain the same.

The following Rectilinear Biphasic Waveform is produced when the M Series with Biphasic option is discharged into a 50 ohm load at the default energy setting of 120 joules. The vertical axis is in amperes; the horizontal axis is in milliseconds. (See following diagram.)







# Appendix


## Overview

This appendix includes:

- Interconnect Diagram for the M Series Monophasic Unit
- Interconnect Diagram for the M Series Biphasic Unit
- Interconnect Diagram for the M Series CCT Biphasic Unit


This appendix also includes the ZOLL M Series Maintenance Tests Checklist. Photocopy the checklist and use the copy to record the results of the maintenance tests performed on the M Series equipment.

# Interconnect Diagram for the M Series Monophasic Unit





# Interconnect Diagram for the M Series CCT Biphasic Unit



# Index

## Numerics

200J MAX BIPHASIC message 58  
2PFE SpO2 Simulator 4  
50J MAX message 58

## A

AC Charger  
  photograph 78  
AC or DC Power Supply 117  
AC/DC Charger module 127  
ADJUST ECG message 58  
Advisory message test 30  
ANALYSIS 58  
ANALYSIS HALTED message 58  
ANALYSIS RESTARTED message 58  
AUDIO FAULT 13 message 59  
AUDIO FAULT 136 59  
AUDIO FAULT 136 message 59  
Audio I/O 117  
AUDIO NOT RECORDING message 59  
AUDIO QUEUE FULL message 59

## B

BATT 59  
BATT HIGH CURRENT message 59  
BATT LOW CURRENT message 59  
BATT LOW VOLTAGE 59  
BATT LOW VOLTAGE message 59  
BATT OVERCHARGE message 59  
Battery Interconnect Board  
  removing and reinstalling 90  
Battery Interconnect Board Assembly 78  
Battery Pack 117  
Biphasic Bridge Assembly 79  
Biphasic Module 122  
BRIDGE TEST FAILED message 60

## C

CABLE FAULT 60  
CABLE FAULT message 60  
CANNOT CHARGE 60  
CANNOT CHARGE message 60  
CHECK 61  
CHECK PADS/POOR PAD CONTACT message 55  
CHECK PULSE 61  
CHECK PULSE message 61  
CHECK RECORDER message 61  
CHECK SPO2 SENSOR message 62  
CHECK SPO2 SITE message 62  
CLOCK FAULT 11 message 62  
CLOCK FAULT 12 message 62  
CLOCK FAULT 13 62  
CLOCK FAULT 13 message 62  
CO2 COMM ERROR message 62  
CO2 SENSOR WARM UP 62  
CO2 SENSOR WARM UP Message 62  
CO2 SENSOR WARM UP message 62  
CONFIRM MANUAL MODE message 62  
Contacting Technical Service vii  
Control Board Panel 84, 85  
Conventions vii

## D

DEFIB 55  
Defib Capacitor 122  
DEFIB DISABLED message 62  
DEFIB FAULT 108 message 64  
DEFIB FAULT 109 message 64  
DEFIB FAULT 111 message 64  
DEFIB FAULT 71 message 63  
DEFIB FAULT 72 message 63  
DEFIB FAULT 76 message 63  
DEFIB FAULT 77 message 63  
DEFIB FAULT 78 message 63  
DEFIB FAULT 79 message 63  
DEFIB FAULT 80 message 63  
DEFIB FAULT 81 message 63  
DEFIB FAULT 84 message 63  
DEFIB FAULT 85 message 63  
DEFIB FAULT 86 message 64  
DEFIB FAULT 87 message 64  
DEFIB FAULT 94 message 64

DEFIB FAULT 95 message 64  
DEFIB FAULT 96 message 64  
DEFIB NOT CHARGED message 64  
DEFIB PAD SHORT 55  
DEFIB PAD SHORT message 64  
Defibrillator  
  charging and discharging 123  
  energy levels 123  
  Defibrillator Choke 122  
  Defibrillator self test 23  
DISABLE SYNC message 65  
DISCHARGE FAULT message 65  
Discharging 124  
Display 83  
Disassembly procedures  
  modules to disassemble 78

## E

ECG FAULT 4 message 65  
ECG FAULT 5 message 65  
ECG Front End 121  
ECG Front End Signal Acquisition 117  
ECG LEAD OFF message 65  
ECG LEADS OFF message 54  
ECG TOO LARGE message 65  
ECG TOO SMALL message 65  
ECG V LEADS OFF message 65  
ECG V2 LEAD OFF message 65  
ECG V3 LEAD OFF message 65  
ECG V4 LEAD OFF message 66  
ECG V5 LEAD OFF message 66  
ECG V6 LEAD OFF message 66  
ECU CRC FAULT message 66  
ENTER ACCESS CODE message 66  
ERASING REPORT message 66  
Error Message  
  50J MAX 58  
  ANALYSIS HALTED 58  
  ANALYSIS RESTARTED 58  
  AUDIO FAULT 136 59  
  AUDIO NOT RECORDING 59  
  BATT HIGH CURRENT 59  
  BATT OVERCHARGE 59  
  BRIDGE TEST FAILED 60  
  CABLE FAULT 60  
  CANNOT CHARGE 60

CHECK RECORDER 61  
CLOCK FAULT 12 62  
CLOCK FAULT 13 62  
CONFIRM MANUAL MODE 62  
DEFIB DISABLED 62  
DEFIB FAULT 108 64  
DEFIB FAULT 109 64  
DEFIB FAULT 111 64  
DEFIB FAULT 71 63  
DEFIB FAULT 72 63  
DEFIB FAULT 76 63  
DEFIB FAULT 77 63  
DEFIB FAULT 78 63  
DEFIB FAULT 79 63  
DEFIB FAULT 80 63  
DEFIB FAULT 81 63  
DEFIB FAULT 84 63  
DEFIB FAULT 85 63  
DEFIB FAULT 86 64  
DEFIB FAULT 87 64  
DEFIB FAULT 94 64  
DEFIB FAULT 95 64  
DEFIB FAULT 96 64  
DEFIB NOT CHARGED 64  
DEFIB PAD SHORT 64  
DISABLE SYNC 65  
DISCHARGE FAULT 65  
ECG FAULT 5 65  
ECG LEAD OFF 65  
ECG TOO LARGE 65  
ECG V2 LEAD OFF 65  
ECG V3 LEAD OFF 65  
ECG V4 LEAD OFF 66  
ECG V6 LEAD OFF 66  
ECU CRC FAULT 66  
ENTER ACCESS CODE 66  
ERASING REPORT 66  
EtCO2 COM ERROR 66  
FAX DIALING 66  
FAX DONE 66  
FAX PREPARING 66  
FAX SENDING 66  
INSERT CARD 66  
LOW BATTERY 67  
NO QRS DETECT 67  
RETRY ANALYSIS 70

EtCO2  
Isolated Power Supply Assembly 129  
EtCO2 COM ERROR message 66  
EtCO2 Isolated Power Supply Assembly 117

### F

FAX DIALING message 66  
FAX DONE message 66  
FAX PREPARING message 66  
FAX SENDING message 66  
Front End Protection Circuitry  
High Voltage Module 122

### G

Gel droop 56

### H

High Voltage Capacitor monitor 124  
High Voltage Module 122  
Reinstalling 91, 94, 96  
High Voltage Module Assembly 122  
removing 94  
High Voltage/Charger Assembly  
removing 91  
HV 122

### I

INSERT CARD message 66  
Interference  
radio and cell phone vi  
Internal Discharge Resistor Module 127  
Components of Charge 127  
Isolated Power Supply  
EtCO2 module 78  
function 129

### L

Logic Control Signal 126  
LOW BATTERY message 67

### M

M Series 78

description 115  
main components 116  
Main Selector Knob 76  
Main System Board 117  
Maintenance Tests  
equipment 4  
Messages 55  
Module types 78  
Monitoring  
User interface module 128  
Monophasic  
High Voltage Module 122

### N

NO QRS DETECT message 67  
NO SHOCK AD message 67  
NOISY ECG message 67  
Notch filter test 20

### O

Ordering Replacement Parts  
part numbers 112

### P

PACER 125  
Pacer Defibrillator function 125  
PACER OUTPUT knob  
pacing 125  
Pacer/Output/Rate knob 76  
Paddle Release Latch  
Reinstalling 105  
Removing 105  
PCMCIA card slot  
removing and reinstalling 104  
PCMCIA Interface 117  
Power Supply  
MCU 120  
Power supply  
function 120  
Print Head Assembly 103  
Printer/Recorder Motor  
removing and reinstalling 100

## R

- Radio frequency emissions vi
- Radio Frequency Interference vi
- reinstall 91, 94, 96, 105
- Reinstalling
  - Control Board 85
  - High Voltage Module 91, 94, 96
  - Paddle Release Latch 105
  - System Interconnect Board 98
  - Upper Housing Assembly 86
  - ZIF Keeper 82
  - remove 92
  - Removing 96, 98
- Battery Interconnect Board 90
- High Voltage Capacitor Assembly 96
- High Voltage Module Assembly 94
- High Voltage/Charger Assembly 91
- Paddle Release Latch 105
- Printer/Recorder Motor 100
- System Interconnect Board 98
- RETRY ANALYSIS message 70

## S

- Safety Consideration v
- Safety Precautions 77
- Safety Relay 122
- Safety Warnings v
- SELECT 30J FOR TEST message 70
- SELECT DEFIB MODE message 70
- SELECT LIMB LEADS message 70
- SELECT PADS message 71
- Service
  - Address viii
  - Service Policy Warranty vii
  - United States
  - Outside of the United States vii
- SET CLOCK message 71

- SET PACE RATE message 71
- SHOCK ADVISED message 71
- Solid State Patient Relay
- High Voltage Module 122
- SpO2 and EtCO2 modules with bracket photograph 78
- SpO2 Isolated Power Supply Assembly 117
- SpO2 Module (without bracket) photograph 79
- STAND 72
- STAND CLEAR 72
- Storing pads 56
- Strip Chart Recorder 117
- System Board Assembly 87
  - photograph 78
- SYSTEM FAULT 2 message 72
- SYSTEM FAULT 5 message 72
- SYSTEM FAULT 6 message 72
- SYSTEM FAULT 7 message 72
- System Interconnect Board 98
  - photograph 78
  - removing and reinstalling 98

## T

- Technical Service Department
  - Address 107
- Test
  - 3, 5, and 12 Leads 10
  - Advisory Message 30
  - Calibrating Pulses on Strip Chart 19
  - Check Summary Report 29
  - Defibrillator Self Test 23
  - EtCO2 Monitor 36, 51
  - Front Panel Button 7
  - Heart Rate Alarm 21
  - Heart Rate Display 18
  - Leakage Current 15

- Notch Filter 20
- Pacer 31
- Paddles 16
- Power Supply 11
- Shock 26
- SpO2 Monitor 33
- Synchronized Cardioversion 25
- Troubleshooting
  - Displayed HR not accurate 57
  - Electronic interference 57
  - Wandering baseline 57

## U

- User Interface Module 128

## V

- V LEADS OFF Messages 55
- Voltage
  - Power supply 120

## W

- Warnings v

## Z

- ZERO CO2 ADAPTER message 73
- ZERO CO2 SENSOR message 73
- Zoll M Series Text Messages 54
- ZOLL Medical Corporation
  - Address
  - Ordering parts for M Series 107
  - Service request number viii
- ZOLL Technical Service Department vii
  - Fax number vii
  - Telephone number vii





# ZOLL M Series Maintenance Test Checklist

Serial No. _____	Location _____	Date _____
Tester _____	Signature _____	

Use this checklist to record the results of the M Series maintenance tests, and keep it for your records.

**Result of Check:**

- No action required
- Minor problems corrected
- Disposable supplies replaced
- Major problems identified (unit out of service)

**Additional Remarks**

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

**1.0 Physical Inspection**

	Pass	Fail	N/A
1.1	0	0	0
1.2	0	0	0
1.3	0	0	0
1.4	0	0	0
1.5	0	0	0
1.6	0	0	0
1.7	0	0	0
1.8	0	0	0
1.9	0	0	0
1.10	0	0	0
1.11	0	0	0
1.12	0	0	0
1.13	0	0	0
1.14	0	0	0
1.15	0	0	0

**2.0 Front Panel Button Test**

	Pass	Fail	N/A
2.1	0	0	0
2.2	0	0	0
2.3s	0	0	0
2.4	0	0	0
2.5	0	0	0
2.6	0	0	0
2.7	0	0	0
2.8	0	0	0
2.9	0	0	0
2.10	0	0	0
2.11	0	0	0
2.12	0	0	0
2.13	0	0	0
2.14	0	0	0
2.15	0	0	0
2.16	0	0	0
2.17	0	0	0
2.18	0	0	0
2.19	0	0	0
2.20	0	0	0
2.21	0	0	0
2.22	0	0	0
2.23	0	0	0
2.24	0	0	0

**3.0 Leads Test.**

	Pass	Fail	N/A
3.1	0	0	0
3.2	0	0	0
3.3	0	0	0
3.4	0	0	0

**4.0 Power Supply Test**

	Pass	Fail	N/A
4.1	0	0	0
4.3	0	0	0
4.4	0	0	0
4.5	0	0	0
4.6	0	0	0
4.8	0	0	0
4.9	0	0	0
4.10	0	0	0

**5.0 Leakage Current Test**

	Pass	Fail	N/A
5.1	0	0	0

**6.0 Paddles Test**

	Pass	Fail	N/A
6.1	0	0	0
6.2	0	0	0
6.3	0	0	0
6.4	0	0	0
6.5	0	0	0
6.6	0	0	0
6.7	0	0	0
6.8	0	0	0

**7.0 Heart Rate Display Test**

	Pass	Fail	N/A
7.1	0	0	0

**8.0 Calibrating Pulses Test**

	Pass	Fail	N/A
8.2	0	0	0

**9.0 Notch Filter Test**

	Pass	Fail	N/A
9.3	0	0	0

**10.0 Heart Rate Alarm Test**

	Pass	Fail	N/A
10.1	0	0	0
10.2	0	0	0
10.3	0	0	0
10.4	0	0	0
10.5	0	0	0
10.6	0	0	0
10.7	0	0	0
10.8	0	0	0
10.9	0	0	0
10.10	0	0	0
10.11	0	0	0
10.12	0	0	0
10.13	0	0	0
10.14	0	0	0

**11.0 Defibrillator Self Test**

	Pass	Fail	N/A
11.1	0	0	0
11.2	0	0	0
11.3	0	0	0
11.4	0	0	0
11.5	0	0	0
11.6	0	0	0
11.7	0	0	0

**12.0 Synchronized Cardioversion Test**

	Pass	Fail	N/A
12.2	0	0	0
12.3	0	0	0

# ZOLL M Series Maintenance Test Checklist

Serial No. _____	Location _____	Date _____
Tester _____	Signature _____	

### 13.0 Shock Test

	Pass	Fail	N/A
13.2	0	0	0
13.3	0	0	0
13.4	0	0	0
13.5	0	0	0
13.6	0	0	0
13.7	0	0	0
13.8	0	0	0
13.9	0	0	0
13.10	0	0	0
13.11	0	0	0
13.12	0	0	0
13.13	0	0	0
13.14	0	0	0
13.15	0	0	0
13.16	0	0	0
13.17	0	0	0
13.18	0	0	0
13.19	0	0	0
13.20	0	0	0
13.21	0	0	0
13.22	0	0	0
13.23	0	0	0
13.24	0	0	0

### 14.0 Summary Report Test

	Pass	Fail	N/A
14.1	0	0	0
14.2	0	0	0
14.3	0	0	0
14.4	0	0	0

### 15.0 Advisory Message Test

	Pass	Fail	N/A
15.2	0	0	0
15.3	0	0	0
15.4	0	0	0

### 16.0 Pacer Test

	Pass	Fail	N/A
16.1	0	0	0
16.2	0	0	0
16.3	0	0	0
16.4	0	0	0
16.5	0	0	0
16.6	0	0	0
16.7	0	0	0
16.8	0	0	0
16.9	0	0	0
16.10	0	0	0
16.11	0	0	0
16.12	0	0	0
16.13	0	0	0

### 17.0 SpO<sub>2</sub> Monitor Test

	Pass	Fail	N/A
17.1	0	0	0
17.2	0	0	0
17.3	0	0	0
17.4	0	0	0
17.5a	0	0	0
17.5b	0	0	0
17.6a	0	0	0
17.6b	0	0	0
17.7	0	0	0
17.8	0	0	0
17.9	0	0	0

### 18.0 EtCO<sub>2</sub> Monitor Test

	Pass	Fail	N/A
18.1	0	0	0
18.2	0	0	0
18.3	0	0	0
18.4	0	0	0
18.5	0	0	0
18.6	0	0	0

### 19.0 Temperature Test

	Pass	Fail	N/A
19.3	0	0	0
19.4	0	0	0
19.5	0	0	0
19.6	0	0	0

### 20.0 IBP Monitoring Test

	Pass	Fail	N/A
20.1	0	0	0
20.2	0	0	0
20.3	0	0	0
20.4	0	0	0
20.5	0	0	0
20.6	0	0	0
20.7	0	0	0
20.8a	0	0	0
20.8b	0	0	0
20.10a	0	0	0
20.10b	0	0	0
20.11	0	0	0
20.13	0	0	0
20.14	0	0	0
20.15	0	0	0
20.16	0	0	0
20.17	0	0	0
20.18	0	0	0
20.19	0	0	0
20.20	0	0	0
20.21	0	0	0
20.22	0	0	0
20.23	0	0	0
20.24	0	0	0
20.25	0	0	0
20.27	0	0	0
20.29	0	0	0
20.30	0	0	0

20.32	0	0	0
20.33	0	0	0
20.35	0	0	0
20.36	0	0	0
20.38	0	0	0
20.39	0	0	0

### 21.0 NIBP Transducer Calibration

	Pass	Fail	N/A
21.1	0	0	0
21.2	0	0	0
21.3	0	0	0
21.4	0	0	0
21.5	0	0	0
21.6	0	0	0
21.7	0	0	0
21.8	0	0	0
21.9	0	0	0

### 22.0 NIBP Monitor Test

	Pass	Fail	N/A
22.2	0	0	0
22.3	0	0	0
22.4	0	0	0

### 23.0 NIBP Leak Test

	Pass	Fail	N/A
23.1	0	0	0
23.2	0	0	0
23.3	0	0	0
23.4	0	0	0
23.5	0	0	0
23.6	0	0	0
23.7	0	0	0

### 24.0 Bluetooth Test

	Pass	Fail	N/A
24.1	0	0	0
24.3	0	0	0
24.9a	0	0	0
24.9b	0	0	0
24.9c	0	0	0