

METALLOGRAPHY LABORATORY MANUAL

Er. Tushar Dasgattanaik

Lecturer in Metallurgy

DIPLOMA

(3rd YEAR-5th SEMESTER)

**Department of Metallurgical
Engineering**



ORISSA SCHOOL OF MINING ENGINEERING, KEONJHAR

A Government of Odisha institution with National Repute Established in the Year 1956 (Approved by AICTE, New Delhi & Affiliated to SCTE&VT, Odisha,BBSR)

EXPERIMENT -1

Aim: - To study the complete working operation of Metallurgical Microscope

Objectives: -

1. Familiarization with the different components of metallurgical components.
2. Familiarization with compound optical microscope and metallographic.
3. To study magnification system and how to increase in magnification.

Introduction: -

The metallurgical microscope is the most important tool of the metallurgist. It consists an objective and an eye-piece. Its primary function is to reveal the details of the object. The clarity and the extent to which the details are revealed depends on the degree to which these optical systems are created.

Principle: - A horizontal beam of light from the light source is reflected by means of a plane glass reflector downwards through the microscope objective on the surface of the specimen some of these incident light reflected from the specimen surface will be magnified and passing through the plane glass reflector and magnified again by upper lens system of the eye-piece.

Constructional Details: - The table type microscopes are consisting of

1. Stage: - A flat movable table supporting specimen. This can be moved up or down by knobs.
2. Tubes: - The vertically movable tube containing eyepiece, objective and plane reflector. The tube length varies from 160 mm to 250 mm.
3. Rough & fine focus Adjustments - The limbs of microscope carry the coarse & fine adjustments.
4. Objective – The body tube carries revolving nose piece carrying the three objectives. This enables quick change of the objective which helps for a quick resolving the structure of metal, the magnification of lenses is enlarged on focal length of the lens used

The important properties of an objective are-

- 1) Magnifying Power
- 2) Resolving Power.

It is the property by which an objective shows distinctly represented two small adjacent bonds in the structure of the object. This is usually expressed as number of lines per mm that can be separated which depends on the numerical operator, the wavelength of the light used. Resolution is particularly important during the microscopy of the micro constituents of metals consisting of fine lamination with core resolution which appears as one uniform area, where as an objective with higher numerical appearance reveals deeper nature of the structure.

Total magnification of microscope may be calculated as

$$M = L * E / F$$

Where, L- The distance from back of objective to eyepiece.

E- Magnification of Eye piece. & F- The focal length of objective.

Eyepiece - It is named, as it is near to the eye. It is made up of various Powers such as _ 5, _ 10, _ 15 etc.

Uses: -The metallurgical microscope is useful in quality control department in Industries to observe & study

- 1) Differential phases
- 2) Porosity or defects.

All these have a great effect on mechanical properties of material.

List of Modern Microscopes –

- i) Watson Royal Microscope.
- ii) Van Lanes Hock Microscope.
- iii) Glass led Microscope.
- iv) Baker series Microscope.
- v) Leitr Microscope.

Conclusion: -

Proper selection of magnification is very much essential because –

- a. To travel inclusions in steel low magnification is used.
- b. To travel grain size, grain structure, twin boundaries etc. medium magnification is used.
- c. To travel particular e.g. coarse perlite, fine martensite, bainite, etc. high magnification is necessary.

Therefore, proper study of metallurgical microscope is necessary.



Fig- Metallurgical Microscope

EXPERIMENT-2

Aim- Preparation of metallic specimen for metallographic study by grinding, polishing and etching.

Objectives-

1. Familiarization with the procedure for preparation of material specimen for microscopic examination.
2. Familiarization with the compound optical microscopic and metallographic.
3. Examination of surface characteristics of engineering materials.
4. Grain size determination of metals.
5. To study structural characteristics or constitution of metal or alloy in relation to its physical and mechanical properties.

Introduction-

There are two examination methods-

1. Macroscopy
2. Microscopy

In macroscopy the examination of structural characteristic or chemical characteristic of metal or an alloy is done by unaided eye or with aid of low power microscope or binocular usually under 10X.

In microscopy similar examination is done with prepared metal specimens, employing magnifications with optical microscope of from 100X to as high as 200X. Apart from observation of micro structural details in a metal or alloy, other defects such as grain boundaries, twins, precipitates can be observed readily in microscopic examination.

The preparation of metallurgical specimen generally can be divided into series of stages sectioning/sampling, mounting, grinding, polishing & etching.

Procedure of Test

1] Sectioning / Sampling-

The choice of sample for microscopic study may be very important. If a failure is to be investigated the sampling should be chosen as close as possible to the area of the failure & should be compared with one taken from the normal section. If the material is soft, such as nonferrous metals or alloy & non heat treated steels, the section is obtained by manual hack sawing /power saw. If the material is hard, the section may be obtained by use of an abrasive cut off wheels. This wheel is thin disk of suitable cutting abrasive rotating at high speed. The specimen should be kept cool during the cutting operation.

2] Rough Grinding –

Whenever possible the specimen should be of a size & shape that is convenient to handle. A soft sample may be made flat by slowly moving it up to & back across the surface of a flat smooth file. The soft hard may be rough ground on a belt sander with specimen kept cool by frequent dipping in water during the grinding operation. In all grinding & polishing operation, the specimen should be moved perpendicular to the existing scratches this will facilitate, recognition of stage when the deeper scratches are replaced by shallower one characteristic of the finer abrasives. The rough grinding is continued until the surface is flat & free from wire brushes & all scratches due to hacksaw or cut-off wheel are no longer visible.

3] Intermediate Polishing –

After the previous processes the specimen is polishing on a series of emery paper containing successively finer abrasive (Si-C). The first paper is usually no. 1 than 1/0, 2/0, 3/0, & finally 4/0. The intermediate polishing operation using emery paper is usually done dry. However, in certain case such as preparation of soft material, Silicon Carbide has greater removal rate & as it is resin bonded, can be used with a lubricant, which prevents overshooting of the sample, minimizes shearing of soft metals & also provides a rising action to flush away surface removal product so the paper won't be clogged.

4] Fine polishing –

The time consumed & the success of fine polishing depends largely on the case that we exercised during the previous polishing processes. The final approximation to the flat, scratch free surface is obtained by the use of a wet rotating wheel covered with a special cloth that is charged by carefully sized abrasive particles. A wide range of abrasive is available for final polishing, while many will do a satisfying job, these appear to be presence of gamma form of aluminium-oxides (Al_2O_3), for ferrous & copper based materials & Cerium oxide for Aluminium, Magnesium & their alloys, other final polishing abrasives often used are diamond, chromium oxide & magnesium oxide etc. A choice of proper polishing cloth depends upon the particular material being polished & the purpose of metallographic study. Many cloths are available of varying lap or pile, from those having no pile, such as silk, to those of intermediate pile such as broad cloth, billiard cloth, canvas cloth & finally to a deep pile such as velvet synthetic clothes are also available for general purpose of which two under the trival names of gamal & micro cloth are most widely used.

5] Etching –

The purpose of etching is to make the many structural characteristics of the metal or alloy visible. The process should be such that the various parts of the microstructure may be clearly differentiated. This is to subject the polished surface to chemical action. In the alloys composed of two or more shapes. The competent are revealed during etching by a preferential attack of one or more of the constituents by the reagent because of difference in chemical composition of the phases. In uniform single phase alloy contact is obtained and the grain boundaries are made visible because of difference in the rate at which various grains are attacked by the reagent This difference in the rate of attack by reagent which is mainly associated with angle of the different grain structure section to the plane of the polished surface. Because of chemical attack of the chemical reagent the grain boundary appears as valleys in the polished surface light from the microscope hitting the side of these valleys will be reflected but of the microscope making the grain boundaries appears dark lines. The section of the appropriate etching reagent is determined by metal or alloys & the specific structure desired for viewing.

Etching Reagent	Composition	Use
Nitric Acid	White nitric acid '1 to 5 ml ethyl or methyl alcohol' 100 ml(95%	Carbon steels and cast iron
Picric acid	Picric acid '49 ethyl or methyl alcohol'	Carbon steels
Ferric chloride and hydrochloric	Ferric chloride and hydrochloric acid 50 ml water 100 ml	Structure of austenitic nickels stainless steel
Ammonium hydroxide and hydrogen	Ammonium hydroxide 90 parts water 5 parts hydrogen peroxide 2 to 5	Copper and its alloys
Palmerton reagent	Chromic oxide 200 gm sodium sulphate 15 gm water	Zinc and its alloys
Ammonium molybdenum	Molybdic acid 100 gm ammonium hydro 140ml water 240 ml, filters add to nitric acid S.P.gr. 60 ml	Lead and its alloys
Hydrofluoric acid	Hydrofluoric acid 0.5 ml	Aluminum and its

Result and Conclusion: -

Correct sample techniques if followed properly, will only reveal the microstructure of a metallic specimen. Otherwise it will be most difficult to distinguish between different phase in the microstructure, the fine scratches on meta surface and grain boundaries etc.

EXPERIMENT-3


1.Aim-: Preparation and study of the Micro Structure of pure metals Mild Steel, Low Carbon steel and High Carbon Steel.

2.Objectives:

- a.) To Learn the preparation of specimen for microscopic observation.
- b) To understand what Microscopy is, and how it can be used to observe Microstructure of metals.

3. APPARATUS / EQUIPMENT

3.1 Abrasive Cat-Off Wheel Machine




- 3 HP Motor, 1440 rpm
- Manual model
- 10” Diameter-wheel
- Fully enclosed
- Cutting Capacity 50mm Diameter & 75mm
- Height X 200mm Length.


This Equipment is used to cut the sample pieces from the component of all kind of materials. These sample pieces can be used for microscopic observation.

3.2 Specimen Mounting Press

- Manual Operation
- 1” Mould Capacity
- Digital Temperature/ Time control unit
- Pressure gauge for reading pressure
- Alarm indicator for cycle completion.




This Equipment is used to mount specimen in bakelite. Particularly is useful in hot mounting process.



3.3 Belt Finishing Machine.


This Equipment is used for rough Grinding of the Specimen 3.4 Polishing Machine Machine.



- Variable speeds 100rpm - 800rpm
- High torque PMDC motor with D. C. Drive.
- Water In and Drain
- Corrosion resistant fiber reinforced cover.
- 8" Disc

This Equipment is used to polish the specimen after rough grinding with the help of polishing clothe. It also useful for rough grinding with grinding disc.

3.5 Metallurgical Microscope.



Binocular model

Fig 5 Fe-Fe₃C Phase Diagram

Figure 5 shows the equilibrium diagram for combinations of carbon in a solid solution of iron. The diagram shows iron and carbons combined to form Fe-Fe₃C at the 6.67%C end of the diagram. The left side of the diagram is pure iron combined with carbon, resulting in steel alloys. Three significant regions can be made relative to the steel portion of the diagram. They are the

eutectoid E, the hypoeutectoid A, and the hypereutectoid B. The right side of the pure iron line is carbon in combination with various forms of iron called alpha iron (ferrite), gamma iron (austenite), and delta iron. The black dots mark clickable sections of the diagram.

Allotropic changes take place when there is a change in crystal lattice structure. From 2802°-2552°F the delta iron has a body-centered cubic lattice structure. At 2552°F, the lattice changes from a body-centered cubic to a face-centered cubic lattice type. At 1400°F, the curve shows a plateau but this does not signify an allotropic change. It is called the Curie temperature, where the metal changes its magnetic properties. Two very important phase changes take place at 0.83%C and at 4.3% C. At 0.83%C, the transformation is eutectoid, called pearlite.

gamma (austenite) —> alpha + Fe₃C (cementite)

At 4.3% C and 2066°F, the transformation is eutectic, called ledeburite.

L(liquid) --> gamma (austenite) + + Fe₃C (cementite)

Definitions

Eutectoid: A eutectoid system occurs when a single-phase solid transforms directly to a two- phase solid.

Hypereutectoid: Hypereutectoid systems exist below the eutectoid temperature.

Hypoeutectoid: Hypoeutectoid systems exist above the eutectoid temperature.

Ferrite: Body-centered cubic iron or an iron alloy based on this structure.

Austenite: Face-centered cubic iron or an iron alloy based on this structure.

Delta iron: The body-centered cubic phase which results when austenite is no longer the most Stable form of iron. Exists between 2802 and 2552 degrees F. has BCC lattice structure and is magnetic.

Body-centered: A structure in which every atom is surrounded by eight adjacent atoms, whether the atom is located at a corner or at the center of a unit cell.

Face-centered: A structure in which there is an atom at the corner of each unit cell and one in the center of each face, but no atom in the center of the cube.

Pearlite: A lamellar mixture of ferrite and carbide formed by decomposing austenite of Eutectoid composition.

Cementite: The second phase formed when carbon is in excess of the solubility limit.

Ledeburite: Eutectic of cast iron. It exists when the carbon content is greater than 2 percent. It contains 4.3 percent carbon in combination with iron.

4.2 Classification of Steel: The steels are classified by various methods and each method is based on a definite criterion as follows,

- i) Amount of carbon
 - a) Low carbon steels (0.008 - 0.3% C)
 - b) Medium carbon steels (0.30 - 0.60% C)
 - c) High carbon steels (0.60 - 2.00% C) ii) Amount of alloying elements and carbon
- iii) Amount of deoxidation
- iv) Method of manufacture
- v) Form and use.

5. PROCEDURE:

The preparation of metallic or other materials for microscopic examination and micro structural characterization is in principal very simple. There are four basic processes that you will need to become familiar with: sample cutting and sectioning, metallographic mounting, Surface grinding and surface polishing.

5.1 Sample Cutting and Sectioning

5.1.1. Sectioning

Sectioning means removal of convenient size specimen from large sample with minimal damage to microstructure with the help of abrasive cut off machine. Abrasive cutting wheel/saw is attached to cutting machine and for work piece holding proper vice is provided on machine. The primary concern in this process is to minimize the heating of the sample due to the cutting. For this reason, the cut-off saws that is equipped with either water-cooling system.

5.2 Mounting

If sample is large enough (about 25 mm square or larger) than do not need to mount it, as it will be able to control the sample during polishing without a mount. For smaller samples there are two basic mounting techniques used in this laboratory. Mounting facilitates handling during preparation and handling. It also avoids damage to polishing wheels during polishing. The most common uses a thermosetting plastic compound (Bakelite) to encapsulate the specimen known as hot-mounting process, and the second uses a room temperature curing epoxy known as cold- mounting process. The Bakelite mounting is by far the most common and easiest. The room temperature curing epoxy mount should only be used for samples that are extremely sensitive to heat. The Bakelite process uses a sample mounting press that applies a pressure to the Bakelite/sample system during the cure to remove voids and gaps and to fully fill the sample spaces. Bakelite comes in a variety of colours, which can be combined to produce easy sample identification.

5.2.2 The procedure for Hot mounting process

The following process is used to encapsulate your specimen in Bakelite:

- 1) Place sample face down on the small piston inside the press, and lower the piston into the cylinder by opening slightly the valve on the front of the press.
- 2) Approximately three tablespoons of Bakelite is poured over the sample, and the top of the press gently screwed into place. **DO NOT TIGHTEN THE TOP OF THE PRESS.** It is only necessary to engage all of the screw threads; you do not have to tightly secure the top.
- 3) The cylindrical heater is plugged in and turned on (the red light should turn on). Place the heater around the mold. The heater is thermostatically controlled and will heat the mold to about 135-150 °C. Close the valve, and pump up the cylinder using the hand lever. As the Bakelite heats, it will begin to flow to fill the void spaces, and the pressure will drop. Maintain the constant.
- 4) When the pressure stops rapidly dropping, the whole mold has reached 150°C. Begin timing for 5-7 minutes to fully cure the Bakelite. Maintain the pressure during the entire heating and cooling cycle.
- 5) At the end of the heating cycle, remove the heater and place the cooling collar on the mold for an additional 6-8 minutes.

6) Crack the valve to release the pressure, and unscrew the top of the mold. When the mold top is fully unscrewed (it may not come out due to adhesion with the Bakelite), close the valve and slowly pump up the cylinder to push the sample fully out of the press. Mark the sample on the back

7) Clean any residual Bakelite off of the mold surfaces.

5.2.2 The procedure for Cold Mounting process.

The following process is used to room temperature curing epoxy process:

- 1) Apply mold release agent to mold. Place specimen in the mold.
- 2) Mix epoxy powder and bonding liquid in 1:2 ration in a cup.
- 3) Pour into mo Id... wait for 10 minutes.
- 4) Eject the mold.

5.3 Sample Surface Polishing

The goal of the surface polishing is to end up with a planar cross section of sample free from scratches or disturbed metal introduced by the cutting and sectioning. This process is a step-wise process that can be broken into three loosely separate parts: grinding, coarse polishing, and final polishing.

5.3.1. Grinding

The first step in preparing your sample is to ensure that you have a flat surface to begin with. A water-cooled abrasive grinder is available to form a flat initial surface from which to begin. After getting a flat sample on the belt grinder, WASH sample thoroughly. The hand lapping station has four graded abrasive papers to produce a sequentially finer surface finish. Be sure the water is turned on and flowing uniformly over the abrasives. Start with the coarsest grit (240) and, using a firm and uniform pressure, slowly move the specimen forward and back across the abrasive. This will produce parallel scratches of uniform size. Continue this step until the entire surface of your sample is flat and contains only scratches of the size of 240 grit abrasive. When the sample is flat and the only scratches remaining are those due to the 240 grit abrasive, WASH your sample and your hands thoroughly, and move to the 320 grit abrasive. Repeat this procedure for the 400 grit and the 600 grit abrasive, checking after each step to be sure that only those scratches remain that are due to the smallest grit.

5.3.2 Rough Grinding

Before proceeding to the first polishing wheel (leftmost wheel), wash sample with water.

- 1) First, apply a small amount of water to the wheel, turn on the motor, and gently clean off the wheel with your fingers.
- 2) Apply a small amount of abrasive slurry to the wheel. This wheel uses an Al₂O₃ abrasive in a water suspension. The abrasive particles are 5 micrometres in diameter.
- 3) Carefully place your sample on the wheel while gripping it tightly. Slowly move the sample in a circular motion against the rotation of the wheel. Use a moderate and even pressure. It is important to ensure that you keep the sample flat on the wheel so that the final surface will be completely planar.
- 4) After several minutes on the wheel, hold the sample in one place for a moment. This will provide lots of parallel scratches that you can use to determine if you have removed the damage from the grinding steps.
- 5) Examine the sample under the microscope to determine if all the scratches are the same size.
- 6) Repeat steps 1-5 on the middle polishing wheel. This wheel uses a 0.3 micrometre Al₂O₃ abrasive in a water suspension.

5.3.3. Final Polishing

- 1) Repeat steps 1-5 above on the right polishing wheel. This wheel uses a 0.05 micrometre Al₂O₃ abrasive in a water suspension. At this point, the sample will be very smooth to the eye and even the oils and dirt on your fingers will scratch it with larger scratches than the abrasive. **DO NOT TOUCH THE SAMPLE SURFACE FROM THIS POINT ON.**
- 2) The last step in the process is to etch the sample to bring out the microstructure.
- 3) Use a cotton swab and a petri dish for the etching. Gently swab the surface of your sample with the etchant. Roughly spreading the etchant will scratch your surface. Let the etchant stand for 15 seconds or so and rinse the sample with water to stop the etching, and rinse again with methanol. Rinse the swab with water and throw into the trash bin.
- 4) Examine specimen under the microscope. You may require several etching steps to bring out the microstructure.

5) If the sample is over-etched, repeat the final polishing step and re-etch for a shorter time.

Samples to be examined at high magnification generally require shorter etching times than those to be viewed at lower magnifications.

After last polishing stage the sample looks mirror like.

5.4 Etching

Grains cannot be seen without etching. Cracks, pores and defects are observed without etching. Etchant reacts with atoms and dissolves them. Atoms at grain boundaries dissolve quickly. Dissolved grain boundaries appear dark.

Steps:

- 1) Apply etchant to polished surface for some time
- 2) Rinse with distilled water

Etchants:

- 1) 50/50 HCl: equal parts hydrochloric acid (HCl) mixed with water.
- 2) Alcoholic Ferric Chloride: 5 grams FeCl_3 ; 2 ml concentrated HCl acid; 95 ml methyl alcohol.
- 3) Aqueous Ferric Chloride: 10 grams FeCl_3 ; 20 ml concentrated HCl acid; 80 ml water.
- 4) Ammonia/Hydrogen peroxide: 1-part strong ammonia; 1-part hydrogen peroxide; 2 parts water; FRESHLY MADE.
- 5) Mixed acids: 95 ml water; 1.5 ml concentrated HCl acid; 2.5 ml concentrated nitric (HNO_3) acid; 0.5 ml hydrofluoric (HF) acid.
- 6) 2% Nital: 2 ml concentrated HNO_3 ; 98 ml methyl alcohol.

5.5. Metallographic Observation

Observe microstructure, Place specimen on metallography and adjust magnification, focus and position s adjust micro High magnification - to study phases and Low magnification -to study grain size.

5.5.1 Microphotography

In this laboratory, you will report the microstructures of prepared samples in specific formats. You will be expected to sketch the microstructure that you see under the

microscope by hand. In sketching the microstructure there are several things to keep in mind. First, the magnification that you use depends upon the scale of the microstructure you are looking for. It is **IMPORTANT** to know in advance of the lab class what the expected microstructure for your samples are and at what scale they should appear. In sketching the microstructure, you should indicate only the important features of the structure that you observe-don't make a photographic reproduction of the microstructure. Simple sketches show that you know what the important structures are and have identified them in the cross section. An example of what is considered to be a good laboratory report sketch of the microstructures is included in the appendix. Your sketch **MUST INCLUDE:**

- 1) The sample name and composition,
- 2) The metallurgical history of the sample,
- 3) A simple sketch of the important microstructure indicating a) The magnification used (i.e. "MAG= 100 X") Important phases and features noted
- 4) Etchant used.

EXPERIMENT-4

1.Aim: Study of the Microstructures of Cast Iron.

2.Objectives:

- a.) To understand what Microscopy is, and how it can be used to observe Microstructure of metals like cast iron.
- b) To conduct typical engineering microscopic observation and be able to recognize the microstructure of typical type of cast iron metal.
- c) To learn the differences in microstructure of different type of cast irons.

3. APPERATUS:

Metallurgical Microscope, Standard Specimens of Cast Iron materials.

3.1 Metallurgical Microscope.



Binocular model

- Eye pieces: 10X,15X & 20X
- Magnification: 10X-500X.

This Microscope is used for -microscopic observation of specimen with the adjustable magnification.


4. THEORY:

4.1 Cast Iron

Cast irons typically contain 2-4 wt.% of carbon with a high silicon concentrations and a greater concentration of impurities than steels. The carbon equivalent (CE) of a cast iron helps to distinguish the grey irons which cool into a microstructure containing graphite and the white irons where the carbon is present mainly as cementite. The term cast iron, like the term steel, identifies a large family of ferrous alloys. Cast irons are multi component ferrous alloys. They contain major (iron, carbon, silicon), minor ($<0.01\%$), and often alloying ($>0.01\%$) elements. Cast iron has higher carbon and silicon contents than steel. Because of the higher carbon content, the structure of cast iron, as opposed to that of steel, exhibits a rich carbon phase. Depending primarily on composition, cooling rate and melt treatment, cast iron can solidify according to the thermodynamically metastable Fe-Fe₃C system or the stable Fe- Gr system. When the metastable path is followed, the rich carbon phase in the eutectic is the iron carbide; when the stable solidification path is followed, the rich carbon phase is graphite. Referring only to the binary Fe-Fe₃C or Fe-Gr system, cast iron can be defined as an iron-carbon alloy with more than 2% C. Important notice is that silicon and other alloying elements may considerably change the maximum solubility of carbon in austenite (g). Therefore, in exceptional cases, alloys with less than 2% C can solidify with a eutectic structure and therefore still belong to the family of cast iron.

4.2 Types of Cast Iron

On the basis of Microstructure cast irons are classified as follows.



4.2.1 White Cast iron

Composition of the iron is appropriate or the cooling rate of the metal is sufficiently rapid during solidification, the metal will solidify with the C combined with iron as iron carbide. This compound, also called cementite, is hard and brittle and dominates the microstructure of white iron. Thus, white iron is hard and brittle and has a white crystalline fracture because it is essentially free of graphite.

White iron has a high compressive strength and excellent wear resistance, and it retains its hardness for limited periods even up to a red heat. It can be produced in selected areas of a casting such as. on the periphery of a cam—by causing localized rapid solidification of the iron. White iron at the surface of a casting is called chill. It is produced by making that portion of the mold where the white iron is desired—of a material that can extract heat very rapidly, such as iron or graphite. White iron does not have the easy castability of other irons because its solidification temperature is generally higher, and it solidifies with C in its combined form as iron carbide. Application includes rollers of rolling mills, dies of metal extrusion and where high wear resistance is necessary.





Fig. 2 Microstructure of white cast iron Fig. 3. Microstructure of Gray cast iron

4.2.1 Gray Cast Iron

When the composition of the molten iron and its cooling rate are appropriate, the C in the iron separates during solidification and forms separate graphite flakes that are interconnected within each eutectic cell. The graphite grows edgewise into the liquid and forms the characteristic flake shape. When gray iron is broken, most of the fracture occurs along the graphite, thereby accounting for the characteristic gray color of the fractured surface. Because the large majority of the iron castings produced are of gray iron, the generic term, cast iron, is often improperly used to mean gray iron specifically. The properties of gray iron are influenced by the size, amount and

distribution of the graphite flakes, and by the relative hardness of the matrix metal around the graphite. These factors are controlled mainly by the C and Si contents of the metal and the cooling rate of the casting. Slower cooling and higher C and Si contents tend to produce more and larger graphite flakes, a softer matrix structure and lower strength. The flake graphite provides gray iron with unique properties such as excellent machinability at hardness levels that produce superior wear-resisting characteristics, the ability to resist galling and excellent vibration damping. The amount of graphite present, as well as its size and distribution, are important to the properties of the iron. Whenever possible, it is preferable to specify the desired properties rather than the factors that influence them.

4.2.2 Nodular cast Iron (Ductile iron, S. G Iron)

Ductile iron also referred to as nodular iron or spheroidal graphite iron, was patented in 1948. After a decade of intensive development work in the 1950s, ductile iron had a phenomenal nine- fold increase in use as an engineering material during the 1960s, and the rapid increase in commercial application continues today. An unusual combination of properties is obtained in ductile iron because the graphite occurs as spheroids rather than as individual flakes as in gray iron. This mode of solidification is obtained by adding a very small, but specific, amount of Mg to molten iron of a proper composition. The base iron is severely restricted in the allowable contents of certain minor elements that can interfere with the graphite spheroid formation. The added Mg reacts with the sulfur and oxygen in the molten iron and changes the way the graphite is formed. Control procedures have been developed to make the processing of ductile iron dependable. The high C and Si content of ductile iron provide the casting process advantages, but the graphite spheroids have only a nominal influence on the mechanical properties of the metal. Ductile iron, like malleable iron, exhibits a linear stress-strain relation, a considerable range of yield strengths and, as its name implies, ductility. Castings are made in a wide range of sizes with sections that can be either very thin or very thick. The different grades are produced by controlling the matrix structure around the graphite either as- cast or by subsequent heat treatment. Only minor compositional differences exist among the regular grades, and these adjustments are made to promote the desired matrix microstructures. Alloy additions may be made to ductile iron to assist in controlling the matrix structure as-cast or to provide response to heat treatment. Special analysis ductile irons and high-alloy ductile irons provide unusual properties for special applications.

4.2.3 Malleable Cast Iron

This type of iron is characterized by having the majority of its C content occur in the

microstructure as irregularly shaped nodules of graphite. This form of graphite is called temper carbon because it is formed in the solid state during heat treatment. The iron is cast as a white iron of a suitable chemical composition. After the castings are removed from the mold, they are given an extended heat treatment starting at a temperature above 1650°F (900°C). This causes the iron carbide to dissociate and the free carbon precipitates in the solid iron as graphite. The rapid solidification rate that is necessary to form the white iron limits the metal thickness in the casting that is practical for the malleable iron process. A wide range of mechanical properties can be obtained in malleable iron by controlling the matrix structure around the graphite. Pearlitic and martensitic matrices are obtained both by rapid cooling through the critical temperature and with alloy additions. Malleable irons containing some combined carbon in the matrix often are referred to as pearlitic malleable, although the microstructure may be martensitic or a spheroidized pearlite.




Fig. 4 Microstructure Nodular Cast Iron




Fig. 5. Microstructure Malleable Cast Iron




Fig. 6 Iron-Carbon diagram (cast iron portion)

4. PROCEDURE:

In this laboratory, you will report the microstructures of prepared samples in specific formats. You will be expected to sketch the microstructure that you see under the microscope by hand. In sketching the microstructure there are several things to keep in mind. First, the magnification that you use depends upon the scale of the microstructure you are looking for. It is IMPORTANT to know in advance of the lab class what the expected microstructure for your samples are and at what scale they should appear. In sketching the microstructure, you should indicate only the important features of the structure that you observe-don't make a photographic reproduction of the' microstructure. Simple sketches show that you know what the important structures are and have identified them in the cross section.

EXPERIMENT-5

1.AIM: Study of the Microstructures of Non Ferrous Metals.

2. Objectives:

- a.) To understand what Microscopy is, and how it can be used to observe Microstructure of non-ferrous metals.
- b) To conduct typical engineering microscopic observation and be able to recognize the microstructure of typical type of copper, aluminium and its alloys.
- c) To learn the differences in microstructure of different type of cast irons.

3. APPERATUS:

Metallurgical Microscope, Specimens of Non Ferrous Metals.

3.1 Metallurgical Microscope.



Binocular model

- Eye pieces: 10X, 15X & 20X
- Magnification: 10X - 500X. This Microscope is used for microscopic observation of specimen with the adjustable magnification.

4.THEORY:

4.1 COPPER AND COPPER ALLOYS

4.1.1 COPPER -: Copper is a face centered cubic (FCC metal) with very high ductility, electrical and thermal conductivity. This is very soft too. It is used in electrical applications, industrial machinery, hearth for arc furnaces. In the alloyed form, it also finds enormous importance. One such application is duralumin (Al-4.5% Cu), which is common material for airplane body material. When it is alloyed with Zn, it is called brass. 70:30 brass is very common in making pipe lines for water.

Copper properties: -

- Its melting pt is 1083C
- Specific gravity is 8.9.
- Has high electrical and thermal conductivity.
- Its good ductility and malleability r due to its FCC structure.
- Cu containing 0.3% As is called as ARSENICAL Cu
- Cu containing 0.6% Al is called as FREE CUTTING Cu.

4.1.2 COPPER ALLOYS.

There are two types of Cu alloys popularly known: - BRASS AND BRONZE.

BRASS: - $\text{Cu} + \text{Zn} = \text{Brass}$.

- Brass are again divided as alpha brass and alpha + beta brass.

// TYPES OF ALPHA BRASS //

1) CAP BRASS: -

- Consists of 2-5% Zn'
- Very ductile alloy.
- Zn is used as deoxidiser for deoxidation of Cu.
- Used for cap of detonation.

2) GUILDING METAL: -

- Zn-(5-15%).
- Has colour of gold.
- Uses-artificial jewellery, condenser tubes, coils, needles, etc.

3) CARTRIDGE BRASS: -

- Known as (70-30) brass
- Has max ductility and malleability of all brasses.
- Applications=radiator fins, lamp fixtures, rivets, springs, etc.

4) ADMIRALITY BRASS: -

- 1% Sn is added to cartridge brass to improve corrosion resistance.
- Uses in condenser tubes and heat exchangers.
- Containing 22%Zn,2%Al,0.04%Sn is used in marine applications.

// (ALPHA +BETA) BRASS//: -

- Contains up to (32-40) %Zn.
- Hard, strong and less corrosion resistance compared to alpha brasses.

1) MUNTZ METAL: -

- Contains 40%Zn.
- Alloy becomes single phase at about 700C
- Also called as (60-40) brass.
- Uses in utensils, shafts, nuts, bolts, condenser tubes, etc.

2) NAVAL BRASS: -

- ADDING 1%Sn to muntz metal inc corrosion resistance to marine environment
- Also called as TOBIN BRONZE.

3) LEADED BRASS: -.

- Lead added in small amount (1-3) % to improve machinability.
- Appears as globules in microstructure.

// BRONZE // (Cu-Sn alloy)

1) AL BRONZE: -

- An alloy of Cu and Ni.
- Max solubility of Cu in Ni is 9.4%
- Has good ductility, strength, toughness, corrosion resistance and fatigue resistance.
- Called as IMITATION GOLD.
- Applications “-components of valves, pump castings, park plug body.

2) TIN BRONZE: -

- ALLOY OF Cu and Sn.
- Tin has affinity towards oxygen. the tin oxide formed reduces ductility and malleability. Zn or P are added to deoxidize these.

COINAGE BRONZE

Contains 5% Sn and 1% Zn.


- It is a ductile metal, having better formability and strength.
- Application=manufacturing of coins.

1) PHOSPHOROUS BRONZE: -

- Alloy of Cu and Zn with P.
- Phosphorous is a strong deoxidizer and helps in inc fluidity.
- Phosphorous content is max upto 1% depending on application,
- (2.5-8) % Zn, (0.1 -0.35) % P known as WROUGHT PHOSPHOROUS BRONZE. • (5-13) % Zn and (0.3-1) % P is known as CAST PHOSPHOROUS BRONZE.

2) MONEL METAL: -

- ALLOY OF Cu and Ni.
- Cupronickel metal alloy consists of 15-30%Ni.
- Cu and Ni are completely soluble into each other in both solid in liquid solution.
- Application—blades of turbine, coins, bullet envelopes.



Phase diagram of Cu-Zn system

ALUMINIUM Properties:

- Light in wt.
- Cast ability and formability are better.
- Corrosion resistance of Al is excellent.
- Powerful deoxidiser.
- Carries more electricity than Cu.
- Costlier metal in the family of light metal alloys.
- Has good electrical and thermal conductivity.

Al ALLOYS

1) Al-Si-Cu: -

- LM2, LM6, LM8, LM13.

- Used for production of castings, due to their excellent fluidity and casting characteristics.

- Higher Si content better are the mechanical properties.

2) Al-Mg (MAGNALIUM ALLOY):-

- LM5NDLM10.

- Used in marine environments. Due to its high strength and resistant to corrosion.

- Has good surface finish.

3) Al-Cu (DURALIUM): -

- LM11 (Cu-4.5%)

- Susceptible to hot tearing. Applications are casted components used in aircrafts.

4) Y ALLOY (HIGH STRENGTH Al ALLOY):-

- LM14.

- CONTAINS 4% Cu, 2% Ni, and 1.5% Mg.

- Application=heavy-duty petrol engine, piston block, etc.

5) RR350 (HINDUMINIUM):-

- Contains 5% Cu, 1.5% Ni with small amounts of Mn, Ti, Sb, Co, Zr.

- Used in aero engines and other continuous elevated temp service applications up to 300°C.

4. PROCEDURE:

In this laboratory, you will report the microstructures of prepared samples in specific formats. You will be expected to sketch the microstructure that you see under the microscope by hand. In sketching the microstructure there are several things to keep in mind. First, the magnification that you use depends upon the scale of the microstructure you are looking for. It is IMPORTANT to know in advance of the lab class what the expected microstructure for your samples are and at what scale they should appear. In sketching the microstructure, you

should indicate only the important features of the structure that you observe- don't make a photographic reproduction of the microstructure. Simple sketches show that you know what the important structures are and have identified them in the cross section.

An example of what is considered to be a good laboratory report sketch of the microstructures is included in the appendix. Your sketch **MUST INCLUDE**:

- 1) The sample name and composition,
- 2) The metallurgical history of the sample,
- 3) A simple sketch of the important microstructure indicating
 - a) The magnification used (i.e. "MAG= 100 X"),
 - b) Important phases and features noted,
- 4) Etchant used.

EXPERIMENT-6


1. AIM: Study of the Microstructures of Heat Treated Steels.

2. OBJECTIVE:

- a. To observe microstructure of heat traded steels.
- b. To study the changes in macrostructure due to heat Treatment.
- c. To study the one of the heat treatment process

3. APPERATUS:

- I) Specimens of High Carbon Steel subjected
- II) Metallurgical Microscope.



Binocular model

- Eye pieces: 10X, 15X & 20X
- Magnification: 10X-500X.

This Microscope is used for microscopic observation of Specimen with the adjustable magnification

4. THEORY:

4.1 Heat Treatment is the controlled heating and cooling of metals to alter their physical and mechanical properties without changing the product shape. Heat treatment is sometimes done inadvertently due to manufacturing processes that either heat or cool the metal such as welding or forming.

Heat Treatment is often associated with increasing the strength of material, but it can also be used to alter certain manufacturability objectives such as improve machining, improve formability, and restore ductility after a cold working operation. Thus it is a very enabling manufacturing process that can not only help other manufacturing process, but can also improve product performance by increasing strength or other desirable characteristics. Steels are particularly suitable for heat treatment, since they respond well to heat treatment and the commercial use of steels exceeds that of any other material. Steels are heat treated for one of the following reasons: Softening, Hardening, Material Modification

4.2 Heat Treatment Processes:

4.2.1. Annealing

The process consists of heating the steel to above A3 temperature for hypo-eutectoid steels and above A1 temperature for hypereutectoid steels by 30-50°C, holding at this temperature for a definite period and slow cooling to below A1 or to room temperature in the furnace. Due to slow cooling, eutectoid phase transformation occurs very nearly in accordance with conditions represented by Fe-C phase diagram.

Purpose of Annealing -:

- To relieve the internal stresses induced due to cold working, welding.
- To reduce hardness and to increase ductility.
- To refine the grain size.
- To make the material homogeneous in respect of chemical composition.
- To increase machinability
- To increase the uniformity of phase distribution and to make the material isotropic in respect to mechanical properties.

Types of Annealing -: i.

Full Annealing

ii.Spheroidized Annealing

iii.Bright Annealing

iv. Box annealing

v. Isothermal Annealing

4.2.2. Normalizing

The process consists of heating the steels above the upper critical temperature (A_3) for hypoeutectoid steels, and above A_{cm} for hypoeutectoid steels by 30 to 50°C, holding long enough at this temperature for homogeneous austenitization and cooling to room temperature in air. Due to air cooling which is slightly fast as compared to furnace cooling employed to annealing, normalized components show slightly different structure and properties than annealed components. After normalizing, microstructure shows more pearlite than observed in annealed components.

Purpose -:

- This process is used to eliminate the cementite network which formed due to slow cooling in the temperature range from A_{cm} to A_1 .
- To relieve internal stresses.
- To make the material homogeneous in respect of chemical composition.
- To increase machinability
- To increase the uniformity of phase distribution and to make the material isotropic in respect to mechanical properties.

4.2.3 Hardening:

The hardening process consists of heating the steel to above A₃ temperature for hypoeutectoid steels and above A₁ temperature.

Draw hardening cycle.

4.2.4 Tempering:

Tempering is followed to reduce the effect of retained austenite in which the hardened steel is reheated below the lower critical temperature and cooled slowly.

4.3 Time — Temperature — Transformation Diagram.

Draw the Time — Temperature — Transformation Diagram.

5. PROCEDURE:

In this laboratory, you will report the microstructures of prepared samples in specific formats. You will be expected to sketch the microstructure that you see under the microscope by hand, in sketching the microstructure there are several things to keep in mind. First, the magnification that you use depends upon the scale of the microstructure you are looking for. It is **IMPORTANT** to know in advance of the lab class what the expected microstructure for your samples are and at what scale they should appear. In sketching the microstructure, you should indicate only the important features of the structure that you observe—don't make a photographic reproduction of the microstructure. Simple sketches show that you know what the important structures are and have identified them in the cross section.

An example of what is considered to be a good laboratory report sketch of the microstructures is included in the appendix. Your sketch **MUST INCLUDE**:

- 1) The sample name and composition,
- 2) The metallurgical history of the sample,
- 3) A simple sketch of the important microstructure indicating a) The magnification used (i.e. "MAG= 100 X"),
 - b) Important phases and features noted,
- 4) Etchant used.