ARC InstaTemp MDTM

Non-Touch Thermometer

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Do not cover the sensor aperture wholly or partially, or insert any objects or material into the aperture.

Do not attempt to replace the internal batteries. The internal batteries are non-removable, non-replaceable and non-serviceable.

Do not open the thermometer casing. If the case is not correctly closed, do not use the device.

Avoid Strong magnetic fields.
Avoid X-ray.
Avoid oxygen rich environment.

DISPOSAL


The ARC InstaTemp MD™ contains batteries: do not dispose in fire, or domestic waste. Dispose in accordance with local regulations at local authority collection points

For more information about disposal please contact ARC Devices customer service (see page 35).

ICON IDENTIFICATION


This instruction manual belongs to this device. It contains important information about operating the ARC InstaTemp MD™ non-touch thermometer. Read the instruction manual carefully as non-observance of these instructions can result in inaccurate temperature readings or damage to the device

Manufacturer

Caution; consult section titled 'Understanding the Temperature Range Indicator System' and 'Taking a Temperature'.

Degree of protection against electrical shock; there is no conductive connection to the person; Type BF Protection

ICON IDENTIFICATION

Device to be disposed of in segregated waste (WEEE Directive 2012/19/EU)

Protect from liquids

Device conforms with RoHS directive (2002/95/EC)

The CE Mark of this device refers to the EC Directive MDD 93/42/EEC

TECHNICAL SPECIFICATIONS

<i>Item</i>	<i>Value</i>
Model	ARC InstaTemp MD™
Display System	Digital segmented display (LCD) with backlight
Reported measurement range	89.6 - 109.1°F (32 - 42.8°C)
Operating temperature / humidity range	59 - 95 °F (15 - 35°C) <i>Operational up to 95% humidity non-condensing</i>
Storage/ Transport Temperature range	45°F - 90°F (7°C - 30°C)
Temperature measurement accuracy	+/-0.4°F and +/-0.2°C in the measurement range 89.6 - 109.1°F (32 - 42.8°C) whilst in an ambient temperature range of 59-95°F (15-35°C) (In accordance with ASTM E1965-98 and ISO 80601-2-56 pt 2-56)
Temperature display	Units °F or °C, user selectable

TECHNICAL SPECIFICATIONS

<i>Item</i>	<i>Value</i>
Power supply	Contains two AAAA alkaline batteries, non-replaceable
Battery life	Approximately 13,000 measurements. Low battery warning displays when approximately 500 measurements left.
Thermometer Dimensions	1.37" x 0.55" x 5.08" (3.5cm x 1.4cm x 12.9cm)
Weight	1.94 oz (55 gm)
Wall Mount & Retractable Tether	<p>Wall Mount:</p> <ul style="list-style-type: none"> • 1.24" x 2.125" x 0.42" (31.5mm x 54mm x 10.7mm) <p>Tether:</p> <ul style="list-style-type: none"> • 1.4" x 2.12" x 0.56" (35.3mm x 53.8mm x 14.2mm) • Retractable cord length: 27" (69cm)
Compliance with standards	<p>Electromagnetic compatibility:</p> <p>IEC60601-1-2 Clinical thermometer:</p> <p>ASTME1965-98(Reapproved 2009), ISO 80601-2-56</p>

LIMITED TWO YEAR WARRANTY

ARC Devices warrants the product to be free of defects in material and workmanship and to perform in accordance with manufacturer's specifications for the period of two years from the date of purchase from ARC Devices, its authorized distributors or agents. The warranty period shall start on the date of purchase. The date of purchase is: 1) the invoiced ship date if the device was purchased directly from ARC Devices, 2) the date specified during product registration, 3) the date of purchase of the product from an ARC Devices authorized distributor as documented from a receipt from said distributor. This warranty does not cover damage caused by: 1) handling during shipping, 2) use or maintenance contrary to labeled instructions, 3) alteration or repair by anyone not authorized by ARC Devices, and 4) accidents.

Shipping cost to return a device to ARC Devices or an ARC Devices Service center is not included.

A service notification number must be obtained from ARC Devices prior to returning any products or accessories to ARC Devices' designated service centers for repair. To obtain a service notification number, contact ARC Devices Technical Support.

THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. ARC DEVICES'S OBLIGATION UNDER THIS WARRANTY IS LIMITED TO REPAIR OR REPLACEMENT OF PRODUCTS CONTAINING A DEFECT. ARC DEVICES IS NOT RESPONSIBLE FOR ANY INDIRECT OR CONSEQUENTIAL DAMAGES RESULTING FROM A PRODUCT DEFECT COVERED BY THE WARRANTY

LIMITED TWO YEAR WARRANTY

How to return your thermometer

Should service be required during or after the warranty period, you must call ARC InstaTemp MD™ Customer Service Department to obtain a return authorization number and shipping address. Repackage the thermometer carefully in its original packaging or securely pack to avoid damage during shipping. Include the original sales slip indicating the date of purchase, a note describing the problem, your return authorization number, and your return address.

Send the ARC InstaTemp MD™ thermometer prepaid and insured.

For Customer Service Department contact:

USA and Canada: +1 561 287 5487

Monday - Friday 8:00 AM to 5:00 PM ET

EMEA and RoW: +353 91 886 002

Monday - Friday 8:00 AM to 5:00 PM GMT

CONTACT ARC DEVICES

Designed, manufactured and distributed in all territories except USA by:

ARC Devices Ltd
2 Deerpark Industrial Estate
Oranmore
Co Galway
Ireland

Phone: +353 (91) 398 191

Email: clinical@arcdevices.com

Web: www.arcdevices.com

Distributed in USA by:

ARC Devices USA, Inc.
1200 N. Federal Hwy, Ste 207,
Boca Raton, FL 33432
USA

Phone: 561.282.6074

CE
0123

PT-100153-IMD Rev 00

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Copyright ARC Devices Ltd. 2015