



Instructions for Use


---

REBOUND<sup>®</sup> DIABETIC WALKER

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	11
ES	Instrucciones para el uso	14
IT	Istruzioni per l'uso	17
NO	Bruksanvisning	20
DA	Brugsanvisning	23
SV	Bruksanvisning	26
EL	Οδηγίες χρήσης	29
FI	Käyttöohjeet	32
NL	Gebruiksaanwijzing	35
PT	Instruções de Utilização	38
PL	Instrukcja użytkownika	41
CZ	Návod k použití	44
TR	Kullanım Talimatları	47
RU	Инструкция по использованию	50
JA	取扱説明書	53
ZH	中文说明书	56
KO	사용 설명서	58

## INDICATION FOR USE

Promotes healing of plantar ulcers by offloading and by protection and immobilization of the foot.

## CAUTIONS/WARNINGS

- **This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.**
- Read the following instructions carefully before using the device.
- Rx Only.
- Modifications to this product are to be made only by a physician's order.
- Weight limit: User's maximum weight is 136 kg (300 lbs).
- Hexplugs should only be removed or altered by a healthcare provider as instructed by a physician.
- Wounds under treatment should be monitored closely for maintenance of pressure relief and healing progress.
- In case of emergency the compliance strap can be cut with scissors for removal of the walker.
- The insole should not be used in conjunction with any other walker or footwear.
- Before reapplying the device inspect:
  - The foot for any rash, sores, new ulcers or other indications of skin irritation.
  - The boot, bladders and liners for any damage.
- User should be cautious when walking on slippery or wet surfaces to avoid injury.
- If skin irritation of any kind develops or any problems occur with the use of this product, immediately contact your healthcare provider.
- Do not remove the walker unless instructed by a physician or for emergencies.
- If the walker has become loose or uncomfortable to wear, readjust the straps as necessary and contact your healthcare provider if discomfort is not relieved.
- This product has been designed and tested for single patient usage and is not recommended for multiple patient usage.
- This product is not intended for use on patients with chronic Charcot feet with a rocker bottom deformity.
- This product is only intended for use in treating plantar ulcers of the foot.

## INSTRUCTIONS FOR FITTING

1. Unfasten the straps, remove the front shell and open the liner.
2. Cover the wound with an appropriate dressing.
3. To heat mold the insole:
  - Apply the supplied sock over the patient's foot/dressing, prior to placing foot on heated insole.
  - Do not remove hexplugs prior to heat molding.
  - Pre-heat the oven to 110°C (230°F).
  - Heat the insole for three minutes at 110°C (230°F).
  - Have the patient stand, wearing the sock, on the pre-heated insole for one minute (**Figure 1**).
  - Have the patient remove foot from the insole.

4. Mark the wound location on the dressing with a transferable mark using iodine, lipstick, indelible ink pencil, etc (**Figure 2**).
5. Reposition the insole within the liner, with the heel properly located at the back of the walker.
6. While seated, place the foot inside walker. Ensure proper location of patient's foot before transferring mark to the insole (**Figure 3**).
7. Wrap liner securely around foot and leg. Have user bear weight to ensure transfer of relief zone mark.
8. With patient seated, remove walker. Remove liner and insole. Lightly press on the transferred mark and remove elevated hexplugs (**Figure 4**). **Caution:** Do not relieve too large an area that may allow ulcer to contact the bottom of the walker while weight bearing.
9. Replace insole in the liner with the hexplugs facing away from the patient's foot. Replace the liner, with the insole in the walker. Place foot inside walker and wrap liner securely around foot and leg.
10. Replace the front shell and fasten ankle (A), foot (B), and calf (C) straps in this order until snug and comfortable (**Figure 5**).  
**Note:** Use of the product without the front shell may result in pressure points and is not recommended.
11. Ensure proper location using the holes in the back of the walker to confirm the ankle is centered and that the heel is not pressed against the posterior of the walker (**Figure 6**).
12. Inflate air chambers by pressing "+" on the air bulb (**Figure 7**), until the walker is snug. **Note:** The Air System is designed with a feature to prevent over-inflation.
13. As needed, deflate the air chambers by pressing "-" on the air bulb (**Figure 8**). **Note:** Always deflate the Air System prior to removing and replacing the foot in the walker.
14. To customize the fit of the walker, the circumference around the calf can be widened by cutting the Flex Edge™ (**Figure 9**).
15. A disposable compliance strap may be placed around the ankle of the walker to inhibit removal. Place the strap through the strap guide on the front shell before securing the compliance strap (**Figure 10**).
  - The compliance strap is for single use only and should be replaced by the clinician when walker is removed.
  - Cut the compliance strap with scissors to remove the walker.
  - Use of the compliance strap is recommended.

## DONNING AND DOFFING INSTRUCTIONS

- To remove the walker, deflate air chambers, unfasten the straps, remove the front shell and open the liner.
- To replace the walker, place the leg inside walker and wrap liner securely around foot and leg. Place front shell over anterior calf and dorsum of foot. Fasten ankle, foot, and calf straps in this order. Tighten straps until snug and comfortable. Inflate air chambers by pressing "+" on the air bulb (**Figure 7 & 8**), until snug.

## CLEANING INSTRUCTIONS

### Insole

- **Do not use a washing machine or dryer.**
- **Do not use solvents or alcohol to clean.**
- Remove the insole from the liner.
- Rinse the insole with warm or cold water [ $\leq 105^{\circ}\text{F}$  ( $40^{\circ}\text{C}$ )].
- The insole should be air dried. To expedite the drying process, a hair dryer may be used at the lowest temperature setting.
- Once the insole is fully dried it can be reinserted into the liner.

- Always follow the fitting instructions to replace the insole and liner into the walker.

### **Walker Frame**

- To clean the walker frame, use a damp cloth, mild soap and warm or cold water [ $\leq 105^{\circ}\text{F}$  ( $40^{\circ}\text{C}$ )]. Air dry or dry with a towel.

### **Liner**

- To hand wash the liner, use warm or cold water [ $\leq 105^{\circ}\text{F}$  ( $40^{\circ}\text{C}$ )] and mild soap. **Do Not Use Bleach.**
- To machine wash the liner, use warm or cold setting [ $\leq 85^{\circ}\text{F}$  ( $30^{\circ}\text{C}$ )] and mild soap. Tumble dry on a low heat setting.

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**EN – Caution:** This product has been designed and tested as a single use product, and is not intended for multiple use. In addition to increasing the risk of infection, used product components may contain small unnoticeable defects that may lead to non-optimal performance, and even mechanical failure. If any problems occur with the use of this product, immediately contact your medical professional. This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

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