



The Right Time.
The Right Fit.

CASEY

Paralegal
NATRELLE® Style 15

YOUR SURGERY PLANNER

BREAST AUGMENTATION WITH *NATRELLE*®
SILICONE-FILLED BREAST IMPLANTS

Natrelle®

 ALLERGAN

Dear Patient,

Allergan has developed this *AUGMENTATION SURGERY WITH NATRELLE® SILICONE GEL-FILLED BREAST IMPLANTS PATIENT PLANNER* to function as a resource for all aspects of your surgery. Please give yourself at least 1 to 2 weeks time to consider this information before deciding to proceed with surgery unless an earlier surgery is deemed medically necessary by your surgeon.

This patient planner should serve primarily as your source of information on the risks and benefits of surgery with **NATRELLE®** Silicone-Filled Breast Implants but also as a convenient place where everything necessary for planning, follow-up and record keeping can be securely stored.

The information contained in Section I is intended to provide you with an understanding of the risks and benefits of surgery with silicone gel-filled breast implants, as well as provide an overview of the experience of patients in the Allergan Core Clinical Study.

Please thoroughly review this information. Following your review, complete the Patient Self Assessment. This assessment will help determine your understanding of the information presented and help your surgeon ensure that your preoperative consultation is effective and comprehensive. Make notes about issues that you would like to further discuss with your plastic surgeon, and ask questions. Give yourself time to consider your choices and proceed with surgery only after you are satisfied you understand the risks and follow-up recommendations associated with silicone gel-filled breast implants, and that the decision to proceed is the right decision for you.

You should become familiar with and use the following components provided in this planner:

SECTION I

- Important Information for Women About Breast Augmentation with **NATRELLE®** Silicone-Filled Breast Implants

SECTION II – FORMS

- Preoperative and Postoperative Checklists and Instructions
- Patient Self Assessment
- Acceptance of Risk and Surgery Consent
- Patient Surgery Record
- Device Tracking Enrollment
- Optional *ConfidencePlus®* Premier Warranty Enrollment Form
- Mammography Information
- To the Primary Care Physician

KATELYN

Mother of Four
NATRELLE® Style 20

BREAST AUGMENTATION WITH NATRELLE®
SILICONE-FILLED BREAST IMPLANTS



Table of Contents

Glossary	1
1. Considering Silicone Gel-Filled Breast Implant Surgery	9
1.1 What gives the breast its shape?	10
1.2 What is a silicone filled breast implant?	10
1.3 Are silicone gel-filled breast implants right for you?	11
1.4 What important factors should you consider in choosing silicone gel-filled implants?	12
2. Breast Implant Complications	19
2.1 What are the potential complications?	19
2.2 What are other reported conditions?	29
3. Allergan’s Clinical Study Results	35
3.1 What are the overview findings of Allergan’s Core Study?	35
3.2 What are the 7-year follow-up rates?	36
3.3 What are the benefits?	36
3.4 What are the 7-year complication rates?	38
3.5 What are the main reasons for reoperation?	41
3.6 What are the main reasons for implant removal?	42
3.7 What are other clinical data findings?	44
4. Surgical Considerations for Breast Augmentation	47
4.1 What are the alternatives to breast augmentation with silicone gel-filled breast implants?	47
4.2 What are questions to consider when choosing a surgeon?	47
4.3 What are other choices and options associated with the surgery?	48
4.4 What follow up examinations are important?	53
5. Additional Information	57
5.1 What types of NATRELLE [®] silicone-filled breast implants are available from Allergan?	57
5.2 What if I experience a problem?	58
5.3 What is device tracking?	59
5.4 What are the <i>ConfidencePlus</i> [®] warranties?	60
5.5 How can I receive more information?	60
For Further Reading and Information	63
Index	68

Glossary

Note: A glossary word appears in **bold** the first time it occurs in the text of this brochure.

Areola	The pigmented or darker colored area of skin surrounding the nipple of the breast.
Asymmetry	Lack of proportion of shape, size, and/or position between the two breasts.
Atrophy	Thinning or diminishing of tissues or muscle.
Autoimmune disease	A disease in which the body mounts an “attack” response to its own tissues or cell types. Normally, the body’s immune mechanism is able to distinguish clearly between what is a normal substance and what is foreign. In autoimmune diseases, this system becomes defective and mounts an attack against normal parts of the body, causing tissue injury. Certain diseases such as rheumatoid arthritis, lupus, and scleroderma are considered to be autoimmune diseases.
Axillary	Pertaining to the armpit area.
Biocompatible	The condition of being compatible with living tissues or systems without being toxic.
Biopsy	The removal and examination of tissues, cells, or fluid from the body.
Body Esteem Scale	A questionnaire which asks about a person’s body image.
Breast augmentation	A surgical procedure to increase breast size. For this document, it refers to placement of a breast implant. The first time a breast implant is placed to increase breast size, it is called primary augmentation. All subsequent times the implant is replaced, it is called revision-augmentation.
Breast implant	An internal artificial device or implant intended to replace the breast.

Breast mass	A lump in the breast.
Breast reconstruction	A surgical procedure to replace breast tissue that has been removed due to cancer or trauma or that has failed to develop properly due to a severe breast abnormality.
Calcification	Process of hardening by calcium salts.
Capsular contracture	<p>A tightening of the tissue capsule surrounding an implant, resulting in firmness or hardening of the breast and in squeezing of the implant if severe. Capsular contracture is classified by Baker Grades. Baker Grades III or IV are the most severe. Baker Grade III often results in the need for additional surgery (reoperation) because of pain and possibly abnormal appearance. Baker Grade IV usually results in the need for additional surgery (reoperation) because of pain and unacceptable appearance. Capsular contracture Baker Grade II may also result in the need for additional surgery. Capsular contracture is a risk for implant rupture. Below is a description of each Baker Grade.</p> <ul style="list-style-type: none"> • Baker Grade I – Normally soft and natural appearance • Baker Grade II – A little firm, but breast looks normal • Baker Grade III – More firm than normal, and looks abnormal (change in shape) • Baker Grade IV – Hard, obvious distortion, and tenderness with pain
Capsule	Scar tissue which forms around the breast implant. Sometimes this capsule squeezes the implant, resulting in capsular contracture.
Capsulectomy	Surgical removal of the scar tissue capsule around the implant.
Capsulotomy (closed)	An attempt to break the scar tissue capsule around the implant by pressing or pushing on the outside of the breast. This method does not require surgery but is a known risk for rupture of the implant and is contraindicated.

Capsulotomy (open)	An attempt to break the scar tissue capsule around the implant by surgical incision into the capsule.
Congenital abnormality	An abnormal development in part of the body, present in some form since birth.
Connective tissue disease/disorder (CTD)	A disease, group of diseases, or conditions affecting connective tissue, such as muscles, ligaments, skin, etc., and/or the immune system. Connective Tissue Diseases (“CTDs”) that involve the immune system include autoimmune diseases such as rheumatoid arthritis, lupus, and scleroderma.
Contraindication	A use that is improper and should not be followed. Failure to follow contraindications identified in the labeling could cause serious harm.
Contralateral	Opposite side.
Core Study	The primary clinical study of augmentation, reconstruction, and revision (revision-augmentation and revision-reconstruction) patients that supported the approval of the premarket approval (PMA) application. Safety and effectiveness data are collected yearly through 10 years, with the follow-up from years 5 through 10 being performed as part of a postapproval Core Study.
Delayed wound healing	Delayed progress in the healing of an opened wound.
Diffusion	Movement from one location to another.
Displacement	Movement of the implant from the usual or proper place.
Dysmorphic disorder	A psychological condition characterized by an obsession with a minor or an imagined physical flaw to the point that it can interfere with normal functioning.
Epidemiological	Relating to the science of explaining the relationships of factors that determine disease frequency and distribution.

Extracapsular rupture	A type of rupture in which the silicone gel is outside of the scar tissue capsule surrounding the implant.
Extrusion	Skin breakdown with the pressing out of the implant through the surgical wound or skin.
Fibromyalgia	A disorder characterized by chronic pain in the muscles and soft tissues surrounding joints, with tenderness at specific sites in the body. It is often accompanied by fatigue.
Fibrous tissues	Connective tissues composed mostly of fibers.
Granuloma	A lump or mass made of inflammatory cells surrounding a foreign substance due to longstanding inflammation.
Hematoma	A collection of blood within a space.
Hypertrophic scarring	An enlarged scar remaining after the healing of a wound.
Immune response	A bodily response to the presence of a foreign substance.
Infection	Invasion with microorganisms (for example, bacteria, viruses). An infection usually results in fever, swelling, redness, and/or pain.
Inflammation	The response of the body to infection or injury that is characterized by redness, swelling, warmth, pain, and/or loss of function.
Inframammary	Below the breast.
Inframammary fold	The crease at the base of the breast.
Inframammary incision	An incision made in the fold below the breast.
Inpatient surgery	A surgical procedure in which the patient is required to stay overnight in the hospital.
Intracapsular rupture	A type of rupture in which the silicone gel remains inside the scar tissue capsule surrounding the implant.

Lactation	The production and secretion of milk by the breast glands.
Low molecular weight silicones	Small silicone molecules that might leak out of the implant.
Lymphadenopathy	Enlargement of the lymph node(s).
Lymphedema	Swelling of the lymph nodes.
MRI	Magnetic resonance imaging. A radiographic examination that currently has the best ability to detect rupture of silicone gel-filled breast implants.
Malposition	Implant malposition or displacement is when the implant is not in the correct spot in the breast. This could have been due to incorrect placement of the implant during the surgery or due to shifting of the implant position over time.
Mammary	Pertaining to the breast.
Mammography	A type of X-ray examination of the breasts used for detection of cancer. Screening mammography – X-ray examination of the breast that is performed on women with no complaints or symptoms of breast cancer; the goal is to detect breast cancer when it is still too small to be felt by a physician or the patient. Diagnostic mammography – X-ray examination in order to evaluate a breast complaint or abnormality detected by physical exam or screening mammography; additional views of the breast are usually taken.
Mammoplasty	Plastic surgery of the breast.
Mastitis	Inflammation of the breast.
Mastopexy	Plastic surgery to lift sagging breasts higher.
Metastatic Disease	Spreading of cancer cells from the original site to other parts of the body.

Migration	Movement of silicone materials outside the breast implant.
Necrosis	Death of cells or tissues.
Outpatient surgery	A surgical procedure in which the patient is not required to stay in the hospital overnight.
Palpability	The ability to feel the implant.
Palpate/palpable	To feel with the hand.
Paresthesia	The feelings of pins and needles in a particular area of your body (particularly the arms and legs).
Pectoralis	Major muscle of the chest.
Periareolar	Around the darkened or pigmented area surrounding the nipple of the breast.
Plastic surgery	Surgery intended for the improvement of appearance of the body.
Postoperatively	After surgery.
Primary breast augmentation	The first time a breast implant is placed for the purpose of breast augmentation.
Ptoxis	Breast sagging that is usually the result of normal aging, pregnancy, or weight loss.
Reoperation	An additional surgery after your first breast implantation.
Revision-augmentation	Refers to the correction or improvement of a primary augmentation. In the context of this document, it refers to surgical removal and replacement of breast implants that were placed originally for primary breast augmentation.
Rheumatological disease/disorder	A variety of diseases involving connective tissue structures of the body, especially the joints and fibrous tissue. These diseases are often associated with pain, inflammation, stiffness, and/or limitation of motion of the affected parts. Can include autoimmune diseases. Fibromyalgia is a rheumatological disorder.

Rosenberg Self-Esteem Scale	A 10-item questionnaire intended to measure overall self-esteem through statements related to feelings of self-worth and self-acceptance.
Rupture	A tear or hole in the implant shell. Silicone implant ruptures may be with or without symptoms. Ruptures can be intracapsular or extracapsular.
Saline	A solution that is made up of water and a small amount of salt.
Scar revision	A surgical procedure to improve the appearance of a scar.
Seroma	A build-up of the watery portion of the blood in a tissue location.
SF-36 Scale	A 36-item questionnaire intended to measure patient health in areas such as vitality, physical functioning, bodily pain, general health, social and emotional functioning, and mental health.
Silent rupture	A breast implant rupture without symptoms and which is not apparent except through appropriate imaging techniques such as MRI. Most silicone gel-filled breast implant ruptures are silent (see symptomatic rupture next page).
Silicone elastomer	A type of silicone that has elastic properties similar to rubber.
Subglandular placement	Placement of a breast implant underneath and within the breast glands but on top of the chest muscle.
Submuscular placement	Placement of a breast implant wholly or partially underneath the chest muscle.
Surgical incision	A cut made to body tissue during surgery.
Symptom	Any perceptible change in the body or its functions that indicates disease or a phase of a disease.
Symptomatic	Any evidence or sign of disease or disorder reported by the patient.

**Symptomatic
rupture**

A breast implant rupture that is associated with symptoms (such as lumps, persistent pain, swelling, hardening, or change in implant shape). Some silicone breast implant ruptures are symptomatic, but most are silent.

Systemic

Pertaining to or affecting the body as a whole.

**Tennessee
Self-Concept Scale**

A questionnaire intended to measure the patient's view of her body and state of health, as well as her attitude about appearance, skills, and sexuality. The questionnaire administered in the Core Study consisted of 18 items.

**Toxic shock
syndrome**

Infection from staphylococci, occurring most often in the vagina of menstruating women using superabsorbent tampons but can also occur in other soft tissue infections. Symptoms include high fever, vomiting, diarrhea, rash, decreased blood pressure and shock, which can result in death.

Transaxillary

Through the axilla (armpit); an incision made under the arm.



1. CONSIDERING SILICONE GEL-FILLED BREAST IMPLANT SURGERY

You may be considering **breast implant** surgery to increase the size of your breasts. This is referred to as **breast augmentation**. Or you may need revision of a previous breast augmentation, which is called **revision-augmentation**. Allergan has prepared this information to help you better understand the **breast implant** procedure and assist you in making an informed decision about breast augmentation or revision-augmentation surgery. It will help to answer some of the questions you may have about the surgery and about breast implants in general. It will also provide you with specific information about the risks and benefits of **NATRELLE®** Silicone-Filled Breast Implants.


This information cannot and should not replace discussing your surgery with your plastic surgeon. Your decision whether or not to get breast implants should be based on realistic expectations of the outcome. There is no guarantee that your results will match those of other women. Your results will depend on many individual factors, such as your overall health (including age), chest structure, breast/nipple shape and position, skin texture, healing capabilities (which may be slowed by radiation and chemotherapy treatment, smoking, alcohol, and various medications), tendency to bleed, prior breast surgery, surgical team's skill and experience, type of surgical procedure, and type and size of implant. Make sure you speak with your surgeon about your expectations of the results, as well as what you can expect regarding the length of the surgery, your recovery, and any risks and potential complications of the surgery. Ask questions.

As part of your decision, both you and your surgeon will be required to sign Allergan's consent to surgery form that confirms your understanding of what you have read. This Allergan consent document will be provided to you by your surgeon.

You should wait at least 1-2 weeks after reviewing and considering this information before deciding whether to have **primary breast augmentation** surgery. In the case of a revision-augmentation, however, your surgeon may find it medically necessary to perform surgery sooner.

1.1 What Gives the Breast Its Shape?

The breast consists of milk ducts and glands, surrounded by fatty tissue that provides its shape and feel. Beneath the breast is the chest muscle (**pectoralis** major muscle).



Factors such as pregnancy (when milk glands are temporarily enlarged), rapid weight loss, and the effects of gravity as you age combine to stretch the skin, which may cause the breast to droop or sag. However, it is important to realize that implants are used to make the breast larger. The implants alone may not adequately lift the breast, or correct the effects of pregnancy, weight loss, or skin stretching. Your surgeon may suggest additional procedures at the time of the breast augmentation, such as **mastopexy**, to help achieve improved breast lift.

1.2 What Is a Silicone Filled Breast Implant?

A silicone gel-filled breast implant is a sac (implant shell) of **silicone elastomer** (rubber) filled with silicone gel. It is surgically implanted either under your breast tissue or under your chest muscle.



1.3 Are Silicone Gel-Filled Breast Implants Right for You?

NATRELLE® Silicone-Filled Breast Implants are indicated for females for the following uses (procedures):

- **Breast augmentation for women at least 22 years old.** Breast augmentation includes primary breast augmentation to increase the breast size, as well as revision surgery to correct or improve the result of a primary breast augmentation surgery.
- **Breast reconstruction.** Breast reconstruction includes primary reconstruction to replace breast tissue that has been removed due to cancer or trauma or that has failed to develop properly due to a severe breast abnormality. Breast reconstruction also includes revision surgery to correct or improve the result of a primary breast reconstruction surgery. (A separate patient brochure is available for those women considering breast reconstruction surgery and should be read prior to reaching a decision to undergo breast reconstruction).

CONTRAINDICATIONS

Breast implant surgery should not be performed in:

- Women with active **infection** anywhere in their body.
- Women with existing cancer or pre-cancer of their breast who have not received adequate treatment for those conditions.
- Women who are currently pregnant or nursing.

PRECAUTIONS

Safety and effectiveness have not been established in patients with the following:

- **Autoimmune diseases** (for example, lupus and scleroderma).
- A weakened immune system (for example, currently taking drugs that weaken the body's natural resistance to disease).
- Conditions that interfere with wound healing and blood clotting.
- Reduced blood supply to breast tissue.
- Radiation to the breast following implantation.
- Clinical diagnosis of depression or other mental health disorders, including body **dysmorphic disorder** and eating disorders. Please discuss any history of mental health disorders with your surgeon prior to surgery. Patients with a diagnosis of depression, or other mental health disorders, should wait until resolution or stabilization of these conditions prior to undergoing breast implantation surgery.

1.4 What Important Factors Should You Consider in Choosing Silicone Gel-Filled Implants?

- **Breast implants are not lifetime devices, and breast implantation is likely not a one-time surgery.** You will likely need additional unplanned surgeries on your breasts because of complications or unacceptable cosmetic outcomes. These additional surgeries can include implant removal with or without replacement, or they can include other surgical procedures. When you have your implants replaced (revision-augmentation), your risk of future complications increases compared to first time (primary) augmentation surgery, so you should also review the complication rates for revision-augmentation patients to see what future risk rates you may experience.



ERICA

Bookkeeper
NATRELLE® Style 20

- Many of the changes to your breast following implantation are irreversible (cannot be undone). If you later choose to have your implant(s) removed and not replaced, you may experience unacceptable dimpling, puckering, wrinkling, or other cosmetic changes of the breast, which can be permanent.
- Breast implants may affect your ability to breastfeed, either by reducing or eliminating milk production.
- **Rupture** of a silicone gel-filled breast implant is most often without **symptoms** (silent). This means that most of the time neither you nor your surgeon will know that your implants have a rupture. In fact, the ability of a physical examination by a plastic surgeon who is familiar with breast implants to detect silicone breast implant rupture is 30%³ compared to 89% for **MRI**.⁴ You will need regular screening MRI examinations over your lifetime in order to determine if **silent rupture** is present. You should have your first MRI at 3 years after your initial implant surgery and then every 2 years, thereafter. The cost of MRI screenings may exceed the cost of your initial surgery over your lifetime. This cost, which may not be covered by your insurance, should be considered in making your decision.
- If implant rupture is noted on MRI, you should have the implant removed, with or without replacement.
- With breast implants, routine screening **mammography** for breast cancer will be more difficult. If you are of the proper age for mammography screening, you should continue to undergo routine mammography screening as recommended by your primary care physician. The implant may interfere with finding breast cancer during mammography. Because the breast and implant are squeezed during mammography, an implant may rupture during the procedure. More x-ray views are necessary for women with breast implants; therefore, you will receive more exposure to radiation. However, the benefit of having the mammogram to find cancer outweighs the risk of the additional X-rays. Be sure to inform the mammography technologist that you have implants.

- You should perform an examination of your breasts every month for cancer screening; however, this may be more difficult with implants. You should ask your surgeon to help you distinguish the implant from your breast tissue.
- You should perform an examination of your breasts for the presence of lumps, swelling, hardening, or change in implant shape, which may be signs of **symptomatic rupture** of the implant. Report any of these symptoms or persistent pain to your surgeon. Your surgeon may recommend an evaluation via MRI to screen for rupture.
- After undergoing breast augmentation surgery (either primary or revision), your health insurance premiums may increase, your insurance coverage may be dropped, and/or future coverage may be denied. Treatment of complications may not be covered as well. You should discuss the complete extent of your insurance coverage with your insurance company before undergoing surgery.
- You should inform any other doctor who treats you of the presence of your implants to minimize the risk of damage to the implants.
- Allergan will continue its ongoing clinical **Core Study** through 10 years to further evaluate the long-term safety and effectiveness of these products. In addition, Allergan has initiated a separate 10-year postapproval study (the Breast Implant Follow-Up Study, or BIFS) to address specific issues for which the Allergan Core Study was not designed to fully answer, as well as to provide a real-world assessment of some endpoints. The endpoints in the large postapproval study include long-term local complications, **connective tissue disease (CTD)**, CTD signs and symptoms, neurological disease, neurological signs and symptoms, offspring issues, reproductive issues, **lactation** issues, cancer, suicide, mammography issues, and MRI compliance and results. Allergan will update its labeling on a regular basis with the results of these two studies. You should also ask your surgeon for any available updated Allergan clinical information.



MINELA

Student
NATRELLE® Style 15

- It is important that you read this entire brochure because you need to understand the risks and benefits and to have realistic expectations of the outcome of your surgery.

2. BREAST IMPLANT COMPLICATIONS

Undergoing any type of surgical procedure involves risks (some serious) such as the effects of anesthesia, infection, swelling, redness, bleeding, pain, and even death, which need to be balanced against the benefits of the surgery itself. There are potential complications specific to breast implant surgery and breast implants, as described below. Located at the end of this brochure is a list of the study reports used to gather the information discussed in the sections below. These may be helpful to you if you wish to learn more about a specific complication or condition. The reference list is not complete because studies are being conducted all the time; your physician may have other resources for further reading as well. It should be noted that the references include augmentation and/or reconstruction patients, as well as implants of different types and from a variety of manufacturers.

2.1 What Are the Potential Complications?

• Rupture

Breast implants are not lifetime devices. Breast implants rupture when the shell develops a tear or hole. Ruptures can occur at any time after implantation, but they are more likely to occur the longer the implant is implanted. The following things may cause your implant to rupture: damage by surgical instruments; stressing the implant during implantation which may weaken it; folding or wrinkling of the implant shell; excessive force to the chest (for example, during closed capsulotomy, which is contraindicated); trauma; compression during mammographic imaging; and severe **capsular contracture**. Breast implants may also simply wear out over time. Laboratory studies have identified some of the types of rupture for Allergan's product; however, it is not known whether these tests have identified all causes of rupture. These laboratory studies will continue post-approval.

Silicone gel-filled breast implant ruptures are most often silent. (MRI examination is currently the best method to screen for silent rupture). This means that most of the time neither you nor your plastic surgeon will know if the implant has a tear or hole in the shell. This is why MRI is recommended at 3 years and then every 2 years, thereafter, to screen for rupture. However, sometimes there are symptoms associated with gel implant rupture. These symptoms include hard knots or lumps surrounding the implant or in the armpit, change or loss of size or shape of the breast or implant, pain, tingling, swelling, numbness, burning, or hardening of the breast.

When MRI signs of rupture are found, or if your surgeon determines you have signs or symptoms of rupture, you should have the implant and any gel removed, with or without replacement of the implant. It also may be necessary to remove the tissue **capsule** as well as the implant, which will involve additional surgery, with associated costs. If you have symptoms such as breast hardness, a change in breast shape or size, and/or breast pain, you should have an MRI to determine whether rupture is present.^{1,5}

There are also consequences of rupture. If rupture occurs, silicone gel may either remain within the scar tissue capsule surrounding the implant (**intracapsular rupture**), move outside the capsule (**extracapsular rupture**), or gel may move beyond the breast (**gel migration**). There is also a possibility that rupture may progress from intracapsular to extracapsular and beyond. There have also been health consequences reported in the literature. See below for details.

Rupture Information on Allergan Implants

In Allergan's Core Study, rupture was assessed for patients who had scheduled MRIs to screen for silent rupture (i.e., part of the MRI cohort) and those who were not assessed for rupture by MRI (i.e., part of the non-MRI cohort). For primary augmentation patients in the MRI cohort, the by-patient rupture rate was 8.6% and the by-implant rupture rate was 5.1% through 7 years. For revision-augmentation patients in the MRI cohort, the rupture rate was 0% through 7 years. This means that through 7 years, approximately 9 of every 100 primary augmentation women had at least one ruptured breast implant. For the non-MRI cohort, there were 3 primary augmentation and 2 revision-augmentation patients who had reported rupture through 7 years.

The rupture rate for the whole MRI cohort in the Core Study (including augmentation, revision-augmentation, reconstruction, and revision-reconstruction patients) through 7 years was 7.3% for patients and 4.5% for implants. Across all patients in the Core Study, all ruptures were intracapsular with 1 case of extracapsular gel (one rupture progressed to extracapsular gel following exploratory surgery to confirm the rupture and then implant replacement was delayed). There were no cases of migrated gel.

Further rupture rate information on **NATRELLE**[®] Silicone-Filled Breast Implants is provided from a published European study known as the International MRI Study.² Silent rupture data were collected via a single MRI on 77 augmentation, 11 reconstruction, and 18 revision patients implanted with smooth and textured **NATRELLE**[®] implants by five surgeons. The average age of the implants was approximately 11 years. Silent rupture was found in approximately 15% of the combined group of augmentation, reconstruction, and revision patients and 8% of the implants. There was one possible case of extracapsular rupture with the remainder classified as intracapsular ruptures. No cases of gel migration were found.

Additional information on rupture will be collected through Allergan's postapproval studies: the continuing Core Study and Breast Implant Follow-Up Study (BIFS).

Additional Information on Consequences of Rupture from Literature

Since silicone implants were not available in the United States for many years but were used in Europe during that time, some information on rupture rates comes from studies conducted in Europe. Studies of Danish women evaluated with MRI involving a variety of manufacturers and implant models showed that about three-fourths of implant ruptures are intracapsular and the remaining one-fourth are extracapsular.⁶ Additional studies of Danish women indicate that over a 2-year period, about 10% of the implants with intracapsular rupture progressed to extracapsular rupture as detected by MRI.⁵ This means that for women with silicone gel rupture within the scar tissue capsule detected via MRI after 2 years, 1 in 10 of these women had progression of the gel outside the scar tissue capsule. In about half of these cases of progression from intracapsular to extracapsular rupture, the women had experienced trauma or mammography.



QUYNH

Lawyer, Mother of Two
NATRELLE® Style 15

In the other half, no cause was given. In the women with extracapsular rupture, after 2 years, the amount of silicone outside the scar tissue capsule increased for about 14% of these women. This means that for 100 women with silicone gel rupture outside the scar tissue capsule, the amount of gel outside the scar tissue capsule increased for 14 women 2 years later. This type of information pertains to a variety of silicone implants from a variety of manufacturers and implant models, and it is not specific to Allergan implants.

Below is a summary of information related to the health consequences of implant rupture, which have not been fully established. These reports were in women who had implants from a variety of manufacturers and implant models.

- Local breast complications reported in the published literature which were associated with rupture include breast hardness, a change in breast shape or size, and breast pain.⁵ These symptoms are not specific to rupture, as they also are experienced by women who have capsular contracture.
- There have been rare reports of gel movement to nearby tissues such as the chest wall, armpit, or upper abdominal wall, and to more distant locations down the arm or into the groin. This has led to nerve damage, **granuloma** formation and/or breakdown of tissues in direct contact with the gel in a few cases. There have been reports of silicone presence in the liver of patients with silicone breast implants. Movement of silicone gel material to lymph nodes in the axilla also has been reported, even in women without evidence of rupture, leading to **lymphadenopathy**.⁴³
- Concerns have been raised over whether ruptured implants are associated with the development of connective tissue or rheumatic diseases and/or symptoms such as fatigue and **fibromyalgia**.^{11,12,18,19} A number of epidemiology studies have evaluated large populations of women with breast implants. These studies do not, taken together, support a significant association of breast implants with a typical, diagnosed rheumatic disease. Other than one small study,¹² these studies do not distinguish whether the women had ruptured or intact implants.

• Capsular Contracture

The scar tissue (capsule) that normally forms around the implant may tighten over time and compress the implant, making it feel firm and leading to what is called capsular contracture. Capsular contracture may be more common following infection, **hematoma**, and **seroma**, and the chance of it happening may increase over time. Capsular contracture occurs more commonly in revision-augmentation than in primary augmentation. Because you may have your initial implants replaced, you should be aware that your risk of capsular contracture increases with revision-augmentation. Capsular contracture is a risk factor for implant rupture, and it is the most common reason for reoperation.

Symptoms of capsular contracture range from mild firmness and mild discomfort to severe pain, distorted shape of the implant, and **palpability** (ability to feel the implant). Capsular contracture is graded into 4 levels depending on its severity.* Baker Grades III and IV are considered severe, and often additional surgery is needed to correct these grades.

- Baker Grade I:** the breast is normally soft and looks natural
- Baker Grade II:** the breast is a little firm but looks normal
- Baker Grade III:** the breast is firm and looks abnormal
- Baker Grade IV:** the breast is hard, painful, and looks abnormal

In Allergan's Core Study, for women receiving augmentation implants for the first time, the risk of severe capsular contracture was 16% through 7 years. This means that 16 out of every 100 women who received Allergan implants for primary breast augmentation had severe capsular contracture at least once during the first 7 years after receiving the implants.

For women receiving revision-augmentation implants, the risk of severe capsular contracture was 20% through 7 years.

Additional surgery may be needed in cases where pain and/or firmness are severe. This surgery ranges from removal of the implant capsule tissue to removal and possible replacement of the implant itself. This surgery may result in loss of your breast tissue. Capsular contracture may happen again after these additional surgeries. Capsular contracture may increase the risk of rupture.¹

* Baker, J.L. Augmentation mammoplasty. In: Owsley, J.Q. and Peterson, R., Eds. *Symposium on aesthetic surgery of the breast*. St. Louis, MO: Mosby; 1978:256-263.

• **Additional Surgeries (Reoperations)**

You should assume that you will need to have additional surgeries (**reoperations**). In Allergan's Core Study, the reoperation rate was 30% for primary augmentation patients and 41% for revision-augmentation patients, which means that 30 out of every 100 women who received Allergan implants for primary augmentation and 41 out of every 100 women who received Allergan implants for revision-augmentation had a reoperation during the first 7 years after receiving the implants.

Patients may decide to change the size or type of their implants, requiring additional surgery. In addition, problems such as rupture, capsular contracture, **hypertrophic scarring** (irregular, raised scar), **asymmetry**, infection, and shifting can require additional surgery. Tables 3 and 4, which summarize the main reasons for performing reoperations in the Core Study, are located in section 3.5. For women receiving primary augmentation implants, the three most common reasons for reoperations were capsular contracture, implant **malposition**, and **ptosis** (sagging). For women receiving revision-augmentation implants, the three most common reasons for additional surgery were capsular contracture, hematoma/seroma, and implant malposition.

• **Implant Removal**

Because these are not lifetime devices, the longer you have your implants the more likely it will be for you to have them removed for any reason, either because of dissatisfaction, an unacceptable cosmetic result, or a complication such as capsular contracture. Having your implants removed and replaced increases your chances of getting future complications.

For women receiving primary augmentation implants in Allergan's Core Study, 14% had their implants removed at least once through 7 years. Capsular contracture and patient request for style/size change were the most common reasons for implant removal. For women receiving revision-augmentation implants in Allergan's Core Study, 24% had their implants removed at least once through 7 years. The most common reasons for implant removal were capsular contracture, patient request for style/size change, and implant malposition.

Most women who have their implants removed have them replaced with new implants, but some women do not. If you choose not to replace your implants, you may have cosmetically unacceptable dimpling, puckering, wrinkling, and/or other potentially permanent cosmetic changes of the breast following removal of the implant. Even if you have your implants replaced, implant removal may result in loss of your breast tissue. Also, implant replacement increases your risks of future complications. For example, the risks of capsular contracture and reoperation increase for patients with implant replacement compared to first time placement. You should consider the possibility of having your implants replaced and its consequences when making your decision to have implants.

- **Unsatisfactory Results**

Unsatisfactory results such as wrinkling, asymmetry, implant **displacement** (shifting), incorrect size, unanticipated shape, implant palpability, scar deformity, and/or hypertrophic scarring, may occur. Some of these results may cause discomfort. Pre-existing asymmetry may not be entirely correctable by implant surgery. Revision surgery may be recommended to maintain patient satisfaction, but carries additional considerations and risks. Selecting an experienced plastic surgeon may minimize, but not necessarily prevent, unsatisfactory results.

- **Pain**

Pain of varying intensity and length of time may occur and persist following breast implant surgery. In addition, improper size, placement, surgical technique, or capsular contracture may result in pain. You should tell your surgeon about significant pain or if your pain persists.

• **Changes in Nipple and Breast Sensation**

Feeling in the nipple and breast can increase or decrease after implant surgery. The range of changes varies from intense sensitivity to no feeling in the nipple or breast following surgery. While some of these changes can be temporary, they can also be permanent, and may affect your sexual response or your ability to nurse a baby. (See the paragraph on breastfeeding below.)

• **Infection**

Infection can occur with any surgery or implant. Most infections resulting from surgery appear within a few days to weeks after the operation. However, infection is possible at any time after surgery. In addition, breast and nipple piercing procedures may increase the possibility of infection. Infections in tissue with an implant present are harder to treat than infections in tissue without an implant. If an infection does not respond to antibiotics, the implant may have to be removed, and another implant may be placed after the infection is resolved (cleared up). As with many other surgical procedures, in rare instances, **toxic shock syndrome** has been noted in women after breast implant surgery, and it is a life-threatening condition. Symptoms include sudden fever, vomiting, diarrhea, fainting, dizziness, and/or sunburn-like rash. You should contact a doctor immediately for diagnosis and treatment if you have these symptoms.

• **Hematoma/Seroma**

Hematoma is a collection of blood within the space around the implant, and a seroma is a build-up of fluid around the implant. Having a hematoma and/or seroma following surgery may result in infection and/or capsular contracture later on. Symptoms from a hematoma or seroma may include swelling, pain, and bruising. If a hematoma or seroma occurs, it will usually be soon after surgery. However, this can also occur at any time after injury to the breast. While the body absorbs small hematomas and seromas, some will require surgery, typically involving draining and potentially placing a surgical drain in the wound temporarily for proper healing. A small scar can result from surgical draining. Implant rupture also can occur from surgical draining if there is damage to the implant during the draining procedure.

- **Breastfeeding**

Breastfeeding difficulties have been reported following breast surgery, including breast reduction and breast augmentation. If your surgeon uses a periareolar surgical approach (an incision around the colored portion surrounding the nipple), it may further increase the chance of breastfeeding difficulties.

- **Calcium Deposits in the Tissue Around the Implant**

Calcium deposits can form in the tissue capsule surrounding the implant. Symptoms may include pain and firmness. Deposits of calcium can be seen on mammograms and can be mistaken for possible cancer, resulting in additional surgery for **biopsy** and/or removal of the implant to distinguish calcium deposits from cancer. If additional surgery is necessary to examine and/or remove **calcifications**, this may cause damage to the implants. Calcium deposits also occur in women who undergo breast reduction procedures, in patients who have had hematoma formation, and even in the breasts of women who have not undergone any breast surgery. The occurrence of calcium deposits increases significantly with age.

- **Extrusion**

Extrusion is when the breast implant comes through your skin. This may occur, for example, when your wound has not closed or when breast tissue covering your implants weakens. Radiation therapy has been reported to increase the likelihood of extrusion. Extrusion requires additional surgery and possible removal of the implant, which may result in additional scarring and/or loss of your breast tissue.

- **Necrosis**

Necrosis is the death of cells or tissues. This may prevent or delay wound healing and require surgical correction, which may result in additional scarring and/or loss of your breast tissue. Implant removal may also be necessary. Factors associated with increased necrosis include infection, use of steroids, smoking, chemotherapy/radiation, and excessive heat or cold therapy.

- **Delayed Wound Healing**

Some patients may experience a prolonged wound healing time. **Delayed wound healing** may increase the risk of infection, extrusion, and necrosis. Depending on the type

of surgery or the incision, wound healing times may vary. Smoking may interfere with the healing process. You should contact your surgeon immediately if your wound does not heal within the period of time he/she has discussed with you.

- **Breast Tissue Atrophy/Chest Wall Deformity**

The pressure of the breast implant may cause breast tissue thinning (with increased implant visibility and palpability) and chest wall deformity. This can occur while implants are still in place or following implant removal without replacement. Either of these conditions may result in additional surgeries and/or unacceptable dimpling/puckering of the breast.

- **Lymphadenopathy**

Lymphadenopathy is a chronic enlargement of the lymph nodes. A lymph node is a round mass of tissue which makes cells as part of your immune system. The lymph nodes in the armpit (axilla) drain the breast area of fluid. Some patients with breast implants report having enlarged lymph nodes in the armpit(s). Sometimes the enlarged lymph nodes are painful. If they become too large or painful, the lymph node(s) may need to be surgically removed. Painful and/or enlarged lymph nodes should be reported to your doctor.

Literature reports associate lymphadenopathy with both intact and ruptured silicone breast implants. One study reported that armpit lymph nodes from women with both intact and ruptured silicone gel implants had abnormal tissue reactions, granulomas, and the presence of silicone.⁴³ These reports were in women who had implants from a variety of manufacturers and implant models.

2.2 What Are Other Reported Conditions?

There have been reports in the literature of other conditions in women with silicone gel-filled breast implants. Many of these conditions have been studied to evaluate their potential association with breast implants. Although no cause-and-effect relationship has been established between breast implants and the conditions listed below, you should be aware of these reports. Furthermore, there is the possibility of risks, yet unknown, which in the future could be determined to be associated with breast implants.

• **Connective Tissue Disease (CTD)**

Connective tissue diseases include diseases such as lupus, scleroderma, and rheumatoid arthritis. **Fibromyalgia** is a disorder characterized by chronic pain in the muscles and soft tissues surrounding joints, with tenderness at specific sites in the body. It is often accompanied by fatigue. There have been a number of published **epidemiological** studies which have looked at whether having a breast implant is associated with having a typical or defined connective tissue disease.^{1,10-18}

The published studies overall show that breast implants are not significantly associated with a risk of developing a typical or defined connective tissue disease.^{1,13-15} However, the study size needed to conclusively rule out a small risk of connective tissue disease among women with silicone gel-filled implants would need to be very large. These studies do not distinguish between women with intact and ruptured implants. Only one study evaluated specific connective tissue disease diagnoses and symptoms in women with silent ruptured versus intact implants, but the study was too small to rule out a small risk.¹²

• **CTD Signs and Symptoms**

Literature reports have also been made associating silicone breast implants with various rheumatological signs and symptoms such as fatigue, exhaustion, joint pain and swelling, muscle pain and cramping, tingling, numbness, weakness, and skin rashes. Scientific expert panels and literature reports have found no evidence of a consistent pattern of signs and symptoms in women with silicone breast implants.^{1,19-22} Having these rheumatological signs and symptoms does not necessarily mean you have a connective tissue disease; however, you should be aware that you may experience these signs and symptoms after undergoing breast implantation. If you notice an increase in these signs or symptoms, you should consider seeing a rheumatologist to determine whether these signs or symptoms are due to a connective tissue disorder or autoimmune disease.

• **Cancer**

Breast Cancer – Reports in the medical literature indicate that patients with breast implants are not at a greater risk than those without breast implants for developing breast cancer.^{23,25,27,30,35} Some reports have suggested that breast

implants may interfere with or delay breast cancer detection by mammography and/or biopsy; however, other reports in the published medical literature indicate that breast implants neither significantly delay breast cancer detection nor adversely affect cancer survival of women with breast implants.^{23,28,31,34,35}

Brain Cancer – One recent study has reported an increased incidence of brain cancer in women with breast implants as compared to the general population.²⁴ The incidence of brain cancer, however, was not significantly increased in women with breast implants when compared to women who had other plastic surgeries. Another recently published review of four large studies of women with cosmetic implants concluded that the evidence does not support an association between brain cancer and breast implants.³³

Respiratory/Lung Cancer – One study has reported an increased incidence of respiratory/lung cancer in women with breast implants.²⁴ Other studies of women in Sweden and Denmark have found that women who get breast implants are more likely to be current smokers than women who get breast reduction surgery or other types of cosmetic surgery.^{26,29,32}

Cervical/Vulvar Cancer – One study has reported an increased incidence of cervical/vulvar cancer in women with breast implants.²⁴ The cause of this increase is unknown.

Other Cancers – One study has reported an increased incidence of stomach cancer and leukemia in women with breast implants compared to the general population.²⁴ This increase was not significant when compared to women who had other types of plastic surgeries.

• **Neurological Disease, Signs, and Symptoms**

Some women with breast implants have complained of neurological symptoms (such as difficulties with vision, sensation, muscle strength, walking, balance, thinking or remembering things) or diseases (such as multiple sclerosis), which they believe are related to their implants. A scientific expert panel report found that the evidence for a neurological disease or syndrome caused by or associated with breast implants is insufficient or flawed.¹

• Suicide

In several studies, a higher incidence of suicide was observed in women with breast implants.³⁶⁻³⁹ The reason for the observed increase is unknown, but it was found that women with breast implants had higher rates of hospital admission due to psychiatric causes prior to surgery, as compared with women who had breast reduction or in the general population of Danish women.³⁷

• Effects on Children

At this time, it is not known if a small amount of silicone may pass through from the breast implant silicone shell into breast milk during breastfeeding. Although there are no current established methods for accurately detecting silicone levels in breast milk, a study measuring silicon (one component in silicone) levels did not indicate higher levels in breast milk from women with silicone gel-filled implants when compared to women without implants.

In addition, concerns have been raised regarding potential damaging effects on children born to mothers with implants. Two studies in humans have found that the risk of birth defects overall is not increased in children born after breast implant surgery.^{41,42} Although low birth weight was reported in a third study, other factors (for example, lower pre-pregnancy weight) may explain this finding.⁴⁰ The author recommended further research on infant health.

• Potential Health Consequences of Gel Bleed

Small quantities of **low molecular weight (LMW) silicone** compounds, as well as platinum (in zero oxidation state), have been found to diffuse (bleed) through an intact implant shell.^{1,45} The evidence is mixed as to whether there are any clinical consequences associated with gel bleed. For instance, studies on implants implanted for a long duration have suggested that such bleed may be a contributing factor in the development of capsular contracture¹ and lymphadenopathy.⁴³ However, evidence against gel bleed being a significant contributing factor to capsular contracture and other local complications is provided by the fact that there are similar or lower complication rates for silicone gel-filled breast implants than for **saline**-filled breast implants. Saline-filled breast implants do not contain silicone gel and, therefore, gel bleed is not an issue for those products. Furthermore, toxicology testing has indicated that the silicone material used in

Allergan's implants does not cause toxic reactions when large amounts are administered to test animals. It should also be noted that studies reported in the literature have demonstrated that the low concentration of platinum contained in breast implants is in the zero oxidation (most **biocompatible**) state.^{44,46,47,49}

Allergan performed a laboratory test to analyze the silicones and platinum (used in the manufacturing process), which may bleed out of intact implants into the body. Over 99% of the LMW silicones and platinum stayed in the implant. The overall body of available evidence supports that the extremely low level of gel bleed is of no clinical consequence.

3. ALLERGAN'S CLINICAL STUDY RESULTS

This section of the brochure summarizes the results of the Allergan Core Study conducted on **NATRELLE**[®] Silicone-Filled Breast Implants for primary augmentation and revision-augmentation. The Allergan Core Study is the primary clinical study for this product. The results of the Core Study give you useful information on the experience of other women with **NATRELLE**[®] Silicone-Filled Breast Implants. While the results cannot be used to predict your individual outcome, they can be used as a general guide of what you may expect. Your own complications and benefits depend on many individual factors.

As a note, supplemental safety information was also obtained from another Allergan clinical study (the Adjunct Study), the Danish Breast Implant Registry, an international clinical MRI study, and the literature to help assess long-term rupture rate and the consequences of rupture for this product. The literature, which had the most available information on the consequences of rupture, was also used to assess other potential complications associated with silicone gel-filled breast implants. The key literature information was discussed throughout the **Breast Implant Complications** section above, and the references can be found at the end of this brochure.

3.1 What Are the Overview Findings of Allergan's Core Study?

The Allergan Core Study is a 10-year study to assess safety and effectiveness in augmentation, reconstruction, and revision (revision-augmentation and revision-reconstruction) patients. Patient follow-up is at 0-4 weeks, 6 months, 12 months, 24 months, and annually through 10 years. Safety is assessed by complications, such as implant rupture, capsular contracture, and reoperation. Benefit (effectiveness) is assessed by breast size change, patient satisfaction and measures of quality of life (QoL).

The Allergan Core Study consists of 715 patients. This includes 455 primary augmentation patients, 147 revision-augmentation patients, 98 primary reconstruction patients, and 15 revision-reconstruction patients. Of these patients, 158 primary augmentation patients, 50 revision-augmentation

patients, 51 primary reconstruction patients, and 5 revision-reconstruction patients are in the MRI cohort, which means that they are assessed for silent rupture by MRI at years 1, 3, 5, 7 and 9. The study is currently ongoing, with the results through 7 years reported in this brochure. Allergan will periodically update this brochure as more information becomes available. You should also ask your surgeon for any available updated Allergan clinical information.

Allergan's Core Study results indicate that the risk of at least one occurrence of any complication (including reoperation) at some point through 7 years after implant surgery is 45% for primary augmentation patients and 57% for revision-augmentation patients. The information below provides more details about the complications and benefits you may experience. More detailed data tables are found in the Appendix of this brochure. Please refer to the glossary for the definition of any complication you may not understand.

3.2 What Are the 7-Year Follow-Up Rates?

Follow-up rates from a clinical study show you how many women continue to provide information on their experience with breast implants.

The Allergan Core Study enrolled 455 augmentation patients. Of the women expected to be seen at the 7-year follow-up visit, 74% were seen.

The Allergan Core Study enrolled 147 revision-augmentation patients. Of the women expected to be seen at the 7-year follow-up visit, 72% were seen.

3.3 What Are the Benefits?


The benefits of **NATRELLE**[®] Silicone-Filled Breast Implants were assessed by a variety of outcomes, including bra cup size change and assessments of patient satisfaction and quality of life. Data were collected before implantation and at scheduled follow-up visits at 1, 2, 4 and 6 years postimplant for those patients who still had their original implants and who came back for these visits.

Breast Measurement: For primary augmentation patients, 396 (87%) of the original 455 patients had a breast measurement within 18 months of surgery. Of these 396 patients, 41% increased by 1 cup size; 45% increased by 2 cup sizes; 8% increased by more

than 2 cup sizes; and 6% had no increase or decrease. See Figure 1 below.

Revision-augmentation patients did not undergo a measurement of breast cup size change because they were undergoing replacement of an existing breast implant.


Figure 1
Cup Size Changes in Primary Augmentation Patients



Patient Satisfaction: Allergan’s patient satisfaction was based on a 5-point scale assessment of satisfaction with their implants at the time of the follow-up visits. Of the original 455 primary augmentation patients, 317 (70%) provided a satisfaction rating at 7 years after implantation, with 301 (95%) of these patients indicating that they were satisfied with their breast implants.

Of the original 147 revision-augmentation patients, 91 (62%) provided a satisfaction rating at 7 years. Of these 81 patients, 81 (89%) indicated that they were satisfied with their breast implants. See Figure 2 below.

Figure 2
Primary Augmentation and Revision-Augmentation Patient Satisfaction Through 7 Years



Quality of Life Assessments: Quality of life assessments were obtained prior to implantation and at 1, 2, 4, and 6 years post-surgery. The 6-year data is provided here. For primary augmentation patients, the **SF-36**, which is a collection of scales that measure mental and physical health, showed an improvement in one scale and a worsening in six scales after 6 years compared to before breast implantation, although all scales remained higher than the general U.S. female population. For patient responses to questions regarding overall self-concept/self-esteem, there was a decrease in self-concept on the **Tennessee Self-Concept Scale** and no change in overall self-esteem on the **Rosenberg Self-Esteem Scale** 6 years after receiving implants. Patient responses to questions on the **Body Esteem Scale** regarding overall body image showed no changes, but a decrease in weight concern and physical condition, and increase with regard to sexual attractiveness were shown.

For revision-augmentation patients, the SF-36 showed no significant changes in all of the scales but one, which showed a decrease after 6 years. Patient responses to questions on the Tennessee Self-Concept Scale and Rosenberg Self-Esteem Scale showed no changes 6 years after receiving implants. Patient responses to questions on the Body Esteem Scale regarding overall body image showed no changes, but a decrease with regard to physical condition was shown.

3.4 What Are the 7-Year Complication Rates?

The complications observed in primary augmentation and revision-augmentation women through 7 years are presented in Table 1 and Table 2, respectively. The rates reflect the percentage of patients who experienced the listed complication at least once within the first 7 years after their implantation. Some complications occurred more than once for some patients. Please refer to the Glossary at the front of this brochure for the definition of any complication you may not understand.

The most common complications experienced within the first 7 years of implantation for primary augmentation patients were reoperation (30% or approximately 30 patients out of 100) and capsular contracture (16% or 16 patients out of 100). The most common complications experienced within the first 7 years of

implantation for revision- augmentation patients were also reoperation (41%) and capsular contracture (20%).

Table 1
7-Year Complication Rates for Primary Augmentation Patients N = 455 Patients

Key Complications ¹	%
Reoperation	30.1%
Capsular Contracture Baker Grade III/IV	15.5%
Implant Removal with Replacement	11.0%
Implant Rupture (MRI cohort)	8.6%
Implant Removal Without Replacement	3.1%
Other Complications Occurring in ≥ 1% of patients ^{2,3}	%
Breast Pain	11.4%
Swelling	9.2%
Nipple Complications	6.7%
Implant Malposition	5.2%
Scarring/Hypertrophic Scarring	3.7%
Asymmetry	3.3%
Ptosis	2.2%
Breast/Skin Sensation Changes	1.6%
Hematoma	1.6%
Implant Palpability/Visibility	1.6%
Seroma/Fluid Accumulation	1.6%
Wrinkling/Rippling	1.2%
Delayed Wound Healing	1.1%

¹ Most events were assessed with severity ratings, and the rates shown in the table include only complications rated moderate, severe or very severe (excludes mild and very mild ratings). All occurrences of reoperation, implant removal, implant rupture, implant extrusion and pneumothorax are included.

² The following complications were reported at a rate less than 1%: bruising, infection, implant extrusion, lymphedema, other complications, redness, skin rash, and tissue/skin necrosis.

³ The following complications were reported at a rate of 0%: capsule calcification, gel migration, irritation, lymphadenopathy, and pneumothorax.

Table 2
7-Year Complication Rates for Revision-Augmentation Patients N = 147 Patients

Key Complications ¹	%
Reoperation	40.5%
Capsular Contracture Baker Grade III/IV	20.4%
Implant Removal with Replacement	20.9%
Implant Removal without Replacement	4.3%
Implant Rupture (MRI cohort)	0%
Other Complications Occurring in ≥ 1% of patients ^{2,3}	%
Breast Pain	10.6%
Swelling	8.4%
Implant Palpability/Visibility	6.8%
Implant Malposition	6.1%
Seroma/Fluid Accumulation	6.1%
Scarring/Hypertrophic Scarring	6.0%
Ptosis	4.8%
Wrinkling/Rippling	4.6%
Asymmetry	3.7%
Bruising	3.0%
Breast/Skin Sensation Changes	2.2%
Hematoma	2.1%
Infection	1.4%
<p>¹ Most events were assessed with severity ratings, and the rates shown in the table include only complications rated moderate, severe or very severe (excludes mild and very mild ratings). All occurrences of reoperation, implant removal, implant rupture, implant extrusion and pneumothorax are included.</p> <p>² The following complications were reported at a rate less than 1%: delayed wound healing, irritation, nipple complications, other complications, redness, and skin rash.</p> <p>³ The following complications were reported at a rate of 0%: capsule calcification, gel migration, implant extrusion, lymphadenopathy, lymphedema, pneumothorax, and tissue/skin necrosis.</p>	

3.5 What Are the Main Reasons for Reoperation?

The reasons for reoperation observed in primary augmentation and revision-augmentation women through 7 years are presented in Table 3 and Table 4, respectively. There may be one or more reasons identified for having a reoperation (additional surgery after the primary or revision breast augmentation). Furthermore, there may be multiple surgical procedures (for example, implant removal with or without replacement, capsule procedures, incision and drainage, implant reposition, **scar revision**, etc.) performed during a reoperation.

Table 3 Main Reasons for Reoperation in Primary Augmentation Patients Through 7 Years	
Reason for Reoperation	n
Capsular Contracture	47
Implant Malposition	24
Ptoxis	21
Need for Biopsy	19
Suspected Rupture	13
Hematoma/Seroma	12
Patient Request for Size/Style Change	9
Scarring/Hypertrophic Scarring	8
Asymmetry	6
Breast Cancer	3
Breast Pain	3
Delayed Wound Healing	3
Implant Palpability/Visibility	2
Infection	2
Wrinkling	2
Implant Extrusion	1
Necrosis	1
Nipple Complications	1
TOTAL	177

The most common reason for reoperation through 7 years in primary augmentation patients was because of capsular contracture (47 of 177 reoperations). In Allergan’s Core Study, there were 455 surgical procedures performed during 177 reoperations involving 130 primary augmentation patients.

Table 4
Main Reasons for Reoperation in Revision-Augmentation Patients Through 7 Years

Reason for Reoperation	n
Capsular Contracture	20
Hematoma/Seroma	13
Implant Malposition	11
Ptosis	9
Need for Biopsy	8
Scarring/Hypertrophic Scarring	7
Suspected Rupture	5
Asymmetry	3
Breast Cancer	3
Infection	3
Patient Request for Style/Size Change	3
Nipple Complications	3
Delayed Wound Healing	2
Wrinkling	2
Breast Pain, Breast Tissue Contour Deformity, Device Injury, Implant Palpability/Visibility, Implant Extrusion	1 each
TOTAL	97


The most common reason for reoperation through 7 years in revision-augmentation patients was also because of capsular contracture (20 of 97 reoperations). In Allergan’s Core Study, there were 267 surgical procedures performed during 97 reoperations involving 56 revision-augmentation patients.

3.6 What Are the Main Reasons for Implant Removal?

The main reasons for implant removal observed in primary augmentation and revision-augmentation women through 7 years are presented in Figures 3 and 4, respectively.


For primary augmentation, there were 107 implants removed in 58 patients. Of these 107 implants, 85 were replaced. The most common reason for implant removal was capsular contracture (35 of the 107 implants removed).

Figure 3. Main Reasons for Implant Removal Through 7 Years Primary Augmentation (n=107)



For revision-augmentation there were 59 implants removed in 32 patients. Of these 59 implants, 51 were replaced. The most common reason for implant removal was also due to capsular contracture (18 of the 59 implants removed).

Figure 4. Main Reasons for Implant Removal Through 7 Years Revision-Augmentation (n=59)



3.7 What Are Other Clinical Data Findings?

Below is a summary of clinical findings from the Allergan Core Study with regard to connective tissue disease (CTD), CTD signs and symptoms, cancer, lactation complications, reproduction complications, and suicide. These issues, along with others, are being further evaluated as part of an Allergan postapproval study of a large number of patients followed through 10 years (Breast Implant Follow-Up Study, or BIFS).

CTD DIAGNOSES

Four primary augmentation patients (0.9%) were reported to have a new diagnosis of CTD; 2 with rheumatoid arthritis at 7 months and at 3 years, respectively, after implantation in the Allergan Core Study and 2 patients with fibromyalgia at 4.5 years after implantation for both. One revision-augmentation patient (0.7%) was reported to have a new diagnosis of fibromyalgia at 10 months after implantation. It cannot be concluded that these CTD diagnoses were caused by the implants because there was no comparison group of similar women without implants.

CTD SIGNS AND SYMPTOMS

In Allergan's Core Study, numerous signs and symptoms were collected at 2, 4, and 6 years post-implant. For primary augmentation patients at 6 years after implantation, statistically significant increases were found for the symptom category of Joint (includes joint pain, stiff in morning, swelling in other joints, and swelling of hands), Muscular (includes back pain, muscle pain, aches, or cramps, muscle weakness, neck pain, paralysis of arms or legs), Gastrointestinal (includes constipation, diarrhea, gastrointestinal pain, heartburn, loss of appetite, stomach pain or cramps, and vomiting), Neurological (includes headaches, loss of balance, numbness/tingle of arms or legs, problems with memory, problems with thinking, and ringing in ears), Urinary (includes problems with urination and urinating too often), and Fibromyalgia (includes back pain, fatigue, neck pain, pain, and pain in the chest). No significant increases were found in the categories of General, Global, Pain, Skin, Fatigue, and Other symptoms. For revision-augmentation patients at 6 years after implantation, no statistically significant increases were found in any of the symptom categories.

The Core Study was not designed to evaluate cause-and-effect associations because there is no comparison group of women without implants, and because other contributing factors, such as medications and lifestyle/exercise, were not studied. Therefore, it cannot be determined whether this increase was due to the implants or not, based on the Core Study. However, you should be aware that you may experience an increase in these symptoms after receiving breast implants.

CANCER

There was 1 primary augmentation patient with a new diagnosis of breast cancer through 7 years in the Allergan Core Study. There was a 13% benign breast disease rate and a 1% malignant breast disease rate through 7 years. In primary augmentation patients there was 1 report of thyroid cancer and 1 report of brain cancer.

For revision-augmentation patients, there was 1 patient with a new diagnosis of breast cancer. There was a 15% benign breast disease rate and a 1% malignant breast disease rate through 7 years. There were no reports of other cancers, such as respiratory or cervical/vulvar, in revision-augmentation patients.

LACTATION COMPLICATIONS

Sixteen (21%) of the 75 primary augmentation patients who attempted to breastfeed following breast implantation in the Allergan Core Study through 7 years experienced difficulty with breastfeeding. The most common difficulty was inadequate milk production. For the 19 revision-augmentation patients who attempted to breastfeed after receiving breast implants, 6 (32%) had difficulty breastfeeding, 5 due to inadequate milk production and 1 due to pain.

REPRODUCTION COMPLICATIONS

Twenty-nine (6%) of the primary augmentation patients in the Allergan Core Study reported a reproduction problem through 7 years, most commonly miscarriage. For the 5 (3%) revision-augmentation patients who experienced a reproduction problem through 7 years, the most common problem was infertility.

SUICIDE

There was 1 report of suicide in the primary augmentation patients and 2 reports of suicide in the revision-augmentation patients in the Allergan Core Study through 7 years.

4. SURGICAL CONSIDERATIONS FOR BREAST AUGMENTATION

4.1 What Are the Alternatives to Breast Augmentation with Silicone Gel-Filled Breast Implants?

For primary augmentation patients, alternatives may include:

- Accepting your breasts as they are and having no surgery.
- Wearing a padded bra or external prostheses.
- Having mastopexy surgery (breast lift) without an implant.
- Having surgery with saline implants.

For revision-augmentation patients, alternatives may include:

- No revision.
- Removal with or without replacement.

4.2 What Are Questions to Consider When Choosing a Surgeon?

When choosing a surgeon who is experienced with breast augmentation, you should find out the answers to the following types of questions:

- How many breast augmentation implantation procedures does he/she perform per year?
- How many years has he/she performed breast augmentation procedures?
- Has he/she completed Allergan's Physician Certification Program for the use of its silicone-filled breast implants?
- Is he/she board certified, and if so, with which board?
- In which state(s) is he/she licensed to practice surgery? (Note that some states provide information on disciplinary action and malpractice claims/settlements to prospective patients, either by request or on the internet.)
- What is the most common complication he/she encounters with breast augmentation?

- What is his/her reoperation rate with breast augmentation, and what is the most common type of reoperation he/she performs?
- Can he/she perform this surgery in a hospital, as well as in the surgeon's independent surgery center? (Note that hospitals require evidence of appropriate training in specific procedures before allowing surgeons to operate in their facilities.)

4.3 What are Other Choices and Options Associated with the Surgery?

There are two approved types of breast implant fillers, saline and silicone, which gives more options to you in terms of the type of implant to achieve the effect you desire. Your surgeon can discuss these options with you and may make recommendations to you based upon the physical contours of your body. The focus of this brochure is silicone-filled breast implants; a separate brochure is available for saline-filled implants. Carefully review the section on complications and the section on Allergan's clinical study so that you may make an informed choice. Be sure to ask your surgeon to see and touch samples of both silicone and saline breast implants.

IMPLANT SHAPE AND SIZE

Depending on the desired shape you wish to achieve, you and your surgeon have implants with different profiles, or styles, from which to choose. Generally, the larger you want your cup size, the larger the breast implant the surgeon will consider (measured in cubic centimeters, or cc's, not in cup sizes, because cup size depends on the size and shape of the individual woman's chest).

Your surgeon will also evaluate your existing breast and skin tissue to determine if you have enough tissue to cover the breast implant you are considering, or, in some cases such as after pregnancy, too much extra skin. If you desire a breast implant size that is too large for your tissue, the surgeon may warn you that breast implant edges may be visible or **palpable postoperatively**. Also, excessively large breast implants may speed up the effects of gravity on the breast, and can result in droop or sag at an earlier age. A recent report indicates that larger-sized implants (greater than 350 cc) may be too large for many women, increasing the risk of developing complications

such as implant extrusion, hematoma, infection, palpable implant folds, and visible skin wrinkling requiring surgical intervention to correct these complications.⁷

SURFACE TEXTURING


Some studies suggest that surface texturing reduces the chance of severe capsular contracture,⁹ while other studies do not.^{7,8} Allergan's Core Study did not show a difference in the likelihood of developing capsular contracture with textured implants compared to smooth implants.

A textured implant may require a larger incision because the rougher textured surface may make it harder to place into the pocket without undue stress, which might damage the implant or decrease its durability.


IMPLANT PLACEMENT

The breast implant can be placed either partially under the pectoralis major muscle (**submuscular**) or on top of the muscle and under the breast gland (**subglandular**). You should discuss with your surgeon the advantages and disadvantages of the implant placement selected for you, as described in the table below.


BREAST BEFORE & AFTER IMPLANTATION



Breast before augmentation



Breast after subglandular augmentation



Breast after submuscular augmentation


Comparison Between Submuscular Versus Subglandular Placement

Submuscular Placement	Subglandular Placement
Surgery may be longer	Surgery may be shorter
Recovery may be longer	Recovery may be shorter
May be more painful	May be less painful
Reoperation may be more difficult	May provide easier access for reoperation
Less visible and palpable implants	More visible and palpable implants
Less likelihood of capsular contracture	Greater likelihood of capsular contracture
Easier imaging during mammography exam	More difficult imaging during mammography exam
May be preferable if you have thin or weakened breast tissue	May not be recommended if you have thin or weakened breast tissue

INCISION SITES

You should discuss with your surgeon the pros and cons for the incision site specifically recommended for you.

The incision size will be larger than for a saline breast augmentation. There are 3 common incision sites: under the arm (**axillary**), around the nipple (**periareolar**), or within the breast fold (**inframammary**).



- **Periareolar** - This incision is typically more concealed, but since it also involves cutting through the breast tissue, it is associated with a higher likelihood of breastfeeding difficulties, as compared to the other incision sites. Cutting through the tissue may make a change in sensation or infection more of a concern.
- **Inframammary** - This incision is generally less concealed than periareolar and associated with fewer breastfeeding difficulties than the periareolar incision site. It is also the most commonly used incision site at the present time, and is felt to give the best access to and control of the breast implant pocket.
- **Transaxillary** - This incision is less concealed than periareolar and associated with fewer breastfeeding difficulties than the periareolar incision site. If the incision is made under the arm, the surgeon may use a probe fitted with a miniature camera, along with minimally invasive (very small) instruments, to create a “pocket” for the breast implant. This approach is more difficult, and may increase the risk of damage to, and unexpected location of, the implant.
- **Umbilical (belly button)** - This incision site has not been studied in Allergan’s Core Study and should not be used for a wide variety of reasons, including potential damage to the implant shell.

ADDITIONAL PROCEDURES AT THE TIME OF BREAST AUGMENTATION

Your surgeon will examine your breasts and help you make decisions to obtain the best result in your individual situation. In some cases, particularly after pregnancy or significant weight loss, implants alone may not address all of the issues, such as sagging or extra skin, affecting your breasts. This is particularly true when there is extra skin remaining from when the breasts were engorged with milk, or when you might have been carrying more weight.

In these situations, your surgeon may recommend a breast lift (mastopexy) to remove some of the extra skin, or to lift the breasts, at the time of implant placement. Mastopexy involves removing a strip of skin from under the breast or around the

nipple to lift the nipple and breast location, and tighten the skin over the breast. Your surgeon will discuss the potential risks, and the location of the additional scars which might be required to lift your breasts or to remove the extra skin.

IMPLANT PALPABILITY

Implants may be more palpable or noticeable if there is an insufficient amount of skin/tissue available to cover the implant and/or when the implant is placed subglandularly.

SURGICAL SETTING AND ANESTHESIA

Augmentation surgery is usually performed on an outpatient basis, in a specialized operating room which may be located in a hospital, a surgery center, or surgical suite in the surgeon's office. General anesthesia is commonly used, and local anesthesia with sedation is also an option. You should be sure to check with your surgeon and with the facility where the surgery will take place, to become aware of the tests, presurgical examinations, and length of time you need to be without food or your routine medications prior to the surgical procedure.

POSTOPERATIVE CARE

You will probably feel somewhat tired and sore for several days following the operation, and your breasts may remain swollen and sensitive to physical contact for a month or longer. You may also experience a feeling of tightness in the breast area as your skin adjusts to your new breast size. The feeling in the breasts and nipple area also may be diminished during this time of swelling and immediate post surgery recovery. Other possible complications are described in the Breast Implant Complications section.

Postoperative care depends on each patient's situation, may involve the use of a special postoperative bra, compression bandage, or jog bra for extra support and positioning while you heal. Some surgeons may not want you to wear a bra at all for a period of time following the surgery. At your surgeon's recommendation, you will most likely be able to return to work within a few days, although for at least a couple of weeks you should avoid any strenuous activities that could raise your pulse and blood pressure, or require strenuous use of your arms and chest. Your surgeon may also recommend breast massage exercises.

Note: If you experience fever, do not feel well, or see noticeable swelling and/or redness or drainage in your implanted breast(s), you should contact your surgeon immediately.

OTHER FACTORS TO CONSIDER IN REVISION-AUGMENTATION SURGERY

Some revision surgeries require removal of an intact implant (for example, **capsulotomy** and pocket adjustments), while others do not require removal of the implant. Any device that has been removed during revision surgery should not be reimplanted. Allergan breast implants are “for single use only.”

4.4 What Follow-Up Examinations are Important?

BREAST SELF-EXAMINATIONS

Following breast augmentation, you should continue to perform a breast self-examination monthly. This may be more difficult with a breast implant in place. To continue to perform a monthly breast self examination efficiently, you should ask your surgeon to help you identify the difference between the implant and your breast tissue. Being able to identify the implant from breast tissue will decrease the necessity of excessive squeezing of the implant during examination. Any new lumps should be evaluated with a biopsy, as appropriate. If a biopsy is performed, be sure to inform the medical professional performing the biopsy that you have breast implants so that care will be taken to avoid injuring the implant.

SCREENING FOR SILENT RUPTURE

Because most ruptures of silicone-filled breast implants are silent, in most cases neither you nor your surgeon will be able to find evidence of rupture. Therefore, evaluation of your implants is needed to screen for implant rupture. The best method of screening is currently MRI at a center with a breast coil, with a magnet of at least 1.5 Tesla. The MRI should be read by a radiologist who is familiar with looking for implant rupture.

It is recommended that your first MRI evaluation take place starting at 3 years after implant surgery and then every 2 years, thereafter, even if you are experiencing no problems with your implant. If signs of rupture are seen on MRI, then you should have your implant removed, with or without replacement. Your doctor should assist you in locating a radiology/screening center, as well as a radiologist who is familiar with the technique and equipment for proper MRI screening for silent rupture of your breast implant.

Symptomatic Rupture

Symptoms associated with rupture may include hard knots or lumps surrounding the implant or in the armpit, loss of size of the breast or implant, pain, tingling, swelling, numbness, burning, or hardening of the breast. If you notice any of these changes, see your plastic surgeon so that he or she can examine the implants and determine whether you need to have an MRI examination to find out if your symptoms are due to rupture of the implant. If rupture has occurred, you should have your implant removed. Consult with your doctor regarding this and any other medical decisions related to your implants. More information on rupture is provided in section 3 of this brochure.

MAMMOGRAPHY

The current recommendations for getting screening/preoperative mammograms are no different for women with breast implants than for those without implants. Mammography exams should be interpreted by radiologists experienced in the evaluation of women with breast implants. It is essential that you tell your mammography technologist before the procedure that you have an implant. You should request a diagnostic mammogram rather than a screening mammogram, because more pictures are taken with diagnostic mammography. The technologist can use special techniques to reduce the possibility of rupture and to get the best possible views of the breast tissue.

5. ADDITIONAL INFORMATION

5.1 What Types of **NATRELLE**[®] Silicone-Filled Breast Implants Are Available from Allergan?

NATRELLE[®] Silicone-Filled Breast Implants come in a variety of profiles and sizes with either a textured shell or smooth surface shell. Your plastic surgeon will discuss with you the implant design that will best help you achieve the result that is right for you.

Examples of **NATRELLE**[®] Smooth and Textured Implant Styles




Smooth Implant



Textured Implant

The following diagram may help you to understand the projections of implants as your surgeon discusses the various options with you.

SILICONE-FILLED BREAST IMPLANT MATRIX



A=Width, B=Projection



Approved Allergan Implant Styles

Style Number	Breast Implant Description	Size Range
Style 10	Smooth shell surface, moderate profile	120cc–800cc
Style 15	Smooth shell surface, moderate-plus profile	155cc–752cc
Style 20	Smooth shell surface, high profile	120cc–800cc
Style 40	Smooth shell surface, moderate profile	80cc–560cc
Style 45	Smooth shell surface, extra-high profile	120cc–800cc
Style 110	BIOCELL® textured shell surface, moderate profile	90cc–510cc
Style 115	BIOCELL® textured shell surface, moderate-plus profile	150cc–716cc
Style 120	BIOCELL® textured shell surface, high profile	180cc–650cc

5.2 What If I Experience a Problem?

You will be given a device identification card with the style and serial number of your breast implant(s). This card is for your permanent record and should be kept in a safe place. In the event you have a concern or problem with your implant, you can use this card to describe the implant to your health care provider or to Allergan.

You should immediately report any problems that you notice with your implants to your plastic surgeon. If you believe that you have experienced a serious problem(s) related to your breast implants, you should have your health professional report the problem(s) to the Food and Drug Administration (FDA) and/or to Allergan. You may also report any serious problem directly through the FDA's MedWatch voluntary reporting system. An adverse event is considered serious and should be reported when it results in a hospitalization, disability, congenital problem with your child, or other medical or surgical intervention. The information reported to MedWatch is entered into databases to be used to follow safety trends (patterns) of a device and to determine whether further follow-up of any potential safety issues related to the device is needed.

To report, use MedWatch form 3500, which may be obtained through FDA's website at <http://www.fda.gov/medwatch/index.html>. You may also call 1.888.463.INFO.FDA (1.888.463.6332), from 10am-4pm Eastern Time, Monday through Friday, to receive an additional FDA MedWatch Package. Keep a copy of the MedWatch form completed by your surgeon for your records.

5.3 What Is Device Tracking?

Silicone gel-filled breast implants are subject to Device Tracking by Federal regulation. This means that your physician will be required to submit to Allergan the serial number of the implant(s) you receive, the date of surgery, information relating to the physician's practice and information on the patient receiving the implant(s). Your surgeon will write this information on the **Device Tracking Form** supplied by Allergan with each silicone-filled breast implant. Your surgeon will return the top portion of the form to Allergan following surgery. The bottom portion of the form will be provided to you following surgery. You have the right to remove your personal information from Allergan's Device Tracking program. If you choose NOT to participate in Device Tracking, please check the appropriate box on the Device Tracking form and return to Allergan. You also have the right to have your personal information withheld from disclosure to third parties who may request information from Allergan, such as the FDA. If you choose to participate in the Device Tracking program but do NOT want your personal information to be released to third parties, please also check the appropriate box.

Allergan strongly recommends that all patients receiving **NATRELLE**[®] Silicone-Filled Breast Implants participate in Allergan's Device Tracking program. This will help ensure that Allergan has a record of each patient's contact information so that all patients, including you, can be contacted in the case of a recall or other problems with your implants that you should be made aware of.

You are encouraged to complete the Device Tracking Form you received following surgery and return it to Allergan in the postage paid business reply envelope provided. Please inform Allergan whenever your contact information changes.

ASSESSMENT OF INFORMATION EFFECTIVENESS

The “Required Information” section of the Device Tracking Form also has a question designed to assess the effectiveness of the *Breast Augmentation with NATRELLE® Silicone-Filled Breast Implants* patient planner provided prior to your surgery. This question asks you to verify that you received and had adequate time to review this patient labeling information. Please check either yes or no. When the Required Information section is complete, return it to Allergan in the postage-paid business reply envelope provided.

5.4 What are the ConfidencePlus® Limited Warranties?

The ConfidencePlus® Limited Warranties provide lifetime replacement and limited financial reimbursement in the event of shell leakage or breakage resulting in implant rupture, subject to certain conditions as fully discussed in the ConfidencePlus® literature. Allergan offers two levels of coverage under its warranty program. Our standard ConfidencePlus® Limited Warranty program applies automatically to every Allergan breast implant recipient subject to the conditions discussed in the ConfidencePlus® literature. The optional ConfidencePlus® Premier Limited Warranty program is available for a low enrollment fee and increases the financial benefit in the event of implant rupture, subject to the conditions discussed in the ConfidencePlus® literature. For more information, please visit www.cppwarranty.com or contact Allergan’s Product Support Department at 1.800.362.4426.

5.5 How Can I Receive More Information?

Upon request, you will be provided with a copy of the package insert (Directions for Use; NATRELLE® Silicone-Filled Breast Implants document). You can request a copy from your surgeon or from Allergan. It can also be found on www.NATRELLE.com. The package insert has many undefined medical and technical terms because it contains information directed only to the surgeon.

For more detailed information on the preclinical and clinical studies conducted by Allergan, you are referred to the Summary of Safety and Effectiveness Data (SSED) for this product which may be accessed at <http://www.fda.gov/cdrh/breastimplants>.

If, after reading this information, you have additional questions about breast implants or breast implant surgery, there are a number of resources available to you.

TOLL-FREE NUMBER

If you are a patient or a prospective patient and wish to speak to an Allergan Breast Implant Support Specialist to inquire about breast implants, discuss any concerns, or request a copy of the patient labeling or package insert (Directions for Use), call toll free at 1.800.362.4426 (7 am to 5 pm Pacific Time).

ADDITIONAL RESOURCES

Allergan

1.800.624.4261

www.NATRELLE.com

www.allergan.com

www.breastimplantanswers.com

Institute of Medicine Report on the Safety of Silicone Implants

www.nap.edu/catalog/9618.html

Food and Drug Administration

1.888.INFO.FDA or 1.240.276.3103

www.fda.gov/cdrh/breastimplants

FOR FURTHER READING AND INFORMATION

Overall Safety Assessment

1. Bondurant, S., Ernster, V., and Herdman, R., Eds. 2000. Safety of silicone breast implants. Committee on the Safety of Silicone Breast Implants, Division of Health Promotion and Disease Prevention, Institute of Medicine. Washington, D.C.: *National Academy Press*.

Implant Rupture

2. Hedén, P., et al. 2006. Prevalence of rupture in Inamed silicone breast implants. *Plast. Reconstr. Surg.* 118:303-8.
3. Hölmich, L.R., et al. 2005. The diagnosis of silicone breast implant rupture. Clinical findings compared to findings at MRI. *Ann. Plast. Surg.* 54(6):583-9.
4. Hölmich, L.R., et al. 2005. The diagnosis of breast implant rupture: MRI findings compared to findings at explantation. 2005. *Eur. J. Radiol.* 53:213-25.
5. Hölmich, L.R., et al. 2004. Untreated silicone breast implant rupture. *Plast. Reconstr. Surg.* 114:204-14.
6. Hölmich, L.R., et al. 2001. Prevalence of silicone breast implant rupture among Danish women. *Plast. Reconstr. Surg.* 108(4):848-58.

Capsular Contracture

7. Henriksen, T.F., et al. 2005. Surgical intervention and capsular contracture after breast augmentation: a prospective study of risk factors. *Ann. Plast. Surg.* 54(4):343-51.
8. Kulmala, I., et al. 2004. Local complications after cosmetic breast implant surgery in Finland. *Ann. Plast. Surg.* 53(5):413-9.
9. Seify, H., et al. 2005. Preliminary (3 years) experience with smooth wall silicone gel implants for primary breast augmentation. *Ann. Plast. Surg.* 54(3):231-5.

Connective Tissue Disease (CTD)

10. Brinton, L.A., et al. 2004. Risk of connective tissue disorders among breast implant patients. *Am. J. Epidemiol.* 160(7):619-27.

11. Brown, S.L., et al. 2001. Silicone gel breast implant rupture, extra capsular silicone, and health status in a population of women. *J. Rheumatol.* 28:996-1003.
12. Hölmich, L.R., et al. 2003. Self-reported diseases and symptoms by rupture status among unselected Danish women with cosmetic silicone breast implants. *Plast. Reconstr. Surg.* 111:723-32.
13. Janowsky, E.C., et al. 2000. Meta-analyses of the relation between silicone breast implants and the risk of connective-tissue diseases. *N. Engl. J. Med.* 342(11):781-90.
14. Lipworth, L., et al. 2004. Silicone breast implants and connective tissue disease: An updated review of the epidemiologic evidence. *Ann. Plast. Surg.* 52:598-601.
15. Tugwell, P., et al. 2001. Do silicone breast implants cause rheumatologic disorders? A systematic review for a court-appointed national science panel. *Arthritis Rheum.* 44(11):2477-84.
16. Weisman, M.H., et al. 1988. Connective-tissue disease following breast augmentation: A preliminary test of the human adjuvant tissue hypothesis. *Plast. Reconstr. Surg.* 82(4):626-30.
17. Williams, H.J., et al. 1997. Breast implants in patients with differentiated and undifferentiated connective tissue disease. *Arthritis Rheum.* 40(3):437-40.
18. Wolfe, F. and Anderson, J. 1999. Silicone filled breast implants and the risk of fibromyalgia and rheumatoid arthritis. *J. Rheumatol.* 26:2025-28.

CTD Signs and Symptoms

19. Berner, I., et al. 2002. Comparative examination of complaints of patients with breast-cancer with and without silicone implants. *Eur. J. Obstet. Gynecol. Reprod. Biol.* 102:61-6.
20. Breiting, V.B., et al. 2004. Long-term health status of Danish women with silicone breast implants. *Plast. Reconstr. Surg.* 114:217-26.
21. Fryzek, J.P., et al. 2001. Self-reported symptoms among women after cosmetic breast implant and breast reduction surgery. *Plast. Reconstr. Surg.* 107:206-13.
22. Kjølner, K., et al. 2004. Self-reported musculoskeletal symptoms among Danish women with cosmetic breast implants. *Ann. Plast. Surg.* 52(1):1-7.

Cancer

23. Brinton, L.A., et al. 2000. Breast cancer following augmentation mammoplasty (United States). *Cancer Causes Control*. 11(9):819-27.
24. Brinton, L.A., et al. 2001. Cancer risk at sites other than the breast following augmentation mammoplasty. *Ann. Epidemiol.* 11:248-56.
25. Bryant, H., and Brasher, P. 1995. Breast implants and breast cancer—reanalysis of a linkage study. *N. Engl. J. Med.* 332(23):1535-9.
26. Cook, L.S. 1997. Characteristics of women with and without breast augmentation. *JAMA*. 20:1612-7.
27. Deapen, D.M., et al. 1997. Are breast implants anticarcinogenic? A 14-year follow-up of the Los Angeles Study. *Plast. Reconstr. Surg.* 99:1346-53.
28. Deapen, D., et al. 2000. Breast cancer stage at diagnosis and survival among patients with prior breast implants. *Plast. Reconstr. Surg.* 105:535-40.
29. Fryzek, J.P., et al. 2000. Characteristics of women with cosmetic breast augmentation surgery compared with breast reduction surgery patients and women in the general population of Sweden. *Ann. Plast. Surg.* 45(4):349-56.
30. Herdman, R.C., et al. 2001. Silicone breast implants and cancer. *Cancer Invest.* 19(8):821-32.
31. Jakubietz, M.G., et al. 2004. Breast augmentation: Cancer concerns and mammography—A literature review. *Plast. Reconstr. Surg.* 113:117e-22e.
32. Kjølner K., et al. 2003. Characteristics of women with cosmetic breast implants compared with women with other types of cosmetic surgery and population-based controls in Denmark. *Ann. Plast. Surg.* 50(1):6-12.
33. McLaughlin, J.K. and Lipworth, L. 2004. Brain cancer and cosmetic breast implants: A review of the epidemiological evidence. *Ann. Plast. Surg.* 52(2):15-17.
34. Miglioretti, D.L., et al. 2004. Effect of breast augmentation on the accuracy of mammography and cancer characteristics. *JAMA*. 291(4):442-50.

35. Pukkala, E., et al. 2002. Incidence of breast and other cancers among Finnish women with cosmetic breast implants, 1970-1999. *J. Long Term Eff. Med. Implants.* 12(4):271-9.

Suicide

36. Brinton, L.A., et al. 2001a. Mortality among augmentation mammo-
plasty patients. *Epidemiol.* 12(3):321-6.
37. Jacobsen, P.H., et al. 2004. Mortality and suicide among Danish
women with cosmetic breast implants. *Arch. Int. Med.*
164(22):2450-5.
38. Koot, V., et al. 2003. Total and cause specific mortality among
Swedish women with cosmetic breast implants: prospective study.
BMJ. 326(7388):527-8.
39. Pukkala, E., et al. 2003. Causes of death among Finnish women
with cosmetic breast implants, 1971-2001. *Ann. Plast. Surg.*
51(4):339-42.

Effects on Breastfeeding/Children

40. Hemminki, E., et al. 2004. Births and perinatal health of infants
among women who have had silicone breast implantation in Finland,
1967-2000. *Acta Obstet. Gynecol. Scand.* 83(12):1135-40.
41. Kjølner, K., et al. 2002. Health outcomes in offspring of Danish
mothers with cosmetic breast implants. *Ann. Plast. Surg.* 48:238-45.
42. Signorello, L.B., et al. 2001. Offspring health risk after cosmetic
breast implantation in Sweden. *Ann. Plast. Surg.* 46:279-86.

Silicone Gel Migration

43. Katzin, W.E., et al. 2005. Pathology of lymph nodes from patients
with breast implants: a histologic and spectroscopic evaluation.
Am. J. Surg. Pathol. 29(4):506-11.

Gel Bleed

44. Chandra, G., et al. 1987. A convenient and novel route to
bis(alkyne)platinum(0) and other platinum(0) complexes from Speier's
hydrosilylation catalyst. *Organometallics.* 6:191-2.
45. Flassbeck, D.B., et al. 2003. Determination of siloxanes, silicon,
and platinum in tissues of women with silicone gel-filled implants.
Anal. Bioanal. Chem. 375(3):356-62 (for example, data from
Patients B & C).

46. Lappert, M.F. and Scott, F.P.A. 1995. The reaction pathway from Speier's to Karstedt's hydrosilylation catalyst. *J. Organomet. Chem.* 492(2):C11-C13.
47. Lewis, L.N., et al. 1995. Mechanism of formation of platinum(0) complexes containing silicon-vinyl ligands. *Organometallics.* 14:2202-13.
48. Lugowski, S.J., et al. 2000. Analysis of silicon in human tissues with special reference to silicone breast implants. *J. Trace Elem. Med. Biol.* 14(1):31-42.
49. Stein, J., et al. 1999. In situ determination of the active catalyst in hydrosilylation reactions using highly reactive Pt(0) catalyst precursors. *J. Am. Chem. Soc.* 121(15):3693-703.



INDEX

Alternatives.....	47
Anesthesia	51
Breast Implant Follow-Up Study, or BIFS	15
Breast measurement	36
Breast reconstruction	11
Breast self-examination	53
Breast tissue atrophy	29
Breastfeeding	28, 32, 45
Calcium deposits.....	28
Cancer.....	30, 45
Capsular contracture.....	24
Chest wall deformity	29
Complications.....	19, 38
<i>ConfidencePlus</i> ®.....	60
Connective tissue disease (CTD).....	30
Contraindications	11
Core Study.....	35
CTD.....	44
CTD signs and symptoms	30, 44
Delayed wound healing	28
Device identification card.....	58
Device tracking	59
Effects on children	32
Extrusion	28
Food and Drug Administration (FDA)	58
Gel diffusion.....	32
Hematoma/seroma	27
Implant palpability	52
Implant placement	49
Implant removal.....	25, 42
Implant Shape and Size.....	48
Incision Sites.....	50
Indications	11

Infection.....	27
Inframammary.....	50
Lactation complications.....	45
Lymphadenopathy.....	29
Mammography.....	14, 55
Mastopexy.....	51
MedWatch.....	58
MRI.....	54
Necrosis.....	28
Neurological disease.....	31
Nipple and breast sensation.....	27
Pain.....	26
Patient satisfaction.....	37
Periareolar.....	50
Postoperative care.....	52
Precautions.....	12
Quality of life assessments.....	38
Reoperation.....	41
Reoperations.....	25
Reproduction complications.....	45
Rupture.....	14, 19
Rupture from literature.....	21
Rupture Information on Allergan Implants.....	20
Screening.....	54
Silent Rupture.....	54
Suicide.....	32, 46
Summary of Safety and Effectiveness Data (SSED).....	61
Surgeon.....	47
Surgical setting.....	52
Symptomatic Rupture.....	54
Transaxillary.....	51
Umbilical.....	51
Unsatisfactory Results.....	26

Preoperative Checklist

❑ PATIENT SELF-ASSESSMENT

Completed by the patient prior to surgery for discussion with the physician.

❑ ACCEPTANCE OF RISK AND SURGERY CONSENT

Reviewed and initialed by the patient and physician and retained in patient medical file.

❑ PATIENT SURGERY RECORD

Important pre- and post-operative appointments and related information recorded by the patient.

Postoperative Checklist

- ❑ **ALLERGAN DEVICE IDENTIFICATION CARD(S)**
Supplied following surgery and retained by patient in designated area of the Breast Augmentation Surgery Planner.
- ❑ **DEVICE TRACKING ENROLLMENT FORM**
Completed and returned by the patient to Allergan in the business reply envelope provided.
- ❑ **OPTIONAL CONFIDENCEPLUS® PREMIER WARRANTY ENROLLMENT FORM**
Completed and returned by the patient to Allergan in the business reply envelope provided.
- ❑ **INFORMATION FOR YOUR HEALTHCARE PROVIDERS**
Completed by the patient to give to her mammography center and primary care physician.



Patient Self-Assessment

Following your review of Section 1, Important Information for Women About Breast Augmentation with **NATRELLE**[®] Silicone-Filled Breast Implants, use this Patient Self-Assessment to evaluate your understanding of the information presented. Be sure to bring this breast surgery planner with the completed Patient Self-Assessment to your consultation with your doctor. He or she will review the assessment and use it to help guide additional discussion about the risks and benefits of surgery. There is additional space at the end of the Self-Assessment to make notes about the information or record specific questions that you would like to discuss with your surgeon.

Each of the following statements is clearly true or false. Indicate your answers by checking true or false. Your surgeon will review your answers with you.

If signs of rupture are seen on an MRI, you should have your implant removed.

- TRUE FALSE

Additional surgery to your breast and/or implant will be likely over the course of your life.

- TRUE FALSE

Your implants are not considered lifetime devices and you will likely undergo implant removal, with or without replacement, during your life.

- TRUE FALSE

You should inform your mammographers about the presence of your implants.

- TRUE FALSE

Your breast implants may interfere with your ability to successfully breastfeed.

- TRUE FALSE

You should perform breast self-examinations monthly and should make sure you know how to distinguish the implant from your breast tissue.

- TRUE FALSE

Silicone gel-filled breast implants have not been clinically tested in women with autoimmune diseases like lupus or scleroderma.

- TRUE FALSE

If you have serious health problems or conditions such as a weakened immune system or compromised blood supply to the breast, you should discuss with your surgeon whether breast augmentation surgery is appropriate for you.

- TRUE FALSE

To detect possible silent rupture, an MRI is recommended 3 years following initial surgery and every 2 years thereafter.

- TRUE FALSE

Although rare, there have been reports in the scientific literature providing evidence that the silicone gel fill may move beyond the fibrous capsule and into the breast tissue or away from the breast (gel migration), particularly if the scar capsule is ruptured, causing local complications such as pain and neuropathy.

- TRUE FALSE

Capsular contracture or hardening of the tissue surrounding the breast implant may result in the need for additional surgery.

- TRUE FALSE

Acceptance of Risk and Surgery Consent

Surgeon and patient initial each

	SURGEON	PATIENT
If signs of rupture are seen on an MRI, then you should have your implant removed.	_____	_____
Additional surgery to your breast and/or implant will be likely over the course of your life.	_____	_____
Your implants are not considered lifetime devices and you will likely undergo implant removal, with or without replacement, during your life.	_____	_____
You should inform your mammography technologist about the presence of your implants.	_____	_____
Your breast implants may interfere with your ability to successfully breastfeed.	_____	_____
You should perform breast self-examinations monthly and should make sure you know how to distinguish the implant from your breast tissue.	_____	_____
To monitor your breast implants for silent rupture, an MRI is recommended three (3) years following surgery and then every two (2) years after that.	_____	_____
The scar tissue or capsule that normally forms around the implant may tighten (contracture) and squeeze the implant, making your breast feel firmer and sometimes painful.	_____	_____
Allergan maintains a breast implant device tracking database and your participation in this database is strongly recommended.	_____	_____

Consent to Surgery

My surgeon has provided me with the AUGMENTATION SURGERY WITH SILICONE GEL-FILLED BREAST IMPLANTS PATIENT PLANNER to inform me prior to my surgery.

I have had adequate time to review and understand the information presented in the AUGMENTATION SURGERY WITH SILICONE GEL-FILLED BREAST IMPLANTS PATIENT PLANNER. My concerns and questions have been addressed by my doctor. I have considered alternatives to augmentation surgery, including use of external prostheses or surgery with saline-filled breast implants.

I am choosing to proceed with silicone gel-filled breast implant surgery.

Patient Name (Printed): _____

Patient Signature: _____

Date: _____

Surgeon Name (Printed): _____

Surgeon Signature: _____

Date: _____



Patient Surgery Record

Use this section to record important dates and contact information related to your breast surgery.

Pre-operative mammogram
baseline (if necessary): _____

Pre-operative appointment date: _____

Surgery date: _____

Surgery location: _____

Contact person at
surgery location: _____

Contact phone number: _____

First post-operative
appointment date: _____

Subsequent post-operative
appointment dates: _____

MRIs are recommended at 3 years following initial surgery and every 2 years thereafter to monitor implant integrity.

Month and year of my first scheduled MRI: _____

Post-surgery mammogram (6 months to 1 year following surgery): _____

Allergan Device Identification Card(s) Information

Record information from your Allergan device identification card(s) below and then place your card(s) in the pocket on the front cover of this planner for a permanent record.

Catalog Number: Left _____

Catalog Number: Right _____

Serial Number: Left _____

Serial Number: Right _____



Device Tracking Instructions

NATRELLE® Silicone-Filled Breast Implants are subject to device tracking per federal regulation. Because your implant is a tracked device, your healthcare provider is required to report certain information to Allergan following surgery, including implant-specific information like serial number and catalog number as well as the date of surgery, implanting surgeon's name and contact information and patient-specific information.



The device tracking form is a 2-part document and the top portion has been forwarded to Allergan by your surgeon. Upon receipt by Allergan of the healthcare provider portion of the device tracking form, you are entered in the device tracking database. The patient portion of the form was supplied to you following surgery in order to complete enrollment in Allergan's device tracking database. Your device tracking form should have stickers with information about your implant(s).

If the stickers are not already attached, use your Patient/Device ID card supplied following surgery to enter the serial number (SN) and catalog number (REF) of your implant(s). Allergan strongly encourages you to participate in device tracking. As information regarding the long-term safety of silicone gel-filled breast implants becomes available, your participation will allow us to provide this information to you. If you do *not* wish to participate in device tracking, check the indicated box and Allergan will remove your personal information from the database.

If you do wish to participate in the device tracking program but do *not* want Allergan to release your personal information to any third parties, such as the FDA, please check the appropriate box.

Finally, please indicate if you received this patient planner and had adequate time to review the information and consider your decision to proceed with breast surgery. Your answer here will help us evaluate the effectiveness of the communication program that Allergan has developed and make improvements if needed.


Place your completed device tracking form in the envelope provided and return it to Allergan.

I. Complete Upon Implant		
Device and Surgery Information		
DATE OF IMPLANTATION mm _____ /dd _____ /yy _____		
Affix Left Device Tracking Label OR Fill In The Device Data		
	(Left) REF _____ (Left) SN _____ <input type="checkbox"/> Reconstruction <input type="checkbox"/> Augmentation <input type="checkbox"/> Revision	
Affix Right Device Tracking Label OR Fill In The Device Data		
	(Right) REF _____ (Right) SN _____ <input type="checkbox"/> Reconstruction <input type="checkbox"/> Augmentation <input type="checkbox"/> Revision	
IMPLANTING/EXPLANTING PHYSICIAN INFORMATION		
LAST NAME _____	FIRST NAME _____	
ADDRESS _____	CITY, STATE/PROVINCE _____	ZIP/POSTAL CODE _____
E-MAIL _____	TELEPHONE _____	FAX _____
ATTENDING/FOLLOWING PHYSICIAN INFORMATION (if different from above)		
LAST NAME _____	FIRST NAME _____	
ADDRESS _____	CITY, STATE/PROVINCE _____	ZIP/POSTAL CODE _____
E-MAIL _____	TELEPHONE _____	FAX _____
PATIENT INFORMATION		
LAST NAME _____	FIRST NAME _____	
ADDRESS _____	CITY, STATE/PROVINCE _____	ZIP/POSTAL CODE _____
DATE OF BIRTH _____	SOCIAL SECURITY NUMBER _____	TELEPHONE _____
II. Complete if Device Discarded or Destroyed <input type="checkbox"/> N/A		
<input type="checkbox"/> Device(s) on this form were discarded/destroyed prior to completion of this procedure		
Date: mm _____ /dd _____ /yy _____ Reason/Comments: _____		
III. Complete if NATRELLE™ Silicone-Filled Breast Implants Were Removed <input type="checkbox"/> N/A		
Explanted Device Information		
Date of explant mm _____ /dd _____ /yy _____		
(Left) Serial # _____ <input type="checkbox"/> Unknown	(Right) Serial # _____ <input type="checkbox"/> Unknown	
(Left) Ref # _____ <input type="checkbox"/> Unknown	(Right) Ref # _____ <input type="checkbox"/> Unknown	
Reason for removal _____	Reason for removal _____	
Original implant date: mm _____ /dd _____ /yy _____ <input type="checkbox"/> Unknown	Original implant date: mm _____ /dd _____ /yy _____ <input type="checkbox"/> Unknown	
Original implanting physician _____ <input type="checkbox"/> Unknown	Original implanting physician _____ <input type="checkbox"/> Unknown	
COMPLETE AND RETURN THIS PAGE TO ALLERGAN IN THE ATTACHED ENVELOPE OR FAX TO 800.432.8803		


L030-06

For reference only





NO POSTAGE
NECESSARY
IF MAILED
IN THE
UNITED STATES



BUSINESS REPLY MAIL

FIRST-CLASS MAIL PERMIT NO. 1639 SANTA BARBARA, CA

POSTAGE WILL BE PAID BY ADDRESSEE

DEVICE TRACKING
ALLERGAN
71 S LOS CARNEROS RD
SANTA BARBARA CA 93117-9817



Optional *ConfidencePlus*[®] Premier Warranty Form

Confidence...it's more than a sense of well-being.

It's the peace of mind that comes with the knowledge your breast implants are covered by an industry leading warranty program. Allergan *ConfidencePlus*[®] breast implant limited warranty programs offer you coverage in the event of implant rupture, including product replacement and financial assistance to cover expenses not reimbursed by your insurance carrier.

Our standard *ConfidencePlus*[®] applies automatically to every **NATRELLE**[®] breast implant recipient and includes lifetime product replacement and up to \$1200 in financial assistance subject to the conditions discussed in the *ConfidencePlus*[®] literature.

The optional *ConfidencePlus*[®] Premier breast implant limited warranty provides all the peace of mind included with our standard *ConfidencePlus*[®] program, but increases the financial assistance to \$2400 and offers free contralateral implant replacement. For the low enrollment fee of \$100, you have access to lifetime product replacement, 10 years of coverage, the freedom to change styles or size as part of your replacement surgery, free contralateral implant replacement and up to \$2400 in financial assistance.¹

That's peace of mind... That's *ConfidencePlus*[®] Premier!

To enroll in our optional *ConfidencePlus*[®] Premier breast implant limited warranty program, use the information contained on your Allergan Device Identification Card(s) supplied to you after surgery or complete the purchase form that follows this page. Once complete, detach the form from this breast surgery planner and return to Allergan in the envelope provided. You may also fax your completed enrollment form with credit card information to 888.647.4029

Your purchase form and \$100 must be received or postmarked within 45 days of surgery and must accompany a certified check, money order, or valid credit card number in order to process your purchase. *Do not send a personal check.* *ConfidencePlus*[®] Premier enrollment forms that accompany a personal check will *not* be processed.

¹A charge may apply on product with a higher list price. The optional *ConfidencePlus*[®] Premier warranty is non-transferable and non-refundable. For complete program details see the *ConfidencePlus*[®] warranty program and terms at www.allergan.com or call Allergan at 1.800.624.4261.

CONFIDENCEPLUS® PREMIER LIMITED WARRANTY OPTION

Mail or fax this completed enrollment form along with your payment to:

Mail: Allergan Warranty Processing

71 S. Los Carneros Road, Goleta, CA 93117

Fax: 888.647.4029

Do not send a personal check. Enrollment forms that include a personal check **will not** be processed.

TO PAY BY CREDIT CARD:

Credit Card Type: Visa MasterCard American Express

Credit Card Number: _____

Expiration Date: _____

Cardholder Name (if other than Patient): _____

Cardholder Signature: _____

TO PAY BY CERTIFIED CHECK OR MONEY ORDER:

Please make payable to:

Allergan *ConfidencePlus*® Premier Limited Warranty

Name: _____

Address: _____

City: _____

State/Zip: _____

Social Security Number: _____

Date of Birth: _____

Home Phone: _____

Implanting Physician Name: _____

Physician Telephone: _____

Date of Surgery: _____

Implant Serial Number(s): Left Side _____


Implant Serial Number(s): Right Side _____

(The serial number follows the letters SN on your device identification card provided by your surgeon)


You will receive a certificate verifying your enrollment in the *ConfidencePlus*® Premier Warranty. Please allow 8–10 weeks for processing.

Note: Enrollment in the *ConfidencePlus*® Premier Warranty may be voided if the information provided is incorrect. The *ConfidencePlus*® Premier Warranty is non-refundable and non-transferable.

*For full details on the *ConfidencePlus*® Warranty and *ConfidencePlus*® Premier Warranty, refer to the *ConfidencePlus*® Warranty Program & Terms, on www.allergan.com, or available from Allergan at 1.800.624.4261.



NO POSTAGE
NECESSARY
IF MAILED
IN THE
UNITED STATES



BUSINESS REPLY MAIL

FIRST-CLASS MAIL

PERMIT NO. 1639

SANTA BARBARA, CA

POSTAGE WILL BE PAID BY ADDRESSEE

WARRANTY PROCESSING
ALLERGAN

71 S LOS CARNEROS RD
SANTA BARBARA CA 93117-9817



Information for the Mammography Center

Please update my patient file to reflect the presence of **NATRELLE®** Silicone-Filled Breast Implants. Since the examination of breasts augmented with breast implants requires more time, please allow additional time when scheduling my next mammogram and alert the physician and technologists performing the exam about the presence of my implants.

You may be aware that the FDA has approved Allergan's **NATRELLE®** Silicone-Filled Breast Implants for use in augmentation, reconstruction, and revision surgery. As part of a women's healthcare network, it is important that you are aware of the latest information on the safety of silicone gel-filled breast implants. For additional information, please consider the following resources:

Institute of Medicine Report on the Safety of Silicone Breast Implants
www.nap.edu/catalog/9618.html

Food and Drug Administration
www.fda.gov/cdrh/breastimplants

Breast Implant Safety
www.breastimplantsafety.org

PATIENT INSTRUCTIONS

Please record the catalog and serial numbers exactly as they appear on your Allergan Device Identification Card(s) before giving this page to your Mammography Center.

Location of implants
(submuscular or subglandular): _____

Catalog Number: Left _____

Catalog Number: Right _____

Serial Number: Left _____

Serial Number: Right _____



Information for Your Primary Care Physician

Your patient has been implanted with **NATRELLE®** Silicone-Filled Breast Implants. It is important that you include this information in her chart because while silicone-filled breast implants have been proven safe in thousands of patients worldwide, they can present additional challenges for attending physicians. So to ensure your patient receives the care she needs, when appropriate, please alert other physicians about the presence of her implants.

You may be aware that the FDA has approved Allergan's **NATRELLE®** Silicone-Filled Breast Implants for use in augmentation, reconstruction, and revision surgery. As part of a women's healthcare network, it is important that you are aware of the latest information on the safety of silicone gel-filled breast implants. For additional information, please consider the following resources:

Institute of Medicine Report on the Safety of Silicone Breast Implants
www.nap.edu/catalog/9618.html

Food and Drug Administration
www.fda.gov/cdrh/breastimplants

Breast Implant Safety
www.breastimplantsafety.org

PATIENT INSTRUCTIONS

Please record the catalog and serial numbers exactly as they appear on your Allergan device identification card(s) before giving this page to your primary care physician. If you have multiple primary care physicians, make copies of this form before providing to your physician.

Location of implants
(submuscular or subglandular): _____

Catalog Number: Left _____

Catalog Number: Right _____

Serial Number: Left _____

Serial Number: Right _____



Natrelle.com[®]

ALLERGAN
71 South Los Carneros Road, Santa Barbara, CA 93117

Tel: 1.800.624.4261
www.allergan.com

v021



ALLERGAN

© 2012 Allergan, Inc. ® and ™ marks owned by Allergan, Inc.

LDOC-03992 Rev.02