



ULTRASONIC SCALING

SWERV [®] Magnetostrictive Scaler	C1
Ultrasonic Inserts	C3
Magneto Reference Section	C8
Symmetry Piezo Tips	C11
Piezo Reference Section	C14





ULTRASONIC SET-UPS


POWER SCALING SET-UP 1

DESCRIPTION	SUGGESTIONS
IMS ULTRASONIC CASSETTE*	IM30HU_
100 Thin SWIVEL, 25k	UI25SS100
Swivel XT, 30k	UI30SSXT

POWER SCALING SET-UP 2

DESCRIPTION	SUGGESTIONS
Handpiece IMS Cassette	IM60HO
AF Swivel Direct Flow Set, 25K	UI25SDSET or UI30SD7

*See Section A for rails color coding options of IMS cassette.





SWERV³® MAGNETOSTRICTIVE SCALER

The SWERV³ magnetostrictive power scaler combines power and handling in a unique ergonomic design that ensures the comfort of dental professionals, and their patients. It integrates finely tuned electronics that improve control and scaling efficiency in all clinical indications, allowing dental professionals to provide optimal care for every patient, every day.

SWERV³ Series – 30K

UM3-30 UK

UM3-30 EU

SWERV³ Series – 25K

UM3-25 UK

UM3-25 EU




POINTS OF PERFORMANCE

- Finely tuned electronics provide superior patient and clinician comfort
- Switching power settings is effortless with the color-coded dual power modes
- The smooth surface makes disinfection simple and effective
- The color-coded, lighted display is easy to identify and adjust with the touch of a button
- The auto-purge function helps remove air from waterlines, saving the clinician time between patients



SWERV³ HANDPIECE

- A lightweight, well-balanced design reduces clinician fatigue.
- Large diameter handpiece reduces finger pinching for increased clinical comfort.
- Enables use of a relaxed grasp that eliminates undesired pressure points.
- Reduces drag and minimises strain on clinicians' hands.



REPLACEMENT PART CODES

SWERV ³ Barrier Sleeves	UMBS
SWERV ³ Foot Control, On-Off	UMFC
SWERV ³ Power Cord UK	UMCORD-UK
SWERV ³ Power Cord EU	UMCORD-EU
SWERV ³ Water Filters, 10-pack	UMWF
SWERV ³ Water Supply Connector	UMWSC
SWERV ³ Water Supply Set	UMWSS

SWERV³ 30K
Handpiece Assembly
UMHP30

SWERV³ 25K
Handpiece Assembly
UMHP25

MAGNETOSTRICTIVE INSERTS DESIGNS

(All Shown in 30 kHz size.)



SWIVEL DIRECT FLOW INSERTS



SWIVEL INSERTS



STREAMLINE INSERTS



ULTRASONIC INSERTS

SWIVEL DIRECT FLOW INSERTS



The Swivel Direct Flow inserts with a gray cone have all of the features of the revolutionary Swivel inserts plus a through-tip water delivery. This design focuses the water flow directly to the tip, reducing excess spray and increasing visibility.



25 kHz Part Code:
30 kHz Part Code:

#10 Universal
UI25SD10
UI30SD10

#100 Thin
UI25SD100
UI30SD100

#1000 Triple Bend
UI25SD1000
UI30SD1000

POWER LEVEL

LOW—HIGH

LOW—MEDIUM

LOW—HIGH

- Light to heavy deposits
- Orthodontic cement

- Light to moderate deposits even in periodontal pockets

- Light to heavy deposits
- Beveled edges for tenacious and heavy deposits
- Orthodontic cement

Hu-Friedy ultrasonic inserts are compatible with SWERV³ and all Cavitron* style magnetostrictive ultrasonic scalers.

*Cavitron is a registered trademark of Dentsply International Corporation.



SWIVEL DIRECT FLOW AFTER FIVE INSERTS



Ideal for:

- Light to moderate deposits even in periodontal pockets
- Furcation areas
- Developmental concavities



Left
UI25SDL
UI30SDL



Straight
UI25SDS
UI30SDS



Right
UI25SDR
UI30SDR

25 kHz Part Code:
 30 kHz Part Code:

POWER LEVEL

LOW—MEDIUM

LOW—MEDIUM

LOW—MEDIUM





ULTRASONIC INSERTS

SWIVEL INSERTS

In the Swivel inserts with a black cone, water flows from the base allowing thorough irrigation and lavage.



#10 Universal

UI25SS10
UI30SS10

#100 Thin

UI25SS100
UI30SS100

#1000 Triple Bend

UI25SS1000
UI30SS1000

Swivel XT

UI25SSXT
UI30SSXT

25 kHz Part Code:
30 kHz Part Code:

POWER LEVEL

LOW—HIGH

LOW—MEDIUM

LOW—MEDIUM

25 kHz: LOW
30 kHz: LOW—MEDIUM

- Light to heavy deposits
- Tenacious deposits
- Orthodontic cement

- Light to moderate deposits even in periodontal pockets

- Light to heavy deposits
- Beveled edges for tenacious and heavy deposits
- Orthodontic cement


- Tenacious deposits
- Furcation areas

EXPERIENCE FREEDOM AT YOUR FINGERTIPS!

The maneuverability and adaptability you want in ultrasonic scaling, with the exclusive Swivel inserts. With a slight rolling motion of your fingertips, rotate the insert 360° to eliminate the old “start, stop, remove, and readjust” routine required with traditional inserts. Swivel inserts also reduce the torque from the handpiece cord, increasing overall balance for hassle-free scaling.

With an 18% larger grip,* finger pinching is reduced, resulting in less hand fatigue** and enhance tactile sensitivity. The lightly textured silicone grips provide a secure, comfortable grasp and are color-coded for easy tip identification.

360° SWIVEL TECHNOLOGY



Swivel XT is thin enough to provide access into hard-to-reach areas, yet strong enough to effectively remove tenacious deposits even when used on low to medium power.

* Compared to Hu-Friedy Streamline Inserts.

** Source: Gerwatowski, L.J., McFall, D.B., Stach, D.: Carpal Tunnel Syndrome; Risk Factors and Preventive Strategies for the Dental Hygienist. Journal of Dental Hygiene, February 1992.



STREAMLINE INSERTS

- Routine scaling is more efficient with our Streamline inserts.
- Inserts have internal water flow delivery – water comes directly from the base
- The color-coded resin handles help in instantly identifying the tip style



#10 Universal
UI1025K
UI1030K

25 kHz Part Code:
 30 kHz Part Code:



#100 Thin
UI25K100S
UI30K100S

POWER LEVEL

LOW—HIGH

LOW—MEDIUM

- Light to heavy and tenacious deposits
- Orthodontic cement

- Light to moderate deposits even in periodontal pockets



#3 Beavertail
UI325K
UI330K

25 kHz Part Code:
 30 kHz Part Code:



#1000 Triple Bend
UI100025K
UI100030K

POWER LEVEL

MEDIUM—HIGH


LOW—HIGH

- Heavy deposits
- Orthodontic cement

- Light to heavy deposits
- Beveled edges for tenacious and heavy deposits
- Orthodontic cement



REFERENCE SECTION



MAGNETOSTRICTIVE INSERT CLEANING AND CARE

CLEANING

To clean, rinse thoroughly or, if you choose, fully immerse in a mild ultrasonic cleaning solution (pH 7-10.5) like Hu-Friedy's IMS Enzymax or Enzymax Earth cleaning solutions. Place in ultrasonic cleaning unit for 7-10 minutes or 16-20 minutes if the insert is in a cassette. Thoroughly rinse the insert with tap water and dry completely prior to sterilisation.

CARE AND MAINTENANCE

To increase the o-ring life of Hu-Friedy ultrasonic inserts and to prevent water flow from the tip after release of the foot control, take the following steps before placing the insert into the handpiece: Hold the handpiece upright over the sink, activate the foot control and flush the water line for 2 minutes. Fill entire handpiece with water. Gently twist the insert into the handpiece until completely seated.

Swivel and Swivel Direct Flow Inserts:

To maintain a smooth swivel, activate the foot switch and rotate insert 2-3 times.

All Thin Inserts (#100 Swivel, #100 Swivel Direct Flow, Swivel XT 30k and After Five): Use on LOW TO MEDIUM power settings only. Swivel XT 25 kHz should be used on low power only.

WARNING: Do not expose resin or silicone handled ultrasonic inserts to solutions containing phenols, gluteraldehydes, iodophors or to dry heat sterilisation.

Use steam sterilisation (autoclave) only for Swivel and Swivel Direct Flow inserts. Detailed care and maintenance instructions are available upon request.

STERILISATION

Sterilise Hu-Friedy ultrasonic inserts in an instrument cassette (i.e. IMS instrument cassette) or paper pouch in a steam autoclave. Chemical vapor sterilisation is unacceptable for Swivel inserts. Rapid heat transfer sterilisation and cold liquid/chemical sterilants are not recommended and may void warranty. Avoid contact with corrosive chemicals such as chlorides, sulfates and surface disinfectants. Do not expose to solutions containing phenols.

OPTIMUM EFFICIENCY


Scaling efficiency can significantly diminish with worn insert tips. Inferior performance, patient discomfort, and poor water delivery can result from worn, damaged, bent, or altered tips.

For optimal scaling efficiency, we recommend checking inserts monthly for signs of wear.

Replacement of Worn Inserts Actual Size Guide



All thin inserts (Swivel #100, Swivel XT, Swivel Direct Flow #100) should be used on LOW TO MEDIUM power settings only.




IMPORTANT NOTE: Hu-Friedy and Dentsply are competitors. This chart is not intended to suggest that the respective Dentsply and Hu-Friedy ultrasonic inserts are equivalent or comparable. Instead, this chart suggests an alternative Hu-Friedy brand ultrasonic insert based on the experience and recommendations of clinicians who have previously used the listed Dentsply inserts. Cavitron and Dentsply are registered trademarks of Dentsply International.

HU-FRIEDY ULTRASONIC INSERTS

Tip Style	Handle Color	25 kHz	30 kHz
SWIVEL			
#10 Universal	lavender	UI25SS10	UI30SS10
#100 Thin	black	UI25SS100	UI30SS100
#1000 Triple Bend	orange	UI25SS1000	UI30SS1000
Swivel XT	purple	UI25SSXT	UI30SSXT
SWIVEL DIRECT FLOW			
#10 Universal	lavender	UI25SD10	UI30SD10
#100 Thin	black	UI25SD100	UI30SD100
#1000 Triple Bend	orange	UI25SD1000	UI30SD1000
After Five Straight	dark blue	UI25SDS	UI30SDS
After Five Right	red	UI25SDR	UI30SDR
After Five Left	teal	UI25SDL	UI30SDL
STREAMLINE			
#10 Universal	lavender	UI1025K	UI1030K
#100 Thin	black	UI25K100S	UI30K100S
#1000 Triple Bend	orange	UI100025K	UI100030K
#3 Beavertail	yellow	UI325K	UI330K

DENTSPLY CAVITRON INSERTS**

Tip Style	25 kHz	30 kHz
TFI #10 Straight		
	grey	black
TFI #1000 Straight Triple		
	grey	black
CTI-10S THINsert		
		purple
FSI-10		
	dark blue	dark blue
FSI-SLI 10S		
	green	green
FSI-1000		
	dark blue	dark blue
FSI-SLI 10S		
	green	green
FSI-SLI 10R		
	green	green
FSI-SLI 10L		
	green	green
TFI #10 Straight		
	grey	black
TFI #1000 Straight Triple		
	grey	black
Bend TFI #3 Thin Beavertail		
	grey	black



Clinicians who use FSI® and Focused Spray® and FSI-SLI Focused Spray Slimline inserts may consider Swivel Direct Flow as an alternative.**



REFERENCE SECTION

POWER DISPERSION

Most Powerful

#1 point of tip
#2 concave surface
#3 convex surface
#4 lateral surfaces

Least Powerful

Insert Surfaces

After Five inserts are 40% thinner than the Original Prophy designs

After Five inserts

#10 Universal

MAGNETOSTRICTIVE INSERT APPLICATION CHART

	#10 Universal	#100 Thin	#1000 Triple Bend	Swivel XT	#3 Beaver-tail	After Five Inserts (Left, Straight, Right)	Original Prophy
Light to moderate deposits	•	•	•				•
Heavy deposits Orthodontic cement	•		•		•		•
Light to moderate deposits in periodontal pockets		•		•		•	
Tenacious subgingival deposits	•		•	•			•
Furcation areas				•		•	
Developmental concavities						•	



By varying the surface of the magnetostrictive insert tip on the tooth you can increase deposit removal efficiency (concave and convex surfaces) or increase patient comfort (lateral surfaces).



SYMMETRY PIEZO TIPS:

E-SERIES TIPS, S-SERIES TIPS

Hu-Friedy Symmetry tips are compatible with: E-Series tips are compatible with all EMS generators. S-Series tips are compatible with Symmetry IQ®, Satelec/Acteon and NSK piezo generators.

SCALING TIPS



S10 Universal
UES10
USS10

Standard diameter tip used to remove heavy supragingival deposits.

POWER LEVEL MED-HIGH



P Universal
UEP

S-Series: not available
Standard diameter tip used to remove moderate to heavy deposits located supragingivally and in shallow pockets.

MED-HIGH



P10 Universal
E-Series: not available
USP10

Standard diameter tip used to remove light deposits and plaque. Tip features smooth edges for safe removal of soft and hard deposits.

POWER LEVEL MED



100 Universal
UE100
US100

Probe-sized tip used to remove light deposits in shallow to moderate pockets and for periodontal maintenance.

MED-HIGH

AFTER FIVE

These probe sized tips are designed to enhance access into periodontal pockets and improve adaptation in interproximal areas, allowing the clinician to more thoroughly debride difficult to reach areas to achieve better clinical results.



S10 Universal
USAFS

POWER LEVEL MED



Right
USAFR

MED



Left
USAFL

MED





PERIO



Universal Perio

E-Series: not available

USU

Thin tip used to remove light deposits and plaque. Tip features smooth edges for safe removal of soft and hard deposits.

POWER LEVEL MED




Left Perio

E-Series: not available

US4L

Standard diameter contra-angle tip used to remove calculus on all surfaces of premolars and molars, especially difficult to reach interproximal surfaces.

POWER LEVEL LOW



Right Perio

E-Series: not available

US4R

Standard diameter contra-angle tip used to remove calculus on all surfaces of premolars and molars, especially difficult to reach interproximal surfaces.

LOW



3 Perio

UEP3

USP3

Thin-bladed curette-like tip used to remove heavy deposits and necrotic cementum on anteriors and premolars.

LOW

DIAMOND-COATED



Left Diamond

UE2LD

US2LD

Thin, contra-angle diamond-coated finishing instrument used for removal of residual debris and designed to enhance access to deep periodontal pockets and furcations.

POWER LEVEL LOW



Right Diamond

UE2RD

US2RD

Thin, contra-angle diamond-coated finishing instrument used for removal of residual debris and designed to enhance access to deep periodontal pockets and furcations.

LOW



Straight Diamond

UE1D

US1D

Thin, straight diamond-coated finishing instrument used for removal of residual debris and designed to enhance access to deep periodontal pockets primarily in anterior areas.

LOW

*Satelec and NSK are trademarks of Acteon, Inc., and Nakanishi Inc. respectively.



SYMMETRY HANDPIECE AND ACCESSORIES

GUARDIAN WRENCH/TIP CARRIER

Symmetry Piezo Tips come in the Guardian Wrench Tip Carrier for safe and easy connection to the handpiece.

Wrenches also hold the tip during cleaning and sterilisation to prevent loss of or damage to the tip.



Scaling Wrench/
Tip Carrier



Scaling Wrench/
Tip Carrier



Perio Wrench/
Tip Carrier

POWER LEVEL	MEDIUM—HIGH	LOW—MEDIUM	LOW
	Use primarily on medium to high power (in orange mode or scaling mode if applicable).	Use primarily on low to medium power (in blue mode or scaling mode if applicable).	Use primarily on low power (in green mode or perio mode if applicable).

CASSETTES

HANDPIECE CASSETTE
DOUBLE SHOWN




Handpiece Cassette Single
IM30IQ8

Symmetry handpiece not included.
Dimensions: 1.5" x 8" x 1.25"
(38 mm x 203 mm x 32 mm)

Handpiece Cassette Double
IM60IQ8

Symmetry handpiece not included.
Dimensions: 3" x 8" x 1.25"
(76 mm x 203 mm x 32 mm)



Wrench Cassette (Blue)
IMI38

Accommodates three Symmetry Piezo
Tips in Guardian Wrenches.
Exterior Dimensions: 5.63" x 1.95" x 1.77"
(143 mm x 49.5 mm x 45 mm)

REPLACEMENT PART

Black Replacement Wrench	USRTW
Handpiece Lux	UPHPL
Handpiece Sheath	UPHPLS
Foot Control	UPFC
Water Tube (2 mm)	UPWT
Water Filter Set (2 mm)	UPWTS
Water Filter	UPWF

Water Filter Case	UPWFC
Irrigation Connector (120 V only)	UPIC
Bulb (set of 3 replacements)	UPB
Spanner Wrench 5 x 8 inch (2 pcs)	UPSW
Power Cord (AC 230 V) EU	UPCORD-EU
Power Cord (AC 230 V) UK	UPCORD-UK



PIEZO REFERENCE

CHARACTERISTICS OF TIP MOVEMENT

The smooth linear movement generated by piezoelectric systems creates a stroke pattern that is similar to hand instrumentation. Adapting the lateral surface of the piezo tip against the tooth surface (similar to a hand instrument) will increase clinical efficiency.



ESSENTIALS OF POWER DISPERSION



- Most of the energy is dispersed at the point of tip (1).
- The lateral surfaces of the tip (2) are moving in a linear pattern, similar to the stroke of a bladed hand instrument.
- Adaptation of this lateral surface provides the most efficient deposit removal.
- Adapt the terminal 2-3 mm of tip for optimal performance.
- When concave front or convex back surfaces (3) are adapted to the tooth surface, the instrument will chatter against the surface rather than effectively remove deposit.

PIEZO TIP ADAPTATION

- Adapt the lateral surface to the tooth for optimal performance.
- Maintain tip angulation near 0° – parallel to tooth surface.
- To minimise the risk of root damage, never adapt tip at a 90° angle.
- Use little or no lateral pressure and let the tip do the work for you.
- When concave front or convex back surfaces are adapted to the tooth surface, the instrument will chatter against the surface rather than effectively instrument.

SENSORY OBSERVATIONS

- Smooth vibration enhances tactile sense for instrumentation.
- The smooth motion of the tip is an excellent option for patients with root sensitivity.
- Piezo tip vibration is very smooth and quiet when adapted correctly on a lateral surface.
- Deposits are pulverised, resulting in a slurry of debris that is visible in the water at the back of the mouth.
- Calculus is generally not removed in big chunks as commonly seen during magnetostrictive instrumentation.



Lateral surface of tip in use

To achieve the greatest efficiency and patient comfort, it is best to adapt the lateral sides of your piezoelectric instrument to the surface of the tooth.



PIEZO TIP CLEANING AND CARE

CLEANING & STERILISATION

- Piezo tips remain in the Wrench during the complete cleaning, disinfection and sterilisation reprocessing cycle.
- Clean using ultrasonic cleaner or thermodisinfectant.
- Steam sterilise for at least 4 minutes at 270°F/132°C or 30 minutes at 250°F/121°C. Do not heat above 275°F/135°C.
- Recommend at least 30 minute dry time after sterilisation cycle.
- Clean, disinfect, inspect and steam sterilise before each use following the respective equipment manufacturer's users instructions.
- Do not expose to phenols, iodophors or dry heat sterilisation.
- Detailed care and maintenance instructions are available upon request or go to Hu-Friedy.eu

TIP INSTALLATION

1. Your Hu-Friedy Piezo tip comes with its own wrench.
2. Attach tip to handpiece using the Wrench by aligning the threaded end of the tip to the handpiece.
3. Rotate the wrench clockwise. Once the instrument is screwed all the way in, give an extra quarter of a turn to obtain the required torque. Incomplete tightening may cause the tip to become loose or lack vibration.
4. When using the Wrench without a tip previously inserted into the tip carrier, take the following step: place the Wrench over the tip by slipping tip through the small hole in the bottom, closed end of the tool. This will position the tip within the Wrench for safety and convenience.

TIP REMOVAL

1. Hold the handpiece vertically. Insert the tip through the hole in the bottom of the Wrench.
2. Depress wrench slightly while rotating to fully secure tip to wrench.
3. Turn the wrench counter-clockwise to loosen and remove tip from handpiece.
4. Once removed from handpiece, do not remove the tip from the wrench carrier.


OPTIMUM EFFICIENCY


Scaling efficiency can significantly diminish with worn tips. Inferior performance and poor water delivery can result from worn, damaged, bent, or altered tips. For optimal scaling efficiency, we recommend checking tips regularly for signs of wear.

APPLICATION CHART

	S10	U	100	P10	P3	4R, 4L	1D	2RD, 2LD	P	AFS, AFL, AFR
Light to moderate supragingival deposits	●	●	●		●				●	
Heavy supragingival deposits	●			●					●	
Light to moderate deposits in shallow pockets		●	●		●		●		●	●
Light to moderate deposits in deep periodontal pockets					●		●	●		●
Moderate to heavy deposits in deep periodontal pockets						●				
Furcation areas								●		●

Replacement of Worn Tips Actual Size Guide





E-SERIES & S-SERIES POWER SETTINGS CROSS REFERENCE

HU-FRIEDY	APPLICATION				POWER P5	POWER SETTING			
	HU-FRIEDY	SATELEC	NSK			SATELEC P5	SATELEC SP NEUTRON	NSK VARIOS	EMS MINI-PIEZON
USS10	SCALING	SCALING	SCALING		MED	7-9	8-10	Min-3	
UES10	SCALING				MED				1-6
USP10	SCALING	SCALING	SCALING		MED	7-8	8-10	Min-3	
UEP	SCALING				MED				1-6
US100	SCALING	SCALING	SCALING		MED	4-6	6-8	Min-1	
UE100	PERIO				LOW				1-3
UE100M	PERIO				LOW				1-3
USU	PERIO	SCALING	PERIO		LOW	7-8	1-6	Min-4	
USP3	PERIO	PERIO	PERIO		LOW	1-3	1-4	Min-4	
UEP3	PERIO				LOW				1-3
US4R	PERIO	PERIO	PERIO		LOW	1-3	1-4	Min-4	
US4L	PERIO	PERIO	PERIO		LOW	1-3	1-4	Min-4	
US1D	PERIO	PERIO	PERIO		LOW	1-3	1-4	Min-4	
UE1D	PERIO				LOW				1-3
US2RD	PERIO	PERIO	PERIO		LOW	1-3	1-4	Min-4	
UE2RD	PERIO				LOW				1-3
US2LD	PERIO	PERIO	PERIO		LOW	1-3	1-4	Min-4	
UE2LD	PERIO				LOW				1-3

HU-FRIEDY E-SERIES & S-SERIES ITEM CODE CROSS REFERENCE

HU-FRIEDY	USS10	UES10	US100	UE100	UE100M	USP10	UEP	USU	USP3	UEP3	US4R	US4L	US1D	UE1D	US2RD	UE2RD	US2LD	UE2LD
SATELEC*	S1		10Z			10P		PFU	H3		H4R	H4L	H1		H2R		H2L	
NSK VARIOS*	G1		G6			G4		P20	P10						PD3		PD2	
EMS*		A		PS	PL3		P											