



Webinar #8 2020

Orthotic Management of Ulnar and/ or Median Nerve Dysfunction

Deborah A. Schwartz, OTD, OTR/L, CHT

Product and Educational Specialist, Physical Rehabilitation

Orfit Industries America

Debby.schwartz@orfit.com

1



Learning Objectives

At the conclusion of this session, participants will be able to:

- 1.** Recognize peripheral nerve dysfunction of the hand and fingers that benefit from orthotic management.
- 2.** Learn tips and tricks for working with low temperature thermoplastic materials that benefit specific orthotic fabrication.
- 3.** Identify the steps of fabrication for 2-4 custom orthoses for the thumb and fingers to address the above conditions.
- 4.** Understand the current levels of evidence to support these orthoses as therapeutic interventions.

© Orfit Industries
2018

2

Peripheral Nerve Dysfunction

Causes:

1. Infection or disease— polio, leprosy
2. Neurologic – Charcot-Marie-Tooth, spinal muscular atrophy
3. Congenital – absence of thenar muscles
4. Trauma – cervical spine, brachial plexus, lacerations
5. Compression or entrapment

© Orfit Industries
2018


3



The Ulnar Nerve

The Ulnar Nerve

- Muscles innervated
 - Flexor carpi ulnaris, flexor digitorum profundus, adductor pollicis, small digital muscles
- Motor functions
 - Finger adduction and abduction other than thumb; thumb adduction, flexion of digits 4 & 5; wrist flexion and adduction
- Sensory
 - Skin over medial surface of the hand through the superficial branch




Google Images

© Orfit Industries
2018

4

Ulnar Nerve Dysfunction



Google Images


Causes:

- Cubital tunnel syndrome
- Impact to the ulnar nerve at the medial epicondyle
- Excessive valgus stress at the elbow (throwing athletes)
- Compression by flexor carpi ulnaris
- Bony spurs at the olecranon and medial epicondyle
- Carpal bone dislocation
- Colles fracture or humeral fracture

© Orfit Industries
2018

5

Ulnar Nerve Dysfunction



Google Images

Low Nerve Injury

Loss of flexion of the proximal phalanges - paralysis of the interossei and other intrinsic muscles.

Clawing results from the extrinsic muscles hyperextending the proximal phalanges and from the pull FDP muscle, which contributes to poor grasp.

© Orfit Industries
2018

6


Ulnar Nerve Dysfunction

High Nerve Injury

FDP muscle is also without innervation
DIP joints are no longer flexed in
digits 4 and 5.

Milder appearing hand deformity.

Over time, both types (high and low) have
deformities that become fixed.



ANATOMY OF THE ULNAR NERVE

Google Images

© Orfit Industries
2018


7

Ulnar Nerve Dysfunction

Flattening of the normal arches of
the hand

Hyper-extension of MCP and
flexion in PIP and DIP of 4, 5th

Unable to ABD and ADD fingers



Froment's Sign:
Compensatory thumb
MP hyperextension and
flexion by FPL during
pinch



Claw hand deformity:
MCP joint hyper
flexion and PIP joint
flexion is caused by
loss of intrinsic
muscles to combat
force of extrinsic
flexors- imbalance
between extrinsic and
intrinsic muscle forces

Google Images

© Orfit Industries
2018

8

Ulnar Nerve Dysfunction

Functional Loss

Grip and pinch are affected
 Loss of finger abduction and adduction
 Inability to flex the 4th & 5th MCP joints while simultaneously extending the IP joints
 Decrease in grip strength
 Loss of fine prehension

Sensory Loss

Loss of sensation to the volar surface of the ulnar aspect of the palm distally and the volar surface of the small and ulnar half of the ring finger

© Orfit Industries
2018

9

Orthotic Management for Ulnar Clawing

Anti- Claw Orthosis

Goal :

Functional positioning of the digits in MCP joint flexion
 (loss of intrinsic muscle function)

Duration: as needed



© Orfit Industries
2018

10

Fabrication of the Ulnar Claw Orthosis

<https://www.youtube.com/watch?v=LhOvT3NjNr4>



© Orfit Industries
2018

11


Materials:
Orfit Strips 3.2mm x 12"
or
Orficast 6 cm- folded lengthwise
(either in half or tripled)



© Orfit Industries
2018

12

The Median Nerve



Google Images

© Orfit Industries
2018

13

The Median Nerve

Four small muscles of the hand are supplied the median nerve

Remember "LOAF":


- lateral two lumbricals
- opponens pollicis
- abductor pollicis brevis
- flexor pollicis brevis

Median Nerve

motor innervation

The median nerve innervates the LOAF muscles:

L = lateral two Lumbricals
O = Opponens pollicis
A = Abductor pollicis brevis
F = Flexor pollicis brevis



Google Images

© Orfit Industries
2018

14

The Median Nerve


Sites of compression:

Pronator syndrome – between the 2 heads of pronator teres

(tenderness in forearm, pain with repetitive pronation)

AIN compression – motor branch off median nerve (unable to make an ok sign)

Carpal tunnel syndrome (nocturnal symptoms, numbness and tingling, atrophy)



Google Images

© Orfit Industries
2018

15

Median Nerve Dysfunction

Differential Diagnosis:

The clinical evaluation of CTS vs Pronator Teres syndrome differs in the following ways:

Tinel sign is typically absent at the wrist, but may be positive over the proximal anterior forearm

Phalen's test is usually negative in pronator syndrome


Palpation demonstrates tenderness over the pronator teres and likely over the medial epicondyle

© Orfit Industries
2018

16

Median Nerve Dysfunction


High Nerve Injury
Above origin of AIN



Google Images

Deborah A. Schwartz

Low Nerve Injury
Thenar intrinsic muscles
paralyzed




© Orfit Industries
2018

17

Median Nerve Dysfunction

Symptoms include

- Weakness and /or loss of functional grip and pinch
- Atrophy of thenar muscles
- Sensory Loss- thumb, index finger, long finger, and the radial aspect of the ring finger



Google Images

© Orfit Industries
2018

18

Median Nerve Dysfunction

ADLs are affected including all fine motor tasks

- Difficulty grasping, pinching
- Problems with opening containers
- Pain with holding objects
- Decreased power grip
- Lumbrical muscles of index and middle finger are weak



Google Images

© Orfit Industries
2018

19

Orthotic Management

Recommended Orthoses:

Hand based thumb spica with MCP joint included

Hand based CMC joint orthosis w/out MCP joint



Goal of the Orthosis:

Positioning for maintaining the first web space to provide support, and improve ADL function

Duration: as needed

© Orfit Industries
2018

20



Orthotic Options for the Thumb



© Orfit Industries
2018


21

Fabrication of a Thumb Orthosis

Materials:

Orficast More 15 cm / 6"
Orfit Precuts

Sheet Materials:
Classic
Orflight
Orfit Colors NS



© Orfit Industries
2018

22

Fabrication of a Thumb Orthosis

Check out Orfit
videos on
[www.youtube.com!](http://www.youtube.com)

<https://www.youtube.com/watch?v=OvF18kafb1Y>

https://www.youtube.com/watch?v=7GrJJ2jVa_M

© Orfit Industries
2018

23

Additional Components of Orthotic Management

First Web Spacer-
when there is limited
opposition and
abduction



Functional Web
Spacer- when there is
limited opposition and
abduction



AIN-Pinch Assist
With limited FDP
and FPL

© Orfit Industries
2018

24

Combined Median Nerve/ Ulnar Nerve Dysfunction



Characteristics:

- Clawing of all four digits
- Adduction of thumb
- Inability to oppose and abduct the thumb
- Inability to extend PIP joints

© Orfit Industries
2018

25



Indications:

- Spinal Cord Injuries
- Charcot Marie Tooth
- ALS
- Trauma

© Orfit Industries
2018

26

Create a Functional Orthosis



Without thumb



With thumb

© Orfit Industries
2018

27

Create a Functional Orthosis




Dynamic thumb attachment




© Orfit Industries
2018

28



29



30

For Show

Thumb and Index
Figure of 8 Orthoses for
enhanced pinch
Orfit Strips 1/12"

Thumb Abduction
Strap
Orficast or Orficast More 2.5"

First Web Spacer
Orficast 2.5"

© Orfit Industries
2018

31

Evidence

Chan, R. K. (2002). Splinting for peripheral nerve injury in upper limb. *Hand Surgery*, 7(02), 251-259.

Choi, J. S., Mun, J. H., Lee, J. Y., Jeon, J. H., Jung, Y. J., Seo, C. H., & Jang, K. U. (2011). Effects of modified dynamic metacarpophalangeal joint flexion orthoses after hand burn. *Annals of rehabilitation medicine*, 35(6), 880.

Colditz, J. C. (2002). Splinting the hand with a peripheral nerve injury. In *Rehabilitation of the hand and upper extremity* (pp. 622-634). Mosby, Inc, St. Louis, MO.

Dauzère, Florence & Delclaux, S. & Pham, T.T. & Rongières, Michel & Mansat, Pierre. (2018). Combined Median and Ulnar Nerve Palsy Complicating Distal Radius Fractures. *Orthopaedics & Traumatology: Surgery & Research*. 104. 10.1016/j.otsr.2018.04.026.

Dell, P. C., & Sforzo, C. R. (2005). Ulnar intrinsic anatomy and dysfunction. *Journal of Hand Therapy*, 18(2), 198-207.

© Orfit Industries
2018

32

Evidence

Gajiwala, K. J., Sams, S. B., Pandya, N., & Wagh, A. (1991). A new dynamic lumbrical simulating splint for claw hand deformity. *Plastic and reconstructive surgery*, 87(1), 170-173.

McKee, P., & Rivard, A. (2004). Orthoses as enablers of occupation: client-centered splinting for better outcomes. *Canadian Journal of Occupational Therapy*, 71(5), 306-314.

Ruijs, A. C., Jaquet, J. B., Kalmijn, S., Giele, H., & Hovius, S. E. (2005). Median and ulnar nerve injuries: a meta-analysis of predictors of motor and sensory recovery after modern microsurgical nerve repair. *Plastic and reconstructive surgery*, 116(2), 484-494.

Seu, M., & Pasqualetto, M. (2012). Hand therapy for dysfunction of the intrinsic muscles. *Hand clinics*, 28(1), 87-100.

Sousa, G. G., & de Macêdo, M. P. (2015). Effects of a dynamic orthosis in an individual with claw deformity. *Journal of Hand Therapy*, 28(4), 425-428.

Watanabe, H., Ogata, K., Okabe, T., & Amano, T. (1978). Hand orthosis for various finger impairments—the KU finger splint. *Prosthetics and orthotics international*, 2(2), 95-100.

© Orfit Industries
2018


33

Tips for Increasing Your Client's Compliance with Orthotic Wear

- Make sure client understand purpose of orthosis
- Make sure client understands wearing schedule
- Make sure client has some say in final design (choice of color of material, straps)
- Have client keep a log of orthotic use
- Assess functional status with and without the orthosis
- ** Measure active and passive range of motion and/or perform a functional assessment prior to orthotic intervention.

© Orfit Industries
2018

34



Thank you for your attention!

Debby.schwartz@orfit.com

© Orfit Industries
2018

35



Join Orfit on Social Media!

Facebook	www.facebook.com/orfit.industries
Orfit Splinting & Rehabilitation Group	https://www.facebook.com/groups/Orfit.splinting
Twitter	@Orfit
LinkedIn	www.linkedin.com/companies/orfit-industries
YouTube	www.youtube.com/orfitindustries
Instagram	https://www.instagram.com/orfitindustries/
Blog	Blog.orfit.com

© Orfit Industries
2018

36



www.orfit.com

