**WHO Prequalification of Diagnostics Programme
PUBLIC REPORT**

**Product: Determine HIV-1/2¹
Number: PQDx 0033-013-00**

Determine HIV-1/2 with product codes **7D2342**, **7D2343**, **7D2343SET**, and **7D2343SETS** manufactured by **Abbott Diagnostics Medical Co., Ltd²**, **Rest-of-World regulatory version** (non-CE-marked regulatory version) was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 25 November 2011.

Summary of prequalification status for Determine HIV-1/2

	Date	Outcome
Prequalification listing	25-Nov-2011	listed
Dossier review	22-Nov-2011	MR
Site inspection (s) of quality management system	24-May-2019	MR
Product performance evaluation	11-Nov-2011	MR

MR: Meets Requirements

Report amendments and/or product changes

This public report has since been amended. Amendments may have arisen because of changes to the prequalified product for which WHO has been notified and has undertaken a review. Amendments to the report are summarized in the following table, and details of each amendment are provided below.

Version	Summary of amendment	Date of report amendment
2.0	Change of shelf life, extended from 14 months to 18 months	5-Jul-2016
3.0	Change of physical address of the Chiba Logistics Center to 483-2 Matsuhidai, Matsudo-shi, Chiba, 270-2214 Japan	18-Oct-2016
4.0	Change of both the HIV-1 and HIV-2 antibodies, and their reagent formulation, in the Control bar solution in order to address customer complaints regarding the weakness of control line intensity.	15-Feb-2017

¹ Product name was changed from Alere Determine HIV-1/2 to Determine HIV 1/2.

² Manufacturer name changed from Alere Medical Co.Ltd to Abbott Diagnostics Medical Co.Ltd.

5.0	Introduction of safety blood lancets as a component of whole blood assay SET. Addition of new configurations of its IVDs, where lancets currently provided with the IVDs are replaced with contact-activated safety lancets. Addition of a product code that corresponds with the new safety lancet (7D2343SETS (Alere Determine HIV-1/2)).	1-Mar-2017
6.0	Shelf life extension of the Chase Buffer (7D2243) from 18 months to 24 months. The Chase Buffer is included in the following SET products: PQDx 0033-013-00: Alere Determine HIV-1/2 SET (7D2343SET) (100 tests/kit set for whole blood assay). Alere Determine HIV-1/2 SETS (7D2343SETS) (100 tests/kit set for whole blood assay with safety lancet).	30-Jul-2018
7.0	Change of the name of the manufacturer from Alere Medical Co., Ltd. to Abbott Diagnostics Medical Co., Ltd. Along with this, product names were changed. Product name changed from Alere Determine HIV-1/2 to Determine HIV 1/2. Products themselves are completely the same.	29-Jan-2020

Intended use:³

According to the claim of Abbott Diagnostics Medical Co., Ltd, “*the DetermineHIV-1/2 is an In Vitro, visually read, qualitative immunoassay for the detection of antibodies to HIV-1 and HIV-2 in human serum, plasma or whole blood. The test is intended as an aid to detect antibodies to HIV-1/HIV-2 from infected individuals. The test is for professional use only*”.

Assay principle:

According to the claim of Abbott Diagnostics Medical Co., Ltd, “*Determine HIV-1/2 is an immunochromatographic test for the qualitative detection of antibodies to HIV-1 and HIV-2. Sample is added to the sample pad. As the sample migrates through the conjugate pad, it reconstitutes and mixes with the selenium colloid-antigen conjugate. This mixture continues to migrate through the solid phase to the immobilized recombinant antigens and synthetic peptides at the patient window site.*

If antibodies to HIV-1 and/or HIV-2 are present in the sample, the antibodies bind to the antigen-selenium colloid and to the antigen at the patient window, forming a red line at the patient window site.

³ This product is one that uses Protein A to detect human IgG antibodies. Protein A is also able to detect other classes of human antibody (IgA, IgD, IgE and IgM) but not as reliably as it does IgG. This product has been prequalified with respect to its ability to detect human IgG antibodies. Any claim to detect other types of antibodies on this kind of product has not been validated based on WHO prequalification requirements.

If antibodies to HIV-1 and/or HIV-2 are absent, the antigen-selenium colloid flows past the patient window, and no red line is formed at the patient window site.

To insure assay validity, a procedural control bar is incorporated in the assay device"

Test kit contents:

Component	20 tests (7D2342)	100 tests (7D2343)	100 tests (7D2343SET)	100 tests (7D2343SETS)
Test cards	2 cards of 10 tests/cards	10 cards of 10 tests/card	10 cards of 10 tests/card	10 cards of 10 tests/card
Chase buffer (7D2243)	N/A	N/A	1 bottle x 2.5ml	1 bottle x 2.5ml
EDTA capillary tubes (7D2222)	N/A	N/A	100	100
Blood lancets (7D2233)	N/A	N/A	100	N/A
Blood lancets, sterile (safety)	N/A	N/A	N/A	100

Items required but not provided within test kit:

Storage:

The test kit should be stored at 2-30 °C.

Shelf life upon manufacture:

18 months.

24 months for Chase Buffer.

Prioritization for prequalification

Based on the established eligibility criteria, Determine HIV-1/2 was given priority for WHO prequalification assessment.

Product dossier assessment

Abbott Diagnostics Medical Co., Ltd (formerly called Alere Medical Co.Ltd) submitted a product dossier for Determine HIV-1/2 as per the "Instructions for compilation of a product dossier" (PQDx_018 v1). The information submitted in the product dossier was reviewed by WHO staff and external experts (assessors) appointed by WHO.

Commitments for prequalification:

1. analytical performance studies
2. clinical performance studies
3. stability studies
4. a new version of the instructions for use.

WHO will follow up implementation of the commitments on the next re-inspection.

Based on the product dossier screening and assessment findings, the product dossier for Determine HIV-1/2 meets WHO prequalification requirements.

Manufacturing site inspection

A comprehensive inspection was performed at the site of manufacture of Determine HIV-1/2 manufactured by Abbott Diagnostics Medical Co., Ltd., (formerly called Alere Medical Co. Ltd) at 357 Matsuhidai Matsudo-shi, Chiba-ken 270-2214, Japan on 24 May 2019⁴. The inspection procedure is described in "Information for manufacturers on WHO prequalification inspection procedures for the sites of manufacture of diagnostics" (PQDx_014 v4).

The inspection found that Abbott Diagnostics Medical Co.Ltd had an established quality management system and manufacturing practices in place that should ensure the manufacture of a product of consistent quality.

Commitments for prequalification:

Abbott Diagnostics Medical Co., Ltd commit to provide WHO with the evidence of the improvement of its data capturing as follows:

- Abbott Diagnostics Medical Co., Ltd will change the HIV-1/2 scoring from current score (-, ±/- ~ 4/+) to new scoring (0-10) by 30 November 2019.
- Abbott Diagnostics Medical Co., Ltd will investigate digital data recording by 20 December 2019.

The manufacturer's responses to the nonconformities found at the time of the inspection were accepted on 24 May 2019. Implementation of commitments is under review.

Based on the site inspection and corrective action plan review, the quality management system for Determine HIV-1/2 meets WHO prequalification requirements.

⁴ Initial inspection took place on 28 September to 1 October 2010

Product performance evaluation

Determine HIV-1/2 was evaluated by WHO in the third quarter of 2011 at the Institute of Tropical Medicine, Belgium – a WHO Collaborating Centre for HIV/AIDS Diagnostics and Laboratory Support. The laboratory evaluation was conducted according to the “WHO Protocol for the laboratory evaluation of HIV serology assays” (PQDx_030 V1.0), and drew the following conclusions:

Determine HIV-1/2 is an immunochromatographic rapid diagnostic test for the detection of antibodies to HIV-1/2 in human serum, plasma, and whole blood. A volume of 50 μ l of serum, plasma or venous/capillary whole blood is required to perform the test procedure. This type of assay requires no sophisticated equipment and can therefore be performed in laboratories with limited facilities. Reading of the results is performed visually i.e. subjective reading.


In this limited performance evaluation using a panel of 1079 biological specimen, we observed an initial sensitivity (95% CI) of 100% (99.1% - 100%) and an initial specificity (95% CI) of 97.87% (96.4% - 98.8%) compared to the reference assays. The final sensitivity (95% CI) was 100% (99.1% - 100%) and the final specificity (95% CI) was 98.93% (97.8% - 99.6%) compared to the reference assays. In this study, 0.3% of the overall results were recorded as indeterminate. Results were interpreted independently by three technicians; the inter-reader variability was 1.4%. The invalid rate was 0.3% for initial testing and 0.1% for repeat testing.

Labelling

- 1. Labels**
- 2. Instructions for use**


1. Labels

1.1 7D2442 Label




Determine HIV-1/2 ROW 20 Test Pouch Label PN: 241563/R8	Pouch Size: 544mm(w) x 160mm(h) Label Artwork Size: 204mm(w) x 132mm(h)	PMS 303 PMS 2269 White
--	--	---

1.2 7D2343 Label


Determine HIV-1/2 ROW 100 Test Pouch Label	Pouch Size: 544mm(w) x 160mm(h)	PMS 303
PN: 241564/R8	Label Artwork Size: 204mm(w) x 132mm(h)	PMS 2269 White

1.3 7D2443SET


1.4 Chase buffer

Chase Buffer	 Black	PN: 1135904901
Box		
Size: (W)22mm x (L)44mm x (H)56mm		

1.5 EDTA Capillary tubes

EDTA Capillary Tubes	 Black
7D2222 Label	PN: 241609/R10
Size: 3.5" x 3"	Date of Last Revision: 2019/03/14

2.0 Instructions for use⁵

⁵ English version of the IFU was the one that was assessed by WHO. It is the responsibility of the manufacturer to ensure correct translation into other languages.

DETERMINE™ HIV-1/2

REF 7D2342, 7D2343

August 2019
241574/R12

EN

Key to symbols used	
REF Catalogue Number	Σ20 Contains Sufficient for 20 tests
IVD In Vitro Diagnostic Medical Device	Σ100 Contains Sufficient for 100 tests
2°C-30°C Store at 2-30°C	Do not reuse
	Keep away from sunlight
	Do not use if package is damaged

This package insert must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of assay results cannot be guaranteed if there are deviations from the instructions in this package insert.

NAME AND INTENDED USE

The Determine™ HIV-1/2 is an *In Vitro*, visually read, qualitative immunoassay for the detection of antibodies to HIV-1 and HIV-2 in human serum, plasma or whole blood. The test is intended as an aid to detect antibodies to HIV-1/HIV-2 from infected individuals. The test is for professional use only.

SUMMARY AND EXPLANATION OF THE TEST

AIDS (Acquired Immunodeficiency Syndrome) is characterized by changes in the population of T-cell lymphocytes. In an infected individual, the virus causes depletion of helper T-cells, which leaves the person susceptible to opportunistic infections and some malignancies. The virus that causes AIDS exists as two related types known as HIV-1 and HIV-2.

The presence of the AIDS virus elicits the production of specific antibodies to either HIV-1 or HIV-2.^{1,2,3}

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

Determine™ HIV-1/2 is an immunochromatographic test for the qualitative detection of antibodies to HIV-1 and HIV-2.

Sample is added to the sample pad. As the sample migrates through the conjugate pad, it reconstitutes and mixes with the selenium colloid-antigen conjugate. This mixture continues to migrate through the solid phase to the immobilized recombinant antigens and synthetic peptides at the patient window site.

If antibodies to HIV-1 and/or HIV-2 are present in the sample, the antibodies bind to the antigen-selenium colloid and to the antigen at the patient window, forming a red line at the patient window site. If antibodies to HIV-1 and/or HIV-2 are absent, the antigen-selenium colloid flows past the patient window, and no red line is formed at the patient window site.

To insure assay validity, a procedural control bar is incorporated in the assay device.

CONTENTS

Determine™ HIV-1/2 Serum/Plasma Assay, 20 Tests (7D2342) or 100 Tests (7D2343)

- Determine™ HIV-1/2 Test Card, 2 cards or 10 cards (10 tests/card), HIV-1/2 recombinant antigen and synthetic peptide coated.

ACCESSORIES (required but not provided)

For testing Whole Blood samples

- 1 Bottle (2.5 mL) Chase Buffer (7D2243) prepared in phosphate buffer. Preservatives: Antimicrobial Agents.

Whole Blood (fingerstick assay)

- EDTA Capillary Tubes (7D2222)

Materials Required But Not Provided

Disposable gloves, timing device

Micropipette capable of delivering 50 µL (other than fingerstick)

Alcohol swab, gauze pad, Lancet (for fingerstick)

WARNINGS AND PRECAUTIONS

For In Vitro Diagnostic Use.

CAUTION:

Appropriate biosafety practices^{4,5} should be used when handling specimens and reagents. These precautions include, but are not limited to the following:

- Wear gloves.
- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect all spills of specimens or reagents using a suitable disinfectant, such as 0.5% sodium hypochlorite.^{6,7}
- Decontaminate and dispose of all specimens, reagents, and other potentially contaminated materials in accordance with local regulations.^{8,9}

STORAGE

The Determine™ HIV-1/2 Test Cards and Chase Buffer must be stored at 2-30°C until expiration date.

Kit components are stable until expiration date when handled and stored as directed. Do not use kit components beyond expiration date.

Immediately reseal all unused tests in the foil pouch containing the desiccant by pressing seal from end to end to close.

SPECIMEN COLLECTION

Serum, Plasma, and Whole Blood Collection by Venipuncture

Human serum, plasma, and whole blood collected by venipuncture should be collected aseptically in such a way as to avoid hemolysis.

NOTE: For whole blood and plasma specimens, EDTA collection tubes must be used.

Whole Blood Collection by Fingerstick¹⁰

Before collecting a fingerstick specimen, place an EDTA capillary tube on a clean dry surface.

- Choose the fingertip of the middle, ring, or index finger (whichever is the least calloused) for adults and children older than one year. Warm the hand as needed with a warm, moist towel or warm water to increase blood flow.
- Clean fingertip with alcohol; allow to air dry. Position the hand palm-side up.
- Use a new lancet for each person. Place the lancet off-center on the fingertip. Firmly press the lancet against the finger and puncture the skin. Dispose of the lancet in an appropriate biohazard sharps container.
- Wipe away the first drop of blood with a sterile gauze pad.
- Hold the finger lower than the elbow and apply gentle, intermittent pressure to the base of the punctured finger several times. Touch the tip of the EDTA Capillary Tube to the drop of blood*. Avoid air bubbles.

*If EDTA Capillary Tubes (7D2222) will be used, fill the tube with blood between the 2 marked lines (50 µL).

SPECIMEN STORAGE

- Serum and plasma specimens should be stored at 2-8°C if the test is to be run within 7 days of collection. If testing is delayed more than 7 days, the specimen should be frozen (-20°C or colder).
- Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 7 days of collection. Do not freeze whole blood specimens.
- Whole blood collected by fingerstick should be tested immediately.

TEST PROCEDURE

The desired number of test units from the 10-test card can be removed by bending and tearing at the perforation.

NOTES:

- Removal of the test units should start from the right side of the test card to preserve the lot number which appears on the left side of the test card.
- Assay should be initiated within 2 hours after removing the protective foil cover from each test.

1. Remove the protective foil cover from each test.

2. For serum or plasma samples:

- Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
- Wait a minimum of 15 minutes (up to 60 minutes) and read result.

3. For whole blood (venipuncture) samples:

- Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
- Wait one minute, then apply one drop of Chase Buffer to the sample pad.
- Wait a minimum of 15 minutes (up to 60 minutes) and read result.

4. For whole blood (fingerstick) samples:

- Apply 50 µL of sample (by EDTA capillary tube) to the sample pad (marked by the arrow symbol).
- Wait until blood is absorbed into the sample pad, then apply one drop of Chase Buffer to the sample pad.
- Wait a minimum of 15 minutes (up to 60 minutes) and read result.


QUALITY CONTROL

To insure assay validity, a procedural control is incorporated in the device and is labeled "Control". If the control bar does not turn red by assay completion, the test result is invalid and the sample should be retested.

INTERPRETATION OF RESULTS


POSITIVE (Two Bars)

Red bars appear in both the control window (labeled "Control") and the patient window (labeled "Patient") of the strip. Any visible red bar in the patient window should be interpreted as positive.


NEGATIVE (One Bar)

One red bar appears in the control window of the strip (labeled "Control"), and no red bar appears in the patient window of the strip (labeled "Patient").

INVALID (No Bar)

If there is no red bar in the control window of the strip, and even if a red bar appears in the patient window of the strip, the result is invalid and should be repeated.

NOTES:

- The test result is positive even if the patient bar appears lighter or darker than the control bar.
- The control bar may exhibit a weak intensity for some patient samples, particularly those with high titer HIV.
- Upon appearance of a red bar in the control window, no matter how faint, the test result is considered valid.
- If an invalid test result occurs repeatedly, or for technical assistance, contact your local distributor or call Technical Support.

LIMITATIONS OF THE PROCEDURE

- The Determine™ HIV-1/2 test is designed to detect antibodies to HIV-1 and HIV-2 in human serum, plasma, and whole blood. Other body fluids or pooled specimens may not give accurate results.
- The intensity of the patient bar does not necessarily correlate to the titer of antibody in the specimen.
- A negative result with Determine™ HIV-1/2 does not exclude the possibility of infection with HIV. A false negative result can occur in the following circumstances:
 - low levels of antibody (e.g., early seroconversion specimens) are below the detection limit of the test
 - infection with a variant of the virus that is less detectable by the Determine™ HIV-1/2 assay configuration
 - HIV antibodies in the patient that do not react with specific antigens utilized in the assay configuration
 - specimen handling conditions which result in loss of HIV antibody multivalency
 - HIV-infected persons taking antiretroviral medication^{11,12,13}

For these reasons care should be taken in interpreting negative results. Other clinical data (e.g., symptoms or risk factors) should be used in conjunction with the test results.

- Positive specimens should be retested using another method and the results should be evaluated in light of the overall clinical evaluation before a diagnosis is made.
- Whole blood or plasma specimens containing anticoagulants other than EDTA may give incorrect results.
- Neonates of HIV-infected mothers may carry maternal antibodies to HIV for up to around eighteen months, which may not necessarily indicate the true infection status of the newborn.

PERFORMANCE CHARACTERISTICS

SPECIFICITY

A total of 1594 serum and plasma specimens from Asia, West Africa, and North America were tested by Determine™ HIV-1/2 and a commercially available test (Table I).

Table I
Specificity of Determine™ HIV-1/2

Population	Number of Specimens Tested	Negative by Determine™ HIV-1/2	Negative by a Commercially Available Test***
Seronegative			
Serum	908	907/908 (99.89%)	908/908 (100.00%)
Plasma	403	403/403 (100.00%)	403/403 (100.00%)
Pregnant Females	58*	57/57 (100.00%)	57/57 (100.00%)
West Africans	49	48/49 (97.96%)	48/49 (97.96%)
Disease States Other than HIV and Potentially Interfering Substances	176*	173/175 (98.86%)	174/175 (99.45%)
Total	1594**	1588/1592 (99.75%)	1590/1592 (99.87%)

* One specimen from a pregnant female and an HCV positive patient were positive by both Determine™ HIV-1/2 and the commercially available test. Both specimens confirmed positive by HIV-1 Western Blot.

** 456 specimens were from North America, 1089 specimens were from Asia, and 49 specimens were from Africa.

***The reference method of a commercially available test is particle agglutination.



A total of 3663 seronegative serum and plasma specimens from North America, Asia, and Africa were tested by Determine™ HIV-1/2 and commercially available tests (Table II). The specimens from North America, Asia, and 49 of 2118 specimens from Africa (referred to as "West Africans" in Table I) were included in Table I. Discordant specimens were confirmed negative by either Western blot or HIV-1 PCR assays.

Table II
A Comparison of Determine™ HIV-1/2 Specificity by Geographic Area

Area	Number of Specimens Tested	Negative by Determine™ HIV-1/2	Negative by Commercially Available Tests*
North America	456	451/454 (99.34%)	453/454 (99.78%)
Asia	1089	1089/1089 (100.00%)	1089/1089 (100.00%)
Africa	2118	2079/2118 (98.16%)	2100/2118 (99.15%)

* The reference methods of commercially available tests are particle agglutination, enzyme immunoassay and chemiluminescent immunoassay.

Whole Blood

A total of 368 seronegative whole blood specimens from Thailand were tested with paired serum and plasma by Determine™ HIV-1/2. Thirty-nine of the whole blood specimens were collected by both venipuncture and fingerstick (Table III).

Table III
A Comparison of Determine™ HIV-1/2 Specificity in Seronegative Whole Blood and Paired Serum and Plasma Specimens

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