remel

BactiDrop™

Voges-Proskauer B EN

REF R21562

INTENDED USE

BactiDrop Voges-Proskauer B is a reagent for use as part of a two-step qualitative procedure for the detection of the ability of a microorganism to produce acetylmethylcarbinol (acetoin) from glucose fermentation. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having bacterial infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

SUMMARY AND EXPLANATION

Voges and Proskauer found that certain bacteria, when grown in an appropriate medium, produced acetylmethylcarbinol.¹ The medium used for this determination is MR-VP Broth, a modification of Clark and Lubs medium.² Barritt recommended employing 5% alpha naphthol and 40% potassium hydroxide as the reagents used to detect the presence of acetylmethylcarbinol.³

PRINCIPLE

The Voges-Proskauer test is based on the detection of acetoin (acetylmethylcarbinol), a precursor to 2,3-butanediol production. One molecule of acetoin is formed by the decarboxylation of two molecules of pyruvic acid. In the presence of oxygen and alkali, acetoin and 2,3-butanediol are oxidized to diacetyl, which is the reactant for the pink color produced in the Voges-Proskauer test. 5% alpha naphthol acts as a color intensifier, increasing the sensitivity of the reaction. 40% potassium hydroxide aids in the absorption of CO₂ to hasten the oxidation of acetoin to diacetyl.

REAGENTS (CLASSICAL FORMULA)*

Potassium hydroxide (CAS 1310-58-3)..... 400.0 g
Demineralized Water (CAS 7732-18-5)..... 1000.0 ml

*Adjusted as required to meet performance standards.

PRECAUTIONS

DANGER! Causes eye and skin burns. Causes digestive and respiratory tract burns. Harmful if swallowed.

This product is for *In Vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after use. Directions should be read and followed carefully. Refer to Safety Data Sheet on company website for detailed information on reagent chemicals.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction, do not use the device.

STORAGE

This product is ready for use and no further preparation is necessary. Store product in its original container at 20-25°C until used. Do not freeze or overheat. Protect from light.

PRODUCT DETERIORATION

This product should not be used if (1) there is evidence of dehydration, (2) the color has changed, (3) the expiration date has passed, or (4) there are other signs of deterioration. The expiration date applies to the product in its intact container when stored as directed.

SPECIMEN COLLECTION, STORAGE, TRANSPORT

Specimens should be collected and handled following recommended guidelines.^{6,7}

MATERIALS PROVIDED

50 BactiDrop Voges-Proskauer B – glass ampule enclosed in a protective plastic sheath with a dropper dispenser tip (0.75/ampule)
1 Ampule crusher (reusable)

MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) MR-VP Broth (R061430), (7) BactiDrop Voges-Proskauer A (R21560).

PROCEDURE

Place dropper in the assembled, reusable ampule crusher provided. Hold the dropper/crusher in an upright position and lightly tap the bottom to dislodge any bubbles that may have formed. Grasp the middle of the dropper/crusher with the thumb and forefinger, and with the tip pointing away, press gently to crush the ampule. Invert dropper and squeeze slightly to dispense in a dropwise fashion. Discard remaining portion of partially used ampule at end of workday.

Test Procedure (Barritt Method³):

1. Dispense 18 drops (0.6 ml) of BactiDrop Voges-Proskauer A into 1 ml of inoculated MR-VP Broth that has been incubated at 35-37°C for 24-72 hours. Shake well.
2. Add 6 drops (0.2 ml) of BactiDrop Voges-Proskauer B. Shake well.
3. Allow tube to stand for 15 minutes. Observe for the development of a pink color.

INTERPRETATION

Positive Test - Pink color development
Negative Test - No color development

QUALITY CONTROL

All lot numbers of BactiDrop Voges-Proskauer B have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

CONTROL

Enterobacter cloacae ATCC® 13047
Escherichia coli ATCC® 25922

RESULTS

Positive*
Negative*

*After 24 hours incubation


LIMITATIONS

1. Occasionally a known positive organism may fail to give a positive test; it may be necessary to gently heat the culture containing the Voges-Proskauer (VP) reagents to effect a pink color development.⁵
2. The Methyl Red and Voges-Proskauer tests should not be used alone to identify the *Enterobacteriaceae* due to the possibility that acetoin may have been destroyed. Prolonged incubation should be avoided.⁵
3. The order of adding reagents must be followed exactly; reversal of the order may cause a weak positive or negative test.⁵
4. The exact amount of Voges-Proskauer B must be added; an excess of this reagent may mask a weak positive VP reaction by the development of a copper color. The development of a copper-like color is not considered a positive reaction.⁵

BIBLIOGRAPHY

1. Voges, O. and B. Proskauer. 1898. Zeit. Hyg. 28:20-32.
2. Clark, W.M. and H.A. Lubs. 1915. J. Infect. Dis. 17:160-173.
3. Barritt, M.M. 1936. J. Pathol. Bacteriol. 42:441-454.
4. Blazevic, D.J. and G.M. Ederer. 1975. Principles of Biochemical Tests in Diagnostic Microbiology. John Wiley & Sons, New York, NY.
5. MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams & Wilkins, Philadelphia, PA.
6. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, and D.W. Warnock. 2011. Manual of Clinical Microbiology. 10th ed. ASM Press, Washington, D.C.
7. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12th ed. Mosby Elsevier, St. Louis, MO.

Symbol Legend

	Catalog Number
	In Vitro Diagnostic Medical Device
	Batch Code (Lot Number)
	Consult Instructions for Use (IFU)
	Temperature Limitation (Storage Temp.)
	Use By (Expiration Date)
	Contains sufficient for <n> tests
	Do not use if packaging damaged
	Authorized European Representative
	Manufacturer
	UK Conformity Assessment
	European Conformity Assessment

Remel Inc
12076 Santa Fe Trail Drive
Lenexa, KS 66215, USA
www.thermofisher.com/microbiology
(800) 255-6730
International: (913) 888-0939

For technical information contact your local distributor.
©2022 Thermo Fisher Scientific Inc. All rights reserved.
BactiDrop is a trademark of Thermo Fisher Scientific and its subsidiaries.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
CAS (Chemical Abstracts Service Registry No.)

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Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA
©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.