



Aspira®
DRAINAGE SYSTEM

A NEW KIND OF FREEDOM

The Aspira® Drainage System is a compassionate home treatment option for end-stage cancer patients with malignant pleural effusions (MPE) and malignant ascites. Aspira allows patients to manage their symptoms at home, eliminating the need for frequent hospital visits.

PATIENT COMFORT

Gravity-based design produces a consistent pressure flow for minimized discomfort.

EASE-OF-USE

After connecting the drainage line to the catheter valve, one gentle squeeze will initiate drainage.

With the Aspira Drainage System, there is no risk of accidentally discharging vacuum. Vacuum is activated with one squeeze of the manual pump.


FASTER DRAINAGE

Faster drainage flow rate means more downtime for the patient.



DISTANCE LENGTH

The Aspira Drainage Catheter is 72 cm long. The Peritoneal Catheter has a space of 15 cm between the cuff to the proximal fenestration. The Pleural Catheter has a space of 11 cm between the cuff to the proximal fenestration.



With three placement method options (antegrade, retrograde and over-the-wire), and multiple kit options designed to accommodate customer and patient needs, the Aspira Drainage System gives you more placement flexibility than the PleurX[®] Drainage System.

ENHANCED TUNNELER DESIGN


- Malleable and pre-bent tunneler
- Tunneler is 37% longer than the PleurX tunneler for placement flexibility*



Aspira Drainage System Tunneler



Dilator




Splittable Sheath and Introducer



18G Introducer Needle

VALVED INTRODUCER

- Peel-away introducer sheath helps reduce fluid leakage and/or air aspiration during placement



REPLACEABLE VALVE

- Valve can be replaced without the need for catheter exchange



PRELOADED STIFFENING STYLET

- Stiffens catheter to facilitate insertion
- Stylet with Hydro-Glide™ coating designed to decrease fluid leakage during insertion



15.5F ASPIRA DRAINAGE CATHETER

- Distally trimmable to accommodate smaller cavities
- Proximally trimmable to enable valve replacement

0.038" x 80 cm
J-tip Guide Wire



Slide Clamp



Attachable
Suture Wing



Scalpel



Luer Adapter



Universal Tubing Adapter

*Data on file.

More consistent pressure during drainage.


The Aspira Drainage System is designed to drain quickly without the use of high pressure. Low, consistent vacuum pressure is preferable to high vacuum pressure that may cause patient discomfort.

Simulated testing shows the Aspira Drainage System:

- Has a 68% higher average flow rate than the PleurX Drainage system.^{1,2,4}
- Demonstrates consistent pressure and flow rate, unlike the PleurX Drainage System, which has a high initial vacuum that drops quickly.³
- Has more fenestrations⁴ and a larger valve¹ than the PleurX Drainage System.



Aspira Drainage System



72%
larger valve¹

Simulated Average Flow Rate Comparison²

Aspira Drainage System	PleurX Drainage System
7.46 mL/sec	4.45 mL/sec

68%
higher
average flow
rate²

PleurX Drainage System

Aspira Drainage System

67%
more
fenestrations⁴

1. Data on file. Internal diameters measured at the narrowest point. N=10 per device.
2. Data on file. Simulated water drainage. N=9 for Aspira Drainage Bag, N=8 for PleurX Drainage Bottle. Bench data may not be representative of clinical outcomes. Different test methods may yield different results.
3. Data on file. Simulated water drainage under similar conditions. N=40 for Aspira Drainage Bag, N=8 for PleurX Drainage Bottle. Bench data may not be representative of clinical outcomes. Different test methods may yield different results.
4. Data on file.

Aspira kit components



INSERTION TRAY COMPONENTS

- 1 – 15.5F Peritoneal or Pleural Drainage Catheter with Stylet
- 1 - #11 Futura® Safety Scalpel
- 1 – 18G Introducer Needle
- 1 – Guide Wire 0.038" x 80 cm J-tip
- 1 – Dilator 8F
- 1 – Dilator 12F
- 1 – 16.5F Splittable Sheath Introducer
- 1 – Tunneler
- 1 – Slide Clamp
- 1 – Aspira Valve Assembly
- 1 – Aspira Valve Cap
- 1 – Suture Wing
- 1 – Luer Adapter
- 1 – Universal Tubing Adapter
- 2 – Split Gauze 4" x 4" (10 cm x 10 cm)
- 6 – Gauze 4" x 4" (10 cm x 10 cm)
- 1 – Surgical Tape
- 1 – Adhesive Dressing



DRESSING KIT COMPONENTS

- 1 – Adhesive Dressing
- 2 – Split Gauze 4" x 4" (10 cm x 10 cm)
- 4 – Gauze 4" x 4" (10 cm x 10 cm)
- 1 – Gloves (1 Pair)
- 1 – Tape Strips (3 Per)
- 3 – Alcohol Wipe
- 1 – CSR Wrap



DRAINAGE KIT COMPONENTS

- 1 – Drainage Bag (1,000 mL)
- 2 – Alcohol Wipe
- 1 – Slide Clamp
- 1 – Tape Strips (3 Per)
- 1 – Aspira Valve Cap



LUER ADAPTER COMPONENTS

- 1 – Luer Adapter
- 1 – Aspira Valve Cap



UNIVERSAL TUBING ADAPTER COMPONENTS

- 1 – Universal Tubing Adapter
- 1 – Aspira Valve Cap



REPAIR KIT COMPONENTS

- 1 – Slide Clamp
- 1 – Aspira Valve Assembly
- 1 – Aspira Valve Cap

SIMPLE TO USE

A simple and easy-to-use design makes the Aspira Drainage System an excellent choice for caregivers and patients.



EASY OCCLUSION MANAGEMENT

A replaceable valve allows for catheter repair, rather than catheter exchange.



Aspira Drainage System

Peritoneal Product Code	Pleural Product Code	Description
4992206	4992506	Catheter Kit 1 Insertion Tray with Catheter
4992207	4992507	Complete Kit 1 Insertion Tray with Catheter, 5 Drainage Kits, 5 Dressing Kits, Voucher for 20 Drainage Kits and 20 Dressing Kits
4992208	4992508	Essential Kit 1 Insertion Tray with Catheter, 20 Drainage Kits, 20 Dressing Kits
4992209	4992509	Starter Kit 1 Insertion Tray with Catheter, 5 Drainage Kits, 5 Dressing Kits


Aspira Drainage Accessories**

Peritoneal Product Code	Description
4991503	Dressing Kit
4992301	Drainage Kit
4992305	Luer Adapter
4992306	Repair Kit
4992307	Universal Tubing Adapter

** Packaged 5 units per box

68%
less storage space*

Twenty Aspira drainage bags and dressings require approximately the same amount of space as one box of 4 PleurX bottles and dressings, reducing storage space and waste.



Visit myaspira.com for more information.

* Data on file. PleurX packaging measurements represent five boxes of PleurX Starter Kits containing four 1,000 mL drainage bottles and dressings per box. Aspira packaging measurements represent one box of an Aspira Essential Kit containing twenty 1,000 mL drainage bags and dressings per box.

All trademarks are the property of their respective owners.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

merit.com | myaspira.com

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE
Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733