



Transitions™
XTRActive®

TransitionsPRO.com

Photochromic performance and polarization are influenced
by temperature, UV exposure and lens material
©2020 Transitions Optical Ltd.

Frames by icl berlin

2020 TRANSITIONS® PRODUCT AVAILABILITY GUIDE



Reinforce the value of your patients' purchase with the Transitions® lenses Certificate of Authenticity!

	<i>New Technology, New Frontier of Performance</i>		<i>Our Best Lenses for Extra Protection, Even in the Car</i>	<i>Our Best Lenses for Glare Reduction Outdoors</i>	<i>Adaptive Sunglass Lenses for Behind-the-Wheel</i>
ICONIC LENS COLORS	<p><i>Gray, Brown or Graphite Green</i></p>	<p><i>Gray, Brown or Graphite Green</i></p>	<p><i>Gray, Brown or Graphite Green</i></p>	<p><i>Gray</i></p>	<p><i>Olive Green to Copper to Dark Red-Brown</i></p>
NEW COLORS	<p><i>Style Colors: Sapphire, Amethyst, Amber and Emerald</i></p>	<p><i>Style Colors: Sapphire, Amethyst, Amber and Emerald</i></p>	<p><i>Style Mirrors: Gold, Silver Shadow, Pink, Red, Green and Blue</i></p>		
INDOOR CLARITY	FULLY CLEAR	FULLY CLEAR	HINT OF TINT FOR PROTECTION	HINT OF TINT FOR PROTECTION	OLIVE GREEN <i>Not Recommended for Indoor Use</i>
OUTDOOR DARKNESS	DARK	DARK	DARKEST <i>even in high temperatures</i>	DARK	DARK
BEHIND-THE-WHEEL ACTIVATION	NO	NO	YES	NO	YES
POLARIZATION	NO	NO	NO	VARIABLE POLARIZATION	YES POLARIZED

Transitions, Transitions Signature, XTRActive and Vantage are registered trademarks and the Transitions logo and Transitions Light Intelligent Lenses are trademarks of Transitions Optical Inc., used under license by Transitions Optical Ltd. Light Under Control is a trademark of Transitions Optical Ltd. GEN 8 is a trademark of Transitions Optical Ltd. ©2020 Transitions Optical Ltd. Drivewear is a registered trademark of Younger Mfg. Co.

Photochromic performance and polarization are influenced by temperature, UV exposure and lens material.



2020 TRANSITIONS® PRODUCT AVAILABILITY GUIDE

Table of Contents (Click on category to go to page)

FSV 1.50	4	Signature® Bifocals 1.50	9
SFSV 1.50	4	Blended Bifocals 1.50 - 1.67	10
1.53	5	Progressives 1.50	11-13
1.59	6	1.53	14-15
1.60 - 1.67	7	1.59	16-20
1.67	8	1.60	20
		1.67	21-24

FSV 1.50

1.50 Standard Index

ESSILOR LENSES

Essilor Lens FSV Spheric

Transitions® XTRActive®

Gray

Brown

Graphite Green

Sphere

Cylinder to

●

●

●

GR

-4.00 to +4.00

-2.00

SFSV 1.50

1.50 Standard Index Single Vision

ESSILOR LENSES

Essilor Lens SFSV

Essilor Stylistic® Single Vision

Eyezen™+

Eyezen™ Start

HOYA VISION CARE

HOYA SFSV 1.50

HOYA SV iQ 1.50

NIKON OPTICAL

Nikon Infinite Single Vision

Nikon Smart Boost

SHAMIR INSIGHTS INC.

Shamir Attitude III Single Vision

Shamir Autograph® III Single Vision

Shamir Relax™ Single Vision

SIGNET ARMORLITE

Signet Armorlite SFSV

Signetek Processed SV

KODAK Digital SV

X-CEL OPTICAL

X-Cel SFSV

●

●

●

GR

BR

GG

-10.50 to +10.00

-4.00

GR

BR

GG

-3.00 to +3.00

-1.00

GR

BR

-10.00 to +6.00

-6.00

BR

-10.00 to +6.00

-6.00

GR

-8.00 to +9.00

-4.00

GR

-8.00 to +9.00

-4.00

GR

BR

-10.00 to + 8.00

-6.00

GR

BR

-10.00 to + 8.00

-6.00

GR

-6.50 to +6.00

-6.00

GR

-12.00 to +7.25

-6.00

GR

-12.25 to +6.00

-6.00

GR

-6.00 to +10.00

-4.00

GR

-8.00 to +6.00

-6.00

GR

-8.00 to +6.00

-6.00

GR

BR

-9.00 to +6.00

-4.00

TRANSITIONS® LIGHT INTELLIGENT LENSES™

		Transitions® XTRActive®			Sphere	Cylinder to
		Gray	Brown	Graphite Green		
1.53 TRIVEX® Single Vision		●	●	●		
ESSILOR LENSES						
	Essilor Lens SFSV TREXA®	GR	BR	GG	-10.00 to +6.00	-4.00
	Essilor Eyecode® Single Vision TREXA®	GR	BR	GG	-10.00 to +6.00	-6.00
	Essilor FIT® Single Vision TREXA®	GR	BR	GG	-10.00 to +6.00	-6.00
	Essilor Azio® Single Vision TREXA®	GR	BR	GG	-10.00 to +6.00	-6.00
	Essilor 360® Single Vision TREXA®	GR	BR	GG	-10.00 to +6.00	-6.00
HOYA VISION CARE						
	HOYA SFSV Phoenix®	GR	BR		-10.00 to +6.00	-4.00
	HOYA SV iQ Phoenix®	GR	BR		-10.00 to +6.00	-4.00
SEIKO OPTICAL PRODUCTS						
	Seiko PL 1.53/Trivex® Transitions®	GR			-10.00 to +7.00	-5.00
SHAMIR INSIGHTS INC.						
	Shamir Attitude III Single Vision	GR			-6.50 to +6.00	-6.50
	Shamir Autograph® III Single Vision	GR			-14.00 to +8.50	-6.50
	Shamir Relax® Single Vision	GR			-13.25 to +6.00	-6.50
SIGNET ARMORLITE						
	Signetek Processed SV	GR			-10.00 to +6.00	-6.00
	KODAK Digital SV	GR			-10.00 to +6.00	-6.00

	Transitions® XTRActive®			Sphere	Cylinder to
	Gray	Brown	Graphite Green		
1.59 Polycarbonate Single Vision					
CARL ZEISS VISION					
synchrony Aspheric SFSV	GR			-10.00 to +6.00	-4.00
synchrony Access HD (75, 125 & 175)	GR			-10.00 to +8.00	-6.00
synchrony SV HDC	GR			-10.00 to +8.00	-6.00
ZEISS 3D Single Vision	GR			-10.00 to +8.00	-4.00
ZEISS Individual® SV	GR			-10.00 to +8.00	-4.00
ZEISS Individual® SV Wrap	GR			-4.00 to +4.00	-4.00
ESSILOR LENSES					
Essilor SFSV Spherical Airwear®	GR •	BR •	GG	-10.00 to +8.00	-4.00
Essilor SFSV Aspheric Airwear®	GR	BR	GG	-10.00 to +8.00	-4.00
Essilor Eyecode® Aspheric SV Airwear®	GR •	BR •	GG	-10.00 to +8.00	-6.00
Essilor FIT® Aspheric SV Airwear®	GR •	BR •	GG	-10.00 to +8.00	-6.00
Essilor AZIO® Aspheric SV Airwear®	GR •	BR •	GG	-10.00 to +8.00	-6.00
Essilor SV 360® Aspheric Airwear®	GR •	BR •	GG	-10.00 to +8.00	-6.00
Essilor Stylistic® Single Vision Airwear®	GR •	BR •	GG	-4.00 to +4.00	-4.00
Eyezen®+ Airwear®	GR •	BR •	GG	-10.00 to +8.00	-6.00
Eyezen®+ Start Airwear®	GR	BR		-10.00 to +6.00	-6.00
HOYA VISION CARE					
HOYA SFSV Polycarbonate	GR			-6.00 to +6.00	-4.00
HOYA SV iQ Polycarbonate	GR			-6.00 to +6.00	-4.00
Sync III	GR			-13.00 to +5.25	-4.00
NIKON OPTICAL					
Nikon Infinite Single Vision	GR	BR		-10.00 to +8.00	-6.00
Nikon Smart Boost	GR	BR		-10.00 to +8.00	-6.00
SHAMIR INSIGHTS INC.					
Shamir Attitude III Single Vision	GR			-6.50 to +6.00	-7.00
Shamir Autograph® III Single Vision	GR			-14.00 to +8.75	-7.00
Shamir Relax® SFSV	GR			-14.00 to +8.00	-7.00
SIGNET ARMORLITE					
Signet Armorlite SFSV	GR			-8.00 to +10.00	-4.00
Signetek Processed SV	GR			-10.00 to +6.50	-6.00
KODAK Digital SV	GR			-10.00 to +6.50	-6.00



• Shows availability of *Transitions® XTRActive®* style mirrors in any of the six colors: Gold, Silver Shadow, Pink, Red, Green and Blue.
 All colors can be applied to either gray or brown *Transitions XTRActive* lenses.

TRANSITIONS® LIGHT INTELLIGENT LENSES™

		Transitions® XTRActive®			Sphere	Cylinder to
		Gray	Brown	Graphite Green		
1.60 High Index Single Vision		●	●	●		
HOYA VISION CARE						
	Sync III	GR			-10.00 to +6.00	-4.00
SHAMIR INSIGHTS INC.						
	Shamir Attitude III Single Vision	GR			-6.50 to +6.00	-7.00
	Shamir Autograph® III Single Vision	GR			-14.00 to +8.75	-7.00
	Shamir Relax™ SFSV	GR			-14.00 to +8.00	-7.00
1.67 Super High Index Single Vision						
CARL ZEISS VISION						
	synchrony 1.67 Aspheric SFSV	GR			-12.00 to +8.00	-6.00
	synchrony Access HD (75, 125 & 175)	GR			-12.00 to +8.00	-6.00
	synchrony SV HDC	GR			-12.00 to +8.00	-6.00
	ZEISS 3D Single Vision	GR			-12.00 to +8.00	-4.00
	ZEISS Individual® SV	GR			-12.00 to +8.00	-4.00
ESSILOR LENSES						
	Essilor SFSV Spherical Thin&Lite® 1.67	GR	BR	GG	-13.00 to +9.00	-4.00
	Essilor Eyecode® SV Thin&Lite® 1.67	GR	BR	GG	-12.00 to +9.00	-6.00
	Essilor FIT® SV Thin&Lite® 1.67	GR	BR	GG	-12.00 to +9.00	-6.00
	Essilor Azio® SV Thin&Lite® 1.67	GR	BR	GG	-12.00 to +9.00	-6.00
	Essilor SV 360° Spherical Thin&Lite® 1.67	GR	BR	GG	-12.00 to +9.00	-6.00
	Essilor Stylistic® Single Vision Thin&Lite® 1.67	GR	BR	GG	-5.00 to +5.00	-3.00
	Eyezen®+ Thin&Lite® 1.67	GR	BR		-14.00 to +9.00	-6.00
	Eyezen®+ Start Thin&Lite® 1.67	GR	BR		-14.00 to +9.00	-6.00



- Shows availability of Transitions® XTRActive® style mirrors in any of the six colors: Gold, Silver Shadow, Pink, Red, Green and Blue. All colors can be applied to either gray or brown Transitions XTRActive lenses.

		Transitions® XTRActive®			Sphere	Cylinder to
		Gray	Brown	Graphite Green		
1.67 Super High Index Single Vision <small>(continued)</small>						
HOYA VISION CARE						
	HOYA SV iQ 1.67	GR			-10.00 to +8.00	-4.00
	Sync III	GR			-13.00 to +7.50	-4.00
NIKON OPTICAL						
	Nikon Infinite Single Vision 1.67	GR	BR		-12.00 to + 8.00	-6.00
	Nikon Smart Boost	GR	BR		-12.00 to + 8.00	-6.00
SEIKO OPTICAL PRODUCTS						
	Seiko PL 1.67 Transitions®	GR			-12.00 to +10.00	-5.00
SHAMIR INSIGHTS INC.						
	Shamir Attitude III Single Vision	GR			-6.50 to +6.00	-8.00
	Shamir Autograph® III Single Vision	GR			-16.00 to +16.00	-8.00
	Shamir Relax®	GR			-16.00 to +8.25	-8.00



- Shows availability of *Transitions® XTRActive®* style mirrors in any of the six colors: Gold, Silver Shadow, Pink, Red, Green and Blue. All colors can be applied to either gray or brown *Transitions XTRActive* lenses.

TRANSITIONS® LIGHT INTELLIGENT LENSES™

SIGNATURE BIFOCALS 1.50

		Transitions XTRActive					
		Gray	Brown	Graphite Green	Sphere	Add	Cylinder to
1.50 Standard Index SFST-28		●	●	●			
ESSILOR LENSES							
	Bifocal FT28	GR			-10.00 to +6.00	+1.00 to +3.00	-4.00
HOYA VISION CARE							
	HOYA Vision Care SFST 28	GR			-11.00 to +6.00	+1.00 to +3.00	-4.00
	HOYA Vision Care SFST 28 iQ	GR			-8.00 to +6.00	+1.00 to +3.00	-4.00
X-CEL OPTICAL							
	X-Cel D-28 78mm	GR			-4.00 to +4.25	+0.75 to +4.00	-4.00
	X-Cel D-28 78mm		BR		-4.00 to +4.25	+3.25 to +4.00	-4.00

BLENDED BIFOCAL

SHAMIR INSIGHT INC.

		Transitions XTRActive						
		Gray	Brown	Graphite Green	Sphere	Add	Cylinder to	MFH
1.50 Standard Index Blended Bifocal <small>(digital)</small>		●	●	●				Minimum Fitting Height
	Shamir Duo	GR			-9.00 to +6.00	+0.75 to +4.00	-6.00	15
1.53 Trivex® Blended Bifocal <small>(digital)</small>								
	Shamir Duo	GR			-13.25 to +6.00	+0.75 to +4.00	-6.50	15
1.59 Polycarbonate Blended Bifocal <small>(digital)</small>								
	Shamir Duo	GR			-14.75 to +6.00	+0.75 to +4.00	-7.00	15
1.60 High Index Blended Bifocal <small>(digital)</small>								
	Shamir Duo	GR			-15.00 to +6.00	+0.75 to +4.00	-7.00	15
1.67 Super High Index Blended Bifocal <small>(digital)</small>								
	Shamir Duo	GR			-16.75 to +8.00	+0.75 to +4.00	-8.00	15

TRANSITIONS® LIGHT INTELLIGENT LENSES™

PROGRESSIVES 1.50

	Transitions XTRActive			Sphere	Add	Cylinder to	MFH
	Gray	Brown	Graphite Green				
1.50 Standard Index Progressives	●	●	●		For Varilux S, at fitting height less than 17mm, ADD power is up to +3.50		Minimum Fitting Height
ESSILOR LENSES							
Varilux® X 4D ⁻	GR	BR		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT ⁻	GR	BR		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design ⁻	GR	BR		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux S 4D ⁻	GR	BR		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux S FIT ⁻	GR	BR		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux S Design ⁻	GR	BR		-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode	GR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT	GR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+	GR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® W2+ Eyecode	GR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+ FIT	GR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+	GR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx ⁻	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx ⁻	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Varilux Comfort Enhanced FIT ⁻	GR			-10.00 to +7.00	+0.75 to +3.50	-6.00	14
Varilux Comfort Enhanced ⁻	GR			-10.00 to +7.00	+0.75 to +3.50	-6.00	14
Varilux Comfort DRx ⁻	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx ⁻	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Varilux Comfort®	GR			-10.00 to +7.00	+0.75 to +3.50	-4.00	17
Varilux Comfort Short ⁻	GR			-10.00 to +7.00	+0.75 to +3.00	-4.00	14
Varilux® Stylistic ⁻ Wrap	GR	BR	GG	-3.00 to +3.00	+0.75 to +3.50	-1.00	17

PROGRESSIVES 1.50

	Transitions XTRActive			Sphere	Add	Cylinder to	MFH
	Gray	Brown	Graphite Green				
1.50 Standard Index Progressives <small>(continued)</small>							
ESSILOR LENSES	●	●	●				Minimum Fitting Height
Essilor Ideal Advanced*	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced* Wrap	GR	BR	GG	-4.00 to +3.00	+0.75 to +3.50	-4.00	17
Essilor Ideal Advanced* FIT	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Essilor Ideal* FIT	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Essilor Ideal*	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Essilor Ideal* Sport	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	18
Essilor Ideal* Computer	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	13
Essilor Ideal* Driving	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Essilor Natural* Digital	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Adaptar* Digital	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Adaptar* Short Digital	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Ovation* Digital	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Essilor SmallFit* Digital	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.00	-6.00	14
HOYA VISION CARE							
HOYALUX Summit CD* 1.50	GR			-10.00 to +6.00	+1.00 to +3.00	-4.00	14
HOYALUX Summit CD* iQ 1.50	GR			-10.00 to +6.00	+1.00 to +3.00	-6.00	14
HOYALUX Summit ECP* 1.50	GR			-10.00 to +6.00	+1.00 to +3.00	-4.00	18
HOYALUX Summit ECP* iQ 1.50	GR			-10.00 to +6.00	+1.00 to +3.00	-6.00	18
iD LifeStyle*2 Clarity 1.50	GR			-8.00 to +6.00	+1.00 to +3.50	-4.00	18
iD LifeStyle*2 Harmony 1.50	GR			-8.00 to +6.00	+1.00 to +3.50	-4.00	18
iD LifeStyle*2 Clarity CD 1.50	GR			-8.00 to +6.00	+1.00 to +3.50	-4.00	14
iD LifeStyle*2 Harmony CD 1.50	GR			-8.00 to +6.00	+1.00 to +3.50	-4.00	14
HOYALUX ARRAY* 11, 13, 15, 17, VL 1.50	GR			-10.00 to +3.75	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY* Wrap 11, 13, 15, 17, VL 1.50	GR			-5.00 to +4.00	+0.75 to +3.50	-4.00	14, 17, 19, 21
Amplitude BKS 1.50	GR			-10.00 to +3.75	+0.75 to +3.50	-4.00	19
Amplitude Mini BKS 1.50	GR			-10.00 to +3.75	+0.75 to +3.50	-4.00	14

PROGRESSIVES 1.50

		Transitions XTRActive			Sphere	Add	Cylinder to	MFH
		Gray	Brown	Graphite Green				
1.50 Standard Index Progressives <small>(continued)</small>								
NIKON OPTICAL		●	●	●				Minimum Fitting Height
	Nikon SeeMax Power		BR		-10.00 to +8.00	+0.75 to +4.00	-6.00	13, 15, 17
	Nikon Presio Power		BR		-10.00 to +8.00	+0.75 to +4.00	-6.00	13, 15, 17
	Nikon DigiLife		BR		-10.00 to +8.00	+0.75 to +4.00	-6.00	13, 15, 17
	Nikon Presio i Digital	GR	BR		-7.00 to +4.00	+0.75 to +4.00	-6.00	15, 17
	Nikon Move	GR	BR		-7.00 to +4.00	+0.75 to +3.50	-4.00	16, 18
SHAMIR INSIGHTS INC.								
	Shamir Autograph® III	GR			-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
	Shamir Autograph® II+	GR			-12.25 to +6.00	+0.75 to +4.00	-6.00	11, 13, 15, 18
	Shamir Intouch®	GR			-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
	Shamir Spectrum+®	GR			-12.25 to +6.00	+0.75 to +3.50	-6.00	14, 16, 18
	Shamir Computer®	GR			-7.00 to +4.00	+0.75 to +2.00	-6.00	14 Lower, 10 Upper
	Shamir Workspace®	GR			-6.00 to +6.00	+0.75 to +2.50	-6.00	14 Lower, 10 Upper
	Shamir Attitude III Sport	GR			-6.50 to +5.25	+0.75 to +3.50	-6.00	18
	Shamir Attitude III Fashion	GR			-6.50 to +5.25	+0.75 to +3.50	-6.00	15, 18
SIGNET ARMORLITE								
	KODAK Unique® 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
	KODAK Unique® HD 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
	KODAK Unique DRO® 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	13
	KODAK Unique DRO® HD 1.50	GR			-8.00 to +4.00	+0.75 to +3.50	-6.00	13

PROGRESSIVES 1.53

	Transitions XTRActive			Sphere	Add	Cylinder to	MFH
	Gray	Brown	Graphite Green				
1.53 TRIVEX Progressives	●	●	●		For Varilux S, at fitting height less than 17mm, ADD power is up to +3.50		Minimum Fitting Height
ESSILOR LENSES							
Varilux® Physio® W3+ Eyecode TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® W2+ Eyecode	GR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+ FIT	GR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+	GR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Physio Enhanced Eyecode® TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio Enhanced FIT® TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio Enhanced® TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio DRx® TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Physio Short DRx® TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Varilux® Physio® TREXA®	GR	BR	GG	-10.00 to +6.00	+1.00 to +3.50 & +4.00	-4.00	17
Varilux Comfort DRx® TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx® TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Essilor Ideal Advanced® TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Essilor Ideal® TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	18
Essilor Ideal® Computer TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	13
Essilor Ideal® Driving TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Essilor Natural® Digital TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom® TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-4.00	15
Essilor Accolade® TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-4.00	17
Adaptar® Digital TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Ovation® Digital TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Essilor SmallFit® Digital TREXA®	GR	BR	GG	-10.00 to +6.00	+0.75 to +3.00	-6.00	14

PROGRESSIVES 1.53

	Transitions XTRActive							
	Gray	Brown	Graphite Green	Sphere	Add	Cylinder to	MFH	
	●	●	●					Minimum Fitting Height
1.53 TRIVEX Progressives <small>(continued)</small>								
HOYA VISION CARE								
HOYALUX Summit CD Phoenix [®]	GR			-10.00 to +6.00	+1.00 to +3.00	-4.00		14
HOYALUX Summit CD iQ Phoenix [®]	GR			-10.00 to +6.00	+1.00 to +3.00	-6.00		14
HOYALUX Summit ECP Phoenix [®]	GR			-10.00 to +6.00	+1.00 to +3.00	-4.00		18
iD LifeStyle CD Phoenix [®]	GR			-8.00 to +6.00	+0.75 to +3.50	-4.00		14
iD LifeStyle2 Clarity Phoenix [®]	GR			-8.00 to +6.00	+0.75 to +3.50	-4.00		18
iD LifeStyle2 Harmony Phoenix [®]	GR			-8.00 to +6.00	+0.75 to +3.50	-4.00		18
iD LifeStyle2 Clarity CD Phoenix [®]	GR			-8.00 to +6.00	+0.75 to +3.50	-4.00		14
iD LifeStyle2 Harmony CD Phoenix [®]	GR			-8.00 to +6.00	+0.75 to +3.50	-4.00		14
HOYALUX ARRAY [®] 11, 13, 15, 17, VL Phoenix [®]	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-4.00		14, 17, 19, 21
HOYALUX ARRAY [®] Wrap 11, 13, 15, 17, VL Phoenix [®]	GR	BR		-5.00 to +4.00	+0.75 to +3.50	-4.00		14, 17, 19, 21
Amplitude BKS Phoenix [®]	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-4.00		19
Amplitude Mini BKS Phoenix [®]	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-4.00		14
SHAMIR INSIGHTS INC.								
Shamir Autograph [®] III	GR			-13.25 to +7.00	+0.75 to +4.00	-6.50		11, 13, 15, 18
Shamir Autograph [®] II+	GR			-13.25 to +6.00	+0.75 to +4.00	-6.50		11, 13, 15, 18
Shamir Intouch [®]	GR			-13.25 to +6.00	+0.75 to +3.50	-6.50		14, 16, 18
Shamir Golf [®]	GR			-5.00 to +4.00	+1.00 to +3.50	-6.50		19
Shamir Spectrum+ [®]	GR			-13.25 to +6.00	+0.75 to +3.50	-6.50		14, 16, 18
Shamir Element [®]	GR			-13.25 to +6.00	+0.75 to +3.50	-6.50		16, 19
Shamir Computer [®]	GR			-6.00 to +8.00	+0.75 to +2.00	-6.50		14 Lower, 10 Upper
Shamir Workspace [®]	GR			-6.00 to +8.00	+0.75 to +2.50	-6.50		14 Lower, 10 Upper
Shamir Attitude III Sport	GR			-6.50 to +6.00	+0.75 to +3.50	-6.50		18
Shamir Attitude III Fashion	GR			-6.50 to +6.00	+0.75 to +3.50	-6.50		15, 18
SIGNET ARMORLITE								
KODAK Precise [®] PB Trivex [®]	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00		17
KODAK Precise [®] Short PB Trivex [®]	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00		13
KODAK Unique Trivex [®]	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00		13
KODAK Unique HD Trivex [®]	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00		13
KODAK Unique DRO [®] Trivex [®]	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00		13
KODAK Unique DRO [®] HD Trivex [®]	GR			-10.00 to +4.50	+0.75 to +3.50	-6.00		13

PROGRESSIVES 1.59

	Transitions XTRActive						
	Gray	Brown	Graphite Green	Sphere	Add	Cylinder to	MFH
1.59 Polycarbonate Progressives	●	●	●				Minimum Fitting Height
CARL ZEISS VISION							
synchrony Easy Adapt HD	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	18
synchrony Easy M HD	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	15
synchrony Easy S HD	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Easy View HD	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	17
synchrony Easy Wear HD	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	18
synchrony Performance HD13	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Performance HD15	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	15
synchrony Performance HD17	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	17
synchrony Performance HD19	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	19
synchrony Performance HDV	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Starter HD	GR			-9.00 to +6.00	+0.50 to +2.00	-6.00	14
synchrony Ultra HD 13	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Ultra HD 15	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	15
synchrony Ultra HD 17	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	17
synchrony Ultra HD 19	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	19
synchrony Ultra HD 21	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	21
synchrony Ultra HDC	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Ultra HDV	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	13
synchrony Work and Go HD	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	17
synchrony Work and Office HD	GR			-9.00 to +6.00	+0.75 to +3.50	-6.00	16
synchrony Work and Read HD	GR			-9.00 to +6.00	+0.75 to +3.00	-6.00	16
ZEISS Choice 13	GR			-10.00 to +6.00	+0.75 to +3.50	-4.00	13
ZEISS Choice 15	GR			-10.00 to +6.00	+0.75 to +3.50	-4.00	15
ZEISS Choice 17	GR			-10.00 to +6.00	+0.75 to +3.50	-4.00	17
ZEISS Choice 19	GR			-10.00 to +6.00	+0.75 to +3.50	-4.00	19
ZEISS Digital 1000	GR			-10.00 to +6.00	+1.00	-4.00	13
ZEISS Digital 1250	GR			-10.00 to +6.00	+1.25	-4.00	13
ZEISS Digital 500	GR			-10.00 to +6.00	+0.50	-4.00	13
ZEISS Digital 750	GR			-10.00 to +6.00	+0.75	-4.00	13
ZEISS Digital Lens	GR			-10.00 to +6.00	+0.50 to +1.25	-4.00	13

PROGRESSIVES 1.59

	Transitions XTRActive			Sphere	Add	Cylinder to	MFH
	Gray	Brown	Graphite Green				
1.59 Polycarbonate Progressives <small>(continued)</small>	●	●	●				Minimum Fitting Height
CARL ZEISS VISION					<small>For Varilux S, at fitting height less than 17mm, ADD power is up to +3.50</small>		
ZEISS Digital Sport	GR			-10.00 to +6.00	+0.50 to +1.25	-4.00	13
ZEISS Digital Sport 1000	GR			-10.00 to +6.00	+1.00	-4.00	13
ZEISS Digital Sport 1250	GR			-10.00 to +6.00	+1.25	-4.00	13
ZEISS Digital Sport 500	GR			-10.00 to +6.00	+0.50	-4.00	13
ZEISS Digital Sport 750	GR			-10.00 to +6.00	+0.75	-4.00	13
ZEISS Individual* 2	GR			-10.00 to +6.00	+0.75 to +3.50	-4.00	13
ZEISS Individual* 2I	GR			-10.00 to +6.00	+0.75 to +3.50	-4.00	13
ZEISS Individual* 2N	GR			-10.00 to +6.00	+0.75 to +3.50	-4.00	13
ZEISS Precision Plus	GR			-10.00 to +6.00	+0.75 to +3.50	-4.00	13
ZEISS Precision Pure 13	GR			-10.00 to +6.00	+0.75 to +3.50	-4.00	13
ZEISS Precision Pure 15	GR			-10.00 to +6.00	+0.75 to +3.50	-4.00	15
ZEISS Precision Pure 17	GR			-10.00 to +6.00	+0.75 to +3.50	-4.00	17
ZEISS Precision Pure 19	GR			-10.00 to +6.00	+0.75 to +3.50	-4.00	19
ZEISS Precision Pure 21	GR			-10.00 to +6.00	+0.75 to +3.50	-4.00	21
ZEISS Precision Superb	GR			-10.00 to +6.00	+0.75 to +3.50	-4.00	13
ESSILOR LENSES							
Varilux® X 4D™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux S 4D™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux S FIT™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux S Design™ Airwear®	GR	BR	GG	-10.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode™ Airwear®	GR	BR		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Airwear®	GR	BR		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Airwear®	GR	BR		-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® W2+ Eyecode™ Airwear®	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+ FIT Airwear®	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+ Airwear®	GR	BR		-10.00 to +6.00	+0.75 to +3.50	-6.00	14

PROGRESSIVES 1.59

	Transitions® XTRActive®			Sphere	Add	Cylinder to	MFH
	Gray	Brown	Graphite Green				
1.59 Polycarbonate Progressives <small>(continued)</small>	●	●	●				For Varilux S, at fitting height less than 17mm, ADD power is up to +3.50
ESSILOR LENSES							Minimum Fitting Height
Varilux Physio Enhanced Eyecode® Airwear®	GR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio Enhanced FIT® Airwear®	GR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio Enhanced® Airwear®	GR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Physio DRx® Airwear®	GR • BR • GG			-10.00 to +8.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx® Airwear®	GR • BR • GG			-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Varilux® Physio® Airwear®	GR			-10.00 to +6.00	+0.75 to +3.50 & +4.00	-4.00	17
Varilux Comfort Enhanced FIT® Airwear®	GR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort Enhanced® Airwear®	GR			-10.00 to +6.00	+0.75 to +3.50	-6.00	14
Varilux Comfort DRx® Airwear®	GR • BR • GG			-10.00 to +6.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx® Airwear®	GR • BR • GG			-10.00 to +6.00	+0.75 to +3.00	-6.00	14
Varilux Comfort® Airwear®	GR			-10.00 to +6.00	+0.75 to +3.50	-4.00	17
Varilux Comfort Short® Airwear®	GR			-10.00 to +6.00	+0.75 to +3.00	-4.00	14
Varilux® Stylistic® Wrap Airwear®	GR • BR • GG			-4.00 to +4.00	+0.75 to +3.50	-2.00	17
Essilor Ideal Advanced® Airwear®	GR • BR • GG			-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced Wrap® Airwear®	GR • BR • GG			-4.00 to +3.00	+0.75 to +3.00	-4.00	17
Essilor Ideal Advanced® FIT	GR • BR • GG			-10.00 to +8.00	+0.75 to +3.00	-6.00	17
Essilor Ideal® FIT Airwear®	GR • BR • GG			-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Airwear®	GR • BR • GG			-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport Airwear®	GR • BR • GG			-10.00 to +8.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Computer Airwear®	GR • BR • GG			-10.00 to +8.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Driving Airwear®	GR • BR • GG			-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Airwear®	GR • BR • GG			-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Accolade Freedom® Airwear®	GR			-10.00 to +6.00	+0.75 to +3.50	-4.00	15
Essilor Accolade® Airwear®	GR			-10.00 to +6.00	+0.75 to +3.50	-4.00	17
Adaptar® Digital Airwear®	GR BR GG			-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Airwear®	GR BR GG			-10.00 to +8.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Airwear®	GR BR GG			-10.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor SmallFit® Digital Airwear®	GR BR GG			-10.00 to +8.00	+0.75 to +3.00	-6.00	14



- Shows availability of *Transitions® XTRActive®* style mirrors in any of the six colors: Gold, Silver Shadow, Pink, Red, Green and Blue. All colors can be applied to either gray or brown *Transitions XTRActive* lenses.

TRANSITIONS® LIGHT INTELLIGENT LENSES™

PROGRESSIVES 1.59

	Transitions XTRActive						
	Gray	Brown	Graphite Green	Sphere	Add	Cylinder to	MFH
1.59 Polycarbonate Progressives <small>(continued)</small>	●	●	●				Minimum Fitting Height
HOYA VISION CARE							
HOYALUX ARRAY® 11, 13, 15, 17, VL Poly	GR			-13.00 to +5.25	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY® Wrap 11, 13, 15, 17, VL Poly	GR			-6.00 to +4.00	+0.75 to +3.50	-4.00	14, 17, 19, 21
Amplitude BKS Poly	GR			-13.00 to +5.25	+0.75 to +3.50	-4.00	19
Amplitude Mini BKS Poly	GR			-13.00 to +5.25	+0.75 to +3.50	-4.00	14
NIKON OPTICAL							
Nikon SeeMax Power	GR	BR		-10.00 to +8.00	+0.75 to +4.00	-6.00	13, 15, 17
Nikon Presio Power	GR	BR		-10.00 to +8.00	+0.75 to +4.00	-6.00	13, 15, 17
Nikon DigiLife	GR	BR		-10.00 to +8.00	+0.75 to +4.00	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR		-10.00 to +5.00	+0.75 to +4.00	-4.00	15, 17
Nikon Move	GR	BR		-10.00 to +5.00	+0.75 to +4.00	-4.00	16, 18
SIGNET ARMORLITE							
KODAK Unique DRO® Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Precise® PB Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	17
KODAK Precise® Short PB Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO® HD Poly	GR			-10.00 to +5.00	+0.75 to +3.50	-6.00	13

PROGRESSIVES 1.59 - 1.60

	Transitions XTRActive							
	Gray	Brown	Graphite Green	Sphere	Add	Cylinder to	MFH	
1.59 Polycarbonate Progressives <small>(continued)</small>	●	●	●					Minimum Fitting Height
SHAMIR INSIGHTS INC.								
Shamir Autograph® III	GR			-14.75 to +8.00	+0.75 to +4.00	-7.00	11, 13, 15, 18	
Shamir Autograph® II+	GR			-14.75 to +6.50	+0.75 to +4.00	-7.00	11, 13, 15, 18	
Shamir Intouch™	GR			-14.75 to +6.00	+0.75 to +3.50	-7.00	15, 18	
Shamir Golf™	GR			-6.00 to +5.25	+1.00 to +3.50	-7.00	19	
Shamir Spectrum+™	GR			-14.75 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18	
Shamir Element™	GR			-10.00 to +6.00	+0.75 to +3.50	-7.00	16, 19	
Shamir Computer™	GR			-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper	
Shamir Workspace™	GR			-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper	
Shamir Attitude III Sport	GR			-6.50 to +6.00	+0.75 to +3.50	-6.50	18	
Shamir Attitude III Fashion	GR			-6.50 to +6.00	+0.75 to +3.50	-6.50	15, 18	
1.60 High Index Progressives								
SHAMIR INSIGHTS INC.								
Shamir Autograph Intelligence™	GR			-14.75 to +8.75	+0.75 to +4.00	-7.00	11, 13, 15, 18	
Shamir Autograph® III	GR			-15.00 to +8.00	+0.75 to +3.50	-7.00	15, 18	
Shamir Autograph® II+	GR			-15.00 to +6.00	+0.75 to +4.00	-7.00	11, 13, 15, 18	
Shamir Golf™ SuperLite™ 1.6	GR			-6.00 to +5.25	+1.00 to +3.50	-7.00	19	
Shamir Spectrum+™	GR			+15.00 to +6.00	+0.75 to +3.50	-7.00	14, 16, 18	
Shamir Computer™	GR			-6.00 to +8.00	+0.75 to +2.00	-7.00	14 Lower, 10 Upper	
Shamir Workspace™	GR			-6.00 to +8.00	+0.75 to +2.50	-7.00	14 Lower, 10 Upper	
Shamir InTouch™	GR			-15.00 to +6.00	+0.75 to +3.50	-7.00	15, 18	
Shamir Attitude III Sport	GR			-6.50 to +6.00	+0.75 to +3.50	-7.00	18	
Shamir Attitude III Fashion	GR			-6.50 to +6.00	+0.75 to +3.50	-7.00	15, 18	

PROGRESSIVES 1.67

	Transitions XTRActive						
	Gray	Brown	Graphite Green	Sphere	Add	Cylinder to	MFH
1.67 Super High Index Progressives	●	●	●				Minimum Fitting Height
CARL ZEISS VISION							
synchrony Easy Adapt HD	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	18
synchrony Easy M HD	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	15
synchrony Easy S HD	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	13
synchrony Easy View HD	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	17
synchrony Easy Wear HD	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	18
synchrony Performance HD13	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	13
synchrony Performance HD15	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	15
synchrony Performance HD17	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	17
synchrony Performance HD19	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	19
synchrony Performance HDV	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	13
synchrony Starter HD	GR			-12.00 to +8.00	+0.50 to +2.00	-6.00	14
synchrony Ultra HD 13	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	13
synchrony Ultra HD 15	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	15
synchrony Ultra HD 17	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	17
synchrony Ultra HD 19	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	19
synchrony Ultra HD 21	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	21
synchrony Ultra HDC	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	13
synchrony Ultra HDV	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	13
synchrony Work and Go HD	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	17
synchrony Work and Office HD	GR			-12.00 to +8.00	+0.75 to +3.50	-6.00	16
synchrony Work and Read HD	GR			-12.00 to +8.00	+0.75 to +3.00	-6.00	16
ZEISS Choice 13	GR			-12.00 to +8.00	+0.75 to +3.50	-4.00	13
ZEISS Choice 15	GR			-12.00 to +8.00	+0.75 to +3.50	-4.00	15
ZEISS Choice 17	GR			-12.00 to +8.00	+0.75 to +3.50	-4.00	17
ZEISS Choice 19	GR			-12.00 to +8.00	+0.75 to +3.50	-4.00	19
ZEISS Digital 1000	GR			-12.00 to +8.00	+1.00	-4.00	13
ZEISS Digital 1250	GR			-12.00 to +8.00	+1.25	-4.00	13
ZEISS Digital 500	GR			-12.00 to +8.00	+0.50	-4.00	13
ZEISS Digital 750	GR			-12.00 to +8.00	+0.75	-4.00	13
ZEISS Digital Lens	GR			-12.00 to +8.00	+0.50 to +1.25	-4.00	13

PROGRESSIVES 1.67

	Transitions® XTRActive®			Sphere	Add	Cylinder to	MFH
	Gray	Brown	Graphite Green				
1.67 Super High Index Progressives (continued)							
For Varilux S, at fitting height less than 17mm, ADD power is up to +3.50							
Minimum Fitting Height							
CARL ZEISS VISION (continued)							
ZEISS Individual® 2	GR			-12.00 to +8.00	+0.75 to +3.50	-4.00	13
ZEISS Individual® 2I	GR			-12.00 to +8.00	+0.75 to +3.50	-4.00	13
ZEISS Individual® 2N	GR			-12.00 to +8.00	+0.75 to +3.50	-4.00	13
ZEISS Precision Plus	GR			-12.00 to +8.00	+0.75 to +3.50	-4.00	13
ZEISS Precision Pure 13	GR			-12.00 to +8.00	+0.75 to +3.50	-4.00	13
ZEISS Precision Pure 15	GR			-12.00 to +8.00	+0.75 to +3.50	-4.00	15
ZEISS Precision Pure 17	GR			-12.00 to +8.00	+0.75 to +3.50	-4.00	17
ZEISS Precision Pure 19	GR			-12.00 to +8.00	+0.75 to +3.50	-4.00	19
ZEISS Precision Pure 21	GR			-12.00 to +8.00	+0.75 to +3.50	-4.00	21
ZEISS Precision Superb	GR			-12.00 to +8.00	+0.75 to +3.50	-4.00	13
ESSILOR LENSES							
Varilux® X 4D® Thin&Lite® 1.67	GR	BR	GG	-12.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X FIT® Thin&Lite® 1.67	GR	BR	GG	-12.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® X Design® Thin&Lite® 1.67	GR	BR	GG	-12.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux S 4D® Thin&Lite® 1.67	GR	BR	GG	-12.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux S FIT® Thin&Lite® 1.67	GR	BR	GG	-12.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux S Design® Thin&Lite® 1.67	GR	BR	GG	-12.00 to +6.00	+0.75 to +4.00	-6.00	14
Varilux® Physio® W3+ Eyecode Thin&Lite® 1.67	GR			-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ FIT Thin&Lite® 1.67	GR			-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux® Physio® W3+ Thin&Lite® 1.67	GR			-14.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	14
Varilux Comfort® W2+ Eyecode Thin&Lite® 1.67	GR			-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+ FIT Thin&Lite® 1.67	GR			-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Comfort® W2+ Thin&Lite® 1.67	GR			-14.00 to +9.00	+0.75 to +3.50	-6.00	14
Varilux Physio DRx® Thin&Lite® 1.67	GR•	BR•	GG	-12.00 to +9.00	+0.75 to +3.50 & +4.00	-6.00	17
Varilux Physio Short DRx® Thin&Lite® 1.67	GR•	BR•	GG	-12.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux Comfort DRx® Thin&Lite® 1.67	GR•	BR•	GG	-12.00 to +9.00	+0.75 to +3.50	-6.00	17
Varilux Comfort Short DRx® Thin&Lite® 1.67	GR•	BR•	GG	-12.00 to +9.00	+0.75 to +3.00	-6.00	14
Varilux® Stylistic® Wrap Thin&Lite® 1.67	GR•	BR•	GG	-5.00 to +5.00	+0.75 to +3.50	-3.00	17



- Shows availability of *Transitions® XTRActive®* style mirrors in any of the six colors: Gold, Silver Shadow, Pink, Red, Green and Blue. All colors can be applied to either gray or brown *Transitions XTRActive* lenses.

TRANSITIONS® LIGHT INTELLIGENT LENSES™

PROGRESSIVES 1.67

	Transitions® XTRActive®			Sphere	Add	Cylinder to	MFH
	Gray	Brown	Graphite Green				
1.67 Super High Index Progressives (continued)							
ESSILOR LENSES (continued)							
Essilor Ideal Advanced® Thin&Lite® 1.67	GR•BR•	GG		-12.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced Wrap® Thin&Lite® 1.67	GR•BR•	GG		-5.00 to +3.00	+0.75 to +3.50	-6.00	17
Essilor Ideal Advanced® FIT Thin&Lite® 1.67	GR•BR•	GG		-14.00 to +8.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® FIT Thin&Lite® 1.67	GR•BR•	GG		-12.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Thin&Lite® 1.67	GR•BR•	GG		-12.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor Ideal® Sport Thin&Lite® 1.67	GR•BR•	GG		-12.00 to +9.00	+0.75 to +3.00	-6.00	18
Essilor Ideal® Computer Thin&Lite® 1.67	GR•BR•	GG		-12.00 to +9.00	+0.75 to +3.00	-6.00	13
Essilor Ideal® Driving Thin&Lite® 1.67	GR•BR•	GG		-12.00 to +9.00	+0.75 to +3.00	-6.00	14
Essilor Natural® Digital Thin&Lite® 1.67	GR•BR•	GG		-12.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Digital Thin&Lite® 1.67	GR•BR•	GG		-12.00 to +9.00	+0.75 to +3.50	-6.00	17
Adaptar® Short Digital Thin&Lite® 1.67	GR•BR•	GG		-12.00 to +9.00	+0.75 to +3.00	-6.00	14
Ovation® Digital Thin&Lite® 1.67	GR•BR•	GG		-12.00 to +9.00	+0.75 to +3.50	-6.00	17
Essilor SmallFit® Digital Thin&Lite® 1.67	GR•BR•	GG		-12.00 to +9.00	+0.75 to +3.00	-6.00	14
HOYA VISION CARE							
HOYALUX ARRAY 11, 13,15,17,VL 1.67	GR			-13.00 to +7.50	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY Wrap 11,13,15,17,VL 1.67	GR			-6.00 to +4.00	+0.75 to +3.50	-4.00	14, 17, 19, 21
HOYALUX ARRAY 2 11, 13,15,17,VL 1.67	GR			-13.00 to +7.50	+0.75 to +3.50	-4.00	14, 17, 19, 21
Amplitude BKS	GR			-13.00 to +7.50	+0.75 to +3.50	-4.00	19
Amplitude Mini BKS	GR			-13.00 to +7.50	+0.75 to +3.50	-4.00	14
NIKON OPTICAL							
Nikon SeeMax Power	GR	BR		-10.00 to +8.00	+0.75 to +3.50	-6.00	13, 15, 17
Nikon Presio Power	GR	BR		-14.00 to +8.00	+0.75 to +3.50	-6.00	13, 15, 17
Nikon DigiLife	GR	BR		-14.00 to +8.00	+0.75 to +3.50	-6.00	13, 15, 17
Nikon Presio i Digital	GR	BR		-12.00 to +8.00	+0.75 to +3.50	-6.00	15, 17
Nikon Move	GR	BR		-12.00 to +5.00	+0.75 to +3.50	-4.00	16, 18



- Shows availability of *Transitions® XTRActive®* style mirrors in any of the six colors: Gold, Silver Shadow, Pink, Red, Green and Blue. All colors can be applied to either gray or brown *Transitions XTRActive* lenses.

TRANSITIONS® LIGHT INTELLIGENT LENSES™

PROGRESSIVES 1.67

	Transitions XTRActive			Sphere	Add	Cylinder to	MFH
	Gray	Brown	Graphite Green				
1.67 Super High Index Progressives <small>(continued)</small>							
SHAMIR INSIGHTS INC.							
Shamir Autograph Intelligence™	GR			-16.75 to +10.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph™ III	GR			-16.75 to +9.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Autograph™ II+	GR			-16.75 to +8.25	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Spectrum+™	GR			-16.75 to +8.00	+0.75 to +4.00	-8.00	11, 13, 15, 18
Shamir Element™ SuperLite™ 1.67	GR			-12.50 to +8.00	+0.75 to +3.50	-8.00	16, 19
Shamir InTouch™	GR			-16.75 to +8.00	+0.75 to +3.50	-8.00	15, 18
Shamir Attitude III Sport	GR			-6.50 to +6.00	+0.75 to +3.50	-8.00	18
Shamir Attitude III Fashion	GR			-6.50 to +6.00	+0.75 to +3.50	-8.00	15, 18
SIGNET ARMORLITE							
KODAK Unique	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique HD	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ 1.67	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13
KODAK Unique DRO™ HD 1.67	GR			-11.00 to +6.00	+0.75 to +3.50	-6.00	13