



Compression
Garment
Collection

CAROMED®
recover beautifully

Quality

As the first post-surgical garment manufacturer, Caromed has defined industry standards for excellence and integrity. Today, our garments are available in new, elegantly designed styles that allow you to recover in comfort—and recover beautifully.

Throughout our history, we've earned a proven reputation for producing superior medical-grade compression garments with unmatched dedication and exacting attention to detail. We have evolved from offering a few select garments, to being a market leader with an extensive array of elite products for cosmetic and reconstructive patients around the world.

All Caromed garments are produced in North Carolina using world-class fabrics from the finest mills throughout the U.S. Each is selected only if it meets our stringent guidelines for medical-grade compressive properties that satisfy scrupulous physician and patient needs. We use only the best available zippers, hook & eye closures and elastics which is why we meticulously cut, sew and produce each garment with accuracy and care.

Caromed's dedication comes from putting our hearts into making sure every garment measures up to a level of excellence medical professionals and patients have come to expect. More than anything, every Caromed garment is designed to assure you can recover in style, so you can feel beautiful about the new you.

Table of Contents

HEAD AND NECK COLLECTION

1-8006	Chin-Neck Garment	2
1-8007	Full Face Garment	3
1-8011	Ear Garment	4
1-8018	Therma-Jaw Garment with Gel Packs	5
1-8026	Face Lift Garment	6
1-8201	Head-Neck Garment	7
	Custom Head and Neck Garment Measuring Instructions	9

ARMS COLLECTION

5-8030	Vest Style Arm Garment	12
5-8035	Bolero Style Arm Garment	14
	Custom Arm Garment Measuring Instructions	16

BREAST AND CHEST COLLECTION

2-8002	Shaped Compression Vest	19
2-8008L	Bra Vest with Cups, Long Style	20
2-8008S	Bra Vest with Cups, Short Style	21
2-8014	Compression Vest	23
2-8015	The Bandeau – Standard Design (Double Bandeau)	25
2-8016	The Bandeau – Single Strap Design	26
2-8020	The Bandeau – One Size Fits All (3” Wide)	28
2-8024	Post-Augmentation Bra/Bandeau	30
2-8120	SofShape™ Breast Augmentation Bra	32
2-8124	Post-Augmentation Bra	34
2-8125	Semi-Shaped Cotton Bra	36
2-8126	Semi-Shaped Satin Bra with Velcro® Straps	38
2-8126EL	Reconstruction Sports Bra with Drain Pouches	40
2-8127	Semi-Shaped White Cotton Bra with Velcro® Straps	41
2-8224	Post-Augmentation Bra with Velcro® Straps (Plus Size)	43
2-100	Caromed Augmentation Sizer Bra	44
2-200	First Impression Bra (by AAdvantage)	45
2-500	First Choice Bra (by AAdvantage)	46
	Custom Breast and Chest Garment Measuring Instructions	48

UNISEX WAISTBINDER COLLECTION

3-8004	Unisex One-Piece Waistbinder	51
3-8527	Unisex Multi-Flap Waistbinder	53
	Custom Waistbinder Measuring Instructions	55

GIRDLE AND COMPRESSION SUIT COLLECTION

4-4000	Powerknit Suspended Girdle	58
4-4500	Satin Suspended Girdle	60
4-4700	Powerknit Suspended Brazilian Lift Garment	62
4-4750	Satin Suspended Brazilian Lift Garment	64
4-8000	Powerknit Girdle	66
3-8005	Female Abdominoplasty Garment	68
3-8025	Adjustable Abdominoplasty Garment	69
4-8200	Abdomino/Lipoplasty Combination Girdle	71
4-8500	Satin Girdle	73
4-8550	Full Body Compression Suit	75
	Custom Girdles and Compression Suit Garment Measuring Instructions	77

MALE GARMENT COLLECTION

2-8003	Male Compression Vest	80
3-9000	Male Abdominoplasty Garment, without Zippers	81
3-9005	Male Abdominoplasty Garment, with Zippers	81
4-8100	Male Compression Girdle	83

PEDIATRIC COLLECTION

6-8004	Pediatric Compression Vest	86
6-8006	Pediatric Chin-Neck Garment	88
6-8035	Pediatric Bolero Style Arm Garment	90

ADDITIONAL INFORMATION

Custom Garment Measuring Instructions	By Collection
Garment Care Instructions	92
Custom Garment Order Form	93
Caromed Fabric Guide	95

CAROMED®

Head and Neck
Garment Collection



Chin-Neck Garment

1-8006

Head & Neck



Features

Provides compression for the chin and upper neck.

- Soft, lightweight, breathable fabric
- 100% latex-free

Ideal for

- Facial surgery
- Neck surgery
- Submental region surgery
- Mentoplasty and/or chin implants
- Ear or periaural procedures
- Oral and maxillofacial surgery

Head/Neck Sizing Chart

Size	Inches		Centimeters	
	Neck	Head	Neck	Head
XS	13	18	33	45.7
S	14	20	35.6	50.8
M	15	22	38.1	55.9
L	16	24	40.6	61
XL	17	26	43.2	66
XXL	18	28	45.7	71.1

Custom sizes available.

For sizing guide, go to page 9.

Colors

Beige

Fabrics

Tricot

Full Face Garment

1-8007

Head & Neck



Features

Provides comfortable compression for all areas of the face and neck, but not the ears.

- Lightweight, soft, breathable fabric
- 100% latex-free

Ideal for

- Facial surgery
- Rhytidectomy and dermabrasion procedures
- Malar implant surgery
- Pediatric facial surgery

Head/Neck Sizing Chart

Size	Inches		Centimeters	
	Neck	Head	Neck	Head
XS	13	18	33	45.7
S	14	20	35.6	50.8
M	15	22	38.1	55.9
L	16	24	40.6	61
XL	17	26	43.2	66
XXL	18	28	45.7	71.1

Custom sizes available.

For sizing guide, go to page 9.

Colors

Beige

Fabrics

Tricot

Ear Garment

1-8011

Head & Neck



Features

The 8011 Ear Garment is made specifically to protect the ear after surgery, holding it comfortably yet firmly in place during healing.

- Lightweight, soft, breathable fabric
- Includes an under-chin strap to hold it in place, if needed
- 100% latex-free

Ideal for

- Ear reconstruction
- Skin or cartilage grafts
- Otoplasty
- Pediatric facial surgery

Head/Neck Sizing Chart

Size	Inches		Centimeters	
	Neck	Head	Neck	Head
XS	13	18	33	45.7
S	14	20	35.6	50.8
M	15	22	38.1	55.9
L	16	24	40.6	61
XL	17	26	43.2	66
XXL	18	28	45.7	71.1

Custom sizes available.

For sizing guide, go to page 9.

Colors

Beige

Fabrics

Tricot

Therma-Jaw Garment

1-8018

Head & Neck



Order Options

- 1-8017 Two Gel Packs (no Garment)
- 1-8018 Garment *with* Gel Packs
- 1-8019 Garment *without* Gel Packs



Features

All of our gel packs allow for immediate cold therapy, and can also be conveniently heated for warm therapy.

- Lightweight, soft, breathable fabric
- Bilateral pockets of soft fabric securely hold the gel packs in place
- 100% latex-free

Ideal for

- Cold/Warm therapy
- Temporomandibular joint surgery or disease
- Craniofacial surgery
- Orthognathic procedures
- Rhytidectomy, parotidectomy surgery
- Oral and maxillofacial surgery
- Pediatric facial surgery

Head/Neck Sizing Chart

Size	Inches		Centimeters	
	Neck	Head	Neck	Head
XS	13	18	33	45.7
S	14	20	35.6	50.8
M	15	22	38.1	55.9
L	16	24	40.6	61
XL	17	26	43.2	66
XXL	18	28	45.7	71.1

Custom sizes available.

For sizing guide, go to page 9.

Colors

Beige

Fabrics

Tricot

Face Lift Garment

1-8026

Head & Neck



Features

Provides compression for neck, jaw, and ears.

- Lightweight, soft, breathable fabric
- Ear slits can be cut for added comfort without the fabric unraveling
- 100% latex-free

Ideal for

- Facial surgery
- Rhytidectomy surgery
- Neck surgery
- Submental region surgery
- Mentoplasty and/or chin implants
- Ear or periaural procedures
- Oral and maxillofacial surgery
- Pediatric facial surgery

Head/Neck Sizing Chart

Size	Inches		Centimeters	
	Neck	Head	Neck	Head
XS	13	18	33	45.7
S	14	20	35.6	50.8
M	15	22	38.1	55.9
L	16	24	40.6	61
XL	17	26	43.2	66
XXL	18	28	45.7	71.1

Custom sizes available.

For sizing guide, go to page 9.

Colors

Beige

Fabrics

Tricot

Head-Neck Garment

1-8201

Head & Neck



Features

Upper band circles the head for added comfort and stability.

- Lightweight, soft, breathable elastic
- One size fits all
- 100% latex-free

One size fits all

Ideal for

- Liposuction on submental region
- Rhytidectomy on submental region
- Scar revisions on submental region
- Pediatric facial surgery

Colors

White

Fabrics

Elastic



Take Two

**One to wear
One to wash**

Because of occasional bleeding after surgery, purchasing two garments allows one to be worn while the other is washed.

Custom Garment Measurement


Head and Neck

The precise fit of a compression garment requires specific measurements. The following guidelines provide you with needed measurements and how-to-measure graphs. All custom orders must be faxed to 888-311-7596 or emailed to sales@caromed.us to ensure accuracy.

Custom Head and Neck Garments

Record the following measurements:

- (A) horizontal head circumference
- (B) vertical head circumference
- (C) neck circumference



CAROMED®

Arm Garment Collection





Vest Style Arm Garment

5-8030



Black Satin

Arms



Features

Includes graded compression from shoulder to wrist.

- Cotton fabric axilla coverage
- Front zipper closure for easy application
- Hook & eye closure behind zipper has a soft under padding
- 100% latex-free

5-8030 Arm Sizing Chart

Size	Upper Arm	Elbow	Wrist	Nominal Length
XS	9.25-9.75	9.25-9.5	5.25-5.375	21.5
S	10-10.5	9.75-10	5.5-5.675	22
M	10.75-11.5	10.25-10.5	5.75-5.875	22.5
L	11.75-12.5	10.75-11	6-6.125	23
XL	12.75-13.5	11.25-11.5	6.1875-6.25	23.5
XXL	13.75-14.5	11.75-12	6.375-6.5	23.75

In inches. Custom sizes available.

Ideal for

- Dermolipectomy
- Lipoplasty of the upper arms
- Lipoplasty of anterior and lateral axillary folds or rolls

For sizing guide, go to page 16.

Colors

Black

Champagne

Pink

Fabrics

Satin



Bolero Style Arm Garment

5-8035



Axilla Detail

Arms



Features

Focuses compression on arms only.

- Cotton axillary and outside seams for greater comfort
- Full back for greater stability and comfort
- Attaches with hook and eye closures across the chest, with a soft under padding
- Graded compression from shoulders to wrists
- 100% latex-free

5-8035 Arm Sizing Chart

Size	Shoulder to Wrist	Mid-Upper Arm	Wrist
XS	17-18	8-9	4.5
S	19-20	9-10	5
M	21-22	10-11	5.5
L	22.5-23	11-12	6
XL	23.5-24	12-13	6.5
XXL	24.5-25	13-14	7

In inches. Custom sizes available.

Colors

Black

Champagne

Fabrics

Satin

Ideal for

- Dermolipectomy
- Lipoplasty of the upper arms
- Lipoplasty of anterior and lateral axillary folds/rolls

For sizing guide, go to page 16.



Take Two

**One to wear
One to wash**

Because of occasional bleeding after surgery, purchasing two garments allows one to be worn while the other is washed.

Custom Garment Measurement

Arm Garments

The precise fit of a compression garment requires specific measurements. The following guidelines provide you with needed measurements and how-to-measure graphs. All custom orders must be faxed to 888-311-7596 or emailed to sales@caromed.us to ensure accuracy.



Custom Arm Garments

Record the following measurements:

Arm Length Measurements

- (V1.0) shoulder to elbow
- (v1.1) shoulder to wrist

Arm Circumferential Measurements

- (W1.0) upper arm deltoid (next to shoulder)
- (W1.1) mid-upper arm
- (W1.2) elbow
- (W1.3) mid-forearm
- (W1.4) wrist
- (X) shoulder to shoulder
- (Y) bra and cup size

CAROMED®

Breast & Chest Garment Collection





Shaped Compression Vest

2-8002



Pink

Breast & Chest



Features

Caromed's shaped compression vests are made from a soft, satin fabric that closely resembles the look and style you enjoy with your everyday undergarments.

- Slightly contoured cup area
- Front detachable straps
- Front zipper closure for easy application
- Hook & eye closures behind zipper have a soft under padding
- 100% latex-free

Ideal for


- Breast reduction
- Breast reconstruction

Chest Sizing Chart

Inches	Centimeters
28	71.1
30	76.2
32	81.3
34	86.4
36	91.4
38	96.5
40	101.6
42	106.7
44	111.8
46	116.8
48	121.9
50	127
52	132.1

Custom sizes available.

Colors



Fabrics



For sizing guide, go to page 48.

Bra Vest with Cups, Long Style

2-8008L



Champagne

Breast & Chest



Features

Provides support for breasts and torso.

- Front detachable straps
- Front zipper closure for easy application
- Hook & eye closures behind zipper have a soft under padding
- Unilateral cup design is available for post-operative mastectomy patients
- 100% latex-free

Ideal for

- Breast reduction
- Breast reconstruction
- Lipoplasty of the back, midriff and flanks
- Mastectomy

Chest Sizing Chart

Inches	Centimeters
28	71.1
30	76.2
32	81.3
34	86.4
36	91.4
38	96.5
40	101.6
42	106.7
44	111.8
46	116.8
48	121.9
50	127
52	132.1

Custom sizes available.

Cup Size

Available in cup sizes A to DD

Colors

Champagne

Fabrics

Satin

For sizing guide, go to page 48.

Bra Vest with Cups, Short Style

2-8008S



Champagne

Breast & Chest



Features

Provides support and stability for recovery.

- Front detachable straps
- Front zipper closure for easy application
- Hook & eye closures behind zipper have a soft under padding
- Unilateral cup design is available for post-operative mastectomy patients
- 100% latex-free

Ideal for

- Breast reduction
- Breast reconstruction
- Lipoplasty of the back, midriff and flanks
- Mastectomy

Chest Sizing Chart

Inches	Centimeters
28	71.1
30	76.2
32	81.3
34	86.4
36	91.4
38	96.5
40	101.6
42	106.7
44	111.8
46	116.8
48	121.9
50	127
52	132.1

Custom sizes available.

Cup Size

Available in cup sizes A to DD

Colors

Champagne

Fabrics

Satin

For sizing guide, go to page 48.



Compression Vest

2-8014



Pink

Breast & Chest



Features

Caromed's compression vests are expertly crafted from soft, breathable Powerknit fabric so you can recover beautifully.

- Soft, breathable Powerknit fabric
- Front detachable straps
- Front zipper closure for easy application
- Hook & eye closures behind zipper have a soft under padding
- 100% latex-free

Ideal for

- Submuscular or subglandular breast augmentation
- Mastectomy

Chest Sizing Chart

Inches	Centimeters
28	71.1
30	76.2
32	81.3
34	86.4
36	91.4
38	96.5
40	101.6
42	106.7
44	111.8
46	116.8
48	121.9
50	127
52	132.1

Custom sizes available.

Colors



Fabrics



For sizing guide, go to page 48.



The Bandeau - Double Strap Design

2-8015

Available in Standard and Lace styles.



Champagne, Lace Style



Champagne, Standard Style

Breast & Chest



Features

Provides downward compression on top to keep implants in place while providing support for the breasts with the lower band.

- Velcro® closures on both ends of the top strap so you can adjust angle and compression
- Soft padding strip to protect inframammary crease
- 100% latex-free

Ideal for

- Submuscular breast augmentation
- Subglandular breast augmentation

Chest Sizing Chart

Inches	Centimeters
28	71.1
30	76.2
32	81.3
34	86.4
36	91.4
38	96.5
40	101.6
42	106.7
44	111.8
46	116.8
48	121.9
50	127
52	132.1

Custom sizes available.

Styles

Available in Standard and Lace styles

Colors

- Black
- Champagne

Fabrics

- Satin

For sizing guide, go to page 48.

The Bandeau - Single Strap Design

2-8016



Champagne

Breast & Chest



Features

Provides downward compression on implants to keep them properly positioned.

- Velcro® closures on both ends of the strap so you can adjust angle and compression
- 100% latex-free

Ideal for

- Submuscular breast augmentation
- Subglandular breast augmentation

Chest Sizing Chart

Inches	Centimeters
28	71.1
30	76.2
32	81.3
34	86.4
36	91.4
38	96.5
40	101.6
42	106.7
44	111.8
46	116.8
48	121.9
50	127
52	132.1

Custom sizes available.

Colors

Black

Champagne

Fabrics

Satin

For sizing guide, go to page 48.



The Bandeau - One Size Fits All (3")

2-8020



Champagne



Breast & Chest



Features

Provides downward compression on the top of breasts to keep implants in place during healing.

- Velcro® closures on both ends of the strap so you can adjust angle and compression easily
- 100% latex-free

Ideal for

- Submuscular breast augmentation
- Subglandular breast augmentation

One size fits all

Colors

Black

Champagne

Fabrics

Cotton




Post-Augmentation Bra/Bandeau

2-8024



Breast & Chest



Features

Combines the 2-8124 post-augmentation bra with a bandeau to stabilize implants.

- Break-away adjustable straps
- Front zipper closure for easy application
- Formed cups
- Simulated underwire support
- Bandeau strap can be adjusted for angle and compression
- 100% latex-free

Ideal for


- Breast augmentation or reduction
- Breast reconstruction
- Submuscular breast augmentation
- Subglandular breast augmentation

Chest Sizing Chart

Inches	Centimeters
28	71.1
30	76.2
32	81.3
34	86.4
36	91.4
38	96.5
40	101.6
42	106.7
44	111.8
46	116.8
48	121.9
50	127
52	132.1

Custom sizes available.

Colors



Fabrics



For sizing guide, go to page 48.



SofShape™ Breast Augmentation Bra

New and Popular

2-8120



Breast support detail

Breast & Chest



Features

The SofShape™ Bra is the only bra with an internal support system designed to maintain contour and protect the incision site with comfort beyond comparison.

- Break-away adjustable straps
- Front zipper closure for easy application
- Formed cups
- Shaped, comfortable internal support system assures a natural breast shape after healing
- 100% latex-free

Ideal for

- Breast augmentation, reduction
- Breast reconstruction

Chest Sizing Chart

Inches	Centimeters
28	71.1
30	76.2
32	81.3
34	86.4
36	91.4
38	96.5
40	101.6
42	106.7
44	111.8
46	116.8
48	121.9
50	127
52	132.1

Custom sizes available.

Cup Size

Available in cup sizes A to DD

Colors

Black

Champagne

Fabrics

Satin

For sizing guide, go to page 48.



Post-Augmentation Bra

2-8124



Black

Breast & Chest



Features

Caromed post augmentation bras feature medical grade compression fabrics that help minimize swelling and bruising, while providing a slightly contoured design for support and comfort

- Front zipper closure for easy application
- Formed cups
- Simulated underwire support

Ideal for

- Breast augmentation
- Breast reduction
- Breast reconstruction


Chest Sizing Chart

Inches	Centimeters
28	71.1
30	76.2
32	81.3
34	86.4
36	91.4
38	96.5
40	101.6

For sizes above 40, see 2-8224 on page 55.
Custom sizes also available.

For sizing guide, go to page 48.

Colors



Fabrics





Semi-Shaped Cotton Bra

2-8125



Clasp Detail

Breast & Chest



Features

Slightly contoured for comfort and support.

- Front hook & eye closures have a soft under padding
- 100% latex-free

Ideal for

- Breast augmentation
- Breast reduction
- Breast reconstruction

Chest Sizing Chart

Inches	Centimeters
28	71.1
30	76.2
32	81.3
34	86.4
36	91.4
38	96.5
40	101.6
42	106.7
44	111.8
46	116.8
48	121.9
50	127
52	132.1

Custom sizes available.

Colors

Black

White

Fabrics

Cotton

For sizing guide, go to page 48.



Semi-Shaped Bra w/ Velcro® Straps

2-8126



Clasp Detail

Breast & Chest



Features

Slightly contoured for comfort and support.

- Detachable Velcro® straps
- Front hook & eye closures have a soft under padding
- 100% latex-free

Ideal for

- Breast augmentation
- Breast reduction
- Breast reconstruction

Chest Sizing Chart

Inches	Centimeters
28	71.1
30	76.2
32	81.3
34	86.4
36	91.4
38	96.5
40	101.6
42	106.7
44	111.8
46	116.8
48	121.9
50	127
52	132.1

Custom sizes available.

Colors

Black

Champagne

Fabrics

Satin

For sizing guide, go to page 48.



Reconstruction Sports Bra

2-8126EL



Pink, with Drain Pouches

Breast & Chest



Features

Provides active support after surgery.

- Satin fabric that is internally lined with soft, breathable cotton
- Adjustable Velcro® front closure
- Adjustable Velcro® straps
- Two drain pouches with Velcro® attachments on both sides can be easily removed
- Additional drain pouches can also be ordered
- 100% latex-free

Ideal for

- Breast reduction
- Breast reconstruction

Chest Sizing Chart

Inches	Centimeters
28	71.1
30	76.2
32	81.3
34	86.4
36	91.4
38	96.5
40	101.6
42	106.7
44	111.8
46	116.8
48	121.9
50	127
52	132.1

Custom sizes available.

Cup Size

Available in cup sizes A to DD

Pouches

Drain pouches can be easily removed.

Colors

Champagne

Pink

Fabrics

Satin

For sizing guide, go to page 48.

Semi-Shaped Bra with Velcro® Straps

2-8127



White

Breast & Chest



Features

Slightly contoured for comfort and support.

- Detachable Velcro® straps
- Front hook & eye closures have a soft under padding
- 100% latex-free

Ideal for

- Breast augmentation
- Breast reduction
- Breast reconstruction

Chest Sizing Chart

Inches	Centimeters
28	71.1
30	76.2
32	81.3
34	86.4
36	91.4
38	96.5
40	101.6
42	106.7
44	111.8
46	116.8
48	121.9
50	127
52	132.1

Custom sizes available.

Colors

Black

White

Fabrics

Cotton

For sizing guide, go to page 48.



Post-Augmentation Bra with Velcro® Strap (Plus Size)

2-8224



Breast & Chest



Features

Caromed post augmentation bras feature medical grade compression fabrics that help minimize swelling and bruising, while providing a slightly contoured design for support and comfort

- Break-away adjustable straps
- Front zipper closure for easy application
- Formed cups
- Heavy straps for greater comfort
- Hook & eye closures behind zipper have a soft under padding

Ideal for

- Breast augmentation
- Breast reduction
- Breast reconstruction

Chest Sizing Chart

Inches	Centimeters
42	106.7
44	111.8
46	116.8
48	121.9
50	127
52	132.1

For sizes below 42, see 2-8024 on page 45.
Custom sizes available.

For sizing guide, go to page 48.

Colors

Black
Champagne
Pink

Fabrics

Satin

Caromed Augmentation Sizer Bra

2-100



White Sizer Bra

Breast & Chest



Features

Provides a way to support the use of breast prosthesis and is used to demonstrate the effects of different size implants.

- Break-away adjustable straps
- Front zipper closure for easy application
- Simulated underwire support
- Provides spaces to insert implants to demonstrate size changes

Chest Sizing Chart

Size	Inches
S	32-34
M	34.5-36
L	36.5-38
XL	38.5-40
XXL	40.5+

Colors

White

Ideal for

- Breast augmentation
- Breast reduction
- Breast reconstruction

For sizing guide, go to page 48.

Fabrics

Cotton

First Impression (by AAdvantage)

2-200



White Nylon/Lycra®, with Pouch

Breast & Chest



Features

- Adjustable straps with Velcro® closure
- Zipper front closure
- Seamless one-piece design provides coverage and protection of incision area
- Single knit nylon and Lycra® Spandex® provides maximum support and security with a superior fit
- Removable and adjustable pouch with Velcro® strap easily accommodates drain or pain pump

Chest Sizing Chart

Size	Inches
S	32-34
M	34.5-36
L	36.5-38
XL	38.5-40
XXL	40.5+

For sizing guide, go to page 48.

Colors

White

Fabrics

Nylon/Lycra®

Ideal for

- Breast Augmentation or Reduction
- Breast Reconstruction
- Biopsy
- Lumpectomy
- Mastectomy

First Choice (by AAdvantage)

2-500



White Nylon/Lycra[®], with Pouch

Breast & Chest



Features

- Adjustable straps with Velcro[®] closures
- Hook-and-Eye front closure
- Wide elastic band provides maximum support and security
- Crisscross back for posture support
- 3-piece cups with flat seams
- Removable and adjustable pouch with Velcro[®] strap easily accommodates drain or pain pump

Chest Sizing Chart

Size	Inches
S	32-34
M	34.5-36
L	36.5-38
XL	38.5-40
2XL	40.5-42
3XL	42.5-44
4XL	44.5-46

Colors

White

Fabrics

Nylon/Lycra[®]

Ideal for

- Breast Augmentation
- Breast Reconstruction
- Breast Reduction
- Biopsy
- Lumpectomy
- Mastectomy

For sizing guide, go to page 48.

The background is a solid light green color with intricate, light-colored floral and scrollwork patterns. A central dark green rectangular box with a thin white border contains the main text. Below the box, there is a paragraph of text. At the bottom of the page, there is a large, light green wavy shape that resembles a stylized hill or a decorative wave.

Take Two

**One to wear
One to wash**

Because of occasional bleeding after surgery, purchasing two garments allows one to be worn while the other is washed.

Custom Garment Measurement


Breast and Chest Garments

The precise fit of a compression garment requires specific measurements. The following guidelines provide you with needed measurements and how-to-measure graphs. All custom orders must be faxed to 888-311-7596 or emailed to sales@caromed.us to ensure accuracy.

Custom Breast and Chest Garments

Record the following measurements:

- (D) shoulder width
- (E) shoulder to natural waistline
- (F) bust
- (G) inframammary fold
- (H) waist



CAROMED®

Waistbinder Collection





Unisex One-Flap Waistbinder

(white) 3-8004
(black) 3-8504



White Waistbinder



Black Waistbinder

Waistbinders



Features

- Soft cotton liner wicks moisture away from the body
- X-shape elastic across the back ensures a “hug the body” fit
- Reinforced back panel provides extra support
- Breathable fabric
- 100% latex-free

Ideal for

- Abdominoplasty
- Lipoplasty of the abdomen, flanks and midriff

Waist Sizing Chart

Size	Waist	
	Inches	Centimeters
XS	22-24	55.9-61
S	26-28	66-71.1
M	30-36	76.2-91.4
L	38-42	96.5-106.7
XL	44-46	111.8-116.8
XXL	47-48	119.4-121.9
XXXL	50-51	127-129.5

Custom sizes available.

Widths/Colors

White:

3-8004-9 9"
3-8004-12 12"
3-8004-15 15"

Black:

3-8504-9 9"
3-8504-12 12"
3-8504-15 15"

Colors

Black

White

Fabrics

Cotton

For sizing guide, go to page 55.



Unisex Multi-Flap Waistbinder

(white) 3-8027
(black) 3-8527



Black Waistbinder



White Waistbinder

Waistbinders



Features

- Soft, cotton liner wicks away moisture from the body
- X-shape elastic across the back ensures a “hug the body” fit
- Reinforced back panel provides extra support
- Breathable fabric
- 100% latex-free

Ideal for

- Abdominoplasty
- Lipoplasty of the abdomen, flanks and midriff

Waist Sizing Chart

Size	Waist	
	Inches	Centimeters
XS	22-24	55.9-61
S	26-28	66-71.1
M	30-36	76.2-91.4
L	38-42	96.5-106.7
XL	44-46	111.8-116.8
XXL	47-48	119.4-121.9
XXXL	50-51	127-129.5

Custom sizes available.

Widths/Colors

White:	
3-8027	9"
3-8037	12"
3-8047	15"
Black:	
3-8527	9"
3-8537	12"
3-8547	15"

Colors

Black

White

Fabrics

Cotton

For sizing guide, go to page 55.



Take Two

**One to wear
One to wash**

Because of occasional bleeding after surgery, purchasing two garments allows one to be worn while the other is washed.

Custom Garment Measurement


Waistbinders

The precise fit of a compression garment requires specific measurements. The following guidelines provide you with needed measurements and how-to-measure graphs. All custom orders must be faxed to 888-311-7596 or emailed to sales@caromed.us to ensure accuracy.

Custom Waistbinders

Record the following measurements:

- (E) shoulder to natural waistline
- (H) waist
- (I) waist to crotch
- (K) hip



CAROMED®

Girdle Collection





Powerknit Suspendered Girdle

4-4000



Black, Legless, Back



Champagne, Capri

Order by Length:

- 4-4010 Legless
- 4-4011 Mid-Thigh
- 4-4012 Above-the-Knee
- 4-4000 Capri
- 4-4013 Ankle



Girdles/Suits



Features

- Lace and hidden seaming on a reinforced panel for a smooth appearance
- Zippers with hook & eye closures for easy application in the OR and at home
- Front detachable suspender strap
- Available without zippers as a second stage garment
- Open crotch design
- **Can be special ordered with a high back**
- 100% latex-free

Ideal for

- Lipoplasty of the trochanteric area, abdomen, thighs, hips, buttocks and flanks

Waist Sizing Chart

Size	Waist	
	Inches	Centimeters
XS	22-24	55.9-61
S	26-28	66-71.1
M	30-36	76.2-91.4
L	38-42	96.5-106.7
XL	44-46	111.8-116.8
XXL	47-48	119.4-121.9
XXXL	50-51	127-129.5

Custom sizes available.

Colors



Fabrics

Powerknit

For sizing guide, go to page 77.



Satin Suspended Girdle

4-4500



Pink, Capri, **High Back**



Mid-Thigh Length

Order by Length:

- 4-4510 Legless
- 4-4511 Mid-Thigh
- 4-4512 Above-the-Knee
- 4-4500 Capri
- 4-4513 Ankle



Girdles/Suits



Features

- Lace and hidden seaming on a reinforced panel for a smooth appearance
- Zippers with hook & eye closures for easy application in the OR and at home
- Front detachable suspender strap
- Available without zippers as a second stage garment
- Open crotch design
- **Can be special ordered with a high back**
- 100% latex-free

Waist Sizing Chart

Size	Waist	
	Inches	Centimeters
XS	22-24	55.9-61
S	26-28	66-71.1
M	30-36	76.2-91.4
L	38-42	96.5-106.7
XL	44-46	111.8-116.8
XXL	47-48	119.4-121.9
XXXL	50-51	127-129.5

Custom sizes available.

Colors



Fabrics



Ideal for

- Lipoplasty of the trochanteric area, abdomen, thighs, hips, buttocks and flanks

For sizing guide, go to page 77.



Powerknit Suspended Brazilian Lift

Garment

4-4700



Black Satin, Front



Black Satin, Side

Note:
Images shown are for the 4-4750 in satin. 4-4700 is the same pattern, but in Powerknit fabric.

Order by Length:

- 4-4701 Mid-Thigh
- 4-4702 Above-the-Knee
- 4-4700 Capri
- 4-4703 Ankle



Girdles/Suits



Features

- Designed with gluteal cutouts and a soft liner to cover buttocks
- Side zippers with hook & eye closures for easy application in the operating room and at home
- Open crotch design for added comfort and convenience
- Suspended to keep the girdle in place
- 100% latex-free

Abdomen & Waist Sizing Chart

Size	Inches		Centimeters	
	Hip	Waist	Hip	Waist
XS	33-34	23-24	83.8-86.4	58.4-61
S	35-36	25-26	88.9-91.4	63.5-66
M	37-38	27-28	94-96.5	68.6-71.1
L	40-43	30-33	101.6-109.2	76.2-83.8
XL	44-46	34-36	111.8-116.8	86.4-91.4
XXL	48-50	38-40	121.9-127	96.5-101.6

Custom sizes available.

Colors

Black

Fabrics

Powerknit

Ideal for

- Brazilian Buttocks Lifts
- Buttocks implants
- Natural fat transfers to buttocks

For sizing guide, go to page 77.



Satin Suspended Brazilian Lift Garment

4-4750



Black Satin, Front



Black Satin, Side

Order by Length:

- 4-4751 Mid-Thigh
- 4-4752 Above-the-Knee
- 4-4750 Capri
- 4-4753 Ankle



Girdles/Suits



Features

- Designed with gluteal cutouts and a soft liner to cover buttocks
- Side zippers with hook & eye closures for easy application in the operating room and at home
- Open crotch design for added comfort and convenience
- Suspended to keep the girdle in place
- 100% latex-free

Abdomen & Waist Sizing Chart

Size	Inches		Centimeters	
	Hip	Waist	Hip	Waist
XS	33-34	23-24	83.8-86.4	58.4-61
S	35-36	25-26	88.9-91.4	63.5-66
M	37-38	27-28	94-96.5	68.6-71.1
L	40-43	30-33	101.6-109.2	76.2-83.8
XL	44-46	34-36	111.8-116.8	86.4-91.4
XXL	48-50	38-40	121.9-127	96.5-101.6

Custom sizes available.

Ideal for

- Brazilian Buttocks Lifts
- Buttocks implants
- Natural fat transfers to buttocks

For sizing guide, go to page 77.

Colors

Black

Fabrics

Satin



Powerknit Girdle

4-8000



Pink, Above-the-Knee



Champagne, Capri

Order by Length:

- 4-8011 Mid-Thigh
- 4-8012 Above-the-Knee
- 4-8000 Capri
- 4-8013 Ankle



Girdles/Suits



Features

- Lace and hidden seaming on a reinforced panel for a smooth appearance
- Standard 4" waistband
- 2" and 6" waistbands available as custom option
- Zippers with hook & eye closures for easy application in the OR and at home
- Open crotch design
- 100% latex-free

Abdomen & Waist Sizing Chart

Size	Inches		Centimeters	
	Hip	Waist	Hip	Waist
XS	33-34	23-24	83.8-86.4	58.4-61
S	35-36	25-26	88.9-91.4	63.5-66
M	37-38	27-28	94-96.5	68.6-71.1
L	40-43	30-33	101.6-109.2	76.2-83.8
XL	44-46	34-36	111.8-116.8	86.4-91.4
XXL	48-50	38-40	121.9-127	96.5-101.6


Custom sizes available.

For sizing guide, go to page 77.

Ideal for

- Lipoplasty of the trochanteric area, abdomen, thighs, hips, buttocks, flanks and knees

Colors



Fabrics





Female Abdominoplasty Garment

(without Drain Pouch) 3-8005
(Lined with Pouch) 3-8005-LP



Champagne

Girdles/Suits



Features

- Made with lightweight, breathable Powerknit fabric
- Standard 4" waistband
- 2" and 6" waistbands available as custom option
- Slit opening in the flap to accommodate urinary catheters
- Additional drain pouches can also be ordered
- Padded, easy-to-use zippers
- Velcro® flap crotch closure
- 100% latex-free

Ideal for

- Lipectomy of the abdomen
- Traditional abdominoplasty

Abdomen & Waist Sizing Chart

Size	Inches		Centimeters	
	Hip	Waist	Hip	Waist
XS	33-34	23-24	83.8-86.4	58.4-61
S	35-36	25-26	88.9-91.4	63.5-66
M	37-38	27-28	94-96.5	68.6-71.1
L	40-43	30-33	101.6-109.2	76.2-83.8
XL	44-46	34-36	111.8-116.8	86.4-91.4
XXL	48-50	38-40	121.9-127	96.5-101.6

Custom sizes available.

Drain Pouches

Drain pouches can be added to this garment by ordering 3-8005-LP

Colors

Black

Champagne

Pink

Fabrics

Powerknit

For sizing guide, go to page 77.


Adjustable Abdominoplasty Garment

3-8025



Champagne

Girdles/Suits



Features

- Adjustable Velcro® openings at the waist, abdomen and crotch
- Powerknit back and bottom combined with cotton waistbinder
- Can be ordered for men with scrotal support pouch
- Slit opening in the flap to accommodate urinary catheters
- 100% latex-free

Abdomen & Waist Sizing Chart

Size	Inches		Centimeters	
	Hip	Waist	Hip	Waist
XS	33-34	23-24	83.8-86.4	58.4-61
S	35-36	25-26	88.9-91.4	63.5-66
M	37-38	27-28	94-96.5	68.6-71.1
L	40-43	30-33	101.6-109.2	76.2-83.8
XL	44-46	34-36	111.8-116.8	86.4-91.4
XXL	48-50	38-40	121.9-127	96.5-101.6

Custom sizes available.

Colors



Fabrics

Powerknit

Ideal for

- Lipoplasty of the trochanteric area, abdomen and flanks

For sizing guide, go to page 77.



Abdomino/Lipo Combination Girdle

4-8200



Champagne, Capri Length



Black, Capri Length

Order by Length:

- 4-8212 Above-the-Knee
- 4-8200 Capri
- 4-8213 Ankle



Girdles/Suits



Features


- Lightweight, breathable Powerknit fabric
- Hidden seaming on a reinforced panel for a smooth appearance
- Standard 6" waistband
- 2" and 4" waistbands available as custom option
- Velcro® flap closure for precise abdominal compression
- Zippers with hook & eye closures for easy application in the operating room and at home
- 100% latex-free

Abdomen & Waist Sizing Chart

Size	Inches		Centimeters	
	Hip	Waist	Hip	Waist
XS	33-34	23-24	83.8-86.4	58.4-61
S	35-36	25-26	88.9-91.4	63.5-66
M	37-38	27-28	94-96.5	68.6-71.1
L	40-43	30-33	101.6-109.2	76.2-83.8
XL	44-46	34-36	111.8-116.8	86.4-91.4
XXL	48-50	38-40	121.9-127	96.5-101.6
Chest sizes 32-44			Chest sizes 81.3-111.8	

Custom sizes available.

Colors



Fabrics



Ideal for

- Abdominoplasty plus lipoplasty of the trochanteric area, abdomen, thighs, hips, flanks and knees

For sizing guide, go to page 77.



Satin Girdle

4-8500



Champagne, Above Knee



Black, Above Knee

Order by Length:

- 4-8511 Mid-Thigh
- 4-8512 Above-the-Knee
- 4-8500 Capri
- 4-8513 Ankle



Girdles/Suits



Features

- Lace and hidden seaming on a reinforced panel for a smooth appearance
- Standard 4" waistband
- 2" and 6" waistbands available as custom option
- Zippers with hook & eye closures for easy application in the OR and at home
- Open crotch design
- 100% latex-free

Abdomen & Waist Sizing Chart

Size	Inches		Centimeters	
	Hip	Waist	Hip	Waist
XS	33-34	23-24	83.8-86.4	58.4-61
S	35-36	25-26	88.9-91.4	63.5-66
M	37-38	27-28	94-96.5	68.6-71.1
L	40-43	30-33	101.6-109.2	76.2-83.8
XL	44-46	34-36	111.8-116.8	86.4-91.4
XXL	48-50	38-40	121.9-127	96.5-101.6

Custom sizes available.

Ideal for

- Lipoplasty of the trochanteric area, abdomen, thighs, hips, buttocks, flanks and knees

For sizing guide, go to page 77.

Colors



Fabrics





Full Body Compression Suit

(Capri) 4-8550
(Above-the-Knee) 4-8552



Black, Capri

Girdles/Suits



Features

- Combination of lightweight, breathable Powerknit body and a soft satin bra top
- Semi-shaped bra with adjustable straps
- Two full-length zippers
- Available as a second stage bodysuit with or without zippers and without an open crotch
- Open crotch design
- 100% latex-free

Ideal for

- Procedures combining breast reduction, reconstruction or augmentation with multiple area lipoplasty

Abdomen & Waist Sizing Chart

Size	Inches		Centimeters	
	Hip	Waist	Hip	Waist
XS	33-34	23-24	83.8-86.4	58.4-61
S	35-36	25-26	88.9-91.4	63.5-66
M	37-38	27-28	94-96.5	68.6-71.1
L	40-43	30-33	101.6-109.2	76.2-83.8
XL	44-46	34-36	111.8-116.8	86.4-91.4
XXL	48-50	38-40	121.9-127	96.5-101.6
Chest sizes 32-44			Chest sizes 81.3-111.8	

Custom sizes available.

For sizing guide, go to page 77.

Colors



Fabrics





Take Two

**One to wear
One to wash**

Because of occasional bleeding after surgery, purchasing two garments allows one to be worn while the other is washed.

Custom Garment Measurement

Girdles and Compression Suits

The precise fit of a compression garment requires specific measurements. The following guidelines provide you with needed measurements and how-to-measure graphs. All custom orders must be faxed to 888-311-7596 or emailed to sales@caromed.us to ensure accuracy.

Custom Girdles and Compression Suits


Record the following measurements:

Circumferential Measurements

- (F) bust
- (G) inframammary fold
- (H) waist
- (K) hip
- (N) upper thigh
- (O) mid thigh
- (P) 2" above knee
- (Q) knee
- (R) 2" below knee
- (S) upper calf
- (T) mid calf
- (U) ankle

Lengths


- (I) waist to crotch
- (J) crotch
- (L) waist to knee
- (M) inseam



CAROMED®

Male Garment Collection





Male Compression Vest

2-8003



Black

Men's



Features

- Lightweight, breathable Powerknit fabric
- Double-reinforced front panels
- Comfortable Velcro® shoulder straps
- Easy-to-use zipper and hook & eye closures have a soft under padding
- 100% latex free

Ideal for

- Gynecomastia procedures
- Lipoplasty of the back, midriff and flanks

Male Chest Sizing Chart

Size	Inches		Centimeters	
	Chest	Waist	Chest	Waist
XS	30-32	34	76.2-81.3	86.4
S	34-36	36	86.4-91.4	91.4
M	38-40	38	96.5-101.6	96.5
L	42-44	40	106.7-111.8	101.6
XL	46-48	42	116.8-121.9	106.7
XXL	50-52	44	127-132.1	111.8
XXXL	54-56	46	137.1-142.3	116.9

Custom sizes available.

For sizing guide, go to page 48.

Colors

Black

Beige

Fabrics

Powerknit

Male Abdominoplasty Garment



Beige 3-9005, with Zippers

(no zippers) 3-9000
(with zippers) 3-9005

Note:

Image shown is the 3-9005.
The 3-9000 does not have zippers.

Men's



Features

- Lightweight, breathable Powerknit fabric
- Velcro® flap closure
- Reinforced, double-zipper front closure (3-9005 only)
- Scrotal support pouch
- 4" waistband
- 100% latex-free

Ideal for


- Traditional abdominoplasty
- Lipectomy of the abdomen

Male Waist Sizing Chart

Size	Waist	
	Inches	Centimeters
XS	26-27	66-68.6
S	28-30	71.1-76.2
M	32-34	81.3-86.4
L	36-38	91.4-96.5
XL	40-42	101.6-106.7
XXL	44-46	111.8-116.8
XXXL	48-50	122-126.9

Custom sizes available.

Colors



Fabrics



For sizing guide, go to page 55.



Male Compression Girdle

4-8100



Beige



Black

Order by Length:

- 4-8100 Above-the-Knee
- 4-8112 Below-the-Knee
- 4-8113 Ankle



Men's



Features

- Lightweight, breathable Powerknit fabric
- Convenient open crotch
- Double-reinforced front panels
- Optional scrotal support (available upon request)
- 4" waistband
- 100% latex-free

Ideal for


- Lipoplasty of the trochanteric area, the abdomen or thighs, hips, buttocks, flanks and knees

Male Waist Sizing Chart

Size	Waist	
	Inches	Centimeters
XS	26-27	66-68.6
S	28-30	71.1-76.2
M	32-34	81.3-86.4
L	36-38	91.4-96.5
XL	40-42	101.6-106.7
XXL	44-46	111.8-116.8
XXXL	48-50	122-126.9

Custom sizes available.

Colors



Fabrics



For sizing guide, go to page 55.

CAROMED®

Pediatric Garment Collection





Pediatric Compression Vest

6-8004



Black

Pediatric



Features

Caromed's compression vests are expertly crafted from soft, breathable Powerknit fabric.

- Soft, breathable Powerknit fabric
- Front detachable straps
- Front zipper closure for easy application
- Hook & eye closures behind zipper have a soft under padding
- 100% latex-free

Ideal for

- Burn skin grafting and scar revision
- Pediatric reconstruction
- Congenital abnormalities
- Soft tissue tumors

Chest/Vest Sizing Chart

Pediatric 6-8004

Size	Inches		Centimeters	
	Chest	Waist	Chest	Waist
Toddler				
S	17-19.25	17-19	43.2-48.9	43.2-48.3
M	19.25-21.25	19-20.5	48.9-54	48.3-52.1
L	21.25-23.25	21-23	54-59.1	53.3-58.4
XL	23.25-25.25	23-25	59.1-64.1	58.4-63.5
Child				
S	23.5-25	20.5-22	59.7-63.5	52.1-55.9
M	25-27	22-24	63.5-68.6	55.9-61
L	27-29	24-26	68.6-73.7	61-66
XL	29-31	26-28	73.7-78.7	66-71.1

Custom sizes available.

Colors

Black

Beige

Fabrics

Powerknit

For sizing guide, go to page 48.



Pediatric Chin-Neck Garment

6-8006



Beige, Back

Pediatric



Features

Provides support and compression for the chin, upper neck, and cheeks.

- Soft, lightweight, breathable fabric
- 100% latex-free

Ideal for

- Facial surgery
- Neck surgery
- Submental region surgery
- Mentoplasty and/or chin implants
- Ear or periaural procedures
- Oral and maxillofacial surgery
- Pediatric facial surgery
- Scar revision

Head/Neck Sizing Chart

Pediatric 6-8006 & 6-8026

Size	Inches		Centimeters	
	Neck	Head	Neck	Head
Toddler				
S	9-10	19-20	22.9-25.4	48.3-50.8
M	10-11	20.5-21.5	25.4-27.9	52.1-54.6
L	11-12	22-24.5	27.9-30.5	55.9-62.2
XL	12-13	25-27	30.5-33	63.5-68.6
Child				
S	10-11	20-21	25.4-27.9	50.8-53.3
M	11.5-12.5	21.5-22.5	29.2-31.8	54.6-57.2
L	13-14	23-25.5	33-35.6	58.4-64.8
XL	14.5-15.5	26-28	36.8-39.4	66-71.1

Custom sizes available.

Colors

Beige

Fabrics

Tricot

For sizing guide, go to page 9.




Pediatric Bolero Style Arm Garment

6-8035



Pediatric



Features

Focuses compression on arms only.

- Outside seams for greater comfort
- Attach with Velcro® straps across the chest
- Graded compression from shoulders to wrists
- 100% latex-free

Ideal for

- Burn skin grafting and scar revision
- Pediatric reconstruction
- Congenital anomalies
- Tumor removal

Arm Sizing Chart (inches) Pediatric 6-8035

Size	Shoulder to Wrist	Mid-Upper Arm	Wrist
Toddler			
S	17.75	8.5-9.5	5.125-5.25
M	19.75	9.5-10.5	5.25-5.375
L	21.75	10.5-11.5	5.375-5.5
XL	23.75	11.5-12.5	5.5-5.625
Child			
S	19.75	9.5-10.5	5.375-5.5
M	21.75	10.5-11.5	5.5-5.625
L	23.75	11.5-12.5	5.625-5.75
XL	25.75	12.5-13.5	5.75-5.875

In inches. Custom sizes available.

Colors

Black

Beige

Fabrics

Satin

For sizing guide, go to page 16.

CAROMED®

Additional Information

Garment Care Instructions

When applying and caring for your Caromed compression garment at home, please follow these instructions carefully. Your garment is a medical device, designed and manufactured to help provide you with optimal cosmetic results.

Caromed recommends that someone assist you with your at-home garment application. Again, follow the instructions one step at a time, and take your time.

Application of your Caromed garment is a step-by-step process. Caromed garments are designed for ease of application; however, since this garment is a medical device, we ask that you follow all instructions carefully for best results.

HELPFUL TIPS

To apply your Caromed garment at home: If you have just showered, make certain that you are completely dry before applying your Caromed garment. Assume a seated position on a bed or comfortable chair. If at any time during the application of your girdle, you become fatigued, please take a minute to rest.

GIRDLES

While seated, fasten the girdle up to the knee area. Secure both the hook & eyes closures and zippers.

Remain seated, step into the girdle and position the girdle properly, gently pulling the girdle into place. While standing propped near your chair or bed, pull the girdle up making sure to position the crotch opening as far up as possible.

Fasten the remaining hook & eye closures on one side only. Carefully pull up the zipper. Now fasten the hook & eye closures on the remaining side. This side will be somewhat more difficult to secure and we recommend that someone be present to provide assistance.

After the hook & eye closures are secured, carefully zip up this side.

PLEASE USE EXTRA CAUTION TO AVOID CATCHING THE ZIPPER IN THE INNER FACING.

Important: If the zipper becomes snagged – **STOP IMMEDIATELY. DO NOT** attempt to force the zipper in any direction. Remove the entire garment completely. This alleviates the tension placed on the zipper. You now will be able to gently work with the zipper until it is freed from the facing.

Urination: Caromed recommends that all female patients purchase a small, flexible plastic funnel to position close to the perineum when voiding. This will prevent accidental splashing or spraying of the garment. Voiding in the standing position will probably be more comfortable the first few days following surgery.

Elimination: We recommend that you unzip one side, drop the flap, hold securely to one side and elimination can be accomplished quite easily without soiling the garment.

Laundering: Caromed recommends that you hand wash your Caromed garment in lukewarm water with a mild soap or detergent, then air dry. Please do not dry your garment in an electric or gas dryer as this may ruin the garment.

MOST PATIENTS REQUIRE TWO GARMENTS IN ORDER TO HAVE A FRESHLY LAUNDERED ONE AVAILABLE

To remove bloodstains, simply soak your garment in a basin of cool water to which 1/4 cup of peroxide has been added. Then launder and let air-dry.

CAROMED®

Custom Garment Order Form

Fax this form to:
888.311.7596

For sales, call **800-948-1938**

Physician Name _____ Date _____
 Phone Number _____ Fax Number _____
 Account # _____ Customer Order # _____
 Surgery Date _____ Date Garment(s) Needed _____
 Custom Items Ordered _____
 Procedure _____
 Patient Height _____ Patient Weight _____
 Patient Name _____

MEASUREMENTS

Refer to diagrams at the end of each section in order to assure that the areas measured are accurate, needed for your specific garment and correspond directly to each letter.

Measurement	Description	Custom Measurements	
A	Circumference of head measured horizontally		
B	Circumference of head measured vertically		
C	Circumference of neck		
D	Shoulder width		
E	Shoulder to natural waistline		
F	Bust (fullest part)		
G	Circumference of chest at inframammary fold		
H	Waistline		
I	Waist to crotch		
J	Crotch (measure from natural waistline through crotch area to back of waistline)		
K	Hip (fullest part)		
L	Natural waistline to knee		
M	Inseam (measure from upper inner thigh to ankle)	Left:	Right:
N	Circumference of upper thigh	Left:	Right:
O	Circumference of mid-thigh	Left:	Right:
P	Circumference of leg (2" above the knee)	Left:	Right:
Q	Circumference of knee	Left:	Right:
R	Circumference of leg (2" below the knee)	Left:	Right:
S	Circumference of upper calf	Left:	Right:
T	Circumference of mid-calf	Left:	Right:
U	Circumference of ankle	Left:	Right:
V 1.0*	Arm length measurement – shoulder to elbow		
V 1.1*	Arm length measurement – shoulder to wrist		
W 1.0*	Arm circumferential measurement – upper arm deltoid (next to shoulder)		
W 1.1*	Arm circumferential measurement – mid-upper arm		
W 1.2*	Arm circumferential measurement – elbow		
W 1.3*	Arm circumferential measurement – mid-forearm		
W 1.4*	Arm circumferential measurement – wrist		
X	Arm circumferential measurement – Shoulder to shoulder		
Y	Male chest size		
Y	Bra and cup size		

* To order custom arm bandage #5-8030, all measurements under sections V and W are required.



Take Two

One to wear
One to wash

Because of occasional bleeding after surgery, purchasing two garments allows one to be worn while the other is washed.

Caromed Fabrics




All Caromed garments are latex-free to prevent allergic reactions.

Caromed uses premium fabrics to provide the best results. Specially designed for Caromed for their compression properties, these fabrics are soft and luxurious to feel and see.

Powerknit

Powerknit is Caromed's most breathable fabric, with a mesh weave and sheer appearance.



Black

Champagne

Pink

Satin

Satin has a soft, elegant feel and provides a shimmer and more coverage than Powerknit.




Black

Champagne

Pink

Cotton

Cotton is soft and natural, but provides less compression and stain fighting than the others.



White


Black

Tricot

Tricot is our softest fabric (although slightly less compression), and is used exclusively on our face garments.



Beige



Made in USA

Caromed is proud to manufacture all of our products in our facility in Raleigh, North Carolina, USA. We started the industry of post-surgical compression garments and we assure that American ingenuity guides our continued success and service to our customers.

CAROMED[®]
recover beautifully