





Innovative. Intuitive. Intelligent. IntelliCart

Zimmer Biomet is committed to providing advanced, easy-to-use surgical equipment solutions for your most challenging OR needs. The IntelliCart System delivers on this commitment with significant advancements in day-to-day OR fluid and smoke waste management technology.

Innovative

- Redesigned vacuum pump with greater usable vacuum flow at a sound level of only 47.5dB (3.87 sones)
- Dual 17-liter lighted reservoirs with ConSeal™ Tint Technology selectively hide surgical fluids
- Clog-free manifold with integrated high-flow filter protects against vacuum system moisture ingress





Intuitive

- Simple and interactive, high-definition touch screen display
- Flexible docking system simplifies Cart processing in any dirty utility room
- PlumeVac[™] Smoke Evacuator available on or off the cart



Intelligent

- Reservoir doors close
 when manifolds are removed
- Reusable reservoirs and recyclable manifolds reduce OR red bag waste by up to 70%
- Advanced cleaning technology uses hot water, enzymatic, and an EPA approved disinfectant
- Backward compatible with UltrafleX Fluid Cart Manifolds and Ultra Evac Station



Duo Fluid Cart

With market leading 34-liter fluid capacity, extra quiet vacuum pump, clog-free suction manifolds and portable smoke evacuation, the IntelliCart System makes managing surgical fluid and smoke waste hazards simple and convenient.

- Control Panel Display
 Touch operation with
 310° rotation
- 2 IV Pole Manual pole standard, power IV pole shown
- Manifolds
 Clog-free and recyclable
- Reservoir Windows
 Equipped with proprietary
 ConSeal Tint Technology
- Redesigned Vacuum Pump
 Extra quiet with HEPA air
 filtration, increased flow
 capabilities and odor control
- 6 Brake Wheels
 Stabilize cart



Reservoir windows with lights on/tint off



Evac Station Connected to Duo Fluid Cart

The IntelliCart System takes infection control seriously. That is why the Evac Station uses an EPA approved disinfectant when it cleans and processes the Duo Fluid Cart. The simple to install station is placed in a surgical decontamination room or utility closet.

- Control Panel Display
 Customizable Evac operation
- 2 External Equipment Surfaces
 Withstands a wide variety
 of disinfectant wipes
- Manifold Storage
 Easy access to new manifolds
- 4 Evac Coupler
 Secure and orientation-independent docking
- PlumeVac Smoke Evacuator
 Optional and portable
- Evac Cleaning Solutions
 EPA approved disinfectant and
 Ultra Enzyme Cleaner

Evac Station 00514010900



Duo Fluid Cart with Smoke Evacuation 00514010200

PlumeVac[™] Smoke Evacuator

The portable PlumeVac Smoke Evacuator gives you the flexibility to put smoke plume evacuation where you need it – on or off the cart.



IntelliCart System Equipment

Duo Fluid Cart	00514010100
Duo Fluid Cart with Smoke Evacuation	00514010200

IntelliCart System Accessories

Power IV Pole	00514012000
Single-Use Manifold-24PK	00514020100
Specimen Trap-10PK	ULST100
Ultra Clot Buster-4oz Bottle, 12PK	ULCB100
Ultra Enzyme-1.0 gal, 4PK	ULEZ200
Verification Log Kit	TPVL200



Evac Station

00514010900

Smoke Evacuation & Accessories			
	PlumeVac Smoke Evacuator	AEVV120	
	PlumeGuard ULPA Filter-2PK	AEVS35302	
	Automatic Activation Device	AAD01	
	Smoke Plume Filter-10PK	ULSF100	



Single Use Manifold 00514020100



Do you have CSTAR^{\otimes} Protection on your Zimmer Biomet Fluid Waste capital equipment? Learn more at service.zimmer.com.



The Zimmer Biomet Clinical Excellence Team consists of PeriOperative Nurses who provide product education, training and support for the IntelliCart System. For more information about the Clinical Excellence Team please contact SurgicalSales-ClinicalExcellence@zimmerbiomet.com.

Contact your Zimmer Biomet representative or visit us at www.zimmerbiomet.com



This documentation is intended exclusively for physicians and is not intended for laypersons. Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendations. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice in whole or in part. Please refer to the package inserts for important product information, including, but not limited to, indications, contraindications, warning, precautions, and adverse effects.

©2016 Zimmer Biomet