

Osteoarthritis Braces - Single Uprights

- Bledsoe Brace™ 210
- Breg® 211
- DonJoy (DJO) 211
- Ossur Americas® 212-213
- Trulife 213
- VQ Orthocare® 213-214

Osteoarthritis Braces - Double Uprights

- Bledsoe Brace™ 214
- Breg® 215-216
- DonJoy (DJO) 216-217
- Medi USA 217
- Ossur Americas® 218

Hyperextension

- Allard USA 218-219
- DeRoyal® 220
- Trulife 220

Ligament Bracing

- Bledsoe Brace™ 221
- Breg® 222-223
- DonJoy (DJO) 223-226
- Medi USA 226
- Ossur Americas® 227-228
- Professional Products, Inc 229
- Trulife 229

ROM - Post Operative/Rehab

- Bledsoe Brace™ 230-232
- Breg® 233-234
- DonJoy (DJO) 235-237
- Ossur Americas® 237-238
- ProCare (DJO) 238-239
- RCAI 240

Knee Immobilizers

- CoreLINE™ 241
- AliMed® 241
- Bledsoe Brace™ 241
- Breg® 242
- Corflex® 242
- DeRoyal® 243
- DonJoy (DJO) 244
- Hely & Weber™ 244-245
- Med Spec 245
- Ossur Americas® 245
- ProCare (DJO) 246-247
- Professional Products, Inc 247
- RCAI® 247-248
- Trulife 248

Patellofemoral

- Bauerfeind 249
- Bledsoe Brace™ 249

- Breg® 249-250
- Corflex® 251
- DeRoyal® 252
- DonJoy (DJO) 253-256
- Hely & Weber™ 257-259
- Medi USA 260
- Med Spec 260
- New Options Sports 260-262
- Ossur Americas® 263
- Palumbo 263-264
- ProCare (DJO) 265-267
- Professional Products, Inc 267
- Truform® 267-268
- Trulife 268

Hinged Knee Braces

- Bauerfeind 269
- Bell-Horn® 269
- Bledsoe Brace™ 269
- Breg® 270
- Corflex® 271
- DeRoyal® 271
- DonJoy (DJO) 272-274
- Hely & Weber™ 274-277
- Med Spec 277
- New Options Sports 278
- Ossur Americas® 278-280
- ProCare (DJO) 280
- U.S. Orthotics 281

Knee Supports

- Alps 281
- Bauerfeind 281
- Bell-Horn® 282
- DeRoyal® 282
- DonJoy (DJO) 282-284
- Ossur Americas® 284
- ProCare (DJO) 284
- Truform® 285
- Trulife 285

Positional/Stretching Knee Orthoses

- AliMed® 286
- Anatomical Concepts 286
- DeRoyal® 286-287
- Lenjoy 287
- OCSI/Neuroflex 288
- RCAI® 288

Miscellaneous

- Breg® 288, 290
- DonJoy (DJO) 289, 291-292
- Knit-Rite® 289-290
- Cho-Pat® 290-291
- Hely & Weber™ 291
- ProCare (DJO) 292

SINGLE UPRIGHTS

Bledsoe Brace - Aligner ESE

- Applies the principles of the Thruster in a lower profile design
- Ideal for patients with mild to moderate osteoarthritis
- Patient adjustable hinge

Indications: Medial or lateral unicompartmental osteoarthritis



Size	Thigh Circum. (6" above mid-patella)	Calf Circum. (6" below mid-patella)	Medial		Lateral	
			Left	Right	Left	Right
X-Small	15"-17 ¹ / ₂ "	12"-13 ¹ / ₂ "	AR079101	AR079201	AR078101	AR078201
Small	17 ¹ / ₂ "-20"	13 ¹ / ₂ "-15"	AR079103	AR079203	AR078103	AR078203
Medium	20"-22 ¹ / ₂ "	15"-16 ¹ / ₂ "	AR079105	AR079205	AR078105	AR078205
Large	22 ¹ / ₂ "-25"	16 ¹ / ₂ "-18"	AR079107	AR079207	AR078107	AR078207
X-Large	25"-27 ¹ / ₂ "	18"-19 ¹ / ₂ "	AR079109	AR079209	AR078109	AR078209

Thruster RLF

- Correction increases as the leg goes into extension, then decreases as leg goes into flexion
- Pull versus push design eliminates pressure points and increases patient compliance
- Ideal for severe osteoarthritis
- Patient adjustable hinge

Indications: Medial or lateral unicompartmental osteoarthritis



Size	Thigh Circum. (6" above mid-patella)	Calf Circum. (6" below mid-patella)	Medial		Lateral	
			Left	Right	Left	Right
Standard						
X-Small	15"-17 ¹ / ₂ "	12"-13 ¹ / ₂ "	TA009101	TA009201	TA008101	TA008201
Small	17 ¹ / ₂ "-20"	13 ¹ / ₂ "-15"	TA009103	TA009203	TA008103	TA008203
Medium	20"-22 ¹ / ₂ "	15"-16 ¹ / ₂ "	TA009105	TA009205	TA008105	TA008205
Large	22 ¹ / ₂ "-25"	16 ¹ / ₂ "-18"	TA009107	TA009207	TA009107	TA008207
X-Large	25"-27 ¹ / ₂ "	18"-19 ¹ / ₂ "	TA009109	TA009209	TA008109	TA008209
Short						
X-Small	15"-17 ¹ / ₂ "	12"-13 ¹ / ₂ "	TA019101	TA019201	TA018101	TA018201
Small	17 ¹ / ₂ "-20"	13 ¹ / ₂ "-15"	TA019103	TA019203	TA018103	TA018203
Medium	20"-22 ¹ / ₂ "	15"-16 ¹ / ₂ "	TA019105	TA019205	TA018105	TA018205
Large	22 ¹ / ₂ "-25"	16 ¹ / ₂ "-18"	TA019107	TA019207	TA018107	TA018207
X-Large	25"-27 ¹ / ₂ "	18"-19 ¹ / ₂ "	TA019109	TA019209	TA018109	TA018209
Athletic						
X-Small	15"-17 ¹ / ₂ "	10 ¹ / ₂ "-12"	TA029101	TA029201	TA028101	TA028201
Small	17 ¹ / ₂ "-20"	12"-13 ¹ / ₂ "	TA029103	TA029203	TA028103	TA028203
Medium	20"-22 ¹ / ₂ "	13 ¹ / ₂ "-15"	TA029105	TA029205	TA028105	TA028205
Large	22 ¹ / ₂ "-25"	15"-16 ¹ / ₂ "	TA029107	TA029207	TA028107	TA028207
X-Large	25"-27 ¹ / ₂ "	16 ¹ / ₂ "-18"	TA029109	TA029209	TA028109	TA028209

Breg - Solus® Knee Brace

- ProForm technology provides a contoured fit and protection
- PivotPoint strap allow for comfort during motion
- Contoured frame padding that conforms to the leg
- Lightweight, high strength 2024 tempered aluminum

Indications: Unicompartmental OA medial; Chondral defect procedures; Meniscal procedures



Size	Thigh Circumference		Product No.	
	6" Above Mid-Patella	6" Below Mid-Patella	Left	Right
X-Small	13"-15"	10"-12"	15410	15610
Small	15"-18"	12"-14"	15420	15620
Medium	18"-19 ¹ / ₂ "	14"-15"	15430	15630
Medium +	19 ¹ / ₂ "-21"	15"-16"	15435	15635
Large	21"-24"	16"-18"	15440	15640
X-Large	24"-27"	18"-20"	15450	15650
2X-Large	27"-30"	20"-23"	15460	15660

DonJoy (DJO) - OA Lite™

- Single upright is malleable to allow further load adjustments with thigh and calf paddles to mitigate hinge rotation
- Pre-loaded hinge (up to 6° of correction) incorporates flexion and extension control
- Breathe-o-Prene® fabric eliminates skin irritation and provides comfort by wicking away moisture
- Brace hinge is on opposite side of affected compartment



Size	Thigh Circum. (6" above mid-patella)	Medial Compartment		Lateral Compartment	
		Left	Right	Left	Right
X-Small	13"-15 ¹ / ₂ "	11-1532-1	11-1531-1	11-1534-1	11-1533-1
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-1532-2	11-1531-2	11-1534-2	11-1533-2
Medium	18 ¹ / ₂ "-21"	11-1532-3	11-1531-3	11-1534-3	11-1533-3
Large	21"-23 ¹ / ₂ "	11-1532-4	11-1531-4	11-1534-4	11-1533-4
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-1532-5	11-1531-5	11-1534-5	11-1533-5
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-1532-6	11-1531-6	11-1534-6	11-1533-6
3X-Large	29 ¹ / ₂ "-32"	11-1532-7	11-1531-7	11-1534-7	11-1533-7


Osteoarthritis
Braces

SINGLE UPRIGHTS

Ossur Americas - Unloader® One

- Single upright with flexible upright that goes on same side as affected compartment
- For 5'1" and under, use short length
- Flexible shells with breathable liners coated with Sensil® silicone for improved comfort and reduced migration

Indications: Mild to severe unicompartmental osteoarthritis; Unicompartmental knee conditions requiring load reduction, such as articular cartilage defect repair, meniscal cartilage repair, avascular necrosis or tibial plateau fracture; Proportionate height-to-weight, moderately proportionate thigh-to-calf girth and moderate varus/valgus angulation; Low to moderate activities of daily living



Right Medial Unloader One with 3-Point Leverage System

Size	Leg Circum. (6" below mid-patella)	Medial Hinge		Lateral Hinge	
		Left	Right	Left	Right
Standard Length (15¹/₂")					
Small	12 ¹ / ₅ "-14"	B-240619712	B-240519712	B-240629712	B-240529712
Medium	14"-16 ³ / ₄ "	B-240619713	B-240519713	B-240629713	B-240529713
Large	16 ³ / ₄ "-20"	B-240619714	B-240519714	B-240629714	B-240529714
X-Large	20"-24"	B-240619715	B-240519715	B-240629715	B-240529715
Short Length (14¹/₂")					
X-Small	9 ⁴ / ₅ "-12 ¹ / ₅ "	B-241619711	B-241519711	B-241629711	B-241529711
Small	12 ¹ / ₅ "-14"	B-241619712	B-241519712	B-241629712	B-241529712
Medium	14"-16 ³ / ₄ "	B-241619713	B-241519713	B-241629713	B-241529713



Medial View of Right Medial Unloader One

Unloader® Spirit

- Recommended for mild to severe unicompartmental osteoarthritis
- Adjustable Dynamic Joint (ADJ) can accommodate a wide range of varus/valgus alignments
- Single upright with joint that goes on the same side as the affected compartment

Indications: Mild to severe unicompartmental osteoarthritis; Unicompartmental knee conditions requiring load reduction, such as articular cartilage defect repair, meniscal cartilage repair, avascular necrosis or tibial plateau fracture; Proportionate height-to-weight, moderately proportionate thigh-to-calf girth and moderate varus/valgus angulation; Low to moderate activities of daily living

Size	Thigh Circum. (6" above mid-patella)	Medial Hinge		Lateral Hinge	
		Left	Right	Left	Right
X-Small	13 ¹ / ₂ "-16"	B-217610101	B-217510101	B-217620101	B-217520101
Small	16"-18"	B-217610102	B-217510102	B-217620102	B-217520102
Medium	18"-20"	B-217610103	B-217510103	B-217620103	B-217520103
Large	20"-23 ¹ / ₂ "	B-217610104	B-217510104	B-217620104	B-217520104
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	B-217610105	B-217510105	B-217620105	B-217520105



Unloader Express

- Recommended for mild to severe unicompartmental osteoarthritis
- For low to moderate activities of daily living
- Dynamic Force Strap (DFS) applies a counteracting force to the adduction moment, considered to be the root cause of knee OA



Size	Thigh Circum. (6" above mid-patella)	Medial Hinge		Lateral Hinge	
		Left	Right	Left	Right
X-Small	13 ¹ / ₂ "-15 ¹ / ₂ "	B-215610001	B-215510001	B-215620001	B-215520001
Small	15 ¹ / ₂ "-17 ¹ / ₂ "	B-215610002	B-215510002	B-215620002	B-215520002
Medium	17 ¹ / ₂ "-20 ¹ / ₂ "	B-215610003	B-215510003	B-215620003	B-215520003
Large	20 ¹ / ₂ "-23 ¹ / ₂ "	B-215610004	B-215510004	B-215620004	B-215520004
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	B-215610005	B-215510005	B-215620005	B-215520005

Trulife - Lenox Hill Cruiser O.A.

- Silver, fiberglass frame
- Provides correction during extension and relief during flexion
- Unique stainless steel joint allows 10° of varus/valgus adjustment for customized fit
- Available with medial or lateral hinge (for OA on medial side, order lateral brace; for OA on lateral side, order medial brace)

Indications: Osteoarthritis of the knee; Varus/valgus correction



Size	Thigh Circum. (7" above knee)	Calf Circum. (7" below knee)	Medial Hinge		Lateral Hinge	
			Left	Right	Left	Right
X-Small	13"-15 ¹ / ₂ "	10 ¹ / ₂ "-12"	OALM1	OARM1	OALL1	OARL1
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	12"-13 ¹ / ₂ "	OALM2	OARM2	OALL2	OARL2
Medium	18 ¹ / ₂ "-21"	13 ¹ / ₂ "-15"	OALM3	OARM3	OALL3	OARL3
Large	21"-23 ¹ / ₂ "	15"-16 ¹ / ₂ "	OALM4	OARM4	OALL4	OARL4
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	16 ¹ / ₂ "-19"	OALM5	OARM5	OALL5	OARL5

VQ Orthocare - BioniCare® Knee System

- Clinically proven signal is delivered through advanced technology electrodes that avoid the need for conductive gel
- 4 ounce battery powered device which generates the proprietary electrical signal
- Low-profile, lightweight brace that provides 3 points of leverage to reduce pressure inside the joint
- 2 electrodes: One covers the knee and the other rests on the thigh to deliver the signal carried via wires from signal generator

Size	Knee Circumference	Product No.	
		Left	Right
Small	13"-15"	55-92003-201	55-92004-201
Medium	15"-18"	55-92003-301	55-92004-301
Large	18"-26"	55-92003-401	55-92004-401



Osteoarthritis
Braces

SINGLE UPRIGHTS

OActive Knee Brace

- Low-profile, lightweight design enhances comfort and compliance
- Unique design provides necessary leverage for medial or lateral compartment unloading
- Lateral application avoids contra-lateral limb and tibial crest pressure issues
- Quick release buckles for easy removal

Size	Patient Height	Product No.	
		Left	Right
Standard	up to 6'1"	08-12006-001	08-12008-001
Tall	6'1" +	08-12026-001	08-12028-001

Description	Product No.
Miscellaneous Items	
Night-Wrap System	55-90000-001
Night-Wrap Accessory Kit	55-92000-001
Electrodes for Night-Wrap (1 mo.)	55-94000-002
BioniCare Accessories	
Electrodes for OActive (Small) (1 mo.)	55-94000-201
Electrodes for OActive (Medium/Large) (1 mo.)	55-94000-351



DOUBLE UPRIGHTS

Bledsoe Brace - Z-12 OA

- Provides mild to moderate OA pain relief while stabilizing the knee and protecting the ACL, PCL, MCL and LCL
- OTS version has 3° of correction capability

Indications: Medial or lateral unicompartmental osteoarthritis while providing support of the ACL, PCL, MCL or LCL and menisci

Size	Thigh Circum. (6" above mid-patella)	Calf Circum. (6" below mid-patella)	Medial		Lateral	
			Left	Right	Left	Right
X-Small	13 ¹ / ₂ "-16"	12 ¹ / ₂ "-14"	AZ294101	AZ294201	AZ284101	AZ284201
Small	16"-18 ³ / ₄ "	14"-15 ¹ / ₂ "	AZ294103	AZ294203	AZ284103	AZ284203
Medium	18 ³ / ₄ "-21 ¹ / ₂ "	15 ¹ / ₂ "-17"	AZ294105	AZ294205	AZ284105	AZ284205
Large	21 ¹ / ₂ "-24 ¹ / ₄ "	17"-18 ¹ / ₂ "	AZ294107	AZ294207	AZ284107	AZ284207
X-Large	24 ¹ / ₄ "-27"	18 ¹ / ₂ "-20"	AZ294109	AZ294209	AZ284109	AZ284209
2X-Large	27"-29 ¹ / ₂ "	20"-21 ³ / ₈ "	AZ294111	AZ294211	AZ284111	AZ284211



Breg - Fusion® Lateral OA

- Innovative thumbwheel dial makes off-loading adjustments easy
- PivotPoint straps provide comfort during motion
- ProForm technology offers contoured fit and maximum protection
- Sleek frame provides unhindered mobility

Indications: Unicompartmental OA medial; Chondral defect procedures; Meniscal procedures

Size	Thigh Circumference		Product No.	
	6" Above Mid-Patella	6" Below Mid-Patella	Left	Right
X-Small	13"-15"	10"-12"	03710	03810
Small	15"-18"	12"-14"	03720	03820
Medium	18"-19½"	14"-15"	03730	03830
Medium +	19½"-21"	15"-16"	03735	03835
Large	21"-24"	16"-18"	03740	03840
X-Large	24"-27"	18"-20"	03750	03850
2X-Large	27"-30"	20"-23"	03760	03860



Fusion® Women's OA

- Unhindered mobility and enhanced comfort
- Ligament protection with sleek contoured frame design

Indications: Unicompartmental OA medial; Chondral defect procedures; Meniscal procedures

Size	Thigh Circumference		Product No.	
	6" Above Mid-Patella	6" Below Mid-Patella	Left	Right
X-Small	13"-15"	10"-12"	10810	10910
Small	15"-18"	12"-14"	10820	10920
Medium	18"-19½"	14"-15"	10830	10930
Medium +	19½"-21"	15"-16"	10835	10935
Large	21"-24"	16"-18"	10840	10940
X-Large	24"-27"	18"-20"	10850	10950
2X-Large	27"-30"	20"-23"	10860	10960



Fusion OA Knee Brace

- Low-profile brace provides pain relief
- Proform technology provides a contoured fit and protection
- Convenient thumbwheel dial hinge for easy adjustments

Indications: Unicompartmental OA medial; Chondral defect procedures; Meniscal procedures

Size	Thigh Circumference		Product No.	
	6" Above Mid-Patella	6" Below Mid-Patella	Left	Right
X-Small	13"-15"	10"-12"	12010	12110
Small	15"-18"	12"-14"	12020	12120
Medium	18"-19½"	14"-15"	12030	12130
Medium +	19½"-21"	15"-16"	12035	12135
Large	21"-24"	16"-18"	12040	12140
X-Large	24"-27"	18"-20"	12050	12150
2X-Large	27"-30"	20"-23"	12060	12160



DOUBLE UPRIGHTS

Compact X2K-OA (OTS)

- Enhanced adjustable hinge – 12 mm of travel (twice the loading capabilities as the standard adjustable hinge)
- Gel comfort condyle pad
- Siliconized strap padding with over-sized tibial frame pad

Indications: Unicompartamental OA; Chondral defect procedures; Meniscal procedures

Size	Thigh Circum. (6" above mid-patella)	Knee Joint Circum.	Compact		Women's	
			Left	Right	Left	Right
X-Small	13"-15 ¹ / ₂ "	11"-12 ³ / ₄ "	25310	25410	25510	25610
Small	15 ¹ / ₂ "-18"	12 ³ / ₄ "-14 ¹ / ₂ "	25320	25420	25520	25620
Medium	18"-19 ¹ / ₂ "	14 ¹ / ₂ "-15 ¹ / ₄ "	25330	25430	25530	25630
Medium+	19 ¹ / ₂ "-21"	15 ¹ / ₄ "-16 ¹ / ₄ "	25335	25435	25535	25635
Large	21"-24"	16 ¹ / ₄ "-18"	25340	25440	25540	25640
X-Large	24"-27"	18"-19 ³ / ₄ "	25350	25450	25550	25650
2X-Large	27"-31"	19 ³ / ₄ "-22"	25360	25460	25560	25660



Compact



Women's

DonJoy (DJO) - OA Adjuster™

- 4-Points-of-Leverage™ system
- Lock and Off-Load Technology™ shifts compressive knee joint forces from the degenerative area to the healthy compartment
- Universal OA key enables patient-controlled load adjustment to meet daily needs

Size	Thigh Circumference (6" above mid-patella)	Medial Compartment		Lateral Compartment	
		Left	Right	Left	Right
X-Small	13"-15 ¹ / ₂ "	11-0871-1	11-0870-1	11-0873-1	11-0872-1
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-0871-2	11-0870-2	11-0873-2	11-0872-2
Medium	18 ¹ / ₂ "-21"	11-0871-3	11-0870-3	11-0873-3	11-0872-3
Large	21"-23 ¹ / ₂ "	11-0871-4	11-0870-4	11-0873-4	11-0872-4
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-0871-5	11-0870-5	11-0873-5	11-0872-5
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-0871-6	11-0870-6	11-0873-6	11-0872-6
3X-Large	29 ¹ / ₂ "-32"	11-0871-7	11-0870-7	11-0873-7	11-0872-7



OA Everyday™

- 3 point load shifts compressive knee joint forces from the degenerative area to the healthy compartment
- Telescoping condyle pad (up to 6° of correction) enables patients to adjust load according to activity level
- Swiveling, quick release buckles facilitate ease of application/removal and accommodates leg movement
- Short brace length (12½") accommodates short inseams



Size	Thigh Circum. (6" above mid-patella)	Medial Compartment		Lateral Compartment	
		Left	Right	Left	Right
X-Small	13"-15½"	11-1340-1	11-1338-1	11-1341-1	11-1339-1
Small	15½"-18½"	11-1340-2	11-1338-2	11-1341-2	11-1339-2
Medium	18½"-21"	11-1340-3	11-1338-3	11-1341-3	11-1339-3
Large	21"-23½"	11-1340-4	11-1338-4	11-1341-4	11-1339-4
X-Large	23½"-26½"	11-1340-5	11-1338-5	11-1341-5	11-1339-5
2X-Large	26½"-29½"	11-1340-6	11-1338-6	11-1341-6	11-1339-6
3X-Large	29½"-32"	11-1340-7	11-1338-7	11-1341-7	11-1339-7

Sili-Grip Strap Pad Replacement Kit

- For OA Adjuster and OA Everyday

Product No. 11-1415

Medi USA - M4S OA Varus (Graphite)

- OA angle adjusters positioned close to the hinge provide easy pressure adjustment
- Flexible frame opposite the loading mechanism allows for greater patient comfort while seated and effective transfer of pressure while active
- Physioglide® hinge design closely mimics the roll and glide motion of the knee
- Lightweight Avional® alloy non-slip frame and flexible floating buckle system ensure high patient comfort and confidence with minimal migration



Size	Thigh Circumference (6" above mid-patella)	Product No.	
		Left	Right
X-Small	12¼"-14½"	62601	62501
Small	14½"-17¾"	62602	62502
Medium	17¾"-20"	62603	62503
Large	20"-22½"	62604	62504
X-Large	22½"-25½"	62605	62505
2X-Large	25½"-29"	62606	62506

Osteoarthritis
Braces

DOUBLE UPRIGHTS

Ossur Americas - OAsys™

- Double upright knee orthosis
- Acculign™ slide bar can be adjusted to provide medial or lateral compartment unloading
- Quick release buckles
- Provides functional support with or without varus/valgus alignment
- Orthosis comes with 10, 20, 30 and 40 degree extension stops, additional flexion stop kit may be purchased separately

Indications: Post-op, rehab and functional support plus unloading for protocols requiring reduced medial or lateral joint loads during recovery; ACL, MCL, PCL, LCL, rotary or combined instabilities; Low to moderate contact/impact activities



Acculign™ Slide Bar

Size	Caliper*	Circumference (just below distal border of the patella)	Product No.	
			Left	Right
Silver				
Small	3 ³ / ₈ "-3 ⁷ / ₈ "	11"-12 ¹ / ₄ "	OA3SL	OA3SR
Medium	3 ⁷ / ₈ "-4 ³ / ₈ "	12 ¹ / ₄ "-13 ⁷ / ₈ "	OA3ML	OA3MR
Large	4 ³ / ₈ "-5"	13 ⁷ / ₈ "-15 ⁷ / ₈ "	OA3LL	OA3LR
X-Large	5"-5 ¹ / ₂ "	15 ⁷ / ₈ "-17 ³ / ₈ "	OA3XLL	OA3XLR
Clear Carbon				
Small	3 ³ / ₈ "-3 ⁷ / ₈ "	11"-12 ¹ / ₄ "	OA3CSL	OA3CSR
Medium	3 ⁷ / ₈ "-4 ³ / ₈ "	12 ¹ / ₄ "-13 ⁷ / ₈ "	OA3CML	OA3CMR
Large	4 ³ / ₈ "-5"	13 ⁷ / ₈ "-15 ⁷ / ₈ "	OA3CLL	OA3CLR
X-Large	5"-5 ¹ / ₂ "	15 ⁷ / ₈ "-17 ³ / ₈ "	OA3CXLL	OA3CXLR

*M/L at joint width in a standard, weight-bearing position



Medial View of OAsys

HYPEREXTENSION

Allard USA - C.H.E.C.K. Comfortable Hyperextension Control Knee

- Dual axis joints allow for full knee flexion and allow each of the four uprights to move independently
- Low-profile hinges with slim uprights that are easily shaped to fit the contours of the leg
- One size fits knee widths 3¹/₂"-5¹/₄"

Description	Weight	Color	Product No.
Aluminum	16 oz	Gray	286500000
Stainless Steel	22 oz	Black	286520000
Replacement Pads	—	Gray	286510000




COMBO

- Low-profile and functional orthotic solution for management of genu recurvatum combined with foot drop
- Uprights easily shaped and adjusted to optimum height for patient
- Modular components allow you to accommodate variances in thigh and calf circumferences
- Weight Limit: 330 lbs; Height Limit: 6' 6"

Five Components Required to Build a COMBO:

1. Joints with Uprights - universal size
2. Tibial Connector - generally same size as a ToeOFF or BlueRocker
3. Knee Control with Popliteal Interface - determined by knee width at the knee
4. Thigh Cuff - Take circumference of thigh at location where cuff to be worn. Cannot be less than 4" to Tibia Plateau
5. Thigh Cuff SoftKIT - same as Thigh Cuff



Size	Product No.
1. Joints w/ Uprights	
Universal	287800000
2. Tibial Connector	
Small	287840011
Medium	287840012
Large/X-Large	287840013

Size	Knee Width	Product No.
3. Knee Control		
Small	2 ³ / ₄ "-4 ¹ / ₄ "	287820011
Medium	3 ³ / ₄ "-5 ¹ / ₄ "	287820012
Large	5"-6 ¹ / ₄ "	287820013

Size	Thigh Circum.	Medial/ Lateral	Product No.	
			Left	Right
4. Thigh Cuff				
Small	11 ³ / ₄ "-15 ³ / ₄ "	4 ¹ / ₄ "-6"	287811011	287812011
Medium	15"-19 ³ / ₄ "	5 ¹ / ₄ "-6"	287811012	287812012
Large	19"-23 ³ / ₄ "	6"-7 ¹ / ₂ "	287811013	287812013
X-Large	22 ³ / ₄ "-27 ¹ / ₂ "	6 ³ / ₄ "-8 ¹ / ₄ "	287811014	287812014

Size	Product No.
5. Thigh Cuff SoftKIT	
Small	287830011
Medium	287830012
Large	287830013
X-Large	287830014



Rear View

Hyperextension

DeRoyal - Hypercontrol® Knee Brace

- Tri-tex construction serves as a breathable alternative to neoprene
- Medial and lateral hinges help prevent hyperextension of the knee
- Anatomical design features unique strapping system at the back of the knee to assist in preventing hyperextension

Indications: Grade I or II collateral ligament sprain; Mild ACL/PCL sprain; Grade I varus and valgus instability

Size	Measurement	Product No.
Small	15 ¹ / ₂ "-18"	14910005
Medium	18"-20 ¹ / ₂ "	14910006
Large	20 ¹ / ₂ "-23"	14910007
X-Large	23"-25 ¹ / ₂ "	14910008
2X-Large	25 ¹ / ₂ "-28"	14910009
3X-Large	28"-30 ¹ / ₂ "	14910010



Trulife - Lenox Hill Swedish Knee Cage

- Lightweight vinyl coated aluminum frame
- Adjustable posterior pressure pad
- Permits knee flexion, restricts extension

Indications: Knee hyperextension, genu recurvatum

Size	Knee Circumference (4" above knee center)	Product No.
Small	10"-12"	0866802
Medium	13"-14"	0866803
Large	15"-16"	0866804
X-Large	17"-18"	0866805



Lenox Hill Hinged Swedish Knee Cage

- Lightweight, aluminum frame with padded straps
- Prevents hyperextension of knee joint
- Circumference taken 4" above knee center

Indications: Knee hyperextension, genu recurvatum

Size	Knee Circumference	Product No.
Small	11"-15"	A24700
Regular	15"-21"	A24701

Replacement Parts for Lenox Hill Hinged Swedish Knee Cage	
Description	Product No.
Top Strap, Regular (with ring attached)	011342
Anterior Upper Strap, Regular	011343
Lower Strap, Regular (loop and ring)	011345



Bledsoe Brace - Z-12

- Low-profile in a 13" length
- Quick-release buckles for easy application
- Pivoting gastroc strapping system addresses migration
- Strong and durable magnesium material

Indications: Provides stability and support following injury or reconstruction of the ACL, PCL, MCL or LCL and menisci



Size	Thigh Circum. (6" above mid-patella)	Calf Circum. (6" below mid-patella)	Product No.	
			Left	Right
X-Small	13 ¹ / ₂ "-16"	12 ¹ / ₂ "-14"	AZ114101	AZ114201
Small	16"-18 ³ / ₄ "	14"-15 ¹ / ₂ "	AZ114103	AZ114203
Medium	18 ³ / ₄ "-21 ¹ / ₂ "	15 ¹ / ₂ "-17"	AZ114105	AZ114205
Large	21 ¹ / ₂ "-24 ¹ / ₄ "	17"-18 ¹ / ₂ "	AZ114107	AZ114207
X-Large	24 ¹ / ₄ "-27"	18 ¹ / ₂ "-20"	AZ114109	AZ114209
2X-Large	27"-29 ¹ / ₂ "	20"-21 ³ / ₈ "	AZ114111	AZ114211

Axiom

- Pivoting gastroc strapping system addresses migration
- Aluminum, impact resistant
- Measure thigh circumference 6" above mid-patella

Indications: Provides stability and support following injury or reconstruction of the ACL, PCL, MCL, LCL and menisci; Also recommended for prophylactic use



Size	Thigh Circum	Calf Circum.	Product No.	
			Left	Right
Standard				
X-Small	13 ¹ / ₂ "-16"	12 ¹ / ₂ "-14"	PK114101	PK114201
Small	16"-18 ³ / ₄ "	14"-15 ¹ / ₂ "	PK114103	PK114203
Medium	18 ³ / ₄ "-21 ¹ / ₂ "	15 ¹ / ₂ "-17"	PK114105	PK114205
Large	21 ¹ / ₂ "-24 ¹ / ₄ "	17"-18 ¹ / ₂ "	PK114107	PK114207
X-Large	24 ¹ / ₄ "-27"	18 ¹ / ₂ "-20"	PK114109	PK114209
2X-Large	27"-29 ¹ / ₂ "	20"-21 ³ / ₈ "	PK114111	PK114211
3X-Large	29 ¹ / ₂ "-31"	21 ³ / ₈ "-22 ³ / ₄ "	PK114113	PK114213
Athletic				
X-Small	13 ¹ / ₂ "-16"	11"-12 ¹ / ₂ "	PK119101	PK119201
Small	16"-18 ³ / ₄ "	12 ¹ / ₂ "-14"	PK119103	PK119203
Medium	18 ³ / ₄ "-21 ¹ / ₂ "	14"-15 ¹ / ₂ "	PK119105	PK119205
Large	21 ¹ / ₂ "-24 ¹ / ₄ "	15 ¹ / ₂ "-17"	PK119107	PK119207
X-Large	24 ¹ / ₄ "-27"	17"-18 ¹ / ₂ "	PK119109	PK119209
2X-Large	27"-29 ¹ / ₂ "	18 ¹ / ₂ "-20"	PK119111	PK119211
3X-Large	29 ¹ / ₂ "-31"	20"-21 ³ / ₈ "	PK119113	PK119213

Breg - Fusion® Knee Brace

- AirTech™ frame pads have innovative windows and channels which permit cooling air flow while allowing moisture and heat to move away from the skin
- Pivot Point Strap tabs allow strap to contour to the leg for superior fit, comfort and suspension
- ProForm Technology - medial side structure designed for superior contoured fit and protection

Indications: ACL; PCL; Collateral ligament protection

Size	Thigh Circumference		Fusion		Fusion Women's	
	6" Above Mid-Patella	6" Below Mid-Patella	Left	Right	Left	Right
X-Small	13"-15"	10"-12"	00510	00710	00910	01110
Small	15"-18"	12"-14"	00520	00720	00920	01120
Medium	18"-19 ¹ / ₂ "	14"-15"	00530	00730	00930	01130
Medium +	19 ¹ / ₂ "-21"	15"-16"	00535	00735	00935	01135
Large	21"-24"	16"-18"	00540	00740	00940	01140
X-Large	24"-27"	18"-20"	00550	00750	00950	01150
2X-Large	27"-30"	20"-23"	00560	00760	00960	01160



Fusion Knee Brace



Fusion Women's

Fusion® XT

- AirTech™ frame pads have innovative windows and channels which permit cooling air flow while allowing moisture and heat to move away from the skin
- Pivot Point Strap tabs allow straps to contour to the leg for superior fit, comfort and suspension

Indications: ACL; PCL; Collateral ligament protection

Size	Thigh Circumference		Product No.	
	6" Above Mid-Patella	6" Below Mid-Patella	Left	Right
Small	15"-18"	12"-14"	00620	00820
Medium	18"-19 ¹ / ₂ "	14"-15"	00630	00830
Medium +	19 ¹ / ₂ "-21"	15"-16"	00635	00835
Large	21"-24"	16"-18"	00640	00840
X-Large	24"-27"	18"-20"	00650	00850
2X-Large	27"-30"	20"-23"	00660	00860



X2K[®] Knee Brace

- Diamond design for optimal varus and valgus stiffness
- Lightweight frame for reduced migration
- Adjustable hinges allow for optimal loading across the joint
- Internally mounted straps for ideal brace suspension

Indications: ACL; PCL; Collateral ligament protection



Size	Thigh Circum. (6" above mid-patella)	Knee Joint Circum.	Product No.	
			Left	Right
X-Small	13"-15 ¹ / ₂ "	11"-12 ³ / ₄ "	20110	20210
Small	15 ¹ / ₂ "-18"	12 ³ / ₄ "-14 ¹ / ₂ "	20120	20220
Medium	18"-19 ¹ / ₂ "	14 ¹ / ₂ "-15 ¹ / ₄ "	20130	20230
Medium +	19 ¹ / ₂ "-21"	15 ¹ / ₄ "-16 ¹ / ₄ "	20135	20235
Large	21"-24"	16 ¹ / ₄ "-18"	20140	20240
X-Large	24"-27"	18"-19 ³ / ₄ "	20150	20250
2X-Large	27"-31"	19 ³ / ₄ "-22"	20160	20260

DonJoy (DJO) - FULLFORCE[™]

- 4-Points-Of-Leverage system provides ligament stability and reduces ACL strain
- Clinically-proven FourcePoint hinge technology
- Internally mounted swiveling strap tabs help accommodate leg movement
- Swooping thigh cuff for inner thigh clearance and bilateral brace use

Indications: Provides ligament stability and significantly reduces ACL strain; Increases flexion angles and reduces anterior shear forces on the tibia



Size	Thigh Circumference (6" above mid-patella)	Product No.			
		ACL Left	ACL Right	CI Left	CI Right
Short Calf Length					
X-Small	13"-15 ¹ / ₂ "	11-3221-1	11-3220-1	11-3223-1	11-3222-1
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-3221-2	11-3220-2	11-3223-2	11-3222-2
Medium	18 ¹ / ₂ "-21"	11-3221-3	11-3220-3	11-3223-3	11-3222-3
Large	21"-23 ¹ / ₂ "	11-3221-4	11-3220-4	11-3223-4	11-3222-4
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-3221-5	11-3220-5	11-3223-5	11-3222-5
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-3221-6	11-3220-6	11-3223-6	11-3222-6
3X-Large	29 ¹ / ₂ "-32"	11-3221-7	11-3220-7	11-3223-7	11-3222-7
Standard Calf Length					
X-Small	13"-15 ¹ / ₂ "	11-0259-1	11-0258-1	11-0265-1	11-0264-1
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-0259-2	11-0258-2	11-0265-2	11-0264-2
Medium	18 ¹ / ₂ "-21"	11-0259-3	11-0258-3	11-0265-3	11-0264-3
Large	21"-23 ¹ / ₂ "	11-0259-4	11-0258-4	11-0265-4	11-0264-4
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-0259-5	11-0258-5	11-0265-5	11-0264-5
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-0259-6	11-0258-6	11-0265-6	11-0264-6
3X-Large	29 ¹ / ₂ "-32"	11-0259-7	11-0258-7	11-0265-7	11-0264-7

ACL Everyday

- Ergonomic, swooping thigh and calf cuffs provide clearance for the contralateral leg
- Contoured, flexible cuff shells conform to patient anatomy
- 10° soft stops to reduce impact noise
- Short brace (13") accommodates short seams

Indications: Mild ACL/MCL/LCL instabilities; Support for the activities of daily living

Size	Thigh Circumference (6" above mid-patella)	Product No.	
		Left	Right
X-Small	13"-15 ¹ / ₂ "	11-1600-1	11-1601-1
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-1600-2	11-1601-2
Medium	18 ¹ / ₂ "-21"	11-1600-3	11-1601-3
Large	21"-23 ¹ / ₂ "	11-1600-4	11-1601-4
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-1600-5	11-1601-5
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-1600-6	11-1601-6
3X-Large	29 ¹ / ₂ "-32"	11-1600-7	11-1601-7



Legend™

- 4 points-of-leverage system
 - Low-profile (0.080 - 0.125" thick), lightweight aircraft grade aluminum frame with durable, non-chip Kraton coating
 - Extra set of suede frame liners included
 - Available in Black
 - Replace ? with size (1-7) from sizing chart
- Indications:** Moderate to severe ACL, PCL, MCL, and LCL instabilities

Description	Product No.	
	Left	Right
ACL	11-0423-?-06060	11-0422-?-06060
PCL	11-0443-?-06060	11-0442-?-06060
CI	11-0453-?-06060	11-0452-?-06060



Legend

4Titide®

- FourcePoint technology to protect the ACL
- Malleable thigh and calf cuffs for an improved fit
- Low-profile (0.080-0.125" thick), lightweight aircraft grade aluminum frame
- Replace ? with size (1-7) from sizing chart

Indications: Mild to severe ACL, PCL, MCL, and LCL instabilities; Solid support for the activities of daily living or non-contact sports

Description	Product No.	
	Left	Right
ACL	11-0831-?-06060	11-0830-?-06060
PCL	11-0833-?-06060	11-0832-?-06060
CI	11-0835-?-06060	11-0834-?-06060
ACL Short	11-0837-?-06060	11-0836-?-06060
PCL Short	11-0839-?-06060	11-0838-?-06060
CI Short	11-0841-?-06060	11-0840-?-06060



4Titide

Legend and 4Titide Sizing Chart

Size	Thigh Circumference (6" above mid-patella)
1 (X-Small)	13"-15 ¹ / ₂ "
2 (Small)	15 ¹ / ₂ "-18 ¹ / ₂ "
3 (Medium)	18 ¹ / ₂ "-21"
4 (Large)	21"-23 ¹ / ₂ "
5 (X-Large)	23 ¹ / ₂ "-26 ¹ / ₂ "
6 (2X-Large)	26 ¹ / ₂ "-29 ¹ / ₂ "
7 (3X-Large)	29 ¹ / ₂ "-32"

Armor

- Clinically proven 4 Points-of-Leverage system
- Swooping thigh cuff for inner thigh clearance and bilateral brace use
- Strong aircraft grade aluminum frame (0.125" thick)

Indications: Mild to severe ACL instabilities; Prophylactic/bilateral use; MCL/LCL instabilities; ACL reconstructions



Size	Thigh Circum. (6" above mid-patella)	Standard Calf		Short Calf	
		Left	Right	Left	Right
FourcePoint Hinge					
X-Small	13"-15 ¹ / ₂ "	11-1443-1	11-1442-1	11-1441-1	11-1440-1
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-1443-2	11-1442-2	11-1441-2	11-1440-2
Medium	18 ¹ / ₂ "-21"	11-1443-3	11-1442-3	11-1441-3	11-1440-3
Large	21"-23 ¹ / ₂ "	11-1443-4	11-1442-4	11-1441-4	11-1440-4
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-1443-5	11-1442-5	11-1441-5	11-1440-5
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-1443-6	11-1442-6	11-1441-6	11-1440-6
3X-Large	29 ¹ / ₂ "-32"	11-1443-7	11-1442-7	11-1441-7	11-1440-7

Size	Thigh Circum. (6" above mid-patella)	Standard Calf		Short Calf	
		Left	Right	Left	Right
Standard Hinge					
X-Small	13"-15 ¹ / ₂ "	11-1445-1-26000	11-1444-1-26000	11-1005-1-26000	11-1004-1-26000
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-1445-2-26000	11-1444-2-26000	11-1005-2-26000	11-1004-2-26000
Medium	18 ¹ / ₂ "-21"	11-1445-3-26000	11-1444-3-26000	11-1005-3-26000	11-1004-3-26000
Large	21"-23 ¹ / ₂ "	11-1445-4-26000	11-1444-4-26000	11-1005-4-26000	11-1004-4-26000
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-1445-5-26000	11-1444-5-26000	11-1005-5-26000	11-1004-5-26000
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-1445-6-26000	11-1444-6-26000	11-1005-6-26000	11-1004-6-26000
3X-Large	29 ¹ / ₂ "-32"	11-1445-7-26000	11-1444-7-26000	11-1005-7-26000	11-1004-7-26000

FourcePoint

- Clinically proven 4 Points-of-Leverage system
- Patented FourcePoint technology is clinically proven to protect ACL
- Moldable thigh and calf cuffs for an improved fit

Indications: Ideal for patients returning to sports specific activity after ACL reconstruction, an ACL deficient knee or for prophylaxis



Size	Thigh Circumference (6" above mid-patella)	Standard Calf		Short Calf	
		Left	Right	Left	Right
X-Small	13"-15 ¹ / ₂ "	11-1301-1	11-1300-1	11-1303-1	11-1302-1
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-1301-2	11-1300-2	11-1303-2	11-1302-2
Medium	18 ¹ / ₂ "-21"	11-1301-3	11-1300-3	11-1303-3	11-1302-3
Large	21"-23 ¹ / ₂ "	11-1301-4	11-1300-4	11-1303-4	11-1302-4
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-1301-5	11-1300-5	11-1303-5	11-1302-5
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-1301-6	11-1300-6	11-1303-6	11-1302-6
3X-Large	29 ¹ / ₂ "-32"	11-1301-7	11-1300-7	11-1303-7	11-1302-7

Female Force™ ACL

- Lightweight, low-profile brace designed specifically for females
- ForcePoint hinges react during leg extension and offer an increasing resistance in the last 35°
- Reduces extension shock and increases posterior tibial loads
- Medial swooping thigh frame for increased medial clearance

Indications: Addresses biomechanical factors that have been implicated in the high incidence of ACL tears in females

Size	Thigh Circumference (6" above mid-patella)	Left	Right	Anti Migration Calf Strap
X-Small	13"-15 ¹ / ₂ "	11-1431-1	11-1430-1	11-2005-1
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-1431-2	11-1430-2	11-2005-2
Medium	18 ¹ / ₂ "-21"	11-1431-3	11-1430-3	11-2005-3
Large	21"-23 ¹ / ₂ "	11-1431-4	11-1430-4	11-2005-4
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-1431-5	11-1430-5	11-2005-5
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-1431-6	11-1430-6	11-2005-6
3X-Large	29 ¹ / ₂ "-32"	11-1431-7	11-1430-7	11-2005-7



Tru-Pull Advanced Patellofemoral Attachment

- Can be utilized on DonJoy's rigid ligament and OA braces for additional support
- Elastomeric pull device provides a strong dynamic pull on the patella as the knee moves into extension
- Choice of 1/4" and 3/8" trimmable buttresses for versatility and comfort

Indications: Symptomatic pain relief of patellofemoral dysfunction including lateral maltracking; Subluxation/dislocation; Patellar tendonitis

Size	Product No.
X-Small/Small	11-1525-2
Medium/Large/X-Large	11-1525-3
2X-Large/3X-Large	11-1525-4



Medi USA - M4S Functional Knee Brace

- S-shape follows anatomy of the femur and tibia; along with the Physioglide® hinge, it closely matches knee anatomy/movement
- Extremely lightweight Avional® alloy frame and flexible floating buckle system ensure high patient comfort and confidence with minimal migration
- Non-slip padding grips without compromising skin condition
- Measure thigh circumference 6" above mid-patella

Size	Product No.	
	Left	Right
X-Small	62451	62351
Small	62452	62352
Medium	62453	62353
Large	62454	62354
X-Large	62455	62355
2X-Large	62456	62356



Ossur Americas - Paradigm® OTS

- Flared and contoured medial thigh shell with flexible subshell
- Calf shell with adjustable strap that captures the tibia
- “Cage” frame to provide total ligament support through multiple points of leverage

Indications: ACL, PCL, MCL, LCL, rotary and combined instabilities; Protection and stabilization of surgical grafts; Recommended for low to medium contact/impact

Size	Thigh Circumference (6" above mid-patella)	Product No.	
		Left Leg	Right Leg
X-Small	13 ¹ / ₂ "-15 ¹ / ₂ "	B-210600001	B-210500001
Small	15 ¹ / ₂ "-17 ¹ / ₂ "	B-210600002	B-210500002
Medium	17 ¹ / ₂ "-20 ¹ / ₂ "	B-210600003	B-210500003
Large	20 ¹ / ₂ "-23 ¹ / ₂ "	B-210600004	B-210500004
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	B-210600005	B-210500005
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	B-210600006	B-210500006



Medial View of Paradigm

Flex®

- Gel condylar pads improve suspension and comfort
- Sport model shaped for lean, athletic legs
- Adjustable Tibial Strap (ATS) captures the tibia and allows adjustments in hinge depth for a more customized fit
- Comes with 0°, 10°, 20°, 30° and 40° stops, additional flexible stop kit may be purchased

Indications: ACL, PCL, MCL, LCL, rotary and combined instabilities; Long term use for post-op rehab and functional support; Recommended for low to medium contact/impact

Size	Caliper*	Leg Circum.**	Flex		Flex Sport	
			Left Leg	Right Leg	Left Leg	Right Leg
Small	3 ¹ / ₂ "-4"	11"-12 ⁵ / ₈ "	FLEX-SL	FLEX-SR	FLEXSP-SL	FLEXSP-SR
Medium	4"-4 ⁵ / ₈ "	12 ⁵ / ₈ "-14 ¹ / ₂ "	FLEX-ML	FLEX-MR	FLEXSP-ML	FLEXSP-MR
Large	4 ⁵ / ₈ "-5 ¹ / ₈ "	14 ¹ / ₂ "-16"	FLEX-LL	FLEX-LR	FLEXSP-LL	FLEXSP-LR
X-Large	5 ¹ / ₈ "-5 ⁵ / ₈ "	16"-17 ⁵ / ₈ "	FLEX-XLL	FLEX-XLR	FLEXSP-XLL	FLEXSP-XLR

*M/L at joint width in a standard, weight-bearing position

**Just below the distal border of the patella



Ligament
Bracing

LIGAMENT BRACING

CTi® OTS

- Carbon frame provides a rigid exoskeleton that stabilizes the knee joint
- Anatomically correct Accutrac® hinges with extension stops
- Breathable liners coated with bio-inert Sensil™ silicone for increased suspension

Indications: ACL, PCL, MCL, LCL, rotary and combined instabilities; Standard model for medium to high impact levels; Pro Sport model for highest impact level



Silicone padded liners



Size	Caliper*	Leg Circum.**	Product No.			
			PCL Version		Non-PCL Version	
			Left	Right	Left	Right
Small	3 ¹ / ₂ "-4"	11"-12 ⁵ / ₈ "	B-238640112	B-238560112	B-238600112	B-238500112
Medium	4"-4 ¹ / ₂ "	12 ⁵ / ₈ "-14 ¹ / ₄ "	B-238640113	B-238560113	B-238600113	B-238500113
Large	4 ¹ / ₂ "-4 ⁷ / ₈ "	14 ¹ / ₄ "-15 ³ / ₈ "	B-238640114	B-238560114	B-238600114	B-238500114
X-Large	4 ⁷ / ₈ "-5 ¹ / ₄ "	15 ³ / ₈ "-16 ⁵ / ₈ "	B-238640115	B-238560115	B-238600115	B-238500115
2X-Large	5 ¹ / ₄ "-5 ³ / ₄ "	16 ⁵ / ₈ "-18 ¹ / ₈ "	B-238640116	B-238560116	B-238600116	B-238500116

*M/L at joint width in a standard, weight-bearing position **Just below the distal border of the patella

CTi OTS Pro Sport

- Same features as CTi OTS Standard model, except designed for highest impact levels

Size	Caliper*	Leg Circum.**	Product No.			
			PCL Version		Non-PCL Version	
			Left	Right	Left	Right
Small	3 ¹ / ₂ "-4"	11"-12 ⁵ / ₈ "	B-239640112	B-239560112	B-239600112	B-239500112
Medium	4"-4 ¹ / ₂ "	12 ⁵ / ₈ "-14 ¹ / ₄ "	B-239640113	B-239560113	B-239600113	B-239500113
Large	4 ¹ / ₂ "-4 ⁷ / ₈ "	14 ¹ / ₄ "-15 ³ / ₈ "	B-239640114	B-239560114	B-239600114	B-239500114
X-Large	4 ⁷ / ₈ "-5 ¹ / ₄ "	15 ³ / ₈ "-16 ⁵ / ₈ "	B-239640115	B-239560115	B-239600115	B-239500115
2X-Large	5 ¹ / ₄ "-5 ³ / ₄ "	16 ⁵ / ₈ "-18 ¹ / ₈ "	B-239640116	B-239560116	B-239600116	B-239500116

*M/L at joint width in a standard, weight-bearing position **Just below the distal border of the patella

Size	Width	Product No.
CTi Strap Pads (Black)		
Small	3 ¹ / ₂ "	B-304010057
Medium/Large	5"	B-304010058
X-Large/2X-Large	7"	B-304010059
CTi Accessories		
Patella Cup		B-711130000

Size	Product No.	
	Left	Right
CTi Silicone Liners		
Small/Medium	B-411600112	B-411500112
Large	B-411600114	B-411500114
X-Large	B-411630114	B-411530114

Professional Products - Wrap-Around Single Axis Knee Supports

- Wrap around design for ease of application
- Black Lycra covered 1/8" neoprene with black nylon inner lining
- Removable Single Axis Hinge with extension stops at 0°, 15°, 30°, and 45°
- Open patella with open popliteal

Indications: Mild medial/lateral collateral ligament instability

Product No. **05427-T8-C2**

Size	Above Knee Measurement	Below Knee Measurement
X-Small	11 ¹ / ₂ "-13 ¹ / ₂ "	10"-12"
Small	13 ¹ / ₂ "-15 ¹ / ₂ "	12"-13 ¹ / ₂ "
Medium	15 ¹ / ₂ "-17 ¹ / ₂ "	13 ¹ / ₂ "-15"
Large	17 ¹ / ₂ "-20 ¹ / ₂ "	15"-17"
X-Large	20 ¹ / ₂ "-23 ¹ / ₂ "	17"-19"
2X-Large	23 ¹ / ₂ "-27 ¹ / ₂ "	19"-21 ¹ / ₂ "

Measure circumference from patella midline, 3" above and 3" below from a standing position against the skin



Trulife - Lenox Hill Precision Pro

- Low-profile, custom-fit aluminum brace
- 4-point knee fixation
- Polycentric joint simulates natural knee movement
- Numbered straps for easy donning and doffing
- If measurements are borderline between sizes, choose the smaller size

Indications: Mild medial/lateral collateral ligament instabilities; After ligament reconstruction; After meniscus reconstruction; After meniscus resection

Size	Thigh Circum.	Calf Circum.	Product No.	
			Left	Right
Small	15"-18 ¹ / ₂ "	12"-13"	PRECPRO-L1	PRECPRO-R1
Medium	18 ¹ / ₂ "-21"	13"-15"	PRECPRO-L2	PRECPRO-R2
Large	21"-23 ¹ / ₂ "	15"-16 ¹ / ₂ "	PRECPRO-L3	PRECPRO-R3
X-Large	23 ¹ / ₂ "-26"	16 ¹ / ₂ "-18 ¹ / ₂ "	PRECPRO-L4	PRECPRO-R4
2X-Large	26"-30"	18 ¹ / ₂ "-21"	PRECPRO-L5	PRECPRO-R5

Measure circumference 6" above and 6" below the knee center



Ligament
Bracing

Bledsoe Brace - Crossover ROM

- Easily adjustable ROM hinges for ideal anatomical placement
- Front closure wraparound design for easy application

Indications: Application to knees for locked or limited motion control during rehabilitation after operative procedures or injury

Size	Thigh Circumference (6" above mid-patella)	Product No.	
		Tritech	3D Neoprene
X-Small	13 ¹ / ₂ "-16"	RK420001	RK320001
Small	16"-18 ³ / ₄ "	RK420003	RK320003
Medium	18 ³ / ₄ "-21 ¹ / ₂ "	RK420005	RK320005
Large	21 ¹ / ₂ "-24 ¹ / ₄ "	RK420007	RK320007
X-Large	24 ¹ / ₄ "-27"	RK420009	RK090009
2X-Large	27"-29 ¹ / ₂ "	RK420011	RK320011



Original Post-Op

- Optional posterior and anterior fracture plates
- Easily converts to KAFO and HKAFO with Bledsoe's Hi-Top Boot and/or Simple Pelvic Braces

Indications: For immediate application to knees for locked or limited motion control of knee during rehabilitation after operative procedures or injury to knee ligaments, patellar ligament, articular cartilage or meniscus; Stable or internally fixed fractures of tibia plateau, condyles, proximal tibia or distal femur; Stable, non-displaced or internally fixed fractures of the middle to distal femur and proximal to middle tibia or fibula (when using the optional fracture bracing plates and shoe insert)

Leg Length	Brace Size	Product No.
17"-19"	18"	EK020018
19"-21"	20"	EK020020
21"-23"	22"	EK020022
23"-25"	24"	EK020024
25"-27"	26"	EK020026
27"-29"	28"	EK020028
29"-31"	30"	EK020030
31"-33"	32"	EK020032
33"-35"	34"	EK020034

Measure leg length from 1" below crotch to ankle



Basic Brace

- Straps are attached to hinge arms for quick connection
- Anti-migration strap retards slippage of calf cuff
- Range-of-motion can be set without tools
- Available with simple hinge or QuikLok lever to rigidly lock the brace at full extension

Indications: For immediate application to knee during rehabilitation after operative procedures or injury to knee ligaments, cartilage, or stable or internally fixed fractures

Size	Length	Product No.	
		Simple Hinge	QuikLok Hinge
Short	16"	AK061001B	AK561001B
Regular	20"	AK061005B	AK561005B
Long	26"	AK061007B	AK561007B



Original Bledsoe Brace

- Offers cast-like support with a weight under 4 pounds; the brace gives support in all planes while allowing muscles to maintain their natural shape
- Rigid medial and lateral stays and molded posterior supports provide a thorough support system
- Fully adjustable for atrophy or edema
- Can easily be removed and reapplied in one piece

Indications: For immediate application to knees after operative procedures or injury to knee ligaments, cartilage, or stable or internally fixed fractures



Leg Length	Brace Size Length	Multi-Centric Hinge	Product No.	
			Simple Hinge	QuikLok Hinge
17"-19"	18"	AK010018BB	AK020018BB	AK520018BB
19"-21"	20"	AK010020BB	AK020020BB	AK520020BB
21"-23"	22"	AK010022BB	AK020022BB	AK520022BB
23"-25"	24"	AK010024BB	AK020024BB	AK520024BB
25"-27"	26"	AK010026BB	AK020026BB	AK520026BB
27"-29"	28"	AK010028BB	AK020028BB	AK520028BB
29"-31"	30"	AK010030BB	AK020030BB	AK520030BB
31"-33"	32"	AK010032BB	AK020032BB	AK520032BB
33"-35"	34"	AK010034BB	AK020034BB	AK520034BB

* For leg length, measure 1" below crotch to ankle

Optional Shoe Insert

- Complete with absorbent sock and free ankle hinge

Size	Shoe Sizes		Product No.	
	Men's	Women's	Left	Right
Small	up to 4	up to 5	AM005101B	AM005201B
Medium	5-9	6-10	AM005105B	AM005205B
Large	10-14	11-15	AM005109B	AM005209B

Extender Plus

- Versatile and quick application with push button length adjustment from 22"-34"
- Accommodates circumference up to 34"
- Double-padded ankle plates for maximum suspension
- Quick release buckles

Indications: Locked or limited motion control of knee during rehabilitation, after operative procedures or injury to the knee ligaments, patellar ligament, articular cartilage or meniscus; Stable or internally fixed fractures of tibial plateau, condyles, proximal tibia or distal femur

Description	Product No.
Extender Plus	EK087000
Extender Plus Short	EK009000



Merit OR

- Vinyl hinge upright pockets with hook-to-pile closures for quick attachment to wraps
- Post Op Hinge provides fast and easy adjustment with flexion limits of -10° to 110° and extension limits of -10° to 110° in 10° increments; no tools to adjust
- Hinge can be locked in position from -10° to 40°

Indications: For immediate application to knees for locked or limited motion control of knee during rehabilitation after operative procedures or injury to knee ligaments, patellar ligament, articular cartilage or meniscus; Stable or internally fixed fractures of tibial plateau, condyles, proximal tibia or distal femur

Size	Leg Length	Product No.
Short	16"-22"	EK080001
Regular	22"-26"	EK080005
Long	26"-30"	EK080007
X-Long	30"+	EK080009



OR Brace

- Vinyl hinge arm pockets with hook-and-pile closures for quick attachment to foam cuffs
- Can be removed in one piece for easy attachment
- Simple hinge is quick setting and limits extension from 0° to 90°; flexion 0° to 135°; and locks 0° to 90° in 15° increments
- Quiklok will enable the knee hinge to be locked rigidly in full extension for ambulation or sleeping

Indications: For motion control of knee during rehabilitation after operative procedures or injury to knee ligaments, cartilage, or stable or internally fixed fractures

Size	Length	Product No.	
		Simple Hinge	Quiklok Hinge
Regular	20"	AK080005B	AK580005B
Long	24"	AK080007B	AK580007B
X-Long	28"	AK080009B	AK580009B



Revolution 3

- Economy post-operative brace
- Patient-friendly adjustment tabs
- Full range-of-motion hinge adjustable from -10° to 110°

Indications: For immediate application to knee for locked or limited motion control during rehabilitation after operative procedures or injury to knee ligaments, patellar ligament, articular cartilage or meniscus; Stable or internally fixed fractures of tibial plateau, condyles, proximal tibia or distal femur

Size	Leg Length	Product No.
Short	18"-22"	EK019003
Long	24"-28"	EK019005



Breg - Post-Op

- Full foam
- Spring loaded flexion/extension stops for easy range-of-motion from 0° to 120°
- Ability to lock out in full extension
- X-Large: 27"-35" thigh circumference
- Scored anodized bars

Indications: ACL, PCL, MCL, LCL injuries; Tibial plateau fractures; Osteochondral repairs; Meniscal repairs; Patella tendon repairs; Condylar fractures; HTOs (high tibial osteotomy); Acute sprains/strains of the knee

Length	Measurement (6" above mid-patella)	Product No.	
		Regular	X-Large
Short	Less than 32" inseam	00112	00113
Long	Greater than 32" inseam	00114	00115



Post-Op Rehab

- Anodized hinge bars
- Spring loaded flexion/extension stops for easy range-of-motion from 0° to 120°
- Full foam
- Latex free

Indications: ACL, PCL, MCL, LCL injuries; Tibial plateau fractures; Osteochondral repairs; Meniscal repairs; Patella tendon repairs; Condylar fractures; HTOs (high tibial osteotomy); Acute sprains/strains of the knee

Product No. **00123**



Post-Op Lite

- Reduced foam for warmer climates
- Spring-loaded flexion/extension stops for easy range-of-motion from 0°-120°
- Ability to lock out in full extension
- X-Large: 27"-35" thigh circumference
- Scored anodized bars

Indications: ACL, PCL, MCL, LCL injuries; Tibial plateau fractures; Osteochondral repairs; Meniscal repairs; Patella tendon repairs; Condylar fractures; HTOs (high tibial osteotomy); Acute sprains/strains of the knee

Length	Measurement (6" above mid-patella)	Product No.	
		Regular	X-Large
Short	Less than 32" inseam	00116	00117
Long	Greater than 32" inseam	00118	00119



ROM -
Post Op/Rehab

Post-Op with Shells

- Spring-loaded flexion/extension stops for easy range-of-motion from 0°-120°
- Full circumferential shells
- Ability to lock out in full extension
- Latex free

Indications: ACL, PCL, MCL, LCL injuries; Tibial plateau fractures; Osteochondral repairs; Meniscal repairs; Patella tendon repairs; Condylar fractures; HTOs (high tibial osteotomy); Acute sprains/strains of the knee

Length	Length	Product No.
Short	23 ¹ / ₂ "	00121
Long	26 ¹ / ₂ "	00122



T-Scope® Post Op Knee Brace

- Telescoping hinge bars extendable from 17" to 27" accommodate a larger range of patient leg lengths
- Quick lock offers ability to knock out brace from -10° to 30° in increments of 10°

Indications: ACL, PCL, MCL, LCL injuries; Tibial plateau fractures; Osteochondral repairs; Meniscal repairs; Patella tendon repairs; Condylar fractures; HTOs (high tibial osteotomy)

Description	Thigh Circumference	Product No.
T-Scope Universal	up to 30 ¹ / ₂ "	07714
T-Scope X-Large	27"-35"	07715
T-Scope Full Foam Universal	up to 30 ¹ / ₂ "	07716
T-Scope Full Foam X-Large	27"-35"	07717



Quick Fit EPO

- Easy-to-move set pins for range-of-motion
- Ability to lock-out in full extension
- Full foam wrap for comfort and compression
- Easy-to-use hook and loop straps

Indications: ACL, PCL, MCL, LCL injuries; Tibial plateau fractures; Osteochondral repairs; Meniscal repairs; Patella tendon repairs; Condylar fractures; Sprains/strains of the knee; Total knees; HTOs (high tibial osteotomy)

Size	Thigh Circumference (6" above mid-patella)	Product No.	
		Short	Long
Small	14"-18"	09342	09352
Medium	18"-22"	09343	09353
Large	22"-27"	09344	09354
X-Large	27"-31"	09345	09355



DonJoy (DJO) - Cool IROM®

- Formable hinge bars allow for varus/valgus positioning
- 10° incremental settings
- 4 straps and abbreviated foam at key control points
- Integrated strap/hinge design with hinge bar breakdown capability

Indications: Immobilization and precise range-of-motion control for ACL, PCL, MCL and LCL surgeries in warmer climates

Size	Product No.	
	Gray	Black
Short (up to 32")	11-0175-2-13066	11-0175-2-06000
Long (32"+)	11-0175-4-13066	11-0175-4-06000



Telescoping Cool ELS®

- Slide lock mechanism that can lock at full extension or can have a full range-of-motion
- Brace can extend to ankle to help prevent migration
- Abbreviated foam

Indications: Post-op use when full extension weight-bearing protocols are used

Product No. **11-0967-9-13066**



Telescoping ELS Knee Support with Full Foam

- Slide lock mechanism that can lock at full extension or can have a full range-of-motion
- Brace can extend to ankle to help prevent migration

Indications: Post-op use when full extension weight-bearing protocols are used

Product No. **11-0965-9-13066**



ROM - POST OPERATIVE/REHAB

Fastfit and Rehab TROM™

- Sized, ready to fit brace with a simple, one piece application
- Easy-to-use TROM hinge provides controlled range-of-motion and quick locking capabilities
- FastFit TROM is 24" in length
- Rehab TROM cuffs designed for extended rehabilitation and is 17" in length

Indications: Ligament surgeries, meniscal repairs and patella realignments

Size	Measurements (6" above mid-patella)	Product No.	
		FastFit TROM	Rehab TROM
Small	14"-18"	11-0842-2-06000	11-0295-2-06000
Medium	18"-22"	11-0842-3-06000	11-0295-3-06000
Large	22"-27"	11-0842-4-06000	11-0295-4-06000
X-Large	27"-32"	11-0842-5-06000	11-0295-5-06000



FastFit TROM

Rehab TROM

Telescoping TROM™ Advance

- Tele-Fit allows independent strap movement proximal and distal of knee, allowing an overall better fit
- Telescoping adjustability from 17³/₄" to 23¹/₄"
- Range-of-motion control from 10° to 90°
- Quick release buckles

Indications: ACL, PCL, MCL & LCL surgeries; Meniscal repairs; Patella realignment; Regenerative chondroplasty; Stable femoral fractures; Total knee replacements; High tibial osteotomy

Description	Thigh Circumference (6" above mid-patella)	Product No.
Cool	Less than 26"	11-9114-9
Cool, X-Large	Greater than 26"	11-9116-9
Full Foam	Less than 26"	11-9115-9



Telescoping TROM™

- Allows flexion and extension adjustments at the touch of a button and a quick lock mechanism from 0° to 90° in 10° increments
- Telescoping from 17³/₄" to 23¹/₄"
- Buckles allow for fast and easy strapping

Indications: Major ligament surgery; Meniscal repair; Patella realignment; Regenerative chondroplasty; Stable femoral fractures

Description	Thigh Circumference (6" above mid-patella)	Product No.
Regular Full Foam	Less than 26"	11-0961-9-06000
X-Large Full Foam	26"-32"	11-0964-9-06000
Cool Universal	Less than 26"	11-0963-9-06000

Competitor™

- Full and limited range-of-motion settings at 0°, 15°, 30° and 60° and 90°
- Dial in adjustment of flexion and extension settings
- Swing-cap dual axis hinge
- Integrated strap/hinge design for circumferential compression
- Breakdown capability on long version; fits up to 30" circumference

Description	Length	Product No.
Short (up to 28" inseam)	18"	11-0523-2-06060
Long (greater than 28" inseam)	23"	11-0523-4-06060

Competitor



Telescoping IROM® Splint

- Precise ROM in 10° increments
- Telescoping bars allow the brace to fit a variety of patients
- Extends to ankle to help prevent migration

Indications: ACL, PCL, MCL and LCL surgeries

Description	Thigh Circumference	Product No.
Universal	—	11-0968-9-06000
X-Large Foam	26-32"	11-0971-9-06000

IROM Splint



Telescoping Cool IROM

- For warmer climates, this brace provides immobilization and precise range-of-motion control for ACL, PCL, MCL and LCL surgeries

Product No. [11-0970-9-06000](#)

Cool IROM



Ossur Americas - Innovator DLX®

- One touch™ hinge for ease of adjustment
- Reusable ROM locks for easy, accurate and quick fitting that ensures patient compliance
- Quick release buckles for easy donning and doffing
- Hyperextension stop at 10°
- Slide-to-size struts allow for easy length adjustments from 17"-28"

Indications: Post-op/post-injury ROM control or immobilization

Size	Thigh Circumference (6" above mid-patella)	Product No.	
		Cool	Full
Regular	Up to 27"	423100I	423200I
X-Large	27"-35"	423101I	423201I

Cool Foam

Full Foam



Innovator®

- Locks quickly at 0°, 10°, 20°, 30°, 45°, 60°, 75°, and 90°
- Simply press, rotate and release flexion/extension settings
- Low-profile hinge with no exposed metal
- Breakaway struts for easy, tool-free length adjustment

Indications: Post-op/post-injury ROM control or immobilization

Size	Thigh Circumference (6" above mid-patella)	Length	Product No.	
			Cool	Full
Regular	Up to 27"	16"	—	403300I
X-Large	27"-35"	16"	—	403301I
Regular	Up to 27"	24"	403100I	403200I
X-Large	27"-35"	24"	403102I	403202I
Regular	Up to 27"	28"	403101I	403201I
X-Large	27"-35"	28"	403103I	403203I



Rehab Universal

- Push-button system for easy length adjustment
- ROM control in 15° increments: Flexion 0°-105°, Extension 15°-90°
- Length adjustable from 17"-28"
- Universal design will fit either leg

Indications: Post-op/post-injury ROM control or immobilization

Size	Drop Lock	Product No.	
		Air Light	Full Foam
Regular	No	B-222000000	B-223000000
Regular	Yes	B-236000000	B-235000000



ProCare (DJO) - KneeRANGER® Sized

- 24" length brace offers post-operative immobilization
- Protected range-of-motion with 7.5° incremental settings between 0°-120° of flexion and 0°-70° of extension
- One-piece easy on/off design

Indications: Immobilization and protected range-of-motion associated with ACL, PCL, LCL and MCL surgeries/repairs

Size	Thigh Circumference (4" above mid-patella)	Product No.
Small	14"-16"	79-94003
Medium	16"-18"	79-94005
Large	18"-20"	79-94007
X-Large	20"-22"	79-94008
2X-Large	22"-24"	79-94009



KneeRANGER® Lite

- Narrow band circumferential cuffs replace full cut leg wraps
- Provides immobilization and protected range-of-motion between 0°-120° of flexion and 0°-70° of extension in 7.5° increments
- Universal length: 24"
- Maximum circumferences: thigh, 31"; calf, 22"

Indications: Immobilization and protected range-of-motion associated with ACL, PCL, LCL and MCL surgeries/repairs

Product No. **79-94040**



KneeRANGER® II - Universal

- Dual axis, polycentric hinge allows for dial-in adjustment
- Flexion and extension stops at 0°, 15°, 30°, 60° and 90°
- Trimmable foam leg wraps for customized fit
- Medial/lateral plastic shells riveted to hinge bars provide additional support following post-operative procedures

Indications: Immobilization and protected range-of-motion associated with ACL, PCL, LCL and MCL surgeries/repairs

Size	Length	Max. Thigh Circumference	Product No.
Universal/Short	17"	36"	79-94230
Universal/Long	26"	36"	79-94240



KneeRANGER® - Universal

- Post-operative immobilization and protected range-of-motion
- 7.5° incremental settings from 0° to 120° of flexion and 0° to 70° of extension
- Hinge bar breakdown allows for up to 5" reduction in length to convert from long leg to short for long-term rehabilitation
- Trimmable medial/lateral plastic shells

Indications: Patients recovering from injury or surgery providing full and limited range-of-motion setting

Size	Length	Product No.	
		Each	6 Pack
X-Short	16" overall length	79-94130	—
Short*	32" inseam or less	79-94000*	79-94100*
Long*	32" plus	79-94001*	79-94101*

*Breakdown bars. Maximum Circumference: Thigh 31", Calf 22"



HD Post Operative Pin® Knee Brace

- Flexion and extension control settings at 5° increments from 0° to 120°
- Allows for free ROM and can be locked in full extension
- Up to 27" thigh and up to 19" calf circumference
- Heavy duty POP with lightweight DL hinge

Indications: ACL; PCL; MCL and LCL surgeries

Size	Length	Breakaway Length	Product No.
Long Full	24 ¹ / ₂ "	18"	73HDP-LX-F
Long Cool	24 ¹ / ₂ "	18"	73HDP-LX-C
Short Full	18"	—	73HDP-SX-F
Short Cool	18"	—	73HDP-SX-C



RCAI - Ratchet POP Knee Brace

- Flexion/extension ROM control in 5° increments
- Angle of limb can be extended in 10° increments from 0°-60° without using tools
- Knee cap provides 3-point pressure leverage

Indications: Post-surgical ROM limitations; Ligament repairs; Instability and soft tissue contractures

Size	Thigh Circumference (4" above knee center)	Product No.
Pediatric	9"-12"	67RKPOP-P
Small	12"-15"	67RKPOP-S
Medium	15"-20"	67RKPOP-M
Large	20"-24"	67RKPOP-L



Universal Leg Brace

- Infinite degree of adjustment settings offer precise flexion and extension control of the knee
- 7° tibial offset follows the tibia in correct anatomical alignment
- Maximum thigh circumference: 27"; Calf circumference: 19"

Indications: Strains; Sprains; Stable tibia or femoral fractures; Ligament repairs; Total knee replacement and reconstruction

Size	Leg Length (1" below groin to medial malleolus)	Regular
Short	18" - 23"	87ULB-SX
Average	24" - 29"	87ULB-AV
Long	30" - 35"	87ULB-LX



Original Post Operative Pin® Knee Brace

- Lightweight, drop pin hinge controls flexion and extension in 7.5° increments
- Allows for free range-of-motion and can be locked in full extension
- Maximum thigh circumference: 27"; Calf circumference: 19"

Indications: Post operative immobilization and protected range-of-motion associated with ACL, PCL, MCL, and LCL surgeries

Size	Length	Breakaway Length	Product No.
Long Full	22 ¹ / ₂ "	16"	77POP-LX-F
Long Cool	22 ¹ / ₂ "	16"	77POP-LX-C
Short Full	16"	—	77POP-SX-F
Short Cool	16"	—	77POP-SX-C

Size	Length	Maximum Thigh	Maximum Calf Circum.	Product No.
Pediatric	13"	18"	13 ¹ / ₂ "	77POP-P



CoreLINE - Tri-Panel Knee Immobilizer

- Adjustable three piece design accommodates multiple circumferences
- Lightweight aluminum posterior, medial and lateral stays are malleable to conform to the contours of the knee
- Side panels can be repositioned to optimally locate the medial and lateral stays as needed for improved fit
- Loop lock closure with plastic D rings

Indications: Post-operative knee support and immobilization; Grade 1 & grade 2 collateral ligament sprains; Patellar subluxations

Length	Max Thigh Circum.	Max Calf Circum.	Product No.
12"	20 ¹ / ₂ "	16 ¹ / ₂ "	CL0131
14"	20 ¹ / ₂ "	16 ¹ / ₂ "	CL0132
16"	22"	18"	CL0133
18"	23"	19"	CL0134
20"	24"	20"	CL0135
22"	24"	20"	CL0136
24"	24"	21"	CL0137



AliMed - Universal Knee Immobilizer

- Adjustable stays and straps facilitate correct medial and lateral positioning
- 3/8" thick laminated foam with alternating hook-and-loop closures
- Dual aluminum stays on lateral and medial stay/strap assemblies
- Fits left or right

Length	Product No.
12"	60141
16"	60142
18"	60143
20"	60144
24"	60145



Bledsoe Brace - Lever Lock

- Easily locks in full extension to immobilize the leg and unlocks to allow motion for rehab, passive motion, or sitting
- Ideal for arthroscopic procedures
- Available in full cuff design only

Indications: For application following arthroscopic reconstruction of the ACL or other procedures or injuries which require locking of the knee in full extension with intermediate periods of full range-of-motion

Size	Leg Length	Product No.
Regular	21"-25"	AK070005B
Long	25"-29"	AK070007B
X-Long	29"+	AK070009B



Breg - Tri-Panel Knee Immobilizer

- Adjustable side panels with posterior rigid stays for additional support
- One size fits all (up to 25" thigh) with easy to use hook and loop

Indications: ACL, PCL, MCL and LCL injuries; Sprains/strains of the knee; Soft tissue injuries

Universal Length	Product No.
16"	96500
20"	96510
24"	96520



Corflex - Universal Knee Immobilizer

- Breathable foam laminate with aluminum stays
- Contoured posterior stays and movable medial/lateral stay assemblies
- 4" elastic patella straps with hook and loop closure
- Measure 4" above mid-patella

Indications: Post-op knee immobilization; Knee injuries where full extension of the knee is desired

Length	Thigh Circumference		Product No.	
	Universal	X-Wide	Universal	X-Wide
6"	11"	—	51-1506	—
9"	14"	—	51-1509	—
12"	20"	25"	51-1512	92-5112
14"	20"	25"	51-1514	92-5114
16"	24"	29"	51-1516	92-5116
18"	24"	29"	51-1518	92-5118
20"	25"	30"	51-1520	92-5120
22"	25"	30"	51-1522	92-5122
24"	27"	32"	51-1524	92-5124



Compression Knee Immobilizer

- Sewn-in contoured posterior stays, movable medial/lateral stays and elastic compression straps for optimal immobilization
- Breathable foam laminate, aluminum stays, premium quality elastic straps and contact closure; No pockets

Indications: Post-op immobilization; Mild to severe knee injuries where full extension of the leg is desired

Length	Thigh Circumference*	Product No.
6"	11"	75-0651
9"	14"	75-0951
10 ¹ / ₂ "	14"	75-1051
12"	24"	75-1251
14"	28"	75-1451
16"	28"	75-1651
18"	28"	75-1851
20"	29"	75-2051
22"	29"	75-2251
24"	31"	75-2451

*Measurement taken from leg circumference 4" above mid-patella



DeRoyal - Cutaway Knee Immobilizer with Patella Strap

- T-bar stays provide optimal strength
- Contoured posterior stays and adjustable medial/lateral stays
- Elastic patella strap helps reduce knee edema

Indications: Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

Universal Length	Thigh Circumference	Product No.
14"	12"-24"	7066-14
16"	12"-24"	7066-16
18"	12"-24"	7066-18
20"	12"-24"	7066-20
22"	12"-24"	7066-22
24"	12"-24"	7066-24
26"	12"-24"	7066-26
28"	12"-24"	7066-28



Knee Immobilizer with Elastic Straps

- T-bar stays provide optimal strength
- Adjustable medial/lateral stays allow for universal sizing
- Elastic straps provide compression, helping to reduce edema

Indications: Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

Universal Length	Thigh Circumference	Product No.
13"	12"-24"	7071-11
16"	12"-24"	7071-12
19"	12"-24"	7071-13
23"	12"-24"	7071-15
28"	12"-24"	7071-17



Deluxe Tri-Panel Knee Immobilizer

- Foam body with durable blue canvas panels
- T-bar stays provide optimal strength
- Contoured posterior stays for correct anatomical fit
- Adjustable medial and lateral panels allow for universal sizing

Indications: Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

Universal Length	Product No.	
	Canvas	Foam
12"	1120917	BF1120917
14"	1140917	—
16"	1160917	BF1160917
19"	1190917	BF1160917
22"	1220917	BF1220917
24"	1240917	BF1240917



KNEE IMMOBILIZERS

DonJoy (DJO) - Basic Knee Immobilizer

- Wide elastic strapping to fit all legs with minimal migration and twisting
- Contoured posterior aluminum stay
- Trimmable open-cell foam

Length	Product No.
12"	11-0389-1-13150
17"	11-0389-2-13150
20"	11-0389-3-13150
24"	11-0389-4-13150
27"	11-0389-5-13150



Deluxe Knee Immobilizer

- Provides lightweight immobilization following injury or surgery
- 6-strap fastening system with extra wide hook and loop
- Anatomically contoured aluminum stays - medial, lateral and posterior

Size	Circumference (4" above mid-patella)	Length				
		12"	16"	20"	23"	26"
Small	16"	11-0226-2	11-0227-2	11-0228-2	11-0229-2	11-0230-2
Medium	19"	11-0226-3	11-0227-3	11-0228-3	11-0229-3	11-0230-3
Large	22"	11-0226-4	11-0227-4	11-0228-4	11-0229-4	11-0230-4
X-Large	26"	11-0226-5	11-0227-5	11-0228-4	11-0229-5	11-0230-5



Hely & Weber - Knee Immobilizer

- Reinforced, moldable posterior stays provide additional rigidity
- Medial and lateral movable stay pods offer extra support
- $\frac{3}{8}$ " pile foam covered with alternating elastic straps make the brace easy to apply and adjust
- Fits mid-patella circumferences 12"-23"

Length	Product No.
12"	201
14"	202
16"	203
18"	204
20"	205
22"	206
24"	207
26"	208



Three Panel Knee Immobilizer

- 3/8" pile foam laminate with black flannel interior for patient comfort
- Two detachable panels for circumferential adjustments for custom fitting
- Posterior and medial/lateral stays are malleable
- Optional patella strap provides extra knee control

Length	Product No.
12"	211
14"	212
16"	213
18"	214
20"	215
22"	216
24"	217
26"	218



Med Spec - Universal Foam Knee Immobilizer

- Lightweight foam knee immobilizer with a wide range of fittings
- Five elastic circumferential straps hold immobilizer firmly in place and reduce migration
- Malleable medial and lateral stays provide customization for additional support
- Sturdy posterior stays keep leg in full extension

Length	Thigh Circumference	Product No.
12"	up to 20"	223690
16"	up to 24 1/2"	223660
18"	up to 26 1/2"	223665
20"	up to 29"	223670
20"	up to 33"	223677
22"	up to 31 1/2"	223675
24"	up to 33 1/4"	223680
24"	up to 37"	223687



Ossur Americas - Exoform® Knee Immobilizer

- Slide-to-size adjustment allows variety of lengths and widths
- Dual-cuff design
- Popliteal pad provides superior fit and helps limit migration
- Perforated padded liner for breathability and comfort
- Cool foam version can be worn directly against the skin for greater comfort on warm days

Description	Leg Circum.	Length	Product No.
Cool	up to 28"	18"-24"	222100
Full	up to 28"	18"-24"	222000
XXXL-Full	up to 38"	18"-24"	222050



Cool Foam

Full Foam

Knee Immobilizers

ProCare (DJO) - Clinic 3-Panel Knee Splint

- Universal 3-piece design for ease of adjustment
- Foam posterior panel may be trimmed for a customized fit
- Foam/mesh side panels with I-beam stays and loop-lock straps

Indications: Post-operative immobilization of the knee joint; Grade 1 and 2 collateral ligament sprains; Patellar subluxations

Length	Max. Thigh Circumference (4" above patella)	Product No.
12"	32"	79-80410
16"	32"	79-80420
20"	32"	79-80170
24"	32"	79-80180
26"	32"	79-80290



3-Panel Knee Splint

- Universal size, 3-piece nylon/fiber laminate construction
- Adjustable side panels with I-beam stays and patella strap
- Three posterior aluminum stays for proper immobilization

Indications: Post-operative immobilization of the knee joint; Grade 1 and 2 collateral ligament sprains; Patellar subluxations

Max. Length	Size	Max. Thigh Circumference (4" above patella)	Product No.
12"	Standard	29"	79-80000
14"	Standard	29"	79-80010
16"	Standard	29"	79-80110
20"	Standard	29"	79-80020
20"	X-Large	36"	79-80021
22"	Standard	29"	79-80040
22"	X-Large	36"	79-80041
24"	Standard	29"	79-80030
24"	X-Large	36"	79-80031



Quick-Fit® Basic Knee Splint (Universal)

- Foam laminate construction with 3" elastic contact closure straps provides compressive support and immobilization
- Dual contoured I-beam posterior stays and moveable medial/lateral double stay casings for correct fit

Length	Max. Thigh Circumference (4" above patella)	Product No.
10"	29"	79-96011
12"	29"	79-96012
14"	29"	79-96014
16"	29"	79-96016
18"	36"	79-96018
20"	36"	79-96019
22"	36"	79-96021
24"	36"	79-96024
26"	36"	79-96026



Universal Foam Knee Splint

- Constructed of open cell foam and detachable medial/lateral aluminum I-beam stay panels for easy adjustment and custom fit
- Three double-bend posterior aluminum stays

Indications: Post-operative immobilization of the knee joint; Grade 1 and 2 collateral ligament sprain; Patella subluxations

Length	Max. Thigh Circumference	Product No.
8"	29"	79-80600
10"	29"	79-96870
12"	29"	79-96871
14"	29"	79-96914
16"	29"	79-96880
18"	29"	79-96881
20"	29"	79-96882
22"	29"	79-96890
24"	29"	79-96891
26"	29"	79-96900



Professional Products - Deluxe Foam Lined Knee Immobilizer

- 3" wide elastic compression straps with hook and loop closure (12" and 14" have 2" wide straps)
- Three aircraft aluminum contoured posterior stays (12" and 14" have two stays) and two aircraft aluminum medial/lateral stays

Length	Product No.
12"	01259-C-B-12
14"	01259-C-B-14
16"	01259-C-B-16
18"	01259-C-B-18
20"	01259-C-B-20
22"	01259-C-B-22
24"	01259-C-B-24
26"	01259-C-B-26
28"	01259-C-B-28
32"	01259-C-B-32
34"	01259-C-B-34



RCAI - Bariatric Knee Immobilizer

- Aluminum stays located at the posterior, medial and lateral sections of the brace provide stability
- Material wicks moisture, keeping patient cool, dry and comfortable

Length	Max. Thigh Circumference (4" above patella)	Product No.
15"	32"	420-BKI-15
18"	34"	420-BKI-18
21"	36"	420-BKI-21
24"	39"	420-BKI-24



Single Panel Knee Immobilizer

- Utilizes elastic closer straps to better grip the leg
- Closed patella area
- Breathable foam material “wicks” away moisture from the skin

Length	Max. Thigh Circumference (4" above patella)	Product No.
6"	10"	42SPKI-06
9"	13"	42SPKI-09
12"	21"	42SPKI-12
14"	26"	42SPKI-14
16"	26"	42SPKI-16
18"	27"	42SPKI-18
20"	28"	42SPKI-20
22"	28"	42SPKI-22
24"	29"	42SPKI-24



Tri-Panel Knee Immobilizer

- Movable tri-panel design accommodates multiple circumferences
- The posterior aspect of the splint is uniquely tailored for a more custom anatomical fit
- The splint features lightweight aluminum posterior, medial and lateral stays

Length	Min. Thigh Circumference	Max. Thigh Circumference	Product No.
6"	6"	10 ¹ / ₄ "	42KI-06*
9"	11"	16"	42KI-09
12"	14"	20"	42KI-12
15"	16"	22"	42KI-15
18"	18"	28"	42KI-18
21"	19"	30"	42KI-21
24"	19"	32"	42KI-24

* 6" Length is single panel; Measure thigh circumference 4" above the knee



Trulife - Lenox Hill Three Panel Knee Immobilizer

- Removable steel stays offer increased support
- Adjusts to leg circumference
- Open patella design
- Hook and loop closure for easy donning and doffing

Indications: Immediate immobilization post-operative or post-injury

Size	Length	Product No.
X-Small	16"	KN300-01
Small	18"	KN300-02
Medium	20"	KN300-03
Large	22"	KN300-04
X-Large	24"	KN300-05



Bauerfeind - GenuTrain P3 Knee Support

- Silicone insert with plastic guide, surrounds kneecap and counteracts patellar shifting through corrective tensioning
- Intermittent compression promotes rapid reabsorption of edema, providing pain relief and improving knee function
- Replace ? with 1 (Nature/Beige), 7 (Black), or 8 (Titanium)

Size	Circumference		Product No.	
	6" below knee	6 ³ / ₄ " above knee	Left	Right
1	11"-12 ¹ / ₄ "	15 ³ / ₈ "-16 ¹ / ₂ "	110414020?0701	110414020?0601
2	12 ¹ / ₄ "-13 ³ / ₈ "	16 ¹ / ₂ "-17 ³ / ₄ "	110414020?0702	110414020?0602
3	13 ³ / ₈ "-14 ¹ / ₂ "	17 ³ / ₄ "-18 ¹ / ₈ "	110414020?0703	110414020?0603
4	14 ¹ / ₂ "-15 ³ / ₄ "	18 ¹ / ₈ "-20 ¹ / ₈ "	110414020?0704	110414020?0604
5	15 ³ / ₄ "-17"	20 ¹ / ₈ "-21 ¹ / ₄ "	110414020?0705	110414020?0605
6	17"-18 ¹ / ₈ "	21 ¹ / ₄ "-22 ¹ / ₂ "	110414020?0706	110414020?0606



Bledsoe Brace - 20.50 Patellofemoral

- Dynamic control in flexion to control patellar subluxation, dislocation and improve patellofemoral tracking
- Lightweight magnesium and short 12" length
- Provides additional MCL and LCL support

Indications: Patellofemoral pain syndrome; Patella Alta and Baja; Patellofemoral malalignment, subluxation and dislocation; Patellar tendonitis; Post-operative effusion; Jumper's Knee; Tight iliotibial band or quadriceps laxity

Size	Thigh Circum. (6" above mid-patella)	Calf Circum. (6" below mid-patella)	Product No.	
			Left	Right
X-Small	13 ¹ / ₂ "-16"	11"-12 ¹ / ₂ "	PF010101	PF010201
Small	16"-18 ³ / ₄ "	12 ¹ / ₂ "-14"	PF010103	PF010203
Medium	18 ³ / ₄ "-21 ¹ / ₂ "	14"-15 ¹ / ₂ "	PF010105	PF010205
Large	21 ¹ / ₂ "-24 ¹ / ₄ "	15 ¹ / ₂ "-17"	PF010107	PF010207
X-Large	24 ¹ / ₄ "-27"	17"-18 ¹ / ₂ "	PF010109	PF010209
2X-Large	27"-29 ¹ / ₂ "	18 ¹ / ₂ "-20"	PF010111	PF010211



Patellofemoral

Breg - Lateral Stabilizer with Hinge

- Available in airmesh and neoprene
- Designed specifically for a comfortable ultra low-profile fit
- J-shaped buttress that is ideal for lateral stabilization of the patella

Indications: Lateral patella subluxation; Chronic patellar maltracking; Mild chondromalacia

Size	Thigh Circumference (6" above mid-patella)	Airmesh Product No.		Neoprene Product No.	
		Left	Right	Left	Right
X-Small	12"-15"	20121	20141	20151	20171
Small	15"-18"	20122	20142	20152	20172
Medium	18"-21"	20123	20143	20153	20173
Large	21"-24"	20124	20144	20154	20174
X-Large	24"-27"	20125	20145	20155	20175
2X-Large	27"-30"	20126	20146	20156	20176



PTO®

- Available with aluminum high performance hinge bars
- Dynamic function for normal gait
- Cawley Tension Hinge

Indications: Mild chondromalacia; Patella subluxation; Patella dislocation; Patella tendonitis; Chronic patellar maltracking; Patellofemoral pain syndrome; Postlateral release; Post-patellar realignment

Size	Thigh Circum.	Closed Back				Open Back			
		Airmesh		Neoprene		Airmesh		Neoprene	
		Left	Right	Left	Right	Left	Right	Left	Right
X-Small	12"-15"	14191	14181	14171	14161	14231	14221	14211	14201
Small	15"-18"	14192	14182	14172	14162	14232	14222	14212	14202
Medium	18"-21"	14193	14183	14173	14163	14233	14223	14213	14203
Large	21"-24"	14194	14184	14174	14164	14234	14224	14214	14204
X-Large	24"-27"	14195	14185	14175	14165	14235	14225	14215	14205
2X-Large	27"-30"	14196	14186	14176	14166	14236	14226	14216	14206

Circumference taken 6" above mid patella



PTO Airmesh



PTO Neoprene

Lateral Stabilizer

- J-shaped buttress that is ideal for lateral stabilization of the patella
- Lateral straps for additional support
- Medial and lateral stays help prevent migration

Indications: Lateral patella subluxation; Chronic patellar maltracking; Mild chondromalacia

Size	Thigh Circum. (6" above mid-patella)	Product No.			
		Neoprene		Airmesh	
		Left	Right	Left	Right
X-Small	12"-15"	07081	07071	11041	11051
Small	15"-18"	07082	07072	11042	11052
Medium	18"-21"	07083	07073	11043	11053
Large	21"-24"	07084	07074	11044	11054
X-Large	24"-27"	07085	07075	11045	11055
2X-Large	27"-30"	07086	07076	11046	11056



Lateral Stabilizer

Buttress Support

- Straps above and below knee with adjustable buttress
- Medial and lateral stays help prevent migration
- Aluminum stays to prevent migration

Indications: Patella subluxation; Chronic patellar maltracking; Mild chondromalacia

Size	Thigh Circumference (6" above mid-patella)	Product No.
X-Small	12"-15"	07211
Small	15"-18"	07212
Medium	18"-21"	07213
Large	21"-24"	07214
X-Large	24"-27"	07215
2X-Large	27"-30"	07216



Corflex - U-Shaped Patella Stabilizer

- Sewn-in U-shaped solid tubular buttress to provide optimal contact and patella stability with spiral stays to prevent rollover; some styles available with superior buttress placement
- Manufactured from premium grade 1/8" latex-free neoprene

Indications: Patellar subluxation; Patello-femoral tracking syndrome; Treatment of patella alta/baja; As part of a rehabilitation program for patella tendon injuries



Size	(6" above mid-patella)	(mid-patella)	(6" below mid-patella)	Product No.
X-Small	16"-18"	12"-13"	11"-12"	88-5092
Small	17"-19"	13"-14"	12"-13"	88-5093
Medium	18"-20"	14"-15"	13"-14"	88-5094
Large	19"-21"	15"-16"	14"-15"	88-5095
X-Large	20"-22"	16"-18"	15"-16"	88-5096
2X-Large	21"-23"	18"-20"	16"-17"	88-5097
3X-Large	24"-26"	20"-22"	17"-18"	88-5098
4X-Large	27"-29"	22"-24"	18"-19"	88-5099

Knee-O-Trakker with Hinge

- Textured C-shaped buttress for maximizing patella control and versatile strapping system that provides dynamic treatment with adjustable tension control
- Manufactured from premium grade latex-free 3/16" neoprene and compression molded textured neoprene buttress

Indications: Patellar subluxation; Patello-femoral syndrome; As part of a rehabilitation program for patella tendon injuries where medial/lateral support is also required



Size	(6" above mid-patella)	(mid-patella)	(6" below mid-patella)	Product No.
X-Small	16"-18"	12"-13"	11"-12"	88-7272
Small	17"-19"	13"-14"	12"-13"	88-7273
Medium	18"-20"	14"-15"	13"-14"	88-7274
Large	19"-21"	15"-16"	14"-15"	88-7275
X-Large	20"-22"	16"-18"	15"-16"	88-7276
2X-Large	21"-23"	18"-20"	16"-17"	88-7277
3X-Large	24"-26"	20"-22"	17"-18"	88-7278
4X-Large	27"-29"	22"-24"	18"-19"	88-7279

Patellofemoral

DeRoyal - Full Length Patella Stabilizer

- Breathable Tri-tex material
- Full-length design for maximum support and compression
- Horseshoe-shaped pad can be placed on medial or lateral side of patella
- Medial and lateral hinges for additional ligament support

Indications: Medial or lateral patellar subluxation; Chondromalacia or patellofemoral tracking dysfunction

Size	Thigh Circumference (6" above mid-patella)	Product No.	
		Left	Right
X-Small	13 ¹ / ₂ "-15 ¹ / ₂ "	14840104	14840204
Small	15 ¹ / ₂ "-18"	14840105	14840205
Medium	18"-20 ¹ / ₂ "	14840106	14840206
Large	20 ¹ / ₂ "-23"	14840107	14840207
X-Large	23"-25 ¹ / ₂ "	14840108	14840208
2X-Large	25 ¹ / ₂ "-28"	14840109	14840209



Concise Patella Stabilizer

- Breathable Tri-tex material
- Low-profile design allows for increased range-of-motion
- Unique strapping system helps pull patella into desired position
- Universal left or right

Indications: Medial or lateral patellar subluxation; Chondromalacia; Patellofemoral tracking dysfunction

Size	Thigh Circumference (6" above mid-patella)	Product No.
X-Small	13 ¹ / ₂ "-15 ¹ / ₂ "	14750004
Small	15 ¹ / ₂ "-18"	14750005
Medium	18"-20 ¹ / ₂ "	14750006
Large	20 ¹ / ₂ "-23"	14750007
X-Large	23"-25 ¹ / ₂ "	14750008
2X-Large	25 ¹ / ₂ "-28"	14750009



Tru-Pull Wraparound

- Independent straps, which provide a consistent pull on the patella
- Choice of two adjustable buttress thicknesses (1/4" and 3/8")
- Universal left or right
- Wrap-style design for easy application and fit

Indications: Patellar subluxation, dislocation, malalignment and instabilities

Size	Thigh Circumference (6" above mid-patella)	Product No.
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-1111-2-06000
Medium	18 ¹ / ₂ "-21"	11-1111-3-06000
Large	21"-23 ¹ / ₂ "	11-1111-4-06000
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-1111-5-06000
Buttress Kit (universal)	—	11-1112-0-06000



DonJoy (DJO) - Tru-Pull® Lite™

- Bifurcated strap dynamically pulls the patella into proper alignment, reducing pain associated with patellofemoral dysfunction
- Dual Durometer Buttress System interfaces directly with the patella, applying a consistent corrective force to reduce pain
- Anti-microbial wicking fabric is hypoallergenic, breathable and lightweight
- Set of removable plastic hinges included with each brace to provide medial/lateral stability

Indications: Acute mild to moderate patellar subluxation or dislocation and maltracking



Size	Thigh Circumference (6" above mid-patella)	Product No.	
		Left	Right
X-Small	13"-15 ¹ / ₂ "	11-0261-1	11-0260-1
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-0261-2	11-0260-2
Medium	18 ¹ / ₂ "-21"	11-0261-3	11-0260-3
Large	21"-23 ¹ / ₂ "	11-0261-4	11-0260-4
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-0261-5	11-0260-5
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-0261-6	11-0260-6
3X-Large	29 ¹ / ₂ "-32"	11-0261-7	11-0260-7

Tru-Pull Advanced System

- Bifurcated, elastomeric strap dynamically pulls the patella as the knee extends
- Two buttresses (1/4" and 3/8") included, trimmable to provide superior, inferior, and lateral patellar stabilization
- Independent anchors prevent brace rotation

Indications: Chronic patellar subluxation or dislocation



Size	Thigh Circumference (6" above mid-patella)	Product No.			
		Standard		with Popliteal Cutout	
		Left	Right	Left	Right
X-Small	13"-15 ¹ / ₂ "	11-1290-1	11-1289-1	11-1293-1	11-1292-1
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-1290-2	11-1289-2	11-1293-2	11-1292-2
Medium	18 ¹ / ₂ "-21"	11-1290-3	11-1289-3	11-1293-3	11-1292-3
Large	21"-23 ¹ / ₂ "	11-1290-4	11-1289-4	11-1293-4	11-1292-4
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-1290-5	11-1289-5	11-1293-5	11-1292-5
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-1290-6	11-1289-6	11-1293-6	11-1292-6
3X-Large	29 ¹ / ₂ "-32"	11-1290-7	11-1289-7	11-1293-7	11-1292-7

Patellofemoral

Tru-Pull Advanced System with Hinges

- Bifurcated, elastomeric strap dynamically pulls the patella as the knee extends
- Two buttresses ($1/4$ " and $3/8$ ") included, trimmable to provide superior, inferior, and lateral patellar stabilization
- Independent anchors prevent brace rotation
- Offers medial/lateral support and flexion/extension control

Indications: Chronic patellar subluxation or dislocation



Size	Thigh Circumference (6" above mid-patella)	Product No.			
		Standard Left	Standard Right	with Popliteal Cutout Left	with Popliteal Cutout Right
X-Small	13"-15 $1/2$ "	11-2290-1	11-2289-1	11-2293-1	11-2292-1
Small	15 $1/2$ "-18 $1/2$ "	11-2290-2	11-2289-2	11-2293-2	11-2292-2
Medium	18 $1/2$ "-21"	11-2290-3	11-2289-3	11-2293-3	11-2292-3
Large	21"-23 $1/2$ "	11-2290-4	11-2289-4	11-2293-4	11-2292-4
X-Large	23 $1/2$ "-26 $1/2$ "	11-2290-5	11-2289-5	11-2293-5	11-2292-5
2X-Large	26 $1/2$ "-29 $1/2$ "	11-2290-6	11-2289-6	11-2293-6	11-2292-6
3X-Large	29 $1/2$ "-32"	11-2290-7	11-2289-7	11-2293-7	11-2292-7

Lateral "J"

- Sewn-in tubular "J" buttress provides inferior and lateral patellar stabilization
- Lateral to medial stabilization strap
- Dual spiral stays offer medial and lateral support



Size	Thigh Circumference (6" above mid-patella)	Product No.	
		Left	Right
Drytex Construction			
X-Small	13"-15 $1/2$ "	11-0660-1-06060	11-0659-1-06060
Small	15 $1/2$ "-18 $1/2$ "	11-0660-2-06060	11-0659-2-06060
Medium	18 $1/2$ "-21"	11-0660-3-06060	11-0659-3-06060
Large	21"-23 $1/2$ "	11-0660-4-06060	11-0659-4-06060
X-Large	23 $1/2$ "-26 $1/2$ "	11-0660-5-06060	11-0659-5-06060
2X-Large	26 $1/2$ "-29 $1/2$ "	11-0660-6-06060	11-0659-6-06060
Neoprene Construction			
X-Small	13"-15 $1/2$ "	11-0321-1-06060	11-0320-1-06060
Small	15 $1/2$ "-18 $1/2$ "	11-0321-2-06060	11-0320-2-06060
Medium	18 $1/2$ "-21"	11-0321-3-06060	11-0320-3-06060
Large	21"-23 $1/2$ "	11-0321-4-06060	11-0320-4-06060
X-Large	23 $1/2$ "-26 $1/2$ "	11-0321-5-06060	11-0320-5-06060
2X-Large	26 $1/2$ "-29 $1/2$ "	11-0321-6-06060	11-0320-6-06060

Hinged Lateral J

- Drytex construction
- Enhanced knee stability with malleable hinges
- Tubular lateral “J” buttress provides stabilization of the patella
- Lateral to medial strapping allows an intimate fit of the “J” buttress along the lateral aspect of the patella

Size	Thigh Circumference (6” above mid-patella)	Drytex		Neoprene	
		Left	Right	Left	Right
X-Small	13”-15 ¹ / ₂ ”	11-0778-1	11-0777-1	11-0776-1	11-0775-1
Small	15 ¹ / ₂ ”-18 ¹ / ₂ ”	11-0778-2	11-0777-2	11-0776-2	11-0775-2
Medium	18 ¹ / ₂ ”-21”	11-0778-3	11-0777-3	11-0776-3	11-0775-3
Large	21”-23 ¹ / ₂ ”	11-0778-4	11-0777-4	11-0776-4	11-0775-4
X-Large	23 ¹ / ₂ ”-26 ¹ / ₂ ”	11-0778-5	11-0777-5	11-0776-5	11-0775-5
2X-Large	26 ¹ / ₂ ”-29 ¹ / ₂ ”	11-0778-6	11-0777-6	11-0776-6	11-0775-6



“H” Buttress

- Sewn in “H” buttress provides inferior, medial, and lateral stabilization
- Superior/inferior straps for compression and support
- Dual spiral stays offer medial and lateral support
- Popliteal cutout

Size	Thigh Circumference (6” above mid-patella)	Product No.
X-Small	13”-15 ¹ / ₂ ”	11-0323-1-06060
Small	15 ¹ / ₂ ”-18 ¹ / ₂ ”	11-0323-2-06060
Medium	18 ¹ / ₂ ”-21”	11-0323-3-06060
Large	21”-23 ¹ / ₂ ”	11-0323-4-06060
X-Large	23 ¹ / ₂ ”-26 ¹ / ₂ ”	11-0323-5-06060
2X-Large	26 ¹ / ₂ ”-29 ¹ / ₂ ”	11-0323-6-06060



Hinged “H” Buttress

- Sewn in “H” buttress provides inferior, medial, and lateral stabilization
- Polycentric hinges offer medial/lateral support and flexion/extension control
- Anterior cover for open patella and hinges
- Superior/inferior straps for compression and support

Size	Thigh Circumference (6” above mid-patella)	Product No.
X-Small	13”-15 ¹ / ₂ ”	11-0828-1-06060
Small	15 ¹ / ₂ ”-18 ¹ / ₂ ”	11-0828-2-06060
Medium	18 ¹ / ₂ ”-21”	11-0828-3-06060
Large	21”-23 ¹ / ₂ ”	11-0828-4-06060
X-Large	23 ¹ / ₂ ”-26 ¹ / ₂ ”	11-0828-5-06060
2X-Large	26 ¹ / ₂ ”-29 ¹ / ₂ ”	11-0828-6-06060
3X-Large	29 ¹ / ₂ ”-32”	11-0828-7-06060



Deluxe Donut

- Sewn-in, trimmable, tubular donut provides inferior, superior, medial, and lateral patellar stabilization
- Superior/inferior straps for compression and support
- Dual spiral stays offer medial and lateral support

Size	Thigh Circum. (6" above mid-patella)	Product No.	
		Deluxe Donut	With Popliteal Cutout
X-Small	13 ¹ / ₂ "-15 ¹ / ₂ "	11-0327-1-06060	11-0380-1-06060
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-0327-2-06060	11-0380-2-06060
Medium	18 ¹ / ₂ "-21"	11-0327-3-06060	11-0380-3-06060
Large	21"-23 ¹ / ₂ "	11-0327-4-06060	11-0380-4-06060
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-0327-5-06060	11-0380-5-06060
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-0327-6-06060	11-0380-6-06060



Horseshoe Patella

- Adjustable, horseshoe-shaped buttress for patellar stabilization
- Superior/inferior straps for compression and support
- Dual spiral stays offer medial and lateral support
- Warm and compressive neoprene construction

Size	Thigh Circum. (6" above mid-patella)	Product No.	
		Horseshoe Patella	With Popliteal Cutout
X-Small	13 ¹ / ₂ "-15 ¹ / ₂ "	11-0334-1-06060	11-0344-1-06060
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-0334-2-06060	11-0344-2-06060
Medium	18 ¹ / ₂ "-21"	11-0334-3-06060	11-0344-3-06060
Large	21"-23 ¹ / ₂ "	11-0334-4-06060	11-0344-4-06060
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-0334-5-06060	11-0344-5-06060
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-0334-6-06060	11-0344-6-06060



Hinged Air DonJoy®

- Inflatable buttresses for medial, lateral or full patellar stabilization
- Polycentric hinges offer medial/lateral support and flexion/extension control
- Wraparound thigh and sleeve calf design facilitates ease of application/removal
- Cool and breathable Drytex construction

Size	Thigh Circumference (6" above mid-patella)	Product No.
X-Small	13"-15 ¹ / ₂ "	11-0855-1-06000
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-0855-2-06000
Medium	18 ¹ / ₂ "-21"	11-0855-3-06000
Large	21"-23 ¹ / ₂ "	11-0855-4-06000
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-0855-5-06000
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-0855-6-06000



Hely & Weber - Patella Stabilizer - Universal Tubular Buttress

- 12" nylon covered neoprene base offers comfort, heat and compression
- Medial and lateral spiral stays limit bunching while allowing full range-of-motion
- Universal tubular buttress is movable and removable for individualized patient care
- Measure at mid-patella with leg fully extended

Size	Knee Circum. (Open and Closed Popliteal)	Product No.		Knee Circum. (Posterior Closure)	Product No. Posterior Closure
		Open Popliteal	Closed Popliteal		
X-Small	12"-13"	3630XS	3730XS	10"-13"	3630PCXS
Small	13"-14"	3630S	3730S	12"-15"	3630PCS
Medium	14"-15"	3630M	3730M	14"-17"	3630PCM
Large	15"-17"	3630L	3730L	16"-19"	3630PCL
X-Large	16"-18"	3630XL	3730XL	18"-22"	3630PCXL
2X-Large	18"-20"	3630XXL	3730XXL	20"-24"	3630PCXXL



Posterior Closure

Hinged Patella Stabilizer - Universal Tubular Buttress

- 12" nylon covered neoprene base offers comfort, heat and compression
- Wide bi-axial hinges provide medial/lateral ligamentous support
- Measure at mid-patella with leg fully extended

Size	Knee Circum. (Open and Closed Popliteal)	Product No.		Knee Circum. (Posterior Closure)	Product No. Posterior Closure
		Open Popliteal	Closed Popliteal		
X-Small	12"-13"	3630HXS	3730HXS	10"-13"	3630HPCXS
Small	13"-14"	3630HS	3730HS	12"-15"	3630HPCS
Medium	14"-15"	3630HM	3730HM	14"-17"	3630HPCM
Large	15"-17"	3630HL	3730HL	16"-19"	3630HPCL
X-Large	16"-18"	3630HXL	3730HXL	18"-22"	3630HPCXL
2X-Large	18"-20"	3630HXXL	3730HXXL	20"-24"	3630HPCXXL



Posterior Closure

Patellofemoral

Patella Stabilizer with Medial/Lateral Buttress

- 13" nylon covered neoprene offers heat, comfort and compression
- Medial and lateral spiral stays limit bunching, while allowing full range-of-motion
- Open popliteal
- Measure at mid-patella with leg fully extended

Size	Knee Circumference	Product No.
X-Small	12"-13"	3671XS
Small	13"-14"	3671S
Medium	14"-15"	3671M
Large	15"-17"	3671L
X-Large	16"-18"	3671XL
2X-Large	18"-20"	3671XXL



Patella Stabilizer with Inferior U Tubular Buttress

- 13" nylon covered neoprene offers comfort, heat and compression
- Medial and lateral spiral stays limit bunching, while allowing full range-of-motion
- Inferior tubular buttress provides patella control
- Measure at mid-patella with leg fully extended

Size	Knee Circum.	Product No.	
		Open Popliteal	Closed Popliteal
X-Small	12"-13"	3641XS	3741XS
Small	13"-14"	3641S	3741S
Medium	14"-15"	3641M	3741M
Large	15"-17"	3641L	3741L
X-Large	16"-18"	3641XL	3741XL
2X-Large	18"-20"	3641XXL	3741XXL



Hinged Patella Stabilizer with Inferior U Tubular Buttress

- 12" nylon covered neoprene offers comfort, heat and compression
- Wide bi-axial hinges provide medial/lateral ligamentous support
- Neoprene hinge covers allow full sports participation
- Measure at mid-patella with leg fully extended

Size	Knee Circum.	Product No.	
		Open Popliteal	Closed Popliteal
X-Small	12"-13"	3641HXS	3741HXS
Small	13"-14"	3641HS	3741HS
Medium	14"-15"	3641HM	3741HM
Large	15"-17"	3641HL	3741HL
X-Large	16"-18"	3641HXL	3741HXL
2X-Large	18"-20"	3641HXXL	3741HXXL

Hinged Patella Stabilizer with Medial-Lateral Buttress with Open Popliteal

- 12" nylon covered neoprene offers heat, comfort and compression
- Wide bi-axial hinges provide medial/lateral ligamentous support
- Neoprene hinge covers allow full sports participation
- Measure at mid-patella with leg fully extended

Size	Knee Circumference	Product No.
X-Small	12"-13"	3671HXS
Small	13"-14"	3671HS
Medium	14"-15"	3671HM
Large	15"-17"	3671HL
X-Large	16"-18"	3671HXL
2X-Large	18"-20"	3671HXXL



Lateral J Patella Stabilizer

- Medial and lateral spiral stays limit bunching and allow a full range-of-motion
- Measure at mid-patella with leg fully extended
- Replace ? with Left (L) or Right (R)

Size	Knee Circumference	Product No.	
		Open Popliteal	Closed Popliteal
X-Small	12"-13"	5690?XS	5790?XS
Small	13"-14"	5690?S	5790?S
Medium	14"-15"	5690?M	5790?M
Large	15"-17"	5690?L	5790?L
X-Large	16"-18"	5690?XL	5790?XL
2X-Large	18"-20"	5690?XXL	5790?XXL



Shields™ Brace

- Kuhl perforated neoprene construction provides comfortable compression
- Calf and thigh strap application adjusts for most leg contours
- Controls malalignments of glide, medial/lateral tilt and anterior/posterior tilt
- Measure at mid-patella with leg extended

Size	Knee Circumference	Product No.
X-Small	10"-12"	5674XS
Small	12"-14"	5674S
Medium	14"-16"	5674M
Large	16"-18"	5674L
X-Large	18"-20"	5674XL



Hinged Shields Brace

- Independent calf and thigh strap application adjusts for most leg contours
- Controls malalignments of glide, medial/lateral tilt and anterior/posterior tilt
- Wide bi-axial hinges enhance medial/lateral stability and prevent hyperextension
- Measure at mid-patella with leg extended

Size	Knee Circumference	Product No.
X-Small	10"-12"	5675XS
Small	12"-14"	5675S
Medium	14"-16"	5675M
Large	16"-18"	5675L
X-Large	18"-20"	5675XL



Medi USA - Genumedi II Knee Support

- Comfortable two-way stretch support
- Lightweight support for weak, sore or swollen knees
- Available in Sand

Indications: Knee sprains and strains, weak, sore or swollen knees and arthritis

Size	Thigh Circum. (6" above center of patella)	Calf Circum. (2" below center of patella)	Product No.
I	14 ¹ / ₂ "-16 ¹ / ₄ "	9 ³ / ₄ "-11"	61171
II	16 ¹ / ₄ "-17"	11"-12 ¹ / ₄ "	61172
III	17"-18 ¹ / ₄ "	12 ¹ / ₄ "-13 ¹ / ₂ "	61173
IV	18 ¹ / ₄ "-19 ¹ / ₄ "	13 ¹ / ₂ "-14 ¹ / ₂ "	61174
V	19 ¹ / ₄ "-20 ¹ / ₂ "	14 ¹ / ₂ "-15 ³ / ₄ "	61175
VI	20 ¹ / ₂ "-22"	15 ³ / ₄ "-17"	61176
VII	22"-23"	17"-18 ¹ / ₄ "	61177



Med Spec - Dynatrack™ Plus Patella Stabilizer with Neoprene Material

- Contains medial and lateral plastic hinges that enhance the brace's ability to conform to the knee during flexion
- Anterior closure design for easy application and optimum fit
- Large popliteal opening with bilateral design

Size	Knee Circumference (taken at mid-patella)	Product No.
X-Small	10"-12"	117001
Small	12"-14"	117002
Medium	14"-16"	117004
Large	16"-18"	117005
X-Large	18"-20"	117006
2X-Large	20"-22"	117007
3X-Large	over 22"	117008



New Options Sports - K2-U Padded Hinge Knee Brace

- Covered multi-position exterior hinges and universal felt buttress
- 3/16" nylon, 2 sides neoprene; upper and lower circumferential strap

Indications: ACL/PCL sprains and tears; Patella pain; Medial/lateral instability; Any joint instability related to sprains or strains associated with athletic activity

Size	Knee Circumference	Product No.
X-Small	12"-13"	K2-U XS
Small	13"-14"	K2-U S
Medium	14"-15"	K2-U M
Large	15"-16"	K2-U L
X-Large	16"-18"	K2-U XL
2X-Large	18"-20"	K2-U XXL
3X-Large	20"-22"	K2-U XXXL
4X-Large	22"+	K2-U XXXXL



K1-U Covered Hinge Knee Brace with Universal Felt Buttress and Criss-Cross Fan Patella Control

- 3/16" nylon, 2 sides neoprene
- Patella opening for positioning and comfort
- Upper and lower circumferential straps with additional straps above and below patella for support

Indications: Medial/lateral instability; Mild to moderate anterior/posterior cruciate instabilities; Meniscus injuries; Support of the knee joints following an injury or after surgery

Size	Knee Circumference	Product No.
X-Small	12"-13"	K1-U XS
Small	13"-14"	K1-U S
Medium	14"-15"	K1-U M
Large	15"-16"	K1-U L
X-Large	16"-18"	K1-U XL
2X-Large	18"-20"	K1-U XXL
3X-Large	20"-22"	K1-U XXXL
4X-Large	22"+	K1-U XXXXL



K4-U Patella Knee Stabilizer with Universal Felt Buttress

- 3/16" patella knee; medial and lateral spiral stays

Size	Knee Circumference	Product No.
X-Small	12"-13"	K4-U XS
Small	13"-14"	K4-U S
Medium	14"-15"	K4-U M
Large	15"-16"	K4-U L
X-Large	16"-18"	K4-U XL
2X-Large	18"-20"	K4-U XXL
3X-Large	20"-22"	K4-U XXXL
4X-Large	22"+	K4-U XXXXL



K14-D Patella Knee Stabilizer with Adjustable Tubular Buttress Straps and Spiral Stays

- 1/8" nylon, 2 sides neoprene
 - Patella opening for positioning and comfort
- Indications:** Medial and lateral support; Hook and loop strapping and internal support stays provide a custom fit

Size	Knee Circumference	Product No.
X-Small	12"-13"	K14-D XS
Small	13"-14"	K14-D S
Medium	14"-15"	K14-D M
Large	15"-16"	K14-D L
X-Large	16"-18"	K14-D XL
2X-Large	18"-20"	K14-D XXL
3X-Large	20"-22"	K14-D XXXL
4X-Large	22"+	K14-D XXXXL



Patellofemoral

K15 Dynamic Patella Stabilizer with “J” Buttress

- 1/8” nylon, 2 sides neoprene
- Patella pull panel with adjustable straps; upper/lower circumferential straps; Medial and lateral spiral stays

Indications: Support and compression for patients suffering from a wide range of patella/femoral disorders; Stabilization for the kneecap during demanding activities

Size	Knee Circumference	Product No.	
		Left	Right
X-Small	12”-13”	K15-L XS	K15-R XS
Small	13”-14”	K15-L S	K15-R S
Medium	14”-15”	K15-L M	K15-R M
Large	15”-16”	K15-L L	K15-R L
X-Large	16”-18”	K15-L XL	K15-R XL
2X-Large	18”-20”	K15-L XXL	K15-R XXL
3X-Large	20”-22”	K15-L XXXL	K15-R XXXL
4X-Large	22”+	K15-L XXXXL	K15-R XXXXL



K18 Osgood Schlater Support with Shock Absorbing V-Shaped Compression Pad

- 1/8” nylon 2 sides neoprene; Upper and lower circumferential straps

Indications: Osgood Schlater syndrome; Patellar tendonitis

Size	Knee Circumference	Product No.
X-Small	12”-13”	K18 XS
Small	13”-14”	K18 S
Medium	14”-15”	K18 M
Large	15”-16”	K18 L
X-Large	16”-18”	K18 XL
2X-Large	18”-20”	K18 XXL
3X-Large	20”-22”	K18 XXXL
4X-Large	22”+	K18 XXXXL



K194 Universal Wraparound Knee with Dynamic Pull Buttress

- 1/8” UBL one side, one side neoprene
- Patella opening for positioning and comfort
- Fits up to 26” thighs
- Universal right or left; one size fits most

Indications: Provides knee support while providing medial and lateral stability; Hook and loop strapping provides user support for normal daily activities and exercise

Product No. **K194**



Ossur Americas - Form Fit Neoprene Support with Stabilized Patella

- Accordion molding in the popliteal area reduces bunching
- Lycra lined

Size	Knee Circum. (at mid-patella)	Universal Buttress	Product No.	
			Sewn In Donut Buttress	Sewn In U- Shaped Buttress
Small	11 ¹ / ₂ "-13"	302563	303563	304563
Medium	13"-14 ¹ / ₂ "	302565	303565	304565
Large	14 ¹ / ₂ "-16"	302567	303567	304567
X-Large	16"-17 ¹ / ₂ "	302568	303568	304568
2X-Large	17 ¹ / ₂ "-19"	302569	303569	304569



Neoprene Wraparound Knee with Patella Support

- Wraparound design
- Universal sizing: one size fits all, right or left
- Internal U-shaped buttress strapping system

Product No. 202100



Palumbo - Patella Stabilizing Brace

- "Y" strap creates a dynamic, medially displacing force to the patella
- Prevents stretching of the medial retinaculum, dislocation and subluxation of the patella
- Reduces arthritic changes and chondromalacia by dissipating abnormal forces
- 1/8" neoprene/nylon on two sides

Size	Circumference			Product No.	
	3" above mid-patella	3" below mid-patella	3" above thigh	Blue	Tan
X-Small	10"-11 ³ / ₄ "	8 ¹ / ₂ "-10"	—	110XSBLUE	110220
Small	11 ³ / ₄ "-13 ¹ / ₂ "	10"-11 ¹ / ₂ "	11 ³ / ₄ "-18 ¹ / ₂ "	110SBLUE	110230
Medium/Small	13 ¹ / ₂ "-15 ¹ / ₄ "	11 ¹ / ₂ "-13"	13 ¹ / ₂ "-20 ¹ / ₄ "	110MSBLUE	110240
Medium	15 ¹ / ₄ "-17"	13"-14 ¹ / ₂ "	15 ¹ / ₄ "-22"	110MBLUE	110250
Medium/Large	17"-19"	14 ¹ / ₂ "-16"	17"-24"	110MLBLUE	110260
Large	19"-21"	16"-17 ¹ / ₂ "	19"-26 ¹ / ₂ "	110LBLUE	110270
X-Large	21"-23 ¹ / ₂ "	17 ¹ / ₂ "-19"	21"-30"	110XLBLUE	110280



Patellofemoral

PATELLOFEMORAL

Palumbo with Dyno Sleeve

- 12" sleeve supports the soft tissues above and below the knee, controlling edema and promote healthy function
- Prevents subluxation and dislocation of the patella
- Stretches contracting patellofemoral ligament and lateral retinaculum and prevents stretching of the medial retinaculum

Indications: Chondromalacia

Size	Circumference		Product No.
	3" above mid-patella	3" below mid-patella	
Small	11 ³ / ₄ "-13 ¹ / ₂ "	10"-11 ¹ / ₂ "	113SWHITEKNIT
Medium	13 ¹ / ₂ "-17"	11 ¹ / ₂ "-14 ¹ / ₂ "	113MWHITEKNIT
Large	17"-21"	14 ¹ / ₂ "-17 ¹ / ₄ "	113LWHITEKNIT



Palumbo Patella Tracker

- Laterally positioned adjustable buttress "floats" maintaining effective seating of the patella both horizontally and vertically during full ROM

Indications: Chondromalacia, patella subluxation, dysplasia, vastus medialis rehab in growing children (postponing or making surgery unnecessary)

Size	Circumference		Product No.	
	3" above patella	3" below patella	Blue	Beige
X-Small	10"-11 ³ / ₄ "	8 ¹ / ₂ "-10"	108XSBLUE	108XSBEIGE
Small	11 ³ / ₄ "-13 ¹ / ₂ "	10"-11 ¹ / ₂ "	108SBLUE	108SBEIGE
Medium/Small	13 ¹ / ₂ "-15 ¹ / ₄ "	11 ¹ / ₂ "-13"	108MSBLUE	108MSBEIGE
Medium	15 ¹ / ₄ "-17"	13"-14 ¹ / ₂ "	108MBLUE	108MBEIGE
Medium/Large	17"-19"	14 ¹ / ₂ "-16"	108MLBLUE	108MLBEIGE
Large	19"-21"	16"-17 ¹ / ₂ "	108LBLUE	108LBEIGE
X-Large	21"-23 ¹ / ₂ "	17 ¹ / ₂ "-19"	108XLBLUE	108XLBEIGE



Premium Palumbo Tracker

- Same as The Palumbo Patella Tracker, but with blue Lycra covered ³/₁₆" neoprene with shorter sleeve for the serious athlete

Size	Circumference		Product No.
	3" above patella	3" below patella	
X-Small	10"-11 ³ / ₄ "	8 ¹ / ₂ "-10"	—
Small	11 ³ / ₄ "-13 ¹ / ₂ "	10"-11 ¹ / ₂ "	109SBLUE
Medium/Small	13 ¹ / ₂ "-15 ¹ / ₄ "	11 ¹ / ₂ "-13"	109MSBLUE
Medium	15 ¹ / ₄ "-17"	13"-14 ¹ / ₂ "	109MBLUE
Medium/Large	17"-19"	14 ¹ / ₂ "-16"	109MLBLUE
Large	19"-21"	16"-17 ¹ / ₂ "	109LBLUE
X-Large	21"-23 ¹ / ₂ "	17 ¹ / ₂ "-19"	109XLBLUE



ProCare (DJO) - Universal Patella Knee with Buttress

- 1/8" neoprene wrap with open patella and spiral stays
- Accommodates buttress pad for additional patellar support
- Measure circumference of knee fully extended

Indications: Mild to moderate compression and patellar control

Size	Knee Circumference	Product No.
Universal	up to 21"	79-82460
X-Large	up to 24"	79-82461



Reinforced Patella Stabilizer

- 1/8" neoprene with sewn-in horseshoe tubular buttress
- Adjustable superior/inferior straps with contact closure for additional compression and stability

Indications: Patella femoral tracking dysfunction; Chondromalacia; Grade 1 collateral ligament sprains

Size	Knee Circumference (6" above mid-patella)	Product No.
Small	15 1/2"-18"	79-94433
Medium	18"-20 1/2"	79-94435
Large	20 1/2"-23"	79-94437
X-Large	23"-25 1/2"	79-94438
2X-Large	25 1/2"-28"	79-94439



Hinged Patella Stabilizer

- 1/8" neoprene with sewn-in tubular buttress and removable dual-axis polycentric hinges for additional medial/lateral support
- Adjustable superior/inferior straps

Indications: Medial/lateral instabilities with patella tendonitis; Chondromalacia; Subluxation; Grade 1 collateral ligament sprains

Size	Knee Circumference (6" above mid-patella)	Product No.
Small	15 1/2"-18"	79-94483
Medium	18"-20 1/2"	79-94485
Large	20 1/2"-23"	79-94487
X-Large	23"-25 1/2"	79-94488
2X-Large	25 1/2"-28"	79-94489



Patellofemoral

Lateral Patella Stabilizer

- 1/8" neoprene with sewn-in tubular "J" buttress and dynamic side-pull compression strap

Indications: Patella misalignment; Dislocation; Subluxation

Size	Knee Circumference	Product No.	
		Left	Right
Small	15 ¹ / ₂ "-18"	79-94473	79-94463
Medium	18"-20 ¹ / ₂ "	79-94475	79-94465
Large	20 ¹ / ₂ "-23"	79-94477	79-94467
X-Large	23"-25 ¹ / ₂ "	79-94478	79-94468
2X-Large	25 ¹ / ₂ "-28"	79-94479	79-94469
3X-Large	28"-31"	79-94479-10	79-94469-10
4X-Large	31"-34"	79-94479-11	79-94469-11

Measure circumference of knee fully extended



Stabilized Knee Support

- 1/8" neoprene with medial/lateral spiral stay for added support
- Compression straps help prevent migration and provide support
- Universal neoprene buttress may be adjusted to allow for proper patella management and stabilization

Indications: Mild patella tendonitis; Chondromalacia; Subluxations

Size	Thigh Circum.	Product No.	
		Closed Pop	Open Pop
X-Small	13 ¹ / ₂ "-15 ¹ / ₂ "	79-82722	79-82752
Small	15 ¹ / ₂ "-18"	79-82723	79-82753
Medium	18"-20 ¹ / ₂ "	79-82725	79-82755
Large	20 ¹ / ₂ "-23"	79-82727	79-82757
X-Large	23"-25 ¹ / ₂ "	79-82728	79-82758
2X-Large	25 ¹ / ₂ "-28"	79-82729	79-82759
3X-Large	28"-31"	79-82729-10	—
4X-Large	31"	79-82729-11	—



Dynamic Pull Universal Knee Wrap

- 1/8" neoprene provides compressive support and warmth
- Adjustable medial/lateral patella pull strap to help improve patella tracking
- 24" circumference

Indications: Patellofemoral tracking dysfunction; Chondromalacia

Description	Product No.
Nylon 1 side	79-82150
Nylon 2 sides	79-82160



Patella Stabilizer with Buttress

- 1/8" neoprene support with adjustable felt horseshoe-buttress or sewn-in donut buttress to accommodate
- Medial/lateral spiral stays and thigh/calf compression straps provide additional stability
- Open popliteal available with donut buttress only

Indications: Mild patella tendonitis; Chondromalacia; Subluxations

Size	Thigh Circumference (6" above mid-patella)	Product No.	
		Horseshoe Buttress	Donut Buttress with Open Pop
X-Small	13 1/2"-15 1/2"	79-92852	79-92842
Small	15 1/2"-18"	79-92853	79-92843
Medium	18"-20 1/2"	79-92855	79-92845
Large	20 1/2"-23"	79-92857	79-92847
X-Large	23"-25 1/2"	79-92858	79-92848
2X-Large	25 1/2"-28"	79-92859	79-92849
3X-Large	28"-31"	79-92859-10	—
4X-Large	31"-34"	79-92859-11	—



Professional Products - Lateral "J" Patella Support

- 1/8" neoprene that provides compression and support
- Open patella and open popliteal with lateral "J" buttress and anterior pad
- Medial/lateral spiral stays for additional stability
- Hook and loop closure

Indications: Subluxation of the patella; Chondromalacia; Tendonitis

Description	Product No.	
	Left	Right
Lateral J Patella Support	05483-T8-C2	05482-T8-C2

Lateral "J" Patella



Truform - Orthotex Knee Stabilizer with Condyle Pad, Spiral Stays

- Adjustable multi-function buttress surrounds patella and provides compression and stabilization
- Spiral stays provide moderate medial/lateral stability
- Opening over back of knee minimizes binding, provides additional wearing comfort
- Encircling straps maintain stabilizer in correct position

Product No. **2541**

Orthotex 2541



Orthotex Knee Stabilizer with Condyle Pad, Hinged Bars

- Polycentric hinged bars with flexion and extension stops provide excellent medial-lateral support and stability
- Adjustable multi-function buttress surrounds patella and provides compression and stabilization
- Opening over back of knee minimizes binding, provides additional wearing comfort
- Encircling straps maintain stabilizer in correct position

Product No. **2543**

Orthotex 2543



Orthotex Knee Support with Stabilizer Pad

- Adjustable multi-function buttress surrounds patella and provides compression and stabilization
- Four-way stretch material provides even compression
- Plush lined interior absorbs moisture, wicks away perspiration
- Sheer elastic opening over back of knee

Product No. **2546**



Trulife - Platinum Open Patella Knee Support

- Open patella knee support with repositionable buttress stabilizes patella
- Medial and lateral boning for added support
- Antibacterial TruTex promotes freshness
- Breathable stretch foam

Indications: Knee strain; Mild arthritis; Runners knee; Osgood-Schlatter disease

Size	Knee Circumference	Product No.
Small	14"-15 ¹ / ₂ "	PT12202
Medium	15 ¹ / ₂ "-18"	PT12203
Large	18"-19 ¹ / ₂ "	PT12204
X-Large	19 ¹ / ₂ "-22"	PT12205



Lenox Hill Neoprene Patella Stabilizer

- Stabilizes patella without compression
- Open popliteal area prevents supports from bunching or binding
- Hook and loop closure
- May be worn under other knee supports

Indications: Patella stabilization in conjunction with a rehabilitation program involving all non-surgical patella dislocation or mis-alignments, including chondromalacia patellae (CMP); Dislocation of the patella (DPL) and Subluxation of the patella (SLP)

Size	Knee Circumference	Product No.
Small - Medium	12"-16"	0837003
Large X-Large	16+"	0837005



Lenox Hill Knee Support with Medial/Lateral Stays

- Double spring steel medial and lateral stays offer added support
- Open popliteal area prevents support from bunching or binding
- Hook and loop closure
- One size fits 7"-20" knee circumference

Indications: Patella stabilization in conjunction with a rehabilitation program involving all non-surgical patella dislocation or mis-alignments, including chondromalacia patellae (CMP), dislocation of the patella (DPL) and subluxation of the patella (SLP)

Product No. **08364**



Bauerfeind - GenuTrain S Knee Support

- Lightweight, breathable, knitted fabric with light, high-tech plastic lateral sidebars
- Silicone insert redistributes pressure from patella to surrounding tissues
- Intermittent compression during movement results in improved metabolism, and reduction in pain and swelling
- Replace ? with 1 (Nature/Beige), 7 (Black), or 8 (Titanium)



Size	Circumference		Product No.	
	4 ³ / ₄ " below knee	5 ¹ / ₂ " above knee	Left	Right
1	11"-12 ¹ / ₄ "	15"-16 ¹ / ₈ "	110413040?0701	110413040?0601
2	12 ¹ / ₄ "-13 ³ / ₈ "	16 ¹ / ₈ "-17 ¹ / ₄ "	110413040?0702	110413040?0602
3	13 ³ / ₈ "-14 ¹ / ₂ "	17 ¹ / ₄ "-18 ¹ / ₂ "	110413040?0703	110413040?0603
4	14 ¹ / ₂ "-15 ³ / ₄ "	18 ¹ / ₂ "-19 ³ / ₄ "	110413040?0704	110413040?0604
5	15 ³ / ₄ "-17"	19 ³ / ₄ "-20 ⁷ / ₈ "	110413040?0705	110413040?0605
6	17"-18 ¹ / ₈ "	20 ⁷ / ₈ "-22"	110413040?0706	110413040?0606

Bell-Horn - Drop Lock Elastic Brace

- Encircling straps
- Beige elastic one way stretch; 12¹/₄"
- Stainless steel hinges with drop locks



Size	Knee Circumference (at knee center)	Product No.
Small	11 ¹ / ₂ "-13"	112 S
Medium	13"-14 ³ / ₄ "	112 M
Large	15"-16 ¹ / ₂ "	112 L
X-Large	16 ³ / ₄ "-18 ¹ / ₄ "	112 XL
2X-Large	18 ¹ / ₄ "-20 ¹ / ₂ "	112 XXL

Bledsoe Brace - Hinged Knee FC (Front Closure)

- Front closure wraparound design for easy application
- Removable polycentric hinges
- Open popliteal for patient comfort

Indications: For sore, strained or sprained knee ligaments



Size	Thigh Circumference (6" above mid-patella)	Product No.
X-Small	15"-17 ¹ / ₂ "	RK171301
Small	17 ¹ / ₂ "-20"	RK171303
Medium	20"-22 ¹ / ₂ "	RK171305
Large	22 ¹ / ₂ "-25"	RK171307
X-Large	25"-27 ¹ / ₂ "	RK171309
2X-Large	27 ¹ / ₂ "-29"	RK171311

HINGED KNEE BRACES

Breg - Economy Hinged Knee

- Full circumferential straps
- Covered hinges are ideal for contact activities
- Length: 13"
- Replace ? with size (1-6) from sizing chart

Description	Product No.	
	Airmesh	Neoprene
Pull-On	2170?	0670?
Open Back, Pull-On	2171?	0672?
Open Back, Pull-On, Wraparound	1123?	—
Open Back w/Adjustable Horseshoe, Pull-On	—	0673?
Open Back, Wraparound	—	1121?



Economy Hinged Knee Sizing Chart

Size	Thigh Circumference (6" above mid-patella)
1 (X-Small)	12"-15"
2 (Small)	15"-18"
3 (Medium)	18"-21"
4 (Large)	21"-24"
5 (X-Large)	24"-27"
6 (2X-Large)	27"-30"

RoadRunner/ShortRunner

- Available in Airmesh and neoprene (ShortRunner available in Airmesh only)
- Low-profile with full circumferential straps
- Polycentric adjustable hinges provide extension control from 10°-40° in 10° increments (10° installed), and flexion control from 45°-90° in 15° increments

Indications: Most ACL, PCL, MCL, LCL injuries and instabilities; Meniscal injuries; Mild osteoarthritis

RoadRunner					
Size	Thigh Circum. (6" above mid-patella)	Closed Back	Open Back	Open Back	
		Pull-On Airmesh	Pull-On Airmesh	Neoprene	Wrap-Around Airmesh
X-Small	12"-15"	14111	14141	07001	14151
Small	15"-18"	14112	14142	07002	14152
Medium	18"-21"	14113	14143	07003	14153
Large	21"-24"	14114	14144	07004	14154
X-Large	24"-27"	14115	14145	07005	14155
2X-Large	27"-30"	14116	14146	07006	14156



RoadRunner

ShortRunner				
Size	Thigh Circum. (6" above mid-patella)	Closed Back	Open Back	
		Pull On	Pull On	Wrap Around
X-Small	12"-15"	21721	21731	21751
Small	15"-18"	21722	21732	21752
Medium	18"-21"	21723	21733	21753
Large	21"-24"	21724	21734	21754
X-Large	24"-27"	21725	21735	21755
2X-Large	27"-30"	21726	21736	21756



ShortRunner

Corflex - Knee Sleeve with Hinge

- Sleeve design with heavy-duty hinges for optimal medial/lateral support
- Manufactured from premium grade $\frac{1}{8}$ " latex-free neoprene

Indications: Medial/lateral instability; Muscle/ligament weakness or injury; Degenerative joint disease

Size	Circumference			Product No.
	6" above mid-patella	mid-patella	6" below mid-patella	
X-Small	16"-18"	12"-13"	11"-12"	88-0302
Small	17"-19"	13"-14"	12"-13"	88-0303
Medium	18"-20"	14"-15"	13"-14"	88-0304
Large	19"-21"	15"-16"	14"-15"	88-0305
X-Large	20"-22"	16"-18"	15"-16"	88-0306
2X-Large	21"-23"	18"-20"	16"-17"	88-0307
3X-Large	24"-26"	20"-22"	17"-18"	88-0308
4X-Large	27"-29"	22"-24"	18"-19"	88-0309

**DeRoyal - Deluxe Hinged Knee Support**

- $\frac{1}{4}$ " neoprene construction for warmth, compression and support
- Medial and lateral hinges with range-of-motion at 0°, 15°, 30° and 45° to limit extension
- Covered hinges allow brace to be worn during contact sports
- Elastic straps above and below knee help prevent migration

Indications: Grade I collateral ligament sprain; Mild ACL/PCL sprain; Grade I varus and valgus instability; Patellar subluxation; Osteoarthritis

Size	Knee Circumference (taken below patella)	Product No.
Small	13"-14"	NE7722-72
Medium	14"-15"	NE7722-73
Large	15"-17"	NE7722-74
X-Large	17"-18"	NE7722-75
2X-Large	18"-20"	NE7722-76
3X-Large	20"-22"	NE7722-77



HINGED KNEE BRACES

DonJoy (DJO) - Economy Hinged Knee

- 4 mm neoprene construction in sleeve style or cool and breathable Drytex construction
- Polycentric hinge with 10³/₄" aluminum uprights provides light support for medial/lateral ligaments
- Uprights are removable for washing
- Available in Black



Neoprene Economy Hinged Knee				
Size	Thigh Circumference (6" above mid-patella)	Product No.		
		Standard	w/Popliteal	w/Inferior Buttress
X-Small	13"-15 ¹ / ₂ "	11-0160-1	11-0372-1	11-0378-1
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-0160-2	11-0372-2	11-0378-2
Medium	18 ¹ / ₂ "-21"	11-0160-3	11-0372-3	11-0378-3
Large	21"-23 ¹ / ₂ "	11-0160-4	11-0372-4	11-0378-4
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-0160-5	11-0372-5	11-0378-5
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-0160-6	11-0372-6	11-0378-6

Drytex Economy Hinged Knee					
Size	Thigh Circumference (6" above mid-patella)	Product No.			
		Standard Sleeve	Sleeve w/Popliteal	Standard Wraparound	Wraparound w/Popliteal
X-Small	13"-15 ¹ / ₂ "	11-0670-1	11-0671-1	11-0672-1	11-0673-1
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-0670-2	11-0671-2	11-0672-2	11-0673-2
Medium	18 ¹ / ₂ "-21"	11-0670-3	11-0671-3	11-0672-3	11-0673-3
Large	21"-23 ¹ / ₂ "	11-0670-4	11-0671-4	11-0672-4	11-0673-4
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-0670-5	11-0671-5	11-0672-5	11-0673-5
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-0670-6	11-0671-6	11-0672-6	11-0673-6
3X-Large	29 ¹ / ₂ "-32"	—	—	11-0672-7	11-0673-7

Hinged Knee Support

- Hinged support combines polycentric hinges and condyle pads to provide medial and lateral stability in a 1¹/₄" neoprene support
- Sewn in leather pull-up loops allow for ease of application
- Extension Stops: 12°, 20°, 30°, and 40° (5° installed)

Size	Thigh Circumference	Product No.
X-Small	13"-15 ¹ / ₂ "	11-0006-1
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-0006-2
Medium	18 ¹ / ₂ "-21"	11-0006-3
Large	21"-23 ¹ / ₂ "	11-0006-4
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-0006-5
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-0006-6
3X-Large	29 ¹ / ₂ "-32 ¹ / ₂ "	11-0006-7



Hinged Knee

- Designed for mild to moderate medial lateral control
- Neoprene or Drytex construction
- Sleeve or Wraparound styles for optimal fit and comfort
- Available with or without popliteal cutout



Drytex Hinged Knee, Standard Sleeve					
Size	Thigh Circum. (6" above mid-patella)	Standard	Popliteal Cutout	Horseshoe Buttress	Popliteal Cutout & Horse- shoe Buttress
X-Small	13"-15 ¹ / ₂ "	11-0555-1	11-0556-1	11-2010-1	11-2011-1
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-0555-2	11-0556-2	11-2010-2	11-2011-2
Medium	18 ¹ / ₂ "-21"	11-0555-3	11-0556-3	11-2010-3	11-2011-3
Large	21"-23 ¹ / ₂ "	11-0555-4	11-0556-4	11-2010-4	11-2011-4
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-0555-5	11-0556-5	11-2010-5	11-2011-5
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-0555-6	11-0556-6	11-2010-6	11-2011-6
3X-Large	29 ¹ / ₂ "-32"	—	—	11-2010-7	11-2011-7

Drytex Hinged Knee, Wraparound			
Size	Thigh Circum. (6" above mid-patella)	Product No.	
		Wraparound	w/Popliteal Cutout
X-Small	13"-15 ¹ / ₂ "	11-2012-1	11-2013-1
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-2012-2	11-2013-2
Medium	18 ¹ / ₂ "-21"	11-2012-3	11-2013-3
Large	21"-23 ¹ / ₂ "	11-2012-4	11-2013-4
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-2012-5	11-2013-5
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-2012-6	11-2013-6
3X-Large	29 ¹ / ₂ "-32"	11-2012-7	11-2013-7

Neoprene Hinged Knee					
Size	Thigh Circum. (6" above mid-patella)	Horseshoe Buttress	Popliteal Cutout & Horse- shoe Buttress	Wraparound	Wraparound w/Popliteal
X-Small	13"-15 ¹ / ₂ "	11-2006-1	11-2007-1	11-2008-1	11-2009-1
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-2006-2	11-2007-2	11-2008-2	11-2009-2
Medium	18 ¹ / ₂ "-21"	11-2006-3	11-2007-3	11-2008-3	11-2009-3
Large	21"-23 ¹ / ₂ "	11-2006-4	11-2007-4	11-2008-4	11-2009-4
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-2006-5	11-2007-5	11-2008-5	11-2009-5
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-2006-6	11-2007-6	11-2008-6	11-2009-6
3X-Large	29 ¹ / ₂ "-32"	11-2006-7	11-2007-7	11-2008-7	11-2009-7

HINGED KNEE BRACES

Playmaker™

- Sleeve or wraparound styles for optimal fit and comfort
- Popliteal cutout for improved fit and comfort
- Versatile strapping allows ACL, PCL or CI configurations
- Standard hinge
- Replace ? with size (1-7) from sizing chart

Description	Product No.	
	Neoprene	Drytex
Sleeve style	11-0860-?	11-0557-?
Sleeve style with popliteal cutout	11-0861-?	11-0558-?
Sleeve style with popliteal cutout & patella donut	11-0863-?	11-1395-?
Wraparound style	11-0864-?	11-0666-?
Wraparound with popliteal cutout	11-0865-?	11-0759-?
Droplock with popliteal cutout	—	11-1482-?
Sleeve Style With Donut	11-0862-?	11-1396-?



Playmaker Sizing Chart

Size	Thigh Circumference (6" above mid-patella)
1 (X-Small)	13"-15½"
2 (Small)	15½"-18½"
3 (Medium)	18½"-21"
4 (Large)	21"-23½"
5 (X-Large)	23½"-26½"
6 (2X-Large)	26½"-29½"
7 (3X-Large)	29½"-32"

Hely & Weber - Velocity with Patella Stabilization

- Kuhl™ perforated neoprene which breathes well and offers warmth and maximum compression
- Range-of-motion hinges feature a polycentric adjustable design to control extension and flexion with comfortable shock absorbing condular pads and shells which ensure patient comfort and compliance
- Extension stops allow 10°, 20°, 30°, 40° and 50° settings, flexion stops at 10°, 50°, 60°, 70°, 80°, 90° and 100° increments
- Offers effective support for ligament/meniscus tears while having the patented "u-shaped" patella locator buttress and adjustable patella control of the Shields™ design

Size	Thigh Circumference (6" above patella)	Product No.
X-Small	13"-15"	5645PS-XS
Small	14"-17"	5645PS-S
Medium	15"-18"	5645PS-M
Large	17"-20"	5645PS-L
X-Large	19"-22"	5645PS-XL
2X-Large	21"-23"	5645PS-XXL
3X-Large	22"-24"	5645PS-XXXL



12" Velocity Hinged Knee

- Kuhl™ perforated neoprene to maximize patient comfort/compliance
- Range-of-motion hinges with stops included; brace can function as a knee immobilizer, full ROM brace or anywhere in between
- Extension stops allow 20°, 30°, and 40° settings, flexion stops at 45°, 60°, 75°, 90°
- Can be locked at full extension

Indications: Treatment of ACL, PCL, MCL and LCL sprains; Postoperative and non-operative meniscus injuries

Size	Thigh Circumference (6" above patella)	Product No.
Pull-Up Knee		
X-Small	13"-16"	5645XS
Small	15"-18"	5645S
Medium	16"-19"	5645M
Large	18"-20"	5645L
X-Large	19"-23"	5645XL
2X-Large	21"-24"	5645XXL
3X-Large	22"-25"	5645XXXL
Anterior Closure Knee		
X-Small	12"-16"	5646XS
Small	14"-18"	5646S
Medium	16"-20"	5646M
Large	17"-22"	5646L
X-Large	19"-25"	5646XL
2X-Large	23"-28"	5646XXL
3X-Large	25"-29"	5646XXXL



12" Knapp Hinged Knee Orthosis (NoSwet®)

- NoSwet non-neoprene construction
- Elasticized band in thigh strap ensures comfort
- Hinges are fully covered for sports activity
- Measure circumference at mid-patella with leg fully extended

Indications: Pre/post-operative care of MCL, LCL and meniscal repairs

Size	Knee Circumference	Product No.
X-Small	10"-12"	7656XS
Small	12"-14"	7656S
Medium	14"-16"	7656M
Large	16"-18"	7656L
X-Large	18"-20"	7656XL
2X-Large	20"-22"	7656XXL



HINGED KNEE BRACES

Lateral “J” Stabilizer with Hinge and Condylar Pads

- J buttress, in conjunction with anterior flap, produces dynamic patella control
- Biaxial hinges with condyle pads provide medial and lateral stability
- Kuhl perforated neoprene allows moisture to escape and the skin to remain dry



Size	Circumference	Product No.	
		Left	Right
Open Popliteal			
X-Small	13"-15"	5694-LT-XS	5694-RT-XS
Small	14"-17"	5694-LT-S	5694-RT-S
Medium	15"-18"	5694-LT-M	5694-RT-M
Large	17"-20"	5694-LT-L	5694-RT-L
X-Large	19"-22"	5694-LT-XL	5694-RT-XL
2X-Large	21"-23"	5694-LT-XXL	5694-RT-XXL
3X-Large	22"-24"	5694-LT-XXXL	5694-RT-XXXL
Closed Popliteal			
X-Small	13"-15"	5794-LT-XS	5794-RT-XS
Small	14"-17"	5794-LT-S	5794-RT-S
Medium	15"-18"	5794-LT-M	5794-RT-M
Large	17"-20"	5794-LT-L	5794-RT-L
X-Large	19"-22"	5794-LT-XL	5794-RT-XL
2X-Large	21"-23"	5794-LT-XXL	5794-RT-XXL
3X-Large	22"-24"	5794-LT-XXXL	5794-RT-XXXL

Axis Hinged Knee Sleeve with Open Popliteal

- 1/8" neoprene with 2-sided nylon outer covering
- Hinges feature hyperextension stops
- Covered hinges allow full sports participation

Size	Knee Circumference	Product No.
X-Small	12"-13"	3655XS
Small	13"-14"	3655S
Medium	14"-15"	3655M
Large	15"-17"	3655L
X-Large	16"-18"	3655XL
2X-Large	18"-20"	3655XXL



3655

Knapp™ Hinged Knee Orthosis

- Kuhl perforated neoprene for cool, comfortable compression
- Easily accessible hinge allows interchangeability of hinge styles
- Hinges are fully covered for sports activity
- Measure circumference at mid-patella with leg fully extended

Indications: Pre/post-operative care of MCL, LCL and meniscal repairs; 16" length can be used for mild to moderate ACL/PCL sprains and strains or as an adjunct to patella fracture

Size	Knee Circumference	Product No.	
		12" Length	16" Length
Standard Closure			
X-Small	10"-12"	5656XS	5658XS
Small	12"-14"	5656S	5658S
Medium	14"-16"	5656M	5658M
Large	16"-18"	5656L	5658L
X-Large	18"-20"	5656XL	5658XL
2X-Large	20"-22"	5656XXL	5658XXL
Heavy Hinges			
X-Small	10"-12"	5656HHXS	5658HHXS
Small	12"-14"	5656HHS	5658HHS
Medium	14"-16"	5656HHM	5658HHM
Large	16"-18"	5656HHL	5658HHL
X-Large	18"-20"	5656HHXL	5658HHXL
2X-Large	20"-22"	5656HHXXL	5658HHXXL
Range-of-Motion			
X-Small	10"-12"	5656ROM	5658ROM
Small	12"-14"	5656ROM	5658ROM
Medium	14"-16"	5656ROM	5658ROM
Large	16"-18"	5656ROM	5658ROM
X-Large	18"-20"	5656ROM	5658ROM
2X-Large	20"-22"	5656ROM	5658ROM



5656



5658

Med Spec - Gripper™ 12" Hinged Knee

- Provides medial and lateral stability and resists distal migration
- Internal SkinLoc™ pad works in concert with outer closure strap to grip proximal end of the calf muscle
- Hinge pockets reinforced with ballistic nylon to enhance durability
- CoolFlex material is breathable and hypoallergenic

Size	Knee Circumference (taken at mid-patella)	Product No.
X-Small	10"-12"	117031
Small	12"-14"	117032
Medium	14"-16"	117034
Large	16"-18"	117035
X-Large	18"-20"	117036
2X-Large	20"-22"	117037
3X-Large	22"-24"	117038
4X-Large	24"+	117039



HINGED KNEE BRACES

New Options Sports - K64 “Knee Mate™” Hinged

- 3/16” nylon, 2 sides neoprene; 1/8” Koolflex version
- Condyle pads behind hinges
- Upper/lower circumferential straps with finger loops and pull tabs
- Covered stainless steel light-duty hinges with hyperextension stop

Size	Knee Circumference	Product No.	
		Neoprene	KoolFlex™
X-Small	12”-13”	K64 XS	KC64 XS
Small	13”-14”	K64 S	KC64 S
Medium	14”-15”	K64 M	KC64 M
Large	15”-17”	K64 L	KC64 L
X-Large	16”-18”	K64 XL	KC64 XL
2X-Large	18”-20”	K64 XXL	KC64 XXL
3X-Large	20”-22”	K64 XXXL	KC64 XXXL
4X-Large	22”+	K64 XXXXL	KC64 XXXXL



K67 “The Hybrid™” Half Sleeve-Half Wrap Hinged Brace Pull Up Series

- 3/16” nylon, 2 sides neoprene
- Anterior wrap around closure on top, closed bottom around calf area; better fit for the cone shaped leg
- Upper/lower circumferential straps with posterior adjustable elastic segment
- Adjustable covered stainless steel hinges with hyperextension stop

Size	Knee Circumference	Product No.
X-Small	12”-13”	K67 XS
Small	13”-14”	K67 S
Medium	14”-15”	K67 M
Large	15”-17”	K67 L
X-Large	16”-18”	K67 XL
2X-Large	18”-20”	K67 XXL
3X-Large	20”-22”	K67 XXXL
4X-Large	22”+	K67 XXXXL



Ossur Americas - Form Fit® Neoprene Hinged Knee Support

- Dual action hinges are molded into the brace to maintain a low-profile
- Accordion molding in the popliteal area reduces bunching behind the knee
- Lycra lined

Size	Knee Circumference (at mid-patella)	Product No.	
		Black	Blue
Small	11 1/2”-13”	302153BLKC	302153BLU
Medium	13”-14 1/2”	302155BLKC	302155BLU
Large	14 1/2”-16”	302157BLKC	302157BLU
X-Large	16”-17 1/2”	302158BLKC	302158BLU
2X-Large	17 1/2”-19”	302159BLKC	302159BLU



Rebound® Knee Brace

- Breathable 3-layer Cooltech fabric reduces heat retention and dries ultrafast
- Stretches in extension and relaxes to original shape in flexion, minimizing popliteal bunching
- Removable, universal patella buttress included
- Dual lower cruciate straps enable fine-tuned control of anterior/posterior hinge placement and cinch above the gastroc to reduce migration
- Replace ? with 1 (long) or 0 (short)

Size	Circumference (6" above mid-patella)	ROM Hinge		Non-ROM Hinge	
		Sleeve	Wrap	Sleeve	Wrap
X-Small	13"-15 ¹ / ₂ "	708?52	707?52	704?52	703?52
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	708?53	707?53	704?53	703?53
Medium	18 ¹ / ₂ "-21"	708?55	707?55	704?55	703?55
Large	21"-23 ¹ / ₂ "	708?57	707?57	704?57	703?57
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	708?58	707?58	704?58	703?58
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	708?59	707?59	704?59	703?59
3X-Large	29 ¹ / ₂ "-32"	—	707?60	—	703?60



Ultrawrap™

- Accutrac hinges with semi-rigid uprights and interchangeable extension stop
- Textured interior grips skin surface to reduce migration and stimulate neuroreceptors
- Removable, semi-rigid uprights provide medial/lateral support
- Non-binding, breathable popliteal cutout improves comfort

Indications: Mild to moderate strains or ACL, MCL, LCL and PCL sprains; UltraWrap ROM model designed for support with flexion and extension control

Size	Knee Circumference (around center of knee cap)	Product No.
X-Small	12"-13"	UWXS
Small	13"-14"	UWS
Medium	14"-15"	UWM
Large	15"-16"	UWL
X-Large	16"-18"	UWXL
2X-Large	18"-20"	UWXXL
3X-Large	20"-22"	UW3XL
4X-Large	22"-24"	UW4XL



HINGED KNEE BRACES

Neoprene Wraparound Hinged Knee Support

- Open popliteal area reduces bunching
- Front hook and loop closures
- Wraparound style with removable dual-action hinges
- 1/8" neoprene

Size	Knee Circumference (at mid-patella)	Product No.
Small	13"-14"	212003
Medium	14"-15"	212005
Large	15"-16"	212007
X-Large	16"-17"	212008
2X-Large	17"-18"	212009
3X-Large	19"-24"	212010
4X-Large	25"-26"	212011



ProCare (DJO) - Hinged Knee Support

- Compressive neoprene, available with and without buttress
- Universal buttress may be adjusted to allow for medial, lateral, inferior, or superior patella stabilization
- Sturdy removable dual axis hinges and encircling contact closure straps provide mild support to medial/ lateral ligaments

Size	Thigh Circum. (6" above mid-patella)	Product No.		
		3/16" with Open Pop.	1/8" with Reinforced Pad	1/8" with Univ. Buttress
X-Small	13 ¹ / ₂ "-15 ¹ / ₂ "	79-82152	79-82162	79-82732
Small	15 ¹ / ₂ "-18"	79-82153	79-82163	79-82733
Medium	18"-20 ¹ / ₂ "	79-82155	79-82165	79-82735
Large	20 ¹ / ₂ "-23"	79-82157	79-82167	79-82737
X-Large	23"-25 ¹ / ₂ "	79-82158	79-82168	79-82738
2X-Large	25 ¹ / ₂ "-28"	79-82159	79-82169	79-82739
3X-Large	28"-30 ¹ / ₂ "	79-82159-10	—	79-82739-10
4X-Large	30 ¹ / ₂ "-33"	—	—	79-82739-11



Reddie® Brace

- 1/8" neoprene wrap-around configuration and strap tabs
- Posterior strap adjustment allows for one-time adjustment and proper positioning of hinges
- Removable dual axis polycentric hinges provide medial/ lateral support

Size	Thigh Circumference (6" above mid-patella)	Product No.
X-Small	13 ¹ / ₂ "-15 ¹ / ₂ "	79-82392
Small	15 ¹ / ₂ "-18"	79-82393
Medium	18"-20 ¹ / ₂ "	79-82395
Large	20 ¹ / ₂ "-23"	79-82397
X-Large	23"-25 ¹ / ₂ "	79-82398
2X-Large	25 ¹ / ₂ "-28"	79-82399
3X-Large	28"-30 ¹ / ₂ "	79-82399-10



Only available with Open Popliteal

U.S. Orthotics - Anterior Opening Knee Support With Hinges and Open Popliteal

- Wrap-around neoprene with polycentric side hinges
- 2" adjustable hook and loop straps
- Excellent for patients who cannot pull on a traditional support

Size	Knee Circumference	Product No.
X-Small	12"	KW6-XS
Small	13"	KW6-SM
Medium	14"	KW6-MD
Large	15"	KW6-LG
X-Large	16"	KW6-XL
2X-Large	17"	KW6-XXL
3X-Large	18"	KW6-3XL



KNEE SUPPORTS

Alps - Knee Brace

- Fabric airslots in the popliteal region for less bunching during knee flexion and edema control

Size	Measurements (4" above knee cap while standing)	Hinged Knee		Without Hinges	
		Pull-On	Wrap-Around	Pull-On	Wrap-Around
X-Small	11"-13"	KBHP-XS	KBHW-XS	KBP-XS	KBW-XS
Small	13 ¹ / ₂ "-15"	KBHP-S	KBHW-S	KBP-S	KBW-S
Medium	15 ¹ / ₂ "-17"	KBHP-M	KBHW-M	KBP-M	KBW-M
Large	17 ¹ / ₂ "-22"	KBHP-L	KBHW-L	KBP-L	KBW-L
X-Large	22 ¹ / ₂ "-24"	KBHP-XL	KBHW-XL	KBP-XL	KBW-XL
2X-Large	24 ¹ / ₂ "-27"	KBHP-XXL	KBHW-XXL	KBP-XXL	KBW-XXL
3X-Large	27 ¹ / ₂ "-30"	KBHP-XXXL	KBHW-XXXL	KBP-XXXL	KBW-XXXL



Knee Supports

Bauerfeind - GenuTrain® Knee Support

- Lightweight, breathable, knitted fabric to enhance proprioception
- Circular silicone insert redistributes pressure from patella to surrounding soft tissues
- Insert provides intermittent compression resulting in improved metabolism
- Replace ? with 1 (Natural/Beige), 7 (Black), or 8 (Titanium)

Size	Circumference		Product No.
	4 ³ / ₄ " below knee	5 ¹ / ₂ " above knee	
1	11"-12 ¹ / ₄ "	15"-16 ¹ / ₈ "	110412050?0001
2	12 ¹ / ₄ "-13 ³ / ₈ "	16 ¹ / ₈ "-17 ¹ / ₄ "	110412050?0002
3	13 ³ / ₈ "-14 ¹ / ₂ "	17 ¹ / ₄ "-18 ¹ / ₂ "	110412050?0003
4	14 ¹ / ₂ "-15 ³ / ₄ "	18 ¹ / ₂ "-19 ³ / ₄ "	110412050?0004
5	15 ³ / ₄ "-17"	19 ³ / ₄ "-20 ⁷ / ₈ "	110412050?0005
6	17"-18 ¹ / ₂ "	20 ⁷ / ₈ "-22"	110412050?0006
7	18 ¹ / ₈ "-19 ¹ / ₄ "	22"-23 ¹ / ₄ "	110412050?0007



Bell-Horn - 12" Elastic Knee Sleeve

- Recommended for stabilization of the knee joint when ligaments or muscles have been damaged or weakened
- Made of one-way stretch elastic
- Can be used as a preventative support for athletes

Size	Knee Circumference (at knee center)	Product No.
Small	11 ¹ / ₂ "-13"	114 S
Medium	13"-14 ³ / ₄ "	114 M
Large	15"-16 ¹ / ₂ "	114 L
X-Large	16 ³ / ₄ "-18 ¹ / ₄ "	114 XL



DeRoyal - Open Patella Knee Support

- ¹/₈" neoprene construction
- Available with oval pad for additional patella control
- Open patella for patient comfort

Indications: Grade I collateral ligament sprain; Knee sprains and strains; Osteoarthritis; Mild joint effusion; Baker's cyst; Bursitis; Tendonitis

Size	Thigh Circumference (6" above mid-patella)	Product No.	
		Without Pad	With Pad
X-Small	13 ¹ / ₂ "-15 ¹ / ₂ "	NE7702-71	NE7702-01
Small	15 ¹ / ₂ "-18"	NE7702-72	NE7702-02
Medium	18"-20 ¹ / ₂ "	NE7702-73	NE7702-03
Large	20 ¹ / ₂ "-23"	NE7702-74	NE7702-04
X-Large	23"-25 ¹ / ₂ "	NE7702-75	NE7702-05
2X-Large	25 ¹ / ₂ "-28"	NE7702-76	NE7702-06
3X-Large	28"-30 ¹ / ₂ "	NE7702-77	NE7702-07
4X-Large	30 ¹ / ₂ "-33"	—	NE7702-08



DonJoy (DJO) - Adjustable Patella Donut

- Drytex construction
- Sewn-in, trimmable tubular donut provides inferior, superior, medial and lateral patellar stabilization
- Popliteal cutout

Size	Thigh Circumference (6" above mid-patella)	Product No.
X-Small	13"-15 ¹ / ₂ "	11-0554-1-06000
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-0554-2-06000
Medium	18 ¹ / ₂ "-21"	11-0554-3-06000
Large	21"-23 ¹ / ₂ "	11-0554-4-06000
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-0554-5-06000
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-0554-6-06000



Knee Support

- High-quality neoprene construction
- 4 mm black neoprene with double sided nylon

Size	Thigh Circum. (6" above mid-patella)	Deluxe Knee	With Patella Cut-Out	With Patella and Popliteal Cut-Out
X-Small	13"-15 ¹ / ₂ "	11-0008-1-06000	11-0358-1-06000	11-0371-1-06000
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-0008-2-06000	11-0358-2-06000	11-0371-2-06000
Medium	18 ¹ / ₂ "-21"	11-0008-3-06000	11-0358-3-06000	11-0371-3-06000
Large	21"-23 ¹ / ₂ "	11-0008-4-06000	11-0358-4-06000	11-0371-4-06000
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-0008-5-06000	11-0358-5-06000	11-0371-5-06000
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-0008-6-06000	11-0358-6-06000	11-0371-6-06000



Drytex Knee Support

- A neoprene alternative with a nylon core and polyester Lycra fabric that allows for improved airflow and breathability
- Maintains warmth and provides compression and support
- Popliteal cut-out
- Available in Black

Size	Thigh Circumference (6" above mid-patella)	Product No.
X-Small	13"-15 ¹ / ₂ "	11-0548-1-06000
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-0548-2-06000
Medium	18 ¹ / ₂ "-21"	11-0548-3-06000
Large	21"-23 ¹ / ₂ "	11-0548-4-06000
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-0548-5-06000
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-0548-6-06000
3X-Large	29 ¹ / ₂ "-32 ¹ / ₂ "	11-0548-7-06000



Spiral Elastic Knee Support

- Seamless two-way stretch cotton elastic
- Open patella
- Even compression without heat
- Available in White

Size	Thigh Circumference (6" above mid-patella)	Product No.
Small	15 ¹ / ₂ " -18 ¹ / ₂ "	11-0240-2
Medium	18 ¹ / ₂ " -21"	11-0240-3
Large	21"-23 ¹ / ₂ "	11-0240-4
X-Large	23 ¹ / ₂ " -26 ¹ / ₂ "	11-0240-5



Universal Basic Knee

- Recommended when a sleeve design cannot be applied due to swelling, wound dressings, or loss of motion
- Double sided nylon; Open patella
- Available in Black
- Universal size up to 26¹/₂"
- Measure 6" above mid-patella

Product No. **11-0267-5-06000**



Ossur Americas - Form Fit Neoprene Knee Sleeve

- Accordion molding in popliteal area reduces bunching
- Available in Black or Blue with open or closed patella

Size	Knee Circum. (at mid-patella)	Blue Product No.		Black Product No.	
		Open Patella	Closed Patella	Open Patella	Closed Patella
1/8" Neoprene Knee Sleeve					
Small	11 ¹ / ₂ "-13"	302013BLU	312013BLU	302013BLK	312013BLK
Medium	13"-14 ¹ / ₂ "	302015BLU	312015BLU	302015BLK	312015BLK
Large	14 ¹ / ₂ "-16"	302017BLU	312017BLU	302017BLK	312017BLK
X-Large	16"-17 ¹ / ₂ "	302018BLU	312018BLU	302018BLK	312018BLK
2X-Large	17 ¹ / ₂ "-19"	302019BLU	312019BLU	302019BLK	312019BLK
1/4" Neoprene Knee Sleeve					
Small	11 ¹ / ₂ "-13"	402013BLU	412013BLU	402013BLK	412013BLK
Medium	13"-14 ¹ / ₂ "	402015BLU	412015BLU	402015BLK	412015BLK
Large	14 ¹ / ₂ "-16"	402017BLU	412017BLU	402017BLK	412017BLK
X-Large	16"-17 ¹ / ₂ "	402018BLU	412018BLU	402018BLK	412018BLK
2X-Large	17 ¹ / ₂ "-19"	402019BLU	412019BLU	402019BLK	412019BLK



ProCare (DJO) - Universal Knee Wrap with Closed Patella

- 1/8" neoprene wraparound style
 - Medial and lateral spiral stays for compressive support
- Indications:** Mild to moderate compression and patellar control following strains, sprains or other minor indications

Description	Circumference (6" above mid-patella)	Product No.
Nylon one side	21"	79-82431
Nylon two sides	21"	79-82441



Truform - Knee Support with Condyle Pads

- Completely open front with opposing contact closure straps permits easy application over tender or painfully swollen knees
- Heavy duty elastic material delivers firm compression to help reduce swelling and minimize pain
- Bilateral pressure pads cradle the kneecap for added stability and help maintain proper joint alignment

Size	Circumference (around bend of the knee)	Product No.
X-Small	9 ¹ / ₂ " - 11 ¹ / ₂ "	2545-XS
Small	11 ¹ / ₂ " - 13"	2545-S
Medium	13" - 14 ³ / ₄ "	2545-M
Large	14 ³ / ₄ " - 16 ¹ / ₂ "	2545-L
X-Large	16 ¹ / ₂ " - 18 ¹ / ₂ "	2545-XL
2X-Large	18 ¹ / ₂ " - 21"	2545-2L
3X-Large	21" - 23 ¹ / ₂ "	2545-3L



Orthotex Stabilizer Wrap with Spiral Stays

- Opens flat for easy application over injured, swollen or tender knees
- Spiral stays provide moderate medial/lateral stability, yet flex and extend naturally with the joint
- Opening over back of knee minimizes binding, provides additional wearing comfort

Size	Circumference (around bend of the knee)	Product No.
X-Small	12" - 13"	2542-XS
Small	13" - 14"	2542-S
Medium	14" - 15"	2542-M
Large	15" - 16"	2542-L
X-Large	16" - 17"	2542-XL
2X-Large	17" - 18"	2542-2L
3X-Large	18" - 19"	2542-3L
4X-Large	19" - 20"	2542-4L
5X-Large	20" - 22"	2542-5L



Trulife - Lenox Hill Neoprene Knee Sleeve

- 1/16" neoprene knee sleeve with patella opening
- Sleeve can be worn individually or as an undersleeve

Indications: Soft tissue injuries; Strains; Sprains

Size	Knee Circumference	Product No.
X-Small	12" - 13"	USPRBK-000
Small	13" - 14"	USPRBK-010
Medium	14" - 15"	USPRBK-020
Large	15" - 16"	USPRBK-030
X-Large	16" - 18"	USPRBK-040
2X-Large	18" - 20"	USPRBK-050
3X-Large	20" - 22"	USPRBK-060
4X-Large	22" - 24"	USPRBK-070



AliMed - Turnbuckle Knee Orthosis

- Steel rod and turnbuckle mechanism provide a low-load stretch to the joint to gradually reduce contracture and progress the range-of-motion incrementally
- Kydex® splint shell is hinged at the knee
- Use for knee contractures from 45° flexion to 10° hyperextension
- Pre-set at 30° flexion and 25° hyperextension

Size	Circumference (around knee cap)	Product No.
Small	13"-15"	66425/NA/B
Medium	15"-17"	66425/NA/C
Large	17"-18"	66425/NA/D
X-Large	18"-20"	66425/NA/E



Anatomical Concepts - Knee Management Orthosis

- For an adjustable range of stabilization of the knee in the sagittal plane as the patient's conditions warrant
- Pre- and post-operative management of bony and soft tissue injuries
- For immediate prophylactic protection at end range until treatment plans are finalized
- Flexion up to 110°

Indications: Application following hamstring releases; Ligament repairs; Total knee arthroplasty; Chronic systemic conditions

Description	Height Adjustment	KMO Product No.	KMO with Varum/Valgum Adjustment
Pediatric	12 ³ / ₄ "-17"	850	860
Adult	17"-26"	950	960
Stabilizing Shells		Replacement Liners	
Pediatric	P85	L85	
Adult	P95	L95	

Size	Length	Calf Circum.	Thigh Circum.
Sizing Chart for KMO			
Pediatric	12 ¹ / ₂ "-16 ¹ / ₂ "	8"-15"	12"-20"
Adult	17"-26"	11"-21"	15"-28 ¹ / ₂ "
Sizing Chart for KMO with Varum/Valgum Adjustment			
Pediatric	13 ¹ / ₂ "-17 ¹ / ₂ "	8"-15"	12"-20"
Adult	18"-27"	11"-21"	15"-28 ¹ / ₂ "



Positional/
Stretching KO

DeRoyal - Pucci® Inflation Knee Orthosis with Hinge

- Innovative therapeutic orthosis
- Provides continuous low stretch therapy to the knee
- Inflatable double upright air bladders

Size	Thigh Circumference	Calf Circumference	Product No.
2X-Small	10 ¹ / ₂ "-12 ¹ / ₂ "	8 ¹ / ₂ "-10"	4270AA
X-Small	13"-15"	10"-11 ¹ / ₂ "	4270A
Small	15 ¹ / ₂ "-18"	11 ¹ / ₂ "-14"	4270B
Medium	18 ¹ / ₂ "-22"	14 ¹ / ₂ "-16 ¹ / ₂ "	4270C
Large	22 ¹ / ₂ "-25"	17"-18 ¹ / ₂ "	4270D



Pro-Glide® Knee

- Dynamic splint offers flexion and extension assist in one unit
- Aircraft quality aluminum frame is strong yet lightweight
- Tension lockout at 0°, 45°, and 90°
- Universal left or right

Indications: Joint stiffness; ROM limitations; Soft tissue contractures

Size	Thigh Circum.	Calf Circum.	Replacement Softgoods	Product No.
Ped. Medium	up to 14"	8"-10 ¹ / ₂ "	4107CK	4107C
Ped. Large	14"+	10 ¹ / ₂ "-14"	4107DK	4107D
Small	14"-16 ¹ / ₂ "	10"-13 ¹ / ₂ "	4105BK	4105B
Medium	16 ¹ / ₂ "-21 ¹ / ₂ "	13 ¹ / ₂ "-16 ¹ / ₂ "	4105CK	4105C
Large	21 ¹ / ₂ "-27"	16 ¹ / ₂ "-20"	4105DK	4105D



Lenjoy - Comfy Goniometer Knee Orthosis - "Set ROM"

- Treats non-fixed contractures of the knee
- Allows for progressive extension in 10° increments
- Thigh Circumference: 19"-23"; Calf Circumference: 16"-20"
- Length: 14"

Product No. **GK-101**



Comfy Spring Loaded Goniometer Knee Orthosis - "Set ROM w/Extension"

- Dynamic splint works to help increase knee extension following surgery, CVA, injury or neurological complications
- Allows for progressive extension in 10° increments
- Thigh Circumference: 19"-23"; Calf Circumference: 16"-20"
- Length: 14"

Product No. **SGK-101**



Comfy Knee Orthosis - "Bend to Fit"

- Excellent for the flaccid or weak extremity or to support/immobilize a painful extremity
- Splint spine can bend to the desired ROM; cuffs are also bendable to allow for optimal fit
- Thigh Circumference: 19"-23"; Calf Circumference: 16"-20"
- Length: 14"

Product No. **K-101**



Positional/
Stretching KO

OCSI/Neuroflex - DynaPro™ Flex Knee

- Designed to treat mild/moderate flexion contractures
- Provides dynamic low load prolonged stretch
- Negative non-linear “flex” properties accommodate involuntary muscle contractions
- Proprietary gel padding in cuff and gel knee pad included

Size	Thigh Circum.	Calf Circum.	Product No.
Small	11”-15”	10”-13”	DPFK20725-S
Medium	15”-19”	12”-15”	DPFK20725-M
Large	19”-21”	14”-18”	DPFK20725-L



RCAI - Premier Knee Orthosis

- Progressive static lock settings permit gradual increases in knee joint extension and positioning
- Lightweight nylon hinges can accommodate patients with 0 to moderate tone and mild spasticity; thermoplastic cuffs
- Length: Small, Medium and Large, 16”; Child and Youth, 12¹/₂”

Size	Measure Mid Thigh Circum.	Product No.	
		Premier KCO	Kodel
Small	14”-17”	27KCO-S	23KCO-S
Medium	17”-20”	27KCO-M	23KCO-M
Large	20”-23”	27KCO-L	23KCO-L
Child	10”-12”	27KCO-C	23KCO-C
Youth	12”-14”	27KCO-Y	23KCO-Y
Knee Cap (for use with Premier Knee Orthosis, Kodel)			
Youth/Child		27KC-Y/C	
Small/Medium		27KC-S/M	
Large		27KC-L	
3-point pressure leverage			



Knee Cap

Positional/
Stretching KO

MISCELLANEOUS

Breg - Undersleeve

- Available in cotton or 1/8” neoprene,
- Fits comfortably under Breg’s functional braces and are ideal as a barrier for sensitive skin, as well as for soft tissue containment
- Cotton undersleeve has closed patella; Cotton/Lycra blend - no elastic
- Neoprene undersleeve has open patella

Size	Length		Product No.	
	Cotton	Neoprene	Cotton	Neoprene
X-Small	19”	17”	09851	07351
Small	19”	18”	09852	07352
Medium	19”	19”	09853	07353
Large	19”	20”	09854	07354
X-Large	19”	21”	09855	07355
2X-Large	19”	22”	09856	—



DonJoy (DJO) - Defiance Sports Cover

- Provides a cover and padding over a knee brace to protect the brace and others from hitting the brace
- Polyester with neoprene padding around the frame and hinges
- Compatible with Defiance, Defiance III, Armor, Legend, Female Fource, 4titude, ACL Everyday, OA Defiance, and OA Adjuster
- Available in Black

Size	Circumference (6" above mid-patella)	Product No.	
		Standard (16" long)	Short (14" long)
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-0016-2-06000	11-0103-2-06000
Medium	18 ¹ / ₂ "-21"	11-0016-3-06000	11-0103-3-06000
Large	21"-23 ¹ / ₂ "	11-0016-4-06000	11-0103-4-06000
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-0016-5-06000	11-0103-5-06000
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-0016-6-06000	11-0103-6-06000



Neoprene/Lycra Undersleeves

- These undersleeves provide skin protection offered in either a thin Cotton/Lycra material or perforated neoprene
- Reduces skin abrasions and prevents rashes
- Compatible with Defiance, Defiance III, Armor, Legend, Female Fource, 4titude, ACL Everyday, OA Defiance, OA Adjuster, OA Everyday
- Available in Black

Size	Circumference (6" above mid-patella)	Product No.	
		Neoprene	Lycra
X-Small	13"-15 ¹ / ₂ "	11-0122-1	—
Small	15 ¹ / ₂ "-18 ¹ / ₂ "	11-0122-2	11-0075-2-00000
Medium	18 ¹ / ₂ "-21"	11-0122-3	11-0075-3-00000
Large	21"-23 ¹ / ₂ "	11-0122-4	11-0075-4-00000
X-Large	23 ¹ / ₂ "-26 ¹ / ₂ "	11-0122-5	11-0075-5-00000
2X-Large	26 ¹ / ₂ "-29 ¹ / ₂ "	11-0122-6	11-0075-6-00000



Knit-Rite - KAFO Interface Socks

- Sold each

Description	Width/Length	White	Pink	Purple
Child Long	2 ¹ / ₂ " x 18"	1SP1AFCL	1S21AFCL	1S31AFCL
Adult Narrow	4" x 2 ¹ / ₂ " x 24"	1SP1NAXL	—	—
Adult Regular	6" x 3 ¹ / ₂ " x 28"	1SP1RGXL	—	—

KAFO Interface Socks with X-Static - The Silver Fiber

- Sold each

Description	Width/Length	Navy	Gray	Black
Child Long	2 ¹ / ₂ " x 18"	1NX1AFCL	1SX1AFCL	1BX1AFCL
Adult Narrow	4" x 2 ¹ / ₂ " x 24"	—	—	1BX1NAXL
Adult Regular	6" x 3 ¹ / ₂ " x 28"	—	—	1BX1RGXL



Knee Interface

- Knee support garment to be worn as an interface with knee orthoses or alone for compressive support
- Moisture wicking CoolMax® or X-STATIC® fibers keep skin cool and dry
- Less bulky than neoprene garment; better air circulation
- Measure thigh circumference at largest part of leg to be covered

Flat Width Size	Length	Thigh Circum.	White CoolMax/Lycra/ Anti-Microbial	Black Polyester/ X-STATIC	Gray Polyester/ X-STATIC
Straight Sleeves					
4"	18"	11"-14"	2PCS0418	2BCX0418	2PCX0418
5"	18"	14"-17"	2PCS0518	2BCX0518	2PCX0518
6"	18"	17"-20"	2PCS0618	2BCX0618	2PCX0618
7"	18"	20"-23"	2PCS0718	2BCX0718	2PCX0718
7"	20"	20"-23"	—	2BCX0720	—
8"	20"	23"-29"	2PCS0820	2BCX0820	2PCX0820
9"	20"	26"-34"	2PCS0920	2BCX0920	2PCX0920
Tapered Sleeves					
5" x 4"	18"	14"-17"	2PCST518	2BCXT518	2PCXT518
6" x 4"	18"	17"-20"	2PCST618	2BCXT618	2PCXT618
7" x 5"	18"	20"-23"	2PCST718	2BCXT718	2PCXT718
8" x 6"	18"	23"-29"	2PCST818	2BCXT818	2PCXT818
9" x 7"	18"	26"-34"	2PCST918	2BCXT918	2PCXT918



Breg - Patella Tendon Strap

- Strapped distal to the knee cap, the rubber buttress can ease the tension of chronic inflammation at the insertion points

Indications: Patellar tendonitis (Jumper's knee)

Size	Knee Circumference (at joint line distal to patella)	Product No.
X-Small	12"-13"	96531
Small	13"-14"	96532
Medium	14"-15"	96533
Large	15"-16"	96534
X-Large	16"-18"	96535



Cho-Pat® - Knee Strap

- Symptomatic relief of chondromalacia patella
- Puts pressure on the tendon below the kneecap to guide it into the groove; Strap elevates the kneecap slightly
- Available in Tan with hook and loop closure

Size	Circumference (below patella)	Product No.
X-Small	up to 10"	CP13
Small	10"-12 ¹ / ₂ "	CP15
Medium	12 ¹ / ₂ "-14 ¹ / ₂ "	CP17
Large	14 ¹ / ₂ "-16 ¹ / ₂ "	CP19
X-Large	16 ¹ / ₂ "-18"	CP21
2X-Large	18"-21"	CP22



Dual Action Knee Strap

- Strengthens kneecap mechanism by applying pressure to tendons above and below kneecap
- Reduces forces of the quadriceps on the tendon attachments and erosion of kneecap under surface
- Reduces development of iliotibial band and overuse syndrome

Size	Knee Circumference (at mid-patella)	Product No.
X-Small	10"-12"	CP-01
Small	12"-14"	CP-02
Medium	14"-16"	CP-03
Large	16"-18"	CP-04
X-Large	18"-20"	CP-05
2X-Large	20"-22"	CP-07



Hely & Weber - Universal Matt Strap

- Works for arm, knee, upper arm and upper leg tendonitis
- Trimmable strap allows customization from 7"-20" circumference
- Movable loop pod allows for sizing and custom fitting
- One size fits all (Measure around affected area and trim to fit)

Product No. **3709**



Matt Strap™ Shock Absorption

- Useful for patella tendonitis, "Jumpers Knee," and minor knee pain
- Shock absorption buttress
- Loop & lock design provides focused compression over the shock absorbing buttress
- For sizing measure circumference of affected area

Size	Circumference	Product No.
X-Small	7"-8"	3704XS
Small	8"-9"	3704S
Medium	9"-10 ¹ / ₂ "	3704M
Large	10 ¹ / ₂ "-12"	3704L
X-Large	12"-14"	3704XL
2X-Large	13"-16"	3704XXL



DonJoy (DJO) - McGuire Knee Wrap with Open Patella

- Designed with a hook engageable material and stretch band to easily wrap around the affected area
- Non-sterile
- Includes both cold pad and wrap
- Size: 10.85" x 11.34"

Product No. **11-1268-9-0000**



IceMan® Cold Therapy Unit Model 1100

- Two-position locking top; Water tight seal and hose stirrup
- Inlet temperature control and monitoring
- Compatibility with a variety of pads
- Universal one size; Available in White/Blue

Product No. **11-0494**



Ice-Man Cold Pads

- Even distribution throughout the pad
- Sterile and non-sterile pads
- Universal Sizes: 12" x 11.34"; Ankle Size: 12.26" x 8.28"
- Available in White

Description	Product No.
Universal Cold Pad, X-Large (Sterile)	11-0804-9-00000
Universal Cold Pad, X-Large (Non-Sterile)	11-0682-9-00000
Ankle Cold Pad (Sterile)	11-0815-9-00000



Dura Soft™ Knee Wrap

- Packs provide 2 hours of cold with one application
- Removable inserts can be exchanged for extended cryotherapy
- The ice packs are malleable to provide comfortable cold therapy
- Ultra-wick fabric with antimicrobial protection
- Universal; Includes 2 ice inserts

Product No. **11-0913-0-02000**



ProCare (DJO) - Surround® Patella Strap

- Helps relieve anterior knee pain caused by patellar tendonitis, chondromalacia and Osgood Schlatter
- Primary body of strap is soft, latex-free nylon/foam laminate with contact closures
- Unique cylindrical Floam™ bladder inside the body of the strap may be adjusted for added compression and customized fit

Size	Calf Circumference (just below patella)	Product No.
X-Small	less than 10"	79-80322
Small	10 ¹ / ₂ "-12 ¹ / ₂ "	79-80323
Medium	13"-15"	79-80325
Large	15 ¹ / ₂ "-17 ¹ / ₂ "	79-80327
X-Large	18"+	79-80328

