


# limbs&things





The medical simulation market, which in 1990 was in an embryonic state, is now a bustling adolescent with a very loud voice.

Our fast moving and talented team is continually focused on maintaining first place in the design, manufacture and sale of part-task trainers.

We have enabled safer training for healthcare professionals, and a more comfortable journey for the patient. A big contribution.

A huge thank you to all our dedicated people who have worked and continue to work at Limbs & Things, ensuring that we serve the changing requirements of medical training and reach all corners of the world.








OUR MISSION

IMPROVE PATIENT OUTCOMES THROUGH REALISTIC LEARNING EXPERIENCE





- CATHETERISATION TRAINERS
- EPIDURAL & LUMBAR PUNCTURE RANGE
- PROMPT FLEX DARK SKIN
- PROMPT FLEX CERVICAL
   CERCLAGE MODULE
- PERINEAL REPAIR & EPISIOTOMY











# CONTENTS



8 Female examination

12 male examination

16 abdominal examination

# PROCEDURAL SKILLS

20 venepuncture

24 INJECTION TECHNIQUE

24 central line management

25 paracentesis

26 CATHETERISATION

28 CHEST TRAUMA MANAGEMENT

30 DIAGNOSTIC PROCEDURES



34 BASIC SURGICAL SKILLS

**36** MINOR SURGERY

38 OBSTETRICS & MIDWIFERY

48 gynaecology

# RANGES

50 undergraduate medicine

51 NURSE TRAINING


52 paramedic training


53 osce training

54 physycian associate

55 obstetrics & midwifery

# PHYSICAL EXAMINATION SKILLS











# Clinical Female Pelvic Trainer (CFPT) Mk 3 - Standard & Advanced

#### 60900

STANDARD



The CFPT Mk 3 is an anatomically accurate and tactile representation of the female pelvis. It is the ideal platform for hands-on examination as well as diagnosis of female conditions and minor pathologies. It can be used for many levels of training from undergraduate onwards, as well as in family health.



Validated in the International Journal of Gynaecology and Obstetrics. Volume 124, pp. 270-3, 2014.

#### SKILLS

- Recognition of perineal and pelvic anatomy including bony landmarks
- Digital vaginal examination
- Bi-manual examination
- Cervical smear procedure (including use of speculum)
- Digital rectal examination
- Professional-to-patient communication







 Bi-manual examination
 Use of speculum and cervical broom



Abdominal wall lifts off...
... allowing uterine module to be lifted out

#### QUALITIES

#### ANATOMY

- Abdomen, pelvis and genitalia
- Vagina, cervix, anus and lower bowel
- Advanced model comes with interchangeable uterine modules with different complications
- Modules:
  - Normal Nulliparous Cervix
  - Large Fibroid Nulliparous Ectropion Cervix
  - Small Fibroid Nulliparous Polyp Cervix
  - Ovarian Cyst Multiparous Cervix
  - Retroverted Multiparous Cervix
  - 10-12 Weeks Pregnant
  - 14-16 Weeks Pregnant

#### REALISM

- Abdominal wall with fat layer makes palpation realistic
- Labia can be parted realistically
- Each uterus is presented at the correct anatomical angle

#### **KEY FEATURES**

- Partial thighs aid anatomical orientation
- Soft and strong perineum and labia

#### VERSATILITY

- Abdominal wall can be removed quickly and easily
- Modules can be changed quickly and easily

#### SAFETY

• Latex free

#### **CLEANING & CARE**

• Skin surface is washable using soap and water





60921 MODULE 2 Large Fibroid - Nulliparous Ectropion Cervix



60924 MODULE 5 Retroverted - Multiparous Cervix



60922 MODULE 3 Small Fibroid - Nulliparous Polyp Cervix



MODULE 6 10-12 Weeks Pregnant



60923 MODULE 4 Ovarian Cyst - Multiparous Cervix



60926 MODULE 7 14-16 Weeks Pregnant

#### WHAT'S INCLUDED 60900 60905 60910 CFPT Abdominal Insert 60920 CFPT Module 1: Normal - Nulliparous Cervix $\checkmark$ 60921 CFPT Module 2: Large Fibroid - Nulliparous Ectropion Cervix × 1 **60922** CFPT Module 3: Small Fibroid - Nulliparous Polyp Cervix × $\checkmark$ 60923 CFPT Module 4: Ovarian Cyst - Multiparous Cervix × $\checkmark$ 60924 CFPT Module 5: Retroverted - Multiparous Cervix × $\checkmark$ 60925 CFPT Module 6: 10-12 Weeks Pregnant × $\checkmark$ × $\checkmark$ 60926 CFPT Module 7: 14-16 Weeks Pregnant V. $\checkmark$ 2 Lower torso $\checkmark$ $\checkmark$ -Base $\checkmark$ $\checkmark$ Tube of aqueous gel 4 $\checkmark$ ✓(x2) -Carry case

# Breast Examination Trainer - Standard, Advanced & Basic







Our Breast Examination Trainer provides a highly realistic learning platform for acquiring the skills required to perform Clinical Breast Examination (CBE).

Featuring 6 readily interchangeable and multi-positional pathologies, providing healthcare professionals with the tools to identify various complications and pathologies, including carcinomas, cysts, fibrocystic disease and fibroadenoma.

Both *Simulated Patient* and *bench top* training can be used in alignment with healthcare professionals promoting best practice techniques to trainees and patients, or for any undergraduate programs running **OSCE**s.



This product was developed in close collaboration with leading clinicians:

- Mr Peter Donnelly, Torbay Hospital, UK
- Clinical Skills Resources Centre, University of Liverpool, UK
- Clinical Practice Centre, St James's University Hospital, Leeds, UK

#### QUALITIES

#### ANATOMY

- Realistic soft tissue breast anatomy
- Pathologies supplied: carcinomas: 2cm, 3cm, 5cm, cyst, fibrocystic disease, fibroadenoma

#### REALISM

- Soft tissue breasts look and feel realistic
- Clavicular and axilla pads for accurate lymph node placement

#### **KEY FEATURES**

- Pathologies can be placed in various predetermined location points and are easily changeable
- Dual purpose product: bench top and Simulated Patient
- Hard torso available for bench top use

#### VERSATILITY

• Comfortable to wear for long periods of time during OSCEs and assessements

#### SAFETY

• Latex free

#### CLEANING & CARE

• Skin surface is washable using soap and water



- Soft tissue breasts look and feel realistic
- Pathologies fit securely into supports
- Clavicular and axilla pads for accurate lymph node placement
- Breasts worn by a simulated patient
- Self Breast Examination training











#### SKILLS

- Clinical Breast Examination (CBE)
- Self Breast Examination (SBE)
- Identification of anatomical landmarks
- Identification of lymph nodes (axillary, supra & infraclavicular)
- Location and diagnosis of pathologies
- Professional-to-patient communication

WHAT'S I	NCLUDED	40201	40200	40202
_	Wearable Examination Breasts	$\checkmark$	$\checkmark$	$\checkmark$
40204	Pathologies	$\checkmark$	$\checkmark$	$\checkmark$
40205	Breast Examination Inserts	$\checkmark$	$\checkmark$	$\checkmark$
40206	Breast Pathologies Supports	$\checkmark$	$\checkmark$	$\checkmark$
40207	Breast Back Plates	$\checkmark$	$\checkmark$	$\checkmark$
40208	Breast Examination Torso	$\checkmark$	$\checkmark$	×
40203	Lymph Node Pads and Lymph Nodes	$\checkmark$	×	×

#### OPTIONAL EXTRAS

40209 Medium Carry Case

# Male Rectal Examination Trainer - Standard & Advanced



The Rectal Examination Trainer offers realistic, repeatable training in the core skill of digital rectal or PR examination.

Both models include a contractible anal sphincter, allowing for anal tone assessment and a faecal impaction clip in module.



The **Standard** model focuses on learning core rectal nursing procedures and excludes the pathologies.

The **Advanced** model comes with a range of prostates and 2 perineums, allowing trainees to become familiar with normal and abnormal findings and acquire key diagnostic skills.

External ano-rectal illustrations are ideal for incorporation into scenarios and **OSCE** assessments.

#### SKILLS

- Digital examination of the anus, rectum and prostate
- Assessment of anal tone
- Identification of faecal matter in the rectum
- Identification of external ano-rectal conditions\*
- Professional-to-patient communication
- Insertion of rectal medication (micro enemas, suppositories)

#### QUALITIES

\* Advanced model

#### ΑΝΑΤΟΜΥ

- Buttocks, anus, rectum, prostate and perineum
- Pathological perineum with polyp and rectal cancer\*
- 4 abnormal prostates\*:
- Benign unilateral enlarged Benign bilateral enlarged
- Unilateral carcinoma Bilateral carcinoma

#### REALISM

- Soft, partable buttocks
- Realistic anus with resting tone and ability to contract

#### **KEY FEATURES**

- Simulation of sphincter contraction allows trainees to assess anal tone
- Interchangeable prostates can be quickly and easily inserted, out of sight of the trainee\*
- Left lateral positioning
- Addition of impacted faecal matter allows trainees to recognise and distinguish this common finding
- High quality illustration pack of 9 external anal conditions

#### VERSATILITY

• Trainer can also be presented in 'semi-standing' position using the optional **Standing Position Stand** (60172)

#### SAFETY

• Latex free

#### **CLEANING & CARE**

• Skin surface is washable using soap and water



WHAT'S	WHAT'S INCLUDED		60171	
60177	Base Unit	$\checkmark$	$\checkmark$	
60174	Left Lateral Stand	$\checkmark$	$\checkmark$	
60175	Normal Perineum with Anal Tone	$\checkmark$	$\checkmark$	
60176	Pathological Perineum Anal Tone	×	$\checkmark$	
60179	Normal Prostate	$\checkmark$	$\checkmark$	
60180	Benign Unilateral Enlarged Prostate	×	$\checkmark$	
60181	Benign Bilateral Enlarged Prostate	×	$\checkmark$	
60182	Unilateral Carcinoma Prostate	×	$\checkmark$	
60183	Bilateral Carcinoma Prostate	×	$\checkmark$	
60184	Impacted Faecal Matter	$\checkmark$	$\checkmark$	
60173	External Conditions Illustration Pack	×	$\checkmark$	
00293	Water-based Lubricant	$\checkmark$	$\checkmark$	

#### OPTIONAL EXTRAS

60172	Standing Position Stand
40209	Medium Carry Case

# Clinical Male Pelvic Trainer (CMPT) Mk 2 - Standard & Advanced

#### 60950

STANDARD



#### SKILLS

- Correct examination procedure of male pelvis
- Examination and evaluation of:
  - Normal anatomy
  - Testicular abnormalities
  - Abdominal and pelvic pain
- Dry catheterisation
- Professional-to-patient communication





*Realistic training experiencePenile cancer* 





Pelvic examination, supine position
 Diagnosis of pathologies

An ideal platform for teaching and learning hands-on male pelvic examination and diagnosis.

The user-friendly design, combined with clear anatomic landmarks and a range of male conditions and pathologies creates a realistic training experience for students across all levels of healthcare education from undergraduate level upwards.



Validated in the *Journal of Surgical Education* Volume 74, Issue 3, May–June 2017.



Developed in collaboration with:

- The Clinical Skills Resource Centre, University of Liverpool, UK
- Imperial College, London, UK
- Guy's Hospital, London, UK
- Southmead Hospital, Bristol, UK
- University of Southampton, UK

#### QUALITIES

#### ΑΝΑΤΟΜΥ

• Abdomen, pelvis and genitalia

Anatomy of the groin with clear anatomical landmarks:

- Anterior Superior Iliac Spine (ASIS)
- Pubic symphysis and tubercles

Genitalia includes:

- Penis (both circumcised and non-circumcised)
- Scrotum (containing testicles with epididymis and vas deferens)

Pathologies (Advanced model only):

- Varicocele
- Testicular tumour and penile cancer
- Epididymal cyst
- Hydrocele
- Orchitis/Epididymo-orchitis
- Indirect inguinal hernia

#### REALISM

• Soft tissue genitals and abdominal wall

#### VERSATILITY

- Soft tissue inserts can be removed and replaced
- Can be used in 2 positions (standing and supine)

#### SAFETY

• Modules 5 & 7 contain latex

#### **CLEANING & CARE**

• Skin surface is washable using soap and water





60953

MODULE 2 Varicocele



60956 MODULE 5 Hydrocele



60954

MODULE 3 Testicular Tumours & Penile Cancer



60957 MODULE 6 Orchitis/Epididymo-orchitis



60955

MODULE 4 Epididymal Cysts



60958 MODULE 7 Indirect Inguinal Hernial

60951

60950

#### WHAT'S INCLUDED

60952	CMPT Module 1: Normal	$\checkmark$	$\checkmark$
60953	CMPT Module 2: Varicocele	×	$\checkmark$
60954	CMPT Module 3: Tumours	×	$\checkmark$
60955	CMPT Module 4: Epididymal Cysts	×	$\checkmark$
60956	CMPT Module 5: Hydrocele	×	$\checkmark$
60957	CMPT Module 6: Orchitis/Epididymo-orchitis	×	$\checkmark$
60958	CMPT Module 7: Indirect Inguinal Hernia	×	$\checkmark$
60959	CMPT Abdominal Insert	$\checkmark$	$\checkmark$
60960	CMPT Pubic Bone	$\checkmark$	✓(x2)
60961	CMPT Foreskins Pack (each pack contains 2 foreskins & 3 lubricant sachets)	$\checkmark$	✓(x2)
-	Lower torso	$\checkmark$	$\checkmark$
	Support base	$\checkmark$	$\checkmark$
60963	LED pen torch	×	$\checkmark$
	Carry case	$\checkmark$	✓(x2)

Patents granted: • US6336812 • EP099027 • DE69804793 • CA 2293585 • ZL98806331 • HK1023832 • 755575 Australia

# Abdominal Examination Trainer

#### 60000

An anatomically accurate adult male torso, used to teach and practise the palpation, auscultation and percussion elements of abdominal or gastrointestinal (GI) examination.

With interchangeable organs of varying sizes, and an integral MP3 player providing realistic sounds, the trainer provides the ability to recognise and differentiate a range of abnormal organs and pathologies.



#### SKILLS

- Familiarity with the abdominal regions and underlying anatomy
- Practice in abdominal palpation, auscultation and percussion
- Ability to differentiate normal from abnormal anatomy
- Identification of ascites, specifically shifting dullness and fluid thrill
- Identification of gaseous distension and bowel obstruction
- Ballottement of kidneys
- Professional-to-patient communication











#### QUALITIES

#### REALISM

- Organs feel realistic on palpation and respond appropriately to percussion
- Lifelike quality of abdominal skin accommodates stretching for gaseous distension and ascites simulation
- Realistic ballottement of enlarged kidney

#### VERSATILITY

- Abdominal skin is simple to remove, allowing quick and easy changeover of organs
- Suitable for both bench top use and hybrid use with Simulated Patient
- Model can be rolled onto side for examination of ascites and shifting dullness

#### **KEY FEATURES**

- Interchangeable organs and pathologies of varying sizes can be quickly and easily inserted into the abdomen
- Trainer or Simulated Patient can vary respiratory movement of liver and spleen by the turn of a wheel
- Integral amplifier and surface mounted MP3 player allow training in auscultation of normal and high-pitched or obstructed bowel sounds, renal and aortic bruits in variable locations
- Volume of bowel sounds can be adjusted and MP3 player allows for additional sounds to be loaded
- Distension set allows for the following:
  - Realistic checking for ascites using percussion, shifting dullness or fluid thrill technique
  - Simulation of bowel obstruction on percussion and auscultation
- Pulse bulb allows simulation of normal and aneurysmal aortic pulse

#### SAFETY

• Latex free

#### **CLEANING & CARE**

• Skin surface is washable using soap and water

#### ANATOMY

- Torso featuring abdomen, pelvis and lower part of thorax
- Bony landmarks include ribs, costal margin, xiphisternum, pubic crest and anterior superior iliac spines
- 3 Livers: slightly enlarged, enlarged with smooth edge and enlarged with irregular edge
- 2 Spleens: slightly enlarged and markedly enlarged
- 2 Enlarged Kidneys
- Distended Bladder
- 2 Aortas: normal and aneurysmal
- Set of 6 Abdominal Pathologies includes:
  - 4 Smooth masses: 2 sizes, each with hard and soft version - 2 Irregular hard masses
- Distension Set includes:
  - Ascites bag
  - Gaseous distension bag, pump and foam insert
- Simplified representation of lower thoracic & lumbar spine

#### WHAT'S INCLUDED

60005	Abdominal Examination Base Unit
60006	Abdominal Examination Skin
60007	Abdominal Examination Rectus Bag
60008	Abdominal Examination Organ Locator
60009	Abdominal Wall
60010	Abdominal Examination Livers (Set of 3)
60011	Abdominal Examination Spleens (Set of 2)
60012	Abdominal Examination Kidneys (x2)
60013	Abdominal Examination Bladder
60014	Abdominal Examination Pathologies (Set of 6)
60015	Abdominal Examination Aortas (Set of 2)
60016	Abdominal Examination - Distension Set
60053	International 5V DC Power Supply
-	AA Batteries (x4)
-	Carry case (x2)



You may also be interested in: 60002 Abdominal Examination Trainer Value Set Visit <u>limbsandthings.com</u>

# SEE HOW A PRODUCT WORKS

SHORT & INFORMATIVE PRODUCT VIDEOS

5 LANGUAGE OPTIONS


HANDS ON DEMONSTRATIONS

LATEST PRODUCTS





# PROCEDURAL SKILLS



# Standard Venepuncture Arm

#### 00330 00331

The Standard Venepuncture Arm provides the ideal platform for practising the skills related to venepuncture and IV cannulation. Comes with a bag and stand.



# M rans M rans M rans

The bag & stand mock blood supply is perfect for a clinical skills training environment where large numbers of trainees are practising on a regular basis, or for self-directed learning time in a skills laboratory

# Bag & Stand Advanced Venepuncture Arm



The Bag & Stand Advanced Venepuncture Arm contains a more complex vein structure in combination with the bag, stand and mock blood supply.



## Advanced Venepuncture Arm



The Advanced Venepuncture Arm features a complex vein system, which is connected to a closed blood supply, for easy adjustment of pressure.





#### SKILLS

- Venepuncture
- IV Cannulation
- IV Infusion
- Professional-to-patient communication







- Visual inspection and palpation
- Vacuum blood collection systems to be used
- Skin to be stretched to 'anchor' veins for easier cannula insertion

#### QUALITIES

#### ANATOMY

- Soft, flexible skin
- Underlying palpable veins:
  - All models: cephalic, basilic and dorsal metacarpal
  - Advanced model only: median cubital

#### KEY FEATURES

- Realistic blood flashback
- Arm allows for digital pressure to stem blood flow
- Infusion tube allows for volume IV fluids to be administered

#### VERSATILITY

- Skin is durable
- Veins last for up to 500 insertions of 21g needle
- Blood has been reformulated to improve leak resistance of veins (concentrated liquid to save you money and space)
- Can be used with vacuum blood collection systems, needle and syringe and IV cannulas
- Can be used for hybrid simulation and for professionalto-patient communication when using the optional Venepuncture Arm Harness

#### SAFETY

• Veins contain latex

#### **CLEANING & CARE**

- Veins are easy to replace
- Easy to clean, service and maintain

		00330	00440	00290
AT'S I	NCLUDED	00331	00441	00298
00292	Advanced Venepuncture Arm Veins	×	✓	~
00295	Advanced Venepuncture Vein Module	×	$\checkmark$	$\checkmark$
00333	Standard Venepuncture Arm Veins	$\checkmark$	×	×
00334	Standard Venepuncture Vein Module	$\checkmark$	×	×
00335	Venepuncture Arm Infusion Tube	$\checkmark$	$\checkmark$	$\checkmark$
00336	Bag & Stand Venepuncture Arm Mock Blood Supply	$\checkmark$	$\checkmark$	×
00337	Pressurised Venepuncture Arm Mock Blood Supply	×	×	$\checkmark$
00020	Concentrated Venous Blood - Starter Pack	$\checkmark$	$\checkmark$	$\checkmark$
00302	Venepuncture Arm Rest	$\checkmark$	$\checkmark$	$\checkmark$
00303	Refill Bottle	×	×	$\checkmark$
-	Venepuncture arm skin	$\checkmark$	$\checkmark$	$\checkmark$
-	Water-based lubricant pack (x3)	$\checkmark$	$\checkmark$	$\checkmark$
-	Arm shell	$\checkmark$	$\checkmark$	$\checkmark$

#### OPTIONAL EXTRAS



# ACF Pad - Venepuncture





This soft tissue strap-on pad can be used to train in venepuncture and represents the antecubital fossa of the right arm.

#### SKILLS

- Professional-to-patient communication
- Recognition of the vein pattern through palpation
- Insertion of a needle and cannula
- Management of blood flow

### Three Vein Pad - Venepuncture





This soft tissue strap-on pad features 3 straight veins, one of which is paediatric size.

#### SKILLS

- Professional-to-patient communication
- Insertion of a needle
- Insertion of cannula
- Management of blood flow

# QUALITIES

#### **KEY FEATURES**

- Ideal for professional-to-patient communication and scenario based training using Simulated Patient
- Strap onto the antecubital fossa of the arm
- Veins are: self-sealing for repeated use, rechargeable via a one-way valve, supplied with mock blood under pressure and replaceable

#### VERSATILITY

- The pad can be connected to the Mock Blood Giving Set (60651) to simulate blood flow
- Only use Mock Blood Venous (250ml) (00223) to refill these veins

#### SAFETY

• Tough backing plate to prevent needlestick injuries

Extended ACF Pad - Venepuncture

This soft tissue pad is supplied on a base for bench top training in venepuncture and insertion of cannulae. It represents the antecubital fossa of the right arm.

• Epidermis and veins contain latex

#### **CLEANING & CARE**

- Replaceable epidermis
- Washable using soap and water


#### WHAT'S INCLUDED

00180	Vein for ACF Pad
00170	Epidermis for ACF Pad
-	Pad base with strap

#### WHAT'S INCLUDED

00213 Vein for Three Vein Pad

**00200** Epidermis for Three Vein Pad Pad base with strap



Improved, unidirectional skin texture.

\$ 20

CE CEOT

# Cutdown Pad



A multi-layered pad containing filled veins for practising venous cutdown procedure.

#### SKILLS

- Palpation to locate vein position
- Incising tissue layers
- Inserting a cannula
- Vein injection both distally and proximally
- Tying off veins

#### QUALITIES

#### ANATOMY

- Improved, unidirectional skin texture
- Layers simulate the epidermis and subdermal layers
- 3 straight veins (each 150mm long) filled with mock blood

#### **KEY FEATURES**

 The pad can be connected to the Mock Blood Giving Set (60651) to simulate blood flow. Only use Mock Blood -Venous (250ml) (00223)

#### VERSATILITY

• The veins can be used by multiple trainees if tied off after each use

#### SAFETY

• Contains latex

#### WHAT'S INCLUDED

#### 00280 Cutdown Pad

## Mock Blood - Venous

#### 00223

250ml of synthetic solution made to the correct colour and viscosity of venous blood.



# Mock Blood Giving Set



A 1500ml blood bag with luer lock connector tube and stand.

This product can be used with Limbs & Things' products requiring a mock blood supply.

#### QUALITIES

#### KEY FEATURES

- Quick and easy to assemble
- Easy storage collapsible stand
- Bag holds up to 1500ml
- Fluid flow shut-off clamp

#### VERSATILITY

- This product can be used with Limbs & Things products requiring a mock blood supply including:
  - Cutdown Pad (**00280**)
  - ACF Pad (00140) and Extended ACF Pad (00211)
  - Three Vein Pad (00161)
  - Vermiform Appendix (50122, 50123 and 50124)
  - Gall Bladders (50127, 50128, 50129 and 50132)

#### WHAT'S INCLUDED

60652	Mock	Blood/Fluid	Bag

- Luer lock tube adapter
- Metal stand with base

# Mock Blood - Arterial

#### 60653

250ml of synthetic solution made to the correct colour and viscosity of arterial blood.



# Injection Trainer

#### 00310



Our soft tissue injection pad is designed for practising intradermal, subcutaneous and intramuscular tissue injection techniques.

The Injection Trainer has multiple tissue layers representing the epidermis, dermis, fat and muscle layer, and can easily attach to an arm or thigh to help teach professional-to-patient communication.

#### SKILLS

- Subcutaneous injection
- Intradermal and intramuscular injection
- Management of tissue
- Professional-to-patient communication

#### QUALITIES

#### ANATOMY

- Epidermis, dermis, fat and muscle
- REALISM

liquid

#### SAFETY

#### **CLEANING & CARE**

• Durable replaceable epidermis

#### WHAT'S INCLUDED

00311	Skin	Pad	&	Muscle	Block

**00320** Epidermis for Injection Trainer (x2) Strap-on base plate

Soft tissue tray with liner

- Tissues are soft and warm to the touch

#### **KEY FEATURES**

• Epidermis layer peels back to quickly release subcuticular

#### VERSATILITY

• Strap for hybrid simulation

• Latex free

## Ultrasound Vascular Access Trainer

#### 60530



An echolucent and contoured slab with representation of the jugular vein for teaching needle insertion under ultrasound.

#### SKILLS

- Use of ultrasound technology
- Handling ultrasound transducers
- Insertion of needles under ultrasound guidance
- Handling and management of guidelines and catheters under ultrasound
- Recognition of vessels using ultrasound

#### QUALITIES

#### REALISM

• Echolucent material provides realistic ultrasound image

#### **KEY FEATURES**

- Vein is self-sealing for repeated use and rechargeable with water via a one-way valve
- Designed to be used on the bench top

#### VERSATILITY

• Fluid volume and pressure can be changed in the vessel to represent variability in patient circulation, providing various levels of difficulty for needle insertion

#### SAFETY

• Veins contain latex

#### WHAT'S INCLUDED

60531	Ultrasound Vascular Access Slab
60532	Ultrasound Vascular Access Vein
-	Base
-	10ml syringe

# Paracentesis Trainer

#### 60100



This trainer fulfills the fundamental component of core medical training in diagnostic and therapeutic techniques of Paracentesis. Both landmark and ultrasound guided techniques can be practised.

Trainees are able to identify the echogenic anatomy enabling them to learn how to safely insert a needle or catheter into the peritoneal cavity.

#### QUALITIES

#### ANATOMY

• Torso featuring bony landmarks (pubis symphysis, iliac crest) and umbilicus

- Bowel

- Floating Bowel

- Internal anatomy includes:
  - Liver
- Spleen
- REALISM
- Realistic tissue and needle response

#### **KEY FEATURES**

- Palpation or ultrasound techniques can be practiced (side by side)
- Internal echogenic anatomy to allow recognition of landmarks under ultrasound
- Two 3.5 litre chambers can be filled with water for practise drainage
- Cost effective consumable pads

#### VERSATILITY

- Self-sealing pads, withstanding up to 200 needle or up to 100 rocket catheter insertions
- Ability to insert and remove a drain
- Positions both supine and lying on side

#### SAFETY

• Latex free

#### **CLEANING & CARE**

• Skin surface is washable using soap and water









- Trainer with skin removed, showing internal anatomy
- Catheter insertion for therapeutic drainage
- Needle insertion for diagnostic purposes
- Ultrasound guided technique

#### SKILLS

- Familiarity with the abdominal regions and underlying anatomy
- Palpation of anatomical landmarks
- Identification of excess fluid
- Through ultrasound, identification of safe needle/catheter insertion sites
- Insertion of needle into peritoneal cavity for therapeutic or diagnostic purposes
- · Professional-to-patient communication

#### WHAT'S INCLUDED

60101	Paracentesis Skin
60101	Paracentesis skin
60102	Paracentesis Pads (x2)
60103	Paracentesis Liver
60104	Paracentesis Spleen
60105	Paracentesis Bowel Insert
60106	Paracentesis Floating Bowel
60107	Paracentesis Clamping Frame
60108	Paracentesis Large Volume Refill System
-	Base unit
-	Carry case

# CATHETERISATION TRAINER RANGE



# Male Catheterisation Trainer

60850



The male trainer has soft and compliant genitalia to enhance the realistic experience of catheter insertion.

Replaceable and removable foreskin.

The Catheterisation Range facilitates the teaching of urethral and suprapubic catheterisation across all levels of training.

Interchangeable modules allow practice of the fundamental procedures of male and female urethral catheterisation, as well as more advance skills of suprapubic management and ultrasound guided insertion.

The anatomically accurate trainers include a supple urethra and resistant sphincter providing a closer to life training experience.



# Female Catheterisation Trainer

60851



The female trainer has soft labia allowing trainees to practice parting using the aseptic no touch technique, as well as management and care of a preinserted suprapubic catheter via the pre-ported bung.



# Suprapubic Insertion Unit

#### 60852



Optional module, designed for advanced skills training, allowing trainees to master ultrasound guided suprapubic catheterisation insertion.

Fits into **60850** & **60851**.



NEW

#### SKILLS

- Correct handling of male and female anatomy
- Aseptic catheterisation technique
- Withdrawal of catheter
- Urinary catheterisation
- Suprapubic catheter management Female module
- Suprapubic catheter insertion (blind and ultrasound guided)\*
- Ultrasound bladder scanning on a full bladder\*
- Fluid management

#### \* 60852 Only









- Aseptic catheterisation technique
- Ability to view catheter path when module used without pelvic shell
- Suprapubic management
- Suprapubic catheter insertion

## QUALITIES

#### REALISM

- Supple urethra & resistant sphincter providing realistic response
- Flaccid penis and removable foreskin
- Vaginal opening and realistic urethral meatus
- Anatomically accurate bladder

#### **KEY FEATURES**

- Ability to assess the student's progress
- Ultrasound compatible\*
- Easy to set-up
- Non-drip valve
- Modular

#### VERSATILITY

- Easily interchangeable male, female and suprapubic catheterisation modules
- Removable foreskin

#### SAFETY


• Trainer is latex free (catheters provided contain latex)

# Standard Catheterisation Trainer Set

#### 60853



# Advanced Catheterisation Trainer Set



WHAT'S INCLUDED		60850	60851	60852	60853	60854
60855	Male Catheterisation Module	$\checkmark$	×	×	$\checkmark$	~
60856	Female Catheterisation Module	×	$\checkmark$	×	$\checkmark$	$\checkmark$
60857	Suprapubic (Ultrasound) Catheterisation Module	×	×	$\checkmark$	×	×
60858	Suprapubic Insertion Bungs (x2)	×	×	$\checkmark$	×	×
60859	Catheterisation Foreskins (x3)	$\checkmark$	×	×	$\checkmark$	$\checkmark$
60862	Aseptic Catheterisation Sleeves (x2)	$\checkmark$	$\checkmark$	×	$\checkmark$	✓(x2)
60863	Water System	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	✓(x2)
60864	Pelvic Shell & Stand	$\checkmark$	$\checkmark$	×	$\checkmark$	✓(x2)
60867	Water Extension Tubing	$\checkmark$	$\checkmark$	×	$\checkmark$	✓(x2)
50183	Mock LA Gel with Syringe	$\checkmark$	$\checkmark$	×	$\checkmark$	✓(x2)
50185	Antifungal Agent (Milton Sterilising Fluid)	$\checkmark$	$\checkmark$	$\checkmark$	$\checkmark$	✓(x2)
-	14 French Foley Catheter	$\checkmark$	$\checkmark$	×	$\checkmark$	✓(x2)
-	Plastic Kidney Dish	$\checkmark$	$\checkmark$	×	$\checkmark$	✓(x2)

#### OPTIONAL EXTRAS


# Chest Drain & Needle Decompression Trainer

#### 60230

Our Chest Drain and Needle Decompression Trainer has been designed to meet the specific requirements of healthcare professionals training in surgical or guidewireassisted thoracostomy and thoracentesis.

This product allows for a variety of chest drain insertion techniques to be performed including ultrasound-guided techniques.

Ideal for ATLS V10 courses



#### QUALITIES

#### ANATOMY

- Representation of adult male thorax with arms raised
- Bony and soft tissue landmarks: manubriosternal joint, clavicles, ribs, pectoralis major and latissimus dorsi
- Internal ultrasound anatomy: diaphragmatic structures and collapsed lung

#### REALISM

- Can give the impression of breathing under ultrasound when using the advanced pad
- Needle decompression air reservoirs provide realistic release of air on insertion of needle

#### **KEY FEATURES**

- Reservoirs can be filled with fluid or mock blood to represent pleural effusion
- Bilateral chest drain and needle decompression pads

#### VERSATILITY

- Suitable for supine, sitting or leaning forwards positions
- Works with thoracic seals when using the standard pad
- Affordable replaceable pads


#### SAFETY

• Latex free

#### SKILLS

- Needle decompression of a tension pneumothorax (at both the 2nd and 5th intercostal space)
- Open, or cut-down chest drain insertion: recognition of correct position, surgical incision, blunt dissection through chest wall, perforation of pleura and finger sweep
- Suture of tube to chest wall
- Ultrasound-guided chest drain insertion (Seldinger-type), including insertion of needle under direct vision and ultrasonic recognition of chest structures
- Management of pleural effusion





- Blunt dissection
- Chest tube insertion
- \_ Seldinger technique
- Internal anatomy visible under ultrasound Realistic needle decompression \_
- \_

## PAD FEATURE COMPARISON

ADVANCED STANDARD

Ultrasoundable	$\checkmark$	×
<ul> <li>Use with liquids</li> </ul>	$\checkmark$	×
Effusion	$\checkmark$	×
<ul> <li>Haemothorax simulation</li> </ul>	$\checkmark$	×
Pleural layer	$\checkmark$	$\checkmark$
<ul> <li>Blunt dissection</li> </ul>	$\checkmark$	$\checkmark$
<ul> <li>Seldinger technique</li> </ul>	$\checkmark$	×
<ul> <li>Can be sutured</li> </ul>	$\checkmark$	$\checkmark$
<ul> <li>Needle insertion</li> </ul>	$\checkmark$	$\checkmark$
<ul> <li>Guidewire insertion</li> </ul>	$\checkmark$	$\checkmark$
<ul> <li>Dilator insertion</li> </ul>	$\checkmark$	$\checkmark$
<ul> <li>Drain tube insertion</li> </ul>	$\checkmark$	$\checkmark$
<ul> <li>Use with adhesive dressings</li> </ul>	$\checkmark$	$\checkmark$
<ul> <li>Shows respiratory swing</li> </ul>	$\checkmark$	×
Latex free	$\checkmark$	$\checkmark$

#### WHAT'S INCLUDED

60231	Needle Decompression Pads (x2)
60232	Advanced Chest Drain Pads (x2)
60234	Standard Chest Drain Pads (x2)
60235	Chest Drain Rib Frame (Left)
60237	Chest Drain Rib Frame (Right)
-	Chest Drain Lungs (x2)
-	Chest Drain Diaphragms (x2)
-	Chest Drain Ribs (x12)
-	Chest Drain Pump: Respiratory Swing
-	Chest Drain Pump: Needle Decompression
-	Chest Drain Slotted Retaining Knobs
-	Key with keyring
-	Torso
-	Carry case

# LUMBAR PUNCTURE & EPIDURAL RANGE



This range offers realistic, repeatable and cost effective training in the core skills of Lumbar Puncture and Epidural Anaesthesia for both palpation and ultrasound guided techniques.

All three models share a common base but offer different levels of complexity to support trainees from early through to advance learning.

Designed in conjunction with clinicians, the models feature realistic feel and response while offering a safe environment for learning and skill maintenance.

# Lumbar Puncture Model

#### 61000

This model is designed for early years training. It features replaceable tissue inserts, durable spine and dura that accurately recreate the feel of puncturing the dura and the collection and measurement of CSF.

Its separate spine and replaceable tissue means the model is a cost-effective way to train large numbers of students.



### Advanced Epidural & Lumbar Puncture Model

#### 61001

The Advanced model allows training in both lumbar puncture and epidural administration. Featuring two epidural inserts: standard and an advanced version with steeper spine, the model is ideal for both early and advanced training for more complex procedures.

This model also features optional fat pads to fit below the skin to increase realism and difficulty for trainees.

Ultrasound Epidural & Lumbar Puncture Model

#### 61002

Ideal for those wishing to teach lumbar puncture or epidurals using the ultrasound guided technique.

Tested to work on all major ultrasound brands, this model includes visible facet joints.





#### QUALITIES

#### ANATOMY

• Palpable sacrum, iliac crests and vertebrae L2-5

#### REALISM

- Can feel when the dura has been punctured
- Able to practice loss of resistance technique\*
- Obese skin options\*
- Optional steeper spine and smaller vertebral spaces\*

#### **KEY FEATURES**

• Position in left lateral and sitting

#### VERSATILITY

- All inserts compatible with common base
- Rotatable skin to ensure cost-effective training

\* 61001 Only








- Epidural inserts
- Realistic feel and response
- Lumbar puncture training


SKILLS	61000	61001	61002
Injection of local anaesthetics	$\checkmark$	$\checkmark$	$\checkmark$
Lumbar puncture	$\checkmark$	$\checkmark$	$\checkmark$
<ul> <li>Collect and measure CSF</li> </ul>	$\checkmark$	$\checkmark$	$\checkmark$
<ul> <li>Epidural administration</li> </ul>	×	$\checkmark$	$\checkmark$
<ul> <li>Ultrasound-guided technique</li> </ul>	×	×	$\checkmark$
Obese option	×	$\checkmark$	×

WHAT'S I	NCLUDED	61000	61001	61002
-	Epidural & Lumbar Puncture Torso with Torso Skin	$\checkmark$	$\checkmark$	$\checkmark$
-	Epidural & Lumbar Puncture Stand	$\checkmark$	$\checkmark$	$\checkmark$
61010	Lumbar Puncture Insert - Bone	$\checkmark$	$\checkmark$	×
61011	Lumbar Puncture Insert - Tissue (x2)	$\checkmark$	$\checkmark$	×
61012	Epidural Insert (x2)	×	$\checkmark$	×
61013	Advanced Epidural Insert (x2)	×	$\checkmark$	×
61014	Epidural & Lumbar Puncture Ultrasound Insert (x2)	×	×	$\checkmark$
61015	Epidural & Lumbar Puncture Skin Pad	$\checkmark$	$\checkmark$	×
61016	Epidural & Lumbar Puncture Fat Pad (x3)	×	$\checkmark$	×
61017	Epidural & Lumbar Puncture Standard Dura (x2)	$\checkmark$	$\checkmark$	$\checkmark$
61018	Epidural & Lumbar Puncture Advanced Dura (x4)	×	$\checkmark$	$\checkmark$
61019	Epidural & Lumbar Puncture Iliac Crest Insert	$\checkmark$	$\checkmark$	$\checkmark$
61020	Epidural & Lumbar Puncture Insert Cartridge	$\checkmark$	$\checkmark$	$\checkmark$
61021	Epidural & Lumbar Puncture Sacral Insert	$\checkmark$	$\checkmark$	×
61022	Water Reservoir	$\checkmark$	$\checkmark$	$\checkmark$
-	Marker pen	×	×	$\checkmark$

# LIMBS & THINGS WEBSITE



# limbsandthings.com



# BSS Day 1 Kit



Basic surgical skills kit facilitating the practice of a range of knot-tying, suturing and tissue handling techniques.

#### SKILLS

- Knot-tying: one-handed reef knot, instrument tie, surgeon's knot, tying at depth
- Suturing techniques: holding/manipulation of needles, interrupted, simple and mattress, continuous, subcuticular
- Skin lesions and LA techniques: excising a skin lesion, excising a sebaceous cyst
- Haemostasis: clip tie, continuity tie, pedicle transfixion
- Tissue handling bowel: end-to-end interrupted sutures
- Fine tissue handling: tendon repair

#### WHAT'S INCLUDED

50061	Double Layer Bowel 30mm
50074	BSS Lesion & Cyst
50051	BSS Surgical Skills Board
50056	BSS Tendon Repair Trainer
-	BSS knot tying trainer
	Skin pad jig with skin pad
-	BSS soft tissue retaining set
	Haemostasis pad

# BSS Day 1 Refill Pack



#### WHAT'S INCLUDED

- Skin pad
- Bowel
- Haemostasis pad
- BSS lesion & cyst
- Tendon

## Knot Tying Trainer



A comprehensive trainer for teaching all surgical knot tying techniques, including those specified in the GB and Ireland Intercollegiate Basic Surgical Skills Course.

#### SKILLS

- One-handed reef knot technique
- Instrument tie
- Surgeon's knot
- Slip knot
- Tying in a small opening
- Tying at depth vertically in a large opening
- Tying at depth, at an angle, in a large opening

#### QUALITIES

#### ANATOMY

- 2 perioperative openings represented by:
  - Small, shallow fixed cylinder for tying in a small opening
  - Large, deep removable cylinder, reversible for angled abdominal and gynaecological depth tying

#### REALISM

- Unique magnetic system to represent tissue strength
- Parallel knotting tubes are elastic for a realistic tissue response

#### **KEY FEATURES**

- Light and compact
- Cylinders are transparent to allow the trainer to observe and assess trainee competence

#### WHAT'S INCLUDED

50164	Two Coloured Cord (x3)
50158	Knotting Tube (x3)
	Base with cylinders
-	Knot tying tension hooks (x3)

50072

# Double Layer Bowel

# 50061 50065 30mm OUTSIDE DIAMETER

# Skin Pad Jig Mk 3

#### 00550



Realistic 2 layer bowel simulation for training in anastomosis techniques. Both models are approximately 200mm in length.

#### SKILLS

- Tissue handling
- Single layer extra mucosal suturing
- Continuous suturing technique
- Stapling
- End-to-end extra mucosal anastomosis
- End-to-side anastomosis

#### QUALITIES

#### REALISM

- Realistic tissue response
- Will withstand a fluid flush test to demonstrate integrity of the anastomosis

#### VERSATILITY

- Light and compact
- Mobile when held in Soft Tissue Retaining Set (50151)

#### SAFETY

- Latex free
- Non-biological material: no hygiene issues or problems

# Skin Pad Jig Strap (x3)

#### 00510

This low cost addition for the Skin Pad Jig Mk 3 (00550) allows the jig to be strapped onto the limbs of a Simulated Patient.

This jig supports the entire range of Limbs & Things synthetic soft tissue pads.

Can be strapped onto the limbs or body of a Simulated Patient when used with the optional strap (**00510**).

#### QUALITIES

#### REALISM

- Presents pads on a curved lifelike profile allowing incisions to gape, akin to real life
- Flexible plastic allows for realistic 'give' of skin pad when pressure is applied

#### **KEY FEATURES**

- Sucker feet work effectively on smooth work surfaces
- Stable and non-slip

#### VERSATILITY

- Optional jig strap allows use with Simulated Patient
- Use with:
  - Professional Skin Pad Mk 2
- Wound Closure Pad
- Lipoma Pad
- Sebaceous Cyst Pad, Curettage, Snip & Shave Excision Pad
- Minor Skin Procedures Pads





# Wound Closure Pad



A low cost, 1 layer pad, ideal for demonstrating, training and practising core skills associated with wound closure.

#### SKILLS

- Making an incision
- Application of skin adhesives
- Interrupted suturing and knotting
- Stapling

An advanced 3 layer skin pad for demonstrating and practising a variety of incisions and a wide range of suturing techniques.

#### SKILLS

- Incisions: linear, ellipse, flaps, shaped
- Subcuticular undermining
- Simple and advanced interrupted suturing techniques
- Subcuticular suturing

00094

00093

SMALL (x4) LARGE (x2)

- Continuous suturing
- Use of adhesive strips



#### QUALITIES

#### REALISM

• Works with all major topical skin adhesives

#### **KEY FEATURES**


- An excellent starter pad for wound closure training
- Pad dimensions: large 145mm x 125mm. small 125mm x 72mm

#### VERSATILITY

• Use with the Skin Pad Jig Mk 3 (00550)

#### SAFETY

• Latex free



#### QUALITIES

#### REALISM

- Improved, unidirectional skin texture
- Realistic tissue response with a similar drag and strength to human skin
- All layers have realistic retention of sutures

#### **KEY FEATURES**

• Pad dimensions: large 145mm x 125mm, small 125mm x 72mm

#### VERSATILITY

- Serves the needs of both specialist and beginner
- Use with the Skin Pad Jig Mk 3 (00550)

#### SAFETY

• Latex free

# ..... IMPROVED


Professional Skin Pad Mk 2

00092



## Lipoma Pad


### 00081



For practising the full lipoma removal procedure. 3 layer skin: epidermis, dermis and subdermal layer. Pad contains 2 lipomas.

#### SKILLS

- Designing the incision
- Marking out the lipoma
- Incising the skin
- Blunt dissection
- Removal
- Closure



For practising the full sebaceous cyst removal procedure. A 3 layer skin pad with epidermis, dermis and subdermal layer. Contains 2 cysts.

#### SKILLS

- Marking out the cyst
- Planning and marking the ellipse
- Incising the skin
- Sharp dissection
- Blunt dissection
- Management of a burst cyst
- Removal
- Closure





### QUALITIES

#### REALISM

- Improved, unidirectional skin texture
- Lipoma is consistent with that found in a patient
- Lipoma sited within the fat layer
- Precise replication of in-vivo adhesion
- All layers of tissue can be identified

#### **KEY FEATURES**

• Pad dimensions: 125mm x 72mm

#### VERSATILITY

• Use with the Skin Pad Jig Mk 3 (00550)

#### SAFETY

• Latex free





### QUALITIES

#### ΑΝΑΤΟΜΥ

- Improved, unidirectional skin texture
- 3 layer skin: epidermis, dermis and subdermal layer

## Contains 2 Cysts

- REALISM
- Cyst is consistent with that found in a patient
- Cyst exhibits appropriate adhesion to other surrounding tissue
- Cyst sited immediately below skin layer and is attached firmly to its underside
- Fragile epidermoid wall filled with simulated pus

#### KEY FEATURES

• Pad dimensions: 125mm x 72mm

#### VERSATILITY

• Use with the Skin Pad Jig Mk 3 (00550)

#### SAFETY

• Latex free

## PROMPT FLEX RANGE




The **PROMPT Flex Range** is an integral part of many training courses around the world, suitable for all levels of training and skill maintenance.

Now also available in a dark skin tone, this range offers the flexibility to train on routine and a number of difficult birthing management scenarios.



80200	PROMPT Flex - Standard - Dark	
80201	PROMPT Flex - PPH Module - Dark	NE
80202	PROMPT Flex - Cervical Dilatation & Effacement Module - Dark	
80203	PROMPT Flex - Caesarean Section Module - Dark	
80206	PROMPT Flex - Advanced - Dark	
80208	PROMPT Flex - Advanced - with Tablet - Dark	
80212	PROMPT Flex - PostPartum Haemorrhage Kit - Dark	
80221	PROMPT Flex - Standard Baby - Dark	
80280	PROMPT Flex - Cervical Cerclage Module - Dark	



### PROMPT Flex - Standard

80100 80200

A modular design birthing simulator allowing for numerous training scenarios.

Used either for stand-alone training or integrated into hybrid simulation, it is an ideal training solution for all skills relating to routine and difficult deliveries.


## **PROMPT Flex - Advanced**

80106 80206

Incorporates Wireless Force Monitoring where force applied to the baby's head is measured during a shoulder dystocia drill. Includes an app with inbuilt scenarios allowing monitoring and recording of force, interventions made and time taken to deliver the baby.



Tablet not included Download app at <u>limbsandthings.com</u>



80108 PROMPT Flex Advanced - with Tablet 80208 PROMPT Flex Advanced - with Tablet - Dark

#### QUALITIES

#### ΑΝΑΤΟΜΥ

- Birth canal and cervix
- Ischial spines and pubic bone
- Gynaecoid pelvis
- Articulating thighs
- IM injecting
- Fully articulated baby with detachable umbilical cord and placenta
- Baby includes clavicles, fontanelles and flexible head

#### **KEY FEATURES**

- Improved and more realistic pelvic floor
- Articulating thighs for McRobert's procedure
- Stretchable perineum
- Soft flexible birthing canal
- Additional modules to extend learning
- Optional lower legs for extended realism
- Optional wireless force monitoring for training in shoulder dystocia management

#### SAFETY

- Latex free
- CLEANING & CARE
- Skin is washable using soap and water





S	ΚI	LL	S	

#### 80100 80106

• Communication & team work skills (SP)	$\checkmark$	$\checkmark$
• Normal	$\checkmark$	$\checkmark$
• Breech	$\checkmark$	$\checkmark$
<ul> <li>Shoulder dystocia management</li> </ul>	$\checkmark$	$\checkmark$
• Instrumental (forceps & vacuum devices)	$\checkmark$	$\checkmark$
<ul> <li>Delivery of placenta</li> </ul>	$\checkmark$	$\checkmark$
Cord prolapse	$\checkmark$	$\checkmark$
<ul> <li>Urinary catheter placement</li> </ul>	$\checkmark$	$\checkmark$
<ul> <li>IM injection pads</li> </ul>	$\checkmark$	$\checkmark$
<ul> <li>Force feedback of baby's head</li> </ul>	×	$\checkmark$
App with comparing & force magnituring	~	/

App with scenarios & force monitoring ×

#### WHAT'S INCLUDED

80100 80106

80120	Birthing Mother (+Upper Legs)	$\checkmark$	$\checkmark$
80121	Standard Baby	$\checkmark$	×
80122FM	Wireless Force Monitoring Baby	×	$\checkmark$
80123	Placenta	$\checkmark$	$\checkmark$
80124	Abdomen for PROMPT Flex	$\checkmark$	$\checkmark$
80125	Perineum and Birth Canal	$\checkmark$	$\checkmark$
80130	Bed Straps	$\checkmark$	$\checkmark$
10193	Lubricant	$\checkmark$	$\checkmark$
01000	Wheeled Carry Case	$\checkmark$	$\checkmark$

## PROMPT Flex Standard Baby

80121 80221

The PROMPT Flex Standard Baby has close to life range of articulation in the shoulder and hip joints, allowing an accurate experience with both delivering and handling the baby. The addition of flexion at the thoracolumbar region provides accurate training experience and enhances realism in vaginal deliveries and c-section manoeuvres.



### QUALITIES

#### REALISM

- Fully articulated baby with umbilical cord
- Full term infant size

#### **KEY FEATURES**

- Improved and more realistic limb articulation
- Correct anatomical landmarks for delivery
- Versatile use for wide range of delivery types
- Flexion allows for life-like manoeuvring


#### VERSATILITY

• Can be quickly and easily used in conjunction with the Birthing Mother and C-Section module

#### SAFETY

• Latex free







## PROMPT Flex - Cervical Dilatation & Effacement Module

80102 80202

Allows for training in assessment of both the latent and active first stages of labour.

This versatile module is an optional add-on to the PROMPT Flex Standard or Advanced.















Fetal station adjustment (-3 to +3)

Cervical Dilation (1-10cm)

#### QUALITIES

#### REALISM

• Realistic representation of cervixes, including anterior lip, and presenting parts in soft birth canal, with palpable ischial spines

#### **KEY FEATURES**

- Inserts to represent early labour cervixes effacement, dilation and ripeness in line with Bishop's scoring
- Positioning mechanism allows adjustment of station, dilation and tilt without removal from the simulator
- Dynamic positioning mechanism allows adjustment of dilation in active labour
- Markers allow tutor to read positioning in situ
- Numerous presenting part inserts including
  - flexed
  - deflexed
  - brow
  - face
  - breech
  - caput and moulding

#### VERSATILITY

• Presenting parts and cervix can be set up outside the model for demonstration, then placed inside for examination

#### SAFETY

• Latex free

#### SKILLS

- Assessment and Bishop's scoring of:
  - Cervical dilation (1-10cm)
  - Cervical effacement (0-100%)
  - Cervical ripeness/consistency (soft, medium, hard)
  - Cervical position (anterior, mid, posterior)
  - Fetal station (-3 to +3)
- Assessment of and artificial rupture of membranes
- Assessment of presenting part flexed, deflexed, brow, face, breech, caput and moulding and caput formation

#### WHAT'S INCLUDED

80138	Perineum
80139	Pelvic Ring
80140	Static Rig
80141	Dynamic Rig
80142	Set of 7 Latent Labour Cervixes
80143	Dynamic Cervix
80144	Set of 5 Presenting Parts
80145	Set of 2 Difficult Labour Heads
80146	Amniotic Membrane Set
80162	Set of 3 Caputs (x2)

## PROMPT Flex - Caesarean Section Module

80103 80203

Ideal for practising the necessary skills required when performing a caesarean section. Trainees can have exposure to delivering a difficult baby such as a deeply impacted head, or a transverse lie.

It includes 2 surgical pads for incision and suturing, along with a pre-incised abdomen for repeated delivery practice to help improve confidence and familiarity of manoeuvres.



#### QUALITIES

#### REALISM

- Abdominal skin has realistic appearance and soft feel of atterm abdomen
- Realistic feel of uterus interior during delivery
- Head presents at realistic level once incision is made

#### **KEY FEATURES**

- Uterus holds baby firmly in position for normal and transverse lies, breech and cephalic presentations
- Uterus support positions baby and directs it towards the pelvis/abdominal incision when fundal pressure is applied
- Pre-incised skin option allows easy repeat delivery practice
- Surgical pad allows realistic dissection and closure and is reversible to allow 2 uses
- Single/double in situ uterine wall closure
- Suturing of rectus sheath, fat and subcuticular or cutaneous skin closure
- Accommodates retraction using metal, plastic or rubber retractors

#### VERSATILITY

• Optional insert to allow presentation of head just below ischial spines for caesarean section at full dilation

#### SAFETY

• Latex free



The module can be used to practise any of the following

methods of delivering a deeply impacted fetal head:

• Assist from a vaginal hand (push technique)

• Reverse breech extraction (pull technique)

• Routine manoeuvres

• Use of a fetal pillow

#### SKILLS

- Abdominal and uterine wall opening and closure
- Cephalic and breech delivery
- C-section at full dilation:
  - Head disimpaction using vaginal push methods
  - Use of balloon disimpacting device
  - Reverse breech extraction
- Transverse lie
- Instrument assisted C-section delivery
- Delivery of the placenta
- Teamwork and human factors training in emergency C-section

#### WHAT'S INCLUDED

80104	Caesarean Section Surgical Pad (x2)
80138	Perineum
80170	Uterus for Caesarean Section Module
80171	Pre-incised Delivery Only Skin
80172	Surgical Skin
80173	Cervix for Caesarean Section Module
80175	Foam Uterus Support
80176	Full Dilatation Pelvic Insert

## PROMPT Flex - PPH Module

#### 80101 80201

Allowing training in the management of postpartum bleeding, it offers a realistic simulation of an atonic and contracting uterus.

This module is an optional add-on for PROMPT Flex Standard or Advanced.

#### SKILLS

- Management of postpartum bleeding:
  - Fundal massage
  - Bi-manual compression
  - Haemostatic balloon insertion and management
- Communication skills with mother when used with Simulated Patient
- Intramuscular injection
- Estimating blood loss

#### WHAT'S INCLUDED

80127 PROMPT Flex PPH Uterus Assembly

**00020** Concentrated Venous Blood - Starter Pack (x2)

PROMPT Flex PPH blood container

#### QUALITIES

#### REALISM

• Atonic uterus can simulate contraction on fundal massage

#### **KEY FEATURES**

- Blood loss of up to 2 litres can be simulated
- Air bulbs are used to stimulate blood flow and control atonic state of the uterus

#### VERSATILITY

• Can be quickly and easily fitted into the birthing mother to allow simulation of situations where bleeding occurs directly after delivery

## PROMPT Flex - Compression Suture Uterus

#### 80111

This module allows for the training in performing a compression suture and ligating the uterine artery.

Used in conjunction with the **80171** Pre-incised Delivery Skin, it provides the ideal training platform for repeat practise of uterine compression suturing.




This product is based on concepts and collaboration with:

- Dr Carol L. Simmons, MD SimTunes, LLC.
- Dr Mona Sharma, Consultant Obstetrician and Gynaecologist at Queen's Hospital, Romford BHR NHS Trust.

#### SKILLS

- All common types of Compression Suture procedures (including B-Lynch, Hayman and Pereira)
- Dissection of broad ligament
- Ligation of uterine arteries



#### QUALITIES

#### ANATOMY

• A representation of post-partum uterus including uterine arteries and broad ligament

#### **KEY FEATURES**

- Uterine arteries can be fluid filled
- Ligation of arteries is possible
- Multiple uses when conducting the compression suture

#### VERSATILITY

• Can be quickly and easily fitted into the birthing mother

## PROMPT Flex - PostPartum Haemorrhage Kit

#### 80112 80212

This kit consists of the PPH module, Compression Suture Uterus and a pre-incised skin. It offers everything required for the management of postpartum bleeding and care including performing a compression suture and ligating the uterine artery.

To be used with the PROMPT Flex Standard or Advanced.

#### SKILLS

- Management of postpartum bleeding:
  - Fundal massage
  - Bi-manual compression
  - Haemostatic balloon insertion and management
- Estimating blood loss
- All common types of Compression Suture procedures including:
  - B-Lynch
  - Hayman
  - Pereira
- Dissection of broad ligament
- Ligation of uterine arteries
- Professional-to-patient communication skills



### WHAT'S INCLUDED

80101	PPH Module
80111	Compression
80171	Pre-incised D

Compression Suture Uterus (Pack of 2) Pre-incised Delivery Skin for PROMPT Flex

- PROMPT Flex Cervical Cerclage Module
  80180 80280 New
  Ideal for practising the necessary skills required for performing cervical cerclage. Trainees can have exposure to performing both elective and emergency cervical cerclage procedures.
  SKILLS
  Management of cervical incompetence
  Elective cerclage
  McDonald cerclage technique
  QUALITIES
  REALISM
- Realistic representation of a dilating cervix with presenting amniotic sac for emergency procedures
- Restrictive space within the vagina to perform the procedure as in life

#### **KEY FEATURES**

- Detachable cervix to allow for repeated training in both elective and emergency cerclage procedures
- The detachable OS provides instant feedback to the trainee
- Fluid filled amniotic sac replicates bulging membranes for emergency procedures

#### WHAT'S INCLUDED

80181	Cervical Cerclage Frame
80182	Cervical Cerclage Vagina
80183	Cervical Cerclage Elective OS (x5)
80184	Cervical Cerclage Emergency OS (x5)
80185	Cervical Cerclage Amniotic Sacs (x5)
80186	Cervical Cerclage Fluid System

## PERINEAL REPAIR AND EPISIOTOMY TRAINER RANGE

This multifaceted trainer allows for the development and maintenance of the skills of episiotomy and perineal repair. A common stand accommodates a choice of two different perineal repair blocks, one anal sphincter repair block and one unincised block. The addition of the episiotomy module gives a cost-effective and realistic training platform for this fundamental skill.

## Midline Perineal Repair Trainer

#### 61100 61120

The module reflects the current teaching practise in the USA, supporting the learning of 2nd degree midline perineal repair.



QUALITIES (Common in all Perineal Repair Range)

#### ANATOMY

• Hymenal remnants, vaginal mucosa, perineal skin, perineal muscles, external anal sphincter, internal anal sphincter, anal mucosa and pararectal fat

#### REALISM

- All essential layers of tissue can accept repeated sutures
- Soft skin with a similar drag and strength to human skin

#### **KEY FEATURES**

- Cost-effective as all layers can be sutured repeatedly
- Durable, long-lasting mediolateral perineal repair block and anal sphincter repair block
- Illustration pack representing each degree of perineal trauma to aid identification of key landmarks and internal anatomy

#### VERSATILITY

- Can be used alongside PROMPT Flex by placing in front of it to add further realism to the simulation
- Quick and easy to set up
- Clamps provided to secure trainer to benchtop
- Compact and portable

SKILLS (Common in both Perineal Repair Trainers)

- 2nd degree perineal repair
  - Tissue layer recognition and handling
  - Suturing in two planes
  - Vaginal mucosal suturing
  - Deep muscle suturing
  - Subcuticular suturing
  - Identification and management of perineal tears
  - Enables digital rectal examination before and after repair

## Mediolateral Perineal Repair Trainer




The module reflects the current teaching practise in the UK, supporting the learning of 2nd degree mediolateral perineal repair.



## OPTIONAL EXTRA

For those looking to extend their training, the following perineal repair blocks are available, and also fit into the **Perineal Repair Stand** (61109).



### Unincised Perineal Repair Block (x2)

#### 108 61128

Ideal to create specific scenarios, this trainer can alter the severity or placement of the tear to match the teaching objectives.



## Anal Sphincter Repair Trainer

#### 61102 61122

Ideal for practising repair of 4th degree anal sphincter injuries (OASI), its realistic anatomy allows for repair of both the internal and external anal sphincter and anal mucosa.



#### SKILLS

- 4th degree repair
- Recognition of internal anatomy
- Tissue Handling
- End-to-end repair of the external anal sphincter
- Overlap repair of the external anal sphincter
- Repair of the internal anal sphincter and anal mucosa

#### WHAT'S INCLUDED

#### 61100 61101 61102

61105	2nd Degree Midline Perineal Repair Block (x2)	✓	×	×
61106	2nd Degree Mediolateral			
	Perineal Repair Block (x2)	×	$\checkmark$	×
61107	Anal Sphincter Repair Block (x2)	×	×	$\checkmark$
61109	Perineal Repair Stand with Clamps	$\checkmark$	$\checkmark$	$\checkmark$
61110	Perineal Repair Block loading frame	$\checkmark$	$\checkmark$	$\checkmark$
61112	Perineal Repair Illustration Pack - Light	$\checkmark$	$\checkmark$	$\checkmark$
00293	Water Based Lubricant	$\checkmark$	$\checkmark$	$\checkmark$



## OPTIONAL EXTRA

This cost-effective, simplified model has been created to use in classroom settings to ensure the acquisition of confidence in the training of episiotomy. The crowning of the head gives trainees a true-to-life experience.

## Episiotomy Trainer Module



#### SKILLS

- Identification of fontanelles
- Handling of a stretched perineum
- Reducing the risk of damaging the baby or mother
- Infiltration of perineum prior to episiotomy
- Performance of mediolateral (or midline) episiotomy
- Use of instruments

#### QUALITIES

#### REALISM

• Realistic representation of tissue to represent a stretched perineum

#### **KEY FEATURES**

- Baby head has realistic look and feel, providing landmarks and appropriate tactile response when inserting fingers to guard the baby's head
- Repair tape increases the longevity of the episiotomy pad, allowing for multiple uses per pad
- Illustration pack to help trainees get familiarised with key landmarks and where to place the cut

#### WHAT'S INCLUDED\*

61113	Episiotomy Incision Pad (x2)
61114	Episiotomy Repair Tape (x24)
61115	Episiotomy Baby Head
61116	Episiotomy Illustration Pack - Light

\* Perineal Repair Stand (61109) to be purchased separately

## Surgical Female Pelvic Trainer (SFPT) Mk 2

#### 60251

A simulator for training in a wide range of laparoscopic gynaecological surgical techniques and procedures. The uterus component can be easily removed and replaced after use.



#### QUALITIES

#### ANATOMY

- Includes:
  - Uterus
  - Cervix
  - Bladder
  - Ovaries
  - Ureters
  - Ovarian artery

- Fibroid

- Dunn's pouch
- Ovarian cyst
- Ectopic pregnancy
- Uterosacral ligament
- Cardinal ligament
- Abdominal skin with realistic tissue properties for trocar insertion
- Pelvic floor with sympathetic features

#### REALISM

• Insufflation of abdomen is represented

#### **KEY FEATURES**

- Provides practice and role playing for gynaecological perioperative teams
- Includes a realistic vaginal canal and cervix for inspection, in which a uterine manipulator can be used
- Realistic bleeding when arteries are cut

#### VERSATILITY

• Compatible for use with harmonic scalpels

#### SAFETY

• Capsule in uterus contains latex

#### SKILLS

- All skills relevant to performing the following procedures:
  - Salpingectomy
  - Myomectomy
  - Hysterectomy
- Oophorectomy
- Dissection down to and location of ureters
- Insertion and use of uterine manipulator
- Working with various gynaecology instruments
- Identification of anatomical landmarks

#### WHAT'S INCLUDED

-	SFPT Mk 2 Abdomen
60252	SFPT Mk 2 Surgical Uterus
60253	SFPT Mk 2 Abdominal Skin Pad
60254	SFPT Mk 2 Perineum
60255	SFPT Mk 2 Bowel
00020	Venous Blood Starter Kit
-	SFPT Mock Blood Kit (includes reservoir and pump)
-	Blunt dispensing needle and syringe

#### CONSUMABLES

60256	SFPT	Uterus	(x5)
			( ( ( )

- 60257 SFPT Uterus (x10) 00021 Concentrated Venous Blood



- Salpingostomy


- Cystectomy

- Vagina - External introitus
  - Fallopian tube

  - Broad ligament
  - Uterine artery
  - Pouch of Douglas

## - Bowel and sigmoid



















## UNDERGRADUATE MEDICINE

### Designed to meet the core training needs and essential skills required at undergraduate level.

ABDOMINAL EXAMINATION

MALE PELVIC EXAMINATION

MALE RECTAL EXAMINATION



FEMALE PELVIC EXAMINATION



LUMBAR PUNCTURE



VENEPUNCTURE








BREAST EXAMINATION



CATHETERISATION



BASIC SURGICAL SKILLS



# NURSE TRAINING

Nursing is at the heart of healthcare, and these products reflect the core simulation training needs required for undergraduate and postgraduate nursing.





FEMALE PELVIC EXAMINATION



MALE PELVIC EXAMINATION



CATHETERISATION



ABDOMINAL EXAMINATION



BREAST EXAMINATION



VENEPUNCTURE



# PARAMEDIC TRAINING

From patient assessment and basic care in years 1 & 2, to more advanced skills in year 3, this range covers essential training elements required for training paramedics.



ABDOMINAL EXAMINATION

OBSTETRICS



CHEST TRAUMA MANAGEMENT



#### VENEPUNCTURE ARMS




#### INJECTIONS



# OSCE TRAINING

It's important students feel confident taking the Objective Structured Clinical Examination. These products are designed to train and test clinical skill performance for the OSCE.



LUMBAR PUCNTURE & EPIDURAL

ABDOMINAL EXAMINATION

CATHETERISATION



MALE RECTAL EXAMINATION



#### PELVIC EXAMINATION



#### BREAST EXAMINATION



VENEPUNCTURE



# PHYSICIAN ASSOCIATE

These products cover the broad set of skills required from Physician Associates including: physical exams, diagnosing illnesses, care coordination as well as a variety of procedures.



BREAST EXAMINATION



ABDOMINAL EXAMINATION

MALE RECTAL EXAMINATION



CATHETERISATION



#### VENEPUNCTUR



MALE PELVIC EXAMINATION



FEMALE PELVIC EXAMINATION



## OBSTETRICS & MIDWIFERY

The PROMPT Flex Range has become an integral part of many birthing training courses around the world, suitable for all levels of teaching, learning and skill maintenance.



C-SECTION

COMPRESSION SUTURE UTERUS





CERVICAL DILATATION & EFFACEMENT



PERINEAL REPAIR TRAINER



PPH











LIMBS & THINGS PRODUCTS DISTRIBUTED BY:



#### LAERDAL GERMANY

Laerdal Medical GmbH 89 86 49 54 0 info@laerdal.de

#### LAERDAL SPAIN

Laerdal España 902 29 11 10 laerdal.spain@laerdal.no

#### LAERDAL POLAND

Laerdal Medical Poland 32 493 70 20 Imp@laerdal.com

#### LAERDAL FRANCE

Laerdal Médical France 4 72 52 02 52 info@laerdal.fr

#### LAERDAL KOREA

Laerdal Medical Korea 2 6714 8800 kcsi@laerdal.com

#### LAERDAL BELGIUM

Laerdal Benelux nv/sa 02 253 36 96 info@laerdal.be

#### LAERDAL NORWAY

Laerdal Medical Norway 51 51 76 66 kundeservice@laerdal.no

#### LAERDAL BRAZIL

Laerdal do Brasil 11 4193 8007 comercial@laerdal.com

#### LAERDAL CHINA

Laerdal China Ltd 3168 2518 sales@laerdal.com.hk

#### LAERDAL ITALY

Laerdal Italia Srl 051 355587 laerdal.italy@laerdal.com

#### LAERDAL INDIA

Laerdal Medical India Pvt. Ltd 44 4261 4773 customer.service@laerdal.co.in

#### LAERDAL MALAYSIA

Laerdal Malaysia Sdn Bhd 603 7782 7002 mycustomerservice@laerdal.com

#### LAERDAL SINGAPORE

Laerdal Singapore Pte Ltd 6282 1912 sg.customer.service@laerdal.com



limbsandthings.com