

LAMPIRAN INAEQAS
KODE REAGENSIA KIMIA KLINIK 1 DAN 2
SIKLUS 1-2022

MEREK/PRODUSEN	NO KATALOG	KODE KIMIA 1	KODE KIMIA 2	METODE
Parameter : Protein Total 1 / 2				
Abbott	7D73	01A101	18A101	Colorimetric (Biuret) Cu-protein complex
ABX Pentra	A11A01669	01A201	18A201	Colorimetric (Biuret) Cu-protein complex
Advia Siemens	01803350	01A301	18A301	Colorimetric (Biuret) Cu-protein complex
AIM	C5T2E0400R	01A401	18A401	Colorimetric (Biuret) Cu-protein complex
AMP Diagnostic	BR 5202	01A501	18A501	Colorimetric (Biuret) Cu-protein complex
Agappe	51013002	01A601	18A601	Colorimetric (Biuret) Cu-protein complex
Bio Meriux	61602	01B201	18B201	Colorimetric (Biuret) Cu-protein complex
Biocon Diagnostics	9102 / 9104 / 9106	01B101	18B101	Colorimetric (Biuret) Cu-protein complex
Biolabo	LP807016	01B301	18B301	Biuret
Cypress Diagnostic	HB019	01C101	18C101	Colorimetric (Biuret) Cu-protein complex
Dialab	D00686 / D95680	01D301	18D301	Colorimetric (Biuret) Cu-protein complex
Dimension Siemens	DF73	01D101	18D101	Colorimetric (Biuret) Cu-protein complex
DSI/DiaSys	1 2311 99 10 021 / 1 2311 99 10 026	01D401	18D401	Colorimetric (Biuret) Cu-protein complex
DSI/Linear	1 2311 99 83 021 / 1 2311 99 83 026	01D501	18D501	Colorimetric (Biuret) Cu-protein complex
Dialine Diagnostic	2341052	01D601	18D601	Biuret with sample blank
Elitech	PRTB-0600	01E101	18E101	Colorimetric (Biuret) Cu-protein complex
Fast	2-221/2-240/2-236/2-237/2-295/4-236/4-436/4-536	01F101	18F101	Colorimetric (Biuret) Cu-protein complex
Human Diagnostic	157004/10570/10570150/10570300/10570600	01H101	18H101	Colorimetric (Biuret) Cu-protein complex
LABTEST	99-100 / 99-250	01L101	18L101	Colorimetric (Biuret) Cu-protein complex
Mindray	TP0102/TP0103/TP0105	01M201	18M201	Biuret Method, End Point
Olympus/ Beckman	OSR6132* / OSR6232*	01O101	18O101	Colorimetric (Biuret) Cu-protein complex
Ortho (Johnson & Johnson)	8392292	01J101	18J101	Colorimetric (Biuret) Cu-protein complex
Pacific / Indec Diagnostic	UA-K032	01P101	18P101	Colorimetric (Biuret) Cu-protein complex
Proline	1 2311 99 10 022/1 2311 99 10 191	01P201	18P201	Colorimetric (Biuret) Cu-protein complex
Randox	TP 245	01R101	18R101	Colorimetric (Biuret) Cu-protein complex
Roche	3183734190 / 4657586190 / 11553836316 / 11929917216	01R202	18R202	Colorimetric (Biuret) Cu-protein complex
Rajawali Nusindo	157092	01R301	18R301	Colorimetric (Biuret) Cu-protein complex
Sysmex	76396	01S201	18S201	Colorimetric (Biuret) Cu-protein complex
Sysmex	TP Reagent	01S401	18S401	Biuret Method, End Point
Sclavo	B74182707/B74182731/B74182750	01S601	18S601	Biuret
Thermo Scientific	981826/981827	01T101	18T101	Colorimetric (Biuret) Cu-protein complex
Lain-lain	-	01XX99	18XX99	Lain-lain
Parameter : Albumin 1/ 2				
Abbott	7D53	02A101	19A101	Colorimetric (BCG) ,Monoreagent
ABX Pentra	A11A01664	02A201	19A201	Colorimetric (BCG) ,Monoreagent
Advia Siemens	07622536	02A301	19A301	Colorimetric (BCG) ,Monoreagent
AIM	C5A2E0400R	02A401	19A401	Colorimetric (BCG) ,Monoreagent
AMP Diagnostic	BR 5002	02A501	19A501	Colorimetric (BCG) ,Monoreagent
AXIS SHIELD	1100296	02A601	19A601	Colorimetric (BCG) ,Monoreagent
Agappe	51001002	02A701	19A701	BCG + Modifikasi IFCC
Biocon Diagnostics	9132 / 9134 / 9136	02B101	19B101	Colorimetric (BCG) ,Monoreagent

MEREK/PRODUSEN	NO KATALOG	KODE KIMIA 1	KODE KIMIA 2	METODE
Bio Meriux	61051 / 61075 / 61142	02B201	19B201	Colorimetric (BCG) ,Monoreagent
Biolabo	80002	02B301	19B301	Bromcresol Green
Cypress Diagnostic	HB0010	02C101	19C101	Colorimetric (BCG) ,Monoreagent
Dimension Siemens	DF13	02D101	19D101	Colorimetric (BCP)
Dialab	D00205 / D97203	02D301	19D301	Colorimetric (BCG) ,Monoreagent
DSI/DiaSys	1 0220 99 10 021 / 1 0220 99 10 026	02D401	19D401	Colorimetric (BCG) ,Monoreagent
DSI/Linear	1 0220 99 83 021 / 1 0220 99 83 026	02D501	19D501	Colorimetric (BCG) ,Monoreagent
Dialine Diagnostic	2115022	02D601	19D601	BCG + Modifikasi IFCC
Elitech	ALBU-0600	02E101	19E101	Colorimetric (BCG) ,Monoreagent
Erba	XSYS001 XSYS001	02E201	19E201	Colorimetric (BCG) ,Monoreagent
Fast	2-213/2-239/2-238/2-294/4-238/4-438/4-538	02F101	19F101	Colorimetric (BCG) ,Monoreagent
Human Diagnostic	10560/156004/10560150/10560300/10560600	02H101	19H101	Bromcresol Green
Ortho (Johnson & Johnson)	1988211 / 8196057	02J101	19J101	Colorimetric (BCG) ,Monoreagent
LABTEST	19-250 / 19-500	02L101	19L101	Colorimetric (BCG) ,Monoreagent
Mindray	ALB0102/ALB0103/ALB0105	02M201	19M201	Bromcresol Green
Olympus/ Beckman	OSR6102 / OSR6202	02O101	19O101	Colorimetric (BCG) ,Monoreagent
Pacific / Indec Diagnostic	UA-K010	02P101	19P101	Colorimetric (BCG) ,Monoreagent
Proline	1 0220 99 10 022/1 0220 99 10 025/1 0220 99 10 029/1 0220 99 10 192	02P201	19P201	Bromcresol Green
Randox	AB 362	02R101	19R101	Colorimetric (BCG) ,Monoreagent
Roche	3183688122 / 4657357190 / 11970569216/11970909216	02R201	19R201	Colorimetric (BCG) ,Monoreagent
Roche	11970569216	02R204	19R204	Colorimetric (BCG) ,Monoreagent
Rajawali Nusindo	156092	02R301	19R301	Colorimetric (BCG) ,Monoreagent
Sysmex	76400	02S201	19S201	Colorimetric (BCG) ,Monoreagent
Sysmex	ALB Reagent	02S401	19S401	Bromcresol Green
Sclavo Diagnostics International S.r.l.	B74182708	02S601	19S601	Bromcresol Green
Thermo Scientific	981357 / 981766	02T101	19T101	Colorimetric (BCP)
Thermo Scientific	981767/981358	02T104	19T104	Bromcresol Green
Lain-lain	-	02XX99	19XX99	Lain-lain
Parameter : AST/SGOT 1/2				
Abbott	7D81, 37° C	03A102	20A102	IFCC without P5P, 37° C (& JSCC)
ABX Pentra	A11A01629, 30° C	03A201	20A201	IFCC without P5P, 30° C (& JSCC)
ABX Pentra	A11A01629, 37° C	03A202	20A202	IFCC without P5P, 37° C (& JSCC)
Advia Siemens	00952166, 30° C / 03039631, 30° C / 07499718, 30° C without P5P & JSCC	03A301	20A301	IFCC without P5P, 30° C (& JSCC)
Advia Siemens	00952166, 30° C / 03039631, 30° C / 07499718, 30° C with P5P	03A302	20A302	IFCC with P5P, 30° C
Advia Siemens	00952166, 37° C / 03039631, 37° C / 07499718, 37° C without P5P & JSCC	03A303	20A320	IFCC without P5P, 37° C (& JSCC)
Advia Siemens	00952166, 37° C / 03039631, 37° C / 07499718, 37° C with P5P	03A304	20A304	IFCC with P5P, 37° C
AIM	C5A6K0120R, 30° C / C5A6K0480R, 30° C	03A401	20A401	IFCC without P5P, 30° C (& JSCC)
AIM	C5A6K0120R, 37° C / C5A6K0480R, 37° C	03A402	20A402	IFCC without P5P, 37° C (& JSCC)
AMP Diagnostic	BR 0612, 30° C	03A501	20A501	IFCC without P5P, 30° C (& JSCC)
AMP Diagnostic	BR 0612, 37° C	03A502	20A502	IFCC without P5P, 37° C (& JSCC)
Agappe	51408003	03A601	20A601	IFCC
Biocon Diagnostics	1171, 30° C / 1176, 30° C / 1177, 30° C / 1178, 30° C / H7001, 30° C	03B101	20B101	IFCC without P5P, 30° C (& JSCC)
Biocon Diagnostics	1171, 37° C / 1176, 37° C / 1177, 37° C / 1178, 37° C / H7001, 37° C	03B102	20B102	IFCC without P5P, 37° C (& JSCC)
Bio Meriux	63212, 30° C	03B201	20B201	IFCC without P5P, 30° C (& JSCC)
Bio Meriux	63212, 37° C	03B202	20B202	IFCC without P5P, 37° C (& JSCC)
Biolabo	80025/80125/80225/80325	03B301	20B301	IFCC
Cypress Diagnostic	HBE06/HBEL01	03C101	20C101	IFCC without P5P, 37° C (& JSCC)

MEREK/PRODUSEN	NO KATALOG	KODE KIMIA 1	KODE KIMIA 2	METODE
Dimension Siemens	DF41A, 30° C	03D101	20D101	IFCC with P5P, 30° C
Dimension Siemens	DF41A, 37° C	03D102	20D102	IFCC with P5P, 37° C
Dialine Diagnostic	2251052	03D201	20D201	IFCC without P5P, 37° C (& JSCC)
Dialab	D00678, 30° C / D94610, 30° C / D98616, 30° C	03D301	20D301	IFCC without P5P, 30° C (& JSCC)
Dialab	D00678, 37° C / D94610, 37° C / D98616, 37° C	03D302	20D302	IFCC without P5P, 37° C (& JSCC)
DSI/DiaSys	1 2601 99 10 021, 30° C / 1 2601 99 10 026, 30° C	03D401	20D401	IFCC without P5P, 30° C (& JSCC)
DSI/DiaSys	1 2601 99 10 021, 37° C / 1 2601 99 10 026, 37° C	03D402	20D402	IFCC without P5P, 37° C (& JSCC)
DSI/Linear	1 2601 99 83 021, 30° C / 1 2601 99 83 026, 30° C	03D501	20D501	IFCC without P5P, 30° C (& JSCC)
DSI/Linear	1 2601 99 83 021, 37° C / 1 2601 99 83 026, 37° C	03D502	20D502	IFCC without P5P, 37° C (& JSCC)
Elitech	ASSL-0410, 30° C	03E101	20E101	IFCC without P5P, 30° C (& JSCC)
Elitech	ASSL-0410, 37° C	03E102	20E102	IFCC without P5P, 37° C (& JSCC)
Erba	XSYS016	03E201	20E201	IFCC, Modified
Fast	1-290/1-222/1-214/1-215/1-313/4-214/4-414/4-514	03F101	20F101	IFCC without P5P, 37° C (& JSCC)
Human Diagnostic	12211/12011/12021/12031/12021150/12021300/12021600	03H101	20H101	IFCC without P5P, 37° C (& JSCC)
Ortho (Johnson & Johnson)	6801453, 30° C	03J101	20J101	IFCC without P5P, 30° C (& JSCC)
Ortho (Johnson & Johnson)	6801453, 37° C	03J102	20J102	IFCC without P5P, 37° C (& JSCC)
Ortho (Johnson & Johnson)	8433815, 30° C	03J103	20J103	IFCC with P5P, 30° C
Ortho (Johnson & Johnson)	8433815, 37° C	03J104	20J104	IFCC with P5P, 37° C
LABTEST	75-4/30, 30° C	03L101	20L101	IFCC without P5P, 30° C (& JSCC)
LABTEST	75-4/30, 37° C	03L102	20L102	IFCC without P5P, 37° C (& JSCC)
Mindray	AST0102/AST0103/AST0105	03M201	20M201	IFCC Method, Kinetic
Olympus/ Beckman	OE60100, 30° C / OE60106, 30° C / QSR6009, 30° C / QSR6109, 30° C /	03O101	20O101	IFCC with P5P, 30° C
Olympus/ Beckman	OE60100, 37° C / OE60106, 37° C / QSR6009, 37° C / QSR6109, 37° C /	03O102	20O102	IFCC with P5P, 37° C
Pacific / Indec Diagnostic	UA-K023, 30° C	03P101	20P101	IFCC without P5P, 30° C (& JSCC)
Pacific / Indec Diagnostic	UA-K023, 37° C	03P102	20P102	IFCC without P5P, 37° C (& JSCC)
Proline	1 2601 99 10 021/1 2601 99 10 025/1 2601 99 10 029/1 2601 99 10 191	03P201	20P201	IFCC
Randox	AS 1204, 30° C / AS 1267, 30° C / AS 2359, 30° C	03R101	20R101	IFCC without P5P, 30° C (& JSCC)
Randox	AS 1204, 37° C / AS 1267, 37° C / AS 2359, 37° C	03R102	20R102	IFCC without P5P, 37° C (& JSCC)
Roche	10851124216, 30° C / 11876848216, 30° C / 20764949322, 30° C / 4657543190, 30° C	03R201	20R201	IFCC without P5P, 30° C (& JSCC)
Roche	10851124216, 37° C / 11876848216, 37° C / 20764949322, 37° C / 4657543190, 37° C	03R202	20R202	IFCC without P5P, 37° C (& JSCC)
Rajawali Nusindo	100191	03R301	20R301	IFCC without P5P, 30° C (& JSCC)
Rajawali Nusindo	100193	03R302	20R302	IFCC without P5P, 37° C (& JSCC)
Sysmex	89796, 30° C / 89806, 30° C	03S201	20S201	IFCC without P5P, 30° C (& JSCC)
Sysmex	89796, 37° C / 89806, 37° C	03S202	20S202	IFCC without P5P, 37° C (& JSCC)
Sysmex	AST.L-R1/AST.L.R2	03S301	20S301	UV Method, IFCC Method
Seimitsu Diagnostics	11S0011	03S501	20S501	IFCC without P5P, 37° C
Sclavo Diagnostics International S.r.l.	B74182714/B74182734/B74182755	03S601	20S601	Kinetic Enzymatic
Thermo Scientific	981771/981363	03T101	20T101	IFCC without P5P, 37° C (& JSCC)
Lain-lain	-	03XX99	20XX99	Lain-lain
Parameter : ALT/SGPT 1/ 2				
Abbott	7D56, 37° C	04A102	21A102	IFCC without P5P, 37° C (& JSCC)
ABX Pentra	A11A01627, 30° C	04A201	21A201	IFCC without P5P, 30° C (& JSCC)
ABX Pentra	A11A01627, 37° C	04A202	21A202	IFCC without P5P, 37° C (& JSCC)
Advia Siemens	02760612, 30° C / 03036926, 30° C / 07501976, 30° C without P5P & JSCC	04A301	21A301	IFCC without P5P, 30° C (& JSCC)
Advia Siemens	02760612, 30° C / 03036926, 30° C / 07501976, 30° C with P5P	04A302	21A302	IFCC with P5P, 30° C
Advia Siemens	02760612, 37° C / 03036926, 37° C / 07501976, 37° C without P5P & JSCC	04A303	21A303	IFCC without P5P, 37° C (& JSCC)
Advia Siemens	02760612, 37° C / 03036926, 37° C / 07501976, 37° C with P5P	04A304	21A304	IFCC with P5P, 37° C

MEREK/PRODUSEN	NO KATALOG	KODE KIMIA 1	KODE KIMIA 2	METODE
AIM	C5A4K0120R, 30° C / C5A4K0480R, 30° C	04A401	21A401	IFCC without P5P, 30° C (& JSCC)
AIM	C5A4K0120R, 37° C / C5A4K0480R, 37° C	04A402	21A402	IFCC without P5P, 37° C (& JSCC)
AMP Diagnostic	BR 0412, 30° C	04A501	21A501	IFCC without P5P, 30° C (& JSCC)
AMP Diagnostic	BR 0412, 37° C	04A502	21A502	IFCC without P5P, 37° C (& JSCC)
Agappe	51409003	04A601	21A601	IFCC
Biocon Diagnostics	1181, 30° C / 1186, 30° C / 1187, 30° C / 1188, 30° C / H6501, 30° C	04B101	21B101	IFCC without P5P, 30° C (& JSCC)
Biocon Diagnostics	1181, 37° C / 1186, 37° C / 1187, 37° C / 1188, 37° C / H6501, 37° C	04B102	21B102	IFCC without P5P, 37° C (& JSCC)
Bio Meriux	63312, 30° C	04B201	21B201	IFCC without P5P, 30° C (& JSCC)
Bio Meriux	63312, 37° C	04B202	21B202	IFCC without P5P, 37° C (& JSCC)
Biolabo	80027/80127/80227/80327	04B301	21B301	IFCC
Cypress Diagnostic	HBE07/HBEL02	04C101	21C101	IFCC without P5P, 37° C (& JSCC)
Dimension Siemens	DF143, 30° C	04D101	21D101	IFCC with P5P, 30° C
Dimension Siemens	DF143, 37° C	04D102	21D102	IFCC with P5P, 37° C
Dialab	D00640, 30° C / D94602, 30° C / D98624, 30° C	04D301	21D301	IFCC without P5P, 30° C (& JSCC)
Dialab	D00640, 37° C / D94602, 37° C / D98624, 37° C	04D302	21D302	IFCC without P5P, 37° C (& JSCC)
DSI/DiaSys	1 2701 99 10 021, 30° C / 1 2701 99 10 026, 30° C	04D401	21D401	IFCC without P5P, 30° C (& JSCC)
DSI/DiaSys	1 2701 99 10 021, 37° C / 1 2701 99 10 026, 37° C	04D402	21D402	IFCC without P5P, 37° C (& JSCC)
DSI/Linear	1 2701 99 83 021, 30° C / 1 2701 99 83 026, 30° C	04D501	21D501	IFCC without P5P, 30° C (& JSCC)
DSI/Linear	1 2701 99 83 021, 37° C / 1 2701 99 83 026, 37° C	04D502	21D502	IFCC without P5P, 37° C (& JSCC)
Dialine Diagnostic	2261052	04D602	21D602	IFCC without P5P, 37° C (& JSCC)
Elitech	ALSL-0410, 30° C	04E101	21E101	IFCC without P5P, 30° C (& JSCC)
Elitech	ALSL-0410, 37° C	04E102	21E102	IFCC without P5P, 37° C (& JSCC)
Erba	XSYS017	04E201	21E201	IFCC, Modified
Fast	1-289/1-221/1-216/1-217/1-312/4-216/4-416/4-516	04F101	21F101	IFCC without P5P, 37° C (& JSCC)
Human Diagnostic	12212/12012/12022/12032/12022150/12022300/12022600	04H101	21H101	IFCC with P5P, 37° C
Ortho (Johnson & Johnson)	1655281, 30° C	04J101	21J101	IFCC with P5P, 30° C
Ortho (Johnson & Johnson)	1655281, 37° C	04J102	21J102	IFCC with P5P, 37° C
Ortho (Johnson & Johnson)	6801454, 30° C	04J103	21J103	IFCC without P5P, 30° C (& JSCC)
Ortho (Johnson & Johnson)	6801454, 37° C	04J104	21J121	IFCC without P5P, 37° C (& JSCC)
LABTEST	74-4/30, 30° C	04L101	21L101	IFCC without P5P, 30° C (& JSCC)
LABTEST	74-4/30, 37° C	04L102	21L102	IFCC without P5P, 37° C (& JSCC)
Mindray	ALT0102/ALT0103/ALT0105	04M301	21M301	IFCC Method, Kinetic
Olympus/ Beckman	OE60100, 30° C / OE60106, 30° C / OSR6007, 30° C / OSR6507, 30° C	04O101	21O101	IFCC with P5P, 30° C
Olympus/ Beckman	OE60100, 37° C / OE60106, 37° C / OSR6007, 37° C / OSR6507, 37° C	04O102	21O102	IFCC with P5P, 37° C
Olympus/ Beckman	OSR6107, 30° C	04O107	21O107	IFCC without P5P, 30° C (& JSCC)
Olympus/ Beckman	OSR6107, 37° C	04O108	21O108	IFCC without P5P, 37° C (& JSCC)
Pacific / Indec Diagnostic	UA-K024, 30° C	04P101	21P101	IFCC without P5P, 30° C (& JSCC)
Pacific / Indec Diagnostic	UA-K024, 37° C	04P102	21P102	IFCC without P5P, 37° C (& JSCC)
Proline	1 2701 99 10 021/ 1 2701 99 10 025/ 1 2701 99 10 029/ 1 2701 99 10 191	04P201	21P201	IFCC
Randox	AL 1205, 30° C / AL 1268, 30° C / AL 2360, 30° C	04R101	21R101	IFCC without P5P, 30° C (& JSCC)
Randox	AL 1205, 37° C / AL 1268, 37° C / AL 2360, 37° C	04R102	21R102	IFCC without P5P, 37° C (& JSCC)
Roche	10851132216, 30° C / 11876805216, 30° C / 20764957322, 30° C / 4718569190, 30° C	04R201	21R201	IFCC without P5P, 30° C (& JSCC)
Roche	10851132216, 37° C / 11876805216, 37° C / 20764957322, 37° C / 4718569190, 37° C	04R202	21R202	IFCC without P5P, 37° C (& JSCC)
Rajawali Nusindo	100291	04R301	21R301	IFCC without P5P, 30° C (& JSCC)
Rajawali Nusindo	100293	04R302	21R302	IFCC without P5P, 37° C (& JSCC)
Sysmex	89850, 30° C / 89860, 30° C	04S201	21S201	IFCC without P5P, 30° C (& JSCC)

MEREK/PRODUSEN	NO KATALOG	KODE KIMIA 1	KODE KIMIA 2	METODE
Systemex	89850, 37° C / 89860, 37° C	04S202	21S202	IFCC without P5P, 37° C (& JSCC)
Systemex	ALT.L-R1/ALT.L.R2	04S401	21S401	UV Method, IFCC Method
Seimitsu Diagnostics	11L0011	04S501	21S501	IFCC without P5P, 37° C
Sclavo Diagnostics International S.r.l.	B74182713/B74182733/B74182754	04S601	21S601	Kinetic Enzymatic
Thermo Scientific	981361 / 981769	04T101	21T101	IFCC
Thermo Scientific	981769 / 981361	04T102	21T102	IFCC without Pxxodoxal Phospat
Lain-lain	-	04XX99	21XX99	Lain-lain
Parameter : Gamma GT 1 / 2				
Abbott	7D65-21	05A102	22A102	IFCC
ABX Pentra	A11A01630, 30° C	05A201	22A201	Szasz, 30° C
ABX Pentra	A11A01630, 37° C	05A202	22A202	Szasz, 37° C
Advia Siemens	02011954, 30° C / 07498649, 30° C	05A301	22A301	Szasz - IFCC recommended method, 30° C
Advia Siemens	02011954, 37° C / 07498649, 37° C	05A302	22A302	Szasz - IFCC recommended method, 37° C
AIM	C5G1K0120R, 30° C / C5G1K0480R, 30° C	05A401	22A401	Szasz, 30° C
AIM	C5G1K0120R, 37° C / C5G1K0480R, 37° C	05A402	22A402	Szasz, 37° C
AMP Diagnostic	BR 1202, 30° C	05A501	22A501	Szasz, 30° C
AMP Diagnostic	BR 1202, 37° C	05A503	22A503	Szasz, 37° C
Agappe	51416002	05A601	22A601	Szasz - IFCC recommended method
Biocon Diagnostics	2621, 30° C / 2625, 30° C / 2628, 30° C / H7501, 30° C	05B101	22B101	Szasz, 30° C
Biocon Diagnostics	2621, 37° C / 2625, 37° C / 2628, 37° C / H7501, 37° C	05B102	22B102	Szasz, 37° C
Bio Meriux	63712, 30° C	05B201	22B201	Szasz, 30° C
Bio Meriux	63712, 37° C	05B202	22B202	Szasz, 37° C
Biolabo	80110/80210/80310	05B301	22B301	Recomended IFCC
Cypress Diagnostic	HBE05	05C101	22C101	Szasz - IFCC recommended method
Dimension Siemens	DF45A, 30° C	05D101	22D101	Szasz - IFCC recommended method, 30° C
Dimension Siemens	DF45A, 37° C	05D102	22D102	Szasz - IFCC recommended method, 37° C
Dialab	D95604, 30° C / D98631, 30° C	05D301	22D301	Szasz - IFCC recommended method, 30° C
Dialab	D95604, 37° C / D98631, 37° C	05D302	22D302	Szasz - IFCC recommended method, 37° C
DSI/DiaSys	1 2801 99 10 021, 30° C / 1 2801 99 10 026, 30° C	05D401	22D401	Szasz - IFCC recommended method, 30° C
DSI/DiaSys	1 2801 99 10 021, 37° C / 1 2801 99 10 026, 37° C	05D402	22D402	Szasz - IFCC recommended method, 37° C
DSI/Linear	1 2801 99 83 021, 30° C / 1 2801 99 83 026, 30° C	05D501	22D501	Szasz - IFCC recommended method, 30° C
DSI/Linear	1 2801 99 83 021, 37° C / 1 2801 99 83 026, 37° C	05D502	22D502	Szasz - IFCC recommended method, 37° C
Dialine Diagnostic	2275052	05D601	22D601	Szasz - IFCC recommended method, 37° C
Elitech	GASL-0400, 30° C	05E101	22E101	Szasz - IFCC recommended method, 30° C
Elitech	GASL-0400, 37° C	05E102	22E102	Szasz - IFCC recommended method, 37° C
Erba	XSYS011	05E201	22E201	Szasz, 37° C
Fast	1-291/1-226/1-224/1-318/4-224/4-424/4-524	05F101	22F101	Szasz - IFCC recommended method
Human Diagnostic	12213/12013/12023/12033	05H101	22H101	Szasz - IFCC recommended method
Ortho (Johnson & Johnson)	8257289, 30° C	05J101	22J101	Szasz, 30° C
Ortho (Johnson & Johnson)	8257289, 37° C	05J102	22J102	Szasz, 37° C
LABTEST	105-2/30, 30° C	05L101	22L101	Szasz - IFCC recommended method, 30° C
LABTEST	105-2/30, 37° C	05L102	22L102	Szasz - IFCC recommended method, 37° C
Mindray	GGT0102/GGT0103/GGT0105	05M201	22M201	IFCC Method, Kinetic
Olympus/ Beckman	OSR6020, 30° C / OSR6120, 30° C / OSR6520, 30° C	05O101	22O101	Szasz - IFCC recommended method, 30° C
Olympus/ Beckman	OSR6020, 37° C / OSR6120, 37° C / OSR6520, 37° C	05O102	22O102	Szasz - IFCC recommended method, 37° C
Olympus/ Beckman	OSR6119, 30° C	05O103	22O103	Szasz, 30° C
Olympus/ Beckman	OSR6119, 37° C	05O104	22O104	Szasz, 37° C

MERЕК/PRODUSEN	NO KATALOG	KODE KIMIA 1	KODE KIMIA 2	METODE
Pacific / Indec Diagnostic	UA-K021, 30° C	05P101	22P101	Szasz, 30° C
Pacific / Indec Diagnostic	UA-K021, 37° C	05P102	22P102	Szasz, 37° C
Proline	1 2801 99 10 021/1 2801 99 10 191	05P201	22P201	Szasz - IFCC recommended method
Randox	GT 2750, 30° C	05R101	22R101	Szasz, 30° C
Randox	GT 2750, 37° C	05R102	22R102	Szasz, 37° C
Roche	12016788216, 30° C / 12016958216, 30° C / 3002721122, 30° C / 4657519190, 30° C	05R201	22R201	Szasz, 30° C
Roche	12016788216, 37° C / 12016958216, 37° C / 3002721122, 37° C / 4657519190, 37° C	05R202	22R202	Szasz, 37° C
Rajawali Nusindo	100391	05R301	22R301	Szasz, 30° C
Rajawali Nusindo	100392	05R302	22R302	Szasz, 37° C
Sysmex	80466, 30° C / 83536, 30° C	05S201	22S201	Szasz, 30° C
Sysmex	80466, 37° C / 83536, 37° C	05S202	22S202	Szasz, 37° C
Sysmex	Gamma.GTP(3C)-L.R1/Gamma.GTP(3C)-L.R2	05S401	22S401	Gamma glutamyltranspeptidase method, IFCC method
Seimitsu Diagnostics	11G0021	05S501	22S501	IFCC
Sclavo Diagnostics International S.r.l.	B74182723/B74182761	05S601	22S601	Enzymatic Kinetic Szasz
Thermo Scientific	981377 / 981778	05T101	22T101	IFCC
Lain-lain	-	05XX99	22XX99	Lain-lain
Parameter : Alkaline Phosphatase 1 / 2				
Abbott	7D55-22	06A102	23A102	AMP,IFCC
ABX Pentra	A11A01626, 30° C	06A201	23A201	IFCC, 30° C
ABX Pentra	A11A01626, 37° C	06A202	23A202	IFCC, 37° C
Advia Siemens	03035814, 30° C / 06258679, 30° C	06A301	23A301	IFCC, modified, 30° C
Advia Siemens	03035814, 37° C / 06258679, 37° C	06A302	23A302	IFCC, modified, 37° C
AIM	C5A3K0120R, 30° C / C5A3K0480R, 30° C	06A401	23A401	DGKC, 30° C
AIM	C5A3K0120R, 37° C / C5A3K0480R, 37° C	06A402	23A402	DGKC, 37° C
AMP Diagnostic	BR 0202, 30° C	06A501	23A501	DGKC, 30° C
AMP Diagnostic	BR 0202, 37° C	06A502	23A502	DGKC, 37° C
Agappe	51401003	06A601	23A601	DGKC
Biocon Diagnostics	1621, 30° C / 1622, 30° C / 1625, 30° C / H1501, 30° C	06B101	23B101	DGKC, 30° C
Biocon Diagnostics	1621, 37° C / 1622, 37° C / 1625, 37° C / H1501, 37° C	06B102	23B102	DGKC, 37° C
Biocon Diagnostics	1721, 30° C / 1722, 30° C / 1725, 30° C / H1601, 30° C / H1602, 30° C	06B107	23B107	IFCC, 30° C
Biocon Diagnostics	1721, 37° C / 1722, 37° C / 1725, 37° C / H1601, 37° C / H1602, 37° C	06B108	23B108	IFCC, 37° C
Bio Meriux	63609, 30° C	06B201	23B201	DGKC, 30° C
Bio Meriux	63609, 37° C	06B202	23B202	DGKC, 37° C
Bio Meriux	63657, 30° C	06B203	23B203	DGKC, 30° C
Bio Meriux	63657, 37° C	06B204	23B204	DGKC, 37° C
Biolabo	92214/92314	06B301	23B301	DGKC
Cypress Diagnostic	HBE02	06C101	23C101	DGKC, 37° C
Dimension Siemens	DF15A, 30° C	06D101	23D101	IFCC, modified, 30° C
Dimension Siemens	DF15A, 37° C	06D102	23D102	IFCC, modified, 37° C
Dialab	D00523, 30° C / D95565, 30° C	06D301	23D301	IFCC, modified, 30° C
Dialab	D00523, 37° C / D95565, 37° C	06D302	23D302	IFCC, modified, 37° C
DSI/DiaSys	1 0401 99 10 021, 30° C / 1 0401 99 10 026, 30° C	06D401	23D401	DGKC, 30° C
DSI/DiaSys	1 0401 99 10 021, 37° C / 1 0401 99 10 026, 37° C	06D402	23D402	DGKC, 37° C
DSI/DiaSys	1 0441 99 10 021, 30° C / 1 0441 99 10 026, 30° C	06D405	23D405	IFCC, 30° C
DSI/DiaSys	1 0441 99 10 021, 37° C / 1 0441 99 10 026, 37° C	06D406	23D423	IFCC, 37° C
DSI/Linear	1 0401 99 83 021, 30° C / 1 0401 99 83 026, 30° C	06D501	23D501	DGKC, 30° C
DSI/Linear	1 0401 99 83 021, 37° C / 1 0401 99 83 026, 37° C	06D502	23D502	DGKC, 37° C

MEREK/PRODUSEN	NO KATALOG	KODE KIMIA 1	KODE KIMIA 2	METODE
DSI/Linear	1 0441 99 83 021, 30° C / 1 0441 99 83 026, 30° C	06D505	23D505	IFCC, 30° C
DSI/Linear	1 0441 99 83 021, 37° C / 1 0441 99 83 026, 37° C	06D506	23D506	IFCC, 37° C
Dialine Diagnostic	2126052	06D601	23D601	IFCC, modified, 37° C
Elitech	PASL-0400, 30° C	06E101	23E101	DGKC, 30° C
Elitech	PASL-0400, 37° C	06E102	23E102	DGKC, 37° C
Erba	XSYS002	06E201	23E201	AMP,IFCC
Fast	2-291/1-218/1-212/1-317/4-212/4-412/4-512	06F101	23F101	DGKC, 37° C
Human Diagnostic	12017, 30° C/12017, 37° C/12217, 30° C/12217, 37° C/12027, 30° C/12027, 37° C	06H101	23H101	DGKC
Human Diagnostic	12117/12117600	06H102	23H102	IFCC
Human Diagnostic	12117, 30° C	06H401	23H401	IFCC
Human Diagnostic	12117, 37° C	06H402	23H402	IFCC
Ortho (Johnson & Johnson)	1053180, 30° C	06J101	23J101	IFCC, modified, 30° C
Ortho (Johnson & Johnson)	1053180, 37° C	06J102	23J102	IFCC, modified, 37° C
LABTEST	79-4/30, 30° C	06L101	23L101	IFCC, modified, 30° C
LABTEST	79-4/30, 37° C	06L102	23L102	IFCC, modified, 37° C
Mindray	ALP0102/ALP0103/ALP0105	06M201	23M201	AMP Buffer Method, Kinetic
Olympus/ Beckman	OSR6004, 30° C / OSR6104, 30° C / OSR6204, 30° C	06O101	23O101	IFCC, 30° C
Olympus/ Beckman	OSR6004, 37° C / OSR6104, 37° C / OSR6204, 37° C	06O102	23O102	IFCC, 37° C
Olympus/ Beckman	OSR6103, 30° C	06O103	23O103	EAE Buffer, JSCC, 30° C
Olympus/ Beckman	OSR6103, 37° C	06O104	23O104	EAE Buffer, JSCC, 37° C
Pacific / Indec Diagnostic	UA-K011, 30° C	06P101	23P101	DGKC, 30° C
Pacific / Indec Diagnostic	UA-K011, 37° C	06P102	23P102	DGKC, 37° C
Proline	1 0441 99 10 190	06P301	23P301	IFCC
Randox	AP 132, 30° C / AP 311, 30° C	06R101	23R101	EAE Buffer, JSCC, 30° C
Randox	AP 132, 37° C / AP 311, 37° C	06R102	23R102	EAE Buffer, JSCC, 37° C
Randox	AP 307, 30° C / AP 542, 30° C	06R103	23R103	DGKC, 30° C
Randox	AP 307, 37° C / AP 542, 37° C	06R104	23R104	DGKC, 37° C
Roche	10816388216, 30° C / 11877348216, 30° C	06R201	23R201	DGKC, 30° C
Roche	10816388216, 37° C / 11877348216, 37° C	06R202	23R202	DGKC, 37° C
Roche	12172933122, 30° C / 12173107122, 30° C / 3333701190, 30° C / 3333752190, 30° C /	06R205	23R205	IFCC, 30° C
Roche	12172933122, 37° C / 12173107122, 37° C / 3333701190, 37° C / 3333752190, 37° C /	06R206	23R206	IFCC, 37° C
Rajawali Nusindo	107091	06R301	23R301	DGKC, 37° C
Rajawali Nusindo	107092	06R302	23R302	DGKC, 30° C
Sysmex	80420, 30° C / 82820, 30° C	06S201	23S201	EAE Buffer, JSCC, 30° C
Sysmex	80420, 37° C / 82820, 37° C	06S202	23S202	EAE Buffer, JSCC, 37° C
Sysmex	ALP-L.R1/ALP-L.R2	06S301	23S301	pNPP Method, EAE Buffer Method
Seimitsu Diagnostics	11A0011	06S501	23S501	IFCC
Sclavo Diagnostics International S.r.l.	B74182722/B74182760	06S601	23S601	IFCC Method, Kinetic
Thermo Scentific	981832/981833	06T101	23T101	IFCC Plus
Lain-lain	-	06XX99	23XX99	Lain-lain
Parameter : Bilirubin Total 1 / 2				
Abbott	6L45	07A101	24A101	Colorimetric/DPD, Bilirubin Direct
ABX Pentra	A11A01639	07A201	24A201	Colorimetric/DPD, Bilirubin Direct
Advia Siemens	02188064	07A301	24A301	Jendrassik-Groff
AIM	C5B1T0200R / C5B1T0400R	07A401	24A401	Colorimetric/DPD, Bilirubin Direct
AMP Diagnostic	BR 2202	07A501	24A501	Colorimetric DCA
Agappe	52002005	07A601	24A601	TOOS Enzymatic Double End Point

MEREK/PRODUSEN	NO KATALOG	KODE KIMIA 1	KODE KIMIA 2	METODE
Biocon Diagnostics	804	07B101	24B101	Colorimetric DCA
Biocon Diagnostics	1965 / 1969	07B102	24B102	Jendrassik-Groff
Bio Meriux	61037	07B201	24B201	Colorimetric DCA
Biolabo	80403/80443/97443	07B301	24B301	Colorimetric DCA
Cypress Diagnostic	HB027/HB026	07C101	24C101	Colorimetric/DPD, Bilirubin Direct
Dimension Siemens	DF67A	07D101	24D101	Colorimetric DCA
Dialab	102011 / 102012 / 102021 / 102022	07D301	24D301	Jendrassik-Groff
Dialab	D96531 / D96543	07D305	24D305	Colorimetric/DPD, Bilirubin Direct
DSI/DiaSys	1 0811 99 10 021 / 1 0811 99 10 026	07D401	24D401	Colorimetric/DPD, Bilirubin Direct
DSI/DiaSys	1 0849 99 90 336	07D403	24D403	Jendrassik-Groff
DSI/DiaSys	108119910962	07D404	24D404	Colorimetric DCA
DSI/Linear	1 0811 99 83 021 / 1 0811 99 83 026	07D501	24D501	Colorimetric/DPD, Bilirubin Direct
DSI/Linear	1 0849 99 83 336	07D503	24D503	Jendrassik-Groff
Dialine Diagnostic	2141052	07D601	24D601	Colorimetric DCA
Elitech	BILI-0600	07E101	24E101	Colorimetric DCA
Elitech	BITO-0600	07E102	24E102	malloy evelyn modified-end point
Erba	XSYS023	07E201	24E201	DIAZO + Hexa Fluoro Phosphat (VITROS)
Fast	2-214/2-245/2-246/2-296/4-245/4-445/4-545	07F101	24F101	Jendrassik-Groff
Human Diagnostic	10012	07H101	24H101	Colorimetric DCA
Human Diagnostic	10740	07H102	24H102	Jendrassik-Groff
Human Diagnostic	10741/10742/10741150/10741300/10741600/10742150/10742300/1742600	07H103	24H103	Colorimetric/DPD/DIAZO method
Ortho (Johnson & Johnson)	8159931	07J101	24J101	DIAZO + Hexa Fluoro Phosphat (VITROS)
LABTEST	94	07L102	24L102	Colorimetric DCA
Mindray	TBI0102/TBI0103/TBI0105	07M201	24M201	DSA Method, End Point
Mindray	TBI0202/TBI0203/TBI0205	07M401	24M401	VOX Method, End Point
Olympus/ Beckman	OSR6112 / OSR6212	07O101	24O101	Colorimetric/DPD, Bilirubin Total, DPD + Caffeine
Pacific / Indec Diagnostic	UA-K014	07P101	24P101	Colorimetric/DPD, Bilirubin Total, DPD + Caffeine
Proline	1 0811 99 10 021/1 0811 99 0 191	07P201	24P201	Colorimetric DCA
Randox	BR 2361	07R101	24R101	Jendrassik-Groff
Roche	3146022122 / 3261638190 / 4774256190 / 11489194216 / 11552414216 /	07R201	24R201	Colorimetric/DPD, Bilirubin Total, DPD + Caffeine
Rajawali Nusindo	101291	07R301	24R301	Jendrassik-Groff
Rajawali Nusindo	101292	07R302	24R302	Jendrassik-Groff
Rajawali Nusindo	101293	07R303	24R303	Jendrassik-Groff
Rajawali Nusindo	101294	07R304	24R304	Jendrassik-Groff
Sysmex	70876 / 76346	07S201	24S201	Colorimetric/DPD, Bilirubin Direct
Sysmex	T.BIL Buffer (R1)/T.Bil Diazo Reagent (R2)	07S401	24S401	Stabilized diazo method
Seimitsu Diagnostics	11B0031	07S501	24S501	DPD
Sclavo Diagnostics International S.r.l.	B74182701/B74182746	07S601	24S601	Winsten-cehelik
Thermo Scientific	981364	07T101	24T101	Colorimetric (DIAZO method) + Sulfanilic Acid
Thermo Scientific	981793	07T103	24T103	NBD
Lain-lain	-	07XX99	24XX99	Lain-lain
Parameter : Glukosa 1 / 2				
Abbott	3L82	08A101	25A101	Hexokinase
ABX Pentra	A11A01667	08A201	25A201	Hexokinase
ABX Pentra	A11A01668	08A202	25A202	GOD-PAP
Advia Siemens	04005536	08A301	25A301	GOD-PAP
Advia Siemens	04903429	08A302	25A302	Hexokinase

MEREK/PRODUSEN	NO KATALOG	KODE KIMIA 1	KODE KIMIA 2	METODE
AIM	C5G2O0400R / C5G2O1000R	08A401	25A401	GOD-PAP
AMP Diagnostic	BR 2906	08A501	25A501	GOD-PAP
AMP Diagnostic	BR 3006	08A502	25A502	Hexokinase
Agappe	51406001	08A601	25A601	GOD-PAP
Biocon Diagnostics	4341 / 4342 / 4347 / H5701	08B101	25B101	GOD-PAP
Bio Meriux	61269 / 61270	08B201	25B201	GOD-PAP
Biolabo	LP80209/87109/LP80809/16GL8	08B301	25B301	GOD-PAP
Cypress Diagnostic	HB009/HB010/HBL04/HBL05	08C101	25C101	GPO-PAP
Dimension Siemens	DF39A	08D101	25D101	Hexokinase
Dialab	D00221 / D00223 / D95218B	08D301	25D301	GOD-PAP
Dialab	D06227	08D303	25D303	Hexokinase
DSI/DiaSys	1 2500 99 10 021 / 1 2500 99 10 023 / 1 2500 99 10 026	08D401	25D401	GOD-PAP
DSI/DiaSys	1 2511 99 10 021 / 1 2511 99 10 026	08D404	25D404	Hexokinase
DSI/Linear	1 2500 99 83 021 / 1 2500 99 83 023 / 1 2500 99 83 026	08D501	25D501	GOD-PAP
DSI/Linear	1 2511 99 83 021 / 1 2511 99 83 026	08D504	25D504	Hexokinase
Dialine Diagnostic	2245018	08D601	25D601	GOD-PAP
Dialine Diagnostic	2249052	08D602	25D602	Hexokinase
Elitech	GHSL-0600	08E101	25E101	Hexokinase
Elitech	GPSL-0490	08E102	25E102	GOD-PAP
Erba	XSYS012	08E201	25E201	GOD POD
Fast	2-218/2-219/2-201/2-203/4-201/4-401/4-501	08F101	25F101	GOD-PAP
Human Diagnostic	10260/10121/10260150/10260300/10260600	08H101	25H101	GOD-PAP
Human Diagnostic	10786	08H102	25H102	Hexokinase
INDOREAGEN	G0520845 / G0530845	08I102	25I102	GOD-PAP
Ortho (Johnson & Johnson)	1707801	08J101	25J101	GOD-PAP
LABTEST	84-1/500 / 84-2/500 / 84-3/500	08L101	25L101	GOD-PAP
LABTEST	85-4/50	08L104	25L104	Hexokinase
Mindray	GLU0102/GLU0103/GLU0105	08M201	25M201	GPO-POD Method, End Point
Mindray	GLU0202/GLU0203/GLU0205	08M301	25M301	HK Method, End Point
Olympus/ Beckman	OSR6121 / OSR6140 / OSR6221 / OSR6521	08O101	25O101	Hexokinase
Pacific / Indec Diagnostic	UA-K022	08P101	25P101	GOD-PAP
Proline	1 2511 99 10 022/1 2511 99 10 025/1 2511 99 10 029/1 2511 99 10 0191	08P201	25P201	Hexokinase
Randox	GL 2614 / GL 2623 / GL 364 / GL 366	08R101	25R101	GOD-PAP
Roche	4404483190 / 4657527190 / 11876899216 / 20767131322/11447513216	08R201	25R201	Hexokinase
Roche	11448668216 / 11491253216	08R205	25R205	GOD-PAP
Rajawali Nusindo	112191	08R301	25R301	GOD-PAP
Rajawali Nusindo	112192	08R302	25R302	GOD-PAP
Sysmex	80586 / 82946	08S201	25S201	Hexokinase
Sysmex	Glu-L.R1/Glu-L.R2	08S401	25S401	Enzyme method, Hexokinase UV method
Seimitsu Diagnostics	11G0011	08S501	25S501	Hexokinase
Sclavo Diagnostics International S.r.l.	B74182718/B74182738/B74182757	08S601	25S601	Enzymatic Oxidase
Thermo Scientific	981379 / 981780	08T101	25T101	GOD-PAP
Thermo Scientific	981779/981304	08T102	25T102	Hexokinase
Lain-lain	-	08XX99	25XX99	Lain-lain
Parameter : Ureum 1 / 2				
Abbott	7D75	09A101	26A101	Urease - GLDH b
ABX Pentra	A11A01641	09A201	26A201	Urease - GLDH b

MEREK/PRODUSEN	NO KATALOG	KODE KIMIA 1	KODE KIMIA 2	METODE
Advia Siemens	01574998 / 03040257	09A301	26A301	Urease - GLDH b
AIM	C5B2K0120R / C5B2K0480R	09A401	26A401	Urease - GLDH b
AMP Diagnostic	BR 4002	09A501	26A501	Urease - GLDH b
Agappe	51412002	09A601	26A601	Urease - GLDH
Biocon Diagnostics	4701 / 4705 / 4708	09B101	26B101	Urease - GLDH b
Biocon Diagnostics	5111 / 5113	09B104	26B104	Berthelot (Urease-Colorimetric) , R1 + standar
Bio Meriux	61912 / 61913	09B201	26B201	Berthelot (Urease-Colorimetric) , R1 + standar
Bio Meriux	61974	09B203	26B203	Urease - GLDH b
Biolabo	80221	09B301	26B301	Colorimetric Berthelot
Biolabo	80221	09B302	26B302	Kinetic
Cypress Diagnostic	HB022/HB029	09C101	26C101	Berthelot (Urease-Colorimetric) , R1 + standar
Dimension Siemens	DF21	09D101	26D101	Urease - GLDH b
Dialab	402999	09D301	26D301	Berthelot (Urease-Colorimetric) , R1 + standar
Dialab	D00716 / D98707	09D302	26D302	Urease - GLDH b
DSI/DiaSys	1 3101 99 10 021 / 1 3101 99 10 026	09D401	26D401	Urease - GLDH b
DSI/DiaSys	1 3115 99 10 021 / 1 3115 99 10 026	09D403	26D403	Berthelot (Urease-Colorimetric) , R1 + standar
DSI/Linear	1 3101 99 83 021 / 1 3101 99 83 026	09D501	26D501	Urease - GLDH b
DSI/Linear	1 3115 99 83 021 / 1 3115 99 83 026	09D503	26D503	Berthelot (Urease-Colorimetric) , R1 + standar
Dialine Diagnostic	2371052	09D601	26D601	Urease - GLDH b
Elitech	URCO-0605	09E101	26E101	Berthelot (Urease-Colorimetric) , R1 + standar
Elitech	URSL-0407	09E102	26E102	Urease - GLDH b
Erba	XSYS020	09E201	26E201	Urease - GLDH
Fast	2-223/2-261/2-206/2-207/2-310/4-206/4-406/4-506	09F101	26F101	Urease - GLDH b
Human Diagnostic	10521/10521150/10521300/10521600	09H102	26H102	Urease - GLDH b
Ortho (Johnson & Johnson)	8102204	09J101	26J101	Berthelot (Urease-Colorimetric) , R1 + standar
LABTEST	27	09L101	26L101	Berthelot (Urease-Colorimetric) , R1 + standar
LABTEST	104-4/250	09L102	26L102	Urease - GLDH b
Mindray	URE0102/URE0103/URE0105	09M201	26M201	Urease-GLDH, UV Method, Fixed Time
Olympus/ Beckman	OSR6134 / OSR6141 / OSR6234 / OSR6534	09O101	26O101	Urease - GLDH b
Pacific / Indec Diagnostic	UA-K034	09P101	26P101	Urease - GLDH
Proline	1 3101 99 10 022/1 3101 99 10 025/1 3101 99 10 029/1 3101 99 10 191	09P201	26P201	Urease - GLDH
Randox	UR 220 / UR 221	09R101	26R101	Urease - GLDH b
Randox	UR107	09R103	26R103	Berthelot (Urease-Colorimetric) , R1 + standar
Roche	4460715190 / 4657616190 / 11489364216 / 11729691216	09R201	26R201	Urease - GLDH b
Rajawali Nusindo	110491	09R301	26R301	Urease - GLDH
Rajawali Nusindo	110493	09R302	26R302	Urease - GLDH
Rajawali Nusindo	110494	09R303	26R303	Urease - GLDH
Rajawali Nusindo	110495	09R304	26R304	Urease - GLDH
Systemex	84896 / 84897	09S201	26S201	Urease - GLDH b
Seimitsu Diagnostics	11U0031	09S501	26S501	Urease - GLDH
Sclavo Diagnostics International S.r.l.	B74182721/B74182741/B74182743	09S601	26S601	Enzymatic (Urease) UV
Thermo Scientific	981820/981818	09T101	26T101	Urease - UV
Lain-lain	-	09XX99	26XX99	Lain-lain
Parameter : Kreatinine 1 / 2				
Abbott	3L81-22	10A101	27A101	Jaffe, Modified
Abbott	8L24-31	10A103	27A103	Enzymatic Photometric (TOOS)
ABX Pentra	A11A01666	10A201	27A201	Jaffe method - Kinetic for automation

MEREK/PRODUSEN	NO KATALOG	KODE KIMIA 1	KODE KIMIA 2	METODE
ABX Pentra	A11A01868	10A202	27A202	Jaffe method - Kinetic for automation
ABX Pentra	A11A01907	10A203	27A203	Enzymatic Colorimetric
ABX Pentra	A11A01933	10A204	27A204	Jaffe method - Kinetic for automation
Advia Siemens	0309070 / 01437451	10A301	27A301	Jaffe method, rate-blanked compensated
Advia Siemens	04992596	10A303	27A303	Enzymatic Colorimetric
AIM	C5C8K0200R	10A401	27A401	Jaffe method - Kinetic for automation
AMP Diagnostic	BR 2802	10A501	27A501	Jaffe method - Kinetic for automation
Agappe	51009001	10A601	27A601	Jaffe method - Kinetic
Biocon Diagnostics	426	10B101	27B101	Enzymatic Colorimetric
Biocon Diagnostics	447 / 448	10B102	27B102	Jaffe method - Kinetic for automation
Bio Meriux	61162	10B201	27B201	Jaffe method - Kinetic for automation
Biolabo	80107	10B301	27B301	Kinetic
Cypress Diagnostic	HB008	10C101	27C101	Jaffe method - Kinetic
Dimension Siemens	DF33A	10D101	27D101	Jaffe method - Kinetic for automation
Dialab	D00617 / D95595	10D301	27D301	Jaffe method - Kinetic for automation
DSI/DiaSys	1 1711 99 10 021 / 1 1711 99 10 026	10D401	27D401	Jaffe method - Kinetic for automation
DSI/DiaSys	117599910963	10D403	27D403	Enzymatic Colorimetric PAP
DSI/Linear	1 1711 99 83 021 / 1 1711 99 83 026	10D501	27D501	Jaffe method - Kinetic for automation
Dialine Diagnostic	2235052	10D701	27D701	Enzymatic Colorimetric
Dialine Diagnostic	2230052	10D703	27D703	Jaffee without deproteinization, with compensation
Elitech	CRCO-0600	10E101	27E101	Jaffe method - Kinetic for automation
Elitech	CRSL-0630	10E102	27E102	Enzymatic Colorimetric
Erba	XSYS024	10E201	27E201	Alkaline Picrate
Fast	2-220/2-232/2-233/2-299/4-233/4-433/4-533	10F201	27F201	Jaffe method - Kinetic
Human Diagnostic	10051/10052/10052150/10052300/10052600	10H101	27H101	Jaffe method - Kinetic
Ortho (Johnson & Johnson)	8141947	10J101	27J101	Enzymatic Colorimetric
LABTEST	96-300	10L101	27L101	Jaffe method - Kinetic for automation
Mindray	CRE0102/CRE0103/CRE0105	10M201	27M201	Jaffe, Modified
Mindray	CRE0202/CRE0203/CRE0205	10M301	27M301	Sarcosine Oxidase Method, End Point
Olympus/ Beckman	OSR61204	10O101	27O101	Enzymatic Colorimetric
Olympus/ Beckman	OSR6178* / OSR6578*	10O102	27O102	Jaffe method - Kinetic for automation
Pacific / Indec Diagnostic	UA-K020	10P101	27P101	Jaffe method, rate-blanked compensated
Proline	1 1759 99 10 022/ 1 1759 99 10 191	10P201	27P201	Enzymatic Colorimetric
Proline	1 1711 99 10 022/1 1711 99 10 025/1 1711 99 10 029/1 1711 99 10 191	10P202	27P202	Jaffe method - Kinetic for automation
Randox	CR 2336 / CR 2337 / CR 510 / CR 524	10R101	27R101	Jaffe method - Kinetic for automation
Roche	3263991190 / 4657489190 / 11775642216 / 11775685216/5401470190	10R201	27R201	Enzymatic Colorimetric
Roche	4810716190 / 4810716190 / 11875418216/11489291216	10R203	27R203	Jaffe method, rate compensated
Roche	11489291216	10R208	27R208	Jaffe method, rate-blanked compensated
Rajawali Nusindo	105191	10R301	27R301	Jaffe method - Kinetic for automation
Sysmex	83446 / 83466	10S201	27S201	Enzymatic Colorimetric
Sysmex	CRE-L.R1/CRE-L.R2	10S301	27S301	Enzyme Method, SOD-PAP colorimetric method
Seimitsu Diagnostics	11C0021	10S502	27S502	Enzymatic Colorimetric
Sclavo Diagnostics International S.r.l.	B74182712/B74182732/B74182753	10S601	27S601	Colorimetric Method Jaffe
Thermo Scientific	981845	10T101	27T101	Enzymatic Colorimetric
Thermo Scientific	981811/981810	10T102	27T102	Jaffe, Modified
Lain-lain	-	10XX99	27XX99	Lain-lain

Parameter : Asam Urat 1 / 2

MEREK/PRODUSEN	NO KATALOG	KODE KIMIA 1	KODE KIMIA 2	METODE
Abbott	3P39	11A101	28A101	Uricase-Peroxidase
ABX Pentra	A11A01670	11A201	28A201	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
Advia Siemens	02871651 / 03051305 / 07497014	11A301	28A301	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
AIM	C5U1E0120R / C5U1E0480R	11A401	28A401	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
AMP Diagnostic	BR 4203	11A501	28A501	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
Agappe	51413002	11A601	28A601	Uricase-Peroxidase
Biocon Diagnostics	4841 / 4842 / 4845 / 4848	11B101	28B101	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
Bio Meriux	61923	11B201	28B201	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
Biolabo	80351/80001/87601	11B301	28B301	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
Cypress Diagnostic	HB023/HB024/HBL02	11C101	28C101	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
Dimension Siemens	DF77	11D101	28D101	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
Dialab	D95459 / D98461	11D301	28D301	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
DSI/DiaSys	1 3021 99 10 021 / 1 3021 99 10 026	11D401	28D401	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
DSI/DiaSys	130019910962	11D403	28D403	Enzymatic Photometric (TOOS)
DSI/Linear	1 3021 99 83 021 / 1 3021 99 83 026	11D501	28D501	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
DSI/Linear	1 3001 99 83 021	11D502	28D502	Uricase-Peroxidase
Dialine Diagnostic	2385052	11D601	28D601	Uricase-Peroxidase
Elitech	AUML-0427	11E101	28E101	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
Fast	2-224/2-235/2-208/2-209/2-308/4-208/4-408/4-508	11F101	28F101	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
Human Diagnostic	10690/10691/10694/10694150/10694300/10694600	11H101	28H101	Uricase-Peroxidase
Ortho (Johnson & Johnson)	1943927	11J101	28J101	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
LABTEST	73-2/100 / 73-4/30	11L101	28L101	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
Mindray	UA0102/UA0103/UA0105	11M201	28M201	Uricase-Peroxidase
Olympus/ Beckman	OSR6098 / OSR6198 / OSR6298	11O101	28O101	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
Pacific / Indec Diagnostic	UA-K035	11P101	28P101	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
Proline	1 3021 99 10 022/1 3021 99 10 025/1 3021 99 10 029/1 3021 99 10 191	11P301	28P301	Fotometri Enzimatik menggunakan TBHA
Randox	UA 230 / UA 233	11R101	28R101	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
Roche	3183807190 / 4657608190 / 11661850216 / 11875426216	11R201	28R201	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
Rajawali Nusindo	103691	11R301	28R301	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
Rajawali Nusindo	103692	11R302	28R302	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
Sysmex	80606 / 89236	11S201	28S201	Uricase-Peroxidase, Mono reagent PAP with LCF (Lipid Clearing
Sysmex	UA-L.R1/UA-L.R2	11S301	28S301	Enzyme method, Uricase-PAP colorimetric method
Seimitsu Diagnostics	11U0011	11S501	28S501	TOOS Enzymatic Double End Point
Sclavo Diagnostics International S.r.l.	B74182764/B74182765/B74182766	11S601	28S601	Enzymatic Method (Uricase) With the Detection System Emerson-
Thermo Scientific	981391/981788	11T101	28T101	Uricase, Trinder
Thermo Scientific	981816/981817	11T102	28T102	Uricase, Trinder
Lain-lain	-	11XX99	28XX99	Lain-lain
Parameter : Kolesterol Total 1 / 2				
Abbott	7D62	12A101	29A101	CHOD-PAP Enzymatic photometric test
ABX Pentra	A11A01634	12A201	29A201	CHOD-PAP Enzymatic photometric test
Advia Siemens	01482198	12A301	29A301	CHOD-PAP Enzymatic photometric test
AIM	C5C4E0100R / C5C4E0400R	12A401	29A401	CHOD-PAP Enzymatic photometric test
AMP Diagnostic	BR 2606	12A501	29A501	CHOD-PAP Enzymatic photometric test
Agappe	51403003	12A601	29A601	CHOD-PAP Enzymatic photometric test
Biocon Diagnostics	4241 / 4242 / 4245 / 4248 / H4001	12B101	29B101	CHOD-PAP Enzymatic photometric test
Bio Meriux	61218 / 61219	12B201	29B201	CHOD-PAP Enzymatic photometric test
Biolabo	80106/87656/87356/LP80106	12B301	29B301	CHOD-PAP Enzymatic photometric test

MEREK/PRODUSEN	NO KATALOG	KODE KIMIA 1	KODE KIMIA 2	METODE
Cypress Diagnostic	HB006	12C101	29C101	CHOD-PAP Enzymatic photometric test
Dimension Siemens	DF27	12D101	29D101	CHOD-PAP Enzymatic photometric test
Sekisui	241673 / 241680	12D201	29D201	CHOD-PAP Enzymatic photometric test
Dialab	D00119 / D95116 / D96112B	12D301	29D301	CHOD-PAP Enzymatic photometric test
DSI/DiaSys	1 1300 99 10 021 / 1 1300 99 10 023 / 1 1300 99 10 026	12D401	29D401	CHOD-PAP Enzymatic photometric test
DSI/Linear	1 1300 99 83 021 / 1 1300 99 83 023 / 1 1300 99 83 026	12D501	29D501	CHOD-PAP Enzymatic photometric test
Dialine Diagnostic	2185022	12D601	29D601	CHOD-PAP Enzymatic photometric test
Elitech	CHSL-0490	12E101	29E101	CHOD-PAP Enzymatic photometric test
Elitech	CHSL-5505	12E103	29E102	CHOD-PAP Enzymatic photometric test
Erba	XSYS009	12E201	29E201	CHOD-PAP Enzymatic photometric test
Fast	2-212/2-211/2-204/2-205/2-298/4-236/4-436/4-536	12F101	29F101	CHOD-PAP Enzymatic photometric test
Human Diagnostic	100177/10019/10028/10028150/10028300/10028600	12H101	29H101	CHOD-PAP Enzymatic photometric test
Ortho (Johnson & Johnson)	1669829	12J101	29J101	CHOD-PAP Enzymatic photometric test
LABTEST	76-2/100 / 76-2/250	12L101	29L101	CHOD-PAP Enzymatic photometric test
Mindray	TC0102/TC0103/TC0105	12M201	29M201	CHOD-POD Method, End Point
Olympus/ Beckman	OSR6116 / OSR6216 / OSR6516	12O101	29O101	CHOD-PAP Enzymatic photometric test
Pacific / Indec Diagnostic	UA-K016	12P101	29P101	CHOD-PAP Enzymatic photometric test
Proline	1 1300 99 10 022/1 1300 99 10 025/1 1300 99 10 029/1 1300 99 10 192	12P201	29P201	CHOD-PAP Enzymatic photometric test
Randox	CH 200 / CH 201 / CH 202	12R101	29R101	CHOD-PAP Enzymatic photometric test
Roche	3039773190 / 4718917190 / 11489232216 / 11491458216	12R201	29R201	CHOD-PAP Enzymatic photometric test
Rajawali Nusindo	101592	12R301	29R301	CHOD-PAP Enzymatic photometric test
Sysmex	89036 / 89056	12S201	29S201	CHOD-PAP Enzymatic photometric test
Sysmex	T.CHO-LA.R1/T.CHO-LA.R2	12S301	29S301	Enzyme method, COD-PAP colorimetric method
Seimitsu Diagnostics	11T0031	12S502	29S502	Enzymatic Colorimetric
Sclavo Diagnostics International S.r.l.	B74182716/B74182736/B74182756	12S601	29S601	Enzymatic Methode (esterase)
Thermo Scientific	981813/981812	12T101	29T101	CHOD-PAP Enzymatic photometric test
Lain-lain	-	12XX99	29XX99	Lain-lain
Parameter : Trigliserida 1 / 2				
Abbott	7D74	13A101	30A101	GPO-PAP
ABX Pentra	A11A01640	13A201	30A201	GPO-PAP
Advia Siemens	09580156	13A301	30A301	GPO-PAP
AIM	C5T3E0100R / C5T3E0400R	13A401	30A401	GPO-PAP
AMP Diagnostic	BR 4503	13A501	30A501	GPO-PAP
Agappe	51410002	13A601	30A601	GPO-PAP
Biocon Diagnostics	5741 / 5742 / 5745 / 5748 / H9001	13B101	30B101	GPO-PAP
Bio Meriux	61236	13B201	30B201	GPO-PAP
Biolabo	80019/87319	13B301	30B301	GPO
Cypress Diagnostic	HB021	13C101	30C101	GPO-PAP
Dimension Siemens	DF69A	13D101	30D101	GPO-PAP
Sekisui	225796 / 225802	13D201	30D201	GPO-PAP - with Glycerol blanking
Dialab	D00389 / D96338 / D98386B	13D301	30D301	GPO-PAP
DSI/DiaSys	1 5710 99 10 021 / 1 5710 99 10 023 / 1 5710 99 10 026	13D401	30D401	GPO-PAP
DSI/Linear	1 5710 99 83 021 / 1 5710 99 83 023 / 1 5710 99 83 026	13D501	30D501	GPO-PAP
Dialine Diagnostic	2365022	13D601	30D601	GPO-PAP
Elitech	TGML-0427	13E101	30E101	GPO-PAP
Erba	XSYS041	13E201	30E201	GPO
Fast	2-222/2-262/2-253/2-254/2-311/4-253/4-453/4-553	13F101	30F101	GPO-PAP

MEREK/PRODUSEN	NO KATALOG	KODE KIMIA 1	KODE KIMIA 2	METODE
Human Diagnostic	10720P/10724/10725/10724150/10724300/ 10724600	13H101	30H101	GPO-PAP
Ortho (Johnson & Johnson)	1336544 / 8329930	13J101	30J101	GPO-PAP
LABTEST	87-2/100 / 87-2/250	13L101	30L101	GPO-PAP
Mindray	TG0102/TG0103/TG0105	13M201	30M201	GPO-POD Method, End Point
Olympus/ Beckman	OSR60118 / OSR61118 / OSR6133 / OSR62118	13O101	30O101	GPO-PAP
Pacific / Indec Diagnostic	UA-K033	13P101	30P101	GPO-PAP
Proline	1 5710 99 10 022/1 5710 99 10 025/1 5710 99 10 029/1 5710 99 10 192	13P301	30P301	GPO-PAP
Randox	TR 210 / TR 212 / TR 213	13R101	30R101	GPO-PAP
Roche	4657594190 / 11488872216 / 11730711216 / 20767107322	13R201	30R201	GPO-PAP
Rajawali Nusindo	116392	13R301	30R301	GPO-PAP
Sysmex	80943 / 80944	13S201	30S201	GPO-PAP - with Glycerol blanking
Systemx	TG-L.R1/TG-L.R2	13S301	30S301	Enzyme colorimetric method of eliminating glycerol
Seimitsu Diagnostics	11T0041	13S501	30S501	Enzymatic with Glycerol Blanking
Sclavo Diagnostics International S.r.l.	B74182719/B74182739/B74182758	13S601	30S601	lipoprotein-lipase (LPL)
Thermo Scientific	981786/981301	13T101	30T101	GPO-PAP
Lain-lain	-	13XX99	30XX99	Lain-lain
Parameter : Kalsium Total (Ca) 1 / 2				
Abbott	3L79	14A101	31A101	Arsenazo III
ABX Pentra	A11A01633	14A201	31A201	OCPC
ABX Pentra	A11A01894	14A202	31A202	Arsenazo III
ABX Pentra	A11A01954	14A204	31A204	OCPC
Advia Siemens	02794363 / 07506552	14A301	31A301	OCPC
AIM	ECAE0200	14A401	31A401	OCPC
AMP Diagnostic	BR 7002	14A501	31A501	OCPC
Agappe	51006001	14A601	31A601	Arsenazo III
Biocon Diagnostics	1903	14B101	31B101	OCPC
Biocon Diagnostics	2003	14B102	31B102	Arsenazo III
Bio Meriux	61041	14B203	31B203	Methyl Thymol Blue
Cypress Diagnostic	HB003	14C101	31C101	Arsenazo III
Cypress Diagnostic	HB004	14C102	31C102	OCPC
Cornley	CORNLEY K-LITE SERIES	14C201	31C201	ISE - Direct
Dimension Siemens	DF23A	14D101	31D101	OCPC
Dialab	D01375 / D01376	14D301	31D301	Arsenazo III
Dialab	D95098	14D303	31D303	OCPC
DSI/DiaSys	1 1121 99 10 021	14D401	31D401	OCPC
DSI/DiaSys	1 1130 99 10 021	14D402	31D402	Arsenazo III
DSI/DiaSys	1 1181 99 10 962	14D403	31D403	Phosphonazo III
DSI/Linear	1 1121 99 83 021	14D501	31D501	OCPC
DSI/Linear	1 1130 99 83 021	14D502	31D502	Arsenazo III
Elitech	CALA-0600	14E101	31E101	Arsenazo III
Elitech	CALO-0600	14E102	31E102	OCPC
Erba	XSYS007	14E201	31E201	Arsenazo III
Fast	3-249/3-250/3-251/3-252/3-324/4-251/4- 451/4551	14F101	31F101	OCPC
Fast	3-253/3-254/3-255/3-256/3-325/4-247/4- 447/4-547	14F102	31F102	Arsenazo III
Human Diagnostic	10011/10011150/10011300/10011600	14H101	31H101	OCPC
Ortho (Johnson & Johnson)	1450261	14J101	31J101	Arsenazo III
LABTEST	95-2/50	14L101	31L101	Arsenazo III

MEREK/PRODUSEN	NO KATALOG	KODE KIMIA 1	KODE KIMIA 2	METODE
Mindray	CA0102/CA0103/CA0105	14M201	31M201	Arsenazo III
Olympus/ Beckman	OSR60117 / OSR61117 / OSR6176 / OSR6276	14O101	31O101	Arsenazo III
Olympus/ Beckman	OSR6113 / OSR6213	14O103	31O103	OCPC
Pacific / Indec Diagnostic	UA-K015	14P101	31P101	Arsenazo III
Phoenix Diagnostic	3-800/3-801/3-802/3-223/3-224	14P201	31P201	ISE - Direct
Randox	CA 590	14R101	31R101	OCPC
Roche	4718933190 / 11489216216 / 11730240216 / 20763128322	14R201	31R201	OCPC
Roche	5061482190/05061334190	14R202	31R202	NM-BAPTA
Rajawali Nusindo	101191	14R301	31R301	OCPC
Sysmex	76416 / 76426	14S201	31S201	OCPC
Sysmex	Ca Buffer (R1)/ Ca Color Solution (R2)	14S401	31S401	OCPC
Sclavo Diagnostics International S.r.l.	B74182715/B74182735	14S601	31S601	Arsenazo III
Thermo Alko	A478-855/A478-856/A478-861/A478-858	14T101	31T101	ISE - Direct
Thermo Scientific	981367/981772	14T201	31T201	Arsenazo III
Lain-lain	-	14XX99	31XX99	Lain-lain
Parameter : Natrium (Na) 1 / 2				
Abbott	2P32	15A101	32A101	ISE - Indirect
ABX Pentra	-	15A201	32A201	ISE - Direct
Advia Siemens	-	15A301	32A301	ISE - Indirect
AIM	ENAE100	15A401	32A401	Mg-uranylacetate / Colorimetric
Biocare	R0010 / 1030108	15B101	32B101	ISE - Direct
Cypress Diagnostic	HB016	15C101	32C101	Mg-uranylacetate / Colorimetric
Cornley	CORNLEY K-LITE SERIES	15C201	32C201	ISE - Direct
Dimension Siemens	S500A	15D101	32D101	IMT - Direct
Dimension Siemens	S600	15D102	32D102	IMT - Indirect
Human Diagnostic	573351	15H101	32H101	ISE - Direct
Human Diagnostic	17410	15H102	32H102	ISE - Direct
Ortho (Johnson & Johnson)	8379034	15J101	32J101	ISE - Direct
Jokoh	JOKOH EX-Z/D	15J301	32J301	ISE - Direct
MEDICA	2102	15M101	32M101	ISE - Direct
Nova Biomedical	15646/15648/15196/14608	15N101	32N101	ISE - Direct
Olympus/ Beckman	OE66313 ISE Direct / OE66314 ISE Direct / OE66315 ISE Direct / OE66316 ISE Direct /	15O101	32O101	ISE - Direct
Olympus/ Beckman	OE66313 ISE In-direct / OE66314 ISE In-direct / OE66315 ISE In-direct / OE66316	15O102	32O102	ISE - Indirect
Phoenix Diagnostic	14-103	15P201	32P201	ISE - Direct
Phoenix Diagnostic	3-800/3-801/3-802/3-223/3-224	15P202	32P202	ISE - Direct
Randox	NA 7167	15R101	32R101	Mg-uranylacetate / Colorimetric
Roche	3110419180/ 3144046001/3112349180	15R201	32R201	ISE - Direct
Roche	10825468001 / 21029371001	15R202	32R202	ISE - Indirect
SFRI	R0401 / R0404	15S101	32S101	ISE - Direct
Audicom		15S601	32S601	Ion Selective Elektrode (ISE)
Thermo Alko	A478-855/A478-856/A478-861/A478-858	15T201	32T201	ISE - Direct
Lain-lain	-	15XX99	32XX99	Lain-lain
Parameter : Kalium(K) 1 / 2				
Abbott	2P32	16A101	33A101	ISE - Indirect
ABX Pentra	-	16A201	33A201	ISE - Direct
Advia Siemens	-	16A301	33A301	ISE - Indirect
AIM	EKAE0100	16A401	33A401	Sodium tetraphenylboron (turbidimetric/nephelometric)

MEREK/PRODUSEN	NO KATALOG	KODE KIMIA 1	KODE KIMIA 2	METODE
Audicom	-	16A601	33A601	Ion Selective Elektrode (ISE)
Biocare	R0010 / 1030108	16B101	33B101	ISE - Direct
Cypress Diagnostic	HB015	16C101	33C101	Sodium tetraphenylboron (turbidimetric/nephelometric)
Cornley	CORNLEY K-LITE SERIES	16C201	33C201	ISE - Direct
Dimension Siemens	S500A	16D101	33D101	IMT - Direct
Dimension Siemens	S600	16D102	33D102	IMT - Indirect
Human Diagnostic	10118	16H101	33H101	Sodium tetraphenylboron (turbidimetric/nephelometric)
Human Diagnostic	17410	16H102	33H102	ISE - Direct
Ortho (Johnson & Johnson)	8157596	16J101	33J101	ISE - Direct
Jokoh	JOKOH EX-Z/D	16J301	33J301	ISE - Direct
MEDICA	2101	16M101	33M101	ISE - Direct
Nova Biomedical	15646/15648/15196/14608	16N101	33N101	ISE - Direct
Olympus/ Beckman	OE66313 ISE Direct / OE66314 ISE Direct / OE66315 ISE Direct / OE66316 ISE Direct /	16O101	33O101	ISE - Direct
Olympus/ Beckman	OE66313 ISE In-direct / OE66314 ISE In-direct / OE66315 ISE In-direct / OE66316	16O102	33O102	ISE - Indirect
Phoenix Diagnostic	14-103	16P201	33P201	ISE - Direct
Phoenix Diagnostic	3-800/3-801/3-802/3-223/3-224	16P202	33P202	ISE - Direct
Randox	PT 1600	16R101	33R101	Sodium tetraphenylboron (turbidimetric/nephelometric)
Roche	3110338180/ 3144046001/3112349180	16R201	33R201	ISE - Direct
Roche	10825441001 / 21029355001	16R202	33R202	ISE - Indirect
SFRI	R0401 / R0404	16S101	33S101	ISE - Direct
Thermo Alko	A478-855/A478-856/A478-861/A478-858	16T101	33T101	ISE - Direct
Lain-lain	-	16XX99	33XX99	Lain-lain
Parameter : Klorida (Cl) 1 / 2				
Abbott	2P32	17A101	34A101	ISE - Indirect
ABX Pentra	-	17A201	34A201	ISE - Direct
Advia Siemens	-	17A301	34A301	ISE - Indirect
AIM	ECLE0100	17A401	34A401	Colorimetric / Photometric thvocyanat
AMP Diagnostic	BR 7302	17A501	34A501	Colorimetric / Photometric thvocyanat
Audicom	-	17A601	34A601	Ion Selective Elektrode (ISE)
Biocon Diagnostics	619	17B101	34B101	Colorimetric / Photometric thvocyanat
Biocare	R0010 / 1030108	17B102	34B102	ISE - Direct
Biolabo	80005	17B103	34B103	Colorimetric (BCP)
Cypress Diagnostic	HB005	17C101	34C101	Colorimetric / Photometric thvocyanat
Cornley	CORNLEY K-LITE SERIES	17C201	34C201	ISE - Direct
Dimension Siemens	S500A	17D101	34D101	IMT - Direct
Dimension Siemens	S600	17D102	34D102	IMT - Indirect
Dialab	D01223	17D301	34D301	Colorimetric / Photometric thvocyanat
DSI/DiaSys	1 1200 99 10 021	17D401	34D401	Colorimetric / Photometric thvocyanat
DSI/Linear	1 1200 99 83 021	17D501	34D501	Colorimetric / Photometric thvocyanat
Elitech	CHLO-0600	17E101	34E101	Colorimetric / Photometric thvocyanat
Erba	XSYS008	17E201	34E201	Mercuric Thiocyanate
Human Diagnostic	10115	17H101	34H101	Colorimetric / Photometric thvocyanat
Human Diagnostic	17410	17H102	34H102	ISE - Direct
Ortho (Johnson & Johnson)	8445207	17J101	34J101	ISE - Direct
Jokoh	JOKOH EX-Z/D	17J301	34J301	ISE - Direct
LABTEST	49	17L101	34L101	Colorimetric / Photometric thvocyanat
MEDICA	2113	17M101	34M101	ISE - Direct

MERKEK/PRODUSEN	NO KATALOG	KODE KIMIA 1	KODE KIMIA 2	METODE
Nova Biomedical	15646/15648/15196/14608	17N101	34N101	ISE - Direct
Olympus/ Beckman	OE66313 ISE Direct / OE66314 ISE Direct / OE66315 ISE Direct / OE66316 ISE Direct / OE66317 ISE Direct / OE66318 ISE Direct / OE66319 ISE Direct / OE66320 ISE Direct	17O101	34O101	ISE - Direct
Olympus/ Beckman	OE66313 ISE In-direct / OE66314 ISE In- direct / OE66315 ISE In-direct / OE66316 ISE In-direct / OE66317 ISE In-direct / OE66318 ISE In-direct / OE66319 ISE In- direct / OE66320 ISE In-direct	17O102	34O102	ISE - Indirect
Pacific / Indec Diagnostic	UA-K045	17P101	34P101	Modified cyan methemoglobin
Phoenix Diagnostic	14-103	17P201	34P201	ISE - Direct
Phoenix Diagnostic	3-800/3-801/3-802/3-223/3-224	17P202	34P202	ISE - Direct
Randox	CL 1645	17R101	34R101	Colorimetric / Photometric thiocyanat
Roche	3003523001 / 3246353001	17R201	34R201	ISE - Indirect
Roche	3110451180/ 3144046001/3112349180	17R202	34R202	ISE - Direct
SFRI	R0401 / R0404	17S101	34S101	ISE - Direct
Thermo Alko	A478-855/A478-856/A478-861/A478-858	17T101	34T101	ISE - Direct
Lain-lain	-	17XX99	34XX99	Lain-lain