

Multidimensional NMR Experiments

Chem 4010/5326:
Organic Spectroscopic Analysis

2D NMR Spectroscopy

General Information

- More complicated experiments to set up than ^1H and ^{13}C
- Changes in pulses (#, length, angles, mixing times, etc.)
- Observe effects based on relationship of nuclei
 - **Can be homonuclear (same nuclei) (e.g. H–H) or heteronuclear (different nuclei) (e.g. H–C, H–P, etc.)**
- Will only go over the what the experiments tell you and how to interpret, and only for the most common and widely used for solving organic structures
 - DEPT, H–H COSY, HMQC (HETCOR), HMBC, INADEQUATE
- There is a whole alphabet soup of other experiments (both 1D and 2D)
 - EXSY, TOCSY, HOHAHA, INEPT, WATERGATE, and many more

2D NMR Spectroscopy

Number of Protons on Carbon

DEPT (*Distortionless Enhancement by Polarisation Transfer*)

- Used to be known as APT (Attached Proton Test)
 - DEPT is ^1H -detected; APT is ^{13}C -detected
- Tells you how many protons are attached to a particular carbon
 - negative peaks = CH_2
 - positive peaks = CH and CH_3 (distinguishable with further processing)
 - “missing” peaks = carbons w/o protons
- With a little help from IR and chemical shift of ^1H and ^{13}C , can get a rough idea of molecular weight

2D NMR Spectroscopy


Number of Protons on Carbon

DEPT (Distortionless Enhancement by Polarisation Transfer)

1 = Carbon spectrum

2 = ↓ CH₂
↑ CH, CH₃

3 = ↑ CH




2D NMR Spectroscopy


Who is Talking to Who?

^1H - ^1H COSY (Correlation Spectroscopy)

- Tells you which protons are coupled to one another
- Very useful when peaks are overlapping in ^1H NMR and you are unable to calculate coupling constants, or when there are a lot of similar coupling constants
- Cross peaks are coupled to each other



1 bond H-H coupling




- Newer method is DQF (Double Quantum Filtered)-COSY
 - same information, but looks “cleaner”

2D NMR Spectroscopy

Who is Talking to Who?

^1H - ^1H COSY (Correlation Spectroscopy)

- Overlapping protons and a lot of similar coupling constants




Raffinose—a trisaccharide

2D NMR Spectroscopy

Who is Talking to Who?

^1H - ^1H COSY (Correlation Spectroscopy)




Raffinose—a trisaccharide

2D NMR Spectroscopy

Who is Talking to Who?

^1H - ^1H COSY (Correlation Spectroscopy)





2D NMR Spectroscopy

Who is Talking to Who?

^1H - ^1H COSY (Correlation Spectroscopy)

– DQF-COSY: *Double Quantum Filtered COSY* – cleans up the spectrum by reducing noncoupled systems (e.g. CH_3 singlets)





2D NMR Spectroscopy

Who is Talking to Who?

^1H - ^1H COSY (Correlation Spectroscopy)



– DQF-COSY: *Double Quantum Filtered COSY* – cleans up the spectrum by reducing noncoupled systems (e.g. CH_3 singlets)



2D NMR Spectroscopy

Who is Talking to Who?


^1H - ^1H COSY (Correlation Spectroscopy)



2D NMR Spectroscopy

Who is Talking to Who?


^1H - ^1H COSY (Correlation Spectroscopy)




Recall...

^1H - ^1H COSY

- Tells you which protons are coupled to one another
- Very useful when peaks are overlapping in ^1H NMR and you are unable to calculate coupling constants, or when there are a lot of similar coupling constants
- Cross peaks are coupled to each other



2 & 3 bond H-H coupling



Recall...


^1H - ^1H COSY

This can still cause ambiguities!


– Tells you which protons are coupled to one another

– Very useful when peaks are overlapping in ^1H NMR and you are unable to calculate coupling constants, or when there are a lot of similar coupling constants

– Cross peaks are coupled to each other



2 & 3 bond H-H coupling



2D NMR Spectroscopy

Who is Talking to Who?

^1H - ^1H TOCSY (Total Correlation Spectroscopy)

- Tells you which protons are in the same spin system (a continuous chain of spin-spin coupled protons)
- Magnetization from H_A is transferred to H_B , which then transfers to H_C and on down the line




- Can be run as either a 1D or 2D experiment
- Running as a 1D experiment greatly simplifies spectra with severe signal overlap
- Related experiment (HOHAHA) gives essentially the same information
- Can also observe one-bond ^1H - ^{13}C couplings (HMQC-TOCSY, hetero-TOCSY, HEHAHA) - will not discuss

1D TOCSY

Spin System Identification



^1H - ^1H TOCSY (Total Correlation Spectroscopy)

- Irradiate at frequency of proton of interest
- Wait a period of time (typically 20 to 400 ms) and then collect 1D spectrum
- Only observe protons to which magnetization has been transferred
- Short periods of time (20 ms) will only give rise to “single step” transfers; longer periods of time allow magnetization to propagate further





1D TOCSY



Spin System Identification



1D TOCSY



1D TOCSY






2D TOCSY

Spin System Identification


^1H - ^1H TOCSY (Total Correlation Spectroscopy)

– Cross peaks are in the same spin system




2D TOCSY

– Cross peaks are in the same spin system




Intermediate in biomimetic synthesis of manzamine




2D TOCSY

– Cross peaks are in the same spin system





Intermediate in biomimetic synthesis of manzamine



2D TOCSY

– Cross peaks are in the same spin system




New Developments

Deconvolution via MDEC


MDEC (Multi Frequency Homonuclear Decoupling)

- Conceptually the same as the homonuclear decoupling experiment mentioned earlier, but allows multiple frequencies to be decoupled *at the same time*



- (a) ¹H NMR
- (b) irr @ H-1
- (c) irr @ H-1/H-7
- (d) irr @ H-1/H-3a/H-7



J. Am. Chem. Soc. **2009**, *131*,
15994–15995



New Developments

Deconvolution via MDEC

1D-TOCSY-MDEC




- (a) ^1H NMR
- (b) 1D-TOCSY
irr @ H-3
- (c) 1D-TOCSY
irr @ H-3 with
MDEC @ H-1b,
H-2b, H-3

2D NMR Spectroscopy

Who is Talking to Who?

^1H - ^{13}C COSY

- HETCOR (Heteronuclear Correlation)
 - older experiment; ^{13}C -detected
- HMQC (Heteronuclear Multiple Quantum Correlation) and HSQC (Heteronuclear Single Quantum Correlation)
 - newer experiments; ^1H -detected; largely replaced HETCOR
- Both give same information, experimentally very different
- Peaks have one-bond coupling (i.e. attached directly)
- Compliments DEPT
- Particularly useful for diastereotopic protons





1 bond H-C coupling

2D NMR Spectroscopy

Who is Talking to Who?

HMQC




1 bond H-C coupling


2D NMR Spectroscopy

Who is Talking to Who?

HMQC



Diastereotopic
Protons




1 bond H-C coupling


2D NMR Spectroscopy

Who is Talking to Who?

HMQC



Can “see into”
multiplets




Caryophyllene Oxide



2D NMR Spectroscopy

Who is Talking to Who?


HMQC




2D NMR Spectroscopy



2D NMR Spectroscopy




2D NMR Spectroscopy



2D NMR Spectroscopy

Who is Talking to Who?

- COLOC (Correlated spectroscopy for Long range Couplings)
 - older experiment; ^{13}C -detected
- HMBC (Heteronuclear Multiple Bond Coherence)
 - newer experiment; ^1H -detected; completely replaced COLOC
- Both give same information, experimentally very different
- Peaks have two- or three-bond coupling
- “Sees through” heteroatoms and quaternary carbons
- Can be very complicated, but is very powerful




2 and 3 bond H-C couplings

2D NMR Spectroscopy

Who is Talking to Who?

^1H - ^{13}C COSY (Long Range)

- HETCOR (Heteronuclear Correlation)
 - older experiment; ^{13}C -detected
- HMQC (Heteronuclear Multiple Quantum Correlation) and HSQC (Heteronuclear Single Quantum Correlation)
 - newer experiments; ^1H -detected; largely replaced HETCOR
- Both give same information, experimentally very different
- Peaks have one-bond coupling (i.e. attached directly)
- Compliments DEPT
- Particularly useful for diastereotopic protons





1 bond H-C coupling

2D NMR Spectroscopy

Who is Talking to Who?

HMQC





1 bond H-C coupling

2D NMR Spectroscopy

Who is Talking to Who?



HMBC




2 and 3 bond H-C couplings

2D NMR Spectroscopy

Who is Talking to Who?





or





2D NMR Spectroscopy

Who is Talking to Who?





2D NMR Spectroscopy

Who is Talking to Who?





2D NMR Spectroscopy

Who is Talking to Who?



2D NMR Spectroscopy

Who is Talking to Who?



2D NMR Spectroscopy

Who is Talking to Who?

^{13}C – ^{13}C COSY


- INADEQUATE (Incredible Natural Abundance Double Quantum Transfer Experiment)
- tells what carbons are attached to each other
- if you know what type of carbon it is (C, C=O, CH, CH₂, CH₃, etc.) from DEPT, you can almost write down the entire gross structure by running two NMR experiments

– **BUT** it is ^{13}C – ^{13}C coupling

– probability of one ^{13}C is 0.01

– two next to each other $0.01 \times 0.01 = 0.0001$
(~1 molecule in 10,000)

- Need lots of sample and instrument time to overcome
- In our facility: 80% v/v, overnight, 500 MHz = nothing




1 bond C–C coupling

2D NMR Spectroscopy

Who is Talking to Who?

INADEQUATE




2D NMR Spectroscopy


Who is Talking to Who?

INADEQUATE



- Cross peaks show up as doublets = J_{CC}
- Diagonal is midway between the two doublets



2D NMR Spectroscopy




2D NMR Spectroscopy




2D NMR Spectroscopy


MASS




$^{13}\text{C}/\text{DEPT}$ NMR 150.9 MHz




2D NMR Spectroscopy




2D NMR Spectroscopy



2D NMR Spectroscopy



2D NMR Spectroscopy



2D NMR Spectroscopy

