



Π

TE

TK

T

T

T



countrycaregroup.com.au

National Distribution, Local Service



Country Care Group is a national group of service providers that carry the full range of products required to meet the clinical needs of its clients. The group has a proven track record in delivering equipment solutions within tight timeframes to the DVA, Hospitals, **Occupational Therapists, Health Care** Centres and other service providers.



CONTENTS

Mobility

Showering Furniture

Toileting Equipment

Patient Chairs

Transfer Aids

Pressure Cushions

Country Care Group consists of privately owned and operated businesses where the owners work within the business providing you with the highest level of service possible. Currently the group have over 80 showrooms and warehouses in various locations across Australia, and operate in excess of 300 delivery vehicles.

The groups experienced team members have a vast wealth of knowledge within the rehabilitation products field, and understand fully the needs of their clients.

CONTACT CCG

FREE CALL 1800 727 382

FACSIMILE

1800 329 382

EMAIL contracts@countrycaregroup.com.au

MEMBER LOCATIONS See page 218

countrycaregroup.com.au







8	
80	
95	
95	T
107	
136	



COUNTRY CARE GROUP SERVICES

HOME MODIFICATIONS

Country Care Group provide a home modifications service which specialises in falls prevention and home access solutions. We work alongside Occupational Therapists and Architects whose clients are elderly or mobility impaired in their homes to make them safe and accessible.



Country Care Group has a large network of qualified tradesmen who have been trained to specialise in all areas of the home accessibility and falls prevention. Whatever the project, whether bathroom modifications, ramp constructions, platform steps, bannisters, grab rails, lever taps or any minor modifications, we are happy to help. We currently provide services to the Department of Veterans Affairs, TAC, SWEP, Aged Care Facilities and the general public in your local area.

Call us to discuss your home modification requirements. Onsite appointments available to ensure the job is performed to your requirements.

COUNTRY CARE GROUP SERVICES

SERVICE AND MAINTENANCE

Country Care Group offers a repairs and maintenance service through our highly qualified technicians. Country Care Group service and repair all mobility products, electric lift chairs, bedding, electric lifting aids including ceiling hoist and stair lifts.

- Qualified technicians
- Test and tag service
- Repairs to all equipment
- General maintenance
- Supply and installation of parts
- Call out service









COUNTRY CARE GROUP SERVICES

HIRE DEPARTMENT

At Country Care Group we offer a wide variety of hire equipment for everything you need for home care and recovery.

Bathroom & Toilet

Bath Boards, Over Toilet Aids, Shower Chairs & Stools, Commodes, Transfer Bench.

Bedroom

COUNTRY CARE GROUP

Electric Beds & Recliner Chairs, Bed Stick, Alternating Air Mattress, Bed Cradles.

Manual Handling

Lifters Hoists and Slings, Transfer Boards.

Mobility Aids

Walking Frames, Rollators, Crutches, Gutter Frames, Transit Wheelchairs, Self Propelling Wheelchairs, Power Chairs, Scooters, Portable Ramps.

Bariatric Equipment

Shower Chairs, Wheelchair, Bedside Commodes, Electric Beds, Pick up Frames, Walking Frames.

A STATE OF A

COUNTRY CARE GROUP SERVICES

DELIVERY

National Distribution, Local Service

Country Care Group supports a total delivery system, which ensures that your products are transported by caring professionals who will deliver, assemble and demonstrate correct operations and procedures of all our items and equipment.



81 LOCATIONS AUSTRALIA WIDE

MOBILITY



KARMA FLEXX SELF PROPELLED FOLDING WHEELCHAIR

Flexible, Robust and Stylish. Sturdy double crossbrace. Adjustable center of gravity for greater stability. Heavy duty functions to assist plus-size clients in need of a mobility device. Designed with key adjustments to create a suitable sitting environment and fit individual needs. Folding lightweight frame. Removable, height adjustable, swing in/out footrests. Angle adjustable backrest. Swing away, height adjustable arm rests with ergonomic pad. Quick release wheels. Different sizes and accessories available.



Country Care

GROUP National Distribution, Local Service



Don't forget to check out the entire range of CCG Mobility products at countrycaregroup.com.au

Wheelchair

3.01 Aluminium Self Propelling Folding Ultra Lightweight

Karma Ergo Lite Deluxe Self-Propel



· Adjusts and stabilises seating pose Double cross braces provide strength and durability

- Available in deluxe style with swing-away removable footrests
- Folding frame
- Push to lock brakes, fixed seat height and folding back

KM2512F20-16 Karma Ergo Lite Deluxe Self Propelled Wheelchair 16"

Backrest Height	420 mm	Floor To Seat Height	480 mm	Folded Width	310 mm	Max User Weight	100 kg
Overall Height	910 mm	Overall Length	930 mm	Overall Width	630 mm	Product Weight	10.2 kg
Seat Depth	405 mm	Seat Width	405 mm				

KM2512F20-18 Karma Ergo Lite Deluxe Self Propelled Wheelchair 18"

Backrest Height	420 mm	Floor To Seat Height	480 mm	Folded Width	310 mm	Max User Weight	100 kg
Overall Height	910 mm	Overall Length	930 mm	Overall Width	650 mm	Product Weight	10.3 kg
Seat Depth	405 mm	Seat Width	455 mm				

Wheelchair

Karma S-Ergo 125 Light Weight Self-Propelled Wheelchair



- Uses Karma's S-Ergo patented seating system · Promotes pressure relief
- Helps prevent the user slipping off the seat
- Adjusts and stabilises seating pose
- Treated with aegie microbe shield reducing stains and odours
- Various accessories/customisations available

KMS125-16 Karma S-Ergo 125 Lightweight Self Propelled Wheelchair 16" Bac Ove Sea

Backrest Height	400 mm	Floor To Seat Height	510 mm	Folded Width	300 mm	Max User Weight	115 kg
Overall Height	910 mm	Overall Length	1000 mm	Overall Width	590 mm	Product Weight	14.2 kg
Geat Depth	430 mm	Seat Width	405 mm				

KMS125-18 Karma S-Ergo 125 Lightweight Self Propelled Wheelchair 18"

Backrest Height	400 mm	Floor To Seat Height	510 mm	Folded Width	300 mm	Max User Weight	115 kg
Overall Height	910 mm	Overall Length	1000 mm	Overall Width	640 mm	Product Weight	14.5 kg
Seat Depth	430 mm	Seat Width	455 mm				

• Uses Karma's S-Ergo patented seating system · Promotes pressure relief and helps prevent the user slipping off the seat

• The deluxe is the lightest self propelling wheelchair from Karma

3.05 Aluminium Self Propelling Folding Lightweight

· Streamline oval tubes add strength to this lightweight wheelchair

KMS125-20 Karma S-Ergo 125 Lightweight Self Propelled Wheelchair 20"

Backrest Height	400 mm	Floor To Seat Height	510 mm	Folded Width	300 mm	Max User Weight	115 kg
Overall Height	910 mm	Overall Length	1000 mm	Overall Width	690 mm	Product Weight	14.8 kg
Seat Depth	430 mm	Seat Width	505 mm				

Breezy BasiX 2 Wheelchair



• The Breezy Basix 2 is a lightweight aluminium, folding wheelchair, suitable for a wide range of users

- Adjustable axle plate Adjustable seat depth
- Adjustable seat height
- 14 to 20 Inch seat width
- Height adjustable armrest
- Also available in transit option Tension adjustable backrest
- Push-to-lock brakes & adjustable wheel positions (vertical and horizontal)

Backrest Height	435 mm	Max User Weight	125 kg	Product Weight	14.9 kg	Max Seat Depth	460 mm
Min Seat Depth	410 mm	Seat Width	410 mm				
0 74101-006 Bre	ezy BasiX 2	Self Propelled Whee	elchair 5	10mm Silver with Dr	um Brak	es	
Backrest Height	435 mm	Max User Weight	125 kg	Product Weight	14.9 kg	Max Seat Depth	460 mm
Min Seat Depth	410 mm	Seat Width	510 mm				
D74101-003 Bree Backrest Height	435 mm	Max User Weight	125 kg	10mm Silver with Dr Product Weight	14.9 kg	Max Seat Depth	460 mm
Dealment Unight	425 mm	Max lloar Weight	105 kg	Droduct Woight	14.0 kg	May Cast Donth	400 mm
Min Seat Depth	410 mm	Seat Width	410 mm				
		Colf Dropollod W/bc	olohoir	ICOmm Dod			
		Self Propelled Whe Max User Weight	elchair 4 125 kg	160mm Red Product Weight	14.9 kg	Max Seat Depth	460 mm
0 74100-023 Bre	ezy BasiX 2				14.9 kg	Max Seat Depth	460 mm
074100-023 Bre Backrest Height Min Seat Depth 074100-007 Bre	ezy BasiX 2 435 mm 410 mm	Max User Weight Seat Width Self Propelled Whe	125 kg 460 mm elchair 3	Product Weight S50mm Silver		Max Seat Depth	
074100-023 Bre Backrest Height Min Seat Depth	ezy BasiX 2 435 mm 410 mm	Max User Weight Seat Width	125 kg 460 mm	Product Weight	14.9 kg 14.9 kg	Max Seat Depth Max Seat Depth	460 mm
074100-023 Bre Backrest Height Min Seat Depth 074100-007 Bre	ezy BasiX 2 435 mm 410 mm ezy BasiX 2	Max User Weight Seat Width Self Propelled Whe	125 kg 460 mm elchair 3	Product Weight S50mm Silver			
074100-023 Bre Backrest Height Min Seat Depth 074100-007 Bre Backrest Height Min Seat Depth	ezy BasiX 2 435 mm 410 mm ezy BasiX 2 435 mm 410 mm	Max User Weight Seat Width Self Propelled Whe Max User Weight	125 kg 460 mm elchair 3 125 kg 355 mm	Product Weight SOmm Silver Product Weight			
074100-023 Bre Backrest Height Min Seat Depth 074100-007 Bre Backrest Height Min Seat Depth	ezy BasiX 2 435 mm 410 mm ezy BasiX 2 435 mm 410 mm	Max User Weight Seat Width Self Propelled Whe Max User Weight Seat Width	125 kg 460 mm elchair 3 125 kg 355 mm	Product Weight SOmm Silver Product Weight			

074100-011 Breezy BasiX 2 Self Propelled Wheelchair 460mm Silver

Backrest Height	435 mm	Max User Weight	125 kg	Product Weight	14.9 kg	Max Seat Depth	460 mm
Min Seat Depth	410 mm	Seat Width	460 mm				


074100-012 Breezy BasiX 2 Self Propelled Wheelchair 510mm Silver

Backrest Height	435 mm	Max User Weight	125 kg	Product Weight	14.9 kg	Max Seat Depth	460 mm
Min Seat Depth	410 mm	Seat Width	510 mm				

Quickie QXi Lightweight Folding Wheelchair



• The Quickie QXi has a durable, low-maintenance frame · Simple serviceability, adjustments are quick and intuitive



EIF61616 Quickie QXi Folding Wheelchair 16" x 16"

Max Backrest Height	510 mm	Max User Weight	120 kg	Min Backrest Height	330 mm	Product Weight	13.2 kg
Seat Depth	405 mm	Seat Width	405 mm				

EIF61816 Quickie QXi Folding Wheelchair 18" x 16"

Max Backrest Height	510 mm	Max User Weight	120 kg	Min Backrest Height	330 mm	Product Weight	13.2 kg
Seat Depth	405 mm	Seat Width	455 mm				

EIF61818 Quickie QXi Folding Wheelchair 18" x 18"

Max Backrest Height	510 mm	Max User Weight	120 kg	Min Backrest Height	330 mm	Product Weight	13.2 kg
Seat Depth	455 mm	Seat Width	455 mm				

EIF62018 Quickie QXi Folding Wheelchair 20" x 18"

Max Backrest Height	510 mm	Max User Weight	120 kg	Min Backrest Height	330 mm	Product Weight	13.2 kg
Seat Depth	455 mm	Seat Width	505 mm				

EIF62020 Quickie QXi Folding Wheelchair 20" x 20"

Max Backrest Height	510 mm	Max User Weight	120 kg	Min Backrest Height	330 mm	Product Weight	13.2 kg
Seat Depth	505 mm	Seat Width	505 mm				

Wheelchair

Quickie QRi Rigid Wheelchair



[•] The Quickie QRi is an ultra lightweight rigid wheelchair • With a strong open frame design perfect for your on-the-go lifestyle

[•] Designed for form and function • A hybrid of durability and low-maintenance style at an unbelievable value



EIR51414 Quickie Q	Ri Rigid \	Wheelchair 14" x 14"		
Max User Weight	120 kg	Product Weight	10 kg	

10



• Perfect for the person looking for simplicity, longevity, and Quickie's most popular options at the right price

Min Backrest Height	330 mm	Product Weight	13.2 kg

3.07 Aluminium Self Propelling Rigid Lightweight

Seat Depth	355 mm	Seat Width	355 mm

EIR51616 Quickie QRi Rigid Wheelchair 16" x 16"									
Max User Weight	120 kg	Product Weight	10 kg	Seat Depth	405 mm	Seat Width	405 mm		
EIR51816 Quickie QRi Rigid Wheelchair 18" x 16"									
Max User Weight	120 kg	Product Weight	10 kg	Seat Depth	405 mm	Seat Width	455 mm		
EIR51818 Quickie QI	Ri Rigid \	Wheelchair 18" x 18"							
Max User Weight	120 kg	Product Weight	10 kg	Seat Depth	455 mm	Seat Width	455 mm		
EIR52020 Quickie QRi Rigid Wheelchair 20" x 20"									
Aax User Weight 120 kg Product Weight 10 kg Seat Depth 505 mm Seat Width 505 mm									

Wheelchair

3.13 Self Propelling Amputee Setting

Karma Eagle Heavy Duty Self Propelled Wheelchair



Mobility

- Lightweight, durable, folding and reliable manual self propelled wheelchair
- Slung and padded seat and backrest Swing away armrest
- Swing away removable footrests
- Adjustable rear-wheel and front castor positioning
- Available with fixed or folding backrest which is height adjustable
- Quick release wheels
- Different sizes and accessories available

KM80200-20 Karma Eagle Heavy Duty Self Propelled Wheelchair 20" x 18"

•	3	3 3 1					
Folded Width	415 mm	Max Backrest Height	475 mm	Max User Weight	160 kg	Min Backrest Height	400 mm
Overall Height	930 mm	Overall Length	1080 mm	Overall Width	725 mm	Seat Depth	455 mm
Seat Width	505 mm						

KM8020Q-22 Karma Eagle Heavy Duty Self Propelled Wheelchair 22" x 18"

Folded Width	425 mm	Max Backrest Height	475 mm	Max User Weight	160 kg	Min Backrest Height	400 mm
Overall Height	930 mm	Overall Length	1080 mm	Overall Width	770 mm	Seat Depth	455 mm
Seat Width	560 mm						

Karma Flexx Heavy Duty Self Propelled Folding Wheelchair



 Flexible, robust and stylish
 Sturdy double crossbrace
Adjustable center of gravity for greater stability
 How we during the provident place alignets in

- · Heavy duty functions to assist plus-size clients in need of a mobility device
- Designed with key adjustments to create a suitable sitting environment and fit individual needs
- Folding lightweight frame • Removable, height adjustable, swing in/out footrests
- Angle adjustable backrest
- · Swing away, height adjustable arm rests with ergonomic pad

KM8022-20 Karma Flexx HD Self Propelled Folding Wheelchair 20"

Floor To Seat Height	520 mm	Folded Width	350 mm	Max Backrest Height	530 mm	Max User Weight	170 kg
Min Backrest Height	380 mm	Overall Height	930 mm	Overall Length	1130 mm	Product Weight	17.2 kg
Seat Depth	460 mm	Seat Width	505 mm				

KM8022-22 Karma Flexx HD Self Propelled Folding Wheelchair 22"											
Floor To Seat Height	520 mm	Folded Width	350 mm	Max Backrest Height	530 mm	Max User Weight	170 kg				
Min Backrest Height	380 mm	Overall Height	930 mm	Overall Length	1130 mm	Product Weight	18.3 kg				
Seat Depth	460 mm	Seat Width	560 mm								

KM8022-24 Karma Flexx HD Self Propelled Folding Wheelchair 24"

			-				
Floor To Seat Height	520 mm	Folded Width	350 mm	Max Backrest Height	530 mm	Max User Weight	170 kg
Min Backrest Height	380 mm	Overall Height	930 mm	Overall Length	1130 mm	Product Weight	19.4 kg
Seat Depth	460 mm	Seat Width	610 mm				

Wheelchair

• Uses Karma's S-Ergo patented seating system

- Promotes pressure relief • Helps prevent the user slipping off the seat
 - Adjusts and stabilises seating pose
 - · Double cross braces provide strength and durability
- Folding, lightweight frame
- Height adjustable foot rest
- Fixed ergonomic arm rest pad

KM2501F14-16 Karma Ergo Lite Transit Wheelchair 16"

Karma Ergo Lite Transit Wheelchair

Backrest Height	420 mm	Floor To Seat Height	480 mm	Folded Width	275 mm	Max User Weight	100 kg
Overall Height	910 mm	Overall Length	940 mm	Overall Width	570 mm	Product Weight	8.1 kg
Seat Depth	410 mm	Seat Width	405 mm				

KM2501F14-18 Karma Ergo Lite Transit Wheelchair 18"

Backrest Height	420 mm	Floor To Seat Height	480 mm	Folded Width	275 mm	Max User Weight	100 kg
Overall Height	910 mm	Overall Length	940 mm	Overall Width	620 mm	Product Weight	8.3 kg
Seat Depth	410 mm	Seat Width	455 mm				

Karma Ergo Lite Deluxe Transit Wheelchair



• Helps prevent the user slipping off the seat • Adjusts and stabilises seating pose Double cross braces provide strength and durability Aegis microbe shield permanently bonded antimicrobial barrier Comfortable and durable upholstery

Promotes pressure relief

KM2512F14-16 Karma Ergo Lite Deluxe Transit Wheelchair

Backrest Height	420 mm	Floor To Seat Height	480 mm	Folded Width	270 mm	Max User Weight	100 kg
Overall Height	910 mm	Overall Length	780 mm	Overall Width	620 mm	Product Weight	8.7 kg
Seat Depth	405 mm	Seat Width	405 mm				



3.14 Aluminium Transit Folding Ultra Lightweight

• Uses Karma's S-Ergo patented seating system

-	1	C	п
	I	ο	

KM2512F14-18 Karma Ergo Lite Deluxe Transit Wheelchair 18" 420 mm Floor To Seat Height 480 mm 310 mm Backrest Height Folded Width Max User Weight 100 kg 910 mm **Overall Length** 780 mm **Overall Width** 650 mm 8.8 kg **Overall Height** Product Weight Seat Depth 405 mm Seat Width 455 mm

Wheelchair

3.18 Aluminium Transit Folding

Karma Transit 2 Heavy Duty Wheelchair



- Folding lightweight, aluminium frame Ideal for those who travel with assistance
- Aircraft grade materials provide strength
- Detachable and height adjustable armrests and footrests Foldable backrest allows for ease of transfers
- Pull to lock brakes

KM8520F12-20 Karma Transit 2 Heavy Duty Wheelchair 20"

Floor To Seat Height	550 mm	Max Backrest Height	475 mm	Max User Weight	160 kg	Min Backrest Height	400 mm		
Overall Height	920 mm	Overall Length	1050 mm	Overall Width	690 mm	Product Weight	15.2 kg		
Seat Depth	455 mm	Seat Width	505 mm						

KM8520F12-22 Karma Transit 2 Heavy Duty Wheelchair 22"

Floor To Seat Height	550 mm	Max Backrest Height	475 mm	Max User Weight	160 kg	Min Backrest Height	400 mm
Overall Height	920 mm	Overall Length	1050 mm	Overall Width	740 mm	Product Weight	15.8 kg
Seat Depth	455 mm	Seat Width	560 mm				

Wheelchair

3.23 Self-Propelling Tilt In Space

Karma VIP Tilt in Space Self Propelled Wheelchair



- Lightweight, compact, folding tilt in space wheelchair
- Designed for excellent pressure distribution Uses exclusive perfectly-matched gas struts for an intelligent centre-of-gravity weight shifting seat
- Provides wheelchair users a wide choice of tilt-in-space angles, 90-160°
- Front wheels are 7" solid while the rear wheels are 14" solid
- Intelligent 0 35 degree centre-of-gravity weight shifting tilt -in-space seat, shifts pressure to reduce bedsores
- Aegis microbe shield protects from odour, staining and deterioration from bacteria or fungus

KMVIP515-16 Karma VIP Tilt in Space Self Propelled Wheelchair 16"

Max Backrest Height	880 mm	Max User Weight	115 kg	Min Backrest Height	460 mm	Overall Height	1400 mm
Overall Length	1070 mm	Overall Width	630 mm	Product Weight	17.7 kg	Seat Depth	410 mm
Seat Width	410 mm						

KMVIP515-18 Karma VIP Tilt in Space Self Propelled Wheelchair 18"											
Max Backrest Height	880 mm	Max User Weight	115 kg	Min Backrest Height	460 mm	Overall Height	1400 mm				
Overall Length	1070 mm	Overall Width	680 mm	Product Weight	17.9 kg	Seat Depth	410 mm				
Seat Width	460 mm										

Wheelchair

Karma VIP 515 Tilt in Space Transit Wheelchair



- Wide choice of tilt-in-space angles Unprecedented pressure releasing effects Reduces the risk of pressure sores • Swing away removable footrests • Swing away, height adjustable ergonomic armrest pad
- Push to lock brakes
- Removeable headrest

KMVIP515T-16 Karma VIP 515 Tilt in Space Transit Wheelchair 16"

Floor To Seat Height	500 mm	Folded Width	310 mm	Max Backrest Height	880 mm	Max User Weight	115 kg
Min Backrest Height	640 mm	Overall Height	140 mm	Overall Length	1070 mm	Overall Width	630 mm
Product Weight	17.5 kg	Seat Depth	430 mm	Seat Width	405 mm		

KMVIP515T-18 Karma VIP 515 Tilt in Space Transit Wheelchair 18"

Floor To Seat Height	500 mm	Folded Width	340 mm	Max Backrest Height	880 mm	Max User Weight	115 kg
Min Backrest Height	640 mm	Overall Height	140 mm	Overall Length	1070 mm	Overall Width	680 mm
Product Weight	17.5 kg	Seat Depth	455 mm	Seat Width	455 mm		

Wheelchair

Quickie SR45 Tilt In Space Wheelchair



- Proven durability and functionality Pre-configured for ease of ordering Foot tilt and composite foot plates with transit tie downs • Dual post height adjustable flip back with full length arm pads
- 70 degrees swing in/out hangers, solid seat pan and rear anti tips

EIZ9-1414 Quickie SR45 Tilt In Space Wheelchair 14" x 14" Self Propelled

Max Backrest Height	610 mm	Max User Weight	120 kg	Min Backrest Height	380 mm	Overall Height	1372 mm
Overall Length	1195 mm	Overall Width	1015 mm	Product Weight	18 kg	Seat Depth	355 mm
Seat Width	355 mm						

14



3.24 Transit Tilt In Space

Aluminium lightweight folding wheelchair

Perfectly matched gas struts provide intelligent centre-of-gravity weight shifting

3.25 Tilt In Space Paediatric

• A well appointed economical option for tilt in space

It features 45 degrees tilt in space range with rotation-in-space technology

- Comes standard with non folding back canes, height adjustable back and stroller handles

EIZ9-1414T Quickie SR45 Tilt In Space Wheelchair 14" x 14" Transit

Max Backrest Height	610 mm	Max User Weight	120 kg	Min Backrest Height	380 mm	Overall Height	1372 mm
Overall Length	1195 mm	Overall Width	1015 mm	Product Weight	18 kg	Seat Depth	355 mm
Seat Width	355 mm						

EIZ9-1616T Quickie SR45 Tilt In Space Wheelchair 16" x 16" Transit

Max Backrest Height	610 mm	Max User Weight	120 kg	Min Backrest Height	380 mm	Overall Height	1372 mm
Overall Length	1195 mm	Overall Width	1015 mm	Product Weight	18 kg	Seat Depth	405 mm
Seat Width	405 mm						

EIZ9-1616 Quickie SR45 Tilt In Space Wheelchair 16" x 16" Self Propelled

Max Backrest Heigh	1	610 mm	Max User Weight	120 kg	Min Backrest Height	380 mm	Overall Height	1372 mm
Overall Length		1195 mm	Overall Width	1015 mm	Product Weight	18 kg	Seat Depth	405 mm
Seat Width		405 mm						

EIZ9-1818 Quickie SR45 Tilt In Space Wheelchair 18" x 18" Self Propelled

Max Backrest Height	610 mm	Max User Weight	120 kg	Min Backrest Height	380 mm	Overall Height	1372 mm
Overall Length	1195 mm	Overall Width	1015 mm	Product Weight	18 kg	Seat Depth	455 mm
Seat Width	455 mm						

EIZ9-2018 Quickie SR45 Tilt In Space Wheelchair 20" x 18" Self Propelled

		· · · · · · · · · · · · · · · · · · ·					
Max Backrest Height	610 mm	Max User Weight	120 kg	Min Backrest Height	380 mm	Overall Height	1372 mm
Overall Length	1195 mm	Overall Width	1015 mm	Product Weight	18 kg	Seat Depth	455 mm
Seat Width	505 mm						

Wheelchair

3.26 Manual Wheelchair Accessories

Body Point 2-Point Padded Push Button Hip Belt Center Pull

- Unique Bodypoint design
- · Ffeatures abound in this simple-to-use belt with 2 points of attachment
- Designed with user comfort in mind
- A seamless transition from the creep-resistant webbing to the contoured pads, with a rounded edge binding where they contact the user
- Knowing the challenge of keeping belts clean, materials are chosen and sewing techniques are used that resist trapping dirt

HB235XS-A1 Body F	oint 2-P	oint Padded Push B	utton Hi	p Belt Center Pull XS	mall wit	h Flat Mount	
Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
Package Type	Single						

HB235XS-A2 Body Point 2-Point Padded Push Button Hip Belt Center Pull XSmall with Cinch Mount

Max Seated Hip Width	1	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
Package Type		Single						

Body Point 2-Point Side Release Padded Hip Belt Center Pull



Designed with user comfort in mind

Unique Bodypoint design

- contact the user

ax Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
ackage Type	Single						
3215XS-A1 Body	y Point 2-	Point Padded Side	Release	Center Pull Hip Be	elt XSmall wi	th Flat Mount	
lax Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
ackage Type	Single						
B215S-A1 Body	Point 2-P	oint Padded Side R	lelease C	enter Pull Hip Bel	t Small with	Flat Mount	
lax Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
ackage Type	Single						
-		oint Padded Side F		-			
lax Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
ackage Type	Single						
B215M-A1 Body	Point 2-P	oint Padded Side F	Release (Center Pull Hip Be	lt Medium w	ith Flat Mount	
lax Seated Hip Width	380	Min Seated Hip Width	000	A H H			
iax Seated hip width	300	Mill Seated hip width	230	Overall Length	1520 mm	Package Quantity	1
· ·	Single	min Seated nip Width	230	Uverali Length	1520 mm	Package Quantity	
ackage Type	Single	Point Padded Side					
ackage Type B215M-A2 Body	Single						
hackage Type B215M-A2 Body Aax Seated Hip Width	Single	Point Padded Side	Release	Center Pull Hip Be	elt Medium w	ith Cinch Mount	
ackage Type B215M-A2 Body Iax Seated Hip Width Iackage Type	Single / Point 2-F 380 Single	Point Padded Side	Release	Center Pull Hip Be	elt Medium w	rith Cinch Mount Package Quantity	
B215M-A2 Body Iax Seated Hip Width Iackage Type B215L-A2 Body	Single / Point 2-F 380 Single	Point Padded Side	Release	Center Pull Hip Be	elt Medium w	rith Cinch Mount Package Quantity	
ackage Type B215M-A2 Body Iax Seated Hip Width ackage Type B215L-A2 Body Iax Seated Hip Width	/ Point 2-F 380 Single Point 2-P	Point Padded Side	Release C	Center Pull Hip Be	elt Medium w 1520 mm t Large with	rith Cinch Mount Package Quantity Cinch Mount	1
Package Type B215M-A2 Body Max Seated Hip Width Package Type B215L-A2 Body Max Seated Hip Width Package Type	/ Point 2-F 380 Single Point 2-P 480 Single	Point Padded Side	Release (230 Release C 330	Center Pull Hip Be	elt Medium w 1520 mm t Large with 1650 mm	rith Cinch Mount Package Quantity Cinch Mount Package Quantity	1
ackage Type B215M-A2 Body Iax Seated Hip Width ackage Type B215L-A2 Body Iax Seated Hip Width ackage Type B215L-A1 Body	/ Point 2-F 380 Single Point 2-P 480 Single	Point Padded Side Min Seated Hip Width oint Padded Side R Min Seated Hip Width	Release (230 Release C 330	Center Pull Hip Be	elt Medium w 1520 mm t Large with 1650 mm	rith Cinch Mount Package Quantity Cinch Mount Package Quantity	1
ackage Type B215M-A2 Body Iax Seated Hip Width ackage Type B215L-A2 Body Iax Seated Hip Width ackage Type B215L-A1 Body Iax Seated Hip Width	/ Point 2-F 380 Single Point 2-P 480 Single Point 2-Pc	Point Padded Side Min Seated Hip Width oint Padded Side R Min Seated Hip Width Dint Padded Side R	Release C 230 Release C 330 elease C	Center Pull Hip Be Overall Length enter Pull Hip Bel Overall Length	elt Medium w 1520 mm It Large with 1650 mm It Large with	rith Cinch Mount Package Quantity Cinch Mount Package Quantity Flat Mount	
ackage Type B215M-A2 Body Iax Seated Hip Width ackage Type B215L-A2 Body Iax Seated Hip Width ackage Type B215L-A1 Body Iax Seated Hip Width ackage Type	Single / Point 2-F 380 Single Point 2-P 480 Single Point 2-P 480 Single Point 2-Pin 480 Single	Point Padded Side Min Seated Hip Width oint Padded Side R Min Seated Hip Width Dint Padded Side R Min Seated Hip Width	Release C 330 elease C 330	Center Pull Hip Bei Overall Length enter Pull Hip Bei Overall Length enter Pull Hip Bei Overall Length	elt Medium w 1520 mm It Large with 1650 mm It Large with	rith Cinch Mount Package Quantity Cinch Mount Package Quantity Flat Mount	
Aackage Type B215M-A2 Body Aax Seated Hip Width Package Type B215L-A2 Body Aax Seated Hip Width Package Type B215L-A1 Body Aax Seated Hip Width Package Type	Single / Point 2-F 380 Single Point 2-P 480 Single Point 2-P 480 Single Point 2-Pin 480 Single	Point Padded Side Min Seated Hip Width oint Padded Side R Min Seated Hip Width Dint Padded Side R	Release C 330 elease C 330	Center Pull Hip Bei Overall Length enter Pull Hip Bei Overall Length enter Pull Hip Bei Overall Length	elt Medium w 1520 mm It Large with 1650 mm It Large with	rith Cinch Mount Package Quantity Cinch Mount Package Quantity Flat Mount	



contact the user



· Ffeatures abound in this simple-to-use belt with 2 points of attachment

• A seamless transition from the creep-resistant webbing to the contoured pads, with a rounded edge binding where they

• Knowing the challenge of keeping belts clean, materials are chosen and sewing techniques are used that resist trapping dirt

• A seamless transition from the creep-resistant webbing to the contoured pads, with a rounded edge binding where they

• Knowing the challenge of keeping belts clean, materials are chosen and sewing techniques are used that resist trapping dirt

Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
Package Type	Single			-			
						- -	
B219XS-A2 Body	y Point 2-I	Point Side Release	Padded I	Hip Belt Rear Pull	XSmall with	Cinch Mount	
Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
	Single	pint Side Delegee D		n Polt Door Dull S	moll with El	at Mount	
B219S-A1 Body		pint Side Release P	Padded Hi	p Belt Rear Pull S Overall Length	mall with Fla	at Mount Package Quantity	1
IB219S-A1 Body	Point 2-Po						1
B219S-A1 Body	Point 2-Pc						1
B219S-A1 Body Max Seated Hip Width Package Type	Point 2-Po		180	Overall Length	1400 mm	Package Quantity	1
Max Seated Hip Width Package Type	Point 2-Po	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1

HB219M-A1 Body Point 2-Point Side Release Padded Hip Belt Rear Pull Medium with Flat Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB219M-A2 Body Point 2-Point Side Release Padded Hip Belt Rear Pull Medium with Cinch Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

Body Point 2-Point Padded Push Button Hip Belt



- Unique Bodypoint design · Ffeatures abound in this simple-to-use belt with 2 points of attachment
- Designed with user comfort in mind
- A seamless transition from the creep-resistant webbing to the contoured pads, with a rounded edge binding where they contact the user
- Knowing the challenge of keeping belts clean, materials are chosen and sewing techniques are used that resist trapping dirt

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						
B205S-A2 Body	y Point 2-P	oint Padded Push E	utton Hi				
	000						
Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Max Seated Hip Width Package Type	280 Single	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single	Min Seated Hip Width Point Padded Push E Min Seated Hip Width					1

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB205L-A1 Body Point 2-Point Padded Push Button Hip Belt Large with Flat Mount

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

HB205L-A2 Body Po	oint 2-Po	int Padded Push Bu	tton Hip	В
Max Seated Hip Width	480	Min Seated Hip Width	330	0
Package Type	Single			
				_

Body Point 2-Point Rehab Latch Padded Hip Belts



	orinque be	a point	~	.0.9.		
•	Ffeatures	abound	in	this	simple	e-1

- Designed with user comfort in mind
- contact the user

HB225	S-A1 Body Po	oint 2-Po	int Rehab Latch Pac	lded Hip	Belts Small with Fla	at Mount		
Max Seate	ted Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package 1	Туре	Single						

HB225S-A2 Body Point 2-Point Rehab Latch Padded Hip Belts Small with Cinch Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

HB225M-A1 Body Point 2-Point Rehab Latch Padded Hip Belts Medium with Flat Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB225M-A2 Body Point 2-Point Rehab Latch Padded Hip Belts Medium with Cinch Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB225L-A1	Body Point 2-Point Rehab Latch Padded Hip B
-----------	---

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

HB225L-A2 Body Point 2-Point Rehab Latch Padded Hip Belts Large with Cinch Mount

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

Body Point 4-Point Side Release Hip Belt, Hook and Loop Pad Dual Pull



aligned, stable position and stays securely against the seat and back supports

- anterior pelvic tilt, pelvic rotation and thrusting

LB416XS-A1 Body Po	oint 4-Po	oint Side Release Hip	o Belt, Ho	ook and Loop Pad Di	ual Pull)	Small with Flat Mou	int
Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
Package Type	Single						



Belt Large with Cinch Mount

Overall Length	1650 mm	Package Quantity	1

e-to-use belt with 2 points of attachment

• A seamless transition from the creep-resistant webbing to the contoured pads, with a rounded edge binding where they

• Knowing the challenge of keeping belts clean, materials are chosen and sewing techniques are used that resist trapping dirt

Belts Large with Flat Mount

• 4-Point Padded Hip Belt Keeping the belt firmly in position

• The 4-Point Padded Hip Belt attaches to the wheelchair in four places to assure that the pelvis does not move from an

• The primary straps position the belt like the 2-point belt and the secondary straps anchor it into position

These extra points prevent the belt from riding up into the abdomen or twisting and moving out of position making it ideal for

LB416XS-A2 Body Point 4-Point Side Release Hip Belt, Hook and Loop Pad Dual Pull XSmall with Cinch Mount

Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
Package Type	Single						

Body Point 2-Point Push Button Hip Belt, Hook and Loop Pad Center Pull

- Unique Bodypoint design
- Ffeatures abound in this simple-to-use belt with 2 points of attachment
- Designed with user comfort in mind • A seamless transition from the creep-resistant webbing to the contoured pads, with a rounded edge binding where they
- contact the user
- . Knowing the challenge of keeping belts clean, materials are chosen and sewing techniques are used that resist trapping dirt

LB205M-A1 Body Point 2-Point Push Button Hip Belt, Hook and Loop Pad Center Pull Medium with Flat Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

LB205M-A2 Body Point 2-Point Push Button Hip Belt, Hook and Loop Pad Center Pull Medium with Cinch Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

Body Point 2-Point Padded Push Button Hip Belt Center Pull



- Unique Bodypoint design
- · Ffeatures abound in this simple-to-use belt with 2 points of attachment • Designed with user comfort in mind
- A seamless transition from the creep-resistant webbing to the contoured pads, with a rounded edge binding where they contact the user
- Knowing the challenge of keeping belts clean, materials are chosen and sewing techniques are used that resist trapping dirt

1	Package Quantity	1400 mm	Overall Length	180	Min Seated Hip Width	280	Max Seated Hip Width
						Single	Package Type
y	Package Quantit	1400 mm	Overall Length	180	Min Seated Hip Width		

Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
Package Type	Single						

Body Point 2-Point Push Button Padded Hip Belt Rear Pull

- Unique Bodypoint design
- · Ffeatures abound in this simple-to-use belt with 2 points of attachment
- Designed with user comfort in mind
- A seamless transition from the creep-resistant webbing to the contoured pads, with a rounded edge binding where they contact the user
- Knowing the challenge of keeping belts clean, materials are chosen and sewing techniques are used that resist trapping dirt

HB239XS-A1 Body Point 2-Point Push Button Padded Hip Belt Rear Pull XSmall with Flat Mount Max Seated Hip Width 230 Min Seated Hip Width 130 Package Type Single

Body Point 2-Point Padded Rehab Latch Hip Belt Rear Pull

Н



- Designed with user comfort in mind
- contact the user

HB229M-A1 Body Point 2-Point Padded Rehab Latch Hip Belt Rear Pull Medium with Flat Mount							
Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						
HR229S-M1 Rody Point 2-Point Padded Pohab Latch Hip Polt Poar Pull Small with Flat Mount							

HB229S-A1 Body Point 2-Point Padded Rehab Latch Hip Belt Rear Pull Small with Flat Mount Max Seated Hip Width 280 Min Seated Hip Width 180 Package Type Single

IB229M-A2 Body Point 2-Point Padded Rehab Latch Hip Belt Rear Pull Medium with Cinch Mount									
Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1		
Package Type	Single								

HB229S-A2 Body Point 2-Point Padded Rehab Latch Hip Belt Rear Pull Small with Cinch Mount

1	Nax Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
F	Package Type	Single						

Body Point 4-Point Padded Push Button Hip Belt Center Pull



- aligned, stable position and stays securely against the seat and back supports
- anterior pelvic tilt, pelvic rotation and thrusting

ŀ	HB435XS-A1 Body Point 4-Point Padded Push Button Hip Belt Center Pull XSmall with Flat Mount								
	Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1	
	Package Type	Single							
			· · · · · · · · · · · · · · · · · · ·						

HB435XS-A2 Body Point 4-Point Padded Push Button Hip Belt Center Pull XSmall with Cinch Mount								
Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1	
Package Type	Single							



20



Overall Length	1270 mm	Package Quantity	1

Ffeatures abound in this simple-to-use belt with 2 points of attachment

• A seamless transition from the creep-resistant webbing to the contoured pads, with a rounded edge binding where they

. Knowing the challenge of keeping belts clean, materials are chosen and sewing techniques are used that resist trapping dirt

Overall Length	1400 mm	Package Quantity	1

• The 4-Point Padded Hip Belt attaches to the wheelchair in four places to assure that the pelvis does not move from an

. The primary straps position the belt like the 2-point belt and the secondary straps anchor it into position

. These extra points prevent the belt from riding up into the abdomen or twisting and moving out of position making it ideal for

Body Point 2-Point Padded Rehab Latch Hip Belt Center Pull

- Unique Bodypoint design
- · Ffeatures abound in this simple-to-use belt with 2 points of attachment
- Designed with user comfort in mind
- A seamless transition from the creep-resistant webbing to the contoured pads, with a rounded edge binding where they contact the user
- Knowing the challenge of keeping belts clean, materials are chosen and sewing techniques are used that resist trapping dirt

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						
B225M Body Poin	nt 2-Point	t Padded Rehab Lato	ch Hip Be	elt Center Pull Mediu	ım Pivot	Mount	
SB225M Body Poin Max Seated Hip Width	nt 2-Point 380	t Padded Rehab Lato	ch Hip Be	elt Center Pull Mediu	IM Pivot	Mount Package Quantity	1

Package Type

SB225XL Body Point 2-Point Padded Rehab Latch Hip Belt Center Pull XLarge Pivot Mount

Package Quantity

Mobility

Body Point 2-Point Rehab Latch Hip Belt, Hook and Loop Pad center Pull



- Unique Bodypoint design · Ffeatures abound in this simple-to-use belt with 2 points of attachment
- Designed with user comfort in mind
- A seamless transition from the creep-resistant webbing to the contoured pads, with a rounded edge binding where they contact the user
- . Knowing the challenge of keeping belts clean, materials are chosen and sewing techniques are used that resist trapping dirt

Single

LB225L-A1 Body Point 2-Point Rehab Latch Hip Belt, Hook and Loop Pad center Pull Large with Flat Mount

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

LB225M-A1 Body Point 2-Point Rehab Latch Hip Belt, Hook and Loop Pad center Pull Medium with Flat Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

LB225L-A2 Body Point 2-Point Rehab Latch Hip Belt, Hook and Loop Pad center Pull Large with Cinch Mount

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

LB225M-A2 Body Point 2-Point Rehab Latch Hip Belt, Hook and Loop Pad center Pull Medium with Cinch Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

Body Point Aeromesh Rapid Dry Bath Belt Two Piece



- the bath belt release knob
- Padding at the middle of the belt allows for greater user comfort
- Helps keep the user in the chair during bathing

Package Quantity 1 BB101L-1 Body Point Aeromesh Rapid Dry Ba	Package Type	Single
BB101L-1 Body Point Aeromesh Rapid Dry Ba	ath Polt Two Diago Lorgo	
· · · · · · · · · · · · · · · · · · ·		
Package Quantity 1	Package Type	Single

Body Point 4-Point Padded Compact Push Button Hip Belt Center Pull



- The 4-Point Padded Hip Belt attaches to the wheelchair in four places to assure that the pelvis does not move from an
- aligned, stable position and stays securely against the seat and back supports
- The primary straps position the belt like the 2-point belt and the secondary straps anchor it into position
- anterior pelvic tilt, pelvic rotation and thrusting
- . These extra points prevent the belt from riding up into the abdomen or twisting and moving out of position making it ideal for

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Max Seated hip width	400	will Seated hip width	330	overall celigui	1030 11111	Fackage Quantity	1
Package Type	Single						
B445M-A1 Body F	oint 4-P	pint Padded Compac	ct Push I	Button Hip Belt Cent	er Pull M	ledium with Flat Mor	unt
Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						
IB445S-A1 Body P	oint 4-Pc	int Padded Compac	t Push B	Button Hip Belt Cent	er Pull Si	nall with Flat Mount	
	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Max Seated Hip Width	200	with ocated hip width					

HB445L-A2 Body Point 4-Point Padded Compact Push Button Hip Belt Center Pull Large with Cinch Mount

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

HB445M-A2 Body Point 4-Point Padded Compact F	Push	B
---	------	---

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB445S-A2 Body Point 4-Point Padded Compact Push Button Hip Belt Center Pull Small with Cinch Mount

[Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
	Package Type	Single						



• The Aeromesh Rapid Dry Bath Belt has the great design and features you're used to . The original two piece bath belt has grommets to enable the belt to be mounted directly to the bath chair using a screw or

• 4-Point Padded Hip Belt Keeping the belt firmly in position

HB4451-A1 Body Point 4-Point Padded Compact Push Button Hin Belt Center Pull Large with Flat Mount

utton Hip Belt Center Pull Medium with Cinch Mount

Body Point 4-Point Padded Compact Push Button Hip Belt Rear Pull



- 4-Point Padded Hip Belt Keeping the belt firmly in position
- The 4-Point Padded Hip Belt attaches to the wheelchair in four places to assure that the pelvis does not move from an aligned, stable position and stays securely against the seat and back supports
- The primary straps position the belt like the 2-point belt and the secondary straps anchor it into position
- These extra points prevent the belt from riding up into the abdomen or twisting and moving out of position making it ideal for anterior pelvic tilt, pelvic rotation and thrusting



Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
ackage Type	Single						
B216XS-A2 Bod	y Point 2-	Point Padded Pusl	h Button I	Dual Pull Hip Belt	XSmall with	Cinch Mount	
Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
Package Type	Single						
B216XS-A1 Body	y Point 2-	Point Padded Push	n Button [)ual Pull Hip Belt	XSmall with	Flat Mount	
Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
Package Type	Single						
B216S-A1 Body	Point 2-P	oint Padded Push	Button Di	ual Pull Hip Belt S	small with Fla	it Mount	
Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						
B216M-A1 Body	Point 2-F	Point Padded Push	Button D	ual Pull Hip Belt N	Medium with	Flat Mount	
Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						
B216M-A2 Body	/ Point 2-F	Point Padded Push	Button D	Overall Length	Medium with	Cinch Mount	1
Package Type	Single	win seated np width	230		1320 11111	Fackage quantity	-
B216L-A1 Body I	Point 2-Po	oint Padded Push I	Button Du	al Pull Hip Belt La	arge with Fla	t Mount	
Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						
	Point 2-P	oint Padded Push	Button Di	ual Pull Hip Belt L	arge with Cir	nch Mount	
B216L-A2 Body	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
-							
				-	-		1
1B216L-A2 Body Max Seated Hip Width Package Type	Single						
Max Seated Hip Width Package Type		Rapid Dry Bath	Belt On	e Piece			

HB449M-A1 Body P	oint 4-Po	pint Padded Compac	ct Push I	Button Hip Belt Rear	Pull Med	dium with Flat Moun	it
Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB449S-A1 Body Point 4-Point Padded Compact Push Button Hip Belt Rear Pull Small with Flat Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

HB449M-A2 Body Point 4-Point Padded Compact Push Button Hip Belt Rear Pull Medium with Cinch Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB449S-A2 Body Point 4-Point Padded Compact Push Button Hip Belt Rear Pull Small with Cinch Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

Body Point 4-Point Rehab Latch Hip Belt, Hook and Loop Pad Rear Pull

• 4-Point Padded Hip Belt Keeping the belt firmly in position

• The 4-Point Padded Hip Belt attaches to the wheelchair in four places to assure that the pelvis does not move from an aligned, stable position and stays securely against the seat and back supports

• The primary straps position the belt like the 2-point belt and the secondary straps anchor it into position



. These extra points prevent the belt from riding up into the abdomen or twisting and moving out of position making it ideal for anterior pelvic tilt, pelvic rotation and thrusting

LB429M-A1 Body Po	oint 4-Po	int Rehab Latch Hip	Belt, Ho	ok and Loop Pad Re	ar Pull M	edium with Flat Mo	unt
Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

LB429M-A2 Body Point 4-Point Rehab Latch Hip Belt, Hook and Loop Pad Rear Pull Medium with Cinch Mount

I	Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
I	Package Type	Single						



Body Point 2-Point Padded Push Button Hip Belt Dual Pull

• Ffeatures abound in this simple-to-use belt with 2 points of attachment

• A seamless transition from the creep-resistant webbing to the contoured pads, with a rounded edge binding where they

• Knowing the challenge of keeping belts clean, materials are chosen and sewing techniques are used that resist trapping dirt

National Distribution, Local Service

Package Quantity 1 Package Type Single											
BB112L-1 Body Point Aeromesh Rapid Dry Bath Belt One Piece Large with Plastic Slides											
Package Quantity	1	Package Type	Single								
Package quality			onigic								
	Aeromesh Rapid Dry Ba	th Belt One Piece XLarge	Single								
BillixL-1 Body Point / Package Quantity	1	th Belt One Piece XLarge	Single								

Body Point 2-Point Padded Push Button Hip Belt Rear Pull



- Unique Bodypoint design
- Ffeatures abound in this simple-to-use belt with 2 points of attachment
- Designed with user comfort in mind A seamless transition from the creep-resistant webbing to the contoured pads, with a rounded edge binding where they contact the user
- Knowing the challenge of keeping belts clean, materials are chosen and sewing techniques are used that resist trapping dirt



HB239XS-A2 Body I	Point 2-F	Point Push Button Pa	added Hi	p Belt Rear Pull XSm	all with	Cinch Mount	
Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
Package Type	Single						

HB209M-A1 Body Point 2-Point Padded Push Button Hip Belt Rear Pull Medium with Flat Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB209S-A1 Body Point 2-Point Padded Push Button Hip Belt Rear Pull Small with Flat Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

HB209M-A2 Body Point 2-Point Padded Push Button Hip Belt Rear Pull Medium with Cinch Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB209S-A2 Body Point 2-Point Padded Push Button Hip Belt Rear Pull Small with Cinch Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

Body Point 2-Point Padded Push Button Hip Belt Rear Pull

- Unique Bodypoint design
- Ffeatures abound in this simple-to-use belt with 2 points of attachment
- Designed with user comfort in mind
- A seamless transition from the creep-resistant webbing to the contoured pads, with a rounded edge binding where they contact the user
- Knowing the challenge of keeping belts clean, materials are chosen and sewing techniques are used that resist trapping dirt

HB249M-A1 Body Point 2-Point Padded Push Button Hip Belt Rear Pull Medium with Flat Mount

		•						
Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1	
Package Type	Single							

HB249S-A1 Body Point 2-Point Padded Push Button Hip Belt Rear Pull Small with Flat Mount

-			•				
Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

HB249M-A2 Body Point 2-Point Padded Push Button Hip Belt Rear Pull Medium with Cinch Mount

Max Seate	ed Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package T	Гуре	Single						

HB249S-A2 Body Point 2-Point Padded Push Button Hip Belt Rear Pull Small with Cinch Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

Body Point 2-Point Padded Push Button Hip Belt Center Pull





- contact the user

HB245L-A1 Body Point 2-Point Padded Push Button Hip Belt Center Pull Large with Flat Mount						lat Mount		
	Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
	Package Type	Single						

HB245M-A1 Body Point 2-Point Padded Push Button Hip Belt Center Pull Medium with Flat Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB245S-A1 Body Point 2-Point Padded Push Button Hip Belt Center Pull Small with Flat Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

HB245L-A2 Body Point 2-Point Padded Push Button Hip Belt Center Pull Large with Cinch Mount

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

HB245M-A2 Body Point 2-Point Padded Push Button Hip Belt Center Pull Medium with Cinch Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB245S-A2 Body Point 2-Point Padded Push Button Hip Belt Center Pull Small with Cinch Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						




Country Care GROUF

• A seamless transition from the creep-resistant webbing to the contoured pads, with a rounded edge binding where they

• Knowing the challenge of keeping belts clean, materials are chosen and sewing techniques are used that resist trapping dirt

Body Point 2-Point Rehab Latch Hip Belt, Hook and Loop Pad Rear Pull

- Unique Bodypoint design
- · Ffeatures abound in this simple-to-use belt with 2 points of attachment
- Designed with user comfort in mind



Knowing the challenge of keeping belts clean, materials are chosen and sewing techniques are used that resist trapping dirt

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

LB229M-A2 Body Point 2-Point Rehab Latch Hip Belt, Hook and Loop Pad Rear Pull Medium with Cinch Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

Body Point 4-Point Push Button Hip Belt, Hook and Loop Pad Center Pull



- The 4-Point Padded Hip Belt attaches to the wheelchair in four places to assure that the pelvis does not move from an aligned, stable position and stays securely against the seat and back supports
- The primary straps position the belt like the 2-point belt and the secondary straps anchor it into position
- These extra points prevent the belt from riding up into the abdomen or twisting and moving out of position making it ideal for anterior pelvic tilt, pelvic rotation and thrusting

LB435XS-A1 Body Point 4-Point Push Button Hip Belt, Hook and Loop Pad Center Pull XSmall with Flat Mount

Package Type Single	Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
	Package Type	Single						


LB435XS-A2 Body Point 4-Point Push Button Hip Belt, Hook and Loop Pad Center Pull XSmall with Cinch Mount

Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
Package Type	Single						

Body Point 4-Point Side Release Padded Hip Belt Center Pull

• 4-Point Padded Hip Belt Keeping the belt firmly in position

- The 4-Point Padded Hip Belt attaches to the wheelchair in four places to assure that the pelvis does not move from an aligned, stable position and stays securely against the seat and back supports
- The primary straps position the belt like the 2-point belt and the secondary straps anchor it into position



. These extra points prevent the belt from riding up into the abdomen or twisting and moving out of position making it ideal for anterior pelvic tilt, pelvic rotation and thrusting

HB415S-A2 Body Point 4-Point Side Release Padded Hip Belt Cer	nter Pull Small with Cinch Mount
---	----------------------------------

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

HB415XS-A2 Body Point 4-Point Side Release Padded Hip Belt Center Pull XSmall with Cinch Mount							
Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
Package Type	Single						

HB415XS-A1 Body Point 4-Point Side Release Padded Hip Belt Center Pull XSmall with Flat Mount

Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
Package Type	Single						

HB415S-A1 Body Point 4-Point Side Release Padded Hip Belt Center Pull Small with Flat Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

HB415M-A1 Body Point 4-Point Side Release Padded Hip Belt Center Pull Medium with Flat Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB415M-A2 Body Point 4-Point Side Release Padded Hip Belt Center Pull Medium with Cinch Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB415L-A2 Body Point 4-Point Side Release Padded Hip Belt Center Pull Large with Cinch Mount

Max Seated Hip Width 480	0	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type Singl	ngle						

HB415L-A1 Body Point 4-Point Side Release Padded Hip Belt Center Pull Large with Flat Mount

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

Body Point 4-Point Side Release Padded Hip Belt Dual Pull



P

- 4-Point Padded Hip Belt Keeping the belt firmly in position
- aligned, stable position and stays securely against the seat and back supports
- anterior pelvic tilt, pelvic rotation and thrusting

HB416S-A2 Body Point 4-Point Side Release Padded Hip Belt Dual Pull Small with Cinch Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

HB416XS-A2 Body Point 4-Point Side Release Padded Hip Belt Dual Pull XSmall with Cinch Mount

Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
Package Type	Single						

HB416XS-A1 Body Point 4-Point Side R	Release Padded Hip
--------------------------------------	--------------------

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1270 mm	Package Quantity	1			
Package Type	Single									

HB416S-A1	Body Point 4-Point Side Release Padded Hip B
-----------	--

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1400 mm	Package Quantity	1			
Package Type	Single									

Country Care GROUF

• The 4-Point Padded Hip Belt attaches to the wheelchair in four places to assure that the pelvis does not move from an

• The primary straps position the belt like the 2-point belt and the secondary straps anchor it into position

• These extra points prevent the belt from riding up into the abdomen or twisting and moving out of position making it ideal for

Belt Dual Pull XSmall with Flat Mount

Belt Dual Pull Small with Flat Mount

HB416M-A1 Body Point 4-Point Side Release Padded Hip Belt Dual Pull Medium with Flat Mount

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB416M-A2 Body Point 4-Point Side Release Padded Hip Belt Dual Pull Medium with Cinch Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB416L-A1 Body Point 4-Point Side Release Padded Hip Belt Dual Pull Large with Flat Mount

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

HB416L-A2 Body Point 4-Point Side Release Padded Hip Belt Dual Pull Large with Cinch Mount

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

Body Point 4-Point Padded Rehab Latch Hip Belt Center Pull

- 4-Point Padded Hip Belt Keeping the belt firmly in position
- The 4-Point Padded Hip Belt attaches to the wheelchair in four places to assure that the pelvis does not move from an aligned, stable position and stays securely against the seat and back supports
- The primary straps position the belt like the 2-point belt and the secondary straps anchor it into position



• These extra points prevent the belt from riding up into the abdomen or twisting and moving out of position making it ideal for anterior pelvic tilt, pelvic rotation and thrusting

HB425L-A1 Body Point 4-Point Padded Rehab Latch Hip Belt Center Pull Large with Flat Mount

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

HB425M-A1 Body Point 4-Point Padded Rehab Latch Hip Belt Center Pull Medium with Flat Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB425S-A1 Body Point 4-Point Padded Rehab Latch Hip Belt Center Pull Small with Flat Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

HB425L-A2 Body Point 4-Point Padded Rehab Latch Hip Belt Center Pull Large with Cinch Mount

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

HB425M-A2 Body Point 4-Point Padded Rehab Latch Hip Belt Center Pull Medium with Cinch Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB425S-A2 Body Point 4-Point Padded Rehab Latch Hip Belt Center Pull Small with Cinch Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

Body Point 4-Point Padded Push Button Hip Belt Rear Pull



 4-Point Padded Hip Belt Keeping the belt firmly in position aligned, stable position and stays securely against the seat and back supports

anterior pelvic tilt, pelvic rotation and thrusting

HB439XS-A1 Body F	B439XS-A1 Body Point 4-Point Padded Push Button Hip Belt Rear Pull XSmall with Flat Mount										
Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1				
Package Type	Single										

HB439XS-A2 Body Point 4-Point Padded Push Button Hip Belt Rear Pull XSmall with Cinch Mount

Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
Package Type	Single						

Body Point 4-Point Push Button Hip Belt, Hook and Loop Pad Center Pull



aligned, stable position and stays securely against the seat and back supports • The primary straps position the belt like the 2-point belt and the secondary straps anchor it into position

anterior pelvic tilt, pelvic rotation and thrusting



		ner den Batten nip	bent, mot			Large main lat mou	
Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

LB405M-A1 Body Point 4-Point Push Button Hip Belt, Hook and Loop Pad Center Pull Medium with Flat Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

LB405L-A2 Body Point 4-Point Push Button Hip Belt, Hook and Loop Pad Center Pull Large with Cinch Mount

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

LB405M-A2 Body Point 4-Point Push Button Hip Belt, Hook and Loop Pad Center Pull Medium with Cinch Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB405S-A1 Body Point 4-Point Padded Push Button Hip Belt Small with Flat Mount F

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

HB405S-A2 Body Point 4-Point Padded Push Button Hip Belt Small with Cinch Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						



- The 4-Point Padded Hip Belt attaches to the wheelchair in four places to assure that the pelvis does not move from an
- The primary straps position the belt like the 2-point belt and the secondary straps anchor it into position
- These extra points prevent the belt from riding up into the abdomen or twisting and moving out of position making it ideal for

• 4-Point Padded Hip Belt Keeping the belt firmly in position

- The 4-Point Padded Hip Belt attaches to the wheelchair in four places to assure that the pelvis does not move from an
- These extra points prevent the belt from riding up into the abdomen or twisting and moving out of position making it ideal for

LB405L-A1 Body Point 4-Point Push Button Hip Belt. Hook and Loop Pad Center Pull Large with Flat Mount

HE4USIN-AI BODY	Point 4-P	oint Padded Push Bi	utton Hip	b Belt Medium with I	lat Mou	nt					
Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1				
Package Type	Single										
HB405M-A2 Body Point 4-Point Padded Push Button Hip Belt Medium with Cinch Mount											
Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1				
Package Type	Single										
HB405L-A1 Body Point 4-Point Padded Push Button Hip Belt Large with Flat Mount											
May Seated Hin Width	/180	Min Sected Hin Width	330	Overall Length	1650 mm	Package Quantity	1				

HB405L-A2 Body Point 4-Point Padded Push Button Hip Belt Large with Cinch Mount

UDAOEM A1 Dedy Deint 4 Deint Dedded Duch Dutten Lin Delt Medium with Elet Meunt

·			•				
Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

Body Point 4-Point Rehab Latch Hip Belt, Hook and Loop Pad Center Pull

- 4-Point Padded Hip Belt Keeping the belt firmly in position
- The 4-Point Padded Hip Belt attaches to the wheelchair in four places to assure that the pelvis does not move from an aligned, stable position and stays securely against the seat and back supports
- The primary straps position the belt like the 2-point belt and the secondary straps anchor it into position



Package Type

Single

• These extra points prevent the belt from riding up into the abdomen or twisting and moving out of position making it ideal for anterior pelvic tilt, pelvic rotation and thrusting

LB425L-A1 Body Point 4-Point Rehab Latch Hip Belt, Hook and Loop Pad Center Pull Large with Flat Mount Framesaver

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

LB425L-B1 Body Point 4-Point Rehab Latch Hip Belt, Hook and Loop Pad Center Pull Large with Flat Mount

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

LB425M-A1 Body Point 4-Point Rehab Latch Hip Belt, Hook and Loop Pad Center Pull Medium with Flat Mount Framesaver

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

LB425M-B1 Body Point 4-Point Rehab Latch Hip Belt, Hook and Loop Pad Center Pull Medium with Flat Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

LB425L-A2 Body Point 4-Point Rehab Latch Hip Belt, Hook and Loop Pad Center Pull Large with Cinch Mount Framesaver

М	ax Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Pa	ackage Type	Single						

LB425L-B2 Body Point 4-Point Rehab Latch Hip Belt, Hook and Loop Pad Center Pull Large with Cinch Mount

Max Seated Hip Width	480	Min Seated Hip Width	330	Overall Length	1650 mm	Package Quantity	1
Package Type	Single						

LB425M-A2 Body Point 4-Point Rehab Latch Hip Belt, Hook and Loop Pad Center Pull Medium with Cinch Mount Framesaver

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

LB425M-B2 Body Point 4-Point Rehab Latch Hip Belt, Hook and Loop Pad Center Pull Medium with Cinch Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

Body Point 4-Point Side Release Padded Hip Belt Rear Pull



aligned, stable position and stays secu
 The primary straps position the belt like
TI 1 1 1 1 1 1 1 1 1

anterior pelvic tilt, pelvic rotation and thrusting

HB419M-A1 Body Point 4-Point Side Release Padded Hip Belt Rear Pull Medium with Flat Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB419S-A1 Body Point 4-Point Side Release Padded Hip Belt Rear Pull Small with Flat Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

HB419XS-A2 Body Point 4-Point Side Release Padded Hip Belt Rear Pull XSmall with Cinch Mount

Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
Package Type	Single						

HB419M-A2 Body Point 4-Point Side Release Padded Hip Belt Rear Pull Medium with Cinch Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB419S-A2 Body Point 4-Point Side Release Padded Hip Belt Rear Pull Small with Cinch Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

HB419XS-A1 Body Point 4-Point Side Release Padded Hip Belt Rear Pull XSmall with Flat Mount

Max Seated Hip Width	230	Min Seated Hip Width	130	Overall Length	1270 mm	Package Quantity	1
Package Type	Single						

Body Point 4-Point Padded Push Button Hip Belt Rear Pull



P

anterior pelvic tilt, pelvic rotation and thrusting

Mobility



• 4-Point Padded Hip Belt Keeping the belt firmly in position

- The 4-Point Padded Hip Belt attaches to the wheelchair in four places to assure that the pelvis does not move from an
 - urely against the seat and back supports
 - e the 2-point belt and the secondary straps anchor it into position
- . These extra points prevent the belt from riding up into the abdomen or twisting and moving out of position making it ideal for

• 4-Point Padded Hip Belt Keeping the belt firmly in position

- The 4-Point Padded Hip Belt attaches to the wheelchair in four places to assure that the pelvis does not move from an aligned, stable position and stays securely against the seat and back supports
- The primary straps position the belt like the 2-point belt and the secondary straps anchor it into position
- . These extra points prevent the belt from riding up into the abdomen or twisting and moving out of position making it ideal for

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1	
Package Type	Single							
HB409S-A1 Body Point 4-Point Padded Push Button Hip Belt Rear Pull Small Flat Mount								
Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1	
Package Type	Single							
		Point Padded Push B	utton Hi	p Belt Rear Pull Med	ium Cinc	h Mount Package Quantity	1	

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

Body Point 4-Point Padded Rehab Latch Hip Belt Rear Pull

- 4-Point Padded Hip Belt Keeping the belt firmly in position
- The 4-Point Padded Hip Belt attaches to the wheelchair in four places to assure that the pelvis does not move from an aligned, stable position and stays securely against the seat and back supports
- The primary straps position the belt like the 2-point belt and the secondary straps anchor it into position
- These extra points prevent the belt from riding up into the abdomen or twisting and moving out of position making it ideal for anterior pelvic tilt, pelvic rotation and thrusting

HB429M-A1 Body Point 4-Point Padded Rehab Latch Hip Belt Rear Pull Medium with Flat Mount Max Seated Hip Width 380 Min Seated Hip Width **Overall Length** 1520 mm Package Quantity 230 Package Type Single

HB429S-A1 Body Point 4-Point Padded Rehab Latch Hip Belt Rear Pull Small with Flat Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

HB429M-A2 Body Point 4-Point Padded Rehab Latch Hip Belt Rear Pull Medium with Cinch Mount

Max Seated Hip Width	380	Min Seated Hip Width	230	Overall Length	1520 mm	Package Quantity	1
Package Type	Single						

HB429S-A2 Body Point 4-Point Padded Rehab Latch Hip Belt Rear Pull Small with Cinch Mount

Max Seated Hip Width	280	Min Seated Hip Width	180	Overall Length	1400 mm	Package Quantity	1
Package Type	Single						

Karma Wheelchair Accessories

Accessories for Karma wheelchairs

Huge range

• To help users take a standard wheelchair and customise it

Can make the wheelchair more comfortable and usable

2-KATBC	Karma Anti Tip Bar Cover
2-KCAS5X1S	Karma Castor 5"x1" Solid
2-KCAS6X1S	Karma Castor 6"x1" Solid
2-KCAS7X1S	Karma Castor 7"x1" Solid
2-KCAS8X1S	Karma Castor 8"x1" Solid
2-KSEBOP	Karma Seat Belt Plastic 1 Pi
2-KCAS6X1P	Karma Castor 6"x1" Pneuma
2-KRCB	Karma Rear Carry Bag
2-KFK-KM8-S	Karma Fork Short KM8
2-KCAS8X1P	Karma Castor 8"x1" Pneuma
2-KFK-125-W	Karma Fork Wide S125 - 7" x
2-KFK-KM8-W	Karma Fork Wide KM8 - 8" x
2-KFP-AADJ-16	Karma Footplate 16" Angle A
2-KFP-AADJ-18	Karma Footplate 18" Angle A
2-KFP-ALUM	Karma Footplate Aluminium
2-KATBKP25	Karma Anti Tip Bars KP25
2-KB-EXT-L	Karma Brake Extension Left
2-KB-EXT-R	Karma Brake Extension Righ
2-KCAS8X2MAC	Karma Castor 8" x 2" MCP
2-KSEXT2	Karma Seat Extension 2"
2-KSEBAUTO	Karma Seat Belt (Auto Style
2-KCAS7X134	Karma Castor 7" x 1.75"
2-KUP-EV-20-B	Karma Vinyl Back 20" Black
2-KUP-EV-1616-S	Karma Vinyl Seat 16" x 16" B
2-KUP-EV-16-BF	Karma Vinyl Back Fold 16"
2-KUP-EV-1816-S	Karma Vinyl Seat 18" x 16"
2-KUP-EV-1818-S	Karma Vinyl Seat 18" x 18"
2-KUP-EV-2018-S-RED	Karma Vinyl Seat 20" x 18" F
2-KWH24175	Karma Wheel 24" x 1.75" QR
2-KWH24175-F	Karma Wheel 24" x 1.75" Fix
2-KBR-CE-KM8	Karma Backrest KM8 Extens
2-KLRE-F16-L	Karma Legrest Elevating for
2-KATB-EC	Karma Anti Tip Bars for Ecor
2-KATB-ERGO	Karma Anti Tip Bars Ergo Se
2-KUP-EV-16-B	Karma Vinyl Back 16"
2-KUP-EV-18-B	Karma Vinyl Back 18"
2-KUP-EV-18-BF	Karma Vinyl Back Folding 18



ece	
atic	
atic	
(1.25"	
2"	
Adjustable	
Adjustable	
a & L Tube	
1	
ht	
Buckle)	
(
Black	
Red	
Solid Tyre	
red	
sion 22"	
r Flexx 16" Left	
nomy 800	
ries	
3" Black	

2-KUP-EV-2018-S	Karma Vinyl Seat 20" Black
2-KUP-EV-20-B-RED	Karma Vinyl Back 20" Red
2-KATBESW	Karma Anti Tip Bars Eagle Swivel
2-KSG	Karma Spoke Guards
2-KBR-RB-16	Karma Backrest Rigidiser Bar 16"
2-KBR-RB-18	Karma Backrest Rigidiser Bar 18"
2-KBR-RB-20	Karma Backrest Rigidiser Bar 20"
2-KBR-RB-22"	Karma Backrest Rigidiser Bar 22"
2-KB-WS-S125	Karma Brake Wire System S125
2-KH-E-EL	Karma Handle Ext Ergo Lite 7cm
2-KH-E-S125	Karma Handle Ext S125 7cm
2-KWH20138	Karma Wheel 20" x 1" 3/8
2-KB-WS-S125-T	Karma Brake Wire System S125 Transit
2-KWH24138S	Karma Wheel 24" x 1"3/8
2-KWH22138	Karma Wheel 22" x 1" 3/8
2-KWH24138M	Karma Wheel 24" x 1" 3/8 Mag
2-KLRE-125/EC-L	Karma Leg Rest Elevating S125/Econ Left
2-KLRE-125/EC-R	Karma Leg Rest Elevating S125/Econ Right
2-KLRE-KM8-L	Karma Legrest Elevating KM8 Left
2-KLRE-KM8-R	Karma Legrest Elevating KM8 Right
2-KUP-TA-16	Karma Tension Adjust 16"
2-KUP-TA-18	Karma Tension Adjust 18"
2-KTB	Karma Tie Down Brackets
2-KFP-KM8-HD-L	Karma Footplate Heavy Duty KM8 Left
2-KFP-KM8-HD-R	Karma Footplate Heavy Duty KM8 Right
2-KHRF	Karma Headrest Folding Large
2-KSS-F-L	Karma Stump Support Left for Flexx
2-KSS-F-R	Karma Stump Support Right for Flexx
2-KSSU	Karma Stump Support Ergo 125 /Econ
2-KWCB	Karma Wheelchair Travel Bag
2-KSEBP	Karma Car Buckle Seat Belt
2-KATB-E	Karma Anti-Tippers for Karma Eagle Wheelchairs
2-KTRAYB	Karma Tray for Karma Manual Wheelchairs
2-KUP-TA-20	Karma Tension Adjustable Backrest Upholstery for Karma 8020
2-KLRE-F16-R	Karma Legrest Elevating for Flexx 16" Right
2-KLRE-F18-L	Karma Legrest Elevating for Flexx 18" Left
2-KLRE-F18-R	Karma Legrest Elevating for Flexx 18" Right
2-KB-WS-EA	Karma Brake Wire System for Eagle
	Karma Brake Wire System for Eagle Karma Brake Wire System for Transit 2

	Karma Anti Tinnara far Karm
2-KATB-F	Karma Anti-Tippers for Karm
2-KSSL	Karma Stump Support for Ka
2-KSSR	Karma Stump Support for Ka
2-KUP-TA-22	Karma Tension Adjustable B
2-KIVPOLE	Karma IV Pole plus Bracket

Quickie Wheelchair Accessories

Accessories for Quickie Wheelchairs

• This variety of accessories aids the customisation of the standard wheelchair to suit the users and carers needs

184BU16	Quickie Reinforced Upholstery
184A120	Quickie Positioning Belt Slide Clar
184A2	Quickie Backpack - Black - Kids
184A1	Quickie Backpack - Black - Adult
184A117	Quickie Belt 1.5" Auto Buckle with
184CU1	Quickie 2" Cushion
184CS5	Quickie Castor Pin Locks
184A500	Quickie Gel Cover Hanger/Frame
184A994	Quickie Gel Cover Swing Away Arm
184BU23	Quickie Depth Adjustable Solid Ba
184SZ1	Quickie Stabiliser Bar
184BP33	Quickie DAB - Swing Away Stroller
184A118	Quickie Belt 2" Aircraft Buckle Pac
184SZ2	Quickie Folding Stabiliser Bar
184SG29	Quickie Plastic Low 6.5" High
184SG33	Quickie Plastic Standard 8.5" High
184A117	Quickie Belt 2" Aircraft Buckle wit
184CF53	Quickie Frog Legs Suspension 0-4
184CF54	Quickie Frog Legs Suspension 45-
184CF55	Quickie Frog Legs Suspension 86-
184A23	Quickie Spoke Guards for Quickie
184A302	Quickie Spoke Guards Clear
184AL16	Quickie One Arm Drive Adjustable
184AT6	Quickie Anti-tip Bars for Quickie 2
184A8	Quickie Under Seat Pouch for Quic
184BU8	Quickie Tension Adjustable Backre
184BP3	Quickie Push Handle 8"
184EH7	Quickie Extension Handle 6" Pair

na Flexx

arma Wheelchairs Left

arma Wheelchairs Right

Backrest Upholstery for Karma 8022

tery
de Clamps
Kids
Adult
le with Clamps
rame
vay Arm
olid Back
Stroller Handles
kle Pad with Clamp
Bar
gh
.5" High
kle with Clamps
ion 0-45kg
ion 45-86kg
ion 86-120kg
uickie 2 Pair
r
stable Axle Plate
ickie 2 Lite Pair
for Quickie 2 Lite
Backrest Quickie 2 Lite



Wheelchair

3.27 Manual Wheelchair Power Conversion Kits

Pride Power Assist - Wheelchair Conversion



- A wheelchair conversion kit for manual wheelchairs
- A kit with a drive unit which attaches to a manual wheelchair · Converts it to an attendant operated electric wheelchair
- Suitable for most manual wheelchairs including folding frames

PMW-PSS Pride Power Assist - Wheelchair Conversion

Battery x2	17 Ah	Charger Capacity	3 A	Max Speed	5 km/h	Max User Weight	120 kg
Overall Height	285 mm	Overall Length	280 mm	Overall Width	220 mm	Range	10 km
Total Weight Without Battery	14 kg						

PMW-PSHD Pride Power Assist - Wheelchair Conversion HD

Battery x2	12 Ah	Charger Capacity	2 A	Max Speed	6.4 km/h	Max User Weight	182 kg
Overall Height	380 mm	Overall Length	380 mm	Overall Width	250 mm	Range	16 km
Total Weight Without Battery	19.2 kg						

Merits Power Pack for Manual Wheelchair



- Twin-wheel design with much greater traction on grass and gravel
- · Easy attachment: attach to any wheelchair
- Let the power-pack to do all the work
- Full attendant controls, making pushing any wheelchair simple and easy • With speeds up to 5.5kph and a range of up to 16km
- Convenient forward and reverse direction switch
- Portable: lightweight for easy transport, packs neatly into a travel bag and enables extra mobility anywhere

F	POOID Merits Power Pack for Manual Wheelchair							
	Package Quantity	1	Package Type	Single				

Wheelchair

ADI AL Positioning Back Support, Active Contour & Low Height



• Ergonomically designed for proper posture • Eases the strain on the lower back during long active days in a chair

• It equally distributes forces through the pelvic base of support Results in less pain

• Designed for active individuals in lightweight wheelchairs

Solid seating using the lightest and strongest components available

• The Al Series back supports are designed to give as much or as little trunk support needed through two different contours

ALBR10LA2 ADI Aluminium Back Support 10" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	254 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBRIOLAR ADI Aluminium Back Support 10" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	254 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR11LA2 ADI Aluminium Back Support 11" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	280 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR11LAR ADI Aluminium Back Support 11" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	280 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR12LA2 ADI Aluminium Back Support 12" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	305 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR12LAR ADI Aluminium Back Support 12" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	305 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR13LA2 ADI Aluminium Back Support 13" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	330 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR13LAR ADI Aluminium Back Support 13" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	330 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR14LA2 ADI Aluminium Back Support 14" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		



3.28 Wheelchair Backrests Standard and Heavy Duty

ALBR14LAR ADI Aluminium Back Support 14" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR15LA2 ADI Aluminium Back Support 15" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR15LAR ADI Aluminium Back Support 15" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR16LA2 ADI Aluminium Back Support 16" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR16LAR ADI Aluminium Back Support 16" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR17LA2 ADI Aluminium Back Support 17" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR17LAR ADI Aluminium Back Support 17" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR18LA2 ADI Aluminium Back Support 18" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR18LAR ADI Aluminium Back Support 18" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR19LA2 ADI Aluminium Back Support 19" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR19LAR ADI Aluminium Back Support 19" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ADI AL Positioning Back Support, Extra Deep Contour & Tall Height



Ergonomically designed for proper posture

- It equally distributes forces through the pelvic base of support Results in less pain
- Designed for active individuals in lightweight wheelchairs
- ALBR16T62 ADI Aluminium Back Support 16" Wide with 2 Point Hardware 150 mm Max User Weight 110 kg Contour Depth Package Quantity Package Type Single ALBR17T62 ADI Aluminium Back Support 17" Wide with 2 Point Hardware Contour Depth 150 mm Max User Weight 110 kg Single Package Quantity Package Type

ALBR18T62 ADI Aluminium Back Support 18" Wide with 2 Point Hardware

Contour Depth	150 mm	Max User Weight	110 kg	Overall Height	457 mm	Overall Width	457 mm		
Package Quantity	1	Package Type	Single	Product Weight	1 kg				
ALBR19T62 ADI Alur									

Contour Depth	150 mm	Max User Weight	110 kg	Overall Height	457 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR20T62 ADI Aluminium Back Support 20" Wide with 2 Point Hardware

Contour Depth	150 mm	Max User Weight	110 kg	Overall Height	457 mm	Overall Width	510 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ADI AL Positioning Back Support, Deep Contour & XTal Height



· Eases the strain on the lower back during long active days in a chair • It equally distributes forces through the pelvic base of support

· Results in less pain

• Designed for active individuals in lightweight wheelchairs

ALBR16XDR ADI Alur	LBR16XDR ADI Aluminium Back Support 16" Wide with Quick Release Hardware							
Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	510 mm	Overall Width	406 mm	
Package Quantity	1	Package Type	Single	Product Weight	1 kg			

ALBR17XDR	ADI Alur	ninium B	ack Support 17" Wid	de with Q)u
Contour Depth		100 mm	Max User Weight	110 kg	0
Package Quantity		1	Package Type	Single	P

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	510 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR18XDR	ADI Aluminium	Back Support	18" Wide with Q

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	510 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		



• Eases the strain on the lower back during long active days in a chair

• Solid seating using the lightest and strongest components available

• The AI Series back supports are designed to give as much or as little trunk support needed through two different contours

Overall Height	457 mm	Overall Width	406 mm
Product Weight	1 kg		

Overall Height	457 mm	Overall Width	430 mm
Product Weight	1 kg		

· Ergonomically designed for proper posture

· Solid seating using the lightest and strongest components available

• The AI Series back supports are designed to give as much or as little trunk support needed through two different contours

uick Release Hardware

uick Release Hardware

ALBR19XDR ADI Aluminium Back Support 19" Wide with Quick Release Hardware									
Contour Depth		100 mm	Max User Weight	110 kg	Overall Height	510 mm	Overall Width	482 mm	
Package Quantity		1	Package Type	Single	Product Weight	1 kg			
				-	-	-			
ALBR20XDR #	ADI Alur	ninium I	Back Support 20" W	ide with	Quick Release Hard	ware			

Single Product Weight

ADI AL Positioning Back Support, Deep Contour & Low Height

Package Type



Package Quantity

- Ergonomically designed for proper posture
- Eases the strain on the lower back during long active days in a chair
- It equally distributes forces through the pelvic base of support Results in less pain
- Designed for active individuals in lightweight wheelchairs
- Solid seating using the lightest and strongest components available
- The AI Series back supports are designed to give as much or as little trunk support needed through two different contours

1 kg

ALBR10LD2 ADI Aluminium Back Support 10" Wide with 2 Point Hardware								
	Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	254 mm
	Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBRIOLDR ADI Aluminium Back Support 10" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	254 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR11LD2 ADI Aluminium Back Support 11" Wide with 2 Point Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	280 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR11LDR ADI Aluminium Back Support 11" Wide with Quick Release Hardware

[Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	280 mm
	Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR12LD2 ADI Aluminium Back Support 12" Wide with 2 Point Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	305 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR12LDR ADI Aluminium Back Support 12" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	305 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR13LD2 ADI Aluminium Back Support 13" Wide with 2 Point Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	330 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR13LDR ADI Aluminium Back Support 13" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	330 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR14LD2 ADI Aluminium Back Support 14" Wide with 2 Point Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR14LDR ADI Aluminium Back Support 14" Wide wit	h Ç
---	-----

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR15LD2 ADI Aluminium Back Support 15" Wide with 2 Point Hardware

(F

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR15LDR ADI Aluminium Back Support 15" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR16LD2 ADI Aluminium Back Support 16" Wide with 2 Point Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBRIGLDR ADI Aluminium Back Support 16" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR17LD2 ADI Aluminium Back Support 17" Wide with 2 Point Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR17LDR ADI Aluminium Back Support 17" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR18LD2 ADI Aluminium Back Support 18" Wide with 2 Point Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR18LDR ADI Aluminium Back Support 18" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR19LD2 ADI Aluminium Back Support 19" Wide with 2 Point Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR19LDR ADI Aluminium Back Support 19" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		



Ouick Release Hardware

ADI AL Positioning Back Support, Extra Deep Contour & Standard Height



Ergonomically designed for proper posture

- Eases the strain on the lower back during long active days in a chair
- It equally distributes forces through the pelvic base of support
- Results in less pain Designed for active individuals in lightweight wheelchairs
- Solid seating using the lightest and strongest components available
- The AI Series back supports are designed to give as much or as little trunk support needed through two different contours

ALBR16S62 ADI Aluminium Back Support 16" Wide with 2 Point Hardware

Contour Depth	150 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	406 mm				
Package Quantity	1	Package Type	Single	Product Weight	1 kg						

ALBR17S62 ADI Aluminium Back Support 17" Wide with 2 Point Hardware

Contour Depth	150 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR18S62 ADI Aluminium Back Support 18" Wide with 2 Point Hardware

Contour Depth	150 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR19S62 ADI Aluminium Back Support 19" Wide with 2 Point Hardware

Contour Depth	150 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR20S62 ADI Aluminium Back Support 20" Wide with 2 Point Hardware

Contour Depth	150 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	510 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ADI AL Positioning Back Support, Active Contour & Medium Height



- Ergonomically designed for proper posture
- Eases the strain on the lower back during long active days in a chair
- It equally distributes forces through the pelvic base of support
- · Results in less pain
- Designed for active individuals in lightweight wheelchairs
- Solid seating using the lightest and strongest components available
- The AI Series back supports are designed to give as much or as little trunk support needed through two different contours

ALBR12MA2 ADI Aluminium Back Support 12" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	305 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR12MAR ADI Aluminium Back Support 12" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	305 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR13MA2 ADI Aluminium Back Support 13" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	330 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR13MAR ADI Aluminium Back Support 13" Wide with Quick Release Hardware Contour Depth 63.5 mm Max User Weight 110 kg Package Quantity Package Type Single ALBR14MA2 ADI Aluminium Back Support 14" Wide with 2 Point Hardware Contour Depth 63.5 mm Max User Weight 110 kg Package Quantity Single Package Type ALBR14MAR ADI Aluminium Back Support 14" Wide with Quick Release Hardware Contour Depth 63.5 mm Max User Weight 110 kg Package Quantity Package Type Single ALBR15MA2 ADI Aluminium Back Support 15" Wide with 2 Point Hardware Contour Depth 63.5 mm Max User Weight 110 kg Package Quantity Package Type Single ALBR15MAR ADI Aluminium Back Support 15" Wide with Quick Release Hardware Contour Depth 63.5 mm Max User Weight 110 kg Single Package Quantity Package Type ALBR16MA2 ADI Aluminium Back Support 16" Wide with 2 Point Hardware 63.5 mm Max User Weight Contour Depth 110 kg Package Quantity Package Type Single ALREIGMAD ADI Aluminium Back Support 16" Wide with Ouick Pelease Ha

ALDRIGIVIAR ADI AIU												
Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	406 mm					
Package Quantity	1	Package Type	Single	Product Weight	1 kg							

ALBR17MA2	ADI Alu	minium I	Back Support 17" Wi	de with	2
Contour Depth		63.5 mm	Max User Weight	110 kg	0

Package Type

Package Quantity

_ _ _ _ _ _

ALBR17MAR ADI Aluminium Back Support 17" Wide with Quick Release Hardware

				•			
Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

Single

ALBR18MA2 ADI Aluminium B	lack Support 18" Wide with 2
---------------------------	------------------------------

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR18MAR ADI Aluminium Back Support 18" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR19MA2 ADI Aluminium Back Support 19" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR19MAR ADI Aluminium Back Support 19" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		



Overall Height	330 mm	Overall Width	330 mm			
Product Weight	1 kg					

Overall Height	330 mm	Overall Width	355 mm
Product Weight	1 kg		

Overall Height	330 mm	Overall Width	355 mm
Product Weight	1 kg		

Overall Height	330 mm	Overall Width	380 mm
Product Weight	1 kg		

Overall Height	330 mm	Overall Width	380 mm
Product Weight	1 kg		

Overall Height	330 mm	Overall Width	406 mm
Product Weight	1 kg		

Point Hardware

Overall Height	330 mm	Overall Width	430 mm
Product Weight	1 kg		

2 Point Hardware

ADI AL Positioning Back Support, Active Contour & Standard Height



Ergonomically designed for proper posture

- Eases the strain on the lower back during long active days in a chair • It equally distributes forces through the pelvic base of support
- Results in less pain
- Designed for active individuals in lightweight wheelchairs
- Solid seating using the lightest and strongest components available
- The AI Series back supports are designed to give as much or as little trunk support needed through two different contours

ALBR14SA2 ADI Aluminium Back Support 14" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR15SA2 ADI Aluminium Back Support 15" Wide with 2 Point Hardware

Co	ontour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	380 mm
Ра	ackage Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR16SA2 ADI Aluminium Back Support 16" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR17SA2 ADI Aluminium Back Support 17" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR18SA2 ADI Aluminium Back Support 18" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR19SA2 ADI Aluminium Back Support 19" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR20SA2 ADI Aluminium Back Support 20" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	510 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR14SAR ADI Aluminium Back Support 14" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR15SAR ADI Aluminium Back Support 15" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR16SAR ADI Aluminium Back Support 16" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR17SAR ADI Aluminium Back Support 17" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR18SAR ADI Aluminium Back Support 18" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBRI9SAR ADI Aluminium Back Support 19" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR20SAR ADI Aluminium Back Support 20" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	510 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ADI AL Positioning Back Support, Deep Contour & Medium Height



F

Ergonomically designed for proper posture

• Eases the strain on the lower back during long active days in a chair • It equally distributes forces through the pelvic base of support

- Results in less pain
- Designed for active individuals in lightweight wheelchairs
- · Solid seating using the lightest and strongest components available

ALBR12MD2 ADI Aluminium Back Support 12" Wide with 2 Point Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	305 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR12MDR ADI Aluminium Back Support 12" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	305 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR13MD2 ADI Aluminium Back Support 13" Wide with 2 Point Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	330 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR13MDR ADI Aluminium Back Support 13" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	330 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR14MD2 ADI Aluminium Back Support 14" Wide with 2 Point Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR14MDR ADI Aluminium Back Support 14" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR15MD2 ADI Aluminium Back Support 15" Wide with 2 Point Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

The AI Series back supports are designed to give as much or as little trunk support needed through two different contours



ALBR15MDR ADI Aluminium Back Support 15" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR16MD2 ADI Aluminium Back Support 16" Wide with 2 Point Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR16MDR ADI Aluminium Back Support 16" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR17MD2 ADI Aluminium Back Support 17" Wide with 2 Point Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR18MD2 ADI Aluminium Back Support 18" Wide with 2 Point Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR18MDR ADI Aluminium Back Support 18" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR19MD2 ADI Aluminium Back Support 19" Wide with 2 Point Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR19MDR ADI Aluminium Back Support 19" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR17MDR ADI Aluminium Back Support 17" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ADI AL Positioning Back Support, Deep Contour & Standard Height



•	Ergonomicall	y d	lesigned	for	proper	posture	
---	--------------	-----	----------	-----	--------	---------	--

- · Eases the strain on the lower back during long active days in a chair
- It equally distributes forces through the pelvic base of support
- Results in less pain
- Designed for active individuals in lightweight wheelchairs Solid seating using the lightest and strongest components available
- . The AI Series back supports are designed to give as much or as little trunk support needed through two different contours

ALBR14SD2 ADI Aluminium Back Support 14" Wide with 2 Point Hardware

[Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	355 mm
	Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR14SDR ADI Aluminium Back Support 14" Wide with Quick Release Hardware 100 mm Max User Weight 110 kg Contour Depth Package Quantity Package Type Single ALBR15SD2 ADI Aluminium Back Support 15" Wide with 2 Point Hardware Contour Depth 100 mm Max User Weight 110 kg Package Quantity Package Type Single ALBR15SDR ADI Aluminium Back Support 15" Wide with Quick Release Hardware Contour Depth 100 mm Max User Weight 110 kg Single Package Quantity Package Type ALBR16SD2 ADI Aluminium Back Support 16" Wide with 2 Point Hardware 100 mm Max User Weight Contour Depth 110 kg Single Package Quantity Package Type ALBR16SDR ADI Aluminium Back Support 16" Wide with Quick Release Hardware 100 mm Max User Weight Contour Depth 110 kg Package Quantity Package Type Single ALBR17SD2 ADI Aluminium Back Support 17" Wide with 2 Point Hardware 100 mm Max User Weight Contour Depth 110 kg Single Package Quantity Package Type ALBR17SDR ADI Aluminium Back Support 17" Wide with Quick Release Hardware Contour Depth 100 mm Max User Weight 110 kg Package Quantity Package Type Single ALBR18SD2 ADI Aluminium Back Support 18" Wide with 2 Point Hardware Contour Depth 100 mm Max User Weight 110 kg Package Quantity Single Package Type ALBR18SDR ADI Aluminium Back Support 18" Wide with Quick Release Hardware Contour Depth 100 mm Max User Weight 110 kg Package Quantity Package Type Single ALBR19SD2 ADI Aluminium Back Support 19" Wide with 2 Point Hardware Contour Depth 100 mm Max User Weight 110 kg Package Quantity Package Type Single

ALBRI9SDR ADI Aluminium Back Support 19" Wide with Quick Release Hardware									
Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	482 mm		
Package Quantity	1	Package Type	Single	Product Weight	1 kg				

ALBR20SD2 ADI Aluminium Back Support 20" Wide with 2 Point Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	510 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR2OSDR ADI Aluminium Back Support 20" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	510 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		



Overall Height	406 mm	Overall Width	355 mm					
Product Weight	1 kg							

Overall Height	406 mm	Overall Width	380 mm
Product Weight	1 kg		

Overall Height	406 mm	Overall Width	380 mm
Product Weight	1 kg		

Overall Height	406 mm	Overall Width	406 mm
Product Weight	1 kg		

Overall Height	406 mm	Overall Width	406 mm
Product Weight	1 kg		

Overall Height	406 mm	Overall Width	430 mm
Product Weight	1 kg		

Overall Height	406 mm	Overall Width	430 mm
Product Weight	1 kg		

Overall Height	406 mm	Overall Width	457 mm
Product Weight	1 kg		

Overall Height	406 mm	Overall Width	457 mm
Product Weight	1 kg		

Overall Height	406 mm	Overall Width	482 mm
Product Weight	1 kg		

ADI AL Positioning Back Support, Deep Contour & Tal Height



- Ergonomically designed for proper posture
- Eases the strain on the lower back during long active days in a chair
- It equally distributes forces through the pelvic base of support
- Results in less pain
- Designed for active individuals in lightweight wheelchairs • Solid seating using the lightest and strongest components available
- The AI Series back supports are designed to give as much or as little trunk support needed through two different contours

ALBRIGTDR ADI Aluminium Back Support 16" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	457 mm	Overall Width	406 mm	
Package Quantity	1	Package Type	Single	Product Weight	1 kg			

ALBR17TDR ADI Aluminium Back Support 20" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	457 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR18TDR ADI Aluminium Back Support 17" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	457 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR19TDR ADI Aluminium Back Support 18" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	457 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR20TDR ADI Aluminium Back Support 19" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	457 mm	Overall Width	510 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ADI AL Positioning Back Support, Active Contour & Tall Height



- Ergonomically designed for proper posture
- Eases the strain on the lower back during long active days in a chair
- It equally distributes forces through the pelvic base of support
- · Results in less pain
- Designed for active individuals in lightweight wheelchairs
- Solid seating using the lightest and strongest components available
- The AI Series back supports are designed to give as much or as little trunk support needed through two different contours

ALBR16TAR ADI Aluminium Back Support 16" Wide with Quick Release Hardware											
Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	457 mm	Overall Width	406 mm				
Package Quantity	1	Package Type	Single	Product Weight	1 kg						

ALBR17TAR ADI Aluminium Back Support 17" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	457 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR18TAR ADI Aluminium Back Support 18" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	457 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	457 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR20TAR ADI Aluminium Back Support 20" Wide with Quick Release Hardware

				•			
Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	457 mm	Overall Width	510 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ADI AL Positioning Back Support, Active Contour & XTall Height



- Results in less pain Designed for active individuals in lightweight wheelchairs
- Solid seating using the lightest and strongest components available

ALBR16XAR ADI Aluminium Back Support 16" Wide with Quick Release Hardware											
Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	510 mm	Overall Width	406 mm				
Package Quantity	1	Package Type	Single	Product Weight	1 kg						

ALBR17XAR ADI Aluminium Back Support 17" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	510 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR18XAR ADI Aluminium Back Support 18" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	510 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR19XAR ADI Aluminium Back Support 19" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	510 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR20XAR ADI Aluminium Back Support 20" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	510 mm	Overall Width	510 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ADI AL Positioning Back Support, Extra Deep Contour & Extra Tall Height



[•] It equally distributes forces through the pelvic base of support Results in less pain

- Designed for active individuals in lightweight wheelchairs
- Solid seating using the lightest and strongest components available

ALBR16X62 ADI Aluminium Back Support 16" Wide with 2 Point Hardware

Contour Depth	150 mm	Max User Weight	110 kg	Overall Height	510 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

50



)uick Release Hardware

• Ergonomically designed for proper posture

• Eases the strain on the lower back during long active days in a chair • It equally distributes forces through the pelvic base of support

• The AI Series back supports are designed to give as much or as little trunk support needed through two different contours

· Ergonomically designed for proper posture

• Eases the strain on the lower back during long active days in a chair

• The AI Series back supports are designed to give as much or as little trunk support needed through two different contours

ALBR17X62 ADI Aluminium Back Support 17" Wide with 2 Point Hardware

Contour Depth	150 mm	Max User Weight	110 kg	Overall Height	510 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR18X62 ADI Aluminium Back Support 18" Wide with 2 Point Hardware

Contour Depth	150 mm	Max User Weight	110 kg	Overall Height	510 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR19X62 ADI Aluminium Back Support 19" Wide with 2 Point Hardware

[Contour Depth	150 mm	Max User Weight	110 kg	Overall Height	510 mm	Overall Width	482 mm
	Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ALBR20X62 ADI Aluminium Back Support 20" Wide with 2 Point Hardware

Contour Depth	150 mm	Max User Weight	110 kg	Overall Height	510 mm	Overall Width	510 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ADI CF Positioning Back Support, Active Contour & Low Height



- Ergonomically designed for proper posture
- Eases the strain on the lower back during long active days in a chair • It equally distributes forces through the pelvic base of support
- Results in less pain
- Designed for active individuals in lightweight wheelchairs
- Solid seating using the lightest and strongest components available
- The CF Series Back Supports are an ultra lightweight seat back
- · Using the lightest and strongest components available

CFBR13LA2 ADI Carbon Back Support 13" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	330 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR13LAF ADI Carbon Back Support 13" Wide with Fixed Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	330 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR14LA2 ADI Carbon Back Support 14" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR14LAF ADI Carbon Back Support 14" Wide with Fixed Hardware

Contour Depth	63.5 n	mm Max User Weight	110 kg	Overall Height	250 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR15LA2 ADI Carbon Back Support 15" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR15LAF ADI Carbon Back Support 15" Wide with Fixed Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR16LA2 ADI Carbon Back Support 16" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

	CFBR16LAF ADI Ca	rbon Bac	k Support 16" Wide v	with Fixe	d Hardware			
Packane Quantity 1 Packane Type Single Product Weight 1 kg	Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	406 mm
radiage quantity i radiage type ongre Product weight rive	Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR17LA2 ADI Carbon Back Support 17" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR17LAF ADI Carbon Back Support 17" Wide with Fixed Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR18LA2 ADI Carbon Back Support 18" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR18LAF ADI Carbon Back Support 18" Wide with Fixed Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR19LA2 ADI Carbon Back Support 19" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR19LAF ADI Carbon Back Support 19" Wide with Fixed Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR13LAR ADI Carbon Back Support 13" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	330 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR14LAR ADI Carbon Back Support 14" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR15LAR ADI Carbon Back Support 15" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR16LAR ADI Carbon Back Support 16" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR17LAR ADI Carbon Back Support 17" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

Mobility

CFBRIBLAR ADI Ca	bon Bac	k Support 18" Wide v	with Quio	k Release Hardwar	е		
Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		
CFBR19LAR ADI Ca	bon Bac	k Support 19" Wide v	with Quid	ck Release Hardward	9		
CFBR19LAR ADI Car Contour Depth	bon Bac	k Support 19" Wide v Max User Weight	vith Quid	ck Release Hardward	e 250 mm	Overall Width	482 mm

ADI CF Positioning Back Support, Active Contour & Medium Height



- Ergonomically designed for proper posture
- Eases the strain on the lower back during long active days in a chair
- It equally distributes forces through the pelvic base of support Results in less pain
- Designed for active individuals in lightweight wheelchairs
- Solid seating using the lightest and strongest components available
- The CF Series Back Supports are an ultra lightweight seat back
- Using the lightest and strongest components available

CFBR13MAF ADJ Carbon Back Support 13" Wide with Fixed Hardware

CI DITI SIMAI ADI CAI	DOIT Dac	k Support 15 Wide		Eu Haluwale			
Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	330 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR13MAR ADI Carbon Back Support 13" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	330 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR15MA2 ADI Carbon Back Support 15" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR18MAR ADI Carbon Back Support 18" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR15MAF ADI Carbon Back Support 15" Wide with Fixed Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR19MA2 ADI Carbon Back Support 19" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR16MAR ADI Carbon Back Support 16" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR17MAR ADI Carbon Back Support 17" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR13MA2 ADI Carbon Back Support 13" Wide with 2 Point Hardware

Package Quantity 1 Package Type Circle Package Weight 1/2	Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	330 mm
Package Quantity I Package Type Single Product weight I kg	Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR14MA2 ADI Carbon Back Support 14" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR14MAF ADI Carbon Back Support 14" Wide with Fixed Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR15MAR ADI Carbon Back Support 15" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR16MA2 ADI Carbon Back Support 16" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR19MAF ADI Carbon Back Support 19" Wide with Fixed Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR19MAR ADI Carbon Back Support 19" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR17MA2 ADI Carbon Back Support 17" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR14MAR ADI Carbon Back Support 14" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR18MA2 ADI Carbon Back Support 18" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR18MAF ADI Carbon Back Support 18" Wide with Fixed Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR16MAF ADI Carbon Back Support 16" Wide with Fixed Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR17MAF ADI Carbon Back Support 17" Wide with Fixed Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

Mobility



ADI CF Positioning Back Support, Deep Contour & Low Height



Ergonomically designed for proper posture

- Eases the strain on the lower back during long active days in a chair
- It equally distributes forces through the pelvic base of support
- Results in less pain Designed for active individuals in lightweight wheelchairs
- Solid seating using the lightest and strongest components available
- The CF Series Back Supports are an ultra lightweight seat back
- · Using the lightest and strongest components available

C	ontour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	482 mm
Р	ackage Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR14LDR ADI Carbon Back Support 14" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR17LDR ADI Carbon Back Support 17" Wide with Quick Release Hardware

[Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	430 mm
	Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR15LDR ADI Carbon Back Support 15" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR18LDR ADI Carbon Back Support 18" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR13LDR ADI Carbon Back Support 13" Wide with Quick Release Hardware

Contour Dept	th	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	330 mm
Package Qua	antity	1	Package Type	Single	Product Weight	1 kg		

CFBR16LDR ADI Carbon Back Support 16" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	250 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ADI CF Positioning Back Support, Deep Contour & Medium Height



P

- Ergonomically designed for proper posture • Eases the strain on the lower back during long active days in a chair
- It equally distributes forces through the pelvic base of support Results in less pain
- Designed for active individuals in lightweight wheelchairs
- · Solid seating using the lightest and strongest components available
- The CF Series Back Supports are an ultra lightweight seat back
- Using the lightest and strongest components available

Contour Deoth 100 mm Max User Weight 110 kg Overall Height 330 mm Overall Width 33	
	330 mm
Package Quantity 1 Package Type Single Product Weight 1 kg	

CFBR14MDR ADI Carbon Back Support 14" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR15MDF ADI Carbon Back Support 15" Wide with Fixed Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		
		·					

CFBR18MDF ADI Carbon Back Support 18" Wide with Fixed Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR19MDF ADI Carbon Back Support 19" Wide with Fixed Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR19MDR ADI Carbon Back Support 19" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR17MDR ADI Carbon Back Support 17" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR13MDR ADI Carbon Back Support 13" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	330 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR14MDF ADI Carbon Back Support 14" Wide with Fixed Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR15MDR ADI Carbon Back Support 15" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		



CFBR16MDF ADI Ca	rbon Bac	k Support 16" Wide	with Fixe	ed Hardware			
Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		
CFBR16MDR ADI Ca	irbon Bad	k Support 16" Wide	with Qui	ck Release Hardwa	re		
Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		
CFBR17MDF ADI Ca	rbon Bac	k Support 17" Wide	with Fixe	ed Hardware			
Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	430 mm

Single

CFBR18MDR ADI Carbon Back Support 18" Wide with Quick Release Hardware

Package Type

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	330 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

Product Weight

1 kg

ADI CF Positioning Back Support, Active Contour & Standard Height



Package Quantity

- · Ergonomically designed for proper posture
- Eases the strain on the lower back during long active days in a chair • It equally distributes forces through the pelvic base of support
- Results in less pain
- Designed for active individuals in lightweight wheelchairs
- Solid seating using the lightest and strongest components available
- The CF Series Back Supports are an ultra lightweight seat back
- Using the lightest and strongest components available

CFBR13SAF ADI Carbon Back Support 13" Wide with Fixed Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	330 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR13SAR ADI Carbon Back Support 13" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	330 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR14SA2 ADI Carbon Back Support 14" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR15SAR ADI Carbon Back Support 15" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR16SAF ADI Carbon Back Support 16" Wide with Fixed Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR16SAR ADI Carbon Back Support 16" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

Contour Depth	63.5 mm	Max User Weight
Package Quantity	1	Package Type
	Carbon Bac	k Support 18" Wide
	Carbon Bac	

CFBR13SA2 ADI Carbon Back Support 13" Wide with 2 Point Hardware

Package Quantity 1 Package Type Single Product Weight 1 kg	Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	330 mm
Toologe quantity in Toologe type onight Trouble treight the	Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR14SAF ADI Carbon Back Support 14" Wide with Fixed Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR14SAR ADI Carbon Back Support 14" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR15SA2 ADI Carbon Back Support 15" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR15SAF ADI Carbon Back Support 15" Wide with Fixed Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR16SA2 ADI Carbon Back Support 16" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR17SAR ADI Carbon Back Support 17" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR19SA2 ADI Carbon Back Support 19" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR19SAR ADI Carbon Back Support 19" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR17SA2 ADI Carbon Back Support 17" Wide with 2 Point Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR17SAF ADI Carbon Back Support 17" Wide with Fixed Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		



2 Point Hardware

Overall Height	406 mm	Overall Width	457 mm
Product Weight	1 kg		

Fixed Hardware

Overall Height	406 mm	Overall Width	457 mm
Product Weight	1 kg		

CFBR18SAR ADI Carbon Back Support 18" Wide with Quick Release Hardware

Contour Depth	63.5 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

ADI CF Positioning Back Support, Deep Contour & Standard Height



- Ergonomically designed for proper posture
- Eases the strain on the lower back during long active days in a chair
- It equally distributes forces through the pelvic base of support Results in less pain
- Designed for active individuals in lightweight wheelchairs
- Solid seating using the lightest and strongest components available
- The CF Series Back Supports are an ultra lightweight seat back
- Using the lightest and strongest components available

CFBR14SDR ADI Carbon Back Support 14" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	355 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR15SDR ADI Carbon Back Support 15" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	380 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR17SDR ADI Carbon Back Support 17" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	430 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR13SDR ADI Carbon Back Support 13" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	330 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR16SDR ADI Carbon Back Support 16" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	406 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		


CFBR18SDR ADI Carbon Back Support 18" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	406 mm	Overall Width	457 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

CFBR19SDR ADI Carbon Back Support 19" Wide with Quick Release Hardware

Contour Depth	100 mm	Max User Weight	110 kg	Overall Height	409 mm	Overall Width	482 mm
Package Quantity	1	Package Type	Single	Product Weight	1 kg		

Jay 3 Posterior Mid Thoracic Back Support with Standard Hardware



• Designed for moderate to high postural needs

J3PAMTS12SH Jay 3 PMT Back Support with Standard Hardware Short 12"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.1 kg						

J3PAMTS14SH Jay 3 PMT Back Support with Standard Hardware Short 14"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.1 kg						

J3PAMTS16SH Jay 3 PMT Back Support with Standard Hardware Short 16"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.1 kg						

J3PAMTS18SH Jay 3 PMT Back Support with Standard Hardware Short 18"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.1 kg						

J3PAMTS20SH Jay 3 PMT Back Support with Standard Hardware Short 20"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.1 kg						

J3PAMTM12SH Jay 3 PMT Back Support with Standard Hardware Medium 12"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.1 kg						

J3PAMTM14SH Jay 3 PMT Back Support with Standard Hardware Medium 14"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.1 kg						

J3PAMTM16SH Jay 3 PMT Back Support with Standard Hardware Medium 16"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.1 kg						

Mobility



• The comprehensive range meets functional needs of wheelchair users featuring many sizes and easy attachment hardware

J3PAMTM18SH Jay 3 PMT Back Support with Standard Hardware Medium 18"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.1 kg						

J3PAMTM20SH Jay 3 PMT Back Support with Standard Hardware Medium 20"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.1 kg						

J3PAMTT12SH Jay 3 PMT Back Support with Standard Hardware Tall 12"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.1 kg						

J3PAMTT14SH Jay 3 PMT Back Support with Standard Hardware Tall 14"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.1 kg						

J3PAMTT16SH Jay 3 PMT Back Support with Standard Hardware Tall 16"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.1 kg						

J3PAMTT18SH Jay 3 PMT Back Support with Standard Hardware Tall 18"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.1 kg						

J3PAMTT20SH Jay 3 PMT Back Support with Standard Hardware Tall 20"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.1 kg						

Jay 3 Posterior Upper Thoracic Back Support with 4-Point Hardware



C N P

C

• Jay 3 Posterior Upper Thoracic Back Support with 4 Point Hardware • Support the back through the Upper Thoracic support position



- - Designed for moderate to high postural needs

J3PAUTS12FS Jay 3 PUT Back Support with 4-Point Hardware Short 12"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.3 kg						

J3PAUTS14FS Jay 3 PUT Back Support with 4-Point Hardware Short 14"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.3 kg						

J3PAUTS16FS Jay 3 PUT Back Support with 4-Point Hardware Short 16"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.3 kg						

J3PAUTS18FS Jay 3 PUT Back Support with 4-Point Hardware Short 18"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.3 kg						

J3PAUTS20FS Jay 3 PUT Back Support with 4-Point Hardware Short 20"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.3 kg						

J3PAUTM14FS Jay 3 PUT Back Support with 4-Point Hardware Medium 14"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.3 kg						

J3PAUTM16FS Jay 3 PUT Back Support with 4-Point Hardware Medium 16"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.3 kg						

J3PAUTM18FS Jay 3 PUT Back Support with 4-Point Hardware Medium 18"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.3 kg						



• Available in a variety of sizes, levels of support, contour depths and mounting hardware

• Comprehensive range to meet functional needs of wheelchair users featuring many sizes and easy attachment hardware

· Provides moderate contouring for clients that would benefit for more lateral support

J3PAUTM20FS Jay 3 PUT Back Support with 4-Point Hardware Medium 20"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.3 kg						

J3PAUTM22FS Jay 3 PUT Back Support with 4-Point Hardware Medium 22"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	225 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.3 kg						

J3PAUTM24FS Jay 3 PUT Back Support with 4-Point Hardware Medium 24"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	225 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.3 kg						

J3PAUTT14FS Jay 3 PUT Back Support with 4-Point Hardware Tall 14"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.3 kg						

J3PAUTT16FS Jay 3 PUT Back Support with 4-Point Hardware Tall 16"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.3 kg						

J3PAUTT18FS Jay 3 PUT Back Support with 4-Point Hardware Tall 18"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.3 kg						


J3PAUTT20FS Jay 3 PUT Back Support with 4-Point Hardware Tall 20"

Contour Depth	55 mm	Max Overall Height	500 mm	Max Overall Width	510 mm	Max User Weight	136 kg
Min Overall Height	420 mm	Min Overall Width	320 mm	Package Quantity	1	Package Type	Single
Product Weight	2.3 kg						

Jay 3 Posterior Lateral Mid Thoracic Back Support with Standard Hardware



- The J3 is designed for moderate to high postural needs



505 mm Package Quantity

Overall Width



• Posterior and Lateral support through the Upper Thoracic support position • Available in a variety of sizes, levels of support, contour depths and mounting hardware • This J3 comprehensive back range meets functional needs of wheelchair users

Max User Weight	136 kg	Min Overall Height	305 mm
Package Type	Single	Product Weight	2.1 kg

Max User Weight	136 kg	Min Overall Height	305 mm
Package Type	Single	Product Weight	2.1 kg

Max User Weight	136 kg	Min Overall Height	305 mm
Package Type	Single	Product Weight	2.1 kg

Max User Weight	136 kg	Min Overall Height	305 mm
Package Type	Single	Product Weight	2.1 kg

Max User Weight	136 kg	Min Overall Height	305 mm
Package Type	Single	Product Weight	2.1 kg

Max User Weight	136 kg	Min Overall Height	305 mm
Package Type	Single	Product Weight	2.1 kg

Max User Weight	136 kg	Min Overall Height	305 mm
Package Type	Single	Product Weight	2.1 kg

Max User Weight	136 kg	Min Overall Height	305 mm
Package Type	Single	Product Weight	2.1 kg

Max User Weight	136 kg	Min Overall Height	305 mm
Package Type	Single	Product Weight	2.1 kg

Max User Weight	136 kg	Min Overall Height	305 mm
Package Type	Single	Product Weight	2.1 kg

1001 MTT10011 -

Overall Width

J3PLM1112SH	Jay 3	PL UT Ba	ack Support with St	andard H	lardware Tall 12"			
Contour Depth		76 mm	Max Overall Height	610 mm	Max User Weight	136 kg	Min Overall Height	305 mm
Overall Width		305 mm	Package Quantity	1	Package Type	Single	Product Weight	2.1 kg
J3PLMTT14SH	Jay 3	PL UT Ba	ack Support with St	andard H	lardware Tall 14"			
Contour Depth		76 mm	Max Overall Height	610 mm	Max User Weight	136 ka	Min Overall Height	305 mm

. . .

... ...

J3PLMTT16SH Jay 3 PL UT Back Support with Standard Hardware Tall 16"

Package Quantity

355 mm

Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	136 kg	Min Overall Height	305 mm
Overall Width	405 mm	Package Quantity	1	Package Type	Single	Product Weight	2.1 kg

Package Type

Single

Product Weight

2.1 kg

J3PLMTT18SH Jay 3 PL UT Back Support with Standard Hardware Tall 18"

Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	136 kg	Min Overall Height	305 mm
Overall Width	455 mm	Package Quantity	1	Package Type	Single	Product Weight	2.1 kg

J3PLMTT20SH Jay 3 PL UT Back Support with Standard Hardware Tall 20"

[Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	136 kg	Min Overall Height	305 mm
	Overall Width	505 mm	Package Quantity	1	Package Type	Single	Product Weight	2.1 kg

Jay 3 Posterior Lateral Shoulder Height Back Support, 4-Point Hardware



- Available in a variety of sizes, levels of support, contour depths and mounting hardware
- · Comprehensive back range to meet functional needs of wheelchair users featuring many sizes and easy attachment hardware
- Designed for moderate to high postural needs
- Provides moderate contouring for clients that would benefit for more lateral support

J3PLSHS14FS Jay 3 PL SH Back Support 4-Point Hardware Short 14"

Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm
Overall Width	355 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PLSHS16FS Jay 3 PL SH Back Support 4-Point Hardware Short 16"

Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm
Overall Width	405 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PLSHS18FS Jay 3 PL SH Back Support 4-Point Hardware Short 18"

Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm
Overall Width	455 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PLSHS20FS Jay 3 PL SH Back Support 4-Point Hardware Short 20"

Conto	ur Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm
Overal	ll Width	505 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PLSHM16FS Jay 3 PL SH Back Support 4-Point Hardware Medium 16"

	-							
Contour Depth		76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm
Overall Width		405 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PLSHM18FS Jay 3	PL SH B	ack Support 4-Point	t Hardwa	are Medium 18"			
Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm
Overall Width	455 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg
J3PLSHM2OFS Jay :	3 PL SH E	Back Support 4-Poin	t Hardw	are Medium 20"			
Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm
					0. 1	Dreduct Weight	0.0 1/2
Overall Width	505 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg
		Package Quantity	Hardwa		Single	Product weight	2.3 Ky
			Hardwa		Single 225 kg	Min Overall Height	2.3 kg 305 mm
J3PLSHT18FS Jay 3	PL SH Ba	ack Support 4-Point		re Tall 18"			
J3PLSHT18FS Jay 3 Contour Depth Overall Width	PL SH Ba 76 mm 455 mm	ack Support 4-Point Max Overall Height	610 mm 1	re Tall 18" Max User Weight Package Type	225 kg	Min Overall Height	305 mm
J3PLSHT18FS Jay 3 Contour Depth Overall Width	PL SH Ba 76 mm 455 mm	ack Support 4-Point Max Overall Height Package Quantity	610 mm 1	re Tall 18" Max User Weight Package Type	225 kg	Min Overall Height	305 mm

Jay 3 Posterior Deep Lateral Shoulder Height Back Support, 4-Point Hardware



Contour Overall V

• Posterior and Lateral support through the Upper Thoracic support position

- hardware • The Jay J3 is designed for moderate to high postural needs

J3PDLSHS14FS	Jay 3 PDL SH Back Support 4-Point Hardw
--------------	---

ır Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
l Width	355 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PDLSHS16FS Jay 3 PDL SH Back Support 4-Point Hardware Short 16"

Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
Overall Width	405 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
Overall Width	455 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PDLSHS20FS Jay 3 PDL SH Back Support 4-Point Hardware Short 20"

Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
Overall Width	505 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PDLSHM16FS Jay 3 PDL SH Back Support 4-Point Hardware Medium 16"

Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
Overall Width	405 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PDLSHM20FS Jay 3 PDL SH Back Support 4-Point Hardware Medium 20"

Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
Overall Width	505 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg



Available in a variety of sizes, levels of support, contour depths and mounting hardware

• This comprehensive back range meets functional needs of wheelchair users featuring many sizes and easy attachment

• The posterior lateral back provides moderate contouring for clients that would benefit for more lateral support

are Short 14"

ware Short 18"

J3PDLSHT18FS Jay 3 PDL SH Back Support 4-Point Hardware Tall 18"										
contour Depth 15	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm			
verall Width 45	455 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg			

J3PDLSHT20FS Jay 3 PDL SH Back Support 4-Point Hardware Tall 20"

Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
Overall Width	505 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

Jay 3 Posterior Deep Lateral Upper Thoracic Back Support, 4-Point Hardware



- Posterior and Lateral support designed to support the back through the Upper Thoracic support position
- Available in a variety of sizes, levels of support, contour depths and mounting hardware
- Comprehensive back range to meet functional needs of wheelchair users featuring many sizes and easy attachment
- hardware
- Designed for moderate to high postural needs Provides moderate contouring for clients that would benefit for more lateral support

J3PDLUTS14FS Jay 3 PDL UT Back Support with 4-Point Hardware Short 14"

Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
Overall Width	355 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PDLUTS16FS Jay 3 PDL UT Back Support with 4-Point Hardware Short 16"

Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
Overall Width	405 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PDLUTS18FS Jay 3 PDL UT Back Support with 4-Point Hardware Short 18"

Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
Overall Width	455 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PDLUTS20FS Jay 3 PDL UT Back Support with 4-Point Hardware Short 20"

Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
Overall Width	505 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PDLUTM14FS Jay 3 PDL UT Back Support with 4-Point Hardware Medium 14"

Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
Overall Width	355 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PDLUTM16FS Jay 3 PDL UT Back Support with 4-Point Hardware Medium 16"

Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
Overall Width	405 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PDLUTM18FS Jay 3 PDL UT Back Support with 4-Point Hardware Medium 18"

Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
Overall Width	455 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PDLUTM20FS Jay 3 PDL UT Back Support with 4-Point Hardware Medium 20"

Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
Overall Width	505 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PDLUTT14FS	Jay 3	B PDL UT	Back Support with 4	4-Point H
Contour Denth		152 mm	Max Overall Height	610 mm

Overall Width 355 mm Package Quantity 1 Package Type Single Product Weight 23 kg	Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
	Overall Width	355 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PDLUTT16FS Jay 3 PDL UT Back Support with 4-Point Hardware Tall 16"

Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
Overall Width	405 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PDLUTT18FS Jay 3 PDL UT Back Support with 4-Point Hardware Tall 18"

Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
Overall Width	455 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PDLUTT20FS Jay 3 PDL UT Back Support with 4-Point Hardware Tall 20"

Contour Depth	152 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	420 mm
Overall Width	505 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

Jay 3 Posterior Lateral Upper Thoracic Back Support, 4-Point Hardware



- Designed to support the back through the Upper Thoracic support position
- The J3 comprehensive back range meets functional needs of wheelchair users
 - The J3 is designed for moderate to high postural needs

J3PLUTS12FS Jay 3	J3PLUTS12FS Jay 3 PL UP Back Support with 4-Point Hardware Short 12"												
Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm						
Overall Width	305 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg						
J3PLUTS14FS Jay 3 PL UP Back Support with 4-Point Hardware Short 14"													
Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm						
Overall Width	355 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg						

J3PLUTS16FS Jay 3 PL UP Back Support with 4-Point Hardware Short 16" Contour Depth 76 mm Max Overall Height 610 mm Overall Width 405 mm Package Quantity

J3PLUTS18FS Jay 3 PL UP Back Support with 4-Point Hardware Short 18"

Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm
Overall Width	455 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PLUTS20FS Jay 3 PL UP Back Support with 4-Point Hardware Short 20"

Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm
Overall Width	505 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PLUTM14FS Jay 3 PL UP Back Support with 4-Point Hardware Medium 14"

Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm
Overall Width	355 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg



Hardware Tall 14"

• Available in a variety of sizes, levels of support, contour depths and mounting hardware

Max User Weight	225 kg	Min Overall Height	305 mm
Package Type	Single	Product Weight	2.3 kg

J3PLUTM16FS Jay 3 PL UP Back Support with 4-Point Hardware Medium 16"

Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm
Overall Width	405 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PLUTM18FS Jay 3 PL UP Back Support with 4-Point Hardware Medium 18"

Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm
Overall Width	455 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PLUTM20FS Jay 3 PL UP Back Support with 4-Point Hardware Medium 20"

Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm
Overall Width	505 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PLUTT14FS Jay 3 PL UP Back Support with 4-Point Hardware Tall 14"

Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm
Overall Width	355 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PLUTT16FS Jay 3 PL UP Back Support with 4-Point Hardware Tall 16"

Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm
Overall Width	405 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PLUTT18FS Jay 3 PL UP Back Support with 4-Point Hardware Tall 18"

Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm
Overall Width	455 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

J3PLUTT20FS Jay 3 PL UP Back Support with 4-Point Hardware Tall 20"

Contour Depth	76 mm	Max Overall Height	610 mm	Max User Weight	225 kg	Min Overall Height	305 mm
Overall Width	505 mm	Package Quantity	1	Package Type	Single	Product Weight	2.3 kg

Jay 3 Backrest - Posterior Lower Thoracic



- Provides outstanding stability and optimal pressure distribution for the client at high risk of skin breakdown
- Offers the choice of Jay flow[™] fluid or Roho Dry Floatation® air inserts, which conform to each individual's shape Adjust to sitting positions, and ensure proper fluid placement beneath bony prominences
- The J3 features a contoured foam base constructed from closed cell foam
- The J3's customizable pelvic loading area sizes are selected separately from the wheelchair cushion width
- This provides ideal forward and lateral stability, optimal immersion, and weight distribution that accommodate a broad range of client needs

Contour Depth	55 mm	Max User Weight	135 kg	Overall Height	205 mm	Overall Width	305 mm
I 3PALTM14SH Jay 3	Backre	st - Posterior Lower	Thoracic	Medium 14"			
Contour Depth	55 mm	Max User Weight	135 kg	Overall Height	205 mm	Overall Width	355 mm
,	r						
,	r	st - Posterior Lower		Medium 16"			
Contour Depth	55 mm	Max User Weight	135 kg	Overall Height	205 mm	Overall Width	405 mm
		Max User Weight			205 mm	Overall Width	405 mm
					205 mm 205 mm	Overall Width Overall Width	405 mm 455 mm
J3PALTM18SH Jay 3 Contour Depth	Backres	st - Posterior Lower	Thoracic	Medium 18" Overall Height			

	Juyo	Backres	t - Posterior Lower T	horacic	Short 16"			
Contour Depth		55 mm	Max User Weight	135 kg	Overall Height	165 mm	Overall Width	405 mm
I3PALTS18SH	Jay 3	Backres	t - Posterior Lower T	horacic	Short 18"			
Contour Depth		55 mm	Max User Weight	135 kg	Overall Height	165 mm	Overall Width	455 mm
I3PALTT12SH	Jay 3	Backres	t - Posterior Lower T	horacic	Tall 12"			
Contour Depth		55 mm	Max User Weight	135 kg	Overall Height	240 mm	Overall Width	305 mm
	Jay 3		t - Posterior Lower T					
Contour Depth		55 mm	Max User Weight	135 kg	Overall Height	240 mm	Overall Width	355 mm
I3PALTT16SH	Jay 3	Backres	t - Posterior Lower T	horacic	Tall 16"			
Contour Depth		55 mm	Max User Weight	135 kg	Overall Height	240 mm	Overall Width	405 mm
	Jay 3 I		Max User Weight t - Posterior Lower T			240 mm	Overall Width	405 mm
	Jay 3 I					240 mm 240 mm	Overall Width Overall Width	405 mm 455 mm
13PALTT18SH Contour Depth		Backres	t - Posterior Lower T	horacic	Tall 18" Overall Height			

Jay 3 Backrest - Posterior Shoulder Height



- of client needs

Contour Depth	55 mm	Max User Weight	135 kg	Overall Height	570 mm	Overall Width	405 mm
						<u></u>	
3PASHM18FS Jay	/ 3 Backre	est - Posterior Should	der Heigl	nt Medium 18"			
Contour Depth	55 mm	Max User Weight	135 kg	Overall Height	570 mm	Overall Width	455 mm
I 3PASHS14FS Jay	3 Backre	st - Posterior Should	ler Heigh	it Short 14"			
Contour Depth	55 mm	Max User Weight	135 kg	Overall Height	480 mm	Overall Width	355 mm
I 3PASHS16FS Jay	3 Backre	st - Posterior Should	ler Heigh	it Short 16"			
I3PASHS16FS Jay	3 Backre	st - Posterior Should	ler Heigh	nt Short 16"			
Contour Depth	55 mm	Max User Weight	135 kg	Overall Height	480 mm	Overall Width	405 mm
I3PASHS18FS Jay	3 Backre	st - Posterior Should	ler Heigh	it Short 18"			
I3PASHS18FS Jay Contour Depth	3 Backre	st - Posterior Should	ler Heigh 135 kg	t Short 18" Overall Height	480 mm	Overall Width	455 mm
Contour Depth	55 mm		135 kg	Overall Height	480 mm	Overall Width	455 mm
Contour Depth	55 mm	Max User Weight	135 kg	Overall Height	480 mm 610 mm	Overall Width Overall Width	455 mm 455 mm
Contour Depth I3PASHT18FS Jay Contour Depth	55 mm 3 Backre 55 mm	Max User Weight st - Posterior Should	135 kg ler Heigh 135 kg	Overall Height t Tall 18" Overall Height			

Mobility

• Provides outstanding stability and optimal pressure distribution for the client at high risk of skin breakdown

• Offers the choice of Jay flow[™] fluid or Roho Dry Floatation® air inserts, which conform to each individual's shape

Adjust to sitting positions, and ensure proper fluid placement beneath bony prominences

• The J3 features a contoured foam base constructed from closed cell foam

• The J3's customizable pelvic loading area sizes are selected separately from the wheelchair cushion width

• This provides ideal forward and lateral stability, optimal immersion, and weight distribution that accommodate a broad range

J3PASHS20FS Jay 3 Backrest - Posterior Shoulder Height Short 20"									
Contour Depth	55 mm	Max User Weight	135 kg	Overall Height	480 mm	Overall Width	505 mm		
J3PASHT2OFS Jay 3	Backres	st - Posterior Should	ler Heigh	it Tall 20"					
Contour Depth	55 mm	Max User Weight	135 kg	Overall Height	610 mm	Overall Width	505 mm		

Jay 3 Backrest - Posterior & Deep Contour Lower Thoracic

	 Provides outstanding stability and optimal pressure distribution for the client at high risk of skin breakdown Offers the choice of Jay flow[™] fluid or Roho Dry Floatation® air inserts, which conform to each individual's shape Adjust to sitting positions, and ensure proper fluid placement beneath bony prominences The J3 features a contoured foam base constructed from closed cell foam The J3's customizable pelvic loading area sizes are selected separately from the wheelchair cushion width This provides ideal forward and lateral stability, optimal immersion, and weight distribution that accommodate a broad range of client needs
13000	154 Jay 2 Postariar, Doop Contour Lower Therasis Medium 14"

J3PUCLIMI4SH Jay	3 Postei	rior, Deep Contour Lo	ower Tho	pracic Medium 14"			
Contour Depth	150 mm	Max User Weight	135 kg	Overall Height	205 mm	Overall Width	355 mm
3PDCLTM16SH Jay	3 Postei	rior, Deep Contour Lo	ower Tho	oracic Medium 16"			
Contour Depth	150 mm	Max User Weight	135 kg	Overall Height	205 mm	Overall Width	405 mm
I 3PDCLTM18SH Jay	3 Postei	rior, Deep Contour Lo	ower Tho	oracic Medium 18"			
Contour Depth	150 mm	Max User Weight	135 kg	Overall Height	205 mm	Overall Width	455 mm
3PDCLTT14SH Jay : Contour Depth	B Posteri	or, Deep Contour Lo Max User Weight	Wer Thoi 135 kg	racic Tall 14" Overall Height	240 mm	Overall Width	355 mm
I 3PDCLTT16SH Jay 3	3 Posteri	or, Deep Contour Lo	wer Tho	racic Tall 16"			
Contour Depth	150 mm	Max User Weight	135 kg	Overall Height	240 mm	Overall Width	405 mm
3PDCLTT18SH Jay :	3 Posteri	or, Deep Contour Lo	wer Tho	racic Tall 18"			
Contour Depth	150 mm	Max User Weight	135 kg	Overall Height	240 mm	Overall Width	455 mm
	3 Postor	ior, Deep Contour Lo	wer Tho	racic Tall 20"			
Japucliizuan Jay	5103101						

Jay 3 Backrest - Posterior & Deep Contour Mid Thoracic



- Adjust to sitting positions, and ensure proper fluid placement beneath bony prominences
- The J3 features a contoured foam base constructed from closed cell foam
- of client needs

J3PDCMTS14SH	Jay 3	Poste	rior, Deep	Contour N	lid Thora	acic Short 14"			
Contour Depth	15	0 mm	Max User Weig	ht	135 kg	Overall Height	305 mm	Overall Width	355 mm
J3PDCMTS16SH	Jay 3	Poste	rior, Deep	Contour N	lid Thora	acic Short 16″			
Contour Depth	15	0 mm	Max User Weig	ht	135 kg	Overall Height	305 mm	Overall Width	405 mm
J3PDCMTS18SH	Jay 3	Poste	rior, Deep	Contour N	lid Thora	acic Short 18″			
Contour Depth	15	0 mm	Max User Weig	ht	135 kg	Overall Height	305 mm	Overall Width	455 mm
J3PDCMTT14SH	Jay 3	Postei	rior, Deep	Contour M	lid Thora	acic Tall 14"			
Contour Depth	15	0 mm	Max User Weig	ht	135 kg	Overall Height	380 mm	Overall Width	355 mm
J3PDCMTT16SH	Jay 3	Poster	ior, Deep	Contour M	lid Thora	acic Tall 16"			
Contour Depth	15	0 mm	Max User Weig	ht	135 kg	Overall Height	380 mm	Overall Width	405 mm
J3PDCMTT18SH	Jay 3	Poster	ior, Deep	Contour M	lid Thora	acic Tall 18"			
Contour Depth	15	0 mm	Max User Weig	ht	135 kg	Overall Height	380 mm	Overall Width	455 mm
J3PDCMTTXXSH	Jay 3	Poste	rior, Deep	Contour N	lid Thora	acic Tall 14"-20"			
Contour Depth		150 r	nm	Max User Weight	t	135 kg	Overall Height		380 mm
J3PDCMTS20SH	Jay 3	Poste	rior, Deep	Contour N	Aid Thor	acic Short 20"			
Contour Depth	15	0 mm	Max User Weig	ht	135 kg	Overall Height	305 mm	Overall Width	505 mm
J3PDCMTT20SH	Jay 3	Poste	rior, Deep	Contour N	/lid Thora	acic Tall 20"			
Contour Depth	15	0 mm	Max User Weig		135 kg	Overall Height	380 mm	Overall Width	505 mm



- Provides outstanding stability and optimal pressure distribution for the client at high risk of skin breakdown
- Offers the choice of Jay flow[™] fluid or Roho Dry Floatation® air inserts, which conform to each individual's shape
- The J3's customizable pelvic loading area sizes are selected separately from the wheelchair cushion width
 This provides ideal forward and lateral stability, optimal immersion, and weight distribution that accommodate a broad range



Jay 3 Backrest - Posterior & Deep Contour Shoulder Height

Jay 3 Backrest - Posterior & Deep Contour Upper Thoracic



- Provides outstanding stability and optimal pressure distribution for the client at high risk of skin breakdown
- Offers the choice of Jay flow[™] fluid or Roho Dry Floatation® air inserts, which conform to each individual's shape
- Adjust to sitting positions, and ensure proper fluid placement beneath bony prominences The J3 features a contoured foam base constructed from closed cell foam
- The J3's customizable pelvic loading area sizes are selected separately from the wheelchair cushion width
- This provides ideal forward and lateral stability, optimal immersion, and weight distribution that accommodate a broad range of client needs

Contour Depth		150 mm	Max User Weight	135 kg	Overall Height	480 mm	Overall Width	405 mm		
J3PDCSHS18FS Jay 3 Posterior, Deep Contour Shoulder Height Short 18"										
Contour Depth		150 mm	Max User Weight	135 ka	Overall Height	480 mm	Overall Width	455 mm		
3PDCSHSXXFS	lav							100 1111		
3PDCSHSXXFS	Jay		rior, Deep Contour S				Min Overall Width	355 mm		



- Provides outstanding stability and optimal pressure distribution for the client at high risk of skin breakdown
- Offers the choice of Jay flow[™] fluid or Roho Dry Floatation® air inserts, which conform to each individual's shape
- Adjust to sitting positions, and ensure proper fluid placement beneath bony prominences
- The J3 features a contoured foam base constructed from closed cell foam
- The J3's customizable pelvic loading area sizes are selected separately from the wheelchair cushion width
- This provides ideal forward and lateral stability, optimal immersion, and weight distribution that accommodate a broad range of client needs

Contour Depth	150 mm	Max User Weight	135 kg	Overall Height	455 mm	Overall Width	355 mm
					·		
I3PDCUTM16FS Ja	y 3 Poste	rior, Deep Contour L	Ipper Tho	pracic Medium 16"			
Contour Depth	150 mm	Max User Weight	135 kg	Overall Height	455 mm	Overall Width	405 mm
3PDCUTM18FS Ja	y 3 Poste	rior, Deep Contour L	Ipper Tho	oracic Medium 18"			
J 3PDCUTM18FS Ja	y 3 Poste	rior, Deep Contour L	Ipper Tho	oracic Medium 18"			
Contour Depth	150 mm	Max User Weight	135 kg	Overall Height	455 mm	Overall Width	455 mm
3PDCUTMXXFS Ja	y 3 Poste	erior, Deep Contour l	Jpper The	pracic Medium 14"-2	0"		
		Mary Oscaral Middle	505 mm	Max User Weight	135 kg	Min Overall Width	355 mm
Contour Depth	150 mm	Max Overall Width	000 11111	······			
	150 mm 455 mm	Max overall width					
Contour Depth		Max uverall width					
Contour Depth Overall Height	455 mm	erior, Deep Contour					

Jay 3 Backrest - Posterior Mid Thoracic



- · Adjust to sitting positions, and ensure proper fluid placement beneath bony prominences

- of client needs

J3PAMTS22FS	Jay 3	Backre	st - Posterior Mid Th	oracic S	h
Contour Depth		55 mm	Max User Weight	225 kg	0
J3PAMTS24FS	Jay 3	Backres	st - Posterior Mid Th	oracic S	h
Contour Depth		55 mm	Max User Weight	225 kg	C
J3PAMTS26FS	Jay 3	Backre	st - Posterior Mid Th	oracic S	h
Contour Depth		55 mm	Max User Weight	225 kg	0
J3PAMTT22FS	Jay 3	Backres	st - Posterior Mid Th	oracic Ta	all
Contour Depth		55 mm	Max User Weight	225 kg	0
J3PAMTT24FS	Jay 3	Backres	st - Posterior Mid Th	oracic Ta	all
Contour Depth		55 mm	Max User Weight	225 kg	0
J3PAMTT26FS	Jay 3	Backres	st - Posterior Mid Th	oracic Ta	all
Contour Depth		55 mm	Max User Weight	225 kg	0

Roho Agility Cushion Mid Contour with Air



- Features Roho air floatation technology · Allows individuals to adjust air volume to their desired comfort level
- Supports postural changes



R-AG3-1410AL-1	Roh	Roho Agility Cushion with Air 14" x 10"							
Overall Depth		255 mm	Overall Width	355 mm	Pa				
R-AG3-1413AL-1	Roh	o Agility	Cushion with Air 14	" x 13"					
Overall Depth		330 mm	Overall Width	355 mm	Pa				
R-AG3-1416AL-1	Roh	o Agility	Cushion with Air 14	" x 16"					
Overall Depth		405 mm	Overall Width	355 mm	Pa				
R-AG3-1418AL-1	Roh	o Agility	Cushion with Air 14	" x 18"					
Overall Depth		455 mm	Overall Width	355 mm	Pa				



• Provides outstanding stability and optimal pressure distribution for the client at high risk of skin breakdown

• Offers the choice of Jay flow™ fluid or Roho Dry Floatation® air inserts, which conform to each individual's shape

The J3 features a contoured foam base constructed from closed cell foam

• The J3's customizable pelvic loading area sizes are selected separately from the wheelchair cushion width

• This provides ideal forward and lateral stability, optimal immersion, and weight distribution that accommodate a broad range

ort 22"					
Overall Height	305 mm	Overall Width	560 mm		
ort 24"					
Overall Height	305 mm	Overall Width	610 mm		
			· · · · · · · · · · · · · · · · · · ·		
ort 26"					
Overall Height	305 mm	Overall Width	660 mm		
ll 22"					
Overall Height	380 mm	Overall Width	560 mm		
l 24"					
Overall Height	380 mm	Overall Width	610 mm		
		· · · · · · · · · · · · · · · · · · ·			3
ll 26″					
Overall Height	380 mm	Overall Width	560 mm		

• The perfect blend of innovative features that give wheelchair users comfortable support without sacrificing function

• The back shell can be flipped to accommodate wider torsos and broader shoulders

• Wide range of tilt angle, depth, width and vertical adjustment capabilities to address specific positioning needs

• Full range of contours for varying postural support and positioning needs

Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single

R-AG3-1610AL-1	Roho Ag	gility Cushion w	ith Air 16" x 10"				
Overall Depth	255 m	om Overall Width	405 mm	Package Quantity	1	Package Type	Single
R-AG3-1613AL-1	Roho Ag	jility Cushion w	ith Air 16" x 13"				
Overall Depth	330 m	owerall Width	405 mm	Package Quantity	1	Package Type	Single
R-AG3-1616AL-1	Roho Ag	jility Cushion w	ith Air 16" x 16"				
Overall Depth	405 m	overall Width	405 mm	Package Quantity	1	Package Type	Single
R-AG3-1618AL-1	Roho Ag	jility Cushion w	ith Air 16" x 18"				
Overall Depth	455 m	overall Width	405 mm	Package Quantity	1	Package Type	Single
R-AG3-1810AL-1	Roho Ag	jility Cushion w	ith Air 18" x 10"				
Overall Depth	255 m	owerall Width	455 mm	Package Quantity	1	Package Type	Single
R-AG3-1813AL-1	Roho Ag	ility Cushion w	ith Air 18" x 13"				
Overall Depth	330 m	overall Width	455 mm	Package Quantity	1	Package Type	Single
R-AG3-1816AL-1	Roho Ag	jility Cushion w	ith Air 18" x 16"				
Overall Depth	405 m	m Overall Width	455 mm	Package Quantity	1	Package Type	Single
R-AG3-1818AL-1	Roho Ag	jility Cushion w	ith Air 18" x 18"				
Overall Depth	455 m	m Overall Width	455 mm	Package Quantity	1	Package Type	Single
R-AG3-2010AL-1	Roho A	gility Cushion w	vith Air 20" x 10"				
Overall Depth	255 m	om Overall Width	505 mm	Package Quantity	1	Package Type	Single
			uith Air 20" y 12"				
R-AG3-2013AL-1	Roho Ag	gility Cushion w					
R-AG3-2013AL-1 Overall Depth	Roho Ag		505 mm	Package Quantity	1	Package Type	Single
Overall Depth	330 m	overall Width		Package Quantity	1	Package Type	Single
Overall Depth	330 m	overall Width	505 mm	Package Quantity Package Quantity	1	Package Type Package Type	Single Single
Overall Depth R-AG3-2016AL-1 Overall Depth	330 m Roho Ag 405 m Roho Ag	Overall Width gility Cushion w Overall Width gility Cushion w	505 mm rith Air 20" x 16" 505 mm rith Air 20" x 18"	Package Quantity			Single
Overall Depth R-AG3-2016AL-1 Overall Depth	330 m Roho Ag 405 m Roho Ag	Overall Width gility Cushion w m Overall Width	505 mm rith Air 20" x 16" 505 mm				
R-AG3-2016AL-1 Overall Depth R-AG3-2018AL-1 Overall Depth	330 m Roho Ag 405 m Roho Ag 455 m	Overall Width gility Cushion w um Overall Width gility Cushion w um Overall Width	505 mm rith Air 20" x 16" 505 mm rith Air 20" x 18"	Package Quantity		Package Type	Single
Overall Depth R-AG3-2016AL-1 Overall Depth R-AG3-2018AL-1 Overall Depth	330 m Roho Ag 405 m Roho Ag 455 m	Overall Width gility Cushion w Overall Width gility Cushion w Overall Width gility Cushion w gility Cushion w	505 mm with Air 20" x 16" 505 mm with Air 20" x 18" 505 mm	Package Quantity		Package Type	Single
Overall Depth R-AG3-2016AL-1 Overall Depth R-AG3-2018AL-1 Overall Depth R-AG3-1420AL-1 Overall Depth	330 m Roho Ag 405 m Roho Ag 455 m Roho Ag 505 m	Overall Width gility Cushion w overall Width gility Cushion w overall Width gility Cushion w gility Cushion w overall Width	505 mm rith Air 20" x 16" 505 mm rith Air 20" x 18" 505 mm rith Air 14" x 20"	Package Quantity Package Quantity	1	Package Type Package Type	Single Single
Overall Depth R-AG3-2016AL-1 Overall Depth R-AG3-2018AL-1 Overall Depth R-AG3-1420AL-1 Overall Depth	330 m Roho Ag 405 m Roho Ag 455 m Roho Ag 505 m	Overall Width gility Cushion w m Overall Width	505 mm rith Air 20" x 16" 505 mm rith Air 20" x 18" 505 mm rith Air 14" x 20" 355 mm	Package Quantity Package Quantity	1	Package Type Package Type	Single Single
Overall Depth R-AG3-2016AL-1 Overall Depth R-AG3-2018AL-1 Overall Depth R-AG3-1420AL-1 Overall Depth R-AG3-1620AL-1 Overall Depth	330 m Roho Ag 405 m Roho Ag 505 m Roho Ag 505 m Roho Ag	Overall Width gility Cushion w m Overall Width gility Cushion w m Overall Width gility Cushion w m Overall Width gility Cushion w gility Cushion w m Overall Width gility Cushion w m Overall Width gility Cushion w overall Width	505 mm rith Air 20" x 16" 505 mm rith Air 20" x 18" 505 mm rith Air 14" x 20" 355 mm rith Air 16" x 20"	Package Quantity Package Quantity Package Quantity		Package Type Package Type Package Type	Single Single Single
Overall Depth R-AG3-2016AL-1 Overall Depth R-AG3-2018AL-1 Overall Depth R-AG3-1420AL-1 Overall Depth R-AG3-1620AL-1 Overall Depth	330 m Roho Ag 405 m Roho Ag 505 m Roho Ag 505 m Roho Ag	Overall Width gility Cushion w overall Width	505 mm rith Air 20" x 16" 505 mm rith Air 20" x 18" 505 mm rith Air 14" x 20" 355 mm rith Air 16" x 20" 405 mm	Package Quantity Package Quantity Package Quantity		Package Type Package Type Package Type	Single Single Single
Overall Depth R-AG3-2016AL-1 Overall Depth R-AG3-2018AL-1 Overall Depth R-AG3-1420AL-1 Overall Depth R-AG3-1620AL-1 Overall Depth R-AG3-1820AL-1 Overall Depth	330 m Roho Ag 405 m Roho Ag 455 m Roho Ag 505 m Soft m	Overall Width gility Cushion w m Overall Width gility Cushion w overall Width	505 mm vith Air 20" x 16" 505 mm vith Air 20" x 18" 505 mm vith Air 14" x 20" 355 mm vith Air 16" x 20" 405 mm vith Air 18" x 20"	Package Quantity Package Quantity Package Quantity Package Quantity Package Quantity Package Quantity		Package Type Package Type Package Type Package Type Package Type	Single Single Single Single
Overall Depth R-AG3-2016AL-1 Overall Depth R-AG3-2018AL-1 Overall Depth R-AG3-1420AL-1 Overall Depth R-AG3-1620AL-1 Overall Depth R-AG3-1820AL-1 Overall Depth	330 m Roho Ag 405 m Roho Ag 455 m Roho Ag 505 m Soft m	Overall Width gility Cushion w overall Width	505 mm rith Air 20" x 16" 505 mm rith Air 20" x 18" 505 mm rith Air 14" x 20" 355 mm rith Air 16" x 20" 405 mm rith Air 18" x 20" 455 mm	Package Quantity Package Quantity Package Quantity Package Quantity Package Quantity Package Quantity		Package Type Package Type Package Type Package Type Package Type	Single Single Single Single

Roho Agility Cushion Max Contour with Air



- Allows individuals to adjust air volume to their desired comfort level

R-AG6-1416AL-1	Roh	o Agility	Cushion with Air 14	" x 16"	
Overall Depth		405 mm	Overall Width	355 mm	
R-AG6-1418AL-1	Roh	o Agility	Cushion with Air 14	" x 18"	
Overall Depth		455 mm	Overall Width	355 mm	1
R-AG6-1420AL-1	Rol	no Agility	Cushion with Air 14	l" x 20"	
Overall Depth		505 mm	Overall Width	355 mm	
R-AG6-1616AL-1	Roh	o Agility	Cushion with Air 16	" x 16"	
Overall Depth		405 mm	Overall Width	405 mm	
R-AG6-1618AL-1	Roh	o Agility	Cushion with Air 16	" x 18"	
Overall Depth		455 mm	Overall Width	405 mm	
R-AG6-1620AL-1	Rol	no Agility	/ Cushion with Air 16	5" x 20"	
Overall Depth		505 mm	Overall Width	405 mm	
R-AG6-1816AL-1	Roh	o Agility	Cushion with Air 18	" x 16"	
Overall Depth		405 mm	Overall Width	455 mm	
R-AG6-1818AL-1	Roh	o Agility	Cushion with Air 18	" x 18"	
Overall Depth		455 mm	Overall Width	455 mm	
R-AG6-1820AL-1	Roł	no Agility	Cushion with Air 18	8" x 20"	
Overall Depth		505 mm	Overall Width	455 mm	
R-AG6-2016AL-1	Roł	no Agility	Cushion with Air 20)" x 16"	
Overall Depth		405 mm	Overall Width	505 mm	
R-AG6-2018AL-1	Rol	no Agility	Cushion with Air 20)" x 18"	
Overall Depth		455 mm	Overall Width	505 mm	
R-AG6-2020AL-1	Ro	ho Agilit	y Cushion with Air 2	0" x 20"	
Overall Depth		505 mm	Overall Width	505 mm	



Features Roho air floatation technology
The perfect blend of innovative features that give wheelchair users comfortable support without sacrificing function

Allows individuals to adjust an volume to their desired connectence.
Supports postural changes
The back shell can be flipped to accomodate wider torsos and broader shoulders
Wide range of tilt angle, depth, width and vertical adjustment capabilities to address specific positioning needs
Full range of contours for varying postural support and positioning needs

Package Quantity	1	Package Type	Single	
Package Quantity	1	Package Type	Single	
Package Quantity	1	Package Type	Single	
Package Quantity	1	Package Type	Single	
Package Quantity	1	Package Type	Single	
Package Quantity	1	Package Type	Single	
Package Quantity	1	Package Type	Single	
Package Quantity	1	Package Type	Single	
i donage quantity		r uokuge type	olingic	
Package Quantity	1	Package Type	Single	
Package Quantity	1	Package Type	Single	
Package Quantity	1	Package Type	Single	
Package Quantity	1	Package Type	Single	
r annage éngulary		i uonage ijpe	Unigic	

Roho Agility Cushion Min Contour with Air



- Features Roho air floatation technology
 The perfect blend of innovative features that give wheelchair users comfortable support without sacrificing function
 Allows individuals to adjust air volume to their desired comfort level
- Supports postural changes
- The back shell can be flipped to accomodate wider torsos and broader shoulders
 Wide range of tilt angle, depth, width and vertical adjustment capabilities to address specific positioning needs
- Full range of contours for varying postural support and positioning needs

	ROU	o Agility	Cushion with Air 14	" x 10"				
Overall Depth		255 mm	Overall Width	355 mm	Package Quantity	1	Package Type	Single
R-AG2-1413AL-1	Roh	o Agility	Cushion with Air 14	" x 13"				
Overall Depth		330 mm	Overall Width	355 mm	Package Quantity	1	Package Type	Single
R-AG2-1416AL-1	Roh	o Agility	Cushion with Air 14	" x 16"				
Overall Depth		405 mm	Overall Width	355 mm	Package Quantity	1	Package Type	Single
R-AG2-1418AL-1	Roh	o Agility	Cushion with Air 14	" x 18"				
Overall Depth		455 mm	Overall Width	355 mm	Package Quantity	1	Package Type	Single
R-AG2-1610AL-1	Roh	o Agility	Cushion with Air 16	" x 10"				
Overall Depth		255 mm	Overall Width	405 mm	Package Quantity	1	Package Type	Single
R-AG2-1613AL-1	Roh	o Agility	Cushion with Air 16	" x 13"				
Overall Depth		330 mm	Overall Width	405 mm	Package Quantity	1	Package Type	Single
R-AG2-1616AL-1	Roh	o Agility	Cushion with Air 16	" x 16"				
Overall Depth		405 mm	Overall Width	405 mm	Package Quantity	1	Package Type	Single
R-AG2-1618AL-1	Roh	o Agility	Cushion with Air 16	" x 18"				
Overall Depth		455 mm	Overall Width	405 mm	Package Quantity	1	Package Type	Single
R-AG2-1810AL-1	Roh	o Agility	Cushion with Air 18	" x 10"				
Overall Depth		255 mm	Overall Width	455 mm	Package Quantity	1	Package Type	Single
R-AG2-1813AL-1	Roh	o Agility	Cushion with Air 18	" x 13"				
Overall Depth		330 mm	Overall Width	455 mm	Package Quantity	1	Package Type	Single
R-AG2-1816AL-1	Roh	o Agility	Cushion with Air 18	" x 16"				
Overall Depth		405 mm	Overall Width	455 mm	Package Quantity	1	Package Type	Single
	Roh	o Agility	Cushion with Air 18	" x 18"				
R-AG2-1818AL-1				455	Package Quantity	1	Deelesse Ture	Single
R-AG2-1818AL-1 Overall Depth		455 mm	Overall Width	455 mm	Fackage quantity		Package Type	Jillyle
R-AG2-1818AL-1 Overall Depth R-AG2-2010AL-1	Roh		Cushion with Air 20		rackaye quantiy		Раскаде туре	Jingie

R-AG2-2013AL-1	Rol	no Agility	/ Cushion with Air 20)" x 13"				
Overall Depth		330 mm	Overall Width	505 mm				
R-AG2-2016AL-1	AG2-2016AL-1 Roho Agility Cushion with Air 20" x 10							
Overall Depth		405 mm	Overall Width	505 mm				

455 mm Overall Width

505 mm

Overall Depth



Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single

National Distribution, Local Service

SHOWERING FURNITURE



The Premium Padded Shower Stool has been Australian Standard Approved to 250kg. It is manufactured out of 32mm aluminium tube, meeting the Australian standards for handrail tubing. Larger tubing makes it easier to grip and stand yourself up. Aluminium frame is heat treated for extra strength.

PREMIUM PADDED SHOWER STOOL

- 250kg weight capacity
 - Australian Standards Approved PU padded seat

TUBALCO ADULT TILTING TRANSIT

MOBILE SHOWER COMMODE

Australian made

face operation

anxiety

Swingaway footrest

Adjustable footrest hangers

Headrest

Manufactured 32mm aluminium • Height adjustable legs

Loop locking armrests with side guards

Remote tilt control allowing face to

Smooth tilting action reduces user

Corrosion resistant components

tube

FOR FULL

PRODUCT DETAILS

87



National Distribution, Local Service

Shower Furniture

R & R Healthcare Equipment Economy Shower Chair



 Manufactured 22mm steel tube Zinc dipped for rust proofing

- Height adjustable legs
 - Australian Standards tested to type B corrosion resistance

12028S R & R Healt							
Max Overall Height	730 mm	Max Seat Height	550 mm	Max User Weight	110 kg	Min Overall Height	610 mm
Min Seat Height	430 mm	Seat Depth	480 mm	Seat Width	460 mm		

Shower Furniture

R & R Healthcare Equipment Premium Padded Shower Stool



- Height adjustable legs • Enhanced with a PU padded seat for extra comfort Manufactured from 32mm aluminium tube · Powder coated white in colour
 - Meets the Australian Standards for handrail tubing · Larger tubing increases users ablility to grip and stand up easily
 - Available with various seat widths: 45cm, 52cm & 60cm
- 12325A45 R & R Healthcare Equipment Premium Padded Shower Stool 450mm

Max	x Overall Height	810 mm	Max Seat Height	640 mm	Max User Weight	250 kg	Min Overall Height	670 mm		
Min	ı Seat Height	480 mm	Seat Depth	300 mm	Seat Width	450 mm				

12325A52 R & R Healthcare Equipment Premium Padded Shower Stool 520mm

Max Overall Height	810 mm	Max Seat Height	640 mm	Max User Weight	250 kg	Min Overall Height	670 mm
Min Seat Height	480 mm	Seat Depth	300 mm	Seat Width	520 mm		

12325A60 R & R Healthcare Equipment Premium Padded Shower Stool 600mm

Max Overall Height	810 mm	Max Seat Height	640 mm	Max User Weight	250 kg	Min Overall Height	670 mm
Min Seat Height	480 mm	Seat Depth	300 mm	Seat Width	600 mm		



Don't forget to check out the entire range of CCG Shower Furniture products at countrycaregroup.com.au

Powder coated white



FOR FULL

PRODUCT DETAILS

PAGE 85

ADULT TRANSIT MOBILE SHOWER COMMODE

- Australian made shower commode. Available with standard features or
- custom made to suit the individual. Corrosion resistant, heavy duty stainless
- steel frame. 3 Flexible PVC backrest straps. Height adjustable swing away footrests,
- slide out footrest or weightbearing platform.
- Swing up arms, swing out arms (locking), loop arms or pan rails.
- 2 swivel castors at rear and 2 total locking castors at front.

Country Care Group | countrycaregroup.com.au | FREE CALL 1800 727 382



4.01 Shower Chair Static Adjustable

· Powder coated white small drain holes in seat (no catching areas)

4.04 Shower Stool Static

· Designed to assist those who have difficulty showering



4.12 Mobile Shower Commode Standard Self Propelled

Aquatec Ocean Mobile Shower Commode



- Stable stainless steel frame • Tool-less height adjustable seat
- Machine washable backrest
- Space saving push bar • Ergonomic swing-away armrests
- Height adjustable swing-away footrests
- Fold-away texture footplates
- Available in transit or self propelled

1564775 Aquatec Ocean Mobile Shower Commode Transit

Armrest Type	Flat	Footrest Type	Swingaway	Max Seat Height	600 mm	Max User Weight	150 kg
Min Seat Height	475 mm	Overall Length	900 mm	Overall Width	560 mm	Product Weight	14.5 kg
Seat Depth	450 mm	Seat Width	480 mm				

PKG-7.12.910 Aquatec Ocean Mobile Shower Commode Self Propelled

Armrest Type	Flat	Footrest Type	Swingaway	Max Overall Height	1070 mm	Max Seat Height	600 mm
Max User Weight	150 kg	Min Overall Height	940 mm	Min Seat Height	475 mm	Overall Length	1060 mm
Overall Width	675 mm	Product Weight	19.5 kg	Seat Depth	450 mm	Seat Width	480 mm

Showerbuddy Roll in Shower Chair



- Uncomplicated chair with as little structure as possible
- Allows easier care and cleaning process for the carer and user Consideration for comfort and form
- Comes standard with a new design commode bucket
- Soft foam hand grips extend around and down the push handle
- Autolocking, easy release foldaway armrests
- Fully height adjustable
- One piece fold away tray for the footrest
- Available with transit or self propelled

SB6-W Showerbuddy Roll in Self Propelled Shower Chair

Armrest Type	Flat	Backrest Height	350 mm	Max Seat Height	595 mm	Max User Weight	136 kg
Min Seat Height	570 mm	Overall Length	1060 mm	Overall Width	705 mm	Seat Depth	430 mm
Seat Width	460 mm						

SB6-C Showerbuddy Roll in Transit Shower Chair


Armrest Type	Flat	Backrest Height	350 mm	Footrest Type	Foldaway	Max Seat Height	575 mm
Max User Weight	136 kg	Min Seat Height	500 mm	Overall Length	900 mm	Overall Width	555 mm
Seat Depth	430 mm	Seat Width	460 mm				

Shower Furniture

Tubalco Transit Mobile Shower Commode



A



Australian made

- 3 flexible PVC backrest straps
- 2 swivel castors at rear and 2 total locking castors at front

00-920A Tubalco Adult Transit Mobile Shower Commode Swingaway Footrest with Loop Arm 460mm

				5 5		•	
Armrest Type	Loop	Floor To Seat Height	605 mm	Footrest Type	Swingaway	Max User Weight	150 kg
Overall Height	1020 mm	Overall Length	880 mm	Overall Width	550 mm	Product Weight	11 kg
Seat Depth	410 mm	Seat Width	460 mm				

00-930A Tubalco Adult Transit Mobile Shower Commode Slideout Footrest and Loop Arm 460mm

Armrest Type	Loop	Floor To Seat Height	605 mm	Footrest Type	Slideout	Max User Weight	150 kg
Overall Height	1020 mm	Overall Length	880 mm	Overall Width	550 mm	Product Weight	11 kg
Seat Depth	410 mm	Seat Width	460 mm				

00-930B Tubalco Adult Transit Mobile Shower Commode Slideout Weightbearing Platform and Flat Arm 460mm

Armrest Type	Flat Lock- ing	Floor To Seat Height	605 mm	Footrest Type	Weight- bearing	Max User Weight	200 kg
Overall Height	1020 mm	Overall Length	805 mm	Overall Width	545 mm	Product Weight	11.5 kg
Seat Depth	410 mm	Seat Width	460 mm				

00-930C Tubalco Adult Transit Mobile Shower Commode Slideout Footrest and Flat Arm 460mm

Armrest Type	Flat Lock- ing	Floor To Seat Height	605 mm	Footrest Type	Slideout	Max User Weight	125 kg
Overall Height	1020 mm	Overall Length	930 mm	Overall Width	550 mm	Product Weight	11 kg
Seat Depth	410 mm	Seat Width	460 mm				



4.13 Mobile Shower Commode Standard **Transit Fixed Height**

Height adjustable swing away, slide out or weightbearing footrest



4.14 Mobile Shower Commode Standard Self **Propelled Fixed Height**

Tubalco Self Propelled Mobile Shower Commode



- Australian made • Custom options available
- Corrosion resistant componenets
- Swing away footrest
- Loop arms
- 3 flexible PVC backrest straps • 22" Self propelling wheels
- 00-950 Tubalco Adult Self Propelled Mobile Shower Commode 460mm

Armrest Type	Loop	Floor To Seat Height	605 mm	Footrest Type	Swingaway	Max User Weight	125 kg
Overall Height	1020 mm	Overall Length	1000 mm	Overall Width	650 mm	Product Weight	14.5 kg
Seat Depth	410 mm	Seat Width	460 mm				

00-950A Tubalco Adult Self Propelled Mobile Shower Commode Powdercoated 460mm

Armrest Type	Loop	Floor To Seat Height	605 mm	Footrest Type	Swingaway	Max User Weight	125 kg
Overall Height	1020 mm	Overall Length	1000 mm	Overall Width	650 mm	Product Weight	14.5 kg
Seat Depth	410 mm	Seat Width	460 mm				

Tubalco Paediatric Self Propelled Shower Commode

1	5
Y DY	VI
3	j Of

- Durable stainless steel frame
- Very stable base with large footprint Special foot rest hangers for shorter legs available
- 3 flexible PVC backrest straps
- 22" self propelling wheels
- Flat teen or flat locking armrests
- Swing away footrest

00-966 Tubalco Paediatric Self Propelled Mobile Shower Commode 320mm

Armrest Type	Flat Lock- ing	Floor To Seat Height	600 mm	Footrest Type	Swingaway	Max User Weight	100 kg
Overall Height	1010 mm	Overall Length	940 mm	Overall Width	625 mm	Paediatric ?	у
Product Weight	13.5 kg	Seat Depth	320 mm	Seat Width	320 mm		


00-966P Tubalco Child Self Propelled Shower Commode 360mm

Armrest Type	Flat Lock- ing	Floor To Seat Height	600 mm	Footrest Type	Swingaway	Max User Weight	100 kg
Overall Height	1010 mm	Overall Length	940 mm	Overall Width	655 mm	Paediatric ?	у
Product Weight	13.5 kg	Seat Depth	320 mm	Seat Width	360 mm		

Shower Furniture

Tubalco Transit Heavy Duty Mobile Shower Commode





Australian made

• Custom options available

• 3 flexible PVC backrest straps • 2 swivel castors at rear and 2 total locking castors at front

00-940-50 Tubalco Heavy Duty Transit Mobile Shower Commode Swingaway Footrest with Loop Arm 500mm

Armrest Type	Loop	Floor To Seat Height	615 mm	Footrest Type	Swingaway	Max User Weight	200 kg
Overall Height	1020 mm	Overall Length	880 mm	Overall Width	560 mm	Product Weight	11 kg
Seat Depth	410 mm	Seat Width	500 mm				

Tubalco Bariatric Transit Mobile Shower Commode



Australian made

- Custom options available
- Corrosion resistant components
- 4 total locking castors
- Loop locking padded arms
- Swingaway or slide out footrest

00-975-50 Tubalco Bariatric Transit Mobile Shower Commode Slideout Footrest 500mm

Armrest Type	Loop Locking Padded	Floor To Seat Height	620 mm	Footrest Type	Slideout	Max User Weight	180 kg
Overall Height	1025 mm	Overall Length	895 mm	Overall Width	750 mm	Product Weight	15 kg
Seat Width	500 mm						

00-970-50 Tubalco Bariatric Transit Mobile Shower Commode Fixed Height 590mm

Armrest Type	Loop Locking Padded	Footrest Type	Swingaway	Max User Weight	180 kg	Overall Height	1025 mm
Overall Length	855 mm	Overall Width	700 mm	Product Weight	15 kg	Seat Width	500 mm



4.15 Mobile Shower Commode Bariatric **Transit and Self Propelled**

4.17 Mobile Shower Commode Standard **Transit Tilting**

Rifton Hygiene & Toileting System

- The Rifton Hygiene and Toileting System (HTS) is designed to improve the health of a child and the life of a caregiver
- Good toileting is a necessity for childrens health
- Hygienic, simple, affordable the Rifton HTS promotes the forward positioning that is ideal and natural for effective toileting . Its versitility and simplicity for use on, over and off the toilet with no tools required



Z110 Rifton Hygien	e & Toilet	ting System Small					
Backrest Height	420 mm	Max Seat Height	508 mm	Max User Weight	34 kg	Min Seat Height	393 mm
Overall Width	540 mm	Paediatric ?	Y	Product Weight	6.8 kg	Max Seat Depth	305 mm
Min Seat Depth	203 mm	Max Seat Width	305 mm	Min Seat Width	203 mm		

Z120 Rifton Hygiene & Toileting System Medium

Backrest Height	495 mm	Max Seat Height	508 mm	Max User Weight	68 kg	Min Seat Height	393 mm
Overall Width	540 mm	Paediatric ?	Y	Product Weight	12.5 kg	Max Seat Depth	368 mm
Min Seat Depth	292 mm	Max Seat Width	355 mm	Min Seat Width	254 mm		

Showerbuddy Roll-In Buddy Tilt



- Revolutionary tilt option allows for pressure release • Gives a more thorough bathing experience for the patient
- Provides easier reach for the caregiver
- Incorporates safety features and corrosion resistant materials as standard
- Adjustable tilt angle from 0 30 degrees
- Rolls over standard and elongated toilets
- Low maintenance
- Removable, lockable & foldaway arms Height adjustable
- Fold up footrest on rolling and bath base

SB3T Showerbuddy Roll-In Buddy Tilt

		· · ·					
Armrest Type	Flat	Backrest Height	440 mm	Footrest Type	Foldaway	Max Seat Height	610 mm
Max User Weight	135 kg	Min Seat Height	510 mm	Overall Length	1100 mm	Overall Width	490 mm
Seat Depth	420 mm	Seat Width	450 mm				

Shower Furniture

Tubalco Adult Tilting Transit Mobile Shower Commode



- Australian made Headrest
- Loop locking armrests with side guards Swingaway footrest
- Adjustable footrest hangers
 Remote tilt control allowing face to face operation
 - Smooth tilting action reduces user anxiety
 - Corrosion resistant components

00-960 Tubalco Adult Tilting Transit Mobile Shower Commode 460mm

Armrest Type	Loop Locking	Floor To Seat Height	610 mm	Footrest Type	Swingaway	Max User Weight	160 kg
Overall Height	1510 mm	Overall Length	1120 mm	Overall Width	575 mm	Product Weight	19.5 kg
Seat Depth	410 mm	Seat Width	460 mm				

Tubalco Paediatric-Teen Tilting Transit Mobile Shower Commode



- Height adjustable swing away footrestsDurable stainless steel frame Very stable base Large footprint Special foot rest hangers for shorter legs available
 - 35 Degree tilt
- Remote attendant tilt release trigger Swing up locking flat arms

00-961P Tubalco Child Transit Tilting Mobile Shower Commode 360mm

Armrest Type	Flat Lock- ing	Floor To Seat Height	600 mm	Footrest Type	Swingaway	Max User Weight	125 kg
Overall Height	1370 mm	Overall Length	1045 mm	Overall Width	510 mm	Product Weight	15.5 kg
Seat Depth	320 mm	Seat Width	360 mm				

00-961 Tubalco Teen Transit Tilting Mobile Shower Commode 400mm

Armrest Type	Flat Lock- ing	Floor To Seat Height	600 mm	Footrest Type	Swingaway	Max User Weight	100 kg
Overall Height	1370 mm	Overall Length	1045 mm	Overall Width	510 mm	Product Weight	15.5 kg
Seat Depth	360 mm	Seat Width	400 mm				



4.18 Mobile Shower Commode Standard **Transit Tilting Fixed**



4.2 Mobile Shower Commode Accessories

Aquatec Accessories for Orca and Beluga Bath Lifters

Accessories and replacement parts for Aquatec Orca bath lifts

1528702 Aquatec Safety Belt for	Aquatec Orca Bath Lift	:	
Package Quantity	1	Package Type	Single

Aquatec Accessories/Replacement Parts for Mobile Shower Commodes

• Accessories and replacement parts for Aquatec mobile shower commodes

Aquatec Sanitary Pan and Lid Aquatec Ocean, VIP, Dual VIP & E-VIP
Aquatec Sanitary Pan Holder for Aquatec Ocean Shower Commode
Aquatec Soft Seat Key Cutout for Aquatec Ocean and E-VIP Shower Commodes
Aquatec Calf Support Height and Angle Adjustable Left or Right Ocean Shower Com- modes
Aquatec Sanitary Pan Holder for Aquatec Ocean VIP, XL VIP & Dual VIP
Aquatec Armrest Lock Pair for Aquatec Ocean VIP, VIP XL, Dual VIP & E-VIP
Aquatec Safety Bar XL for Aquatec Ocean XL and VIP Shower Commodes
Aquatec Chest Belt for Aquatec Ocean Shower Commodes
Aquatec Pelvis Belt for Aquatec Ocean Shower Commodes
Aquatec Amputation Support Above Knee for Aquatec Ocean Shower Commodes Left
Aquatec Amputation Support Below Knee for Aquatec Ocean Shower Commodes Left
Aquatec Armrest for Aquatec Ocean and E-VIP XL Shower Commodes Size XL
Aquatec Amputation Support Above Knee for Aquatec Ocean Shower Commodes Right
Aquatec Amputation Support Below Knee for Aquatec Ocean Shower Commodes Right
Aquatec Anti-Tippers for Aquatec Ocean Shower Commodes
Aquatec Universal Soft Seat for Aquatec Ocean Shower Commodes
Aquatec Safety Belt for Aquatec Orca Bath Lift
Aquatec Ergonomic Hygiene Recess Seat for Aquatec Ocean Shower Commodes
Aquatec Safety Bar for Aquatec Ocean Shower Commodes
Aquatec Leg Extensions (Tall)
Aquatec XL Facelift Kit for Ocean Commodes
Aquatec Plastic Back Rest for Ocean Commodes
Aquatec Amputation Support One Piece Aquatec Ocean Shower Commodes
Quickie Anti-tip Bars for Quickie 2 Lite Pair
Quickie Under Seat Pouch for Quickie 2 Lite
Quickie Tension Adjustable Backrest Quickie 2 Lite
Quickie Push Handle 8"
Quickie Extension Handle 6" Pair


Tubalco Mobile Shower Commode Open Seats



00-908B Tubalco O	pen Sea	t Blue 460mm		
Overall Depth	410 mm	Overall Width	460 mm	P
00-908W Tubalco ()pen Sea	at White 460mm		
Overall Depth	410 mm	Overall Width	460 mm	P
00-902 Tubalco Op	en Seat	White 500mm		
Overall Depth	400 mm	Overall Width	500 mm	P
00-990 Tubalco Op	en Seat	White 550mm		
Overall Depth	430 mm	Overall Width	550 mm	P
00-997 Tubalco Op	en Seat	White 600mm		
Overall Depth	460 mm	Overall Width	600 mm	P
00-996 Tubalco Op	en Seat	White 650mm		
Overall Depth	460 mm	Overall Width	650 mm	P

Tubalco Mobile Shower Commode Seats





- With stainless steel staples and glued joints Various sizes available
- Also available as an open front seat

00-909B Tubalco C	losed Se	eat Blue 460mm					
Overall Depth	460 mm	Overall Width	460 mm	Package Quantity	1	Package Type	Single
00-909W Tubalco (Closed S	eat White 460mm					
Overall Depth	410 mm	Overall Width	460 mm	Package Quantity	1	Package Type	Single
00-903 Tubalco Clo Overall Depth	osed Sea	at White 500mm Overall Width	500 mm	Package Quantity	1	Package Type	Single
00-991 Tubalco Clo	sed Sea	t White 550mm		1		1	
Overall Depth	430 mm	Overall Width	550 mm	Package Quantity	1	Package Type	Single
00-995 Tubalco Clo	osed Sea	at White 600mm					
Overall Depth	460 mm	Overall Width	600 mm	Package Quantity	1	Package Type	Single



Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single



Overall Depth	460 mm	Overall Width	650 mm	Package Quantity	1	Package Type	Single
ubalaa Shau		de Pressure Relie	ving Or	on Soato			
		ue riessuie keile	wing of	Jell Seals			
		For Tubalco commodes					
		White or blue leather loo	,				
		 5cm thick foam padding With stainless steel stap 		d ininte			
		Various sizes available	lies allu giue	u joints			
		 Also available as a close 	d front seat				
		Policying Cost Port		latery Onen Front FO	0		
		_		lstery Open Front 50	Omm		
	co Pressure	Relieving Seat Dart	ex Upho	Istery Open Front 50 Package Quantity	1 1	Package Type	Single
Overall Depth	400 mm	Overall Width	500 mm	Package Quantity	1	Package Type	Single
Overall Depth	400 mm	Overall Width	500 mm		1	Package Type	Single
Overall Depth 10-990P Tubal	400 mm	Overall Width	500 mm	Package Quantity	1	Package Type Package Type	Single Single
Overall Depth 1 0-990P Tubal	400 mm	Overall Width Relieving Seat Dart	500 mm	Package Quantity	1 Omm		
Overall Depth IO-990P Tubal Overall Depth	400 mm co Pressure 430 mm	Overall Width Relieving Seat Dart Overall Width	500 mm ex Upho 550 mm	Package Quantity	1 0mm 1		
Overall Depth O-990P Tubal Overall Depth O-994P Tubal	400 mm co Pressure 430 mm	Overall Width Relieving Seat Dart Overall Width	500 mm ex Upho 550 mm	Package Quantity Istery Open Front 55 Package Quantity	1 0mm 1		
Overall Depth O-990P Tubal Overall Depth O-994P Tubal	400 mm co Pressure 430 mm co Pressure	Overall Width Relieving Seat Dart Overall Width Relieving Seat Dart	500 mm ex Upho 550 mm ex Upho	Package Quantity Istery Open Front 55 Package Quantity Istery Open Front 60	1 0mm 1	Package Type	Single
Overall Depth D0-990P Tubal Overall Depth D0-994P Tubal Overall Depth	400 mm co Pressure 430 mm co Pressure 460 mm	Overall Width Relieving Seat Dart Overall Width Relieving Seat Dart Overall Width	500 mm ex Upho 550 mm ex Upho 600 mm	Package Quantity Istery Open Front 55 Package Quantity Istery Open Front 60 Package Quantity	1 0mm 1 00mm 1	Package Type	Single
Overall Depth DO-990P Tubal Overall Depth DO-994P Tubal Overall Depth	400 mm co Pressure 430 mm co Pressure 460 mm	Overall Width Relieving Seat Dart Overall Width Relieving Seat Dart Overall Width	500 mm ex Upho 550 mm ex Upho 600 mm	Package Quantity Istery Open Front 55 Package Quantity Istery Open Front 60	1 0mm 1 00mm 1	Package Type	Single

Tubalco Shower Commode Pressure Relieving Closed Seats

- For Tubalco commodes White or blue leather look vinyl 5cm thick foam padding
 - With stainless steel staples and glued joints
 - Various sizes available Also available in an open front seat

00-903P Tubalco Pressure Relieving Seat Dartex Upholstery Closed Front 500mm										
Overall Depth	400 mm	Overall Width	500 mm	Package Quantity	1	Package Type	Single			
00-991P Tubalco Pressure Relieving Seat Dartex Upholstery Closed Front 550mm										
Overall Depth	430 mm	Overall Width	550 mm	Package Quantity	1	Package Type	Single			
•										
00-997P Tubalco F	ressure	Relieving Seat Darte	ex Uphol	stery Closed Front 6	50mm					
00-997P Tubalco F	460 mm	Relieving Seat Darte	ex Uphol	stery Closed Front 6 Package Quantity	50mm	Package Type	Single			
Overall Depth	460 mm	Overall Width	650 mm	-	1	Package Type	Single			

Tubalco Shower Commode Accessories/Options

• Replacement parts and accessories for Tubalco Shower Commodes • Help to customise the shower commode to suit the users needs

722	Tubalco Height Adjustable Swing Away Foot Rest
793	Tubalco Swivel Castor 125mm x 35mm
794	Tubalco Total Locking Castor 125mm x 35mm
26	Tubalco Safety Flip Bars Pair
782	Tubalco Deluxe Pan Carrier
785	Tubalco Standard Bed Pan
00-1000	Tubalco Articulated Head Rest Complete with Pad
00-994	Tubalco Teen Pan
00-998	Tubalco Elevated Leg Rest
737B	Tubalco Padded Armrests Blue On Chair Welded Se
737W	Tubalco Padded Armrests White On Chair Welded
737PB	Tubalco Padded Back 600mm Wide

Showerbuddy Accessories/Replacement Parts for Shower Buddy

Accessories and replacement parts for Shower Buddy shower commodes

117	Showerbuddy Safety Restraints for Shower Buddy
NRA	Showerbuddy Neckrest Assembly for Shower Buddy
PBC	Showerbuddy Padded Foam Backrest for Shower Buddy
SCCF	Showerbuddy Seat Cushion Closed Front for Shower Buddy
SCLO	Showerbuddy Seat Cushion Left Opening for Shower Buddy
SCOC	Showerbuddy Seat Cushion Overlay for Shower Buddy
SCRO	Showerbuddy Seat Cushion Right Opening for Shower Buddy
WW/22	Showerbuddy Wheelchair Wheel Kit for Shower Buddy
737PB	Tubalco Padded Back 600mm Wide



Seam Pads d Seam Pads



Body Point Aeromesh Rapid Dry Bath Belt Two Piece



- The Aeromesh Rapid Dry Bath Belt has the great design and features you're used to • The original two piece bath belt has grommets to enable the belt to be mounted directly to the bath chair using a screw or
- the bath belt release knob • Padding at the middle of the belt allows for greater user comfort Helps keep the user in the chair during bathing

BB101M-1 Body Point Aeromesh Rapid Dry Bath Belt Two Piece Medium

BB101L-1 Body Point Aeromesh Rapid Dry Bath Belt Two Piece Large

Body Point Aeromesh Rapid Dry Bath Belt One Piece



Showering Furniture

- The Aeromesh Rapid Dry Bath Belt has the great design and features you're used to
- . The new one piece bath belt comes with either grommets to enable the belt to be mounted directly to the bath chair using a screw or the bath belt release knob
- It is also available with plastic slides so the belt can be mounted directly to the bath chair frame by wrapping the webbing end around the tubing and securing with the slide
- Padding at the middle of the belt allows for greater user comfort
- Helps keep the user in the chair during bathing

Package Quantity 1 Package Type Single Paediatric ? Y										
BB112L-1 Body Point Aeromesh Rapid Dry Bath Belt One Piece Large with Plastic Slides										
BBIIZE-I BOUY POINT AEI	omesn kapiu	Diy Batil Belt Olle Piece	Large with P	lastic Silues						
Package Quantity 1 Package Type Single Paediatric ? Y										
r acrage quantity	·	i dokage iype	Single							
	romesh Rapid	d Dry Bath Belt One Piece	5							
	romesh Rapio		5	Paediatric ?	Ŷ					
BB111XL-1 Body Point Ae	romesh Rapio	d Dry Bath Belt One Piece	e XLarge		Y					
BB111XL-1 Body Point Ae Package Quantity	1	d Dry Bath Belt One Piece	2 XLarge Single	Paediatric ?	Ŷ					





TOILETING EQUIPMENT

Don't forget to check out the entire range of CCG Toileting Equipment products at countrycaregroup.com.au

National Distribution, Local Service

Toileting Equipment

4.21 Bedside Commode Standard

R & R Healthcare Equipment Premium Deluxe Bedside Commode



- Hinged seat Lightweight yet strong Padded vinyl contoured backrest Toilet seat and backrest add comfort when toileting Height adjustable legs
- Front and side skirt to hide the bucket underneath Comes in various widths and two colours

12087A45 R & R Healthcare Equipment Premium Deluxe Bedside Commode 450mm

N	lax Seat Height	620 mm	Max User Weight	250 kg	Min Seat Height	470 mm	Seat Depth	460 mm
S	eat Width	450 mm						

Toileting Equipment

Toileting Equipment

4.22 Bedside Commode Heavy Duty

R & R Healthcare Equipment Premium Deluxe Bedside Commode



- Hinged seat Lightweight yet strong Padded vinyl contoured backrest
- Toilet seat and backrest add comfort when toileting
- Height adjustable legs • Front and side skirt to hide the bucket underneath
- Comes in various widths and two colours

12087A52 R & R Healthcare Equipment Premium Deluxe Bedside Commode Wide Fawn 520mm

Max Seat Height	620 mm	Max User Weight	250 kg	Min Seat Height	470 mm	Seat Depth	460 mm
Seat Width	520 mm						

12087A60 R & R Healthcare Equipment Premium Deluxe Bedside Commode Extra Wide Fawn 600mm

Max Seat Height	620 mm	Max User Weight	250 kg	Min Seat Height	470 mm	Seat Depth	460 mm
Seat Width	600 mm						

PATENT



EVOLUTION SUPREN

- Made from visco elastic foam in a carflex cover which provides multiple stretch
- The user's weight is distributed over a large area relieving pressure
- Doesn't compromise postural support and allows natural alignment of the spine
- increased lumber support, drop down positions
- Delivers a comfortable long term sitting



NINGHAM LIFT CHAIR

- Simple operation and space saver .
- Lift assist function .
- .
- Can be custom made
- . Max user weight 130kg
- Emergency lower battery back up





environment and helps protect against pressure sores

- This Chair also allows for the semi fowler position which increases blood circulation and helps the user to breathe more easily
- With fully adjustable head support, This chair has an extra high leg position to reduce fluid retention
- sides and a wide variety of seating This chair is designed to be low maintenance as well as user friendly for both the user and the carer

Triple roll back comfort and can be modified to cater for back conditions Choice of fabrics, sizes and optional extras available

Don't forget to check out the entire range of CCG Patient Chair products at countrycaregroup.com.au

Patient Chairs

5.06 Accessories for Power Lift Chairs

Oscar Furniture Over Chair Activity Table



• Suitable for eating or other activities that are made easier with a table Slides up to recliner or static chair

Wooden Over Chair Activity Table

- Convenient and easy to use
- Cleverly designed to be used with your Oscar electric lift recliner
- Four castors for easy manoeuvrability

31054 Oscar Furniture Over Chair Activity Table										
Overall Height 680 mm Overall Length 510 mm Overall Width 860 mm										

Patient Chairs

5.06 Power Lift Chairs

Oscar Furniture Menningham Lift Chair Single Motor



- Simple operation and space saver Lift assist function
- Triple roll back comfort and can be modified to cater for back conditions
- · Choice of fabrics, sizes and optional extras available
- Can be custom made
- Max user weight 130kg
- Emergency lower battery back up

CWSEL-ME-M-SO Oscar Furniture Menningham Lift Chair Single Motor Mini Olive Suede

Armrest Height	190 mm	Extended Footrest	840 mm	Max Backrest Height	725 mm	Max Seat Height	485 mm
Max User Weight	130 kg	Min Backrest Height	575 mm	Min Seat Height	435 mm	Max Seat Depth	510 mm
Min Seat Depth	410 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-ME-M-SM Oscar Furniture Menningham Lift Chair Single Motor Mini Moss Suede

Armrest Height	190 mm	Extended Footrest	840 mm	Max Backrest Height	725 mm	Max Seat Height	485 mm
Max User Weight	130 kg	Min Backrest Height	575 mm	Min Seat Height	435 mm	Max Seat Depth	510 mm
Min Seat Depth	410 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-ME-M-SN Oscar Furniture Menningham Lift Chair Single Motor Mini Nutmeg Suede

Armrest Height	190 mm	Extended Footrest	840 mm	Max Backrest Height	725 mm	Max Seat Height	485 mm
Max User Weight	130 kg	Min Backrest Height	575 mm	Min Seat Height	435 mm	Max Seat Depth	510 mm
Min Seat Depth	410 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-ME-A-SO Oscar Furniture Menningham Lift Chair Single Motor A Olive Suede

Armrest Height	190 mm	Extended Footrest	960 mm	Max Backrest Height	750 mm	Max Seat Height	515 mm
Max User Weight	130 kg	Min Backrest Height	600 mm	Min Seat Height	465 mm	Max Seat Depth	550 mm
Min Seat Depth	450 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-ME-A-SM)scar Furr	niture Menningham I	lift Chair	Single Motor A Mos	s Suede		
Armrest Height	190 mm	Extended Footrest	960 mm	Max Backrest Height	750 mm	Max Seat Height	515 mm
Max User Weight	130 kg	Min Backrest Height	600 mm	Min Seat Height	465 mm	Max Seat Depth	550 mm
Min Seat Depth	450 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

Armrest Height	190 mm	Extended Footrest	960 mm	Max Backrest Height	750 mm	Max Seat Height	515 mm
Max User Weight	130 kg	Min Backrest Height	600 mm	Min Seat Height	465 mm	Max Seat Depth	550 mm
Min Seat Depth	450 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEK-ME-A-SN Oscar Furniture Menningham Lift Chair Single Motor A Nutmeg Suede

Armres	st Height	190 mm	Extended Footrest	960 mm	Max Backrest Height	750 mm	Max Seat Height	515 mm
Max Us	ser Weight	130 kg	Min Backrest Height	600 mm	Min Seat Height	465 mm	Max Seat Depth	550 mm
Min Sea	at Depth	450 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-ME-B-SO Oscar Furniture Menningham Lift Chair Single Motor B Olive Suede

Armrest Height	170 mm	Extended Footrest	1015 mm	Max Backrest Height	775 mm	Max Seat Height	535 mm
Max User Weight	130 kg	Min Backrest Height	625 mm	Min Seat Height	485 mm	Max Seat Depth	575 mm
Min Seat Depth	500 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-ME-B-SN Oscar Furniture Menningham Lift Chair Single Motor B Moss Suede

Armrest Height	170 mm	Extended Footrest	1015 mm	Max Backrest Height	775 mm	Max Seat Height	535 mm
Max User Weight	130 kg	Min Backrest Height	625 mm	Min Seat Height	485 mm	Max Seat Depth	575 mm
Min Seat Depth	500 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-ME-B-SM Oscar Furniture Menningham Lift Chair Single Motor B Nutmeg Suede

Armrest Height	170 mm	Extended Footrest	1015 mm	Max Backrest Height	775 mm	Max Seat Height	535 mm
Max User Weight	130 kg	Min Backrest Height	625 mm	Min Seat Height	485 mm	Max Seat Depth	575 mm
Min Seat Depth	500 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-ME-C-SO Oscar Furniture Menningham Lift Chair Single Motor C Olive Suede

Armrest Height	170 mm	Extended Footrest	1110 mm	Max Backrest Height	840 mm	Max Seat Height	575 mm
Max User Weight	130 kg	Min Backrest Height	690 mm	Min Seat Height	525 mm	Max Seat Depth	575 mm
Min Seat Depth	525 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-ME-C-SM Oscar Furniture Menningham Lift Chair Single Motor C Moss Suede

[Armrest Height	170 mm	Extended Footrest	1110 mm	Max Backrest Height	840 mm	Max Seat Height	575 mm
	Max User Weight	130 kg	Min Backrest Height	690 mm	Min Seat Height	525 mm	Max Seat Depth	575 mm
	Min Seat Depth	525 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-ME-C-SN Oscar Furniture Menningham Lift Chair Single Motor C Nutmeg Suede

Armrest Height	170 mm	Extended Footrest	1110 mm	Max Backrest Height	840 mm	Max Seat Height	575 mm
Max User Weight	130 kg	Min Backrest Height	690 mm	Min Seat Height	525 mm	Max Seat Depth	575 mm
Min Seat Depth	525 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		





Oscar Furniture Menningham Lift Chair Dual Motor



- Simple operation allowing the chair to slowly rise to the required position
- Triple roll back comfort, can be modified to cater for back conditions • The dual motor allows independent footrest/recline function
- Choice of fabric and options including custom sizing
- Ultra fresh antimicrobial dunlop foam
- Max user weight 130kg
- Controls come on right but can be moved to left

CDALR-ME-M-SO Oscar Furniture Menningham Lift Chair Dual Motor Mini Olive Suede

Armrest Height	190 mm	Extended Footrest	840 mm	Max Backrest Height	725 mm	Max Seat Height	485 mm
Max User Weight	130 kg	Min Backrest Height	575 mm	Min Seat Height	435 mm	Max Seat Depth	510 mm
Min Seat Depth	410 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-ME-M-SM Oscar Furniture Menningham Lift Chair Dual Motor Mini Moss Suede

Armrest Height	190 mm	Extended Footrest	840 mm	Max Backrest Height	725 mm	Max Seat Height	485 mm
Max User Weight	130 kg	Min Backrest Height	575 mm	Min Seat Height	435 mm	Max Seat Depth	510 mm
Min Seat Depth	410 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-ME-M-SN Oscar Furniture Menningham Lift Chair Dual Motor Mini Nutmeg Suede

Armrest Height	190 mm	Extended Footrest	840 mm	Max Backrest Height	725 mm	Max Seat Height	485 mm
Max User Weight	130 kg	Min Backrest Height	575 mm	Min Seat Height	435 mm	Max Seat Depth	510 mm
Min Seat Depth	410 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-ME-A-SO Oscar Furniture Menningham Lift Chair Dual Motor A Olive Suede

Armrest Height	190 mm	Extended Footrest	960 mm	Max Backrest Height	750 mm	Max Seat Height	515 mm
Max User Weight	130 kg	Min Backrest Height	600 mm	Min Seat Height	465 mm	Max Seat Depth	550 mm
Min Seat Depth	450 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-ME-A-SM Oscar Furniture Menningham Lift Chair Dual Motor A Moss Suede

Armrest Height	190 mm	Extended Footrest	960 mm	Max Backrest Height	750 mm	Max Seat Height	515 mm
Max User Weight	130 kg	Min Backrest Height	600 mm	Min Seat Height	465 mm	Max Seat Depth	550 mm
Min Seat Depth	450 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-ME-A-SN Oscar Furniture Menningham Lift Chair Dual Motor A Nutmeg Suede

Armrest Height	190 mm	Extended Footrest	960 mm	Max Backrest Height	750 mm	Max Seat Height	515 mm
Max User Weight	130 kg	Min Backrest Height	600 mm	Min Seat Height	465 mm	Max Seat Depth	550 mm
Min Seat Depth	450 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-ME-B-SO Oscar Furniture Menningham Lift Chair Dual Motor B Olive Suede

Armrest Height	170 mm	Extended Footrest	1015 mm	Max Backrest Height	775 mm	Max Seat Height	535 mm
Max User Weight	130 kg	Min Backrest Height	625 mm	Min Seat Height	485 mm	Max Seat Depth	575 mm
Min Seat Depth	500 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-ME-B-SM Oscar Furniture Menningham Lift Chair Dual Motor B Moss Suede

Armrest Height	170 mm	Extended Footrest	1015 mm	Max Backrest Height	775 mm	Max Seat Height	535 mm
Max User Weight	130 kg	Min Backrest Height	625 mm	Min Seat Height	485 mm	Max Seat Depth	575 mm
Min Seat Depth	500 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-ME-B-SN Oscar Furniture Menningham Lift Chair Dual Motor B Nutmeg Suede

Armrest Height	170 mm	Extended Footrest	1015 mm	Max Backrest Height	775 mm	Max Seat Height	535 mm
Max User Weight	130 kg	Min Backrest Height	625 mm	Min Seat Height	485 mm	Max Seat Depth	575 mm
Min Seat Depth	500 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-ME-C-SO Oscar Furniture Menningham Lift Chair Dual Motor C Olive Suede

Armrest Height	170 mm	Extended Footrest	1110 mm	Max Backrest Height	840 mm	Max Seat Height	575 mm
Max User Weight	130 kg	Min Backrest Height	690 mm	Min Seat Height	525 mm	Max Seat Depth	575 mm
Min Seat Depth	525 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-ME-C-SM Oscar Furniture Menningham Lift Chair Dual Motor C Moss Suede

Armrest Height	170 mm	Extended Footrest	1110 mm	Max Backrest Height	840 mm	Max Seat Height	575 mm
Max User Weight	130 kg	Min Backrest Height	690 mm	Min Seat Height	525 mm	Max Seat Depth	575 mm
Min Seat Depth	525 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-ME-C-SN Oscar Furniture Menningham Lift Chair Dual Motor C Nutmeg Suede

Armre	rest Height	170 mm	Extended Footrest	1110 mm	Max Backrest Height	840 mm	Max Seat Height	575 mm
Max L	User Weight	130 kg	Min Backrest Height	690 mm	Min Seat Height	525 mm	Max Seat Depth	575 mm
Min S	Seat Depth	525 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

Oscar Furniture Brumby Lift Chair Single Motor



- Wall saver action, one motor, saves space Dual roll backrest comfort
- Lift assist for help getting in and out of chair
- A range of fabrics, sizes and options available including custom sizing

- Emergency lower battery back up

CWSEL-BR-C-SM Oscar Furniture Brumby Lift Chair Single Motor C Moss Suede

Armrest Height	170 mm	Extended Footrest	1050 mm	Max Backrest Height	820 mm	Max Seat Height	575 mm
Max User Weight	130 kg	Min Backrest Height	670 mm	Min Seat Height	525 mm	Max Seat Depth	540 mm
Min Seat Depth	465 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-BR-C-SN Oscar Furniture Brumby Lift Chair Single Motor C Nutmeg Suede

				· · · · · · · · · · · · · · · · · · ·			
Armrest Height	170 mm	Extended Footrest	1050 mm	Max Backrest Height	820 mm	Max Seat Height	575 mm
Max User Weight	130 kg	Min Backrest Height	670 mm	Min Seat Height	525 mm	Max Seat Depth	540 mm
Min Seat Depth	465 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-BR-C-SO Oscar Furniture Brumby Lift Chair Single Motor C Olive Suede

Armrest Height	170 mm	Extended Footrest	1050 mm	Max Backrest Height	820 mm	Max Seat Height	575 mm
Max User Weight	130 kg	Min Backrest Height	670 mm	Min Seat Height	525 mm	Max Seat Depth	540 mm
Min Seat Depth	465 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-BR-A-SM Oscar Furniture Brumby Lift Chair Single Motor A Moss Suede

Armrest Height	190 mm	Extended Footrest	930 mm	Max Backrest Height	730 mm	Max Seat Height	515 mm
Max User Weight	130 kg	Min Backrest Height	580 mm	Min Seat Height	465 mm	Max Seat Depth	520 mm
Min Seat Depth	420 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		



· Simple operation with controls on the right, can be moved to the left

CWSEL-BR-A-SN Oscar Furniture Brumby Lift Chair Single Motor A Nutmeg Suede

Armrest Height	190 mm	Extended Footrest	930 mm	Max Backrest Height	730 mm	Max Seat Height	515 mm
Max User Weight	130 kg	Min Backrest Height	580 mm	Min Seat Height	465 mm	Max Seat Depth	520 mm
Min Seat Depth	420 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-BR-A-SO Oscar Furniture Brumby Lift Chair Single Motor A Olive Suede

Armrest Height	190 mm	Extended Footrest	930 mm	Max Backrest Height	730 mm	Max Seat Height	515 mm
Max User Weight	130 kg	Min Backrest Height	580 mm	Min Seat Height	465 mm	Max Seat Depth	520 mm
Min Seat Depth	420 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-BR-B-SM Oscar Furniture Brumby Lift Chair Single Motor B Moss Suede

Armrest Height	170 mm	Extended Footrest	980 mm	Max Backrest Height	775 mm	Max Seat Height	535 mm
Max User Weight	130 kg	Min Backrest Height	625 mm	Min Seat Height	485 mm	Max Seat Depth	540 mm
Min Seat Depth	465 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-BR-B-SN Oscar Furniture Brumby Lift Chair Single Motor B Nutmeg Suede

Armrest Height	170 mm	Extended Footrest	980 mm	Max Backrest Height	775 mm	Max Seat Height	535 mm
Max User Weight	130 kg	Min Backrest Height	625 mm	Min Seat Height	485 mm	Max Seat Depth	540 mm
Min Seat Depth	465 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-BR-B-SO Oscar Furniture Brumby Lift Chair Single Motor B Olive Suede

Armrest Height	170 mm	Extended Footrest	980 mm	Max Backrest Height	775 mm	Max Seat Height	535 mm
Max User Weight	130 kg	Min Backrest Height	625 mm	Min Seat Height	485 mm	Max Seat Depth	540 mm
Min Seat Depth	465 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-BR-M-SM Oscar Furniture Brumby Lift Chair Single Motor Mini Moss Suede

Armrest Height	190 mm	Extended Footrest	840 mm	Max Backrest Height	710 mm	Max Seat Height	485 mm
Max User Weight	130 kg	Min Backrest Height	560 mm	Min Seat Height	460 mm	Max Seat Depth	510 mm
Min Seat Depth	410 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-BR-M-SN Oscar Furniture Brumby Lift Chair Single Motor Mini Nutmeg Suede

Armrest Height	190 mm	Extended Footrest	840 mm	Max Backrest Height	710 mm	Max Seat Height	485 mm
Max User Weight	130 kg	Min Backrest Height	560 mm	Min Seat Height	460 mm	Max Seat Depth	510 mm
Min Seat Depth	410 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CWSEL-BR-M-SO Oscar Furniture Brumby Lift Chair Single Motor Mini Olive Suede

Armrest Height	190 mm	Extended Footrest	840 mm	Max Backrest Height	710 mm	Max Seat Height	485 mm
Max User Weight	130 kg	Min Backrest Height	560 mm	Min Seat Height	460 mm	Max Seat Depth	510 mm
Min Seat Depth	410 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

Oscar Furniture Brumby Lift Chair Dual Motor



A

A

• Dual motor chair allows independent footrest/recline function • Dual roll backrest comfort can be modified to cater for back conditions • Simple operation, controls come on right, can be moved to left

• Lift assist for getting in and out of chair Choice of fabrics, options, sizes and custom sizing available

• Ultra fresh antimicrobial dunlop foam with walnut show wood frames Emergency lower battery back up system

CDALR-BR-A-SM Oscar Furniture Brumby Lift Chair Dual Motor A Moss Suede

		-					
Armrest Height	190 mm	Extended Footrest	930 mm	Max Backrest Height	730 mm	Max Seat Height	515 mm
Max User Weight	130 kg	Min Backrest Height	580 mm	Min Seat Height	465 mm	Max Seat Depth	520 mm
Min Seat Depth	420 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-BR-A-SN Oscar Furniture Brumby Lift Chair Dual Motor A Nutmeg Suede

Armrest Height	190 mm	Extended Footrest	930 mm	Max Backrest Height	730 mm	Max Seat Height	515 mm
Max User Weight	130 kg	Min Backrest Height	580 mm	Min Seat Height	465 mm	Max Seat Depth	520 mm
Min Seat Depth	420 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-BR-A-SO Oscar Furniture Brumby Lift Chair Dual Motor A Olive Suede

Armrest Height	190 mm	Extended Footrest	930 mm	Max Backrest Height	730 mm	Max Seat Height	515 mm
Max User Weight	130 kg	Min Backrest Height	580 mm	Min Seat Height	465 mm	Max Seat Depth	520 mm
Min Seat Depth	420 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-BR-B-SM Oscar Furniture Brumby Lift Chair Dual Motor B Moss Suede

Armrest Height	170 mm	Extended Footrest	980 mm	Max Backrest Height	775 mm	Max Seat Height	535 mm
Max User Weight	130 kg	Min Backrest Height	625 mm	Min Seat Height	485 mm	Max Seat Depth	540 mm
Min Seat Depth	465 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-BR-B-SN Oscar Furniture Brumby Lift Chair Dual Motor B Nutmeg Suede

Armrest Height	170 mm	Extended Footrest	980 mm	Max Backrest Height	775 mm	Max Seat Height	535 mm
Max User Weight	130 kg	Min Backrest Height	625 mm	Min Seat Height	485 mm	Max Seat Depth	540 mm
Min Seat Depth	465 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-BR-B-SO Oscar Furniture Brumby Lift Chair Dual Motor B Olive Suede

Armrest Height	170 mm	Extended Footrest	980 mm	Max Backrest Height	775 mm	Max Seat Height	535 mm
Max User Weight	130 kg	Min Backrest Height	625 mm	Min Seat Height	485 mm	Max Seat Depth	540 mm
Min Seat Depth	465 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-BR-C-SM Oscar Furniture Brumby Lift Chair Dual Motor C Moss Suede

Armrest Height	170 mm	Extended Footrest	1050 mm	Max Backrest Height	820 mm	Max Seat Height	575 mm
Max User Weight	130 kg	Min Backrest Height	670 mm	Min Seat Height	525 mm	Max Seat Depth	540 mm
Min Seat Depth	465 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-BR-C-SN Oscar Furniture Brumby Lift Chair Dual Motor C Nutmeg Suede

Armrest Height	170 mm	Extended Footrest	1050 mm	Max Backrest Height	820 mm	Max Seat Height	575 mm
Max User Weight	130 kg	Min Backrest Height	670 mm	Min Seat Height	525 mm	Max Seat Depth	540 mm
Min Seat Depth	465 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		



CDALR-BR-C-SO Oscar Furniture Brumby Lift Chair Dual Motor C Olive Suede

Armrest Height	170 mm	Extended Footrest	1050 mm	Max Backrest Height	820 mm	Max Seat Height	575 mm
Max User Weight	130 kg	Min Backrest Height	670 mm	Min Seat Height	525 mm	Max Seat Depth	540 mm
Min Seat Depth	465 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-BR-M-SM Oscar Furniture Brumby Lift Chair Dual Motor Mini Moss Suede

Armrest Height	190 mm	Extended Footrest	840 mm	Max Backrest Height	710 mm	Max Seat Height	485 mm
Max User Weight	130 kg	Min Backrest Height	560 mm	Min Seat Height	460 mm	Max Seat Depth	510 mm
Min Seat Depth	410 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		


CDALR-BR-M-SN Oscar Furniture Brumby Lift Chair Dual Motor Mini Nutmeg Suede

Armrest Height	190 mm	Extended Footrest	840 mm	Max Backrest Height	710 mm	Max Seat Height	485 mm
Max User Weight	130 kg	Min Backrest Height	560 mm	Min Seat Height	460 mm	Max Seat Depth	510 mm
Min Seat Depth	410 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

CDALR-BR-M-SO Oscar Furniture Brumby Lift Chair Dual Motor Mini Olive Suede

Armrest Height	190 mm	Extended Footrest	840 mm	Max Backrest Height	710 mm	Max Seat Height	485 mm
Max User Weight	130 kg	Min Backrest Height	560 mm	Min Seat Height	460 mm	Max Seat Depth	510 mm
Min Seat Depth	410 mm	Max Seat Width	590 mm	Min Seat Width	490 mm		

Oscar Furniture B24 Mobile Lift Chair



- Easy movement and operation
- Powered by the 24V rechargeable battery system • Large range of options available
- Triple roll back
- Control is standard on the right hand side, (can be fitted to the left hand side)

CHCEL-B24-FB Oscar Furniture B24 Battery Lift Chair Basil Fusion

Armrest Height	190 mm	Backrest Height	630 mm	Floor To Seat Height	480 mm	Extended Footrest	985 mm
Max User Weight	130 kg	Max Seat Depth	550 mm	Min Seat Depth	450 mm	Max Seat Width	590 mm
Min Seat Width	490 mm						

CHCEL-B24-FS Oscar Furniture B24 Battery Lift Chair Sage Fusion

Armrest Height	190 mm	Backrest Height	630 mm	Floor To Seat Height	480 mm	Extended Footrest	985 mm
Max User Weight	130 kg	Max Seat Depth	550 mm	Min Seat Depth	450 mm	Max Seat Width	590 mm
Min Seat Width	490 mm						

CHCEL-B24-FSB Oscar Furniture B24 Battery Lift Chair Steel Blue Fusion

Armrest Height	190 mm	Backrest Height	630 mm	Floor To Seat Height	480 mm	Extended Footrest	985 mm
Max User Weight	130 kg	Max Seat Depth	550 mm	Min Seat Depth	450 mm	Max Seat Width	590 mm
Min Seat Width	490 mm						

Oscar Furniture M5 Bariatric Lift Chair



- Bariatric dual action lift recline chair · Lift assist makes it easy to get into and out of the chair
- · Style, performance, fit and comfort
 - Pressure care visco foam in 4 way stretch fabric
 - Extra wide wood show armrests
 - 75mm swivel locking castors at rear
 - Controls on right side as standard (can be moved to left)

CHCEL-M5-650-SM Oscar Furniture M5 Bariatric Lift Chair 650mm Wide Moss Micro Suede

Backrest Height	690 mm	Floor To Seat Height	500 mm	Max User Weight	240 kg	Seat Depth	530 mm
Seat Width	650 mm						

CHCEL-M5-650-SN Oscar Furniture M5 Bariatric Lift Chair 650mm Wide Nutmeg Micro Suede

Backrest Height	690 mm	Floor To Seat Height	500 mm	Max User Weight	240 kg	Seat Depth	530 mm
Seat Width	650 mm						

CHCEL-M5-650-S0 Oscar Furniture M5 Bariatric Lift Chair 650mm Wide Olive Micro Suede

Backrest Height	690 mm	Floor To Seat Height	500 mm	Max User Weight	240 kg	Seat Depth	530 mm
Seat Width	650 mm						

CHCEL-M5-850-SM Oscar Furniture M5 Bariatric Lift Chair 850mm Wide Moss Micro Suede

Backrest Height	690 mm	Floor To Seat Height	500 mm	Max User Weight	320 kg	Seat Depth	530 mm
Seat Width	850 mm						

CHCEL-M5-850-SN Oscar Furniture M5 Bariatric Lift Chair 850mm Wide Nutmeg Micro Suede

Backrest Height	690 mm	Floor To Seat Height	500 mm	Max User Weight	320 kg	Seat Depth	530 mm
Seat Width	850 mm						

CHCEL-M5-850-S0 Oscar Furniture M5 Bariatric Lift Chair 850mm Wide Olive Micro Suede

Backrest Height	690 mm	Floor To Seat Height	500 mm	Max User Weight	320 kg	Seat Depth	530 mm
Seat Width	850 mm						

Pride LL-805 Lift Chair



- 2-position wall hugger Stylish look, unique design with sewn split back design
- Requires only 100mm space between wall and lounge
- Colour options: Cocoa, Fern, Sandal
- Single motor
- Suits person 5'4" 6'

LL805C Pride TMR-805 Cocoa									
Backrest Height	685 mm	Floor To Seat Height	480 mm	Max User Weight	136 kg	Overall Width	900 mm		
Seat Depth	510 mm	Seat Width	510 mm						



• Medium to large wall hugger chaise lounger

LL805F Pride TMR-805 Fern

Backrest Height	685 mm	Floor To Seat Height	480 mm	Max User Weight	136 kg	Overall Width	900 mm
Seat Depth	510 mm	Seat Width	510 mm				

LL805S Pride TMR-805 Sandal

[Backrest Height	685 mm	Floor To Seat Height	480 mm	Max User Weight	136 kg	Overall Width	900 mm
	Seat Depth	510 mm	Seat Width	510 mm				

Pride L560 Lift Chair



- Small to medium chaise lounger • 3 position, single motor Pillow back design
 - Ultimate in lift chair comfort, style and performance
 - Emergency battery back up system
 - Colour options: ocean , burgundy , spruce
 - Suits person 5'4" and below

L560BF Pride L560 Lift Chair Burgund

LS6UBF Pride LS60 Lift Chair Burgundy										
Backrest Height	635 mm	Floor To Seat Height	480 mm	Max User Weight	170 kg	Seat Depth	480 mm			
Seat Width	480 mm									

L5600F Pride L560 Lift Chair Ocean

Backrest Height	635 mm	Floor To Seat Height	480 mm	Max User Weight	170 kg	Seat Depth	480 mm
Seat Width	480 mm						

L560SF Pride L560 Lift Chair Spruce

Backrest Height	635 mm	Floor To Seat Height	480 mm	Max User Weight	170 kg	Seat Depth	480 mm
Seat Width	480 mm						

Pride Serta Lift Chair



- Features superior comfort, performance and style • A zero gravity position, unique back style, hand controls
 - Super soft Serta perfect fabric selections
 - The ultimate in lift chair comfort
 - Available in a single or dual motor
 - The single motor provides some positioning where as the dual motor provides infinite positioning

SR-358M Pride Serta Lift Chair Single Motor

Backrest Height	620 mm	Floor To Seat Height	510 mm	Max User Weight	170 kg	Seat Depth	510 mm
Seat Width	510 mm						

SR525M Pride Serta Lift Chair Dual Motor

Backre	est Height	670 mm	Floor To Seat Height	510 mm	Max User Weight	170 kg	Seat Depth	510 mm
Seat W	Vidth	510 mm						

Pride LL-806 Lift Chair



- The Pride LL806 Lift Chair is exclusive to Peak Care • The LL-806 lift chair is a wall hugger making it suitable for use in a small room
- - The LL-806 has a stylish look and a unique design
 - With a pillow back design
 - Suitable for clients 5'4" to 6'0" tall. Available in 3 colours

LL806 Pride LL-806 Lift Chair

Backrest Height	685 mm	Floor To Seat Height	480 mm	Max User Weight	136 kg	Seat Depth	510 mm	
Seat Width	510 mm							

Pride T3 Trendelenburg Lift Chair



Provides infinite positioning

- A medium to large lift chair
- Features trendelenburg position
- Elevates legs above head for greater relaxation and a stand up function
- Zero gravity position

T3BL Pride T3 Trendelenburg Lift Chair

		9 0					
Backrest Height	685 mm	Floor To Seat Height	480 mm	Max User Weight	170 kg	Seat Depth	510 mm
Seat Width	510 mm						

Pride LC-107 Lift Chair



- The LC-107 is one of Pride's most versatile lift chairs · Features dual motor, which means infinite positions for maximum comfort
- Independent backrest and footrest movement
- An electrically operated wooden framed, fully upholstered lounge chair Available in burgundy , pecan and ocean colours

LC107 Pride LC-107 Lift Chair

		Line onioi	•					
В	ackrest Height	700 mm	Floor To Seat Height	480 mm	Max User Weight	147 kg	Seat Depth	480 mm
S	eat Width	510 mm						



• Requires only 10cm between the back of the chair and the wall • A 2 position chaise lounger that saves space without sacrificing comfort

Patient Chairs

5.07 Care Chair

R & R Healthcare Equipment Evolution Supreme Chair



- · Made from visco elastic foam in a carflex cover which provides multiple stretch
- The user's weight is distributed over a large area relieving pressure
- Doesn't compromise postural support and allows natural alignment of the spine • With fully adjustable head support, increased lumber support, drop down sides and a wide variety of seating positions
- Delivers a comfortable long term sitting environment and helps protect against pressure sores
- This Chair also allows for the semi fowler position which increases blood circulation and helps the user to breathe more easily
- This chair has an extra high leg position to reduce fluid retention
- This chair is designed to be low maintenance as well as user friendly for both the user and the carer

20009 R & R Healthcare Equipment Evolution Supreme Chair Wide Ice Blue

Backrest Height	800 mm	Floor To Seat Height	620 mm	Max Overall Depth	1810 mm	Max Overall Height	1320 mm
Max User Weight	190 kg	Min Overall Depth	1020 mm	Min Overall Height	930 mm	Overall Width	830 mm
Product Weight	50 kg	Seat Depth	585 mm	Seat Width	585 mm		

20016 R & R Healthcare Equipment Evolution Supreme Chair Small Dot Forest

Backrest Height	800 mm	Floor To Seat Height	620 mm	Max Overall Depth	1810 mm	Max Overall Height	1320 mm
Max User Weight	145 kg	Min Overall Depth	1020 mm	Min Overall Height	930 mm	Overall Width	700 mm
Product Weight	50 kg	Seat Depth	485 mm	Seat Width	455 mm		

20017 R & R Healthcare Equipment Evolution Supreme Chair Small Dot Coffee

Backrest Height	800 mm	Floor To Seat Height	620 mm	Max Overall Depth	1810 mm	Max Overall Height	1320 mm
Max User Weight	145 kg	Min Overall Depth	1020 mm	Min Overall Height	930 mm	Overall Width	700 mm
Product Weight	50 kg	Seat Depth	485 mm	Seat Width	455 mm		

20018 R & R Healthcare Equipment Evolution Supreme Chair Standard Dot Forest

Backrest Height	800 mm	Floor To Seat Height	620 mm	Max Overall Depth	1810 mm	Max Overall Height	1320 mm
Max User Weight	145 kg	Min Overall Depth	1020 mm	Min Overall Height	930 mm	Overall Width	730 mm
Product Weight	50 kg	Seat Depth	585 mm	Seat Width	485 mm		

20019 R & R Healthcare Equipment Evolution Supreme Chair Standard Dot Coffee

Backrest Height	800 mm	Floor To Seat Height	620 mm	Max Overall Depth	1810 mm	Max Overall Height	1320 mm
Max User Weight	145 kg	Min Overall Depth	1020 mm	Min Overall Height	930 mm	Overall Width	730 mm
Product Weight	50 kg	Seat Depth	585 mm	Seat Width	485 mm		

20007 R & R Healthcare Equipment Evolution Supreme Chair Small Ice Blue

Backrest Height	800 mm	Floor To Seat Height	620 mm	Max Overall Depth	1810 mm	Max Overall Height	1320 mm
Max User Weight	145 kg	Min Overall Depth	1020 mm	Min Overall Height	930 mm	Overall Width	700 mm
Product Weight	50 kg	Seat Depth	485 mm	Seat Width	455 mm		

20008 R & R Healthcare Equipment Evolution Supreme Chair Standard Ice Blue

Backrest Height	800 mm	Floor To Seat Height	620 mm	Max Overall Depth	1810 mm	Max Overall Height	1320 mm
Max User Weight	145 kg	Min Overall Depth	1020 mm	Min Overall Height	930 mm	Overall Width	730 mm
Product Weight	50 kg	Seat Depth	585 mm	Seat Width	485 mm		

FOR FULL PRODUCT DETAILS

6. HERCULES **150 PATIENT LIFT**

When you need a lifter that supports multiple users within a home care environment then the Hercules 150 is for you.

The Hercules 150 patient lifter is designed for home use. Ideal for use in areas where there is limited room for larger hoists and yet maintains a lift range from 400mm – 1,600mm.

Hoist patients from the floor or seated positions with ease up to a user weight of 150kg.

TRANSFER

The magnetic controller allows the user to place the controller anywhere on the lifter.

The manual leg operation function is simple to operate with a foot pedal at the rear of the mast, alleviating any strain on the carer.

Multiple push handle heights to accommodate for a range of carer heights.

Comes with a four point yoke as standard. Optional pivot frame head available.

Folds for easy storage and transport.



Don't forget to check out the entire range of CCG Transfer Aid products at countrycaregroup.com.au









HUMAN CARE A TARCEI HOISI

Is a lightweight portable lift with a discreet and modern desian

Its portability ensures that the same lift easily can be used in multiple rooms

Advanced technology platform that allows tracking of lift usage and visualization of collected data

Features a unique telescoping hanger bar that offers different widths in one product

Fits any patient and resulting in a comfortable sitting position

Can be used with either a 2, 3 or 4-point suspension



Transfer Aids

6.01 Mobile Lifting Hoists

R & R Healthcare Equipment Hercules 150 Lifter



Manual leg operation Magnetic controller allows you to place the controller anywhere on the lifter Four point yoke as standard Pivot frame head available

- Multiple push positions on push handle
- 150kg weight capacity
- Australian Standards approved
- Easy roll dual castors
- Folding mast

H150ML R & R Healthcare Equipment Hercules 150 Manual Leg Institution Linak Lifter

Min Lift Height	lin Lift Height 400 mm Max Lift Height		1600 mm	Min Overall Height	1400 mm	Max Overall Height	1850 mm
Overall Length	1220 mm	Overall Width	670 mm	Min Overall Width	670 mm	Max Overall Width	1190 mm
Base Width Open	1200 mm	Base Width Closed	625 mm	Lifter Base Clearance	125 mm	Max User Weight	150 kg

HI50MT R & R Healthcare Equipment Hercules 150 Manual Leg Home Care T Motion Lifter

Min Lift Height	400 mm	Max Lift Height	1600 mm	Min Overall Height	1400 mm	Max Overall Height	1850 mm
Overall Length	1220 mm	Overall Width	670 mm	Min Overall Width	670 mm	Max Overall Width	1190 mm
Base Width Open	1200 mm	Base Width Closed	625 mm	Lifter Base Clearance	125 mm	Max User Weight	150 kg

Oxford Calibre



- Ergonomically designed patient lift for the safe lifting of bariatric patients
- · Sturdy and robust, yet incredibly versatile
- 6 point positioning cradle
- Integrated weigh-scale Impressive safe working load of 385 kg
- Meets the unique demands of bariatric care in today's modern healthcare environment
- Powered leg positioning and twin pivoted front castors

OXF-CALIBRE Oxford Calibre Bariatric Lifter

Min Lift Height	755 mm	Max Lift Height	1930 mm	Min Overall Height	1590 mm	Max Overall Height	2250 mm
Overall Length	1700 mm	Base Width Open	1400 mm	Base Width Closed	845 mm	Turning Radius	1940 mm
Product Weight	112 kg	Max User Weight	385 kg				

Molift Mover 180



- Easily manoeuvred • Large hoisting interval of 27-163 cm • Allows hoisting from high surfaces as well as from the floor • Option of 2 or 4 point suspension Molift's unique push bar provides multiple grip possibilities for caregivers

18000 Molift Mover 180

Min Lift Height	265 mm	Max Lift Height	1630 mm	Overall Height	1940 mm	Overall Length	1250 mm	
Base Width Open	1100 mm	Base Width Closed	695 mm	Product Weight	32 kg	Max User Weight	180 kg	

Molift Partner 255



- Excellent stability
- With a total weight of 42.5 kg and a safe working load of 255kg
- Equipped with environmentally friendly NiMh-battery

13003 Molift Partner 255 Standard Base

Min Lift Height	390 mm	Max Lift Height	1680 mm	Overall Height	1965 mm	Overall Length	1285 mm
Base Width Open	1065 mm	Base Width Closed	690 mm	Product Weight	42.5 kg	Max User Weight	255 kg

13100 Molift Partner 255 with Support Arms

Min Lift Height	390 mm	Max Lift Height	1680 mm	Overall Height	1965 mm	Overall Length	1285 mm
Base Width Open	1065 mm	Base Width Closed	690 mm	Product Weight	50 kg	Max User Weight	255 kg

Molift Smart 150



• Foldable in one piece without using any tools Hoisting range of 27-168 cm Standard 4-point sling suspension

• Easily wheeled and transported

Comfortable and spacious hoisting position for the user

920008 Molift Smart 150

Min Lift Height	270 mm	Max Lift Height	1665 mm	Overall Length	1200 mm	Base Width Open	660 mm
Base Width Closed	940 mm	Turning Radius	650 mm	Product Weight	26 kg	Max User Weight	150 kg



· Lightweight and compact aluminium construction

· Equipped with service notification system that notifies when a service is required

• Ideal solution for hoisting with high complexity situations in acute and critical care areas in both hospitals and care facilities • Excellent hoisting range 37-179 cm allows easy hoisting from the floor as well as high surfaces

• A straight vertical column in combination with the fixed suspension maintains the user at a good distance from the column

· Combined with excellent construction design and wheel solution it becomes very easy to manoeuvre

• Designed for home care, travelling and easy storage

· Small footprint and low weight distribution make manoeuvering in small spaces easy



Molift Mover 300



- With its high hoisting capacity of 300 kg
- Ideal solution for hospitals and care facilities with bariatric clients who have frequent and everyday moving and handling needs
- With a total weight of only 39.3 kg the Molift Mover 300 is the lightest in its category
- The unique construction and wheel positioning provides excellent manoeuvrability
- Molift Mover 300 hoists easily from the floor, reaching down to 37 cm at the lowest
- Excellent maximum hoisting height at 177 cm provides good free hoisting space and possibility for gait training
- The hoist is also equipped with an environmentally friendly NiMh-battery

MM300C Molift Mover 300 Standard

Min Lift Height	375 mm	Max Lift Height	1775 mm	Overall Height	2100 mm	Overall Length	1285 mm				
Base Width Open	1065 mm	Base Width Closed	690 mm	Product Weight	39.3 kg	Max User Weight	300 kg				

Molift Mover 205



Transfer Aids

- Flexible all-round mobile passive hoist Suitable for nursing homes, institutions and hospital settings
- Excellent hoisting height of 36-180 cm
- Low total weight of 39.3 kg facilitating wheeling and manoeuvrability
- Standard 4-point suspension that used together with a Molift Sling provides a comfortable and spacious position for the user
- Quick release function provides option of other suspensions

14000 Molift Mover 205 Patient Lifter with Standard Base

Min Lift Height	355 mm	Max Lift Height	1805 mm	Overall Height	2010 mm	Overall Length	1320 mm		
Base Width Open	1070 mm	Base Width Closed	700 mm	Product Weight	39.3 kg	Max User Weight	205 kg		

14050 Molift Mover 205 Patient Lifter with Ambulating Arms

			5				
Min Lift Height	355 mm	Max Lift Height	1805 mm	Overall Height	2010 mm	Overall Length	1320 mm
Base Width Open	1070 mm	Base Width Closed	700 mm	Product Weight	39.3 kg	Max User Weight	205 kg

Transfer Aids

6.04 Ceiling Hoist - Fixed

Human Care HeliQ Ceiling Hoist



- A modern, discreet look that blends into any environment
- Its small size ensures the lift can easily be used in either institutional or home applications
- Advanced technology platform makes it possible to measure the use and performance of the lift
- Allows a wide adjustment range with 2-point and 4-point hanger bars in different lengths to suit virtually any patient

55120-AU	Human C	are Heli() Ceiling Hoist 300		
Overall Height		200 mm	Quoroll Longth	225 mm	Overall Width

Overall Height	300 mm	Overall Length	235 mm	Overall Width	140 mm	Product Weight	7 kg
Max User Weight	300 kg						


55110-AU	Human Care HeliQ Ceiling Hoist 220
----------	------------------------------------

Overall Height	300 mm	Overall Length	235 mm	Overall Width	140 mm	Product Weight	7 kg
Max User Weight	220 kg						
55100-AU Human C	are Heli	Q Ceiling Hoist 150					
55100-AU Human C Overall Height	are Heli(Q Ceiling Hoist 150 Overall Length	235 mm	Overall Width	140 mm	Product Weight	7 kg

verall Height	300 mm	Overall Length	235 mm	Overall Width	140 mm	Product Weight	7 kg
ax User Weight 220 kg							
100-AU Human Care HeliQ Ceiling Hoist 150							
5 100-AU Human C	are Heli(Q Ceiling Hoist 150					
5 100-AU Human C verall Height	are Heli(300 mm	Q Ceiling Hoist 150 Overall Length	235 mm	Overall Width	140 mm	Product Weight	7 kg
			235 mm	Overall Width	140 mm	Product Weight	7 kg

Molift Air







Transfer Aids

Human Care Altair Ceiling Hoist



- Is a lightweight portable lift with a discreet and modern design • Its portability ensures that the same lift easily can be used in multiple rooms Advanced technology platform that allows tracking of lift usage and visualization of collected data Features a unique telescoping hanger bar that offers different widths in one product
- atient and resulting in a comfortable sitting position ed with either a 2, 3 or 4-point suspension

and the second se	i cuturco u t
	 Fits any pat
	 Can be used
14	

55220-AU Human Care Altair Ceiling Hoist 300							
Overall Height	400 mm	Overall Length	300 mm	Overall Width	140 mm	Product Weight	10.9 kg
Max User Weight	300 kg						
55200-AU Human C	are Alta	ir Ceiling Hoist 150	1			1	

Overall Height	400 mm	Overall Length	300 mm	Overall Width	140 mm	Product Weight	10.9 kg
Max User Weight	150 kg						

55210-AU Human Care Altair Ceiling Hoist 220

Overall Height	400 mm	Overall Length	300 mm	Overall Width	140 mm	Product Weight	10.9 kg
Max User Weight	220 kg						



• An innovative, functional and well-designed hoisting aid for disabled users

• Combined with the Molift Rail System, it solves all hoisting situations in health care and other patient handling environments • Light and flexible, suitable for sitting and horizontal transfers as well as standing and gait training situations

. Easy to handle during installation and servicing, due to a very low weight and the clever quick release function · Flexible solutions such as remote control and possibilities for room-to-room transfer

7.5 kg Max User Weight 300 kg	7.5 kg	Max User Weight	205 kg
7.5 kg Max User Weight 300 kg			
7.5 kg Max User Weight 300 kg			

6.06 Ceiling Hoist - Portable

Transfer Aids



Human Care Roomer S Overhead Lifter **Molift Nomad** • Room-to-room transfer without modifications to doorways Integrated carry handles for easy handling A portable overhead lift offering the possibility to transfer from one room to another in one single motion • Together with the Human Care sling range the Roomer S offers optimum comfort for the patient • Suitable for rooms with sliding doors or fire doors as no modificaiton to the doorway is required • Eliminates the need for turntables and rail switches Can be used with all passive Molift Slings 6 emergency stop 55300-AU Human Care Roomer S Overhead Lifter 22160 Molift Nomad 160 Ceiling Hoist Overall Length 160 mm Overall Width 400 mm Product Weight Max User Weight 220 kg 7 kg Max Lift Height 2000 mm Overall Height Product Weight 6.8 kg Max User Weight **Oxford Voyager Ceiling Hoist Transfer Aids** · Will satisfy the majority of manual handling requirements with a maximum user weight of 190 kgs 4-way model (powered vertically and horizontally) • A full range of track layout options is available • Emergency lower and emergency stop One of the most trusted overhead systems available **Oxford Portable Gantry** • Suitable for institutional or homecare use Quick and easy to use • A maximum user weight of 200 kgs • A full range of track layout options is available VOY-ALB4-4FB-UK Oxford Voyager 420 4-Way Ceiling Hoist Emergency lower and emergency stop One of the most trusted overhead systems available 12.7 kg Product Weight Max User Weight 190 kg Suitable for institutional or homecare use Quick and easy to use **Oxford Voyager Portable Ceiling Hoist** V0Y-9720010 Oxford Portable Gantry • Will satisfy the majority of manual handling requirements Min Overall Height 1830 mm Max Overall Height A maximum user weight of 200 kgs A portable ceiling hoist which is lightweight and removable Max Overall Width 3230 mm Product Weight • With dual controls and 24hour of battery use • A full range of track layout options is available Emergency lower and emergency stop **Molift Rail System** • One of the most trusted overhead systems available Suitable for institutional or homecare use Quick and easy to use

VO	Y-9800005.04.UK Oxford Voyac	ger Portable Ceiling Ho	ist	
Pro	oduct Weight	5.5 kg	Max User Weight	200 kg

160kg to 300kg

MLRAILSYS3M Molift Rail System

MLRAILSYSN3M Molift Rail System



522 mm

160 kg

2406 mm

36 kg

- In combination with a trolley it can climb from one rail system to another without a transition rail . With the extension arm it can be easily connected and disconnected by a single caregiver
- Curved design and integrated 4-point suspension provides a comfortable and spacious position
- · Safety features electronic overload protection, electrical emergency lowering, automatic service notification and
- The weight of the Nomad is only 6.8 kg and has a hoisting capacity up to 255 kg • Can be used to climb from one rail system to another giving easy transfer from room to room

Overall Length	330 mm	Overall Width	200 mm

6.07 Ceiling Hoist - Portable - Associated Accessories

• Will satisfy the majority of manual handling requirements

• A portable gantery which is lightweight and removable

Overall Length	1075 mm	Min Overall Width	2210 mm
Max User Weight	200 kg		

• Has rail profiles with innovative solutions for ceiling integration

- It also includes switches and other transition options for added flexibility
- The rail system has few parts, making it one of the easiest systems to install on the market
- These features combine to make the Molift Rail System extremely competitive in terms of return on investment
- Depending on ceiling hoist motor, Molift Air or Molift Nomad, the system can be designed to offer hoisting capacity from



Human Care Portable Stand Niklas



 Tooless height and width adjustment
Easily managable by one caregiver
 No strain on patient or carer
 Easily transfer patient from the floor, bed, shower, wheelchair

• Lifts and moves heavy patients easily and safely

Human Care Portable Stand Niklas								
	Overall Length	1145 mm	Base Width Open	2450 mm	Base Width Closed	1445 mm	Product Weight	30 kg
	Max User Weight	220 kg						

Transfer Aids

6.08 Hoist Slings

Human Care Patient Lifting Full Body Net Sling



- The Full Body Net Sling is a comfortable sling without reinforcements, handles or padding in the back
- It offers good support for the entire body and can be used in most lifting situations, fits most patients
- The Full Body Net Sling can be used in washing and bathing situations • Its polyester net fabric so can be left behind the patient once the transfer is completed

25005-S Human Care Patient Lifting Full Body Net Sling Small		
Max User Weight 250 kg		
25005-M Human Care Patient Lifting Full Body Net Sling Medium		
Max User Weight 250 kg		
25005-L Human Care Patient Lifting Full Body Net Sling Large		
Max User Weight 250 kg		
25005-XL Human Care Patient Lifting Full Body Net Sling XLarge		
25005-XL Human Care Patient Lifting Full Body Net Sling XLar	ge	

Human Care Patient Lifting Full Body Sling



• It offers good support for the entire back and extra support for the head • This sling can be used in most lifting situations, fits most patients • The reinforced head support makes it suitable for lifting patients from flat surfaces

Max User Weight	
25010-M Human Care Patient Liftin	ag Full Rody Sling Mod
	ing Full Body Sling Med
Max User Weight	
25010-L Human Care Patient Liftin	g Full Body Sling Large
Max User Weight	

25010-S Human Care Patient Lifting Full Body Sling Sma

25010-XL Human Care Patient Lifting Full Body Sling XLa Max User Weight

25010-XXL Human Care Patient Lifting Full Body Sling XX Max User Weight

25010-XXXL Human Care Patient Lifting Full Body Sling X Max User Weight

Human Care Patient Lifting Soft Sling High



 The Silhouette Sling is a comfortable sling without reinforcements, handles or padding • The sling adapts to the patients body and offers good support for the entire back

- seats or seating shells

25020-JS Human Care Patient Lifting Soft Sling High Juni Max User Weight

25020-JM Human Care Patient Lifting Soft Sling High Jun Max User Weight

25020-S Human Care Patient Lifting Soft Sling High Sma Max User Weight

25020-M Human Care Patient Lifting Soft Sling High Med Max User Weight



• The Full Body Sling is a comfortable sling without reinforcements, handles or padding in the back

• The low friction polyester fabric makes the sling application fast and simple

II		
	250 kg	
lium		
	250 kg	
2		
	250 kg	
rge		
	250 kg	Tra
Large	2	Transfer Aids
	250 kg	I A
XXLar	rge	ids
	250 kg	

• This sling can be used in most lifting situations but is especially suitable for transferring to and from wheelchairs, molded

• The thin polyester net material allows moisture and heat to transfer efficiently from the patient . This allows the sling to be left in the wheelchair or bed after the transfer is completed

ior Si	mall
	250 kg
ior N	ledium
	250 kg
II	
	250 kg
lium	
	250 kg

NSW	Health
GOVERNMENT	Health

25020-L Human Care Patient Lifting Soft Sling High Large		
Max User Weight	250 kg	
25020-XL Human Care Patient Lifting Soft Sling High XLarge		
Max User Weight	250 ka	

Human Care Patient Lifting Net High



- The Silhouette Sling is a comfortable sling without reinforcements, handles or padding • The sling adapts to the patients body and offers good support for the entire back
- This sling can be used in most lifting situations but is especially suitable for transferring to and from wheelchairs, molded seats or seating shells
- The thin polyester net material allows moisture and heat to transfer efficiently from the patient
- This allows the sling to be left in the wheelchair or bed after the transfer is completed

25015-JS Human Care Patient Lifting Net High Junior Small		
Max User Weight	250 kg	
25015-JM Human Care Patient Lifting Net High Junior Medium		
Max User Weight	250 kg	
25015-S Human Care Patient Lifting Net High Small		
Max User Weight	250 kg	
25015-M Human Care Patient Lifting Net High Medium		
25015-M Human Care Patient Lifting Net High Medium		
25015-M Human Care Patient Lifting Net High Medium Max User Weight	250 kg	
Max User Weight	250 kg	
	250 kg	
Max User Weight 25015-L Human Care Patient Lifting Net High Large		

Human Care Patient Lifting Sit Sling



- The Sit Sling is a shoulder high model with divided leg parts and reinforced back and legs
- It is designed to give the patient extra back support and is a good solution for patients with poor trunk stability
- This sling can be used in most lifting situations and fits most patients
- The reinforeced back and low friction polyester fabric makes the sling application fast and simple

25030-JS Human Care Patient Lifting Sit Sling Junior Small		
Max User Weight 250 kg		
25030-JM Human Care Patient Lifting Sit Sling Jumior Medium		
Max User Weight 250 kg		

Max User Weight	250 kg
25030-M Human Care	Patient Lifting Sit Sling Medium
Max User Weight	
25030-L Human Care	Patient Lifting Sit Sling Large
Max User Weight	250 kg
25030-XL Human Care	e Patient Lifting Sit Sling XLarge
Max User Weight	250 kg
Human Care Patie	nt Lifting Bath Sling
	 It is designed to give the patient extra back support making the Bath Sling a good solution for patients with poor tr stability This sling can be used in most lifting situations and fits most patients The Bath Sling is made with fast drying polyester net making it suitable for use during bathing and showering
	Lifting Bath Sling Junior Medium
Max User Weight	250 kg
Max User Weight 25035-S Human Care	250 kg Patient Lifting Bath Sling Small
Max User Weight	250 kg
Max User Weight 25035-S Human Care Max User Weight	250 kg Patient Lifting Bath Sling Small
Max User Weight 25035-S Human Care Max User Weight	250 kg Patient Lifting Bath Sling Small 250 kg
Max User Weight 25035-S Human Care Max User Weight 25035-M Human Care Max User Weight	250 kg Patient Lifting Bath Sling Small 250 kg Patient Lifting Bath Sling Medium 250 kg 250 kg 250 kg
Max User Weight 25035-S Human Care Max User Weight 25035-M Human Care Max User Weight 25035-L Human Care	250 kg Patient Lifting Bath Sling Small 250 kg Patient Lifting Bath Sling Medium 250 kg Patient Lifting Bath Sling Large
Max User Weight 25035-S Human Care Max User Weight 25035-M Human Care Max User Weight	250 kg Patient Lifting Bath Sling Small 250 kg Patient Lifting Bath Sling Medium 250 kg 250 kg 250 kg
Max User Weight 25035-S Human Care Max User Weight 25035-M Human Care Max User Weight 25035-L Human Care Max User Weight	250 kg Patient Lifting Bath Sling Small 250 kg Patient Lifting Bath Sling Medium 250 kg Patient Lifting Bath Sling Large
Max User Weight 25035-S Human Care Max User Weight 25035-M Human Care Max User Weight 25035-L Human Care Max User Weight	250 kg Patient Lifting Bath Sling Small 250 kg Patient Lifting Bath Sling Medium 250 kg Patient Lifting Bath Sling Large 250 kg 250 kg
Max User Weight 25035-S Human Care Max User Weight 25035-M Human Care Max User Weight 25035-L Human Care Max User Weight 25035-XL Human Care	250 kg Patient Lifting Bath Sling Small 250 kg Patient Lifting Bath Sling Medium 250 kg Patient Lifting Bath Sling Large 250 kg
Max User Weight 25035-S Human Care Max User Weight 25035-M Human Care Max User Weight 25035-L Human Care Max User Weight 25035-XL Human Care Max User Weight	250 kg Patient Lifting Bath Sling Small 250 kg Patient Lifting Bath Sling Medium 250 kg Patient Lifting Bath Sling Large 250 kg
Max User Weight 25035-S Human Care Max User Weight 25035-M Human Care Max User Weight 25035-L Human Care Max User Weight 25035-XL Human Care Max User Weight	250 kg Patient Lifting Bath Sling Small 250 kg Patient Lifting Bath Sling Medium 250 kg Patient Lifting Bath Sling Large 250 kg

25040-JM Human Care Patient Lifting Amputee Sling Ju Max User Weight

Country Care
GROUP
National Distribution, Local Service

• The sling is made in low friction polyester fabric making sling application fast and simple

unior M	ledium
	250 kg

Transfer Aids


Health
Health

25040-S Human Care Patient Lifting Amputee Sling Small			
Max User Weight	250 kg		
25040-M Human Care Patient Lifting Amputee Sling Medium			
Max User Weight	250 kg		
25040-L Human Care Patient Lifting Amputee Sling Large			
Max User Weight	250 kg		
25040-XL Human Care Patient Lifting Amputee Sling XLarge			
Max User Weight	250 kg		

Human Care Patient Lifting Net Sling



- The Basic Sling is a comfortable sling without reinforcements, handles or padding
- The sling adapts to the patients body and offers good support for the entire back
- This sling can be used in most lifting situations, fits most patients and is a good choice for basic, all-round usage • The thin polyester net material allows moisture and heat to transfer efficiently from the patient
- This allows the sling to be left in the wheelchair or bed after the transfer is completed



25045-JS Human Care Patient Lifting Net Sling Junior Small		
Max User Weight	250 kg	
25045-JM Human Care Patient Lifting Net Sling Junior Mediu	m	
Max User Weight	250 kg	
25045-S Human Care Patient Lifting Net Sling Small		
Max User Weight	250 kg	
25045-M Human Care Patient Lifting Net Sling Medium		
Max User Weight	250 kg	
25045-L Human Care Patient Lifting Net Sling Large		
Max User Weight	250 kg	
25045-XL Human Care Patient Lifting Net Sling XLarge		
Max User Weight	250 kg	

Human Care Patient Lifting Comfort Sling Low



- The Comfort Sling is a soft sling without reinforcements, handles or padding The sling adapts easily to the patients body
- The sling is designed to give the patient comfort and a sense of safety
- The Comfort Sling is a sit on sling without divided leg parts that distributes weight evenly in order to reduce pressure
- The trevira fabric allows it to be left in the chair or wheelchair after the transfer is completed
- The sling is applied while the patient is lying down

25060-JM Human Care Patient Lifting Comfort Sling Low Junio	or I
Max User Weight	25
25060-S Human Care Patient Lifting Comfort Sling Low Small	
Max User Weight	25
25060-M Human Care Patient Lifting Comfort Sling Low Medic	1 11 25
Max User Weight	25
25060-L Human Care Patient Lifting Comfort Sling Low Large	
Max User Weight	25
25060-XL Human Care Patient Lifting Comfort Sling Low XLarg	e

Human Care Patient Lifting Comfort Sling High



Max User Weight

• The sling adapts easily to the patients body • The sling is designed to give the patient comfort and a sense of safety

25065-S Human Care Patient Lifting Comfort Sling High Small Max User Weight

25065-M Human Care Patient Lifting Comfort Sling High Medium Max User Weight

25065-L Human Care Patient Lifting Comfort Sling High Large Max User Weight

25065-XL Human Care Patient Lifting Comfort Sling High XLarge Max User Weight

Oxford Comfort Access Sling Padded with Head Support



• Padded with integrated head support for additional comfort • Provides access to allow the removal of clothing

• It is not a general purpose sling and will only suit around 25% of clients

It is an easy to fit sling

- MX-SL25576 Oxford Comfort Access Sling Padded with Head Support Small Red Max User Weight 68 kg

Country Care
GROUP
National Distribution, Local Service

v Junior Medium

	250 kg
Small	

250 kg

Medium

250 kg

Large

250 kg

250 kg

• The Comfort Sling High is a soft sling without reinforcements, handles or padding

- The Comfort Sling is a sit on sling without divided leg parts that distributes weight evenly in order to reduce pressure
- The trevira fabric allows it to be left in the chair or wheelchair after the transfer is completed
- It's suitable for patients with need for extra head support and the sling is applied while the patient is lying down

250 kg 250 kg 250 kg 250 kg

Designed specifically to facilitate the toileting process

• The sling utilises the securi3 clip attachment system in conjunction with the 4-point cradle

Min User Weight

34 kg



MX-SL25577 Oxford Comfort Access Sling Padded with Head Support Medium Yellow				
Max User Weight	91 kg	Min User Weight	57 kg	
MX-SL25578 Oxford Comfort Access Sling Padded with Head Support Large Green				
Max User Weight	136 kg	Min User Weight	80 kg	

Human Care Patient Lifting Net Sling



- The Basic Sling is a comfortable sling without reinforcements, handles or padding
- The sling adapts to the patients body and offers good support for the entire back • This sling can be used in most lifting situations, fits most patients and is a good choice for basic, all-round usage
- The thin polyester net material allows moisture and heat to transfer efficiently from the patient
- This allows the sling to be left in the wheelchair or bed after the transfer is completed

25050-JS Human Care Patient Lifting Soft Sling Junior Small	
Max User Weight	250 kg
25050-JM Human Care Patient Lifting Soft Sling Junior Mediu	ım
Max User Weight	250 kg
25050-S Human Care Patient Lifting Soft Sling Small	
Max User Weight	250 kg
25050-M Human Care Patient Lifting Soft Sling Medium	
Max User Weight	250 kg
25050-L Human Care Patient Lifting Soft Sling Large	
Max User Weight	250 kg
25050-XL Human Care Patient Lifting Soft Sling XLarge	
Max User Weight	250 kg

Oxford Full Back Sling



- An easy to fit, contoured sling Designed to fit 85-90% of clients
- Offers greater comfort and negates the need for side suspenders
- It has integrated, boned head support and leg padding as standard
- FIM Scores: 1 and 2 and any other person from the floor
- Available in 3 material types: polyester padded, net (for wet environments) and silkfit

MX-SL1068	Oxford Full Back Sling	Small Red		
Max User Weight		68 kg	Min User Weight	34 kg
MX-SL1069 Oxford Full Back Sling Medium Yellow				
Max User Weight		91 kg	Min User Weight	57 kg
MX-SL1070 Oxford Full Back Sling Large Green				
Max User Weight		136 kg	Min User Weight	80 kg

Oxford Quick Fit Padded Sling Polyester



• The Quickfit Sling is an easy-to-fit universal sling Designed to suit approximately 85% of patients It is available with head support and can be fleece lined for additional comfort · Available in polyester, padded polyester and net (for wet environments)

 Special sling adaptions also available on request • Also available with side suspenders

MX-SL1000P	Oxford Quick Fit Padded Sling Poly Small Red		
Max User Weight		68 kg	Min User Weight
MX-SL1003P	Oxford Quick Fit Pad	ded Sling Poly Medium	Yellow
MX-SL1003P Max User Weight	Oxford Quick Fit Pad	ded Sling Poly Medium	Yellow Min User Weight

MX-SLIUU6P	Oxford Quick Fit Pad	ded Sling Poly Large Gr	1
Max User Weight		136 kg	

Oxford Quick Fit Sling Net



Yellow Min User Weight 57 kg MY-SI 1006D Oxford Quick Eit Doddod Sling Doly L reen Min User Weight 80 kg • The Quickfit Sling is an easy-to-fit universal sling Designed to suit approximately 85% of patients • It is available with head support and can be fleece lined for additional comfort · Available in polyester, padded polyester and net (for wet environments) Special sling adaptations also available on request Also available with side suspenders 68 kg Min User Weight 34 kg 91 kg Min User Weight 57 kg 136 kg Min User Weight 80 kg

- MX-SL1001 Oxford Quick Fit Net Sling Small Red Max User Weight MX-SL1004 Oxford Quick Fit Net Sling Medium Yellow Max User Weight MX-SL1007 Oxford Quick Fit Net Sling Large Green Max User Weight

Oxford Access Padded Toileting Sling



It is an easy to fit sling

esidents should be in a sitting position when using this sling

Re				
			X	XXX
	9		P	Į
		2	C	1

SL1601 Oxford Access Toileting Sling Padded Small Red Overall Length 670 mm Overall Width 700 mm

120

Transfer Aids



• The Access Sling has been designed specifically to facilitate the toileting procedure

· Padded for additional comfort which allows access to and removal of clothing

	Max User Weight	68 kg	Min User Weight	34 kg
--	-----------------	-------	-----------------	-------

34 kg

NSW	Health
GOVERNMENT	rieaith

	890 mm	Overall Width	850 mm	Max User Weight	91 kg	Min User Weight	57 kg
11000 0(A To !!!						
		eting Sling Pado	-		100 km		00 hr
Overall Length	1030 mm	Overall Width	1060 mm	Max User Weight	136 kg	Min User Weight	80 kg
Oxford Comf	ort Poly Sl	ing Padded					
		-					
		The Comfort Sline		ontoured sling			
X II III	6	Accommodates 8 Used in conjuncti	5-90% of clients ion with a 4-point cr	adle			
AL	1			ther a seated or reclined	position		
				patented Securi3 clip sys			
		Available in other	materials, the net n	naterial is suitable for we	t environments		
T	T						
	-						
	-						
//X-SL25507 Ox	cford Comfo	rt Poly Sling Pac	dded Small R	ed			
/X-SL25507 OX Max User Weight	cford Comfo	rt Poly Sling Pac	lded Small R	ed Min User Weight	_	34 kg	
	cford Comfo		dded Small R		-	34 kg	
Max User Weight				Min User Weight	-	34 kg	
Max User Weight		68 kg		Min User Weight		34 kg 57 kg	
Max User Weight		68 kg		Min User Weight			
Max User Weight L25508 Oxform Max User Weight	d Comfort P	68 kg	d Medium Ye	Min User Weight Ilow Min User Weight	_		
Max User Weight L25508 Oxform Max User Weight	d Comfort P	68 kg oly Sling Padde 91 kg	d Medium Ye	Min User Weight Ilow Min User Weight			
Max User Weight L25508 Oxford Max User Weight L25578 Oxford	d Comfort P	68 kg oly Sling Padde 91 kg oly Sling Padded	d Medium Ye	Min User Weight Ilow Min User Weight		57 kg	
Max User Weight GL25508 Oxford Max User Weight GL25578 Oxford Max User Weight	d Comfort P d Comfort Po	68 kg oly Sling Padde 91 kg oly Sling Padded 136 kg	d Medium Ye d Large Gree	Min User Weight IIOW Min User Weight N Min User Weight Min User Weight		57 kg	
Max User Weight L25508 Oxford Max User Weight L25578 Oxford Max User Weight	d Comfort P d Comfort Po	68 kg oly Sling Padde 91 kg oly Sling Padded	d Medium Ye d Large Gree	Min User Weight IIOW Min User Weight N Min User Weight Min User Weight		57 kg	

Transfer Aids

- An easy fitting sling with adjustable straps for position control • An inbuilt head support for users with less control
- This sling suits the majority of pivot frame lifters on the market
- Strong but light and breathable
- The sling is made from a quick drying mesh fabric
- Has a 3mm quilted cushioning foam sewn into the leg and head support sections

• Can be used for general lifting from the floor, bed or chair in a lying or sitting position

SL146191 Novis Pivot Sling General Purpose with Head Support Small Max Shoulder Width 640 mm Min Shoulder Width 600 mm 700 mm Max User Weight Min Coccyx To Top of Head 205 kg Max Coccyx To Top of Head 780 mm SL146291 Novis Pivot Sling General Purpose with Head Support Medium

Min Shoulder Width	700 mm	Max Shoulder Width	810 mm	Max User Weight	205 kg	Min Coccyx To Top of Head	750 mm
Max Coccyx To Top of Head	900 mm						

SL146391 Novis Pivot Sling General Purpose with Head Support Large

Min Shoulder Width	760 mm	Max Shoulder Width	870 mm	Max User Weight	300 kg	Min Coccyx To Top of Head	830 mm
Max Coccyx To Top of Head	1000 mm						

Novis Pivot Frame Hygiene Access Sling SL247291 Novis Pivot Frame Sling Hygiene Access Media 530 mm Min Shoulder Width Max Shoulder Width SL247191 Novis Pivot Frame Sling Hygiene Access Small Min Shoulder Width 500 mm Max Shoulder Width

SL247391 Novis Pivot Frame Sling Hygiene Access Large Min Shoulder Width 560 mm Max Shoulder Width

Novis General Purpose Sling Mesh with Head Support



• Versatile sling, easy fitting and ideal for users with less control Suits general lifting from the floor, bed and chair

- The head support to provide a safer and more secure transfer for patients • Made from high grade mesh which can be washed and is quick drying
 - Strong but light and breathable
 - Availbale in yoke or pivot frame models
 - Padded head and legs

SL146181 Novis General Purpose Sling Mesh with Head Support Small

Leg Width	230 mm	Leg Length	540 mm	Overall Length	1070 mm	Head Strap	300 mm
Overall Width	830 mm	Leg Strap	450 mm	Depth	790 mm	Max User Weight	205 kg

SL146281 Novis General Purpose Sling Mesh with Head Support Medium							
Leg Width	240 mm	Leg Length	550 mm	Overall Length	1220 mm	Head Strap	300 mm
Overall Width	950 mm	Leg Strap	450 mm	Depth	920 mm	Max User Weight	205 kg

Leg Width	260 mm	Leg Length	700 mm	Overall Length	1350 mm	Head Strap	300 mm
Overall Width	1155 mm	Leg Strap	450 mm	Depth	1000 mm	Max User Weight	300 kg

SL146481 Novis General Purpose Sling Mesh with Head Support XLarge

Overall Width 1350 mm Leg Strap 450 mm Depth 1150 mm Max User Weight 300 kg	Leg Width	260 mm	Leg Length	750 mm	Overall Length	1500 mm	Head Strap	300 mm
	Overall Width	1350 mm	Leg Strap	450 mm	Depth	1150 mm	Max User Weight	300 kg

SL146581 Novis General Purpose Sling Mesh with Head Support XXLarge

[Leg Width	270 mm	Leg Length	800 mm	Overall Length	1600 mm	Head Strap	300 mm
	Overall Width	1450 mm	Leg Strap	450 mm	Depth	1200 mm	Max User Weight	300 kg



• Is specifically designed for toilet transfer with an open area from mid back to mid thigh

· Suits lifting from commodes or profiled beds with its increased positioning control for correct seating

• The single attachment, reinforced nylon resin pivot clips and buckles makes it suitable for the majority of pivot frame lifters

• Its quick drying mesh fabric is machine washable up to 85°C for microbiological control

• This heavy duty, sling has 3mm of quilted cushioning foam sewn into the leg and head support sections

Manufactured from the highest quality material, in compliance with industry standards

um			
	630 mm	Max User Weight	205 kg
	600 mm	Max User Weight	205 kg
5			
	700 mm	Max User Weight	300 kg

Support Large

123



Novis General Purpose Sling Mesh



- Versatile, easy fitting and provides good support for users with limited upper body control Suits general lifting from the floor, bed or chair
- The 3 loop support strap system controls positioning of head and leg supports
- Can be used with the majority of yoke frame lifters on the market
- Brightly coloured handles and adjusters are secured to the sling for easy access and visibility
- Made from a quick drying mesh fabric
- Has 3mm quilted cushioning foam in the leg sections for greater comfort

SL147181 Novis General Purpose Mesh Yoke Sling Small

			5					
Leg Width	230 mm	Leg Length	670 mm	Overall Length	1160 mm	Head Strap	450 mm	
Overall Width	840 mm	Leg Strap	450 mm	Depth	490 mm	Max User Weight	205 kg	

SL147281 Novis General Purpose Mesh Yoke Sling Medium

Leg Width	240 mm	Leg Length	655 mm	Overall Length	1200 mm	Head Strap	450 mm
Overall Width	890 mm	Leg Strap	450 mm	Depth	545 mm	Max User Weight	205 kg

SL147381 Novis General Purpose Mesh Yoke Sling Large

Leg Width	275 mm	Leg Length	685 mm	Overall Length	1325 mm	Head Strap	450 mm
Overall Width	1135 mm	Leg Strap	440 mm	Depth	640 mm	Max User Weight	300 kg

Novis Hygiene Access Sling with Head Support



- · Specifically designed for toilet transfers, for users with less control who require additional head support
- The open area is from mid back to mid thigh for hygiene access
- There is a buckle waist strap for added security with generous padding under the arms
- The 3 loop support strap system controls positioning of head and leg supports with its brightly coloured handles and adjusters, secured to sling, for easy access and visibility
- The quick drying mesh fabric is machine washable
- This sling suits the majority of yoke frame lifters on the market

SL246281 Novis Hygiene Access with Head Support Medium

Leg Width	190 mm	Leg Length	580 mm	Overall Length	1010 mm	Head Strap	290 mm
Overall Width	1030 mm	Leg Strap	440 mm	Depth	690 mm	Torso Band	980 mm
Max User Weight	205 kg						

SL246181 Novis Hygiene Access with Head Support Small

Leg Width	190 mm	Leg Length	570 mm	Overall Length	995 mm	Head Strap	290 mm
Overall Width	900 mm	Leg Strap	440 mm	Depth	620 mm	Torso Band	850 mm
Max User Weight	205 kg						

SL246381 Novis Hygiene Access with Head Support Large

Leg Width	190 mm	Leg Length	630 mm	Overall Length	1100 mm	Head Strap	290 mm
Overall Width	1100 mm	Leg Strap	440 mm	Depth	800 mm	Torso Band	1050 mm
Max User Weight	300 kg						

SL246481 Novis Hygiene Access with Head Support XLarge

Leg Width	200 mm	Leg Length	700 mm	Overall Length	1230 mm	Head Strap	300 mm
Overall Width	1300 mm	Leg Strap	450 mm	Depth	880 mm	Torso Band	1200 mm
Max User Weight	300 kg						

SL246581	Novis Hygiene Access	with Head Support XXL
----------	----------------------	-----------------------

Leg Width	210 mm	Leg Length	750 mm	Overall Length	1370 mm	Head Strap	300 mm
Overall Width	1450 mm	Leg Strap	450 mm	Depth	980 mm	Torso Band	1300 mm
Max User Weight	300 kg						

Novis Yoke Hygiene Sling



 A buckle waist strap provides added security Strong but light and breathable Recommended for general toilet transfers

- Hook and loop strap for added security • Padded body, legs, under arms and mid section

Hygiene fabric sling

Available for yoke frames

SL267181 Novis Yoke Hygiene Sling Small

		-					
Leg Width	175 mm	Leg Length	625 mm	Overall Length	820 mm	Head Strap	320 mm
Overall Width	940 mm	Leg Strap	645 mm	Depth	195 mm	Torso Band	890 mm
Max User Weight	205 kg						

SL267281 Novis Yoke Hygiene Sling Medium

Leg Width	175 mm	Leg Length	695 mm	Overall Length	890 mm	Head Strap	320 mm
Overall Width	1040 mm	Leg Strap	645 mm	Depth	195 mm	Torso Band	990 mm
Max User Weight	205 kg						

SL267381 Novis Yoke Hygiene Sling Large

Leg Width	175 mm	Leg Length	740 mm	Overall Length	935 mm	Head Strap	320 mm
Overall Width	1170 mm	Leg Strap	655 mm	Depth	195 mm	Torso Band	1120 mm
Max User Weight	300 kg						

Molift RgoSling Ambulating Vest



• Also available are groin straps which not only provide additional support but prevent pinching when more weight is put on the vest

. The front buckle system provides increased possibilities for the user to more independently perform gait training or toileting

1721300 Molift RgoSling Ambulating Vest XXSmall					
Max User Weight		17 kg	I		
1721310 Molift R	RgoSling Ambulati	ing Vest XSmall			
Max User Weight		25 kg	I		
1721320 Molift F	RgoSling Ambulat	ing Vest Small			
Max User Weight		50 kg	I		
1721330 Molift F	RgoSling Ambulat	ing Vest Medium			
Max User Weight		95 kg			

124



Large

• The sling includes an open area from mid back to mid thigh to facilitate easy clothing removal and hygiene access

Developed for rehabilitation exercises such as gait training, standing and balance exercises

Double buckled wide padded belt distributes pressure more evenly

Adopted features from professional climbing equipment to provide the best possible comfort and safety

• The self-tightening suspension loops prevent users from slipping or falling

Min User Weight	12 kg
Min User Weight	17 kg
Min User Weight	25 kg
Min User Weight	45 kg

NSW GOVERNMENT	Health

1721340 Molift RgoSling Ambulat	ing Vest Large		
Max User Weight	160 kg	Min User Weight	90 kg
1721350 Molift RgoSling Ambulat	ing Vest XLarge		
Max User Weight	240 kg	Min User Weight	160 kg
Molift RgoSling Amputee Hi	nhBack		
• Fc • Ar • Le • Ar • Ar • Ar • Ar	ditional straps prevent the legs from	support throughout the entire lift ength for application under both legs n sliding out of position k for users with reduced trunk tone and head stability	
1721200 Molift RgoSling Ampute	e HighBack XXSmall		
Max User Weight	17 kg	Min User Weight	12 kg
1721210 Molift RgoSling Amputee	e HighBack XSmall		
Max User Weight	25 kg	Min User Weight	17 kg
1721220 Molift RgoSling Ampute	e HighBack Small		
Max User Weight	50 kg	Min User Weight	25 kg
1721230 Molift RgoSling Ampute	e HighBack Medium		
Max User Weight	95 kg	Min User Weight	45 kg
1721240 Molift RgoSling Ampute	e HighBack Large		
Max User Weight	160 kg	Min User Weight	90 kg
1721250 Molift RgoSling Ampute	e HighBack XLarge		
Max User Weight	240 kg	Min User Weight	160 kg
1721260 Molift RgoSling Ampute	e HighBack XXLarge		
Max User Weight	300 kg	Min User Weight	230 kg

Molift RgoSling Amputee MediumBack



- Provides additional support necessary for hoisting users with impaired muscle tone or amputated limbs
 Anatomical cut and padding provides support throughout the entire lift
 Leg parts are divided and have extra length for application under both legs

- Additional straps prevent the user's legs from sliding out of position
- Available in mediumback and highback for users with reduced trunk tone and head stability Can be used for hoisting from horizontal and sitting positions

1721100 Molift RgoSling Amputee MediumBack XXSmall			
Max User Weight	17 kg	Min User Weight	12 kg

1721110 Molift RgoSling Amputee MediumBack XSmall								
Max User We	ight	25 kg	Min User Weight					
1721120	1721120 Molift RgoSling Amputee MediumBack Small							
Max User We	ight	50 kg	Min User Weight					
1721130	1721130 Molift RgoSling Amputee MediumBack Medium							
Max User We	ight	95 kg	Min User Weight					
1721140	Molift RgoSling Amputee	MediumBack Large						
Max User We	ight	160 kg	Min User Weight					
1721150	Molift RgoSling Amputee	MediumBack XLarge						
Max User We	ight	240 kg	Min User Weight					
1721160	Molift RgoSling Amputee	MediumBack XXLarge						
Max User We	ight	300 kg	Min User Weight					

Molift RgoSling HighBack Net



All-round sling that fits most users andProvides support for the body and head Allows the user to adopt a slightly inclin
Suitable for hoisting from the floor • The sling is especially suitable for show Net-mesh material that breathes well as • Hoisting straps in several different colo

1720410 Molift RgoSling HighBack Net XSmall Max User Weight 25 kg 1720420 Molift RgoSling HighBack Net Small Max User Weight 50 kg 1720430 Molift RgoSling HighBack Net Medium					
Max User Weight	17 kg	N			
1720410 Molift RgoSling HighBac	k Net XSmall				
Max User Weight	25 kg	N			
1720420 Molift RgoSling HighBack Net Small					
Max User Weight	50 kg	N			
1720430 Molift RgoSling HighBack Net Medium					
Max User Weight	95 kg	N			
1720440 Molift RgoSling HighBac	k Net Large				
Max User Weight	160 kg	N			
1720450 Molift RgoSling HighBac	k Net XLarge				
Max User Weight	240 kg	N			
1720460 Molift RgoSling HighBac	k Net XXLarge				
Max User Weight	300 kg	N			

Transfer Aids

Transfer Aids

Min User Weight	45 kg	
Min User Weight	90 kg	
Min User Weight	160 kg	
Min User Weight	230 kg	
hoisting situations d ned position wer and bathing situations and dries fast pur-coded loops in order to find the right fit		_
Min User Weight	12 kg	
Min User Weight	17 kg	

17 kg

25 kg

in oser weight	2.3 Ky
in User Weight	45 kg
in User Weight	90 kg
in User Weight	160 kg
in User Weight	230 kg

25 kg

Country Care GROUP

Molift RgoSling Toilet HighBack



- A flexible sling with optimal functions for toilet situations Provides exceptional trunk support
- Secured by double belt buckles
- Flexible leg support allows different positioning for user
- Several alternative loops adjust the seating angle according to the user's needs
- Available in low back and high back for users with little or no head stability

 Get weight and 	l pressure in optimal	position whilst	maintaining good	d stability and	l safety

720600 Molift RgoSling Toilet HighBack XXSmall						
Max User Weight	17 kg	Min User Weight	12 kg			
720610 Molift RgoSling Toilet HighBack XSmall						
Max User Weight	25 kg	Min User Weight	17 kg			
720630 Molift RgoSling Toilet HighBack Medium						
Max User Weight	95 kg	Min User Weight	45 kg			
1720620 Molift RgoSling Toilet HighBack Small						
Max User Weight	50 kg	Min User Weight	25 kg			
1720640 Molift RgoSling Toilet Hig	ghBack Large					
Max User Weight	160 kg	Min User Weight	90 kg			
1720650 Molift RgoSling Toilet HighBack XLarge						
Max User Weight	240 kg	Min User Weight	160 kg			
1720660 Molift RgoSling Toilet HighBack XXLarge						
Max User Weight	300 kg	Min User Weight	230 kg			

Molift RgoSling Toilet LowBack



A flexible sling with optimal functions for toilet situations
 Provides exceptional trunk support
 Secured by double belt buckles
Get weight and proceure in optimal position whilet maintainin

- whilst maintaining good stability and safety Get weight and pressure in optimal po
- Flexible leg supports facilitates individual fitting Several alternative loops adjust the seating angle according to the user's needs
- Available in low back and high back for users with little or no head stability

· ·	Availa	11044	Dack	anu	myn	Dack	101	useis	VVILII	iittie	iicau	อเฉมเม

720500 Molift RgoSling Toilet LowBack XXSmall									
Max User Weight	17 kg	Min User Weight	12 kg						
1720510 Molift RgoSling Toilet LowBack XSmall									
Max User Weight	25 kg	Min User Weight	17 kg						
1720520 Molift RgoSling Toilet LowBack Small									
Max User Weight	50 kg	Min User Weight	25 kg						

1720530 Molift RgoSling Toilet LowBack Medium								
Max User Weight	95 kg	1						
1720540 Molift RgoSling Toilet Lo	wBack Large							
Max User Weight	160 kg	1						
1720550 Molift RgoSling Toilet Lo	wBack XLarge							
1720550 Molift RgoSling Toilet Lo Max User Weight	240 kg							
	240 kg							
Max User Weight	240 kg							

Molift RgoSling Comfort HighBack



 Adapts to body shape • Support for the entire back; full body support sling Made of a mesh polyester material that breathes and dries quickly

• Suitable when the sling is to be left underneath the user in a chair or wheelchair, and for users with musclar or skeletal pain A high back model which puts the user into a more posteriorly tilted position
The shoulder straps have the unique sliding loop, which provides smooth transition from lying to sitting position and vice

versa • Must be applied with the user in a lying position

1721420 Molift RgoSling Comfort	HighBack Small
Max User Weight	50 kg
1721430 Molift RgoSling Comfort	HighBack Medium
1721430 Molift RgoSling Comfort Max User Weight	HighBack Medium

1721440 Molift RgoSling Comfort HighBack Large Max User Weight 160 kg

Molift RgoSling MediumBack Net



 Allows the user to adpot a slightly inclined position Suitable for hoisting from the floor

Net-mesh material that breathes well and dries fast

172340 Molift RgoSling MediumBack Net Large								
Max User Weight	160 kg							
172330 Molift RgoSling MediumB	ack Net Medium							
Max User Weight	95 kg							
172320 Molift RgoSling MediumBa	ack Net Small							
Max User Weight	50 kg							

Transfer Aids

Min User Weight	45 kg	
Min User Weight	90 kg	
Min User Weight	160 kg	
Min User Weight	230 kg	

Country Care GROUP ition. Loca

Min User Weight	25 kg
Min User Weight	45 kg
Min User Weight	90 kg

 All-round sling that fits most users and hoisting situations · Provides support for the back and body and is suitable for a patient with headcontrol

• The sling is especially suitable for shower and bathing situations

• Hoisting straps in several different colour-coded loops in order to find the right fit

Min User Weight 90 kg Min User Weight 45 kg Min User Weight 25 kg **Transfer Aids**



172310 Molift RgoSling MediumBack Net XSmall										
Max User Weight 25 kg Min User Weight 17 kg										
172300 Molift RgoSling MediumBack Net XXSmall										
Max User Weight	17 kg	Min User Weight	12 kg							
172350 Molift RgoSling MediumB	ack Net XLarge									
Max User Weight	240 kg	Min User Weight	160 kg							
172360 Molift RgoSling MediumBack Net XXLarge										
Max User Weight	300 kg	Min User Weight	230 kg							

Human Care Patient Lifting Net Sling



- The Basic Sling is a comfortable sling without reinforcements, handles or padding
- The sling adapts to the patients body and offers good support for the entire back • This sling can be used in most lifting situations, fits most patients and is a good choice for basic, all-round usage
- The thin polyester net material allows moisture and heat to transfer efficiently from the patient
- This allows the sling to be left in the wheelchair or bed after the transfer is completed



Human Care Basic Sling Junior Small	
Max User Weight	250 kg
Human Care Basic Sling Junior Medium	
Max User Weight	250 kg
Human Care Basic Sling Small	
Max User Weight	250 kg
Human Care Basic Sling Medium	
Max User Weight	250 kg
Human Care Basic Sling Large	
Max User Weight	250 kg
Human Care Basic Sling XLarge	
Max User Weight	250 kg

Transfer Aids Novis Disposable General Purpose Sling with Head Support



- Is suitable for the majority of yoke frame lifters on the market

SL156281 Novis Disposable General Purpose with Head Support Medium

ſ	Leg Width	240 mm	Leg Length	410 mm	Overall Length	1020 mm	Head Strap	300 mm
	Overall Width	980 mm	Leg Strap	450 mm	Depth	880 mm	Max User Weight	205 kg

SL156381 Novis Disposable General Purpose with Head Support Large

Leg Width	260 mm	Leg Length	700 mm	Overall Length	1350 mm	Head Strap	300 mm
Overall Width	1155 mm	Leg Strap	450 mm	Depth	1000 mm	Max User Weight	205 kg

Transfer Aids

Molift Quick Raiser



• Near perfect weight distribution with a small turning radius Imitates the natural pattern of movement with some degree of weight bearing capacity • Easily moves around toilets, beds, chairs and wheelchairs Compact and easy; fits almost everywhere Low weight, stable and easy to manoeuvre even in narrow spaces • Molift Quick Raiser 1 has a safe working load of 160 kg, extra low base 7 cm and narrow front 34 cm Molift Quick Raiser 2 has a safe working load of 160 kg and a low base 10 cm • Molift Quick Raiser 2+ has been adapted for higher safe working load, up to 200 kg; it has a reinforced column and suspension with a low base 10 cm

4100 Molift Quick Raiser 1

Overall Height 1280 mm Overall Length

30 kg

Max User Weight

4400 Molift Ouick Raiser 2

Product Weight

Overall Height	1730 mm	Overall Length	1150 mm	Overall Width	610 mm	Turning Radius	550 mm
Product Weight	36 kg	Max User Weight	160 kg				

1070 mm

160 kg

4450 Molift Ouick Raiser 2+

Overall Height	1280 mm	Overall Length	1265 mm	Overall Width	610 mm	Turning Radius	725 mm
Product Weight	38.5 kg	Max User Weight	200 kg				



6.09 Single Use Slings for Lifting Hoists

• A versatile general purpose sling with inbuilt head control for users with less controlled body movements

Designed for a single user only for maximum infection control management

• Suits general lifting from the floor, bed or chair and general patient transfers, toileting or hygiene care

• The 3 loop support strap system controls positioning of head and leg supports with its brightly coloured handles and This sling is manufactured from non-woven, spun bond polypropylene for increased comfort against bare skin

• It is heavy duty with a 3 mm quilted cushioning foam in leg sections for greater comfort

6.10 Mobile Standing Hoists

Overall Width	610 mm	Turning Radius	550 mm



Transfer Aids		6.11 Mobile Star	nding Hoist Slings				
Oxford Deluxe Sit to Stand	Sling						
 A padded sling covered in a synthetic fabric, which fits behind the back Features a click together fastening strap, which fastens around the users trunk Webbing loop straps attach to the sling with a choice of three positions 							
MX-SL13498N Oxford Deluxe Sit t	o Stand Sling Large						
Max User Weight	136 kg	Min User Weight	80 kg				
MX-SL13499N Oxford Deluxe Sit t	o Stand Sling Small						
Max User Weight	68 kg	Min User Weight	34 kg				
MX-SL13500N Oxford Deluxe Sit t	o Stand Sling Medium						
Max User Weight	91 kg	Min User Weight	57 kg				
Oxford Deluxe Transport Sli	Oxford Deluxe Transport Sling						
 Easy to fit general purpose sling giving full support to clients This sling is designed for use with stand-aid type lifts only It will not suit double amputees or lifting from a horizontal position FIM Scores: 2, 3 and 4 Special sling adaptions also available on request 							

MX-SL 13501N Oxford Deluxe Transport Sling Medium							
Max User Weight 91 kg Min User Weight 57 kg							
MX-SL13502N Oxford Deluxe Transport Sling Small							
Max User Weight	68 kg	Min User Weight	34 kg				
MX-SL 13503N Oxford Deluxe Transport Sling Large							
Max User Weight	136 kg	Min User Weight	80 kg				

Molift RgoSling StandUp with Support



 Padded for active hoisting • For use with the Molift Quick Raiser

- Perfect lumbar support
 - Ideal sling for moving between different seating positions
 - Ideal for use in toilet situations
 - Anti-slip lining to prevent sliding up
 - Wide-back support distributes weight evenly for a more comfortable lift
 - Waist belt provides additional support

1720810 Molift RgoSling Stand Up with Support XSmall							
Min Waist Size	550 mm	Max Waist Size	600 mm				
1720820 Molift RgoSling Stand Up with Support Small							
Min Waist Size	600 mm	Max Waist Size	750 mm				
1720830 Molift RgoSling Stand Up with Support Medium							
Min Waist Size	750 mm	Max Waist Size	1000 mm				
1720840 Molift RgoSling Stand Up with Support Large							
Min Waist Size	1000 mm	Max Waist Size	1500 mm				
1720850 Molift RgoSling Stand Up with Support XLarge							
Min Waist Size	1500 mm	Max Waist Size	1750 mm				

Molift RgoSling StandUp



- Padded sling for active hoisting
 For use with the Molift Quick Raiser
- Perfect lumbar support
- Ideal hoist for moving between different seating positions • Quickly and easily placed around the user
- Anti-Slip lining to prevent sliding up
- Wide-back support distributes weight evenly for a more comfortable lift
- Waist belt provides additional support
- 1720710 Molift RgoSling StandUp XSmall Min Waist Size 550 mm Max Waist Size 600 mm 1720720 Molift RgoSling StandUp Small 600 mm Max Waist Size Min Waist Size 750 mm 1720730 Molift RgoSling StandUp Medium Min Waist Size 750 mm Max Waist Size 1000 mm 1720740 Molift RgoSling StandUp Large Min Waist Size 1000 mm Max Waist Size 1500 mm 1720750 Molift RgoSling StandUp XLarge Min Waist Size 1500 mm Max Waist Size 1750 mm



Product Weight	2 kg	Max User Weight	300 kg
Product Weight	2 kg	Max User Weight	300 kg
1			
Product Weight	2 kg	Max User Weight	300 kg
Product Weight	2 kg	Max User Weight	300 kg
Product Weight	2 kg	Max User Weight	300 kg

Product Weight	2 kg	Max User Weight	300 kg
Product Weight	2 kg	Max User Weight	300 kg
Product Weight	2 kg	Max User Weight	300 kg
Product Weight	2 kg	Max User Weight	300 kg
Product Weight	2 kg	Max User Weight	300 kg



	G R O U P National Distribution, Local Service
vith Scale	
350 kg	
e Power Cradle with Scale	
350 kg	
e secrui3 safety clip system dle handle either down or up to obtain an upright (available	or reclined position

rg	
227 kg	
227 kg	

· By attaching the sling via various loops, the patient's position can be adjusted from reclined to upright, or the flexion of the

227 kg	
227 kg	
227 kg	

PRESSURE **CUSHIONS**





ROHO OUADTRO SELECT

- Allows for quick and simple position selection for every user while seated
- Push the knob to the open position allows adjustment then push the . knob closed to maintain this position
- Great for maximum function and posture
- Stability is enhanced by locking air into place in each of the four quadrants
- Minimises side to side or front to back motion
- Lightweight, portable, repairable and washable Air is able to circulate around the cells, reaching the surface of the skin and providing a cooling environment



Don't forget to check out the entire range of CCG Pressure Cushion products at countrycaregroup.com.au

Pressure Cushions Non Powered

Skil-Care Bariatric Tri-Foam Contoured Bottom Cushion



 Maximum comfort for bariatric users Three part construction prevents bottoming out • Top layer of visco-elastic foam reduces pressure Middle layer of high-density foam adds comfort • A firm foundation controls sling-seat hammocking • Fluid proof, wipe clean low shear II cover prevents skin damaging friction

• Safety strap keeps cushion in place

SK754821B	Skil-Ca	re Bariat	ric Tri-Foam Contou	red Botto	on
Max User Weight		292 kg	Overall Depth	505 mm	0

Skil-Care Bariatric Tri-Foam Flat Bottom Cushic



 Maximum comfort for bariatric users · Three part construction prevents bottom

- Top layer of visco-elastic foam reduces • Middle layer of high-density foam adds
- A firm foundation controls sling-seat ha
- · Fluid proof, wipe clean low shear II cover prevents skin damaging friction
- Safety strap keeps cushion in place
- Available with contoured base for sling-seat chairs and flat base for solid-seat chairs

SK754833B	Skil-Care Bariatric Tri-Foam Flat Bottom Cush				
Max User Weight		292 kg	Overall Depth	405 mm	0
SK754835B	Skil-Ca	re Baria	tric Tri-Foam Flat Bo	ttom Cu	sł
Max User Weight		292 kg	Overall Depth	455 mm	0
SK754837B	Skil-Ca	re Bariat	ric Tri-Foam Flat Bo	ttom Cus	sh
Max User Weight		292 kg	Overall Depth	405 mm	0
SK754839B	Skil-Ca	re Baria	tric Tri-Foam Flat Bo	ttom Cu	sł
Max User Weight		292 kg	Overall Depth	455 mm	(
Max User Weight SK754841B	Skil-Ca		overall Depth		
	Skil-Ca				
SK754841B		r e Bariat	ric Tri-Foam Flat Bot	ttom Cus	sh
SK754841B Max User Weight		r e Bariat	ric Tri-Foam Flat Bo	ttom Cus	sh
SK754841B Max User Weight SK754843B	Skil-Ca	re Bariat 292 kg re Bariat 292 kg	ric Tri-Foam Flat Bo Overall Depth tric Tri-Foam Flat Bo	455 mm ttom Cus	sh
SK754841B Max User Weight SK754843B Max User Weight	Skil-Ca	re Bariat 292 kg re Bariat 292 kg	ric Tri-Foam Flat Bo Overall Depth tric Tri-Foam Flat Bo Overall Depth	455 mm ttom Cus	sh
SK754841B Max User Weight SK754843B Max User Weight SK754845B	Skil-Ca Skil-Ca	re Bariat 292 kg re Bariat 292 kg re Bariat 292 kg	ric Tri-Foam Flat Bo Overall Depth tric Tri-Foam Flat Bo Overall Depth tric Tri-Foam Flat Bo	ttom Cus 455 mm ttom Cus 505 mm ttom Cus 455 mm	sh sh
SK754841B Max User Weight SK754843B Max User Weight SK754845B Max User Weight	Skil-Ca Skil-Ca	re Bariat 292 kg re Bariat 292 kg re Bariat 292 kg	ric Tri-Foam Flat Bo Overall Depth tric Tri-Foam Flat Bo Overall Depth tric Tri-Foam Flat Bo Overall Depth	ttom Cus 455 mm ttom Cus 505 mm ttom Cus 455 mm	sh c sh



7.01 Pressure Cushions Foam

Available with contoured base for sling-seat chairs and flat base for solid-seat chairs

m Cushion 24" x 20"				
Overall Width	610 mm	Thickness	100 mm	
on				
ning out				
pressure comfort				
immocking				

hion 20" x 16"			
Overall Width	505 mm	Thickness	100 mm
hion 20" x 18"			
Overall Width	505 mm	Thickness	100 mm
hion 22" x 16"			
Overall Width	560 mm	Thickness	100 mm
hion 22" x 18"			
Overall Width	560 mm	Thickness	100 mm
nion 24" x 18"			
Overall Width	610 mm	Thickness	100 mm
hion 24" x 20"			
Overall Width	610 mm	Thickness	100 mm
hion 26" x 18"			
Overall Width	660 mm	Thickness	100 mm
hion 26" x 20"			
Overall Width	660 mm	Thickness	100 mm



Max User Weight	292 kg Overall Dep	th 455 mm	Overall Width	710 mm	Thickness	100 mm
		oam Flat Bottom C		710		100
Max User Weight	292 kg Overall Dep	505 mm	Overall Width	710 mm	Thickness	100 mm
	 Multivary The F Great 	Sure Check Cus Purpose Cushions for the ma with each user Pressure Check Cushion is a lig t for any kind of chair and can Pressure Check Cus Overall Width	nagement of pressure ulcers, ht/heavy foam cushion with be attached for safety			ning and alignment all
ay Ion Cushio	n					
	For o The t Distr The i Its ac	on cushion features a dual-lay ptimized immersion the foam' op layer of soft, visco memory jbutes weight over a larger are nner cover is water-resistant : qua-guard zipper and anti-wic nachine-washable, stretch ou	foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam	shion width and th idual maximum comfor n base from fluids	nerefore user weight t	t
	For o The t Distr The i Its a The r addit	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant t qua-guard zipper and anti-wic nachine-washable, stretch ou ional skin protection	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam	shion width and th idual maximum comfor n base from fluids	nerefore user weight t	t
,	• For o • The t • Distr • The i • Its au • The r • addit	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant a qua-guard zipper and anti-wic nachine-washable, stretch ou ional skin protection	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam ter cover reduces surface ten	shion width and th idual maximum comfor n base from fluids sion and promote	t s airflow, dissipates	t s heat and moisture for
ION1823 Jay Ion Max User Weight	For o The t Distr The i Its a The r addit	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant t qua-guard zipper and anti-wic nachine-washable, stretch ou ional skin protection	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam	shion width and th idual maximum comfor n base from fluids	t s airflow, dissipates	t
Max User Weight	• For o • The t • Distr • The i • Its au • The r • addit	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant : qua-guard zipper and anti-wic nachine-washable, stretch ou tional skin protection	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam ter cover reduces surface ten	shion width and th idual maximum comfor n base from fluids sion and promote	t s airflow, dissipates	t s heat and moisture for
Max User Weight	For o The t Distr The i Its ac The r addit Cushion 18" x 23" 136 kg	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant : qua-guard zipper and anti-wic nachine-washable, stretch ou tional skin protection	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam ter cover reduces surface ten	shion width and th idual maximum comfor n base from fluids sion and promote	nerefore user weight t s airflow, dissipates	t s heat and moisture for
Max User Weight	Cushion 18" x 24"	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant t qua-guard zipper and anti-wic nachine-washable, stretch ou ional skin protection Overall Depth	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam ter cover reduces surface ten 585 mm	shion width and th idual maximum comfor n base from fluids sion and promote	nerefore user weight t s airflow, dissipates	t s heat and moisture for 455 mm
Max User Weight	Cushion 18" x 24"	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant t qua-guard zipper and anti-wic nachine-washable, stretch ou ional skin protection Overall Depth	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam ter cover reduces surface ten 585 mm	shion width and th idual maximum comfor n base from fluids sion and promote	t s airflow, dissipates	t s heat and moisture for 455 mm
Max User Weight IION1824 Jay Ion Max User Weight IION1915 Jay Ion (Max User Weight	• For o • The t • Distr • The i • Its au • The r • The r • Its au • Its	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant t qua-guard zipper and anti-wic machine-washable, stretch ou cional skin protection Overall Depth	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam ter cover reduces surface ten 585 mm 610 mm	shion width and th idual maximum comfor n base from fluids sion and promote Overall Width	t s airflow, dissipates	t s heat and moisture for 455 mm 455 mm
Max User Weight ION1824 Jay Ion Max User Weight ION1915 Jay Ion (Max User Weight ION1916 Jay Ion (For o The t Distr The i Distr The i Its ar The r addil Cushion 18" x 23" 136 kg Cushion 18" x 24" 136 kg Cushion 19" x 15" 136 kg Cushion 19" x 15"	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant t ua-guard zipper and anti-wic machine-washable, stretch ou cional skin protection Overall Depth Overall Depth	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam ter cover reduces surface ten 585 mm 610 mm 380 mm	shion width and th idual maximum comfor n base from fluids sion and promote Overall Width Overall Width	t s airflow, dissipates	t s heat and moisture for 455 mm 455 mm 480 mm
Max User Weight ION1824 Jay Ion Max User Weight ION1915 Jay Ion (Max User Weight	• For o • The t • Distr • The i • Its au • The r • The r • Its au • Its	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant t qua-guard zipper and anti-wic machine-washable, stretch ou cional skin protection Overall Depth	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam ter cover reduces surface ten 585 mm 610 mm	shion width and th idual maximum comfor n base from fluids sion and promote Overall Width	t s airflow, dissipates	t s heat and moisture for 455 mm 455 mm
Max User Weight IION1824 Jay Ion Max User Weight IION1915 Jay Ion (Max User Weight IION1916 Jay Ion (Max User Weight	For o The t Distr The i Distr The i Its ar The r addit Cushion 18" x 23" 136 kg Cushion 18" x 24" 136 kg Cushion 19" x 15" 136 kg Cushion 19" x 16" 136 kg	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant t ua-guard zipper and anti-wic machine-washable, stretch ou cional skin protection Overall Depth Overall Depth	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam ter cover reduces surface ten 585 mm 610 mm 380 mm	shion width and th idual maximum comfor n base from fluids sion and promote Overall Width Overall Width	t s airflow, dissipates	t s heat and moisture for 455 mm 455 mm 480 mm
Max User Weight IION1824 Jay Ion Max User Weight IION1915 Jay Ion (Max User Weight IION1916 Jay Ion (Max User Weight	For o The t Distr The i Distr The i Its ar The r addit Cushion 18" x 23" 136 kg Cushion 18" x 24" 136 kg Cushion 19" x 15" 136 kg Cushion 19" x 16" 136 kg	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant t ua-guard zipper and anti-wic machine-washable, stretch ou cional skin protection Overall Depth Overall Depth	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam ter cover reduces surface ten 585 mm 610 mm 380 mm	shion width and th idual maximum comfor n base from fluids sion and promote Overall Width Overall Width	erefore user weight t s airflow, dissipates	t s heat and moisture for 455 mm 455 mm 480 mm
Max User Weight IION1824 Jay Ion Max User Weight IION1915 Jay Ion (Max User Weight IION1916 Jay Ion (Max User Weight IION1917 Jay Ion (Max User Weight	For o The t Distr Tis k Distr Dis Dis Dis Dis Dis	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant i qua-guard zipper and anti-wic nachine-washable, stretch ou ional skin protection Overall Depth Overall Depth Overall Depth	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam ter cover reduces surface ten 585 mm 610 mm 380 mm 405 mm	shion width and th idual maximum comfor n base from fluids sion and promote Overall Width Overall Width	erefore user weight t s airflow, dissipates	t s heat and moisture for 455 mm 455 mm 455 mm 480 mm 480 mm
Max User Weight IION1824 Jay Ion Max User Weight IION1915 Jay Ion (Max User Weight IION1916 Jay Ion (Max User Weight IION1917 Jay Ion (Max User Weight IION1918 Jay Ion (Max User Weight IION1918 Jay Ion (Max User Weight	For o The t Distr Distr The t Distr The t Distr The t Distr Distret Distret Distret Distret Distret Distr Dis Dis Distr Distr Dist	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant a qua-guard zipper and anti-wic nachine-washable, stretch ou ional skin protection Overall Depth Overall Depth Overall Depth Overall Depth	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam ter cover reduces surface ten 585 mm 610 mm 380 mm 405 mm 430 mm	shion width and th idual maximum comfor n base from fluids sion and promote Overall Width Overall Width Overall Width	erefore user weight t s airflow, dissipates	t s heat and moisture for 455 mm 455 mm 455 mm 480 mm 480 mm 480 mm
Max User Weight IION1824 Jay Ion Max User Weight IION1915 Jay Ion (Max User Weight IION1916 Jay Ion (Max User Weight IION1917 Jay Ion (Max User Weight	For o The t Distr Tis k Distr Dis Dis Dis Dis Dis	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant i qua-guard zipper and anti-wic nachine-washable, stretch ou ional skin protection Overall Depth Overall Depth Overall Depth	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam ter cover reduces surface ten 585 mm 610 mm 380 mm 405 mm	shion width and th idual maximum comfor n base from fluids sion and promote Overall Width Overall Width	erefore user weight t s airflow, dissipates	t s heat and moisture for 455 mm 455 mm 455 mm 480 mm 480 mm
Max User Weight IION1824 Jay Ion Max User Weight IION1915 Jay Ion (Max User Weight IION1916 Jay Ion (Max User Weight IION1917 Jay Ion (Max User Weight IION1918 Jay Ion (Max User Weight IION1918 Jay Ion (Max User Weight	• For o • The t • Distr • Distr • The t • Distr • The t • The t • The r • 136 kg Cushion 19" x 17" • 136 kg Cushion 19" x 18" • 136 kg	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant a qua-guard zipper and anti-wic nachine-washable, stretch ou ional skin protection Overall Depth Overall Depth Overall Depth Overall Depth	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam ter cover reduces surface ten 585 mm 610 mm 380 mm 405 mm 430 mm	shion width and th idual maximum comfor n base from fluids sion and promote Overall Width Overall Width Overall Width	erefore user weight t s airflow, dissipates	t s heat and moisture for 455 mm 455 mm 455 mm 480 mm 480 mm 480 mm
Max User Weight IION1824 Jay Ion Max User Weight IION1915 Jay Ion (Max User Weight IION1916 Jay Ion (Max User Weight IION1917 Jay Ion (Max User Weight IION1918 Jay Ion (Max User Weight	• For o • The t • Distr • Distr • The t • Distr • The t • The t • The r • 136 kg Cushion 19" x 17" • 136 kg Cushion 19" x 18" • 136 kg	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant a qua-guard zipper and anti-wic nachine-washable, stretch ou ional skin protection Overall Depth Overall Depth Overall Depth Overall Depth	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam ter cover reduces surface ten 585 mm 610 mm 380 mm 405 mm 430 mm	shion width and th idual maximum comfor n base from fluids sion and promote Overall Width Overall Width Overall Width	erefore user weight t s airflow, dissipates	t s heat and moisture for 455 mm 455 mm 455 mm 480 mm 480 mm 480 mm
Max User Weight IION1824 Jay Ion Max User Weight IION1915 Jay Ion (Max User Weight IION1916 Jay Ion (Max User Weight IION1917 Jay Ion (Max User Weight IION1918 Jay Ion (Max User Weight IION1918 Jay Ion (Max User Weight IION1919 Jay Ion (Max User Weight	• For o • The t • Distr • The t • Distr • The t • The t <td< td=""><td>ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant a qua-guard zipper and anti-wic machine-washable, stretch ou ional skin protection Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth</td><td>s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam ter cover reduces surface ten 610 mm 610 mm 405 mm 405 mm 455 mm</td><td>shion width and th idual maximum comfor n base from fluids sion and promote Overall Width Overall Width Overall Width</td><td>erefore user weight t s airflow, dissipates</td><td>t s heat and moisture for 455 mm 455 mm 455 mm 480 mm 480 mm 480 mm 480 mm</td></td<>	ptimized immersion the foam op layer of soft, visco memory ibutes weight over a larger are nner cover is water-resistant a qua-guard zipper and anti-wic machine-washable, stretch ou ional skin protection Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth	s firmness increases with cus foam conforms to each indiv a for pressure reduction and and easy to wipe off and clea king thread protect the foam ter cover reduces surface ten 610 mm 610 mm 405 mm 405 mm 455 mm	shion width and th idual maximum comfor n base from fluids sion and promote Overall Width Overall Width Overall Width	erefore user weight t s airflow, dissipates	t s heat and moisture for 455 mm 455 mm 455 mm 480 mm 480 mm 480 mm 480 mm

JION1921 Jay Ion Cushio	n 19" x 21"	
Max User Weight	136 kg	Overall Depth
JION1922 Jay Ion Cushic	n 10" v 22"	
Max User Weight	136 kg	Overall Depth
	lookg	
JION1923 Jay Ion Cushic	on 19" x 23"	
Max User Weight	136 kg	Overall Depth
JION1924 Jay Ion Cushic	on 19" x 24"	
Max User Weight	136 kg	Overall Depth
1000015		
JION2015 Jay Ion Cushic		Querell Denth
Max User Weight	136 kg	Overall Depth
JION2016 Jay Ion Cushic	on 20" x 16"	
Max User Weight	136 kg	Overall Depth
JION2017 Jay Ion Cushic	n 20″ v 17″	
Max User Weight	136 kg	Overall Depth
JION2018 Jay Ion Cushic	on 20" x 18"	
Max User Weight	136 kg	Overall Depth
JION2019 Jay Ion Cushic	on 20" x 19"	
Max User Weight	136 kg	Overall Depth
· · · · · · · · · · · · · · · · · · ·	130 Kg	orenan bepan
JION2021 Jay Ion Cushic	on 20" x 21"	
		Overall Depth
JION2021 Jay Ion Cushic	on 20" x 21"	
JION2021 Jay Ion Cushic Max User Weight	on 20" x 21"	
JION2021 Jay Ion Cushic Max User Weight JION2023 Jay Ion Cushic Max User Weight	on 20" x 21" 136 kg on 20" x 23" 136 kg	Overall Depth
JION2021 Jay Ion Cushic Max User Weight JION2023 Jay Ion Cushic	on 20" x 21" 136 kg on 20" x 23" 136 kg	Overall Depth
JION2021 Jay Ion Cushic Max User Weight JION2023 Jay Ion Cushic Max User Weight JION2024 Jay Ion Cushic Max User Weight	on 20" x 21" 136 kg on 20" x 23" 136 kg on 20" x 24" 136 kg	Overall Depth Overall Depth
JION2021 Jay Ion Cushid Max User Weight JION2023 Jay Ion Cushid Max User Weight JION2024 Jay Ion Cushid Max User Weight JION2117 Jay Ion Cushio	on 20" x 21" 136 kg on 20" x 23" 136 kg on 20" x 24" 136 kg n 21" x 17"	Overall Depth Overall Depth Overall Depth
JION2021 Jay Ion Cushic Max User Weight JION2023 Jay Ion Cushic Max User Weight JION2024 Jay Ion Cushic Max User Weight	on 20" x 21" 136 kg on 20" x 23" 136 kg on 20" x 24" 136 kg	Overall Depth Overall Depth
JION2021 Jay Ion Cushid Max User Weight JION2023 Jay Ion Cushid Max User Weight JION2024 Jay Ion Cushid Max User Weight JION2117 Jay Ion Cushio	on 20" x 21" 136 kg on 20" x 23" 136 kg on 20" x 24" 136 kg n 21" x 17" 136 kg	Overall Depth Overall Depth Overall Depth
JION2021 Jay Ion Cushid Max User Weight JION2023 Jay Ion Cushid Max User Weight JION2024 Jay Ion Cushid Max User Weight JION2117 Jay Ion Cushio Max User Weight	on 20" x 21" 136 kg on 20" x 23" 136 kg on 20" x 24" 136 kg n 21" x 17" 136 kg	Overall Depth Overall Depth Overall Depth
JION2021 Jay Ion Cushid Max User Weight JION2023 Jay Ion Cushid Max User Weight JION2024 Jay Ion Cushid Max User Weight JION2117 Jay Ion Cushio Max User Weight JION2118 Jay Ion Cushio Max User Weight	on 20" x 21" 136 kg on 20" x 23" 136 kg on 20" x 24" 136 kg n 21" x 17" 136 kg n 21" x 18" 136 kg	Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth
JION2021 Jay Ion Cushic Max User Weight JION2023 Jay Ion Cushic Max User Weight JION2024 Jay Ion Cushic Max User Weight JION2117 Jay Ion Cushic Max User Weight JION2118 Jay Ion Cushic Max User Weight JION2119 Jay Ion Cushic	on 20" x 21" 136 kg on 20" x 23" 136 kg on 20" x 24" 136 kg n 21" x 17" 136 kg n 21" x 18" 136 kg n 21" x 18" 136 kg	Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth
JION2021 Jay Ion Cushid Max User Weight JION2023 Jay Ion Cushid Max User Weight JION2024 Jay Ion Cushid Max User Weight JION2117 Jay Ion Cushio Max User Weight JION2118 Jay Ion Cushio Max User Weight	on 20" x 21" 136 kg on 20" x 23" 136 kg on 20" x 24" 136 kg n 21" x 17" 136 kg n 21" x 18" 136 kg	Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth
JION2021 Jay Ion Cushic Max User Weight JION2023 Jay Ion Cushic Max User Weight JION2024 Jay Ion Cushic Max User Weight JION2117 Jay Ion Cushic Max User Weight JION2118 Jay Ion Cushic Max User Weight JION2119 Jay Ion Cushic	on 20" x 21" 136 kg on 20" x 23" 136 kg on 20" x 24" 136 kg n 21" x 17" 136 kg n 21" x 18" 136 kg n 21" x 18" 136 kg n 21" x 19" 136 kg	Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth
JION2021 Jay Ion Cushic Max User Weight JION2023 Jay Ion Cushic Max User Weight JION2024 Jay Ion Cushic Max User Weight JION2117 Jay Ion Cushic Max User Weight JION2118 Jay Ion Cushic Max User Weight JION2119 Jay Ion Cushic Max User Weight	on 20" x 21" 136 kg on 20" x 23" 136 kg on 20" x 24" 136 kg n 21" x 17" 136 kg n 21" x 18" 136 kg n 21" x 18" 136 kg n 21" x 19" 136 kg	Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth
JION2021 Jay Ion Cushic Max User Weight JION2023 Jay Ion Cushic Max User Weight JION2024 Jay Ion Cushic Max User Weight JION2117 Jay Ion Cushic Max User Weight JION2118 Jay Ion Cushic Max User Weight JION2119 Jay Ion Cushic Max User Weight JION2120 Jay Ion Cushic Max User Weight	on 20" x 21" 136 kg on 20" x 23" 136 kg on 20" x 24" 136 kg n 21" x 17" 136 kg n 21" x 18" 136 kg n 21" x 19" 136 kg on 21" x 20" 136 kg	Overall Depth
JION2021 Jay Ion Cushic Max User Weight JION2023 Jay Ion Cushic Max User Weight JION2024 Jay Ion Cushic Max User Weight JION2117 Jay Ion Cushic Max User Weight JION2118 Jay Ion Cushic Max User Weight JION2119 Jay Ion Cushic Max User Weight JION2120 Jay Ion Cushic Max User Weight	on 20" x 21" 136 kg on 20" x 23" 136 kg on 20" x 24" 136 kg n 21" x 17" 136 kg n 21" x 18" 136 kg n 21" x 19" 136 kg on 21" x 20" 136 kg	Overall Depth
JION2021 Jay Ion Cushic Max User Weight JION2023 Jay Ion Cushic Max User Weight JION2024 Jay Ion Cushic Max User Weight JION2117 Jay Ion Cushic Max User Weight JION2118 Jay Ion Cushic Max User Weight JION2120 Jay Ion Cushic Max User Weight JION2121 Jay Ion Cushic Max User Weight	on 20" x 21" 136 kg on 20" x 23" 136 kg on 20" x 24" 136 kg n 21" x 17" 136 kg n 21" x 18" 136 kg n 21" x 19" 136 kg on 21" x 20" 136 kg n 21" x 21" 136 kg	Overall Depth
JION2021 Jay Ion Cushic Max User Weight JION2023 Jay Ion Cushic Max User Weight JION2024 Jay Ion Cushic Max User Weight JION2117 Jay Ion Cushic Max User Weight JION2118 Jay Ion Cushic Max User Weight JION2119 Jay Ion Cushic Max User Weight JION2120 Jay Ion Cushic Max User Weight JION2121 Jay Ion Cushic	on 20" x 21" 136 kg on 20" x 23" 136 kg on 20" x 24" 136 kg n 21" x 17" 136 kg n 21" x 18" 136 kg n 21" x 19" 136 kg on 21" x 20" 136 kg n 21" x 21" 136 kg	Overall Depth

535 mm	Overall Width	480 mm
560 mm	Overall Width	480 mm
585 mm	Overall Width	480 mm
610 mm	Overall Width	480 mm
380 mm	Overall Width	505 mm
405 mm	Overall Width	505 mm
430 mm	Overall Width	505 mm
455 mm	Overall Width	505 mm
480 mm	Overall Width	505 mm
535 mm	Overall Width	505 mm
585 mm	Overall Width	505 mm
610 mm	Overall Width	505 mm
430 mm	Overall Width	535 mm
455 mm	Overall Width	535 mm
480 mm	Overall Width	535 mm
505 mm	Overall Width	535 mm
_		
535 mm	Overall Width	535 mm
_		
560 mm	Overall Width	535 mm

Country Care

oution. Local Service

JION2123 Jay Ion Cushi	on 21" x 23"				
Max User Weight	136 kg	Overall Depth	585 mm	Overall Width	535 mm
JION2124 Jay Ion Cushi	on 21" x 24"				
Max User Weight	136 kg	Overall Depth	610 mm	Overall Width	535 mm
JION2217 Jay Ion Cushi	on 22" x 17"				
Max User Weight	136 kg	Overall Depth	430 mm	Overall Width	560 mm
JION2218 Jay Ion Cushi	on 22" x 18"				
Max User Weight	136 kg	Overall Depth	455 mm	Overall Width	560 mm
JION2219 Jay Ion Cushi	on 22" x 19"				
Max User Weight	136 kg	Overall Depth	480 mm	Overall Width	560 mm
JION2220 Jay Ion Cush	on 22" x 20"				
Max User Weight	136 kg	Overall Depth	505 mm	Overall Width	560 mm
JION2221 Jay Ion Cushi	on 22" v 21"				
Max User Weight	136 kg	Overall Depth	535 mm	Overall Width	560 mm
		or and other			
JION2222 Jay Ion Cush	on 22" x 22"				
Max User Weight	136 kg	Overall Depth	560 mm	Overall Width	560 mm
JION2223 Jay Ion Cush	on 22" x 23"				
Max User Weight	136 kg	Overall Depth	585 mm	Overall Width	560 mm
JION2224 Jay Ion Cushi	on 22" x 24"				
Max User Weight	136 kg	Overall Depth	610 mm	Overall Width	560 mm
JION2319 Jay Ion Cushi	on 23" x 19"				
Max User Weight	136 kg	Overall Depth	480 mm	Overall Width	585 mm
JION2320 Jay Ion Cush	ion 23" x 20"				
Max User Weight	136 kg	Overall Depth	505 mm	Overall Width	585 mm
JION2321 Jay Ion Cushi	on 23" x 21"				
Max User Weight	136 kg	Overall Depth	535 mm	Overall Width	585 mm
JION2322 Jay Ion Cushi	on 23" x 22"				
Max User Weight	136 kg	Overall Depth	560 mm	Overall Width	585 mm
JION2323 Jay Ion Cushi	on 23" x 23"				
Max User Weight	136 kg	Overall Depth	585 mm	Overall Width	585 mm
JION2324 Jay Ion Cushi	on 23" x 24"				
-		Overall Depth	610 mm	Overall Width	585 mm
Max User Weight	136 kg				
-	-			·	
Max User Weight JION2419 Jay Ion Cushi Max User Weight	-	Overall Depth	480 mm	Overall Width	610 mm
JION2419 Jay Ion Cushi	on 24" x 19"		480 mm	Overall Width	610 mm

JION2421 Jay Ion Cushi	on 24" x 21"	
Max User Weight	136 kg	Overall Depth
JION2422 Jay Ion Cushi	ion 24" x 22"	
Max User Weight	136 kg	Overall Depth
		1
JION2423 Jay Ion Cushi		
Max User Weight	136 kg	Overall Depth
JION2424 Jay Ion Cushi	ion 24" x 24"	
Max User Weight	136 kg	Overall Depth
JION1414 Jay Ion Cushi	on 14" x 14"	
Max User Weight	136 kg	Overall Depth
JION1415 Jay Ion Cushie	on 14" x 15"	
Max User Weight	136 kg	Overall Depth
	on 1/// v 10//	
JION1416 Jay Ion Cushie Max User Weight	DN 14" X 16" 136 kg	Overall Depth
mux out neight	100 Ng	oreitan beptit
JION1417 Jay Ion Cushid	on 14" x 17"	
Max User Weight	136 kg	Overall Depth
JION1418 Jay Ion Cushi	on 14" x 18"	
Max User Weight	136 kg	Overall Depth
	15% 14%	
JION1514 Jay Ion Cushid		Queroll Denth
Max User Weight	136 kg	Overall Depth
JION1515 Jay Ion Cushie	on 15" x 15"	
Max User Weight	136 kg	Overall Depth
JION1516 Jay Ion Cushid	on 15" x 16"	
Max User Weight	136 kg	Overall Depth
JION1517 Jay Ion Cushie		
Max User Weight	136 kg	Overall Depth
JION1518 Jay Ion Cushie	on 15" x 18"	
Max User Weight	136 kg	Overall Depth
JION1519 Jay Ion Cushid	on 15" y 10"	
JION1519 Jay Ion Cushio	136 kg	Overall Depth
mux user mergin	100 Kg	
JION1520 Jay Ion Cushi	on 15" x 20"	
Max User Weight	136 kg	Overall Depth
JION1521 Jay Ion Cushid	on 15" x 21"	
Max User Weight	136 kg	Overall Depth
JION1614 Jay Ion Cushi		
Max User Weight	136 kg	Overall Depth

	535 mm	Overall Width	610 mm
	560 mm	Overall Width	610 mm
	585 mm	Overall Width	610 mm
	610 mm	Overall Width	610 mm
	355 mm	Overall Width	355 mm
	380 mm	Overall Width	355 mm
	405 mm	Overall Width	355 mm
	430 mm	Overall Width	355 mm
	455 mm	Overall Width	355 mm
	355 mm	Overall Width	380 mm
	380 mm	Overall Width	380 mm
_			
	405 mm	Overall Width	380 mm
_			
	430 mm	Overall Width	380 mm
_			
	455 mm	Overall Width	380 mm
	480 mm	Overall Width	380 mm
	505 mm	Overall Width	380 mm
	535 mm	Overall Width	380 mm
	355 mm	Overall Width	405 mm



JION1615 Jay Ion Cushi					
Max User Weight	136 kg	Overall Depth	380 mm	Overall Width	405 mm
JION1617 Jay Ion Cushi	on 16" x 17"				
Max User Weight	136 kg	Overall Depth	430 mm	Overall Width	405 mm
110N1619 Joy Jon Cuchi	op 16" v 10"				
JION1618 Jay Ion Cushi		Querell Death	455 mm	Overall Width	405
Max User Weight	136 kg	Overall Depth	455 mm	Overall Width	405 mm
JION1619 Jay Ion Cushi	on 16" x 19"				
Max User Weight	136 kg	Overall Depth	480 mm	Overall Width	405 mm
JION1620 Jay Ion Cush	ion 16" x 20"				
Max User Weight	136 kg	Overall Depth	505 mm	Overall Width	405 mm
-					
JION1621 Jay Ion Cushi	on 16" x 21"		_		
Max User Weight	136 kg	Overall Depth	535 mm	Overall Width	405 mm
JION1622 Jay Ion Cush	ion 16" x 22"				
Max User Weight	136 kg	Overall Depth	560 mm	Overall Width	405 mm
-					
JION1714 Jay Ion Cushi	on 17" x 14"		_		
Max User Weight	136 kg	Overall Depth	355 mm	Overall Width	430 mm
JION1715 Jay Ion Cushi	on 17" x 15"				
Max User Weight	136 kg	Overall Depth	380 mm	Overall Width	430 mm
JION1716 Jay Ion Cushi	on 17" x 16"		_		
Max User Weight	136 kg	Overall Depth	405 mm	Overall Width	430 mm
JION1717 Jay Ion Cushi	on 17" x 17"				
Max User Weight	136 kg	Overall Depth	430 mm	Overall Width	430 mm
JION1718 Jay Ion Cushi	on 17" x 18"			1	
Max User Weight	136 kg	Overall Depth	455 mm	Overall Width	430 mm
JION1719 Jay Ion Cushi	on 17" x 19"				
Max User Weight	136 kg	Overall Depth	480 mm	Overall Width	430 mm
				•	
JION1720 Jay Ion Cush					
Max User Weight	136 kg	Overall Depth	505 mm	Overall Width	430 mm
JION1721 Jay Ion Cushi	on 17" x 21"				
Max User Weight	136 kg	Overall Depth	535 mm	Overall Width	430 mm
1001200					
JION1722 Jay Ion Cush			500		100
Max User Weight	136 kg	Overall Depth	560 mm	Overall Width	430 mm
JION1723 Jay Ion Cush	ion 17" x 23"				
Max User Weight	136 kg	Overall Depth	585 mm	Overall Width	430 mm
	10/ 14/				
JION1814 Jay Ion Cushi			055		
Max User Weight	136 kg	Overall Depth	355 mm	Overall Width	455 mm

1011015		
JION1815 Jay Ion Cushio		
Max User Weight	136 kg	Overall Depth
JION1817 Jay Ion Cushio	n 18" x 17"	
Max User Weight	136 kg	Overall Depth
JION1819 Jay Ion Cushio	n 18" x 19"	
Max User Weight	136 kg	Overall Depth
JION1821 Jay Ion Cushio	n 18" x 21"	
Max User Weight	136 kg	Overall Depth
JION1822 Jay Ion Cushic	on 18" x 22"	
Max User Weight	136 kg	Overall Depth
JION1616 Jay Ion Cushio	n 16" x 16"	_
Max User Weight	136 kg	Overall Depth
JION1816 Jay Ion Cushio	n 18" x 16"	
Max User Weight	136 kg	Overall Depth
JION1818 Jay Ion Cushio	n 18" x 18"	
Max User Weight	136 kg	Overall Depth
JION1820 Jay Ion Cushic	on 18" x 20"	
Max User Weight	136 kg	Overall Depth
JION2020 Jay Ion Cushi	on 20" x 20"	
Max User Weight	136 kg	Overall Depth
JION2022 Jay Ion Cushi	on 20" x 22"	
Max User Weight	136 kg	Overall Depth
Jay Lite Cushion		
	Weighin The foar	B Lite Wheelchair Cushion is ext g a mere 0.8 kg, it's the perfect n base is mildly contoured to en



 It provides ventilation and airflow beneau • The optiwell, cutouts beneath the ischia The cover features a wicking material ar

To dissipate heat and moisture, keeping
Available in a variety of sizes 16"-20" x

JL1516 Jay Lite Standard Cushion 15" x 16"					
Max User Weight	250 kg	Overall Depth	405 mm		
Thickness	85 mm				

JL1517 Jay Lite Standard Cushion 15" x 17"

Max User Weight	250 kg	Overall Depth	430 mm	Overall Width	380 mm	Product Weight	0.77 kg
Thickness	85 mm						

558 mm

143

t complement to an ultra lightwu encourage proper positioning of eath the client ial tuberosities and a reticulated and a layer of spacer fabric ig the top surface and the client	rosities and a reticulated foam insert distribute weight over a greater area to					
Overall Width 380 mm Product Weight 0.77 kg						
	300 11111		U.// KY			

508 mm

Overall Width

380 mm	Overall Width	455 mm
430 mm	Overall Width	455 mm
· · · · · · · · · · · · · · · · · · ·		
480 mm	Overall Width	455 mm
535 mm	Overall Width	455 mm
560 mm	Overall Width	455 mm
406 mm	Overall Width	406 mm
406 mm	Overall Width	457 mm
457 mm	Overall Width	457 mm
508 mm	Overall Width	457 mm
508 mm	Overall Width	508 mm



JL1518 Jay Lite Standard Cushion 15" x 18"

Max User Weight	250 kg	Overall Depth	455 mm	Overall Width	380 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1615 Jay Lite Standard Cushion 16" x 15"

Max User Weight	250 kg	Overall Depth	380 mm	Overall Width	405 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1617 Jay Lite Standard Cushion 16" x 17"

Max User Weight	250 kg	Overall Depth	430 mm	Overall Width	405 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1618 Jay Lite Standard Cushion 16" x 18"

Max User Weight	250 kg	Overall Depth	455 mm	Overall Width	405 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1620 Jay Lite Standard Cushion 16" x 20"

Max User Weight	250 kg	Overall Depth	505 mm	Overall Width	405 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1414 Jay Lite Standard Cushion 14" x 14"

Max User Weight	250 kg	Overall Depth	355 mm	Overall Width	355 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1416 Jay Lite Standard Cushion 14" x 16"

Max User Weight	250 kg	Overall Depth	405 mm	Overall Width	355 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1418 Jay Lite Standard Cushion 14" x 18"

Max User Weight	250 kg	Overall Depth	455 mm	Overall Width	355 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1515 Jay Lite Standard Cushion 15" x 15"

Max User Weight	250 kg	Overall Depth	380 mm	Overall Width	380 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1715 Jay Lite Standard Cushion 17" x 15"

Max User Weight	250 kg	Overall Depth	380 mm	Overall Width	430 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1716 Jay Lite Standard Cushion 17" x 16"

Max User Weight	250 kg	Overall Depth	405 mm	Overall Width	430 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1717 Jay Lite Standard Cushion 17" x 17"

Max User Weight	250 kg	Overall Depth	430 mm	Overall Width	430 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1718 Jay Lite Standard Cushion 17" x 18"

Max User Weight	250 kg	Overall Depth	455 mm	Overall Width	430 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1817 Jay Lite Standard Cushion 18" x 17"

Max User Weight	250 kg	Overall Depth	430 mm	Overall Width	455 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1916 Jay Lite Standard Cushion 19" x 16"

Max User Weight	250 kg	Overall Depth	405 mm	Overall Width	480 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1917 Jay Lite Standard Cushion 19" x 17"

Max User Weight	250 kg	Overall Depth	430 mm	Overall Width	480 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1918 Jay Lite Standard Cushion 19" x 18"

,							
Max User Weight	250 kg	Overall Depth	455 mm	Overall Width	480 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1920 Jay Lite Standard Cushion 19" x 20"

Max User Weight	250 kg	Overall Depth	505 mm	Overall Width	480 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL2016 Jay Lite Standard Cushion 20" x 16"

· · · · · · · · · · · · · · · · · · ·							
Max User Weight	250 kg	Overall Depth	405 mm	Overall Width	505 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL2218 Jay Lite Standard Cushion 22" x 18"

Thickness 85 mm	Max User Weight	250 kg	Overall Depth	455 mm
	Thickness	85 mm		

JL2220 Jay Lite Standard Cushion 22" x 20"

Max User Weight	250 kg	Overall Depth	505 mm	Overall Width	560 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL2222 Jay Lite Standard Cushion 22" x 22"

Max User Weight	250 kg	Overall Depth	560 mm	Overall Width	560 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL2418 Jay Lite Standard Cushion 24" x 18"

Max User Weight	250 kg	Overall Depth	455 mm	Overall Width	610 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL2420 Jay Lite Standard Cushion 24" x 20"

Max User Weight	250 kg	Overall Depth	505 mm	Overall Width	610 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL2422 Jay Lite Standard Cushion 24" x 22"

Max User Weight	250 kg	Overall Depth	560 mm
Thickness	85 mm		

JL2424 Jay Lite Standard Cushion 24" x 24"

-							
Max User Weight	250 kg	Overall Depth	610 mm	Overall Width	610 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1616 Jay Lite Standard Cushion 16" x 16"

Max User Weight	250 kg	Overall Depth	405 mm	Overall Width	405 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL1816 Jay Lite Standard Cushion 18" x 16"

Max User Weight	250 kg	Overall Depth	405 mm	Overall Width	455 mm	Product Weight	0.77 kg
Thickness	85 mm						

Pressure Cushions

144



Overall Width	560 mm	Product Weight	0.77 kg

Overall Width	610 mm	Product Weight	0.77 kg

JL1818 Jay Lite Standard Cushion 18" x 18"

Max User We	eight	250 kg	Overall Depth	455 mm	Overall Width	455 mm	Product Weight	0.77 kg
Thickness		85 mm						

JL1820 Jay Lite Standard Cushion 18" x 20"

Max User Weight	250 kg	Overall Depth	505 mm	Overall Width	455 mm	Product Weight	0.77 kg
Thickness	85 mm						

JL2018 Jay Lite Standard Cushion 20" x 18"

Max User Weight	250 kg	Overall Depth	455 mm	Overall Width	505 mm	Product Weight	0.77 ka
Max user weight	200 Ky		400 11111		505 11111		U.// Ky
Thickness	85 mm						

JL2020 Jay Lite Standard Cushion 20" x 20"

•							
Max User Weight	250 kg	Overall Depth	505 mm	Overall Width	505 mm	Product Weight	0.77 kg
Thickness	85 mm						

Jay Soft Combi P Cushion



- A two part seating system comprising a contoured foam seat and a separate contoured lumbar support backrest
 The cushion is moulded from two different kinds of foam, an inner firm foam core for support and an outer softer layer for
- comfort

 Suitable for postural support for people who have a low risk of developing pressure sores
- Available in mild or soft versions and a variety of sizes

Jscombil418 Jay Soft Combi P Cushion 14" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	355 mm	Product Weight	1 kg
Thickness	75 mm						

Jscombi15.516 Jay Soft Combi P Cushion 15.5" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	390 mm	Product Weight	1 kg
Thickness	75 mm						

Jscombi15.518 Jay Soft Combi P Cushion 15.5" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	390 mm	Product Weight	1 kg
Thickness	75 mm						

Jscombil518 Jay Soft Combi P Cushion 15" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	380 mm	Product Weight	1 kg
Thickness	75 mm						

Jscombi1617 Jay Soft Combi P Cushion 16" x 17"

Max User Weight	113 kg	Overall Depth	430 mm	Overall Width	405 mm	Product Weight	1 kg
Thickness	75 mm						

Jscombil620 Jay Soft Combi P Cushion 16" x 20"

Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	405 mm	Product Weight	1 kg
Thickness	75 mm						

Jscombil716 Jay Soft Combi P Cushion 17" x 16"

Ma	ax User Weight	113 kg	Overall Depth	405 mm	Overall Width	430 mm	Product Weight	1 kg
Thi	ickness	75 mm						

Jscombi1717	Jay Soft Combi P Cushion 17" x 17"
-------------	------------------------------------

Max User Weight	113 kg	Overall Depth	430 mm	Overall Width	430 mm	Product Weight	1 kg		
Thickness	75 mm								

Jscombi1718 Jay Soft Combi P Cushion 17" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	430 mm	Product Weight	1 kg
Thickness	75 mm						

Jscombi1816 Jay Soft Combi P Cushion 18" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	455 mm	Product Weight	1 kg
Thickness	75 mm						

Jscombi1817 Jay Soft Combi P Cushion 18" x 17"

Max User Weight	113 kg	Overall Depth	430 mm	Overall Width	455 mm	Product Weight	1 kg
Thickness	75 mm						

Jscombi1818 Jay Soft Combi P Cushion 18" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	455 mm	Product Weight	1 kg	
Thickness	75 mm							

Jscombi1820 Jay Soft Combi P Cushion 18" x 20"

Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	455 mm	Product Weight	1 kg
Thickness	75 mm						

Jscombi1916 Jay Soft Combi P Cushion 19" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	480 mm	Product Weight	1 kg
Thickness	75 mm						

Jscombil917 Jay Soft Combi P Cushion 19" x 17"

Max User Weight	113 kg	Overall Depth	430 mm	Overall Width	480 mm	Product Weight	1 kg
Thickness	75 mm						

Jscombi1918 Jay Soft Combi P Cushion 19" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	480 mm	Product Weight	1 kg
Thickness	75 mm						

Jscombi1920 Jay Soft Combi P Cushion 19" x 20"

Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	480 mm	Product Weight	1 kg
Thickness	75 mm						

Jscombi2016 Jay Soft Combi P Cushion 20" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	505 mm	Product Weight	1 kg
Thickness	75 mm						

Jscombi2018 Jay Soft Combi P Cushion 20" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	505 mm	Product Weight	1 kg
Thickness	75 mm						

Jscombi2020 Jay Soft Combi P Cushion 20" x 20"

Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	505 mm	Product Weight	1 kg
Thickness	75 mm						

Jscombi2218 Jay Soft Combi P Cushion 22" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	560 mm	Product Weight	1 kg
Thickness	75 mm						



Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	560 mm	Product Weight	1
Thickness	75 mm						
scombi2222 Ja	y Soft Con	1bi P Cushion 22	2" x 22"				
Max User Weight	113 kg	Overall Depth	560 mm	Overall Width	560 mm	Product Weight	i
Thickness	75 mm						
scombi2418 Jay	y Soft Com	bi P Cushion 24	" x 18"				
Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	610 mm	Product Weight	
Thickness	75 mm						
scombi2420 Ja Max User Weight Thickness	113 kg	1bi P Cushion 24 Overall Depth	505 mm	Overall Width	610 mm	Product Weight	
Max User Weight Thickness	113 kg 75 mm	Overall Depth	505 mm	Overall Width	610 mm	Product Weight	
Max User Weight Thickness SCOMDi2422 Ja	113 kg 75 mm y Soft Con	Overall Depth	505 mm				
Max User Weight Thickness SCOMDi2422 Ja Max User Weight	113 kg 75 mm y Soft Con 113 kg	Overall Depth	505 mm	Overall Width Overall Width	610 mm	Product Weight Product Weight	
Max User Weight Thickness SCOMDi2422 Ja	113 kg 75 mm y Soft Con	Overall Depth	505 mm				
Max User Weight Thickness SCOMDi2422 Ja Max User Weight Thickness	113 kg 75 mm y Soft Con 113 kg 75 mm	Overall Depth	505 mm 505 mm 560 mm 560 mm				
Max User Weight Thickness SCOMDi2422 Ja Max User Weight Thickness	113 kg 75 mm y Soft Con 113 kg 75 mm	Overall Depth	505 mm 505 mm 560 mm 560 mm				
Max User Weight Thickness SCOMDi2422 Ja Max User Weight Thickness SCOMDi2424 Ja	113 kg 75 mm y Soft Con 113 kg 75 mm 9 Soft Con 113 kg 75 mm	Overall Depth	505 mm	Overall Width	610 mm	Product Weight	
Max User Weight Thickness SCOMDi2422 Ja Max User Weight Thickness SCOMDi2424 Ja Max User Weight Thickness	113 kg 75 mm y Soft Con 113 kg 75 mm y Soft Con 113 kg 75 mm	Overall Depth	505 mm 505 mm 560 mm 560 mm 610 mm 610 mm	Overall Width	610 mm	Product Weight	
Max User Weight Thickness SCOMDi2422 Ja Max User Weight Thickness SCOMDi2424 Ja Max User Weight Thickness	113 kg 75 mm y Soft Con 113 kg 75 mm y Soft Con 113 kg 75 mm	Overall Depth	505 mm 505 mm 560 mm 560 mm 610 mm 610 mm	Overall Width	610 mm	Product Weight	

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	355 mm	Product Weight	1 kg
Thickness	75 mm						

Jay Zip Cushion



- Dual layered contoured foam base
- Reticulated foam in the outer cover promotes air movement
 Soft and stretchable outer cover
- Provides comfort and support for the user

IZIPO808 Jay Zi Max User Weight	75 kg	Overall Depth	205 mm	Overall Width	205 mm	Paediatric ?	у
ZIPO809 Jay Zij	p Cushion &	3" x 9"					
Max User Weight	75 kg	Overall Depth	230 mm	Overall Width	205 mm	Paediatric ?	у
ZIPO810 Jay Zip	p Cushion 8	" x 10"					
IZIPO810 Jay Zip Max User Weight	p Cushion 8	" x 10" Overall Depth	255 mm	Overall Width	205 mm	Paediatric ?	у
	75 kg	Overall Depth	255 mm	Overall Width	205 mm	Paediatric ?	у

JZIP0812 Jay Zip C	Cushion 8	3" x 12"	
Max User Weight	75 kg	Overall Depth	305 mm
JZIPO813 Jay Zip C	Cushion 8	3" x 13"	
Max User Weight	75 kg	Overall Depth	330 mm
JZIP0814 Jay Zip C			255 mm
Max User Weight	75 kg	Overall Depth	355 mm
JZIP0908 Jay Zip (Cushion	9" x 8"	
Max User Weight	75 kg	Overall Depth	205 mm
JZIP0909 Jay Zip (Cushion	9″ x 9″	
Max User Weight	75 kg	Overall Depth	230 mm
JZIPO910 Jav Zip C	Cushion 9	9" x 10"	
Max User Weight	75 kg	Overall Depth	255 mm
JZIPO911 Jay Zip C	uchion 0	" v 11"	
Max User Weight	75 kg	Overall Depth	280 mm
JZIP0912 Jay Zip C			
Max User Weight	75 kg	Overall Depth	305 mm
JZIPO913 Jay Zip C	Cushion 9	9" x 13"	
Max User Weight	75 kg	Overall Depth	330 mm
JZIP0914 Jay Zip C	Cushion 9)" x 14"	
Max User Weight	75 kg	Overall Depth	355 mm
JZIPO915 Jav Zip C	cushion 9)" v 15"	
Max User Weight	75 kg	Overall Depth	380 mm
17101000			
JZIP1008 Jay Zip C Max User Weight		O" X 8" Overall Depth	205 mm
max user weight	73 Kg	overan beptil	203 1111
JZIP1009 Jay Zip C	Cushion 1	0" x 9"	
Max User Weight	75 kg	Overall Depth	230 mm
JZIP1010 Jay Zip C	ushion 1	0" x 10"	
Max User Weight	75 kg	Overall Depth	255 mm
JZIP1011 Jay Zip Ci	ushion 10)" x 11"	
Max User Weight	75 kg	Overall Depth	280 mm
17101010			
JZIP1012 Jay Zip C Max User Weight	75 kg	U" X 12" Overall Depth	305 mm
		oronan bepun	300 1111
JZIP1013 Jay Zip C		D" x 13"	
Max User Weight	75 kg	Overall Depth	330 mm
JZIP1014 Jay Zip C	ushion 1	0" x 14"	

Country Care

Overall Width	205 mm	Paediatric ?	у
Overall Width	205 mm	Paediatric ?	у
Overall Width	205 mm	Paediatric ?	у
Overall Width	230 mm	Paediatric ?	у
		1	
Overall Width	230 mm	Paediatric ?	у
Overall Width	230 mm	Paediatric ?	у
		· · · · · · · · · · · · · · · · · · ·	
Overall Width	230 mm	Paediatric ?	у
			,
Overall Width	230 mm	Paediatric ?	у
Overall Width	230 mm	Paediatric ?	у
Overall Width	230 mm	Paediatric ?	у
Overall Width	230 mm	Paediatric ?	у
		1	
Overall Width	255 mm	Paediatric ?	у
		1	
Overall Width	255 mm	Paediatric ?	у
		1	
Overall Width	255 mm	Paediatric ?	у
Overall Width	255 mm	Paediatric ?	у
Overall Width	255 mm	Paediatric ?	у
Overall Width	255 mm	Paediatric ?	у
Overall Width	255 mm	Paediatric ?	у

JZIP1015 Jay Zip Cu	ushion 10)" x 15"					
Max User Weight	75 kg	Overall Depth	380 mm	Overall Width	255 mm	Paediatric ?	у
JZIP1016 Jay Zip Co	ushion 10)" x 16"					
Max User Weight	75 kg	Overall Depth	405 mm	Overall Width	255 mm	Paediatric ?	у
JZIP1108 Jay Zip Ci	ushion 11	" x 8"					
Max User Weight	75 kg	Overall Depth	205 mm	Overall Width	280 mm	Paediatric ?	у
17/01100 100 700 00		" O"		·		·	
JZIP1109 Jay Zip Cu Max User Weight	75 kg	Overall Depth	230 mm	Overall Width	280 mm	Paediatric ?	y
-							,
JZIP1110 Jay Zip Cu							
Max User Weight	75 kg	Overall Depth	255 mm	Overall Width	280 mm	Paediatric ?	у
JZIP1111 Jay Zip Cu	shion 11"	x 11"					
Max User Weight	75 kg	Overall Depth	280 mm	Overall Width	280 mm	Paediatric ?	у
JZIP1112 Jay Zip Cu	Ishion 11	" x 12"					
Max User Weight	75 kg	Overall Depth	305 mm	Overall Width	280 mm	Paediatric ?	у
						1	
JZIP1113 Jay Zip Cu							
Max User Weight	75 kg	Overall Depth	330 mm	Overall Width	280 mm	Paediatric ?	у
JZIP1114 Jay Zip Cu	shion 11	" x 14"					
Max User Weight	75 kg	Overall Depth	355 mm	Overall Width	280 mm	Paediatric ?	у
JZIP1115 Jay Zip Cu	shion 11'	" x 15"					
Max User Weight	75 kg	Overall Depth	380 mm	Overall Width	280 mm	Paediatric ?	у
JZIP1116 Jay Zip Cu	shion 11	" x 16"					
Max User Weight	75 kg	Overall Depth	405 mm	Overall Width	280 mm	Paediatric ?	у
JZIP1117 Jay Zip Cu	' Ishion 11'	' x 17"					
Max User Weight	75 kg	Overall Depth	430 mm	Overall Width	280 mm	Paediatric ?	у
17ID1200 104 75 0	uphice 1	2" x 0"					
JZIP1209 Jay Zip C Max User Weight		2" X 9" Overall Depth	230 mm	Overall Width	305 mm	Paediatric ?	у
	-						,
JZIP1210 Jay Zip Cu							
Max User Weight	75 kg	Overall Depth	255 mm	Overall Width	305 mm	Paediatric ?	у
JZIP1211 Jay Zip Cu	shion 12	" x 11"					
Max User Weight	75 kg	Overall Depth	280 mm	Overall Width	305 mm	Paediatric ?	у
JZIP1212 Jay Zip Cu	ushion 12	2" x 12"					
Max User Weight	75 kg	Overall Depth	305 mm	Overall Width	305 mm	Paediatric ?	у
JZIP1213 Jay Zip Cu	ushion 12	2" x 13"					
Max User Weight	75 kg	Overall Depth	330 mm	Overall Width	305 mm	Paediatric ?	у
JZIP1214 Jay Zip Cu	ushion 12	2" x 14"					
Max User Weight	75 kg	Overall Depth	355 mm	Overall Width	305 mm	Paediatric ?	у

JZIP1215 Jay Zip	Cushion 1	2" x 15"	
Max User Weight	75 kg	Overall Depth	380 mm
JZIP1216 Jay Zip	Cushion 1	2" x 16"	
Max User Weight	75 kg	Overall Depth	405 mm
JZIP1217 Jav Zip	Cushion 1	0" v 17"	
			430 mm
Max User Weight	75 kg	Overall Depth	430 11111
JZIP1218 Jay Zip	Cushion 1	2" x 18"	
Max User Weight	75 kg	Overall Depth	455 mm
IZIP1310 Jay Zip	Cushion 1	3" v 10"	
Max User Weight	75 kg	Overall Depth	255 mm
JZIP1311 Jay Zip (Cushion 1	3″ x 11″	
Max User Weight	75 kg	Overall Depth	280 mm
IZIP1312 Jay Zip	Cushion 1	3" x 12"	
Max User Weight	75 kg	Overall Depth	305 mm
		I	
JZIP1313 Jay Zip	Cushion 1	3" x 13"	
Max User Weight	75 kg	Overall Depth	330 mm
IZIP1314 Jay Zip	Cushion 1	3" x 14"	
Max User Weight	75 kg	Overall Depth	355 mm
	Cushion 1		
Max User Weight	75 kg	Overall Depth	380 mm
JZIP1316 Jay Zip	Cushion 1	3" x 16"	
Max User Weight	75 kg	Overall Depth	405 mm
JZIP1317 Jay Zip	Cushion 1	0" x 17"	
Max User Weight	75 kg	Overall Depth	430 mm
max over meight	70 Kg	overall bepair	100 1111
JZIP1318 Jay Zip	Cushion 1	3" x 18"	
Max User Weight	75 kg	Overall Depth	455 mm
JZIP1411 Jay Zip (Cushion 14	4" x 11"	
Max User Weight	75 kg	Overall Depth	280 mm
IZIP1412 Jay Zip	Cushion 1	4" x 12"	
Max User Weight	75 kg	Overall Depth	305 mm
JZIP1413 Jav Zip	Cushion 1	4" x 13"	
Max User Weight	75 kg	Overall Depth	330 mm
	Cushion 1		
Max User Weight	75 kg	Overall Depth	355 mm
JZIP1415 Jay Zip	Cushion 1	4" x 15"	
Max User Weight	75 kg	Overall Depth	380 mm

Country Care G R O U P National Distribution, Local Service

Overall Width	305 mm	Paediatric ?	у
Overall Width	305 mm	Paediatric ?	у
Overall Width	305 mm	Paediatric ?	у
Overall Width	305 mm	Paediatric ?	у
		1	
Overall Width	330 mm	Paediatric ?	у
Overall Width	330 mm	Paediatric ?	у
Overall Width	330 mm	Paediatric ?	у
Overall Width	330 mm	Paediatric ?	у
		1	
Overall Width	330 mm	Paediatric ?	у
		1	
Overall Width	330 mm	Paediatric ?	у
Overall Width	330 mm	Paediatric ?	у
		1	
Overall Width	330 mm	Paediatric ?	у
Overall Width	330 mm	Paediatric ?	у
		1	
Overall Width	355 mm	Paediatric ?	у
Overall Width	355 mm	Paediatric ?	у
Overall Width	355 mm	Paediatric ?	у
Overall Width	355 mm	Paediatric ?	у
Overall Width	355 mm	Paediatric ?	у

JZIP1416 Jay Zip Cu	shion 14	" x 16"					
Max User Weight	75 kg	Overall Depth	405 mm	Overall Width	355 mm	Paediatric ?	у
JZIP1417 Jay Zip Cu	shion 14	" x 17"					
Max User Weight	75 kg	Overall Depth	430 mm	Overall Width	355 mm	Paediatric ?	у
JZIP1418 Jay Zip Cu	shion 14	" x 18"					
Max User Weight	75 kg	Overall Depth	455 mm	Overall Width	355 mm	Paediatric ?	у
JZIP1512 Jay Zip Cu	shion 15	" x 12"					
Max User Weight	75 kg	Overall Depth	305 mm	Overall Width	380 mm	Paediatric ?	у
JZIP1513 Jay Zip Cu	shion 15	" x 13"					
Max User Weight	75 kg	Overall Depth	330 mm	Overall Width	380 mm	Paediatric ?	у
JZIP1514 Jay Zip Cu	shion 15	" x 14"					
Max User Weight	75 kg	Overall Depth	355 mm	Overall Width	380 mm	Paediatric ?	у
JZIP1515 Jay Zip Cu	shion 15	" x 15"					
Max User Weight	75 kg	Overall Depth	380 mm	Overall Width	380 mm	Paediatric ?	у
JZIP1516 Jay Zip Cu	shion 15	" x 16"					
Max User Weight	75 kg	Overall Depth	405 mm	Overall Width	380 mm	Paediatric ?	у
JZIP1517 Jay Zip Cu	shion 15	" x 17"					
Max User Weight	75 kg	Overall Depth	430 mm	Overall Width	380 mm	Paediatric ?	у
JZIP1518 Jay Zip Cu	shion 15	" x 18"					
Max User Weight	75 kg	Overall Depth	455 mm	Overall Width	380 mm	Paediatric ?	у
JZIP1613 Jay Zip Cu	shion 16	" x 13"					
Max User Weight	75 kg	Overall Depth	330 mm	Overall Width	405 mm	Paediatric ?	у
JZIP1614 Jay Zip Cu	shion 16	" x 14"					
Max User Weight	75 kg	Overall Depth	355 mm	Overall Width	405 mm	Paediatric ?	у
JZIP1615 Jay Zip Cu	shion 16	" x 15"					
Max User Weight	75 kg	Overall Depth	380 mm	Overall Width	405 mm	Paediatric ?	у
JZIP1616 Jay Zip Cu	shion 16	" x 16"					
Max User Weight	75 kg	Overall Depth	405 mm	Overall Width	405 mm	Paediatric ?	у
JZIP1617 Jay Zip Cu	shion 16	" x 17"					
Max User Weight	75 kg	Overall Depth	430 mm	Overall Width	405 mm	Paediatric ?	у
JZIP1618 Jay Zip Cu	shion 16	" x 18"					

Pressure Cushions Non Powered	7.02 Pressure Cushions Air
Skil-Care Multi Purpose Cushion Air Equalisation	
 Features high/low profile position equalisatio Variable action of individual air cells Vow surface tension reduces friction and skin Pressure adjustable An easy to use air pump is included with all n Low shear II cloth cover is fluid proof for incore Anti-slip bottom for added user safety 	shear nodels
SK751540B Skil-Care Multi Purpose Cushion Air Equalisation	High 16" x 16"
Overall Depth 405 mm Overall Width	405 mm Thickness 100 mm
SK751545B Skil-Care Multi Purpose Cushion Air Equalisation	Low 16" x 16"
Overall Depth 405 mm Overall Width	405 mm Thickness 50 mm
SK751550B Skil-Care Multi Purpose Cushion Air Equalisation	High 18" x 16"
Overall Depth 405 mm Overall Width	455 mm Thickness 100 mm
CV7E1ECOD Old One Multi Durnana Quakian Air Equalization	Leve 10" v: 10"
SK751560B Skil-Care Multi Purpose Cushion Air Equalisation Overall Depth 405 mm Overall Width	455 mm Thickness 50 mm
SK751562B Skil-Care Multi Purpose Cushion Air Equalisation	High 20" x 16"
Overall Depth 405 mm Overall Width	505 mm Thickness 100 mm
SK751564B Skil-Care Multi Purpose Cushion Air Equalisation	Low 20" x 16"
Overall Depth 405 mm Overall Width	505 mm Thickness 50 mm
Roho Single Valve Low Profile Cushion	
 5.5 cm air cell height to maximise stability of Lightweight, adjustable and portable cushion Protects from the potential risk of tissue/skin Can be adjusted to address the exact body sh Provides easier transfers due to low cell heigh The cushion is ventilated allowing air movement Flame resistant neoprene rubber Low Friction and shear 	breakdown Iape, size and seating needs of the user It
R-1R1112LPC Roho Single Valve LP Cushion 20.25" x 22"	
Overall Depth 560 mm Overall Width 515 mm Product	tt Weight 1.5 kg Thickness 55 mm
R-1R118LPC Roho Single Valve LP Cushion 20.25" x 15"	
	tt Weight 1.5 kg Thickness 55 mm



1.5 kg Thickness 55 mm 1.5 kg Thickness 55 mm

R-1R1210LPC Roho	Single V	alve LP Cushion 22	2″ x 18.5″				
Overall Depth	470 mm	Overall Width	560 mm	Product Weight	1.5 kg	Thickness	55 mm
R-1R1211LPC Roho	Single Va	alve LP Cushion 22	" x 20.25"				
Overall Depth	515 mm	Overall Width	560 mm	Product Weight	1.5 kg	Thickness	55 mm
R-1R1212LPC Rohd) Sinale V	alve LP Cushion 22	2" x 22"				
Overall Depth	560 mm	Overall Width	560 mm	Product Weight	1.5 kg	Thickness	55 mm
R-1R129LPC Roho	Single Va	lve LP Cushion 22'	" x 16.75"				
Overall Depth	425 mm	Overall Width	560 mm	Product Weight	1.5 kg	Thickness	55 mm
R-1R128LPC Roho	Single Va	lve LP Cushion 22'	" x 15"				
Overall Depth	380 mm	Overall Width	560 mm	Product Weight	1.5 kg	Thickness	55 mm
R-1R29LPC Roho S	Single Val	ve LP Cushion 4.75	5" x 16.75"				
Overall Depth	425 mm	Overall Width	120 mm	Product Weight	1.5 kg	Thickness	55 mm
	-	ve LP Cushion 7.5"					
Overall Depth	190 mm	Overall Width	190 mm	Product Weight	1.5 kg	Thickness	55 mm
R-1R54LPC Roho S	Single Val	ve LP Cushion 9.25	5″ x 7.5″				
Overall Depth	190 mm	Overall Width	235 mm	Product Weight	1.5 kg	Thickness	55 mm
R-1R58LPC Roho S	Single Val	ve LP Cushion 9.25	5″ x 15″				
Overall Depth	380 mm	Overall Width	235 mm	Product Weight	1.5 kg	Thickness	55 mm
		ve LP Cushion 9.25			151		
Overall Depth	425 mm	Overall Width	235 mm	Product Weight	1.5 kg	Thickness	55 mm
R-1R65LP Roho Si	ngle Valv	e LP Cushion 11.5"	x 9.25″				
Overall Depth	235 mm	Overall Width	290 mm	Product Weight	1.5 kg	Thickness	55 mm
R-1R66LP Roho Si	ngle Valv	e LP Cushion 11.5"	x 11.5″				
Overall Depth	290 mm	Overall Width	290 mm	Product Weight	1.5 kg	Thickness	55 mm
	lingle Vel		" v 10 05"				
R-1R67LPC Roho S	335 mm	ve LP Cushion 11.5'	290 mm	Product Weight	1.5 kg	Thickness	55 mm
	333 1/111		230 1/111	Floudet weight	1.0 Ky	THICKIESS	53 11111
R-1R68LPC Roho S	Single Val	ve LP Cushion 11.5	" x 15"				
Overall Depth	380 mm	Overall Width	290 mm	Product Weight	1.5 kg	Thickness	55 mm
R-1R74LPC Roho S	Single Val	ve LP Cushion 13.2	5″ x 7.5″				
Overall Depth	190 mm	Overall Width	335 mm	Product Weight	1.5 kg	Thickness	55 mm
R-1R69LPC Robo	Single Val	ve LP Cushion 11 5	" x 16.75"				
R-1R69LPC Roho S	Single Val	ve LP Cushion 11.5	" x 16.75"	Product Weight	1.5 kg	Thickness	55 mm
Overall Depth	425 mm	Overall Width	290 mm	Product Weight	1.5 kg	Thickness	55 mm
Overall Depth R-1R76LPC Roho S	425 mm	Overall Width	290 mm				
Overall Depth	425 mm	Overall Width	290 mm	Product Weight Product Weight	1.5 kg	Thickness Thickness	55 mm 55 mm
Overall Depth R-1R76LPC Roho S Overall Depth	425 mm Single Val	Overall Width	290 mm 25" x 11.5" 335 mm	Product Weight			

R-1R78LPC Roho	Single Va	ve LP Cushion 13	8.25″ x 15″
Overall Depth	380 mm	Overall Width	335 mm
		ve LP Cushion 13	
Overall Depth	425 mm	Overall Width	335 mm
R-1R810LPC Roh	o Single Va	alve LP Cushion 1	5" x 18.5"
Overall Depth	470 mm	Overall Width	380 mm
R-1R84LPC Robo		lve LP Cushion 15	" v 7 6"
Overall Depth	190 mm	Overall Width	380 mm
R-1R85LPC Roho	Single Va	lve LP Cushion 15	5" x 9.25"
Overall Depth	235 mm	Overall Width	380 mm
R-1R86LPC Roho	single Va	lve LP Cushion 15	5" x 11.5"
Overall Depth	290 mm	Overall Width	380 mm
-			
R-1R87LPC Roho	Single Va	ve LP Cushion 15	5" x 13.25"
Overall Depth	335 mm	Overall Width	380 mm
R-1R88LPC Roho) Single Va	lve LP Cushion 15	5" x 15"
Overall Depth	380 mm	Overall Width	380 mm
		Ive LP Cushion 15	
Overall Depth	425 mm	Overall Width	380 mm
R-1R910LPC Roh	o Single Va	alve LP Cushion 1	6.75" x 18.5"
Overall Depth	470 mm	Overall Width	425 mm
R-1R911LPC Roh	n Sinalo Va	lve LP Cushion 1	6 75" v 20 25
Overall Depth	515 mm	Overall Width	425 mm
R-1R1010LPC Rol	ho Single \	alve LP Cushion	18" x 18"
Overall Depth	455 mm	Overall Width	455 mm
R-1R97LPC Roho	Single Va	ve LP Cushion 16	6.75" x 13.25"
Overall Depth	335 mm	Overall Width	425 mm
		Ive LP Cushion 16	
Overall Depth	380 mm	Overall Width	425 mm
R-1R109LPC Roh	o Single Va	alve LP Cushion 1	8″ x 16″
Overall Depth	405 mm	Overall Width	455 mm
D 1011101 D0			
		alve LP Cushion	
Overall Depth	455 mm	Overall Width	505 mm
R-1R1011LPC Rot	no Single V	alve LP Cushion	18.5″ x 20.25
Overall Depth	515 mm	Overall Width	470 mm
R-1R108LPC Roh	o Sinalo V	alve LP Cushion 1	8 5" ¥ 15"
K INTOOLEG ROIT			
Overall Depth	380 mm	Overall Width	470 mm

Country Care
GROUP
National Distribution, Local Service

Product Weight	1.5 kg	Thickness	55 mm
Product Weight	1.5 kg	Thickness	55 mm
Product Weight	1.5 kg	Thickness	55 mm
Product Weight	1.5 kg	Thickness	55 mm
· · · · · · · · · · · · · · · · · · ·	illo kg		
Product Weight	15 km	Thiskness	55 mm
Product Weight	1.5 kg	Thickness	55 mm
Product Weight	1.5 kg	Thickness	55 mm
Product Weight	1.5 kg	Thickness	55 mm
Product Weight	1.5 kg	Thickness	55 mm
Product Weight	1.5 kg	Thickness	55 mm
Product Weight	1.5 kg	Thickness	55 mm
,			
Product Weight	1.5 kg	Thickness	55 mm
Product Weight	1.5 kg	Thickness	55 mm
Product Weight	1.5 kg	Thickness	55 mm
Product Weight	1.5 kg	Thickness	55 mm
Product Weight	1.5 kg	Thickness	55 mm
Product Weight	1.5 kg	Thickness	55 mm
	-		
Product Weight	1.5 kg	Thickness	55 mm
Product Weight	1.5 kg	Thickness	55 mm



R-1R99LPC Roho Single Valve LP Cushion 16" x 16"							
Overall Depth	405 mm	Overall Width	405 mm	Product Weight	1.5 kg	Thickness	55 mm
R-1R1111LPC Roho Single Valve LP Cushion 20.25" x 20.25"							
Overall Depth	515 mm	Overall Width	515 mm	Product Weight	1.5 kg	Thickness	55 mm

Roho Dual Valve High Profile Cushion

 An adjustable, low deformation, counter pressure environment
• An inflatable cushion of neoprene (latex free) rubber "cells" with two inflation valves

- Independently adjust two separate compartments of the cushion to suit the needs of the user
- Is light and conforms to the shape of the user



 Interconnected air cells which allows for maximum therapeutic performance and enhanced sitting times Suitable for people at medium risk of developing pressure sores

	515 mm	Overall Width	470 mm	Product Weight	1.5 kg	Thickness	100 mm
Overall Depth	515 11111		4/0 11111	Floudet weight	1.0 Kg	THICKIESS	100 11111
-2R105C Roho D	ual Valve	HP Cushion 18.5" x	9.25"				
Overall Depth	235 mm	Overall Width	470 mm	Product Weight	1.5 kg	Thickness	100 mm
- 2R107C Roho Di	ual Valve	HP Cushion 18.5" x	13.25″				
Overall Depth	335 mm	Overall Width	470 mm	Product Weight	1.5 kg	Thickness	100 mm
2-2R1111C Roho D	ual Valve	HP Cushion 20.25"	x 20.25"				
Overall Depth	515 mm	Overall Width	515 mm	Product Weight	1.5 kg	Thickness	100 mm
2-2R66C Roho Du	al Valve I	HP Cushion 11.5" x 1	1.5″				
Overall Depth	290 mm	Overall Width	290 mm	Product Weight	1.5 kg	Thickness	100 mm
-2R67C Roho Du	al Valve F	IP Cushion 11.5" x 1	3.25″				
Overall Depth	335 mm	Overall Width	290 mm	Product Weight	1.5 kg	Thickness	100 mm
-2R68C Roho Du	al Valve I	HP Cushion 11.5" x 1	5″				
Overall Depth	380 mm	Overall Width	290 mm	Product Weight	1.5 kg	Thickness	100 mm
-2R74C Roho Du	al Valve H	IP Cushion 13.25" x	7.5″				
	al Valve H	IP Cushion 13.25" x	7.5 "	Product Weight	1.5 kg	Thickness	100 mm
Overall Depth	190 mm		335 mm	Product Weight	1.5 kg	Thickness	100 mm
Overall Depth	190 mm	Overall Width	335 mm	Product Weight Product Weight	1.5 kg 1.5 kg	Thickness Thickness	100 mm 100 mm
Overall Depth R-2R77C Roho Du Overall Depth	190 mm al Valve H 335 mm	Overall Width	335 mm (13.25" 335 mm				
Overall Depth R-2R77C Roho Du Overall Depth R-2R78C Roho Du	190 mm al Valve H 335 mm	Overall Width	335 mm (13.25" 335 mm				
Overall Depth 2-2R77C Roho Du Overall Depth 2-2R78C Roho Du Overall Depth	190 mm al Valve H 335 mm al Valve H 335 mm al Valve H 380 mm	Overall Width HP Cushion 13.25" x Overall Width HP Cushion 13.25" x	335 mm (13.25" 335 mm (15" 335 mm	Product Weight	1.5 kg	Thickness	100 mm
Overall Depth R-2R77C Roho Du Overall Depth R-2R78C Roho Du Overall Depth R-2R810C Roho D	190 mm al Valve H 335 mm al Valve H 335 mm al Valve H 380 mm	Overall Width HP Cushion 13.25" x Overall Width HP Cushion 13.25" x Overall Width	335 mm (13.25" 335 mm (15" 335 mm	Product Weight	1.5 kg	Thickness	100 mm
Overall Depth R-2R77C Roho Du Overall Depth R-2R78C Roho Du Overall Depth R-2R810C Roho Di Overall Depth	190 mm al Valve H 335 mm al Valve H 380 mm ual Valve 470 mm	Overall Width HP Cushion 13.25" x Overall Width HP Cushion 13.25" x Overall Width HP Cushion 15" x 18	335 mm (13.25" 335 mm (15" 335 mm 8.5" 380 mm	Product Weight Product Weight	1.5 kg	Thickness Thickness	100 mm

R-2R86C	Roho Dual	Valve H	P Cushion 15" x 11.5	7	
Overall Depth		290 mm	Overall Width	380 mm	I
R-2R87C	Roho Dual	Valve H	P Cushion 15" x 13.2	5"	
Overall Depth		335 mm	Overall Width	380 mm	
R-2R88C	Roho Dual	Valve H	P Cushion 15" x 15"		
Overall Depth		380 mm	Overall Width	380 mm	
R-2R89C	Roho Dual	Valve H	P Cushion 15" x 16.7	5"	
Overall Depth		425 mm	Overall Width	380 mm	
R-2R910C	Roho Dua	al Valve I	HP Cushion 16.75" x	18.5″	
Overall Depth		470 mm	Overall Width	425 mm	
R-2R98C	Roho Dual	Valve H	P Cushion 16.75" x 1	5″	
Overall Depth		380 mm	Overall Width	425 mm	
R-2R99C	Roho Dual	Valve H	P Cushion 16" x 16"		
Overall Depth		405 mm	Overall Width	405 mm	
R-2R10100	C Roho Du	al Valve	HP Cushion 18" x 18	"	
Overall Depth		455 mm	Overall Width	455 mm	
R-2R109C	Roho Dua	al Valve I	HP Cushion 18" x 16"		
Overall Depth		405 mm	Overall Width	455 mm	
R-2R11100	Roho Dua	al Valve	HP Cushion 20" x 18	n	
Overall Depth		455 mm	Overall Width	505 mm	

Roho Dual Valve Low Profile Cushion



 $\ensuremath{\,\bullet\,}$ An adjustable cushion with low deformation and a counter pressure environment • An inflatable cushion of neoprene (latex free) rubber "cells" • Two inflation valves which independently adjust two separate compartments of the cushion to suit the needs of the user • Is lightweight, portable and conforms to the shape of the user • The interconnected air cells allow for maximum therapeutic performance and enhanced sitting times

• The lower cell height can make transferring easier

Suitable for people at medium risk of developing pressure sores

R-2R1011LPC Roho I	Dual Valv	ve LP Cushion 18.5" >	(20.25"	
Overall Depth	515 mm	Overall Width	470 mm	P
R-2R1111LPC Roho D	ual Valv	e LP Cushion 20.25"	x 20.25	"
Overall Depth	515 mm	Overall Width	515 mm	P
R-2R66LPC Roho Du	ual Valve	LP Cushion 11.5" x 1	1.5"	
Overall Depth	290 mm	Overall Width	290 mm	P
R-2R67LPC Roho Du	ial Valve	LP Cushion 11.5" x 1	3.25″	
Overall Depth	335 mm	Overall Width	290 mm	P



Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
	kg		100 1111
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 Kg	THICKIESS	100 11111
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product weight	1.5 ку	INICKNESS	
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
in a second s			

roduct Weight	50 mm
roduct Weight	50 mm
and and Michaels	50
roduct Weight	50 mm
roduct Weight	50 mm
roduct Weight roduct Weight	50 mm

Overall Depth	335 mm	Overall Width	335 mm	Product Weight	1.5 kg	Thickness	50 mm
-2R78LPC Roho	Dual Valvo	e LP Cushion 13.25	" x 15"				
Overall Depth	380 mm	Overall Width	335 mm	Product Weight	1.5 kg	Thickness	50 mm
R-2R810LPC Rohe	o Dual Valv	ve LP Cushion 15" x	18.5″				
Overall Depth	470 mm	Overall Width	380 mm	Product Weight	1.5 kg	Thickness	50 mm
R-2R87LPC Roho	Dual Valvo	e LP Cushion 15" x	13.25″				
Overall Depth	335 mm	Overall Width	380 mm	Product Weight	1.5 kg	Thickness	50 mm
R-2R88LPC Roho	Dual Valv	e LP Cushion 15" x	15″				
Overall Depth	380 mm	Overall Width	380 mm	Product Weight	1.5 kg	Thickness	50 mm
R-2R89LPC Roho	Dual Valv	e LP Cushion 15" x	16.75″				
Overall Depth	425 mm	Overall Width	380 mm	Product Weight	1.5 kg	Thickness	50 mm
		Overall Width		Product Weight	1.5 kg	Thickness	50 mm
				Product Weight Product Weight	1.5 kg	Thickness Thickness	50 mm
R-2R910LPC Rohd	Dual Valv	ve LP Cushion 16.75	5" x 18.5" 425 mm				
R-2R910LPC Rohd	Dual Valv	/e LP Cushion 16.75 Overall Width	5" x 18.5" 425 mm				
R-2R910LPC Rohd Overall Depth R-2R98LPC Roho Overall Depth	Dual Valv 470 mm Dual Valv 380 mm	ve LP Cushion 16.75 Overall Width e LP Cushion 16.75	5" x 18.5" 425 mm * x 15" 425 mm	Product Weight	1.5 kg	Thickness	50 mm
R-2R910LPC Rohd Overall Depth R-2R98LPC Roho Overall Depth	Dual Valv 470 mm Dual Valv 380 mm	/e LP Cushion 16.75 Overall Width e LP Cushion 16.75 Overall Width	5" x 18.5" 425 mm * x 15" 425 mm	Product Weight	1.5 kg	Thickness	50 mm
R-2R910LPC Roho Overall Depth R-2R98LPC Roho Overall Depth R-2R99LPC Roho Overall Depth	Dual Valv 470 mm Dual Valv 380 mm Dual Valv 405 mm	/e LP Cushion 16.75 Overall Width e LP Cushion 16.75 Overall Width e LP Cushion 16" x	5" x 18.5" 425 mm x 15" 425 mm 425 mm 16" 405 mm	Product Weight Product Weight	1.5 kg	Thickness Thickness	50 mm
R-2R910LPC Roho Overall Depth R-2R98LPC Roho Overall Depth R-2R99LPC Roho Overall Depth	Dual Valv 470 mm Dual Valv 380 mm Dual Valv 405 mm	ve LP Cushion 16.75 Overall Width E LP Cushion 16.75 Overall Width E LP Cushion 16" x Overall Width	5" x 18.5" 425 mm x 15" 425 mm 425 mm 16" 405 mm	Product Weight Product Weight	1.5 kg	Thickness Thickness	50 mm
R-2R910LPC Roho Overall Depth R-2R98LPC Roho Overall Depth R-2R99LPC Roho Overall Depth R-2R1010LPC Roh Overall Depth	Dual Valv 470 mm Dual Valv 380 mm Dual Valv 405 mm to Dual Va	/e LP Cushion 16.75 Overall Width e LP Cushion 16.75 Overall Width e LP Cushion 16" x Overall Width	5" x 18.5" 425 mm x 15" 425 mm 425 mm 405 mm x 18" 455 mm	Product Weight Product Weight Product Weight	1.5 kg	Thickness Thickness Thickness	50 mm 50 mm 50 mm
R-2R910LPC Roho Overall Depth R-2R98LPC Roho Overall Depth R-2R99LPC Roho Overall Depth R-2R1010LPC Roh Overall Depth	Dual Valv 470 mm Dual Valv 380 mm Dual Valv 405 mm to Dual Va	/e LP Cushion 16.75 Overall Width e LP Cushion 16.75 Overall Width e LP Cushion 16" x Overall Width Ive LP Cushion 18" Overall Width	5" x 18.5" 425 mm x 15" 425 mm 425 mm 405 mm x 18" 455 mm	Product Weight Product Weight Product Weight	1.5 kg	Thickness Thickness Thickness	50 mm 50 mm 50 mm
R-2R910LPC Roho Overall Depth R-2R98LPC Roho Overall Depth R-2R99LPC Roho Overall Depth R-2R1010LPC Roho Overall Depth R-2R109LPC Roho Overall Depth	Dual Valv 470 mm Dual Valv 380 mm Dual Valv 405 mm 455 mm Dual Valv 455 mm	/e LP Cushion 16.75 Overall Width e LP Cushion 16.75 Overall Width e LP Cushion 16" x Overall Width Ive LP Cushion 18" Overall Width /e LP Cushion 18" x	5" x 18.5" 425 mm 425 mm 425 mm 425 mm 405 mm x 18" 455 mm 455 mm	Product Weight Product Weight Product Weight Product Weight	1.5 kg 1.5 kg 1.5 kg 1.5 kg 1.5 kg	Thickness Thickness Thickness Thickness	50 mm 50 mm 50 mm 50 mm

Roho Single Valve High Profile Cushion



- Appropriate for high risk clients for pressure sores or skin breakdown
- Interconnected air cells maximise stability of the user
 Adjustable, low deformation and counter-pressure environment
- Prevents and assists in healing of pressure sores
- Conforms to body shape, imerses and envelopes
 Maintains blood flow while reducing risk of skin breakdown
- Comfortable for longer sitting tolerance
 Lightweight and portable, low friction and shear

R-1R109C Roho Sir	igle Valve	HP Cushion 18" x 16	5″						
Overall Depth	405 mm	Overall Width	455 mm	Product Weight	1.5 kg	Thickness	100 mm		
Overall Juppin 405 mm Overall Width 455 mm Product Weight 1.5 kg Linckness 100 mm R-1R1110C Roho Single Valve HP Cushion 20" x 18"									
R-IKITIUL RONO SI	ngle Valv	e HP Cushion 20" x I	8″						

D_1D00C	Dobo Sing		HP Cushion 16" x 16	"
Overall Depth	Rono Sing	405 mm	Overall Width	405 mm
	Daha Oʻaal			
R-1R79C	Roho Sing		HP Cushion 13.25" x	
Overall Depth		425 mm	Overall Width	335 mm
R-1R810C	Roho Sing	gle Valve	HP Cushion 15" x 18	3.5″
Overall Depth		470 mm	Overall Width	380 mm
R-1R84C	Roho Sing	le Valve	HP Cushion 15" x 7.5	5″
Overall Depth		190 mm	Overall Width	380 mm
R-1R85C	Roho Sing	le Valve	HP Cushion 15" x 9.3	25″
Overall Depth		235 mm	Overall Width	380 mm
R-1R86C	Roho Sina	le Valve	HP Cushion 15" x 11.	5"
Overall Depth		290 mm	Overall Width	380 mm
			I	
R-1R87C	Roho Sing		HP Cushion 15" x 13	
Overall Depth		335 mm	Overall Width	380 mm
R-1R88C	Roho Sing	le Valve	HP Cushion 15" x 15	"
Overall Depth		380 mm	Overall Width	380 mm
R-1R89C	Roho Sing	le Valve	HP Cushion 15" x 16	.75″
Overall Depth		425 mm	Overall Width	380 mm
R-1R910C	Roho Sino	ole Valve	e HP Cushion 16.75"	x 18.5″
Overall Depth		470 mm	Overall Width	425 mm
R-1R911C	Poho Sinc	ulo Valvo	HP Cushion 16.75" >	(20 25"
Overall Depth	Kono Sing	515 mm	Overall Width	425 mm
R-1R912C	Roho Sini	nle Valve	HP Cushion 16.75"	
Overall Depth		560 mm	Overall Width	425 mm
			I	
R-1R94C	Roho Sing		HP Cushion 16.75" x	7.5″
Overall Depth		190 mm	Overall Width	425 mm
R-1R95C	Roho Sing	le Valve	HP Cushion 16.75" x	9.25″
Overall Depth		235 mm	Overall Width	425 mm
R-1R96C	Roho Sing	le Valve	HP Cushion 16.75" x	11.5"
Overall Depth	_	290 mm	Overall Width	425 mm
R-1R97C	Roho Sina	le Valve	HP Cushion 16.75" x	13.25"
Overall Depth		335 mm	Overall Width	425 mm
Overall Depth				425 mm
Overall Depth			Overall Width HP Cushion 16.75" x Overall Width	425 mm
Overall Depth R-1R98C Overall Depth	Roho Sing	le Valve	HP Cushion 16.75" x Overall Width	425 mm 1 5″ 425 mm
Overall Depth R-1R98C Overall Depth	Roho Sing	le Valve	HP Cushion 16.75" x	425 mm 1 5″ 425 mm



Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm
Product Weight	1.5 kg	Thickness	100 mm

R-1R1012C Roho	Single Valve HP Cus	shion 18.5" x 22"				
Overall Depth	560 mm Overall Widt	h 470 mm	Product Weight	1.5 kg	Thickness	100 mm
-1 R1013 Roho S	Single Valve HP Cust	nion 18.5" x 23.75"				
Overall Depth	605 mm Overall Widtl		Product Weight	1.5 kg	Thickness	100 mm
R-1R104C Roho	Single Valve HP Cus					
Overall Depth	190 mm Overall Widtl	n 470 mm	Product Weight	1.5 kg	Thickness	100 mm
R-1R107C Roho	Single Valve HP Cus	hion 18.5" x 13.25"				
Overall Depth	335 mm Overall Width	h 470 mm	Product Weight	1.5 kg	Thickness	100 mm
R-1R108C Roho	Single Valve HP Cus	hion 18.5" x 15"				
Overall Depth	380 mm Overall Width		Product Weight	1.5 kg	Thickness	100 mm
			•			
	Single Valve HP Cus					
Overall Depth	515 mm Overall Width	n 515 mm	Product Weight	1.5 kg	Thickness	100 mm
R-1R1112C Roho	Single Valve HP Cus	hion 20.25" x 22"				
Overall Depth	560 mm Overall Widtl		Product Weight	1.5 kg	Thickness	100 mm
	Single Valve HP Cus		1			
Overall Depth	605 mm Overall Width	n 515 mm	Product Weight	1.5 kg	Thickness	100 mm
R-1R118C Roho S	Single Valve HP Cust	1ion 20.25" x 15"				
Overall Depth	380 mm Overall Widtl	n 515 mm	Product Weight	1.5 kg	Thickness	100 mm
			_			
	Single Valve HP Cust					
Overall Depth	425 mm Overall Width	n 515 mm	Product Weight	1.5 kg	Thickness	100 mm
R-1R1210C Roho	Single Valve HP Cus	shion 22" x 18.5"				
Overall Depth	470 mm Overall Widtl	n 560 mm	Product Weight	1.5 kg	Thickness	100 mm
0 1012110 Daha		hian 2011 y 20 251				
	Single Valve HP Cus		Droduct Weight	1.5 kg	Thickness	100 mm
Overall Depth		1 560 IIIII	Product Weight	1.5 Ky	mickness	
R-1R1212C Roho	Single Valve HP Cus	shion 22" x 22"				
Overall Depth	560 mm Overall Widtl	n 560 mm	Product Weight	1.5 kg	Thickness	100 mm
101200 Dates	Circula Value UD Qual	hian 00% x 10 75%				
	A25 mm Overall Widtl		Product Weight	1.5.40	Thickness	100 mm
Overall Depth		h 560 mm		1.5 kg	Thickness	100 mm
R-1R1310C Roho	Single Valve HP Cus	shion 23.75" x 18.5'	1			
Overall Depth	470 mm Overall Width	h 600 mm	Product Weight	1.5 kg	Thickness	100 mm
0 1012110 Date		hinn 00 75" 00 0	- "			
	Single Valve HP Cus			15 kg	Thickness	100 mm
Overall Depth	olo IIIII Uverali Widti		Product Weight	1.5 kg	THICKNESS	
R-1R1312C Roho	Single Valve HP Cus	shion 23.75" x 22"				
Overall Depth	560 mm Overall Widtl	h 600 mm	Product Weight	1.5 kg	Thickness	100 mm
0-1012120 Daka			E "			
	Single Valve HP Cus			15 40	Thickness	100
Overall Depth	605 mm Overall Widtl	h 600 mm	Product Weight	1.5 kg	Thickness	100 mm

	jie valve	HP Cushion 23.75"	X 16.75				
Overall Depth	425 mm	Overall Width	600 mm	Product Weight	1.5 kg	Thickness	100 mm
R-1R48C Roho Singl	e Valve	HP Cushion 7.5" x 15	5″				
Overall Depth	380 mm	Overall Width	190 mm	Product Weight	1.5 kg	Thickness	100 mm
R-1R64C Roho Singl	le Valve	HP Cushion 11.5" x 7	.5"				
Overall Depth	190 mm	Overall Width	290 mm	Product Weight	1.5 kg	Thickness	100 mm
R-1R65C Roho Singl	e Valve	HP Cushion 11.5" x 9	.25″				
Overall Depth	235 mm	Overall Width	290 mm	Product Weight	1.5 kg	Thickness	100 mm
R-1R66C Roho Singl	e Valve	HP Cushion 11.5" x 1	1.5″				
Overall Depth	290 mm	Overall Width	290 mm	Product Weight	1.5 kg	Thickness	100 mm
R-1R67C Roho Singl	e Valve	HP Cushion 11.5" x 1	3.25″				
Overall Depth	335 mm	Overall Width	290 mm	Product Weight	1.5 kg	Thickness	100 mm
R-1R68C Roho Singl	e Valve	HP Cushion 11.5" x 1	5″				
Overall Depth	380 mm	Overall Width	290 mm	Product Weight	1.5 kg	Thickness	100 mm
R-1R69C Roho Singl	e Valve	HP Cushion 11.5" x 1	6.75″				
Overall Depth	425 mm	Overall Width	290 mm	Product Weight	1.5 kg	Thickness	100 mm
R-1R710C Roho Sing	jle Valve	HP Cushion 13.25"	x 18.5″				
Overall Depth	470 mm	Overall Width	335 mm	Product Weight	1.5 kg	Thickness	100 mm
R-1R74C Roho Singl	e Valve I	HP Cushion 13.25" x	7.5″				
Overall Depth	190 mm	Overall Width	335 mm	Product Weight	1.5 kg	Thickness	100 mm
		Overall Width HP Cushion 13.25" x	335 mm	Product Weight	1.5 kg	Thickness	100 mm
			335 mm	Product Weight Product Weight	1.5 kg 1.5 kg	Thickness Thickness	100 mm 100 mm
R-1R75C Roho Singl Overall Depth	e Valve	HP Cushion 13.25" x Overall Width	335 mm 9.25″ 335 mm				1
R-1R75C Roho Singl Overall Depth	e Valve	HP Cushion 13.25" x Overall Width	335 mm 9.25″ 335 mm				1
R-1R75C Roho SingI Overall Depth R-1R76C Roho SingI Overall Depth	e Valve 235 mm e Valve 290 mm	HP Cushion 13.25" x Overall Width HP Cushion 13.25" x Overall Width	335 mm 9.25" 335 mm 11.5" 335 mm	Product Weight	1.5 kg	Thickness	100 mm
R-1R75C Roho SingI Overall Depth R-1R76C Roho SingI Overall Depth	e Valve 235 mm e Valve 290 mm	HP Cushion 13.25" x Overall Width HP Cushion 13.25" x Overall Width	335 mm 9.25" 335 mm 11.5" 335 mm	Product Weight	1.5 kg	Thickness	100 mm
R-1R75C Roho SingI Overall Depth R-1R76C Roho SingI Overall Depth R-1R77C Roho SingI Overall Depth	e Valve 235 mm e Valve 290 mm e Valve 335 mm	HP Cushion 13.25" x Overall Width HP Cushion 13.25" x Overall Width HP Cushion 13.25" x	335 mm 9.25" 335 mm 11.5" 335 mm 13.25" 335 mm	Product Weight Product Weight	1.5 kg	Thickness	100 mm

Roho Contour Select Cushion



Provides the skin protection you expect from RohoWith the positioning and stability needed for increased function • The contoured design positions the pelvis back in the wheelchair and centers the user comfortably in the middle of the

- cushion
- The ISOFLO memory control locks air into each of the four quadrants

These features provide the stability that will allow the user to maximise function while minimising tissue deformation

Pressure Cushions

R-CS1010C Roho Conto	our Select Cushion 18.5" x 1	18.5″	
Overall Depth	470 mm	Overall Width	470 mm
R-CS1011C Roho Conto	ur Select Cushion 18.5" x 2	20.25"	
Overall Depth	515 mm	Overall Width	470 mm
R-CS109C Roho Conto	ur Select Cushion 18.5" x 10	6.75"	
Overall Depth	425 mm	Overall Width	470 mm
R-CS1110C Roho Conto	ur Select Cushion 20.25" x	(18.5"	
Overall Depth	470 mm	Overall Width	515 mm
R-CS1111C Roho Contou	ur Select Cushion 20.25" x	20.25"	
Overall Depth	515 mm	Overall Width	515 mm
R-CS119C Roho Contou	r Select Cushion 20.25" x	16.75"	
Overall Depth	425 mm	Overall Width	515 mm
R-CS810C Roho Contou	ur Select Cushion 15" x 18.	5″	
Overall Depth	470 mm	Overall Width	380 mm
R-CS88C Roho Contou	r Select Cushion 15" x 15"		
Overall Depth	380 mm	Overall Width	380 mm
R-CS89C Roho Contou	r Select Cushion 15" x 16.7	5″	
Overall Depth	425 mm	Overall Width	380 mm
R-CS910C Roho Conto	ur Select Cushion 16.75" x	18.5″	
Overall Depth	470 mm	Overall Width	425 mm
R-CS911C Roho Contou	r Select Cushion 16.75" x 2	20.25"	
Overall Depth	515 mm	Overall Width	425 mm
R-CS97C Roho Contou	r Select Cushion 16.75" x 13	3.25″	
Overall Depth	335 mm	Overall Width	425 mm
R-CS98C Roho Contou	r Select Cushion 16.75" x 1	5"	
Overall Depth	380 mm	Overall Width	425 mm
R-CS99C Roho Contou	r Select Cushion 16.75" x 1	6.75″	

Roho Enhancer Cushion



- The Enhancer Dry Floatation Cushion offers you a unique option in wheelchair seating
 The user can benefit from positive midline channeling of the legs along with increased lateral stability
- Inhibits sacral sitting while enhancing tissue viability
 Unique contoured shape is the result of flexible air cells of varying heights
- Onique concorrect snape is the result of nextile an cens of varying neights
 The two manifold system allows you to adjust the inside cells independently form the outside cells
 A cushion that shapes to the user's body
 Provides superior pressure distribution and minimal shear forces

- Offers the added benefit of progressive position

R-ENH910C Roho Enhanc Overall Depth R R-ENH1010C Roho Enhanc Overall Depth R R-ENH109C Roho Enhanc Overall Depth R R-ENH109C Roho Enhanc Overall Depth R R-ENH1110C Roho Enhanc Overall Depth R R-ENH1111C Roho Enhanc Overall Depth R R-ENH1011C Roho Enhanc Overall Depth R R-ENH1012C Roho Enhanc Overall Depth R R-ENH1012C Roho Enhanc Overall Depth R R-ENH1012C Roho Enhanc Overall Depth R R-ENH119C Roho Enhanc Overall Depth R R-ENH1212C Roho Enhanc Overall Depth R R-ENH66C Roho Enhanc Overall Depth R R-ENH67C Roho Enhanc Overall Depth R R-ENH67C Roho Enhanc Overall Depth R	er Cushior 425 mm er Cushior 470 mm er Cushior 470 mm er Cushior 470 mm er Cushior 515 mm er Cushion 515 mm cer Cushion 560 mm 425 mm er Cushion 425 mm 2000 mm 1000	Overall Width I Overall Width Overall Width Overall Width Overall Width I Overall Width Overall Width I I Overall Width I
R-ENH910C Roho Enhanc Overall Depth Image: Constraint of the second	er Cushior 470 mm cer Cushio 470 mm er Cushior 470 mm er Cushior 515 mm cer Cushion 515 mm cer Cushion 560 mm er Cushion 560 mm cer Cushion 560 mm cer Cushion 560 mm cer Cushion 560 mm	I I I I I I I I I I I I I I I I I I I
Overall Depth A R-ENH1010C Roho Enhanc Overall Depth A R-ENH109C Roho Enhanc Overall Depth A R-ENH1110C Roho Enhanc Overall Depth A R-ENH1110C Roho Enhanc Overall Depth A R-ENH1111C Roho Enhanc Overall Depth A R-ENH1011C Roho Enhanc Overall Depth A R-ENH1012C Roho Enhanc Overall Depth A R-ENH1012C Roho Enhanc Overall Depth A R-ENH1012C Roho Enhanc Overall Depth A R-ENH1212C Roho Enhanc Overall Depth A R-ENH66C Roho Enhanc Overall Depth A R-ENH66C Roho Enhanc Overall Depth A R-ENH67C Roho Enhanc Overall Depth A R-ENH77C Roho Enhanc Overall Depth A R-ENH78C Roho Enhanc	470 mm cer Cushio 470 mm er Cushior 425 mm er Cushior 470 mm er Cushior 515 mm cer Cushior 560 mm cer Cushion 560 mm cer Cushion 560 mm cer Cushion 560 mm r Cushion 290 mm	Overall Width Overall Width Overall Width III.5" x 18.5" Overall Width III.5" x 16.75" Overall Width III.5" x 16.75" Overall Width III.5" x 18.5" Overall Width III.5" x 20.25" Overall Width III.5" x 20.25" Overall Width III.5" x 16.75" Overall Width II.5" x 11.5" Overall Width II.5" x 13.25" Overall Width
R-ENH1010C Roho Enhance Overall Depth R R-ENH109C Roho Enhance Overall Depth R R-ENH1110C Roho Enhance Overall Depth R R-ENH1111C Roho Enhance Overall Depth R R-ENH1011C Roho Enhance Overall Depth R R-ENH1011C Roho Enhance Overall Depth R R-ENH1012C Roho Enhance Overall Depth R R-ENH1012C Roho Enhance Overall Depth R R-ENH1012C Roho Enhance Overall Depth R R-ENH1212C Roho Enhance Overall Depth R R-ENH66C Roho Enhance Overall Depth R R-ENH67C Roho Enhance Overall Depth R R-ENH67C Roho Enhance Overall Depth R R-ENH77C Roho Enhance Overall Depth R R-ENH77C Roho Enhance Overall Depth	cer Cushio 470 mm er Cushior 425 mm er Cushior 470 mm er Cushion 515 mm cer Cushion 560 mm er Cushion 560 mm er Cushion 560 mm er Cushion 560 mm	Image: Product with the second sec
Overall Depth A R-ENH109C Roho Enhanc Overall Depth A R-ENH1110C Roho Enhanc Overall Depth A R-ENH1111C Roho Enhanc Overall Depth A R-ENH1111C Roho Enhanc Overall Depth A R-ENH1011C Roho Enhanc Overall Depth A R-ENH1012C Roho Enhanc Overall Depth A R-ENH1012C Roho Enhanc Overall Depth A R-ENH1012C Roho Enhanc Overall Depth A R-ENH1212C Roho Enhanc Overall Depth A R-ENH66C Roho Enhanc Overall Depth A R-ENH67C Roho Enhanc Overall Depth A R-ENH67C Roho Enhanc Overall Depth A R-ENH77C Roho Enhanc Overall Depth A R-ENH77C Roho Enhanc Overall Depth A R-ENH78C Roho Enhanc	470 mm er Cushior 425 mm er Cushior 470 mm er Cushior 515 mm er Cushion 515 mm cer Cushion 560 mm er Cushion 560 mm cer Cushion 560 mm cr Cushion 560 mm er Cushion 560 mm cr	Overall Width I1.5" x 11.5" Overall Width 11.5" x 13.25" Overall Width
R-ENH109C Roho Enhanc Overall Depth R R-ENH1110C Roho Enhanc Overall Depth R R-ENH1111C Roho Enhanc Overall Depth R R-ENH1011C Roho Enhanc Overall Depth R R-ENH1011C Roho Enhanc Overall Depth R R-ENH1012C Roho Enhanc Overall Depth R R-ENH1012C Roho Enhanc Overall Depth R R-ENH1212C Roho Enhanc Overall Depth R R-ENH66C Roho Enhanc Overall Depth R R-ENH66C Roho Enhanc Overall Depth R R-ENH67C Roho Enhanc Overall Depth R R-ENH67C Roho Enhanc Overall Depth R R-ENH77C Roho Enhanc Overall Depth R R-ENH77C Roho Enhanc Overall Depth R R-ENH78C Roho Enhanc	er Cushior 425 mm er Cushior 470 mm er Cushior 515 mm er Cushion 560 mm er Cushion 560 mm cer Cushion 560 mm er Cushion 290 mm	IB.5" x 16.75" Overall Width IIII 20.25" x 18.5" Overall Width IIIIIIIIIIIIIIIIIIIIIIIIIIIIIIIIIIII
Overall Depth R-ENH1110C Roho Enhance Overall Depth Image: Comparison of the second	er Cushion	Overall Width n 20.25" x 18.5" Overall Width n 20.25" x 20.25" Overall Width n 18.5" x 20.25" Overall Width n 18.5" x 20.25" Overall Width 11.5" x 11.5" overall Width 11.5" x 13.25" overall Width
R-ENH1110C Roho Enhance Overall Depth Image: Control of the second s	er Cushion er Cushion er Cushion er Cushion fis mm er Cushion fis	n 20.25" x 18.5" Overall Width n 20.25" x 20.25" Overall Width n 18.5" x 20.25" Overall Width n 18.5" x 22" Overall Width n 20.25" x 16.75" Overall Width n 22" x 22" Overall Width 11.5" x 11.5" Overall Width 11.5" x 13.25" Overall Width
Overall Depth	470 mm er Cushior 515 mm ser Cushio 515 mm cer Cushio 560 mm er Cushion 560 mm r Cushion 290 mm	Overall Width Overall Width Overall Width Ill Overall Width Ill.5" x 11.5" Overall Width Ill.5" x 13.25" Overall Width
R-ENH1111C Roho Enhanc Overall Depth R-ENH1011C Roho Enhanc Overall Depth R-ENH1012C Roho Enhance Overall Depth R-ENH119C Roho Enhance Overall Depth R-ENH1212C Roho Enhance Overall Depth R-ENH66C Roho Enhance Overall Depth R-ENH67C Roho Enhance Overall Depth R-ENH67C Roho Enhance Overall Depth R-ENH77C Roho Enhance Overall Depth	er Cushior 515 mm er Cushion 515 mm cer Cushion 660 mm 425 mm cer Cushion 560 mm r Cushion 290 mm r Cushion	Image: 20.25" x 20.25" Overall Width Image: 18.5" x 20.25" Overall Width Image: 18.5" x 22" Overall Width Image: 18.5" x 11.5" Overall Width Image: 18.5" x 13.25" Overall Width
Overall Depth Image: Control of the	er Cushion cer Cushion	Overall Width n 18.5" x 20.25" Overall Width on 18.5" x 22" Overall Width 11.5" x 11.5" overall Width 11.5" x 13.25" overall Width
R-ENH1011C Roho Enhance Overall Depth R-ENH1012C Roho Enhance Overall Depth R-ENH119C Roho Enhance Overall Depth R-ENH1212C Roho Enhance Overall Depth R-ENH66C Roho Enhance Overall Depth R-ENH67C Roho Enhance Overall Depth R-ENH77C Roho Enhance Overall Depth R-ENH77C Roho Enhance Overall Depth	er Cushion	n 18.5" x 20.25" Overall Width Overall Width
Overall Depth Image: Control of the second state of the seco	515 mm cer Cushion 560 mm 425 mm cer Cushion 560 mm r Cushion 290 mm r Cushion	Overall Width 11.5" x 11.5" Overall Width 0verall Width 0verall Width 0verall Width
Overall Depth Image: Control of the second state of the seco	515 mm cer Cushion 560 mm 425 mm cer Cushion 560 mm r Cushion 290 mm r Cushion	Overall Width 11.5" x 11.5" Overall Width 0verall Width 0verall Width
R-ENH1012C Roho Enhance Overall Depth R-ENH119C Roho Enhance Overall Depth R-ENH1212C Roho Enhance Overall Depth R-ENH66C Roho Enhance Overall Depth R-ENH67C Roho Enhance Overall Depth R-ENH77C Roho Enhance Overall Depth R-ENH77C Roho Enhance Overall Depth	cer Cushio 560 mm er Cushion 425 mm cer Cushio 560 mm er Cushion 290 mm r Cushion	on 18.5" x 22" Overall Width a 20.25" x 16.75" Overall Width on 22" x 22" Overall Width 11.5" x 11.5" Overall Width 11.5" x 13.25" Overall Width
Overall Depth Image: Construction of the second	er Cushion 425 mm cer Cushio 560 mm r Cushion 290 mm r Cushion	Overall Width 1 20.25" x 16.75" Overall Width 0n 22" x 22" Overall Width 11.5" x 11.5" Overall Width 11.5" x 13.25" Overall Width
R-ENH119C Roho Enhance Overall Depth R-ENH66C Roho Enhance Overall Depth R-ENH67C Roho Enhance Overall Depth R-ENH67C Roho Enhance Overall Depth R-ENH77C Roho Enhance Overall Depth R-ENH77C Roho Enhance	er Cushion 425 mm cer Cushio 560 mm r Cushion 290 mm r Cushion	20.25" x 16.75" Overall Width on 22" x 22" Overall Width 11.5" x 11.5" Overall Width 11.5" x 13.25" Overall Width
Overall Depth Image: Control Contron Contro Control Control Control Contron Control Cont	425 mm cer Cushio 560 mm r Cushion 290 mm r Cushion	Overall Width on 22" x 22" Overall Width 11.5" x 11.5" Overall Width 11.5" x 13.25" Overall Width
R-ENH1212C Roho Enhance Overall Depth R-ENH66C Roho Enhance Overall Depth R-ENH67C Roho Enhance Overall Depth R-ENH77C Roho Enhance Overall Depth R-ENH77C Roho Enhance	cer Cushio 560 mm r Cushion 290 mm r Cushion	on 22" x 22" Overall Width 11.5" x 11.5" Overall Width 11.5" x 13.25" Overall Width
Overall Depth Image: Construction of the second	r Cushion	Overall Width 11.5" x 11.5" Overall Width 11.5" x 13.25" Overall Width
Overall Depth Image: Construction of the second	r Cushion	Overall Width 11.5" x 11.5" Overall Width 11.5" x 13.25" Overall Width
R-ENH66C Roho Enhance Overall Depth R-ENH67C Roho Enhance Overall Depth R-ENH77C Roho Enhance Overall Depth R-ENH78C Roho Enhance	r Cushion 290 mm r Cushion	11.5" x 11.5" Overall Width 11.5" x 13.25" Overall Width
Overall Depth Image: Constraint of the second s	290 mm	Overall Width 11.5" x 13.25" Overall Width
R-ENH67C Roho Enhance Overall Depth R-ENH77C Roho Enhance Overall Depth R-ENH78C Roho Enhance	r Cushion	11.5" x 13.25" Overall Width
Overall Depth Image: Constraint of the second sec		Overall Width
R-ENH77C Roho Enhance Overall Depth R-ENH78C Roho Enhance	335 mm	
Overall Depth R-ENH78C Roho Enhance		13.25" x 13.25"
Overall Depth		13.25" x 13.25"
R-ENH78C Roho Enhance		
	335 mm	Overall Width
Overall Depth		
	380 mm	Overall Width
R-ENH87C Roho Enhance	r Cushion	15" x 13.25"
Overall Depth	335 mm	Overall Width
R-ENH88C Roho Enhance	r Cushion	15" x 15"
Overall Depth	380 mm	Overall Width
R-ENH89C Roho Enhance	r Cushion	15" x 16.75"
Overall Depth	425 mm	Overall Width
R-ENH911C Roho Enhance		
	er Cushion	16.75" x 20.25"

425 mm	Thickness	108 mm
425 mm	Thickness	108 mm
470 mm	Thickness	108 mm
470 mm	Thickness	108 mm
515 mm	Thickness	108 mm
515 mm	Thickness	108 mm
470 mm	Thickness	108 mm
470 mm	Thickness	108 mm
515 mm	Thickness	108 mm
560 mm	Thickness	108 mm
290 mm	Thickness	108 mm
290 mm	Thickness	108 mm
335 mm	Thickness	108 mm
335 mm	Thickness	108 mm
380 mm	Thickness	108 mm
380 mm	Thickness	108 mm
380 mm	Thickness	108 mm
425 mm	Thickness	108 mm



R-ENH98C Roho Enhancer Cushion 16.75" x 15"

Overall Depth	380 mm	Overall Width		425 mm	Thickness		108 mm
Roho Quadtro S	elect Low F	rofile Cushior	1				
		Ilows for quick and simple Push the knob to the open p ireat for maximum functior itability is enhanced by loc Ainimises side to side or frr ightweight, portable, repai ir is able to circulate arour	, oosition allows and posture king air into pla ont to back mol rable and wash	adjustment then push ace in each of the four o ion able	he knob closed quadrants		nt
R-QS1010LPC Roho	Quadtro Sele	ct Low Profile Cu	shion 18"	x 18″			
Overall Depth	455 mm Overal	l Width 4	55 mm Pro	duct Weight	1.4 kg	Thickness	50 mm
R-QS109LPC Roho	Quadtro Solor	t Low Profile Cus	tion 10"	v 16"			
Overall Depth	-			duct Weight	1.4 kg	Thickness	50 mm
R-QS1110LPC Roho	Quadtro Sele	ct Low Profile Cu	shion 20"	x 18″	_		
Overall Depth	455 mm Overal	l Width 5	05 mm Pro	duct Weight	1.4 kg	Thickness	50 mm
R-QS99LPC Roho Q	uadtro Select	Low Profile Cus	nion 16" x	16"			
Overall Depth				duct Weight	1.4 kg	Thickness	50 mm
-		ct Low Profile Cu					
Overall Depth	515 mm Overal	I Width 4	70 mm Pro	duct Weight	1.4 kg	Thickness	50 mm
R-QS1012LPC Roho	Quadtro Sele	ct Low Profile Cu	shion 18.	5" X 22"			
Overall Depth	560 mm Overal	l Width 4	70 mm Pro	duct Weight	1.4 kg	Thickness	50 mm
-	-	t Low Profile Cus					
Overall Depth	335 mm Overal	I Width 4	70 mm Pro	duct Weight	1.4 kg	Thickness	50 mm
R-QS108LPC Roho	Quadtro Selec	t Low Profile Cus	shion 18.5	" X 15"			
Overall Depth				duct Weight	1.4 kg	Thickness	50 mm
			hier 00 f				
R-QS1111LPC Roho	_			25" X 20.25"	14 40	Thickness	50 mm
overall bepth	515 IIIII Uveral			lauct weight	1.4 kg	Inickness	50 11111
R-QS1112LPC Roho	Quadtro Sele	ct Low Profile Cu	shion 20.	25" X 22"			
Overall Depth	560 mm Overal	l Width 5	15 mm Pro	duct Weight	1.4 kg	Thickness	50 mm
R-QS119LPC Roho ()undtre Color	t Low Drofile Cue	hion 20 0	E" V 10 75"			
K-USIISLAC KOUOI	Juadtro Selec			5" X 16.75"	1.4 kg	Thickness	50 mm
Averall Denth	/25 mm Overel	Width					
Overall Depth	425 mm Overal	l Width 5			1.4 Kg	THICKNESS	50 11111
		ct Low Profile Cu			1.4 Kg	THEAHCSS	30 1111
	Quadtro Sele	ct Low Profile Cu	shion 22'		1.4 kg	Thickness	50 mm
R-QS1210LPC Roho	Quadtro Sele	ct Low Profile Cu	shion 22'	' X 18.5" Iduct Weight			

	no Quadtr	o Select Low Prof	ile Cushion	22" X 22"			
Overall Depth	560 mm	Overall Width	560 mm	Product Weight	1.4 kg	Thickness	50 mm
R-QS1310LPC Ro	oho Quadtr	o Select Low Prof	file Cushion	23.75" X 18.5"			
Overall Depth	470 mm	Overall Width	605 mm	Product Weight	1.4 kg	Thickness	50 mm
R-OS1311LPC Ro	ho Quadtro	o Select Low Prof	ile Cushion	23.75" X 20.25"			
Overall Depth	515 mm	Overall Width	605 mm	Product Weight	1.4 kg	Thickness	50 mm
-	-	Select Low Profile					
Overall Depth	290 mm	Overall Width	290 mm	Product Weight	1.4 kg	Thickness	50 mm
R-QS67LPC Roh	o Quadtro S	Select Low Profile	e Cushion 11	" X 13.25"			
Overall Depth	335 mm	Overall Width	290 mm	Product Weight	1.4 kg	Thickness	50 mm
R-OS88LPC Roh	o Quadtra	Soloot Low Drofil	o Cuchion 14				
Overall Depth	380 mm	Select Low Profile		Product Weight	1.4 kg	Thickness	50 mm
	000 mm	J COLUMNIA COLUMN	000 1111	gitt			
R-QS89LPC Roh	o Quadtro	Select Low Profile	e Cushion 15	5" X 16.75"			
Overall Depth	425 mm	Overall Width	380 mm	Product Weight	1.4 kg	Thickness	50 mm
R-OS9101 PC Rol	ho Quadtro	Select Low Profi	le Cushion 1	6.75" X 18.5"			
Overall Depth	470 mm	Overall Width	425 mm	Product Weight	1.4 kg	Thickness	50 mm
R-QS911LPC Roh	io Quadtro	Select Low Profil	e Cushion 1	6.75" X 20.25"			
Overall Depth	515 mm	Overall Width	425 mm	Product Weight	1.4 kg	Thickness	50 mm
R-QS97LPC Roh	o Quadtro (Select Low Profile	e Cushion 16	6.75" X 13.25"			
Overall Depth	335 mm	Overall Width	425 mm	Product Weight	1.4 kg	Thickness	50 mm
-		Select Low Profile			1.41	7111	50
Overall Depth	o Quadtro 380 mm	Select Low Profile Overall Width	e Cushion 16 425 mm	6.75" X 15" Product Weight	1.4 kg	Thickness	50 mm
Overall Depth	380 mm		425 mm	Product Weight	1.4 kg	Thickness	50 mm
Overall Depth	380 mm	Overall Width	425 mm	Product Weight	1.4 kg 1.4 kg	Thickness Thickness	50 mm
Overall Depth R-QS68LPC Roh Overall Depth	380 mm	Overall Width Select Low Profile Overall Width	425 mm e Cushion 11 290 mm	Product Weight .5" X 15" Product Weight			
Overall Depth R-QS68LPC Roh Overall Depth R-QS77LPC Roh	380 mm o Quadtro 380 mm o Quadtro S	Overall Width Select Low Profile Overall Width Select Low Profile	425 mm e Cushion 11 290 mm e Cushion 13	Product Weight .5" X 15" Product Weight 3.25" X 13.25"	1.4 kg	Thickness	50 mm
Overall Depth R-QS68LPC Roh Overall Depth	380 mm	Overall Width Select Low Profile Overall Width	425 mm e Cushion 11 290 mm	Product Weight .5" X 15" Product Weight			
Overall Depth R-QS68LPC Roh Overall Depth R-QS77LPC Roh Overall Depth	o Quadtro 380 mm 380 mm o Quadtro 335 mm	Overall Width Select Low Profile Overall Width Select Low Profile	425 mm e Cushion 11 290 mm e Cushion 13 335 mm	Product Weight .5" X 15" Product Weight 3.25" X 13.25" Product Weight	1.4 kg	Thickness	50 mm
Overall Depth R-QS68LPC Roh Overall Depth R-QS77LPC Roh Overall Depth	o Quadtro 380 mm 380 mm o Quadtro 335 mm	Overall Width Select Low Profile Overall Width Select Low Profile Overall Width	425 mm e Cushion 11 290 mm e Cushion 13 335 mm	Product Weight .5" X 15" Product Weight 3.25" X 13.25" Product Weight	1.4 kg	Thickness	50 mm
Overall Depth R-QS68LPC Roh Overall Depth R-QS77LPC Roh Overall Depth R-QS78LPC Roh Overall Depth Overall Depth	380 mm o Quadtro 380 mm o Quadtro 335 mm o Quadtro 380 mm	Overall Width Select Low Profile Overall Width Select Low Profile Overall Width Select Low Profile Overall Width	425 mm e Cushion 11 290 mm e Cushion 13 335 mm e Cushion 13 335 mm	Product Weight .5" X 15" Product Weight 3.25" X 13.25" Product Weight 3.25" X 15" Product Weight	1.4 kg	Thickness Thickness	50 mm
Overall Depth R-QS68LPC Roh Overall Depth R-QS77LPC Roh Overall Depth R-QS78LPC Roh Overall Depth Overall Depth	380 mm o Quadtro 380 mm o Quadtro 335 mm o Quadtro 380 mm	Overall Width Select Low Profile Overall Width Select Low Profile Overall Width Select Low Profile	425 mm e Cushion 11 290 mm e Cushion 13 335 mm e Cushion 13 335 mm	Product Weight .5" X 15" Product Weight 3.25" X 13.25" Product Weight 3.25" X 15" Product Weight	1.4 kg	Thickness Thickness	50 mm
Overall Depth R-QS68LPC Roh Overall Depth R-QS77LPC Roh Overall Depth R-QS78LPC Roh Overall Depth R-QS79LPC Roh R-QS79LPC Roh	o Quadtro 380 mm o Quadtro 335 mm o Quadtro 380 mm o Quadtro	Overall Width Select Low Profile Overall Width	425 mm e Cushion 11 290 mm e Cushion 13 335 mm e Cushion 13 335 mm e Cushion 13	Product Weight .5" X 15" Product Weight 3.25" X 13.25" Product Weight 3.25" X 15" Product Weight 3.25" X 16.75"	1.4 kg	Thickness Thickness Thickness	50 mm 50 mm 50 mm
Overall Depth R-QS68LPC Roh Overall Depth R-QS77LPC Roh Overall Depth R-QS78LPC Roh Overall Depth R-QS79LPC Roh Overall Depth R-QS79LPC Roh Overall Depth	380 mm o Quadtro 380 mm o Quadtro 335 mm o Quadtro 380 mm o Quadtro 380 mm o Quadtro 425 mm	Overall Width Select Low Profile Overall Width	425 mm e Cushion 11 290 mm e Cushion 13 335 mm e Cushion 13 335 mm e Cushion 13 335 mm	Product Weight .5" X 15" Product Weight 3.25" X 13.25" Product Weight 3.25" X 15" Product Weight 3.25" X 16.75" Product Weight	1.4 kg	Thickness Thickness Thickness	50 mm 50 mm 50 mm
Overall Depth R-QS68LPC Roh Overall Depth R-QS77LPC Roh Overall Depth R-QS78LPC Roh Overall Depth R-QS79LPC Roh Overall Depth R-QS79LPC Roh Overall Depth	380 mm o Quadtro 380 mm o Quadtro 335 mm o Quadtro 380 mm o Quadtro 380 mm o Quadtro 425 mm	Overall Width Select Low Profile Overall Width	425 mm e Cushion 11 290 mm e Cushion 13 335 mm e Cushion 13 335 mm e Cushion 13 335 mm	Product Weight .5" X 15" Product Weight 3.25" X 13.25" Product Weight 3.25" X 15" Product Weight 3.25" X 16.75" Product Weight	1.4 kg	Thickness Thickness Thickness	50 mm 50 mm 50 mm
Overall Depth R-QS68LPC Roh Overall Depth R-QS77LPC Roh Overall Depth R-QS78LPC Roh Overall Depth R-QS79LPC Roh Overall Depth R-QS810LPC Roh Overall Depth	380 mm o Quadtro 380 mm o Quadtro 335 mm o Quadtro 380 mm o Quadtro 380 mm o Quadtro 380 mm o Quadtro 425 mm ho Quadtro 470 mm	Overall Width Select Low Profile Overall Width	425 mm e Cushion 11 290 mm e Cushion 13 335 mm e Cushion 13 335 mm e Cushion 13 335 mm ile Cushion 13 335 mm ile Cushion 13 338 mm	Product Weight .5" X 15" Product Weight 3.25" X 13.25" Product Weight 3.25" X 15" Product Weight 3.25" X 16.75" Product Weight 5" X 18.5" Product Weight	1.4 kg 1.4 kg 1.4 kg 1.4 kg 1.4 kg	Thickness Thickness Thickness Thickness Thickness	50 mm 50 mm 50 mm 50 mm

Pressure Cushions

Country Care
GROUP
National Distribution, Local Service



Roho Quadtro Select High Profile Cushion



- Allows for quick and simple position selection for every user while seated • Push the knob to the open position allows adjustment then push the knob closed to maintain this position
- Great for maximum function and posture
- Stability is enhanced by locking air into place in each of the four quadrants
 Minimises side to side or front to back motion

- Lightweight, portable, repairable and washable
 Air is able to circulate around the cells, reaching the surface of the skin and providing a cooling environment

R-QS1010C Roho	Quadtro S	elect High Profile	Cushion 18	8" x 18"			
Overall Depth	455 mm	Overall Width	455 mm	Product Weight	1.4 kg	Thickness	105 mm
R-QS1011C Roho	Quadtro So	elect High Profile	Cushion 18	8.5" x 20.25"			
Overall Depth	515 mm	Overall Width	470 mm	Product Weight	1.4 kg	Thickness	105 mm
R-QS1012C Roho	Quadtro S	elect High Profile	Cushion 18	3.5" x 22"			
Overall Depth	560 mm	Overall Width	470 mm	Product Weight	1.4 kg	Thickness	105 mm
R-QS108C Roho	Quadtro Se	elect High Profile	Cushion 18	.5″ x 15″			
Overall Depth	380 mm	Overall Width	470 mm	Product Weight	1.4 kg	Thickness	105 mm
R-QS1111C Roho	Quadtro Se	lect High Profile	Cushion 20).25" x 20.25"			
Overall Depth	515 mm	Overall Width	515 mm	Product Weight	1.4 kg	Thickness	105 mm
R-QS1112C Roho	Quadtro So	elect High Profile	Cushion 20).25" x 22"			
Overall Depth	560 mm	Overall Width	515 mm	Product Weight	1.4 kg	Thickness	105 mm
R-QS1113C Roho	Quadtro S	elect High Profile	Cushion 20	0.25" x 23.75"			
Overall Depth	605 mm	Overall Width	515 mm	Product Weight	1.4 kg	Thickness	105 mm
R-QS118C Roho (Quadtro Se	lect High Profile (Cushion 20	.25" x 15"			
Overall Depth	380 mm	Overall Width	515 mm	Product Weight	1.4 kg	Thickness	105 mm
R-QS119C Roho (Quadtro Se	lect High Profile (Cushion 20	.25" x 16.75"			
Overall Depth	425 mm	Overall Width	515 mm	Product Weight	1.4 kg	Thickness	105 mm
R-QS1210C Roho	Quadtro S	elect High Profile	Cushion 2	2" x 18.5"			
Overall Depth	470 mm	Overall Width	560 mm	Product Weight	1.4 kg	Thickness	105 mm
R-QS1211C Roho	Quadtro So	elect High Profile	Cushion 22	2" x 20.25"			
Overall Depth	515 mm	Overall Width	560 mm	Product Weight	1.4 kg	Thickness	105 mm
R-QS1212C Roho	Quadtro S	elect High Profile	Cushion 2	2" x 22"			
Overall Depth	560 mm	Overall Width	560 mm	Product Weight	1.4 kg	Thickness	105 mm
R-QS129C Roho	Quadtro Se	lect High Profile (Cushion 22	." x 16.75"			
Overall Depth	425 mm	Overall Width	560 mm	Product Weight	1.4 kg	Thickness	105 mm
R-QS1310C Roho	Quadtro S	elect High Profile	Cushion 2	3.75" x 18.5"			
Overall Depth	470 mm	Overall Width	605 mm	Product Weight	1.4 kg	Thickness	105 mm

-	o Quadtro Select High I				
Overall Depth	515 mm Overall Width	605 mm Produc	t Weight 1.4 kg	Thickness	105 mm
R-QS1312C Rol	no Quadtro Select High	Profile Cushion 23.75"	x 22"		
Overall Depth	560 mm Overall Width	605 mm Produc	t Weight 1.4 kg	Thickness	105 mm
D 0014100 D-1			0.55"		
R-QS1410C ROI Overall Depth	10 Quadtro Select High 470 mm Overall Width		8.55" :t Weight 1.4 kg	Thickness	105 mm
overall beput			t weight 1.4 kg	THICKNESS	103 11111
R-QS1411C Roh	o Quadtro Select High I	Profile Cushion 24" x 2).25″		
Overall Depth	515 mm Overall Width	610 mm Produc	t Weight 1.4 kg	Thickness	105 mm
R-0S66C Roho	Quadtro Select High Pr	rofile Cushion 11.5" x 11	.5"		
Overall Depth	290 mm Overall Width		t Weight 1.4 kg	Thickness	105 mm
				I	
-	Quadtro Select High Pr				
Overall Depth	335 mm Overall Width	290 mm Produc	t Weight 1.4 kg	Thickness	105 mm
R-QS911C Roho	o Quadtro Select High P	rofile Cushion 16.75" x	20.25"		
Overall Depth	515 mm Overall Width	425 mm Produc	t Weight 1.4 kg	Thickness	105 mm
		ofile Quebics 10 75"			
-	Quadtro Select High Pr			Thiskness	105 mm
Overall Depth	335 mm Overall Width	425 mm Produc	t Weight 1.4 kg	Thickness	105 mm
R-QS98C Roho	Quadtro Select High Pr	rofile Cushion 16.75" x	15″		
Overall Depth	380 mm Overall Width	425 mm Produc	t Weight 1.4 kg	Thickness	105 mm
R-0S77C Roho	Quadtro Select High Pr	ofile Cushion 13 25" v	13 25"		
Overall Depth	335 mm Overall Width		t Weight 1.4 kg	Thickness	105 mm
R-QS78C Roho	Quadtro Select High Pr	ofile Cushion 13.25" x	15"		
Overall Depth	380 mm Overall Width	335 mm Produc	t Weight 1.4 kg	Thickness	105 mm
R-QS79C Roho	Quadtro Select High Pr	ofile Cushion 13.25" x	16.75"		
Overall Depth	425 mm Overall Width		t Weight 1.4 kg	Thickness	105 mm
D 000100 - i					
-	o Quadtro Select High P			This has a second	105
Overall Depth	470 mm Overall Width	380 mm Produc	t Weight 1.4 kg	Thickness	105 mm
R-QS811C Rohe	o Quadtro Select High P	rofile Cushion 15" x 20	.25″		
Overall Depth	515 mm Overall Width	380 mm Produc	t Weight 1.4 kg	Thickness	105 mm
R-0S87C Roho	Quadtra Salaat High Dr	ofilo Cuchion 15" y 12 () F "		
Overall Depth	Quadtro Select High Pr 335 mm Overall Width		zt Weight 1.4 kg	Thickness	105 mm
eteran pepti			1.4 Kg	Interior	100 1111
R-QS88C Roho	Quadtro Select High Pr	ofile Cushion 15" x 15"			
Overall Depth	380 mm Overall Width	380 mm Produc	t Weight 1.4 kg	Thickness	105 mm
	Quadtro Select High Pr	rofile Cushion 15" x 16 "	/5"		
R-0589C Robo	Quadro ocicor nigil Fi			Thickness	105 mm
R-QS89C Roho	425 mm Overall Width	380 mm Produc	t Weight 1.4 kg	THICKNESS	
•	425 mm Overall Width	380 mm Produ u	t weight 1.4 kg	THICKNESS	100 1111
Overall Depth	425 mm Overall Width			THICKIESS	

Country Care
GROUP
National Distribution, Local Service

R-QS99C Roho Quadtro Select High Profile Cushion 16" x 16"										
Overall Depth	405 mm	Overall Width	405 mm	Product Weight	1.4 kg	Thickness	105 mm			
R-QS109C Roho Quadtro Select High Profile Cushion 18" x 16"										
Overall Depth	405 mm	Overall Width	455 mm	Product Weight	1.4 kg	Thickness	105 mm			
R-QS1110C Roho Quadtro Select High Profile Cushion 20" x 18"										
Overall Depth	455 mm	Overall Width	505 mm	Product Weight	1.4 kg	Thickness	105 mm			

Roho Mosaic Cushion

- The Mosaic Cushion has been designed to provide a comfortable, functional and stable sitting environment
- Suitable for people who are at risk or have a history of skin breakdown . This adjustable, easy to use cushion features interconnected air cells that are easy to set up, inflate and adjust
- This lightweight cushion is easy to transport and handle, thereby reducing fatigue on the user
- A cost effective, uniquely designed seat cushion which offers comfort, protection and adjustability to suit individual requirements



R-MOSAIC1616HD Roho Mosaic Cushion with Heavy Duty Cover 16" x 16"										
Max User Weight	157 kg	Overall Depth	405 mm	Overall Width	405 mm	Product Weight	0.5 kg			
Thickness	70 mm									

R-MOSAIC1816HD Roho Mosaic Cushion with Heavy Duty Cover 18" x 16"

Max User Weight	157 kg	Overall Depth	405 mm	Overall Width	455 mm	Product Weight	0.5 kg
Thickness	70 mm						

R-MOSAIC1818HD Roho Mosaic Cushion with Heavy Duty Cover 18" x 18"

	Max User Weight	157 kg	Overall Depth	455 mm	Overall Width	455 mm	Product Weight	0.5 kg
	Thickness	70 mm						

R-MOSAIC1616C Roho Mosaic Cushion with Standard Cover 16" x 16"

Max User	r Weight	157 kg	Overall Depth	405 mm	Overall Width	405 mm	Product Weight	0.5 kg
Thicknes	s	70 mm						

R-MOSAIC1816C Roho Mosaic Cushion with Standard Cover 18" x 16"

Max User Weight	157 kg	Overall Depth	455 mm	Overall Width	455 mm	Product Weight	0.5 kg
Thickness	70 mm						

R-MOSAIC1818C Roho Mosaic Cushion with Standard Cover 18" x 18"

Max User Weight	157 kg	Overall Depth	455 mm	Overall Width	455 mm	Product Weight	0.5 kg
Thickness	70 mm						

Roho Quadtro Select Mid Profile Cushion

- Allows for quick and simple position selection for every user while seated
- Push the knob to the open position allows adjustment then push the knob closed to maintain this position
- Great for maximum function and posture
- Stability is enhanced by locking air into place in each of the four quadrants
- Minimises side to side or front to back motion
- Lightweight, portable, repairable and washable
- Air is able to circulate around the cells, reaching the surface of the skin and providing a cooling environment

R-QS99MPC Roho Quadtro Select	Mid Profile Cushion 16	5.7
Overall Depth	425 mm	0
R-QS910MPC Roho Quadtro Selec	t Mid Profile Cushion 1	6.
Overall Depth	470 mm	01
R-QS98MPC Roho Quadtro Select	Mid Profile Cushion 16	5.7
Overall Depth	380 mm	0
R-QS109MPC Roho Quadtro Selec	t Mid Profile Cushion 1	8.
R-QS88MPC Roho Quadtro Select	Mid Profile Cushion 15	5"
Overall Depth	380 mm	0
R-QS89MPC Roho Quadtro Select	Mid Profile Cushion 15	5″
Overall Depth	425 mm	0
R-QS1010MPC Roho Quadtro Sele	ct Mid Profile Cushion	18
		_
Overall Depth	470 mm	0

Roho LTV Seat Cushion Fabric Cover



R-LTV-FABRIC	Roho LTV Seat Cush	nion Fabric Cover	
Package Quantity		1	Ра

Roho LTV Seat Cushion with Ultra Leather Cover



- Customize the seat to your specific comfort level Roho shape fitting technology makes the difference
- Made of individual, interconnected air cells that allow air to slowly transfer from chamber to chamber • Evenly distributes body weight pressure
- Lightweight and easy to fold
- It can easily be transported from your scooter, to the car or wherever you may need a cushion • The padded standard cover or optional synthetic leather cover adds an extra level of comfort

R-LTV-ULTRA	Roho LTV S	eat Cushion	with Ultra Leather Co
Overall Depth		430 mm	Overall Width



168

Pressure Cushions

75" x 16.75"	
Dverall Width	425 mm
75" 10 5"	
.75" x 18.5"	
Overall Width	425 mm
75" x 15"	
Dverall Width	425 mm
.5" x 16.75"	
Overall Width	470 mm
x 15″	
Dverall Width	380 mm
x 16.75"	
Overall Width	380 mm
8.5" x 18.5"	
8.5" x 18.5" Jverall Width	470 mm
	470 mm
	470 mm
	470 mm
Dverall Width	470 mm
	470 mm
Dverall Width	470 mm
Dverall Width laily wear and tear giving them a longer life os the cushion do its job	

- The secret to the extraordinary comfort of this cushion is the valve



Roho Mini-Max Cushion



- The Mini Max Cushion consists of a thin pad of interconnected air cells
- Designed for an active user at minimal risk of skin breakdown and is designed to provide a level of impact absorption
- It features an air valve in one corner
- It is able to be customised to size if required
- Minimizes sitting and shearing pressures Excellent for sports

R-MINI-15-C Roho Mini			
Overall Depth	380 mm	Overall Width	380 mm
R-MINI-16-C Roho Mini	i-Max Cushion 16" x 16"		
Overall Depth	405 mm	Overall Width	405 mm
R-MINI-18-C Roho Mini	i-Max Cushion 18" x 16"		
R-MINI-18-C Roho Mini	i-Max Cushion 18" x 16" 405 mm	Overall Width	455 mm
		Overall Width	455 mm

Jay 2 High Risk Pressure Care Cushion



- Pre-contoured foam wheelchair cushions
- Features a Jay flow™ fluid tripad with up to 3" of loading for superior skin protection
- An easy to modify base with optional positioning components for optimal stability
- Provides exceptional lateral and forward stability for individuals with changing needs or unique postures
- Conforms to each individual's shape, adjusts to sitting positions Ensures proper fluid placement beneath bony prominences to help protect the skin from breakdown
- · Zero-maintenance design delivers superior results throughout the life of the cushion

J2108 Jay 2 High Risk Pressure Care Cushion 14" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	355 mm	Product Weight	2.4 kg
Thickness	75 mm						

J2107 Jay 2 High Risk Pressure Care Cushion 16" x 20"

Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	405 mm	Product Weight	2.4 kg
Thickness	75 mm						

J2116 Jay 2 High Risk Pressure Care Cushion 20" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	505 mm	Product Weight	2.4 kg
Thickness	75 mm						

J2100 Jay 2 High Risk Pressure Care Cushion 18" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	455 mm	Product Weight	2.4 kg
Thickness	75 mm						

J2101 Jay 2 High Risk Pressure Care Cushion 18" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	
Thickness	75 mm			

J2102 Jay 2 High Risk Pressure Care Cushion 20" x 20"

Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	505 mm	Product Weight	2.4 kg
Thickness	75 mm						

J2103 Jay 2 High Risk Pressure Care Cushion 18" x 20"

Jeree Jay Ernginna			, , , ,				
Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	455 mm	Product Weight	2.4 kg
Thickness	75 mm						

J2104 Jay 2 High Risk Pressure Care Cushion 24" x 20"

Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	610 mm	Product Weight	2.4 kg
Thickness	75 mm						

12105 1av 2 High Risk Pressure Care Cushion 16" x 16"

JEICO Juy 2 mgm	151 1 1 2 3						
Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	405 mm	Product Weight	2.4 kg
Thickness	75 mm						

J2106 Jay 2 High Risk Pressure Care Cushion 16" x 18"

· · · · · · · · · · · · · · · · · · ·							
Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	405 mm	Product Weight	2.4 kg
Thickness	75 mm						

J2117 Jay 2 High Risk Pressure Care Cushion 17" x 17"

Max User Weight	113 kg	Overall Depth	430 mm	Overall Width	430 mm	Product Weight	2.4 kg
Thickness	75 mm						

J2118 Jay 2 High Risk Pressure Care Cushion 20" x 18"

, ,							
Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	508 mm	Product Weight	2.4 kg
Thickness	75 mm						

Jay 2 Deep Contour Cushion



N

Т

The JAY® J2™ pre-contoured foam wheelchair cushions feature a Jay flow™ fluid tripad with up to 3" of loading for superior skin protection and an easy to modify base with optional positioning components for optimal stability. The J2 cushions provide exceptional lateral and forward stability for individuals with changing needs or unique postures. Conforms to each individual's shape, adjusts to sitting positions, and ensures proper fluid placement beneath bony prominences to help protect the skin from breakdown. Its zero-maintenance design delivers superior results throughout the

J2406 Jay 2 Deep Contour Cushion 20" x 16"

	Product Weight	2.7 kg
Thickness 115 mm		

J2400 Jay 2 Deep Contour Cushion 20" x 20"

	-
Thickness 115 mm	


J2408 Jay 2 Deep Contour Cushion 20" x 18"

	, on to an e						
Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	505 mm	Product Weight	2.7 kg
Thickness	115 mm						

he life of the cushion.		

Overall Width 455 mm Product Weight 2.4 kg





J2426 Jay 2 Deep Contour Cushion 22" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	560 mm	Product Weight	2.7 kg
Thickness	115 mm						

J2428 Jay 2 Deep Contour Cushion 22" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	560 mm	Product Weight	2.7 kg
Thickness	115 mm						

J2466 Jay 2 Deep Contour Cushion 16" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	405 mm	Product Weight	2.7 kg
Thickness	115 mm						

J2480 Jay 2 Deep Contour Cushion 18" x 20"

Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	455 mm	Product Weight	2.7 kg
Thickness	115 mm						

J2444 Jay 2 Deep Contour Cushion 14" x 14"

Max User Weight	113 kg	Overall Depth	355 mm	Overall Width	355 mm	Product Weight	2.7 kg
Thickness	115 mm						

J2446 Jay 2 Deep Contour Cushion 14" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	355 mm	Product Weight	2.7 kg
Thickness	115 mm						

J2486 Jay 2 Deep Contour Cushion 18" x 16"

Max	x User Weight	113 kg	Overall Depth	405 mm	Overall Width	455 mm	Product Weight	2.7 kg
Thio	ckness	115 mm						

J2488 Jay 2 Deep Contour Cushion 18" x 18"

[Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	455 mm	Product Weight	2.7 kg
	Thickness	115 mm						

J2455 Jay 2 Deep Contour Cushion 15" x 15"

Max User Weight	113 kg	Overall Depth	380 mm	Overall Width	380 mm	Product Weight	2.7 kg
Thickness	115 mm						

J2457 Jay 2 Deep Contour Cushion 15" x 17"

Max User Weight	113 kg	Overall Depth	430 mm	Overall Width	380 mm	Product Weight	2.7 kg
Thickness	115 mm						

J2460 Jay 2 Deep Contour Cushion 16" x 20"

Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	405 mm	Product Weight	2.7 kg
Thickness	115 mm						

J2468 Jay 2 Deep Contour Cushion 16" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	405 mm	Product Weight	2.7 kg
Thickness	115 mm						

J2475 Jay 2 Deep Contour Cushion 17" x 15"

Max User Weight	113 kg	Overall Depth	380 mm	Overall Width	430 mm	Product Weight	2.7 kg
Thickness	115 mm						

J2477 Jay 2 Deep Contour Cushion 17" x 17"

Max User Weight	113 kg	Overall Depth	430 mm	Overall Width	430 mm	Product Weight	2.7 kg
Thickness	115 mm						

Jay 2 Plus Cushion



 $\bullet\,$ A pre contoured foam cushion for bariatric clients Features either foam or a fluid pad Provides optimal skin protection and a modifiable base for optimal stability

J302018 Jay 2 Plus Cus	hion 20" x 18	n
Max User Weight	295 kg	Overall Depth
J302019 Jay 2 Plus Cus	hion 20″ x 19'	"
Max User Weight	295 kg	Overall Depth
	-	-
J302020 Jay 2 Plus Cus	hion 20" x 20)″
Max User Weight	295 kg	Overall Depth
J302021 Jay 2 Plus Cus	hion 20" x 21'	1
Max User Weight	295 kg	Overall Depth
J302022 Jay 2 Plus Cus	hion 20" x 22	2"
Max User Weight	295 kg	Overall Depth
J302118 Jay 2 Plus Cush	nion 21" x 18"	
Max User Weight	295 kg	Overall Depth
J302119 Jay 2 Plus Cush	nion 21" x 19"	
Max User Weight	295 kg	Overall Depth
J302120 Jay 2 Plus Cus	hion 21" x 20'	7
Max User Weight	295 kg	Overall Depth
J302121 Jay 2 Plus Cust	nion 21" x 21"	
Max User Weight	295 kg	Overall Depth
J302122 Jay 2 Plus Cus	hion 21" x 22'	,
Max User Weight	295 kg	Overall Depth
J302218 Jay 2 Plus Cus	hion 22" x 18'	π
Max User Weight	295 kg	Overall Depth
J302219 Jay 2 Plus Cus	hion 22" x 19'	7
Max User Weight	295 kg	Overall Depth
J302220 Jay 2 Plus Cus	hion 22" x 20)″
Max User Weight	295 kg	Overall Depth
J302221 Jay 2 Plus Cus	hion 22" x 21'	1
J302221 Jay 2 Plus Cus Max User Weight	hion 22" x 21' 295 kg	Overall Depth



455 mm	Overall Width	505 mm
480 mm	Overall Width	505 mm
505 mm	Overall Width	505 mm
535 mm	Overall Width	505 mm
560 mm	Overall Width	505 mm
455 mm	Overall Width	535 mm
480 mm	Overall Width	535 mm
505 mm	Overall Width	535 mm
535 mm	Overall Width	535 mm
560 mm	Overall Width	535 mm
455 mm	Overall Width	560 mm
480 mm	Overall Width	560 mm
505 mm	Overall Width	560 mm
535 mm	Overall Width	560 mm

J302222 Jay 2 Plus Cu	shion 22" x 2	2"			
Max User Weight	295 kg	Overall Depth	560 mm	Overall Width	560 mm
J302318 Jay 2 Plus Cu	shion 23" x 18	n			
Max User Weight	295 kg	Overall Depth	455 mm	Overall Width	585 mm
J302319 Jay 2 Plus Cu	bion 22" y 10	π			
Max User Weight	295 kg	Overall Depth	480 mm	Overall Width	585 mm
_					
J302320 Jay 2 Plus Cu			1		
Max User Weight	295 kg	Overall Depth	505 mm	Overall Width	585 mm
J302321 Jay 2 Plus Cu	shion 23" x 21	π			
Max User Weight	295 kg	Overall Depth	535 mm	Overall Width	585 mm
J302322 Jay 2 Plus Cu	shion 23" x 2	2"			
Max User Weight	295 kg	Overall Depth	560 mm	Overall Width	585 mm
1202410 Jose 0 Diver 0	hien 0.4" 10				
J302418 Jay 2 Plus Cus Max User Weight	295 kg	Overall Depth	455 mm	Overall Width	610 mm
max user weight	255 Ky	overan beptit	433 11111		010 IIIII
J302419 Jay 2 Plus Cu	shion 24" x 19	n			
Max User Weight	295 kg	Overall Depth	480 mm	Overall Width	610 mm
J302420 Jay 2 Plus Cu	shion 24" x 20)″			
Max User Weight	295 kg	Overall Depth	505 mm	Overall Width	610 mm
1202/21 101 2 Dius Cur	bion 04" v 01	"			
J302421 Jay 2 Plus Cu			535 mm	Overall Width	610 mm
J302421 Jay 2 Plus Cus Max User Weight	shion 24" x 21 295 kg	" Overall Depth	535 mm	Overall Width	610 mm
_	295 kg	Overall Depth	535 mm	Overall Width	610 mm
Max User Weight	295 kg	Overall Depth	535 mm 560 mm	Overall Width Overall Width	610 mm 610 mm
Max User Weight	295 kg shion 24" x 22 295 kg	Overall Depth 2″ Overall Depth			
Max User Weight J302422 Jay 2 Plus Cu Max User Weight	295 kg shion 24" x 22 295 kg	Overall Depth 2″ Overall Depth			
Max User Weight J302422 Jay 2 Plus Cu Max User Weight J302518 Jay 2 Plus Cu Max User Weight	295 kg shion 24" x 22 295 kg shion 25" x 18 295 kg	Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth	560 mm	Overall Width	610 mm
Max User Weight J302422 Jay 2 Plus Cu Max User Weight J302518 Jay 2 Plus Cu Max User Weight J302519 Jay 2 Plus Cu	295 kg shion 24" x 22 295 kg shion 25" x 18 295 kg shion 25" x 19	Overall Depth 2" Overall Depth " Overall Depth "	560 mm 455 mm	Overall Width Overall Width	610 mm 635 mm
Max User Weight J302422 Jay 2 Plus Cu Max User Weight J302518 Jay 2 Plus Cu Max User Weight	295 kg shion 24" x 22 295 kg shion 25" x 18 295 kg	Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth	560 mm	Overall Width	610 mm
Max User Weight J302422 Jay 2 Plus Cu Max User Weight J302518 Jay 2 Plus Cu Max User Weight J302519 Jay 2 Plus Cu	295 kg shion 24" x 22 295 kg shion 25" x 18 295 kg shion 25" x 19 295 kg 295 kg	Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth	560 mm 455 mm	Overall Width Overall Width	610 mm 635 mm
Max User Weight J302422 Jay 2 Plus Cu Max User Weight J302518 Jay 2 Plus Cu Max User Weight J302519 Jay 2 Plus Cu Max User Weight	295 kg shion 24" x 22 295 kg shion 25" x 18 295 kg shion 25" x 19 295 kg 295 kg	Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth	560 mm 455 mm	Overall Width Overall Width	610 mm 635 mm
Max User Weight J302422 Jay 2 Plus Cu Max User Weight J302518 Jay 2 Plus Cu Max User Weight J302519 Jay 2 Plus Cu Max User Weight J302520 Jay 2 Plus Cu Max User Weight	295 kg shion 24" x 22 295 kg shion 25" x 18 295 kg shion 25" x 19 295 kg shion 25" x 20 295 kg	Overall Depth 2" Overall Depth " Overall Depth " Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth	560 mm 455 mm 480 mm	Overall Width Overall Width Overall Width	610 mm 635 mm 635 mm
Max User Weight J302422 Jay 2 Plus Cu Max User Weight J302518 Jay 2 Plus Cu Max User Weight J302519 Jay 2 Plus Cu Max User Weight J302520 Jay 2 Plus Cu	295 kg shion 24" x 22 295 kg shion 25" x 18 295 kg shion 25" x 19 295 kg shion 25" x 20 295 kg	Overall Depth 2" Overall Depth " Overall Depth " Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth	560 mm 455 mm 480 mm	Overall Width Overall Width Overall Width	610 mm 635 mm 635 mm
Max User Weight J302422 Jay 2 Plus Cu Max User Weight J302518 Jay 2 Plus Cu Max User Weight J302519 Jay 2 Plus Cu Max User Weight J302520 Jay 2 Plus Cu Max User Weight J302521 Jay 2 Plus Cu Max User Weight	295 kg shion 24" x 22 295 kg shion 25" x 18 295 kg shion 25" x 19 295 kg shion 25" x 20 295 kg shion 25" x 21 295 kg	Overall Depth Overall Depth " Overall Depth " Overall Depth	560 mm 455 mm 480 mm 505 mm	Overall Width Overall Width Overall Width Overall Width	610 mm 635 mm 635 mm 635 mm
Max User Weight J302422 Jay 2 Plus Cu Max User Weight J302518 Jay 2 Plus Cu Max User Weight J302519 Jay 2 Plus Cu Max User Weight J302520 Jay 2 Plus Cu Max User Weight J302521 Jay 2 Plus Cu Max User Weight J302521 Jay 2 Plus Cu Max User Weight J302522 Jay 2 Plus Cu	295 kg shion 24" x 22 295 kg shion 25" x 18 295 kg shion 25" x 19 295 kg shion 25" x 20 295 kg shion 25" x 21 295 kg shion 25" x 21 295 kg	Overall Depth 2" Overall Depth " Overall Depth " Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth Overall Depth 2"	560 mm 455 mm 480 mm 505 mm 535 mm	Overall Width Overall Width Overall Width Overall Width Overall Width	610 mm 635 mm 635 mm 635 mm
Max User Weight J302422 Jay 2 Plus Cu Max User Weight J302518 Jay 2 Plus Cu Max User Weight J302519 Jay 2 Plus Cu Max User Weight J302520 Jay 2 Plus Cu Max User Weight J302521 Jay 2 Plus Cu Max User Weight J302521 Jay 2 Plus Cu Max User Weight J302522 Jay 2 Plus Cu Max User Weight J302522 Jay 2 Plus Cu Max User Weight	295 kg shion 24" x 22 295 kg shion 25" x 18 295 kg shion 25" x 19 295 kg shion 25" x 21 295 kg shion 25" x 21 295 kg shion 25" x 22 295 kg	Overall Depth Overall Depth " Overall Depth " Overall Depth	560 mm 455 mm 480 mm 505 mm	Overall Width Overall Width Overall Width Overall Width	610 mm 635 mm 635 mm 635 mm
Max User Weight J302422 Jay 2 Plus Cu Max User Weight J302518 Jay 2 Plus Cu Max User Weight J302519 Jay 2 Plus Cu Max User Weight J302520 Jay 2 Plus Cu Max User Weight J302521 Jay 2 Plus Cu Max User Weight J302522 Jay 2 Plus Cu Max User Weight J302522 Jay 2 Plus Cu Max User Weight J302524 Jay 2 Plus Cu Max User Weight J302524 Jay 2 Plus Cu Max User Weight J302524 Jay 2 Plus Cu Max User Weight J302618 Jay 2 Plus Cu	295 kg shion 24" x 22 295 kg shion 25" x 18 295 kg shion 25" x 19 295 kg shion 25" x 21 295 kg shion 25" x 21 295 kg shion 25" x 22 295 kg shion 25" x 23 295 kg	Overall Depth 2" Overall Depth " Overall Depth " Overall Depth Overall Depth	560 mm 455 mm 480 mm 505 mm 535 mm 560 mm	Overall Width Overall Width Overall Width Overall Width Overall Width Overall Width	610 mm 635 mm 635 mm 635 mm 635 mm 635 mm
Max User Weight J302422 Jay 2 Plus Cu Max User Weight J302518 Jay 2 Plus Cu Max User Weight J302519 Jay 2 Plus Cu Max User Weight J302520 Jay 2 Plus Cu Max User Weight J302521 Jay 2 Plus Cu Max User Weight J302521 Jay 2 Plus Cu Max User Weight J302522 Jay 2 Plus Cu Max User Weight J302522 Jay 2 Plus Cu Max User Weight	295 kg shion 24" x 22 295 kg shion 25" x 18 295 kg shion 25" x 19 295 kg shion 25" x 21 295 kg shion 25" x 21 295 kg shion 25" x 22 295 kg	Overall Depth Overall Depth " Overall Depth " Overall Depth	560 mm 455 mm 480 mm 505 mm 535 mm	Overall Width Overall Width Overall Width Overall Width Overall Width	610 mm 635 mm 635 mm 635 mm
Max User Weight J302422 Jay 2 Plus Cu Max User Weight J302518 Jay 2 Plus Cu Max User Weight J302519 Jay 2 Plus Cu Max User Weight J302520 Jay 2 Plus Cu Max User Weight J302521 Jay 2 Plus Cu Max User Weight J302522 Jay 2 Plus Cu Max User Weight J302522 Jay 2 Plus Cu Max User Weight J302524 Jay 2 Plus Cu Max User Weight J302524 Jay 2 Plus Cu Max User Weight J302524 Jay 2 Plus Cu Max User Weight J302618 Jay 2 Plus Cu	295 kg shion 24" x 22 295 kg shion 25" x 18 295 kg shion 25" x 19 295 kg shion 25" x 21 295 kg shion 25" x 21 295 kg shion 25" x 22 295 kg shion 25" x 18 295 kg	Overall Depth 2" Overall Depth " Overall Depth " Overall Depth 0verall Depth Overall Depth 2" Overall Depth Overall Depth Overall Depth	560 mm 455 mm 480 mm 505 mm 535 mm 560 mm	Overall Width Overall Width Overall Width Overall Width Overall Width Overall Width	610 mm 635 mm 635 mm 635 mm 635 mm 635 mm

Max User Weight	295 kg	Overall Depth	505 mm	Overall Width	660 mm
1302621 Jay 2 Plus	Cushion 26" x	21″			
Max User Weight	295 kg	Overall Depth	535 mm	Overall Width	660 mm
1302622 Jay 2 Plus	Cushion 26" x	22"		Ì	
Max User Weight	295 kg	Overall Depth	560 mm	Overall Width	660 mm
Jay J3 Cushion P	ressure Care	e High Risk			
	 Offe Adju Feat Cust Prov 	ides outstanding stability and op rs the choice of Jay flow fluid or st to sitting positions, and ensur ures a contoured foam base con romizable pelvic loading area siz ides ideal forward and lateral sta tt needs	Roho dry floatation air inse e proper fluid placement ber structed from closed cell foa es are selected separately fr	erts, which conform to each ind neath bony prominences am rom the wheelchair cushion wic	ividual
J3CS1212 Jay J3 Higi	h Risk Pressur	e Care Cushion Deep	12" x 12"		
Max User Weight	227 kg	Overall Depth	305 mm	Overall Width	305 mm
I3CS1214 Jay J3 Hig	h Risk Pressur	e Care Cushion Deep	12" x 14"		
Max User Weight	227 kg	Overall Depth	355 mm	Overall Width	305 mm
13CS1412 Jay J3 Hig	h Risk Pressur	e Care Cushion Deep	14" x 12"		
Max User Weight	227 kg	Overall Depth	305 mm	Overall Width	355 mm
13CS1416 Jav 13 Hig	h Risk Pressur	e Care Cushion Deep	14" x 16"		
Max User Weight	227 kg	Overall Depth	405 mm	Overall Width	355 mm
13CS1418 Jay J3 Hig	h Risk Pressur	e Care Cushion Deep	14" x 18"		
Max User Weight	227 kg	Overall Depth	455 mm	Overall Width	355 mm
J3CS1515 Jay J3 Higl	h Risk Pressur	e Care Cushion Deep	15" x 15"		
Max User Weight	227 kg	Overall Depth	380 mm	Overall Width	380 mm
13CS1516 Jay J3 Hig	h Risk Pressur	e Care Cushion Deep	15" x 16"		
Max User Weight	227 kg	Overall Depth	405 mm	Overall Width	380 mm
13CS1517 Jay J3 Higl	n Risk Pressur	e Care Cushion Deep	15" x 17"		
	227 kg	Overall Depth	430 mm	Overall Width	380 mm
Max User Weight			10" 15"		
	h Risk Pressur	e Care Cushion Deep	16" X 15"		
	h Risk Pressur	e Care Cushion Deep Overall Depth	380 mm	Overall Width	405 mm
13CS1615 Jay J3 Higi Max User Weight	227 kg	-	380 mm	Overall Width	405 mm
J3CS1615 Jay J3 Higi Max User Weight	227 kg	Overall Depth	380 mm	Overall Width Overall Width	405 mm
J3CS1615 Jay J3 Hig Max User Weight J3CS1617 Jay J3 Hig Max User Weight	227 kg h Risk Pressur 227 kg	Overall Depth	380 mm 16" x 17" 430 mm		





3CS1619 Jay J3 Hig	h Risk Pressur	e Care Cushion Deep	J IG X I9		
Max User Weight	227 kg	Overall Depth	480 mm	Overall Width	405 mm
3CS1620 Jay J3 Hig	h Risk Pressu	re Care Cushion Deep	o 16" x 20"		
Max User Weight	227 kg	Overall Depth	505 mm	Overall Width	405 mm
200171E 1au 10 Utal	h Diels Dreeser	o Come Qualtiere Door	17" 15"		
		e Care Cushion Deep	380 mm	Overall Width	430 mm
Max User Weight	227 kg	Overall Depth	380 11111		430 11111
3CS1716 Jay J3 Higl	h Risk Pressur	e Care Cushion Deep	17" x 16"		
Max User Weight	227 kg	Overall Depth	405 mm	Overall Width	430 mm
13CS1717 Jav J3 High	n Risk Pressur	e Care Cushion Deep	17" x 17"		
Max User Weight	227 kg	Overall Depth	430 mm	Overall Width	430 mm
		e Care Cushion Deep			
Max User Weight	227 kg	Overall Depth	455 mm	Overall Width	430 mm
3CS1816 Jay J3 Hig	h Risk Pressur	re Care Cushion Deep) 18" x 16"		
Max User Weight	227 kg	Overall Depth	405 mm	Overall Width	455 mm
19001017 Jaw 10 Ukal	h Diele Dresser	o Oana Quahian Daan	10" + 17"		
	227 kg	e Care Cushion Deep	430 mm	Overall Width	455 mm
Max User Weight	227 Ky	Overall Depth	430 11111		455 11111
3CS1820 Jay J3 Hig	h Risk Pressu	re Care Cushion Deep	p 18" x 20"		
Max User Weight	227 kg	Overall Depth	505 mm	Overall Width	455 mm
13CS1916 Jay J3 Hig	h Diek Drossu	re Care Cushion Deep	19" x 16"		
Max User Weight	227 kg	Overall Depth	405 mm	Overall Width	480 mm
		e Care Cushion Deep			
Max User Weight	227 kg	Overall Depth	430 mm	Overall Width	480 mm
I3CS1918 Jay J3 Higi	h Risk Pressur	re Care Cushion Deep) 19" x 18"		
Max User Weight	227 kg	Overall Depth	455 mm	Overall Width	480 mm
12001020 1av 10 11	h Diele Dresser	a Cara Quahian Daa	. 10" 00"		
, , ,		re Care Cushion Deep		Querall Width	490 mm
J3CS1920 Jay J3 Hig Max User Weight	Jh Risk Pressu 227 kg	re Care Cushion Deep Overall Depth	p 19" x 20"	Overall Width	480 mm
Max User Weight	227 kg	-	505 mm	Overall Width	480 mm
Max User Weight	227 kg	Overall Depth	505 mm	Overall Width Overall Width	480 mm 505 mm
Max User Weight J3CS2016 Jay J3 Hig Max User Weight	227 kg J h Risk Pressu 227 kg	Overall Depth re Care Cushion Deep Overall Depth	505 mm p 20" x 16" 405 mm		
Max User Weight I3CS2016 Jay J3 Hig Max User Weight	227 kg J h Risk Pressu 227 kg	Overall Depth	505 mm p 20" x 16" 405 mm		
Max User Weight I3CS2016 Jay J3 Hig Max User Weight I3CS2018 Jay J3 Hig	227 kg Jh Risk Pressu 227 kg Jh Risk Pressu	Overall Depth re Care Cushion Deep Overall Depth re Care Cushion Deep	505 mm p 20" x 16" 405 mm p 20" x 18"	Overall Width	505 mm
Max User Weight I3CS2016 Jay J3 Hig Max User Weight I3CS2018 Jay J3 Hig Max User Weight	227 kg Jh Risk Pressu 227 kg Jh Risk Pressu 227 kg	Overall Depth re Care Cushion Deep Overall Depth re Care Cushion Deep	505 mm p 20" x 16" 405 mm p 20" x 18" 455 mm	Overall Width	505 mm
Max User Weight I3CS2016 Jay J3 Hig Max User Weight I3CS2018 Jay J3 Hig Max User Weight	227 kg Jh Risk Pressu 227 kg Jh Risk Pressu 227 kg	Overall Depth re Care Cushion Deep Overall Depth re Care Cushion Deep Overall Depth Overall Depth	505 mm p 20" x 16" 405 mm p 20" x 18" 455 mm	Overall Width	505 mm
Max User Weight J3CS2016 Jay J3 Hig Max User Weight J3CS2018 Jay J3 Hig Max User Weight J3CS2118 Jay J3 Hig Max User Weight	227 kg Jh Risk Pressu 227 kg Jh Risk Pressu 227 kg h Risk Pressur 227 kg	Overall Depth re Care Cushion Deep Overall Depth	505 mm p 20" x 16" 405 mm p 20" x 18" 455 mm p 21" x 18" 455 mm	Overall Width Overall Width	505 mm 505 mm
Max User Weight J3CS2016 Jay J3 Hig Max User Weight J3CS2018 Jay J3 Hig Max User Weight J3CS2118 Jay J3 Hig Max User Weight	227 kg Jh Risk Pressu 227 kg Jh Risk Pressu 227 kg h Risk Pressur 227 kg	Overall Depth re Care Cushion Deep Overall Depth re Care Cushion Deep Overall Depth re Care Cushion Deep	505 mm p 20" x 16" 405 mm p 20" x 18" 455 mm p 21" x 18" 455 mm	Overall Width Overall Width	505 mm 505 mm
Max User Weight I3CS2016 Jay J3 Hig Max User Weight I3CS2018 Jay J3 Hig Max User Weight I3CS2118 Jay J3 Hig Max User Weight I3CS2119 Jay J3 Hig Max User Weight	227 kg Jh Risk Pressu 227 kg Jh Risk Pressu 227 kg h Risk Pressur 227 kg h Risk Pressur 227 kg	Overall Depth re Care Cushion Deep Overall Depth	505 mm p 20" x 16" 405 mm p 20" x 18" 455 mm p 21" x 18" 455 mm p 21" x 19" 480 mm	Overall Width Overall Width Overall Width Overall Width	505 mm 505 mm 505 mm 535 mm
Max User Weight J3CS2016 Jay J3 Hig Max User Weight J3CS2018 Jay J3 Hig Max User Weight J3CS2118 Jay J3 Hig Max User Weight J3CS2119 Jay J3 Hig Max User Weight	227 kg Jh Risk Pressu 227 kg Jh Risk Pressu 227 kg h Risk Pressur 227 kg h Risk Pressur 227 kg	Overall Depth re Care Cushion Deep Te Care Cushion Deep	505 mm p 20" x 16" 405 mm p 20" x 18" 455 mm p 21" x 18" 455 mm p 21" x 19" 480 mm	Overall Width Overall Width Overall Width Overall Width	505 mm 505 mm 505 mm 535 mm

lax User Weight	227 kg	Overall Depth	535 mm	Overall Width	535 mm
3CS2122 Jay J3 Hi	gh Risk Pressu	ire Care Cushion Deep	21" x 22"		
Max User Weight	227 kg	Overall Depth	560 mm	Overall Width	535 mm
3CS2218 Jay J3 Hi	gh Risk Pressu	ire Care Cushion Deep) 22" x 18"		
Max User Weight	227 kg	Overall Depth	455 mm	Overall Width	560 mm
3CS2219 1av 13 Hi	ah Risk Pressu	re Care Cushion Deer	22" x 19"		
Max User Weight	227 kg	Overall Depth	480 mm	Overall Width	560 mm
		ure Care Cushion Dee	o 22" x 20"		
Max User Weight	227 kg	Overall Depth	505 mm	Overall Width	560 mm
3CS2221 Jay 13 Hi	gh Risk Pressu	ire Care Cushion Deep) 22" x 21"		
Max User Weight	227 kg	Overall Depth	535 mm	Overall Width	560 mm
1000000					
-	-	ure Care Cushion Deep			
Max User Weight	227 kg	Overall Depth	560 mm	Overall Width	560 mm
13CS2318 Jay J3 Hi	-	ire Care Cushion Deep	23" x 18"		
Max User Weight	227 kg	Overall Depth	455 mm	Overall Width	585 mm
13CS2319 Jay J3 Hi	gh Risk Pressu	ire Care Cushion Deep	23" x 19"		
Max User Weight	227 kg	Overall Depth	480 mm	Overall Width	585 mm
	Ū.		400 11111	overall width	000 1111
13662320 10V 12 U	igh Dick Proces	ure Care Cuchion Door			000 mm
-	-	ure Care Cushion Deep	p 23" x 20"		
J3CS2320 Jay J3 H Max User Weight	igh Risk Press	ure Care Cushion Dee		Overall Width	585 mm
Max User Weight	227 kg	-	p 23" x 20"		
Max User Weight	227 kg	Overall Depth	p 23" x 20"		
Max User Weight J3CS2321 Jay J3 Hi Max User Weight	227 kg gh Risk Pressu 227 kg	Overall Depth	p 23" x 20" 505 mm 0 23" x 21" 535 mm	Overall Width	585 mm
Max User Weight J3CS2321 Jay J3 Hi Max User Weight	227 kg gh Risk Pressu 227 kg	Overall Depth Ire Care Cushion Deep Overall Depth Overall Depth	p 23" x 20" 505 mm 0 23" x 21" 535 mm	Overall Width	585 mm
Max User Weight 13CS2321 Jay J3 Hi Max User Weight 13CS2322 Jay J3 H Max User Weight	227 kg gh Risk Pressu 227 kg igh Risk Pressu 227 kg	Overall Depth Ire Care Cushion Deep Overall Depth Ure Care Cushion Deep	p 23" x 20" 505 mm 23" x 21" 535 mm p 23" x 22" 560 mm	Overall Width Overall Width	585 mm
Max User Weight I3CS2321 Jay J3 Hi Max User Weight I3CS2322 Jay J3 H Max User Weight	227 kg gh Risk Pressu 227 kg igh Risk Pressu 227 kg	Overall Depth Ire Care Cushion Deep Overall Depth Ure Care Cushion Deep Overall Depth	p 23" x 20" 505 mm 23" x 21" 535 mm p 23" x 22" 560 mm	Overall Width Overall Width	585 mm
Max User Weight J3CS2321 Jay J3 Hi Max User Weight J3CS2322 Jay J3 H Max User Weight J3CS2418 Jay J3 Hi Max User Weight	227 kg 227 kg 227 kg 227 kg 227 kg 227 kg 227 kg 227 kg 227 kg 227 kg	Overall Depth Ire Care Cushion Deep Overall Depth ure Care Cushion Deep Overall Depth Ire Care Cushion Deep	p 23" x 20" 505 mm p 23" x 21" 535 mm p 23" x 22" 560 mm p 24" x 18" 455 mm	Overall Width Overall Width Overall Width Overall Width	585 mm 585 mm 585 mm
Max User Weight I3CS2321 Jay J3 Hi Max User Weight I3CS2322 Jay J3 H Max User Weight I3CS2418 Jay J3 Hi Max User Weight	227 kg 227 kg 227 kg 227 kg 227 kg 227 kg 227 kg 227 kg 227 kg 227 kg	Overall Depth Ire Care Cushion Deep Overall Depth Ure Care Cushion Deep Overall Depth Ire Care Cushion Deep Overall Depth	p 23" x 20" 505 mm p 23" x 21" 535 mm p 23" x 22" 560 mm p 24" x 18" 455 mm	Overall Width Overall Width Overall Width Overall Width	585 mm 585 mm 585 mm
Max User Weight J3CS2321 Jay J3 Hi Max User Weight J3CS2322 Jay J3 H Max User Weight J3CS2418 Jay J3 Hi Max User Weight J3CS2419 Jay J3 Hi Max User Weight	227 kg gh Risk Pressu 227 kg igh Risk Pressu 227 kg gh Risk Pressu 227 kg 227 kg 227 kg 227 kg	Overall Depth Ire Care Cushion Deep Overall Depth ure Care Cushion Deep Overall Depth Ire Care Cushion Deep Overall Depth Ire Care Cushion Deep Overall Depth Ire Care Cushion Deep Ire Care Cushion Deep Ire Care Cushion Deep	p 23" x 20" 505 mm 23" x 21" 535 mm p 23" x 22" 560 mm p 24" x 18" 455 mm p 24" x 19" 480 mm	Overall Width Overall Width Overall Width Overall Width Overall Width	585 mm 585 mm 585 mm 585 mm 610 mm
Max User Weight I3CS2321 Jay J3 Hi Max User Weight I3CS2322 Jay J3 H Max User Weight I3CS2418 Jay J3 Hi Max User Weight I3CS2419 Jay J3 Hi Max User Weight	227 kg gh Risk Pressu 227 kg igh Risk Pressu 227 kg gh Risk Pressu 227 kg 227 kg 227 kg 227 kg	Overall Depth Ire Care Cushion Deep Overall Depth ure Care Cushion Deep Overall Depth Ire Care Cushion Deep Overall Depth Ire Care Cushion Deep Overall Depth Ire Care Cushion Deep Overall Depth Overall Depth Overall Depth	p 23" x 20" 505 mm 23" x 21" 535 mm p 23" x 22" 560 mm p 24" x 18" 455 mm p 24" x 19" 480 mm	Overall Width Overall Width Overall Width Overall Width Overall Width	585 mm 585 mm 585 mm 585 mm 610 mm
Max User Weight J3CS2321 Jay J3 Hi Max User Weight J3CS2322 Jay J3 H Max User Weight J3CS2418 Jay J3 Hi Max User Weight J3CS2419 Jay J3 Hi Max User Weight J3CS2420 Jay J3 H Max User Weight	227 kg gh Risk Pressu 227 kg igh Risk Pressu 227 kg gh Risk Pressu 227 kg gh Risk Pressu 227 kg igh Risk Pressu 227 kg igh Risk Pressu 227 kg	Overall Depth Ire Care Cushion Deep Overall Depth Ure Care Cushion Deep Overall Depth Ire Care Cushion Deep Ire Care Cushion Deep Overall Depth	p 23" x 20" 505 mm p 23" x 21" 535 mm p 23" x 22" 560 mm p 24" x 18" 455 mm p 24" x 19" 480 mm p 24" x 20" 505 mm	Overall Width Overall Width Overall Width Overall Width Overall Width	585 mm 585 mm 585 mm 585 mm 610 mm 610 mm
Max User Weight I3CS2321 Jay J3 Hi Max User Weight I3CS2322 Jay J3 H Max User Weight I3CS2418 Jay J3 Hi Max User Weight I3CS2419 Jay J3 Hi Max User Weight I3CS2420 Jay J3 H Max User Weight	227 kg gh Risk Pressu 227 kg igh Risk Pressu 227 kg gh Risk Pressu 227 kg gh Risk Pressu 227 kg igh Risk Pressu 227 kg igh Risk Pressu 227 kg	Overall Depth Ire Care Cushion Deep Overall Depth ure Care Cushion Deep Overall Depth Ire Care Cushion Deep Overall Depth Overall Depth	p 23" x 20" 505 mm p 23" x 21" 535 mm p 23" x 22" 560 mm p 24" x 18" 455 mm p 24" x 19" 480 mm p 24" x 20" 505 mm	Overall Width Overall Width Overall Width Overall Width Overall Width	585 mm 585 mm 585 mm 585 mm 610 mm 610 mm
Max User Weight J3CS2321 Jay J3 Hi Max User Weight J3CS2322 Jay J3 H Max User Weight J3CS2418 Jay J3 Hi Max User Weight J3CS2419 Jay J3 Hi Max User Weight J3CS2420 Jay J3 H Max User Weight J3CS2421 Jay J3 Hi Max User Weight	gh Risk Pressu 227 kg 227 kg	Overall Depth Ire Care Cushion Deep Overall Depth ure Care Cushion Deep Overall Depth Ire Care Cushion Deep Overall Depth ure Care Cushion Deep Ire Care Cushion Deep Ire Care Cushion Deep	p 23" x 20" 505 mm 23" x 21" 535 mm p 23" x 22" 560 mm p 24" x 18" 455 mm p 24" x 19" 480 mm p 24" x 20" 505 mm	Overall Width	585 mm 585 mm 585 mm 610 mm 610 mm 610 mm 610 mm
Max User Weight J3CS2321 Jay J3 Hi Max User Weight J3CS2322 Jay J3 H Max User Weight J3CS2418 Jay J3 Hi Max User Weight J3CS2419 Jay J3 Hi Max User Weight J3CS2420 Jay J3 H Max User Weight J3CS2421 Jay J3 Hi Max User Weight	gh Risk Pressu 227 kg 227 kg	Overall Depth Ire Care Cushion Deep Overall Depth ure Care Cushion Deep Overall Depth Ire Care Cushion Deep Overall Depth	p 23" x 20" 505 mm 23" x 21" 535 mm p 23" x 22" 560 mm p 24" x 18" 455 mm p 24" x 19" 480 mm p 24" x 20" 505 mm	Overall Width	585 mm 585 mm 585 mm 610 mm 610 mm 610 mm 610 mm
Max User Weight J3CS2321 Jay J3 Hi Max User Weight J3CS2322 Jay J3 H Max User Weight J3CS2418 Jay J3 Hi Max User Weight J3CS2419 Jay J3 Hi Max User Weight J3CS2420 Jay J3 H Max User Weight J3CS2421 Jay J3 Hi Max User Weight J3CS2422 Jay J3 H	gh Risk Pressu 227 kg 227 kg	Overall Depth Ire Care Cushion Deep Overall Depth ure Care Cushion Deep Overall Depth Ire Care Cushion Deep Ire Care Cushion Deep Overall Depth	p 23" x 20" 505 mm 23" x 21" 535 mm p 23" x 22" 560 mm p 24" x 18" 455 mm p 24" x 19" 480 mm p 24" x 20" 505 mm p 24" x 21" 535 mm p 24" x 22" 560 mm	Overall Width Overall Width	585 mm 585 mm 585 mm 610 mm 610 mm 610 mm 610 mm 610 mm 610 mm

Country Care G R O U P National Distribution. Local Service

J3CS1414 Jay J3 High Risk Pressure Care Cushion Deep 14" x 14"									
Max User Weight 227 kg Overall Depth 355 mm Overall Width 355 mm									
J3CS1616 Jay J3 High Risk Pressure Care Cushion Deep 16" x 16"									
Max User Weight	227 kg	Overall Depth	405 mm	Overall Width	405 mm				
J3CS1818 Jay J3 High Risk Pressure Care Cushion Deep 18" x 18"									
I3CS1818 Jay J3 High Ris	sk Pressure C	Care Cushion Deep 18" x 1	18″						
J3CS1818 Jay J3 High Ris Max User Weight	sk Pressure C 227 kg	Care Cushion Deep 18" x 1 Overall Depth	8″ 455 mm	Overall Width	455 mm				
Max User Weight	227 kg	-	455 mm	Overall Width	455 mm				

Jay Duo Cushion



- For shear reduction and low seat to floor height for foot propulsion
- Features an integrated, contoured solid seat pan with a soft foam and fluid seat cushion
- The duo provides optimal comfort and skin protection without affecting seat to floor height
- For the client at moderate risk of skin breakdown, who requires mild positioning

J700 Jay Duo Cushion 18" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	455 mm			
J701 Jay Duo Cushion 18	3° X 18°							
Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	455 mm			
J705 Jay Duo Cushion 16" x 16"								
Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	405 mm			
J706 Jay Duo Cushion 1	6" x 18"							
Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	405 mm			
J709 Jay Duo Cushion 2	0" x 18"							
Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	505 mm			
	^ 							
J726 Jay Duo Cushion 2	0" x 16"							
Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	505 mm			
J776 Jay Duo Cushion 17	7″ x 16″							
Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	430 mm			
J778 Jay Duo Cushion 17	7" x 18"							
Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	430 mm			

Jay Easy Cushion



- A skin protection and positioning wheelchair cushion
- Encourages proper positioning of the pelvis and thighs • The Easy features a small flow fluid insert
- prominences to help protect the skin from breakdown
- Zero-maintenance design delivers superior results for the life of the cushion
- A moisture-resistant cover protects the foam • Its no-slip surface keeps the cushion in place during daily use and transfers

JE1414C Jay Easy Cushion Curved 14" x 14" Max User Weight 113 kg Overall Depth 355 mm 65 mm Thickness JE1414F Jay Easy Cushion Flat 14" x 14" Max User Weight 113 kg Overall Depth 355 mm 65 mm Thickness JE1416C Jay Easy Cushion Curved 14" x 16" Max User Weight 113 kg 405 mm **Overall Depth** Thickness 65 mm JE1416F Jay Easy Cushion Flat 14" x 16" Max User Weight 405 mm 113 kg Overall Depth Thickness 65 mm JE1418C Jay Easy Cushion Curved 14" x 18" Max User Weight 113 kg Overall Depth 455 mm 65 mm Thickness

IE1418F Jay Easy Cushion Flat 14" x 18"							
Max User Weight	113 kg	Overall Depth	455 mm				
Thickness	65 mm						

JE1515C Jay Easy Cushion Curved 15" x 15"

Max User Weight	113 kg	Overall Depth	380 mm	Overall Width	380 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1515F Jay Easy Cushion Flat 15" x 15"

Max User Weight	113 kg	Overall Depth	380 mm	Overall Width	380 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1516C Jay Easy Cushion Curved 15" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	380 mm	Product Weight	2.7 kg
Thickness	65 mm						

1F1516F lav Easy Cushion Flat 15" x 16"

JEIJIOF JAY EASY CUSHIOH FIAL IS X IO							
Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	380 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1517F Jay Easy Cushion Curved 15" x 17"

Max User Weight	113 kg	Overall Depth	430 mm	Overall Width	380 mm	Product Weight	2.7 kg
Thickness	65 mm						



- Featuring a hi-resiliency, contoured foam base that accommodates a curved or flat seating surface and Jay flow™ fluid tripad
- Conforms to each individual's shape, adjusts to sitting positions and ensures proper fluid placement beneath bony

Overall Width	355 mm	Product Weight	2.7 kg

Overall Width	355 mm	Product Weight	2.7 kg

Overall Width	355 mm	Product Weight	2.7 kg

Overall Width	355 mm	Product Weight	2.7 kg

Overall Width	355 mm	Product Weight	2.7 kg

Overall Width	355 mm	Product Weight	2.7 kg

JE1615F Jay Easy Cushion Flat 16" x 15"

Max User Weight	113 kg	Overall Depth	380 mm	Overall Width	405 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1617C Jay Easy Cushion Curved 16" x 17"

Max User Weight	113 kg	Overall Depth	430 mm	Overall Width	405 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1617F Jay Easy Cushion Flat 16" x 17"

Max User Weight	113 kg	Overall Depth	430 mm	Overall Width	405 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1618C Jay Easy Cushion Curved 16" x 18"

[Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	405 mm	Product Weight	2.7 kg
	Thickness	65 mm						

JE1618F Jay Easy Cushion Flat 16" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	405 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1716C Jay Easy Cushion Curved 17" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	430 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1716F Jay Easy Cushion Flat 17" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	430 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1717C Jay Easy Cushion Curved 17" x 17"

Max User Weight	113 kg	Overall Depth	430 mm	Overall Width	430 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1717F Jay Easy Cushion Flat 17" x 17"

Max User Weight	113 kg	Overall Depth	430 mm	Overall Width	430 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1718C Jay Easy Cushion Curved 17" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	430 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1718F Jay Easy Cushion Flat 17" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	430 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1817C Jay Easy Cushion Curved 18" x 17"

Max User Weight	113 kg	Overall Depth	430 mm	Overall Width	455 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1817F Jay Easy Cushion Flat 18" x 17"

Max User Weight	113 kg	Overall Depth	430 mm	Overall Width	455 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1916C Jay Easy Cushion Curved 19" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	480 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1917C Jay Easy Cushion Curved 19" x 17"

Max User Weight	113 kg	Overall Depth	430 mm	Overall Width	480 mm	Product Weight	2.7 kg				
Thickness	65 mm										

JE1918C Jay Easy Cushion Curved 19" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	480 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1918F Jay Easy Cushion Flat 19" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	480 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1920C Jay Easy Cushion Curved 19" x 20"

Max User Weight	113 kg	Overall Depth	505 mm
Thickness	65 mm		

JE2016C Jay Easy Cushion Curved 20" x 16"

	Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	505 mm	Product Weight	2.7 kg	
	Thickness	65 mm							

JE2016F Jay Easy Cushion Flat 20" x 16"

Max User Weight	113 kg	Overall Depth	405 mm
Thickness	65 mm		

JE2218C Jay Easy Cushion Curved 22" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	560 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE2218F Jay Easy Cushion Flat 22" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	560 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE2220C Jay Easy Cushion Curved 22" x 20"

lax User Weight	113 kg	Overall Depth	505 mm	Overall Width	560 mm	Product Weight	2.7 kg
hickness	65 mm						

JE2220F Jay Easy Cushion Flat 22" x 20"

Thi

Thickness

Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	560 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE2222C Jay Easy Cushion Curved 22" x 22"

65 mm

Max User Weight	113 kg	Overall Depth	560 mm

JE2222F Jay Easy Cushion Flat 22" x 22"

Max User Weight	113 kg	Overall Depth	560 mm	Overall Width	560 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE2418C Jay Easy Cushion Curved 24" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	610 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE2418F Jay Easy Cushion Flat 24" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	610 mm	Product Weight	2.7 kg
Thickness	65 mm						



Overall Width	480 mm	Product Weight	2.7 kg

Overall Width	505 mm	Product Weight	2.7 kg

Overall Width	560 mm	Product Weight	2.7 kg

JE2420C Jay Easy Cushion Curved 24" x 20"

Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	610 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE2420F Jay Easy Cushion Flat 24" x 20"

Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	610 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE2422C Jay Easy Cushion Curved 24" x 22"

Max User Weight	113 kg	Overall Depth	560 mm	Overall Width	610 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE2422F Jay Easy Cushion Flat 24" x 22"

Max User Weight	113 kg	Overall Depth	560 mm	Overall Width	610 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE2424C Jay Easy Cushion Curved 24" x 24"

Max User Weight	113 kg	Overall Depth	610 mm	Overall Width	610 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE2424F Jay Easy Cushion Flat 24" x 24"

Max User We	eight	113 kg	Overall Depth	610 mm	Overall Width	610 mm	Product Weight	2.7 kg
Thickness		65 mm						

JE1616C Jay Easy Cushion Curved 16" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	405 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1616F Jay Easy Cushion Flat 16" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	405 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1816C Jay Easy Cushion Curved 18" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	455 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1816F Jay Easy Cushion Flat 18" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	455 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1818C Jay Easy Cushion Curved 18" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	455 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1818F Jay Easy Cushion Flat 18" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	455 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1820C Jay Easy Cushion Curved 18" x 20"

Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	455 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE1820F Jay Easy Cushion Flat 18" x 20"

[Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	455 mm	Product Weight	2.7 kg
	Thickness	65 mm						

1E2018C lay Fasy Cushion Curved 20" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	505 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE2018F Jay Easy Cushion Flat 20" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	505 mm	Product Weight	2.7 kg
Thickness	65 mm						

JE2020C Jay Easy Cushion Curved 20" x 20"

	Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	505 mm	Product Weight	2.7 kg
	Thickness	65 mm						

JE2020F Jay Easy Cushion Flat 20" x 20"

Max User Weight	113 kg	Overall Depth	505 mm
Thickness	65 mm		

Jay Fusion Cushion



 An adjustable cushion designed for skin protection • Features Jay Flow Fluid or Roho Dry Floatation air insert options • A package thats comfortable, stable and supportive

JFUSION1414	Jay Fusion	Cushion 14"	x 14″
Max User Weight		136 kg	Overall Depth
JFUSION1415	Jay Fusion	Cushion 14" x	x 15″
Max User Weight		136 kg	Overall Depth
JFUSION1416	Jay Fusion	Cushion 14"	x 16″
Max User Weight		136 kg	Overall Depth
JFUSION1417	Jay Fusion	Cushion 14" >	k 17"
Max User Weight		136 kg	Overall Depth
JFUSION1418	Jay Fusion	Cushion 14" x	x 18″
Max User Weight		136 kg	Overall Depth
JFUSION1514	Jay Fusion	Cushion 15" x	x 14"
Max User Weight		136 kg	Overall Depth
JFUSION1515	Jay Fusion	Cushion 15" >	x 15″
Max User Weight		136 kg	Overall Depth
JFUSION1516	Jay Fusion	Cushion 15"	x 16″
Max User Weight		136 kg	Overall Depth
JFUSION1517	Jay Fusion	Cushion 15" >	k 17"
Max User Weight		136 kg	Overall Depth

Overall Width	505 mm	Product Weight	2.7 kg

• The Jay Fusion Cushion combines the J2's skin protection capabilities with the J3's maximum stability

	355 mm	Overall Width	355 mm
	380 mm	Overall Width	355 mm
	405 mm	Overall Width	355 mm
	430 mm	Overall Width	355 mm
	455 mm	Overall Width	355 mm
	355 mm	Overall Width	380 mm
_			
	380 mm	Overall Width	380 mm
	405 mm	Overall Width	380 mm
_			
	430 mm	Overall Width	380 mm

	Jay Fusion	Cushion 15"	x 18″			
Max User Weight		136 kg	Overall Depth	455 mm	Overall Width	380 mm
JFUSION1519	Jay Fusion	Cushion 15"	x 19″			
Max User Weight	-	136 kg	Overall Depth	480 mm	Overall Width	380 mm
JFUSION1520	lav Fusior	1 Cushion 15"	x 20"			
Max User Weight		136 kg	Overall Depth	505 mm	Overall Width	380 mm
JFUSION1521	lay Fusion	Cushion 15"	v 21"			
Max User Weight	Jayrasion	136 kg	Overall Depth	535 mm	Overall Width	380 mm
1EUIGION161/		Cuchico 16"	× 1 <i>4</i> "			
Max User Weight	Jay Fusion	Cushion 16"	X 14 Overall Depth	355 mm	Overall Width	405 mm
				1		
JFUSION1615 Max User Weight	Jay Fusion	Cushion 16"	x 15" Overall Depth	380 mm	Overall Width	405 mm
Max USEI Weight		130 Kg		300 mm		405 1111
JFUSION1616	Jay Fusion	Cushion 16"				
Max User Weight		136 kg	Overall Depth	405 mm	Overall Width	405 mm
JFUSION1617	Jay Fusion	Cushion 16"	x 17"			
Max User Weight		136 kg	Overall Depth	430 mm	Overall Width	405 mm
JFUSION1618	Jay Fusion	Cushion 16"	x 18″			
Max User Weight		136 kg	Overall Depth	455 mm	Overall Width	405 mm
JFUSION1619	Jay Fusion	Cushion 16"	x 19″			
Max User Weight		136 kg	Overall Depth	480 mm	Overall Width	405 mm
JFUSION1620	Jay Fusion	Cuchion 16"	v 20"			
			X 20			
Max User Weight		136 kg	Overall Depth	505 mm	Overall Width	405 mm
	lay Fusion	136 kg	Overall Depth	505 mm	Overall Width	405 mm
	Jay Fusion	1	Overall Depth	505 mm 535 mm	Overall Width Overall Width	405 mm
JFUSION1621 Max User Weight	-	136 kg	Overall Depth X 21" Overall Depth			
JFUSION1621 Max User Weight	-	136 kg	Overall Depth X 21" Overall Depth			
JFUSION1621 Max User Weight JFUSION1622 Max User Weight	Jay Fusion	136 kg Cushion 16" 136 kg Cushion 16" 136 kg 136 kg	Overall Depth x 21" Overall Depth x 22" Overall Depth Overall Depth	535 mm	Overall Width	405 mm
JFUSION1621 Max User Weight JFUSION1622 Max User Weight JFUSION1714	Jay Fusion	136 kg Cushion 16" 136 kg Cushion 16" 136 kg Cushion 16" 136 kg	Overall Depth x 21" Overall Depth x 22" Overall Depth x 14"	535 mm	Overall Width Overall Width	405 mm
JFUSION1621 Max User Weight JFUSION1622 Max User Weight JFUSION1714 Max User Weight	Jay Fusion Jay Fusion	136 kg Cushion 16" 136 kg Cushion 16" 136 kg Cushion 17" 136 kg	Overall Depth x 21" Overall Depth x 22" Overall Depth x 14" Overall Depth	535 mm 560 mm	Overall Width	405 mm
JFUSION1621 Max User Weight JFUSION1622 Max User Weight JFUSION1714 Max User Weight JFUSION1715	Jay Fusion Jay Fusion	136 kg Cushion 16" 136 kg Cushion 16" 136 kg Cushion 17" 136 kg Cushion 17" Cushion 17"	Overall Depth x 21" Overall Depth x 22" Overall Depth x 14" Overall Depth x 15"	535 mm 560 mm 355 mm	Overall Width Overall Width Overall Width	405 mm 405 mm 430 mm
JFUSION1621 Max User Weight JFUSION1622 Max User Weight JFUSION1714 Max User Weight JFUSION1715 Max User Weight	Jay Fusion Jay Fusion Jay Fusion	136 kg Cushion 16" 136 kg Cushion 16" 136 kg Cushion 17" 136 kg Cushion 17" 136 kg	Overall Depth x 21" Overall Depth x 22" Overall Depth x 14" Overall Depth x 15" Overall Depth	535 mm 560 mm	Overall Width Overall Width	405 mm
JFUSION1621 Max User Weight JFUSION1622 Max User Weight JFUSION1714 Max User Weight JFUSION1715 Max User Weight JFUSION1716	Jay Fusion Jay Fusion Jay Fusion	136 kg Cushion 16" 136 kg Cushion 16" 136 kg Cushion 17" 136 kg Cushion 17" 136 kg Cushion 17" 136 kg Cushion 17" Cushion 17" 136 kg Cushion 17" Cushion 17" Cushion 17"	Overall Depth x 21" Overall Depth x 22" Overall Depth x 14" Overall Depth x 15" Overall Depth x 15" Overall Depth	535 mm 560 mm 355 mm 380 mm	Overall Width Overall Width Overall Width	405 mm 405 mm 430 mm 430 mm
JFUSION1621 Max User Weight JFUSION1622 Max User Weight JFUSION1714 Max User Weight JFUSION1715 Max User Weight	Jay Fusion Jay Fusion Jay Fusion	136 kg Cushion 16" 136 kg Cushion 16" 136 kg Cushion 17" 136 kg Cushion 17" 136 kg	Overall Depth x 21" Overall Depth x 22" Overall Depth x 14" Overall Depth x 15" Overall Depth	535 mm 560 mm 355 mm	Overall Width Overall Width Overall Width	405 mm 405 mm 430 mm
JFUSION1621 Max User Weight JFUSION1622 Max User Weight JFUSION1714 Max User Weight JFUSION1715 Max User Weight JFUSION1716	Jay Fusion Jay Fusion Jay Fusion Jay Fusion	136 kg Cushion 16" 136 kg Cushion 16" 136 kg Cushion 17" 136 kg Cushion 17" 136 kg Cushion 17" 136 kg Cushion 17" Cushion 17" 136 kg Cushion 17" Cushion 17" Cushion 17"	Overall Depth x 21" Overall Depth x 22" Overall Depth x 14" Overall Depth x 15" Overall Depth x 15" Overall Depth x 16" Overall Depth	535 mm 560 mm 355 mm 380 mm	Overall Width Overall Width Overall Width	405 mm 405 mm 430 mm 430 mm
JFUSION1621 Max User Weight JFUSION1622 Max User Weight JFUSION1714 Max User Weight JFUSION1715 Max User Weight JFUSION1716 Max User Weight	Jay Fusion Jay Fusion Jay Fusion Jay Fusion	136 kg Cushion 16" 136 kg Cushion 16" 136 kg Cushion 17" 136 kg	Overall Depth x 21" Overall Depth x 22" Overall Depth x 14" Overall Depth x 15" Overall Depth x 15" Overall Depth x 16" Overall Depth	535 mm 560 mm 355 mm 380 mm	Overall Width Overall Width Overall Width	405 mm 405 mm 430 mm 430 mm
JFUSION1621 Max User Weight JFUSION1622 Max User Weight JFUSION1714 Max User Weight JFUSION1715 Max User Weight JFUSION1716 Max User Weight JFUSION1717	Jay Fusion Jay Fusion Jay Fusion Jay Fusion Jay Fusion	136 kg Cushion 16" 136 kg Cushion 16" 136 kg Cushion 17" Cushion 17" Cushion 17" Cushion 17" Cushion 17"	Overall Depth x 21" Overall Depth x 22" Overall Depth x 14" Overall Depth x 15" Overall Depth x 15" Overall Depth x 16" Overall Depth x 17" Overall Depth	535 mm 560 mm 355 mm 380 mm 405 mm	Overall Width Overall Width Overall Width Overall Width Overall Width	405 mm 405 mm 430 mm 430 mm

JFUSION1719	Jay Fusion	Cushion 17"	x 19″
Max User Weight	-	136 kg	Overall Depth
JFUSION1720		Quebien 17	× 20"
Max User Weight	Jay Fusion	Cushion 17"	X 20" Overall Depth
Max user weight		150 Kg	
JFUSION1721	Jay Fusion	Cushion 17" x	x 21″
Max User Weight		136 kg	Overall Depth
IFUSION1722	lay Eucion	Cushion 17"	v 22"
Max User Weight	Jay Tusion	136 kg	Overall Depth
_			
JFUSION1723	Jay Fusion	Cushion 17"	x 23"
Max User Weight		136 kg	Overall Depth
JFUSION1814	lay Fusion	Cushion 18"	х 1 4 ″
Max User Weight	Juy ruoion	136 kg	Overall Depth
JFUSION1815	Jay Fusion	Cushion 18"	x 15″
Max User Weight		136 kg	Overall Depth
JFUSION1816	lav Fusion	Cushion 18"	x 16″
Max User Weight	Jujiuolon	136 kg	Overall Depth
		-	
JFUSION1817	Jay Fusion	Cushion 18"	x 17"
Max User Weight		136 kg	Overall Depth
JFUSION1818	Jay Fusion	Cushion 18"	x 18″
Max User Weight	-	136 kg	Overall Depth
JFUSION1819	Jay Fusion	Cushion 18"	
Max User Weight		136 kg	Overall Depth
JFUSION1820	Jay Fusion	Cushion 18"	x 20"
Max User Weight		136 kg	Overall Depth
JFUSION1821		Cuchien 10	v 01/
Max User Weight	Jay Fusion	Cushion 18"	X 2 1 Overall Depth
max user weight		100 kg	
JFUSION1822	Jay Fusion	Cushion 18"	x 22"
Max User Weight		136 kg	Overall Depth
JFUSION1823		Cuchion 10"	v 00"
Max User Weight	Jay Fusion	Cushion 18"	X 23 Overall Depth
max osci weight		100 kg	overall beptil
JFUSION1824	Jay Fusion	Cushion 18"	x 24"
Max User Weight		136 kg	Overall Depth
JFUSION1915	lay Euclos	Cushion 19"	v 15"
Max User Weight	Jay rusion	136 kg	Overall Depth
			c.c./uii septi
JFUSION1916	Jay Fusion	Cushion 19"	x 16″
Max User Weight		136 kg	Overall Depth

480 mm	Overall Width	430 mm
505 mm	Overall Width	430 mm
_		
535 mm	Overall Width	430 mm
560 mm	Overall Width	430 mm
585 mm	Overall Width	430 mm
355 mm	Overall Width	455 mm
380 mm	Overall Width	455 mm
405 mm	Overall Width	455 mm
430 mm	Overall Width	455 mm
-		
455 mm	Overall Width	455 mm
480 mm	Overall Width	455 mm
505 mm	Overall Width	455 mm
535 mm	Overall Width	455 mm
560 mm	Overall Width	455 mm
585 mm	Overall Width	455 mm
610 mm	Overall Width	455 mm
380 mm	Overall Width	480 mm
405 mm	Overall Width	480 mm

Country Care

ution. Local Service

JFUSION1917	lay Fusion	Cushion 19"	v 17"			
Max User Weight	Jay Tusion	136 kg	Overall Depth	430 mm	Overall Width	480 mm
	Jay Fusion	Cushion 19"				
Max User Weight		136 kg	Overall Depth	455 mm	Overall Width	480 mm
JFUSION1919	Jay Fusion	Cushion 19"	x 19″			
Max User Weight	-	136 kg	Overall Depth	480 mm	Overall Width	480 mm
151101011020	1	0				
JFUSION1920 Max User Weight	Jay Fusion	Cushion 19"	X 20" Overall Depth	505 mm	Overall Width	480 mm
Max user weight		150 kg		505 1111		400 11111
JFUSION1921	Jay Fusion	Cushion 19"	x 21″			
Max User Weight		136 kg	Overall Depth	535 mm	Overall Width	480 mm
JFUSION1922	lay Fusion	Cushion 19"	x 22"			
Max User Weight	Juy Tubion	136 kg	Overall Depth	560 mm	Overall Width	480 mm
JFUSION1923	Jay Fusion		x 23″			
Max User Weight		136 kg	Overall Depth	585 mm	Overall Width	480 mm
JFUSION1924	Jay Fusion	Cushion 19"	x 24"			
Max User Weight	-	136 kg	Overall Depth	610 mm	Overall Width	480 mm
151101010015						
JFUSION2015	Jay Fusion	Cushion 20"		200	Overall Middle	F0F mm
Max User Weight		136 kg	Overall Depth	380 mm	Overall Width	505 mm
JFUSION2016	Jay Fusion	Cushion 20"	x 16″			
Max User Weight		136 kg	Overall Depth	405 mm	Overall Width	505 mm
JFUSION2017	lay Fusion	Cushion 20"	v 17"			
Max User Weight	Jayrasion	136 kg	Overall Depth	430 mm	Overall Width	505 mm
		-				
JFUSION2018	Jay Fusion	Cushion 20"	x 18″			
Max User Weight		136 kg	Overall Depth	455 mm	Overall Width	505 mm
JFUSION2019	Jay Fusion	Cushion 20"	x 19"			
Max User Weight	-	136 kg	Overall Depth	480 mm	Overall Width	505 mm
151101012020	lev Freier	Quebien 00	" x 00"			
JFUSION2020 Max User Weight	Jay Fusior	136 kg	Verall Depth	505 mm	Overall Width	505 mm
wax user weight		130 KY		505 11111		505 11111
JFUSION2021	Jay Fusion	Cushion 20"	x 21″			
Max User Weight		136 kg	Overall Depth	535 mm	Overall Width	505 mm
JFUSION2022	lay Fusion	n Cushion 20'	' x 22"			
Max User Weight	549 1 43101	136 kg	Overall Depth	560 mm	Overall Width	505 mm
		-				
JFUSION2023	Jay Fusior					
Max User Weight		136 kg	Overall Depth	585 mm	Overall Width	505 mm
JFUSION2024	Jay Fusion	Cushion 20'	' x 24"			
Max User Weight	-	136 kg	Overall Depth	610 mm	Overall Width	505 mm

JFUSION2117	Jav Fusion	Cushion 21"	x 17″
Max User Weight	, , , , , , , , , ,	136 kg	Overall Depth
JFUSION2118	Jay Fusion	Cushion 21"	
Max User Weight	_	136 kg	Overall Depth
JFUSION2119	Jay Fusion	Cushion 21"	x 19″
Max User Weight		136 kg	Overall Depth
1FUSION2120		Cushion 21"	v 20"
Max User Weight	Jay Fusion	136 kg	Overall Depth
JFUSION2121	Jay Fusion	Cushion 21"	x 21"
Max User Weight		136 kg	Overall Depth
JFUSION2122	Jay Fusion	Cushion 21"	x 22"
Max User Weight	-	136 kg	Overall Depth
JFUSION2123	Jay Fusion	Cushion 21"	
Max User Weight		136 kg	Overall Depth
JFUSION2124	Jay Fusion	Cushion 21"	x 24"
Max User Weight	-	136 kg	Overall Depth
151101010017			
JFUSION2217	Jay Fusion	Cushion 22"	
Max User Weight		225 kg	Overall Depth
JFUSION2218	Jay Fusion	Cushion 22"	x 18″
Max User Weight		225 kg	Overall Depth
JFUSION2219	lay Fusion	Cushion 22"	v 19"
Max User Weight	Jay Tusion	225 kg	Overall Depth
		-	
JFUSION2220	Jay Fusior	Cushion 22	" x 20"
Max User Weight		225 kg	Overall Depth
JFUSION2221	Jay Fusion	Cushion 22"	x 21″
Max User Weight		225 kg	Overall Depth
	lev E. J		
JFUSION2222	Jay Fusior	Cushion 22	
Max User Weight		225 kg	Overall Depth
JFUSION2223	Jay Fusior	Cushion 22	" x 23"
Max User Weight		225 kg	Overall Depth
161101012224		Cuchion 20	" v 24"
JFUSION2224 Max User Weight	Jay Fusior	225 kg	0verall Depth
max oser weight		LEO NG	oreitan bepan
JFUSION2319	Jay Fusion	Cushion 23"	x 19"
Max User Weight		225 kg	Overall Depth
JFUSION2320	lay Fusion	1 Cushion 23	" x 20"
Max User Weight	Jug 1 03101	225 kg	Overall Depth
		-	

430 mm	Overall Width	535 mm
455 mm	Overall Width	535 mm
480 mm	Overall Width	535 mm
_		
505 mm	Overall Width	535 mm
535 mm	Overall Width	535 mm
560 mm	Overall Width	535 mm
585 mm	Overall Width	535 mm
610 mm	Overall Width	535 mm
430 mm	Overall Width	560 mm
_		
455 mm	Overall Width	560 mm
480 mm	Overall Width	560 mm
505 mm	Overall Width	560 mm
535 mm	Overall Width	560 mm
560 mm	Overall Width	560 mm
585 mm	Overall Width	560 mm
610 mm	Overall Width	650 mm
480 mm	Overall Width	585 mm
505 mm	Overall Width	585 mm

Country Care

oution. Local Service

JFUSION2321 Jay Fusio	on Cushion 23"	' x 21"			
Max User Weight	225 kg	Overall Depth	535 mm	Overall Width	585 mm
JFUSION2322 Jay Fusi	on Cushion 23	" x 22"			
Max User Weight	225 kg	Overall Depth	560 mm	Overall Width	585 mm
JFUSION2323 Jay Fusi	on Cushion 23	" x 23"	_		
Max User Weight	225 kg	Overall Depth	585 mm	Overall Width	585 mm
JFUSION2324 Jay Fusi	on Cushion 23'	" x 24"			
Max User Weight	225 kg	Overall Depth	610 mm	Overall Width	585 mm
JFUSION2419 Jay Fusio	on Cushion 24"	x 19″			
Max User Weight	225 kg	Overall Depth	480 mm	Overall Width	610 mm
JFUSION2420 Jay Fusi	on Cushion 24'	" x 20"			
Max User Weight	225 kg	Overall Depth	505 mm	Overall Width	610 mm
JFUSION2421 Jay Fusio	on Cushion 24"	x 21"			
Max User Weight	225 kg	Overall Depth	535 mm	Overall Width	610 mm
JFUSION2422 Jay Fusi	on Cushion 24'	' x 22"			-
Max User Weight	225 kg	Overall Depth	560 mm	Overall Width	610 mm
JFUSION2423 Jay Fusi	on Cushion 24'	' x 23"			
Max User Weight	225 kg	Overall Depth	585 mm	Overall Width	610 mm
JFUSION2424 Jay Fusio	on Cushion 24'	' x 24"			
Max User Weight	225 kg	Overall Depth	610 mm	Overall Width	610 mm

Jay X2 Cushion



• A moderate to high pressure relieving cushion

- With a softer foam base and a fluid bladder
- It is lightweight and maintenance free A reduced profile that is designed to provide stability and pressure relief to active clients

······	14" x 14"				
Max User Weight	113 kg	Overall Depth	355 mm	Overall Width	355 mm
JX21416 Jay X2 Cushion	14" x 16"				
Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	355 mm
IV91515 Jaw VO Overhiere	15" 15"				
JX21515 Jay X2 Cushion	15" X 15"				
Max User Weight	15" X 15" 113 kg	Overall Depth	380 mm	Overall Width	380 mm
-	113 kg	Overall Depth	380 mm	Overall Width	380 mm

JX21616 Jay X2 Cushion	16" x 16"	
Max User Weight	113 kg	Overall Depth
JX21617 Jay X2 Cushion	16" x 17"	
Max User Weight	113 kg	Overall Depth
JX21618 Jay X2 Cushion	16" x 18"	
Max User Weight	113 kg	Overall Depth
JX21620 Jay X2 Cushion	16" x 20"	
Max User Weight	113 kg	Overall Depth
JX21716 Jay X2 Cushion	17" x 16"	
Max User Weight	113 kg	Overall Depth
JX21717 Jay X2 Cushion	17" x 17"	
Max User Weight	113 kg	Overall Depth
JX21816 Jay X2 Cushion	18" x 16"	
Max User Weight	113 kg	Overall Depth
JX21818 Jay X2 Cushion	18" x 18"	
Max User Weight	113 kg	Overall Depth
JX21820 Jay X2 Cushion	18" x 20"	
Max User Weight	113 kg	Overall Depth
JX22020 Jay X2 Cushior	1 20" x 20"	
Max User Weight	113 kg	Overall Depth

Jay Union Cushion Standard Sizes



• A versatile and comfortable cushion • Protects the skin and allows patient positioning Composed of a dynamic fluid and foam layering system • A moisture resistant inner cover, and an anti-microbial outer cover • The foam base is moderately contoured to encourage proper positioning of the pelvis and thighs • Features a top layer of soft, visco memory foam, with an extra layer covering the pelvic loading area • The ultimate wheelchair seating surface that distributes weight for pressure reduction and maximum comfort • The fluid insert that lies beneath the visco foam conforms to each individual's unique contours for a reduced risk of skin breakdown

JU1816 Jay Union Cu	ushion S	tandard 18" x 16"		
Max User Weight	136 kg	Overall Depth	405 mm	0
JU1817 Jay Union Cu	ushion S	tandard 18" x 17"		
Max User Weight	136 kg	Overall Depth	430 mm	0
JU1818 Jay Union Cu	ushion S	tandard 18" x 18"		
Max User Weight	136 kg	Overall Depth	455 mm	0
JU1819 Jay Union Cu	ushion S	tandard 18" x 19"		
Max User Weight	136 kg	Overall Depth	480 mm	0

405 mm	Overall Width	405 mm
430 mm	Overall Width	405 mm
155	•	105
455 mm	Overall Width	405 mm
505 mm	Overall Width	405 mm
405 mm	Overall Width	430 mm
430 mm	Overall Width	430 mm
105		155
405 mm	Overall Width	455 mm
455 mm	Overall Width	455 mm
505 mm	Overall Width	455 mm
505 mm	Overall Width	505 mm

Overall Width	455 mm	Product Weight	2 kg
Overall Width	455 mm	Product Weight	2 kg
Overall Width	455 mm	Product Weight	2 kg
			,
Ourse II Wildel	455	Due due to We lie ha	0.1-1
Overall Width	455 mm	Product Weight	2 kg

Country Care GROUP

1111820 lav Unio	n Cushion	Standard 18" x 20"	,				
Max User Weight	136 kg	Overall Depth	505 mm	Overall Width	455 mm	Product Weight	2 kg
U1001 1	0	0					
Max User Weight	n Cushion 136 kg	Standard 18" x 21" Overall Depth	535 mm	Overall Width	455 mm	Product Weight	2 kg
max user weight	150 kg		333 1111		455 1111	Fibuuct weight	2 Ng
U1822 Jay Unio	n Cushion	Standard 18" x 22"					
Max User Weight	136 kg	Overall Depth	560 mm	Overall Width	455 mm	Product Weight	2 kg
U1823 Jay Unio	n Cushion	Standard 18" x 23"	,				
Max User Weight	136 kg	Overall Depth	585 mm	Overall Width	455 mm	Product Weight	2 kg
111824 Jay Unio	n Cuchion	Standard 18" x 24"					
Max User Weight	136 kg	Overall Depth	610 mm	Overall Width	455 mm	Product Weight	2 kg
-		Standard 19" x 15"					
Max User Weight	136 kg	Overall Depth	380 mm	Overall Width	480 mm	Product Weight	2 kg
U1916 Jay Unio	n Cushion	Standard 19" x 16"					
Max User Weight	136 kg	Overall Depth	405 mm	Overall Width	480 mm	Product Weight	2 kg
U1917 Jay Unior	n Cuebion	Standard 19" x 17"					
Max User Weight	136 kg	Overall Depth	430 mm	Overall Width	480 mm	Product Weight	2 kg
-		Standard 19" x 18"					
Max User Weight	136 kg	Overall Depth	455 mm	Overall Width	480 mm	Product Weight	2 kg
U1919 Jay Unio	n Cushion	Standard 19" x 19"					
Max User Weight	136 kg	Overall Depth	480 mm	Overall Width	480 mm	Product Weight	2 kg
U1920 Jav Unio	n Cushion	Standard 19" x 20"	,				
Max User Weight	136 kg	Overall Depth	505 mm	Overall Width	480 mm	Product Weight	2 kg
111021 Joy Unio		Ctonderd 10" v 21"	_				
Max User Weight	136 kg	Standard 19" x 21" Overall Depth	535 mm	Overall Width	480 mm	Product Weight	2 kg
	100 kg	orcian Sepan			100 1111	rioudor incigit	Ling
U1922 Jay Unio	n Cushion	Standard 19" x 22"					
Max User Weight	136 kg	Overall Depth	560 mm	Overall Width	480 mm	Product Weight	2 kg
U1923 Jay Unio	n Cushion	Standard 19" x 23"	,				
Max User Weight	136 kg	Overall Depth	585 mm	Overall Width	480 mm	Product Weight	2 kg
11192/ Joy Unio	n Cuchion	Standard 19" x 24"					
Max User Weight	136 kg	Overall Depth	610 mm	Overall Width	480 mm	Product Weight	2 kg
		Standard 20" x 15"					
Max User Weight	136 kg	Overall Depth	380 mm	Overall Width	505 mm	Product Weight	2 kg
IU2016 Jay Unio	n Cushion	Standard 20" x 16"					
Max User Weight	136 kg	Overall Depth	405 mm	Overall Width	505 mm	Product Weight	2 kg
IU2017 Jav Unio	n Cuehion	Standard 20" x 17"					
Max User Weight	136 kg	Overall Depth	430 mm	Overall Width	505 mm	Product Weight	2 kg
							9

	.	Standard 20" x 18"					
U2018 Jay Unior	n Cushion						
Max User Weight	136 kg	Overall Depth	455 mm	Overall Width	505 mm	Product Weight	2 kg
10010				·			
-		Standard 20" x 19"					
Max User Weight	136 kg	Overall Depth	480 mm	Overall Width	505 mm	Product Weight	2 kg
U2020 Jay Unio	n Cushion	n Standard 20" x 20	"				
Max User Weight	136 kg	Overall Depth	505 mm	Overall Width	505 mm	Product Weight	2 kg
-		Standard 20" x 21"					
Max User Weight	136 kg	Overall Depth	535 mm	Overall Width	505 mm	Product Weight	2 kg
U2022 Jay Unio	n Cushion	Standard 20" x 22	"				
Max User Weight	136 kg	Overall Depth	560 mm	Overall Width	505 mm	Product Weight	2 kg
-		Standard 20" x 23			F05 mm	Droduct Weight	0.44
Max User Weight	136 kg	Overall Depth	585 mm	Overall Width	505 mm	Product Weight	2 kg
U2024 Jay Unio	n Cushion	Standard 20" x 24	n				
Max User Weight	136 kg	Overall Depth	610 mm	Overall Width	505 mm	Product Weight	2 kg
U1510	·	0					
-	Cushion	Standard 15" x 18"	_				
	100 km	Overall Danth	455 mm	Overell Width		Deaduct Mainht	
Max User Weight	136 kg	Overall Depth	455 mm	Overall Width	380 mm	Product Weight	2 kg
		Overall Depth Standard 15" x 19"	455 mm	Overall Width	380 mm	Product Weight	2 Kg
			455 mm	Overall Width Overall Width	380 mm 380 mm	Product Weight Product Weight	2 kg
U1519 Jay Union Max User Weight	136 kg	Standard 15" x 19" Overall Depth	480 mm				
U1519 Jay Union Max User Weight U1520 Jay Unior	n Cushion 136 kg	Standard 15" x 19" Overall Depth Standard 15" x 20"	480 mm	Overall Width	380 mm	Product Weight	2 kg
U1519 Jay Union Max User Weight	136 kg	Standard 15" x 19" Overall Depth	480 mm				
U1519 Jay Union Max User Weight U1520 Jay Unior Max User Weight	136 kg 136 kg 136 kg 136 kg	Standard 15" x 19" Overall Depth Standard 15" x 20"	480 mm	Overall Width	380 mm	Product Weight	2 kg
U1519 Jay Union Max User Weight U1520 Jay Unior Max User Weight	136 kg 136 kg 136 kg 136 kg	Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth	480 mm	Overall Width	380 mm	Product Weight	2 kg
U1519 Jay Union Max User Weight U1520 Jay Unior Max User Weight U1521 Jay Union Max User Weight	n Cushion 136 kg n Cushion 136 kg n Cushion 136 kg	Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth Standard 15" x 21" Overall Depth	480 mm	Overall Width Overall Width	380 mm 380 mm	Product Weight Product Weight	2 kg 2 kg
U1519 Jay Union Max User Weight U1520 Jay Union Max User Weight U1521 Jay Union Max User Weight U1614 Jay Union	n Cushion 136 kg n Cushion 136 kg n Cushion 136 kg	Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth Standard 15" x 21" Overall Depth Standard 16" x 14"	480 mm 505 mm 535 mm	Overall Width Overall Width Overall Width	380 mm 380 mm 380 mm	Product Weight Product Weight Product Weight	2 kg 2 kg 2 kg
U1519 Jay Union Max User Weight U1520 Jay Unior Max User Weight U1521 Jay Union Max User Weight	n Cushion 136 kg 136 kg 136 kg 136 kg 136 kg	Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth Standard 15" x 21" Overall Depth	480 mm	Overall Width Overall Width	380 mm 380 mm	Product Weight Product Weight	2 kg 2 kg
U1519 Jay Union Max User Weight U1520 Jay Union Max User Weight U1521 Jay Union Max User Weight U1614 Jay Union Max User Weight	n Cushion 136 kg n Cushion 136 kg n Cushion 136 kg n Cushion 136 kg	Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth Standard 15" x 21" Overall Depth Standard 16" x 14"	480 mm 505 mm 535 mm	Overall Width Overall Width Overall Width	380 mm 380 mm 380 mm	Product Weight Product Weight Product Weight	2 kg 2 kg 2 kg
U1519 Jay Union Max User Weight U1520 Jay Union Max User Weight U1521 Jay Union Max User Weight U1614 Jay Union Max User Weight	n Cushion 136 kg n Cushion 136 kg n Cushion 136 kg n Cushion 136 kg	Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth Standard 15" x 21" Overall Depth Standard 16" x 14" Overall Depth	480 mm 505 mm 535 mm	Overall Width Overall Width Overall Width	380 mm 380 mm 380 mm	Product Weight Product Weight Product Weight	2 kg 2 kg 2 kg
U1519 Jay Union Max User Weight U1520 Jay Union Max User Weight U1521 Jay Union Max User Weight U1614 Jay Union Max User Weight U1615 Jay Union Max User Weight	Cushion I36 kg	Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth Standard 15" x 21" Overall Depth Standard 16" x 14" Overall Depth Standard 16" x 15" Overall Depth	480 mm 505 mm 535 mm 355 mm	Overall Width Overall Width Overall Width Overall Width	380 mm 380 mm 380 mm 405 mm	Product Weight Product Weight Product Weight Product Weight	2 kg 2 kg 2 kg 2 kg
U1519 Jay Union Max User Weight U1520 Jay Union Max User Weight U1521 Jay Union Max User Weight U1614 Jay Union Max User Weight U1615 Jay Union Max User Weight U1616 Jay Union	Cushion 136 kg 136 kg 136 kg Cushion 136 kg Cushion 136 kg 136 kg 136 kg 136 kg 136 kg 136 kg	Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth Standard 15" x 21" Overall Depth Standard 16" x 14" Overall Depth Standard 16" x 15" Overall Depth Standard 16" x 16"	480 mm 505 mm 535 mm 355 mm 380 mm	Overall Width Overall Width Overall Width Overall Width Overall Width	380 mm 380 mm 380 mm 380 mm 405 mm	Product Weight Product Weight Product Weight Product Weight Product Weight	2 kg 2 kg 2 kg 2 kg 2 kg
U1519 Jay Union Max User Weight U1520 Jay Union Max User Weight U1521 Jay Union Max User Weight U1614 Jay Union Max User Weight U1615 Jay Union Max User Weight	Cushion I36 kg	Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth Standard 15" x 21" Overall Depth Standard 16" x 14" Overall Depth Standard 16" x 15" Overall Depth	480 mm 505 mm 535 mm 355 mm	Overall Width Overall Width Overall Width Overall Width	380 mm 380 mm 380 mm 405 mm	Product Weight Product Weight Product Weight Product Weight	2 kg 2 kg 2 kg 2 kg
U1519 Jay Union Max User Weight U1520 Jay Union Max User Weight U1521 Jay Union Max User Weight U1614 Jay Union Max User Weight U1615 Jay Union Max User Weight U1616 Jay Union Max User Weight	Cushion 136 kg 136 kg 136 kg Cushion 136 kg Cushion 136 kg Cushion 136 kg 136 kg Cushion 136 kg	Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth Standard 15" x 21" Overall Depth Standard 16" x 14" Overall Depth Standard 16" x 15" Overall Depth Standard 16" x 16"	480 mm 505 mm 535 mm 355 mm 380 mm	Overall Width Overall Width Overall Width Overall Width Overall Width	380 mm 380 mm 380 mm 380 mm 405 mm	Product Weight Product Weight Product Weight Product Weight Product Weight	2 kg 2 kg 2 kg 2 kg 2 kg
U1519 Jay Union Max User Weight U1520 Jay Union Max User Weight U1521 Jay Union Max User Weight U1614 Jay Union Max User Weight U1615 Jay Union Max User Weight U1616 Jay Union Max User Weight	Cushion 136 kg 136 kg 136 kg Cushion 136 kg Cushion 136 kg Cushion 136 kg 136 kg Cushion 136 kg	Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth Standard 15" x 21" Overall Depth Standard 16" x 14" Overall Depth Standard 16" x 15" Overall Depth Standard 16" x 16" Overall Depth	480 mm 505 mm 535 mm 355 mm 380 mm	Overall Width Overall Width Overall Width Overall Width Overall Width	380 mm 380 mm 380 mm 380 mm 405 mm	Product Weight Product Weight Product Weight Product Weight Product Weight	2 kg 2 kg 2 kg 2 kg 2 kg
U1519 Jay Union Max User Weight U1520 Jay Union Max User Weight U1521 Jay Union Max User Weight U1614 Jay Union Max User Weight U1615 Jay Union Max User Weight U1616 Jay Union Max User Weight U1617 Jay Union Max User Weight		Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth Standard 15" x 21" Overall Depth Standard 16" x 14" Overall Depth Standard 16" x 15" Overall Depth Standard 16" x 16" Overall Depth Standard 16" x 17" Overall Depth	 480 mm 505 mm 535 mm 3355 mm 380 mm 405 mm 	Overall Width	380 mm 380 mm 380 mm 380 mm 405 mm 405 mm	Product Weight Product Weight Product Weight Product Weight Product Weight Product Weight	2 kg 2 kg 2 kg 2 kg 2 kg 2 kg 2 kg 2 kg
U1519 Jay Union Max User Weight U1520 Jay Union Max User Weight U1521 Jay Union Max User Weight U1614 Jay Union Max User Weight U1615 Jay Union Max User Weight U1616 Jay Union Max User Weight U1617 Jay Union Max User Weight U1618 Jay Union	Cushion 136 kg 136 kg 136 kg Cushion 136 kg	Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth Standard 15" x 21" Overall Depth Standard 16" x 14" Overall Depth Standard 16" x 15" Overall Depth Standard 16" x 16" Overall Depth Standard 16" x 17" Overall Depth Standard 16" x 17" Overall Depth	 480 mm 505 mm 535 mm 3355 mm 355 mm 380 mm 405 mm 430 mm 	Overall Width	380 mm 380 mm 380 mm 380 mm 405 mm 405 mm 405 mm	Product Weight	2 kg 2 kg 2 kg 2 kg 2 kg 2 kg 2 kg 2 kg
U1519 Jay Union Max User Weight U1520 Jay Union Max User Weight U1521 Jay Union Max User Weight U1614 Jay Union Max User Weight U1615 Jay Union Max User Weight U1616 Jay Union Max User Weight U1617 Jay Union Max User Weight	Cushion 136 kg	Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth Standard 15" x 21" Overall Depth Standard 16" x 14" Overall Depth Standard 16" x 15" Overall Depth Standard 16" x 16" Overall Depth Standard 16" x 17" Overall Depth	 480 mm 505 mm 535 mm 3355 mm 380 mm 405 mm 	Overall Width	380 mm 380 mm 380 mm 380 mm 405 mm 405 mm	Product Weight Product Weight Product Weight Product Weight Product Weight Product Weight	2 kg 2 kg 2 kg 2 kg 2 kg 2 kg 2 kg 2 kg
U1519 Jay Union Max User Weight U1520 Jay Union Max User Weight U1521 Jay Union Max User Weight U1614 Jay Union Max User Weight U1615 Jay Union Max User Weight U1616 Jay Union Max User Weight U1617 Jay Union Max User Weight U1618 Jay Union Max User Weight	Cushion 136 kg 136 kg 136 kg 136 kg Cushion 136 kg Cushion 136 kg Cushion 136 kg Cushion 136 kg	Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth Standard 15" x 21" Overall Depth Standard 16" x 14" Overall Depth Standard 16" x 15" Overall Depth Standard 16" x 16" Overall Depth Standard 16" x 17" Overall Depth Standard 16" x 17" Overall Depth	 480 mm 505 mm 535 mm 3355 mm 355 mm 380 mm 405 mm 430 mm 	Overall Width	380 mm 380 mm 380 mm 380 mm 405 mm 405 mm 405 mm	Product Weight	2 kg 2 kg 2 kg 2 kg 2 kg 2 kg 2 kg 2 kg
U1519 Jay Union Max User Weight U1520 Jay Union Max User Weight U1521 Jay Union Max User Weight U1614 Jay Union Max User Weight U1615 Jay Union Max User Weight U1616 Jay Union Max User Weight U1617 Jay Union Max User Weight U1618 Jay Union Max User Weight	Cushion 136 kg 136 kg 136 kg 136 kg Cushion 136 kg Cushion 136 kg Cushion 136 kg Cushion 136 kg	Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth Standard 15" x 21" Overall Depth Standard 16" x 14" Overall Depth Standard 16" x 15" Overall Depth Standard 16" x 16" Overall Depth Standard 16" x 17" Overall Depth Standard 16" x 18" Overall Depth	 480 mm 505 mm 535 mm 3355 mm 355 mm 380 mm 405 mm 430 mm 	Overall Width	380 mm 380 mm 380 mm 380 mm 405 mm 405 mm 405 mm	Product Weight	2 kg 2 kg 2 kg 2 kg 2 kg 2 kg 2 kg 2 kg
U1519 Jay Union Max User Weight U1520 Jay Union Max User Weight U1521 Jay Union Max User Weight U1614 Jay Union Max User Weight U1615 Jay Union Max User Weight U1616 Jay Union Max User Weight U1617 Jay Union Max User Weight U1618 Jay Union Max User Weight U1619 Jay Union Max User Weight	Image: Current of the second of the secon	Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth Standard 15" x 21" Overall Depth Standard 16" x 14" Overall Depth Standard 16" x 15" Overall Depth Standard 16" x 16" Overall Depth Standard 16" x 17" Overall Depth Standard 16" x 18" Overall Depth Standard 16" x 18" Overall Depth Standard 16" x 19" Overall Depth	 480 mm 505 mm 535 mm 3355 mm 355 mm 380 mm 405 mm 430 mm 430 mm 	Overall Width Overall Width	380 mm 380 mm 380 mm 380 mm 405 mm 405 mm 405 mm 405 mm	Product Weight Product Weight	2 kg 2 kg 2 kg 2 kg 2 kg 2 kg 2 kg 2 kg
U1519 Jay Union Max User Weight U1520 Jay Union Max User Weight U1521 Jay Union Max User Weight U1614 Jay Union Max User Weight U1615 Jay Union Max User Weight U1616 Jay Union Max User Weight U1617 Jay Union Max User Weight U1618 Jay Union Max User Weight U1619 Jay Union Max User Weight	Image: Current of the second of the secon	Standard 15" x 19" Overall Depth Standard 15" x 20" Overall Depth Standard 15" x 21" Overall Depth Standard 16" x 14" Overall Depth Standard 16" x 15" Overall Depth Standard 16" x 16" Overall Depth Standard 16" x 17" Overall Depth Standard 16" x 18" Overall Depth Standard 16" x 18" Overall Depth	 480 mm 505 mm 535 mm 3355 mm 355 mm 380 mm 405 mm 430 mm 430 mm 	Overall Width Overall Width	380 mm 380 mm 380 mm 380 mm 405 mm 405 mm 405 mm 405 mm	Product Weight Product Weight	2 kg 2 kg 2 kg 2 kg 2 kg 2 kg 2 kg 2 kg



JU1622 Jay Union Cu Max User Weight JU1714 Jay Union Cus Max User Weight JU1715 Jay Union Cus Max User Weight JU1716 Jay Union Cus Max User Weight JU1717 Jay Union Cus Max User Weight JU1718 Jay Union Cus Max User Weight	136 kg shion S1 136 kg shion S1 136 kg shion S1 136 kg shion S1 136 kg	Overall Depth tandard 17" x 14" Overall Depth tandard 17" x 15" Overall Depth tandard 17" x 16" Overall Depth tandard 17" x 17" Overall Depth	535 mm 560 mm 355 mm 380 mm 405 mm	Overall Width Overall Width Overall Width Overall Width Overall Width	405 mm 405 mm 430 mm 430 mm	Product Weight Product Weight Product Weight Product Weight Product Weight	2 kg 2 kg 2 kg 2 kg 2 kg 2 kg
Max User Weight IU1714 Jay Union Cus Max User Weight IU1715 Jay Union Cus Max User Weight IU1715 Jay Union Cus Max User Weight IU1716 Jay Union Cus Max User Weight IU1717 Jay Union Cus Max User Weight IU1718 Jay Union Cus Max User Weight IU1718 Jay Union Cus Max User Weight IU1719 Jay Union Cus	136 kg shion S1 136 kg shion S1 136 kg shion S1 136 kg shion S1 136 kg Shion S1 136 kg	Overall Depth tandard 17" x 14" Overall Depth tandard 17" x 15" Overall Depth tandard 17" x 16" Overall Depth tandard 17" x 17" Overall Depth	355 mm 380 mm 405 mm	Overall Width Overall Width Overall Width	430 mm 430 mm	Product Weight Product Weight	2 kg 2 kg
U1714 Jay Union Cus Max User Weight U1715 Jay Union Cus Max User Weight U1716 Jay Union Cus Max User Weight U1717 Jay Union Cus Max User Weight U1718 Jay Union Cus Max User Weight U1719 Jay Union Cus	shion S1 136 kg shion S1 136 kg shion S1 136 kg shion S1 136 kg	tandard 17" x 14" Overall Depth tandard 17" x 15" Overall Depth tandard 17" x 16" Overall Depth tandard 17" x 17" Overall Depth	355 mm 380 mm 405 mm	Overall Width Overall Width Overall Width	430 mm 430 mm	Product Weight Product Weight	2 kg 2 kg
Max User Weight U1715 Jay Union Cus Max User Weight U1716 Jay Union Cus Max User Weight U1717 Jay Union Cus Max User Weight U1718 Jay Union Cus Max User Weight U1718 Jay Union Cus Max User Weight U1719 Jay Union Cus	136 kg shion S1 136 kg shion S1 136 kg shion S1 136 kg shion S1	Overall Depth tandard 17" x 15" Overall Depth tandard 17" x 16" Overall Depth tandard 17" x 17" Overall Depth	380 mm 405 mm	Overall Width Overall Width	430 mm	Product Weight	2 kg
U1715 Jay Union Custom Max User Weight Image: State Sta	shion S1 136 kg shion S1 136 kg shion S1 136 kg shion S1	tandard 17" x 15" Overall Depth tandard 17" x 16" Overall Depth tandard 17" x 17" Overall Depth	380 mm 405 mm	Overall Width Overall Width	430 mm	Product Weight	2 kg
Max User Weight	136 kg shion St 136 kg shion St 136 kg shion St	Overall Depth tandard 17" x 16" Overall Depth tandard 17" x 17" Overall Depth	405 mm	Overall Width			
U1716 Jay Union Cus Max User Weight U1717 Jay Union Cus Max User Weight U1718 Jay Union Cus Max User Weight	shion S1 136 kg shion St 136 kg shion S1	tandard 17" x 16" Overall Depth candard 17" x 17" Overall Depth	405 mm	Overall Width			
Max User Weight	136 kg shion St 136 kg shion St	Overall Depth candard 17" x 17" Overall Depth	I		430 mm	Product Weight	2 kg
U1717 Jay Union Cus Max User Weight U1718 Jay Union Cus Max User Weight U1719 Jay Union Cus	shion St ^{136 kg} shion St	andard 17" x 17" Overall Depth	I		430 mm	Product Weight	2 kg
Max User Weight	136 kg shion St	Overall Depth	430 mm				
U1718 Jay Union Cus Max User Weight U1719 Jay Union Cus	shion St		430 mm				
Max User Weight 1		tandard 17" x 18"		Overall Width	430 mm	Product Weight	2 kg
Max User Weight							
-		Overall Depth	455 mm	Overall Width	430 mm	Product Weight	2 kg
Max User Weight	shion St	tandard 17" x 19"				·	·
	136 kg	Overall Depth	480 mm	Overall Width	430 mm	Product Weight	2 kg
U1720 Jay Union Cu	ishion S	tandard 17" x 20"					
Max User Weight	136 kg	Overall Depth	505 mm	Overall Width	430 mm	Product Weight	2 kg
U1721 Jay Union Cu	shion St	tandard 17" x 21"					
Max User Weight	136 kg	Overall Depth	535 mm	Overall Width	430 mm	Product Weight	2 kg
U1814 Jay Union Cu	shion S	tandard 18" x 14"					
Max User Weight	136 kg	Overall Depth	355 mm	Overall Width	455 mm	Product Weight	2 kg
U1815 Jay Union Cu	shion St	tandard 18" x 15"					
Max User Weight	136 kg	Overall Depth	380 mm	Overall Width	455 mm	Product Weight	2 kg
U1414 Jay Union Cu	shion S	tandard 14" x 14"					
Max User Weight	136 kg	Overall Depth	355 mm	Overall Width	355 mm	Product Weight	2 kg
U1415 Jay Union Cu	shion St	tandard 14" x 15"					
Max User Weight	136 kg	Overall Depth	380 mm	Overall Width	355 mm	Product Weight	2 kg
U1416 Jay Union Cu	shion S ¹	tandard 14" x 16"					
Max User Weight	136 kg	Overall Depth	405 mm	Overall Width	355 mm	Product Weight	2 kg
U1417 Jay Union Cu	shion St	tandard 14" x 17"					
Max User Weight	136 kg	Overall Depth	430 mm	Overall Width	355 mm	Product Weight	2 kg
U1418 Jay Union Cu	shion S	tandard 14" x 18"					
Max User Weight	136 kg	Overall Depth	455 mm	Overall Width	355 mm	Product Weight	2 kg
111514 Jay Union Cu	shion St	tandard 15" x 14"					
Jay Union Cu		Overall Depth	355 mm	Overall Width			

JU1515 3	lay Union Ci	ushion S	tandard 15" x 15"		
Max User Wei	ght	136 kg	Overall Depth	380 mm	Over
JU1516	lay Union Ci	ushion S	tandard 15″ x 16″		
Max User Wei	ght	136 kg	Overall Depth	405 mm	Over
JU1517 J	lay Union Cu	ushion S	tandard 15" x 17"		
Max User Wei	ght	136 kg	Overall Depth	430 mm	Over
	on Cushid		 A versatile and comfortat Protects the skin and allo Composed of a dynamic f A moisture resistant inne The foam base is modera Features a top layer of so The tuitimate wheelchair s The fluid insert that lies b breakdown 	ws patient po iluid and foan r cover, and a tely contoure ft, visco men seating surfac	n laye an ant d to e nory fo ce tha
		ushion W	/ide 21" x 17"		
Max User Wei	ght	225 kg	Overall Depth	430 mm	Over
JU2118	lay Union Ci	ushion V	/ide 21" x 18"		

JU2118 Jay Union C	ushion V	Vide 21" x 18"		
Max User Weight	225 kg	Overall Depth	455 mm	(
JU2119 Jay Union C	ushion V	Vide 21" x 19"		
Max User Weight	225 kg	Overall Depth	480 mm	(
JU2120 Jay Union C	ushion \	Wide 21" x 20"		
Max User Weight	225 kg	Overall Depth	505 mm	(
JU2121 Jay Union C	ushion V	Vide 21" x 21"		
Max User Weight	225 kg	Overall Depth	535 mm	(
JU2122 Jay Union C	ushion \	Wide 21" x 22"		
Max User Weight	225 kg	Overall Depth	560 mm	
JU2123 Jay Union C	ushion \	Wide 21" x 23"		
Max User Weight	225 kg	Overall Depth	585 mm	(
JU2124 Jay Union C	ushion \	Vide 21" x 24"		
Max User Weight	225 kg	Overall Depth	610 mm	(
JU2217 Jay Union C	ushion V	Vide 22" x 17"		
Max User Weight	225 kg	Overall Depth	430 mm	
JU2218 Jay Union C	ushion \	Wide 22" x 18"		
Max User Weight	225 kg	Overall Depth	455 mm	(
JU2219 Jay Union C	ushion \	Vide 22" x 19"		
Max User Weight	225 ka	Overall Depth	480 mm	



Overall Width	380 mm	Product Weight	2 kg
Overall Width	380 mm	Product Weight	2 kg
Overall Width	380 mm	Product Weight	2 kg

ioning

ayering system anti-microbial outer cover

o encourage proper positioning of the pelvis and thighs

/ foam, with an extra layer covering the pelvic loading area

hat distributes weight for pressure reduction and maximum comfort

o foam conforms to each individual's unique contours for a reduced risk of skin

JU2117 Jay Union	Cushion \	Nide 21" x 17"					
Max User Weight	225 kg	Overall Depth	430 mm	Overall Width	535 mm	Product Weight	2 kg
JU2118 Jay Union	Cushion \	Wide 21" x 18"					
Max User Weight	225 kg	Overall Depth	455 mm	Overall Width	535 mm	Product Weight	2 kg
JU2119 Jay Union	Cushion \	Wide 21" x 19"					
Max User Weight	225 kg	Overall Depth	480 mm	Overall Width	535 mm	Product Weight	2 kg
JU2120 Jay Union	Cushion	Wide 21" x 20"					
Max User Weight	225 kg	Overall Depth	505 mm	Overall Width	535 mm	Product Weight	2 kg
JU2121 Jay Union	Cushion \	Wide 21" x 21"					
Max User Weight	225 kg	Overall Depth	535 mm	Overall Width	535 mm	Product Weight	2 kg
JU2122 Jay Union	Cushion	Wide 21" x 22"					
Max User Weight	225 kg	Overall Depth	560 mm	Overall Width	535 mm	Product Weight	2 kg
JU2123 Jay Union	Cushion	Wide 21" x 23"					
Max User Weight	225 kg	Overall Depth	585 mm	Overall Width	535 mm	Product Weight	2 kg
JU2124 Jay Union	Cushion	Wide 21" x 24"					
Max User Weight	225 kg	Overall Depth	610 mm	Overall Width	535 mm	Product Weight	2 kg
JU2217 Jay Union	Cushion	Wide 22" x 17"					
Max User Weight	225 kg	Overall Depth	430 mm	Overall Width	560 mm	Product Weight	2 kg
JU2218 Jay Union	Cushion	Wide 22" x 18"					
Max User Weight	225 kg	Overall Depth	455 mm	Overall Width	560 mm	Product Weight	2 kg
JU2219 Jay Union	Cushion	Wide 22" x 19"					
Max User Weight	225 kg	Overall Depth	480 mm	Overall Width	560 mm	Product Weight	2 kg

JU2220 Jay Union (Cushion	Wide 22" x 20"					
Max User Weight	225 kg	Overall Depth	505 mm	Overall Width	560 mm	Product Weight	2 kg
JU2221 Jay Union C	Cushion \	Wide 22" x 21"					
Max User Weight	225 kg	Overall Depth	535 mm	Overall Width	560 mm	Product Weight	2 kg
JU2222 Jay Union (Cushion	Wide 22" x 22"					
Max User Weight	225 kg	Overall Depth	560 mm	Overall Width	560 mm	Product Weight	2 kg
JU2223 Jay Union (Cushion	Wide 22" x 23"					
Max User Weight	225 kg	Overall Depth	585 mm	Overall Width	560 mm	Product Weight	2 kg
JU2224 Jay Union (Cushion	Wide 22" x 24"					
Max User Weight	225 kg	Overall Depth	610 mm	Overall Width	560 mm	Product Weight	2 kg
JU2319 Jay Union C	Cushion	Wide 23" x 19"					
Max User Weight	225 kg	Overall Depth	480 mm	Overall Width	585 mm	Product Weight	2 kg
JU2320 Jay Union (Cushion	Wide 23" x 20"					
Max User Weight	225 kg	Overall Depth	505 mm	Overall Width	585 mm	Product Weight	2 kg
JU2321 Jay Union C	ushion \	Wide 23" x 21"					
Max User Weight	225 kg	Overall Depth	535 mm	Overall Width	585 mm	Product Weight	2 kg
JU2322 Jay Union (Pushion						
JOLOLL JAY UNUIT	JUSIIIOII	Wide 23" x 22"					
Max User Weight	225 kg	VICE 23" x 22" Overall Depth	560 mm	Overall Width	585 mm	Product Weight	2 kg
-	225 kg	Overall Depth	560 mm	Overall Width	585 mm	Product Weight	2 kg
Max User Weight	225 kg	Overall Depth	560 mm 585 mm	Overall Width Overall Width	585 mm 585 mm	Product Weight Product Weight	2 kg 2 kg
Max User Weight	225 kg Cushion 225 kg	Overall Depth Wide 23" x 23" Overall Depth					
Max User Weight JU2323 Jay Union (Max User Weight	225 kg Cushion 225 kg	Overall Depth Wide 23" x 23" Overall Depth					
Max User Weight JU2323 Jay Union (Max User Weight JU2324 Jay Union (225 kg Cushion 225 kg Cushion 225 kg	Overall Depth Wide 23" x 23" Overall Depth Wide 23" x 24" Overall Depth	585 mm	Overall Width	585 mm	Product Weight	2 kg
Max User Weight JU2323 Jay Union (Max User Weight JU2324 Jay Union (Max User Weight JU2419 Jay Union (225 kg Cushion 225 kg Cushion 225 kg	Overall Depth Wide 23" x 23" Overall Depth Wide 23" x 24" Overall Depth	585 mm 610 mm	Overall Width	585 mm	Product Weight	2 kg
Max User Weight JU2323 Jay Union (Max User Weight JU2324 Jay Union (Max User Weight JU2419 Jay Union (225 kg Cushion 225 kg Cushion 225 kg Cushion N 225 kg	Overall Depth Wide 23" x 23" Overall Depth Wide 23" x 24" Overall Depth Wide 24" x 19" Overall Depth	585 mm 610 mm	Overall Width Overall Width	585 mm	Product Weight Product Weight	2 kg
Max User Weight JU2323 Jay Union (Max User Weight JU2324 Jay Union (Max User Weight JU2419 Jay Union (Max User Weight	225 kg Cushion 225 kg Cushion 225 kg Cushion N 225 kg	Overall Depth Wide 23" x 23" Overall Depth Wide 23" x 24" Overall Depth Wide 24" x 19" Overall Depth	585 mm 610 mm	Overall Width Overall Width	585 mm	Product Weight Product Weight	2 kg
Max User Weight JU2323 Jay Union (Max User Weight JU2324 Jay Union (Max User Weight JU2419 Jay Union (Max User Weight JU2420 Jay Union (225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg	Overall Depth Wide 23" x 23" Overall Depth Wide 23" x 24" Overall Depth Wide 24" x 19" Overall Depth Wide 24" x 20" Overall Depth	585 mm 610 mm 480 mm	Overall Width Overall Width Overall Width	585 mm 585 mm 610 mm	Product Weight Product Weight Product Weight	2 kg 2 kg 2 kg
Max User Weight JU2323 Jay Union (Max User Weight JU2324 Jay Union (Max User Weight JU2419 Jay Union (Max User Weight JU2420 Jay Union (Max User Weight	225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg	Overall Depth Wide 23" x 23" Overall Depth Wide 23" x 24" Overall Depth Wide 24" x 19" Overall Depth Wide 24" x 20" Overall Depth	585 mm 610 mm 480 mm	Overall Width Overall Width Overall Width	585 mm 585 mm 610 mm	Product Weight Product Weight Product Weight	2 kg 2 kg 2 kg
Max User Weight JU2323 Jay Union (Max User Weight JU2324 Jay Union (Max User Weight JU2419 Jay Union (Max User Weight JU2420 Jay Union (Max User Weight JU2421 Jay Union (225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg	Overall Depth Wide 23" x 23" Overall Depth Wide 23" x 24" Overall Depth Wide 24" x 19" Overall Depth Wide 24" x 20" Overall Depth Wide 24" x 21" Overall Depth	585 mm 610 mm 480 mm 505 mm	Overall Width Overall Width Overall Width Overall Width	585 mm 585 mm 610 mm 610 mm	Product Weight Product Weight Product Weight Product Weight	2 kg 2 kg 2 kg 2 kg
Max User Weight JU2323 Jay Union (Max User Weight JU2324 Jay Union (Max User Weight JU2419 Jay Union (Max User Weight JU2420 Jay Union (Max User Weight JU2421 Jay Union (Max User Weight	225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg	Overall Depth Wide 23" x 23" Overall Depth Wide 23" x 24" Overall Depth Wide 24" x 19" Overall Depth Wide 24" x 20" Overall Depth Wide 24" x 21" Overall Depth	585 mm 610 mm 480 mm 505 mm	Overall Width Overall Width Overall Width Overall Width	585 mm 585 mm 610 mm 610 mm	Product Weight Product Weight Product Weight Product Weight	2 kg 2 kg 2 kg 2 kg
Max User Weight JU2323 Jay Union (Max User Weight JU2324 Jay Union (Max User Weight JU2419 Jay Union (Max User Weight JU2420 Jay Union (Max User Weight JU2421 Jay Union (Max User Weight JU2422 Jay Union (Max User Weight	225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg Cushion	Overall Depth Wide 23" x 23" Overall Depth Wide 23" x 24" Overall Depth Wide 24" x 19" Overall Depth Wide 24" x 20" Overall Depth Wide 24" x 21" Overall Depth Wide 24" x 22"	585 mm 610 mm 480 mm 505 mm 535 mm	Overall Width Overall Width Overall Width Overall Width Overall Width	585 mm 585 mm 610 mm 610 mm	Product Weight Product Weight Product Weight Product Weight Product Weight	2 kg 2 kg 2 kg 2 kg
Max User Weight JU2323 Jay Union (Max User Weight JU2324 Jay Union (Max User Weight JU2419 Jay Union (Max User Weight JU2420 Jay Union (Max User Weight JU2421 Jay Union (Max User Weight JU2422 Jay Union (Max User Weight	225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg Cushion	Overall Depth Wide 23" x 23" Overall Depth Wide 23" x 24" Overall Depth Wide 24" x 19" Overall Depth Wide 24" x 20" Overall Depth Nide 24" x 21" Overall Depth Wide 24" x 22" Overall Depth	585 mm 610 mm 480 mm 505 mm 535 mm	Overall Width Overall Width Overall Width Overall Width Overall Width	585 mm 585 mm 610 mm 610 mm	Product Weight Product Weight Product Weight Product Weight Product Weight	2 kg 2 kg 2 kg 2 kg
Max User Weight JU2323 Jay Union (Max User Weight JU2324 Jay Union (Max User Weight JU2419 Jay Union (Max User Weight JU2420 Jay Union (Max User Weight JU2420 Jay Union (Max User Weight JU2421 Jay Union (Max User Weight JU2422 Jay Union (Max User Weight JU2423 Jay Union (Max User Weight	225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg Cushion 225 kg	Overall Depth Wide 23" x 23" Overall Depth Wide 23" x 24" Overall Depth Wide 24" x 19" Overall Depth Wide 24" x 20" Overall Depth Wide 24" x 21" Overall Depth Wide 24" x 22" Overall Depth Wide 24" x 22" Overall Depth	585 mm 610 mm 480 mm 505 mm 535 mm 5360 mm	Overall Width Overall Width Overall Width Overall Width Overall Width Overall Width	585 mm 585 mm 610 mm 610 mm 610 mm	Product Weight	2 kg 2 kg 2 kg 2 kg 2 kg 2 kg

Jay GS Cushion Designed for paediatric clients

Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	455 mm	Paediatric ?	у					
Product Weight	3 kg	Thickness	50 mm				İ					
11105.1N Jay GS Cushion 16" x 19"												
Max User Weight	113 kg	Overall Depth	480 mm	Overall Width	405 mm	Paediatric ?	у					
Product Weight	3 kg	Thickness	50 mm									
J1110.1N Jay GS Cushion 14" x 17"												
Max User Weight	113 kg	Overall Depth	430 mm	Overall Width	355 mm	Paediatric ?	у					
	nong											
Product Weight	3 kg	Thickness	50 mm									
	3 kg		50 mm									
Product Weight	3 kg		50 mm 330 mm	Overall Width	255 mm	Paediatric ?	у					
Product Weight	3 kg	x 13"		Overall Width	255 mm	Paediatric ?	У					
Product Weight J11111.1N Jay GS C Max User Weight	3 kg Cushion 10" 113 kg 3 kg	X 13" Overall Depth Thickness	330 mm	Overall Width	255 mm	Paediatric ?	у					
Product Weight J11111.1N Jay GS C Max User Weight Product Weight	3 kg Cushion 10" 113 kg 3 kg	X 13" Overall Depth Thickness	330 mm	Overall Width Overall Width	255 mm 305 mm	Paediatric ? Paediatric ?	у					
Product Weight J1111.1N Jay GS C Max User Weight Product Weight J1115.1N Jay GS C	3 kg Cushion 10" 113 kg 3 kg Cushion 12"	x 13" Overall Depth Thickness x 15"	330 mm 50 mm									
Product Weight J11111.1N Jay GS C Max User Weight Product Weight J1115.1N Jay GS C Max User Weight	3 kg 3 kg 113 kg 3 kg	x 13" Overall Depth Thickness x 15" Overall Depth	330 mm 50 mm 305 mm									

Varilite Meridian Bariatric Wheelchair Cushion



Powered

- Creates an integral pre-ischial ridge and wedge effect Enhanced comfort, support and pressure redistribution Also provides skin protection and postural support • The reticulated foam inside the cover promotes air exchange • Lightweight, comfortable and funcitonal cushion

 - No user Weight limit

7203B Varilite Meridian Bariatric Cushion Mesh Cover 22" x 20"										
Overall Depth	505 mm	Overall Width	560 mm	Product Weight	0.5 kg	Thickness	100 mm			
7204B Varilite Meri	dian Baı	iatric Cushion Incon	itinence	Cover 22" x 20"						
Overall Depth	505 mm	Overall Width	560 mm	Product Weight	0.5 ka	Thickness	100 mm			



• With moderate to aggressive asymmetric posture and high postural support requirements • Lightweight foam base with growth capability via well inserts and a wide range of postural support • For added skin protection there is a pressure relief pad with more fluid

Suitable for clients unable to weight shift, limited postural stability and unable to reposition

7.04 Pressure Cushions Combination

• Features dual chambers , air -foam flotation valves for adjusting air volume

• The cushion is available in a mesh or incontinence cover with breathable two-way stretch material

Overall Depth	455 mm	Overall Width	560 mm	Product Weight	0.5 kg	Thickness	100 mm				
7284B Varilite Meridian Bariatric Cushion Incontinence Cover 22" x 18"											
Overall Depth	455 mm	Overall Width	560 mm	Product Weight	0.5 kg	Thickness	100 mm				
	Invidion Dov	detrie Cushien Meek	Cover	4" 00"							
7403B Varilite Meridian Bariatric Cushion Mesh Cover 24" x 20"											
AUSD variate iv	lendian Bar	natric cushion mesh	i cover 2	4" X 20"							
Overall Depth	505 mm	Overall Width	610 mm	Product Weight	0.5 kg	Thickness	100 mm				
					0.5 kg	Thickness	100 mm				
Overall Depth	505 mm	Overall Width	610 mm	Product Weight	0.5 kg	Thickness	100 mm				
Overall Depth	505 mm		610 mm	Product Weight	0.5 kg	Thickness	100 mm				
Overall Depth	505 mm	Overall Width	610 mm	Product Weight	0.5 kg 0.5 kg	Thickness	100 mm				
Overall Depth 7404B Varilite N	505 mm	Overall Width	610 mm	Product Weight Cover 24" x 20"							
Overall Depth 7404B Varilite M Overall Depth	505 mm Aeridian Bar	Overall Width riatric Cushion Incor Overall Width	610 mm tinence 610 mm	Product Weight Cover 24" x 20" Product Weight							
Overall Depth 7404B Varilite M Overall Depth	505 mm Aeridian Bar	Overall Width	610 mm tinence 610 mm	Product Weight Cover 24" x 20" Product Weight							

Varilite Evolution Bariatric Wheelchair Cushion



- Cushion cover to replace original cover Suitable for use with Evolution PSV Bariatric Wheelchair Cushion
- The Evolution Cushion is for users with a high risk of tissue breakdown
- The Evolution Cushion provides outstanding pressure distribution and postural support
- Available in a four-way-stretch mesh or incontinence cover
- Available in a range of sizes 22"-24" x 18"-20"

Overall Depth	505 mm	Overall Width	560 mm
202B Varilite Evolution B	Bariatric Cushion Inco	ntinence Cover 22" x 20"	
Overall Depth	505 mm	Overall Width	560 mm
281B Varilite Evolution Ba	ariatric Cushion Mesł	n Cover 22" x 18"	
Overall Depth	455 mm	Overall Width	560 mm
282B Varilite Evolution B	ariatric Cushion Inco	ntinence Cover 22" x 18"	
Overall Depth	455 mm	Overall Width	560 mm
			560 mm
401B Varilite Evolution Ba			560 mm 610 mm
401B Varilite Evolution Ba	ariatric Cushion Mesh	D Cover 24" x 20" Overall Width	
Overall Depth 401B Varilite Evolution Ba Overall Depth 402B Varilite Evolution B Overall Depth	ariatric Cushion Mesh	D Cover 24" x 20" Overall Width	
401B Varilite Evolution Ba Overall Depth 402B Varilite Evolution B	ariatric Cushion Mesh 505 mm Cariatric Cushion Inco 505 mm	n Cover 24" x 20" Overall Width ntinence Cover 24" x 20" Overall Width	610 mm
401B Varilite Evolution Ba Overall Depth 402B Varilite Evolution B Overall Depth	ariatric Cushion Mesh 505 mm Cariatric Cushion Inco 505 mm	n Cover 24" x 20" Overall Width ntinence Cover 24" x 20" Overall Width	610 mm
401B Varilite Evolution Ba Overall Depth 402B Varilite Evolution B Overall Depth 481B Varilite Evolution Ba	ariatric Cushion Mesh 505 mm cariatric Cushion Inco 505 mm ariatric Cushion Mesh 455 mm	n Cover 24" x 20" Overall Width ntinence Cover 24" x 20" Overall Width n Cover 24" x 18" Overall Width	610 mm 610 mm

Varilite Meridian Wheelchair Cushion • Features dual chambers , air -foam floatation valves for adjusting air volume Creates an integral pre-ischial ridge and wedge effect • Enhanced comfort, support and pressure redistribution Also provides skin protection and postural support • The reticulated foam inside the cover promotes air exchange Lightweight, comfortable and funcitonal cushion • No user Weight limit

78230 Varilite Meridian Cushion Mesh Cover 18" x 22"

78240 Varilite Meridian Cushion Incontinence Cover 18" x 22"

Varilite Zoid PSV Wheelchair Cushion



- Lightweight and comfortable • Low profile without a medial thigh separator
- Enables day to day activities
- Vibration dampening
 - Its tapered shape makes it ideal for rigid chairs

31110 Varilite Zoid PSV Cushion Mesh Cover 10" x 10"

31111 Varilite Zoid PSV Cushion Incontinence Cover 10" x 10"

38310 Varilite Zoid PSV Cushion Incontinence Cover 18" x 16"

Max User Weight	260 kg	Overall Depth	405 mm	Overall Length	457 mm	Overall Width	455 mm
Product Weight	0.65 kg	Thickness	74 mm				

Varilite ProForm NX Dual Chamber Cushion



- Cost effective, durable, adjustable and easy-to-customize product
- A modifiable contoured base, superior pressure distribution, asymmetric positioning options, comfort and simplicity • An ideal off-the shelf solution when a customizable wheelchair cushion is needed
- Closed-cell foam base will not absorb water and is resistant to bacteria
- Contains two self-inflating, independently-adjustable chambers
- The foam conforms to the user, and the user's weight is distributed over the surface of the cushion · Available with a mesh or incontinence cover with breathable four-way-stretch material
- 78712 Varilite ProForm NX Dual Chamber Cushion Mesh (Overall Depth 431 mm Overall Width

Pressure Cushions



- The cushion is available in a mesh or incontinence cover with breathable two-way stretch material

Performance cushion of choice for active wheelchair users at moderate risk of tissue breakdown

• Its cover is designed to protect the thighs of users with tight hamstrings or who tuck their legs under their chair

• Recommended for individuals with leg length discrepancies, amputations, and other pelvic obliquities

Cover 15" x 17"					
381 mm Thickness 12	20 mm				



Roho AirLite Cushion



- A contoured, anatomical-shaped foam with a pre-set, sealed Air Flotation component • Enables users to sit supported and comfortably for an extended time
- Long lasting, lightweight and durable

posture support and comfort

- . Leg (femur) troughs and an ischial well provide support for lower extremity control, stability and comfort
- Fluid resistant polyurethane film laminated to 2 way stretch, polyester knit, and non skid bottom cover

• Appropriate for clients who have partial to full sensation, are able to perform independent weight shifts but require basic

R-A/AL1917 Roho AirLite Cushion 19" x 17"

Max User Weight	:	113 kg	Overall Depth	430 mm	Overall Width	480 mm	Thickness Back Edge	83 mm
Thickness Front I	Edge	50 mm						

• Aids basic positioning needs to assist in the prevention of skin breakdown

R-A/AL1920 Roho AirLite Cushion 19" x 20"

Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	480 mm	Thickness Back Edge	83 mm
Thickness Front Edge	50 mm						

R-A/ALCUSTOM Roho AirLite Cushion Custom

Max User Weight	113 kg	Thickness Back Edge	83 mm	Thickness Front Edge	50 mm
-----------------	--------	---------------------	-------	----------------------	-------

R-A/AL1414 Roho AirLite Cushion 14" x 14"

Max User Weight	113 kg	Overall Depth	355 mm	Overall Width	355 mm	Thickness Back Edge	83 mm
Thickness Front Edge	50 mm						

R-A/AL1416 Roho AirLite Cushion 14" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	355 mm	Thickness Back Edge	83 mm
Thickness Front Edge	50 mm						

R-A/AL1515 Roho AirLite Cushion 15" x 15"

[Max User Weight	113 kg	Overall Depth	380 mm	Overall Width	380 mm	Thickness Back Edge	83 mm
[Thickness Front Edge	50 mm						

R-A/AL1616 Roho AirLite Cushion 16" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	405 mm	Thickness Back Edge	83 mm
Thickness Front Edge	50 mm						

R-A/AL1618 Roho AirLite Cushion 16" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	405 mm	Thickness Back Edge	83 mm
Thickness Front Edge	50 mm						

R-A/AL1716 Roho AirLite Cushion 17" x 16"

Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	430 mm	Thickness Back Edge	83 mm
Thickness Front Edge	50 mm						

R-A/AL1717 Roho AirLite Cushion 17" x 17"

Max User Weight	113 kg	Overall Depth	430 mm	Overall Width	430 mm	Thickness Back Edge	83 mm
Thickness Front Edge	50 mm						

R-A/AL1718 Roho AirLite Cushion 17" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	430 mm	Thickness Back Edge	83 mm
Thickness Front Edge	50 mm						

R-A/AL1816 Roho AirLite Cushion 18" x 16"

· · · · · · · · · · · · · · · · · · ·								
Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	455 mm	Thickness Back Edge	83 mm	
Thickness Front Edge	50 mm							
R-A/AL1818 Roho AirLite Cushion 18" x 18"								

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	455 mm	Thickness Back Edge	83 mm
Thickness Front Edge	50 mm						

R-A/AL1820 Roho AirLite Cushion 18" x 20"

	Max User Weight	113 kg	Overall Depth	505 mm	Overall Width	455 mm	Thickness Back Edge	83 mm
	Thickness Front Edge	50 mm						

R-A/AL2018 Roho AirLite Cushion 20" x 18"

Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	505 mm	Thickness Back Edge	83 mm
Thickness Front Edge	50 mm						

Roho Harmony Cushion



- The Harmony $\ensuremath{^{\rm T}}$ cushion is a lightweight and stable cushioning solution resistant cover
- maximise their functional potential
- This enhances blood flow and provides longer sitting times • Available in a range of sizes 14"-20" x16"-18"

Max User Weight	110 km										
	113 kg	Overall Depth	405 mm	Overall Width	355 mm	Product Weight	0.7 kg				
R-HARMONY-1616C R	Roho Har	mony Cushion 16" x	16″								
Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	405 mm	Product Weight	0.7 kg				
R-HARMONY-1618C R	Roho Har	Roho Harmony Cushion 16" x 18"									
Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	405 mm	Product Weight	0.7 kg				
R-HARMONY-1716C R	Roho Har	mony Cushion 17" x	16″								
Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	430 mm	Product Weight	0.7 kg				
R-HARMONY-1718C R	Roho Har	mony Cushion 17" x	18″								
Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	430 mm	Product Weight	0.7 kg				
R-HARMONY-1816C R	Roho Har	mony Cushion 18" x	16″								
Max User Weight	113 kg	Overall Depth	405 mm	Overall Width	455 mm	Product Weight	0.7 kg				
R-HARMONY-1818C R	Roho Har	mony Cushion 18" x	18″								
Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	455 mm	Product Weight	0.7 kg				
R-HARMONY-2018C	Roho Ha	rmony Cushion 20"	x 18″								
Max User Weight	113 kg	Overall Depth	455 mm	Overall Width	505 mm	Product Weight	0.7 kg				

199

• A lightweight and comfortable Roho air floatation cushion with a pre-contoured foam base, encased within a moisture

 Can maximize abilities, independence and freedom without adding unnecessary weight to the mobility system • The cushion is built on ananatomically shaped, foam platform that enhances the users seated stability and helps them

Country Care GROUF



Roho Hybrid Elite Cushion Single Valve



- Cushion for wheelchairs
- Combines the superb skin protection and shape matching ability of dry floatation with the stability of jay base technology Allows for improved pelvic alignment and stability Improved position and function
- While maintaining optimal immersion and envelopment of both the ischial tuberosities and greater trochanters
 Provides critical pressure redistribution around the pelvis area and a reduced friction and shear environment
- Cushion is adjustable to fit clients shape

R-HE-2RHYHD1414C	Roho Hybrid B	Elite Cushion Heavy Du	ity Cover 14" x 14"	
Overall Depth		355 mm	Overall Width	355 mm
R-HE-2RHYHD1415C	Roho Hybrid E	Elite Cushion Heavy Du	ity Cover 14" x 15"	
Overall Depth		380 mm	Overall Width	355 mm
R-HE-2RHYHD1416C	Roho Hybrid B	Elite Cushion Heavy Du	ity Cover 14" x 16"	
Overall Depth		405 mm	Overall Width	355 mm
R-HE-2RHYHD1515C	Roho Hybrid E	Elite Cushion Heavy Du	ity Cover 15" x 15"	
Overall Depth		380 mm	Overall Width	380 mm
R-HE-2RHYHD1516C	Roho Hybrid E	Elite Cushion Heavy Du	ty Cover 15" x 16"	
Overall Depth		405 mm	Overall Width	380 mm
R-HE-2RHYHD1517C	Roho Hybrid E	lite Cushion Heavy Du	ty Cover 15" x 17"	
Overall Depth		430 mm	Overall Width	380 mm
R-HE-2RHYHD1615C	Roho Hybrid E	lite Cushion Heavy Du	ty Cover 16" x 15"	
Overall Depth		380 mm	Overall Width	405 mm
R-HE-2RHYHD1616C	Roho Hybrid E	lite Cushion Heavy Du	ty Cover 16" x 16"	
Overall Depth		405 mm	Overall Width	405 mm
R-HE-2RHYHD1617C	Roho Hybrid E	lite Cushion Heavy Du	ty Cover 16" x 17"	
Overall Depth		430 mm	Overall Width	405 mm
R-HE-2RHYHD1618C	Roho Hybrid E	lite Cushion Heavy Du	ty Cover 16" x 18"	
Overall Depth		455 mm	Overall Width	405 mm
R-HE-2RHYHD1620C	Roho Hybrid	Elite Cushion Heavy Du	uty Cover 16" x 20"	
Overall Depth		505 mm	Overall Width	405 mm
R-HE-2RHYHD1716C	Roho Hybrid E	lite Cushion Heavy Du	ty Cover 17" x 16"	
Overall Depth		405 mm	Overall Width	430 mm
R-HE-2RHYHD1717C	Roho Hybrid E	lite Cushion Heavy Du	ty Cover 17" x 17"	
Overall Depth		430 mm	Overall Width	430 mm
R-HE-2RHYHD1718C	Roho Hybrid E	lite Cushion Heavy Du	ty Cover 17" x 18"	
Overall Depth		455 mm	Overall Width	430 mm

R-HE-2RHYHD1816C	Dobo Hybrid E	lite Cushion Heavy Du	+v C
Overall Depth		405 mm	Over
overall beptil		403 mm	UVEI
R-HE-2RHYHD1817C	Roho Hybrid E	lite Cushion Heavy Du	ty C
Overall Depth		430 mm	Over
R-HE-2RHYHD1818C	Roho Hybrid E	lite Cushion Heavy Du	ty C
Overall Depth		455 mm	Overa
R-HE-2RHYHD1820C	Roho Hybrid I	Elite Cushion Heavy Du	ity (
Overall Depth		505 mm	Over
R-HE-2RHYHD2018C	Roho Hybrid I	Elite Cushion Heavy Du	ity (
Overall Depth		455 mm	Over
R-HE-2RHYHD2020C	Roho Hybrid	Elite Cushion Heavy D	uty
Overall Depth		505 mm	Over
R-HE-2RHYHD2218C	Roho Hybrid I	Elite Cushion Heavy Du	ity C
Overall Depth		455 mm	Over
R-HE-2RHYHD2220C	Roho Hybrid	Elite Cushion Heavy D	uty (
Overall Depth		505 mm	Over
R-HE-2RHYHD2418C	Roho Hybrid I	Elite Cushion Heavy Du	ity C
Overall Depth		455 mm	Over
R-HE-2RHYHD2420C	Roho Hybrid	Elite Cushion Heavy D	uty
Overall Depth		505 mm	Over
Roho Hybrid Elit	e Cushion S	ingle Valve	
Se de la compañía de	Co Allu Im Wh Pro	shion for wheelchairs mbines the superb skin protection a ows for improved pelvic alignment a proved position and function ille maintaining optimal immersion ovides critical pressure redistribution shion is adjustable to fit clients sha	and sta and er in arou
R-HE-1RHYHD1414C	Roho Hybrid E	lite Cushion Heavy Du	ty C
Overall Depth		355 mm	Over
R-HE-1RHYHD1415C	Roho Hybrid E	lite Cushion Heavy Du	ty C

R-HE-1RHYHD1414C Roho Hybrid Elite Cushion Heavy Duty Cover 14" x 14"				
Overall Depth	355 mm	Overall Width	355 mm	
R-HE-1RHYHD1415C Roho Hybrid Elite Cushion Heavy Duty Cover 14" x 15"				
Overall Depth	380 mm	Overall Width	355 mm	
R-HE-1RHYHD1416C Roho Hybrid Elite Cushion Heavy Duty Cover 14" x 16"				
Overall Depth	405 mm	Overall Width	355 mm	
R-HE-1RHYHD1515C Roho Hybrid Elite Cushion Heavy Duty Cover 15" x 15"				
Overall Depth	380 mm	Overall Width	380 mm	

Pressure Cushions



y Cover 18" x 16"	
Overall Width	455 mm
y Cover 18" x 17"	
Overall Width	455 mm
y Cover 18" x 18"	
Overall Width	455 mm
ty Cover 18" x 20"	
Overall Width	455 mm
ty Cover 20" x 18"	
Overall Width	505 mm
ty Cover 20" x 20"	
Overall Width	505 mm
ty Cover 22" x 18"	
Overall Width	560 mm
ty Cover 22" x 20"	
Overall Width	560 mm
y Cover 24" x 18"	
Overall Width	610 mm
ty Cover 24" x 20"	
Overall Width	610 mm

shape matching ability of dry floatation with the stability of jay base technology tability

envelopment of both the ischial tuberosities and greater trochanters round the pelvis area and a reduced friction and shear environment

B LIE ABOUT STREET				
R-HE-IRHYHD1516C	Roho Hybrid E	Elite Cushion Heavy Du	ty Cover 15" x 16"	
Overall Depth		405 mm	Overall Width	380 mm
R-HE-IRHYHD1517C	Roho Hybrid E	lite Cushion Heavy Du	ty Cover 15" x 17"	
Overall Depth		430 mm	Overall Width	380 mm
		lite Cuchien Lleew Du	ty Cover 16" x 15"	
		Elite Cushion Heavy Du	-	405
Overall Depth		380 mm	Overall Width	405 mm
R-HE-1RHYHD1616C	Roho Hybrid F	Elite Cushion Heavy Dut	tv Cover 16" x 16"	
Overall Depth		405 mm	Overall Width	405 mm
R-HE-1RHYHD1617C	Roho Hybrid E	Elite Cushion Heavy Dut	ty Cover 16" x 17"	
Overall Depth		430 mm	Overall Width	405 mm
L				,
R-HE-1RHYHD1618C	Roho Hybrid E	Elite Cushion Heavy Du	ty Cover 16" x 18"	
Overall Depth		455 mm	Overall Width	405 mm
K-HE-IRHYHD1620C	Roho Hybrid	Elite Cushion Heavy Du	ty Cover 16" x 20"	
Overall Depth		505 mm	Overall Width	405 mm
		lite Cushion Heavy Du	ty Covor 17" y 16"	
Overall Depth		405 mm	Overall Width	430 mm
R-HE-1RHYHD1717C	Roho Hybrid F	lite Cushion Heavy Dut	ty Cover 17" x 17"	
Overall Depth	Kono nyona E	430 mm	Overall Width	430 mm
overan bepar				
R-HE-1RHYHD1718C	Roho Hybrid E	Elite Cushion Heavy Dut	ty Cover 17" x 18"	
Overall Depth	-	455 mm	Overall Width	430 mm
R-HE-1RHYHD1816C	Roho Hybrid E	Elite Cushion Heavy Du	ty Cover 18" x 16"	
Overall Depth		405 mm	Overall Width	455 mm
R-HE-1RHYHD1817C	Daha Uuhuld F			
Overall Depth		lite Cushion Heavy Du	ty Cover 18" x 17"	
oreitan beptil		lite Cushion Heavy Dut	ty Cover 18" x 17" Overall Width	455 mm
	-	430 mm	Overall Width	455 mm
R-HE-1RHYHD1818C	-	430 mm	Overall Width ty Cover 18" x 18"	
	-	430 mm	Overall Width	455 mm 455 mm
R-HE-1RHYHD1818C Overall Depth	Roho Hybrid E	430 mm Elite Cushion Heavy Dut	Overall Width ty Cover 18" x 18" Overall Width	
R-HE-1RHYHD1818C Overall Depth R-HE-1RHYHD1820C	Roho Hybrid E	430 mm Elite Cushion Heavy Dur 455 mm Elite Cushion Heavy Du	Overall Width ty Cover 18" x 18" Overall Width ity Cover 18" x 20"	455 mm
R-HE-1RHYHD1818C Overall Depth	Roho Hybrid E	430 mm Elite Cushion Heavy Dut	Overall Width ty Cover 18" x 18" Overall Width	
R-HE-1RHYHD1818C Overall Depth R-HE-1RHYHD182OC Overall Depth	Roho Hybrid E Roho Hybrid I	430 mm Elite Cushion Heavy Dur 455 mm Elite Cushion Heavy Du	Overall Width ty Cover 18" x 18" Overall Width ty Cover 18" x 20" Overall Width	455 mm
R-HE-1RHYHD1818C Overall Depth R-HE-1RHYHD182OC Overall Depth	Roho Hybrid E Roho Hybrid I	430 mm Elite Cushion Heavy Dur 455 mm Elite Cushion Heavy Dur 505 mm	Overall Width ty Cover 18" x 18" Overall Width ty Cover 18" x 20" Overall Width	455 mm
R-HE-1RHYHD1818C Overall Depth R-HE-1RHYHD182OC Overall Depth R-HE-1RHYHD2018C	Roho Hybrid E Roho Hybrid I	430 mm Elite Cushion Heavy Dut 455 mm Elite Cushion Heavy Du 505 mm Elite Cushion Heavy Du	Overall Width ty Cover 18" x 18" Overall Width ty Cover 18" x 20" Overall Width ty Cover 20" x 18"	455 mm
R-HE-1RHYHD1818C Overall Depth R-HE-1RHYHD182OC Overall Depth R-HE-1RHYHD2018C	Roho Hybrid E Roho Hybrid I Roho Hybrid I	430 mm Elite Cushion Heavy Dut 455 mm Elite Cushion Heavy Du 505 mm Elite Cushion Heavy Du	Overall Width ty Cover 18" x 18" Overall Width ty Cover 18" x 20" Overall Width ty Cover 20" x 18" Overall Width Overall Width	455 mm
R-HE-1RHYHD1818C Overall Depth R-HE-1RHYHD182OC Overall Depth R-HE-1RHYHD2018C Overall Depth	Roho Hybrid E Roho Hybrid I Roho Hybrid I	430 mm Elite Cushion Heavy Dut 455 mm Elite Cushion Heavy Du 505 mm Elite Cushion Heavy Du 455 mm	Overall Width ty Cover 18" x 18" Overall Width ty Cover 18" x 20" Overall Width ty Cover 20" x 18" Overall Width Overall Width	455 mm
R-HE-1RHYHD1818C Overall Depth R-HE-1RHYHD182OC Overall Depth R-HE-1RHYHD2018C Overall Depth R-HE-1RHYHD2020C Overall Depth	Roho Hybrid E Roho Hybrid I Roho Hybrid I Roho Hybrid I	430 mm Elite Cushion Heavy Dut 455 mm Elite Cushion Heavy Du 505 mm Elite Cushion Heavy Du 455 mm Elite Cushion Heavy Du 505 mm	Overall Width ty Cover 18" x 18" Overall Width ty Cover 18" x 20" Overall Width ty Cover 20" x 18" Overall Width uty Cover 20" x 20" Overall Width	455 mm 455 mm 505 mm
R-HE-1RHYHD1818C Overall Depth R-HE-1RHYHD182OC Overall Depth R-HE-1RHYHD2018C Overall Depth R-HE-1RHYHD2020C Overall Depth	Roho Hybrid E Roho Hybrid I Roho Hybrid I Roho Hybrid I	430 mm Elite Cushion Heavy Dut 455 mm Elite Cushion Heavy Du 505 mm Elite Cushion Heavy Du 455 mm Elite Cushion Heavy Du	Overall Width ty Cover 18" x 18" Overall Width ty Cover 18" x 20" Overall Width ty Cover 20" x 18" Overall Width uty Cover 20" x 20" Overall Width	455 mm 455 mm 505 mm
R-HE-1RHYHD1818C Overall Depth R-HE-1RHYHD182OC Overall Depth R-HE-1RHYHD2018C Overall Depth R-HE-1RHYHD2020C Overall Depth	Roho Hybrid E Roho Hybrid I Roho Hybrid I Roho Hybrid I	430 mm Elite Cushion Heavy Dut 455 mm Elite Cushion Heavy Du 505 mm Elite Cushion Heavy Du 455 mm Elite Cushion Heavy Du 505 mm	Overall Width ty Cover 18" x 18" Overall Width ty Cover 18" x 20" Overall Width ty Cover 20" x 18" Overall Width uty Cover 20" x 20" Overall Width	455 mm 455 mm 505 mm
R-HE-1RHYHD1818C Overall Depth R-HE-1RHYHD182OC Overall Depth R-HE-1RHYHD2018C Overall Depth R-HE-1RHYHD202OC Overall Depth R-HE-1RHYHD2218C Overall Depth	Roho Hybrid E Roho Hybrid I Roho Hybrid I Roho Hybrid I Roho Hybrid I	430 mm Elite Cushion Heavy Dut 455 mm Elite Cushion Heavy Du 505 mm Elite Cushion Heavy Du 455 mm Elite Cushion Heavy Du 505 mm Elite Cushion Heavy Du 455 mm	Overall Width ty Cover 18" x 18" Overall Width ty Cover 18" x 20" Overall Width ty Cover 20" x 18" Overall Width Ity Cover 20" x 20" Overall Width ty Cover 22" x 18" Overall Width ty Cover 22" x 18"	455 mm 455 mm 505 mm
R-HE-1RHYHD1818C Overall Depth R-HE-1RHYHD1820C Overall Depth R-HE-1RHYHD2018C Overall Depth R-HE-1RHYHD2020C Overall Depth R-HE-1RHYHD2218C Overall Depth R-HE-1RHYHD2220C	Roho Hybrid E Roho Hybrid I Roho Hybrid I Roho Hybrid I Roho Hybrid I	430 mm Elite Cushion Heavy Dut 455 mm Elite Cushion Heavy Du 505 mm Elite Cushion Heavy Du 455 mm Elite Cushion Heavy Du 505 mm Elite Cushion Heavy Du 455 mm	Overall Width ty Cover 18" x 18" Overall Width ty Cover 18" x 20" Overall Width ty Cover 18" x 20" Overall Width Uty Cover 20" x 18" Overall Width ty Cover 20" x 20" Overall Width ty Cover 22" x 18" Overall Width	455 mm 455 mm 505 mm 505 mm 560 mm
R-HE-1RHYHD1818C Overall Depth R-HE-1RHYHD182OC Overall Depth R-HE-1RHYHD2018C Overall Depth R-HE-1RHYHD202OC Overall Depth R-HE-1RHYHD2218C Overall Depth	Roho Hybrid E Roho Hybrid I Roho Hybrid I Roho Hybrid I Roho Hybrid I	430 mm Elite Cushion Heavy Dut 455 mm Elite Cushion Heavy Du 505 mm Elite Cushion Heavy Du 455 mm Elite Cushion Heavy Du 505 mm Elite Cushion Heavy Du 455 mm	Overall Width ty Cover 18" x 18" Overall Width ty Cover 18" x 20" Overall Width ty Cover 20" x 18" Overall Width Ity Cover 20" x 20" Overall Width ty Cover 22" x 18" Overall Width ty Cover 22" x 18"	455 mm 455 mm 505 mm

overan bepan			400 11111	
R-HE-1RHYHD242	OC Roho Hy	ybrid I	Elite Cushion Heavy	y Duty
Overall Depth			505 mm	0\
Roho Hybrid I	Elite Cush	ion S	ingle Valve	
Contraction of the second seco		 Co All Im Wf Pro 	shion for wheelchairs mbines the superb skin protect ows for improved pelvic alignn proved position and function ille maintaining optimal immer ovides critical pressure redistr shion is adjustable to fit client	nent and rsion and ibution a
R-HE-1RHY1414C	Roho Hybri	d Elite	e Cushion Standard	l Cove
Overall Depth	-		355 mm	01
R-HE-1RHY1415C	Roho Hybri	d Elite	e Cushion Standard	l Cove
Overall Depth			380 mm	0\
R-HE-1RHY1416C	Roho Hybri	d Elite	e Cushion Standard	l Cove
Overall Depth			405 mm	0
R-HE-1RHY1515C	Roho Hybri	d Elite	e Cushion Standard	l Cove
Overall Depth			380 mm	01
R-HE-1RHY1516C	Roho Hybri	d Elite	Cushion Standard	
Overall Depth			405 mm	0\
R-HE-1RHY1517C	Roho Hybri	d Elite	Cushion Standard	l Cove
Overall Depth			430 mm	01
R-HE-1RHY1615C	Roho Hybri	d Elite	Cushion Standard	l Cove
Overall Depth			380 mm	0\
R-HE-1RHY1616C	Roho Hybri	d Elite	e Cushion Standard	l Cove
Overall Depth			405 mm	0\
R-HE-1RHY1617C	Roho Hybri	d Elite	Cushion Standard	
Overall Depth			430 mm	0\
R-HE-1RHY1618C	Roho Hybri	d Elite	e Cushion Standard	l Cove
Overall Depth			455 mm	01
R-HE-1RHY1620C	Roho Hybr	id Elit	e Cushion Standard	d Cove
Overall Depth			505 mm	01
R-HE-1RHY1716C	Roho Hybri	d Elite	Cushion Standard	
Overall Depth			405 mm	0\

455 mm

Overall Depth



R-HE-1RHYHD2418C Roho Hybrid Elite Cushion Heavy Duty Cover 24" x 18"

nd shape matching ability of dry floatation with the stability of jay base technology Id stability

nd envelopment of both the ischial tuberosities and greater trochanters a around the pelvis area and a reduced friction and shear environment be

ver 14" x 14"	
Overall Width	355 mm
ver 14" x 15"	
Overall Width	355 mm
ver 14" x 16"	
Overall Width	355 mm
ver 15" x 15"	
Overall Width	380 mm
ver 15" x 16"	
Overall Width	380 mm
ver 15" x 17"	
Overall Width	380 mm
ver 16" x 15"	
Overall Width	405 mm
ver 16" x 16"	
Overall Width	405 mm
ver 16" x 17"	
Overall Width	405 mm
ver 16" x 18"	
Overall Width	405 mm
ver 16" x 20"	
Overall Width	405 mm
ver 17" x 16"	
Overall Width	430 mm

R-HE-1RHY1717C Roho Hybrid Eli	te Cushion Standard Co	over 17" x 17"	
Overall Depth	430 mm	Overall Width	430 mm
		17/ 10/	
R-HE-1RHY1718C Roho Hybrid Eli			
Overall Depth	455 mm	Overall Width	430 mm
R-HE-1RHY1816C Roho Hybrid Eli	ite Cushion Standard Co	over 18" x 16"	
Overall Depth	405 mm	Overall Width	455 mm
R-HE-1RHY1817C Roho Hybrid Eli	te Cushion Standard Co	over 18" x 17"	
Overall Depth	430 mm	Overall Width	455 mm
R-HE-1RHY1818C Roho Hybrid El	ite Cushion Standard Co	over 18" x 18"	
Overall Depth	455 mm	Overall Width	455 mm
R-HE-1RHY1820C Roho Hybrid El	ite Cushion Standard C	over 18" x 20"	
Overall Depth	505 mm	Overall Width	455 mm
R-HE-1RHY2018C Roho Hybrid El	ite Cushion Standard C	over 20" x 18"	
Overall Depth	455 mm	Overall Width	505 mm
R-HE-1RHY2020C Roho Hybrid E	lite Cushion Standard C	Cover 20" x 20"	
Overall Depth	505 mm	Overall Width	505 mm
R-HE-1RHY2218C Roho Hybrid El	ite Cushion Standard C	over 22" x 18"	
Overall Depth	455 mm	Overall Width	560 mm
R-HE-1RHY2220C Roho Hybrid E	lite Cushion Standard C	Cover 22" x 20"	
Overall Depth	505 mm	Overall Width	560 mm
R-HE-1RHY2418C Roho Hybrid Elite Cushion Standard Cover 24" x 18"			
Overall Depth	455 mm	Overall Width	610 mm
R-HE-1RHY2420C Roho Hybrid E	lite Cushion Standard C	cover 24" x 20"	
Overall Depth	505 mm	Overall Width	610 mm

Roho Hybrid Elite Cushion Dual Valve



 Cushion for wheelchairs Combines the superb skin protection and shape matching ability of dry floatation with the stability of jay base technology 	
Allows for improved pelvic alignment and stability	
Improved position and function	
While maintaining optimal immersion and envelopment of both the ischial tuberosities and greater trochanters	
Provides critical pressure redistribution around the pelvis area and a reduced friction and shear environment	
 Cushion is adjustable to fit clients shape 	

R-HE-2RHY1414C Roho Hybrid Elite Cushion Standard Cover 14" x 14"				
Overall Depth 355 mm Overall Width 355 mm				
R-HE-2RHY1415C Roho Hybrid Elite Cushion Standard Cover 14" x 15"				
Overall Depth 380 mm Overall Width 355 mm				

Depth	Roho Hybrid	Elite Cushion	Standard Co	over 14" x 1
		405 mm		Overall Width
-2RHY1515C	Roho Hvbrid	Elite Cushion	Standard Co	over 15" x 1
Depth		380 mm		Overall Width
-2RHY1516C	Roho Hybrid	Elite Cushion	Standard Co	over 15" x 1
Depth		405 mm		Overall Width
2RHY1517C	Roho Hybrid	Elite Cushion	Standard Co	ver 15" x 1
Depth		430 mm	otanidara ot	Overall Width
2RHY1615C	Pobo Hybrid	Elite Cushion	Standard Co	over 16" v 1
		380 mm		Overall Width
Depth		300 11111		
-2RHY1616C	Roho Hybrid	Elite Cushion	Standard Co	over 16" x 1
Depth	-	405 mm		Overall Width
-2RHY1617C F	Roho Hybrid	Elite Cushion	Standard Co	over 16" x 1
Depth		430 mm		Overall Width
2RHY1618C		Elito Cuchica	Standard C	Wor 16" v 1
		Elite Cushion	i Stanuard Co	
Depth		455 mm		Overall Width
2RHY1620C	Roho Hybrid	Elite Cushio	n Standard C	over 16" x 2
Depth	,	505 mm		Overall Width
-2RHY1716C F	Roho Hybrid	Elite Cushion	Standard Co	over 17" x 1
Depth		405 mm		Overall Width
2RHY1717C F		Flite Quebier	Chandard Ca	17// 17//
		Elite Cushion	Standard Co	
Depth		430 mm		Overall Width
-2RHY1718C F	Roho Hybrid	Elite Cushion	Standard Co	over 17" x 18
Depth	,	455 mm		Overall Width
2RHY1816C	Roho Hybrid	Elite Cushion	Standard Co	over 18" x 1
Depth		405 mm		Overall Width
-2RHY1817C F		Elite Cushion	Standard Co	wor 10" v 1
		430 mm	Stanuaru CC	Over 18" X 1
Depth		430 mm		overall width
	Roho Hybrid	Elite Cushion	Standard Co	over 18" x 1
-2RHY1818C		455 mm		Overall Width
-2RHY1818C F		400 1111		
		400 mm		
Depth	Roho Hybrid	Elite Cushio	n Standard C	over 18" x 2
Depth	Roho Hybrid		n Standard C	OVEr 18" x 2
Depth •2RHY182OC Depth		Elite Cushion		Overall Width
Depth -2RHY182OC Depth -2RHY2018C		Elite Cushion		Overall Width
Depth •2RHY182OC Depth		Elite Cushion		Overall Width
Depth •2RHY182OC Depth •2RHY2018C Depth	Roho Hybrid	Elite Cushion	n Standard C	Overall Width
		/55 mm		Overall Wi

National Distribution, Local Service

er 15" x 15"	
Iverall Width	380 mm
er 15" x 16"	
lverall Width	380 mm
er 15" x 17"	
Iverall Width	380 mm
10" 15"	
er 16" x 15"	
Iverall Width	405 mm
er 16" x 16"	
Iverall Width	405 mm
	400 11111
er 16" x 17"	
Iverall Width	405 mm
er 16" x 18"	
Iverall Width	405 mm
ver 16" x 20"	
lverall Width	405 mm
er 17" x 16"	
Iverall Width	430 mm
er 17" x 17"	
	400
Iverall Width	430 mm
er 17" x 18"	
Iverall Width	430 mm
er 18" x 16"	
lverall Width	455 mm
er 18" x 17"	
Iverall Width	455 mm
er 18" x 18"	
Iverall Width	455 mm
ver 18" x 20"	
Iverall Width	455 mm
ver 20" x 18"	
	F0F mm
Iverall Width	505 mm
ver 20" x 20"	
Iverall Width	505 mm

355 mm



205

R-HE-2RHY2218C Roho Hybrid Elite Cushion Standard Cover 22" x 18"			
Overall Depth	455 mm	Overall Width	560 mm
R-HE-2RHY2220C Roho Hybrid Elite Cushion Standard Cover 22" x 20"			
Overall Depth	505 mm	Overall Width	560 mm
R-HE-2RHY2418C Roho Hybrid Elite Cushion Standard Cover 24" x 18"			
Overall Depth	455 mm	Overall Width	610 mm
R-HE-2RHY2420C Roho Hybrid Elite Cushion Standard Cover 24" x 20"			
Overall Depth	505 mm		

Roho Nexus Spirit Cushion



- The Roho Nexus Spirit combines the stability of a contoured foam base with the therapeutic care expected from a Roho
 This design allows for increased stability for transfers and for an enhanced sitting posture for the user
- Positioned to address critical pressure points, provides you the highest level of skin protection
- The cushion provides an adjustable low friction and low shear environment to protect the user against wounds
- Lightweight, portable, repairable and washable
 Made from flows resistant paperage rubbs; and polywrethans
- Made from flame resistant neoprene rubber and polyurethane
 A large range of sizes is available from 14"-20" x 14"-18"

R-NEXUS-SPIRIT-1414C	Roho Nexus	Spirit Cushion 14" x 14"			
Overall Depth	355 mm	Overall Width	355 mm	Product Weight	1.4 kg
R-NEXUS-SPIRIT-1416C	Roho Nexus	Spirit Cushion 14" x 16"			
Overall Depth	405 mm	Overall Width	355 mm	Product Weight	1.4 kg
R-NEXUS-SPIRIT-1515C	Roho Nexus	Spirit Cushion 15" x 15"			
Overall Depth	380 mm	Overall Width	380 mm	Product Weight	1.4 kg
R-NEXUS-SPIRIT-1616C	Roho Nexus	Spirit Cushion 16" x 16"			
Overall Depth	405 mm	Overall Width	405 mm	Product Weight	1.4 kg
R-NEXUS-SPIRIT-1618C	Roho Nexus	Spirit Cushion 16" x 18"	_		
Overall Depth	455 mm	Overall Width	405 mm	Product Weight	1.4 kg
R-NEXUS-SPIRIT-1717C	Roho Nexus	Spirit Cushion 17" x 17"			
Overall Depth	430 mm	Overall Width	430 mm	Product Weight	1.4 kg
R-NEXUS-SPIRIT-1816C	Roho Nexus	Spirit Cushion 18" x 16"			_
Overall Depth	405 mm	Overall Width	455 mm	Product Weight	1.4 kg
R-NEXUS-SPIRIT-1818C	Roho Nexus	Spirit Cushion 18" x 18"	_		
Overall Depth	455 mm	Overall Width	455 mm	Product Weight	1.4 kg
R-NEXUS-SPIRIT-20180	Roho Nexu	s Spirit Cushion 20" x 18"			
Overall Depth	455 mm	Overall Width	505 mm	Product Weight	1.4 kg

Roho Toilet Seat Cushio	on	
	 Positions the body secure Made from high quality n Designed for standard siz Maximizes stability of the Fits all standard toilets 	eoprene toilets and por
R-A/SEAT-TOILET Roho Toilet Product Weight	Seat Cushion	
Pressure Cushions No Powered	on	
Varilite ProForm NX Cus	hion Cover	
	 The ProForm NX Cushion i obliquities A modifiable contoured b Available with a mesh or i Available in a range of size 	is recommended ase, superior pre incontinence cov tes 12" x 14" - 20"
1740 Varilite ProForm NX Cus	SNION COVER MESN 1:	380 mm Pa
Overall Depth 380 mm 1741 Varilite ProForm NX Customers		
Overall Depth 430 mm	Overall Width	380 mm Pa
1742 Varilite ProForm NX Cu	shion Cover Mesh 17	7" x 17"
Overall Depth 430 mm	Overall Width	430 mm Pa
1743 Varilite ProForm NX Cu	shion Cover Mesh 2	0" x 16"
Overall Depth 405 mm	Overall Width	505 mm Pa
1744 Varilite ProForm NX Cu	shion Cover Mesh 2	0" x 18"
Overall Depth 455 mm	Overall Width	505 mm Pa
1745 Varilite ProForm NX Cus	shion Cover Incontir	nence 15" x
Overall Depth 380 mm	Overall Width	380 mm Pa
1746 Varilite ProForm NX Cu	shion Cover Incontir	nence 15" :
Overall Depth 430 mm	Overall Width	380 mm Pa
1747 Varilite ProForm NX Cu	shion Cover Incontir	nence 17" x

430 mm Overall Width

430 mm

Overall Depth



oortable bedside commodes otecting from a potential tissue breakdown

0.6 kg

7.05 Pressure Cushion Covers

acement cover for the ProForm NX Cushion ed for individuals with leg length discrepancies, amputations, and other pelvic

pressure distribution, asymmetric positioning options, comfort and simplicity cover with breathable four-way-stretch material 20" x 20"

Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
′ x 15″			
Package Quantity	1	Package Type	Single
' x 17"			
Package Quantity	1	Package Type	Single
x 17"			
Package Quantity	1	Package Type	Single

verall Depth	405 mm	Overall Width	505 mm	Package Quantity	1	Package Type	Single
40	_						
		ushion Cover In					
Overall Depth	455 mm	Overall Width	505 mm	Package Quantity	1	Package Type	Single
750 Varilite Pro	oForm NX C	ushion Cover In	continence 1	2" x 14"			
Overall Depth	355 mm	Overall Width	305 mm	Package Quantity	1	Package Type	Single
7E1 Martilla Dua				AU 1 AU			
Variate Pro	355 mm	UShion Cover Inc	355 mm	4" X 14" Package Quantity	1	Package Type	Single
overan beptit	000 1111	overall width	000 1111	Tackage Quantity		i dokuge type	Jungie
752 Varilite Pro	oForm NX C	ushion Cover Ind	continence 1	4" x 16"			
Overall Depth	405 mm	Overall Width	355 mm	Package Quantity	1	Package Type	Single
753 Varilite Pro	Form NX C	ushion Cover In	continence 1	6" x 16"			
Overall Depth	405 mm	Overall Width	405 mm	Package Quantity	1	Package Type	Single
754 Varilite Pro	oForm NX C	ushion Cover In		6" x 18"			
Overall Depth	455 mm	Overall Width	405 mm	Package Quantity	1	Package Type	Single
755 Varilite Pro	oForm NX C	ushion Cover Inc	continence 1	6" x 20"			
Overall Depth	505 mm	Overall Width	405 mm	Package Quantity	1	Package Type	Single
		ushion Cover In					
Overall Depth	405 mm	Overall Width	455 mm	Package Quantity	1	Package Type	Single
757 Varilite Pro	Form NX C	ushion Cover Inc	continence 1	8" x 18"			
Overall Depth	455 mm	Overall Width	455 mm	Package Quantity	1	Package Type	Single
				Package Quantity	1	Package Type	Single
	oForm NX C	ushion Cover Ind	continence 1	Package Quantity 8" x 20"			
758 Varilite Pro				Package Quantity	1	Package Type Package Type	Single Single
758 Varilite Pro	DForm NX C	ushion Cover Ind	continence 1 455 mm	Package Quantity 8" x 20" Package Quantity			
758 Varilite Pro Overall Depth 759 Varilite Pro	DForm NX C	Ushion Cover Ind	continence 1 455 mm	Package Quantity 8" x 20" Package Quantity			
758 Varilite Pro overall Depth 759 Varilite Pro overall Depth	DForm NX C 505 mm DForm NX C 505 mm	ushion Cover Ind Overall Width ushion Cover Ind Overall Width	continence 1 455 mm continence 2 505 mm	Package Quantity 8" x 20" Package Quantity 20" x 20" Package Quantity		Package Type	Single
758 Varilite Pro overall Depth 759 Varilite Pro overall Depth 760 Varilite Pro	DForm NX C 505 mm DForm NX C 505 mm DForm NX C	ushion Cover Ind Overall Width ushion Cover Ind Overall Width ushion Cover Me	continence 1 455 mm continence 2 505 mm esh 12" x 14"	Package Quantity 8" x 20" Package Quantity 20" x 20" Package Quantity		Package Type Package Type	Single
758 Varilite Pro overall Depth 759 Varilite Pro overall Depth 760 Varilite Pro	DForm NX C 505 mm DForm NX C 505 mm	ushion Cover Ind Overall Width ushion Cover Ind Overall Width	continence 1 455 mm continence 2 505 mm	Package Quantity 8" x 20" Package Quantity 20" x 20" Package Quantity	1	Package Type	Single
758 Varilite Pro overall Depth 759 Varilite Pro overall Depth 760 Varilite Pro Overall Depth	DForm NX C 505 mm DForm NX C 505 mm DForm NX C 355 mm	ushion Cover Ind Overall Width ushion Cover Ind Overall Width ushion Cover Me	continence 1 455 mm continence 2 505 mm esh 12" x 14" 305 mm	Package Quantity 8" x 20" Package Quantity 20" x 20" Package Quantity	1	Package Type Package Type	Single
758 Varilite Pro overall Depth 759 Varilite Pro overall Depth 760 Varilite Pro Overall Depth	DForm NX C 505 mm DForm NX C 505 mm DForm NX C 355 mm	ushion Cover Ind Overall Width ushion Cover Ind Overall Width ushion Cover Mo Overall Width	continence 1 455 mm continence 2 505 mm esh 12" x 14" 305 mm	Package Quantity 8" x 20" Package Quantity 20" x 20" Package Quantity	1	Package Type Package Type	Single
758 Varilite Pro Overall Depth 759 Varilite Pro Overall Depth 760 Varilite Pro Overall Depth 761 Varilite Pro Overall Depth	DForm NX C 505 mm DForm NX C 505 mm DForm NX C 355 mm DForm NX Cu 355 mm	ushion Cover Ind Overall Width ushion Cover Ind Overall Width ushion Cover Me Overall Width ushion Cover Me Overall Width	continence 1 455 mm continence 2 505 mm esh 12" x 14" 305 mm esh 14" x 14" 355 mm	Package Quantity 8" x 20" Package Quantity 20" x 20" Package Quantity Package Quantity Package Quantity Package Quantity	1 1 1	Package Type Package Type Package Type	Single Single Single
758 Varilite Pro overall Depth 759 Varilite Pro Overall Depth 760 Varilite Pro Overall Depth 761 Varilite Pro Overall Depth 762 Varilite Pro	DForm NX C 505 mm DForm NX C 505 mm DForm NX C 355 mm DForm NX Cu 355 mm	ushion Cover Ind Overall Width ushion Cover Ind Overall Width ushion Cover Me Overall Width	continence 1 455 mm continence 2 505 mm esh 12" x 14" 305 mm esh 14" x 14" 355 mm	Package Quantity 8" x 20" Package Quantity 20" x 20" Package Quantity Package Quantity Package Quantity Package Quantity	1 1 1	Package Type Package Type Package Type Package Type Package Type	Single Single Single
758 Varilite Pro Overall Depth 759 Varilite Pro Overall Depth 760 Varilite Pro Overall Depth 761 Varilite Pro Overall Depth 762 Varilite Pro	DForm NX C 505 mm DForm NX C 505 mm DForm NX C 355 mm DForm NX C 355 mm	ushion Cover Ind Overall Width ushion Cover Ind Overall Width ushion Cover Me Overall Width ushion Cover Me Overall Width	continence 1 455 mm continence 2 505 mm esh 12" x 14" 305 mm esh 14" x 14" 355 mm esh 14" x 16"	Package Quantity 8" x 20" Package Quantity 20" x 20" Package Quantity Package Quantity Package Quantity Package Quantity		Package Type Package Type Package Type	Single Single Single Single
758 Varilite Pro Overall Depth 759 Varilite Pro Overall Depth 760 Varilite Pro Overall Depth 761 Varilite Pro Overall Depth 762 Varilite Pro Overall Depth	DForm NX C 505 mm DForm NX C 505 mm DForm NX C 355 mm DForm NX C 355 mm DForm NX C 405 mm	ushion Cover Ind Overall Width ushion Cover Ind Overall Width ushion Cover Me Overall Width ushion Cover Me Overall Width	continence 1 455 mm continence 2 505 mm esh 12" x 14" 305 mm esh 14" x 14" 355 mm esh 14" x 16" 355 mm	Package Quantity 8" x 20" Package Quantity 20" x 20" Package Quantity		Package Type Package Type Package Type Package Type Package Type	Single Single Single Single
758 Varilite Pro overall Depth 759 Varilite Pro Overall Depth 760 Varilite Pro Overall Depth 761 Varilite Pro Overall Depth 762 Varilite Pro Overall Depth 763 Varilite Pro	DForm NX C 505 mm DForm NX C 505 mm DForm NX C 355 mm DForm NX C 355 mm DForm NX C 405 mm	ushion Cover Ind Overall Width ushion Cover Ind Overall Width ushion Cover Me Overall Width ushion Cover Me Overall Width	continence 1 455 mm continence 2 505 mm esh 12" x 14" 305 mm esh 14" x 14" 355 mm esh 14" x 16" 355 mm	Package Quantity 8" x 20" Package Quantity 20" x 20" Package Quantity		Package Type Package Type Package Type Package Type Package Type	Single Single Single Single
758 Varilite Pro Overall Depth 759 Varilite Pro Overall Depth 760 Varilite Pro Overall Depth 761 Varilite Pro Overall Depth 762 Varilite Pro Overall Depth 763 Varilite Pro Overall Depth	DForm NX C 505 mm DForm NX C 505 mm DForm NX C 355 mm DForm NX C 355 mm DForm NX C 405 mm DForm NX C	ushion Cover Ind Overall Width Ushion Cover Ind Overall Width Ushion Cover Me Overall Width Ushion Cover Me Overall Width Ushion Cover Me Overall Width	continence 1 455 mm continence 2 505 mm esh 12" x 14" 305 mm esh 14" x 14" 355 mm esh 14" x 16" 355 mm esh 16" x 16" 405 mm	Package Quantity 8" x 20" Package Quantity 20" x 20" Package Quantity		Package Type Package Type Package Type Package Type Package Type Package Type	Single Single Single Single Single Single
758 Varilite Pro Overall Depth 759 Varilite Pro Overall Depth 760 Varilite Pro Overall Depth 761 Varilite Pro Overall Depth 762 Varilite Pro Overall Depth 763 Varilite Pro Overall Depth 764 Varilite Pro	DForm NX C 505 mm DForm NX C 505 mm DForm NX C 355 mm DForm NX C 355 mm DForm NX C 405 mm DForm NX C	ushion Cover Ind Overall Width ushion Cover Ind Overall Width ushion Cover Me Overall Width ushion Cover Me Overall Width ushion Cover Me	continence 1 455 mm continence 2 505 mm esh 12" x 14" 305 mm esh 14" x 14" 355 mm esh 14" x 16" 355 mm esh 16" x 16" 405 mm	Package Quantity 8" x 20" Package Quantity 20" x 20" Package Quantity		Package Type Package Type Package Type Package Type Package Type Package Type	Single Single Single Single Single Single
758 Varilite Pro Overall Depth 759 Varilite Pro 759 Varilite Pro Overall Depth 760 Varilite Pro 760 Varilite Pro Overall Depth 761 Varilite Pro Overall Depth 762 Varilite Pro Overall Depth 763 Varilite Pro Overall Depth 763 Varilite Pro Overall Depth 763 Varilite Pro 764 Varilite Pro	DForm NX C 505 mm DForm NX C 505 mm DForm NX C 355 mm DForm NX C 355 mm DForm NX C 405 mm DForm NX C	ushion Cover Ind Overall Width ushion Cover Ind Overall Width ushion Cover Me Overall Width ushion Cover Me Overall Width ushion Cover Me Overall Width	continence 1 455 mm continence 2 505 mm esh 12" x 14" 305 mm esh 14" x 14" 355 mm esh 14" x 16" 355 mm esh 16" x 16" 405 mm esh 16" x 18"	Package Quantity 8" x 20" Package Quantity 20" x 20" Package Quantity		Package Type Package Type Package Type Package Type Package Type Package Type	Single Single Single Single Single Single Single Single
758 Varilite Pro overall Depth 759 Varilite Pro overall Depth 760 Varilite Pro overall Depth 761 Varilite Pro overall Depth 762 Varilite Pro overall Depth 763 Varilite Pro overall Depth 764 Varilite Pro overall Depth	DForm NX C 505 mm DForm NX C 505 mm DForm NX C 355 mm DForm NX C 355 mm DForm NX C 405 mm DForm NX C 405 mm DForm NX C	ushion Cover Ind Overall Width ushion Cover Ind Overall Width ushion Cover Me Overall Width ushion Cover Me Overall Width ushion Cover Me Overall Width	continence 1 455 mm continence 2 505 mm esh 12" x 14" 305 mm esh 14" x 14" 355 mm esh 14" x 16" 355 mm esh 14" x 16" 355 mm esh 16" x 16" 405 mm esh 16" x 18" 405 mm	Package Quantity 8" x 20" Package Quantity 20" x 20" Package Quantity Package Quantity		Package Type Package Type Package Type Package Type Package Type Package Type	Single Single Single Single Single Single Single Single

1700				
1766 Varilite ProFor	105 mm	Shion Cover Mesh 1	8" x 16" 455 mm	Pi
overall bepar	100 11111		400 11111	
1767 Varilite ProFor	m NX Cu	shion Cover Mesh 18	8" x 18"	_
Overall Depth	455 mm	Overall Width	455 mm	Pa
1768 Varilite ProFor	m NX Cu	shion Cover Mesh 1	B" x 20"	
Overall Depth	505 mm	Overall Width	455 mm	Pa
1769 Varilite ProFor	rm NX Cu	shion Cover Mesh 2	0" x 20"	
Overall Depth	505 mm	Overall Width	505 mm	Pa
Varilite Evolution	& Evolu	ution PSV Bariatri • Cushion cover to replace • Suitable for use with Evo • The Evolution Cushion pr • Available in a four-way-st • Available in a range of siz	original cover lution PSV Bar for users with ovides outsta rretch mesh ou	r iat n a ind r in
4626 Varilite Evolut Overall Depth 4628 Varilite Evolut	455 mm	volution PSV Cushion Overall Width	610 mm	Pa
Overall Depth	505 mm	Overall Width	610 mm	/IE
			0.0	
Varilite Meridian	Cushi	on Cover		
10 10		 The Meridian Cushion Co The Meridian Cushion fee Available in a mesh or inco Available in a large variet 	tures dual ch continence co	nan ver
2521 Varilite Meridi	ian Cush	ion Cover Mesh 16"	x 18″	
Overall Depth	455 mm	Overall Width	405 mm	P
2522 Varilite Merid	ian Cusł	nion Cover Mesh 16"	x 20″	
Overall Depth	505 mm	Overall Width	405 mm	Pa
2523 Varilite Merid	ian Cush	nion Cover Mesh 17"	x 17″	
Overall Depth	430 mm	Overall Width	430 mm	Pa
2524 Varilite Merid				
LUL- Variate Meria	ian Cusr	ion Cover Mesh 18"	x 16″	
Overall Depth	405 mm	Overall Width	x 16" 455 mm	Pa

2525 Varilite Meridian Cushion Cover Mesh 18" x 18"

455 mm Overall Width

455 mm

Overall Depth

Pressure Cushions



Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single

on Cover

atric Wheelchair Cushion

a high risk of tissue breakdown nding pressure distribution and postural support r incontinence cover

8"-20"

esh 24" x 18"			
Package Quantity	1	Package Type	Single
esh 24" x 20"			
Package Quantity	1	Package Type	Single

ement cover for the Meridian Cushion

nambers , air -foam floatation valves for adjusting air volume ver with breathable two-way stretch material n 14" x 14" - 20" x 20"

Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single

	505 mm	Overall Width	455 mm	Package Quantity	1	Package Type	Single
						,	
527 Varilite Meri	dian Cus	hion Cover Mesh	20" x 16"				
Overall Depth	405 mm	Overall Width	505 mm	Package Quantity	1	Package Type	Single
528 Varilite Meri	idian Cus	hion Cover Mesh	20" x 18"				
Overall Depth	455 mm	Overall Width	505 mm	Package Quantity	1	Package Type	Single
2529 Varilite Meri	_	hion Cover Mesh	20" x 20"				
Overall Depth	505 mm	Overall Width	505 mm	Package Quantity	1	Package Type	Single
2534 Varilite Meri	idian Cus	hion Cover Incont	tinence 12"	x 14"			
Overall Depth	355 mm	Overall Width	305 mm	Package Quantity	1	Package Type	Single
		hion Cover Incont					
Overall Depth	355 mm	Overall Width	355 mm	Package Quantity	1	Package Type	Single
2536 Varilite Meri	idian Cus	hion Cover Incont	tinence 14"	x 16″			
Overall Depth	405 mm	Overall Width	355 mm	Package Quantity	1	Package Type	Single
2537 Varilite Meri					_		
Overall Depth	405 mm	Overall Width	405 mm	Package Quantity	1	Package Type	Single
2538 Varilite Meri	idian Cus	hion Cover Incont	tinence 16"	x 18″			
Overall Depth	455 mm	Overall Width	405 mm	Package Quantity	1	Package Type	Single
2539 Varilite Meri							
Overall Depth	505 mm	Overall Width	405 mm	Package Quantity	1	Package Type	Single
2540 Varilite Meri	idian Cus	hion Cover Incont	tinence 17"	x 17″			
2540 Varilite Meri Overall Depth	idian Cus 430 mm	hion Cover Incont	tinence 17" 430 mm	X 17" Package Quantity	1	Package Type	Single
Overall Depth	430 mm	Overall Width	430 mm	Package Quantity	1	Package Type	Single
Overall Depth 2541 Varilite Meric	430 mm	Overall Width	430 mm	Package Quantity x 16"			
Overall Depth	430 mm	Overall Width	430 mm	Package Quantity	1	Package Type Package Type	Single Single
Overall Depth 2541 Varilite Merio Overall Depth	430 mm dian Cus 405 mm	Overall Width hion Cover Incont Overall Width	430 mm inence 18" 455 mm	Package Quantity X 16" Package Quantity			
Overall Depth 2541 Varilite Meric	430 mm dian Cus 405 mm	Overall Width hion Cover Incont Overall Width	430 mm inence 18" 455 mm	Package Quantity X 16" Package Quantity			
Overall Depth 2541 Varilite Meric Overall Depth 2542 Varilite Meri Overall Depth	430 mm dian Cus 405 mm dian Cus 455 mm	Overall Width hion Cover Incont Overall Width whion Cover Incont Overall Width Overall Width	430 mm inence 18" 455 mm tinence 18" 455 mm	Package Quantity x 16" Package Quantity x 18" Package Quantity	1	Package Type	Single
Overall Depth 2541 Varilite Meric Overall Depth 2542 Varilite Meri Overall Depth 2543 Varilite Meri	430 mm dian Cus 405 mm idian Cus idian Cus	Overall Width hion Cover Incont Overall Width hion Cover Incont Overall Width thion Cover Incont	430 mm inence 18" 455 mm tinence 18" 455 mm tinence 18"	Package Quantity x 16" Package Quantity x 18" Package Quantity x 20"	1	Package Type Package Type	Single
Overall Depth 2541 Varilite Meric Overall Depth 2542 Varilite Meri Overall Depth	430 mm dian Cus 405 mm dian Cus 455 mm	Overall Width hion Cover Incont Overall Width whion Cover Incont Overall Width Overall Width	430 mm inence 18" 455 mm tinence 18" 455 mm	Package Quantity x 16" Package Quantity x 18" Package Quantity	1	Package Type	Single
Overall Depth 2541 Varilite Merid Overall Depth 2542 Varilite Meri Overall Depth 2543 Varilite Meri Overall Depth	430 mm dian Cus 405 mm idian Cus 455 mm idian Cus 505 mm	Overall Width hion Cover Incont Overall Width hion Cover Incont Overall Width Overall Width chion Cover Incont Overall Width Overall Width	430 mm inence 18" 455 mm tinence 18" 455 mm tinence 18" 455 mm	Package Quantity x 16" Package Quantity x 18" Package Quantity x 20" Package Quantity	1	Package Type Package Type	Single
Overall Depth 2541 Varilite Merid Overall Depth 2542 Varilite Meri Overall Depth 2543 Varilite Meri Overall Depth	430 mm dian Cus 405 mm idian Cus 455 mm idian Cus 505 mm	Overall Width hion Cover Incont Overall Width hion Cover Incont Overall Width Overall Width chion Cover Incont Overall Width Overall Width	430 mm inence 18" 455 mm tinence 18" 455 mm tinence 18" 455 mm	Package Quantity x 16" Package Quantity x 18" Package Quantity x 20" Package Quantity	1	Package Type Package Type	Single
Overall Depth 2541 Varilite Merid Overall Depth 2542 Varilite Merid Overall Depth 2543 Varilite Merid Overall Depth 2544 Varilite Merid Overall Depth	430 mm dian Cus 405 mm idian Cus 455 mm idian Cus 505 mm idian Cus 405 mm	Overall Width hion Cover Incont Overall Width chion Cover Incont Overall Width	430 mm inence 18" 455 mm tinence 18" 455 mm tinence 18" 455 mm tinence 20' 505 mm	Package Quantity x 16" Package Quantity x 18" Package Quantity x 20" Package Quantity ' x 16" Package Quantity		Package Type Package Type Package Type	Single Single Single
Overall Depth 2541 Varilite Merin Overall Depth 2542 Varilite Merin Overall Depth 2543 Varilite Merin Overall Depth 2544 Varilite Merin Overall Depth 2545 Varilite Merin	430 mm dian Cus 405 mm dian Cus 505 mm dian Cus 405 mm dian Cus	Overall Width hion Cover Incont Overall Width whion Cover Incont Overall Width hion Cover Incont Overall Width hion Cover Incont Overall Width hion Cover Incont	430 mm inence 18" 455 mm tinence 18" 455 mm tinence 18" 455 mm tinence 20' 505 mm tinence 20'	Package Quantity x 16" Package Quantity x 18" Package Quantity x 20" Package Quantity x 16" Package Quantity ' x 16" Package Quantity		Package Type Package Type Package Type Package Type Package Type	Single Single Single Single Single Single
Overall Depth 2541 Varilite Merid Overall Depth 2542 Varilite Merid Overall Depth 2543 Varilite Merid Overall Depth 2544 Varilite Merid Overall Depth	430 mm dian Cus 405 mm idian Cus 455 mm idian Cus 505 mm idian Cus 405 mm	Overall Width hion Cover Incont Overall Width chion Cover Incont Overall Width	430 mm inence 18" 455 mm tinence 18" 455 mm tinence 18" 455 mm tinence 20' 505 mm	Package Quantity x 16" Package Quantity x 18" Package Quantity x 20" Package Quantity ' x 16" Package Quantity		Package Type Package Type Package Type	Single Single Single
Overall Depth 2541 Varilite Meria Overall Depth 2542 Varilite Meria Overall Depth 2543 Varilite Meria Overall Depth 2544 Varilite Meria Overall Depth 2545 Varilite Meria Overall Depth	430 mm dian Cus 405 mm idian Cus idian Cus 505 mm idian Cus	Overall Width hion Cover Incont Overall Width whion Cover Incont Overall Width hion Cover Incont Overall Width hion Cover Incont Overall Width hion Cover Incont	430 mm inence 18" 455 mm tinence 18" 455 mm tinence 18" 455 mm tinence 20' 505 mm tinence 20' 505 mm	Package Quantity x 16" Package Quantity x 18" Package Quantity x 20" Package Quantity ' x 16" Package Quantity ' x 16" Package Quantity ' x 18" Package Quantity		Package Type Package Type Package Type Package Type Package Type	Single Single Single Single Single Single
Overall Depth 2541 Varilite Meria Overall Depth 2542 Varilite Meria Overall Depth 2543 Varilite Meria Overall Depth 2544 Varilite Meria Overall Depth 2545 Varilite Meria Overall Depth	430 mm dian Cus 405 mm idian Cus idian Cus 505 mm idian Cus	Overall Width hion Cover Incont Overall Width chion Cover Incont Overall Width	430 mm inence 18" 455 mm tinence 18" 455 mm tinence 18" 455 mm tinence 20' 505 mm tinence 20' 505 mm	Package Quantity x 16" Package Quantity x 18" Package Quantity x 20" Package Quantity ' x 16" Package Quantity ' x 16" Package Quantity ' x 18" Package Quantity		Package Type Package Type Package Type Package Type Package Type	Single Single Single Single Single Single
Overall Depth 2541 Varilite Merid Overall Depth 2542 Varilite Merid Overall Depth 2543 Varilite Merid Overall Depth 2544 Varilite Merid Overall Depth 2545 Varilite Merid Overall Depth 2546 Varilite Merid Overall Depth	430 mm dian Cus 405 mm idian Cus idian Cus	Overall Width hion Cover Incont Overall Width chion Cover Incont overall Width	430 mm inence 18" 455 mm tinence 18" 455 mm tinence 18" 455 mm tinence 20' 505 mm tinence 20' 505 mm tinence 20' 505 mm tinence 20' 505 mm	Package Quantity x 16" Package Quantity x 18" Package Quantity x 20" Package Quantity ' x 16" Package Quantity ' x 16" Package Quantity ' x 18" Package Quantity ' x 18" Package Quantity ' x 20"		Package Type Package Type Package Type Package Type Package Type Package Type	Single Single Single Single Single Single Single

2620 Varilite Merid	ian Cush	ion Cover Mesh 14"	x 20″	
Overall Depth	505 mm	Overall Width	355 mm	P
2621 Varilite Meridi	an Cush	ion Cover Mesh 15"	x 15″	
Overall Depth	380 mm	Overall Width	380 mm	P
2622 Varilite Merid	ian Cush	ion Cover Mesh 15"	x 17″	
Overall Depth	430 mm	Overall Width	380 mm	P
2623 Varilite Merid	ian Cush	ion Cover Mesh 16"	x 14″	
Overall Depth	355 mm	Overall Width	405 mm	P
2624 Varilite Merid	ian Cush	ion Cover Mesh 14"	x 18″	
Overall Depth	455 mm	Overall Width	355 mm	P
2625 Varilite Merid	ian Cush	ion Cover Mesh 14"	x 20″	
Overall Depth	505 mm	Overall Width	355 mm	P
2626 Varilite Merid	ian Cush	ion Cover Mesh 15"	x 15″	
				_
Overall Depth	380 mm	Overall Width	380 mm	P
				P
Overall Depth 2627 Varilite Meridi Overall Depth				
2627 Varilite Meridi	ian Cush	ion Cover Mesh 15" Overall Width	x 17" 380 mm	P

Varilite Evolution & Evolution PSV Cushion Cover



• Cushion cover to replace original cover Suitable for use with Evolution PSV Wheelchair Cushion

- Available in a two-way-stretch mesh or incontinence cover
 Available in a range of sizes 12"-14" x 20"-20"

4602 Varilite Evolu	tion & Ev	olution PSV Cushior	n Cover I	ncontinence 16" x 2	0"		
Overall Depth	505 mm	Overall Width	405 mm	Package Quantity	1	Package Type	Single
4605 Varilite Evolu	tion & E\	volution PSV Cushior	n Cover I	Mesh 20" x 16"			
Overall Depth	405 mm	Overall Width	505 mm	Package Quantity	1	Package Type	Single
4599 Varilite Evolu Overall Depth	tion & Ev 505 mm	volution PSV Cushior	505 mm	Mesh 20" x 20" Package Quantity	1	Package Type	Single
4600 Varilite Evolu	tion & Ev	volution PSV Cushio	n Cover I	ncontinence 20" x 2	20″	I	
Overall Depth	505 mm	Overall Width	505 mm	Package Quantity	1	Package Type	Single
4601 Varilite Evolut	ion & Ev	olution PSV Cushior	n Cover N	/lesh 16" x 20"			
Overall Depth	505 mm	Overall Width	405 mm	Package Quantity	1	Package Type	Single



Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Раскауе туре	Sillyle
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Deckage Quentity	1	Dookogo Tupo	Cingle
Package Quantity		Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
rackaye Quantity		гаскаде туре	Jillyle
Package Quantity	1	Package Type	Single

• The Evolution Cushion is for users with a high risk of tissue breakdown

• The Evolution Cushion provides outstanding pressure distribution and postural support

verall Depth	405 mm	Overall Width	505 mm	Package Quantity	1	Package Type	Single
607 Varilite Evo	olution & E	volution PSV Cu	shion Cover	Mesh 20" x 18"			
Overall Depth	455 mm	Overall Width	505 mm	Package Quantity	1	Package Type	Single
608 Varilite Evo	olution & E	volution PSV Cu	shion Cover	Incontinence 20	' x 18″		
Overall Depth	455 mm	Overall Width	505 mm	Package Quantity	1	Package Type	Single
620 Varilite Evo	olution & E	volution PSV Cu	shion Cover	Mesh 22" x 18"			
Overall Depth	455 mm	Overall Width	560 mm	Package Quantity	1	Package Type	Single
622 Varilite Evo	olution & E	volution PSV Cu	shion Cover	Incontinence 22'	′ x 18″		
Overall Depth	455 mm	Overall Width	560 mm	Package Quantity	1	Package Type	Single
623 Varilite Evo	olution & E	volution PSV Cu	shion Cover	Mesh 12" x 14"			
Overall Depth	355 mm	Overall Width	305 mm	Package Quantity	1	Package Type	Single
624 Varilite Evo	olution & F	volution PSV Cu	shion Cover	Incontinence 12"	x 14"		
Overall Depth	355 mm	Overall Width	305 mm	Package Quantity	14	Package Type	Single
	000 mm		300 mm	. annage tauning		. actuage type	
640 Varilite Evo	olution & E	volution PSV Cu	shion Cover	Mesh 17" x 17"			
Overall Depth	430 mm	Overall Width	430 mm	Package Quantity	1	Package Type	Single
641 Varilite Evo					17/		
641 Varilite Evo							
				Incontinence 17"		Deckers Ture	Cinala
	430 mm	Overall Width	430 mm	Package Quantity	X I / "	Package Type	Single
Overall Depth	430 mm		430 mm	Package Quantity		Package Type	Single
Overall Depth	430 mm	Overall Width	430 mm	Package Quantity		Package Type Package Type	Single
Overall Depth	430 mm	Overall Width	430 mm	Package Quantity Mesh 14" x 14"	1		
Overall Depth 1643 Varilite Evo Overall Depth	430 mm	Overall Width	430 mm shion Cover 355 mm	Package Quantity Mesh 14" x 14"	1		
Overall Depth 643 Varilite Evo Overall Depth 644 Varilite Evo	430 mm	Overall Width	430 mm shion Cover 355 mm	Package Quantity Mesh 14" x 14" Package Quantity	1		
Overall Depth 643 Varilite Evo Overall Depth 644 Varilite Evo Overall Depth	430 mm blution & E 355 mm blution & E 355 mm	Overall Width	430 mm shion Cover 355 mm shion Cover 355 mm	Package Quantity Mesh 14" x 14" Package Quantity Incontinence 14" Package Quantity	1 1 x 14"	Package Type	Single
Overall Depth 643 Varilite Evo Overall Depth 644 Varilite Evo Overall Depth 645 Varilite Evo	430 mm blution & E 355 mm blution & E 355 mm	Overall Width	430 mm shion Cover 355 mm shion Cover 355 mm	Package Quantity Mesh 14" x 14" Package Quantity Incontinence 14" Package Quantity	1 1 x 14"	Package Type	Single
Overall Depth 643 Varilite Evo Overall Depth 644 Varilite Evo Overall Depth 645 Varilite Evo	430 mm blution & E 355 mm blution & E blution & E	Overall Width	430 mm shion Cover 355 mm shion Cover 355 mm shion Cover	Package Quantity Mesh 14" x 14" Package Quantity Incontinence 14" Package Quantity Mesh 14" x 16"	1 1 x 14"	Package Type Package Type	Single
Overall Depth 643 Varilite Evo Overall Depth 644 Varilite Evo Overall Depth 645 Varilite Evo Overall Depth	430 mm blution & E 355 mm blution & E 355 mm blution & E 405 mm	Overall Width	430 mm shion Cover 355 mm shion Cover 355 mm shion Cover 355 mm	Package Quantity Mesh 14" x 14" Package Quantity Incontinence 14" Package Quantity Mesh 14" x 16"	1 x14" 1 1	Package Type Package Type	Single
Overall Depth 643 Varilite Evo Overall Depth 644 Varilite Evo Overall Depth 645 Varilite Evo Overall Depth 646 Varilite Evo	430 mm blution & E 355 mm blution & E 355 mm blution & E 405 mm	Overall Width	430 mm shion Cover 355 mm shion Cover 355 mm shion Cover 355 mm	Package Quantity Mesh 14" x 14" Package Quantity Incontinence 14" Package Quantity Mesh 14" x 16" Package Quantity	1 x14" 1 1	Package Type Package Type	Single
Overall Depth 643 Varilite Evo Overall Depth 644 Varilite Evo Overall Depth 645 Varilite Evo Overall Depth 646 Varilite Evo Overall Depth	430 mm blution & E 355 mm blution & E 405 mm blution & E 405 mm	Overall Width Evolution PSV Custom Overall Width	430 mm shion Cover 355 mm shion Cover 355 mm shion Cover 355 mm shion Cover 355 mm	Package Quantity Mesh 14" x 14" Package Quantity Incontinence 14" Package Quantity Mesh 14" x 16" Package Quantity Incontinence 14" Package Quantity Package Quantity	1 x 14" 1 1 x 16"	Package Type Package Type Package Type	Single Single Single
overall Depth 1643 Varilite Evo Overall Depth 1644 Varilite Evo Overall Depth 1645 Varilite Evo Overall Depth 1646 Varilite Evo Overall Depth 1665 Varilite Evo	430 mm blution & E 355 mm blution & E 405 mm blution & E 405 mm blution & E	Overall Width Evolution PSV Custor Overall Width	430 mm shion Cover 355 mm shion Cover 355 mm shion Cover 355 mm shion Cover	Package Quantity Mesh 14" x 14" Package Quantity Incontinence 14" Package Quantity Mesh 14" x 16" Package Quantity Incontinence 14" Package Quantity Mesh 16" x 16"	1 x 14" 1 1 x 16" 1	Package Type Package Type Package Type Package Type	Single Single Single Single
Overall Depth 1643 Varilite Evo Overall Depth 1644 Varilite Evo Overall Depth 1645 Varilite Evo Overall Depth 1646 Varilite Evo Overall Depth 1665 Varilite Evo	430 mm blution & E 355 mm blution & E 405 mm blution & E 405 mm	Overall Width Evolution PSV Custom Overall Width	430 mm shion Cover 355 mm shion Cover 355 mm shion Cover 355 mm shion Cover 355 mm	Package Quantity Mesh 14" x 14" Package Quantity Incontinence 14" Package Quantity Mesh 14" x 16" Package Quantity Incontinence 14" Package Quantity Package Quantity	1 x 14" 1 1 x 16"	Package Type Package Type Package Type	Single Single Single
Overall Depth 643 Varilite Evo Overall Depth 644 Varilite Evo Overall Depth 645 Varilite Evo Overall Depth 646 Varilite Evo Overall Depth 665 Varilite Evo Overall Depth	430 mm blution & E 355 mm blution & E 405 mm blution & E 405 mm blution & E 405 mm	Overall Width Evolution PSV Current Overall Width	430 mm shion Cover 355 mm shion Cover 405 mm	Package Quantity Mesh 14" x 14" Package Quantity Incontinence 14" Package Quantity Mesh 14" x 16" Package Quantity Incontinence 14" Package Quantity Mesh 16" x 16"	1 x 14" 1 1 x 16" 1 1	Package Type Package Type Package Type Package Type	Single Single Single Single
Overall Depth 643 Varilite Evo Overall Depth 644 Varilite Evo Overall Depth 645 Varilite Evo Overall Depth 646 Varilite Evo Overall Depth 645 Varilite Evo Overall Depth 665 Varilite Evo Overall Depth 665 Varilite Evo Overall Depth 665 Varilite Evo Overall Depth	430 mm blution & E 355 mm blution & E 405 mm blution & E 405 mm blution & E 405 mm	Overall Width Evolution PSV Current Overall Width	430 mm shion Cover 355 mm shion Cover 405 mm	Package Quantity Mesh 14" x 14" Package Quantity Incontinence 14" Package Quantity Mesh 14" x 16" Package Quantity Incontinence 14" Package Quantity Mesh 16" x 16" Package Quantity	1 x 14" 1 1 x 16" 1 1	Package Type Package Type Package Type Package Type	Single Single Single Single
Overall Depth 643 Varilite Evo Overall Depth 644 Varilite Evo Overall Depth 645 Varilite Evo Overall Depth 646 Varilite Evo Overall Depth 645 Varilite Evo Overall Depth 665 Varilite Evo Overall Depth 665 Varilite Evo Overall Depth 666 Varilite Evo Overall Depth	430 mm blution & E 355 mm blution & E 355 mm blution & E 405 mm	Overall Width Evolution PSV Custor Overall Width	430 mm shion Cover 355 mm shion Cover 405 mm shion Cover 405 mm shion Cover 405 mm	Package Quantity Mesh 14" x 14" Package Quantity Incontinence 14" Package Quantity Mesh 14" x 16" Package Quantity Incontinence 14" Package Quantity Mesh 16" x 16" Package Quantity Incontinence 16" Package Quantity	1 x 14" 1 x 16" 1 x 16" 1 x 16"	Package Type Package Type Package Type Package Type Package Type Package Type	Single Single Single Single Single Single
Overall Depth 643 Varilite Evo Overall Depth 644 Varilite Evo Overall Depth 645 Varilite Evo Overall Depth 646 Varilite Evo Overall Depth 645 Varilite Evo Overall Depth 665 Varilite Evo Overall Depth 665 Varilite Evo Overall Depth 666 Varilite Evo Overall Depth	430 mm blution & E 355 mm blution & E 355 mm blution & E 405 mm	Overall Width Svolution PSV Custor Overall Width	430 mm shion Cover 355 mm shion Cover 405 mm shion Cover 405 mm shion Cover 405 mm	Package Quantity Mesh 14" x 14" Package Quantity Incontinence 14" Package Quantity Mesh 14" x 16" Package Quantity Incontinence 14" Package Quantity Mesh 16" x 16" Package Quantity Incontinence 16" Package Quantity	1 x 14" 1 x 16" 1 x 16" 1 x 16"	Package Type Package Type Package Type Package Type Package Type Package Type	Single Single Single Single Single Single
Overall Depth 643 Varilite Evo Overall Depth 644 Varilite Evo Overall Depth 645 Varilite Evo Overall Depth 646 Varilite Evo Overall Depth 645 Varilite Evo Overall Depth 665 Varilite Evo Overall Depth 666 Varilite Evo Overall Depth	430 mm blution & E 355 mm blution & E 355 mm blution & E 405 mm	Overall Width Evolution PSV Custor Overall Width	430 mm shion Cover 355 mm shion Cover 405 mm shion Cover 405 mm shion Cover 405 mm	Package Quantity Mesh 14" x 14" Package Quantity Incontinence 14" Package Quantity Mesh 14" x 16" Package Quantity Incontinence 14" Package Quantity Mesh 16" x 16" Package Quantity Incontinence 16" Package Quantity	1 x 14" 1 x 16" 1 x 16" 1 x 16"	Package Type Package Type Package Type Package Type Package Type Package Type	Single Single Single Single Single Single
Overall Depth 643 Varilite Evo Overall Depth 644 Varilite Evo Overall Depth 645 Varilite Evo Overall Depth 646 Varilite Evo Overall Depth 645 Varilite Evo Overall Depth 665 Varilite Evo Overall Depth 665 Varilite Evo Overall Depth 666 Varilite Evo Overall Depth 666 Varilite Evo Overall Depth 666 Varilite Evo Overall Depth 667 Varilite Evo Overall Depth	430 mm blution & E 355 mm blution & E 355 mm blution & E 405 mm	Overall Width	430 mm shion Cover 355 mm shion Cover 405 mm shion Cover 405 mm shion Cover 405 mm shion Cover 405 mm	Package Quantity Mesh 14" x 14" Package Quantity Incontinence 14" Package Quantity Mesh 14" x 16" Package Quantity Mesh 14" x 16" Package Quantity Incontinence 14" Package Quantity Incontinence 14" Package Quantity Incontinence 16" Package Quantity Incontinence 16" Package Quantity Mesh 16" x 18" Package Quantity	1 x 14" 1 x 16" 1 x 16" 1	Package Type Package Type Package Type Package Type Package Type Package Type	Single Single Single Single Single Single Single Single
overall Depth 1643 Varilite Evo Overall Depth 1644 Varilite Evo Overall Depth 1645 Varilite Evo Overall Depth 1665 Varilite Evo Overall Depth 1666 Varilite Evo Overall Depth 1666 Varilite Evo Overall Depth 1666 Varilite Evo Overall Depth 1666 Varilite Evo Overall Depth 1668 Varilite Evo	430 mm blution & E 355 mm blution & E 355 mm blution & E 405 mm	Overall Width	430 mm shion Cover 355 mm shion Cover 405 mm shion Cover 405 mm shion Cover 405 mm shion Cover 405 mm	Package Quantity Mesh 14" x 14" Package Quantity Incontinence 14" Package Quantity Mesh 14" x 16" Package Quantity Incontinence 14" Package Quantity Mesh 16" x 16" Package Quantity Incontinence 16" Package Quantity Mesh 16" x 18" Package Quantity Incontinence 16"	1 x 14" 1 x 16" 1 x 16" 1	Package Type Package Type	Single Single Single Single Single Single Single Single Single
Overall Depth 643 Varilite Evo Overall Depth 644 Varilite Evo Overall Depth 645 Varilite Evo Overall Depth 646 Varilite Evo Overall Depth 6665 Varilite Evo Overall Depth 6666 Varilite Evo Overall Depth 6667 Varilite Evo Overall Depth	430 mm blution & E 355 mm 355 mm 355 mm 355 mm blution & E 405 mm blution & E 405 mm blution & E 405 mm blution & E 405 mm	Overall Width Evolution PSV Current Evolution PSV Current Evolution PSV Current Evolution PSV Current Evolution PSV Current	430 mm shion Cover 355 mm shion Cover 405 mm shion Cover 405 mm shion Cover 405 mm shion Cover 405 mm	Package Quantity Mesh 14" x 14" Package Quantity Incontinence 14" Package Quantity Mesh 14" x 16" Package Quantity Mesh 14" x 16" Package Quantity Incontinence 14" Package Quantity Incontinence 14" Package Quantity Incontinence 16" Package Quantity Incontinence 16" Package Quantity Mesh 16" x 18" Package Quantity	1 x 14" 1 x 16" 1 x 16" 1 x 16" 1 x 16" 1 x 18"	Package Type Package Type Package Type Package Type Package Type Package Type	Single
Overall Depth 643 Varilite Evo Overall Depth 644 Varilite Evo Overall Depth 645 Varilite Evo Overall Depth 646 Varilite Evo Overall Depth 646 Varilite Evo Overall Depth 665 Varilite Evo Overall Depth 666 Varilite Evo Overall Depth 668 Varilite Evo Overall Depth	430 mm blution & E 355 mm blution & E 355 mm blution & E 405 mm blution & E 455 mm blution & E 455 mm blution & E 455 mm	Overall Width Evolution PSV Current Evolution PSV Current Evolution PSV Current Evolution PSV Current Evolution PSV Current	430 mm shion Cover 355 mm shion Cover 405 mm	Package Quantity Mesh 14" x 14" Package Quantity Incontinence 14" Package Quantity Mesh 14" x 16" Package Quantity Incontinence 14" Package Quantity Mesh 16" x 16" Package Quantity Incontinence 16" Package Quantity Mesh 16" x 18" Package Quantity Incontinence 16" Package Quantity Incontinence 16" Package Quantity	1 x 14" 1 x 16" 1 x 16" 1 x 16" 1 x 16" 1 x 18"	Package Type Package Type	Single

verall Depth	505 mm	Overall Width	455 mm	Package Quantity	1	Package Type	Single
681 Varilite Evol	ution & E	volution PSV Cushic	on Cover	Mesh 14" x 18"			
Overall Depth	455 mm	Overall Width	355 mm	Package Quantity	1	Package Type	Single
682 Varilite Evol	lution & E	volution PSV Cushi	on Cover	Mesh 14" x 20"			
Overall Depth	505 mm	Overall Width	355 mm	Package Quantity	1	Package Type	Single
1683 Varilite Evol	lution & E	volution PSV Cushi	on Cover	Mesh 15" x 15"			
Overall Depth	380 mm	Overall Width	380 mm	Package Quantity	1	Package Type	Single
1684 Varilite Evo	lution & E	volution PSV Cushi	on Cover	Mesh 15" x 17"			
Overall Depth	430 mm	Overall Width	380 mm	Package Quantity	1	Package Type	Single
1685 Varilite Evol	lution & E	volution PSV Cushi	on Cover	Mesh 18" x 16"			
Overall Depth	405 mm	Overall Width	455 mm	Package Quantity	1	Package Type	Single
1686 Varilite Evo	lution & E	volution PSV Cushi	on Cover	Incontinence 18"	x 16″		
Overall Depth	405 mm	Overall Width	455 mm	Package Quantity	1	Package Type	Single
1687 Varilite Evol	ution & E	volution PSV Cushi	on Cover	Mesh 18" x 18"			
Overall Depth	455						
Overall Depth	455 mm	Overall Width	455 mm	Package Quantity	1	Package Type	Single
		Overall Width			1 x 18″	Package Type	Single
688 Varilite Evo					1 x 18″	Package Type Package Type	Single Single
1688 Varilite Evol	lution & E	volution PSV Cushi	on Cover	Incontinence 18" Package Quantity	1 x 18"		
1688 Varilite Evo Overall Depth 1689 Varilite Evo	lution & E	Volution PSV Cushi	on Cover	Incontinence 18" Package Quantity	1 x 18" 1		
1688 Varilite Evol Overall Depth 1689 Varilite Evol Overall Depth	lution & E 455 mm lution & E 355 mm	Volution PSV Cushin	on Cover 455 mm on Cover 405 mm	Incontinence 18" Package Quantity Mesh 16" x 14" Package Quantity	1	Package Type	Single
1688 Varilite Evol overall Depth 1689 Varilite Evol overall Depth 1690 Varilite Evol	lution & E 455 mm lution & E 355 mm	Overall Width	on Cover 455 mm on Cover 405 mm	Incontinence 18" Package Quantity Mesh 16" x 14" Package Quantity	1	Package Type	Single
1688 Varilite Evol Overall Depth 1689 Varilite Evol Overall Depth 1690 Varilite Evol Overall Depth	lution & E 455 mm lution & E 355 mm lution & E 455 mm	Volution PSV Cushi Overall Width Volution PSV Cushi Overall Width	on Cover 455 mm 0n Cover 405 mm 0n Cover 355 mm	Incontinence 18" Package Quantity Mesh 16" x 14" Package Quantity Incontinence 14" Package Quantity	1 x 18″ 1	Package Type Package Type	Single
1688 Varilite Evol Overall Depth 1689 Varilite Evol Overall Depth 1690 Varilite Evol Overall Depth	lution & E 455 mm lution & E 355 mm lution & E 455 mm	Volution PSV Cushi Overall Width Volution PSV Cushi Overall Width	on Cover 455 mm 0n Cover 405 mm 0n Cover 355 mm	Incontinence 18" Package Quantity Mesh 16" x 14" Package Quantity Incontinence 14" Package Quantity	1 x 18″ 1	Package Type Package Type	Single
1688 Varilite Evol Overall Depth 1689 Varilite Evol Overall Depth 1690 Varilite Evol Overall Depth 1691 Varilite Evol Overall Depth	Iution & E 455 mm Jution & E 355 mm Iution & E 455 mm Jution & E 505 mm	Volution PSV Cushi Overall Width Volution PSV Cushi Overall Width Volution PSV Cushi Overall Width	on Cover 455 mm 455 mm 405 mm 405 mm 355 mm 355 mm 355 mm	Incontinence 18" Package Quantity Mesh 16" x 14" Package Quantity Incontinence 14" Package Quantity Incontinence 14" > Package Quantity	1 x 18" 1 x 20"	Package Type Package Type Package Type Package Type	Single Single Single
1688 Varilite Evol Overall Depth 1689 Varilite Evol Overall Depth 1690 Varilite Evol Overall Depth 1691 Varilite Evol Overall Depth	Iution & E 455 mm Jution & E 355 mm Iution & E 455 mm Jution & E 505 mm	Evolution PSV Cushi Overall Width Overall Width Overall Width Evolution PSV Cushi Overall Width Overall Width Volution PSV Cushi Overall Width	on Cover 455 mm 455 mm 405 mm 405 mm 355 mm 355 mm 355 mm	Incontinence 18" Package Quantity Mesh 16" x 14" Package Quantity Incontinence 14" Package Quantity Incontinence 14" > Package Quantity	1 x 18" 1 x 20"	Package Type Package Type Package Type Package Type	Single Single Single
1688 Varilite Evol Overall Depth 1689 Varilite Evol Overall Depth 1690 Varilite Evol Overall Depth 1691 Varilite Evol Overall Depth 1692 Varilite Evol Overall Depth	Iution & E 455 mm Iution & E 355 mm Iution & E 455 mm Iution & E 505 mm Iution & E 380 mm	Volution PSV Cushi Overall Width Volution PSV Cushi Overall Width Volution PSV Cushi Overall Width Volution PSV Cushi Overall Width	on Cover 455 mm 455 mm 405 mm 0 Cover 355 mm 0 Cover 355 mm 0 Cover 380 mm	Incontinence 18" Package Quantity Mesh 16" x 14" Package Quantity Incontinence 14" Package Quantity Incontinence 14" Package Quantity Incontinence 15" Package Quantity Incontinence 15"	1 x 18" 1 x 20" 1 x 15" 1	Package Type Package Type Package Type Package Type Package Type	Single Single Single Single
1688 Varilite Evol Overall Depth 1689 Varilite Evol Overall Depth 1690 Varilite Evol Overall Depth 1691 Varilite Evol Overall Depth 1692 Varilite Evol Overall Depth	Iution & E 455 mm Iution & E 355 mm Iution & E 455 mm Iution & E 505 mm Iution & E 380 mm	Volution PSV Cushi Overall Width Volution PSV Cushi	on Cover 455 mm 455 mm 405 mm 0 Cover 355 mm 0 Cover 355 mm 0 Cover 380 mm	Incontinence 18" Package Quantity Mesh 16" x 14" Package Quantity Incontinence 14" Package Quantity Incontinence 14" Package Quantity Incontinence 15" Package Quantity Incontinence 15"	1 x 18" 1 x 20" 1 x 15" 1	Package Type Package Type Package Type Package Type Package Type	Single Single Single Single
1688 Varilite Evol Overall Depth 1689 Varilite Evol Overall Depth 1690 Varilite Evol Overall Depth 1691 Varilite Evol Overall Depth 1692 Varilite Evol Overall Depth 1695 Varilite Evol Overall Depth	Iution & E 455 mm Iution & E 355 mm Iution & E 455 mm Iution & E 505 mm Iution & E 380 mm Iution & E 430 mm	Volution PSV Cushi Overall Width Overall Width Overall Width Volution PSV Cushi Overall Width Volution PSV Cushi Overall Width Volution PSV Cushi Overall Width	on Cover 455 mm on Cover 405 mm on Cover 355 mm on Cover 355 mm on Cover 355 mm on Cover 380 mm on Cover 380 mm	Incontinence 18" Package Quantity Mesh 16" x 14" Package Quantity Incontinence 14" Package Quantity Incontinence 14" Package Quantity Incontinence 15" Package Quantity Incontinence 15" Package Quantity	1 x 18" 1 x 20" 1 x 15" 1 x 17" 1	Package Type Package Type Package Type Package Type Package Type Package Type	Single Single Single Single Single





Varilite Zoid Cushion Cover



- The Zoid Cushion Cover is designed as a replacement cover for the Zoid Cushion
- The Zoid Cushion is performance cushion of choice for active wheelchair users at moderate risk of tissue breakdown
- Its tapered shape makes it ideal for rigid chairs and its cover is designed to protect the thighs of users with tight hamstrings
- or who tuck their legs under their chair • A large range of sizes are available 12" - 20" x 8" - 20"
- With 2 cover options mesh or incontinence

Overall Depth	455 mm	Overall Width	505 mm	Package Quantity	1	Package Type	Single
F00							
		Cover Mesh 20					
Overall Depth	505 mm	Overall Width	505 mm	Package Quantity	1	Package Type	Single
4571 Varilite Zo	oid Cushion	Cover Mesh 18"	x 16″				
Overall Depth	405 mm	Overall Width	455 mm	Package Quantity	1	Package Type	Single
		Cover Mesh 18"					
Overall Depth	505 mm	Overall Width	455 mm	Package Quantity	1	Package Type	Single
4573 Varilite Z	oid Cushion	Cover Mesh 20	" x 16"				
Overall Depth	405 mm	Overall Width	505 mm	Package Quantity	1	Package Type	Single
		Cover Mesh 12"					
Overall Depth	455 mm	Overall Width	305 mm	Package Quantity	1	Package Type	Single
4652 Varilite Z	oid Cushion	Cover Mesh 12'	′ x 12″				
Overall Depth	305 mm	Overall Width	305 mm	Package Quantity	1	Package Type	Single
4653 Varilite Z	oid Cushion	Cover Mesh 12'	' x 14"				
Overall Depth	355 mm	Overall Width	305 mm	Package Quantity	1	Package Type	Single
1651 Varilita 7	oid Cuchion	Cover Mesh 14'	' v 1/I"				
Overall Depth	355 mm	Overall Width	355 mm	Package Quantity	1	Package Type	Single
overall bepair	000 11111			r dokuge quantity		r dokuge type	olligic
4655 Varilite Z	oid Cushion	Cover Mesh 14'	′ x 16″				
Overall Depth	405 mm	Overall Width	355 mm	Package Quantity	1	Package Type	Single
4050							
		Cover Mesh 16'					
Overall Depth	405 mm	Overall Width	405 mm	Package Quantity	1	Package Type	Single
4657 Varilite Z	oid Cushion	Cover Mesh 18"	x 18″				
Overall Depth	455 mm	Overall Width	455 mm	Package Quantity	1	Package Type	Single
4659 Varilite Z	oid Cushion	Cover Mesh 16'	′ x 18″				
Overall Depth	455 mm	Overall Width	405 mm	Package Quantity	1	Package Type	Single
1697 Varilita 7	oid Cuphion	Covor Moch 14	v 10"				
4697 Varilite Z	oid Cushion	Cover Mesh 14"	x 18"	Package Quantity	1	Package Type	Single

4698 Varilite Zoi	d Cushion	Cover Mesh 15" x 1	5″
Overall Depth	380 mm	Overall Width	380 mm
4699 Varilite Zoi	d Cushion	Cover Mesh 15" x 17	7"
Overall Depth	430 mm	Overall Width	380 mm
4702 Varilite Zoi	d Cushion	Cover Incontinence	e 12" x 8"
Overall Depth	205 mm	Overall Width	305 mm
4711 Varilite Zoid	Cushion (Cover Incontinence	12" x 12"
Overall Depth	305 mm	Overall Width	305 mm
4712 Varilite Zoid	d Cushion	Cover Incontinence	12" x 14"
Overall Depth	355 mm	Overall Width	305 mm
4715 Varilite Zoid	l Cushion	Cover Incontinence	14" x 14"
Overall Depth	355 mm	Overall Width	355 mm
4719 Varilite Zoid	d Cushion	Cover Incontinence	14" x 16"
Overall Depth	405 mm	Overall Width	355 mm
4720 Varilite Zoi	d Cushion	Cover Incontinence	e 14″ x 18′
Overall Depth	455 mm	Overall Width	355 mm
4721 Varilite Zoid	l Cushion	Cover Incontinence	15" x 15"
Overall Depth	380 mm	Overall Width	380 mm
4722 Varilite Zoi	d Cushion	Cover Incontinence	e 15″ x 17″
Overall Depth	430 mm	Overall Width	380 mm
4725 Varilite Zoi	d Cushion	Cover Incontinence	e 16" x 16'
Overall Depth	405 mm	Overall Width	405 mm
4726 Varilite Zoi	d Cushion	Cover Incontinence	e 16″ x 18′
Overall Depth	455 mm	Overall Width	
		overall math	405 mm
4727 Varilite Zoio	1 Cushion	Cover Incontinence	
4727 Varilite Zoid	d Cushion		
Overall Depth	505 mm	Cover Incontinence	2 16" x 20
Overall Depth	505 mm	Cover Incontinence Overall Width	2 16" x 20
Overall Depth 4728 Varilite Zoid Overall Depth	505 mm d Cushion 430 mm	Cover Incontinence Overall Width Cover Incontinence	2 16" x 20 405 mm 2 17" x 17" 430 mm
Overall Depth 4728 Varilite Zoid Overall Depth	505 mm d Cushion 430 mm	Cover Incontinence Overall Width Cover Incontinence Overall Width	2 16" x 20 405 mm 2 17" x 17" 430 mm
Overall Depth 4728 Varilite Zoid Overall Depth 4729 Varilite Zoid Overall Depth	505 mm	Cover Incontinence Overall Width Cover Incontinence Overall Width Cover Incontinence	e 16" x 20 405 mm e 17" x 17" 430 mm e 18" x 16' 455 mm
Overall Depth 4728 Varilite Zoid Overall Depth 4729 Varilite Zoid Overall Depth	505 mm	Cover Incontinence Overall Width Cover Incontinence Overall Width Cover Incontinence	e 16" x 20 405 mm e 17" x 17" 430 mm e 18" x 16' 455 mm
Overall Depth 4728 Varilite Zoid Overall Depth 4729 Varilite Zoid Overall Depth 4730 Varilite Zoid Overall Depth	505 mm d Cushion 430 mm d Cushion 405 mm d Cushion 455 mm	Cover Incontinence Overall Width Cover Incontinence Overall Width Cover Incontinence Overall Width Cover Incontinence	e 16" x 20 405 mm e 17" x 17" 430 mm e 18" x 16' 455 mm e 18" x 18' 455 mm
Overall Depth 4728 Varilite Zoid Overall Depth 4729 Varilite Zoid Overall Depth 4730 Varilite Zoid Overall Depth	505 mm d Cushion 430 mm d Cushion 405 mm d Cushion 455 mm	Cover Incontinence Overall Width Cover Incontinence Overall Width Cover Incontinence Overall Width Cover Incontinence Overall Width Cover Incontinence	e 16" x 20 405 mm e 17" x 17" 430 mm e 18" x 16' 455 mm e 18" x 18' 455 mm
Overall Depth 4728 Varilite Zoid Overall Depth 4729 Varilite Zoid Overall Depth 4730 Varilite Zoid Overall Depth 4731 Varilite Zoid Overall Depth	505 mm d Cushion 430 mm d Cushion 405 mm d Cushion 455 mm d Cushion 505 mm	Cover Incontinence Overall Width Cover Incontinence Overall Width Cover Incontinence Overall Width Cover Incontinence Overall Width Cover Incontinence	e 16" x 20 405 mm e 17" x 17" 430 mm e 18" x 16' 455 mm e 18" x 18' 455 mm 18" x 20' 455 mm
Overall Depth 4728 Varilite Zoid Overall Depth 4729 Varilite Zoid Overall Depth 4730 Varilite Zoid Overall Depth 4731 Varilite Zoid Overall Depth	505 mm d Cushion 430 mm d Cushion 405 mm d Cushion 455 mm d Cushion 505 mm	Cover Incontinence Overall Width Cover Incontinence	e 16" x 20 405 mm e 17" x 17" 430 mm e 18" x 16 455 mm e 18" x 18 455 mm 18" x 20 455 mm



Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
	_		
Package Quantity	1	Package Type	Single
	_		
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
	_		
Package Quantity	1	Package Type	Single
	_		
Package Quantity	1	Package Type	Single
	_		
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single
Package Quantity	1	Package Type	Single

4781 Varilite Zoid Cushion Cover Incontinence 20" x 18"									
Overall Depth 455 mm Overall Width 505 mm Package Quantity 1 Package Type Single									
4785 Varilite Zoid Cushion Cover Incontinence 20" x 20"									
Overall Depth	505 mm	Overall Width	505 mm	Package Quantity	1	Package Type	Single		
4847 Varilite Zoid (Cushion	Cover Mesh 16" x 20	n						
Overall Depth	505 mm	Overall Width	405 mm	Package Quantity	1	Package Type	Single		
4888 Varilite Zoid Cushion Cover Mesh 17" x 17"									
4888 Varilite Zoid (Cushion	Cover Mesh 17" x 17"							

Roho Cover Select Series



- Roho Cushion Cover Select Series Easy to remove and replace
- Helps to protect the cushion inside it
- Made from a durable fabricMakes sitting more comfortable
- Cover only

R-COV-SS99 Roho Cover Select Series 9" x 9" HP									
Overall Depth 230 mm Overall Width 230 mm Package Quantity 1						Package Type	Single		
R-COV-SS99LP Rob	o Cover (Select Series 9" x 9"	IP						
Overall Depth	230 mm	Overall Width	230 mm	Package Quantity		Package Type			

Pressure Cushions Non Powered

7.06 Manual Pumps for Cushion Inflation

Roho Air Pump



- The Roho Air Pump is an easy way to inflate your Roho mattress or cushion
 The dual action pump
- Pumps air into the cushion with 1 stroke of the pump (Air is pumped in on both the up and down strokes)

AIR PUMP ROHO Roho Air Pump



National Distribution, Local Service

217

COUNTRY CARE GROUP

NATIONAL MEMBER DIRECTORY

MEMBER NAME	ADDRESS	TOWN	STATE	P/C	PHONE	FAX
ACT						
Mobility Matters	35 Townsville ST	FYSHWICK	ACT	2609	02 6280 7244	02 6280 0602
Total Mobility	Unit 3/21 Raws Cres	HUME	ACT	2620	02 6260 1003	02 6260 1156

NSW Astley Mobility Cnr Mamre Rd & Great Western Hwy ST MARYS 02 9673 2212 02 9673 0955 NSW 2760 Astley Mobility Cnr Pennant Hills Rd & Fisher Ave PENNANT HILLS 03 9484 7070 02 9484 7070 NSW 2120 Country Care Group Northern Beaches 577 Pittwater Rd BROOKVALE NSW 2100 02 9907 2696 02 9907 2191 Country Mobility Trends 28 Holland Rd COOMA NSW 2630 02 6452 5525 Eden Healthcare Solutions 575 Olive St 2640 02 6041 9700 AI BURY NSW 02 6041 9711 Hospital At Home 1/30 Heathcote Rd MOOREBAN NSW 2170 02 9601 7757 02 6901 7879 Hospital At Home 69 Munibung Rd CARDIFF 2285 02 4954 4173 02 4954 4761 NSW 2650 Hospital At Home 352 Edward St WAGGA WAGGA NSW 02 6925 4966 02 6925 4977 2500 02 4227 4315 Independent Home Care Supplies 2/21 Denison St WOLLONGONG NSW 02 4227 4316 Independent Living Specialists 67 Mars Rd LANE COVE NSW 2066 02 9427 4995 02 9427 4338 RANDWICK 2031 02 9398 9263 02 9398 7765 Independent Living Specialists 204 Alison Rd NSW Independent Living Specialists 295-297 Princes Hwy ARNCLIFFE NSW 2205 02 9597 6847 02 9633 1577 Independent Living Specialists Shop 6, 15 Bransgrove St WENTWORTHVILLE NSW 2145 02 9633 1155 02 9633 1577 Independent Living Specialists Unit 4, 7-13 Victoria Ave CASTLE HILL NSW 2154 02 9659 9750 02 9659 8526 Mobility 2 You 2/12 Aston Rd ERINA NSW 2250 02 4367 5751 02 4367 5758 Mobility 2 You 6/2 O'Hart Cl CHARMHAVEN NSW 2263 02 4392 0148 Mobility Matters BFGA 2250 117 Carp St NSW 02 6492 3888 02 6492 3889 Mobility Matters BATEMANS BAY 02 4472 1544 34 Orient St NSW 2536 02 4472 1044 New England Mobility 155 Miller St ARMIDAI F NSW 2350 02 6772 6007 02 6772 6003 Northern Rivers Surgical BALLINA NSW 2478 02 6686 6644 02 6686 9383 18 Endeavour Cl Nth Degree Care TURRAMURRA NSW 2074 02 9983 9100 02 9983 1806 1323 Pacific Hwy Patient Handling Australia 1 Mitchell St MARRICKVILLE NSW 2204 1300 734 501 02 9557 0600 Riverina Scooters and Mobility 8 Blake St WAGGA WAGGA NSW 2650 02 6921 4444 02 6921 1842 South Coast Independent Living Supplies NOWRA 2541 Shop 3/92 North St NSW 02 4421 3383 02 4423 3397 Specialised Wheelchair Company BROOKVALE NSW 2100 02 9905 5333 11 Mitchell Rd 02 9905 2208 Total Mobility Unit 10, 1-11 Burns Rd HEATHCOTE NSW 1300 868 662 02 9520 4345 Vital Living Solutions 2/187 Lake Rd PORT MACQUARIE NSW 2444 02 6581 2400 02 6581 2422 Vital Living Solutions 2/6 Enterprise Crt FORSTER NSW 2428 02 6554 9444 02 6554 9777 2450 02 6651 7900 02 6651 1335 Willaid Healthcare Equipment 1/16 Isles Dr COFFS HARBOUR NSW

NI						
Keep Moving	15 Swan Crescent	WINNELLIE	NT	0820	08 8947 5122	08 8947 2531
Keep Moving	8 Gap Rd	ALICE SPRINGS	NT	0870	08 8952 3599	08 8952 4833
Equip4Living	Lot 4643 Export Dr	BERRIMAH	NT	0820	0405 041 386	08 8362 6499

MEMBER NAME	ADDRESS	TOWN	STATE	P/C	PHONE	FAX
QLD						
Austech Medical Pty Ltd	Unit 1, 88 Lear Jet Dr	CABOOLTURE	QLD	4510	07 5495 8663	07 5495 8273
Central Queensland Mobility	1A 57 Alexandra St	ROCKHAMPTON	QLD	4701	1800 201 070	07 4926 1073
Complete Mobility and Rehab	Shop 1, 92 Anderson St	MANUNDA	QLD	4870	07 4032 1140	07 4053 1104
Independent Living Solutions - Think Mobility	8 Keane St	CURRAJONG	QLD	4812	07 4728 1200	07 4728 1201
Lakeside Mobility	2/1 Metier Linkway	BIRTINYA	QLD	4575	07 5325 1231	07 5309 5350
Medtech Mobility Equipment	18 Prescott St	TOOWOOMBA	QLD	4350	07 4639 3977	07 4659 5953
Queensland Rehab Equipment	10 Devlan St	MANSFIELD	QLD	4122	1300 743 710	1300 802 710
Queensland Rehab Equipment		GOLD COAST	QLD		1300 743 710	1300 802 710
Regional Rehab & Mobility	4A 120 Hanson Rd	GLADSTONE	QLD	4680	07 4972 4449	07 4972 4443
Regional Rehab & Mobility		BUNDABERG	QLD	4670	0418 716 953	07 4926 1073
Think Mobility	2037 Sandgate Rd	VIRGINIA	QLD	4014	1300 881 968	07 3318 7601
Think Mobility Mackay	3/53 Central Park Dr	PAGET	QLD	4740	1300 839 036	1300 839 036
SA Equip4Living	57 Magill Rd	STEPNEY	SA	5069	08 8362 6422	08 8362 6499
	17 Lindsay Ave	EDWARDSTOWN	SA	5039	08 8276 1300	08 8362 6499
Xtracare Contracts	622 – 624 Regency Rd	BROADVIEW	SA	5083	08 8266 7000	08 8266 7111
VicLean Healthcare	24-28 Earl St	LAUNCESTON	TAS	7250	03 6331 7400	03 6331 7140
					0000017100	
VIC	156 York St	SALE	VIC	3850	1300 723 830	1300 723 831
VIC Access Rehabilitation Equipment	156 York St 115A Nicholson St	SALE BAIRNSDALE	VIC			1300 723 831
VIC Access Rehabilitation Equipment	115A Nicholson St	-	-	3850	1300 723 830	1300 723 831
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales	115A Nicholson St 4/3 Torca Tce	BAIRNSDALE	VIC	3850 3875	1300 723 830 03 5152 5188	
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave	BAIRNSDALE MORNINGTON BERWICK	VIC VIC VIC	3850 3875 3931 3806	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955	03 5975 8304 03 9707 4933
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales Country Care Group Bendigo	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave 54 Bridge St	BAIRNSDALE MORNINGTON BERWICK BENDIGO	VIC	3850 3875 3931 3806 3550	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955 03 5442 1131	03 5975 8304 03 9707 4933 03 5442 1135
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales Country Care Group Bendigo Country Care Group Echuca	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave	BAIRNSDALE MORNINGTON BERWICK	VIC VIC VIC VIC	3850 3875 3931 3806	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955	03 5975 8304 03 9707 4933
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales Country Care Group Bendigo Country Care Group Echuca Country Care Group Mildura	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave 54 Bridge St Shop 2/110 Northern Hwy 165 Ninth St	BAIRNSDALE MORNINGTON BERWICK BENDIGO ECHUCA MILDURA	VIC VIC VIC VIC VIC VIC VIC	3850 3875 3931 3806 3550 3564 3500	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955 03 5442 1131 03 5482 1422 03 5022 1680	03 5975 8304 03 9707 4933 03 5442 1135 03 5480 3524 03 5021 5633
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales Country Care Group Bendigo Country Care Group Bendigo Country Care Group Echuca Country Care Group Mildura Country Care Group Sunbury	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave 54 Bridge St Shop 2/110 Northern Hwy 165 Ninth St 18/106-126 Gap Rd	BAIRNSDALE MORNINGTON BERWICK BENDIGO ECHUCA MILDURA SUNBURY	VIC VIC VIC VIC VIC	3850 3875 3931 3806 3550 3564 3500 3500 3429	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955 03 5442 1131 03 5482 1422 03 5022 1680 03 9744 7851	03 5975 8304 03 9707 4933 03 5442 1135 03 5480 3524
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales Country Care Group Bendigo Country Care Group Bendigo Country Care Group Echuca Country Care Group Echuca Country Care Group Sunbury Country Care Horsham	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave 54 Bridge St Shop 2/110 Northern Hwy 165 Ninth St 18/106-126 Gap Rd 6 King Drive	BAIRNSDALE MORNINGTON BERWICK BENDIGO ECHUCA MILDURA SUNBURY HORSHAM	VIC VIC VIC VIC VIC VIC VIC VIC	3850 3875 3931 3806 3550 3564 3500 3429 3400	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955 03 5442 1131 03 5482 1422 03 5022 1680 03 9744 7851 03 5381 1404	03 5975 8304 03 9707 4933 03 5442 1135 03 5480 3524 03 5021 5633 03 9744 7183 03 5382 0766
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales Country Care Group Bendigo Country Care Group Sunbury Country Care Horsham Country Healthcare Shepparton	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave 54 Bridge St Shop 2/110 Northern Hwy 165 Ninth St 18/106-126 Gap Rd 6 King Drive 234-236 Wynham St	BAIRNSDALE MORNINGTON BERWICK BENDIGO ECHUCA MILDURA SUNBURY HORSHAM SHEPPARTON	VIC VIC VIC VIC VIC VIC VIC VIC VIC	3850 3875 3931 3806 3550 3564 3500 3429 3400 3420 3400	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955 03 5442 1131 03 5482 1422 03 5022 1680 03 9744 7851 03 5381 1404 03 5821 9371	03 5975 8304 03 9707 4933 03 5442 1135 03 5480 3524 03 5021 5633 03 9744 7183 03 5382 0766 03 5821 9241
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales Country Care Group Bendigo Country Care Group Bendigo Country Care Group Echuca Country Care Group Echuca Country Care Group Sunbury Country Care Group Sunbury Country Care Horsham Country Healthcare Shepparton Elite Mobility	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave 54 Bridge St Shop 2/110 Northern Hwy 165 Ninth St 18/106-126 Gap Rd 6 King Drive 234-236 Wynham St 2 Dawson St	BAIRNSDALE MORNINGTON BERWICK BENDIGO ECHUCA MILDURA SUNBURY HORSHAM SHEPPARTON BALLARAT	VIC VIC VIC VIC VIC VIC VIC VIC	3850 3875 3931 3806 3550 3564 3500 3429 3400 3630 3630 3350	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955 03 5442 1131 03 5482 1422 03 5022 1680 03 9744 7851 03 5381 1404 03 5821 9371 03 5335 9011	03 5975 8304 03 9707 4933 03 5442 1135 03 5480 3524 03 5021 5633 03 9744 7183 03 5382 0766 03 5821 9241 03 5331 9655
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales Country Care Group Bendigo Country Care Group Bendigo Country Care Group Echuca Country Care Group Echuca Country Care Group Mildura Country Care Group Sunbury Country Care Horsham Country Healthcare Shepparton Elite Mobility Fisher Lane Mobility	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave 54 Bridge St Shop 2/110 Northern Hwy 165 Ninth St 18/106-126 Gap Rd 6 King Drive 234-236 Wynham St	BAIRNSDALE MORNINGTON BERWICK BENDIGO ECHUCA MILDURA SUNBURY HORSHAM SHEPPARTON	VIC VIC VIC VIC VIC VIC VIC VIC VIC VIC	3850 3875 3931 3806 3550 3564 3500 3429 3400 3420 3400	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955 03 5442 1131 03 5482 1422 03 5022 1680 03 9744 7851 03 5381 1404 03 5821 9371	03 5975 8304 03 9707 4933 03 5442 1135 03 5480 3524 03 5021 5633 03 9744 7183 03 5382 0766 03 5821 9241
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales Country Care Group Bendigo Country Care Group Bendigo Country Care Group Echuca Country Care Group Echuca Country Care Group Sunbury Country Care Group Sunbury Country Care Horsham Country Healthcare Shepparton Elite Mobility Fisher Lane Mobility Fisher Lane Mobility	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave 54 Bridge St Shop 2/110 Northern Hwy 165 Ninth St 18/106-126 Gap Rd 6 King Drive 234-236 Wynham St 2 Dawson St 26 Victoria Cres 13/178 Duke St	BAIRNSDALE MORNINGTON BERWICK BENDIGO ECHUCA MILDURA SUNBURY HORSHAM SHEPPARTON BALLARAT ABBOTSFORD BRAYBROOK	21V 21V 21V 21V 21V 21V 21V 21V 21V 21V	3850 3875 3931 3806 3550 3564 3500 3429 3400 3429 3400 3400 3420 3400 3400 3400 3400 3400	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955 03 5442 1131 03 5482 1422 03 5022 1680 03 9744 7851 03 5381 1404 03 5381 1404 03 5821 9371 03 5335 9011 03 9419 2250 03 9312 1467	03 5975 8304 03 9707 4933 03 5442 1135 03 5480 3524 03 5021 5633 03 9744 7183 03 5382 0766 03 5821 9241 03 5331 9655 03 9417 2074 03 9417 2074
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales Country Care Group Bendigo Country Care Group Bendigo Country Care Group Echuca Country Care Group Echuca Country Care Group Sunbury Country Care Group Sunbury Country Care Horsham Country Healthcare Shepparton Elite Mobility Fisher Lane Mobility Fisher Lane Mobility G Mobility	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave 54 Bridge St Shop 2/110 Northern Hwy 165 Ninth St 18/106-126 Gap Rd 6 King Drive 234-236 Wynham St 2 Dawson St 26 Victoria Cres 13/178 Duke St Shop 6 39 Princes Way	BAIRNSDALE MORNINGTON BERWICK BENDIGO ECHUCA MILDURA SUNBURY HORSHAM SHEPPARTON BALLARAT ABBOTSFORD BRAYBROOK WARRAGUL	VIC VIC VIC VIC VIC VIC VIC VIC VIC VIC	3850 3875 3931 3806 3550 3564 3500 3429 3400 3429 3400 3630 3420 3630 3067 3019 3820	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955 03 5442 1131 03 5482 1422 03 5022 1680 03 9744 7851 03 5381 1404 03 5821 9371 03 5335 9011 03 9419 2250 03 9312 1467 03 5623 6966	03 5975 8304 03 9707 4933 03 5442 1135 03 5480 3524 03 5021 5633 03 9744 7183 03 5382 0766 03 5821 9241 03 5331 9655 03 9417 2074 03 9417 2074 03 5623 6966
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales Country Care Group Bendigo Country Care Group Bendigo Country Care Group Echuca Country Care Group Echuca Country Care Group Sunbury Country Care Group Sunbury Country Care Horsham Country Healthcare Shepparton Elite Mobility Fisher Lane Mobility Fisher Lane Mobility Gelong Wheelchairs	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave 54 Bridge St Shop 2/110 Northern Hwy 165 Ninth St 18/106-126 Gap Rd 6 King Drive 234-236 Wynham St 2 Dawson St 26 Victoria Cres 13/178 Duke St Shop 6 39 Princes Way 106 Bailey St	BAIRNSDALE MORNINGTON BERWICK BENDIGO ECHUCA MILDURA SUNBURY HORSHAM SHEPPARTON BALLARAT ABBOTSFORD BRAYBROOK WARRAGUL GROVEDALE	VIC VIC VIC VIC VIC VIC VIC VIC VIC VIC	3850 3875 3931 3806 3550 3564 3500 3429 3400 3630 3630 3630 3630 3067 3019 3820 3216	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955 03 5442 1131 03 5442 1131 03 5022 1680 03 5022 1680 03 9744 7851 03 5381 1404 03 5335 9011 03 9419 2250 03 9312 1467 03 5623 6966 03 5244 0844	03 5975 8304 03 9707 4933 03 5442 1135 03 5480 3524 03 5021 5633 03 9744 7183 03 5382 0766 03 5821 9241 03 5331 9655 03 9417 2074 03 9417 2074 03 5623 6966 03 5244 0850
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales Country Care Group Bendigo Country Care Group Bendigo Country Care Group Echuca Country Care Group Echuca Country Care Group Mildura Country Care Group Sunbury Country Care Horsham Country Healthcare Shepparton Elite Mobility Fisher Lane Mobility Fisher Lane Mobility Gelong Wheelchairs Geelong Wheelchairs	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave 54 Bridge St Shop 2/110 Northern Hwy 165 Ninth St 18/106-126 Gap Rd 6 King Drive 234-236 Wynham St 2 Dawson St 2 Dawson St 26 Victoria Cres 13/178 Duke St Shop 6 39 Princes Way 106 Bailey St 215 Murray St	BAIRNSDALE MORNINGTON BERWICK BENDIGO ECHUCA MILDURA SUNBURY HORSHAM SHEPPARTON BALLARAT ABBOTSFORD BRAYBROOK WARRAGUL GROVEDALE COLAC	VIC VIC VIC VIC VIC VIC VIC VIC VIC VIC	3850 3875 3931 3806 3550 3564 3500 3429 3400 3630 3630 3630 3630 3067 3019 3820 3216 3216 3207	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955 03 5442 1131 03 5442 1131 03 5022 1680 03 5022 1680 03 535 1404 03 5331 1404 03 5335 9011 03 9419 2250 03 9312 1467 03 5623 6966 03 5244 0844 03 5231 6050	03 5975 8304 03 9707 4933 03 5442 1135 03 5480 3524 03 5021 5633 03 9744 7183 03 5382 0766 03 5821 9241 03 5331 9655 03 9417 2074 03 9417 2074 03 5623 6966 03 5244 0850 03 5244 0850
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales Country Care Group Bendigo Country Care Group Bendigo Country Care Group Echuca Country Care Group Echuca Country Care Group Sunbury Country Care Group Sunbury Country Care Horsham Country Healthcare Shepparton Elite Mobility Fisher Lane Mobility Fisher Lane Mobility Gelong Wheelchairs Geelong Wheelchairs Life Mobility	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave 54 Bridge St Shop 2/110 Northern Hwy 165 Ninth St 18/106-126 Gap Rd 6 King Drive 234-236 Wynham St 2 Dawson St 26 Victoria Cres 13/178 Duke St Shop 6 39 Princes Way 106 Bailey St 215 Murray St 24 Corporate Blvd	BAIRNSDALE MORNINGTON BERWICK BENDIGO ECHUCA MILDURA SUNBURY HORSHAM SHEPPARTON BALLARAT ABBOTSFORD BRAYBROOK WARRAGUL GROVEDALE COLAC BAYSWATER	VIC VIC VIC VIC VIC VIC VIC VIC VIC VIC	3850 3875 3931 3806 3550 3564 3500 3429 3400 3630 3630 3630 3630 3067 3019 3820 3019 3820 3216 3707 3153	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955 03 5442 1131 03 5442 1131 03 5022 1680 03 9744 7851 03 5381 1404 03 5335 9011 03 9312 1467 03 5623 6966 03 5244 0844 03 5231 6050 03 5231 6050 03 9729 5552	03 5975 8304 03 9707 4933 03 5442 1135 03 5480 3524 03 5021 5633 03 9744 7183 03 5382 0766 03 5821 9241 03 5331 9655 03 9417 2074 03 9417 2074 03 5623 6966 03 5244 0850 03 5244 0850 03 9729 5112
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales Country Care Group Bendigo Country Care Group Bendigo Country Care Group Echuca Country Care Group Echuca Country Care Group Sunbury Country Care Group Sunbury Country Care Horsham Country Healthcare Shepparton Elite Mobility Fisher Lane Mobility Fisher Lane Mobility Gelong Wheelchairs Geelong Wheelchairs Life Mobility Life Mobility	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave 54 Bridge St 54 Dridge St 54 Dridge St 54 Dridge St 18/106-126 Gap Rd 6 King Drive 234-236 Wynham St 2 Dawson St 2 Dawson St 26 Victoria Cres 13/178 Duke St 5hop 6 39 Princes Way 106 Bailey St 215 Murray St 24 Corporate Blvd 34 Railway Ave	BAIRNSDALE MORNINGTON BERWICK BENDIGO ECHUCA MILDURA SUNBURY HORSHAM SHEPPARTON BALLARAT ABBOTSFORD BRAYBROOK WARRAGUL GROVEDALE COLAC BAYSWATER RINGWOOD EAST	VIC VIC VIC VIC VIC VIC VIC VIC VIC VIC	3850 3875 3931 3806 3550 3564 3500 3429 3400 3630 3630 3630 3630 3077 3019 3820 3216 3707 3153 3153	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955 03 5442 1131 03 5442 1131 03 5022 1680 03 9744 7851 03 5381 1404 03 5382 9371 03 535 9011 03 9312 1467 03 5623 6966 03 5244 0844 03 5231 6050 03 9729 5552 03 9879 8885	03 5975 8304 03 9707 4933 03 5442 1135 03 5480 3524 03 5021 5633 03 9744 7183 03 5382 0766 03 5821 9241 03 5331 9655 03 9417 2074 03 5623 6966 03 5244 0850 03 5244 0850 03 9729 5112 03 9879 3133
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales Country Care Group Bendigo Country Care Group Bendigo Country Care Group Echuca Country Care Group Echuca Country Care Group Sunbury Country Care Group Sunbury Country Care Horsham Country Healthcare Shepparton Elite Mobility Fisher Lane Mobility Fisher Lane Mobility Geelong Wheelchairs Geelong Wheelchairs Life Mobility Life Mobility Life Mobility Mitchell Healthcare	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave 54 Bridge St Shop 2/110 Northern Hwy 165 Ninth St 18/106-126 Gap Rd 6 King Drive 234-236 Wynham St 2 Dawson St 2 Dawson St 26 Victoria Cres 13/178 Duke St Shop 6 39 Princes Way 106 Bailey St 215 Murray St 24 Corporate Blvd 34 Railway Ave 14a Sydney Rd	BAIRNSDALE MORNINGTON BERWICK BENDIGO ECHUCA MILDURA SUNBURY HORSHAM SHEPPARTON BALLARAT ABBOTSFORD BRAYBROOK WARRAGUL GROVEDALE COLAC BAYSWATER RINGWOOD EAST BENALLA	VIC VIC VIC VIC VIC VIC VIC VIC VIC VIC	3850 3875 3931 3806 3550 3564 3500 3429 3400 3630 3630 3350 3067 3019 3820 3216 3820 3216 3707 3153 3153 3135	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955 03 5442 1131 03 5442 1131 03 5022 1680 03 9744 7851 03 5381 1404 03 5382 9371 03 5335 9011 03 5335 9011 03 9312 1467 03 5623 6966 03 5244 0844 03 5231 6050 03 9729 5552 03 9879 8885 03 5762 1125	03 5975 8304 03 9707 4933 03 5442 1135 03 5442 1135 03 5442 1135 03 5480 3524 03 5021 5633 03 9744 7183 03 5382 0766 03 5821 9241 03 5331 9655 03 9417 2074 03 5623 6966 03 5623 6966 03 5244 0850 03 5244 0850 03 9729 5112 03 9879 3133 03 5762 1124
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales Country Care Group Bendigo Country Care Group Bendigo Country Care Group Echuca Country Care Group Echuca Country Care Group Sunbury Country Care Group Sunbury Country Care Horsham Country Healthcare Shepparton Elite Mobility Fisher Lane Mobility Fisher Lane Mobility Gaelong Wheelchairs Geelong Wheelchairs Life Mobility Life Mobility Life Mobility Mitchell Healthcare Mitchell Healthcare	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave 54 Bridge St Shop 2/110 Northern Hwy 165 Ninth St 18/106-126 Gap Rd 6 King Drive 234-236 Wynham St 2 Dawson St 26 Victoria Cres 13/178 Duke St Shop 6 39 Princes Way 106 Bailey St 215 Murray St 24 Corporate Blvd 34 Railway Ave 14a Sydney Rd 95 Murphy St	BAIRNSDALE MORNINGTON BERWICK BENDIGO ECHUCA MILDURA SUNBURY HORSHAM SHEPPARTON BALLARAT ABBOTSFORD BRAYBROOK WARRAGUL GROVEDALE COLAC BAYSWATER RINGWOOD EAST BENALLA WANGARATTA	VIC VIC VIC VIC VIC VIC VIC VIC VIC VIC	3850 3875 3931 3806 3550 3564 3500 3429 3400 3630 3630 3630 3067 3019 3820 3067 3019 3820 3067 3153 3153 3153 3135	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955 03 5442 1131 03 5442 1131 03 5022 1680 03 9744 7851 03 5381 1404 03 5382 19371 03 5335 9011 03 5335 9011 03 5321 9467 03 5623 6966 03 5234 0844 03 5231 6050 03 5231 6050 03 9729 5552 03 9879 8885 03 5721 5660	03 5975 8304 03 9707 4933 03 5442 1135 03 5442 1135 03 5442 1135 03 5442 1135 03 5021 5633 03 9744 7183 03 5382 0766 03 5821 9241 03 5331 9655 03 9417 2074 03 5623 6966 03 5244 0850 03 5244 0850 03 9729 5112 03 9879 3133 03 5762 1124 03 5721 9667
VIC Access Rehabilitation Equipment Bairnsdale Mobility & Hire Coastcare Medical Equipment Hire and Sales Coastcare Medical Equipment Hire and Sales Country Care Group Bendigo Country Care Group Bendigo Country Care Group Echuca Country Care Group Echuca Country Care Group Sunbury Country Care Group Sunbury Country Care Horsham Country Healthcare Shepparton Elite Mobility Fisher Lane Mobility Fisher Lane Mobility Fisher Lane Mobility Geelong Wheelchairs Geelong Wheelchairs Life Mobility Life Mobility Life Mobility Mitchell Healthcare Omni Healthcare	115A Nicholson St 4/3 Torca Tce 1/80-86 Enterprise Ave 54 Bridge St Shop 2/110 Northern Hwy 165 Ninth St 18/106-126 Gap Rd 6 King Drive 234-236 Wynham St 2 Dawson St 2 Dawson St 26 Victoria Cres 13/178 Duke St Shop 6 39 Princes Way 106 Bailey St 215 Murray St 24 Corporate Blvd 34 Railway Ave 14a Sydney Rd	BAIRNSDALE MORNINGTON BERWICK BENDIGO ECHUCA MILDURA SUNBURY HORSHAM SHEPPARTON BALLARAT ABBOTSFORD BRAYBROOK WARRAGUL GROVEDALE COLAC BAYSWATER RINGWOOD EAST BENALLA	VIC VIC VIC VIC VIC VIC VIC VIC VIC VIC	3850 3875 3931 3806 3550 3564 3500 3429 3400 3630 3630 3350 3067 3019 3820 3216 3820 3216 3707 3153 3153 3135	1300 723 830 03 5152 5188 03 5975 6107 03 9707 4955 03 5442 1131 03 5442 1131 03 5022 1680 03 9744 7851 03 5381 1404 03 5382 9371 03 5335 9011 03 5335 9011 03 9312 1467 03 5623 6966 03 5244 0844 03 5231 6050 03 9729 5552 03 9879 8885 03 5762 1125	03 5975 8304 03 9707 4933 03 5442 1135 03 5442 1135 03 5442 1135 03 5480 3524 03 5021 5633 03 9744 7183 03 5382 0766 03 5821 9241 03 5331 9655 03 9417 2074 03 5623 6966 03 5623 6966 03 5244 0850 03 5244 0850 03 9729 5112 03 9879 3133 03 5762 1124

VIC						
Access Rehabilitation Equipment	156 York St	SALE	VIC	3850	1300 723 830	1300 723 831
Bairnsdale Mobility & Hire	115A Nicholson St	BAIRNSDALE	VIC	3875	03 5152 5188	
Coastcare Medical Equipment Hire and Sales	4/3 Torca Tce	MORNINGTON	VIC	3931	03 5975 6107	03 5975 8304
Coastcare Medical Equipment Hire and Sales	1/80-86 Enterprise Ave	BERWICK	VIC	3806	03 9707 4955	03 9707 4933
Country Care Group Bendigo	54 Bridge St	BENDIGO	VIC	3550	03 5442 1131	03 5442 1135
Country Care Group Echuca	Shop 2/110 Northern Hwy	ECHUCA	VIC	3564	03 5482 1422	03 5480 3524
Country Care Group Mildura	165 Ninth St	MILDURA	VIC	3500	03 5022 1680	03 5021 5633
Country Care Group Sunbury	18/106-126 Gap Rd	SUNBURY	VIC	3429	03 9744 7851	03 9744 7183
Country Care Horsham	6 King Drive	HORSHAM	VIC	3400	03 5381 1404	03 5382 0766
Country Healthcare Shepparton	234-236 Wynham St	SHEPPARTON	VIC	3630	03 5821 9371	03 5821 9241
Elite Mobility	2 Dawson St	BALLARAT	VIC	3350	03 5335 9011	03 5331 9655
Fisher Lane Mobility	26 Victoria Cres	ABBOTSFORD	VIC	3067	03 9419 2250	03 9417 2074
Fisher Lane Mobility	13/178 Duke St	BRAYBROOK	VIC	3019	03 9312 1467	03 9417 2074
G Mobility	Shop 6 39 Princes Way	WARRAGUL	VIC	3820	03 5623 6966	03 5623 6966
Geelong Wheelchairs	106 Bailey St	GROVEDALE	VIC	3216	03 5244 0844	03 5244 0850
Geelong Wheelchairs	215 Murray St	COLAC	VIC	3707	03 5231 6050	03 5244 0850
Life Mobility	24 Corporate Blvd	BAYSWATER	VIC	3153	03 9729 5552	03 9729 5112
Life Mobility	34 Railway Ave	RINGWOOD EAST	VIC	3135	03 9879 8885	03 9879 3133
Mitchell Healthcare	14a Sydney Rd	BENALLA	VIC	3672	03 5762 1125	03 5762 1124
Mitchell Healthcare	95 Murphy St	WANGARATTA	VIC	3677	03 5721 5660	03 5721 9667
Omni Healthcare	206 Creswick Road	BALLARAT	VIC	3350	03 5333 4006	03 5333 3825
Omni Healthcare	11 Fairy Street	WARRNAMBOOL	VIC	3280	03 5562 6966	03 5562 8944
Whitecross Healthcare	51 West Fyans St	NEWTON	VIC	3220	03 5222 6666	03 5222 6600

WA						
AC Mobility	5 Boulder Rd	MALAGA	WA	6090	08 9209 1777	08 9209 1999
Australian Medical Supplies	56 Catalano Cct	CANNING VALE	WA	6155	1300 770 921	1300 770 945
Mandurah Mobility	45 Reserve Dr	MANDURAH	WA	6210	08 9535 1411	08 9535 1755

Country Care Group





National Distribution, Local Service

FREE CALL 1800 727 382 Facsimile 1800 329 382 Email contracts@countrycaregroup.com.au

countrycaregroup.com.au



CLINICAL FURNITURE CONTRACT CATALOGUE - 217B