

- Tight economy
- Scarce resources

Pressure Ulcers

Try Levabo's product range of inflatable one-patient cushions

- Hygienic
- Comfortable
- Easy
- Economic



LEVABO[®]
medico cushions



Cushions for pressure relief of heels and legs



| Name | | Size in cm. | Item no. | Quantity | Application | About the cushion |
|------------------|--|--------------|----------|----------|-------------|---|
| Heel up - Short | | 18 x 32 x 25 | 72001 | 12 pcs | | Heel protector for mobile and seated users, typically wheelchair. Can be used for both prevention and treatment of pressure sores. |
| Heel up - Medium | | 18 x 39 x 25 | 72007 | 12 pcs | | Heel protector useful for both prevention and treatment of pressure sores. Typically for users with short legs. Most useful for bed-bound patients. It supports the front foot and minimizes risk of foot drop. |
| Heel up - Long | | 18 x 48 x 25 | 72002 | 12 pcs | | Heel protector useful for both prevention and mitigation of pressure sores. Most useful for bed-bound patients. It supports the front foot and minimizes the risk of foot drop. |
| Heel up - Max | | 18 x 50 x 30 | 72018 | 12 pcs | | Heel protector useful for both prevention and mitigation of pressure sores. Most useful for bed-bound patients and for users with a long foot. It has an extra air chamber under the foot, which provides an optimal support of the front foot. |



| | | | | | | |
|-----------------------|--|--------------|-------|--------|--|---|
| All up - Multi Small | | 50 x 25 x 8 | 72003 | 12 pcs | | Positioning cushion with 2 chambers. Multifunction. Can be used for both supine and sitting positions. (between legs, back, very good support for the arm and between the knees). |
| All up - Multi Medium | | 50 x 27 x 8 | 72004 | 12 pcs | | Positioning cushion with 3 chambers. Multifunction. Especially useful for a supine position (between legs, back), and possibly as support under the feet. |
| All up - Multi Large | | 50 x 50 x 8 | 72005 | 12 pcs | | Positioning cushion 4 chambers. Multifunction. With the possibility of being folded/rolled up. (Under the leg and arm; comes with extra velcro strap). |
| All up - Multi Max | | 50 x 50 x 20 | 72012 | 6 pcs | | Positioning cushion, particularly useful for hip fracture patients and when the lower limb need to be elevated above heart level. |
| All up - Wing | | 50 x 85 x 8 | 72006 | 6 pcs | | Relief cushion for lower limbs, used for bed-bound patients. The cushion is easily fixed to the mattress. |

Accessories

| Name | Size in cm. | Item no. |
|---------------------------|-------------|----------|
| Heel up - Inlay X 10 | 35 x 50 | 76008 |
| All up - Cover Medium X 5 | 50 x 50 | 76012 |
| All up - Cover Large X 5 | 50 x 100 | 76013 |





Awards


Levabo was awarded The Innovation Elder Care Award 2015 (Innovationspreis Altenpflege 2015) for the product series Heel up® on Germany's largest elder care conference.

The panel of judges especially emphasized how easy applicable the products are, as the inflator indicator Stop Up® ensures correct inflation without the user being required to have extensive training on how to inflate the cushions.

Additionally, the products were praised for their compact size not taking up space in the depots, and them being easy to transport to patients, particularly and advantage for the domiciliary care.



Remember the exposed areas



The first sign of a pressure sores can develop within a very short time, 2 – 4 hours !

Pressure sores typically develop on the exposed areas of the body, where the bones are only protected by a thin layer of skin.

With Levabo's flexible one-patient product range it is possible to make a quick and effective positioning of the patient before the exposed areas are affected.

Cushions for pressure relief in the seated position



| Name | Size in cm. | Item no. | Quantity | Application | About the cushion |
|--------------------------------|-------------|----------|----------|-------------|---|
| All up - Seat | 42 x 42 x 5 | 72008 | 12 pcs | | Seat cushion, extra strong. Pressure distribution for the mitigation of pressure sores. Standard wheelchair size. Max 150 kg. |
| All up - Seat Large | 52 x 52 x 5 | 72019 | 12 pcs | | Seat cushion, extra strong. Pressure distribution for the prevention and mitigation of pressure sores. Extra length and width. Max 150 kg. |
| All up - Seat Max (bariatrisk) | 52 x 62 x 5 | 72020 | 12 pcs | | Seat cushion, extra strong. Bariatric cushion for user weights up to 250 kg. Pressure distribution for the prevention and mitigation of pressure sores. Extra length and width. |
| All up - Seat Twin | 48 x 52 x 5 | 72021 | 12 pcs | | Seat cushion, extra strong. Suitable for hemiplegic patients with unstable posture. The cushion is split in two to adjust for individual air pressure. For the prevention and relief of pressure sores. Max 150 kg. |
| All up - Seat Tailbone | 42 x 42 x 5 | 72022 | 12 pcs | | Special sitting cushion, extra strong. To apply for full relief of the rectum, typically for radiation damage, a sensitive tailbone or great discomfort caused by hemorrhoids. Max 150 kg. |
| All up - Back & Seat | 27 x 45 x 5 | 72009 | 12 pcs | | Sitting cushion and back cushion. With advantage, usefull for both areas. Good back support, easy cushion for correct seating. Max 150 kg. |
| All up - Donut | 45 x 45 x 8 | 72017 | 12 pcs | | Seat cushion, extra strong. Annular (ring shaped) with no hole in the middle, but with optimum pressure to the relief and prevention of pressure sores. Max 150 kg. |
| All up - Bath | 50 x 50 x 8 | 72013 | 12 pcs | | Special seat cushion, extra strong. For bath and toileting chairs to offer a better distribution of pressure for prolonged use. Can be fastened to the chair with tape. After bathing, the cushion can be dried and reused. Max 150 kg. |

Prevention and treatment

Cushions for special positioning

| | | | | | |
|-----------------|--------------|-------|--------|--|---|
| All up - Belly | 65 x 75 x 15 | 72023 | 6 pcs | | Special positioning cushions with particular application to the radiology department to lift the chest at special examinations. |
| All up - Kidney | 50 x 20 x 20 | 72014 | 12 pcs | | Special positioning cushion for special for a Radiology Department for seating kidney/pelvic area. |

Levabo Emergency Package

Our products really make a difference.

The innovative positioning products from Levabo[®] are developed in collaboration with leading wound care nurses and are produced in Denmark. The products have a series of advantages compared to conventional positioning products.


With the use of Levabo[®]'s positioning cushions you always have the needed cushion at hand for effective prevention of pressure sores. The products are one-patient, inflatable positioning cushions, which can be disposed after use. The products take up very little space, both useful for the depot, but also for easy transportation of Levabo[®]'s emergency package, developed for that purpose.

These advantages are especially useful for the domiciliary care. As a domiciliary care nurse you can initiate effective prevention of pressures sores at an early stage with Levabo[®]'s products, easier and cheaper than with conventional positioning products. Already at suspicions of the development of pressure sores, the domiciliary nurse can immediately position the patient optimally, and thereby avoid unnecessary discomfort for the patient, as well as save the health care system for potential treatment expenses.

The products are multifunctional and can with advantage be used in several different positions. The one-patient use ensures high hygiene.



| Name | Size in cm. | Item no. |
|-------------------------|-------------|----------|
| All up - Seat Cover X 5 | 42 x 42 x 5 | 76009 |

















Prevention of pressure sores

Accessories



| Name | Size in cm. | Item no. |
|-------------------------------|-------------|----------|
| All up - Seat Cover X 5 | 42 x 42 x 5 | 76009 |
| All up - Seat Large Cover X 5 | 52 x 52 x 5 | 76010 |
| All up - Overlay Medium X 5 | 50 x 50 | 76012 |
| All up - Overlay Large X 5 | 50 x 100 | 76013 |

Cushions for pressure relief of body and arms

| Name | Size in cm. | Item no. | Quantity | Application | About the cushion |
|-------------------------|--------------|----------|----------|---|--|
| All up - Multi Small | 50 x 25 x 8 | 72003 | 12 pcs |  | Positioning cushion with 2 chambers. Multifunction. Can be used for both supine and sitting positions. (between legs, back, very good support for the arm and between the knees). |
| All up - Multi Medium | 50 x 27 x 8 | 72004 | 12 pcs |  | Positioning cushion with 3 chambers. Multifunction. Especially useful for a supine position (between legs, back), and possibly as support under the feet. |
| All up - Multi Large | 50 x 50 x 8 | 72005 | 12 pcs |  | Positioning cushion 4 chambers. Multifunction. With the possibility of being folded/rolled up. (Under the leg and arm; comes with extra velcro strap). |
| All up - Multi Max | 50 x 50 x 20 | 72012 | 6 pcs |  | Positioning cushion, particularly useful for hip fracture patients and when the lower limb need to be elevated above heart level. |
| All up - Universal | 50 x 50 x 8 | 72010 | 12 pcs |  | Positioning cushion, typically for the back, sacrum and tailbone in adjustable bed and the lateral positions between the legs. Can be folded if a higher cushion is desired. |
| All up - Universal Long | 55 x 70 x 8 | 72015 | 12 pcs |  | Positioning cushion, typically for the lateral positions between the legs, longer model which protects the knees and malleolus. Can be folded in half if a higher cushion is desired. |
| All up - Universal Slim | 45 x 75 x 8 | 72016 | 12 pcs |  | Positioning cushion, typically for lateral positions between the legs, longer and narrower model which ensures a distance between the knees in the lateral position. Can be folded if a higher cushion is desired. |
| All up - Head | 28 x 30 x 5 | 72011 | 12 pcs |  | Cushion, particularly useful for the preventing and treating of pressure sores on the ears and neck, and is also useful for coma patients. |
| All up - Wing | 50 x 85 x 8 | 72006 | 6 pcs |  | Relief cushion for lower limbs, used for bed-bound patients. The cushion is easily fixed to the mattress. |
| All up - Wedge 30 deg. | 70 x 40 x 15 | 72024 | 6 pcs |  | Positioning cushion for the lateral position. Ensures a stable position at 30°, and the wedge shape is fixed by a flap that is pushed underneath the patient. |
| All up - Lap | 80 x 60 x 10 | 72025 | 6 pcs |  | Positioning cushion, crescent. Extra height. Advantageously used for wheelchair users for relief and support of arms. |
| All up - Lap Max | 80 x 60 x 20 | 72026 | 6 pcs |  | Positioning cushion, crescent. Advantageously used for wheelchair users for relief and support of arms. |
| All up - Arm Rest | 55 x 25 x 8 | 72027 | 12 pcs |  | Positioning cushion for relief of arms. |
| All up - Elbow | 30 x 20 x 10 | 72028 | 12 pcs |  | Positioning cushion for relief of the elbow. |



Levabo's products are space saving in the depot.



Accessories

| Name | Size in cm. | Item no. |
|---------------------------|-------------|----------|
| All up - Velcro Strap X 1 | 55 x 200 | 76011 |
| All up - Cover Medium X 5 | 50 x 50 | 76012 |
| All up - Cover Large X 5 | 50 x 100 | 76013 |

Levabo Pumps


All Levabo® products are delivered with an inflation adapter set, which is included in the packaging of the individual cushion.

| Name | Size in cm. | Item no. |
|---------------------------------|--------------|----------|
| Levabo Manual Pump | 29 x 6 x 6 | 76004 |
| Levabo Battery Pump | 26 x 5 x 5 | 76001 |
| Levabo 240V Pump | 24 x 13 x 11 | 76002 |
| Levabo inflate adapter set X 10 | | 76005 |



Always correct inflation

With Levabo®'s Stop Up® patent system you are always sure to make a correct inflation. The Stop Up® strip indicates when the right amount of air is in the cushion. Check out the video at: www.levabo.com




Levabo®'s Stop Up® system with easily to understandable pictograms.

- Single-patient use ensures increased hygiene.
- Heel up® is designed for the treatment and prevention of bedsores in the heel region.
- All up® was designed for the prevention and treatment of bedsores on the elbows, shoulders, back, etc.
- Easy & simple to use.
- Good comfort
- Soft, skin-friendly, heat dissipation and moisture-absorbent material.
- Space saving in the storage area.
- Delivered in clean, sealed and flat packaging, including an inflation adapter set.
- Tested and approved according to ISO 10993-10 skin irritation and hypersensitivity.
- Available in several sizes.

Levabo Guide


For your iPhone, iPad & Android



Get the free Levabo app on:

<http://www.levabo.dk/apps.html>

Or scan the QR code:



Dealer



www.levabo.dk
www.facebook.com/LevaboAps



LEVABO®
 medico cushions

Levabo ApS
 Højlundsparken 49
 8355 Solbjerg
 Tel.: +45 20 70 09 30
 info@levabo.dk
 www.levabo.com