

Safety and Efficacy of a Leadless Pacemaker: Results from the LEADLESS II clinical trial

Vivek Y. Reddy, MD, T. Jared Bunch, MD, Daniel J. Cantillon, MD, Rahul Doshi, MD, N.A. Mark Estes, MD, Derek V. Exner, MD, Paul Friedman, MD, Gery Tomassoni, MD, John Ip, MD, Kenneth Plunkitt, MD

Icahn School of Medicine at Mount Sinai, New York, NY, Intermountain Medical Center Heart Institute, Salt Lake City, UT, Cleveland Clinic, Cleveland, OH, USC University Hospital, Los Angeles, CA, Tufts University, Boston, MA, Libin Cardiovascular Institute of Alberta, Calgary, Canada, Mayo Clinic, Rochester, MN, Central Baptist Hospital, Lexington, KY, Sparrow Research, Lansing, MI, Naples Community Hospital, Naples, FL.



- COI: St Jude Medical Inc – Grant support & Consultant
- I will be discussing the use of non-FDA approved devices



DECLARATION OF INTEREST


- Research contracts
- Consulting/Royalties/Owner/ Stockholder of a healthcare company



Today's Leadless Pacemaker System

The Nanostim Device

- **Percutaneous femoral vein delivery**
 - 18F introducer /steerable catheter
 - <30 minute skin-to-skin procedure
- **Self-contained device in ventricle**
 - No lead or surgical pocket
 - Inherently MRI compatible
- **Conventional Features**
 - Temperature-Based Rate Response
 - >10-yr battery life
 - Hysteresis
 - Magnet Mode
- **Flexible replacement options**
 - Catheter-based retrieval
 - Deliver additional leadless pacemakers
 - Revert to conventional pacing lead



Leadless II Clinical Trial

Overview

- Prospective, multicenter, non-randomized, FDA IDE study
- Objective:
 - To evaluate the clinical safety and efficacy of non-surgical implantation of the leadless cardiac pacemaker in patients indicated for a VVI(R) pacemaker.
- Primary Cohort: 1st 300 patients followed for 6mo (June 2015)
- Total Cohort: All patients enrolled by June 2015 (n=526)
- Primary Endpoints (by ITT):
 - Safety: Freedom from Serious Adverse Device Effects at 6 months
 - Efficacy: Acceptable pacing capture threshold (d2.0 V at 0.4 msec) and a therapeutically acceptable sensing amplitude (R wave ≥ 5.0 mV, or a value equal to or greater than the value at implantation) through 6 mo.
- 56 Centers in the US, Canada and Australia
 - 100 Operators (only one had prior experience with leadless pacing)



Leadless II Clinical Trial

Outcomes

- **Device was successfully implanted in ~96% of patients**
- **Primary Safety Endpoint (Intent-to-Treat Analysis)**
 - 280 of the 300 patients achieved endpoint (93.3%; 95% CI = 89.9 to 95.9)
 - This exceeded the performance goal of 86% (P<0.001)
- **Primary Efficacy Endpoint (Intent-to-Treat Analysis)**
 - 270 of the 300 patients achieved endpoint (90.0%; 95% CI = 86.0 to 93.2)
 - This exceeded the performance goal of 85% (P = 0.007)
- **Efficacy (Successful implants)**
 - 289 patients with successful device implant
 - 270 of the 289 patients achieved endpoint (93.4%; 95% CI = 89.9 to 96.0)
 - This exceeded the performance goal of 85% (P <0.001)

- Based on device-use characteristics, the battery longevity is estimated to be **15.0 ± 6.7 yrs** (95% CI, 14.2 to 15.8 yrs)

Leadless II Clinical Trial

Device-Related SAEs

Event	Primary Cohort (N = 300)			Total Cohort (N = 526)		
	No. of Events	No. of Patients	Event Rate	No. of Events	No. of Patients	Event Rate
			%			%
Total	22	20	6.7	40	34	6.5
Cardiac perforation	4	4	1.3	8	8	1.5
Cardiac tamponade with intervention	1	1	0.3	5	5	1.0
Cardiac perforation requiring intervention	1	1	0.3	1	1	0.2
Pericardial effusion with no intervention	2	2	0.7	2	2	0.4
Vascular complication	4	4	1.3	6	6	1.1
Arrhythmia during device implantation	2	2	0.6	3	3	0.6
Cardiopulmonary arrest during implantation procedure	0	0	0	1	1	0.2
Device dislodgement	5	5	1.7	6	6	1.1
Device migration during implantation owing to inadequate fixation	0	0	0	2	2	0.4
Pacing threshold elevation with retrieval and implantation of new device	4	4	1.3	4	4	0.8
Other *	3	3	0.9	10	10	1.9

* Includes: ischemic stroke, angina pectoris, pericarditis, acute confusion & expressive aphasia, dysarthria & lethargy post implant, contrast induced nephropathy, orthostatic hypotension with weakness, left leg weakness during implant, probable pulmonary embolism, ischemic stroke



Leadless II Clinical Trial

Conclusions

- The Leadless Pacemaker was successfully implanted in ~96% of attempted patients.
- The trial met the pre-specified Safety and Efficacy endpoints
 - Complication rate similar to that seen with conventional pacemakers
 - Complication rate likely to improve with operator experience (Rem: in this study, 99 of 100 operators had never implanted a leadless device)
- The device was shown to be retrievable in a subgroup of patients (n=7) who needed a replacement (Time from implant = 160 ± 180 days; Range = 1 to 413 days)
- The estimated device longevity based on the 6-month follow-up duration is encouraging



Leadless II Clinical Trial

Limitations

- An Observational study (not Randomized)
- Mean Follow-Up of only 6 months
- How to manage device after battery depletion?
 - Possible to retrieve after ~1 year, but what about 5, 10, 15 yrs?
 - Retrieval *vs* Abandonment
- Limited device diagnostics (eg, no electrogram data)
- Large venous sheath (18Fr)
 - Now increasingly common used for cardiology procedures
 - Low observed rate of hematomas
- Single-chamber (RV) pacing only
 - Device-to-device communication is in development
 - Would lead to dual-chamber, CRT, etc

