

Shielding for Multislice CT Scanners

Donna M. Stevens, MS
 Imaging Physics Department
 Diagnostic Imaging Division
 The University of Texas M.D. Anderson Cancer Center
 Houston, Texas

AAPM Summer School
 St John's University
 Collegeville, MN
 27 July 2007

Site Planning for CT

- Location, room dimensions
- Equip type
- Electrical
- HVAC
- Water
- Structural loading
- Shielding

... and one or two
 \$\$\$\$\$\$

MDACC Imaging Physics

The Shielding Problem

- An administrator's office is adjacent to a CT scan room. The administrator sits 4m from the scanner. How much shielding is required?

3 mm Pb

If only it were this easy!!!

MDACC Imaging Physics

Project Planning


- Neighboring spaces
 - Adj occup factors if needed
 - Dist of closest approach, $\geq 0.3\text{m}$
- Design goals or limits
 - Public, Controlled, Pregnant worker
 - Adjust if needed
- Workload estimate



MDACC Imaging Physics

Preliminary Information

- Architectural drawings (Plan view) of exam room, floor above, and floor below
 - Elevation sections for floor and ceiling
 - Occupancy factors for floors above and below
 - Two rooms away or across hall (remote areas may be more sensitive than adjacent)
- Composition of walls, ceilings and floors
 - Materials and thickness
- Scanner placement from vendor
 - Distance to protected areas beyond barriers
- Scatter contributions from other rooms/floors

MDACC Imaging Physics





Workload Estimates

- Number procedures per week
 - 3 patients per hour
 - 1 to 5 procedures per patient
- Body% and Head%
- Contrast and non-contrast scans
- 120kVp vs other kVp
- Scan parameters of protocol

MDACC Imaging Physics

Multi-Slice Helical CT Shielding

- Thinner slice protocols may require more dose – create more scatter
 - More photons needed to generate adequate photon statistics per slice (smaller voxels, higher noise)
 - Environmental radiation levels typically increase with increase in beam width
 - However, fewer rotations are needed to produce the scan

MDACC Imaging Physics

Multi-Slice Helical CT Shielding

- Over-scan at ends of volume add scatter

- Worst with widest beams
- Ceiling and floor deserve close scrutiny

MDACC Imaging Physics

Barrier Determination

- NCRP 147, section 5.6, pg 94
- CTDI method
- DLP method
- Scatter plots

MDACC Imaging Physics

CTDI Method

- What we need to know:
 - Pitch = table/beam
 - Beam width: T_b or (nT)
 - Rotation time
 - Peripheral CTDI100
 - measure (can scale by kVp²)
 - Look on ImPACT website

CTDI Method

Unshielded weekly exposure calculation:

$$\text{Secondary exposure per procedure at one meter } K_s = \kappa \times \left[\frac{L}{p} \right] \times \left[\text{mAs/Rotation} \right] \times \left[\text{CTDI}_{100, \text{peripheral}} / \text{mAs} \right] \times \left[\frac{\text{Scan kV}}{\text{CTDI kV}} \right]^2$$

Where:

κ is the scatter fraction at one meter per cm scanned.
L is the length of the scanned volume.
p is pitch.

κ (head)	$9 \times 10^{-5} \text{ cm}^{-1}$
κ (body)	$3 \times 10^{-4} \text{ cm}^{-1}$

CTDI Method

- ImPACT (the UK's CT evaluation center)
- measured axial and peripheral CTDI₁₀₀
- for most scanners on the market
- Excel format

www.impactscan.org

NCRP 147 DLP Method

DLP (Dose-Length Product)

- $\text{CTDI}_W = 1/3 \text{ ctr CTDI}_{100} + 2/3 \text{ Surf CTDI}_{100}$
- $\text{CTDI}_{VOL} = \text{CTDI}_W / \text{Pitch}$
- $\text{DLP} = \text{CTDI}_{VOL} * L$
- L = Scan length for average **series** (cm)
- Units of mGy-cm
- From scanner display ... verify these values!

$$\text{DLP} = [1/3 \text{ CTDI}_{100, \text{Center}} + 2/3 \text{ CTDI}_{100, \text{Surface}}] * L/p$$

NCRP 147 DLP Method

- Weekly Air Kerma at 1m (K_{sec}^1)

$$K_{\text{sec}}^1 (\text{head}) = \kappa_{\text{head}} * \text{DLP}$$

$$K_{\text{sec}}^1 (\text{body}) = 1.2 * \kappa_{\text{body}} * \text{DLP}$$

Factor of 1.2 assumes peripheral CTDI₁₀₀ = 2*Center CTDI₁₀₀ for Body



$$\kappa_{\text{head}} = 9 \times 10^{-5} \text{ cm}^{-1}$$

$$\kappa_{\text{body}} = 3 \times 10^{-4} \text{ cm}^{-1}$$

Use inverse square to find unshielded weekly exposure at barrier from K_{sec}^1

Scatter plot Method

- Assume an isotropic exposure distribution w/ the vendor-supplied scatter distribution plots (max is approx. 45° to the scanner axis).
- Overestimates shielding for gantry shadows
- W, in mA*min per week
- Determine weekly exp at shielded point
- Pay attention to
 - Beam width
 - kVp and mAs
 - phantom



- ### ! Use Caution with Scatter Plots
- Choice of plot (Head vs Body)
 - Normalization of data
 - kVp of plot vs clinical
 - mAs per scan
 - Beam width of plot vs clinical
 - Total mAs per scan
 - Pitch, rotation, total beam-on time
 - Accounts for scan acquisition time for diff beam width

Question

Do I really need to put lead in the ceiling of a 64-slice CT scan room?

- ### Method
- Calculate the unshielded weekly exposure rate at 0.5 m beyond the floor above
 - Find the maximum weekly exposure at 1 m from isocenter and inverse-square this out to the occupied area beyond the barrier.
 - Perform barrier thickness calculations
 - Occupancy, permissible dose, attenuation of concrete, etc.

- ### Example
- 180 Procedures/week
 - 150 Abdomen & Pelvis $L = 60\text{cm} + 0.4 \cdot 30\text{cm}$
 - 30 Head $L = 1.4 \cdot 15\text{cm}$
 - 40% w&w/o contrast
 - 13.0' (4.2 m) ceiling height (finished floor to finished floor) $D_{\text{sec}} = 3.7\text{m}$
 - GE LightSpeed 16
 - Ignores overscan at ends!
 - Effect Worsens with wider beams (64-slice)

Protocols

	kVp	mA	Time (sec)	Pitch	Beam (mm)	Table (mm/rot)
Head	120	240	1.0	1.375	10	13.75
Body	120	265	0.8	0.938	20	18.75

NCRP 147 DLP Method

Procedure	CTDI _{vol} (mGy)	Scan Length (L, cm)	DLP (mGy-cm)
Head	60	20	1200
Body	15	35	525
Abdomen	25	25	625
Pelvis	25	20	500
Body (Chest, Abdomen, or Pelvis)			550

Unshielded Weekly Exposure at Barrier

- Average Air Kerma/procedure at 1m (K_{sec}^1) – 40% w&w/o contrast

$$K_{sec}^1 (\text{head}) = 1.4 * K_{head} * \text{DLP}$$

$$= 1.4 * 9 \times 10^{-5} \text{ cm}^{-1} * 1200 \text{ mGy-cm}$$

$$= 0.15 \text{ mGy}$$

$$K_{sec}^1 (\text{body}) = 1.4 * 1.2 * K_{body} * \text{DLP}$$

$$= 1.4 * 1.2 * 3 \times 10^{-4} \text{ cm}^{-1} * 550 \text{ mGy-cm}$$

$$= 0.28 \text{ mGy}$$

Unshielded Weekly Exposure at Barrier

- Weekly Air Kerma (K_{sec}) at Ceiling:
 - 30 head procedures/wk
 - 150 body procedures/wk
 - $D_{sec} = 4.2 \text{ m} + 0.5 \text{ m} - 1 \text{ m} = 3.7 \text{ m}$

$$K_{sec} (\text{head}) = 30 * 0.15 \text{ mGy} * (1\text{m}/3.7\text{m})^2$$

$$= 0.33 \text{ mGy}$$

$$K_{sec} (\text{body}) = 150 * 0.28 \text{ mGy} * (1\text{m}/3.7\text{m})^2$$

$$= 3.04 \text{ mGy}$$

Unshielded Weekly Exposure at Barrier

- Weekly Air Kerma (K_{sec}) at Ceiling:

$$K_{sec} (\text{Total}) = K_{sec} (\text{head}) + K_{sec} (\text{body})$$

$$K_{sec} (\text{Total}) = 0.33 \text{ mGy} + 3.03 \text{ mGy}$$

$$K_{sec} (\text{Total}) = 3.37 \text{ mGy}$$

Required Transmission (B)


$$B = \frac{P}{K_{sec} * T}$$

P = Maximum permissible weekly exposure
T = Occupancy Factor

$$= \frac{0.02 \text{ mGy}}{3.37 \text{ mGy} * 1} = 5.9 \times 10^{-3}$$

Total Shielding Required

Use Simpkin curve fit equations or look up on published attenuation diagrams (NCRP 147 Fig. A-2)



MDACC Imaging Physics


31

Existing Shielding

- Determine attenuation of existing barriers with Tc-99m source and Na-I detector
- Determine lead-equivalence of barrier
- Floors and ceilings
 - Find lead equivalence from documentation of concrete thickness.
 - If necessary, Find thickness by drilling a test hole and measuring.
 - Always assume light weight concrete, unless proven otherwise (30% less dense than standard density, coefficients used in NCRP 147)

MDACC Imaging Physics

32



Existing Shielding

- Subtract existing lead-equivalence from total required
- Convert to 1/32 inch multiples (round up)
- Total lead to add = (Total required) – (Existing)
 - = 1.37 mm – 0.45 mm
 - = 0.9 mm
 - Round up to 1/16" Pb

MDACC Imaging Physics

34



Comparison of Methods

	DLP NCRP 147		CTDI ₁₀₀		Scatter plot	
	Head	Body	Head	Body	Head	Body
K _{sec} @ 1m	4.5	42	2.0	101	2.9	95
Weekly kerma @ 3.7m	3.37		7.53		7.11	
Total Barrier (mm Lead)	1.4		1.7		1.7	

MDACC Imaging Physics





36

Ceiling Considerations

- Pb mounting in ceiling is manually applied and très cher! (very expensive!)
- Isotropic distribution is conservative, but not so realistic
- Consider % of scans helical w/o gantry tilt (tilted axials usually for Head only)
- Smaller area of ceiling to cover = smaller cost ... THIS time
- Additional cost possibly incurred in future renovation

MDACC Imaging Physics

37





Watcha Gonna Do? #1

- Attended waiting room adjacent to CT room
- New PET-CT to be installed on floor below

CT

New PET-CT

Attended waiting

MDACC Imaging Physics 47

Watcha Gonna Do? #2

- CT room on 3rd floor, exterior wall, standard windows
- Lab area across driveway
- Current kerma in labs OK, but close to limit
- New PET-CT to be added adj to existing CT

Lab space

Driveway

Existing CT
New PET-CT

MDACC Imaging Physics 48

Thank you!

Acknowledgements:

Jeff Shepard

Bud Wendt