

ASCRS ♦ ASOA Symposium & Congress
Technicians & Nurses Program
May 6-10, 2016 – New Orleans



Colour Vision Testing

Craig Simms

BSc, COMT, CDOS, ROUB

Chief Ophthalmic Medical Technologist
Program Coordinator, KOTC
Kingston, Ontario

JCAHPO ACE Meeting 2015
Las Vegas



Objectives

- Basic review of colour vision
- Discuss colour vision testing with pseudoisochromatic plates
 - Ishihara plates
 - Hardy-Rand-Rittler plates
- Discuss colour vision testing with D-15
- Discuss colour vision with FM 100 Hue
- Apply colour vision testing in a case study

Colour Vision

- Dependent upon healthy cones
- Cones are thought to contain 3 different photosensitive pigments
 - Red
 - Green
 - Blue



Why do we see in colour?

- To tell the difference between coloured berries
- To tell the good from rotten food
- To tell ripe from unripe food
- Mark Changizi, The Vision Revolution
 - Subtle and not so subtle changes in skin colour (changes in blood oxygenation)
 - primates with colour vision have evolved to have bare patches on their bodies (face, rump, genitals)
 - 10% of all medical disorders list skin colour as a diagnostic sign



Colour Vision

- Colour depends on
 - Hue
 - The colour of an object
 - An object will have a particular hue because it reflects or transmits light of a certain wavelength
 - Shades
 - Addition of black to a hue



Colour Vision

- Colour depends on
 - Saturation
 - Index of purity of a hue
 - The amount of white (highly saturated has little white)
 - Pink is less saturated than red
 - Intensity
 - Brightness of a colour



Colour Vision

- White objects
 - Reflect most light
- Black objects
 - Absorb most light



Colour Vision Defects

- Believed to arise from a deficiency or absence of one or more visual pigments
- Most common
 - Congenital
 - Red-green deficiency
 - X-linked recessive
 - 8% to 10% males
 - 0.4% females



Colour Vision Defects Inheritance

- Females - XX
- Males - XY
- X chromosome carries the defect (recessive)
- Female inherits it from her colour deficient father
 - Normal colour vision from the normal X
 - Carrier
 - Any of her sons that get the defective X will have the colour deficiency



Acquired Colour Vision Defects

- Secondary to disease or trauma
- Macular problems
 - Anything affecting the cones
- Optic nerve problems
 - TRO
 - Optic neuritis
- Chemical poisoning



Colour Vision Defects

Congenital

- Bilateral
- Affects both eyes equally
- Males > females
- Non-progressive
- No cure

Acquired

- Unilateral or bilateral
- Can affect eyes differently
- Males = females
- Progressive
- May go away when cause is removed



Colour Vision Defects Categories

- Trichromatism
- Dichromatism
- Monochromatism



Trichromatism

- Possess all 3 cone pigments
- Normal colour vision



Anomalous Trichromatism

- Partial deficiency of one of the 3 cone pigments
 - Protanomaly
 - Red deficiency
 - Poor red-green discrimination
 - Poor blue-green discrimination



Anomalous Trichromatism

- Partial deficiency of one of the 3 cone pigments
 - Deuteranomaly
 - Green deficiency
 - Poor green-purple discrimination
 - Poor red-purple discrimination



Anomalous Trichromatism

- Partial deficiency of one of the 3 cone pigments
 - Tritanomaly
 - Blue deficiency
 - Poor blue-green discrimination
 - Poor yellow-green discrimination



Dichromatism

- Complete deficiency of one cone pigment
 - Protanopia
 - Absent red cone pigment



Protanopia

- Absent red cone pigment

Used with
permission



Dichromatism

- Complete deficiency of one cone pigment

- Deuteranopia

- Absent green cone pigment



Deuteranopia

- Absent green cone pigment



Used with
permission



Dichromatism

- Complete deficiency of one cone pigment

- Tritanopia

- Absent blue cone pigment



Tritanopia

- Absent blue cone pigment



Used with
permission



Monochromatism

- Only one cone pigment
- Require at least 2 different cone pigments to see any kind of colour
- Will see only shades of gray




Colour Vision Tests

- Pseudoisochromatic plates
 - Ishihara Plates
 - Hardy-Rand Rittler Plates
- Farnworth Munsell D-15
- FM 100 Hue



Ishihara Plates


- Plates composed of coloured dots, differing in shade and hue, that form numbers and patterns that are not visible with color vision defects
- Provides assessment of congenital defects
- Different versions



Hotel Dieu HOSPITAL

Ishihara Plates

Plate 1




Red-green defect	Total colour deficiency
12	12

Hotel Dieu HOSPITAL

Ishihara Plates

Plate 2




Red-green defect	Total colour deficiency
3	-

Hotel Dieu HOSPITAL

Ishihara Plates

Plate 3




Red-green defect	Total colour deficiency
70	-

Hotel Dieu HOSPITAL

Ishihara Plates

Plate 4




Red-green defect	Total colour deficiency
2	-

Hotel Dieu HOSPITAL

Ishihara Plates

Plate 5






Red-green defect	Total colour deficiency
5	-

Hotel Dieu HOSPITAL

Ishihara Plates


Plate 6



	Red-green defect	Total colour deficiency
	17	-

Ishihara Plates


Plate 7



	Red-green defect	Total colour deficiency
	21	-

Ishihara Plates


Plate 8



	Red-green defect	Total colour deficiency
	-	-

Ishihara Plates


Plate 9



	Red-green defect	Total colour deficiency
	-	-

Ishihara Plates


Plate 10



	Red-green defect	Total colour deficiency
	-	-

Ishihara Plates


Plate 11

	Red-green defect	Total colour deficiency
	-	-

Ishihara Plates

•Plate 12




Red-green defect	Total colour deficiency
-	-

Hotel Dieu HOSPITAL

Ishihara Plates

Plate 13




Red-green defect	Total colour deficiency
-	-

Hotel Dieu HOSPITAL

Ishihara Plates

Plate 14




Red-green defect	Total colour deficiency
5	-

Hotel Dieu HOSPITAL

Ishihara Plates

Plate 15




Red-green defect	Total colour deficiency
45	-

Hotel Dieu HOSPITAL

Ishihara Plates

Plate 16




Protan (red)		Deutan (green)	
strong	mild	strong	mild
6	(2)6	2	2(6)

Hotel Dieu HOSPITAL

Ishihara Plates

Plate 17




Protan (red)		Deutan (green)	
strong	mild	strong	mild
2	(4)2	4	4(2)

Hotel Dieu HOSPITAL

Ishihara Plates

Plate 18




Protan (red)		Deutan (green)	
strong	mild	strong	mild
Purple	Purple (red)	Red	Red (purple)

Hotel Dieu HOSPITAL

Ishihara Plates

Plate 19




Normal	Red – green defect	Total colour deficiency
-	Trace along a line	-

Hotel Dieu HOSPITAL

Ishihara Plates

Plate 20




Normal	Red – green defect	Total colour deficiency
Trace bluish-green line	-	-

Hotel Dieu HOSPITAL

Ishihara Plates

Plate 21




Normal	Red – green defect	Total colour deficiency
Trace orange line	-	-

Hotel Dieu HOSPITAL

Ishihara Plates

Plate 22




Normal	Red – green defect	Total colour deficiency
Trace greenish line	Trace greenish line with the purple	-

Hotel Dieu HOSPITAL

Ishihara Plates

Plate 23




Normal	Red – green defect	Total colour deficiency
Trace purple and orange line	Trace purple and bluish green line	-

Hotel Dieu HOSPITAL

Ishihara Plates

Plate 24



Ishihara Plates

•Testing: Monocular? Binocular?

•Review colour vision defects (briefly and basically)

Congenital

- Bilateral
- Affects both eyes equally
- Males>females
- Non-progressive
- No cure

Acquired

- Unilateral or bilateral
- Can affect eyes differently
- Males=females
- Progressive
- May go away when cause is removed



Ishihara Plates

•Testing: Monocular? Binocular?

•Acquired (TRO, chemical poisoning, etc)

- Monocular

•Congenital

- Binocular



Ishihara Plates Scoring

•Scoring

•Give 1 point for each correct number

•Total number correct out of the total number of plates

•15/15

•9/15



Ishihara Plates Scoring

•Scoring

•Score using plate # 1 or not

•Since everyone should see it

•0/14 is equivalent to 1/15



Ishihara Plates Scoring

•Scoring (using 15 plates in the book)

•9 or less plates are read correctly = colour deficiency

•Track which numbers are missed

- Determine if red/green defect is present



Ishihara Plates

- Adequate lighting
 - Daylight or best artificial light
- Held at 75 cm
- Tilted to 90° to the line of vision
- Should answer within 3 seconds
- Do not touch colour plates
- Keep book closed when not in use




Hardy-Rand-Rittler plates

- No longer manufactured
- Pseudoisochromatic plates
 - Plates for detecting protans
 - Plates for detecting deutans
 - Plates for detecting tritans



Hardy-Rand-Rittler plates

- Neutral gray background
- Coloured circles, crosses, and triangles
 - Presented in higher and lower saturations
 - Detect the degree of colour vision deficiency



Hardy-Rand-Rittler plates

- Colour vision deficiency can be graded
 - Mild, moderate, severe
- Can differentiate Red-Green and Yellow-Blue defects
- Impossible to memorize (unlike Ishihara)
 - Identify not only what shape is seen but also in what quadrant it is in
 - Could rotate the plate to make even more combinations
- Detailed instructions included in the book



Farnsworth-Munsell D-15

- Series of 15 coloured caps
- Differentiate between
 - Protan
 - Deutan
 - Tritan
- Will not pick up subtle defects
 - Anomalous trichromat



Farnsworth-Munsell D-15

- Mix up caps and place them on table
- Have the patient arrange the caps in order for the anchor (both eyes open)
- 2 minute limit
 - Not enforced



Farnsworth-Munsell D-15

- Flip box over and record on the score sheet the order of the caps



Farnsworth-Munsell D-15

- Connect the numbers on the score sheet

TEST DICHOTOMIQUE de FARNSWORTH pour la Cécité des Couleurs - Série D-15

NOM: _____ Age: _____ Date: _____ Examineur: _____

Line: d'examen: _____

ANALYSE DICHOTOMIQUE

Type: Ane de Confusion

PROTANE (ROUGE - Bleu vert) MÉSOTIE

DEUTANE (VERT - rouge bleu-vert)

TRITANE (VIOLET - jaune vert) ÉCHIC


Test: Ordre donné par le sujet: 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15

Retest: Ordre donné par le sujet: _____

TEST RETEST

Farnsworth-Munsell D-15

- Normal vision = numbers in order (1-15)
- If all in order = finished



Farnsworth-Munsell D-15

- If any out of order = perform retest (score both)

TEST DICHOTOMIQUE de FARNSWORTH pour la Cécité des Couleurs - Série D-15

NOM: _____ Age: _____ Date: _____ Examineur: _____

Line: d'examen: _____

ANALYSE DICHOTOMIQUE

Type: Ane de Confusion

PROTANE (ROUGE - Bleu vert) MÉSOTIE

DEUTANE (VERT - rouge bleu-vert)

TRITANE (VIOLET - jaune vert) ÉCHIC

Test: Ordre donné par le sujet: 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15


Retest: Ordre donné par le sujet: _____

TEST RETEST

Farnsworth-Munsell D-15

- Defect corresponds with the line to which it is parallel


Protan



Farnsworth-Munsell D-15

- Defect corresponds with the line to which it is parallel


Deutan



Farnsworth-Munsell D-15

- Defect corresponds with the line to which it is parallel

Tritan



Farnsworth-Munsell 100-Hue

- Tests colour discrimination
- Separate people with normal colour vision into
 - Superior colour discrimination
 - Average colour discrimination
 - Low colour discrimination




Farnsworth-Munsell 100-Hue

- Four wooden cases
- 93 coloured caps
 - 85 removable
 - 2 caps in each box are fixed and repeated



Farnsworth-Munsell 100-Hue

- Perform each eye separately
- Natural daylight or equivalent



Farnsworth-Munsell 100-Hue

Procedure

- The object of the test is to arrange the caps in order according to colour. Transfer the caps in to the panel so they form a regular colour series between the 2 caps. It should take you about 2 minutes per panel. However, accuracy is more important than speed. Arrange them as best as you can...but do not dawdle.



Farnsworth-Munsell 100-Hue

Recording Data

- Record the order of the caps



-- -- -- 18 22 16 17 19 -- --
 13 14 15 16 17 18 19 20 21 22



Farnsworth-Munsell 100-Hue

Scoring

- Essentially counting errors
- Score for each cap is the sum of the difference between the number of the cap and the number of the caps adjacent to it

Farnsworth-Munsell 100-Hue

RECORDED ORDER 13 14 15 18 22 16 17 19






Farnsworth-Munsell 100-Hue

RECORDED ORDER 13 14 15 18 22 16 17 19

DIFFERENCE b/w CAPS

1 1 3 4 6 1 2

Farnsworth-Munsell 100-Hue



RECORDED ORDER 13 14 15 18 22 16 17 19

DIFFERENCE b/w CAPS

SCORE FOR CAP


2 4 7 10 7 3

Therefore, cap number 14 has an error score of 2 (1+1)
2 is the lowest possible score
Cap 22 has an error score of 10 (4+6)






Farnsworth-Munsell 100-Hue

- Inner circle corresponds to the number of caps
- First circle represents a score error of 2
- Mark each error score on the radial line
- Connect the points
- Total error = sum the errors on each radial line counting the inner circle as zero
- Re-test as required





Subject -
Reference -
Date of Birth - 01-28-2003
Date of test - Right Eye
Test Duration - 24 (mins)
Total error score (TES) = 208
Classification - Low discrimination



Farnsworth-Munsell 100-Hue

Results

- Superior colour discrimination
 - Total error score of 0 to 16
 - 0 to 4 transpositions
- Average colour discrimination
 - Total error score of 20 to 100
- Low colour discrimination
 - Total error score of greater than 100




Case Study

Case Study




Females - XX
Males - XY
X linked recessive

- Affected Male
- Affected Female
- Unaffected Male
- Unaffected Female
- Female Carrier

3rd generation offspring from youngest daughter

Male - 9 years old


3rd generation offspring from youngest daughter

Male - 9 years old

correctly identified 4 half numbers



- 55 instead of 57
- 17 instead of 15
- 21 instead of 71
- 27 instead of 97

IHCP 0/16






23 → Suggestive for deuteranopia

42

3rd generation offspring from youngest daughter

Male - 9 years old






3rd generation offspring from youngest daughter

Male - 9 years old

ICHP: deuteranopia

Anecdotal: normal / abnormal




3rd generation offspring from youngest daughter
Male - 9 years old

FARNSWORTH DICHO TOMOUS TEST for Color Blindness—Panel D-15
Name: ALECK Age: 9 Date: _____ File No.: _____
Department: _____ Test: _____

The DEMONSTRATION ANSWERS are:
Type: PROBLEM (RED-green) PASS
DEUTAN (GREEN-yellow) FAIL
TRITAN (VIOLET-purple-blue) FAIL

The Infant's Color: _____
The Infant's Color: _____
The Infant's Color: _____
The Infant's Color: _____

TEST: _____ RETEST: _____

© 1992 HANCOCK & COMPANY, INC. All rights reserved. No reproduction or translation without the written permission of Hancock & Company, Inc. 1000 North Washington Street, Boston, MA 02114

3rd generation offspring from youngest daughter
Male - 9 years old

© 1992 HANCOCK & COMPANY, INC. All rights reserved. No reproduction or translation without the written permission of Hancock & Company, Inc. 1000 North Washington Street, Boston, MA 02114

3rd generation offspring from youngest daughter
Male - 9 years old

FARNSWORTH DICHO TOMOUS TEST for Color Blindness—Panel D-15
Name: ALECK Age: 9 Date: _____ File No.: _____
Department: _____ Test: _____

The DEMONSTRATION ANSWERS are:
Type: PROBLEM (RED-green) PASS
DEUTAN (GREEN-yellow) FAIL
TRITAN (VIOLET-purple-blue) FAIL

The Infant's Color: _____
The Infant's Color: _____
The Infant's Color: _____
The Infant's Color: _____

TEST: _____ RETEST: _____

© 1992 HANCOCK & COMPANY, INC. All rights reserved. No reproduction or translation without the written permission of Hancock & Company, Inc. 1000 North Washington Street, Boston, MA 02114

3rd generation offspring from youngest daughter
Male - 9 years old

© 1992 HANCOCK & COMPANY, INC. All rights reserved. No reproduction or translation without the written permission of Hancock & Company, Inc. 1000 North Washington Street, Boston, MA 02114

3rd generation offspring from youngest daughter
Male - 9 years old

ICHP: deuteranopia
Anecdotal: normal / abnormal
D-15: normal

FARNSWORTH DICHO TOMOUS TEST for Color Blindness—Panel D-15
Name: ALECK Age: 9 Date: _____ File No.: _____
Department: _____ Test: _____

The DEMONSTRATION ANSWERS are:
Type: PROBLEM (RED-green) PASS
DEUTAN (GREEN-yellow) FAIL
TRITAN (VIOLET-purple-blue) FAIL

The Infant's Color: _____
The Infant's Color: _____
The Infant's Color: _____
The Infant's Color: _____

TEST: _____ RETEST: _____

© 1992 HANCOCK & COMPANY, INC. All rights reserved. No reproduction or translation without the written permission of Hancock & Company, Inc. 1000 North Washington Street, Boston, MA 02114

3rd generation offspring from youngest daughter
Male - 9 years old

H-R-R Plates

Could not identify plates 3-6 = R-G defect

More checks under deutan = deutan

No further errors = mild

METHOD OF RECORDING SUBJECT'S RESPONSES

FOUR PLATES, DO NOT READ.


SCREENING SERIES	DIAGNOSTIC SERIES	DEUTAN	TRITAN	SCREENING SERIES
7 []	7 []	[]	[]	7 []
8 []	8 []	[]	[]	8 []
9 []	9 []	[]	[]	9 []
10 []	10 []	[]	[]	10 []
11 []	11 []	[]	[]	11 []
12 []	12 []	[]	[]	12 []
13 []	13 []	[]	[]	13 []
14 []	14 []	[]	[]	14 []
15 []	15 []	[]	[]	15 []
16 []	16 []	[]	[]	16 []
TOTAL	TOTAL	Deutan	Deutan	Deutan
17 []	17 []	[]	[]	17 []
18 []	18 []	[]	[]	18 []
19 []	19 []	[]	[]	19 []
20 []	20 []	[]	[]	20 []
TOTAL	TOTAL	Deutan	Deutan	Deutan

© 1992 HANCOCK & COMPANY, INC. All rights reserved. No reproduction or translation without the written permission of Hancock & Company, Inc. 1000 North Washington Street, Boston, MA 02114

3rd generation offspring from youngest daughter

Male - 9 years old

ICHP: deuteranopia
 Anecdotal: normal / abnormal
 D-15: normal
 H-R-R: mild deutan





3rd generation offspring from youngest daughter

Male - 9 years old

100 Hue



Low colour discrimination


3rd generation offspring from youngest daughter

Male - 9 years old

ICHP: deuteranopia
 Anecdotal: normal / abnormal
 D-15: normal
 H-R-R: mild deutan
 100 Hue : low colour discrimination

"This is REO the cat. He is blue"



Colour Vision Testing

Craig Simms
 BSc, COMT, CDOS, ROUB
 Chief Ophthalmic Medical Technologist
 Program Coordinator, KOTC
 Kingston, Ontario

JCAHPO ACE Meeting 2015
 Las Vegas

