



Using high throughput microtissue culture to study the difference in prostate cancer cell behavior and drug response in 2D and 3D co-cultures

Eman Mosaad^{1,2,3}, Karen Chambers⁴, Kathryn Futrega^{1,2}, Judith Clements² and Michael Robert Doran^{1,2,5,6*}

Abstract

Background: There is increasing appreciation that non-cancer cells within the tumour microenvironment influence cancer progression and anti-cancer drug efficacy. For metastatic prostate cancer (PCa), the bone marrow microenvironment influences metastasis, drug response, and possibly drug resistance.

Methods: Using a novel microwell platform, the *Microwell-mesh*, we manufactured hundreds of 3D co-culture microtissues formed from PCa cells and bone marrow stromal cells. We used luciferase-expressing C42B PCa cells to enable quantification of the number of PCa cells in complex microtissue co-cultures. This strategy enabled us to quantify specific PCa cell growth and death in response to drug treatment, in different co-culture conditions. In parallel, we used Transwell migration assays to characterize PCa cell migration towards different 2D and 3D stromal cell populations.

Results: Our results reveal that PCa cell migration varied depending on the relative aggressiveness of the PCa cell lines, the stromal cell composition, and stromal cell 2D or 3D geometry. We found that C42B cell sensitivity to Docetaxel varied depending on culture geometry, and the presence or absence of different stromal cell populations. By contrast, the C42B cell response to Abiraterone Acetate was dependent on geometry, but not on the presence or absence of stromal cells.

Conclusion: In summary, stromal cell composition and geometry influences PCa cell migration, growth and drug response. The Microwell-mesh and microtissues are powerful tools to study these complex 3D interactions.

Keywords: Prostate cancer, Bone marrow stromal cells, Co-culture, Microwell-mesh platform, 3D culture

Background

Despite significant improvements in the survival of prostate cancer (PCa) patients with localized disease, survival drops significantly if the cancer has metastasized to a distal site [1]. Approximately 90% of metastatic prostate cancer patients suffer bone metastasis [2, 3], making modeling of PCa cell behavior within the bone tissue micro-environment especially relevant.

¹Stem Cell Therapies Laboratory, Queensland University of Technology (QUT), Translational Research Institute (TRI), 37 Kent Street, Brisbane, QLD, Australia ²Australian Prostate Cancer Research Centre – Queensland (APCRC-Q), Translational Research Institute (TRI), Brisbane, Australia Full list of author information is available at the end of the article



© The Author(s). 2018 **Open Access** This article is distributed under the terms of the Creative Commons Attribution 4.0 International License (http://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made. The Creative Commons Public Domain Dedication waiver (http://creativecommons.org/publicdomain/zero/1.0/) applies to the data made available in this article, unless otherwise stated.

^{*} Correspondence: michael.doran@qut.edu.au

Within the bone, there is evidence that the first site of metastasis is the hematopoietic stem cell (HSC) niche [1]. Key HSC niche microenvironmental cell populations include bone marrow mesenchymal stromal cells (BMSC), osteoblasts and adipocytes [4–6]. These cell populations are all thought to influence PCa metastasis and disease progression [7–9]. Dissecting the influence played by each stromal cell population in vivo is challenging, and this is an area where in vitro model experimentation may offer an advantage over more complex animal models. An on-going challenge is the establishment of an in vitro model that mimics the in vivo microenvironment sufficiently to yield clinically relevant results or insights. The most common tissue culture models are 2D cell monolayers grown

on tissue culture polystyrene. Monolayer cultures do not facilitate tissue-like cell-cell interactions [10], and cancer cells cultured in 2D monolayers tend to be hypersensitive to anti-cancer drugs [11]. This has motivated a surge in the development of 3D cancer models that are meant to better recapitulate 3D cellular organization and complex tissue microenvironments [12, 13].

Despite the potential advantages of 3D culture models, their use in PCa drug screening remains limited. Traditional 2D tissue culture plates are inexpensive, the majority of imaging systems/protocols are designed to be compatible with 2D culture plates, and a range of automated fluidic systems are compatible with 2D culture systems. These features have not yet been efficiently integrated into 3D culture systems. For example, hydrogel matrix-based 3D cultures can be costly, they commonly suffer from significant 3D tissue size heterogeneity, and harvest from the gel is necessary for many forms of analysis [14]. Our team previously introduced the Microwell-mesh as a high throughput platform suitable for 3D tissue culture [15, 16]. The Microwell-mesh uses a microwell platform to facilitate the manufacture of hundreds of uniform multicellular 3D microtissues. It differs from previous microwell platforms in that it has a nylon mesh fixed over the microwells, and this enables retention of individual microtissues within discrete microwells even during repeat full medium exchanges. The capacity to exchange medium repeatedly is especially useful in drug testing applications. Additionally, the mesh enables establishment of microtissues from one cell type, and then the addition of a second cell type at a later time point. Because specific numbers of cells can be deposited and retained in discrete microwells, this allows the assembly of co-culture microtissues having specific co-culture cellular composition. This design makes the Microwell-mesh platform ideal for use in the simultaneous manufacture, characterization and study of the drug response of hundreds of microtissues in a high throughput manner.

An additional complexity associated with designing co-culture drug assays is that it is challenging, and potentially expensive, to specifically quantify the number of cancer cells without the co-culture population confounding this measurement. For example, simple Alamar blue metabolic readouts would include both metabolic contributions from the cancer and the stromal co-culture cell population(s), making specific cancer cell responses challenging to delineate. To overcome this barrier, we mimicked McMillin and colleagues who used a luciferase reporter system to enable the indirect estimation of cancer cell numbers in complex co-cultures via bioluminescence [17]. In our studies, the PCa cells were transduced to express a luciferase reporter, allowing us to indirectly quantify PCa cell number in complex co-cultures with stromal cell populations that did not express luciferase.

Herein, PCa cell migration and proliferation in response to bone marrow stromal cell populations cultured in 2D and 3D was contrasted. We used the Microwell-mesh system to form microtissues containing both PCa and bone marrow stromal cells, and used the luciferase reporter system to enable indirect quantification of PCa growth as well as death in response to anti-cancer drugs in complex co-cultures. The response of PCa cells to Docetaxel and Abiraterone Acetate in 2D or 3D, and in the presence or absence of stromal cells was characterized.

Methods

PCa cell line culture

PCa cell lines used were PC3 (purchased from ATCC° Number: CRL-1435), C42B (derived and generously shared by Dr. Chung [18]), and LNCaP cells (purchased from ATCC° Number: CRL-1740). Cell lines were authenticated at the Genomic Research Centre (GRC; Brisbane, Australia) using Short Tandem Repeat (STR) analysis. STR profiles of the cell lines were compared to the ATCC STR Database to verify cell line identity; and all cell lines showed \geq 80% match to the corresponding reference STR profile. C42B were mapped back to the LNCaP STR profile, as C42B were derived from LNCaP [18]. Cells were cultured in low glucose Dulbecco's modified Eagle's medium (DMEM-LG; Thermo Fisher) supplemented with 10% fetal bovine serum (FBS; Thermo Fisher) and 1% penicillin/streptomycin (P/S; Thermo Fisher). For some assays, FBS was replaced with 10% charcoal stripped fetal bovine serum (CSS; Thermo Fisher) to mimic androgen deprivation conditions. Cells were grown in a cell culture incubator at 37 °C and 5% CO2 and 2% O2. All cells were passaged when monolayers reached $\sim 80\%$ confluency using 0.25% Trypsin/EDTA (Thermo Fisher).

Human BMSC isolation, culture, characterization, and differentiation

Human bone marrow aspirates were collected at the Mater Hospital (Brisbane, Australia) from two fully informed and consenting healthy male volunteer donors. In accordance with the Australian National Health and Medical Research Council's Statement on Ethical Conduct in Research Involving Humans, ethical approval was granted through the Mater Health Services Human Research Ethics Committee and Queensland University of Technology Ethics Committee (number: 1000000938). Aspirates were collected from the iliac crest of volunteers. Mononuclear cell isolation was achieved by density gradient centrifugation, using Ficoll-Paque Plus (GE Healthcare), as previously described [19]. Bone marrow samples were diluted 1:2 with phosphate buffered saline (PBS; Thermo Fisher) containing 2 mM EDTA (Ambion). Then the diluted sample was carefully overlayed on the Ficoll-Paque plus layer and centrifuged for 30 min at 400×g. The

mononuclear cells collected from the interface were then washed, resuspended in DMEM-LG supplemented with 10% FBS, and 1% P/S. Cells were then cultured overnight in a humidified incubator containing 5% CO₂ with 20% O₂ atmosphere at 37°C. Tissue culture plastic-adherent cells were enriched by removing the medium containing non-adherent cells, and fresh culture medium added to each flask. Subsequent BMSC expansion was performed in 5% CO₂ and 2% O₂ atmosphere at 37°C. Cells were passaged when the monolayer reached 80% confluency using 0.25% Trypsin/EDTA. All experiments were performed using BMSC between passage 2 and 5.

The isolated cells were characterized for the expression of BMSC surface antigens; CD44, CD90, CD105, CD73, CD146, CD45, CD34 and HLADR; and mesodermal trilineage differentiation capacity and confirmed to be in accordance with the standard criteria of multipotent mesenchymal stromal cells reported previously by Dominici et al. [20].

Osteogenic and adipogenic differentiation were induced by culturing 60×10^3 and 40×10^3 cells/cm² in osteogenic or adipogenic induction medium for 14 days, respectively. Both induction media consisted of high glucose DMEM media (DMEM-HG) containing 100 µM sodium pyruvate, 1X GlutaMax, 10% FBS and 1% P/S (all from Thermo Fisher). Additionally, osteogenic medium contained 50 µM L-ascorbic acid 2-phosphate, 100 nM dexamethasone (Sigma-Aldrich) and 10 mM β-glycerol phosphate, while adipogenic medium contained 10 µg/mL insulin, 200 µM indomethacin, 1 µM dexamethasone (Sigma-Aldrich) and 500 µM 3-isobutyl-1-methyl xanthine (all from Sigma-Aldrich). Culture medium was replaced with fresh media twice per week.

Generation of C42B cell line expressing luciferase-GFP (Luc-GFP)


Firefly luciferase-expressing C42B cells were generated using fresh lentiviral particles produced in-house. Luciferase-GFP (Luc-GFP) insertion constructs contained Bioluminescence Imaging Vectors (BLIV, System Biosciences) with MSCV (MSCV-Luc-GFP) promoters (Additional file 1: Figure S1). Plasmid production was achieved by using Stbl3 chemically competent *E.coli* (Thermo Fisher) as per the manufacturer's instructions. This was followed by a purification step using the NucleoBond Xtra EF plasmid purification kit (Midi EF, Macherey-Nagel) to obtain endotoxin-free plasmid DNA. Plasmid packaging was then performed using TGEN packaging plasmid mix with the transfection reagent, Lipofectamine 2000 (Thermo Fisher). The lentiviral particles were produced by 293FT cells (Thermo Fisher) following the manufacturer's instructions. Viral particlecontaining media was then placed onto cancer cells, with the addition of 8 µg/mL polybrene (Sigma-Aldrich) to enhance transduction efficiency. Positively transduced (Luc-GFP) cells were enriched using two rounds of fluorescence-activated cell sorting (FACS; MoFlo Astrios, Beckman Coulter). This yielded a stable population of C42B cells that expressed Luc-GFP driven by a MSCV promoter. We validated the stability of luciferase gene expression in monolayer and Transwell co-culture conditions using quantitative real time-polymerase chain reaction (qRT-PCR) [15] (Additional file 1: Figure S2) and appropriate PCR primer sets (Additional file 1: Table S1).

3D culture system design and fabrication

An in-house fabricated microwell platform was fabricated from polydimethylsiloxane (PDMS; Slygard). PDMS microwell arrays were fabricated as described previously [11, 15]. Briefly, liquid PDMS (1:10 curing agent to polymer ratio) was permitted to cure over a patterned polystyrene mold having the negative of the microwell pattern for 1 h at 80°C. A sheet of PDMS with the microwell array pattern cast into it (each microwell had dimensions of $800 \times 800 \ \mu m$ square and a depth of 500 µm) was produced and peeled from the molds. Discs (1 cm² or 2 cm²) were punched from the PDMS sheets and then glued into culture plates with silicone glue (Selleys). For drug testing experiments, Microwell-mesh inserts were made by fixing a nylon mesh with 36 μ m x 36 μ m pore openings (Amazon) to the top of the microwells using silicone glue. Once the glue had cured, excess mesh was trimmed from the disc inserts using scissors. Inserts were then anchored into individual wells in 24or 48-well plates by placing a small amount of silicone glue at the bottom of the well, and the insert pressed into the well. Plates with microwell inserts were submerged in 70% ethanol for 1 h for sterilisation, followed by rinsing of each culture well 4 times with PBS (Thermo Fisher). To prevent cell adhesion to the PDMS during culture, the PDMS microwell inserts were soaked in a sterile solution of 5% Pluronic-F127 (Sigma-Aldrich) in PBS for 10 min [21], and then rinsed 3 times with PBS before cells were seeded.

Assembly of microtissues

In this study, we formed microtissues assembled from PCa, BMSC (non-induced), osteoblasts or adipocytes alone, or combination co-cultures of PCa with BMSC, osteoblasts or adipocytes. The Microwell-mesh platform was used to study PCa proliferation and drug response in direct co-cultures where multiple medium and drug exchanges were required. Fig. 1 schematically illustrates how the Microwell-mesh differs from traditional open top microwell platforms, and how centrifugation can be used to evenly distribute the seeded cell suspension into the array of microwells. Each insert had approximately 150 microwells/cm², (equivalent to \sim 150 microwells per well in a 48-well plate). Seeding a different number of cells in suspension over the microwells could control the number of cells per microwell or per microtissue. Following seeding of the cell suspension, plates were centrifuged at



400 × g for 5 min to aggregate the cells uniformly at the bottom of each microwell. The aggregation of cells into microwells was visually confirmed using a microscope (Olympus CKX14), and images captured using a digital camera (Olympus DP26) and software (Olympus cellSens Entry). Plates were then transferred to a cell culture incubator maintained at 37°C and 5% CO₂.

PCa cell Transwell migration assay


We were also interested in determining if PCa cell migration towards stromal cells differed depending on the geometry of the stromal cells. In 3D microtissue co-cultures, PCa cells localized to the outside of the microtissue, but this did not provide insight into how different stromal cells might influence cell migration. To overcome this obstacle we developed a modified Transwell assay. Here, we either cultured the stromal cells as 2D monolayers or as 3D microtissues in open top microwell inserts. To quantify PCa cell migration, PCa cells were placed into Transwell inserts (pore size of 8 µm, Merck Millipore) and positioned either on top of 2D stromal cell monolayers or on top of 3D stromal cell microtissues (see Fig. 2). BMSC were seeded in 24-well tissue culture plates, and cells cultured in osteogenic or adipogenic medium for 14 days or in maintenance medium for 24 h. For 2D monolayers, 10-, 20- and 60×10^3 cells/cm² were seeded and cultured in the corresponding culture media. For 3D microtissues, 600 cells/microtissue were seeded in the microwell inserts anchored in the 24-well tissue culture plate as described above. Transwell inserts were seeded with 36×10^3 PC3, C42B or LNCaP cells and permitted to incubate for 24 h. Inserts were then placed on top of either the 2D or 3D stromal cell populations and incubated for 18 h. At the end of the incubation period, Transwell inserts containing PCa cells were washed and moved to a new tissue culture plate. Adherent cells attached to the top surface of the Transwell insert were removed using cotton buds, while cells that had migrated to the bottom surfaces of the Transwell inserts were fixed using ice cold methanol for 15 min. Fixed Transwell inserts were immersed in crystal violet stain (0.5%, diluted in H₂O) for 15 min. Transwell inserts were washed in running tap water to remove excess stain. Crystal violet stain was extracted from cells into 500 μ l of 10% acetic acid for 10 min. The optical density (OD) of the extract was measured at 595 nm (Multiskan Go microplate spectrophotometer, Thermo Fisher). Optical density of extracts from cell-free Transwell inserts functioned as controls for empty wells. For each cell line, parallel Transwell inserts containing PCa cells not exposed to the stromal co-culture conditions functioned as baseline migration controls.

Confocal imaging of 3D PCa-BMSC co-culture

PCa-BMSC microtissues were established using the microwell culture system. Single cell suspensions of C42B and BMSC were stained with green molecular probe (CellTrace green CFSE) and red molecular probe (CellTracker Red CMTPX; both from Thermo Fisher), respectively. A cell suspension combining the two cell types in a 1:1 ratio was generated, and seeded into 48-well tissue culture plates with microwell inserts to obtain microtissues each containing 600 cells (300 C42B cells and 300 BMSC). Following 24 h incubation at 37°C, 5% CO₂, microtissues were collected and imaged using a Zeiss 510 Meta confocal microscope to characterize 3D cellular organization.

Bioluminescence assay

In vitro bioluminescence of Luciferase-tagged PCa cells was used as an indirect method to estimate viable cancer cells in mono- and co-cultures. For the luciferase activity assay, D-luciferin (Promega) was added to the culture medium at a final concentration of 15 μ g/mL, then incubated at 37 °C for 15 min and the bioluminescence acquired using a PHERAstar FS plate reader (BMG LABTECH).



Cell proliferation and drug testing in direct co-culture system

PCa cell proliferation and responses to anti-cancer drugs were tested in both 2D and 3D co-cultures. In 48-well tissue culture plates, co-cultures were established in 2D monolayers or as 3D microtissue cultures. Two weeks prior to establishing co-cultures, BMSC were assembled into 3D microtissues of 300 cells/microtissue or as 2D monolayers of 60×10^3 and 100×10^3 cells/cm² to permit differentiation to osteogenic and adipogenic lineages, respectively. At day zero, C42B Luc-GFP cells were added as a single cell suspension on top of stromal cell (BMSC, osteoblasts or adipocytes) monolayers or stromal cell microtissues in the Microwell-mesh. C42B cells were seeded at either 10×10^3 cells/cm² on top of stromal cell monolayers, or 300 cells per microtissue on top of established stromal microtissues. For cell proliferation experiments, PCa cells were permitted to grow for 24 and 48 h in 2D and 3D co-cultures as mono- or co-cultures then the bioluminescence was measured as described above. Data is presented as relative bioluminescence (RLU) relative to luciferase-tagged PCa cells in mono-cultures.

Docetaxel and Abiraterone Acetate were used in the drug testing studies. Docetaxel and Abiraterone Acetate were dissolved in dimethyl sulfoxide (DMSO; all from Sigma-Aldrich), and then aliquoted and stored at -80° C. On the day of treatment, an aliquot was thawed and diluted in culture medium to the specified concentrations.

PCa cells were permitted to adhere or aggregate into spheriods for 24 h in co-cultures. For Docetaxel treatments, all cultures were treated with the indicated concentrations starting one day after the initiation of the co-cultures, with drug treatment exposure being continuous for the next 48 h. For Abiraterone Acetate treatments, all cultures were first depleted of androgens for 48 h by replacing the FBSsupplemented culture medium with CSS-supplemented medium. Cultures were then exposed to the specified concentrations of Abiraterone Acetate for 48 h.

At the end of the drug treatment period, epifluorescence and phase contrast microscopy images were captured and bioluminescence signals from each culture were measured, as described above. Bioluminescence data is presented as a percentage of the relative bioluminescence units (RLU) compared to vehicle-treated cultures.

Statistical analysis

Results represent two independent experiments using two BMSC donors. Each of the replicate experiments included four biological replicate cultures (n = 4), unless otherwise indicated. Error bars represent one standard deviation. Statistical significance of data was evaluated using two-way analysis of variance (ANOVA), using Prism software, Version 6.0 (GraphPad). *P*-values for each comparison are represented by asterisks as indicated in figure captions.

Results

Indirect Transwell co-culture of PCa cells with 2D and 3D bone marrow stromal cells

Using the Transwell assay, the migration of PCa cells towards BMSC, osteoblasts and adipocytes cultured in 2D monolayers or 3D microtissues was assessed following 18 h of co-culture (Fig. 2a). LNCaP, C42B and PC3 cells were used to represent or model different stages of PCa disease aggressiveness.

Of the 2D cultures, BMSC monolayers induced the greatest migration rates in all PCa cell lines tested. By contrast, the influence of osteoblasts and adipocytes on PCa migration was PCa cell line dependent. For the less aggressive cell lines, C42B and LNCaP, both osteoblasts and adipocytes had minimal influence on PCa cell migration rates. The highly aggressive bone metastatic PC3 cells demonstrated a significantly elevated migration rate towards osteoblasts and adipocytes cultured in 2D monolayers (Fig. 2b).

Unlike 2D BMSC cultures, which increased the migration of all PCa cells tested, 3D BMSC microtissues only increased the migration of PC3 cells. Indirect co-culture with 3D adipocyte microtissues decreased PC3 cell migration, and had no measurable effect on C42B or LNCaP cell migration rates. Similarly, 3D osteoblast microtissues did not increase the migration rate of any of the PCa cell lines tested (Fig. 2b).

Spatial organization of C42B cells and stromal cells in 3D co-cultures

To characterize the spatial organization of 3D co-culture microtissues, we first labeled each cell type with differently colored fluorescent probes to enable the two cell types to be distinguished from each other. C42B cells were labeled with a green probe, while BMSC were labeled with a red probe. Figure 3a shows microtissue co-cultures formed from C42B and BMSC. Confocal images of 3D co-culture microtissues demonstrated a consistent and structured organization of the two cell types across the diameter of the microtissues. BMSC consistently localized within the core of the microtissue, whereas C42B cells were localized in the outer layer of the microtissue after 24 h of co-culture (Fig. 3b).

C42B cell proliferation in co-cultures with bone marrow stromal cells

Conventional methods of quantifying cell proliferation, such as metabolic activity or cell viability assays, do not allow for quantification of the cell growth of a single cell population within a mixed cell population co-culture. To study C42B cell proliferation in co-cultures, we labeled the PCa cell population with a luciferase reporter system. Relative bioluminescence signal from the PCa cell populations functioned to provide an indirect estimate of the number of viable PCa cells in the different co-cultures.


To study PCa cell proliferation, we used C42B cells stably expressing Luc-GFP. These luciferase-expressing C42B cells were cultured in 2D monolayers or 3D microtissue cultures for 24 or 48 h, as either mono-cultures of PCa cells (control) or co-cultures of PCa cells with stromal cells (BMSC, osteoblasts or adipocytes). The bioluminescence assay was then performed to estimate the number of C42B cells in each culture condition at each time point.

In 2D cultures, the bioluminescence values indicated a significant increase in C42B cell number in all co-culture conditions after 24 h, relative to mono-culture controls. After 48 h of culture, the effect of stromal cells on C42B cell proliferation was less pronounced. However, the overall bioluminescence after 48 h was significantly greater than after 24 h for all cultures (Fig. 4a), indicating continual cell proliferation in all culture conditions.

In 3D cultures, co-culture with adipocytes enhanced C42B cell proliferation after 24 and 48 h of culture, while co-culture with osteoblasts did not influence C42B cell proliferation rate. Despite the slight decrease in bioluminescence of C42B cells co-cultured with BMSC at 24 h, the bioluminescence tended to increase (non-significant increase) after 48 h of co-culture (Fig. 4b). Similar to 2D cultures, an overall increase in the bioluminescence of C42B cells in 3D cultures was observed at 48 h, relative to 24 h-cultures.

C42B cell drug response in co-cultures with bone marrow stromal cells


Next, we evaluated the response of PCa cells in 2D monolayers and 3D microtissues to Docetaxel and Abiraterone Acetate, two drugs used to treat advanced




PCa. Luciferase-tagged C42B cells were used in these experiments, and the bioluminescence measurements provided an indirect estimate of the viable cell number in the cultures after 48 h of drug treatments. Three replicate experiments were also performed using 2D and 3D co-culture of osteoblasts and Luciferase-expressing C42B cells. Over the total co-culture period, C42B cell viability fell dramatically, even in control co-cultures with no drug. These data suggest that long-term stability of 2D and 3D co-culture is stromal cell type dependent. In the subsequent analysis below, we focused on results derived from 2D and 3D cultures of C42B cells alone, or in co-cultured with BMSC or adipocytes.

After 24 h of establishing the mono- and co-cultures, Docetaxel treatment was performed for 48 h. In 2D cultures, there was significantly greater bioluminescence signal from the C42B cells co-cultured with BMSC or adipocytes in all Docetaxel concentrations (0.01–10 nM) relative to the bioluminescence in mono-cultures at the same drug concentration (Fig. 5a). Unexpectedly, 3D BMSC co-cultures showed a significant increase in bioluminescence. By contrast, adipocyte co-cultures behaved similarly to corresponding mono-cultures (Fig. 5b). In general, 3D cultures demonstrated reduced sensitivity to Docetaxel in both mono- and co-cultures with BMSC or adipocytes (Fig. 5).

For anti-androgen treatment, C42B cells were cultured in androgen deprived setting (CSS-supplemented culture media), and then treated with Abiraterone Acetate. Abiraterone Acetate is a first-in-class inhibitor of the CYP17A enzyme to prevent the biosynthesis of androgens intracellularly from their steroidal precursor [22]. The bioluminescence assay was used to assess PCa cell response in 2D and 3D mono- and co-cultures. Fig. 6 shows the bioluminescence measurements as a percentage of the corresponding vehicle control culture.






In 2D cultures, co-cultures with BMSC and adipocytes demonstrated no significant change in the anti-androgen treatment response compared to the mono-cultures of C42B cells (Fig. 6a). Similarly, 3D co-cultures did not result in change in the bioluminescence, except with BMSC co-cultures treated with 10 μ M Abiraterone Acetate, which resulted in a decrease in bioluminescence (Fig. 6b). Generally, 3D mono- and co-cultures were less sensitive to increasing concentrations of Abiraterone Acetate, relative to their corresponding 2D monolayer controls.


Discussion

Metastatic, and particularly castrate-resistant prostate cancer (CRPC), remain challenging to treat [23]. It is thought that the bone marrow microenvironment plays a pivotal role in promoting bone metastasis, possibly facilitating the transition to CRPC forms, and impacting on PCa cell drug response [24–27]. A barrier to understanding these interactions, in both drug development

and testing, is the lack of in vitro models that adequately mimic aspects of the bone marrow microenvironment in a practical and high throughput manner.

Our team previously described the development of a high throughput 3D culture platform we termed the *Microwell-mesh* [15]. This platform uses a microwell insert to facilitate the manufacture of hundreds of uniform 3D multicellular microtissues. It differs from previous microwell platforms in that it has a nylon mesh fixed over the microwells, and this enables retention of individual microtissues within discrete microwells even during repeat full medium exchanges. This design is unique, and especially well suited to the assembly of 3D cultures which mimic aspects of the bone marrow microenvironment, and offers the opportunity to perform complex cultures that involve the differentiation of BMSC into different bone-like tissues, subsequent seeding of cultures with PCa cells, and the multiple medium exchanges required to study the interaction of cells and different drugs in these complex cultures.





Using the Microwell-mesh to perform 3D cultures, and traditional 2D culture controls, we evaluated PCa cell migration and proliferation in response to bone marrow stromal cell populations, as well as PCa cell response to Docetaxel and Abiraterone Acetate. The goal of this study was to better understand the difference 2D and 3D stromal cell populations might have on PCa culture outcomes, and to describe models that could advance the field's capacity to study these differences.

To study the impact of bone marrow stromal cells on the migration potential of PCa cells, we used a modified Transwell assay to quantify the migration of three different PCa cell lines towards different populations of bone marrow stromal cells (see Fig. 2). PCa cell migration rates varied depending on the aggressiveness of the PCa cell lines tested. In cell lines derived from less aggressive disease (LNCaP), relative to aggressive disease (C42B and PC3), there was a corresponding reduction in the rate of cell migration towards the bone marrow stromal cells cultured in 2D monolayers. PC3 cells, which model aggressive disease, demonstrated increased migration rates towards 2D monolayers of undifferentiated BMSC, osteoblasts and adipocytes. By contrast, PC3 cells demonstrated an increased rate of migration towards 3D osteoblasts and a reduced rate of migration towards undifferentiated BMSC or adipocytes, relative to controls. This data highlights the difference in PCa cell response depending on the PCa cell phenotype, the bone marrow stromal cell phenotype, and depending on the 2D or 3D organization of the bone marrow stromal cells. Appreciating that these factors influence outcome is an important first step that can inform our understanding and future experimental design. However, it is equally imporant to appreciate that outcomes can be influenced by the selected assay, and that not all in vitro and in vivo assays will necessarily yield the same outcome. Transwell cultures enable quantification of the influence secreted factors have on PCa cell migration, but do not necessarily provide insight into how stromal cell-specific matrix or bound factors may directly influence PCa cell behavior. Thus, Transwell assay outcomes provide only part of the necessary insight.

Next, we investigated how 2D or 3D culture of different bone marrow stromal cell populations impacted on C42B cell proliferation. C42B cell proliferation was greater when these cells were seeded on 2D monolayers of undifferentiated BMSC, adipocytes or osteoblasts (see Fig. 4a). This result is consistent with the general view that stromal cells can play a supportive role in co-cultures, and especially those that mimic aspects of the support environment found in the bone marrow niche [28, 29]. This result is also consistent with a previous report indicating that BMSC-conditioned media supports PCa cell proliferation [30]. In 3D co-cultures, only adipocytes were found to drive significant increases in C42B cell proliferation (see Fig. 4b). This substantial difference in 2D and 3D culture outcomes is interesting, as it indicates that geometry can significantly impact co-culture outcomes. Future studies might compare the secretion profiles of BMSC, adipocytes or osteoblasts in 2D and 3D, with the hypothesis that culture geometry significantly influences what factors are produced by the stromal cell populations. There are already a number of studies that suggest the secretome of BMSC is more supportive when BMSC are assembled into spheroids [31, 32]. Characterizing precise changes in the gene expression or secretion profile of the mesenchymal and PCa cells assembled into spheroids would require digesting the co-culture microtissues into single cell suspensions, followed by cell sorting and then gene or protein analysis. The considerable processing steps and time would likely confound the results. Thus, within this manuscript we focused on platform development and phenomenological characterization of PCa growth and drug response as influenced by the presence or absence of different mesenchymal cell populations. Equally valuable, would be to compare how 2D and 3D co-cultures influence the proliferation of primary PCa cells. Primary PCa cells are particularly challenging to culture in vitro [33], but their response to co-culture might be more meaningful than the response from an adapted cell line. We see these important investigations as beyond the scope of this manuscript, but obvious opportunities that could be explored using the Microwell-mesh as a tool to facilitate these important next steps.

In our studies, we found that the drug response of C42B cells differed in 2D and 3D co-cultures, and that response varied depending on the stromal cell population used in the co-culture (see Figs. 5 and 6). Collectively, the presence of BMSC or adipocytes in 2D or 3D co-culture reduced C42B cell sensitivity to Docetaxel, a drug commonly used to treat metastatic disease. Other groups have reported similar observations [17, 34, 35], suggesting that bone marrow stromal cells likely do influence PCa cell drug sensitivity. In contrast to tests conducted with Docetaxel, the drug response of C42B cells to Abiraterone Acetate did not appear to be influenced by the presence or absence of bone marrow stromal cells. However, the organization of C42B cells into 3D cultures did reduce these cells sensitivity to Abiraterone Acetate, relative to 2D cultures. This outcome suggests that relative proliferation rates, which are generally reduced in 3D cultures [36, 37], may play a greater role than the presence or absence of stromal cells in influencing the impact of anti-androgen treatment.

Conclusions

Overall, our results indicate that C42B cell behaviour can vary depending on the phenotype and geometry of bone marrow stromal cells included in co-culture. Through this work we have described the development of a 3D co-culture platform, the Microwell-mesh, that enables the assembly of 3D bone stromal cell microtissues, the subsequent introduction of PCa cells, and then evaluation of PCa cell proliferation or drug response. Using this novel 3D culture platform, we show that PCa cell response to drugs varies considerably in 2D and in 3D, and culture outcomes are also stromal cell-type dependent. This 3D culture tool provides more complex in vitro analysis, and will hopefully lead to the more efficient identification of improved PCa treatment strategies.

Additional file

Additional file 1: Figure S1. Restriction map of MSCV-Luc-GFP plasmid. Figure S2. Luciferase gene expression in C42B-MSCV cell lines. Table S1. Primers and annealing temperatures used for gRT-PCR. (PDF 609 kb)

Abbreviations

2D: Two-dimensional; 3D: Three-dimensional; BMSC: Bone marrow mesenchymal stromal cells; CRPC: Castrate-resistant prostate cancer; CS-BLI: Cell-specific bioluminescence assay; CSS: Charcoal-stripped serum; HSC: Haematopoietic stem cell; Luc-GFP: Luciferase-GFP; PCa: Prostate cancer

Acknowledgments

The authors would like to thank the National Health and Medical Research Council (NHMRC) of Australia and the Prostate Cancer Foundation of Australia (PCFA) for supporting this research. The authors would like to thank Mr. Eric Franklin for his assistance with BMSC isolation. The authors would like to thank the Egyptian Ministry of Higher Education for supporting Eman Mosaad with a PhD scholarship.

Funding

M.R.D. is funded by the National Health and Medical Research Council (NHMRC) of Australia CDF Fellowship. J.C. was funded by an NHMRC Principal Resarch Fellowship. The project was funded by the Prostate Cancer Foundation of Australia and the NHMRC. The funding bodies did not influence the design of the study and collection, analysis, and interpretation of data and in writing the manuscript.

Availability of data and materials

The datasets used and/or analysed during the current study are available from the corresponding author on reasonable request.

Authors' contributions

All Authors have read and approved this manuscript. EM designed research, performed research, analysed data and wrote the paper. KC and JC designed research, analysed data and wrote the paper. KF designed research, performed research, analysed data and wrote the paper. MRD designed research, analysed data and wrote the paper.

Ethics approval and consent to participate

None of the prostate cancer cell lines used in this study required institutional human ethics approval. Only the collection of bone marrow aspirate donations required human ethics approval. Ethical approval was granted through the Mater Health Services Human Research Ethics Committee and Queensland University of Technology Ethics Committee (number: 100000938), in accordance with the Australian National Health and Medical Research Council's Statement on Ethical Conduct in Research Involving Humans. Volunteer bone marrow donors were informed, and provided written consent prior to participating.

Competing interests

The authors declare that they have no competing interests.

Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Author details

¹Stem Cell Therapies Laboratory, Queensland University of Technology (QUT), Translational Research Institute (TRI), 37 Kent Street, Brisbane, QLD, Australia. ²Australian Prostate Cancer Research Centre – Queensland (APCRC-Q), Translational Research Institute (TRI), Brisbane, Australia. ³Biochemistry division, Chemistry Department, Faculty of Science, Damietta University, Damietta, Egypt. ⁴Quadram Institute Bioscience, Norwich Research Park, Norwich, UK. ⁵Mater Research Institute – University of Queensland, Translational Research Institute (TRI), Brisbane, Australia. ⁶Australian National Centre for the Public Awareness of Science, Australian National University, Canberra, Australia.

Received: 28 November 2017 Accepted: 1 May 2018 Published online: 24 May 2018

References

- Shiozawa Y, Pedersen EA, Havens AM, Jung Y, Mishra A, Joseph J, Kim JK, Patel LR, Ying C, Ziegler AM, Pienta MJ, Song J, Wang J, Loberg RD, Krebsbach PH, Pienta KJ, Taichman RS. Human prostate cancer metastases target the hematopoietic stem cell niche to establish footholds in mouse bone marrow. J Clin Invest. 2011;121(4):1298–312.
- Bubendorf L, Schöpfer A, Wagner U, Sauter G, Moch H, Willi N, Gasser TC, Mihatsch MJ. Metastatic patterns of prostate cancer: an autopsy study of 1,589 patients. Hum Pathol. 2000;31(5):578–83.
- 3. Jacobs SC. Spread of prostatic cancer to bone. Urology. 1983;21(4):337-44.
- Sakaguchi Y, Sekiya I, Yagishita K, Ichinose S, Shinomiya K, Muneta T. Suspended cells from trabecular bone by collagenase digestion become virtually identical to mesenchymal stem cells obtained from marrow aspirates. Blood. 2004;104:9.
- Bianco P, Riminucci M, Gronthos S, Robey PG. Bone marrow stromal stem cells: nature, biology, and potential applications. Stem Cells. 2001; 19(3):180–92.
- Corselli M, Chin CJ, Parekh C, Sahaghian A, Wang W, Ge S, Evseenko D, Wang X, Montelatici E, Lazzari L, Crooks GM, Péault B. Perivascular support of human hematopoietic stem/progenitor cells. Blood. 2013;121:15.
- Bliss SA, Sinha G, Sandiford O, Williams L, Engelberth DJ, Guiro K, Isenalumhe LL, Greco SJ, Ayer S, Bryan M, Kumar R, Ponzio N, Rameshwar P. Mesenchymal stem cell-derived exosomes stimulates cycling quiescence and early breast cancer dormancy in bone marrow. Cancer Res. 2016;1092:2016.
- Sieh S, Lubik AA, Clements JA, Nelson CC, Hutmacher DW. Interactions between human osteoblasts and prostate cancer cells in a novel 3D in vitro model. Organ. 2010;6(3):181.
- Tokuda Y, Satoh Y, Fujiyama C, Toda S, Sugihara H, Masaki Z. Prostate cancer cell growth is modulated by adipocyte-cancer cell interaction. BJU Int. 2003;91(7):716–20.
- de Barros APDN, Takiya CM, Garzoni LR, Leal-Ferreira ML, Dutra HS, Chiarini LB, Meirelles MN, Borojevic R, Rossi MID. Osteoblasts and bone marrow mesenchymal stromal cells control hematopoietic stem cell migration and proliferation in 3D in vitro model. PLoS One. 2010;5(2):e9093.
- Chambers KF, Mosaad EMO, Russell PJ, Clements JA, Doran MR. 3D cultures of prostate cancer cells cultured in a novel high-throughput culture platform are more resistant to chemotherapeutics compared to cells cultured in monolayer. PLoS One. 2014;9(11):e111029.
- 12. Xu X, Farach-Carson MC, Jia X. Three-dimensional in vitro tumor models for cancer research and drug evaluation. Biotechnol Adv. 2014;32(7):1256–68.
- 13. Kimlin LC, Casagrande G, Virador VM. In vitro three-dimensional (3D) models in cancer research: an update. Mol Carcinog. 2013;52(3):167–82.
- Thoma CR, Zimmermann M, Agarkova I, Kelm JM, Krek W. 3D cell culture systems modeling tumor growth determinants in cancer target discovery. Adv Drug Deliv Rev. 2014;69:29–41.
- Futrega K, Palmer JS, Kinney M, Lott WB, Ungrin MD, Zandstra PW, Doran MR. The microwell-mesh: a novel device and protocol for the high throughput manufacturing of cartilage microtissues. Biomaterials. 2015;62:1–12.
- Mosaad EO, Chambers KF, Futrega K, Clements JA, Doran MR. The microwell-mesh: a high-throughput 3D prostate cancer spheroid and drugtesting platform. Sci Rep. 2018;8(1):253.

- McMillin DW, Delmore J, Weisberg E, Negri JM, Geer DC, Klippel S, Mitsiades N, Schlossman RL, Munshi NC, Kung AL, Griffin JD, Richardson PG, Anderson KC, Mitsiades CS. Tumor cell-specific bioluminescence platform to identify stroma-induced changes to anticancer drug activity. Nat Med. 2010;16(4):483–9.
- Thaimann N, Edmund E, Hopwood VL, Pathak S, Von Eschenbach A, Chung LWK, Thalmann GN, Anezinis PE, Chang SM, Zhau HE, Kim EE, Hopwood VL, Pathak S, von Eschenbach AC, Chung LWK. Androgen-independent Cancer progression and bone metastasis in the LNCaP model of human prostate Cancer. Cancer Res. 1994;54(10):2577–82.
- Markway BD, Tan G-K, Brooke G, Hudson JE, Cooper-White JJ, Doran MR. Enhanced chondrogenic differentiation of human bone marrow-derived mesenchymal stem cells in low oxygen environment micropellet cultures. Cell Transplant. 2010;19(1):29–42.
- Dominici M, Le Blanc K, Mueller I, Slaper-Cortenbach I, Marini F, Krause D, Deans R, Keating A, Prockop D, Horwitz E. Minimal criteria for defining multipotent mesenchymal stromal cells. The International Society for Cellular Therapy position statement. Cytotherapy. 2006;8(4):315–7.
- Ungrin MD, Joshi C, Nica A, Bauwens C, Zandstra PW. Reproducible, ultra high-throughput formation of multicellular organization from single cell suspension-derived human embryonic stem cell aggregates. PLoS One. 2008;3(2):e1565.
- Li Z, Alyamani M, Li J, Rogacki K, Abazeed M, Upadhyay SK, Balk SP, Taplin M-E, Auchus RJ, Sharifi N. Redirecting abiraterone metabolism to fine-tune prostate cancer anti-androgen therapy. Nature. 2016;533(7604):547–51.
- 23. Arcangeli S, Zilli T, De Bari B, Alongi F. Hit the primary: a paradigm shift in the treatment of metastatic prostate cancer? Crit Rev Oncol Hematol. 2016;97:231–7.
- Decker AM, Jung Y, Cackowski F, Taichman RS. The role of hematopoietic stem cell niche in prostate cancer bone metastasis. J Bone Oncol. 2016;5(3): 117–20.
- Mimeault M, Batra SK. New concepts on the critical functions of Cancerand metastasis-initiating cells in treatment resistance and disease relapse: molecular mechanisms, signaling transduction elements and novel targeting therapies. Dordrecht: Springer; 2010. p. 175–207.
- Ren G, Esposito M, Kang Y. Bone metastasis and the metastatic niche. J Mol Med. 2015;93(11):1203–12.
- Mimeault M, Hauke R, Batra S. Recent advances on the molecular mechanisms involved in the drug resistance of Cancer cells and novel targeting therapies. Clin Pharmacol Ther. 2008;83(5):673–91.
- Tuxhorn JA, McAlhany SJ, Dang TD, Ayala GE, Rowley DR. Stromal cells promote angiogenesis and growth of human prostate tumors in a differential reactive stroma (DRS) xenograft model. Cancer Res. 2002;62:11.
- Brennen WN, Denmeade SR, Isaacs JT. Mesenchymal stem cells as a vector for the inflammatory prostate microenvironment. Endocr Relat Cancer. 2013; 20(5):269–90.
- Luo J, Ok Lee S, Liang L, Huang C-K, Li L, Wen S, Chang C. Infiltrating bone marrow mesenchymal stem cells increase prostate cancer stem cell population and metastatic ability via secreting cytokines to suppress androgen receptor signaling. Oncogene. 2014;33(21):2768–78.
- YlÖstalo JH, Bartosh TJ, Coble K, Prockop DJ. Human mesenchymal stem/ stromal cells cultured as spheroids are self-activated to produce prostaglandin E2 that directs stimulated macrophages into an antiinflammatory phenotype. Stem Cells. 2012;30(10):2283–96.
- Pinho S, Lacombe J, Hanoun M, Mizoguchi T, Bruns I, Kunisaki Y, Frenette PS. PDGFRa and CD51 mark human nestin+ sphere-forming mesenchymal stem cells capable of hematopoietic progenitor cell expansion. J Exp Med. 2013;210:7.
- 33. Gao D, Vela I, Sboner A, Iaquinta PJ, Karthaus WR, Gopalan A, Dowling C, Wanjala JN, Undvall EA, Arora VK, Wongvipat J, Kossai M, Ramazanoglu S, Barboza LP, Di W, Cao Z, Zhang QF, Sirota I, Ran L, MacDonald TY, Beltran H, Mosquera J-M, Touijer KA, Scardino PT, Laudone VP, Curtis KR, Rathkopf DE, Morris MJ, Danila DC, Slovin SF, Solomon SB, Eastham JA, Chi P, Carver B, Rubin MA, Scher HI, Clevers H, Sawyers CL, Chen Y. Organoid cultures derived from patients with advanced prostate Cancer. Cell. 2014;159(1):176–87.
- Yang N, Yan T, Zhu H, Liang X, Leiss L, Sakariassen P, Skaftnesmo K, Huang B, Costea D, Enger P, Li X, Wang J. A co-culture model with brain tumor-specific bioluminescence demonstrates astrocytes-induced drug resistance in glioblastoma. J Transl Med. 2014;12(1):278.

- Straussman R, Morikawa T, Shee K, Barzily-Rokni M, Qian ZR, Du J, Davis A, Mongare MM, Gould J, Frederick DT, Cooper ZA, Chapman PB, Solit DB, Ribas A, Lo RS, Flaherty KT, Ogino S, Wargo JA, Golub TR. Tumour microenvironment elicits innate resistance to RAF inhibitors through HGF secretion. Nature. 2012;487(7408):500–4.
- Ambrosi D, Mollica F. The role of stress in the growth of a multicell spheroid. J Math Biol. 2004;48(5):477–99.
- Rhee S, Grinnell F. Fibroblast mechanics in 3D collagen matrices. Adv Drug Deliv Rev. 2007;59(13):1299–305.

Ready to submit your research? Choose BMC and benefit from:

- fast, convenient online submission
- thorough peer review by experienced researchers in your field
- rapid publication on acceptance
- support for research data, including large and complex data types
- gold Open Access which fosters wider collaboration and increased citations
- maximum visibility for your research: over 100M website views per year

At BMC, research is always in progress.

Learn more biomedcentral.com/submissions

