## **Embalming Powder**

## Sanisorb

A must have in any embalming room, Dodge Sanisorb can be used in any situation where a dependable absorptive product with a pleasant scent is required.

#### Sanisorb

#### **Performance Highlights:**

- Absorbs copious amount of moisture.
- Perfect inside plastic garments where leakage or oedema may cause problems.
- Excellent for use in coffins of repatriation cases.
- Use in coffins stored in vaults or mausoleums.
- Cleans up any embalming room spills instantly.
- Absorbs moisture in the cavities of autopsy cases as well as inside bone donation incisions.
- Uses for Sanisorb are limited only by your imagination.



Ordering Data b. Sanisorb Catalogue No: 530450 Unit Size: 500gm c. Catalogue No: 530451 Unit Size: 2kgs

## Q-S Powder



### Ordering Data e. Q-S Powder Catalogue No: 560111

Unit Size: 500gm Bottle

Catalogue No: **560140** Unit Size: 3.5kgs

#### Q-S Powder Performance Highlights

- Forms powerful, positive and dependable closures.
- Contains highly effective deodorising and preserving chemicals.
- Pleasant to use, does not create dust in the air.
- Will not cake or turn gummy, Q-S Powder lasts for years if container is kept closed.
- Excellent for incisions, wounds, open lesions, suture needle holes, body orifices and autopsy cases.
- Works virtually on contact. No waiting around for a seal to develop.

A low cost, high bulking absorbent compound that can be used for drying and filling autopsied cases. Hardening Compound delivers limited preservation and helps to supress offensive odours. Hardening Compound will not cake or solidify in the drum and is simple and economical in application.

### Hardening Compound



Paraformaldehyde Prills



A convenient way of applying Paraformaldehyde to the surfaces or cavities of a body in a relatively dust free manner. Paraformaldehyde Prills can be used in lieu of Hypo work on autopsies and difficult cases, or applied to bodies where decomposition has rendered arterial injection impossible. Paraformaldehyde Prills fume heavily and care should always be taken to apply in a well ventilated area.

#### Ordering Date

a. Paraformaldehyde Prills Catalogue No: 910205 Unit Size: 3.5kg

**Ordering Data** 

Unit Size: 3kgs

d. Hardening Compound Catalogue No: 530030

R

# Embalming Powder

## **Action Powder**



A real must have for any and every embalming room. Dodge Action powder delivers fast dependable preservation to areas not receiving fluid through arterial injection. Action powder may be used freely on autopsy cases for treating the trunk walls or cranial vault. Action may also be used for cases requiring deodorising and preservative treatment which are beyond arterial injection.

#### Action Embalming Powder Performance Highlights:

- Assured performance on cases where arterial injection could be difficult if not impossible.
- Excellent inside caskets, grave vaults, mausoleums and crypts.
- Very useful for treating bed sores, cancerous or putrefactive areas.
- Effective in lieu of hypo work on wet or dry decomposed cases or autopsy cases where trunk walls do not receive adequate fluid during injection.
  Use with plastic undergarments to preserve and dry oedematous limbs.

Ordering Date a. Action Embalming Powder Catalogue No: 536098 Unit Size: 4Kg

- Useful for treating infant cases where traditional treatments would be impracticable.
- Used by many for the treatment of exhumed or disinterred cases.

Combining the properties of a preservative powder, hardening compound and absorptive in one, Viscerock Plus is a true time saver in the embalming room. Great for use by fume sensitive embalmers, Viscerock Plus is anti-gassing, free flowing, it will not cake or congeal while stored, absorbs 370% it's own weight of fluid and adheres well to tissue surfaces.

## VP (Viscerock Plus )

b.

#### **Ordering Data**

**b. VP (Viscerock Plus)** Catalogue No: **535012** Unit Size: 2.25kg

#### VP (Viscerock Plus) Performance Highlights:

- Combines the powerful preservative attributes of Dodge Viscerock with the absorptive qualities of Sanisorb.
- Easy to use and store in handy 2.25kg Tub.
- Adheres better to tissue surfaces.
- Pleasant scent overcomes decomposition odours.
- Great for use of oedematous or leaking limbs with plastic undergarments.

C.

## Viscerock FF (Formaldehyde Free)

New from the Dodge Company, Viscerock FF (Formaldehyde Free). This product is nearly dustless, free flowing and will not cake in any environment. Viscerock FF will dry and firm moist tissues and adhere to all tissue surfaces. It can be used in plastic garments to dry oedematous limbs and inside the cavities of autopsied cases. Pleasant in use, Viscerock FF is great for use by fume sensitive embalmers and will create no unpleasant odours during use.



Ordering Data

c. Viscerock FF (formaldehyde Free) Catalogue No: 535115 Unit Size: 2.25kg

As with all products containing Formaldehyde, the operator should be sure to use them in a well ventilated area. For more difficult situations, please see our range of Formaldehyde respirators. 29

## **Absorbent Products**

#### Dodgesorb Absorbent Sheets Performance Highlights:

- Each sheet consists of three layers: a top fabric layer, a Dodgesorb- impregnated middle layer and a plastic bottom layer.
- Even after exposure to moisture, they retain integrity and can be removed from their area of use.
- White Sheet 25" x 72" (A)
- For use over or under the body.
- Absorbs approximately 4.6 litres of fluid.
- White Sheet 63" x 90" (B)
- For use in wrapping the entire body.
- Absorbs approximately 6.2 litres of fluid.

## **DodgeSorb Absorbent Sheets**



#### **Ordering Data:**

a. Dodgesorb Absorbent Sheets (A) White Sheet 25" x 72" Catalogue No: 038345 Unit Size: Five sheets per pack



Ordering Data: b. Dodgesorb Absorbent Sheets (B) White Sheet 63" x 72" Catalogue No: 038340 Unit Size: Five sheets per pack

### **DodgeSorb Absorbent Liners**



#### Ordering Data:

c. Dodgesorb Absorbent Liners White Sheet 25" x 86" Catalogue No: 038375 Unit Size: Five sheets per pack

#### **Dodgesorb Absorbent Liner Performance Highlights:**

- These liners consist of an exposed layer of Dodgesorb impregnated material, which will turn to loose gel in the presence of moisture.
- Once in place, this sheet must remain there.
- Ideal insurance for oedematous cases or those being repatriated by road, air or sea.
- White Sheet 25" x 86"
- For use inside coffins and caskets being repatriated or entombed.
- Absorbs approximately 5 litres of fluid.

#### DodgeSorb Absorbent Inserts Performance Highlights:

These are smaller absorbent pads, with a DodgeSorb layer inside two fabric layers.



#### Small Insert 3.25"x 10"

- Absorbs approximately 232mls of fluid.
- Used for packing cranial, thoracic and tracheal stomata.

#### **Ordering Data:**

#### d. DodgeSorb Absorbent Inserts

Small Insert 3.25" x 10" Catalogue No: **038330** Unit Size: Ten per case



#### Medium Insert 3.25"x 14"

- Absorbs approximately 478mls of fluid.
- A tub like pad through which a pipe prosthetic can be inserted, or use

#### **Ordering Data:**

e. DodgeSorb Absorbent Inserts Medium Insert 3.25" x 14" Catalogue No: 038325 Unit Size: Ten per case

## **DodgeSorb Absorbent Inserts**



#### Large Insert 3.25"x 37.5"

- Absorbs approximately 1000mls of fluid.
- Also tube like, fits around long bone prosthetics.

#### **Ordering Data:**

**f. DodgeSorb Absorbent Inserts** Large Insert 3.25" x 37.5" Catalogue No: **038320** Unit Size: Ten per case

Α

S

## **Absorbent Products**

### **DodgeSorb Garments**

## **DodgeSorb Garments**

#### **Performance Highlights:**

Fluid impervious white plastic sleeves and stockings have an inner liner that contains a layer of DodgeSorb to provide added protection by capturing and containing fluids.

b.

**Ordering Data:** 

**Stockings** 

b. DodgeSorb Absorbent

Catalogue No: 760885

Unit Size: Two per bag

All openings are elasticated to provide a good fit and reduce leakage.



a. DodgeSorb Absorbent Sleeves Catalogue No: 7608825 Unit Size: Two per bag



**DodgeSorb Cranial Caps** 

## **DodgeSorb Cranial Caps**

#### **Performance Highlights:**


- Cranial caps are designed to absorb seepage due to autopsied or mutilated heads.
- Inside the cap is absorbent, while the outside is fluid impervious.
- Conveniently conforms to the shape of the head.
- Slide under the head and the cap sinks into the coffin or casket pillow.

### **Ordering Data:**

c. DodgeSorb Cranial Caps Catalogue No: 750594 Unit Size: Twelve per bag



### DodgeSorb Kit



#### **DodgeSorb Absorbent Kit Performance Highlights:**

- A convenient way of keeping the range of DodgeSorb items in the embalming room.
- Each kit contains five small inserts, five medium inserts, two large inserts, one white liner and one cranial cap.

#### **Ordering Data:**

d. DodgeSorb Absorbent Kit Catalogue No: 038380

31

## **Absorbent Products**

## Superabsorbent Sachets

a.



Sold in boxes of 50 individual sachets, this absorbent product can be used to seal leaking incisions, pack the nose and throat of the deceased and even to close anal and vaginal orifices. Dodge Superabsorbent is applied in the same way as cotton wool or cellulose wadding. This product will absorb a remarkable amount of fluid and form a leak proof seal.

Ordering Data: a. Dodge Superabsorbent Sachets Catalogue No: 750991 Unit Size: 50 Sachets per box

## Superabsorbent Cranium Caps



#### Sold individually, these caps are made especially to deal with any leakage from the cranium of the deceased due to trauma or post mortem examination before and after embalming. The cap is simply placed between the head of the deceased and the coffin or casket pillow. Containing the same ingredients as the Superabsorbent sachet, any leakage is captured and solidified to prevent staining of the interior of the coffin or casket.

#### **Ordering Data:**

**b. Dodge Superabsorbent Cranium Caps** Catalogue No: **750597** Unit Size: Sold Individually

## **Cranial Caps**



These caps are easy to conceal and can be adjusted by the user to suit individual cases. Designed to divert leakage from the back of the head, these caps help to prevent staining of the deceased clothing and the interior of the coffin or casket.

Ordering Data: c. Dodge Cranial Caps Catalogue No: 750596 Unit Size: Sold Individually

32